 IMMULITE[®]
2000

Digoxin

DPC[®]

IMMULITE® 2000 Digoxin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of digoxin in serum or heparinized plasma, as an aid in monitoring the therapeutic administration of this cardioglycoside, while avoiding toxicity.

Catalog Number: **L2KDI2** (200 tests),
L2KDI6 (600 tests)

Test Code: **DGX** Color: **Light Green**

Summary and Explanation

Immunoassay has proved particularly useful in determining drug overdoses in patients treated with cardiac glycosides. The technique is also useful (a) in clarifying situations where a patient's symptoms might be due either to intrinsic heart disease or to digitalis intoxication; (b) where there is doubt concerning the type of digitalis preparation the patient is taking — in this case, a digitoxin immunoassay is also necessary; (c) for measuring the digoxin ingestion of patients with an inadequate history of previous dosage; (d) in documenting cases of underdigitalization as well as digitalis (digoxin) excess; (e) in monitoring the toxic response in patients with myocardial disease associated with hypokalemia, hypomagnesemia, hypercalcemia, hypoxia and alkalosis, which are particularly sensitive to digitalis; and (f) in preventing overdigitalization, particularly in patients whose renal function is deteriorating or for whom an increased digoxin dosage is contemplated.

The high sensitivity of digoxin immunoassay is especially necessary in view of the small differences and occasional overlap that exist between therapeutic and toxic levels of circulating digoxin.¹³ (Intoxication is defined in terms of arrhythmias and disturbances of cardiac conduction due to the drug's presence.) Smith et al reported serum digoxin concentrations of 0.8–2.4 ng/mL (1.0–3.1 nmol/L) in nontoxic patients, and 2.1–

8.7 ng/mL (2.7–11.1 nmol/L) in toxic patients, based on data taken six hours post-dose in patients with normal renal function.⁸ More recent clinical investigations have confirmed the association of toxicity with serum levels above 2 ng/mL.

Principle of the Procedure

IMMULITE 2000 Digoxin is a solid-phase, enzyme-labeled, competitive chemiluminescent immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

Collect blood in plain tubes without gel barriers, noting the time of collection, just before the next dose — or, in any case, no less than six hours after the administration of this cardioglycoside, to avoid misleadingly high digoxin results.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Digoxin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum or plasma.

Storage: 7 days at 2–8°C, or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Digoxin Bead Pack (L2DI12)

With barcode. 200 beads, coated with monoclonal murine anti-digoxin. Stable at 2–8°C until expiration date.

L2KDI2: 1 pack. **L2KDI6:** 3 packs.

Digoxin Reagent Wedge (L2DIA2)

With barcode. Two reagents: 5.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to digoxin, in buffer, and 11.5 mL of buffer. Stable at 2–8°C until expiration date.

L2KDI2: 1 wedge. **L2KDI6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Digoxin Adjustors (L2DIJ3, L2DIJ4)

Two vials (Low and High), 2 mL each, of digoxin in human serum matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KDI2: 1 set. **L2KDI6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

SDCM: Bi-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of digoxin.

Therapeutic Range

The following results are reported in the literature, based on data taken six hours post-dose in patients with normal renal function.⁸

Group	Digoxin, ng/mL	Digoxin, nmol/L
Nontoxic patients	0.8 – 2.4	1.0 – 3.1
Toxic patients	2.1 – 8.7	2.7 – 11.1

Consider these limits as *guidelines* only.

Limitations

Sole reliance on digoxin concentrations for determining digitalis toxicity is not warranted, and must be supplemented with additional clinical and electrocardiographic information.

Many physiological, pharmacological, pathological and genetic factors may affect the interpretation of results (see Summary and Explanation of the Test section).

The presence of endogenous digoxin-like immunoreactive factors (DLIF) has been reported in the serum of patients with renal and hepatic failure, and in newborns and women in the third trimester of pregnancy.¹⁴ These factors can result in falsely elevated digoxin results in a

number of commercially available immunoassays.

In cases of overdosage, specimens should be obtained *before* administering digoxin immune Fab antidote preparations (e.g. Digibind[®]) since these can be expected to interfere with any digoxin immunoassay procedure. Once such a preparation has been administered, accurate determinations of digoxin can be made only by specialized techniques (not routinely available in most clinical settings) until the Fab fragments have been eliminated from the body: this takes several days for patients with intact renal function, and a week or longer for patients with renal impairment. In addition to *in vitro* effects, these preparations, administered intravenously, cause a rapid rise in pharmacologically inactive, Fab-bound, total circulating digoxin that may be detected immunologically.¹⁵

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
ng/mL × 1.281 → nmol/L

Calibration Range: 0.5 to 8 ng/mL
(0.6 to 10 nmol/L).

Analytical Sensitivity: 0.2 ng/mL
(0.3 nmol/L).

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three digoxin solutions (15, 30 and 60 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for digoxin, though Lanatoside C and Methyl digoxin do exhibit significant crossreactivity. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemoglobin: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Matched serum and heparinized samples were collected from 28 volunteers. Varying amounts of digoxin were added to the matched samples, which were analyzed by the IMMULITE 2000 Digoxin procedure. By linear regression:

(Na Heparin) = 0.96 (Serum) – 0.03 ng/mL
r = 0.988

Means:
2.75 ng/mL (Serum)
2.61 ng/mL (Na Heparin)

The results show that heparinized plasma yields essentially equivalent results to serum.

Method Comparison 1: The assay was compared to IMMULITE Digoxin on 97 samples. (Concentration range: approximately 0.5 to 6.1 ng/mL. See graph 1.) By linear regression:

(IML 2000) = 0.94 (IML) + 0.097 ng/mL
r = 0.977

Means:
1.4 ng/mL (IMMULITE 2000)
1.4 ng/mL (IMMULITE)

Method Comparison 2: The assay was compared to DPC's Coat-A-Count (RIA) Digoxin on 168 samples. (Concentration range: approximately 0.5 to 6.0 ng/mL. See graph 2.) By linear regression:

$$(IML\ 2000) = 0.99 (CAC) + 0.004\ ng/mL$$

$$r = 0.975$$

Means:
1.3 ng/mL (IMMULITE 2000)
1.3 ng/mL (Coat-A-Count)

Method Comparison 3: The assay was compared to another automated immunoassay (Kit A) on 169 samples. (Concentration range: approximately 0.5 to 6.2 ng/mL. See graph 3.) By linear regression:

$$(IML\ 2000) = 0.96 (Kit\ A) - 0.09\ ng/mL$$

$$r = 0.971$$

Means:
1.3 ng/mL (IMMULITE 2000)
1.5 ng/mL (Kit A)

References

1) Butler VP, Jr., Lindenbaum J. Serum digitalis measurements in the assessment of digitalis resistance and sensitivity. *Am J Med* 1975;58:460-9. 2) Doherty JE. How and when to use the digitalis serum levels. *J Am Med Assoc* 1978;239:2594-6. 3) Gault MH, et al. Studies of digoxin dosage, kinetics and serum concentrations in renal failure and review of the literature. *Nephron* 1976;17:161-87. 4) Huffman DH, et al. Association between clinical cardiac status, laboratory parameters, and digoxin usage. *Am Heart J* 1976;91:28-34. 5) Rietbrock N, Alken RG. The therapy of heart failure with digitalis. *Deutsche Med Wochenschr* 1980;105:1622-8. 6) Shapiro B, et al. Pitfalls in the application of digoxin determinations. *Sem Nucl Med* 1975;5:205-20. 7) Shapiro W. Current considerations in digoxin usage. *CRC Crit Rev Clin Lab Sci* 1978;9:321-46. 8) Smith TW, et al. Determination of therapeutic and toxic serum digoxin concentrations by radioimmunoassay. *N Engl J Med* 1969;281:1212-16. 9) Smith TW. Medical intelligence. *N Engl J Med* 1973;288:719-22 and 942-6. 10) Smith TW. Digitalis toxicity: epidemiology and clinical use of serum concentration measurements. *Am J Med* 1975;58:470-6. 11) Weintraub M. Interpretation of the serum digoxin concentration. *Clin Pharmacokinetics* 1977;2:205-19. 12) Wettrell G. Digoxin therapy in infants. *Acta Paediatr Scand* 1976; Suppl 257:7-28. 13) Hallworth MJ. Practical therapeutic drug monitoring: digoxin. *Lab Medica* 1991 Mar/Apr;19-22. 14) Valdes R Jr. Endogenous digoxin-like immunoreactive factors: impact on digoxin measurements and

potential physiological implications. *Clin Chem* 1985;9:1525-32. 15) Digoxin immune Fab. In: McEvoy GK, editor. *AHFS Drug Information* 92. American Society of Hospital Pharmacists, 1992: 1977-81.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	0.6	0.05	8.3%	0.05	8.3%	
2	1.2	0.06	5.0%	0.07	5.8%	
3	2.7	0.12	4.4%	0.15	5.6%	
4	3.3	0.14	4.2%	0.14	4.2%	
5	6.6	0.31	4.7%	0.35	5.3%	

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	4 in 4 ⁵	2.1	—	—
	2 in 4	1.1	1.1	100%
	1 in 4	0.5	0.5	100%
2	8 in 8	4.2	—	—
	4 in 8	2.0	2.1	95%
	2 in 8	1.0	1.1	91%
	1 in 8	0.5	0.5	100%
3	8 in 8	4.2	—	—
	4 in 8	1.8	2.1	86%
	2 in 8	0.9	1.1	82%
	1 in 8	0.5	0.5	100%
4	8 in 8	6.4	—	—
	4 in 8	3.2	3.2	100%
	2 in 8	1.5	1.6	94%
	1 in 8	0.7	0.8	88%

5	8 in 8	7.3	—	—
	4 in 8	3.5	3.7	95%
	2 in 8	1.7	1.8	94%
	1 in 8	0.9	0.9	100%

Recovery (ng/mL)




	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	ND	—	—
	A	0.7	0.8	88%
	B	1.5	1.5	100%
	C	2.9	3.0	97%
2	—	<0.5	—	—
	A	0.8	0.8	100%
	B	1.5	1.5	100%
	C	2.8	3.0	93%
3	—	1.0	—	—
	A	1.7	1.7	100%
	B	2.4	2.5	96%
	C	3.7	4.0	93%
4	—	2.0	—	—
	A	2.8	2.7	104%
	B	3.3	3.4	97%
	C	4.7	4.9	96%
5	—	3.2	—	—
	A	3.3	3.8	87%
	B	4.1	4.5	91%
	C	5.6	6.0	93%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
Diagnostic Products Corporation
DPC Biermann GmbH