
T: 02 8212 4074 F: 02 9423 6992
 www.keydiagnostics.com.au
 E: infoxkeydiagnostics.com.au
 PO Box 1038 Gymea NSW 2227

**REF FSK3 MICROPOLYSPORA FAENI (FARMER'S LUNG)
 IMMUNODIFFUSION SYSTEM**

INTENDED USE

The Microgen Bioproducts *Micropolyspora faeni* (*Faeni rectivirgula*) Immunodiffusion System is a serological double gel diffusion assay intended for the detection of precipitating antibodies to the thermophilic actinomycete *M. faeni*, which is the main sensitising agent found in mouldy hay (or other cereal) and responsible for farmer's lung - a pulmonary hypersensitivity reaction to inhaled mouldy hay dust. The kit is intended for professional laboratory use only.

PRINCIPLE OF THE TEST

The precipitating antibodies detected by the kits are generated as part of the immune response to exposure to the antigens of *Micropolyspora faeni*. These antibodies are predominantly IgG class, although other classes may have a precipitating function. The test relies on the principle of double gel diffusion (Ouchterlony). When soluble antigens and homologous antibodies are placed in adjacent wells cut into suitable diffusion media such as agar or agarose, they diffuse towards one another and produce visible precipitation lines along the interface of optimal relative concentrations. Control sera containing suitable titres of the relevant antibodies are tested in wells adjacent to the antigens being reacted with the patient's serum. The development of precipitin lines in control tests is used to confirm reagent and plate performance. With experience, the controls can be used as markers indicating completion of diffusion. By utilising two antigen concentrations (20mg/mL and 2mg/mL) and careful selection of the method of preparation of antigenic extracts, it is believed that the maximum number of all positive sera will be identified (1,3). Consequently, the diagnosis of systemic infection or of hypersensitivity reactions should be made considerably easier. The presence of antibodies resulting from previous exposure to the causative agents can often be differentiated from current infection since the latter is typically characterised by high or rapidly rising titres over a period of time.

CONT		KIT PRESENTATION	
CONTROL	-	FSK/a Negative Control: Sheep serum	2.5mL
CONTROL	+	FSK3/b Positive Control: Sheep anti- <i>Micropolyspora faeni</i> antiserum	2.5mL
CF AG	20	FSK3/c <i>Micropolyspora faeni</i> culture filtrate antigens (20mg/mL):	0.5mL
CF AG	2	FSK3/d <i>Micropolyspora faeni</i> culture filtrate antigens (2mg/mL):	0.5mL

All reagents contain 0.099% sodium azide as a preservative.

Instructions for Use
 Record Cards (20)

Additional Requirements:

- Double diffusion plates (Packs of 20 ready-to-use plates are available from Microgen Bioproducts - Product Code FSK/DDP)

WARNINGS AND PRECAUTIONS

Safety:

1. The reagents supplied in this kit are for *in vitro* diagnostic use only
2. Sodium azide, which is used as a preservative in the kit reagents can react with lead or copper plumbing to form potentially explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build-up.
3. Appropriate precautions should be taken when handling or disposing of potentially pathogenic samples. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30 minutes. Liquid waste containing acid must be neutralised before treatment.
4. The antigen preparations have been inactivated during the manufacturing process. However, they should be handled as though potentially infectious.

Procedural:

1. The kit should be used according to these instructions.
2. During use, the reagents should be kept at ambient temperature for the minimum time possible.
3. Do not dilute any of the kit reagents
4. Do not intermix reagents from different batches of kits.
5. Do not freeze any of the kit reagents

STORAGE AND SHELF LIFE

The *M. faeni* Immunodiffusion System should be stored firmly stoppered at 2-8°C when not in use. Do not freeze. Any slight turbidity in reagents, particularly controls, will not normally interfere with satisfactory performance. The kit should not be used after the expiry date printed on the carton label.

SPECIMENS

5-10mL of venous blood should be collected and allowed to clot in a glass tube without anti-coagulant. Remove the serum and store at 2-8°C if the test is to be performed within 4-5 days. If longer storage is necessary, freeze the serum at -20°C to -70°C. Repeated freezing and thawing causes loss of precipitating antibody. Gradual decline in antibody titre is also to be expected if sera are stored at -20°C for longer than 6 months.

DOUBLE GEL DIFFUSION PLATES (Product Code FSK/DDP)

Each plastic envelope holds a plastic box that divides into two compartments separated by a foam layer. In each compartment, there is a glass slide coated with agar in which wells have been pre-cut. After use (see Test Procedure below), the stained, dried plates are suitable for filing as permanent records.

Chemical composition

Plates are prepared with 1.5% purified agar (Difco) in McIlvaine's citric acid phosphate buffer, pH 7.0 containing 0.099% sodium azide.

Dimensions

The slide is ca. 51x76mm and the agar gel layer is nominally 2mm thick. The large patient's serum well is 12.5mm in diameter with a nominal capacity of 200 μ L. The outermost antiserum control wells are 8mm diameter and will hold 100 μ L. The small 4mm antigen wells have a capacity of 25 μ L.

Filling the wells

Many micro-pipetting devices of both fixed and adjustable volume are commercially available and are recommended for loading the wells according to the volumes specified above. A fresh pipette or tip is essential for each serum specimen or antigenic extract. The large serum well size compared to the antigen well is to increase the amount of antibody available for reaction and is optimal for detection of weak precipitin reactions using unconcentrated patient's serum.

TEST PROCEDURE

Immunodiffusion

1. Bring the control sera and antigen extracts to room temperature. Gently mix any patient's serum that has been frozen.
2. Fill in a test record card with patient's data, number, wells to be used, contents, etc.
3. Remove the plastic box from the polythene envelope. Unhinge the box by gently twisting to form the two compartments. Remove the polythene foam layer if both compartments are needed, otherwise store the unused compartment (with foam) in the resealable envelope provided.
4. Paint the specimen number on the agar layer using indelible stain (e.g. 0.5% Alcian Blue) applied with a fine brush. The dye must be one that will not diffuse.
5. Fill the smallest wells with antigenic extracts. Both concentrations of the extract must be tested separately. The patient's serum is placed in the large well and appropriate control sera in the outer wells. (Sufficient control sera are included to fill 25 x positive control wells and 25 x negative control wells.)
6. Place each compartment in a moisture chamber on a level surface and incubate at 35-37°C for 3 days or room temperature for 5 days. Plates should be examined daily for signs of precipitin lines. N.B. Resealable plastic envelopes are quite convenient, if moistened, for use as moisture chambers.

Washing and staining slides

7. Cut around the glass slide in each compartment with a scalpel blade or knife and raise up the glass agar plate. Place the slide in a dish or large beaker of saline (1% sodium chloride + 0.099% sodium azide). Wash the plate for at least 24 hours and preferably 48 hours changing the saline wash as frequently as is practicable.
8. Finally wash for 2-3 hours in distilled water.
9. Remove the slide, cover with filter paper and air dry at 37°C to reduce the gel to a thin film. Drying may be accelerated by careful use of a hair dryer. If necessary, moisten the filter paper to remove it. Re-dry the slide briefly afterwards.
10. Stain the slide in a dish of 0.5% Coomassie Blue B.L. (or equivalent) in methanol/acetic acid/water (5:1:4v/v) for 10 minutes.
11. Destain the slide in several changes of the same solvent for 15 minutes. Air-dry, record the results and file the slide.

N.B. Some users have found that a final 15 minute rinse of the slide in a 5-10% glycerol-saline solution enhances long-term preservation.

Reading and recording results

The unstained plates are best examined by viewing against a black background using transmitted light. Look for precipitin reaction lines in the area between the antigen wells and the patient's serum well. The positive control can be used as a marker for precipitin line development. The appearance of lines can be provisionally noted on the record card when diffusion appears complete but maximum sensitivity is obtained after full washing and staining.

INTERPRETATION

Precipitin lines should be visible between the positive control wells and the antigen wells. The presence of precipitin lines between the patient's sample and the antigen wells is indicative of a positive result. No test should be reported as negative until the slide has been inspected following staining and destaining as faint precipitin lines may not be visible until this procedure has been completed. When the serum contains a large amount of antibody, the precipitin line will be close to the antigen well. Where there is very little antibody, the precipitin reaction is close to the patient's serum well. When more than one antibody is present, a series of precipitin lines appear at the 20mg/mL and 2mg/mL antigen levels; these can be used as a semi-quantitative guide to titre. In addition, it is useful to record whether a line is sharp or broad and fuzzy.

The precipitin titre of a positive serum may be determined by performing the double diffusion test with serial doubling dilutions of the sample.

Clinical significance

Precipitins against *M. faeni* are regarded as an index of exposure to this organism in response to inhalation of mouldy hay dust. Such precipitins indicate a potential to respond to subsequent dust inhalation by immune complex production. There is evidence that the degree of severity of pulmonary involvement is linked with a high titre antibody response (1,2).

The precipitins are useful to help distinguish farmer's lung from diseases that sometimes present a similar course, e.g. tuberculosis, influenza, brucellosis, mycoplasma infection (1). Occasionally, serum from cases of *Mycoplasma pneumoniae* cross-react with *M. faeni* antigenic extracts (4).

See Bibliography below for further information on correlation of precipitin line formation with *M. faeni* exposure.

LIMITATIONS OF USE

1. Results of a single test procedure should not be relied upon for diagnostic purposes but should be interpreted by the clinician in the context of all available clinical and laboratory information. Detection of precipitating antibodies alone does not constitute a diagnosis of active disease caused by *M. faeni*.
2. The presence of low levels of precipitins may simply indicate past or present exposure to *M. faeni* with no incidence of disease.
3. The absence of precipitins does not preclude the presence of disease caused by *M. faeni* (1,5)
4. Serum from cases of *Mycoplasma pneumoniae* occasionally cross-react with *M. faeni* antigenic extracts (4).

PERFORMANCE CHARACTERISTICS

There is no reference test procedure for precipitating antibodies to *M. faeni* with which to compare and evaluate FSK3.

251 serum samples have been tested with FSK3 and an enzyme immunoassay for IgG to *M. faeni* but correlation between the two sets of data is limited due to the following factors:

- Precipitating antibodies reacting in FSK3 may include classes other than IgG that will not be detected by the EIA.
- Some IgG antibodies detected by EIA may be directed against antigens not involved in the precipitating reaction. E.g. some low molecular weight antigens do not form precipitin bands.
- Not all IgG antibodies are precipitating.
- The EIA may not detect IgG to some antigens if they do not readily adsorb on to the surface of the microtitration plate.
- The EIA is generally more sensitive and will detect low levels of IgG. However, these levels may not be clinically significant. Higher levels of precipitating antibodies correlate with the presence of active disease.
- Diagnostically important levels of precipitating antibodies are not always detected by EIA. Similarly, relatively high levels of IgG in some sera may not be detected by FSK3.

The results of the comparative study (shown below) should therefore be reviewed and interpreted in the light of these factors. There are a number of discrepant results which may not have clinical significance. This emphasises the "Limitations of Use" listed above which state that diagnosis of disease can only be made after considering the results of a range of tests and clinical symptoms.

		FSK3				Total
		+++	++	+	-	
IgG EIA	+++	2	0	0	1	3
	++	0	4	0	2	6
	+	2	8	1	4	15
	-	0	7	7	213	227
	Total	4	19	8	220	251

Symbols: +++ Strong positive reaction
 ++ Moderate positive reaction
 + Weak positive reaction
 - Negative reaction

Specificity is fairly high (213/227 = 93.8%) but sensitivity of FSK3 in comparison with EIA is relatively low (17/24 = 70.8%). It is clear, however, from the above data that whilst some EIA positives (7) are not detected by FSK1, there are a number of samples (14) which are FSK3 positive, EIA negative.

REPRODUCIBILITY

Intra-batch reproducibility was established by testing the 2 positive controls and 1 negative control in one batch of product on three separate occasions. One batch of FSK DDP diffusion plates was used throughout. Whilst there was some variation in the number of precipitin lines visualised after staining, a clear positive reaction was seen on each occasion with both positive controls. No reactions with the negative control were seen.

Inter-batch reproducibility was examined by testing the 2 positive controls and 1 negative control in three batches of product. A clear positive reaction was seen with both positive controls in all batches although there was some variation in the number of precipitin lines generated in the reaction. No reactions were seen with the negative control.

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