

Fill out circle completely

- fill in the center completely
- if center is left open, the scanner will not read it.

Licensure Number

- leave first column blank
- fill in zero for next three columns
- licensure number from the right and the others with zeros
 - Gen - 1
 - Spec - 2
 - Lab Asst - 3
 - Tech - 4
 - Cyto - 5
 - Phleb - 6

CLP.201753-PHL

I.D. NUMBER
0 0 0 2 0 1 7 5 3

The licensure number in the example to the left doesn't have a letter in front. It is in the 200,000 series. Just fill in the number.

CLP.G01753-Gen

I.D. NUMBER
0 0 0 6 0 1 7 5 3

The example on the right is the has only 5 numbers preceded with a letter. The letter has a corresponding number that can be filled in.

Overview of Autoimmune and Drug Induced Hemolytic Anemias

Objectives

- Summarize the clinical findings associated with anemia.
- Define immune hemolytic anemia and indicate the types of antibodies involved.
- Identify and investigate potential serologic discrepancies caused by autoimmune and drug-induced hemolytic anemias.

Introduction to Anemia

- One of the most common problems encountered in clinical medicine
- Is not a disease
 - Expression of an underlying disorder or disease
 - Must determine what is causing the anemia



Introduction to Anemia

- Defined as a decrease in the competence of blood to carry oxygen to tissues.
 - Causing tissue hypoxia
- Decrease in the normal concentration of Hgb and/or RBCs
 - RBCs are destroyed prematurely
 - Normal life span ~120 days



Introduction to Anemia

- Individual may or may not become anemic
- Presence and severity of anemia depends on:
 - Severity of hemolysis
 - Ability of BM to compensate for RBC loss



Introduction to Anemia

- Anemia develops when:
 - RBC loss or destruction exceeds the maximum capacity of the BM RBC production
 - The BM RBC production is impaired




Introduction to Anemia

- Diagnosing Anemia
 - Most often is discovered through laboratory tests
 - Patient History is always important
 - Family history
 - Description and duration of symptoms
 - Medications
 - Blood loss
 - Organomegaly



Laboratory Investigation

- CBC
 - RBC count, Hgb, Hct, RBC indices, WBC count, platelet count
 - Differential with RBC morphology
- Additional tests
 - Reticulocyte count
 - Bilirubin
 - Urine and stool can be examined for blood



Functional Classification of Anemia

- Three pathophysiologic mechanisms
 - Proliferation defect
 - decreased production
 - Maturation defect
 - **Survival defect**
 - increased destruction




Functional Classification

- Survival defects
 - Results from loss of circulating RBCs
 - Either by hemorrhage or hemolysis
- BM maturation is increased and orderly



Functional Classification

- Survival defects
 - No anemia
 - If BM can \uparrow production at same rate as cells being lost or hemolyzed
 - Compensated hemolytic anemia
 - Can rapidly develop into anemia
 - RBC destruction \uparrow beyond capacity of BM (hemolytic crisis)
 - BM stops producing RBCs (aplastic crisis)



Survival defects

- Classify hemolytic anemia based on cause of shortened RBC survival
 - Extrinsic anemia
 - Antagonist in RBC's environment causes injury to RBC
 - toxic substances
 - RBC trauma in circulation
 - **Immune-mediated destruction**



Survival defects

- Common observation in circulation
 - Poikilocytes form after the RBC leaves the BM
 - Schistocytes
 - Result of intravascular mechanical trauma to RBC
 - Spherocytes
 - Indicate extravascular RBC membrane damage



Survival defects

- Usually patients demonstrate normocytic normochromic RBCs
 - May see macrocytosis if reticulocyte production is \uparrow enough
 - May see microcytosis if there is a predominance of schistocytes and or spherocytes



Survival defects

- Sites of destruction in hemolytic anemia
 - Intravascular hemolysis occurs within circulation due to:
 - Activation of complement on RBC membrane
 - Physical or mechanical trauma to RBCs
 - Presence of soluble toxic substances



Sites of destruction in hemolytic anemia

- Intravascular hemolysis
 - release of cell contents into the plasma
 - destruction of red blood cells in the circulation



Sites of destruction in hemolytic anemia

- Intravascular hemolysis
 - Laboratory findings:
 - Hemoglobinemia
 - Hemoglobinuria
 - Hemosiderinuria
 - the presence of hemosiderin in the urine
 - » "brown urine"
 - ↓ haptoglobin
 - ↑ LD (early indicator)
 - Schistocytes on RBC morphology





Sites of destruction in hemolytic anemia

- Extravascular hemolysis
 - Occurs within macrophages of the spleen, liver, and BM due to phagocytes in the tissues removing RBCs from circulation
 - Hgb NOT released directly in the plasma
 - NO Hemoglobinemia
 - NO Hemoglobinuria



Sites of destruction in hemolytic anemia

- Extravascular hemolysis
 - Measurements of the products of heme catabolism
 - ↑ exhaled CO
 - ↑ carboxyhemoglobin
 - ↑ serum bilirubin
 - ↓ haptoglobin if severe or chronic
 - ↑ urine or fecal urobilinogen
 - Spherocytes on RBC morphology



Immune hemolytic anemia (IHA)

- Underlying mechanism involves the reaction of an antibody and/or complement with RBC antigens
 - With subsequent cell destruction



Classification of IHAs

- Based on stimulus for antibody production
 - Alloimmune hemolytic anemia
 - Autoimmune hemolytic anemia
 - Drug-induced hemolytic anemia
- Important to determine process because each type requires specific treatment



Immunity

- In a normally functioning immune system, some individuals will produce antibodies to foreign antigens.



Alloimmune hemolytic anemia

- Antibody (Ab) development to RBC antigen (Ag) that the individual lacks.
 - Does not react with individual's own RBCs
- HDFN
 - mother makes Abs against Ags on fetal RBCs
 - » Inherited from the father
- Transfusion reactions
 - where recipient makes Abs to Ags on transfused (donor) cells




IAT

- IAT (indirect Coombs)
 - detects free antibody in the serum/plasma
 - Used for antibody screen and antibody ID panels



IAT

Indirect Coombs test / Indirect antiglobulin test



IAT

Donor	D	C	c	E	e	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Patient Results
I	+	0	0	+	+	+	+	+	+	+	+	0	0	0	0	+	+	+	3+
II	+	0	+	+	0	0	+	0	+	0	+	0	0	+	0	0	+	0	3+
III	0	0	+	0	+	0	+	+	+	+	+	+	+	+	+	+	+	+	0

Donor	D	C	c	E	e	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Patient Results	
1	-	+	0	0	+	+	+	+	+	+	+	0	0	0	0	+	+	+	3+	
2	+	+	0	0	+	+	+	+	+	+	+	0	0	+	0	+	+	+	3+	
3	+	+	0	0	+	+	+	+	+	+	+	0	0	0	+	0	+	+	3+	
4	+	+	0	+	0	0	+	+	+	0	0	0	0	+	+	+	+	+	3+	
5	+	+	0	+	0	0	+	+	+	0	0	+	+	+	+	+	+	+	3+	
6	0	+	+	0	+	0	+	+	+	+	+	+	+	+	+	+	+	+	0	
7	0	0	+	+	0	+	+	+	+	+	+	0	0	+	+	+	+	+	0	
8	0	0	+	+	0	+	+	+	+	+	+	0	0	+	+	+	+	+	0	
9	0	0	+	+	0	+	+	+	+	+	+	0	0	+	+	+	+	+	0	
10	0	0	+	+	0	+	+	+	+	+	+	0	0	+	+	+	+	+	0	
																			Autocontrol	0

D = no agglutination
+ = agglutination (with or without strength of reaction noted in front)

DAT

- DAT – (direct Coombs)
 - Means for determining whether immune-mediated destruction of RBCs is contributing to anemia
 - Positive result confirms presence of RBCs coated *in vivo* with antibodies or complement


DAT

- Does not distinguish between autoantibodies and alloantibodies


DAT

- Antihuman globulin (AHG) directed against
 - IgG alone
 - Complement components only
 - Specifically C3d and/or C3b
 - Polyspecific
 - Combination of anti-IgG and anti-Complement

DAT



DAT



DAT

- Applications of the DAT
 - Investigation of HDFN
 - Hemolytic Transfusion Rxn (HTR)
 - First immunologic evidence of HTR
 - alloantibody coating transfused cells
 - **Investigation of autoantibodies**
 - autoantibodies coating your own cells
 - **Drug Induced Antibodies**



Autoimmune hemolytic anemia (AIHA)




Autoimmunity

- Occasionally, a malfunction in the mechanism regulating immune responses occurs and antibodies against “self” (autoantibodies) are produced.



Autoimmune hemolytic anemia (AIHA)

- Shortened RBC survival
 - <120 days
- Caused by production of autoantibodies against RBC antigens
 - Antibody attaches to RBC creating an altered RBC which is removed from circulation
 - Most hemolysis is extravascular via splenic macrophages.



AIHA

- Further classified
 - Warm-antibody autoimmune HA (WAIHA)
 - Cold-antibody AIHA
 - Cold Agglutinin Disease (CAD) or syndrome (CAS)
 - Paroxysmal Cold Hemoglobinuria (PCH)
 - Mixed-Type AIHA



Warm-antibody autoimmune (WAIHA)



Warm-antibody autoimmune (WAIHA)

- Most common form of AIHA
 - ~70% of the cases
 - Optimal reactivity at 37° C



Warm-antibody autoimmune (WAIHA)

- Can occur at any age
 - Incidence increases after age 40
 - Childhood incidence peaks in the first 4 years of life



Warm-antibody autoimmune (WAIHA)

- Culprit
 - Usually IgG antibody
 - rarely IgM, IgA
 - Most Abs react with "Rh protein" complex
 - Do not react with Rh null
 - Occasionally have single specificity within Rh system
 - » anti-e
 - In conjunction with or precipitated by formation of alloantibody



Warm AIHA

- Idiopathic WAIHA
 - Accounts for about 60% of cases of warm AIHA
 - Acute idiopathic WAIHA
 - Severe anemia
 - Developing over two to three days
 - » Hemolysis is self limited
 - Duration
 - » Several weeks → several years duration



Warm AIHA

- Chronic idiopathic WAIHS
 - Hemolysis is unabating




Warm AIHA

- Secondary WAIHA associated with:
 - Lymphoproliferative disease
 - CLL , Hodgkin's
 - Many children with idiopathic WAIHA eventually develop a lymphoproliferative disease.




Warm AIHA

- Secondary WAIHA associated with:
 - Neoplastic Diseases
 - referring to both malignant and benign growths
- Few Examples
 - Leukemia
 - Metastatic Carcinoma
 - Multiple Myeloma
 - Osteosarcoma




Warm AIHA

- Secondary WAIHA associated with:
 - Other Autoimmune disorders
 - SLE, RA, Crohn's disease, etc.




WAIHA Blood Bank Testing

- Warm Autoantibodies
 - ABO type generally not affected
 - Testing is at RT (IS)
- Rh grouping
 - only affected if weak D testing is done
 - Weak D is not recommended




WAIHA Blood Bank Testing

- Warm Autoantibodies
 - DAT is positive (usually)
 - cells coated with IgG
- IAT may be positive or negative




Common Blood Bank findings in AIHA

- Positive IAT if the DAT is so strongly positive that Ab is “spilling” into the serum
 - Everything may be positive
 - Including the autocontrol
- Why it matters:
 - Risk the inability to detect underlying alloantibody; must rule-out or ID alloantibody by warm auto/allo adsorption techniques



	IS	37°c	AHG (IAT)	ENZYME (FICIN IAT)
Screen Cell I	0	0	2+	4+
Screen Cell II	0	0	2+	4+
Screen Cell III	0	0	2+	4+
Autocontrol	0	0	2+	4+

- Autocontrol positive...only an autoantibody?
 - No way to know without further testing
 - The autoantibody must be removed so we can see if there is an alloantibody
 - Adsorption and Elution techniques



- Warm antibodies in the serum may show a relative anti-e specificity
 - Reacting weakest with e negative RBCs
 - May react with all RBCs of normal Rh type

CELL	D	C	E	c	e	SERUM/PLASMA				
						RT	37°C	AHG(IAT)		
R1R1-39	1	+	+	0	0	+	0	0	3+	
R1R1-44	2	+	+	0	0	+	0	0	3+	
R2R2-23	3	+	0	+	+	0	0	0	±	
rr-33	4	0	0	0	+	+	0	0	3+	
rr-26	5	0	0	0	+	+	0	0	3+	
r'r-4	6	0	+	0	+	+	0	0	3+	
r'r-17	7	0	0	+	+	+	0	0	3+	
R0r-13	8	+	0	0	+	+	0	0	3+	
R1r-14	9	+	+	0	+	+	0	0	3+	
R1R2-8	10	+	+	+	+	+	0	0	3+	
Autocontrol Untreated RBCs	11								3+	
Autocontrol EGA treated	12	EGA= EDTA glycine acid treated								3+

Common Blood Bank findings in AIHA

- Reasons for negative IAT
 - All autoantibody bound to red cells
 - no “spill over” into serum
 - Why it matters:
 - Generally better prognosis because total antibody production is lower

Common Blood Bank findings in AIHA

- Positive DAT because the antibody is already present *in vivo*
 - Patients who have not been recently transfused

- DAT and autocontrol are not equivalent

Poly AHG	Anti-IgG	Anti-C3	Control
3+	2+	0	0

- Negative DAT
 - May not detect very small amounts of IgG antibodies on RBC
 - Result negative if complete destruction of antibody or complement-coated RBCs occurred

Laboratory Goals

- Establish that antibody on patients cell is autoantibody
- Determine if there is a specificity to the autoantibody
- Detect and identify any clinically significant "underlying" alloantibodies



Identification Techniques


- Adsorption
 - Technique to remove an Ab from a patient's plasma
 - Removing antibody from plasma by binding to the corresponding antigen on RBCs
 - using optimal incubation conditions



- Adsorption
 - useful in differentiating autoantibody from alloantibody
 - Remove autoantibody reactivity so that alloantibody detection tests and diagnostic tests for differentiating the immune hemolytic anemias can be performed.



Auto/allo antibodies




- Adsorption
 - Autologous RBCs could be coated with antibodies (DAT+)
 - Adsorption is facilitated by dissociating autoantibody from the RBC membrane
 - Uncovering Ag sites that can bind free autoantibody to remove it from the serum
 - » Can be performed by heat elution RBCs placed at 56° C for 3-5 minutes
 - » Subsequent treatment of RBCs with enzymes enhances that adsorption process




- Treatment of RBCs with enzymes enhances that adsorption process by removing membrane structures that otherwise hinder the association between antigen and antibody
 - ZZAP reagent
 - Mixture of proteolytic enzyme and sulfhydryl reagent dithiothreitol (DTT)



- ZZAP removes immunoglobulins and complement from the RBCs and enhances the adsorption process




- Autologous adsorption
 - Use when patient has not been transfused within last 3 months
 - Use patient's own cells for adsorption




Autoadsorption

- Patient has not been transfused within last 3 months
 - All RBCs are self
 - Only auto antibodies will be removed with the use of autologous cell for the autoadsorption because alloantibodies are made to RBC Ags the patient lacks.




Autoadsorption




Allogeneic/Homologous adsorption

- Allogeneic/Homologous adsorption
 - Use when patient has been recently transfused
 - Transfused (donor) cells present in circulation are likely to adsorb the alloantibodies being sought.
 - Use cells with complementary phenotypes to remove the autoantibody and leave the alloantibody in the test system
 - » Alloantibodies that remain can be confirmed by testing against a panel of reagent RBCs



Allogeneic/Homologous adsorption


- Treating the adsorbing cells with enzyme or ZZAP typically enhances the adsorption process.
 - Treated cells may lack the antigens
 - Those antigens may be destroyed by DDT and/or enzymes.



Alloadsorption

- Use when the patient has been recently transfused
 - RBCs are both self and donor
 - should not be used for adsorption
 - Use cells with complementary phenotypes
 - To the patient's phenotype
 - Or common phenotypes of the population
 - R1R1, R2R2, rr, etc.

Patient serum added, 37C incubation, spin and remove serum.



Test the adsorbed serum to look for remaining antibody reactivity and determine specificity.

This testing is done by advanced transfusion service or sent to an IRL.


R1, K-, s-, Jk(a+b-)
 R2, K-, s+, Jk(a-b+)
 rr, K+, s-, Jk(a+b-)

- Elution
 - Whenever the DAT is positive we can perform an elution
 - Or after adsorption techniques
 - Recovery of adsorbed antibody by treating antibody-coated RBC's to break the bonds between antigen and antibody

- Elution
 - Autoantibody in eluate will usually react with all cells tested
 - Pan-agglutination
 - Concentration of antibody removed from the patient's cells may be greater than the quantity of antibody in the serum
 - Differentiate autoantibody from alloantibody to high prevalence antigen
 - » by testing eluate against patient's cells

- Prior to an elution, cells must be thoroughly washed to remove serum containing excess antibody
 - Supernatant fluid from the final wash of RBCs to be eluted should be tested in parallel with the eluate
 - Should be non-reactive
 - Ensures antibody detected in the eluate is on RBC bound antibody and not free antibody from the plasma

- Elution
 - Ag-Ab bond must be broken
 - Releasing Ab from the Ag on the RBC into fluid medium
 - Usually saline



CELL	D	C	E	c	e	ELUATE		
						Eluate IAT	Last Wash IAT	
R1R1-39	1	+	+	0	0	+	4+	0✓
R1R1-44	2	+	+	0	0	+	4+	0✓
R2R2-23	3	+	0	+	+	0	4+	0✓
rr-33	4	0	0	0	+	+	4+	0✓
rr-26	5	0	0	0	+	+	4+	0✓
r'r-4	6	0	+	0	+	+	4+	0✓
r'r-17	7	0	0	+	+	+	4+	0✓
R0r-13	8	+	0	0	+	+	4	0✓
R1r-14	9	+	+	0	+	+	4+	0✓
R1R2-8	10	+	+	+	+	+	4+	0✓
Autocontrol Untreated RBCs	11						4+	3+
Autocontrol EGA treated	12						4+	0✓

cell	Warm Autoantibody with underlying alloanti-E														Patient IAT		Auto-ADS						
	D	C	E	c	e	P1	M	N	S	s	Lea	Leb	K	k	Fya	Fyb	Jka	Jkb	37	IAT	Ficin IAT		
1	++	0	0	+	+	+	+	0	+	0	+	0	+	+	0	+	+	0	+	+	0	3+	0✓
2	++	0	0	+	+	+	+	0	+	0	+	0	+	+	+	+	+	0	0	0	0	3+	0✓
3	+	0	+	+	0	+	0	+	0	+	0	+	0	+	+	+	0	+	0	2+	4+	3+	
4	0	0	0	+	+	0	+	0	+	+	+	0	0	+	0	+	0	+	0	0	0	3+	0✓
5	0	0	0	+	+	+	+	+	+	0	+	+	+	+	+	+	+	0	0	0	0	3+	0✓
6	0	+	0	+	+	0	+	0	+	+	0	0	+	+	0	+	0	+	0	0	0	3+	0✓
7	0	0	+	+	+	+	+	+	0	0	+	0	+	+	0	0	+	1+	4+	2+			
8	+	0	0	+	+	+	+	+	0	0	0	0	0	+	0	0	+	+	0	0	0	3+	0✓
9	++	0	+	+	+	0	+	+	+	+	0	0	+	+	0	+	+	0	0	0	0	3+	0✓
10	++	+	+	0	+	0	+	+	0	+	0	+	0	+	+	+	+	1+	4+	2+			
11 AC																			0	3+	3+		
12 AC treat																			0	3+	0✓		

Treatment of WAIHA

- Steroids
 - suppress autoantibody production
- Splenectomy
 - Decreases antibody production
 - Removes potent site of RBC damage and destruction

Treatment of WAIHA

- Immunosuppressive drugs
- To transfusion or not to transfuse?

Treat with Transfusion?

- One Perspective:
 - Transfusion should be avoided unless life threatening anemia.
 - Suitable units can be challenging and time-consuming to find
 - due to the need to rule out alloantibody/ies

Treat with Transfusion?

- Another Perspective:
 - Do not withhold blood transfusion to patient with a justified need
 - RBC destruction of transfused cells should be no different than currently circulating cells



Cold-antibody AIHA




Benign Cold Reactive Autoantibody Characteristics

- Demonstrable in serum of many normal, healthy adults when testing at 4° C
 - Usually does not present a serological problem because routine tests are not done at 4° C
- Usually IgM, anti-I
 - Low titer (<64 at 4° C)



- May interfere with routine tests performed at RT
 - RBCs may be heavily coated
 - May directly agglutinate
- Causing false positives

	Anti-A	Anti-B	Anti-D	D control
Serum suspended RBCs	1+	1+	1+	1+
Warm-washed, saline suspended RBCs	0	0	1+	N/A



	I.S	37C	IAT	Enzyme IAT
Screen cell I	1+	0	0√	3+
Screen cell II	1+	0	0√	3+
Screen cell III	1+	0	0√	3+
Auto Control	2+	0	0√	3+

Cord cell I	0
Cord cell II	0

Cord cells are I negative, i positive, confirming Anti-I in panel above



Pathologic Cold Autoagglutinins

- Cold Agglutinin Syndrome (CAS) or Disease (CAD)
- Paroxysmal Cold Hemoglobinuria (PCH)



Cold AIHA

- Cold Agglutinin Disease (CAD) or syndrome (CAS)
- Accounts for ~ 16–30% of cases of AIHA
 - Less common in children than adults



Cold Agglutinin Disease (CAD)

- Optimal reactivity < 37° C
- Culprit
 - Usually IgM (rarely IgA or IgG) antibody
 - Capable of binding complement
 - Complement mediated lysis
 - Most with anti-I specificity
 - High Titer ≥ 1000



Cold Agglutinin Disease (CAD)

- Idiopathic chronic CAS
 - Usually chronic, occurring after age 50
 - Anti-I
 - Monoclonal IgM/kappa light (κ) chain



Cold Agglutinin Disease (CAD)

- Secondary CAS
 - Acute, self-limiting usually associated with infectious disease
 - Polyclonal IgM auto-Abs with specificity for Ii antigens
 - anti-I—M. pneumoniae
 - anti-i—Infectious mononucleosis
 - anti-Pr—varicella, rubella



Cold Agglutinin Disease (CAD)

- Secondary CAS
 - Chronic form is associated with Lymphoproliferative disorders
 - Lymphoma
 - Waldenstrom's Macroglobulinemia
 - Monoclonal IgM/k light chain



Cold Agglutinin Disease (CAD)

- Most are not clinically significant
 - But can be with resultant IHA.
 - Severity of disease
 - Related to thermal range of the Ab



Cold Agglutinin Disease (CAD)

- Cold-reacting antibodies with a wide range of thermal reactivity (up to 32° C) can cause problems when the peripheral circulation cools to this temperature.
 - Chronic hemolytic or Episodic hemolysis associated with chilling
 - Areas of the body that cool to the Ab thermal range
 - Sludging of blood flow within capillaries



Cold Agglutinin Disease (CAD)

- Extremities affected
 - Nose, fingers, toes and ears
 - Vascular changes
 - Acrocyanosis
 - » Hands and/or feet turn blue
 - Raynaud's phenomenon
 - » Pain with color change patterns in skin
 - » White: spasm of the vessels
 - » Blue: cyanosis
 - » Red: indicates return of blood flow to the area



Cold Agglutinin Disease (CAD)

- Hemoglobinuria
 - Accompanies acute hemolytic attacks
- Automated blood counts
 - Must warm blood
 - Falsely ↓ RBC
 - Falsely ↑ MCV



Pathologic CAS serology

- May interfere with ABO/Rh, antibody detection and ID, and cross-match
- Enhanced reactivity with enzymes
 - Most have anti-I specificity



Pathologic CAS serology

- Positive DAT
 - Due to complement
 - Polyspecific AHG
 - anti-C3b/C3d
- IgM with titer often >1000 at 4° C
 - May react over a wide thermal range up to 37° C



Cold Agglutinin Disease (CAD)

- Most patients require no treatment and are instructed to avoid the cold, keep warm, or move to a milder climate



Paroxysmal Cold Hemoglobinuria (PCH)

Paroxysmal Cold Hemoglobinuria (PCH)

- Least common type of AIHA
 - Can occur at any age
 - Accounts for 30–40% of all AIHA in children
 - Most frequent < age 5
 - Can be associated with viral and bacterial infections
 - » measles, mumps, influenza, adenovirus, chickenpox, CMV, and Epstein-Barr virus (infectious mononucleosis)
 - » syphilis, *Haemophilus influenzae* and *Mycoplasma pneumoniae*

Paroxysmal Cold Hemoglobinuria (PCH)

- Usually transient disorder
 - Resolves spontaneously with resolution of infectious
 - Transfusions may be needed in severe cases
- Chronic was linked to syphilis

Paroxysmal Cold Hemoglobinuria (PCH)

- Culprit
 - Bi-phasic complement fixing IgG antibody
 - Donath Landsteiner antibody
 - usually anti-P specificity
 - At lower temperatures (< 20° C), IgG binds to RBCs and activates complement
 - At 37° C IgG elutes off leaving only complement coating cells
 - » Cells undergo complement-mediated intravascular lysis

Paroxysmal Cold Hemoglobinuria (PCH)



- Clinical findings
 - Hemoglobinuria, Hemoglobinemia, Jaundice
 - ↓ serum complement
 - Erythrophagocytosis
 - Anemia depends on frequency and severity of attack
 - Hgb drops sharply
 - Can ↓ as low as 5 g/dL

PCH serology

- DAT is positive
 - due to complement
 - Polyspecific
 - Anti-C3b/C3d

PCH serology

- Definitive test is Donath-Landsteiner that detects biphasic properties of the antibody
 - Sample must be collected and maintained during clotting at 37° C
 - Immediately separate serum from red cells to isolate the antibody for processing



Mixed Autoimmune Hemolytic Anemia (MAIHA)

MAIHA

- Simultaneously have two types of autoantibodies causing both WAIHA and CAIHA
 - Rare

MAIHA

- Presence of warm-reacting IgG autoAb and cold-reacting IgM autoAb
 - Both have high titer and ↑ thermal amplitude
 - Reacts at >30° C
 - Mixture of both intravascular (IgM) and extravascular (IgG) hemolysis

MAIHA

- Idiopathic
 - 50% of cases are idiopathic
- Remainder are secondary to:
 - lymphoproliferative disorder
 - SLE
 - HIV



MAIHA

- Treatment may lead to one form of AIHA predominating over the other
 - Most patients respond to corticosteroids without transfusions.



Let's "do" drugs



Drug-Induced HA

- Acquired cause of hemolytic anemia
 - Not all individuals taking the same drug develop HA
- > 125 drugs identified
 - The drug itself does not cause RBC injury
 - Immune response to drug-induced alteration of the RBC



Drug-Induced HA

- Drugs implicated
 - Nonsteroidal anti-inflammatory drugs
 - Diuretics
 - Antineoplastic drugs
 - Antimicrobials
 - 3rd generation cephalosporins
 - Cefotetan and ceftriaxone with the majority of cases
 - Most fatalities
 - » >50% with ceftriaxone
 - Piperacillin



Drug-Induced HA

- Suspected when there is no other explanation for the serologic or hematologic findings and if the patient has a history of taking the drug
- Resolution is withdrawal of the offending drug



Drug-Induced HA

- Mechanisms are controversial and not well defined



Drug-Induced HA

- New "unifying" hypothesis
 - Drug binds to RBC membrane
 - Once bound, antibodies can be produced to react with:
 - Ag specific to the drug
 - Ag that resemble a combination of drug and RBC proteins
 - Ag primarily on the RBC membrane
 - Explains how patients develop more than one type of drug-induced Ab



Drug-Induced HA




Figure 20-1. Proposed unifying theory of drug-induced antibody reactions (based on a cartoon by Habibi as cited by Garratty¹¹). The thicker darker lines represent antigen-binding sites on the F(ab) region of the drug-induced antibody. Drugs (haptens) bind loosely, or firmly, to cell membranes and antibodies may be made to: a) the drug (producing in-vitro reactions typical of a drug adsorption (penicillin-type) reaction); b) membrane components, or mainly membrane components (producing in-vitro reactions typical of autoantibody); or c) part-drug, part-membrane components (producing an in-vitro reaction typical of the so-called immune complex mechanism).^{10,12}



Drug-Induced HA

- Two major types of drug induced antibodies have been described
 - Drug dependent antibodies
 - requires presence of drug during testing
 - Drug independent antibodies
 - reacts without presence of drug



Drug-Induced HA

- Drug dependent antibodies
 - Drug adsorption mechanism
 - Drug binds covalently to the RBC
 - Ab produces against the drug Ag
 - Ab attaches to the drug-coated RBC which is removed from circulation
 - » Destruction is usually extravascular, complement is seldom involved
 - Drugs implicated
 - Piperacillin and cefotetan
 - IV penicillin
 - » high dose



Drug-Induced HA

- Drug dependent antibodies
 - Drug binds weakly to the RBCs
 - Ab is formed to Ag that consist of both drug and erythrocyte membrane
 - Ab fixes complement
 - » Intravascular hemolysis
 - Drug implicated
 - Ceftriaxone



Drug-Induced HA

- Drug independent antibodies
 - Drug binds to the RBC
 - Ab is produced against primarily RBC Ag
 - Can be due to:
 - Alteration of the RBC membrane by the drug
 - Molecular mimicry



Drug-Induced HA

- Drug independent antibodies
 - Serological reactions indistinguishable from warm AIHA
 - Thought that the RBC membrane is altered by the drug, thus the immune system sees the alteration as a foreign antigen.
 - Historically linked to anti-hypersensitive drug Methyl dopa (Aldomet)
 - Today Fludarabine (CLL cancer drug)



Drug-Induced HA

- Drug independent antibodies
 - NIPA, non-immunologic protein adsorption
 - Modification of the RBC membrane that results in nonimmunologically adsorbed IgA, IgM, IgG, C3
 - No drug antibody is involved
 - β -lactamase inhibitors
 - » Clavulanic acid
 - » Sulbactam
 - » Tazobactam
 - Cefotetan
 - platinum-based chemotherapeutic drugs



Drug Autoantibody serologic and clinical observations

- Positive DAT –
 - Usually due to IgG
 - Polyspecific
 - Anti-IgG
- Eluate reacts with all normal cells tested, occasionally showing Rh specificity



Drug Autoantibody serologic and clinical observations

- Autoantibody may “spill over” into serum, causing routine test to be positive
- Alloantibodies are difficult to identify or rule out due to pan-agglutination



Not routine testing

- Most blood banks do not do extensive investigation if patient has positive DAT suspected to be drug related unless there is a hematologic complication
 - Drugs may be continued unless a hemolytic anemia results
 - If the drug is implicated, STOP giving it immediately



ParSCORE™ STUDENT ENROLLMENT SHEET										ParSCORE™ SCORE SHEET									
1-6					11-16					21-26					31-36				
ID NUMBER					PHONE NUMBER					1					2				
Licensure Number					leave blank					3					4				
CITY					STATE					5					6				
ALE					PIN					7					8				
MON					SHR					9					10				
Example: Monroe					Fill in Date					11					12				
										13					14				
										15					16				
										17					18				
										19					20				
										21					22				
										23					24				
										25					26				
										27					28				
										29					30				
										31					32				
										33					34				
										35					36				
										37					38				
										39					40				
										41					42				
										43					44				
										45					46				
										47					48				
										49					50				
										51					52				
										53					54				
										55					56				
										57					58				
										59					60				
										61					62				
										63					64				
										65					66				
										67					68				
										69					70				
										71					72				
										73					74				
										75					76				
										77					78				
										79					80				
										81					82				
										83					84				
										85					86				
										87					88				
										89					90				
										91					92				
										93					94				
										95					96				
										97					98				
										99					100				

1-6 11-16 21-26 31-36

1st speaker 2nd speaker 3rd speaker 4th speaker

Evaluation of Speaker

1. Overall quality of presentation.
2. Knowledge of subject matter.
3. Organization of Presentation

Evaluation of Course

4. Rate the session.
5. Achieved stated objectives.
6. Would recommend this course to others?
A. Yes B. No

Comments on reverse

Speaker - strengths/weaknesses
Course - strengths/weaknesses
General - comments about the program or suggestions for future programs

A - excellent
 B - Good
 C - Average
 D - Fair
 E - Poor