SERVICE MANUAL

SERIES 320 XT SERIES 320 XB SERIES 320 IP65

Order no.

| | | | |
|----------|----------|---------------------|------------------------|
| 350-8100 | Index c4 | Service manual cpl. | Service-Anleitung kpl. |
| 350-7900 | Index c4 | Content | Inhalt |

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Copyright

These service instructions are protected by copyright. All rights reserved. No part of the service instructions may be reproduced, processed, duplicated or published in any form by photocopying, microfilming, reprinting or other process, in particular electronic means, without the written agreement of Precisa Gravimetrics AG.

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Section C: Service twin calibration - - - - - C1

| | |
|---|------------|
| 1. Service tools and equipment..... | C2 |
| 1.1 Tool set (350-8537)..... | C2 |
| 1.2 Equipment..... | C3 |
| 2. Opening up a balance | C4 |
| 2.1 General procedure for opening up a Series 320 XT/XB balance..... | C4 |
| 2.2 General procedure for opening up a Series 320 IP65 balance | C4 |
| 3. Introduction / Preparation of the tools | C5 |
| 3.1 Introduction | C5 |
| 3.2 Preparation of the tools | C5 |
| 4. Dismantling the weighing cell | C6 |
| 4.1 Removing the top of balance | C6 |
| 4.2 Removing the weighing cell from the bottom of balance..... | C6 |
| 4.3 Removing the flexure strap | C6 |
| 4.4 Removing the two flexure holders..... | C6 |
| 4.5 Removing the coil..... | C7 |
| 4.6 Removing the support piece | C7 |
| 4.7 Removing the balance arm | C7 |
| 5. Installing the new flexure sheets | C9 |
| 5.1 Removing the old flexure sheets..... | C9 |
| 5.2 Installing the new flexure sheets | C9 |
| 6. Cleaning | C10 |
| 6.1 Cleaning the pot..... | C10 |
| 6.2 Cleaning the coil | C10 |
| 6.3 Final check of cleaning..... | C10 |
| 7. Assembling the weighing cell | C11 |
| 7.1 Assembling the balance arm..... | C11 |
| 7.2 Assembling the coil, the magnet cover and the transport safety device | C12 |
| 7.3 Assembling the sensor mechanism and the stopper | C12 |
| 7.4 Assembling new flexure sheets..... | C12 |
| 7.5 Assembling the support piece and inner cone | C13 |
| 7.6 Assembling the flexure holders..... | C14 |
| 7.7 Assembling a new flexure strap | C14 |
| 7.8 Dismantling the assembly jigs..... | C15 |
| 7.9 Installing the weighing cell in the bottom of balance | C15 |

Section D: Service single calibration - - - - - D1

| | |
|---|------------|
| 1. Service tools and equipment..... | D2 |
| 1.1 Tool set (350-8537)..... | D2 |
| 1.2 Equipment..... | D3 |
| 2. Opening up a balance | D4 |
| 2.1 General procedure for opening up a Series 320 XT/XB balance..... | D4 |
| 2.2 General procedure for opening up a Series 320 IP65 balance..... | D4 |
| 3. Introduction / Preparation of the tools | D5 |
| 3.1 Introduction | D5 |
| 3.2 Preparation of the tools..... | D5 |
| 4. Dismantling the weighing cell..... | D6 |
| 4.1 Removing the top of balance | D6 |
| 4.2 Removing the calibration unit (XT only)..... | D6 |
| 4.3 Removing the weighing cell from the bottom of balance..... | D6 |
| 4.4 Removing the flexure strap | D6 |
| 4.5 Removing the two flexure holders..... | D7 |
| 4.6 Removing the coil..... | D7 |
| 4.7 Removing the support piece | D7 |
| 4.8 Removing the balance arm | D7 |
| 5. Installing the new flexure sheets | D9 |
| 5.1 Removing the old flexure sheets..... | D9 |
| 5.2 Installing the new flexure sheets..... | D9 |
| 6. Cleaning | D10 |
| 6.1 Cleaning the pot..... | D10 |
| 6.2 Cleaning the coil | D10 |
| 6.3 Final check of cleaning..... | D10 |
| 7. Assembling the weighing cell | D11 |
| 7.1 Assembling the balance arm and the support piece | D11 |
| 7.2 Assembling the coil, the magnet cover and the transport safety device | D12 |
| 7.3 Assembling the sensor mechanism, the inner cone and the stopper..... | D12 |
| 7.4 Assembling the flexure holders..... | D13 |
| 7.5 Assembling new flexure sheets..... | D14 |
| 7.6 Assembling a new flexure strap..... | D15 |
| 7.7 Dismantling the assembly jigs..... | D16 |
| 7.8 Installing the weighing cell in the bottom of balance..... | D16 |
| 7.9 Installing the calibration unit in the bottom of balance (XT only)..... | D16 |

Section E: Adjustments - - - - - E1

| | |
|---|------------|
| 1. Adjustments | E2 |
| 1.1 Activation of the service mode | E2 |
| 1.2 Adjusting of the symmetry | E3 |
| 1.3 Checking the pre-load | E3 |
| 1.4 Adjusting of the corner load | E4 |
| 1.5 Checking the lowering process of the calibration weight (XT only) | E7 |
| 1.6 Final assembly of the balance | E7 |
| 1.7 Checking the corner load | E8 |
| 1.8 Linearisation of a 320XT Series balance | E8 |
| 1.9 Linearisation of a 320XB Series balance | E9 |
| 1.10 Final check of the adjustment | E10 |
| 1.11 S-correction of a balance | E10 |
| 1.12 Adjustment of the internal calibration weight (XT only) | E11 |
| 2. Final control | E13 |
| 3. Error messages | E14 |
| 3.1 Operating error | E14 |
| 3.2 Fatale error | E14 |
| 3.3 Hardware error | E14 |
| 4. Discussion with the manufacturer | E15 |

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In order to facilitate the searching of spare parts, this chapter is reorganized in groups.

The following chart shows, to which group the balance belongs, regarding the series number, the order number and the model type.

Example:

Series 320XT

Order No. 320-9301

Model 620M

Group M1

On all subordinated representations, the desired spare part is to find in group **M1**

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
2.1 Measuring resistor / Messwiderstand A2

| | | | |
|---|--------------|-----------------------|-------------------|
| 6 | 320-7032-004 | Measuring resistor A2 | Messwiderstand A2 |
|---|--------------|-----------------------|-------------------|

2.2 Measuring resistor / Messwiderstand M1

| | | | |
|---|--------------|-----------------------|-------------------|
| 6 | 320-7032-003 | Measuring resistor M1 | Messwiderstand M1 |
|---|--------------|-----------------------|-------------------|

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








1. Main board / Hauptprint 320-7215-020 / 021

1.1 Exchange of the main board 320-7215-020 / 021

The new basic version of the Main Board 320-7215-020 or 320-7215-021 is applicable for all balance with software N5x-xxxx-xxx.

For balances with a resolution smaller than **620'000 Digits**, the mainboard with suffix **-020** is suitable. For balances with higher resolutions the mainboard with suffix **-021** has to be used.


The drawing below shows which hardware components are to be converted by an exchange of a main board.

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







2. Main board / Hauptprint 320-7200-020 / 021

2.1 Exchange of the main board 320-7200-020 / 021

The old basic version of the main board 320-7200-020 or 320-7200-021 is no more available.

The old basic version of the main board 320-7200-020 / 021 is full exchangeable with the new version of the main board 320-7215-020 / 021, see "Main board / Hauptprint 320-7215-020 / 021" on page B2.

In case of special software, please contact our PRECISA customer service dept.

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
Section C: Service twin calibration

1. Service tools and equipment

1.1 Tool set (350-8537)

| Position | Part-No. | Description | Quantity |
|----------|--------------|--------------------------------------|----------|
| 1 | W 41-1688 | Mounting plate | 1 |
| 2 | W 41-1688-5 | Support post, long | 2 |
| 3 | W 41-1688-6 | Support post, short | 2 |
| 4 | 350-6203 | Line-head screwdriver M4 | 1 |
| 5 | 240-7124 | Screwdriver, small, no.1 | 1 |
| 6 | PW 29.43.300 | Pin for flexure holder | 2 |
| 7 | W 41-1700 | Assembly jig with slot | 1 |
| 8 | W 41-1701 | Assembly jig with hole | 1 |
| 9 | W 41-1677 | Assembly pin for assembly jig | 4 |
| 10 | W 41-1691 | Centring spigot for flexures | 2 |
| 11 | W 41-1774 | Centring spigot for the coil EP / ES | 1 |
| 12 | W 41-1675 | Centring spigot for the coil XT / XB | 1 |
| 13 | W 41-1674 | Flexure sheet anchor | 1 |
| 14 | PN 1100-172 | Line-head screw M4x12 | 2 |
| 15 | W 41-1699 | Alignment jig | 1 |
| 16 | W 41-1578-8 | Connecting cable | 1 |
| 17 | W41-1681 | Centring bush | 1 |

1.2 Equipment

1.2.1 Implements

- Soldering iron
- Tweezers
- Lens
- Voltmeter

1.2.2 Standard tools

- Box spanner Nos.4, 5.5 and 7
- Flat-bladed screwdriver Nos.1, 2 and 4
- Phillips screwdriver No.2
- Hexagon socket wrench WS 2.5mm


2. Opening up a balance

2.1 General procedure for opening up a Series 320 XT/XB balance

1. Disconnect the balance from the mains (Pull out mains plug)
2. Remove all windshields, protective rings, weighing pans and pan holders
3. Unscrew 3 cylindrical screws on the top of balance using a No. 4 screwdriver
4. Slightly raise the top of balance
5. Release the cover-foil ribbon cable from the main board

2.2 General procedure for opening up a Series 320 IP65 balance

1. Disconnect the balance from the mains (Pull out mains plug)
2. Remove all windshields, protective rings, balance pans and pan holders
3. Unscrew 3 screws IP56 with a No. 1 screwdriver and remove the clamping ring
4. Unscrew 3 cylindrical screws on the top of balance and 2 cylindrical screws on the bottom of balance using a No. 4 screwdriver
5. Slightly raise the top of balance (Do not damage the sealing)
6. Release the cover-foil ribbon cable from the main board

3. Introduction / Preparation of the tools

3.1 Introduction

In the event of a broken flexure, it is recommended that **all** flexures be replaced. The following assembly instructions provide information on how the flexures can be replaced most efficiently. It is important to employ a structured procedure for dismantling since it is essential to replace the **same parts** (particularly screws) **in the same place** when re-assembling.

The various items of the assembly instructions mean:

(26/B3): component items


[5]: tooling items see "Service tools and equipment" on page C2

Example: 3 nuts (26/B6) means 3 nuts as Item 26 in section B on page 6

3.2 Preparation of the tools

Before dismantling the balance it is necessary to carry out some preparation of the tools. 2 long [2] and 2 short [3] support posts must be screwed into the correct places on the assembly plate

The soldering iron should be plugged in or preheated


4. Dismantling the weighing cell

4.1 Removing the top of balance

- see "Opening up a balance" on page C4

Fig. 1

4.2 Removing the weighing cell from the bottom of balance

1. On the main board, unplug the sensor mechanism plug (13/A11) using a small screwdriver
2. On the calibration, unscrew 4 screws (12/A37) from the counterholder (4/A37) and remove counterholder and calibration weight (XT only)
3. Turn over the bottom of balance and unscrew the 3 nuts (26/A11) on the underside of the housing
4. Hold the weighing cell **only** on the chassis and lift it out of the bottom of balance
5. The chassis can now be placed on the assembly plate for easier handling
6. To protect the weigh cell, the two assembly jigs [7] + [8] are introduced in accordance with Fig. 1 and secured with the 4 assembly pins [9]

4.3 Removing the flexure strap

1. Unscrew 2 screws (25/A11), remove flexure strap (3/A11) incl. bases (9/A11) and spring washers (28/A11) → Do not separate screw components, assemble in same relationship
2. Screw the flexure strap anchor [12] on tightly with 2 screws [13] in place of the flexure strap

4.4 Removing the two flexure holders

1. Using the two pins for flexure holder [6] fix the upper flexure holder (2/A11) into the chassis

2. Unscrew the 4 screws (24/A11) of the upper flexure holder and remove together with the 4 spring washers (28/A11) → 2 on the chassis, 2 on the support piece
3. Remove the two pins with the upper flexure holder incl. 4 spring washers and the flexure sheets, then take out the pins from the flexure holder
4. Fix the lower flexure holder (2/A11) into the chassis with the two pins for flexure holder [6]
5. Unscrew 4 screws (24/A11) from the lower flexure holder and remove together with the 4 spring washers → 2 on the chassis, 2 on the support piece
6. Take out both pins with the lower flexure holder and remove the lower flexure holder together with the 4 spring washers and the flexure sheets

4.5 Removing the coil

1. Carefully unsolder 2 connecting wires (1/A11) from the connecting board (17/A11) using the tweezers
2. Unscrew 3 screws (23/A11) from the magnet cover (8/A11), remove the cover together with transport safety device (4/A11) and spring washer (28/A11)
3. Unscrew 2 brass screws (21/A11) from the coil (12/A11) on the balance arm (19/A11) using a No.2 screwdriver, remove the coil together with the washer (27/A11)
4. Carefully take out the coil from the pot (34/A11)

4.6 Removing the support piece

1. Unscrew 3 screws (23/A11) on the inner cone (13/A11) and remove inner cone
2. Unscrew upper screw (balance arm side) of the flexure strap anchor [12] and pull off the support piece together with the flexure strap anchor over the assembly pins [9]

4.7 Removing the balance arm

1. Unscrew brass screw (30/A11) from the stopper (29/A11), turn the stopper inwards ca. 90° and gently re-tighten the brass screw
2. Unscrew the 4 screws (25/A11) from the flexure sheets (6/A11) and remove both flexure sheets together with the bases (9/A11) and spring washers (28/A11)

Note: Since the lug of the sensor mechanism is firmly attached to the balance arm, the sensor mechanism must be unscrewed to dismantle the balance arm and gently tilted back:

3. Unscrew 2 screws (20/A11) from the sensor mechanism plate (5/A11) and remove the plate
4. Unscrew screw (20/A11) from the connection board (17/A11)
5. Unscrew 3 screws (22/A11) of the sensor mechanism and gently tilt the sensor mechanism to the rear
6. Remove the two inner assembly pins [9] → unscrew from the balance arm
7. Remove the two outer assembly pins [9] → unscrew from the chassis

8. Remove the two assembly jigs [7] + [8]
9. Carefully tilt the balance arm (19/A11) over the pot (34/A11) and the sensor mechanism and take out
10. Further dismantling of the weighing cell is unnecessary

5. Installing the new flexure sheets


5.1 Removing the old flexure sheets

- Take out all screws from the flexure holder and remove the old flexure sheets together with spring washers

5.2 Installing the new flexure sheets

1. Place a flexure holder on the assembly plate [1] as in Fig.2
2. Select 4 new flexure sheets (7/A11) → Place the new flexure sheets carefully on the assembly plate [1] as in Fig. 2 → the slot of the sheet is always on the flexure holder
3. Place 4 spring washers (28/A11) on the flexure sheet and screw these tightly to the flexure holder with 4 screws (22/A11)
4. Take the complete flexure holder carefully off the assembly plate → Do not touch the flexures!
5. Prepare the second flexure holder in the same way

Fig. 2

6. Cleaning

6.1 Cleaning the pot

1. Cut strips of a double-sided adhesive tape (ca. 2x2 cm)
2. Wind the adhesive tape around one end of a toothpick
3. Introduce the toothpick, adhesive tape first, vertically into the annular gap
4. With circular movements move it repeatedly around the magnet. In this way any dirt particles stick to the adhesive tape
5. If necessary, repeat the procedure with a freshly prepared toothpick until the annular gap is free from impurities
6. Wipe over the top of the pot magnet with a piece of adhesive tape

6.2 Cleaning the coil

1. Using a lint-free cloth, remove any contamination from the coil
2. Do not clean the coil with adhesive tape; risk of injury

6.3 Final check of cleaning


- Visually inspect the annular gap and the coil for any contamination and clean if necessary

7. Assembling the weighing cell

7.1 Assembling the balance arm

1. To facilitate the assembly, the chassis can be placed on the assembly plate → Attention to the sensor mechanism!
2. Pre-assemble assembly jigs [7] + [8] with 2 assembly pins [9] on the chassis as in Fig. 3
3. Adjust assembly jigs [7] + [8] using the alignment jig [14] horizontally in accordance with Fig.3 and tighten the two assembly pins [9] with a No. 4 box spanner

Fig. 3

4. Carefully introduce the balance arm (19/A11) over the pot (34/A11) and past it to the sensor mechanism
5. Assemble the balance arm with 2 assembly pins [9] loosely behind the assembly jigs. The balance arm must remain moveable
6. The right inner assembly pin of the assembly jig [8] can now be firmly tightened
7. The left inner assembly pin of the assembly jig [7] can now be firmly tightened
8. The balance arm are now fixed

7.2 Assembling the coil, the magnet cover and the transport safety device


1. Place the coil onto the balance arm and centre with the centring spigot [11]
2. Screw the coil tightly onto the balance arm with 2 brass screws (21/A11) and 2 washers (27/A11)
3. The coil must be centred in such a way that the centring spigot [11] can be removed smoothly
4. Pre-assemble the magnet cover (8/A11) with 2 screws (23/A11) → screws left and right
5. Pre-assemble the transport safety device (4/A11) with the third screw (23/A11) and spring washer (28/A11) → gently tighten
6. Centring bush (16/C2) shall be introduced between the balance arm-pilot and transportation safety device to assure constant annular gap. Tighten screw (23/A11) and remove centring bush.

7.3 Assembling the sensor mechanism and the stopper

1. Press the sensor mechanism carefully into the cutout in the chassis in the direction of the pot and secure to the chassis with 3 screws (22/A11)
2. Firmly screw the connection board (17/A11) onto the chassis with screw (20/A11)
3. Fix the sensor mechanism plate (5/A11) with 2 screws (20/A11) to the chassis
4. Solder the two connection wires (1/A11) in parallel and strain-free to the connection board (17/A11)
5. Turn back the stopper (29/A11) with brass screw (30/A11) by 90° and screw tight
6. The stopper must be free of play and lie centrally in the symmetry screw (32/A11) → if necessary, re-adjust symmetry screw

7.4 Assembling new flexure sheets


Abb. 4

1. Select 2 new flexure sheets (6/A11) see "How to order spare parts / Anleitung zur Bestellung von Ersatzteilen" on page A2
2. Bring one flexure sheet into position as in Fig. 4 → slots are on the support piece side (below)
3. With a centring spigot [10], centre the flexure sheet in the hole on the chassis side (above)
4. With the second centring spigot [10], centre the flexure sheet in the slot on the support piece side (below)
5. Tighten the first (upper) centring spigot and then remove the second (lower) centring spigot, without touching the flexure sheet in doing so
6. Fasten the flexure sheet on the balance arm side with base (9/A11), spring washer (28/A11) and screw (25/A11) → Fig. 4
7. Remove the second (upper) centring spigot and fasten the flexure sheet on the chassis side with base (9/A11), spring washer (28/A11) and screw (25/A11)

7.5 Assembling the support piece and inner cone

Fig.5


1. Introduce the pre-assembled support piece over the two inner assembly posts [9] → Fig. 5
2. It must be possible to push the support piece smoothly over the two inner assembly pins. In case the support piece cannot be pushed smoothly over the two assembly pins, unscrew pin [9] and tighten it again after.
3. Locate the support piece and the balance arm using the flexure strap anchor [12] and the two screws [13] → Fig. 5
4. The support piece and balance arm are now fixed.
5. Assemble the inner cone (12//A11) to the support piece with 3 screws (23/A11)

7.6 Assembling the flexure holders

1. Introduce the two pins for flexure holder [6] in the first pre-assembled flexure holder (→ see “Installing the new flexure sheets” on page C9) in such a way that the shoulder of the pins is opposite the screw heads on the flexure holder
2. Bring the two pins with the flexure holder into position in the chassis → Fig. 5
3. Lower the flexure holder carefully over the two pins on the chassis, until the flexure sheets just touches the chassis and the support piece
4. Secure the flexure holder with 4 screws (24/A11) and 4 spring washers (28/A11)
5. Push through the two pins up to the installed flexure holder
6. On the opposite side introduce the second pre-assembled flexure holder (→ see “Installing the new flexure sheets” on page C9) with the screw heads inwards over the two pins and lower onto the chassis, until the flexure sheets just touches the chassis and the support piece
7. Secure the flexure holder with 4 screws (24/A11) and 4 spring washers (28/A11)
8. The two pins must be capable of being withdrawn smoothly from the flexure holder

Fig. 6


7.7 Assembling a new flexure strap

1. Select a new flexure strap (3/A11) see “How to order spare parts / Anleitung zur Bestellung von Ersatzteilen” on page A2
2. Unscrew 2 screws [13] on the flexure strap anchor [12] and remove it
3. Bring the flexure strap into position as shown in Fig.7 → the slot is on the support piece side
4. Using a centring spigot [10], centre the flexure strap in the hole on the balance arm side
5. With the second centring spigot [10] centre the flexure strap in the slot on the support piece side
6. Tighten the first (upper) centring spigot and then carefully remove the second (lower) centring spi-

got, without touching the flexure strap in doing so

7. Fasten the flexure strap on the support piece side with base (9/A11), spring washer (28/A11) and screw (25/A11) → Fig. 7
8. Remove the second (upper) centring spigot and fasten the flexure strap on the balance arm side with base (9/A11), spring washer (28/A11) and screw (25/A11)

Fig. 7

7.8 Dismantling the assembly jigs

1. Unscrew and remove the two inner assembly pins [9] on the support piece
2. Unscrew and remove the two outer assembly pins [9] on the support piece
3. Take out the two assembly jigs [7] + [8] from the side
4. If the two assembly pins [6] on the flexure holder (2/A11) are still assembled, remove them carefully


7.9 Installing the weighing cell in the bottom of balance

1. Carry out a further movement check on the balance arm → the balance arm must not touch the transport safety device (4/A11). Visual check with the lens
2. Place the weighing cell carefully in the bottom of balance → hold the weighing cell only to the chassis
3. Tilt the bottom of balance to the side and screw the weighing cell securely from below with 3 nuts (26/A11)
4. Carefully connect the sensor mechanism plug to the main board (8/A5)
5. The balance is now ready for the adjustments on the following pages
6. Insert the calibration weights and mount the counterholder (3/A37) with 4 screws (10/A37) (XT only)

Section D: Service single calibration

1. Service tools and equipment

1.1 Tool set (350-8537)

| Position | Part-No. | Description | Quantity |
|----------|--------------|--------------------------------------|----------|
| 1 | W 41-1688 | Mounting plate | 1 |
| 2 | W 41-1688-5 | Support post, long | 2 |
| 3 | W 41-1688-6 | Support post, short | 2 |
| 4 | 350-6203 | Line-head screwdriver M4 | 1 |
| 5 | 240-7124 | Screwdriver, small, no.1 | 1 |
| 6 | PW 29.43.300 | Pin for flexure holder | 2 |
| 7 | W 41-1700 | Assembly jig with slot | 1 |
| 8 | W 41-1701 | Assembly jig with hole | 1 |
| 9 | W 41-1677 | Assembly pin for assembly jig | 4 |
| 10 | W 41-1691 | Centring spigot for flexures | 2 |
| 11 | W 41-1774 | Centring spigot for the coil EP / ES | 1 |
| 12 | W 41-1675 | Centring spigot for the coil XT / XB | 1 |
| 13 | W 41-1674 | Flexure sheet anchor | 1 |
| 14 | PN 1100-172 | Line-head screw M4x12 | 2 |
| 15 | W 41-1699 | Alignment jig | 1 |
| 16 | W 41-1578-8 | Connecting cable | 1 |
| 17 | W41-1681 | Centring bush | 1 |

1.2 Equipment

1.2.1 Implements

- Soldering iron
- Tweezers
- Lens
- Voltmeter

1.2.2 Standard tools

- Box spanner Nos.4, 5.5 and 7
- Flat-bladed screwdriver Nos.1, 2 and 4
- Phillips screwdriver No.2
- Hexagon socket wrench WS 2.5mm


2. Opening up a balance

2.1 General procedure for opening up a Series 320 XT/XB balance

1. Disconnect the balance from the mains (Pull out mains plug)
2. Remove all windshields, protective rings, weighing pans and pan holders
3. Unscrew 3 cylindrical screws on the top of balance using a No. 4 screwdriver
4. Slightly raise the top of balance
5. Release the cover-foil ribbon cable from the main board

2.2 General procedure for opening up a Series 320 IP65 balance

1. Disconnect the balance from the mains (Pull out mains plug)
2. Remove all windshields, protective rings, balance pans and pan holders
3. Unscrew 3 screws IP56 with a No. 1 screwdriver and remove the clamping ring
4. Unscrew 3 cylindrical screws on the top of balance and 2 cylindrical screws on the bottom of balance using a No. 4 screwdriver
5. Slightly raise the top of balance (Do not damage the sealing)
6. Release the cover-foil ribbon cable from the main board

3. Introduction / Preparation of the tools

3.1 Introduction

In the event of a broken flexure, it is recommended that **all** flexures be replaced. The following assembly instructions provide information on how the flexures can be replaced most efficiently. It is important to employ a structured procedure for dismantling since it is essential to replace the **same parts** (particularly screws) **in the same place** when re-assembling.

The various items of the assembly instructions mean:

(26/B3): component items


[5]: tooling items see "Service tools and equipment" on page D2

Example: 3 nuts (26/B6) means 3 nuts as Item 26 in section B on page 6

3.2 Preparation of the tools

Before dismantling the balance it is necessary to carry out some preparation of the tools. 2 long [2] and 2 short [3] support posts must be screwed into the correct places on the assembly plate

The soldering iron should be plugged in or preheated

4. Dismantling the weighing cell


4.1 Removing the top of balance

- see "Opening up a balance" on page D4

4.2 Removing the calibration unit (XT only)

Take out the plug of the servo-motor (4/A37) on the connection board (12/A5) and unscrew 3 screws (8/A37), lift complete the calibration unit, incl. calibration weight, out of the bottom of balance.

Fig. 1

4.3 Removing the weighing cell from the bottom of balance

1. On the main board, unplug the sensor mechanism plug (15/A11) using a small screwdriver
2. Turn over the bottom of balance and unscrew the 3 nuts (26/A11) on the underside of the housing
3. Hold the weighing cell **only** on the chassis and lift it out of the bottom of balance
4. The chassis can now be placed on the assembly plate for easier handling
5. To protect the weigh cell, the two assembly jigs [7] + [8] are introduced in accordance with Fig. 1 and secured with the 4 assembly pins [9]

4.4 Removing the flexure strap

1. Undo 2 screws (25/A11), remove flexure strap (3/A11) incl. bases (9/A11) and spring washers (28/A11) → Do not separate screw components, assemble in same relationship
2. Screw the flexure strap anchor [12] on tightly with 2 screws [13] in place of the flexure strap

4.5 Removing the two flexure holders

1. Using the two pins for flexure holder [6] fix the upper flexure holder (2/A11) into the chassis
2. Unscrew the 4 screws (24/A11) of the upper flexure holder and remove together with the 4 spring washers (28/A11) → 2 on the chassis, 2 on the support piece
3. Remove the two pins with the upper flexure holder incl. 4 spring washers and the flexure sheets, then take out the pins from the flexure holder
4. Fix the lower flexure holder (2/A11) into the chassis with the two pins for flexure holder [6]
5. Unscrew 4 screws (24/A11) from the lower flexure holder and remove together with the 4 spring washers → 2 on the chassis, 2 on the support piece
6. Take out both pins with the lower flexure holder and remove the lower flexure holder together with the 4 spring washers and the flexure sheets

4.6 Removing the coil

1. Carefully unsolder 2 connecting wires (1/A11) from the connecting board (17/A11) using the tweezers
2. Unscrew 3 screws (23/A11) from the magnet cover (8/A11), remove the cover together with transport safety device (4/A11) and spring washer (28/A11)
3. Unscrew 2 brass screws (21/A11) from the coil (12/A11) on the balance arm (19/A11) using a No.2 screwdriver, remove the coil together with the washer (27/A11)
4. Carefully take out the coil from the pot (34/A11)

4.7 Removing the support piece

1. Unscrew 3 screws (23/A11) on the inner cone (13/A11) and remove inner cone
2. Dismantle the weight holder left (14/A11) and right (15/A11) (XT only)
3. Unscrew upper screw (balance arm side) of the flexure strap anchor [12] and pull off the support piece together with the flexure strap anchor over the assembly pins [9]

4.8 Removing the balance arm

1. Unscrew brass screw (30/A11) from the stopper (29/A11), turn the stopper inwards ca. 90° and gently re-tighten the brass screw
2. Unscrew the 4 screws (25/A11) from the flexure sheets (6/A11) and remove both flexure sheets together with the bases (9/A11) and spring washers (28/A11)

Note: Since the lug of the sensor mechanism is firmly attached to the balance arm, the sensor mechanism must be unscrewed to dismantle the balance arm and gently tilted back:

3. Unscrew 2 screws (20/A11) from the sensor mechanism plate (5/A11) and remove the plate
4. Unscrew screw (20/A11) from the connection board (17/A11)
5. Unscrew 2 screws (22/A11) of the sensor mechanism and gently tilt the sensor mechanism to the rear

6. Remove the two inner assembly pins [9] → unscrew from the balance arm
7. Remove the two outer assembly pins [9] → unscrew from the chassis
8. Remove the two assembly jigs [7] + [8]
9. Carefully tilt the balance arm (19/A11) over the pot (34/A11) and the sensor mechanism and take out
10. Further dismantling of the weighing cell is unnecessary

5. Installing the new flexure sheets


5.1 Removing the old flexure sheets

- Take out all screws from the flexure holder and remove the old flexure sheets together with spring washers

5.2 Installing the new flexure sheets

1. Place a flexure holder on the assembly plate [1] as in Fig.2
2. Select 4 new flexure sheets (7/A11) → Place the new flexure sheets carefully on the assembly plate [1] as in Fig. 2 → the slot of the sheet is always on the flexure holder
3. Place 4 spring washers (28/A11) on the flexure sheet and screw these tightly to the flexure holder with 4 screws (22/A11)
4. Take the complete flexure holder carefully off the assembly plate → Do not touch the flexures!
5. Prepare the second flexure holder in the same way

Fig. 2

6. Cleaning

6.1 Cleaning the pot

1. Cut strips of a double-sided adhesive tape (ca. 2x2 cm)
2. Wind the adhesive tape around one end of a toothpick
3. Introduce the toothpick, adhesive tape first, vertically into the annular gap
4. With circular movements move it repeatedly around the magnet. In this way any dirt particles stick to the adhesive tape
5. If necessary, repeat the procedure with a freshly prepared toothpick until the annular gap is free from impurities
6. Wipe over the top of the pot magnet with a piece of adhesive tape

6.2 Cleaning the coil

1. Using a lint-free cloth, remove any contamination from the coil
2. Do not clean the coil with adhesive tape; risk of injury

6.3 Final check of cleaning


- Visually inspect the annular gap and the coil for any contamination and clean if necessary

7. Assembling the weighing cell

7.1 Assembling the balance arm and the support piece

1. To facilitate the assembly, the chassis can be placed on the assembly plate → Attention to the sensor mechanism!
2. Pre-assemble assembly jigs [7] + [8] with 2 assembly pins [9] on the chassis as in Fig. 3
3. Adjust assembly jigs [7] + [8] using the alignment jig [14] horizontally in accordance with Fig.3 and tighten the two assembly pins [9] with a No. 4 box spanner


Fig. 3

4. Carefully introduce the balance arm (19/A11) over the pot (34/A11) and past it to the sensor mechanism
5. Assemble the balance arm with 2 assembly pins [9] loosely behind the assembly jigs. The balance arm must remain moveable
6. Introduce the pre-assemble support piece (10/A11) over the two inner assembly posts [9] → Fig. 4
7. It must be possible to push the support piece smoothly over the two assembly pins
8. Locate the support piece and the balance arm using the flexure strap anchor [12] and the two screws [13] → Fig. 4
9. The right inner assembly pin of the assembly jig [8] can now be firmly tightened
10. The left inner assembly pin of the assembly jig [7] can now be firmly tightened
11. Firmly tighten the two screws [13] of the fixed flexure strap anchor [12]

12. The support piece and the balance arm are now fixed

Fig. 4

7.2 Assembling the coil, the magnet cover and the transport safety device

1. Place the coil onto the balance arm and centre with the centring spigot [11]
2. Screw the coil tightly onto the balance arm with 2 brass screws (21/A11) and 2 washers (27/A11)
3. The coil must be centred in such a way that the centring spigot [11] can be removed smoothly
4. Pre-assemble the magnet cover (8/A11) with 2 screws (23/A11) → screws left and right
5. Pre-assemble the transport safety device (4/A11) with the third screw (23/A11) and spring washer (28/A11) → gently tighten
6. Centring bush (16/C2) shall be introduced between the balance arm pilot and transportation safety device to assure constant annular gap. Tighten screw (23/A11) and remove centring bush.

7.3 Assembling the sensor mechanism, the inner cone and the stopper


1. Press the sensor mechanism carefully into the cutout in the chassis in the direction of the pot and secure to the chassis with 2 screws (22/A11)
2. Firmly screw the connection board (17/A11) onto the chassis with screw (20/A11)
3. Fix the sensor mechanism plate (5/A11) with 2 screws (20/A11) to the chassis
4. Solder the two connection wires (1/A11) in parallel and strain-free to the connection board (17/A11)
5. Assemble the inner cone (11/A11) to the support piece with 3 screws (23/A11)
6. Turn back the stopper (29/A11) with brass screw (30/A11) by 90° and screw tight

7. The stopper must be free of play and lie centrally in the symmetry screw (32/A11) → if necessary, re-adjust symmetry screw

7.4 Assembling the flexure holders

1. Introduce the two pins for flexure holder [6] in the first pre-assembled flexure holder (→ see “Installing the new flexure sheets” on page C9) in such a way that the shoulder of the pins is opposite the screw heads on the flexure holder
 - Bring the two pins with the flexure holder into position in the chassis → Fig. 5
 - Lower the flexure holder carefully over the two pins on the chassis, until the flexure sheets just touches the chassis and the support piece
 - Secure the flexure holder with 4 screws (24/A11) and 4 spring washers (28/A11)
 - Push through the two pins up to the installed flexure holder
 - On the opposite side introduce the second pre-assembled flexure holder (→ see “Installing the new flexure sheets” on page C9) with the screw heads inwards over the two pins and lower onto the chassis, until the flexure sheets just touches the chassis and the support piece
 - Secure the flexure holder with 4 screws (24/A11) and 4 spring washers (28/A11)
 - The two pins must be capable of being withdrawn smoothly from the flexure holder


Fig. 5

7.5 Assembling new flexure sheets

1. Select 2 new flexure sheets (6/A11) see “How to order spare parts / Anleitung zur Bestellung von Ersatzteilen” on page A2
2. Bring one flexure sheet into position as in Fig. 6 → slots are on the support piece side (below)
3. With a centring spigot [10], centre the flexure sheet in the hole on the chassis side (above)
4. With the second centring spigot [10], centre the flexure sheet in the slot on the support piece side (below)
5. Tighten the first (upper) centring spigot and then remove the second (lower) centring spigot, without touching the flexure sheet in doing so
6. Fasten the flexure sheet on the support piece side with base (9/A11), spring washer (28/A11) and screw (25/A11) → Fig. 6
7. Remove the second (upper) centring spigot and fasten the flexure sheet on the chassis side with base (9/A11), spring washer (28/A11) and screw (25/A11)


Fig. 6

7.6 Assembling a new flexure strap

1. Select a new flexure strap (3/A11) see “How to order spare parts / Anleitung zur Bestellung von Ersatzteilen” on page A2
2. Unscrew 2 screws [13] on the flexure strap anchor [12] and remove it
3. Bring the flexure strap into position as shown in Fig.7 → the slot is on the support piece side
4. Using a centring spigot [10], centre the flexure strap in the hole on the balance arm side
5. With the second centring spigot [10] centre the flexure strap in the slot on the support piece side
6. Tighten the first (upper) centring spigot and then carefully remove the second (lower) centring spigot, without touching the flexure strap in doing so
7. Fasten the flexure strap on the support piece side with base (9/A11), spring washer (28/A11) and screw (25/A11) → Fig. 7
8. Remove the second (upper) centring spigot and fasten the flexure strap on the balance arm side with base (9/A11), spring washer (28/A11) and screw (25/A11)

Fig. 7

7.7 Dismantling the assembly jigs

1. Unscrew and remove the two inner assembly pins [9] on the support piece
2. Unscrew and remove the two outer assembly pins [9] on the support piece
3. Take out the two assembly jigs [7] + [8] from the side
4. If the two assembly pins [6] on the flexure holder (2/A11) are still assembled, remove them carefully

7.8 Installing the weighing cell in the bottom of balance

1. Carry out a further movement check on the balance arm → the balance arm must not touch the transport safety device (4/A11). Visual check with the lens
2. Place the weighing cell carefully in the bottom of balance → hold the weighing cell only to the chassis
3. Tilt the bottom of balance to the side and screw the weighing cell securely from below with 3 nuts (26/A11)

7.9 Installing the calibration unit in the bottom of balance (XT only)

1. Introduce the calibration unit incl. calibration weight (page /A40)

Warning: Do not touch the weighing cell with the calibration unit or the weight!

2. Align the calibration unit in the housing's bottom and secure with 3 screws (8//A38)
3. Reconnect the plug of the servo motor (4//A38) to the connection board (9/A5)
4. /The balance is now ready for the adjustments on the following page

Section E: Adjustments

1. Adjustments

1.1 Activation of the service mode

1.1.1 Set the open balance to the service mode

1. Remove the jumper J1 from the main board


Fig. 1

2. Level the balance with the with the levelling screws, connect the balance to the mains with the power supply
3. The balance will start up and will be in the service mode

1.1.2 Set the closed balance to the service mode, XB Series

1. Level the balance with the with the levelling screws, connect the balance to the mains with the power supply
2. With the small screwdriver [5] from the tool set, briefly activate the service switch on the rear of the bottom of balance → the display will briefly flash on
3. Press the MODE-key twice, „CODE 0000“ appears in the display
4. A short tap on the MODE-key will chance the value of the figure, by long pressure chance to the next figure
5. Insert the code „5236“ and confirm with the PRINT-key
6. The balance is now in the service mode

1.1.3 Set the closed balance to the service mode, XT Series

1. Level the balance with the with the levelling screws, connect the balance to the mains with the power supply
2. With the small screwdriver [5] from the tool set, briefly activate the service switch on the rear of the bottom of balance → the display will briefly flash on

3. Push the SHIFT-key , „CODE 0“ appears in the display
4. Press ENTER-key, „CODE 0000“ appears in the display
5. With the arrow keys insert the code „5236“ and confirm with the ENTER-key
6. The balance is now in the service mode


Fig. 2

1.2 Adjusting of the symmetry

1. Connect the connecting cable [15] from the tool set to plug S3 on the main board
2. Connect the voltmeter to the connecting cable
3. The voltmeter must show approximately the **same voltage** in the end positions of the balance arm (once positive, once negative)
4. If the difference is too great, the position of the balance arm can be adjusted with the symmetry screw
5. Detach connecting cable and voltmeter from the balance


1.3 Checking the pre-load

6. The pre-load error ERROR 60 (below minimal converter range) must not appear when the pan is in position **without load**
7. The pre-load error ERROR 61 (above maximal converter range) must not appear when the pan is in position with **full load**
8. If an ERROR is displayed, the symmetry must be re-adjusted and the transport safety device adjusted from the start → see “Assembling the coil, the magnet cover and the transport safety device” page B 12

1.4 Adjusting of the corner load


1. Place the upper part of the housing on the left alongside the lower part of the housing and connect the ribbon cable from the keyboard membrane to the main board
2. Place the pan holder and the weighing pan on the balance
3. Tare the balance with the TARE key
4. Check with a test weight for corner load → see “Final control” page D 13 the deviation between the two diagonals A-C and B-D → see next page. The diagonal with the **larger absolute deviation** should be adjusted first
5. Adjustment of the corner load is effected by turning the screws (1) and (2) in the direction of the


arrow

Example of corner load adjustment::

1. The diagonal A-C has the larger difference
2. Place a test weight at point C of the weighing pan Tare the balance using the T-key
3. Move the weight to point A
4. There is a difference of 12 digits between points A and C (A=12, C=Tare=0)
5. The displayed value for point A is reduced by turning the corner load screw 1 above point A in a clockwise direction → the difference between A-C becomes smaller
6. There is a difference of 7 digits between points B and D (B=-7, D=Tare=0)
7. The displayed value for point B is reduced by turning the corner load screw 2 above point B in an anti-clockwise direction → the difference between B-D becomes smaller
8. The diagonals A-C and B-D should be repeatedly adjusted until the two diagonals have reached the permitted tolerances for corner load → see "Final control" page D 13
9. Place a test weight on point Z of the weighing pan, tare balance using the T-key
10. Check the points A-B-C-D against the permitted tolerances and correct if necessary
11. The corner load adjustment for balances with a circular weighing pan is similar to that for those with a rectangular weighing pan..

1.5 Checking the lowering process of the calibration weight (XT only)

1. The balance must be in service mode → see “Set the open balance to the service mode” page D 2
2. The ribbon cable from the membrane keyboard must be connected to the main board (S4)
3. Press and hold down the shift key  until "TESTPROG 1" appears on the Info-Display, then release
4. The calibration weight is actuated with the two function keys F3 (UP) and F4 (DOWN)
5. Check the lowering process: the calibration weight must run cleanly into the outer carrier of the weighing cell → carry out a full re-alignment of the calibration unit in the bottom of balance

1.6 Final assembly of the balance

1. Disconnect the balance from the mains, pull out the mains plug
2. Remove weighing pan incl. pan holder from the balance
3. Replace the jumper J1 on the main board → see page E2
4. The ribbon cable from the membrane keyboard must be attached to the main board (S4)
5. Place the top of balance on the bottom of balance and fasten together with 3 screws
6. In the case of IP65 balances: Screw the clamping ring carefully to the inner cone with 3 screws IP65 → The membrane must be installed without wrinkles
7. Install the weighing pan incl. pan holder and, if necessary, the protective ring with the windshield
8. The balance must have attained its operating temperature before further adjustment can be carried out

1.7 Checking the corner load

1. The balance must have reached its operating temperature
2. Check the corner load tolerance → see “Final control” on page E13

1.8 Linearisation of a 320XT Series balance

1. The balance must have reached its operating temperature → push the service switch → see “Set the closed balance to the service mode, XT Series” page D 2
2. Have weights of 1/2 full load and full load ready → see “Final control” on page E13
3. Press and hold down the shift key , enter service access code 5236, until "LINEARITY" appears on the Info-Display, then release the key
4. With empty weighing pan press the function key F1 "0" → balance will flash on
5. Place a 1/2 full load weight in the centre of the weighing pan
6. Press the function key F2 "1/2" → display will flash on
7. Place a full load weight in the centre of the weighing pan
8. Press the function key F3 "1/1" → display will flash on
9. Press the function key F4 "SET" → display will flash on and show the value for the full load
10. Check the 1/4 and 3/4 full load. If the indicated value lies beyond the linearity-tolerance → see “Final control” on page E13, a S-correction must be carried out → see “S-correction of a balance” on page E10
11. The internal calibration weight must now be adjusted → see “Adjustment of the internal calibration weight (XT only)” on page E11

1.9 Linearisation of a 320XB Series balance

1. The balance must have reached its operating temperature → push the service switch → see "Set the closed balance to the service mode, XB Series" page D 2
2. Have weights of 1/2 full load and full load ready → see "Final control" on page E13
3. Press and hold down the MODE, enter service access code 5236, key until "LINEARISATION" appears on the Info-Display, then release
4. The F-key assignment will be shown on the Info-Display (F1=0, F2=1/2, F3=1/1, F4=SET)
5. With weighing pan empty, press and hold down the REF key until "FKT-KEY 1" appears on the Info-Display, then release → balance determines the zero point and display flashes on
6. Place a weight with 1/2 the full load in the centre of the weighing pan
7. Press and hold down the REF key until "FKT-KEY 2" appears on the Info-Display, then release → balance determines the 1/2 load and display flashes on
8. Place a weight with the full load in the centre of the weighing pan
9. Press and hold down the REF key until "FKT-KEY 3" appears on the Info-Display, then release → balance determines the full load and flashes on
10. Press and hold down the REF key until "FKT-KEY 4" appears on the Info-Display, then release → balance determines the linearity and display flashes on
11. Using the ON/OFF key switch the balance off and on again. It will leave service mode

1.10 Final check of the adjustment

1. The balance must have reached its operating temperature
2. Calibrate the balance using the T-key (XB external only, XT internal and external)
3. Check the corner load → see “Adjusting of the corner load” on page E4
4. Check the linearity with half weight and full weight
5. Carry out any S-corrections → see “S-correction of a balance” on page E10


1.11 S-correction of a balance

1.11.1 General notes on S-correction

- The balance must have reached its operating temperature and must be in service mode → see page E2
- In every linearisation of a balance the correction factor is deleted and must, if necessary, be redefined

1.11.2 Measuring the S-correction

1. Place a weight of 1/4 full load in the centre of the weighing pan, note the difference between the actual weight and the weight indicated by the balance
2. Place a weight of 3/4 full load in the centre of the weighing pan, note the difference between the actual weight and the weight indicated by the balance

1.11.3 Calculating the correction factor

- For the correction factor the sign of the measured value must be reversed
for example: + 4 digits were measured → - 4 digits must be entered
for example: - 2 digits were measured → + 2 digits must be entered

1.11.4 Entry of the correction factor for the 320XT Series balances

1. Press and hold down the shift key  until "SCORRECTION" appears on the Info-Display, then release the key
2. Using the function key F2 "SECTOR" the correction range Sa or Sb is established
3. Using the function key F3 "UP" and F4 "DOWN" the chosen correction factor is entered and set with the Function key F1 "SET"

1.11.5 Entry of the correction factor for the 320XB Series balances

1. Press and hold down the MODE key until "SCORRECTION" appears on the Info-Display, then release
2. The assignment of the F-keys is shown on the Info-Display (F1=SET, F2=SECTOR, F3=UP, F4=DOWN)
3. Press and hold down the REF key until "FKT-KEY 2" appears on the Info-Display, then release
→ determine the correction range Sa or Sb
4. Press and hold down the REF key until "FKT-KEY 3" appears on the Info-Display, then release
→ The digit will be added to the pre-selected range
5. A brief depression of the REF key now adds a digit to the selected range
6. Press and hold down the REF key until "FKT-KEY 4" appears on the Info-Display, then release
→ The digit will be subtracted from the pre-selected range
7. A brief depression of the REF key now subtracts a digit from the pre-selected range
8. Press and hold down the REF key until "FKT-KEY 1" appears on the Info-Display, then release
→ Correction factors Sa and Sb will be set
9. Press and hold down the MODE key until "BALANCING" appears on the Info-Display, then release
10. Submit the balance to a further final check → see "Final control" on page E13

1.12 Adjustment of the internal calibration weight (XT only)

1. The balance must have reached its operating temperature
2. Disconnect the balance from the power supply and reconnect, allow the balance to restart normally
3. Push the service switch, see page E2, Fig. 2

NOTE: The jumper J1 must be in place!

4. Have weights of 1/2 full load ready → see "Final control" on page E13
5. Press and hold down the T-key until "CALIBRATION" appears on the Info-Display, then release
6. The process of internal calibration weight adjustment will be initialised "INTERNAL CALIBRATION" will appear on the Info-Display and the balance will show "INT"
7. "EXTERNAL CALIBRATION" will appear on the Info-Display and the balance will first show "0000"
8. If the balance displays, e.g. "4000", it requires a weight of this value
9. Place a weight of e.g. 4000g in the centre of the weighing pan, the balance will flash
10. If the balance again shows "0000", remove the weight from the weighing pan
11. The balance terminates the external calibration and "INTERNAL CALIBRATION" will again appear

on the Info-Display and the process will automatically terminate

12. The balance is switched off and on again and service mode is abandoned

2. Final control

| Balance XB or XT | Weighing range | Readability | Reproducibility (+ / -) | Linearity weight (g) | Linearity tolerance (+ / -) | Corner load weight (g) | Corner load tolerance (+ / -) |
|----------------------------|-----------------------|--------------------|-----------------------------------|--------------------------------|---------------------------------------|----------------------------------|---|
| 120 A | 120 g | 0.1 mg | 0.1 mg | 2 x 50 | 0.2 mg | 50 | 0.2 mg |
| 220 A | 220 g | 0.1 mg | 0.1 mg | 2 x 100 | 0.2 mg | 100 | 0.2 mg |
| 220 A-FR | 220 g | 1 mg | 1 mg | 2 x 100 | 1 mg | 50 | 0.2 mg |
| | 80 g | 0.1 mg | 0.1 mg | | 0.2 mg | | |
| 160 M | 160 g | 1 mg | 1 mg | 2 x 80 | 1.5 mg | 50 | 2 mg |
| 220 M | 220 g | 1 mg | 1 mg | 2 x 100 | 1.5 mg | 100 | 2 mg |
| 320 M | 320 g | 1 mg | 1 mg | 2 x 150 | 1.5 mg | 100 | 2 mg |
| 420 M | 420 g | 1 mg | 1 mg | 2 x 200 | 1.5 mg | 200 | 2 mg |
| 620 M | 620 g | 1 mg | 1 mg | 2 x 300 | 1.5 mg | 200 | 2 mg |
| 620 M-FR | 620 g | 10 mg | 10 mg | 2 x 300 | 10 mg | 100 | 2 mg |
| | 120 g | 1 mg | 1 mg | | 1.5 mg | | |
| 920 M | 920 g | 1 mg | 1 mg | 2 x 400 | 1.5 mg | 300 | 2 mg |
| 1220 M | 1220 g | 1 mg | 1 mg | 2 x 600 | 2 mg | 400 | 3 mg |
| 1220 M-FR | 1220 g | 10 mg | 10 mg | 2 x 600 | 10 mg | 200 | 2 mg |
| | 240 g | 1 mg | 1 mg | | 1.5 mg | | |
| 2220 M-DR | 2220 g | 10 mg | 10 mg | 2 x 1000 | 10 mg | 500 | 3 mg |
| | 520 g | 1 mg | 1 mg | | 2 mg | | |
| 320 C | 320 g | 10 mg | 5 mg | 2 x 150 | 10 mg | 100 | 20 mg |
| 620 C | 620 g | 10 mg | 5 mg | 2 x 300 | 10 mg | 200 | 20 mg |
| 1200 C | 1200 g | 10 mg | 10 mg | 2 x 600 | 10 mg | 400 | 20 mg |
| 2200 C | 2200 g | 10 mg | 10 mg | 2 x 1000 | 15 mg | 1000 | 20 mg |
| 3200 C | 2200 g | 10 mg | 10 mg | 2 x 1500 | 15 mg | 1000 | 20 mg |
| 4200 C | 3200 g | 10 mg | 10 mg | 2 x 2000 | 15 mg | 2000 | 20 mg |
| 4200 C-FR | 4200 g | 0.1 g | 0.1 g | 2 x 2000 | 0.1 g | 1000 | 30 mg |
| | 1200 g | 10 mg | 10 mg | | 15 mg | | |
| 6200 C | 6200 g | 10 mg | 10 mg | 2 x 3000 | 15 mg | 2000 | 30 mg |
| 6200 C-FR | 6200 g | 0.1 g | 0.1 g | 2 x 3000 | 0.1 g | 2000 | 30 mg |
| | 2200 g | 10 mg | 10 mg | | 15 mg | | |
| 3200 D | 3200 g | 0.1 g | 50 mg | 2 x 1500 | 0.1 g | 1000 | 0.2 g |
| 6200 D | 6200 g | 0.1 g | 0.1 g | 2 x 3000 | 0.1 g | 2000 | 0.2 g |
| 6200 D-FR | 6200 g | 1 g | 1 g | 2 x 3000 | 1 g | 1500 | 0.2 g |
| | 1600 g | 0.1 g | 0.1 g | | 0.1 g | | |
| 8200 D | 8200 g | 0.1 g | 0.1 g | 2 x 4000 | 0.15 g | 3000 | 0.2 g |
| 10200 D | 10200 g | 0,1 g | 0.1 g | 2 x 5000 | 0,15 g | 4000 | 1 g |
| 10200 D-FR | 10200 g | 1 g | 0.5 g | 2 x 5000 | 1 g | 3000 | 0.2 g |
| | 3200 g | 0.1 g | 0.1 g | | 0.1 g | | |
| 10200 G | 10200 g | 1 g | 0.5 g | 2 x 5000 | 1 g | 4000 | 1 g |

3. Error messages

3.1 Operating error

(Error disappears with correct use or application)

| ERROR-No. | ERROR responsibility | ERROR repair |
|------------------|--|---------------------------|
| ERROR 1 | invalid program | load right program |
| ERROR 3 | weight too small | use higher weight |
| ERROR 4 | wrong anti-theft-code | use right code |
| ERROR 5 | wrong anti-theft-code | use right code |
| ERROR 8 | zero point outside range | |
| ERROR 9 / 0 | calibration factor bad | repeat calibration |
| ERROR 9 / 1 | zero point difference too high | repeat calibration |
| ERROR 9 / 2 | difference of internal weight too high | repeat calibration |
| ERROR 9 / 3 | calibration factor bad | repeat calibration |
| ERROR 30 | no internal reference weight | acknowledge with TARE-key |
| ERROR 31 | no external reference weight | acknowledge with TARE-key |
| ERROR 47 | too much values | use less values |
| ERROR 48 | other unit then first value | use the original unit |
| ERROR 49 | too less values | use more values |
| ERROR 53 | values out of range +/- 50% | use right weight |
| ERROR 60 | minimum transformer area fallen below | |
| ERROR 61 | maximum transformer area exceeded | |

3.2 Fatale error

(Balance stays still, error repair in service mode)

| ERROR-No. | ERROR responsibility | ERROR repair |
|------------------|---|----------------------|
| ERROR 10 | no linearisation values | repeat linearisation |
| ERROR 11 | no calibration factor | repeat calibration |
| ERROR 21 | no SCS temperature compensation values | |
| ERROR 22 | no temperature compensation values | |
| ERROR 23 | no SCS temperature no-linearity compensation values | |

3.3 Hardware error

(Balance stays still)

| ERROR-No. | ERROR responsibility | ERROR repair |
|------------------|---|-----------------------------------|
| ERROR 14 | wrong model code | call service |
| ERROR 16 | internal RAM destroyed | replace CPU |
| ERROR 17 | external RAM destroyed | replace RAM |
| ERROR 18 / 1 | wrong program check sum | load right program, replace FLASH |
| ERROR 18 / 2 | no approved program build in approved balance | load right program |
| ERROR 18 / 8 | FLASH empty | load right program, replace FLASH |
| ERROR 18 / 9 | wrong controller check sum | replace CPU |
| ERROR 19 | wrong code figure | load right code figure |
| ERROR 20 | temperature measurement defective | call service |

4. Discussion with the manufacturer

In discussions with the manufacturer's works the following information should be provided:

- **Balance type** e.g. XT 220A
- **Order number** e.g. 320-9201 (Label on the rear of the balance)
- **Sales number** e.g. N 43210 (Label on the rear of the balance)
- **Series number** e.g. 2400-100 (Label on the rear of the balance)
- **EPROM number** e.g. K30 00 00 (Label on EPROM inside the balance)
- **Program number** e.g. 351,00,00 (appears on the display when switching on the balance)
- **As precise a description of the defect as possible**
- **Weighing results**

PRECISA GRAVIMETRICS AG
Moosmattstrasse 32
Po-Box 352
CH-8953 Dietikon
Tel. +41 44 744 28 28
Fax. +41 44 744 28 38
e-mail: service@precisa.ch
Internet: <http://www.precisa.com>

Subject to technical modification

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