Angle-stabile

Distal Radius Plate

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Titanium

Variable angle-stable

Königsee Implantate

www.koenigsee-implantate.de

Distributed in the UK by

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Nearly all fractures of the distal radius can be managed surgically with the 2.7 mm palmar angle-stable radius plate for internal fixation. The anatomical shape of the plate in combination with the variable angle-stability of the two distal rows of holes provides stable management of all extra-articular and nearly all intra-articular fractures (C1-3 in the AO classification) without additional measures. Particularly in the case of intra-articular comminuted fractures, the variable angle-stability provides excellent conditions for adjusting the screw position to the individual fracture situation. The small holes incorporated in the plate can be used for preliminary introduction of K-wires, which makes the intraoperative procedure considerably easier. Preliminary K-wire fixation allows the position of the plate relative to the radius and the reduction result to be assessed easily under image intensifier control and corrected if necessary. All of the required plate lengths are available, each in two different widths at the distal end (23 mm and 27 mm).

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Unstable extra-articular and intra-articular extension fractures distal radius (23-A3/C1/C2/C3). Flexion fractures of the distal radius (Smith fracture and reversed Barton fracture - 23-B3).

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
Treatment with 2.7 mm variable angle-stable 7/2-hole radius plate.

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
***Attention:**
For the variable angle-stable locking
use the variable drill-guide sleeve!

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Positioning splint, where applicable, easily removable wrist orthosis with Velcro fastening. Start independent exercises without range of motion restrictions from the first postoperative day. Physiotherapy is required in only a few cases.

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The operation is performed in supine position. The arm is placed on an arm table. Closed reduction is carried out first, if necessary under image intensifier control. In A3 fractures and the majority of C2 fractures, anatomical reduction is largely achieved when ligamentotaxis is maintained.

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The operation is performed without a tourniquet and meticulous haemostasis is obtained with bipolar coagulation. Short (approx. 3-5 cm) incision over the flexor surface of the wrist. Exposure of the flexor carpi radialis tendon, which is retracted slightly toward the ulnar side.

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The forearm fascia is divided immediately below or just to the radial side of the tendon. The flexor carpi radialis tendon and the finger flexor tendons / muscle bellies beneath it are held aside (using a self-retaining retractor if necessary). The median nerve is not exposed and the carpal tunnel is not touched (exception: pre-existing CTS). Division of pronator quadratus with diathermy on the radial side, and the muscle belly is pushed away from the flexor surface of the radius.

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The fracture line on the flexor surface is exposed cleanly by sparing removal of periosteum. This is followed by definitive anatomical reduction.

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
A suitable plate is selected (the two-hole plate is usually sufficient, 23 mm or 27 mm wide, depending on the width of the radius).

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The plate is fixed to the articular fragment and to the shaft with one or more Kirschner wires through the holes provided in the plate. This achieves preliminary fixation of the plate to the radius and of the reduction (image intensifier control).

The plate position is corrected if necessary with further preliminary Kirschner wire fixation and subsequent image intensifier control.

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During drilling, the opposite cortex is deliberately drilled cautiously with a finger on the dorsal side of the radius. The length is measured through the angle-stable drill bushing or after using the variable angle-stable drill bushing with corresponding distance sleeve on the length gauge, palpating the tip of the gauge on the dorsal side of the radius beneath the skin so that the length of the screws can be established exactly.

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The screws introduced distally should just reach the dorsal cortex but not push away smaller cortical fragments (which are common when there is a dorsal zone of comminution). It is therefore beneficial, as mentioned above, to drill the thin dorsal cortex carefully first. In this technique, at least 3 or 4 screws are introduced in succession into the articular block and the radial screw in particular should diverge into the radial styloid process depending on the inclination. In the case of intra-articular fractures, a screw position in the fracture gap should be avoided.

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Important: when tightening the distal angle-stable screws, counterpressure should be exerted with a finger from the dorsal side in order to press the bone against the plate. That way, the anatomical shape of the plate is used for exact reduction.

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Final image intensifier control and documentation in two planes.

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After X-ray documentation and irrigation, insertion of a vacuum drain and skin closure.

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*Treated with 5/2-hole variable
angle-stable radius plate*

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*Extension fracture with
impression of the articular
surface, type C3*

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
*Treated with 7/2-hole variable
angle-stable radius plate*

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The postoperative CT shows an anatomically reconstructed articular surface without a step, with the screws in correct, strictly subchondral (extra-articular) position.

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*Range of motion on the 3rd postoperative day, early functional follow-up treatment.
No immobilisation!*

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Extension fracture type C2 with depressed ulnar marginal fragment

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Reduction of the ulnar marginal fragment with Kirschner wire from the dorsal side („joy stick“), at the same time preliminary fixation of this fragment with a Kirschner wire from the palmar aspect through the plate.

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„Avoiding“ the intra-articular fracture gap by exploiting the variable angle stability. Treated with variable angle-stable radius plate 7/3-hole.

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

Juxta-articular forearm fracture (70-year old patient, radius flexion fracture). Treated with 5/3-hole plate on the radius and 2.0 mm distal angle-stable ulna plate.

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
Comminuted fracture (C3) 38-year old patient. Healed appearance 6 months postoperatively, metal removal planned.

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


Juxta-articular forearm fracture with C2 fracture of the radius, preexisting bayonet deformity after old radius fracture healed in deformity (78-year old patient). Corrective angle-stable internal fixation of the radius with 7/3-hole plate and angle-stable internal fixation of the fracture of the ulna head with 2.0mm distal ulna plate.

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
Distal forearm comminuted fracture, high-grade osteoporosis (85-year old patient). Long bridging palmar supporting angle-stable internal fixation of the radius with 5/8-hole plate, postoperative immobilisation in a cast.

**Variabel winkelstabile distale Radiusplatte, schmal***Variable angle-stable distal radius plate, small*


Bestell-Nr. Code N° titanium	Länge Length	Lochanzahl N° of hole	Anzahl im Set Quantity in set	
5.9823.52	rechts / right	43 mm	5/2	-
5.9823.53	rechts / right	51 mm	5/3	1
5.9823.54	rechts / right	60 mm	5/4	1
5.9823.56	rechts / right	77 mm	5/6	-
5.9823.58	rechts / right	94 mm	5/8	-
5.9823.510	rechts / right	111 mm	5/10	-
5.9823.513	rechts / right	136 mm	5/13	-
5.9833.52	links / left	43 mm	5/2	-
5.9833.53	links / left	51 mm	5/3	1
5.9833.54	links / left	60 mm	5/4	1
5.9833.56	links / left	77 mm	5/6	-
5.9833.58	links / left	94 mm	5/8	-
5.9833.510	links / left	111 mm	5/10	-
5.9833.513	links / left	136 mm	5/13	-

Variabel winkelstabile distale Radiusplatte, breit*Variable angle-stable distal radius plate, wide*

Bestell-Nr. Code N° titanium	Länge Length	Lochanzahl N° of hole	Anzahl im Set Quantity in set	
5.9843.72	rechts / right	43 mm	7/2	-
5.9843.73	rechts / right	51 mm	7/3	1
5.9843.74	rechts / right	60 mm	7/4	1
5.9843.76	rechts / right	77 mm	7/6	-
5.9843.78	rechts / right	94 mm	7/8	-
5.9843.710	rechts / right	111 mm	7/10	-
5.9843.713	rechts / right	136 mm	7/13	-
5.9853.72	links / left	43 mm	7/2	-
5.9853.73	links / left	51 mm	7/3	1
5.9853.74	links / left	60 mm	7/4	1
5.9853.76	links / left	77 mm	7/6	-
5.9853.78	links / left	94 mm	7/8	-
5.9853.710	links / left	111 mm	7/10	-
5.9853.713	links / left	136 mm	7/13	-

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Optionally only angle-stable distal radius plates are available in narrow and wider version as well.

Kortikalisschrauben, D 2,7 mm winkelstabil selbstschneidend*Cortical screws, diameter 2.7 mm, angle-stable self-tapping*

Bestell-Nr. Titan Code N° titanium	Länge Length	Anzahl im Set Quantity in set
3.125.12	12 mm	5
3.125.14	14 mm	5
3.125.16	16 mm	5
3.125.18	18 mm	5
3.125.20	20 mm	5
3.125.22	22 mm	5
3.125.24	24 mm	5
3.125.26	26 mm	5
3.125.28	28 mm	5

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
No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying, recording or otherwise without the prior written permission of the publisher.

Supplied by:

ATHRODAX
HEALTHCARE INTERNATIONAL LTD

Hawthorn Business Park, Drybrook,
Gloucestershire GL17 9HP
Telephone: (01594) 544440
Fax: (01594) 545800
Email: sales@athrodax.co.uk
Website: www.athrodax.co.uk

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Ltd. OA Dr. med. Klaus-Jürgen Maier