

Cardiac Rhythm Management // Sentus ProMRI® QP

# Sentus ProMRI®

The Smallest MRI QP Lead



.... ———  
**QP ExCELS**

# Clinical Safety and Performance of the Sentus OTW Quadripolar LV Lead

---

## QP ExCELS Study Objective

Confirm the safety and effectiveness of the BIOTRONIK Sentus QP LV lead.

## Study System

- Sentus OTW QP L, Sentus OTW QP S, and Sentus OTW QP S/49 LV lead
- Sentus OTW QP S/49 LV, Sentus OTW QP L/49 and Sentus OTW QP L/49 LV leads were not implanted or analyzed as part of the pre-market study
- BIOTRONIK CRT-D system with IS4 LV port

## Study Design

- Prospective, non-randomized, multi-center, international study
- Combined Pre-Market Investigational Device Exemption (IDE) and Post-Approval Study regulated by Food and Drug Administration
- Pre-Market Analyses Cohort: All subjects implanted (or implant attempt) on or before January 29, 2016 (N=299)



# 97.1% Sentus Complication-Free Rate Through 6 Months

## Clinical Goal

Evaluate the Sentus QP related complication-free rate through 6 months post-implant.

## Clinical Results

Eight of the 279 subjects included in the primary endpoint 1 analysis experienced endpoint-related events as determined by the CEC resulting in a complication-free rate of 97.1% (271/279),  $p < 0.0001$ , 95% CI: (94.4, 98.8)\*.

### Primary Endpoint 1 Adverse Event Summary

Adverse Event Type	Subjects with AE	% with AE	Number of AEs	Rate (per subject year)
Related to Sentus QP Lead and Meets Primary Endpoint Criteria				
Lead dislodgement	4	1.43%	4	0.028
Extracardiac stimulation	3	1.08%	3	0.021
Lead impedance out of range, high	1	0.36%	1	0.007
Total	8	2.87%	8	0.055

Total number of implanted subjects = 279; subject years since successful implant = 144.7

Minimal complication rate at 3 Months.

\*Compared to performance goal of 90%


# 93.4% of Patients Experienced Acceptable LV Pacing Threshold ( $\leq 2.5V$ @ 0.4ms)

## Clinical Goal

Evaluate the percentage of subjects with acceptable pacing threshold in permanently programmed vector at 3 months.

## Clinical Results

A total of 267 of the 286 subjects included in primary endpoint 2 analysis had an acceptable LV Pacing Threshold resulting in a rate of 93.4% (267/286),  $p=0.002$ , 95% CI: (89.8, 96.0)\*.



Acceptable LV Threshold at 3 Months.

# 95.6% of Subjects Had at Least One Novel Vector with an Acceptable Threshold at 3 Months

## Clinical Goal

Evaluate Sentus QP pacing threshold in novel vectors at 3 months.

## Clinical Results

Number of novel vectors with acceptable threshold at 3 months.

Novel Vectors with Acceptable Threshold	Number of Subjects	% (N=155)
7 of 7 Novel Vectors	9	5.8%
6 of 7 Novel Vectors	22	14.2%
5 of 7 Novel Vectors	28	18.1%
4 of 7 Novel Vectors	16	10.3%
3 of 7 Novel Vectors	17	11.0%
2 of 7 Novel Vectors	19	12.3%
1 of 7 Novel Vectors	32	20.6%
0 of 7 Novel Vectors	12	7.7%

52.1% (147/282) entered 3M visit with a traditional vector. Of these 91.2% had one or more novel vectors with an acceptable threshold

Pacing Vector	Mean Threshold (V)+/- SD	Number of Subjects with Pacing Threshold Tests Conducted at 0.4 ms	No Capture Possible N (%)	Acceptable LV Pacing Threshold N (%)
LV1 tip to LV4 ring	1.4 +/-1.26	221	14 (6.3%)	168 (76.0%)
LV2 ring to LV4 ring	2.1 +/-1.71	197	15 (7.6%)	136 (69.0%)
LV3 ring to LV2 ring	3.0 +/-1.82	184	27 (14.7%)	73 (39.7%)
LV3 ring to LV4 ring	2.8 +/-1.85	189	46 (24.3%)	82 (43.4%)
LV3 ring to RV coil	2.2 +/-1.75	197	31 (15.7%)	121 (61.4%)
LV4 ring to LV2 ring	3.6 +/-2.04	201	37 (18.4%)	58 (28.9%)
LV4 ring to RV coil	3.3 +/-2.45	208	67 (32.2%)	72 (34.6%)

# Sentus QP vs. the Competition

The Sentus QP lead demonstrates the best complication-free rate at 6 months, a very low LV lead dislodgement rate, and comparable rates of successful thresholds  $\leq 2.5V$ .

	Primary Endpoint 1		Primary Endpoint 2
	Complication-Free Rate at 6 Months	LV Lead Dislodgement Rate	Threshold $\leq 2.5V$ at 3 Months
BINC – Sentus QP <sup>1</sup>	97.1%	1.43%	93.4%
BSX – Acuity X4 <sup>2</sup>	96.5%	1.30%	94.0%
MDT – Attain Performa <sup>3</sup>	96.0%	2.05%	93.9% <sup>4</sup>
STJ – Quartet <sup>5</sup>	96.0% <sup>6</sup>	3.50%	90.4%

1. Sentus OTW QP Lead Family Technical Manual.

2. BSX – Acuity X4 straight model: Boston Scientific. Clinical Summary Navigate X4 Study (with the AcuityTM X4 lead family). January 2016.

3. MDT – Attain Performa model 4298: Medtronic, ATTAIN® PERFORMA® MODELS 4298, STRAIGHT 4398, AND S 4598 QUADRIPOlar LEFT VENTRICULAR LEADS Safety and effectiveness information. December 17, 2014.

4. At 6 months.

5. STJ – Quartet model 1458Q: Tomassoni G, et al. Postoperative performance of the Quartet® left ventricular heart lead, J Cardiovasc Electrophysiol, 2013, vol. 24 [pg. 449 -56].

6. Through 3 months.

# QP ExCELS Study Design


## Patient Population

### Key Inclusion Criteria

- Standard CRT-D indication according to clinical routine
- De novo implantation or upgrade from existing ICD or pacemaker implant (with no prior attempt at LV lead placement) utilizing a BIOTRONIK CRT-D system with IS4 LV port and Sentus QP LV lead


### Key Exclusion Criteria

- Contraindication to CRT-D therapy
- Currently implanted with an endocardial or epicardial left ventricular lead or had prior attempt to place a left ventricular lead
- Cardiac surgical procedure, such as coronary artery bypass graft, valve surgery, or ablation that is planned to occur within 6 months after implant (ablations planned to occur prior to or at implant are not exclusionary)



# Key Results

---



**Freedom from Sentus QP complications**  
through 6 months post-implant

---

Data demonstrates  
**safety and effectiveness**  
of Sentus QP lead

---

All pre-market study primary safety and  
performance endpoints were  
**met with statistical significance**