Operating Instructions

DKL
YOUR DENTAL UNIT

Manufacturer:
DKL CHAIRS GmbH
An der Ziegelei 3
D-37124 Rosdorf
+49 (0)551-50060
info@dkl.de
www.dkl.de

L2-SUC

Form no. GEN-L2SUC
Rev. 02 / 01.08.2021
Subject to change without
notice.

Contents

Symbols in the Operating Instructions	3
Symbols inside the Unit	4
Introduction	5
Electromagnetic Compatibility (EMC)	6
Safety Instructions - General	7
Safety Notes – Assistant’s Device and Cuspidor	8
Technical Specifications	9
Media Requirements	10
Product Description	13
Working Properly on the L2	15
Moving the Treatment Chair	16
Putting the Treatment Unit into Operation	18
Function Syringe	19
LED.light Operating Lamp	21
Tray	23
Cuspidor	24
Suction System	25
Cleaning and Disinfection of the Suction Handpieces	27
Cleaning and Disinfection of the Suction System	29
Cleaning and Disinfection of the Surfaces	31
Cleaning and Disinfection of the Instrument Holders	34
Bottle Care System	35
Flushing Function with the Bottle Care System	36
Intensive Flushing with the Bottle Care System	36
Water Separation Unit (WSU)	37
Flushing Function with the Water Separation Unit	39
Intensive Flushing with the Water Separation Unit	39
Filling DK-DOX 150 into the Water Separation Unit	40
Maintenance and Inspection	41
Safety Inspections	41
Waste Disposal	42
EMC - Manufacturer’s Declaration	43

Annex

- > Operating Instructions DÜRR Dental Spittoon valve 3
- > Operating Instructions DÜRR Dental CAS 1 Combi-Separator
- > Operating Instructions DÜRR Dental CS 1 Combi-Sepamatic
- > INSTRUCTION MANUAL FARO DENTAL SYRINGE SYR3
- > Use Manual Luzzani Syringe Minilight
- > INSTRUCTION MANUAL FARO ALYA DENTAL LED LIGHT

Symbols in the Operating Instructions

	WARNING! (risk of injury)		CAUTION! (to prevent damage occurring)		General explanations, without risk to persons or objects
	Thermodisinfectable		Sterilisable up to the specified temperature		Call customer service!

Symbols on the unit

	Observe the operating instructions!		ON / OFF		Do not dispose of with household waste.
	CE-marking with identification number of the notified body		Foot controller		Type B application part
	Manufacturing date		Model (designation)		Serial number
	Manufacturer		Medical device		UDI - Product Identification
V	Electrical voltage	AC	Alternating current	VA	Electrical power consumption
A	Current intensity	Hz	Frequency of alternating current		
	Electrical fuse				

Symbols inside the Unit

Earth conductor
connection - Protective
earth


Functional
earth


N

Connection point for neutral
conductor


Symbols on the packaging

Air humidity, limitation

Permissible temperature range

Transport upright; top

Protect from moisture!

Do not stack!


Fragile

Introduction

For your safety and the safety of your patients


These operating instructions are intended to explain how to use your product. However, we must also warn of possible dangerous situations. Your safety, the safety of your team and, of course, the safety of your patients are very important to us.

Please observe the safety instructions!

Intended purpose

This treatment unit is used for the diagnosis and therapy of children and adults in the field of dentistry.

Improper use can damage the treatment unit and thus pose risks and hazards to the patient, user and third parties.

Qualification of the user

The DKL treatment unit may only be used after medically, professionally and practically trained personnel have been instructed. The development and design of the treatment unit were geared towards the target group of dentists, dental hygienists, qualified dental employees (prophylaxis) and dental assistants.

Production according to EU directive
The medical device complies with the provisions of Directive 93/42/EEC.

Responsibility of the manufacturer

The manufacturer can only be held responsible for the impact on the safety, reliability and performance of the treatment unit if the following instructions are observed:

- > The dental unit must be used in accordance with these operating instructions.
- > If assembly, additions, new settings, changes or repair work is carried out by DKL or trained technicians authorised by DKL or personnel of authorized dealers trained by DKL.
- > The electrical installation of the room must comply with the regulations of the IEC 60364-7-710 standard („Erection of electrical installations in rooms used for medical purposes“) or comply with the regulations applicable in your country.
- > The recommended annual maintenance is carried out and any repair work in this context meets the requirements of EN 62353.
„Repeat tests and pre-commissioning tests of medical electrical equipment and systems – general regulations“are fully complied with.
- > The national legal regulations are observed when using the device, in particular the applicable health and safety regulations and accident prevention measures.

Electromagnetic Compatibility (EMC)

Medical electrical equipment is subject to special precautions with regard to EMC and must be installed and commissioned in accordance with the EMC instructions. DKL guarantees that the dental unit complies with the EMC guidelines only if original DKL accessories and spare parts are used. The use of accessories and spare parts not approved by DKL may lead to an increased emission of electromagnetic interference or to a reduced resistance to electromagnetic interference.

The EMC manufacturer declaration can be found on page 43.

HF communication equipment

Do not use portable and mobile HF-communication equipment (such as mobile telephones) during operation. These can affect medical electrical devices.


Risks due to electromagnetic fields

The functionality of implantable systems, such as cardiac pacemakers and implantable cardioverter defibrillators (ICD), can be influenced by electric, magnetic and electromagnetic fields.

- > Before using the product, ask the patient and user about implanted systems and check the use.
- > Perform a risk-benefit analysis.
- > Do not place the product near implanted systems.
- > Do not place the instruments on the patient's body.
- > Take appropriate emergency precautions and respond immediately to health changes.
- > Symptoms such as increased heart rate, irregular pulse, and dizziness may be signs of problems with a pacemaker or ICD.

Safety Instructions - General

- > Before being put into initial operation, the treatment unit must be kept at room temperature for 24 hours.
- > Before each application, check the treatment unit and the instruments with cables for damage and loose parts.
- > Do not operate the treatment unit if it is damaged.
- > Check the set parameters each time you restart the unit.
- > Carry out a test run before each application.
- > The application and timely shutdown of the system is the user's responsibility.
- > Make sure that in the event of a device or instrument failure, the treatment can be completed safely.
- > Use only original DKL fuses.
- > Never touch the patient and the electrical connection at the treatment unit at the same time.
- > Do not lean on the doctor's device, the assistant's device, the tray or the operating lamp.
- > When moving the treatment chair, the doctor's device, the assistant's device, the tray or the operating lamp, pay attention to the patient and the practice personnel.
- > Always switch off the treatment unit before leaving the practice.

Hygiene and care before using the device

- > Clean and disinfect the device immediately before or after each treatment!
- > Wear protective clothing.

Observe your country-specific guidelines, standards and specifications for cleaning, disinfection and sterilisation.

The treatment unit is classified as an „ordinary device“ (closed device without protection against water ingress).

The treatment unit is not suitable for use in an explosive atmosphere or in explosive mixtures of anaesthetics with oxygen or nitrous oxide.

The treatment unit is not suitable for use in rooms with an oxygen enriched atmosphere.

Safety notes – Patient chair

- > Not suitable for patients who cannot remain in a resting position due to mental or physical disabilities.
- > The patient's arms and legs must rest on the upholstered parts of the chair.
- > Do not exceed the maximum patient weight of 150 kg.
- > Do not sit on the head or foot rest of the horizontally aligned patient chair.
- > Position changes must always be carried out under the surveillance of the person giving treatment.
- > Watch patients while moving the treatment chair.
- > Make sure that there are no objects under the treatment chair.

Safety Notes – Assistant's Device and Cuspidor


- > Before putting the device into initial operation and after downtimes (weekends, (public) holidays etc.), flush the water lines intensively.
- > Flush all instrument connections for 2 minutes before starting work.
- > Press the glass fill button several times before starting work.
- > Rinse used instruments for 20 seconds after each treatment.

Risk of injury or infection caused by instruments that are not in use:
The arrangement of the instruments may cause injury or infection to the hand and forearm when accessing the tray or the display.
Therefore, when accessing the tray or the display, pay attention to the arrangement of the instruments.

Highly immunosuppressed patients or patients with specific lung diseases should not get into contact with the water of the treatment unit.
It is recommended to use sterile solutions.

- > Do not exceed the maximum permissible load of 2 kg on the Assistant's device.

Technical Specifications

The motors of the treatment unit are designed for intermittent operation in accordance with the dental treatment method.

Driving motors for patient chair and backrest: duty cycle (max. 25 s „ON“ / 400 s „OFF“)

Supply voltage	230V AC
Nominal voltage	max. 4 A
Frequency	50/60 Hz
Fuse	T 6.3 A H 250 V primary/ time-delay
Maximum power consumption	800 VA
Device class according to the 93/42/EEC Directive	Ila
Protection class	Device of protection class I
Contamination level	2
Over voltage category	II
Power cable	3x1,5 mm ²
Suction control lines to the suction device	5x1,5 mm ²
Potential equalisation	1x 4 mm ²
Relay control line optional special function	3x1,5 mm ²
Free end electrical cables above floor	500 mm
Fuse for domestic installation	Circuit breaker: 16 A medium-lag Recommendation: circuit breaker type C
Degree of protection against ingress of water	Ordinary device (without protection against water ingress)

Permanently connected device. In order to avoid the risk of electric shock, this device may only be connected to a power supply with an earth conductor.

Weight	
L2-SUC	max. 200 kg

Transport and storage conditions	
Ambient temperature	30 to +70 °C
Relative humidity	10 to 80 %
Atmospheric pressure	500 hPa to 1060 hPa

Operating conditions	
Quality and load-bearing capacity of the floor	The floor must be level and horizontal according to EN 18202. Unevenness of the floor along the total length of the chair base up to 2 mm is acceptable. The minimum load-bearing capacity of the floor must be 0.5 N/cm ² (equivalent to approx. 500 kg/m ²).
Ambient temperature	10 to 35 °C
Relative humidity	15 to 80 %
Atmospheric pressure	700 hPa to 1060 hPa
Installation site	≤ 3,000 m above sea level The treatment unit is not suitable for operation in hazardous areas.

Media Requirements

Media water	
Water hardness	1.5 to 2.14 mmol/l
ph-value	6,5 to 8,5
Water filtration on site	≤ 100 µm
Water inflow	Pipe 10x1mm, angle valve outlet 3/8"
Water connection above floor	min. 40 mm, max. 60 mm
Water inlet pressure	2.0 to max. 6.0 bars
Water quality	Cold water in accordance with local and national drinking water regulations.
Minimum flow rate	3 l/min

- Perform the installation according to the national installation requirements (e.g. EN 1717).
- For the reduction of microorganisms in the water supply pipe, please observe the following when laying this pipe to the treatment unit:
 - Avoid long stub lines to the treatment unit.
 - Select the installation in such a way that other essential consumers (e.g. washbasin) are as far as possible behind the connection of the treatment unit can be supplied from the same pipe.
 - Avoid laying the hot water supply pipes in parallel.
- Recommendation: For the water supply of the treatment unit, install an angle valve with 2 outlets and 2 stop cocks. The second outlet allows easy sampling of water for microbiological analysis.

Connection to the public drinking water supply

The treatment unit with a water separation unit complies with the requirements of EN 1717 (free outlet with separation distance ≥ 20 mm) and the DVGW (German Technical and Scientific Association for Gas and Water). It is intrinsically safe in accordance with worksheet W540 and therefore also meets the requirements of W270 and KTW (guideline for hygienic assessment of organic materials in contact with drinking water).

When the treatment unit is equipped with a cuspidor, the bowl rinser ensures the free outlet with a separation distance ≥ 20 mm. When the treatment unit is equipped with a Bottle Care System, the spray supply of the instruments is separate from the public water supply.

Before the treatment unit is installed, the microbiologically perfect water quality of the domestic water supply should be ensured and documented in the form of a microbial count. Sampling and microbial count should be carried out by a competent laboratory.

Media air	
Air inlet pressure	max. 7 bars
Air consumption	80 NI/min
On-site air filtration	≤ 100 particles size 1 - 5 µm referred to one m ³ of air
Oil content	≤ 0.5mg/m ³ , oil-free compressors; the compressor must suck in hygienically perfect air.
Humidity	Pressure dew point ≤ -20 °C at atmospheric pressure
Compressed air supply	Pipe 10x1 mm, angle valve outlet 3/8"
Air connection above the floor	min. 40 mm, max. 60 mm

Clean air and water pipes before installing the unit.

Chips and other foreign substances could be flushed or blown into the treatment unit.

Metal chips can interfere with the function of pneumatic components. Filters are clogged by foreign substances.

- When assembling, make sure that there are no chips or other foreign substances in the pipes.
- Flush the water pipes.
- Blow out the air ducts.
- Make sure that no further foreign substances get into the pipes and ducts after rinsing or blowing out.

Requirements for the suction system	
Vacuum at supply connection	min. 0,12 bar, max. 0,18 bar
Minimum suction power at supply connection	≥750NL/min
Suction system	Type 1: high flow rate wet or dry suction
Diameter of suction handpieces:	small suction handpiece: 6 mm large suction handpiece: 16 mm
Suction pipe	DN40 HT-PP (polypropylene, inside diameter approx. 36.5 mm)
Water drain	DN40 HT-PP (polypropylene, inside diameter approx. 36.5 mm)
Gradient	Min. 10 mm per metre
Wastewater volume	3 l/min


Filter in the treatment unit		maintenance interval	Article number
Particle filter water inlet	50 µm	Replace annually	210623-WBG
Particle filter compressed air inlet	50 µm	Replace annually	210623-LBG
Solid particle filter for the suction system	Mesh size 1 mm	In case of damage, replace at least annually.	514100

Typical pressure in the suction system


Spray mist suction	Vacuum / mbar
90 NL/min	22,6
150 NL/min	38,2
200 NL/min	60,0
250 NL/min	88,8
300 NL/min	124
316 NL/min	137
Saliva ejector	
50 NL/min	100,0
55 NL/min	120,0
60 NL/min	135,2
67 NL/min	162
80 NL/min	200

Requirements for Supply Connections


①	Air: pipe min. 10x1 mm, angle valve outlet 3/8"
②	Water drain DN40 HT-PP
③	Suction line DN40 HT-PP
④	Power cable 3x1.5 mm ²
④	Equipotential bonding 1x 4 mm ²
④	Control line to suction device 5x1.5 mm ²
⑤	Water: pipe min. 10x1 mm, angle valve outlet 3/8" (2 outlets with 2 stopcocks)


Dimensions in Millimetres

L2-SUC Version with short back rest

L2-SUC Version with long back rest

Treatment chair:
lowest position 570 mm
highest position 800 mm


Product Description

①	Treatment chair, seat
②	Treatment chair, back rest
③	Treatment chair, double-jointed head rest
④	Treatment chair, joystick
⑤	Assistant's device
⑥	Cuspidor
⑦	Tray
⑧	Operating lamp LED.light

Product Description

①	Spray mist suction
②	Saliva ejector
③	Syringe



Connecting the instrument hoses and suction tubes.

The instrument hoses can be connected or disconnected via a plug connection under the assistant's device. Depending on the equipment of the model, the assistant's device is equipped with the following instruments (from the left): the syringe, Spray mist suction and finally the Saliva ejector. The instrument holders are marked on the back. Do not connect the tubes crosswise.

Marking	Type of holder
GS	Spray mist suction
KS	Saliva ejector
S	Syringe

Treating patients in a lying position gives you a clear view of the upper and lower jaw. This optimal working posture is ensured by sufficient legroom. During lower-jaw treatments, the patients are positioned horizontally; for treatments in the upper-jaw region the patient's head is slightly overstretched backwards. The double-jointed headrest offers stable support for the patient's head at the starting point of the neck muscles. This way the patient is relaxed and does not need to constantly correct his or her position. Treatment of patients in a lying position offers you an optimal view and, consequently, optimal results.

Moving the Treatment Chair

Joystick at the chair base

Move the joystick downwards

>Treatment chair moves downwards.

Move the joystick upwards

>Treatment chair moves upwards.

Move the joystick to the left

>Back rest tilts backwards.

Move the joystick to the right

>Back rest moves to an upright position.

Joystick

Tap the Joystick Twice

>Treatment chair moves to the entry /exit position.

Tap the joystick twice downwards

>Treatment chair moves to the treatment position p2.

Tap the Joystick Twice to the Left

>Treatment chair moves to the rinsing position and back to the last position when tapping twice again.

Tap the Joystick Twice to the Right

>Treatment chair moves to the treatment position p1.

Programme keys at the back rest

Press "p0"

>Treatment chair moves to the position for getting on/off.

Press "p1"

>Treatment chair moves to treatment position p1.

Press "p2"

>Treatment chair moves to treatment position p2.

Press "p3"

>Treatment chair moves to treatment position p3.

Briefly press "lp"

>Treatment chair moves to the rinsing position and after renewed pressing back to the "last position".

Saving programme keys p0 – p3

For programming the keys, move to the desired programme position manually and then keep the respective programme key pressed for about 3 seconds until you hear a signal tone. Now you have successfully saved your individual treatment position.

Saving programme key LP

Press the "lp" key to move the chair to the rinsing position. By pressing the key "lp" once again, the chair moves back to the previous programme position or to the manually set position. In order to programme the chair, move it manually to the desired rinsing position and then keep the key "lp" pressed for about 3 seconds until you hear a signal tone. Now you have successfully saved your rinsing position.

Moving the Treatment Chair

Emergency-Stop System

Briefly tap the joystick or any programme key at the back rest to stop the active programme immediately. For this purpose, you can also press the start button for the instruments at the foot controller.

You can find an application film at www.youtube.com/DKL_Germany.

Video: DKL CHAIRS L2 SERIES FUNKTIONEN TREATMENT CHAIR POSITIONS

LINK: <https://youtu.be/qzRuFbGAuA0>

The treatment chair can optionally be equipped with arm rests. The arm rest on the assistant's side is firmly screwed on. The arm rest on the doctor's side can be removed by loosening the two knurled screws on the back of the back rest.

Pull the release lever up to freely adjust the double-jointed head rest. Press the lever down to fix the position of the double-jointed head rest. The release lever must always be pulled up completely to move the head rest.

Manually Extractable Double-Jointed Head Rest

Programme Run

Press the programme keys "0" and "lp" simultaneously for about 3 seconds until you hear a signal tone. The treatment chair then moves the seat and the back rest into the lowest position. The motion sequence has been reset now.


Putting the Treatment Unit into Operation

Putting the Unit into Operation

Before putting your treatment unit into initial operation, an intensive flushing must be carried out [→page 36 or page 39].

Activating / Deactivating the Treatment Unit


The treatment unit is equipped with a power switch ① on the chair base. The power switch connects the treatment unit with the power supply. In the event of longer downtimes, the treatment unit should be disconnected from the power supply. The treatment unit contains a device fuse ②. Switch on the treatment unit at the power switch. The power switch lights up green.

Function FARO Syringe SYR3


See also operating instructions of the FARO company for their SYR3 dental syringe!

Press the green button to release water.
Press the light blue button to release air.
Press both buttons simultaneously to release spray.

To remove the sleeve, press button „A“.
Pull out the sleeve in the direction of the arrow.

To remove the nozzle, screw off ring „B“ and pull out the nozzle in the direction of the arrow.

Before reassembling the syringe, lubricate the 2 seals on the inner body and the metal tube of the nozzle with Vaseline grease.

Cleaning the syringe

Clean the sleeve and nozzle under running water or using alcohol-based solutions.

Sterilisation

- Remove the sleeve including the nozzle from the syringe body.
- Clean the sleeve and the nozzle.
- Sterilize the sleeve including the nozzle in the autoclave.

Cleaning the spouts

If the nozzle from which the water and the air come is clogged, clean it with the supplied probe.



Sleeve and spout sterilizable

Function Luzzani Syringe Minilight


See also operating instructions of the Luzzani company for their Minilight syringe!

To blow water into the operating field, just press the left button on the handpiece, symbol: water. To insufflate air into the operating field, just press the right button on the handpiece, symbol: air. To blow a combination of air and water (spray), press both buttons on the handpiece at the same time.

After each use on a patient, the handpiece and tip of the syringe **MUST** be cleaned and sterilised to guarantee maximum hygiene.

Unscrew the nozzle retainer.


Withdraw the nozzle retainer.

Withdraw the nozzle.

You will find the sterilisation procedure


Press the pawl to release the handpiece.

Withdraw handpiece.


LED.light Operating Lamp

When swivelling and moving the operating lamp, always make sure that no objects or persons are in the swivelling range of the arm system. Otherwise, personal injury or property damage may occur. In order to move the operating lamp only touch its handles!

Switching the Lamp on and off

The operating lamp is switched on or off by moving the microswitch to the left / right. The intensity of the operating lamp can be adjusted from 8,000 to 50,000 LUX. Press and hold the switch to the left or right to adjust the light intensity. If the switch is pressed backwards or forwards, the operating lamp dims to 8,000 LUX. When the switch is pressed again, the light intensity set last is selected.

Mirror

Loosen the plastic cover at the front of the protective cap by lightly pressing on its upper end. On the back of the cover there is a mirror. Attach the mirror to the plastic cover.

Optical Properties:

Light field	: 170 mm x 85 mm
Lux	: 8,000 to 50,000 lux at a distance of 700 mm
Colour temperature	: 5,000 K
CRI (colour rendering index)	: > 85


You can find an application film at www.youtube.com/DKLGermany.
Video: DKL CHAIRS L2-D2 SERIES FUNCTIONS LED.LIGHT OPERATING LAMP

LINK: <https://youtu.be/OXSchFNIFXM>

Cleaning and Disinfection of the LED.light Operating Lamp

Pull the protective cap towards you to remove it for cleaning.

Cleaning and care


The reflectors must be cleaned with cotton wool and ethyl alcohol. Do not use detergents that contain surfactants or water-repellent substances (staining).

Sterilisation of the handles

To remove the handles, turn the handle lock on the handles and pull off the handles. To mount the handles, reattach them and push them to the limit. Then lock the handle lock.

The handles can be sterilized 200 times with standard cycles of 121°C/134°C.


Risk of breakage of the plastic parts.

Do not use any cleaning agents or disinfectants with the following ingredients for cleaning the plastic parts:

- AMMONIUM HYDROXIDE
- SODIUM HYDROXIDE
- METHYLENE CHLORIDE
- METHYL ALCOHOL.

Tray


The tray should be set correctly for your standard equipment and thus remain in its desired position. The maximum load of the tray is 2 kg. Do not lean on the tray.

Pull back the bellows until the adjusting screw (1) is visible. Place your standard equipment on the tray (max. 2 kg) and move the tray into a horizontal position.

Set the adjusting screw (1) so that the arm remains in the horizontal position (with a slight upward tendency).


Turn the adjusting screw in the „ + „ direction = high weight.
Turn the adjusting screw in the „ - „ direction = low weight

Secure the 3 grub screws (2) with safety lacquer (blue).

With the adjusting screw, the tray can be aligned in its radially horizontal position.

Pull back the bellows until the adjusting screw (1) is visible. Move the tray into a horizontal position. Secure the adjusting screw (1) for the inclination with thread locking fluid (blue).

You can find an application film at www.youtube.com/DKL_Germany.
Video: DKL CHAIRS L2-D2-RANGE TRAY SETTINGS

LINK: <https://youtu.be/r-Y14eEYtqI>

Cuspidor

①	Operating status with installed water separation unit
②	Bowl rinser
③	Glass filler

The factory setting for the bowl rinser is 7 seconds.
The glass filler is factory-set to 3 seconds.

Activate the bowl rinser at the cuspidor.

Activate the cup filler at the cuspidor.

Changing the bowl rinsing time.

Keep the bowl rinser key at the cuspidor pressed. After 2 seconds, a short signal is audible. Keep the key pressed until the desired rinsing time has been reached. Successful storage is confirmed with another signal tone. Maximum rinsing time: 25 seconds.

Changing the Filling Time for the Cup

Keep the cup filler key at the cuspidor pressed. After 2 seconds, a short signal is audible. Keep the key pressed until the desired fill level has been reached in the glass. Successful storage is confirmed with another signal tone. Maximum filling time: 10 seconds.

After switching on the treatment unit, the bowl rinsing process starts automatically for the saved duration and rinses the bowl.

You can find an application film at www.youtube.com/DKL_Germany.
Video: DKL CHAIRS L2-D2 SERIES FUNCTIONS CUSPIDOR CUP FILLER & BOWL RINSER

LINK: https://youtu.be/SsFSKDpJI_I

Cleaning the Sieve

Clean the sieve in the cuspidor bowl once a day under running water.

To reduce the risk of infection, liquid-tight gloves must be worn during maintenance work.

Never work without a filter, otherwise there is a risk that parts will settle in the suction system and thus impair its function.


Suction System

Reflux Effect

When the suction cannula adheres to the mucous membrane of the oral cavity or to the tongue, a so-called reflux effect may occur. To prevent cross-infection among patients we recommend using suction cannulas with air-bleed openings. Here a defined bypass airstream is introduced into the suction handpiece via lateral recesses. Even if the cannula adheres to the mucous membrane of the oral cavity or to the tongue and is thus blocked, a sufficient airstream from the patient to the suction system (and not the other way round!) is maintained.

Remove the suction tube from the suction tube retainer. By opening the slider or turning the control knob, the suction power is active.

Cleaning the suction filter

To reduce the risk of infection, liquid-tight gloves must be worn during maintenance work.

Open the lid of the filter drawer on the assistant's device to clean the disposable filter (daily) and replace it if damaged. The filter is designed as a disposable filter and cannot be thermally disinfected.

Never work without a filter; otherwise there is a risk of parts settling in the tube holder and impairing its function.

Suction System

For the intended use, please observe the operating instructions issued by DÜRR Dental:



- > Cuspidor valve, Version 3
- > CAS 1 operating instructions issued by DÜRR Dental
- > CS 1 operating instructions issued by DÜRR Dental, depending on the equipment and design of the suction system.



Open the door of the cuspidor fountain by pulling the handle. Depending on the equipment and design of the suction system, you will find the following in the cuspidor:

- > Cuspidor valve, Version 3 (wet suction)
- > CAS 1 operating instructions issued by DÜRR Dental (dry suction)
- > CS 1 operating instructions issued by DÜRR Dental (dry suction)

Cleaning and Disinfection of the Suction Handpieces


Cleaning and Disinfection of the Suction Handpieces

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For further information for manual and automatically preparation of materials refer to the hose manifold Comfort Installation and Operating Instructions, order no.: 9000-606-18. See also download pages at www.duerr.de

3

4

Cleaning and Disinfection of the Suction System

We recommend using the OroCup system of the DÜRR company for cleaning and disinfection.

Scope of delivery

1. OroCup, order no. 0780-350-00
2. 2 x Ø 16 mm insert (for large suction tube), (1x loose, 1x fixed)
3. 2 x Ø 6 mm insert (for the saliva ejector),
4. 1 x Ø 11 mm insert (other)

Orotol® plus

MD 555 cleaner

Consumables

- Orotol® plus suction unit disinfectant CDS110P6150 liquid concentrate
- MD 555 cleaner, special cleaner for suction units CCS555C6150 foam-free concentrate for dental suction units and discharge lines

Product Description

The OroCup care system is a closed dosing system for easy preparation and aspiration of disinfectants and special cleaning agents for suction system. With the OroCup, the suction systems can be equipped with all the components and the cuspidor can be cleaned and disinfected. The OroCup is suitable for standard suction tubes with different diameters. In the lid of the OroCup, there is one fixed insert for Ø 16 mm tubes. Two further inserts can be selected and used, depending on the diameter of the suction tubes. As needed, 1 - 3 suction tubes can be attached at the same time; unused connections have no influence on the function.

1. Selecting and Attaching Inserts

- Select and attach the insert according to the diameter of the suction tube (Ø 16 mm for the large suction tube, Ø 6 mm for the saliva ejector). Unused connections have no influence on the function.


2. Cleaning and Disinfecting the Suction System

Cleaning and disinfecting take place at the end of the treatment day; at higher utilisation levelstwice per day (e.g. at noon and in the evening or as needed).

- Wear personal protective equipment.

2.1 Pre-Cleaning with Water

- Aspirate 2 litres of water.

Cleaning and Disinfection of the Suction System

2.2 Preparation in the OroCup

- Depending on how much you need, prepare 1 or 2 litres of ready-to-use solution. Observe the manufacturer's instructions.

- Close the lid of the OroCup and shake the OroCup.


2.3 Positioning the OroCup and Aspirating

- Open the lid of the OroCup.
- Attach the disinfected suction handpieces and aspirate the ready-to-use solution for use.

Up to 2 litres of ready-to-use solution can be prepared in the OroCup.

Position the OroCup for aspirating the ready-to-use solution as follows:

- When aspirating the liquid, place the OroCup in a vertical position.

2.4 Cleaning and Disinfecting the Cuspidor Bowl

- Use at least 250 ml of ready-to-use solution per cuspidor bowl.

2.5 Final Rinsing

- After the exposure time has elapsed, aspirate 2 litres of water, see 2.1 (Pre-cleaning with water).

You can find an application film at www.youtube.com/DKL_Germany.

Video: DKL CHAIRS L2-D2 SERIES CLEANING AND DISINFECTION OF THE SUCTION SYSTEM

LINK: <https://youtu.be/39Lo60yeZnw>

Cleaning and Disinfection of the Surfaces

Hygiene and Care of the Stainless-Steel Surfaces

The regular cleaning of stainless-steel surfaces is recommended for hygienic as well as aesthetic reasons and serves to remove grease stains or finger marks. These can be easily removed with commercially available chlorine- and acid-free stainless-steel cleaners. We recommend applying Prestan to the surface in question.

Most stainless-steel care products contain silicone oil. Using these products can make your work a lot easier. They effortlessly remove any finger marks, but do not necessarily prevent new ones. Depending on the intensity of use, the protective layer remains in place for a few days. Microfibre cloths slightly moistened with water have also proved to be very effective.

Never use abrasive agents such as scouring powder, scouring milk or steel wool as these may cause scratches. Brushed surfaces must always be wiped in the direction of the finish. For this purpose, we recommend using a microfibre cloth. After cleaning, we recommend always wiping stainless-steel surfaces dry with a lint-free cloth to remove water stains or residual cleaning agent.

Disinfection of Stainless-Steel Surfaces

Do you put emphasis on a germ-free surface? Here, too, stainless steel proves to be extremely robust. Any commercially available chlorine-free disinfectant can be used.

Tests have shown that stainless steel is considerably easier to disinfect than other materials and even a lot less disinfectant is required in order to meet hygiene requirements.

The Most Important Facts at a Glance:

Effective and generally safe to use on surfaces are

- Soft sponges or microfibre cloths,
- Soapy water (to remove greasy stains),
- Diluted vinegar (to remove lime),
- Sodium bicarbonate (to remove coffee stains),
- Soda (to remove tea stains),
- Alcoholic solvents (to remove glue) and
- Special stainless-steel care products (for cleaning and conservation).

Caution is called for with

- Disinfectants containing chlorine and cleaning agents containing bleach (risk of corrosion).

Never use:

- Scrubbing sponges (scratches and extraneous rust),
- Scouring powder (scratches)
- Silver polish (corrosive).

Cleaning and Disinfection of the Surfaces

DentaClean: cleaning agent for imitation leather and plastic surfaces

Properties

DentaClean gently and easily cleans soiled imitation leather and plastic surfaces.

Use

Test on a hidden area first. We recommend cleaning the chair upholstery at the end of every treatment day. This is particularly important with light colours; any visible soiling must be removed immediately. Use the provided sponge to apply DentaClean in circular motions to the surfaces to be cleaned. Then remove moisture and dirt with a soft, absorbent cloth. For a thorough cleaning, use a cleaning brush instead of the sponge twice a week. Finally, seal with DentaProtect. Upholstery that is treated regularly and properly with DentaProtect is easier to clean!

DentaProtect: Care and protection for your imitation leather upholstery

Properties

DentaProtect is a product for caring and protecting heavy-duty imitation leather upholstery. Sealing the surface, it acts as a micro-binding agent and protects the upholstery from damage due to abrasion, soiling and also discolouration by non-fixed colourants in clothes.

Application

After a thorough cleaning in the evening with a soft sponge or a brush, apply the sealant to the dry upholstery. You need 1 – 2 wipes for the complete upholstery set of your treatment chair. Close the box immediately after taking out the wipes.

proPad: Ready-to-use alcohol-free wipes for rapid disinfection of alcohol-sensitive surfaces; particularly suited for use on treatment chair upholstery.

Area of Application According to the Biocides Regulation:

Ready-to-use wipes for disinfecting alcohol-sensitive surfaces.

Particularly soft, lint-free and tear-resistant fleece. 100% PES. High absorption of the disinfectant solution. Weight: 50 g/m². Large wipe: 28 x 30 cm. Without alcohol or perfume.

Application:

For disinfection, take the wipes out individually, wipe all over the surfaces and objects to be disinfected and allow the disinfectant to take effect. After the specified residence time, wipe the surfaces off with a clean cloth.

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Article number PP300

Cleaning and Care

Care set
1 x DentaClean 200 ml spray foam cleaner (article DC200)
1 x DentaProtect wipe dispenser box (article DP100)
1 x cleaning brush, 3 x cleaning sponges

Article number PSET

DentaClean 1000ml
Refill bottle for spray foam cleaner

Article number DC1000

DentaClean 200ml
Spray foam cleaner incl. 2 x cleaning sponges

Article number DC200

DentaProtect wipe dispenser box
100 wipes in a disposable sealing bag

Article number DP100

DentaProtect wipe dispenser set
6 wipe dispenser boxes with 100 wipes each in a disposable sealing bag

Article number DP600

Cleaning brush set, 4 pieces

Article number RB4

Cleaning sponge set, 8 pieces

Article number RS8

You can find an application film at www.youtube.com/DKLGermany.
Video: MEDICAL UPHOLSTERY – CLEANING DISINFECTION PROTECTION

LINK: <https://youtu.be/kUutWxmJI2E>

Cleaning and Disinfection of the Instrument Holders

- Pull out the instrument holder
- Rinse off any soiling with water
- Remove any residual liquid (absorbent cloth; blow dry with compressed air)
- Disinfection with disinfectants; wipe disinfection is recommended.
- Observe the disinfectant manufacturer's instructions for use
- After manual cleaning and disinfection, a steam sterilization (packaged) in a steam sterilizer class B or S (in accordance with EN 13000) is necessary.
- Remove any soiling on the instrument panel with a damp cloth.
- Remove any residual liquid (absorbent cloth; blow dry with compressed air)
- Disinfection of the instrument panel with disinfectants; wipe disinfection is recommended.
- Observe the disinfectant manufacturer's instructions for use
- Put the instrument holders back into their respective positions.

Marking	Type of holder
GS	Spray mist suction
KS	Saliva ejector
S	Function syringe


The supply tubes are not approved for mechanical cleaning (thermo washer disinfectant) and sterilization. Do not twist or fold the supply tubes! Do not roll up the tubes too tightly!

Bottle Care System

The Bottle Care System is a treatment water system for the self-sufficient water supply for all the instruments and the glass filler. There is also the possibility of intensive decontamination of the service water lines in the treatment unit.

Changing the treatment water bottle on the outside of the cuspidor:

1. Set the flip switch (1) on the bottle holder to „off“.
Turn the treatment water bottle (2) slightly to the left until the ventilation becomes audible.
2. Remove the treatment water bottle (2) from its holder from the left-hand side. Fill the bottle and then turn it clockwise into the bottle holder.
3. Set the flip switch (1) on the bottle holder to „on.“
The water supply is ready for operation.

Only use bottles approved by DKL.

Use the bottles before their expiration date (see bottle). Replace the bottles at the latest if they show visible damage or when they reach their expiration date – otherwise they might burst!

Empty the treatment water bottle at the end of a working day and refill the bottle at the beginning of a working day (after initial rinsing 120 sec. RKI) with fresh operating water and DK-DOX 150 Chlorine Dioxide ready-to-use solution (1 ampoule per liter).

You can find an application film at www.youtube.com/DKL_Germany.

Video: DK-DOX 150 READY-TO-USE CHLORINE DIOXIDE SOLUTION GLASS AMPOULES FOR THE BOTTLE CARE SYSTEM

LINK: https://youtu.be/Mj_y2YXAJKQ

Product number: 590013

<https://dkl.de/en/DK-DOX150-Clorine-dioxide-Ready-to-use-solution/590013>

To maintain the water quality in the Bottle Care System of DKL dental units.

CONTENT

- 30 x 5 ml glass ampoule of chlorine dioxide solution
- 1 x ampoule opener

APPLICATION

Open a 5 ml glass ampoule of chlorine dioxide solution using the ampoule opener and dispense it into the bottle of the Bottle Care System per 1 litre of water.

Bottle Care System

Labelling on the bottle

REF Article number

Expiration date year-month

Bottle Disinfection

For disinfection of the inside of the bottle at regular (weekly) intervals, we recommend BC-San 100. Further product information can be obtained from ALPRO Medical GmbH at www.alpro-medical.com.

Rinsing Function with the Bottle Care System

If the device is equipped with a bottle-care system, make sure that the bottle is freshly filled with water (see Bottle Care System).

Remove the first instrument. Activate the instrument with the setting “Spray active” for 120 seconds. For this purpose, hold the instrument over the cuspidor bowl or a sink. Repeat this procedure with all the instruments. Then activate the cup filler.

Then activate the cup filler.

Carry out the initial rinse before starting work without chlorine dioxide solution.

Intensive Flushing Function with the Bottle Care System

We recommend rehabilitating the waterways after longer periods of inactivity (holidays) or at least once a year.

You can find an application film at www.youtube.com DKL Germany.

Video: BOTTLE CARE SYSTEM WITH TOUCH SCREEN – DISINFECTION OF THE WATER SUPPLY

LINK: <https://youtu.be/1ySjnPDyuBQ>

Test water quality with Bottle Care System or Water Separation Unit (WTE)

You can find an application film at www.youtube.com DKL Germany.

Video: DKL CHAIRS L2-D2 TEST STRIPS FOR THE DK-DOX-150 DETECTION IN THE TREATMENT WATER

LINK: <https://youtu.be/3JUx-PY1xHo>

Product number: 590008 Test strips 0,1-0,4 ppm chlorine dioxide (50 pcs)
<https://dkl.de/en/Test-strips-0-1-0-4-ppm-chlorine-dioxide-50-pcs/590008>

Water Separation Unit (WSU)

The water separation unit meets the requirements of ISO 7494-2 and EN 1717 (free outlet with separation distance ≥ 20 mm). It is intrinsically safe in accordance with DVGW (German Technical and Scientific Association for Gas and Water) worksheet 540 and also meets the requirements of W270 and KTW (guideline for hygienic assessment of organic materials in contact with drinking water). It can be directly connected to the public drinking water supply. The water separation unit is a downstream dosing unit for the maintenance of the quality of the supplied water.

The LED button for the WSU is on the cuspidor. The WSU is started automatically by the control system of the dental unit. A manual start of the WSU is only necessary when the corresponding push-button signal is received.

The push-button on the WSU flashes green = start process with start filling. The push-button of the WSU is constantly green = normal operation.

After switching on the dental unit, the water separation unit goes into normal operation after approximately 50 seconds; water release is enabled.

Intensive flushing is activated by a double click on the push-button. The WSU button lights up permanently blue during the phase of double dosage. When the rinsing is activated, the WSU button flashes green-blue. When the rinsing process is finished, the button lights up permanently green. The return to normal operation is indicated on the display.

When the push-button of the WSU flashes yellow, DK-DOX 150 must be refilled. You can continue to work without restriction and interruption until the treatment unit is restarted.

The push-button of the WSU flashes red-blue after a restart of the treatment unit. DK-DOX 150 must be refilled! Caution! Continued operation without DK-DOX 150 only possible after confirmation by briefly pressing the button. If DK-DOX 150 is not refilled, the button of the WSU indicates this by flashing yellow after the start phase.

Error Messages and Service Mode

WARNING!

Push-button of the water separation unit (WSU) is flashing yellow-red. Fault in the dosing unit. WSU continues to work. Call customer service!

Push-button of the WSU is flashing red. Call customer service immediately!

Push-button of the WSU lights up permanently red. Immediately call customer service! Switch off the treatment unit at the main switch! The overflow sensor has tripped! There is a malfunction of the level sensor and / or the traveling valve.

Push-button of the WSU lights up permanently magenta. Immediately call customer service! Switch off the treatment unit at the main switch!

Push-button of the WSU flashes magenta. The WSU stops! The filling process of the supplied water is too slow. Reset the WSU by keeping the button pressed for 8 seconds and restart it by pressing the button briefly. If the magenta flashes continue, call customer service!

SERVICE MODE!

Keep the push-button of the WSU pressed for 8 seconds.

The push-button of the WSU lights up white.

The mixing tank, the double-piston pump and, if necessary, all the tubing of the dental unit are pumped empty. For this purpose, the instruments and the glass filler must be activated to get the water out of the water lines.

Completion of the emptying process is indicated by a white double flashing. After that, the WSU is in standby mode. The LED-ring of the push-button is switched off.

Flushing Function with the Water Separation Unit (WSU)

If the treatment unit is equipped with a water separation unit, make sure that the LED button on the cuspidor lights up green.


Remove the first instrument. Activate the instrument with the setting “Spray active” for 120 seconds. For this purpose, hold the instrument over the cuspidor bowl or a sink. Repeat this procedure with all the instruments. Then activate the cup filler.

Intensive Flushing with the Water Separation Unit (WSU)

After longer downtimes of the treatment unit, such as weekends or holidays, we recommend flushing the water lines intensively before starting work.

The push-button lights up blue: intensive flushing has started. Double click on the button of the cuspidor to activate intensive flushing. Intensive flushing is carried out in two phases. During the first phase with double dosage, the push-button of the WSU flashes blue. During the second phase, the rinsing phase, the push-button of the water separation unit flashes blue-green up to the end of the intensive flushing process.

Remove the first instrument. Activate the instrument with the setting “Spray active” for 600 seconds. For this purpose, hold the instrument over the cuspidor bowl. Repeat this procedure with all the instruments.

Once all the instruments have been flushed, activate the cup filler. Repeat this procedure until normal operation is indicated (the green push-button lights up) on the water separation unit.

After completion of the intensive flushing, the push-button of the cuspidor lights up green to indicate normal operation.

While intensive flushing is activated (push-button lights up blue) do not work on the patient.


Filling DK-DOX 150 into the Water Separation Unit

DK-DOX 150 is used to maintain the quality of the water supplied to DKL dental units. DK-DOX 150 is a pH-neutral, chlorine-free, aqueous chlorine dioxide solution and is used in the water separation unit with a mixing ratio of 0.75mg/l. The product is not classified as a hazardous substance.

DK-DOX 150 bottle Content 250 ml

Order number: 590007

Push-button flashes red/blue.
DK-DOX 150 must be refilled.

Open the cover of the storage container for the water separation unit at the base of the cuspidor.

You can find an application film at www.youtube.com/watch?v=CObils6gkIQ
DKL Germany. Video: DKL CHAIRS L2-D2 SERIES FILLING DK-DOX 150 WATER SEPARATION UNIT (WSU)

LINK: <https://youtu.be/CObils6gkIQ>

Do not fill any other agents into the storage container of the water separation unit. Only refill with DK-DOX150!

Open the cap of the storage container. Fill in 250 ml DK-DOX 150. Make sure that the glass bowl in the storage container does not overflow. After filling, put the cap back on the storage container and close the cover again. The water separation unit automatically returns to normal operation. The button lights up green.

Maintenance and Inspection

In order to ensure the operational and functional reliability of your treatment unit and to avoid damage due to wear and tear it is necessary to perform maintenance once a year. Maintenance is carried out by an authorised technician of your specialist dealer or a DKL CHAIRS technician.

The work steps to be performed and the parts to be replaced are specified in the document "Maintenance Log". The tasks that were performed have to be entered in the maintenance log, which is part of the medical devices logbook.

Safety Inspections

Dental units are designed in such a way that a first fault does not present a hazard to patients, operators or third parties. Therefore, it is important to detect such faults before a second fault occurs, which may result in a hazard.

For this reason, safety inspections should be carried out every 3 years to detect electrical faults in particular (e.g. faulty electrical insulation). These checks are carried out by an authorised technician of your specialist dealer or a DKL CHAIRS technician. The work steps to be performed are specified in the document "Safety Inspections". The measured values have to be documented.

Safety inspections have to be carried out when putting your treatment unit into initial operation, after expansion or retrofitting activities on your treatment unit and after repair jobs. The safety inspections are carried out in accordance with DIN EN 62353.

The treatment unit may only be operated when the safety checks have been passed.

Warranty Declaration

12 Months Warranty

This DKL medical product has been manufactured with the utmost care by highly qualified specialists. Multifarious checks and inspections ensure faultless performance. Please note that warranty claims will only be accepted if all the instructions in this operating manual have been observed.

DKL as the manufacturer shall be liable for material and manufacturing faults within a warranty period of 12 months from the date of purchase. Accessories and consumables (seals, filters, lamps and suction tubes) are excluded from this warranty. We do not accept liability for damages caused by improper treatment or repair work carried out by third parties that are not authorised by DKL!

Any warranty claims must be filed with the supplier or an authorised DKL service partner and the sales slip must be enclosed. Any performance of this warranty does not extend the warranty period.

Waste Disposal

Make sure that the parts that are being disposed of are not contaminated.

Observe your local and national laws, guidelines, standards and regulations for disposal.

> Medical devices

> Waste electrical and electronic equipment

Further information on disposal can be found at <http://dkl.de>

Disposal and recycling of DKL transport packaging is carried out within the scope of the Dual System via the local waste disposal and recycling companies.

DKL transport packaging returned by customers at their own expense is supplied by DKL to the recycling companies set up for this purpose without further costs and without reimbursement.

EMC - Manufacturer's Declaration for the Model L2-SUC

- **WARNING:** The use of accessories that do not conform to the manufacturer's specifications may result in higher interference levels and/or lower interference immunity.
- Operate the equipment in a location as far away as possible from equipment that emits electrical and magnetic disturbances. If it is necessary to operate the device in the immediate vicinity of other devices, make sure that the system functions correctly.

BASIC SAFETY

BASIC SAFETY is ensured if it meets the safety requirements of the IEC 60601-1 standard, in particular the requirements against: electrical shock, mechanical hazards and hazards due to excessive temperatures.

ESSENTIAL PERFORMANCE

The dental unit has no direct clinical function or essential performance according to IEC 60601-1, IEC 80601-2-60, 201.4.3 ESSENTIAL PERFORMANCE.

Performance limitations are permitted according to the following criteria. This is considered in the risk analysis of the system.

Criterion A

The dental unit will withstand the test without damage or other interference. During and after the test, the device will operate perfectly within the specified limits. Basic safety is guaranteed throughout.

Criterion B

The dental unit will withstand the test without damage or other interference. After the test, the device will operate perfectly within the specified limits. Basic safety is guaranteed throughout.

Criterion C

A temporary malfunction is permitted if the function resets itself or if it can be restored by user intervention. Basic safety is guaranteed throughout.

Intended operating environment

Intended operating environments are typical professional health care facilities and areas of home health care.

Technical description

This dental unit has been tested and developed to meet the EMC behaviour in the specified environment. This includes special EMC-filters to reduce the radiation of electromagnetic waves as defined in IEC 60601-1-2.

Please read and follow all technical documentation to avoid adverse events for the patient or user.

IEC STANDARD 60601-1-2:2014, 4th Edition

This device is approved for use in a specific electromagnetic environment. The customer or user of the device must ensure that it is used in an electromagnetic environment in accordance with the description given below.

Emission Measurement	Agreement	Guidelines Regarding the Electromagnetic Environment
RF-emission according to CISPR 11	Group 1	This device uses RF-energy for internal functions only. RF-emissions are therefore very low, and it is unlikely that other nearby electronic equipment will be disturbed.
RF-emission according to CISPR 11	Class B	The device is suitable for use in all environments, including residential areas, and approved for direct connection to the public low-voltage network for residential areas.
Harmonics according to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker according to IEC 61000-3-3	met	

Interference Immunity Test	IEC 60601-test level	Compliance level	Electromagnetic Environment - Guidelinesf
Electrostatic Discharge (ESD) according to IEC 61000-4-2	± 8 kV contact discharge ±2, 4, 8, 15 kV air discharge	± 8 kV contact discharge ±15 kV air discharge	The floor should be wood, concrete or tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%. Criterion B
Fast transient electrical disturbances/bursts according to IEC 61000-4-4 (only for V 300/600)	± 2 kV for mains 100 kHz repeat rate	± 2 kV for mains 100 kHz repeat rate	The quality of the mains power supply should meet the requirements for a normal commercial or clinical environment. Criterion B
Surge voltages according to IEC 61000-4-5 (only for V 300/600)	± 0,5 kV , ± 1 kV L to N ± 0,5 kV , ± 1 kV ± 2 kV L to GND	± 1 kV L to N ± 2 kV L to GND	The quality of the mains power supply should meet the requirements for a normal commercial or clinical environment. Criterion B
Voltage dips, short-term interruptions and voltage fluctuations of the mains supply lines according to IEC 61000-4-11 (only for V 300/600)	0 % UT 0°,45°,90°,135°,180°,225°,270°,315° 0 % UT 0° 0% 70 % UT 0 % UT 0%	0 % UT for 1/2 Period 1 Period 25 /30 Periods (50/60Hz) 250/300 Periods (50/60Hz) for 5 s	The quality of the mains power supply should meet the requirements for a normal commercial or clinical environment. If the user of the product requires continuous operation even with interruptions of the power supply, the product should be connected to an uninterruptible power supply. Criterion A (max. mains voltage) Criterion B (min. mains voltage) Criterion A (max. mains voltage) Criterion B (min. mains voltage) Criterion A (max. mains voltage) Criterion B (min. mains voltage) Criterion A (max. mains voltage) Criterion B (min. mains voltage)
Magnetic field at the mains frequency (50/60 Hz) according to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields at the mains frequency should have levels typical of an application in a commercial or clinical environment.

Note: UT is the alternating mains voltage prior to the application of the test level.

Specifications for Enclosure Port Immunity

Immunity Test	Test condition	IEC 60601 level of conformity	Electromagnetic Environmental Recommendation
Radiated electromagnetic fields from high-frequency wireless communication devices IEC 61000-4-30-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM 1kHz	10 V/m 80 MHz – 2,7 GHz	The quality of the main power supply should correspond to the one for a professional health care facility and be appropriate regarding the environment in areas of domestic health care. Criterion A
	385MHz (18Hz pulse modulation)	27 V/m	
	450MHz (FM+/-5KHz deviation 1kHz sine or 18Hz pulse modulation)	28 V/m	
	710MHz (217Hz PM)	9 V/m	
	745MHz (217Hz PM)	9 V/m	
	780MHz (217Hz PM)	28 V/m	
	810MHz (18Hz PM)	28 V/m	
	870MHz (18Hz PM)	28 V/m	
	930MHz (18Hz PM)	28 V/m	
	1720MHz (217Hz PM)	28 V/m	
	1845MHz (217Hz PM)	28 V/m	
	1970MHz (217Hz PM)	28 V/m	
	2450MHz (217Hz PM)	28 V/m	
	5240MHz (217Hz PM)	9 V/m	
	5500MHz (217Hz PM)	9 V/m	
5785MHz (217Hz PM)	9 V/m		

Interference immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment - guidelines
Conducted RF-disturbance variables according to IEC 61000-4-6 (only for V 300/600) Radiated RF-disturbance variables and near fields from wireless communication equipment according to IEC 61000-4-3	3 Veff 150 kHz to 80 MHz 10 V/m 80 MHz to 2,7 GHz	10 Veff 10 V/m	The distance between portable or mobile RF-communications equipment and parts of the product, including cables, should not be less than the recommended protective distance calculated with the equation applicable to the transmission frequency. Recommended protective distance: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ for 80 MHz to 800MHz $d = 2,3\sqrt{P}$ for 800 MHz to 2.7 GHz Here P is the maximum nominal power output of the transmitter in watts (W) according to the transmitter manufacturer's specifications, and ^d is the recommended distance in metres (m). The field strength of permanently installed RF-transmitters, which was determined by an electromagnetic location test ^a should not exceed the level ^b permitted in any frequency range. Interference may occur in the immediate vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz respectively, the larger frequency range applies.

Note 2: These guidelines may not apply to all situations. The propagation of electromagnetic waves is affected by the absorption and reflection of structures, objects, people and animals.

^a The field strength of permanently installed transmitters, such as base stations for radio telephony (cordless or mobile phones), mobile radio stations, amateur radio transmitters, AM and FM radio and television transmitters, cannot theoretically be calculated with absolute accuracy. To determine the electromagnetic fields generated by fixed RF-transmitters, an electromagnetic site inspection should be carried out. If the measured field strength at the location where the device is used exceeds the permissible RF-field strength specified above, the instrument should be observed. Additional measures may be necessary, e.g. reorientation or change of location of the device.

^b In the frequency range between 150 kHz and 80 MHz, the field strength should be less than 3 V/m.

Manufacturer's Declaration - Electromagnetic Interference Immunity III

The device is approved for use in a specific electromagnetic environment.

The customer or user of the device must ensure that it is used in an electromagnetic environment as described below

Interference Immunity Test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidelines
Fluctuations in the mains frequency and mains voltage according to IEC 601-1, section 10.2.2. a (only for V 300/600)	Nominal frequency: up to 100 Hz: variations of ± 1 Hz of the nominal frequency; variations of $\pm 10\%$ Hz of the nominal voltage	Nennfrequenz: bis zu 100 Hz: Schwankungen von ± 1 Hz der Nennfrequenz; Schwankungen von $\pm 10\%$ Hz der Nennspannung	The quality of the mains voltage supply should meet the requirements of a normal commercial or clinical environment.

Manufacturer's Declaration - Recommended Protective Distances between Portable or Mobile RF-Communication Equipment and the Device

The device is intended for use in an electromagnetic environment where the radiated RF-disturbance variables are checked. The customer or user of the device can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF-communication equipment (transmitters) and the device in accordance with the following recommendations, which are based on the maximum output power and frequency of the communication device.

Maximum nominal power of the transmitter in watts (W)	Protective distance as a function of the frequency of the transmitter in metres (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum output power not specified above, the recommended safety distance d in meters (m) can be calculated with an equation from the transmitter frequency and the maximum nominal output power P of the transmitter in watts (W) based on the transmitter manufacturer's specifications.

Note 1: At 80 MHz and 800 MHz respectively, the larger frequency range applies.

Note 2: These guidelines may not apply to all situations. The propagation of electromagnetic waves is affected by the absorption and reflection of structures, objects, people and animals.

ATTENTION: The use of this device directly adjacent to or coupled to another unit should be avoided as it may lead to unintentional behaviour. However, if this arrangement is unavoidable, both devices must be observed to verify that they are functioning normally.

CAUTION: Portable RF-communication equipment (including antenna cables or external antennas) should not be closer than 30 cm to the ME-equipment or ME-system, including those cables specified by the manufacturer. Otherwise, a power limitation of the device could be caused.

Minilight

Installation and use manual

LUZZANI DENTAL S.R.L.

Via Torino, 3
20030 Senago (MI)
Italy

www.luzzani.it

E-mail: info@luzzani.it

Tel.: +39 029988433
Fax: +39 0299010379

ED. 5 REV. 0

MADE IN ITALY



1 3 7 0

CONTENTS

0 — SYRINGE LEGEND	8 — CLEANING AND STERILIZATION
1 — WARNINGS	9 — MAINTENANCE
2 — PRODUCT DESCRIPTION	10 — DISPOSAL AND SCRAPPING
2.1 — GENERAL	11 — INFORMATION FOR THE DENTIST
2.2 — GENERAL CHARACTERISTICS	12 — FAULT REPORT FORM
2.3 — MODELS	13 — SYMBOLS
2.4 — CONTROLS	14 — WIRING DIAGRAM
2.5 — CE MARKING	
3 — IDENTIFICATION DATA AND WARRANTY	
3.1 — MANUFACTURING BATCH	
3.2 — WARRANTY	
4 — PACKAGING	
5 — TECHNICAL CHARACTERISTICS	
5.1 — OPERATING CONDITIONS	
5.2 — TRANSPORT AND STORAGE CONDITIONS	
6 — INSTALLATION AND CONNECTIONS	
6.1 — CONNECTION TO ELECTRICAL SYSTEM	
6.2 — CONNECTION TO HYDRAULIC SYSTEM	
6.3 — CONNECTION TO COMPRESSED AIR SYSTEM	
6.4 — CONNECTION OF STRESS ABSORBER CORD	
6.5 — NOTES FOR CORRECT CONNECTION	
7 — NORMAL USE	
7.1 — INSUFFLATION OF COLD WATER	
7.2 — INSUFFLATION OF COLD AIR	
7.3 — COMBINED INSUFFLATION OF COLD WATER AND AIR (SPRAY)	
7.4 — INSUFFLATION OF WARM WATER (present in version 6f and L)	
7.5 — INSUFFLATION OF WARM AIR (present in version 6f and L)	
7.6 — COMBINED INSUFFLATION OF WARM WATER AND AIR (SPRAY) (present in version 6f and L)	
7.7 — LIGHT FUNCTION	
7.8 — FIRST TIME USE AND USE AFTER LONG PERIODS OF NON-USE	

0 — SYRINGE LEGEND

1 — WARNINGS

1.1 — Any unauthorised tampering, modification or improper use immediately terminates the warranty and exonerates our company from all liability for injury or damage to persons, animals or property that may be caused by such interference.

1.2 — To ensure maximum user and patient safety, the transformer used must be rated Safe Extra Low Voltage (SELV) with CE marking for medical use, in compliance with IEC 80.601 and IEC 60.601.

1.3 — Connection to terminals with output voltage exceeding 24 V could irreparably damage the device.

1.4 — The stress absorber cord must be anchored to the dental unit. This cord is designed to absorb any strains, thus preventing abnormal traction on the electrical or air/water line connections.

1.5 — Before use, the electric, water and air lines must be correctly connected. The syringe could be damaged if used when not connected to the water and air supplies. The lines must be connected correctly since inverting them would damage the syringe.

1.6 — To meet the requirements of Directive 93/42/EEC, the company has established a post-marketing surveillance procedure to monitor any problems generated by the use of our products. The attached form allows you to report any faults and suggest improvements which will be considered for subsequent versions of the product.

1.7 — With each syringe, the package also includes a User's manual which also includes a "Fault Report Form". Since this is required by law, the user must receive this User's manual. Therefore, the syringe installer is responsible for delivering this User's manual to the dentist. Directive 93/42/EEC requires product traceability: therefore, our customers are asked to ensure that, in case of emergency, we can identify the end customer to whom the product has been sold.

1.8 — Manufacturers and/or installers of dental units are required to comply with all the prescriptions outlined in this document.

1.9 — Use the Minilight syringe only for the applications described in the instructions for use.

1.10 — This product must only be installed by qualified persons.

1.11 — Never modify the syringe in any way. This is strictly forbidden.

1.12 — Use only original parts produced or approved by the manufacturer. If non-OEM accessories or consumables are used, the company cannot guarantee safe operation and function. No claims can be made for damages resulting from improper use.

1.13 — Disposable tip adapters are not part of the Minilight Syringe and may compromise proper functioning.

1.14 — Do not use the device in close contact with anaesthetic gases or in highly oxygenated environments (with an oxygen content >25%) or in areas where there is a risk of explosion.

1.15 — Do not perform any maintenance procedures not indicated in the manual.

1.16 — Before using the syringe, make certain that the water and air supplies have been activated.

1.17 — National regulations regarding dental unit water and air quality must be met.

1.18 — The air used must be dry, clean and free of oil.

1.19 — The Minilight syringe meets the requirements laid out in European Council Directive 93/42/EEC: Class II A

1.20 — The Minilight syringe must not be used near or set on other devices.

1.21 — Use only accessories, cables, transducers specified or supplied by Luzzani Dental.

1.22 — Do not use mobile RF communication devices within a distance less than 30 cm from the device.

2 — PRODUCT DESCRIPTION

2.1 — GENERAL

The Minilight syringe is a medical device designed to blow air and water (separately or together, at room temperature or at body temperature) to clean and/or dry the oral cavity during any dental procedure. It has been designed for use in dentist offices and dental clinics and is built into dental unit used exclusively by dentists. Product life — under proper maintenance conditions — is 5 years.

2.2 — GENERAL CHARACTERISTICS

- The Minilight syringe is a medical device for dental use (class II a).
- Protection rating against direct contacts: B
- Temporary operation: 10 sec. ON/ 20 sec. OFF.

The Minilight syringe has been designed using the latest ergonomic concepts for easy use and immediate cleaning and sterilisation. Both the tip and handpiece can easily be removed for perfect autoclave disinfection and sterilisation at 134°C (see point 8). Several handpieces of different shape are available. Choose the shape you need: angled or stylet. The air and water can also be heated to body temperature, thus sparing the patient even the slightest discomfort induced by insufflation of products at ambient temperature.

2.3 — MODELS

The models are differentiated into versions based on the

number of functions provided:

- 3F: cold water/air/spray
- 6F: cold water/air/spray - warm water/air/spray
- G3F: cold water/air/spray (air and water inverted)
- G6F: warm water/air/spray (air and water inverted)
- A: cold water or air only
- L: cold/warm water, air and spray + light

Note:

Models 3F, G3F, 6F, G6F and L can be recognized by the printing on the back of the syringe.

Various handpieces, that differ in shape, can be mounted on all versions of the Minilight syringe.

The shape of the handpiece can be:

STYLET

ANGLED

This difference lets the dentist choose the tool ergonomically best suited for the purpose. To guarantee maximum hygiene and atoxicity, the handpieces are made of stainless steel. The devices are produced entirely in our workshop, with a tested, constantly updated work cycle using the most sophisticated machinery compliant with current quality system directives (UNI EN ISO 13485 certified).

2.4 — CONTROLS

Every syringe and all of its parts undergo duly documented, 100% complete functional and safety testing to ensure that the technical and functional design requirements are fully met.

2.5 — CE MARKING

All products bear CE marking both on syringe handpiece and inside. (batch number, autoclave symbol, Luzzani Dental logo, product name, CE marking with Notified Body number). The User's manual supplied with the product also includes details of our company, the main product characteristics and instructions for correct use and maintenance.

3 — IDENTIFICATION DATA AND WARRANTY

3.1 — MANUFACTURING BATCH

A number, marked on the inside of each product, identifies the production batch; the number is printed on the central body of the syringe. This number uniquely identifies the production batch thus always guaranteeing traceability of the product and each of its components, with relative test sheets.

3.2 — WARRANTY

The product is guaranteed by our company for 12 months from the date of the delivery document. The warranty covers any device manufacturing defects (materials) and is limited solely to the replacement of defective parts, performed in our workshop. The product must be sent to our premises at the expense of the customer. For the warranty to be valid, the product must be returned intact, complete and showing with no signs of tampering. The syringe has no functional expiration date; its expected life span is 5 years.

4 — PACKAGING

The product is shipped in suitable packaging to prevent problems during transport. The packaging consists of a plastic bag containing the Minilight syringe. A protective handpiece sheath and tip-saving tube protect the syringe during transport. Several bags are placed in one box. The Minilight syringe comes ready for connection to the dental unit once all packaging has been removed. For the correct use of the syringe, the protective tube must also be removed from the tip.

 **IMPORTANT NOTE:**

With each syringe, the package also includes a User's manual which also includes a "Fault Report Form". Since this is required by law, the user must receive this User's manual. Therefore, the syringe installer is responsible for delivering these forms to the dentist. Directive 93/42/EEC requires product traceability: therefore, our customers are asked to ensure that we can identify the end customer to whom the product has been sold.

5 — TECHNICAL CHARACTERISTICS

FUNCTIONS		6F	3F	L
Supply voltage	V~	24	***	24
LED voltage	Vdc	***	***	3,5
Max. absorbed current	A	4,3	***	4,3
Electrical power	W	103	***	103
Water supply pressure	Kpa	250	250	250
Air supply pressure	Kpa	450	450	450
Water flow rate	Cc/min	110	110	110
Air flow rate	NI/min	10	10	10
International protection		IP40	IP40	IP40

5.1 — OPERATING CONDITIONS

Ambient temperature	10°C / +45°C
Relative air humidity	30% / 85%
Atmospheric pressure	80 Kpa - 106 Kpa

5.2 — TRANSPORT AND STORAGE CONDITIONS

Temperature	-20°C / +60°C
Relative Humidity	30% / 85%
Atmospheric Pressure	50 Kpa - 106 Kpa

6 — INSTALLATION AND CONNECTIONS

6.1 — CONNECTION TO ELECTRICAL SYSTEM

The connection involves hooking up the two electrical terminals to the transformer, arranged to provide an output of 24 V (see wiring diagram 14.1). Model L only: voltage for LED light 3.3 Vdc: black wire (+), white wire (-)

 **NOTE**

- The electric power supply and transformer used must be rated Safe Extra Low Voltage (SELV) for medical use, in compliance with the mandatory regulations.
- Connection to terminals with output voltage exceeding 24 V could irreparably damage the unit and compromise safety.

6.2 — CONNECTION TO HYDRAULIC SYSTEM

The syringe's green line must be hooked up to the water supply.

 **NOTE**

- The operating pressure is 250 kPa. A slight increase in pressure would only increase the power of the jet, but would not create any problems and is not dangerous.
- The water used must be potable water, filtered (<25 µm) and free of bacteria, etc.
- For the syringe to function properly, the water pressure must not be lower than indicated.

6.2 — CONNECTION TO COMPRESSED AIR SYSTEM

The syringe's blue line must be connected to the compressed air system.

 **NOTE**

- The recommended operating pressure is around 450 kPa.
- When using the Minilight syringe, national regulations regarding water and air quality must also be met.
- The air must be medical grade, dry and free of oil and bacteria – a 5µm air filter is recommended.

6.4 — CONNECTION OF STRESS ABSORBER CORD

The stress absorber cord must be anchored to the dental

unit. This cord is designed to absorb any strains, thus preventing abnormal traction on the electrical or air/water line connections. The manufacturer cannot be held liable for malfunctions caused by failure to anchor the stress absorber cord.


6.5 — NOTES FOR CORRECT CONNECTION

- Before carrying out functional tests, both the electric, water and air lines must be correctly connected.
- Using the syringe without connecting the water and air supplies could damage the syringe.
- The lines must be connected carefully since inverting them could damage the syringe.


7 — NORMAL USE


7.1 — INSUFFLATION OF COLD WATER

To blow cold water into the operating field, just press the left button on the handpiece, symbol: 


7.2 — INSUFFLATION OF COLD AIR

To insufflate cold air into the operating field, just press the right button on the handpiece, symbol: 

7.3 — COMBINED INSUFFLATION OF COLD WATER AND AIR (SPRAY)

To blow a combination of cold air and water (spray), press both buttons on the handpiece at the same time:


7.4 — INSUFFLATION OF WARM WATER (present in version 6f and L)

To blow warm water into the operating field, turn the selector switch at the base of the handpiece to the right (the green

LED lights up)

and press the left button on the handpiece:


7.5 — INSUFFLATION OF WARM AIR (present in version 6f and L)

To blow warm air into the operating field, turn the selector switch at the base of the handpiece to the right (the green LED lights up):

and press the right button on the handpiece:


7.6 — COMBINED INSUFFLATION OF WARM WATER AND AIR (SPRAY) (present in version 6f and L)

To blow a combination of warm air and water (spray), turn the selector switch at the base of the handpiece to the right (the green LED lights up)

and press both buttons on the handpiece at the same time:

7.7 — LIGHT FUNCTION (present in version L)

This function is present in the Minilight L version. In this model, the handpiece has an optical fibre that conducts the light generated by a LED located in the body of the syringe. We recommend using a switch to turn the LED on and off. For connections, see wiring diagram 14.2.

NOTE

- The sole function of the selector switch is to preselect warm operations. The water and/or air are only heated at the moment in which they are actually used.
- The cleaning (or line washing) procedures must always be performed in the cold position.

WARNING

Do not use the tip improperly. Remove and sterilise the tip after each patient.

IMPORTANT

Air and water must be able to flow freely from the tip. Do not rest the tip on the tooth or on an object. Do not press the tip against impression materials as they could cause obstruction

7.8 — FIRST TIME USE AND USE AFTER LONG INTERVALS

-  Sterilize the handpiece and all accessories before use.
- After prolonged periods of inactivity, clean, treat and sterilise the handpiece.

BEFORE EACH PATIENT

1. Make certain the handpiece has been sterilised.
2. Adjust the supply of fluids from the dental unit (see table in point 5).
3. Press the air button and make certain that there is a clearly perceptible jet of air.
4. Check the water flow rate.
5. Use only filtered water that is free of oil and microorganisms.

6. Check the tip for any obstructions or deposits. Clean if necessary.

NOTE

- Flush out the syringe at the beginning of each work day (minimum flushing time: 2 minutes) and before each patient (minimum flushing time: 20-30 sec.).
- Immediately upstream of the syringe, install filters able to retain the microorganisms coming from the hydro-pneumatic circuit.

8 — CLEANING AND STERILIZATION


-  After each use on a patient, the handpiece and tip of the syringe **MUST** be cleaned and sterilised to guarantee maximum hygiene.

Sterilisation symbol on the handpiece: 

To do this, proceed as follows:

Disconnect the tip by unscrewing the tip retainer


1. unscrew the tip retainer

2. withdraw the tip retainer


3. withdraw the tip

or remove the entire handpiece by pressing the button on the lower part of the handpiece and pulling upwards.

1. press the pawl to release the handpiece

2. withdraw handpiece

wipe with a damp cloth, removing any stains. Set in a steam autoclave at 134°C for AT LEAST 3 minutes (in compliance with CEI EN 13060).

A — WARNINGS

The syringe should always be sterilised, even before using it for the first time.

- Inappropriate sterilisation is hazardous for the patient and for the operators.
- Do not perform spray disinfection.
- Do not immerse in disinfectant liquids.
- Do not perform cold or hot air chemical sterilisation.
- The personnel performing the task must be skilled and specially trained.
- Use disinfectant according to the specifications on the manufacturer label.
- Do not use chlorine-based liquids.
- When simultaneously sterilising more than one item in an autoclave, check that the load does not exceed the maximum allowed.

B — PREPARATION

Eliminate surface dirt using a disposable paper towel. Clean the inside of the lines by running air and water through the syringe for about 30 seconds. Remove the stainless-steel handpiece by pressing the button on its terminal section. Unscrew the tip retainer and remove the tip.

C — MANUAL CLEANING

Use a disposable paper towel and potable water to remove any impurities or dirt that may be present

D — AUTOMATIC CLEANING

Not envisaged

E — MANUAL DISINFECTION

Disinfect only with a disposable cloth and the permitted disinfectant (following the instructions on the label and product technical data sheet).

Recommended disinfectants:

- Incidin liquid
- FD 322 Durr
- Mikrozid AF liquid

F — AUTOMATIC DISINFECTION

Not envisaged

G — MANUAL DRYING

Dry with disposable paper towelling. Dry with clean, dry, uncontaminated compressed air, inside and out, continuing until completely dry. Do not dry with hot air.

H — AUTOMATIC DRYING

Not envisaged

I — MAINTENANCE AND CONTROL

No special maintenance is necessary. There is no objective period of time that limits the useful life of the handpiece: visually check for damage and signs of wear, and if found, replace the part.

L — PACKAGING

Use heat sealable film-paper sterilisation pouches of appropriate size.

M — STERILISATION

The handpiece and tip can be sterilised.

Sterilise in a class B steam autoclave in compliance with EN 13060 ISO 17665-1.

3-phase sterilisation with fractional vacuum system at 134°C +/- 1°C at a pressure of 2.13 bar, applying a 4-minute delay. Never exceed 134°C. The autoclave must be validated.

N — STORAGE

No particular requirements apart from storage in the sealed, sterilised pouches. Store in a suitable place that is dry, out of direct sunlight and possibly with low bioburden.

9 — MAINTENANCE

The instrument requires no specific maintenance apart from normal cleaning and sterilisation as described in the previous paragraph.

10 — DISPOSAL AND SCRAPPING

-  The product does not contain dangerous or toxic-hazardous components. Separate waste collection is required for electrical equipment. Follow the regulations in force in your country.

11 — INFORMATION FOR THE DENTIST


-  The dental unit manufacturer is required to deliver the Luzzani Dental syringe User's manual to the end user.

12 — FAULT REPORT FORM

To meet the requirements of Directive 93/42/CEE as amended, the company has established a post-marketing surveillance procedure to monitor any problems generated by the use of our products. This commitment includes the requirement that both user and manufacturer inform the competent authorities of any incident caused to patient or user by malfunction or deterioration in the characteristics and/or performance of a device, as well as any inadequacy in the instructions for use. We kindly ask you to inform us of any anomalies by sending us the sheet attached to the last page of this manual.


13 — SYMBOLS

SYMBOLS


- | | | | |
|---|---|---|-----------------------|
|  | Do not overturn |  | Batch number |
|  | Fragile |  | Sterilise |
|  | Keep dry |  | Consult User's manual |
|  | Type B device | | |
|  | Alternating current | | |
|  | General warnings | | |
|  | Separate collection for electrical and electronic equipment | | |
|  | Manufacturer | | |
|  | Double insulation | | |

14 — WIRING DIAGRAM


14.1 – VERSION 6F

14.2 – VERSION L

INSTRUCTIONS FOR CLEANING AND STERILIZATION OF MINILIGHT, MINIMATE, MINIBRIGHT SYRINGES IN ACCORDANCE WITH UNI EN ISO 17665 REQUIREMENTS

<p>Warning</p> 	<ul style="list-style-type: none"> • Sterilization must be performed even when using the syringe for the first time. • Inappropriate sterilization is dangerous for patients and operators. • Do not perform spray sterilization. Do not place in any disinfectant liquids. • Do not perform cold chemical or hot air sterilization. • The appointed staff must be specialized and trained. • Use the disinfectant in accordance with the specifications set by the manufacturer indicated on the label. • Do not use chlorine-based solutions. • When sterilizing more than one piece in one autoclave do not exceed its maximum load capacity.
<p>Preparation</p>	<p>Remove dirt from the surface using a disposable paper tissue. Let air and water flow from the syringe for about 30 seconds in order to clean the internal channels. Press the button located at the bottom of the sleeve and slide the stainless steel sleeve off the syringe body. Unscrew the ferrule and remove the tip.</p>
<p>Manual cleansing</p>	<p>Wipe with a disposable paper tissue and with the aid of drinking water to remove any impurities and dirt</p>
<p>Automatic cleansing</p> 	<p>Not available</p>
<p>Manual disinfection</p>	<p>Perform disinfection only with a disposable tissue and with compatible disinfectants (in accordance with the instructions included in the product label and technical data sheet). Recommended disinfectants:</p> <ul style="list-style-type: none"> • Incidin liquid • FD 322 Durr • Mikrozyd AF Liquid
<p>Automatic disinfection</p> 	<p>Not available</p>
<p>Manual drying</p>	<p>Dry with disposable paper tissue. Dry with dry, clean and uncontaminated compressed air both internally and externally until completely dry. Do not dry with hot air.</p>
<p>Automatic drying</p>	<p>Not available</p>
<p>Maintenance and checking</p>	<p>No particular maintenance is required. There is no objective term limiting the usage life of the sleeve: check to see if there are if any damages or signs of wear and tear, replace the part if necessary.</p>
<p>Packaging</p>	<p>Use appropriately sized sterilization packages made of thermoweldable film.</p>
<p>Sterilization</p> 	<p>The sleeve and tip are autoclavable. Sterilize using EN 13060 ISO 17665-1 compliant class-B steam autoclave. Sterilize with 3 phases fractionated vacuum 134° C +/- 1 °C with 2.13 bar pressure, 4 minutes wait. Never exceed 135° C. The autoclave must be validated.</p>
<p>Preservation</p>	<p>No particular requirements other than keeping them in their sealed and sterilized package. Keep them in an appropriate environment and out of direct sunlight and in a dry place, which should have low bioburden where possible.</p>

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*It complies with the following European Directives and their implementations into national laws
 and subsequent modifications.*

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135°C
}}}

SAFETY RULES**GB**

The SYR3 dental syringe is designed exclusively for use in dentistry, by dental professionals or an assistant under a dentist's supervision, in order to perform specific cleaning and drying operations on the operating area. The device must be installed on a specific supply system (dental unit) and be connected to a dedicated FARO hose. The device must be sterilised in an autoclave before it is used on another patient. The device is not supplied sterile and must therefore be sterilised prior to use. Perform only the operations contained in the following manual. In all other cases request technical assistance.

HOSE AND SYRINGE INSTALLATION

Installation must be performed by specialised personnel only.

Check that the package contains the following components:


- Dental syringe.
- Supply hose.
- O-rings.
- Probe.
- Instruction manual.

Connect the BLUE tube to the AIR intake, the GREEN tube to the H₂O intake and fasten the safety cord.


Connect the hose to the syringe by tightening the ring nut. Check the water and air pressure values and make sure the fluids are at ambient temperature (T_{max} < 41°C).

INSTRUCTIONS FOR USE

- For water, press the green button.
- For air, press the blue button.
- For spray, press both buttons simultaneously.

Disassembling sheath and nozzle

- To remove the sheath, press button "A" and pull off in the direction of the arrow.

- To remove the nozzle, unscrew ring "B" and pull off in the direction of the arrow.

- Before reassembling the syringe, grease the two seals on the inner body and the metal tube of the nozzle with petroleum jelly.

Cleaning the device

Clean the sheath and the nozzle under running water or use alcohol-based solutions.

Sterilisation

- Pull the sheath, complete with the nozzle, from the body of the syringe.
- Clean sheath and nozzle.
- Sterilise the sheath complete with the nozzle in an autoclave. The sheath has a valve that prevents contamination of the body of the syringe.

Cleaning nozzle holes

- If the tiny holes in the nozzle that the water and air come out of are plugged, free the by using the specially provided probe.

TECHNICAL SPECIFICATIONS**Three-function syringe SYR3**

- Air pressure : **4 bar**
- Water pressure : **2 bar**
- Air flow : **10 l/min**
- Water flow : **80-100 ml/min**
- Blue tube for passage of : **AIR**
- Green tube for passage of : **WATER**
- Weight : **100 g**
- WATER-AIR temperature system temperature
- Storage with packaging intact for a maximum of 15 weeks with:
 - **Ambient temperature from -20°C to +70°C**
 - **Relative humidity from 10% to 90%**
 - **Atmospheric pressure from 500 to 1060 mbar.**

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GUARANTEE CERTIFICATE

FARO grants to the final customer a **12-months** guarantee starting from the date of the purchase invoice. Repairs under the guarantee shall be done at the **FARO premises or by an authorized FARO Service** person with transportation costs and risks at the expense of the purchaser. The guarantee certificate found on the endleaf of this booklet is considered valid only if it is entirely completed and bears the reseller's stamp. Repairs under guarantee (during the period of effectiveness) will be done only if the article to be repaired is accompanied by the shipping note or purchase invoice. The guarantee covers failures due to poor material quality or manufacturing defects and in case of legitimate claim the guarantee will allow repair or replacement free of charge. **No compensation for damages and/or interest will be recognized.** The guarantee is not deemed valid in the final judgement of FARO in case of tampering, damage, incorrect use, faulty maintenance or normal wear.

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Spittoon valve 3

EN

Installation and Operating Instructions

CE

7560100003L02

 **DÜRR
DENTAL**

1801V002

Contents

Important information

1	About this document	2
1.1	Warnings and symbols	2
1.2	Copyright information	2
2	Safety	3
2.1	Intended purpose	3
2.2	Intended use	3
2.3	Improper use	3
2.4	General safety information	3
2.5	Systems, connection with other devices	3
2.6	Qualified personnel	3
2.7	Protection from electric shock	3
2.8	Only use genuine parts	3
2.9	Transport	4
2.10	Disposal	4

Product description

3	Overview	5
3.1	Scope of delivery	5
3.2	Special accessories	5
3.3	Wear parts and spare parts	5
4	Technical data	6
4.1	Type plate	7
4.2	Conformity assessment	7
5	Operation	8
5.1	Operating function	8
5.2	Cleaning function	8

Installation

6	Requirements	9
6.1	Setup options	9
6.2	Preparing for the installation	9
6.3	Hose materials	9
6.4	Information about electrical connections	9
7	Installation	10
7.1	Installation overview	10
7.2	Installation of the spittoon valve	11
7.3	Establishing the compressed air connection	11
7.4	Electrical connections	12
7.5	Circuit diagram	13
8	Commissioning and first start-up	14

Operation

9	Disinfection and cleaning	15
9.1	Spittoon valve	15
9.2	Suction system	15
10	Maintenance	16

Troubleshooting

11	Tips for operators and service technicians	17
12	Transporting the unit	18

Appendix

13	Information about EMC in accordance with EN 60601-1-2	19
13.1	General information	19
13.2	Abbreviations	19
13.3	Guidelines and manufacturer's information	19

Important information

1 About this document

These installation and operating instructions form part of the unit.

If the instructions and information in these installation and operating instructions are not followed, Dürr Dental will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.

1.1 Warnings and symbols

Warnings

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The following warning symbols are used:

General warning symbol

Warning – dangerous high voltage

Biohazard warning

The warnings are structured as follows:

SIGNAL WORD

Description of the type and source of danger

Here you will find the possible consequences of ignoring the warning

- > Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

- **DANGER**
Immediate danger of severe injury or death
- **WARNING**
Possible danger of severe injury or death
- **CAUTION**
Risk of minor injuries
- **NOTICE**
Risk of extensive material/property damage

Other symbols

These symbols are used in the document and on or in the unit:

Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.

Comply with the Operating Instructions.

Wear hand protection.

Wear protective goggles.

Use a face mask.

Refer to the accompanying electronic documents.

Cleaning button

Air

Vacuum

Manufacturer

Order number

Serial number

1.2 Copyright information

All names of circuits, processes, names, software programs and units used in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from Dürr Dental.

2 Safety

Dürr Dental has designed and constructed this device so that when used properly and for the intended purpose there is no danger to people or property. Nevertheless, residual risks can remain. You should therefore observe the following notes.

2.1 Intended purpose

The spittoon valve is designed for installation in a treatment unit in dental surgeries or dental clinics.

The installation of the spittoon valve into a treatment unit helps to avoid suction noises emanating from the spittoon.

2.2 Intended use

2.3 Improper use

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damages resulting from this. In these cases the user/operator will bear the sole risk.

2.4 General safety information

- › When operating this device always observe all guidelines, laws, and other rules and regulations that are applicable at the site of operation.
- › Prior to each use, check condition of the device and make sure it is in perfect working order.
- › Do not convert or modify the units.
- › Observe the Installation and Operating Instructions.
- › Make the Installation and Operating Instructions available to the person operating the device at all times.

2.5 Systems, connection with other devices

Additional devices connected with medical electrical devices must be proven to conform with their corresponding IEC or ISO standards. All configurations must continue to comply with the standard requirements for medical systems (see IEC 60601-1-1 or section 16 of the 3rd edition of IEC 60601-1 respectively).

Whoever connects additional devices to medical electrical devices automatically becomes the system configurator and is responsible for ensuring that the system corresponds with the standard requirements for systems. Local laws take priority over the requirements outlined above.

2.6 Qualified personnel

Operation

Persons who operate the units must ensure safe and correct handling based on their training and knowledge.

- › Instruct or have every user instructed in handling the unit.

Installation and repairs

- › Installation, readjustments, alterations, upgrades and repairs must be carried out by Dürr Dental or by qualified personnel specifically approved and authorized by Dürr Dental.

2.7 Protection from electric shock

- › When working on the units observe all the relevant electrical safety regulations.
- › Never touch the patient and unshielded plug connections on the device at the same time.
- › Immediately replace any damaged lines and connections.

Observe the EMC rules concerning medical devices

- › Observe specific precautionary measures relating to electromagnetic compatibility (EMC) for medical devices, see "13 Information about EMC in accordance with EN 60601-1-2".
- › The appliance is designed for the use in health care establishments (in accordance with IEC 60601-1-2). If the appliance is operated in another environment, potential effects on electromagnetic compatibility must be taken into account.

2.8 Only use genuine parts

- › Only use Dürr Dental parts or accessories and special accessories specifically approved by Dürr Dental.
- › Only use only genuine working parts and spare parts.

2.9 Transport

The original packaging provides optimum protection for the device during transport.

If required, original packaging for the unit can be ordered from Dürr Dental.

 Dürr Dental does not accept any responsibility or liability for damage occurring during transport due to the use of incorrect packaging, even where the unit is still under guarantee.

- › Only transport the device in its original packaging.
- › Keep the packing materials out of the reach of children.


2.10 Disposal

 The unit may be contaminated. Instruct the company disposing of the waste to take the relevant safety precautions.

- › Decontaminate potentially contaminated parts before disposing of them.
- › Uncontaminated parts (e.g. electronics, plastic and metal parts etc.) should be disposed of in accordance with the local waste disposal regulations.
- › If you have any questions about the correct disposal of parts, please contact your dental trade supplier.

3 Overview

1 Spittoon valve

3.1 Scope of delivery

The following items are included in the scope of delivery (possible variations due to country-specific requirements and/or import regulations):

- Spittoon valve 3 7560700xxx**
- Spittoon valve

3.2 Special accessories

The following optional items can be used with the device:

- Switch control panel 7560-520-00

3.3 Wear parts and spare parts

The following working parts need to be changed at regular intervals (refer to the "Maintenance" section):

- Protective strainer 0700-702-06E

 Information on spare parts can be found on the website portal for authorised specialist dealers under: www.duerrdental.net.

4 Technical data

Electrical data

Safety low voltage	(AC/DC)	
	V	24
Frequency	Hz	50 - 60
Nominal current	A	0.1
Rated power	W	2.4
Type of protection		IP 21

Electrical data, suction unit relay

Switching voltage	(AC/DC)	
min.	V	5
max.	V	24
Switching current		
min.	mA	10
max.	A	2

Connections

Supply and waste water connection		
DürrConnect	mm	Ø 20
Collection vessel vent connection	mm	Ø 9
Compressed air connection	mm	Ø 4

Media

Compressed air		
min.	bar / MPa	3 / 0.3
max.	bar / MPa	5 / 0.5
Fluid flow rate, max.	l/min	3.5
Fluid temperature, max.	°C	35
Suction system pressure		
max.	mbar/hPa	-200
Absolute	mbar/hPa	800

General data

Duty cycle	%	40
Medical device		Class I
Weight	g	240
Dimensions (H x W x D)	cm	143 x 75 x 110

Ambient conditions during storage and transport


Temperature	°C	-10 to +60
Relative humidity	%	< 95

Ambient conditions during operation

Temperature	°C	+10 to +40
Relative humidity	%	< 70
Air pressure	hPa	700 - 1060

4.1 Type plate

The type plates are located on the side of the fluid collector.

1 Type plate

4.2 Conformity assessment

This device has been subjected to conformity acceptance testing in accordance with the current relevant European Union guidelines. This equipment conforms to all relevant requirements.

5 Operation


Figure 1: Idle phase


Figure 2: Operational phase

- 1 Fluid collector
- 2 Air extraction seal
- 3 Vent
- 4 Protective strainer
- 5 Exhaust air damper
- 6 Solenoid valve
- 7 Compressed air connection
- 8 Shut-off valve
- 9 Float sensor
- 10 Magnet in float sensor
- 11 Float sensor detection
- 12 Magnet in cleaning button
- 13 Cleaning button detection sensor

5.1 Operating function

The waste water from the spittoon flows through the coarse filter into the collector vessel. If enough fluid is present then the magnet in the float sensor is detected by the control electronics. The control electronics start up the suction unit with the suction unit relay and actuate the solenoid valve for the compressed air supply. The inflowing compressed air opens the shut-off valve via a piston. The fluid from the collector vessel is then sucked into the suction pipe. As soon as the fill level in the collector vessel has dropped, this is detected by the control electronics and the solenoid valve is switched off. While waste water continues to flow in from the spittoon the collector vessel refills and the process starts again from the beginning.

5.2 Cleaning function

The cleaning function is activated by permanent pressure on the yellow cleaning button on the spittoon valve or on the cleaning button on the switch control panel (if present). As a result the solenoid valve for the compressed air supply, and therefore the shut-off valve, is opened and the suction unit relay is actuated in order to start up the suction unit.

The cleaning and disinfection solutions can now be aspirated without hindrance through the spittoon valve into the suction pipe and into the suction unit. A suction noise can be heard at the spittoon.

6 Requirements

6.1 Setup options

- Installation in treatment units in dental surgeries or dental clinics.

6.2 Preparing for the installation

Prior to installation of the spittoon valve the following media should be checked and if necessary adjusted; refer also to "4 Technical data":

- Vacuum of the suction system
- Compressed air supply
- Water amount from the spittoon

Do not remove the gold collector or the coarse sieves from the spittoon.

6.3 Hose materials

For waste connections and suction lines only use the following hose types:

- Flexible spiral hoses made of PVC with integrated spiral or equivalent hoses
- Hoses that are resistant to dental disinfectants and chemicals

Plastic hoses will display signs of ageing over time. Therefore, they should be inspected regularly and replaced as necessary.

The following types of hoses must not be used:


- Rubber hoses
- Completely PVC hoses
- Hoses that are not sufficiently flexible

6.4 Information about electrical connections

- › The supply voltage to the device must satisfy the requirements for two patient protection (MOPP) protective measures as set out in IEC 60601-1 in relation to the supply network.
- › The supply voltage must satisfy the following voltage/power requirements:
24 V AC/DC, 50-60 Hz, at least 2.4 VA

7 Installation


7.1 Installation overview

- 1 Spittoon valve
- 2 Station selection valve
- 3 Pressure reducer
- 4 Suction pipe connection
- 5 Rinsing unit
- 6 Auxiliary air nozzle
- 7 Hose manifold
- 8 Switch control panel
- 9 Spittoon outlet


7.2 Installation of the spittoon valve

-  The cleaning function can be activated via the yellow button. For this reason the spittoon valve should be positioned in an easily accessible location. If this is not possible, a separate switch control panel can be used as an optional accessory.
- › Disconnect the treatment unit from the power supply and secure it so that it cannot be switched back on again.
- › Firmly screw the spittoon valve onto a suitable place on the treatment unit.
- › Connect the drain hose from the spittoon to the inlet of the spittoon valve.
- › Connect the outlet of the spittoon valve to the suction pipe.

7.3 Establishing the compressed air connection

- › Disconnect a suitable compressed air line from the treatment unit.
- › Install a T-piece with 4 mm branch in the compressed air line.
- › Connect a compressed air hose to the T-piece.
- › Route the compressed air hose to the spittoon valve, cut it off straight and insert it.



-  To pull off the compressed air hose from the spittoon valve, press the black sleeve on the compressed air connection inwards.

7.4 Electrical connections

 Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).


 The requirements of IEC 60601-1 must be satisfied during installation.

- › Open the cover of the control electronics.
- › Route the power supply and control line to the spittoon valve.
- › Attach the connector to the connection lines.
- › Plug in the connector at the corresponding positions on the control electronics.
- › Secure the connection lines with cable ties to the housing.

- X1 Cleaning button for switch control panel
- X3 Solenoid valve
- X4 Control line for suction unit
- X5 Power supply
- K1 Suction unit relay
- N2 Float sensor detection
- N3 Cleaning button detection sensor

7.5 Circuit diagram

- 1 Hose manifold
- 2 Station selection valve
- 3 Rinsing unit
- 4 Spittoon valve
- X1 Cleaning button for switch control panel
- X3 Solenoid valve
- X4 Control line for suction unit
- X5 Power supply
- K1 Suction unit relay
- N2 Float sensor detection
- N3 Cleaning button detection sensor
- 5 Suction machine relay in the treatment unit

8 Commissioning and first start-up

In many countries technical medical products and electrical devices are subject to regular checks at set intervals. The owner must be instructed accordingly.

- › Turn on the unit power switch or the main surgery switch.
- › Carry out an electrical safety check in accordance with applicable local regulations (e.g. the German Ordinance on the Installation, Operation and Use of Medical Devices / Medizinprodukte-Betreiberverordnung) and record the results as appropriate, e.g. in the technical log book.
- › Carry out a functional inspection of the system and check the connections for leaks.
- › Attach and screw on the covers.

9 Disinfection and cleaning



NOTICE

Device malfunctions or damage due to use of incorrect media

Guarantee claims may become invalid as a result.

- › Do not use any foaming agents, e.g. household cleaning agents or instrument disinfection agents.
- › Do not use abrasive cleaners.
- › Do not use agents containing chlorine.
- › Do not use any solvents like acetone.

9.1 Spittoon valve

- › Switch on the rinsing for the spittoon.
- › Keep pressing the yellow cleaning button of the spittoon valve or the cleaning button on the switch control panel (if present) until rinsing of the spittoon is finished.
- › Pour disinfection solution into the spittoon and at the same time press the yellow cleaning button of the spittoon valve or the cleaning button on the switch control panel (if present) until the disinfection solution has been aspirated.

9.2 Suction system

After every treatment

- › Aspirate a glass of cold water through the large and the small suction hoses. Do this even if only the small suction hose was actually used during treatment.

Suction through the large suction hose causes a large amount of air to be drawn up, thereby considerably increasing the cleaning effect.

Daily after the end of treatment

After higher workloads before the midday break and in the evening

The following are required for disinfection/cleaning:

- Material-compatible, non-foaming disinfection/cleaning agents with Dürr Dental approval, e. g. Orotol plus.
- Unit care system, e.g. OroCup
- › To pre-clean, suck up 2 litres of water with the care system.
- › Aspirate the disinfection/cleaning agent with the care system.

Once or twice a week before the midday break

Under harsher conditions (e.g. hard water or frequent use of prophylaxis powders) 1x daily before the midday break

The following are required for cleaning:

- Material-compatible, non-foaming special cleaning agents that have been approved by Dürr Dental, e.g. MD 555 Cleaner
- Unit care system, e.g. OroCup
- › To pre-clean, suck up 2 litres of water with the care system.
- › Aspirate the cleaning agent with the care system.
- › Rinse with ca. 2 l water after the application time.

10 Maintenance

 All maintenance work must be performed by a qualified expert or by one of our Service Technicians.

 Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Maintenance interval

Maintenance work

Monthly	<ul style="list-style-type: none">› Press the yellow cleaning button to empty the collection vessel.› Clean the yellow coarse filter or replace it if required.
Annually	<ul style="list-style-type: none">› Check compressed air supply. *› Perform a function test.
Every 3 years	<ul style="list-style-type: none">› Clean the float monitor in the collection vessel. *› Clean the seal on the float sensor for ventilation and replace if required. *

* Only by customer services service technicians.

11 Tips for operators and service technicians

Any repairs above and beyond routine maintenance must only be carried out by suitably qualified personnel or by one of our service technicians.

WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

Fault	Probable cause	Solution
Spittoon valve not working	No power supply	› Check power supply and restore. *
	Faulty connections	› Check the plug connections. *
	Relay not switching	› Check the switching function of the relay. *
	No compressed air present	› Check the compressed air supply of the spittoon valve. *
	Sensor defective	› Check the function of the sensor with the aid of the button. › Check the function by manually moving the float sensor.
Suction unit does not start up or runs continuously	Float sensor does not move in its housing	› Clean the housing and float sensor. * › Insert the float sensor correctly. *
Fluid does not drain off	Drain blocked	› Clean the drain line. * › Check whether the filter is blocked, clean if necessary.

* Only by customer services service technicians.

12 Transporting the unit

WARNING

Infection due to contaminated unit


- › Disinfect the unit before transport.
- › Close all media connections.

EN

Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

- › Before disassembly, clean and disinfect the suction unit and the unit using a suitable disinfectant approved by Dürr Dental.
- › Disinfect a defective unit using a suitable surface disinfection agent.
- › Seal all connections with sealing caps.
- › Pack the unit securely in preparation for transport.

- 1 DürrConnect dummy bushing (order no. 0700-700-10E)
- 2 Protective cap (order no. 9000-412-85)
- 3 DürrConnect hose connector socket Ø 20 mm (order no. 0700-700-20E)
- 4 Protective cap (order no. 9000-412-98)
- 5 Sealing cap (order no. 9000-310-002)

13 Information about EMC in accordance with EN 60601-1-2

13.1 General information

The information in this leaflet includes excerpts from the relevant European standards for electrical, medical devices. It must be observed when installing Dürr Dental devices or combining them with products of other manufacturers. If you are uncertain about anything, please refer to the complete standard.

13.2 Abbreviations

EMC	Electromagnetic compatibility
HF	High frequency
U_T	Rated voltage of the device (supply voltage)
V_1, V_2	Compliance level for the test in acc. with IEC 61000-4-6
E_1	Compliance level for the test in acc. with IEC61000-4-3
P	Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer
d	Recommended safety distance in metres (m)

13.3 Guidelines and manufacturer's information

Electromagnetic emissions for all devices and systems

The device is designed for operation in an electromagnetic environment as specified below. The customer or operator of the device should ensure that the device is operated in such an environment.

Interference emission measurements	Compliance	Electromagnetic environment - guidelines
HF emissions in accordance with CISPR 11	Group 1	The spittoon valve uses HF energy only for its internal functions. For this reason, HF transmissions are very low and it is unlikely that any interference will be caused to neighbouring electronic devices.
HF emissions in accordance with CISPR 11	Class B	The spittoon valve is suitable for use in all facilities including those in living areas and areas that are directly connected to the public mains electricity supply that also supplies buildings used for residential purposes.
Harmonics in acc. with IEC 61000-3-2	Not applicable	
Voltage fluctuations/flickers in acc. with IEC 61000-3-3	Not applicable	

Resistance to electromagnetic interference (immunity) for all devices and systems

The device is designed for use in electromagnetic environments specified below. The customer or operator of the device should ensure that the device is operated such an environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) in acc. with IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	±8 kV contact discharge	Floors should be made of wood or cement, or covered with ceramic tiles. If the floor is covered by synthetic material, then the relative humidity must be at least 30%.
Electrical fast transient/burst immunity test in accordance with IEC 61000-4-4	±2 kV for mains cables ±1 kV for input and output cables	±2 kV for mains cables ±1 kV for input and output cables	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Voltage surge in accordance with IEC 61000-4-5	±1 kV voltage outer conductor/outer conductor ±2 kV voltage outer conductor/earth	±1 kV push-pull voltage Not applicable	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Magnetic field for a supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	30 A/m	30 A/m	The magnetic fields at mains frequency should be within the range of typical values encountered in a commercial or hospital environment.

Table 1: Resistance to electromagnetic interference (immunity) for all devices and systems

Electromagnetic interference immunity for devices or systems that are not life-sustaining

Portable and mobile communication devices should not be used any closer to the unit (including cables) than the recommended safety distance, which is calculated based on the applicable formula for the transmission frequency.

Interference immunity tests	IEC 60601 - test level	Compliance level	Recommended safety distance
Conducted HF disturbance variables in accordance with IEC 61000-4-6	$3 V_{\text{eff}}$ 150 kHz to 80 MHz	$[V_1] = 10 \text{ V}$	$d = 0.35 \cdot \sqrt{P}$
Emitted HF disturbance variables in accordance with IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	$[E_1] = 10 \text{ V/m}$	$d = 0.35 \cdot \sqrt{P}$ for 80 MHz to 800 MHz $d = 0.7 \cdot \sqrt{P}$ for 800 MHz to 2.7 GHz

Table 2: Electromagnetic interference immunity for units or systems that are operated in health-care facilities

P Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer

d Recommended safety distance in metres (m)

The field strength of stationary communication devices should be lower than the compliance level for all frequencies based on inspections on site^{a, b}.

Interference is possible in the environment of units that have the following symbols.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic radiation is affected by absorption and reflection on the building, objects and people.

^a The field strength of stationary transmitters, such as the base stations of mobile phones and land mobile radios, amateur radio stations, AM and FM radio and television broadcasters, for example, cannot be accurately predicted. In order to determine the electromagnetic environment with regard to stationary transmitters, a study of electromagnetic phenomena at the site should be considered. If the measured field strength at the location where the unit is used exceeds the compliance levels stated above, the unit should be monitored to verify that it works as intended. If unusual performance characteristics are observed additional measures may be required, such as a changing the orientation of the unit or moving it to a different location.

^b Over the frequency range of 150 kHz to 80 MHz, the field strength should be less than $[V_1]$ V/m.

Recommended safety distance between portable and mobile HF communication devices and the unit

The device is designed for use in the electromagnetic environments specified below, in which the HF disturbance variables are controlled. The customer or the operator of the device can help to prevent electromagnetic interference by maintaining the minimum distances between mobile HF communication equipment (transmitters) and the device as recommended below in accordance with the maximum output line of the communication equipment.

Keep a minimum distance of 30 cm between the device and mobile communication devices.

Rated power of the transmitter (W)	Safety distance based on the transmission frequency (m)		
	150 kHz to 80 MHz $d = 0.35 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \cdot \sqrt{P}$
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.11	1.11	2.21
100	3.5	3.5	7

Table 3: Recommended safety distance between portable and mobile HF communication devices and the unit

For transmitters whose maximum rated power is not specified in the table shown above, the recommended safety distance d in metres (m) can be determined from the formula that belongs to the respective column where P is the maximum rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic waves is affected by absorption and reflection on the building, objects and people.


Hersteller/Manufacturer:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
Fon: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

CAS 1 Combi-Separator

EN

Installation and Operating Instructions

CE

7117100018L30

 DÜRR
DENTAL

1712V002

Contents

Important information

1	About this document	3
1.1	Warnings and symbols	3
1.2	Copyright information	3
2	Safety	4
2.1	Intended purpose	4
2.2	Intended use	4
2.3	Improper use	4
2.4	Systems, connection with other devices	4
2.5	General safety information	5
2.6	Qualified personnel	5
2.7	Protection from electric shock	5
2.8	Only use genuine parts	5
2.9	Transport	5
2.10	Disposal	5

Product description

3	Overview	6
3.1	Scope of delivery	6
3.2	Special accessories	6
3.3	Disposable materials	6
3.4	Wear parts and spare parts	6
4	Technical data	7
4.1	CAS 1 Combi-Separator	7
4.2	Type plate	8
4.3	Conformity assessment	8
4.4	Approvals	8
5	Operation	9
5.1	Operation	10
5.2	Separation	10
5.3	Spittoon connections	10
5.4	Station selection valve / safety valve	10
5.5	Amalgam separation	10
5.6	Sediment level measurement	11
5.7	Operating problems	11
5.8	Service key	11

Installation

6	Requirements	12
6.1	Installation/setup room	12
6.2	Setup options	12
6.3	Hose materials	12
6.4	Installation and routing of hoses and pipes	12
6.5	Information about electrical connections	12
6.6	Information about connecting cables	12
7	Installation	13
7.1	Combining devices safely	13
7.2	Installation of the CAS 1 in treatment units	14
7.3	Electrical connections, controller	15
7.4	Electrical connections	15
8	Commissioning and first start-up	16
9	Service program	17
10	Description of the service program	18
10.1	Service program ON/OFF	18
10.2	Display test	18
10.3	Sediment level measurement	18
10.4	Motor start - motor braking	18
10.5	Input and output signals	18

Operation

11 Display/handling	19
11.1 Ready for operation	19
11.2 Amalgam collector vessel is 95% full	19
11.3 Amalgam collector vessel is 100% full	19
11.4 Amalgam collector vessel not in position	19
11.5 Motor fault	20
12 Disinfection and cleaning	20
12.1 After every treatment	20
12.2 Daily after the end of treatment	20
12.3 Once or twice a week before the midday break	21
13 Replace the amalgam collector vessel	21
13.1 Disposal of the collector vessel	22
14 Maintenance	23
14.1 Tests	24

Troubleshooting

15 Tips for operators and service technicians	25
16 Transporting the unit	28
16.1 Close CAS 1	28

Appendix

17 Information about EMC in accordance with EN 60601-1-2	29
17.1 General information	29
17.2 Abbreviations	29
17.3 Guidelines and manufacturer's information	29

Important information

1 About this document

These installation and operating instructions form part of the unit.

 If the instructions and information in these installation and operating instructions are not followed, Dürr Dental will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.

1.1 Warnings and symbols

Warnings

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The following warning symbols are used:

 General warning symbol

 Biohazard warning

The warnings are structured as follows:

 **SIGNAL WORD**
Description of the type and source of danger
 Here you will find the possible consequences of ignoring the warning
 > Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

- **DANGER**
Immediate danger of severe injury or death
- **WARNING**
Possible danger of severe injury or death
- **CAUTION**
Risk of minor injuries
- **NOTICE**
Risk of extensive material/property damage

Other symbols

These symbols are used in the document and on or in the unit:

 Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.

 Comply with the Operating Instructions.

 Wear hand protection.

 Switch off and de-energise the unit (e.g. unplug from mains).

 Hose manifold connection

 Spittoon connections

 Suction unit connection

 Drain connection

 Unit in operation

 Unit operation interrupted

 Audible signal/melody sounds

 Do not reuse

 CE labelling

 Order number

 Serial number

 Manufacturer

1.2 Copyright information

All names of circuits, processes, names, software programs and units used in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from Dürr Dental.

2 Safety

Dürr Dental has designed and constructed this device so that when used properly and for the intended purpose there is no danger to people or property. Nevertheless, residual risks can remain. You should therefore observe the following notes.

2.1 Intended purpose

The CAS 1 Combi-Separator is designed for continuous separation of liquids and air and for separation of amalgam from the entire waste water from dental treatment units.

2.2 Intended use

Installation in accordance with the requirements of the water authorities in the German Federal States or in accordance with local regulations. In accordance with the installation regulations of DIBT Berlin.

The CAS 1 Combi-Separator is designed for the separation of liquids and air as well as for the separation of amalgam from the waste water from a single treatment unit in a dry system.

Installation in dental treatment units and in practice rooms (housing version). Positioned in the suction line after the spittoon and manifold.

The minimum volume of waste water that can be supplied to the unit is 0.1 l/min, but must not exceed 4.0 l/min. Here, a separation efficiency of at least 95% is maintained with amalgam.

A rinsing unit with fresh water can be installed in the suction line upstream of the Combi-Separator.

Installation, servicing and repairs must only be performed by qualified personnel specifically approved and authorized by Dürr Dental.

The disposable amalgam containers must only be used once.

 For surgical treatments and when the Airflow is being used, the CAS 1 Combi-Separator requires a rinsing unit to be installed, which feeds a small amount of water to the device during aspiration. This thins any liquid (e.g. saliva, blood) that occurs so it can be transported more easily.

The CAS 1 Combi Separator for KaVo treatment units must be set up in a defined installation setup in order to meet the relevant safety standards. For this reason it must only be installed in the treatment units that have been designed and approved for this purpose by KaVo.

KaVo-approved treatment units:

New units delivered from 01/2016 onwards:
E50, E50 Life, E70/E80, E70/E80 Vision, 1058, 1058 Life

Spare parts requirements for old units such as 1078 and 1080 among others.

2.3 Improper use

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damages resulting from this. In these cases the user/operator will bear the sole risk.

The Combi-Separator must only be used to process liquids from the oral cavity. It must not be used for the aspiration of any other substances, such as dust, sludge, plaster or similar.

Only chemicals and disinfectants that will not damage the materials, e.g. Orotol Plus or equivalent, may be used.

The unit is not suitable for installation downstream of 2 or more workplaces. The max. water volume of 4.0 l/min must not be exceeded.

The device must not be run in continuous operation; the braking process after the end of the run is required in order to keep the centrifuge drum clean.

The device must not be installed with the drain higher than the connection piece on the device. Do not use any risers. All pipes must have a downward gradient.

Not suitable for wet rooms! Do not use this device to aspirate flammable or explosive mixtures. Do not use the unit in a potentially explosive environment!

2.4 Systems, connection with other devices

Additional devices connected with medical electrical devices must be proven to conform with their corresponding IEC or ISO standards. All configurations must continue to comply with the standard requirements for medical systems (see IEC 60601-1-1 or section 16 of the 3rd edition of IEC 60601-1 respectively).

Whoever connects additional devices to medical electrical devices automatically becomes the system configurator and is responsible for ensuring that the system corresponds with the standard requirements for systems. Local laws take priority over the requirements outlined above.

2.5 General safety information

- › When operating this device always observe all guidelines, laws, and other rules and regulations that are applicable at the site of operation.
- › Prior to each use, check condition of the device and make sure it is in perfect working order.
- › Do not convert or modify the units.
- › Observe the Installation and Operating Instructions.
- › Make the Installation and Operating Instructions available to the person operating the device at all times.

2.6 Qualified personnel

Operation

Persons who operate the units must ensure safe and correct handling based on their training and knowledge.

- › Instruct or have every user instructed in handling the unit.

Installation and repairs

- › Installation, readjustments, alterations, upgrades and repairs must be carried out by Dürr Dental or by qualified personnel specifically approved and authorized by Dürr Dental.

2.7 Protection from electric shock

- › When working on the units observe all the relevant electrical safety regulations.
- › Never touch the patient and unshielded plug connections on the device at the same time.
- › Immediately replace any damaged lines and connections.

Observe the EMC rules concerning medical devices

- › Observe specific precautionary measures relating to electromagnetic compatibility (EMC) for medical devices, see "17 Information about EMC in accordance with EN 60601-1-2".

2.8 Only use genuine parts

- › Only use Dürr Dental parts or accessories and special accessories specifically approved by Dürr Dental.
- › Only use only genuine working parts and spare parts.

DÜRR MEDICAL accepts no liability for damages or injury resulting from the use of non-approved accessories or special accessories, or from the use of non-genuine working parts or spare parts.

The use of non-approved accessories, special accessories or non-genuine working parts / spare parts (e.g. mains cable) can have a negative effect in terms of electrical safety and EMC.

2.9 Transport

The original packaging provides optimum protection for the device during transport.

If required, original packaging for the unit can be ordered from Dürr Dental.

Dürr Dental does not accept any responsibility or liability for damage occurring during transport due to the use of incorrect packaging, even where the unit is still under guarantee.

- › Only transport the device in its original packaging.
- › Keep the packing materials out of the reach of children.

2.10 Disposal


The unit may be contaminated. Instruct the company disposing of the waste to take the relevant safety precautions.

- › Decontaminate potentially contaminated parts before disposing of them.
- › Untamminated parts (e.g. electronics, plastic and metal parts etc.) should be disposed of in accordance with the local waste disposal regulations.
- › If you have any questions about the correct disposal of parts, please contact your dental trade supplier.

3 Overview

1 CAS 1 Combi-Separator

3.1 Scope of delivery

The scope of delivery can vary slightly depending on the version.

The following items are included in the scope of delivery:

- CAS 1**7117-100-51
- Combi-Separator
- Replacement disposable amalgam container
- Installation and operating instructions
- Operating Handbook

3.2 Special accessories

The following optional items can be used with the device:

- Various installation sets are available on request.
- Display panel7805-116-00E
- Cable for display panel, 1 m9000-119-043
- Cable for display panel, 3 m9000-119-042
- Station selection valve.....7560-500-60
- Station selection valve for CAS 1 / CS 1.....7560-500-80
- Vario rinsing unit.....7100-260-50
- OroCup care system.....0780-350-00
- Test vessel7117-064-00

3.3 Disposable materials

The following materials are consumed during operation of the device and must be ordered separately:

- Disposable amalgam container . . . 7117-033-00
- DürrConnect protective strainer, 5 pieces 0700-700-18E
- DürrConnect protective strainer, 5 pieces 0700-700-28E
- Orotol plus (2.5 litre bottle) CDS110P6150
- MD 550 spittoon bowl cleaner (750 ml bottle). CCS550C4500
- MD 555 cleaner (2.5 litre bottle). CCS555C6150

3.4 Wear parts and spare parts

The following working parts need to be changed at regular intervals (refer to the "Maintenance" section):

- Bellows. 7117-420-25E
- Service kit (3-year interval). 7117-980-32
- Service kit (5-year interval). 7117-980-30

Information on spare parts can be found on the website portal for authorised specialist dealers under: www.duerdental.net.

EN

4 Technical data

4.1 CAS 1 Combi-Separator

Electrical data – centrifuge motor		
Nominal voltage	V	24 AC
Frequency	Hz	50 / 60
Rated power	VA	100
Electrical data – electronics		
Nominal voltage	V	24 AC
Nominal current	A	0.2
Signal input from the hose manifold	V	24 AC/DC
Media		
Air flow volume	l/min	≤ 300
Flow rate		high
The suction system must be suitable for a high flow rate in accordance with EN ISO 10637.		
Max. pressure	hPa/mbar	-160
Min. volume of aspiration fluid	l/min	≥ 0.1
max.	l/min	≤ 1.0
Water supply, spittoon	l/min	≤ 3
Total flow of waste liquids	l/min	≤ 4
Usable volume in amalgam collecting container	ccm	c. 90
Replacement interval		4 - 6 months
General data		
Drive motor nominal speed	rpm	2800
Operating mode		S5 95% DC*
Type of protection		IP 20
Protection class		II
Noise level ** approx.	dB(A)	56
Dimensions (H x W x D)	mm	255 x 151 x 110
Weight, approx.	kg	2.7
Separation rate	%	≥ 95
Medical device (class)		I


* DC = duty cycle

** Noise levels in acc. with EN ISO 1680 "airborne noise emissions"; measured in a sound-proofed room. The levels are average values with a tolerance of ±1.5 dB(A). Higher values may be obtained in rooms with reverberating sound characteristics.

Ambient conditions during storage and transport		
Temperature	°C	-10 to +60
Relative humidity	%	< 95
Ambient conditions during operation		
Temperature	°C	+10 to +40
Relative humidity	%	< 70

4.2 Type plate

The type plates are located on the cover of the motor.

1 Type plate

4.3 Conformity assessment

This device has been subjected to conformity acceptance testing in accordance with the current relevant European Union guidelines. This equipment conforms to all relevant requirements.

4.4 Approvals


Centre of Competence in Civil Engineering, Berlin

Test number	Z-64.1-20
-------------	-----------

Separation method compliant with standard

ISO 11143	Type 1
-----------	--------

5 Operation

- 1 Fluid intake
- 2 Vacuum, to suction unit
- 3 Aspiration input
- 4 Fluid output
- 5 Motor
- 6 Separation
- 7 Separation rotor
- 8 Centrifuge
- 9 Light barriers (3x)
- 10 Sensor enclosure
- 11 Cone pump
- 12 Amalgam collector vessel
- 13 Float sensor
- 14 Fluids
- 15 Amalgam particles

5.1 Operation

CAS 1 Combi-Separator

The task of the CAS 1 combi-separator is to provide continuous separation of secretions and air as well as the amalgam separation of all the waste water from the treatment unit.

The waste water flows through the connection (1) from the spittoon directly into the centrifuge (8) and amalgam separation.

During the suction phase the aspirated secretions are separated from the aspirated air in the separation unit (6). The secretions accumulating in the separation unit are continuously transported to the centrifuge (8), where the amalgam particles are then separated.

Underneath the centrifuge is a replaceable amalgam collector vessel (12), into which the separated amalgam particles (15) are rinsed once the centrifuge (8) is switched off. A float sensor (13) checks the level within the collector vessel and sends a signal to the display panel when it needs replacing. In combination with a light barrier (9c), this float sensor also monitors whether a collector vessel is in use.

The compact size of the CAS 1 Combi-Separator allows it to be installed in dental treatment units. This results in short secretion carrying lines. After the centrifuge is switched off, the braking cycle triggers a self-cleaning process. This self-cleaning process also leads to smooth and silent running, as well as providing a separation efficiency of more than 95%, even under heaviest loads.

5.2 Separation

At the inlet connection (3) of the CAS 1, the aspirated fluid/air mix is accelerated and set into a spiral motion in the separation unit (6). The resulting centrifugal forces sling the aspirated particles against the outer wall. The air is continuously separated from the fluid and escapes via the spinning separation rotor (7) to the suction unit.

The aspirated air is subject to high centrifugal forces by the separation rotor (7), which is driven by the motor (1), which ensures that no fluid or blood foam can be carried into the suction unit.

The spiral motion feeds the separated fluid continuously to the pump wheel, which transports the fluid into the collector vessel. The fluid is transported to the centrifuge (8) via a pump cone (11).

An external station selection valve connects the CAS 1 with the suction unit via the vacuum connection (2).

5.3 Spittoon connections

The waste water from the spittoon flows through a protective strainer on the fluid inlet (1) and into the collector vessel (12). Once sufficient fluid has been collected, the float sensor (13) activates a light barrier (9a) and (9b) via a sensor housing (10) and switches on the motor (1). The fluid is transported to the centrifuge (8) via a pump cone (11).

5.4 Station selection valve / safety valve

The station selection valve has 2 tasks:

1st task:

The station selection valve interrupts the suction flow between the hose manifold and the suction unit. As soon as a suction hose is removed from the hose manifold, a solenoid valve opens the station selection valve and suction flow is enabled.

2nd task:

The station selection valve also acts as a safety valve. If the CAS 1 is over-full or not functioning properly, the system will perform a safety shutdown. This safety shutdown prevents fluids from being drawn into the dry suction pipe.

For single station suction systems, the station selection valve takes over the function of the safety valve.

In various types, a station selection valve is already integrated in the CAS 1. The station selection valve is on the connection (2) of the CAS 1.

5.5 Amalgam separation

The switches in the hose manifold or the light barrier of the sensor system switch on the motor and the associated centrifuge (8).

The fluid containing amalgam particles flows continuously to the collector vessel (12). The fluids ejected by the centrifuge are pumped through the fluid output (4) to the central waste water system.

As soon as no further fluid is fed to the amalgam separator, e.g. when the suction hose is placed back in the hose manifold, the centrifuge drum is switched off after a short delay time. This switch-off brakes the motor, as a result of which the ring of water, which continues to rotate due to inertia, rinses the separated particles out of the centrifuge (8) downwards into the collector vessel.

The separated amalgam particles form a sediment in the replaceable collector vessel. The level of fluid in the collector is regulated by the pump cone so that the risk of fluid escaping when the collector vessel is changed can be avoided.

5.6 Sediment level measurement

The fill level in the collector vessel (12) is checked by a float sensor (13) every time the main power switch is switched on.

The centrifuge motor starts, fluid is transported via the pump cone to the centrifuge drum (8) and provides a constant level of fluid (underside of the cone pump) in the collector vessel. The float sensor sinks. Two light barriers (9a) and (9b) measure the fluid level. Once the level reaches 95% in the collector vessel, this is displayed on the display panel.

5.7 Operating problems

If the unit is not ready for operation due to a fault, this will be indicated on the display panel via illuminated LEDs and an audible signal.

5.8 Service key

On the display panel there is a service key that can be used to switch off the audible signal in the event of a fill level warning or if a fault message is indicated. This button can also be used to start the device manually. To do this, press the button for longer than 2 seconds until the drive motor starts up.

6 Requirements

6.1 Installation/setup room

The room chosen for set up must fulfil the following requirements:

- Closed, dry, well-ventilated room
- Should not be a room made for another purpose (e. g. boiler room or wet cell)

6.2 Setup options

CAS 1 Combi-Separator

- Directly in the treatment unit.
- In a special housing in an extension of the treatment unit.

6.3 Hose materials

For waste connections and suction lines only use the following hose types:

- Flexible spiral hoses made of PVC with integrated spiral or equivalent hoses
- Hoses that are resistant to dental disinfectants and chemicals

Plastic hoses will display signs of ageing over time. Therefore, they should be inspected regularly and replaced as necessary.

The following types of hoses must not be used:

- Rubber hoses
- Completely PVC hoses
- Hoses that are not sufficiently flexible

6.4 Installation and routing of hoses and pipes

- › Execute the on-site pipe installation in accordance with the applicable local regulations and standards.
- › Lay the hose installation of the drains to or from the unit at a sufficient incline.

If incorrectly laid, the hoses can become blocked with sedimentation.

6.5 Information about electrical connections

- › Ensure that electrical connections to the mains power supply are carried out in accordance with current valid national and local regulations and standards governing the installation of low voltage units in medical facilities.
- › Install an all-pole disconnect switch with a contact opening width of at least 3 mm in the electrical connection to the mains power supply.
- › Observe the current consumption of the devices that are to be connected.
- › Install electrical lines without mechanical tension.
- › Make the electrical connection via the main power switch of the treatment unit or via the main power switch of the practice.

6.6 Information about connecting cables

Mains supply cable

Installation type	Line layout (minimum requirements)
Fixed installation	- Plastic sheathed cable (e.g. type NYM-J)
Flexible	- PVC flexible line (e.g. H05 VV-F) or - Rubber connection (e.g. H05 RN-F or H05 RR-F)

Control cable

Installation type	Line layout (minimum requirements)
Fixed installation	– Shielded sheathed cable (e.g. (N)YM (St)-J)
Flexible	– PVC data cable with shielded cable sheathing, as used for telecommunications and IT processing systems (e.g. type LiYCY) or – Lightweight PVC control cable with shielded cable sheathing

Wire cross-section

Unit feed:

– 0.75 mm²

Connection external valves / units:

– 0.5 mm²

7 Installation

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

7.1 Combining devices safely

Take care when connecting units together or to parts of other systems as there is always an element of risk (e.g. due to leakage currents).

- › Only connect units when there can be no question of danger to operator or to patient.
- › Only connect units when it is safe to do so and there is no risk of damage or harm to the surroundings.
- › If it is not 100% clear from the unit data sheet that such connections can be safely made or if you are in any doubt, always get a suitably qualified person (e.g. the manufacturer) to verify that the setup is safe.

A copy of the system manufacturer's declaration in accordance with Article 12 of Directive 93/42/EEC can be found in our download section at www.duerrdentel.com (document no. 9000-461-264).

7.2 Installation of the CAS 1 in treatment units


WARNING

Infection due to contaminated unit

- Clean and disinfect the suction before working on the unit.
- Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Attach the unit vertically at a suitable position in the treatment unit. The unit is mounted on rubber pads and suspended in a metal frame. This mounting arrangement prevents the transmission of any vibrations to the treatment unit while the device is running. Vibrations may occur if the unit is not positioned vertically. A minimum distance of 3 mm must be maintained to the surroundings.

Station selection valve

In various types, the station selection valve is directly mounted on the CAS 1. The station selection valve (for separate installation) should be fitted in the suction pipe in the treatment unit, preferably near the end connection in the floor socket. In some installation setups the station selection valve also functions as a safety valve, so its actuation must be implemented via the CAS 1.


For further information, refer to the station selection valve installation and operating instructions

Inlet and outlet hoses

Connect and attach the inlet and outlet hoses with DürrConnect connectors to the relevant connections on the unit. Route the hoses at an incline.

Recommended diameter of the connection hoses: Ø 25 mm.

The minimum nominal width for the outlet hose is 15 mm.

- 1 Hose manifold
- 2 Spittoon
- 3 Outlet
- 4 Suction unit

Spittoon connections

In some dental units it is possible that noises can be heard at the spittoon, which are amplified by the funnel shape of the spittoon itself. In this case, the outlet between spittoon and CAS 1 should be bled. A corresponding siphon trap with ventilation is available as a special accessory.

Rinsing unit

It is recommended that the suction system is equipped with a rinsing unit, e.g. in the treatment unit. The rinsing unit provides a small amount of water during aspiration. This dilutes the aspirated fluids (blood, saliva, rinsing water, etc.), which can then be transported more effectively.

For further information, refer to the rinsing unit installation and operating instructions


7.3 Electrical connections, controller

Power supply:

– Safety transformer order number:
9000-150-46

or

– Safety transformer 24 V AC with a with an isolator consisting of two means of patient protection (MOPP) between the mains circuit and secondary circuit, min. 100 VA, secondary fuse T 4 AH (or IEC 60127-2/V T 4 AH, 250 V)

- X1 Power supply in accordance with EN 60601-1, 24 V AC
- X2 Signal input 24 V AC/DC
- X3.1 Place selection valve / safety valve (only CAS 1, max. output 8 W)
- X3.2 Rinsing unit (CAS 1 only)
- X4 CAN bus
- X6 Display panel, external (X6A = connection for predecessor model)
- X7 Sensor technology
- X8 Production interface
- X9 Motor
- H1 Motor control display
- H2 Manifold control display
- H3 Place selection valve control display
- H4 Control display, collecting container missing

- A Place selection valve
- B Rinsing unit
- C Suction unit relay (alternative)
- D Display panel, external

7.4 Electrical connections

Station selection valve / safety valve

› Connect the station selection valve / safety valve using a 2-core wire with connector to the X3 connection of the control.

Rinsing unit

› Connect the rinsing unit using a 2-core wire with connector to the X3 connection of the control.

At the connection for the rinsing unit, a suction unit relay, for example, can be connected if there is no isolation present between the suction unit signal and station selection valve in the treatment unit. Note the power consumption of the suction unit relay.

Display panel

The display panel is used to indicate messages acoustically and visually (via LEDs).

A display panel is already integrated in the unit and should be visible/audible at all times.

If the display panel is not visible/audible, fit an additional display panel in an easily visible location. The display panel is connected to the X6 socket (RJ-45 socket). An existing Dürr Dental display panel with a 6-pin connector can be connected to the X6A connector when replacing an older device.

If the installation of the amalgam separator in a neighbouring room or in the basement results in distances of more than 3 m, we recommend installing a shielded network cable with RJ-45 sockets.

8 Commissioning and first start-up


In many countries technical medical products and electrical devices are subject to regular checks at set intervals. The owner must be instructed accordingly.

- › Turn on the unit power switch or the main surgery switch.
- › Carry out an electrical safety check in accordance with applicable local regulations (e.g. the German Ordinance on the Installation, Operation and Use of Medical Devices / Medizinprodukte-Betreiberverordnung) and record the results as appropriate, e.g. in the technical log book.
- › Check the aspiration function.
- › Check the start function via the spittoon.
- › Check the connections, hoses and device for leaks.


EN

9 Service program


1


2


3


4

5

6

EN

10 Description of the service program

Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

The various unit functions can be checked with the aid of the service program.

The individual program steps are:

- Display test
- Sediment level measurement
- Motor start and motor braking with rpm check
- Input and output signals

Function of the service key:

- By pressing the service key twice the next individual program step is called up.
- By pressing the service key once that program step is repeated.

A press of the service key is confirmed by an audible signal.

10.1 Service program ON/OFF

On

- Press the service key and switch on the voltage supply to the unit.
- As soon as a signal melody can be heard, release the service key.

The green, yellow and orange LEDs on the display panel light up (display test) and the service program is activated.

Off

Switch off the main supply to the unit.

10.2 Display test

The display test is activated as soon as the service program is started.

The LEDs on the display panel are checked.

All three LEDs must come on. There is also an audible signal, which can be switched off by pressing the service button.

10.3 Sediment level measurement

While the service program is activated, the safety check for the collector vessel is deactivated.

The sediment level measurement can be used to check the function of the sediment sensor and the function of the LEDs.

Every time the service key is pressed, the sediment level is checked. If a test collector vessel is used for this, the different levels can be scanned and made visible on the display panel.

While changing the collectors (collector vessel - test collector vessel) in the service program the unit remains in the ON state.

10.4 Motor start - motor braking

The drive motor starts and, after approx. 5 seconds, braking is applied. If the service key is pressed during these 5 seconds, the motor will immediately be braked.

This procedure can be repeated by pressing the service key 1x again.

The drive motor starts up.


As a result of the rpm monitoring, the LED will change from orange to green upon start-up and from green to orange upon braking.

10.5 Input and output signals

- After this program item is activated, the yellow LED flashes and a cycled DC voltage (approx. 22-30 V) can be measured at the terminal for the rinsing unit.
- If the suction hose is lifted off the hose manifold the green LED will also come on.
- Removal of the collecting container causes the red LED to illuminate.

11 Display/handling

- 1 Display panel
- 2 RED display
- 3 GREEN LED
- 4 Audible signal/melody
- 5 Reset/service key
- 6 YELLOW LED

11.1 Ready for operation

- Green LED is on

11.2 Amalgam collector vessel is 95% full

- Yellow LED is on
- Green LED is on

Audible signal melody sounds

- At a fill level of 95%, the signal melody can be switched off by pressing the reset button. The device is then ready for operation again.
- The yellow LED comes on as a reminder that the amalgam collector vessel is due to be changed. The level display is repeated every time the unit is switched on at the main power switch.

We recommend changing the amalgam collector vessel when it reaches 95% full.

11.3 Amalgam collector vessel is 100% full

- Yellow LED is on
- Red display flashes
- Audible signal melody sounds

- At a fill level of 100% the signal melody can no longer be switched off by pressing the reset button.
- The collecting container needs to be replaced. Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).
- The separator will not be ready for operation again until the amalgam collecting container has been replaced

11.4 Amalgam collector vessel not in position

- Red display flashes
- Audible signal

- Press the reset button briefly to switch off the audible signal.
- Switch off the device.
- Insert the collecting container.
- Switch on the unit.
- Green LED lights up – "Ready for operation"

If this error message occurs when the collecting container is correctly inserted, this indicates that there is a technical defect – inform your Service Technician.

EN

11.5 Motor fault

-  Red display and
 -  green LED flash alternately
 -  Audible signal
- Press the reset button briefly to switch off the audible signal.
 - If the reset button is pressed for longer than 2 seconds the unit can be restarted.
 - Green LED lights up – "Ready for operation"
-  If, after pressing the reset button repeatedly, the fault report reappears again each time, this indicates a technical defect – inform your Service Technician.

12 Disinfection and cleaning

NOTICE

Device malfunctions or damage due to use of incorrect media

Guarantee claims may become invalid as a result.

- › Do not use any foaming agents, e.g. household cleaning agents or instrument disinfection agents.
- › Do not use abrasive cleaners.
- › Do not use agents containing chlorine.
- › Do not use any solvents like acetone.

12.1 After every treatment

- › Aspirate a glass of cold water through the large and the small suction hoses. Do this even if only the small suction hose was actually used during treatment.

Suction through the large suction hose causes a large amount of air to be drawn up, thereby considerably increasing the cleaning effect.

12.2 Daily after the end of treatment

After higher workloads before the midday break and in the evening

The following are required for disinfection/cleaning:

- Material-compatible, non-foaming disinfection/cleaning agents with Dürr Dental approval, e. g. Orotol plus.
- Unit care system, e.g. OroCup
- › To pre-clean, suck up 2 litres of water with the care system.
- › Aspirate the disinfection/cleaning agent with the care system.

12.3 Once or twice a week before the midday break

Under harsher conditions (e.g. hard water or frequent use of prophylaxis powders) 1x daily before the midday break

The following are required for cleaning:

- Material-compatible, non-foaming special cleaning agents that have been approved by Dürr Dental, e.g. MD 555 Cleaner
- Unit care system, e.g. OroCup
- > To pre-clean, suck up 2 litres of water with the care system.
- > Aspirate the cleaning agent with the care system.
- > Rinse with ca. 2 l water after the application time.

13 Replace the amalgam collector vessel

WARNING

Risk of contamination if the amalgam collector vessel is reused since the collector vessel is not water-tight.

- > Do not use the collector vessel more than once (disposable item).

Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

We strongly recommend that the amalgam collector vessel should only be changed in the morning before the start of work. This will prevent fluid from dripping out of the drum while it is being changed.

- > Disconnect all power from the device.
- > Remove the full amalgam collecting container and from the device.
- > Pour disinfectant for suction units (e.g. Orotol plus, 30 ml) into the full amalgam collector vessel.
- > Close and secure the full amalgam collector vessel using the cap. Observe the markings on the cap and on the collector vessel.
- > Place the securely closed amalgam collector vessel into its original packaging and seal.
- > Insert a new amalgam collector vessel in the unit and clamp it in position. Only use original amalgam collector vessels.
- > Switch on the power supply. The device is ready for operation again.

13.1 Disposal of the collector vessel

Used amalgam collector vessels must not be sent in the post!

Dürr Dental is not a waste management company and is not allowed by law to accept any filled amalgam collector vessels.

- › Arrange to have filled amalgam collector vessels collected from the surgery by a local waste management company.
- › New amalgam collector vessels should be ordered from your specialist dental equipment retailer.
- › Document the replacement and legally compliant disposal of the filled waste amalgam collector vessel in the Operating Handbook.

In some countries the owner is required to keep an operating handbook. This operating handbook must document all maintenance work, service work, checks and amalgam disposal.

14 Maintenance

All maintenance work must be performed by a qualified expert or by one of our Service Technicians.

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Maintenance interval	Maintenance work
Dependent upon the level of usage of the device	<ul style="list-style-type: none"> › Replace the amalgam collecting container when a fill level of 95% or 100% is indicated on the display panel › Clean or replace protective sieves during replacement of the amalgam collecting container. At the latest, do this when the suction or draining power of the device decreases.
Annually	<ul style="list-style-type: none"> › Cleaning of the suction unit in accordance with the operating instructions. › Clean the float. * › Replace the bellows. *
Every 3 years	<ul style="list-style-type: none"> › Replace the rubber grommets on the connections. * › Replace the float. *
Every 5 years	<ul style="list-style-type: none"> › Replace the centrifuge drum and seal. * › Replace all O-rings (from the replacement parts kit) in the device. * › Replace the rubber grommets on the connections. * › Replace the float. *

* Only by customer services service technicians.

14.1 Tests

WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

In some countries the owner is required to keep an operating handbook. This operating handbook must document all maintenance work, service work, checks and amalgam disposal.

Annual inspection

This inspection should only be carried out by suitably trained staff.

For inspection, the following are required:

- Test vessel

Work steps to be performed:

- › General functional check (e.g. aspiration, spittoon inlet)
- › Service program

The following measurement times apply to fill level measurements with a test vessel:

- For a fill level of 95%, the measurement result is displayed after approx. 30 sec, whereby the drive motor is briefly switched off during the measurement.
- At a fill level of 100% the measurement result is displayed after approx. 90 sec continuous running.

Inspection of the general operating condition every 5 years

This inspection must be carried out every 5 years (in accordance with the German Waste Water Regulations, Annex 50, Dental Treatment) by an inspector in accordance with national regulations.

For inspection, the following are required:

- Test vessel
- Measuring beaker

Work steps to be performed:

- › Fill the test vessel with water and insert it into the unit.
- › Start the device and wait until it switches off again.
- › Once the device has switched off, remove the test vessel and measure the remaining amount of water.

The unit is working correctly if:

- there is at minimum content of 140 ml in the **test vessel**.

If there is less fluid, clean the centrifuge drum or check the operation of the unit.

15 Tips for operators and service technicians

Any repairs above and beyond routine maintenance must only be carried out by suitably qualified personnel or by one of our service technicians.


WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

- 1 Display panel
- 2 RED display
- 3 GREEN LED
- 4 Audible signal/melody
- 5 Reset/service key
- 6 YELLOW LED

Fault	Probable cause	Solution
Device not "ready for operation"	The main power switch of the treatment unit or surgery is not switched on	› Main power switch ON
No display on the display panel.	If an external display panel is fitted: cable not correctly connected	› Check cable connections

Fault	Probable cause	Solution
Yellow display is on GREEN LED illuminates Audible signal melody sounds	Amalgam collecting container is 95% full	› Change the amalgam collecting container.
	Float sensor dirty or blocked	› If this display occurs repeatedly even when the collecting container is empty, check that the float sensor can move freely.
Yellow display is on Red display flashes Audible signal melody sounds	Amalgam collecting container is 100% full	› Change the amalgam collecting container. Audible signal can no longer be switched off.
	Float sensor dirty or blocked	› If this display occurs repeatedly even when the collecting container is empty, check that the float sensor can move freely.
	Waster water line/siphon trap dirty	› Clean the waste water line/siphon trap. *
The RED and GREEN displays flash alternately Audible signal	Motor is dirty or defective	› Check motor for smooth running; replace the centrifuge if necessary. * › Replace the device. *
	Contact problems at X9	› Plug in the connector correctly. * › Replace the PCB main board and connector on the motor. *
Orange LED flashes Audible signal		Press the service key briefly to switch off the audible signal
	Amalgam collecting container not correctly in position	› Switch OFF the device. › Insert the amalgam collecting container in the correct position. › Switch ON the device.
	Float sensor missing	› Insert the float sensor. *
Water accumulating in the spittoon	Coarse sieve in the fluid inlet blocked	› Clean the coarse sieve.
	Outlet ineffective or not vented	› Check or retrofit the ventilation. *
Suction power too weak or interrupted	Coarse sieve is blocked on the inlet of the aspiration	› Clean the coarse sieve.
	Place selection valve not or incompletely open	› Check the control voltage. * › Clean the place selection valve. *
Device running continuously	Float sensor blocked in water start position	› Clean the float. * › Free up the float sensor linkage so that it can move freely. *
	Start signal at the signal input (X2)	› Check the control voltage. *
	Waster water line/siphon trap dirty	› Clean the waste water line/siphon trap. *
Noise at the spittoon	Outlet ineffective or not vented	› Check or retrofit the ventilation. *

EN

Fault	Probable cause	Solution
Increased vibration of the device	Pump cone dirty	› Clean or replace the pump cone. *
	Centrifuge dirty	› Clean or replace the centrifuge. *
	Water supply too low	› Introduce water into the suction pipe. › Retrofit the rinsing unit. * › Check the rinsing unit for its correct installation position. * › Check the function of the rinsing unit. *
Water cannot be pumped away or only insufficiently	Centrifuge dirty	› Clean or replace the centrifuge
	Waster water line/siphon trap dirty	› Clean the waste water line/siphon trap

* Only by customer services service technicians.

16 Transporting the unit

WARNING

Infection due to contaminated unit


- › Disinfect the unit before transport.
- › Close all media connections.

Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

- › Before disassembly, clean and disinfect the suction unit and the unit using a suitable disinfectant approved by Dürer Dental.
- › Disinfect a defective unit using a suitable surface disinfection agent.
- › Seal all connections with sealing caps.
- › Pack the unit securely in preparation for transport.

16.1 Close CAS 1

- 1 Dummy bushing
- 2 Ring clamp

17 Information about EMC in accordance with EN 60601-1-2

17.1 General information

The information in this leaflet includes excerpts from the relevant European standards for electrical, medical devices. It must be observed when installing Dürr Dental devices or combining them with products of other manufacturers. If you are uncertain about anything, please refer to the complete standard.

17.2 Abbreviations

EMC	Electromagnetic compatibility
HF	High frequency
U_T	Rated voltage of the device (supply voltage)
V_1, V_2	Compliance level for the test in acc. with IEC 61000-4-6
E_1	Compliance level for the test in acc. with IEC61000-4-3
P	Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer
d	Recommended safety distance in metres (m)

17.3 Guidelines and manufacturer's information

Electromagnetic emissions for all devices and systems

The device is designed for operation in an electromagnetic environment as specified below. The customer or operator of the device should ensure that the device is operated in such an environment.

Interference emission measurements	Compliance	Electromagnetic environment - guidelines
HF emissions in accordance with CISPR 11	Group 1	The separator uses HF energy exclusively for internal functions. For this reason, HF emissions are very low and it is unlikely that any interference will be caused to neighbouring electronic devices.
HF emissions in accordance with CISPR 11	Class B	The separator is suitable for use in all facilities including those in living areas and areas that are directly connected to the public mains electricity supply which also supply buildings used for residential purposes.
Harmonics in acc. with IEC 61000-3-2	Class A	
Voltage fluctuations/flickers in acc. with IEC 61000-3-3	Compliant	

Resistance to electromagnetic interference (immunity) for all devices and systems

The device is designed for use in electromagnetic environments specified below. The customer or operator of the device should ensure that the device is operated such an environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) in acc. with IEC 61000-4-2	±8 kV contact discharge ±2; 4; 8; 15 kV air discharge	±8 kV contact discharge ±2; 4; 8; 15 kV air discharge	Floors should be made of wood or cement, or covered with ceramic tiles. If the floor is covered by synthetic material, then the relative humidity must be at least 30%.
Electrical fast transient/burst immunity test in accordance with IEC 61000-4-4	±2 kV for mains cables ±1 kV for input and output cables	±2 kV for mains cables ±1 kV for input and output cables	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Voltage surge in accordance with IEC 61000-4-5	±0,5; 1 kV voltage outer conductor-outer conductor	±0.5; 1 kV symmetrical voltage	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Voltage drops, short-term interruptions and fluctuations of the supply voltage in accordance with IEC 61000-4-11	< 5% U_T (> 95% drop in U_T) for 1/2 period 70% U_T (30% drop in U_T) for 25 periods < 5% U_T (> 95% drop in U_T) for 5 s	< 5% U_T (> 95% drop in U_T) for 1/2 period 70% U_T (30% drop in U_T) for 25 periods < 5% U_T (> 95% drop in U_T) for 5 s	The quality of the supply voltage should correspond to a typical commercial or hospital environment. If the operator of the device needs the unit to continue working even if the mains power supply is interrupted, we recommend powering the device from an uninterruptible power supply (UPS) or from a battery.
Magnetic field for a supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	30 A/m	30 A/m	The magnetic fields at electrical frequency should be within the range of typical values encountered in a commercial or hospital environment.

Table 1: Resistance to electromagnetic interference (immunity) for all devices and systems

Electromagnetic interference immunity for devices or systems that are not life-sustaining

Portable and mobile communication devices should not be used any closer to the unit (including cables) than the recommended safety distance, which is calculated based on the applicable formula for the transmission frequency.

Interference immunity tests	IEC 60601 - test level	Compliance level	Recommended safety distance
Conducted HF disturbance variables in accordance with IEC 61000-4-6	$3 V_{\text{eff}}$ 150 kHz to 80 MHz	$[V_1] = 3 V$	$d = 1.2 \cdot \sqrt{P}$
Emitted HF disturbance variables in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	$[E_1] = 3 V/m$	$d = 1.2 \cdot \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \cdot \sqrt{P}$ for 800 MHz to 2.5 GHz

P Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer

d Recommended safety distance in metres (m)

The field strength of stationary communication devices should be lower than the compliance level for all frequencies based on inspections on site^{a, b}

Interference is possible in the environment of units that have the following symbols.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic radiation is affected by absorption and reflection on the building, objects and people.

^a The field strength of stationary transmitters, such as the base stations of mobile phones and land mobile radios, amateur radio stations, AM and FM radio and television broadcasters, for example, cannot be accurately predicted. In order to determine the electromagnetic environment with regard to stationary transmitters, a study of electromagnetic phenomena at the site should be considered. If the measured field strength at the location where the unit is used exceeds the compliance levels stated above, the unit should be monitored to verify that it works as intended. If unusual performance characteristics are observed additional measures may be required, such as a changing the orientation of the unit or moving it to a different location.

^b Over the frequency range of 150 kHz to 80 MHz, the field strength should be less than $[V_1]$ V/m.

Recommended safety distance between portable and mobile HF communication devices and the unit

The device is designed for use in the electromagnetic environments specified below, in which the HF disturbance variables are controlled. The customer or the operator of the device can help to prevent electromagnetic interference by maintaining the minimum distances between mobile HF communication equipment (transmitters) and the device as recommended below in accordance with the maximum output line of the communication equipment.

Keep a minimum distance of 30 cm between the device and mobile communication devices.

Rated power of the transmitter (W)	Safety distance based on the transmission frequency (m)		
	150 kHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \cdot \sqrt{P}$
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Table 2: Recommended safety distance between portable and mobile HF communication devices and the unit

For transmitters whose maximum rated power is not specified in the table shown above, the recommended safety distance d in metres (m) can be determined from the formula that belongs to the respective column where P is the maximum rated power of the transmitter in watts (W) in accordance with the specifications of the transmitter manufacturer.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic waves is affected by absorption and reflection on the building, objects and people.


Hersteller/Manufacturer:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
Fon: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

CS 1 Combi-Sepamatic

EN

Installation and Operating Instructions

CE

9000-606-39/30

 DÜRR
DENTAL

1709V004

Contents

Important information

1 About this document	2
1.1 Warnings and symbols	2
1.2 Copyright information	2
2 Safety	3
2.1 Intended purpose	3
2.2 Intended use	3
2.3 Improper use	3
2.4 Systems, connection with other devices	3
2.5 General safety information	4
2.6 Qualified personnel	4
2.7 Protection from electric shock	4
2.8 Only use genuine parts	4
2.9 Transport	4
2.10 Disposal	4

Product description

3 Overview	5
3.1 Scope of delivery	5
3.2 Accessories	5
3.3 Special accessories	5
3.4 Disposable materials	5
3.5 Wear parts and spare parts	5
4 Technical data	6
4.1 CS 1 Combi-Sepamatic	6
4.2 Type plate	7
4.3 Conformity assessment	7
5 Operation	8
5.1 Separation	8
5.2 Station selection valve	8

Installation

6 Requirements	9
6.1 Setup options	9
6.2 Hose materials	9
6.3 Installation and routing of hoses and pipes	9
6.4 Information about electrical connections	9
6.5 Information about connecting cables	9
7 Installation	10
7.1 Installation of the CS 1 in treatment units	10
7.2 Electrical connections, controller	11
7.3 Electrical connections	11
8 Commissioning and first start-up	12

Operation

9 Disinfection and cleaning	13
9.1 After every treatment	13
9.2 Daily after the end of treatment	13
9.3 Once or twice a week before the midday break	13
10 Maintenance	14

Troubleshooting

11 Tips for operators and service technicians	15
11.1 Replacing the fuse	16
12 Transporting the unit	16
12.1 Close the CS 1	16

Appendix

13 Information about EMC in accordance with EN 60601-1-2	17
13.1 General information	17
13.2 Abbreviations	17
13.3 Guidelines and manufacturer's information	17
13.4 Calculation table	21

EN

Important information

1 About this document

These installation and operating instructions form part of the unit.

If the instructions and information in these installation and operating instructions are not followed, Dürr Dental will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.

1.1 Warnings and symbols

Warnings

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The following warning symbols are used:

General warning symbol

Biohazard warning

The warnings are structured as follows:

SIGNAL WORD

Description of the type and source of danger

Here you will find the possible consequences of ignoring the warning

➤ Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

- **DANGER**
Immediate danger of severe injury or death
- **WARNING**
Possible danger of severe injury or death
- **CAUTION**
Risk of minor injuries
- **NOTICE**
Risk of extensive material/property damage

Other symbols

These symbols are used in the document and on or in the unit:

Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.

Comply with the Operating Instructions.

Wear hand protection.

Switch off and de-energise the unit (e.g. unplug from mains).

Hose manifold connection

Suction unit connection

Drain connection

CE labelling

Order number

Serial number

Manufacturer

1.2 Copyright information

All names of circuits, processes, names, software programs and units used in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from Dürr Dental.

2 Safety

Dürr Dental has designed and constructed this device so that when used properly and for the intended purpose there is no danger to people or property. Nevertheless, residual risks can remain. You should therefore observe the following notes.

2.1 Intended purpose

The CS 1 Combi-Sepamatic is designed for the continuous separation of liquids in the suction flow of dental treatment units in dry suction systems.

2.2 Intended use

Installation in accordance with the requirements of the water authorities in the German Federal States or in accordance with local regulations.

Installation in dental treatment units and in practice rooms (housing version).

Positioned in the suction line of a dry suction system after the spittoon and manifold.

The CS 1 Combi-Sepamatic is designed for the continuous separation of liquids in the suction flow of an individual dental treatment unit in a dry suction system.

The CS 1 Combi-Sepamatic can be operated in continuous operation.

The permissible volume of supply water amounts to min. 0.1l/min, but must not exceed 2.0l/min.

A rinsing unit with fresh water can be installed upstream of the CS 1 Combi-Sepamatic.

For surgical treatment and when using the Airflow, the CS 1 Combi-Sepamatic requires a rinsing unit to be installed, which feeds a small amount of water to the device during aspiration. This thins any liquid (e.g. saliva, blood) that occurs so it can be transported more easily.

Installation, servicing and repairs must only be performed by qualified personnel specifically approved and authorized by Dürr Dental.

The CS 1 Combi-Sepamatic is unable to separate amalgam. The treatment of waste water containing amalgam requires the connection of an amalgam separator to the CS 1 Combi-Sepamatic.

2.3 Improper use

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damages resulting from this. In these cases the user/operator will bear the sole risk.

The Combi-Sepamatic must only be used to process liquids emerging from the oral cavity and not any other substances such as dust, sludge, plaster or similar.

Only chemicals and disinfectants that will not damage the materials, e.g. OrotolPlus or equivalent, must be used.

The unit is not suitable for installation downstream of 2 or more workplaces. The max. water volume of 2.0l/min must not be exceeded.

The device must not be installed with the drain higher than the connection piece on the device. Do not use any risers. All pipes must have a downward gradient.

Not suitable for wet rooms. Do not use this device to aspirate flammable or explosive mixtures. Do not use the unit in a potentially explosive environment.

2.4 Systems, connection with other devices

Additional devices connected with medical electrical devices must be proven to conform with their corresponding IEC or ISO standards. All configurations must continue to comply with the standard requirements for medical systems (see IEC 60601-1-1 or section 16 of the 3rd edition of IEC 60601-1 respectively).

Whoever connects additional devices to medical electrical devices automatically becomes the system configurator and is responsible for ensuring that the system corresponds with the standard requirements for systems. Local laws take priority over the requirements outlined above.

2.5 General safety information

- › When operating this device always observe all guidelines, laws, and other rules and regulations that are applicable at the site of operation.
- › Prior to each use, check condition of the device and make sure it is in perfect working order.
- › Do not convert or modify the units.
- › Observe the Installation and Operating Instructions.
- › Make the Installation and Operating Instructions available to the person operating the device at all times.

2.6 Qualified personnel

Operation

Persons who operate the units must ensure safe and correct handling based on their training and knowledge.

- › Instruct or have every user instructed in handling the unit.

Installation and repairs

- › Installation, readjustments, alterations, upgrades and repairs must be carried out by Dürr Dental or by qualified personnel specifically approved and authorized by Dürr Dental.

2.7 Protection from electric shock

- › When working on the units observe all the relevant electrical safety regulations.
- › Never touch the patient and unshielded plug connections on the device at the same time.
- › Immediately replace any damaged lines and connections.

Observe the EMC rules concerning medical devices

- › Observe specific precautionary measures relating to electromagnetic compatibility (EMC) for medical devices, see "13 Information about EMC in accordance with EN 60601-1-2".

2.8 Only use genuine parts

- › Only use Dürr Dental parts or accessories and special accessories specifically approved by Dürr Dental.
- › Only use only genuine working parts and spare parts.

DÜRR MEDICAL accepts no liability for damages or injury resulting from the use of non-approved accessories or special accessories, or from the use of non-genuine working parts or spare parts.

The use of non-approved accessories, special accessories or non-genuine working parts / spare parts (e.g. mains cable) can have a negative effect in terms of electrical safety and EMC.

2.9 Transport

The original packaging provides optimum protection for the device during transport.

If required, original packaging for the unit can be ordered from Dürr Dental.

Dürr Dental does not accept any responsibility or liability for damage occurring during transport due to the use of incorrect packaging, even where the unit is still under guarantee.

- › Only transport the device in its original packaging.
- › Keep the packing materials out of the reach of children.

2.10 Disposal


The unit may be contaminated. Instruct the company disposing of the waste to take the relevant safety precautions.

- › Decontaminate potentially contaminated parts before disposing of them.
- › Uncontaminated parts (e.g. electronics, plastic and metal parts etc.) should be disposed of in accordance with the local waste disposal regulations.
- › If you have any questions about the correct disposal of parts, please contact your dental trade supplier.

3 Overview

1 CS 1 Combi-Sepamatic

3.1 Scope of delivery

The scope of delivery can vary slightly depending on the version.

The following items are included in the scope of delivery:

CS 1.....**7117-100-7x**
or

CS 1.....**7117-100-8x**

- Combi-Sepamatic
- or Combi-Sepamatic inc. station selection valve
- Rinsing unit
- Installation and Operating Instructions

3.2 Accessories

The following articles are necessary for the operation of the unit, depending on the application:

Various installation sets are available on request
 Safety transformer 24 V, 100VA...9000-150-46
 Station selection valve for
 CAS 1 / CS 17560-500-80

3.3 Special accessories

The following optional items can be used with the device:

Station selection valve.....7560-500-60
 Rinsing unit II7100-250-50
 OroCup care system.....0780-350-00

3.4 Disposable materials

The following materials are consumed during operation of the device and must be ordered separately:

DürrConnect protective strainer,
 5 pieces0700-700-18E
 Orotol plus
 4 x 2.5 l bottles/carton CDS110P6150
 MD 550 spittoon bowl cleaner
 6 x 800 ml bottles / cardboard
 box..... CCS550A4750
 MD 555 cleaner
 4 x 2.5 L bottle / carton CCS555C6150

3.5 Wear parts and spare parts

The following working parts need to be changed at regular intervals (refer to the "Maintenance" section):

- Protective strainer
- Rubber grommets
- O-rings

Information on spare parts can be found on the website portal for authorised specialist dealers under:
www.duerrdental.net.

4 Technical data

4.1 CS 1 Combi-Sepamatic

Electrical data – centrifuge motor

Nominal voltage	V	24 AC
Frequency	Hz	50 / 60
Rated power	VA	70

Media

Fluid volume		
min.	l/min	≥ 0.1
max.	l/min	≤ 2.0
Air flow volume	l/min	≤ 350
Flow rate		high
The suction system must be suitable for a high flow rate in accordance with EN ISO 10637.		
Max. pressure	hPa/mbar	-160

General data

Drive motor nominal speed	min ⁻¹	2800
Operating mode		S1 100% DC*
Type of protection		IP 20
Protection class		II
Noise level ** approx.	dB(A)	49
Dimensions (H x W x D)	cm	12.5 x 15 x 12
Weight, approx.	kg	1.4
Medical device (class)		I

* DC = duty cycle

** Noise levels in acc. with EN ISO 1680 "airborne noise emissions"; measured in a sound-proofed room. The levels are average values with a tolerance of ±1.5 dB(A). Higher values may be obtained in rooms with reverberating sound characteristics.

Ambient conditions during storage and transport


Temperature	°C	-10 to +60
Relative humidity	%	< 95

Ambient conditions during operation

Temperature	°C	+10 to +40
Relative humidity	%	< 70

4.2 Type plate

The type plates are on the motor cover and on the motor flange.


1 Type plate

4.3 Conformity assessment

This device has been subjected to conformity acceptance testing in accordance with the current relevant European Union guidelines. This equipment conforms to all relevant requirements.

5 Operation

- 1 Motor
- 2 Vacuum, to suction unit
- 3 Separation
- 4 Aspiration input
- 5 Pump wheel
- 6 Separation rotor
- 7 Fluid output
- 8 Waste valve
- 9 Relief valve

5.1 Separation

Every time the suction hose is taken out of the hose manifold, the CS 1 Combi-Sepamatic and the suction unit are started.

The mixture of liquid and air drawn up is accelerated in the intake connection and then set in spiral motion in the separation. The resulting centrifugal forces sling the aspirated particles against the outer wall. The air is continuously separated from the fluid and escapes to the suction unit via the spinning separation rotor.

The aspirated air is subject to high centrifugal forces by the separation rotor, which ensures that no fluid or blood foam can be carried into the suction unit.

The spiral motion serves to continuously transport the separated liquid to the pump wheel, this then pumps the liquid into the central waste water drainage system via the waste water valve.

The air bleed is carried out via the relief valve. If fluid escapes upwards into the air bleed area following a fault, the relief valve closes automatically.

5.2 Station selection valve

The station selection valve interrupts the suction flow between the hose manifold and the suction unit. As soon as a suction hose has been removed from the hose manifold, the station selection valve is opened and suction flow is enabled.

A station selection valve is already integrated in various versions of the CS 1. An external station selection valve can be electrically controlled via the CS 1.

6 Requirements

6.1 Setup options

CS 1 Combi-Sepamatic

- Directly in the treatment unit.
- In a special housing in an extension of the treatment unit.

6.2 Hose materials

For waste connections and suction lines only use the following hose types:

- Flexible spiral hoses made of PVC with integrated spiral or equivalent hoses
- Hoses that are resistant to dental disinfectants and chemicals

Plastic hoses will display signs of ageing over time. Therefore, they should be inspected regularly and replaced as necessary.

The following types of hoses must not be used:

- Rubber hoses
- Completely PVC hoses
- Hoses that are not sufficiently flexible

6.3 Installation and routing of hoses and pipes

- › Execute the on-site pipe installation in accordance with the applicable local regulations and standards.
- › Lay the hose installation of the drains to or from the unit at a sufficient incline.

If incorrectly laid, the hoses can become blocked with sedimentation.

6.4 Information about electrical connections

- › Ensure that the electrical connections to the mains power supply are established in accordance with current valid national and local regulations and standards governing the installation of low voltage units in medical facilities.

- › Install an all-pole disconnect switch with a contact opening width of at least 3 mm in the electrical connection to the mains power supply. It must be possible to secure the disconnect switch so that it cannot be inadvertently switched back on again.
- › Install electrical lines without mechanical tension.
- › Make the electrical connection via the main power switch of the treatment unit or via the main power switch of the practice.

6.5 Information about connecting cables

Mains supply cable

Installation type	Line layout (minimum requirements)
Fixed installation	– Plastic sheathed cable (e.g. type NYM-J)
Flexible	– PVC flexible line (e.g. H05 VV-F) or – Rubber connection (e.g. H05 RN-F or H05 RR-F)

Control cable

Installation type	Line layout (minimum requirements)
Fixed installation	– Shielded sheathed cable (e.g. (N)YM (St)-J)
Flexible	– PVC data cable with shielded cable sheathing, as used for telecommunications and IT processing systems (e.g. type LiYCY) or – Lightweight PVC control cable with shielded cable sheathing

Wire cross-section

Unit feed:

- 0.75 mm²

Connection external valves / units:

- 0.5 mm²

7 Installation

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

7.1 Installation of the CS 1 in treatment units


WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Attach the unit vertically at a suitable position in the treatment unit. The unit is mounted on rubber pads and suspended in a metal frame. This mounting arrangement prevents the transmission of any vibrations to the treatment unit while the device is running. Vibrations may occur if the unit is not positioned vertically. A minimum distance of 3 mm must be maintained to the surroundings.

Station selection valve

In various types, the place selection valve is directly mounted on the CS 1. The station selection valve (for separate installation) should be fitted in the suction pipe in the treatment unit, preferably near the end connection in the floor socket. The electrical connection should then also be carried out on the CS 1.



For further information, refer to the station selection valve installation and operating instructions

Inlet and outlet hoses

Connect and attach the inlet and outlet hoses with DürrConnect connectors to the relevant connections on the unit. Route the hoses at an incline.

Recommended diameter of the connection hoses: \varnothing 25 mm.

The minimum nominal width for the outlet hose is 15 mm.

- 1 Hose manifold
- 2 Vent
- 3 Outlet
- 4 Suction unit

Rinsing unit

It is recommended that the suction system is equipped with a rinsing unit, e.g. in the treatment unit. The rinsing unit provides a small amount of water during aspiration. This dilutes the aspirated fluids (blood, saliva, rinsing water, etc.), which can then be transported more effectively.

For further information, refer to the rinsing unit installation and operating instructions

Installation sets

Installation sets and detailed documentation for various installation situations are available from the manufacturers.


When installed in a housing, ventilation slits should be provided to avoid heat build-up in the housing.

7.2 Electrical connections, controller


Power supply:

- Safety transformer order number: 9000-150-46
- or
- Safety transformer 24 V AC with a with an isolator consisting of two means of patient protection (MOPP) between the mains circuit and secondary circuit, min. 100 VA, secondary fuse T 4 AH (or IEC 60127-2/V T 4 AH, 250 V)

- X1 Power supply in accordance with EN 60601-1, 24 V AC
- X2 Signal input 24 V AC/DC
- X3 Place selection valve and/or rinsing unit (max output 8 W)
- F1 T 4 AH, 250 V in accordance with IEC 60127-2

7.3 Electrical connections

- X1 Power supply
- X2 Hose manifold start signal
- X3 Outgoing signal station selection valve and/or rinsing unit

- › Remove the motor cover of the CS 1.
- › Attach the connector to the connection lines.

 To open, lift the terminal lever upwards.

- › Plug the connector onto the control.
- › Put the motor cover on.

8 Commissioning and first start-up

 In many countries technical medical products and electrical devices are subject to regular checks at set intervals. The owner must be instructed accordingly.

- › Turn on the unit power switch or the main surgery switch.
- › Carry out an electrical safety check in accordance with applicable local regulations (e.g. the German Ordinance on the Installation, Operation and Use of Medical Devices / Medizinprodukte-Betreiberverordnung) and record the results as appropriate, e.g. in the technical log book.
- › Check the aspiration function.
- › Check the connections, hoses and device for leaks.

9 Disinfection and cleaning

NOTICE

Device malfunctions or damage due to use of incorrect media

Guarantee claims may become invalid as a result.

- › Do not use any foaming agents, e.g. household cleaning agents or instrument disinfection agents.
- › Do not use abrasive cleaners.
- › Do not use agents containing chlorine.
- › Do not use any solvents like acetone.

9.1 After every treatment

- › Aspirate a glass of cold water through the large and the small suction hoses. Do this even if only the small suction hose was actually used during treatment.

Suction through the large suction hose causes a large amount of air to be drawn up, thereby considerably increasing the cleaning effect.

9.2 Daily after the end of treatment

After higher workloads before the midday break and in the evening

The following are required for disinfection/cleaning:

- Material-compatible, non-foaming disinfection/cleaning agents with Dürr Dental approval, e. g. Orotol plus.
- Unit care system, e.g. OroCup
- › To pre-clean, suck up 2 litres of water with the care system.
- › Aspirate the disinfection/cleaning agent with the care system.

9.3 Once or twice a week before the midday break

Under harsher conditions (e.g. hard water or frequent use of prophylaxis powders) 1x daily before the midday break

The following are required for cleaning:

- Material-compatible, non-foaming special cleaning agents that have been approved by Dürr Dental, e.g. MD 555 Cleaner
- Unit care system, e.g. OroCup
- › To pre-clean, suck up 2 litres of water with the care system.
- › Aspirate the cleaning agent with the care system.
- › Rinse with ca. 2 l water after the application time.

10 Maintenance

 All maintenance work must be performed by a qualified expert or by one of our Service Technicians.

 Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

WARNING

Infection due to contaminated unit


- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Maintenance interval

Dependent upon the level of usage of the device

Maintenance work

- › Clean or replace the protective sieves at the aspiration inlet. At the latest, however, when the suction power of the unit diminishes.

Annually	<ul style="list-style-type: none"> › Cleaning of the suction unit in accordance with the operating instructions. › Clean or replace the protective sieves at the aspiration inlet. › Where a rinsing unit is present: Clean the sieve/coarse filter in the water intake. * › Perform a function test. *
Every 3 years	› Replace the rubber grommets on the connections. *
Every 5 years	<ul style="list-style-type: none"> › Replace the rubber grommets on the connections. * › Replace all o-rings in the unit. *

* Only by customer services service technicians.

11 Tips for operators and service technicians

Any repairs above and beyond routine maintenance must only be carried out by suitably qualified personnel or by one of our service technicians.

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).

WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

Fault	Probable cause	Solution
Device does not start	No power supply	<ul style="list-style-type: none"> › Check power supply. * › Check the fuses and replace if necessary. *
	No start signal	› Check the control voltage at the signal input. *
Suction power too weak or interrupted	Coarse sieve is blocked on the inlet of the aspiration	› Clean the coarse sieve.
	Place selection valve not or incompletely open	<ul style="list-style-type: none"> › Check the control voltage. * › Clean the place selection valve. *

* Only by customer services service technicians.


11.1 Replacing the fuse

Prior to working on the device or in case of danger, disconnect it from the mains (e. g. pull the mains plug).


Transformer

- › Unscrew and remove the safety cover.
- › Replace the fuse.

Fuse housing

- › Turn the fuse housing to open it.
- › Replace the fuse.

- 1 Fuse housing
- 2 Fuses

12 Transporting the unit

WARNING

Infection due to contaminated unit


- › Disinfect the unit before transport.
- › Close all media connections.

Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

- › Before disassembly, clean and disinfect the suction unit and the unit using a suitable disinfectant approved by Dürr Dental.
- › Disinfect a defective unit using a suitable surface disinfection agent.
- › Seal all connections with sealing caps.
- › Pack the unit securely in preparation for transport.

12.1 Close the CS 1

- 1 Dummy bushing
- 2 Ring clamp

13 Information about EMC in accordance with EN 60601-1-2

13.1 General information

The information in this leaflet includes excerpts from the relevant European standards for electrical, medical devices. It must be observed when installing Dürr Dental devices or combining them with products of other manufacturers. If you are uncertain about anything, please refer to the complete standard.

13.2 Abbreviations

EMC	Electromagnetic compatibility
HF	High frequency
U_T	Rated voltage of the device (supply voltage)
V_1, V_2	Compliance level for the test in acc. with IEC 61000-4-6
E_1	Compliance level for the test in acc. with IEC61000-4-3
P	Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer
d	Recommended safety distance in metres (m)

13.3 Guidelines and manufacturer's information

Electromagnetic emissions for all devices and systems

Interference emission measurements	Compliance	Electromagnetic environment - guidelines
HF emissions in accordance with CISPR 11	Group 1	The CS 1 Combi-Sepamatic uses HF energy exclusively for internal functions. As a result, HF-transmissions are very low and it is highly unlikely that any interference will be caused to neighbouring electronic devices.
HF emissions in accordance with CISPR 11	Class B	The CS 1 Combi-Sepamatic is suitable for use in all facilities including those in living areas and areas that are directly connected to the public mains electricity supply that also supplies buildings used for residential purposes.
Harmonics in acc. with IEC 61000-3-2	Class A	
Voltage fluctuations/flickers in acc. with IEC 61000-3-3	Compliant	

Resistance to electromagnetic interference (immunity) for all devices and systems

The device is designed for use in electromagnetic environments specified below. The customer or operator of the device should ensure that the device is operated such an environment.

Interference immunity tests	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) in acc. with IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be made of wood or cement, or covered with ceramic tiles. If the floor is covered by synthetic material, then the relative humidity must be at least 30%.
Electrical fast transient/burst immunity test in accordance with IEC 61000-4-4	±2 kV for mains cables ±1 kV for input and output cables	±2 kV for mains cables ±1 kV for input and output cables	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Voltage surge in accordance with IEC 61000-4-5	±1 kV voltage outer conductor/outer conductor ±2 kV voltage outer conductor/earth	±1 kV push-pull voltage ±2 kV common mode voltage	The quality of the supply voltage should correspond to a typical commercial or hospital environment.
Voltage drops, short-term interruptions and fluctuations of the supply voltage in accordance with IEC 61000-4-11	< 5% U_T (> 95% drop in U_T) for 1/2 period 40% U_T (60% drop in U_T) for 5 periods 70% U_T (30% drop in U_T) for 25 periods < 5% U_T (> 95% drop in U_T) for 5 s	< 5% U_T (> 95% drop in U_T) for 1/2 period 40% U_T (60% drop in U_T) for 5 periods 70% U_T (30% drop in U_T) for 25 periods < 5% U_T (> 95% drop in U_T) for 5 s	The quality of the supply voltage should correspond to a typical commercial or hospital environment. If the operator of the device needs the unit to continue working even if the mains power supply is interrupted, we recommend powering the device from an uninterruptible power supply (UPS) or from a battery.
Magnetic field for a supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	The magnetic fields at mains frequency should be within the range of typical values encountered in a commercial or hospital environment.

Table 1: Resistance to electromagnetic interference (immunity) for all devices and systems

Electromagnetic interference immunity for devices or systems that are not life-sustaining

Portable and mobile communication devices should not be used any closer to the unit (including cables) than the recommended safety distance, which is calculated based on the applicable formula for the transmission frequency.

Interference immunity tests	IEC 60601 - test level	Compliance level	Recommended safety distance
Conducted HF disturbance variables in accordance with IEC 61000-4-6	$3 V_{\text{eff}}$ 150 kHz to 80 MHz	$[V_1]$ V	$d = [3.5 / V_1] \cdot \sqrt{P}$ $d = 1.2 \cdot \sqrt{P}$
Emitted HF disturbance variables in accordance with IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	$[E_1]$ V/m	$d = [3.5 / E_1] \cdot \sqrt{P}$ for 80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$ for 80 MHz to 800 MHz $d = [7 / E_1] \cdot \sqrt{P}$ for 800 MHz to 2.7 GHz $d = 2.3 \cdot \sqrt{P}$ for 800 MHz to 2.7 GHz

Table 2: Electromagnetic interference immunity for units or systems that are operated in health-care facilities

P Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer

d Recommended safety distance in metres (m)

The field strength of stationary communication devices should be lower than the compliance level for all frequencies based on inspections on site^a.

Interference is possible in the environment of units that have the following symbols.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic radiation is affected by absorption and reflection on the building, objects and people.

^a The field strength of stationary transmitters, such as the base stations of mobile phones and land mobile radios, amateur radio stations, AM and FM radio and television broadcasters, for example, cannot be accurately predicted. In order to determine the electromagnetic environment with regard to stationary transmitters, a study of electromagnetic phenomena at the site should be considered. If the measured field strength at the location where the unit is used exceeds the compliance levels stated above, the unit should be monitored to verify that it works as intended. If unusual performance characteristics are observed additional measures may be required, such as a changing the orientation of the unit or moving it to a different location.

^b Over the frequency range of 150 kHz to 80 MHz, the field strength should be less than $[V_1]$ V/m.

Recommended safety distance between portable and mobile HF communication devices and the unit

The device is designed for use in the electromagnetic environments specified below, in which the HF disturbance variables are controlled. The customer or the operator of the device can help to prevent electromagnetic interference by maintaining the minimum distances between mobile HF communication equipment (transmitters) and the device as recommended below in accordance with the maximum output line of the communication equipment.

Rated power of the transmitter (W)	Safety distance based on the transmission frequency (m)		
	150 kHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Table 3: Recommended safety distance between portable and mobile HF communication devices and the unit

For transmitters whose maximum rated power is not specified in the table shown above, the recommended safety distance d in metres (m) can be determined from the formula that belongs to the respective column where P is the maximum rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer.

Comment 1 The higher frequency range applies for 80 MHz and 800 MHz.

Comment 2 These guidelines may not apply in all cases. The propagation of electromagnetic waves is affected by absorption and reflection on the building, objects and people.

13.4 Calculation table

If the measured values deviate from the standard, the values are specified in chapter "4 Technical data".

The safety distances can then be calculated in the tables shown below.

P:

V_1 :

E_1 :

P Rated power of the transmitter in Watts (W) in accordance with the specifications of the transmitter manufacturer

V_1 Compliance level for the test in acc. with IEC61000-4-6

E_1 Compliance level for the test in acc. with IEC61000-4-3

Interference immunity tests	IEC 60601 - test level	Compliance level	Recommended safety distances
Conducted HF disturbance variables in accordance with IEC 61000-4-6	3 V_{eff} 150 kHz to 80 MHz	$[V_1]$ V	$d = [3.5 / V_1] \cdot \sqrt{P}$
Emitted HF disturbance variables in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	$[E_1]$ V/m	$d = [3.5 / E_1] \cdot \sqrt{P}$ for 80 MHz to 800 MHz $d = [7 / E_1] \cdot \sqrt{P}$ for 800 MHz to 2.5 GHz

Rated power of the transmitter (W)	Safety distance based on the transmission frequency (m)		
	150 kHz to 80 MHz $d = [3.5/V_1] \cdot \sqrt{P}$	80 MHz to 800 MHz $d = [3.5/E_1] \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d = [7 / E_1] \cdot \sqrt{P}$
0.01			
0.1			
1			
10			
100			

Hersteller/Manufacturer:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
Fon: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

01

ALYA

LAMPADA DENTALE A LED
DENTAL LED LIGHT

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
CONTENT

1	SAFETY REQUIREMENTS	36
1.1	SYMBOLS USED.....	36
1.1.1	Symbols used in the manual	36
1.1.2	Symbols on the labels	36
1.1.3	Symbols on the packaging.....	37
1.2	PLANNED USE.....	37
1.3	PLANNED USER.....	37
1.3.1	Professional qualification:	37
1.3.2	Minimum skills.....	37
1.3.3	Experience	37
1.3.4	Possible user handicaps	37
1.4	GENERAL STANDARDS AND MAIN WARNINGS.....	38
1.5	STORAGE AND USE: ENVIRONMENTAL PROVISIONS.....	39
1.6	REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY	39
2	GENERAL FEATURES.....	41
2.1	DESCRIPTION OF THE PRODUCT	41
2.2	DESCRIPTION OF THE PARTS	42
2.3	DEVICE IDENTIFICATION.....	42
3	DEVICE INSTALLATION.....	44
3.1	DIMENSIONS	44
3.2	COMPLETE DENTAL LIGHT	45
3.2.1	Electrical requirements	45
3.2.2	Complete combined version light assembly	46
3.2.3	Remote cable connection.....	46
3.2.4	Complete dental light connection to the ceiling, wall or floor	47
3.2.5	Installation applications.....	47
3.3	HEAD.....	56
3.3.1	Mechanical requirements.....	56
3.3.2	Electrical requirements	57
4	INSTRUCTIONS FOR USE.....	58
4.1	SWITCH ON AND OFF.....	58
4.1.1	Adjustment:	58
4.1.2	Light / Head WITH PROXIMITY.....	58
4.1.3	Light / Complete light with Theia Tech / “ALYA” head WITH REMOTE CONTROL.....	59
4.1.4	SYNCRONIZATION command.....	59
4.2	JOYSTICK ALYA JOYSTICK LEVER ASSEMBLY	61
5	ORDINARY MAINTENANCE.....	61
6	CLEANING	61
6.1	PARABOLA CLEANING	61
6.2	CLEANING THE HEAD.....	62
6.3	ARM CLEANING	62
7	STERILISING THE HANDLES	62
7.1	Removing the handles	62
7.2	Decontamination and disinfectant	62
7.3	Sterilisation.....	63
8	PERIODIC CHECKS.....	63
9	ACOUSTIC ALARM SIGNALS.....	63
9.1	Acoustic Signals.....	63
9.2	Guide to problems	64
10	Technical features.....	65
10.1	WIRING DIAGRAMS.....	66

1. SAFETY REQUIREMENTS

Dear Client,

FARO hopes you enjoy your work with the new high quality light. For safe work and to take full advantage of the performance of the product, carefully read this manual before using the device. In particular, follow all the warnings and the notes given.

1.1 SYMBOLS USED

1.1.1 Symbols used in the manual

	WARNING
The paragraphs marked with this symbol contain instructions that must be carefully followed to avoid damaging the device, harming the operator or the patient.	
	WARNING
These instructions warn you that you must pay attention to avoid situations that could damage the device.	
	FORBIDDEN
This icon highlights what you should not do to avoid damaging the device.	
	NOTES
This icon supplies information that allows you to use the device more efficiently.	

1.1.2 Symbols on the labels






The data plate is fixed:

- for the complete light: on the rear arm
- for the head: under the heat sink cover

and outlines the following data:

Serial Number (SN): year (YY) / range of origin (LD for dental light - TE for head only) plus the progressive number (NNNNNN) e.g.: SN14LD000001 for the complete light SN 14TE000001 for the head.

The following standardised symbols are also present:

	Read the use instructions. Supply electricity
	Manufacturing symbol pursuant to Directive 93/42/EEC
	The use instructions include safety warnings
	WEEE equipment in compliance with the Directive 1012/19/EC. Dispose of the product in compliance with this directive.
	Double insulation. Class 2 device against electrical risk
	Serial Number
	Can be sterilised with humid heat at 134°C
	Use the device at a temperature between 10°C and 40°C
	Use the device at pressure between 800mbar and 1060mbar
	Use the device at relative humidity between 30 RH and 75RH
	Symbol to adjust light intensity
	Symbol to switch on the light
	[Symbol to switch on/off the light on the rear arm]

1.1.3 Symbols on the packaging

HIGH

FRAGILE

DO NOT WET

DO NOT ROLL

DO NOT USE HOOKS


500 kg max


MAXIMUM STACKABLE WEIGHT

**STORAGE TEMPERATURE
TRANSPORT CONDITIONS**

RELATIVE HUMIDITY

ATMOSPHERIC PRESSURE

RECYCLABLE CARDBOARD

1.2 PLANNED USE

The device is intended for use solely in dental practices by dentists and hygienists to light the operative and intervention site in the treatment of oral cavity and dental treatments.

In normal use, the device is positioned over the patient's body at a distance of 700mm, the distance for which the lighting features were designed.

The treated patients can be of all ages with typical dental pathologies.

1.3 PLANNED USER

The planned user is a dentist or dental hygienist.

1.3.1 Professional qualification:

- Degree in medicine with dentistry specialisation
- Degree in dentistry
- Degree in dental hygiene

1.3.2 Minimum skills

- Those planned for the professional qualification
- Understanding of language: Those acquired for the professional qualification

1.3.3 Experience

- Those outlined to conduct the profession

1.3.4 Possible user handicaps

- For use, complete use is necessary of an upper limb.
- Visual faculty compatible with the profession

1.4 GENERAL STANDARDS AND MAIN WARNINGS

- The device can be applied to the dental unit, but also be installed on specific applications. The device can be powered both by the dental unit and by a power supply unit connected directly to the mains. See the specific installation paragraph.
- The device does not have Essential Performances for which inadequacy of the device performance does not prejudice patient safety.
- The device does not provide life support.
- The device must be clean before use (see Device Cleaning paragraph).
- The packaging of the light is suitable to adequately protect it from penetration of external agents.

	<p>Warnings against electrical danger and fire</p> <p>Do not use the light in the event parts are damaged. Installation of the device must only be carried out by qualified staff. The dental light must be installed on a specific control and power supply device, such as dental units, or with an electrical system that meets standard IEC 60364-1 and “national installation regulations for electrical systems in premises for medical use”. The appliance must be installed with an omnipolar separation device from the mains and compliant with Standard IEC 61058-1. Installation and maintenance of device conformity with the standard IEC 60601-1 is the responsibility of the installation technician or the manufacturer of the combined units. Check the power supply voltage, indicated on the data plate, corresponds to the mains voltage. Do not carry out any maintenance on the light when the power supply is inserted: then disconnect the power supply from the mains before intervening.</p>
	<p>Warning against danger of wear of the mechanical parts and falling suspended weights</p> <p>Do not use detergents containing the following to clean plastic parts: AMMONIUM HYDROXIDE - SODIUM HYDROXIDE - METHYLENE CLORIDE – METHYL ALCOHOL. Non-compliance with this specification could cause: RISK OF WEAR ON PLASTIC PARTS RESULTING IN BREAKAGE. Do not spray any kind of chemical agent directly on the light. In particular, use is forbidden of abrasive substances, acids, containing chlorine.</p>
	<p>Warnings for danger of suspended loads falling</p> <p>Strictly comply with compliance for maximum loads planned. Do not knock against or overload the limit switches on the arms and heads.</p>
	<p>Warnings for biological danger and glare</p> <p>Do not fasten or focus the light strip directly in the patient’s eyes, especially patients at greater risk of eye injury (e.g. children with eye diseases). In this case, always use appropriate guards and precautions. The light is classified as exempt from photo-biological risk in compliance with EN 62471. However, it cannot be excluded that particularly photo-sensitive patients or those who have taken photo-sensitising medicine can develop a rash or allergic reactions to light. In this case, suspend the treatment and use very low lighting levels. The articulated arm and the joints of the light allow correct positioning of the light strip.</p>
	<p>Warning for danger of damaging the electrical parts</p> <p>Do not overcharge the arms and the joints with end of stroke knocks. Rotation of the head and arms as well as the limit switches can damage the conductor insulations.</p>
	<p>Warning for danger of explosion</p> <p>The device is not suitable for installation in environments with the presence of inflammable gas or risks of oxygen.</p>
	<p>Warning for danger of patient-patient cross contamination</p> <p>The dentist must use disposable protection on the handles of the light and sterilise them after each patient. To disinfect the surfaces use water-alcohol mixed disinfectant (see maintenance/cleaning paragraph).</p>
	<p>Warning for danger of wrong maintenance</p> <p>Do not carry out maintenance operations or replacements of parts other than those outlined in the manual. Any intervention not indicated in the manual could compromise the safe appearance of the device. Only carry out the maintenance operations in the manual; in any other case, contact technical support.</p>

The product is covered by WEEE Directive 2012/19/EU
 To scrap and dispose of the materials, comply with the standard in force in your country, if necessary contacting recognised and authorised specialist companies.
 At the end of the life cycle, divide the materials based on their type (ferrous, rubber, plastic). Do not leave small parts of the equipment unguarded or within reach of exposed people (children) because they are a potential sources of danger.
 The company FARO does not allow any changes to the product not expressly authorised in writing, otherwise conformity to the safety standards and the warranty expires.
 Other warnings are outlined in the titles of this manual.

1.5 STORAGE AND USE: ENVIRONMENTAL PROVISIONS


The appliance in the original packaging can be transported or kept in a warehouse for a period of 15 weeks if the following environmental conditions are met:

- Environmental temperature from -20°C to + 70°C
- Relative humidity from 10% to 90%
- Atmospheric pressure from 500 to 1060 mbar


The appliance must be used in the following environmental conditions:

- Temperature from 10° to 40°C
- Max altitude: 2000 m
- Relative humidity from 30% to 75%


ATMOSPHERIC PRESSURE

RELATIVE HUMIDITY

USE TEMPERATURE

1.6 REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This medical device requires particular precautionary measures to ensure electromagnetic compatibility, and must be installed and used in compliance with the information provided in the accompanying documentation.

Manufacturer's guide and declaration – Electromagnetic emissions		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Emission tests	Conformity	Guide – Electromagnetic environment
RF Emission CISPR15	Yes	The ALYA lamp uses RF energy only for internal functioning. Therefore, its RF emissions are very low and probably do not cause any interference with other nearby electronic equipment.
RF Emission CISPR15	Yes	The ALYA lamp is suitable for use in all buildings, including domestic buildings and those connected directly to the low voltage public electricity mains that power domestic buildings.
Harmonic emission	Class C	
RF Emission CISPR11 / EN 55011	Yes	The ALYA lamp is not suitable for interconnections with other devices (ceiling model).

ELECTROMAGNETIC IMMUNITY

Manufacturer's guide and declaration – Electromagnetic immunity		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Immunity test	Conformity	Guide – Electromagnetic environment
Electrostatic discharge (ESD) IEC/EN61000-4-2	± 6kV contact ± 8kV air	Flooring must be in wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC/EN61000-4-4	± 2kV power supply ± 1kV for input/output lines	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Surge IEC/EN61000-4-5	± 1kV differential mode ± 2kV common mode	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Voltage dips, short interruption and voltage variation IEC/EN61000-4-11	< 5% Ut for 0,5 cycle 40% Ut for 05 cycle 70% Ut for 25 cycle <5% Ut for 5 sec.	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment. If the user of the ALYA lamp requires continuative use, also when electrical power is not present, it is advisable to have a UPS (uninterrupted power supply).
Power frequency magnetic field IEC/EN61000-4-8	3A/m	Level of the magnetic field at the typical power main frequency of a commercial or hospital environment.
Conducted immunity IEC/EN61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting equipment)	Portable and mobile RF communication devices should not be used near any part of the dental unit, including cables, except for when they comply with the recommended separation distances calculated by the equation applicable for the transmitter frequency. Recommended separation distances: d = 1,2√P d = 1,2√P from 80 Mhz to 800 MHz d = 2,3√P from 800 MHz to 2,5 GHz Where P is the maximum nominal power emitted from the transmitter in Watts (W), according to the transmitter manufacturer, and d is the recommended separation distance in metres (m). The intensity of the fixed RF transmitter field, as determined in an electromagnetic study of the site a, could be less than the level of conformity of each frequency interval. Interference could be expected in proximity to devices bea-
Conducted immunity IEC/EN61000-4-6	3Vrms 80MHz to 2.5GHz (for non life-supporting equipment)	

Note: Ut is the value of the voltage.
 Note 1: At 80 MHz and 800 Mhz the higher frequency interval is applied.
 Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.
 a) The ISM bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.
 b) The levels of conformity of the ISM bands between 150 kHz and 80 MHz and in the bands 80 MHz to 2.5 GHz are indicative. The probability of a portable transmitter causing interference if accidentally brought in the range of a patient decreases depending on environment conditions. This is why an additional factor of 10/3 has been incorporated into the formula used to calculate the separation distance from transmitters.
 c) Field intensities for fixed transmitters like base stations for radio-telephones (cellular and cordless phones) and terrestrial transmitters, amateur radio stations, AM and FM radio transmitters and TV transmitters cannot be predicted theoretically and precisely. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic investigation of the site should be done. If the intensity of the field measured in the area of use of the dental unit is higher than the applicable level of conformity described above, it is necessary to monitor the lamp to check that it functions properly. If abnormal performance is detected, additional measures could prove necessary, such as a different orientation or position of the lamp.
 d) The field intensity in an interval of frequencies from 150 kHz to 80 MHz should be less than 3 V/m.

Recommended separation distances between portable and mobile radio-communication devices and the dental unit

The ALYA lamp is designed to function in an electromagnetic environment in which disturbances from RF emissions are kept under control. The client or operator of the unit can contribute to preventing electromagnetic interference by assuring a minimum distance between mobile and portable RF communication devices (transmitters) and the dental unit, as recommended below, in relation to the maximum output power of the radio-communication devices.

Maximum nominal output power of the transmitter W	Separation distance at the transmitter frequency m		
	150 kHz to 80 MHz d = 1,2 √P	80 MHz to 800 MHz d = 1,2 √P	800 MHz to 2,5 GHz d = 2,3 √P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal power emission not listed above, the recommended separation distance d in metres (m) can be calculated using the equation applied to the transmitter frequency, where P is the maximum nominal power emission of the transmitter in Watts (W) according to the manufacturer.

Note:
 At 80 MHz and 800 MHz, the highest interval is applied.
 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.

2. GENERAL FEATURES

2.1 DESCRIPTION OF THE PRODUCT

The device is used by the user to light the operative field when treating dental pathologies. The light source on the head is composed of two LED whose light reflects against two parabolas passing through two secondary lenses.

The parabolas enable a regular and uniform spotlight to be obtained at each lighting level and to evenly distribute the light in the operative field, without creating shadows or obstruction by the operator.

Adjustment of light intensity can be carried out with a joystick or proximity sensor. The proximity sensor switches on and off the light without having direct contact, thereby eliminating the possibility of cross infection on the command.

The “automatic switch on” or “Auto-on” function allows the light to switch on automatically each time the light power supply is activated.

The remote cable allows bringing the light commands to the combined unit. Comply with the instructions provided in the installation paragraph.

On the head near the joystick and/or the proximity sensor there is a button which allows for activation of the synchronization function with the environment lamp produced by Faro. The synchronization function allows the Alya lamp to control the level of illumination of the environment lamp in order to guarantee a more uniform level of illumination between the operating area and the surrounding zone, in order to reduce the effect of glare and to improve comfort. Visual.

In the version with light on the rear arm called “Alya with Theia Tech” the light source is composed of a series of LED whose light passes through a diffuser for distribution in the environment below.

Adjustment of the lighting level takes place in synchronous mode with that of the head for which on reducing or increasing lighting produced by the head, that on the arm also adjusts.

The “Theia Tech” light version also has a local command intended only to provide on/off on the fixed arm.

Once the local command is switched on, the light automatically synchronises with the level of head intensity. If the light on the head is switched off, the light on the rear arm switches on to maximum intensity. On switching on the light on the head, the light on the arm will synchronise automatically.


The light on the arm performs a visibility increase function in the pre-operative zone to reduce glare that generates following and/or after viewing the operative field.]

Maintenance is facilitated thanks to application of the new technologies which take into consideration the various needs for safety, ergonomics and hygiene.



The handles can be removed and sterilised. Comply with the specifications defined in the specific section.

For the electrical connections, comply with the instructions supplied in the installation paragraph and the wiring diagrams included in this manual


2.2 DESCRIPTION OF THE PARTS

- 1 – Head
- 2 – Central arm
- 3 – Rear arm without transformer with or without light (Theia Tech)
- 4 – Joints
- 5 – Pin for combined or application connection


- 1 – Joystick
- 2 – LED secondary lens
- 3 – Parabola
- 4 – Sterilisable handles
- 5 – Proximity

- 6 – Electronic board slot
- 7 – Switch on and adjustment symbols
- 8 – Synchronization button

Joystick version

Sensor version

2.3 DEVICE IDENTIFICATION

The variants on the market are differentiated by:

- Type of device (complete light, complete light with Theia Tech or head)
- Switch on and adjustment interface (**Joystick or proximity sensor, for complete lights and head**)
- Combined control mode (on-off function, remote control; for complete light and head)
- Type of mount (combined, ceiling; for complete light only)
- Arms length (for complete light only)
- Power supply (with or without transformer, for complete light only)

Codes are developed as follows:

ALYA – Complete light					
Type	3° digit Mounting and combined control	4° digit – Voltage interface	5° digit – rear arm x central arm (mm)	6° digit Available	7° 8° 9° digits Customisation
5 1	0 Standard combined	0 Joystick 17-24 V AC 22-35 V DC	0 810x550	0	000 (std FARO) JJJ
	1 Ceiling	1 Proximity sensor 17-24 V AC 22-35 V DC	1 960x550		
	2 Combined Auto-on		9 810x 855		
	4 Combined Rem cable				



ALYA – Complete light with Theia Tech					
Type	3° digit Mounting and arm dimensions	4° digit Head type	5° digit interface type	6° digit Voltage	7° 8° 9° digits Customisation
52	1 Combined 550*810	1 Standard	1 joystick	1 24Vac 50/60Hz 24 V dc	000 (std FARO) JJJ
	2 Combined 550*960	2 With extension	2 Joystick auto-on		
	3 Combined 855*810		3 sensor		
	4 Combined 855*960		4 Auto-on sensor		
	5 Ceiling 550*810		5 Joystick + remote cable		
	6 Ceiling 550*960		6 Auto-on joystick + remote cable		
	7 Ceiling 855*810		7 Sensor + remote cable		
	8 Ceiling 855*960		8 Auto-on sensor + remote cable		

ALYA HEAD					
Type	3° digit Temperature col and combined control	4° digit – Power supply and control	5° digit Available	6° digit Colours	7° 8° 9° digits Custom
3 1	5 Optical unit 5000 K	0 Joystick 17-24 V AC 22-35 V DC	0	0 White	000 Std FARO
	6 Optical unit 5000 K On/Off	1 Proximity sensor 17-24 V AC 22-35 V DC		3 Grey	
	8 L.U. 4000 K				

3. DEVICE INSTALLATION

	<p>Warnings for electrical danger, for the light and head:</p> <p>The device must be installed by specialist technicians. On installation, the power supply must always be disconnected. Refer to the wiring diagrams in the manual. Check the plate data before installation</p>
	<p>Note for installation</p> <p>The power supply cable on the complete light is supplied without any connector or terminal to allow connection according to the specifications of the combination or application. The functionality and safety of the light does not depend on the polarity of the power supply current. Therefore inversion of the power supply cables will not pose a risk of malfunctioning.</p>

3.1 DIMENSIONS

	A	B	C	D
mm	550	830	265	510
mm	550	980	265	510
mm	855	830	394	706
mm	855	980	394	706

3.2 COMPLETE DENTAL LIGHT

3.2.1 Electrical requirements

The requirements for correct installation for any application (combined, wall, floor or ceiling) are as follows:

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
Complete light version 17-24 Vac 50/60 Hz		transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T1.6AL 250V Minimum requirements: • Output: 17-24 Vac; • Power: 26 VA; • Class B; • Rigidity over 4000 V. • Thermal protection		
Complete light version 22-33Vdc	2 x 0.5 mm ² 300 V 105°C PVC insulation diameter insulation 1.85 mm Only use certified terminals and connectors with resistance to flame V1 or similar.	Power supply unit complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T630mAL 250V Minimum requirements: • Output: 22-33 Vdc; • Power: 14 VA; • Class B; • Rigidity over 4000 V. Continuous protection from short circuit or overcurrent	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.
Complete light version with Theia Tech 24Vac 50/60Hz		transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T2AL 250V Minimum requirements: • Output: 24Vac; • Power: 40VA • Class B; • Rigidity over 4000 V. Thermal protection		

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
Complete light version with Theia Tech 24Vdc	2 x 0.5 mm ² 300 V 105°C PVC insulation diameter insulation 1.85 mm Only use certified terminals and connectors with resistance to flame V1 or similar.	transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T2AL 250V Minimum requirements: • Output: 24 Vdc; • Power 28 VA • Class B; • Rigidity over 4000 V. Thermal protection	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.

Tab 1 – Requirements for electrical connection and compliance with IEC 60601-1.

Check the packaging contains the following parts:

- Dental light / Head (in the required version)
- Sheet to download for site instructions www.faro.it/download

3.2.2 Safe working load

The dental light EDI Led type suits various installation positions:

DENTAL UNIT, CEILING, WALL, FLOOR.

SAFE WORKING LOAD

	Total load (SAFE WORKING LOAD)	Safety load (MINIMUM BREAKING LOAD)
Mobile arm 855 mm	29.2 N	235 N
Mobile arm 550 mm	25.6 N	205 N

3.2.3 Complete combined version light assembly

With a digital level, ensure the connection element is perfectly parallel to the ground.

Install the light by inserting the light terminal pin in the specific combined compartment.


Connect the power supply cable according to the specifications outlined in Tab. 1.

Check the light stays balanced in all positions. If necessary, use the adjustment system on the spring to balance the light.

Check switch on and adjustment and (if present) the Auto-on command and the remote cable.

3.2.4 Remote cable connection

Connect the cable to the two buttons (A and B) with contact normally open (already supplied) according to the following diagram:

3.2.5 Installation applications



- The applications are not supplied with the light
- The device has to be installed by specialized technicians only

! The light has to be installed with FARO applications only.

⊖ The light is supplied with rotation limit switch between the fix and the mobile arm. **THE LIMIT MUST NOT BE PASSED OVER OR FORCED.**

CEILING MOUNTING INSTALLATION

1. Ceiling flange
2. Expander
3. Screw
4. Washer
5. Cable gland
6. Terminal board
7. Flange
8. Screw
9. Nut
10. Washer
11. Washer
12. Nut
13. Screw
14. Support column
15. Ring
16. Ceiling light support
17. Plug
18. Screw
19. Screw
20. Column coupling bushing
21. Lamp pin
22. Key switch
23. Fixing guide

CEILING MOUNTING

NB1. The device has to be installed by specialized technicians only

NB2. Always cut off the power supply prior to installation

NB3. Before proceeding with the assembly, ensure the ceiling can bear the application.

The authorized construction materials of the ceiling are concrete and natural stone. Use the supplied fittings or equivalent ones.


NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mAL250V~).

INSTALLATION STEPS


A. After having fixed as reference point the chair centre (A), install at a distance of 650mm and 150mm, according to the directions given in the figure.

B. Unfit the flange (7) by removing the nuts (12) and washers (11).

C. By using as guide the flange (1), carry out on the ceiling 4 bores with $\varnothing 14$ drill. Fit the expanders (2) in these bores.


D. Take the flange (1). Pass the power cable through the cable gland (5), then push the flange (1) against the ceiling; do not choke the cable between the flange (1) and the ceiling. Pass the screws (3) and washers (4) through the 4 bores used as reference in drilling the ceiling. Lock with the special wrench (installation accessory) the screws (3)

E. Connect the power cable to the terminal board (6) (see electrical diagrams)


CEILING MOUNTING

F. Match the 2 fixing guides (23) to the screws (8) and fix them by means of the nuts (9) and washers (10)

G. Calculate the proper length of the column (14), according to the formula $L=H-1900\text{mm}$. Cut the exceeding column (14) part on the side where NO LATERAL BORES ARE PRESENT

H. Insert the column (14) in the flange (7) and mark on the column (14) the position of the bores on the flange (7). Pay attention to the orientation of the column in relation to the dental unit. Remove the column and carry out two $\text{Ø}8$ bores at the marked points


I. Fit on the column (14) the ring (15) at about 300mm height (it is only a temporary position to allow the assembling)

J. Insert the ceiling light support (16) on the column (14)


K. Fit the column (14) in the special bore of the column fixing flange (7)

L. Lock the screw (13) and the two screws (18) with hexagonal spanners (installation accessory). Tighten sturdily the screw (13) and make sure the screws (18) have passed through the bores on the column (14)


CEILING MOUNTING

M. Hook the unit freshly assembled (column fixing flange (7) + column (14)) to the fixing guides (23), by matching the 4 bores of the flange (7) to the screws (8) of the flanged on the ceiling (1)


N. Screw (without locking) the nuts (12) and the remaining washers (11) on the screws (8) of the ceiling flange (1)

O. Unscrew the three screws (19) of the column (14) and remove the bushing (20)


P. Insert the bushing (20) on the lamp pin (21)

Q. Insert in the pin (21) groove the key switch (22)



CEILING MOUNTING

- R. Slip - from the top - inside the column (14) a traction cord.
- S. Connect the lamp conductor to the traction cord.
- T. Fit the lamp on the column (14) and fix it with the three screws (19); make sure the bores of the bushing (20) match the screws seat on the column (14) and tighten. Simultaneously pull the traction cord to push out the conductor of the lamp from the column fixing flange (7) of about 200 mm.
- U. Connect the lamp conductor to the terminal board (6) (see electrical diagrams).
- V. Verificare la perpendicolarità della colonna agendo sui dadi (9).
- W. Stringere i dadi (12) e le rondelle (11) per fissare la flangia (7), rendendola indipendente dalle guide di fissaggio (23).
- X. Far aderire la plafoniera (16) al soffitto, spingendovi contro l'anello (15).

WIRING DIAGRAM



CEILING MOUNTING WITHOUT transformer

WALL MOUNTING

LIGHT MOUNTING WALL MOUNTING

1. Plug
2. Screw
3. Wall mounting
4. Expanders
5. Terminal board
6. Terminal board cover
7. Cover
8. Screw
8. Washer
8. Cable gland

NB1. The device has to be installed by specialized technicians only

NB2. Always cut off the power supply at the installation site prior to installation

NB3. Before proceeding with the assembly, ensure the wall can bear the application.

The authorized construction materials of the wall are concrete and natural stone. Use the supplied fittings or equivalent ones.


NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mA/250V~).

The resulting medical system has to be declared by the installer as being compliant with IEC/EN 60601-1

WALL MOUNTING


Always cut off the power supply at the installation site prior to installation

- After having fixed as reference point the chair centre (see fig. **A-B**), carry out on the wall three d.14 bores at the bores of the wall mounting (3); pay attention the perpendicularity of the bore.
- Insert the expanders (12) in the d.14 bores previously carried out and block with the special hex spanner (installation accessory) the 2 screws; do not choke the cable between the application (3) and the wall.
- Apply three plugs (1) on the wall mounting (3) bores.
- Unscrew the screw (8) Remove the cover (7), insert the lamp arm inside the application and grease the pin. Connect the lamp cables to the terminal board (5) (see the electrical diagram below), including the earthing cable. Connect the wires which come out of the wall to the terminal board, if previously walled in. In the absence of such arrangement, the connection as to be carried out with a movable cable, to be passed in the cable gland (10).
- Fit back the cover (7) by means of screws (8).

WIRING DIAGRAM



Wall mounting

FLOOR MOUNTING

FLOOR MOUNTING

1. Pin
2. Bushing
3. Screw
4. Support column
5. plugs
6. Screw
7. Washer
8. Screw
9. Cover
10. Terminal board
11. Cable gland
12. Expanders
13. Floor support
14. Dowels
15. Plug

FLOOR MOUNTING

The device has to be installed by specialized technicians only

NB2. Always cut off the power supply prior to installation

NB3. Before proceeding with the assembly, ensure the floor can bear the application.


The authorized construction materials of the wall are concrete and natural stone. Use the supplied fitting or equivalent ones.

NB4. Maximum applicable load: 70 kg

NB5. Installation is allowed only in spaces provided with electrical plant compliant with the national legislation in force as regards spaces for medical use

NB6. The lamp without transformer must be powered with low-voltage (12-24Vac o 17-33Vdc) using a safety transformer (according to IEC/EN 60601-1) with thermal protection or at least one suitable fuse (T500mAL250V~). The resulting medical system has to be declared by the installer as being compliant with IEC/EN 60601-1


A. After having fixed as reference the centre of the chair “a”, install at a distance 650mm and 150mm according to the directions given in figure“C “

- Always cut off the power supply at the installation site prior any operation
- After having fixed as reference point (a) the chair centre (see fig.C), drill four d.14 bores in the floor, so that they match the bores on the support (13).
- Prepare the floor support (13) by passing through the washer (7)and screw (6), then screw the expanders (12) on the screws (6) a few turns, pass the power cable through the cable gland (11).
- Insert the expanders (12) in the d.14 bores previously carried out and block with the special hex spanner (installation accessory) the 6 screws; do not choke the cable between the support (13) and the floor.
- Apply four plugs (5) on the floor support (13) bores.
- Untighten the screws (8) and remove the cover plate (9) Connect the power supply cable to the terminal board (10).
- Fix the column (4) to the floor support (13); while securing it, ensure the column is perpendicular. Secure with the three screws (3) the bushing (2) to the column (4); ensure the holes of the bushing (2) match the screws on the column (4).
- Connect the light cables to the terminal board (10).
- Fit the cover plate (9) to the floor support (13) with two screws (8).

WIRING DIAGRAM

Floor mounting

3.3 HEAD

3.3.1 Mechanical requirement

For mechanical connection you must have an adequate space for the pin in the head and the G. fastening components.

The support system must be designed to support the following loads, multiplied by the safety factors outlined by IEC 60601-1 or IEC 80601-2-60

Alya Head	Alya Screen
1,80 kg	0,35 kg
2,15 kg	

For mechanical connection, follow the procedure below:

<p>1 - Support the head and insert the washers in the threaded pin in compliance with the sequence in the figure. 2 - Then insert ring nut G in compliance with the sequence indicated in the figure and screw in with adequate equipment. The ring nut must be screwed in to give the right rotational force to the head.</p>	
<p>3 – Screw in the 2 safety screws F.</p>	
	<p>Attention</p> <p>The central arm without the head load tends to rise in a sudden manner with the risk of knocking against parts of the body. During the entire installation, keep the central arm in position and do not release it until head installation is complete.</p>
	<p>Warnings for danger of suspended loads falling</p> <p>Attention - Risk of head falling after installation:</p> <ul style="list-style-type: none"> - only use screws supplied by FARO. - screw in the safety screws together.

Once mechanical connection is complete, complete electrical wiring.

3.3.2 Electrical requirements

The requirements for correct installation of **the head** are as follows:

Power supply	Power cable	Type of power supply and safety requirements	Classification	Compliance with IEC 60601-1
17-24 Vac 50/60 Hz	Power cable: 2 red unipolar cables: UL Style 1061 300 V T 80°C 1x26 AWG VW 1 Ø max 1.02mm Std connector: molex 51021-0300 with 3 poles	Transformer complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T1.6AL 250V Minimum requirements: • Output: 17 - 24 Vac; • Power: 26 VA; • Class B; • Rigidity over 4000 V. • Thermal protection	Component built-in	The medical system must be declared compliant with IEC/EN60601-1 by the installation technician or manufacturer. Note for the installation technician: ensure the combined version on which the light is installed is certified to host the complete light.
22-33Vdc		Power supply unit complies with IEC/EN 60601-1 third edition and IEC/EN 60601-1-2 with thermal protection or protection downstream with at least one appropriate fuse: • T630mAL 250V Minimum requirements: • Output: 22-33 Vdc • Power: 14 VA; • Class B; • Rigidity over 4000 V. • Continuous protection from short circuit or overcurrent		

4. ISTRUCTION FOR USE

Read paragraph 1 carefully for safe use of the device.

The device must be clean before use (see Device Cleaning paragraph).

	Attention
	Simultaneous use of the light with electro-surgical scalpels can cause its malfunctioning.
	Attention
	The control joystick must be handled with care to avoid breakages. Never move the light using the joystick to grip.
	Note
	Each time you switch on the light, the light intensity will be saved on previous switch off.
	Warning - danger of contact with powered parts
	do not use the device if parts or casing are damaged.

4.1 SWITCH ON AND OFF

Refer to §1.1 for the switch on and adjustment symbols.

Complete Light

To switch on and off the light, press and release the joystick lever using the right or left side. The light intensity on switch on will always be the last used before switch off.

Complete light with Theia Tech

Same operations of the complete light plus those on the fixed arm switch on and/or off in sync with that of the head.

The light on the fixed arm can be switched on/off even using the button on the arm. If the light is switch on, the light will be automatically synchronised, if the light on the head is off, it will be regulated to maximum intensity.

4.1.1 Adjustment:

a) to reduce light intensity keep the lever pressed on the joystick on the left side (light rear view) until the desired intensity is reached.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

a) to increase light intensity keep the lever pressed on the joystick on the right side (light rear view) until the desired intensity is reached.

On reaching the maximum intensity, you will hear an acoustic signal (1 beep).

c) to skip to maximum intensity, press and release the joystick lever using the front or rear side. On subsequent pressing on the front or rear side, the light intensity will return to that previously saved.

The light on the fixed arm is adjusted in sync with that of the head, it cannot be adjusted independently.

4.1.2 Light / Head WITH PROXIMITY

Switch on/off

Complete Light: To switch on or off, approach the sensor once up to a maximum distance of 3 cm.

Complete light with Theia Tech: same operations of the complete light plus the light on the fixed arm switch on and/or off in sync with that of the head.

Adjustment

Complete light: To adjust light intensity, you need to stay still near the sensor until the desired intensity is obtained. Adjustment enables passage from the maximum value to the minimum value again to the maximum. On reaching the maximum intensity, you will hear an acoustic signal (2 beeps). On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

	Note
	Each time you switch on the light, the light intensity will be saved on previous switch off.

Complete light with Theia Tech same operations of the complete light plus the light on the fixed arm will adjust in sync with that of the head.

4.1.3 Light / Complete light with Theia Tech / “ALYA” head WITH REMOTE CONTROL

Switch on / Switch off / Adjustment

- To switch on and/off press and release button “A”.

- Adjustment:

a) to reduce the light intensity keep button “A” pressed until reaching the desired intensity.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

b) to increase the light intensity keep button “A” pressed until reaching the desired intensity.

On reaching the maximum intensity, you will hear an acoustic signal (2 beeps).

c) to immediately reach the minimum light intensity, press button “B”.

On reaching the minimum intensity, you will hear an acoustic signal (1 beep).

Subsequent button pressing will bring the light back to the previous light intensity selected.

	Note
	Each time you switch on the light, the light intensity will be saved on previous switch off.

4.1.4 LAMP / COMPLETE LAMP WITH THEIA TECH / “ALYA” HEAD WITH SYNCHRONIZATION COMMAND

Where foreseen, it is possible to connect the Alya lamp to the Faro environment lamp in wireless mode in order to turn them into a synchronized lighting system named “Synchro”.

The “Synchro” mode has been purposefully designed to improve the comfort of the medical dentist/orthodontist, in order to reduce the effect of glare generated when passing from observation of a highly illuminated surface (e.g. oral cavity with a dental lamp) to a poorly lit surface (e.g. dental quiver).

With the method named “Synchro”, activatable using the button situated at the head of the Alya lamp, it is possible to change the illumination produced by the Faro environment lamp automatically on the basis of the value of illumination produced by Alya.

Note: Between the dental lamp and the environmental one, a small delay in synchronization may be experienced. This is due to the communications protocol. This effect is normal and does not signify a defect.

The “Synchro” function, in order to be enabled, requires a “Pairing” process (which should only be carried out once) in order to create a link between the two lamps. Subsequently, the “Synchro” function can be enabled and/or disabled as required by the user, from the button situated on the dental lamp.

“PAIRING” PROCEDURE

NOTE:

- The “Pairing” procedure is only necessary at the first connection, although it can be repeated in the event of changing the Alya lamp or the electronics of one of the two interconnected lamps.

- If there should be more environment lamps in the surgery, ensure that the other lamps are turned off or have been on for more than 60 seconds.

To carry out the **“Pairing”**, proceed as follows:

1. Apply current to the Faro environment lamp you wish to couple.
The environment lamp prepares itself for a Pairing connection, for a maximum time of 60 seconds.
2. Within the 60 seconds, press the “synchro” button on the dental lamp for at least 3 seconds, but for no more than 6 seconds, otherwise the procedure will be cancelled.
On receipt of the **“Pairing”** request by the dental lamp, the environment lamp activates the blue LED on the aluminium body. If the blue LED should not activate, it is possible (within 60 seconds from switching on) to carry out other attempts. Once the time has been exceeded, it is necessary to repeat the procedure from point 1.
3. From the blue LED illuminating on the environment lamp, there are 60 seconds to confirm the **“Pairing”** by pressing the programming button  situated on the remote control of the environment lamp. At this point, the blue LED (of the environment lamp) flashes at double speed and then goes out. If the button  on the remote controls is not pressed within 60 seconds, the blue LED goes out and the procedure must be repeated from point 1.

Having carried out the **“Pairing”**, synchronization between the 2 lamps is now enabled.

To **ACTIVATE THE SYNCHRONIZATION FUNCTION**, proceed as follows:

Press the Synchro button for 2 seconds then release it.

On release, a sound signal (Beep) will be heard and the blue LED on the environment lamp will light to indicate that synchronization has been activated.

To **DE-ACTIVATE THE SYNCHRONIZATION FUNCTION**, proceed as follows:

Press the Synchro button for 2 seconds then release it.


On release, a sound signal (Beep) will be heard and the blue LED on the environment lamp will go out to indicate that synchronization has been disabled.

Note connected with synchronization:

- When the Faro environment lamp is synchronized with the Alya lamp (or it adjusts itself automatically), the blue LED on the body is lit and steady. When the LED is off, the synchronization is not active.
- The remote control is always enabled, so it is therefore possible to change the illumination value. Nevertheless, if the environment lamp should find itself in a synchronized state (with the blue LED lit), as soon as a new adjustment is carried out on the Alya lamp, the illumination value will immediately be updated.
- If the Alya lamp should be turned off, the environment lamp will remain lit at the illumination value in use.
- If the Alya lamp should be turned on, the environment lamp will turn on automatically.

4.2 JOYSTICK ALYA JOYSTICK LEVER ASSEMBLY

- Fully insert lever “A” on the pin of the joystick.
- The hole in lever “A” must be positioned corresponding to surface “B”.
- Fully screw in grub screw “C” with the wrench supplied.

	Attention
	The control joystick must be handled with care to avoid breakages. Never move the light using the joystick to grip.

5. ORDINARY MAINTENANCE

There are no ordinary maintenance operations.

6. CLEANING

 	Warning against danger of wear and corrosion and falling suspended weights
	For all metal or plastic parts of the light, it is strictly forbidden to use substances that are abrasive, acid, substances containing chlorine or chloride ions, detergents with Trilene base, petrol, white spirit or similar. It is forbidden to directly spray any chemical substance on the device.

6.1 PARABOLA CLEANING

Cleaning must be carried out using a soft cloth in cotton or absorbent cotton with ethyl alcohol or the specific PERFLEX detergent. Water-alcohol based disinfectants are suitable with 70% isopropyl alcohol or ethanol.

 	Attention - potential damage or wear on the parabolas
	Never spray detergent directly on the parabolas. Cleaning operations on the parabolas must be carried out wearing gloves, to avoid leaving fingerprints on the surfaces. Never use detergents containing surfactants or water-repellents that depositing can leave streaks. Slight streaking will not prejudice the quality of the light. Products differing from those suggested could damage the parabolas. If in doubt, contact FARO customer care.

6.2 CLEANING THE HEAD

Cleaning must be carried out using a soft cloth in cotton or absorbent cotton with ethyl alcohol or the specific PERFLEX detergent. Water-alcohol based disinfectants are suitable with 70% isopropyl alcohol or ethanol.

 	<p>Warning against danger of wear of the plastic and falling suspended weights</p> <p>Never spray detergent directly on the head. Do not use detergents-disinfectants containing the following to clean plastic parts:</p> <ul style="list-style-type: none"> • AMMONIUM HYDROXIDE • SODIUM HYDROXIDE • METHYLENE CLORIDE • ME THYL ALCOHOL. • ALL KINDS OF ACIDS <p>Faro has tested and suggests the following disinfectants: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive</p>

6.3 ARM CLEANING

Always use a cloth soaked in disinfectant approved to disinfect the surfaces and pass it over.


 	<p>Warning against danger of corrosion and mechanical collapse with falling suspended weights</p> <p>Never spray chemical substances directly on the arms or joints and their openings.</p>
	<p>Warning against danger of wear of the plastic with falling suspended weights</p> <p>Donotusedetergents-disinfectantscontainingthefollowingtocleanplasticparts:</p> <ul style="list-style-type: none"> • AMMONIUM HYDROXIDE • SODIUM HYDROXIDE • METHYLENE CLORIDE • ME THYL ALCOHOL. • ALL KINDS OF ACIDS <p>Faro has tested and suggests the following disinfectants: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive</p>
 	

7. STERILISING THE HANDLES

	<p>Warning - danger of cross contamination</p> <p>The handles are not supplied sterile, they must therefore be sterilised before use. The handles must be sterilised before each patient.</p>

7.1 Removing the handles

To remove the handle, unscrew knob “A” and remove it from the support.

7.2 Decontamination and disinfectant

Before sterilising the handles, they must be decontaminated and disinfected.

To disinfect, Faro has tested the following products for disinfection: Eco Jet-1 (Cattani Group) / Sporekin Plus DS (Ims srl) / Zefirol Quick (Molteni Dental) / Durr FD366 Sensitive

 	<p>Attention - danger of plastic breaking</p> <p>The handles cannot be disinfected by thermo-disinfection.</p>

7.3 Sterilisation

The handles must be packaged in compliance with EN 868-5.

The handles can be sterilised with standard cycles 121°/134° C up to two hundred (200) cycles or however up to loss of the mechanical performance.

The parameters of the sterilisation cycle are as follows:

Ciclo EN 13060	Temperature	Pressure	Minimum Holding Time
B	121°C	207 KPa	15 min
B	134°C	308 KPa	3 min

8. PERIODIC CHECKS

Operation	Frequency	Applicability		Procedure	Enabling
		LD	TE		
Check the absence of range between the arm joints	Annual	x	N/A	Check the gap between the joints 5 and the arm has not changed since first use	User
Check data plate is legible	Annual	x	x		User
Check casing is intact	Twice-yearly				Service Engineer
Check electrical safety EN 62353 1. Rigidity 2. Dispersion	Twice-yearly	x	x	Measure dielectric rigidity and dispersion of the casing. Limits defined in IEC 60601-1	Service Engineer
Light parameter checks.	Twice-yearly	x	x	With a spectroradiometer, measure the values of: • Maximum lighting: >35000 lux • CRI decay: <20%. • Underlying value of the blue light on the spectrum emitted measured in: <100 W/m ²	Service Engineer

Service Engineer: competent person in electro-mechanical equipment maintenance.

9. ACOUSTIC ALARM SIGNALS

9.1 Acoustic Signals

OpL** = Beep 30 seconds

OTP* = Beep 30 seconds

* OTP: LED overtemperature protection.

** OpL: Disconnected led load

9.2 GUIDE TO PROBLEMS

The table below represents a guide to the potential defects of the light.

If you cannot resolve the problems, call technical support

Effect	Cause	Action (Service Engineer - SE)	Resp
The light does not switch on	Power supply not inserted or inserted incorrectly.	Check the power supply is inserted and the combined unit is on.	User
	Interference with electro-surgical scalpels or high energy tools	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
	Command on joystick applied incorrectly.	To switch on and off the light, press and release the joystick lever using the right or left side.	User
The light will flicker	Interference with electro-surgical scalpels or high energy tools.	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
The light does not adjust light intensity	Command on joystick applied incorrectly.	Use the command correctly as described in this manual.	User
	Interference with electro-surgical scalpels or high energy tools.	Switch off the electro-surgical scalpels and check the permanence of the effect.	User
The light intensity is considerably reduced	Parabolas or secondary lens dirty.	Clean the parabolas and the secondary lenses.	User
	Use of wrong procedures.	Check you have maximum adjustment with the command.	User
Stains have appeared on the parabolas or the reflective layer has come away.	Use of non-approved products.	Clean the surfaces with specific "Faro Perflex" product. Clean the surfaces with isopropyl alcohol. To restore the surfaces, you need to get service to replace the parabolas.	User
The light does not stay balanced and tends to lower	Excessive load on head (small mirrors, cameras, etc..).	Remove the excess loads.	User
The lamp does not control	Synchro function off	Activate the function; see 4.1.4	User

10. TECHNICAL FEATURES

Complete light:

Power supply voltage (without transformer):

- 17÷24Vac ±10% - 50/ 60Hz;
- 22÷35Vdc ±10%

Absorbed power:

- 26VA (version 17÷24Vac);
- 14VA (version 22÷35Vdc)

Recommended fuses:

- Version 17÷24Vac: T1.6AL 250V
- Version 22÷35Vdc: T630mAL 250V

Protection against electrical danger:

- Class II device

Classification EN 62471:

- Exempt class

Complete light with Theia Tech

Power supply voltage (without transformer):

- 24Vac ±10% - 50/ 60 Hz
- 24Vdc ±10%

Absorbed power:

- 40VA (version 24Vac)
- 28VA (version 24Vdc)

Recommended fuses:

- T2AL 250V

Protection against electrical danger:

- Class II device

Classification EN 62471:

- Exempt class

Optical features of the light produced by the head in compliance with ISO 9680

Spot light dimensions: 180 mm x 90 mm

Lux: 3.000*-50.000* lux @700mm

Colour temperature: 5000 K*


CRI (colour rendering index): >95*

* Typical values subject to tolerances.


Labelling in compliance with EN 62471: not necessary

10.1 WIRING DIAGRAMS


- **Complete light:**
Electrical circuit in light head

Wiring Diagram – Alya with transformer


- Complete light with Theia Tech:
Electrical circuit in light head

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
GUARANTEE CERTIFICATE

FARO offers the final customer a **24 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser. **Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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FARO S.p.A.
FARO SpA si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.
FARO SpA reserves the right to change the specifications of this equipment without notice.
FARO SpA se reserve le droit de modifier, sans préavis, les caractéristiques dans ce manuel.
FARO SpA behält sich recht vor, jederzeit stillschweigend technische oder bauliche Änderung vorzunehmen.
FARO SpA se reserva el derecho de modificar sin aviso previo la características incluidas en el presente manual de uso.