

- Founded in 2007
- HQ in Vienna, Austria
- Lexogen, Inc., NH, the USA
- Multinational team
- Products & Services

Developing **innovative technologies** that will allow to **resolve the complexity** of the **transcriptome**



100%


RNA Sequencing Tools

www.lexogen.com

Enabling Complete Transcriptome Sequencing



RNA-Seq project pipeline



Challenges in RNA-Seq library prep



- Amount of input material



- RNA quality

- High Quality
- Low Quality (FFPE)



- Ribosomal RNA elimination



- Strand-specificity



Challenges in RNA-Seq library prep



- Reproducibility



All Lexogen Kits

- Time of the library prep




- Library complexity




- Bioinformatics



PRODUCT PORTFOLIO IN THE RNA-SEQ WORKFLOW



Product portfolio in the RNA-Seq workflow



Lexogen's solution for gene expression profiling



**QUANTTM
SEQ**
Sequencing that counts

**High-throughput version
(384 samples, 4x96 barcodes)**



Two major types of RNA-Seq – which strategy?



Whole Transcriptome Sequencing

- mRNA or total RNA-Seq
- *de novo* assembly of transcripts
- Isoform detection
- Expression profiling (RPKM needed)



Assembly and FPKM



Multiple fragments per transcript

Expression Profiling Sequencing

- mRNA
- Counting the number of sequences
- Sequencing one fragment per transcript



Counting



One fragment per transcript


Which RNA-Seq Experiment Strategy?



Whole Transcriptome Sequencing



Assembly and FPKM



Multiple fragments per transcript

Expression Profiling Sequencing



Counting



One fragment per transcript




Only 1 read for 1 transcript is enough → higher multiplexing is possible, less computational effort → cheaper and faster


2 types of QuantSeq



- QuantSeq 3' mRNA-Seq Library Prep Kit
- QuantSeq-Flex Targeted mRNA-Seq Library Prep Kit



QuantSeq Workflow



**From 100 pg total RNA input
4.5 hrs library prep time, 2 hrs hands-on time
barcodes available: 96x96 (all unique)**



Cat.No. 015 (FWD):

Gene expression analysis



- NGS reads will be generated **towards** the **poly(A) tail**
- Read 1 sequencing linker is located in **5' part of the second strand synthesis primer**

Cat.No. 016 (REV):

3'UTR analysis



- NGS reads start **directly at the 3' end** of transcripts
- Read 1 sequencing linker is located on the **5' end of the oligodT primer**
- **Custom Sequencing Primer (CSP)** covers poly(T) stretch and replaces Multiplexing Read 1 Sequencing Primer - Information must be provided to sequencing facility along with the CSP (see UG for more information)!
- Recommended for **PE sequencing**

Strand-specific mapping of 3' ends



- ✓ QuantSeq covers the very 3' end and saves more than 90 % sequencing depth while still determining gene expression accurately.






Figure | Coverage versus normalized transcript length in QuantSeq and standard mRNA-Seq. RSeQC-derived coverage is plotted for all transcript (areas) and the ERCC mix only (lines), for QuantSeq (green) and mRNA-Seq (grey). Numbers give the area under the curve (AUC) values as a measure for sequence coverage.

Required number of sequencing reads




1) Liu Y., et al., RNA-seq differential expression studies: more sequence or more replication? *Bioinformatics* 30(3):301-304 (2014) 2) Liu Y., et al., Evaluating the impact of sequencing depth on transcriptome profiling in human adipose. *Plos One* 8(6):e66883 (2013) 3) Bentley, D. R. et al. Accurate whole human genome sequencing using reversible terminator chemistry. *Nature* 456, 53–59 (2008) 4) Rozowsky, J. et al., PeakSeq enables systematic scoring of ChIP-seq experiments relative to controls. *Nature Biotech.* 27, 65-75 (2009).


QuantSeq quantifications of FFPE and fresh frozen (cryo) samples correlate well



Xenograft of MOLP-8 human tumor cell line, split and processed either as fresh frozen cryo-block (RIN 8.3) or FFPE material (RIN 2.8, DV₂₀₀ of 87 %).



QuantSeq outperforms Microarrays in the detection of low and medium abundant transcripts



Minimum detection value



- **Illumina HT12 array:** quantile normalized signal levels (Yu et al.)
- **QuantSeq:** counts based on total gene reads

Yu, J. et al. (2015). *BMC Genomics*. 16:710



Alignment:

- Read mapping is simplified by skipping the junction detection.



Concentration calculation:

- QuantSeq is streamlined for gene-read counting. No calculation of FPKM/RPKM values are needed.
- High strand-specificity reduces ambiguous counting.

Streamlined data analysis pipeline – FREE for QuantSeq users on Bluebee platform




www.bluebee.com/quantseq/

- **User-friendly data analysis** on the **Bluebee** genomics analysis platform
- For any user, also **without bioinformatics background**
- **Free for QuantSeq customers** – code is included in the kit
- **Data analysis only takes 10 minutes** once set up



Data analysis workflow



Five simple steps to have your data analyzed



Register

on the Bluebee platform using the kit code provided with the QuantSeq kit (FWD or REV for Illumina, single read runs only)

Get connected

Connect your sequencer with the Bluebee platform (using the Bluebee Service Connector)

Upload your data

zipped fastq file = sequencing data file

Run your pipeline



Choose QuantSeq FWD or REV for Illumina pipeline

Get your results

Mapping statistics and read counts
DE analysis (list of DE genes, heat map, PCA plot and MA plot)
FastQC, RSeQC

Remark: For a detailed walk through of each of these steps, please have closer look at the **video-tutorials**.

Bluebee data analysis – overview



For more information and support please contact us at bioinfo@lexogen.com

QuantSeq: Cost Calculations



| | | QuantSeq | mRNA-Seq |
|--------------------------------|----------------|---------------------------|-------------------|
| Libraries / lane | | 96 (3M reads) | 15 (20M reads) |
| Cost / rxn | | USD 19.80 | USD 51.00 |
| Total library prep costs | | USD 1900.80 | USD 765.00 |
| 1 seq. lane (300M SR100 reads) | | USD 1500.00 | USD 1500.00 |
| Data evaluation | | included @ BlueBee | significant |
| RNA to data | total | USD 3400.80 | USD 2265 |
| | per rxn | USD 35.43 | USD 151.00 |

**RNA-seq based expression profiling
including sequencing and data evaluation**

~36 USD

per sample

QuantSeq versus standard mRNA-Seq



| | QuantSeq | mRNA-Seq ¹ |
|---|---------------|-----------------------|
| Overall protocol time | 4 h 30' | 12 h |
| Hands-on time | 1h 45' | 4h 30' |
| Computing Time ² | 35' | 2h 50' |
| Total RNA input | 0.1 – 2000 ng | 100 – 1000 ng |
| Input-output correlation (R^2) ³ | 0.97 | 0.82 |
| Strandedness ³ | > 99.9 % | 93.4 % - 97.8 % |
| Differential expression (AUC) ⁴ | 0.84 – 0.90 | 0.74 – 0.78 |

¹ TruSeq mRNA stranded library preparation (Illumina)

² mRNA-Seq needs a splice-junction mapper such as TopHat2, QuantSeq only a short-read aligner such as Bowtie2


³ Calculated on ERCC-mapping reads

⁴ Area under curve (see slide “accurate gene quantification even at low read depths”)

2 types of QuantSeq




- QuantSeq 3' mRNA-Seq Library prep kit
- QuantSeq-Flex Targeted mRNA-Seq Library prep kit



RNA-Seq experiments tailored to your needs



| Application | RT | SSS |
|--------------------------------------|----------|----------|
| QuantSeq 3' mRNA-Seq | Oligo-dT | Random |
| Targeted 3' mRNA-Seq | Oligo-dT | Targeted |
| Targeted RNA-Seq, novel fusions | Targeted | Random |
| Targeted RNA-Seq, known targets only | Targeted | Targeted |





- Highly efficient and very cost-efficient solution for expression profiling (Alternative to microarrays, qPCR and standard RNA-Seq)
- Fast (only 4.5 hours) and simple all-in one protocol (fewer steps = less chance of error)
- Suitable for degraded RNA samples (e.g., FFPE)
- Wide input range:
 - For Illumina: 100 pg – 2 µg of total RNA
 - For Ion Torrent: 5 ng – 2 µg of total RNA
- High level of multiplexing (up to 384 samples/lane)
- > 99.9 % strand-specificity
- Optimized for shorter reads (SR50, 100). Library Sizes around 335-456 bp with insert sizes of 203-324 bp
- Usage of custom primers for targeted RNA-Seq (QuantSeq-FLEX)
- Automation-friendly (autoQuantSeq for Illumina)
 - For Sciclone NGS and Zephyr liquid handlers of PerkinElmer and the Hamilton Microlab STAR Workstations, Agilent Bravo, Beckman.



The effector AWR5 from the plant pathogen *Ralstonia solanacearum* is an inhibitor of the TOR signalling pathway

Crina Popa, Liang Li, Sergio Gil, Laura Tatjer, Keisuke Hashii, Mitsuaki Tabuchi, Núria S. Coll, Joaquín Ariño & Marc Valls

[Sci. Rep. 6, 27058; doi: 10.1038/srep27058 \(2016\)](#)

Transcriptional profiling and muscle cross-section analysis reveal signs of ischemia reperfusion injury following total knee arthroplasty with tourniquet

Jonathan B. Muyskens, Austin D. Hocker, Douglas W. Turnbull, Steven N. Shah, Brick A. Lantz, Brian A. Jewett, Hans C. Dreyer

[Physiological Reports Published 5 January 2016 Vol. 4 no. e12671 DOI: 10.14814/phy2.12671](#)

Find the latest QuantSeq publications at <https://www.lexogen.com/publications/#quantseqpub>

Quote from paper in Nature: Dysregulation of expression correlates with rare-allele burden and fitness loss in maize

[Karl A. G. Kremling¹, Shu-Yun Chen^{2,3}, Mei-hsiu Su², nicholas K. Lepak⁴, M. Cinta Romay², Kelly L. Swarts^{1,5}, Fei Lu^{2,6}, Anne Lorant⁷, Peter J. Bradbury⁴ & Edward S. Buckler^{1,2,4}](#)

We automated a 3' mRNA sequencing method (**QuantSeq**, Lexogen GmbH), which is **more efficient and accurate than mRNA sequencing** and deals well with paralogues.“

Watch our QuantSeq webinars!



WEBINAR:
Integration of 3' mRNA-Seq and iCLIP to derive high-resolution RNA maps for the regulation of alternative polyadenylation



Dr. Jernej Ule (Francis Crick Institute)
Dr. Gregor Rot (Institute of Molecular Life Sciences University of Zurich)

<https://www.lexogen.com/webinar-expressrna/>

WEBINAR:
Analysis of the transcriptome of carriers of pathological variants in PSEN1, PSEN2 and APP that cause Alzheimer's Disease



Oscar Harari (Department of Psychiatry at Washington University in St Louis)

<https://www.lexogen.com/webinar-analysis-of-the-transcriptome/>

WEBINAR



Free QuantSeq data analysis pipeline on the Bluebee platform

Jekaterina Aleksejeva (Lexogen)

<https://www.lexogen.com/webinar-free-quantseq-data-analysis/>

WEBINAR:
Mapping nuclear-exosome targeted poly(A) tails with 3'-RNA seq



Kevin Roy (Stanford University School of Medicine)

<https://www.lexogen.com/webinar-mapping-nuclear-exosome-targeted-polya-tails-with-3-rna-seq/>

WEBINAR:
Gene Expression Analysis Using 3'-RNA Sequencing



Behnam Abasht (University of Delaware)

<https://www.lexogen.com/webinar-gene-expression-analysis-using-3-rna-sequencing/>



**QUANT™
SEQ**
Sequencing that counts

Automation



QuantSeq has been automated on various liquid handling instruments







- ✓ Perkin Elmer Sciclone / Zephyr
- ✓ Hamilton Microlab STAR
- ✓ Beckmann Biomek FxP
- ✓ Agilent Bravo NGS Workstation (B) – Application Note available
- ✓ Eppendorf EpMotion 5075
 - Script files are available at www.lexogen.com
 - **Dummy Reagents** for QuantSeq Automation (Cat. No. 019.384) – free of charge
 - Set of 3 **mock solutions** for **volume checks** and **consistency**
 - **Purification Module** with Magnetic Beads (Cat. No. 022.24)



Efficient, robust, and reproducible library preparation in a **high-throughput format**, with minimal drop out.


QuantSeq Automation on Agilent Bravo shows high consistency for technical replicates





Globin Block Modules






- ✓ Informative
- ✓ low cost
- ✓ minimally invasive
- ✓ highly accessible
- ✓ biomarker discovery and monitoring

- ⊖ high abundance of globin mRNA (HBA, HBB)
- ⊖ 50-80 % of sequencing reads
- ⊖ limiting gene detection and quantification sensitivity

Seamless gene expression profiling of blood samples



- For **human** (*Homo sapiens*, Cat. No. 070.96) and **pig** (*Sus scrofa*, Cat. No. 071.96) **blood samples**
- Efficient **globin depletion** for **less than USD 4.98** per sample
- Block globin mRNA **during QuantSeq** library prep – **no additional steps** required, simple exchange of RNA Removal Solutions (RS vs. RS-GBs)
- **Input starting from 50 ng of total RNA**
- **Reduce globin mapping reads by up to 91 %**
- **Detect and quantify thousands of additional genes**




Efficiently deplete globin mRNA for less than \$ 4.98 per sample!


Efficiently deplete globin mRNA and detect and quantify thousands of additional genes



Globin mRNAs reads reduced from 50 – 80 %, to as low as 5 %



Gene detection is increased by thousands of additional genes





i5 Unique Dual Indexing Kits

**QUANT™
SEQ** **SENSE™**
Sequencing that counts Making sense of RNA sequencing



i5 Unique Dual Indexing Add-on Kit – Your benefits



Enhanced multiplexing capacity

Benefit from up to 9,216 different i5/i7 index combinations, or up to 96 unique dual indexing barcoding options for all QuantSeq and SENSE libraries



Minimize index mis-assignment

Up to 96 unique dual indexing barcode options enable you to detect and quantify index hopping




Index Balance Checker

Use our online help tool to choose the optimal index combination for your RNA-Seq experiment



Problem: Index hopping/switching

= **mis-assignment** that results in the incorrect assignment of libraries from the expected index to a different index (in the multiplexed tool)




Resource: Illumina Whitepaper 770-2017-004-D (2017) *Effects of Index Misassignment on Multiplexing and Downstream Analysis. Learn why it happens and best practices to reduce the impact of index hopping.*

Reason: Index hopping

- **Left-over primers** causing issues with the ExAmp (exclusion amplification chemistry)
- **Patterned flow cells:** HiSeq3000/4000, NovaSeq
- **Slightly elevated levels: 0.1-2 %** depending on type, quality, and handling library

Solution

- **Thorough purification protocols** to remove free adapter (always **two rounds** of post PCR, optional repurification of lane mix)
- **96 i5 x 96 i7 unique dual indices** 




Unique Molecular Identifiers

**QUANTTM
SEQ**
Sequencing that counts




UMI Second Strand Synthesis Module – Your benefits



De-duplicate your RNA-Seq data

Tag individual transcripts to detect and eliminate amplification bias for unbiased gene expression profiling



Fully integrated within QuantSeq

Save your precious time and add UMIs during QuantSeq library prep with a simple solution exchange, no additional steps required



Data analysis tool available

De-duplicate your QuantSeq data using the *collapse_UMI_bam* command-line analysis tool, available from Lexogen

UMI working principle

UMI = Unique Molecular Identifier, also known as Molecular Barcode




Library prep **without** UMI

Problem: Sequence specific PCR bias:

- Pre-PCR: identical cDNA fragment pools (cDNA1 and cDNA2) with the same size
- Post-PCR: cDNA1 and cDNA2 sizes differ

Inaccurate gene and isoform quantification




Library prep **with** UMI

Solution to sequence specific PCR bias:

- Ligation of random oligonucleotides (UMIs) to DNA fragments before PCR
- Distinguish pre-PCR fragments with identical sequence
- Number of distinct UMIs after PCR is the same as the number of fragment copies before PCR

Unbiased gene expression profiling



UMI Second Strand Synthesis Module for QuantSeq FWD (Illumina, Read 1) (Cat. No. 081.96)



- **Kit components:** UMI Second Strand Synthesis Mix (USS)
- **Compatibility:**
 - only for QuantSeq FWD (Cat. No. 015)
 - **NOT** compatible with QuantSeq REV (Cat. No. 016) or Ion Torrent (Cat. No. 012)
- **Design:** random primers including 6 bp UMIs. UMI sequence is located between the partial P5 adapter and the random priming sequence.
- **Usage:**
 - simple solution exchange
 - UMIs added during the second strand synthesis step of QuantSeq
 - USS replaces Second Strand Synthesis Mix 1 (SS1) from standard QuantSeq FWD kit at step 7

Kit Components and Storage Conditions




Figure 2. Location of kit component.

| UMI Second Strand Synthesis Module for QuantSeq FWD (Illumina, Read 1) (Cat. No. 081.96) | Tube Label | Volume* 96 rxn | Storage |
|--|------------|-------------------|---------|
| UMI Second Strand Synthesis Mix | USS ● | 1056 µl | -20 °C |

*including 10% surplus



Second Strand Synthesis

NOTE: This protocol replaces steps 7 and 8 of the detailed protocol from the QuantSeq 3' mRNA-Seq Library Prep Kit for Illumina User Guide (015UG009, 015UG110). Step 8 has not been changed for UMI libraries and is included here for ease of reference.

Follow steps 1 - 6 as indicated in the QuantSeq 3' mRNA-Seq Library Prep Kit for Illumina User Guides (see above).

7 Add 10 µl of UMI Second Strand Synthesis Mix (USS ●) to the reaction. Mix well by pipetting, and seal the plate. **REMARK:** Use a pipette set to 30 µl for efficient mixing.

8 Incubate the plate for 1 minute at 98 °C in a thermocycler, and slowly cool down to 25 °C by setting the ramp speed to 10% (0.5 °C/second). Incubate the reaction for 30 minutes at 25 °C. Quickly spin down the plate at room temperature before removing the sealing foil.

➔ Proceed to step 9 of the detailed protocol in the QuantSeq 3' mRNA-Seq Library Prep Kit for Illumina User Guides (see above).

User Information will be added in the **QuantSeq User Guide**.




- **Tag individual transcripts during QuantSeq** library preparation by a simple solution exchange
- **Identify PCR amplification bias** and **de-duplicate** your **sequencing data** for unbiased gene expression profiling
- **6 bp UMIs** added **during second strand synthesis** – no additional steps, same streamlined QuantSeq workflow
- **Single-read sequencing** (SR 75 or higher) enough – UMI read-out at the start of Read 1
- Compatible with **QuantSeq 3' mRNA-Seq (FWD)** for Illumina
- Compatible with **Globin Block Module** and **Unique Dual Indexing Module** – use for single and dual-indexed libraries
- **Data Analysis Tool** ***collapse_UMI_bam*** for collapsing reads by UMIs available from Lexogen



Full Service



Fully integrated service-workflow: from RNA to evaluated data



- Fully integrated package: **economical and convenient**
- **Includes RNA QC, library preparations, NGS run, and data analysis**
NGS data QC, mapping, counting, (differential) gene expression analysis
- **Project management and support** by the experts in the field


Pricing depends on the number of samples to be processed, the desired degree of multiplexing, and the sequencing mode. For a quote or an example report, please inquire at services@lexogen.com

Now available: QuantSeq SERVICE



- ✓ **From RNA to evaluated data**
- ✓ Fully integrated package: economical, and convenient
- ✓ **Includes RNA QC, library preparations, NGS run (SR75), and data analysis**
- ✓ NGS data QC, mapping, counting, (differential) gene expression analysis
- ✓ **Pricing - depending on number of samples**
- ✓ **Project management and support** by the experts in the field

Product portfolio in the RNA-Seq workflow





SLAMseq

Metabolic RNA-Seq



The transcriptome is not completely described by measuring RNA steady state-levels




- Conventional RNA-Seq experiments measure RNA steady-state levels
- Metabolic RNA-Seq can measure RNA synthesis and degradation separately




Gain novel insights in gene expression and its controls!

RNA synthesis



RNA degradation





Herzog, V. A. et al. (2017) *Nature Methods* 10.1038/nmeth.4435

SLAMseq METABOLIC RNA LABELING KITS



Identify direct transcriptional targets of any gene





| Kit Type | Module | Application |
|----------------------|---|---|
| SLAMseq Explorer Kit | Cell Viability Titration Module (Cat. No. 059.24) | <ul style="list-style-type: none"> •Assess S4U toxicity in target cell lines •Optimize S4U labeling concentrations |
| | S4U Incorporation Module (Cat. No. 060.24) | <ul style="list-style-type: none"> •Measure S4U uptake and incorporation rates using HPLC analysis |
| SLAMseq Kinetics Kit | Anabolic Kinetics Module (Cat. No. 061.24) | <ul style="list-style-type: none"> •Label newly transcribed RNA with S4U •Measure nascent RNA expression •Analyze RNA synthesis kinetics |
| | Catabolic Kinetics Module (Cat. No. 062.24) | <ul style="list-style-type: none"> •Label existing RNA with S4U •Assess transcript stability •Analyze RNA degradation kinetics |

- **Kit size and format: 24-well cell culture plate** format using **24× 0.5 - 1 ml** growth medium, all **examples** in the **User Guide** are based on 24-well format. **Up- and down-scaling** to suit different culture vessel sizes is **possible**.
- Interested in bigger kits sizes? Please contact info@lexogen.com

SLAMseq Kinetics Kits – measure nascent RNA expression and transcript stability



Anabolic Kinetics Module – RNA syns = S4U-labeling of newly synthesized (nascent) RNA





Figure | Anabolic kinetics labeling experiment time course.

Catabolic Kinetics Module – RNA degradation rates = S4U-labeling of existing RNA







Figure | Catabolic kinetics labeling experiment time course.


QuantSeq 3' mRNA-Seq library preparation is ideal for SLAMseq samples



- **Cost-effective** library preparation
- **Low read depth:** 3 M reads for gene expression, 20-30 M reads for SLAMseq
- **Poly(A)-selective:** functional, fully processed RNA pol II transcript, no requirement for rRNA depletion
- **Facilitated data analysis:** highest strandedness, no length normalization
- **Sensitive:** from 100 pg total RNA input; standard input for SLAMseq is 1 μ g
- **Robust quantification of T > C conversions,** covers inherently U-rich 3' UTRs



Product portfolio in the RNA-Seq workflow





- Streamlined and user-friendly SLAMseq-QuantSeq data analysis integrated on the Bluebee genomics analysis platform <https://www.bluebee.com/lexogen>
- **DUNK: Digital Unmasking of Nucleotide conversion-containing k-mers**
- **Modified alignment algorithm** to report **read counts** for **T>C containing reads**
- For any user – **no specialized bioinformatics knowledge required**
- **Highly secure** cloud-based solution
 - **Five simple steps** to have your **SLAMseq-QuantSeq data analyzed: Register, get connected, upload** compressed fastq files, **select** species and **run** pipeline.
 - **Output: statistics** about **T->C conversion rates** and **alignments** falling into **unique 3' UTR** regions (CPM, T coverage, read counts, multimapping, etc). In addition, SLAMdunk produces further statistical and diagnostic information. Results can be downloaded in batch or individually.
- **Complete metabolic RNA-Seq solution** – from **RNA labeling** to **NGS library preparation** and **data analysis**



Developed by the Zuber (IMP), the von Haesler (MFPL) and the Ameres (IMBA) groups


Herzog, V. A. et al. (2017) Nature Methods 10.1038/nmeth.4435

Neumann, Tobias, Herzog, V.A., Muhar, M., von Haeseler, A., Zuber, J., Ameres, S.L. & Rescheneder, P., unpublished data



- Analyze **transcriptome-wide kinetics of RNA synthesis and turnover**
- **Measure newly synthesized RNA expression and transcript stability**
- Gain **novel insights** into the **control of gene expression**
- Assess **toxicity and incorporation of 4-Thiouridine (S4U) in cultured cells**
- **Only two extra steps** added to a standard RNA-Seq workflow:
 - **Labeling of RNA by adding S4U** to the culture medium
 - **Pre-processing of the total RNA with iodoacetamide** to alkylate 4-thiol group
- **Combine** with [QuantSeq 3' mRNA-Seq](#) or [QuantSeq Flex Targeted RNA-Seq Kits](#) for **time- and cost-efficient, high-throughput metabolic sequencing** from **multiple timepoints** and in **replicates**
- **SLAMdunk: available as user-friendly data analysis pipeline** on Bluebee for **analyzing SLAMseq-QuantSeq sequencing data**
 - **Modified alignment algorithm** to report **read counts for T>C containing reads**
 - Herzog et al., *Thiol-linked alkylation of RNA to assess expression dynamics* (Nature Methods, 2017: [DOI:10.1038/nmeth.4435](https://doi.org/10.1038/nmeth.4435))


Product portfolio in the RNA-Seq workflow






SENSE™

Making sense of RNA sequencing



Strand-specific RNA sequencing is a crucial factor



About 25% of human genes share expressed mRNA regions and up to 10% of genes in the human genome participate in bi-directional promoters* .


* Engström, P.G. *et al.*, 2006, **Complex Loci in Human and Mouse Genomes**, *PLoS Genetics* V.2, p.0564-0577.

Stranded libraries **preserve information** about transcribed strand


- Detect **antisense** transcripts
- **Overlapping** genes
- **Fusion** genes

- Measure **gene expression**
- **Annotate** the structures of all transcribed genes
- Measure the extent of **alternative splicing**

RNA-Seq library prep methods




Lexogen's proprietary strand displacement stop/ligation technology




- 1 Random hybridization of starter/stopper heterodimers** including sequencing platform-compatible linker sequences, to the RNA
- 2 Extension** of the starter to the next hybridized heterodimer by reverse transcription
- 3 Ligation** of the newly synthesized cDNA insert to the stopper
- 4 Hydrolysis of the RNA during second strand synthesis** and conversion of the library to **double-stranded DNA**

Avoiding reverse transcription artifacts affecting strand-specificity



SENSE

- No fragmentation
- No new 5' ends
- No hairpin loop artifacts



SENSE

- Efficient stop of extension (of the Starter at the next Stopper)
- No strand displacement

*Perocchi, F., et al., 2007, Antisense artifacts in transcriptome microarray experiments are resolved by actinomycin D, Nucleic Acids Res. 35(19): e128


Highly reproducible strand-specificity



✓ Mean strand-specificity: 99.922% (ERCCs)

| Experimenter 1 | | Experimenter 2 | | Experimenter 3 | |
|----------------|-------------------|----------------|-----------------|----------------|-------------------|
| Sample | Strandedness, % | Sample | Strandedness, % | Sample | Strandedness, % |
| 12782 ACATTA | 99.927 | 12790 AATGAA | 99.981 | 12798 GATCAC | 99.975 |
| 12783 GGTGAG | 99.975 | 12791 GATTGT | 99.968 | 12799 CGCGGA | 99.948 |
| 12784 CGAAGG | 99.911 | 12792 ATAAGA | 99.947 | 12800 CCTAAG | 99.925 |
| 12785 AAGACA | 99.944 | 12793 GCCACA | 99.958 | 12801 GGCTGC | 99.910 |
| 12786 TAATCG | 99.943 | 12794 GAACCT | 99.893 | 12802 ACCAGT | 99.957 |
| 12787 CGCAAC | 99.841 | 12795 CGGTTA | 99.933 | 12803 GTGCCA | 99.810 |
| 12788 AATAGC | 99.693 min | 12796 AACGCC | 99.886 | 12804 AGATAG | 99.993 max |
| 12789 TTAACT | 99.949 | 12797 CAGATG | 99.864 | 12805 TCGAGG | 99.987 |

>99.9% strand-specificity ensures confident antisense transcript detection



| Strand specificity | Antisense transcripts detected |
|--------------------|--------------------------------|
| 99.99% | 2904 (100%) |
| 99.90% | 2415 (83%) |
| 99.00% | 1506 (52%) |
| 98.00% | 1153 (40%) |



CORALL™

NEW

Cover all bases





Excellent 5' representation

Get full transcript coverage starting at the very 5' end



Fast and easy library prep


Prepare libraries in only 5 hours and save research time



Unique Molecular Identifiers included



Tag individual transcripts to detect and eliminate amplification bias

CORALL WORKFLOW






- ✓ Surpassing gene detection levels, even at high TPM thresholds
- ✓ 96% of all genes detected by two other methods covered and additional 20% detected (TPM>10)
- ✓ 76% of exclusively detected genes are protein-coding genes



Product portfolio in the RNA-Seq workflow






Small RNA-Seq Kit



Lexogen's Small RNA-Seq Kit for Illumina



Exceptional miRNA discovery and high reproducibility




Figure | Total number of miRNAs detected across 4 different protocols. Dilution series (6 pg, 60 pg, and 600 pg) of purified plasma RNA used for library prep from different vendors were performed. The obtained libraries were sequenced at equal molarity at ~ 1.5 – 2M Total Raw Reads per sample. Lexogen’s Small RNA-Seq Kit showed much higher numbers of detected miRNAs at ≥ 5 Raw Reads across all concentrations tested.





Figure | a) Dilution series of purified plasma RNA in Small RNA-Seq shows high inter-replicate correlation Reads per Million (RPM) across concentrations tested (60 pg and 600 pg). **b)** Dilution series of purified plasma RNA in Small RNA-Seq shows high correlation RPM across concentrations tested (60 pg and 600 pg). As an example, only data for 60 pg and 600 pg are shown.



- **Discovery and profiling of small RNA** (≤ 200 nt), including miRNA, piRNA, siRNA, snoRNA, tsRNA, srRNA, U-RNA
- **Gel-free user-friendly workflow** (bundle version with Purification Module with Magnetic Beads available (Cat. No. 058))
- **Wide input range from 50 pg to 1 μ g of RNA**
(1 ng – 1,000 ng of cellular total RNA or 50 pg – 1,000 ng enriched small RNA including plasma, serum, and urine)
- **Optimized for low RNA content samples** such as **liquid biopsies** (e.g., plasma, serum, and urine) including **exosomes**
- **All-in-one** protocol – only **75 minutes hands-on** time
- **High reproducibility** for inter-replicates and across concentrations
(R^2 0.999 to 0.952)
- **Exceptional miRNA discovery** compared to other protocols, **in particular for low RNA inputs**
- **Combination with SPLIT RNA Extraction Kit** offers the complete solution for small RNA analysis

Product portfolio in the RNA-Seq workflow





SIRVs™

Spike-in RNA Variant Control Mixes



Do you use controls in your experiments?





Do you use spike in controls in your **RNA** experiments?




100%

Yes.



~ 100 USD

Do you use spike in controls in your **RNA-Seq** experiments?



Aug 2014, poll by RNA-seqblog.com

52% Yes.

7% Not yet.

41% No.



48 samples

your samples

Library preparation

950 -1'630 USD


125 PE, 300 Mreads

2'500 USD

Bioinformatics

substantial

> 3'450 USD





Validation and **Monitoring** of the RNA-seq workflow assesses

- **technical parameters** such as biases and sensitivity
- **functional aspects** like isoform detection and quantification



Know the **performance** and the **limits** of your RNA-seq workflow and data evaluation pipeline.


Concordance determines whether samples are comparable

- **in experiments** and
- **between experiments**



Know the **technical variability** between **samples** and **experiments** before investigating **biological variance**.

SIRV modules address transcriptome complexity in a defined and condensed manner



- **Transcriptome properties:** high complexity and large concentration differences, several RNA classes with specific properties
- **SIRVs as a family of modules:** each module probes a specific aspect
 - **SIRV isoform module:** transcription and splicing variants -> isoforms
 - **ERCC single isoform module** -> concentrations



Current SIRV genome: 7 SIRV isoform genes and 92 ERCC genes




Figure | SIRVome and SIRV transcripts. The SIRV isoform and ERCC genes are lined up on a SIRVome. Pullout, the Compact Coverage Visualization (CCV) shows intron regions common to all transcripts as short standardized gaps irrespective of the original sequence length, which provides a more comparable overview of the actual SIRV transcripts.

- **Sequences and gene definitions:** multi-FASTA and GTF files for the individual SIRV genes (isoforms and ERCCs) and the entire SIRVome are freely accessible at <https://www.lexogen.com/sirvs>
- **GenBank:** SIRV isoform sequences can be accessed using No KX147759 to KX147765 for SIRV1 to SIRV7



- **7 artificial genes** (derived from human model genes) with 6-18 transcript variants each, **69 *in-vitro* transcript variants** in total
- **Length 191 - 2528 nt, poly(A₃₀)-tail, GC-content 29.5 – 51.2 %**
- **Canonical GT-AG exon-intron junctions** with known exceptions
- **Intron sequences:** random, maintain GC content of upstream and downstream exons
- **No significant similarity** to any known sequences
- Three mixes **E0, E1, and E2** with molar ratios of transcripts at **magnitudes of 0, 1, and 2**, respectively
- Also available in **conjunction with ERCCs**

SIRVs mimic the complexity of the transcriptome



- Alternative splicing
- Promoter and poly(A) site usage
- Overlapping genes
- Antisense transcription






Figure | SIRVs design overview. SIRVs genes mimic human model genes to represent in their entirety all main aspects of alternative transcription in numerous repeats and variations. A5SS and A3SS, alternative 5'/3' splice sites, MXE, mutually exclusive exons.



- Developed by the External RNA Controls Consortium (ERCC)
- **Exemplary studies** by the **FDA Sequencing Quality Control (SEQC)** Consortium and the **Association of Biomolecular Resource Facilities (ABRF)**
- **92 artificial transcripts** with **non-overlapping sequences following the 1 gene, 1 exon, and 1 transcript layout**
- Unique sequence identities makes them well suited for measuring **technical parameters irrespective of isoforms**
- **Standardized assessment of gene expression platforms** such as quantitative **RT-PCR, microarrays, and NGS** technologies
 - **Dynamic range**
 - **Dose response**
 - **Lower limit of detection** and **efficiency**
 - **Fold change response** of RNA sequencing pipelines,within the complexity boundaries of monoexonic, non-overlapping RNA sequences.

Spiking of samples



SIRV
— SUITE —
**Experiment
Designer***


- Spike SIRVs into sample: **during or after RNA extraction**
- SIRVs have **exclusive sequences** – spike into **any sample or organism**
- SIRVs are **polyadenylated** – library prep can start from poly(A) selected fractions as well as from total RNA, depleted RNA, etc.
- **Spike-in RNA amount**: only a **small fraction** (usually 1 %) of all **NGS reads** mapping to the SIRV genome, the “SIRVome”.

* For SIRV isoforms

Spike-in RNA Variant Control Mixes (SIRV set 1)




- 3 predefined SIRV mixes: E0, E1, and E2
- Each mix contains all 69 SIRVs in different concentration ratios



SIRV-Set 2 and SIRV-Set 3 contain the isoform mix E0




- **SIRV-Set 2** (Cat. No. 050.xx) contains mix **E0 only**
- **SIRV-Set 3** (Cat. No. 051.xx) contains mix **E0** in mixture **with the ERCC** mix
- These sets may come in **dried format**.



- **A set of freely available software tools** accompanying SIRV isoforms
- **Streamlined and unified data evaluation process** embedded in the **Galaxy environment** (galaxyproject.org)
- **Evaluator report (pdf)** as **standardized documentation** of SIRV isoforms evaluation

- **Modules:**

- **Designer:** plan experiments
- **Specifier:** enter experiment information
- **Evaluator:** evaluate SIRV isoforms data and obtain a standardized but customizable report and data files for downstream applications
- **Comparator:** Cross-sample comparisons and referencing
- **Data Base:** experiment metadata, SIRV isoform reads, quality metrics



Webinar: <https://www.lexogen.com/controlling-rna-seq-experiments-using-spike-in-rna-variants/>



- **Spike-ins** to RNA-Seq experiments with **exclusive sequences** (platform and sample independent)
- **SIRV isoforms** provide **isoform complexity** through **mimicking transcription** and **alternative splicing** variations
- **ERCCs** provide **abundance complexity** through **high dynamic range** (10^6)
- **Assess** your pipeline for **transcript isoform detection** and **quantification performance**
- **Determine accuracy** and **precision** in concentration and differential expression measurements
- **Validation** of **workflows** and **experiments** under authentic conditions
- **Identify** the sources of **errors** and help **improving experimental workflows** from library preparation, over sequencing up to data evaluations
- Provide an universal marker for **comparing RNA-Seq experiments**
- **SIRV Suite: Streamlined** and **unified data evaluation process** embedded in **Galaxy**
- **SIRV-Set 2 and 3** as **dried format** – **easy delivery, storage** and **handling**

Watch our SIRVs webinar



SIRVs™
Spike-in RNA Variant Control Mixes

The only RNA controls designed
for splice variant detection

**WEBINAR: Controlling RNA-Seq Experiments Using
Spike-In RNA Variants**

Continuing Education Credits: P.A.C.E. CE | Florida CE

REGISTER NOW FOR ON DEMAND




Lukas Paul, PhD
Head of Services, Lexogen

<https://www.lexogen.com/controlling-rna-seq-experiments-using-spike-in-rna-variants/>

- The **concept of complex RNA spike-ins** covering **transcription** and **splicing events**
- **Design of SIRV mixes** to allow for **RNA-Seq pipeline quality control** and **validation**
- Assessment of **differential gene expression** at the transcript level
- The **SIRV Suite**, a Galaxy-based platform for easy and complete spike-in experiment design, data evaluation and comparison
- The **concept of the quality matrix** (coverage CoD, accuracy, precision) and **concordance**
- The SIRV Suite data base to find and research **comparable RNA-Seq experiments** - on the transcript level
- **Upcoming RNA spike-in modules**: from transcript length, poly(A) variants to base modifications

Product portfolio in the RNA-Seq workflow






SPLIT

Fractions for pure RNA sequencing




Sources of variation in RNA sequencing data*



*Hoen, P. A. C 't, et al., 2013, **Reproducibility of high-throughput mRNA and small RNA sequencing across laboratories**, *Nature Biotechnology*

SPLIT RNA Extraction Kit




- 1 Sample homogenization** in a highly chaotropic isolation buffer (guanidinium)
- 2 Acidic phenol-chloroform** extraction aided by phase lock gel tubes
- 3 Silica-column purification and fractioning** of RNA in
 - **total RNA** or
 - **small and large RNA** fractionsby using different volumes of isopropanol

User supplied materials:

Phenol solution pH 4.3, Chloroform, Isopropanol, Ethanol

*Free trial kit (3 extractions including splitting option)



Extraction of high quality RNA



- **No genomic DNA** contamination (no slot-retained band on agarose gel)
- Extraction of **RNA > 10,000 nt**
- **No size bias** and **loss of long RNAs** due to gDNA removal columns
- **No RNA degradation** induced by DNase heat inactivation




Figure A | The SPLIT protocol extracts RNA free of gDNA contamination. In comparison gDNA contamination is present after TRIzol extraction and becomes visible as a slot-retained band on a denaturing agarose gel. Importantly, gDNA contamination cannot be detected with Agilent's Bioanalyzer RNA chip (see Figure B).




Figure B | The SPLIT protocol enables the extraction of high quality RNA. Importantly, no major difference between SPLIT extracted and gDNA contaminated RNA extracted with TRIzol / isopropanol precipitation can be detected on Agilent's Bioanalyzer RNA chip. All samples have a RIN of 8.2 - 8.3.




Figure C | The SPLIT protocol preserves RNA integrity and recovers the complete RNA size range. RNA integrity of an RNA marker (200 – 10,000 nt) after employing different gDNA removal methods was assessed on Agilent's TapeStation.



- RNA for **any downstream application**
- Total RNA from **< 17 nt to > 10,000 nt**
- Efficient recovery of **sRNA**
- **Splitting option** into large and small RNA fractions (cut-off ~150 nt)
- **Highest RNA integrity** and **purity** (RIN > 8; up to RIN 10 for cell culture)
- RNA **free of genomic DNA**
 - without RNA degradation induced by DNase heat inactivation
 - without size bias due to gDNA removal columns
- **High yield** (> 100 µg RNA / column) and **extraction efficiency** (down to 0.5 mg tissue)
- **Species-independent** RNA extraction (e.g., mammals, plants, insects, and others)
- User Guide protocols for **animal / plant tissue, cell culture** and **fluid** samples
- **Convenient** and **fast** protocol (only 30 - 35 min)



- **Mammals** e.g., human, mouse, rat, pig
- **Plants** e.g., Arabidopsis thaliana, Crocus sativus, tomato, Norway spruce (*Picea abies*), Craterostigma plantagineum, ryegrass (*Lolium perenne*), Pak choi (*Brassica rapa* ssp. *Chinensis*), wallflower (*Erysimum Brassicaceae*)
- **Insects** e.g., Drosophila
- **Bacteria** e.g. E.coli, Bacillus subtilis
- **Others** e.g., Sea anemone (*Nematostella*), Jellyfish (*Aurelia*), Frog (*Xenopus laevis*), fungi (*Fusarium*), *P. infestans* sporangia, zoospores



InFusion: Advancing Discovery of Fusion Genes and Chimeric Transcripts from Deep RNA-Sequencing Data

Konstantin Okonechnikov, Aki Imai-Matsushima, Lukas Paul, Alexander Seitz, Thomas F. Meyer, Fernando Garcia-Alcalde

[PLoS One. 2016 Dec 1, doi:10.1371/journal.pone.0167417](#)

Differences in DNA Repair Capacity, Cell Death and Transcriptional Response after Irradiation between a Radiosensitive and a Radioresistant Cell Line

Mireia Borràs-Fresneda, Joan-Francesc Barquinero, Maria Gomolka, Sabine Hornhardt, Ute Rössler, Gemma Armengol & Leonardo Barrios

[Sci. Rep. 6, 27043; doi:10.1038/srep27043 \(2016\)](#)

Alternative Splice Forms Influence Functions of Whirlin in Mechanosensory Hair Cell Stereocilia

Seham Ebrahim, Neil J. Ingham, Morag A. Lewis, Michael J.C. Rogers, Runjia Cui, Bechara Kachar, Johanna C. Pass, Karen P. Steel

[doi:10.1016/j.celrep.2016.03.081](#)

Excessive Osteocytic Fgf23 Secretion Contributes to Pyrophosphate Accumulation and Mineralization Defect in Hyp Mice

Sathish K. Murali, Olena Andrukhova, Erica L. Clinkenbeard, Kenneth E. White, Reinhold G. Erben

[PLoS Biol 14\(4\): e1002427. doi: 10.1371/journal.pbio.1002427](#)



Amphiregulin lacks an essential role for the bone anabolic action of parathyroid hormone.


Freya F. Jay, Mithila Vaidya, Sabrina M. Porada, Olena Andrukhova, Marlon R. Schneider, Reinhold G. Erben
Mol Cell Endocrinol. 2015 Sep 28. pii: S0303-7207(15)30097-6. doi: 10.1016/j.mce.2015.09.031.

Nuclear accumulation of CDH1 mRNA in hepatocellular carcinoma cells.

Ghafoory S, Mehrabi A, Hafezi M, Cheng X, Breitkopf-Heinlein K, Hick M, Huichalaf M, Herbel V, Saffari A, Wölfl S.
Oncogenesis. 2015 Jun 1;4:e152. doi: 10.1038/oncsis.2015.11.

Find the latest SPLIT publications at <https://www.lexogen.com/publications/#splitpub>

Product portfolio in the RNA-Seq workflow






Poly(A) RNA Selection Kit



Poly(A) RNA Selection Kit




- 1 Highly specific binding of poly(A) RNA to oligo(dT)30 beads
- 2 Enrichment of poly(A) RNA through magnetic bead-based purification



- **Highly specific poly(A) enrichment** from total RNA (ct rRNA mapping of only 0.0003% of all reads)
- **Straightforward** protocol - only **one binding step** required due to high specificity
- **Recovery of 1-3%** poly(A) RNA of total RNA (depending on the input material)
- **Up- and down-scaling** of the protocol (500 ng - 100 µg of total RNA)
- **Rapid** turnaround – only 20 minutes hands-on time
- **Automation-friendly**
- Poly(A) RNA as **pure eluate** or for **direct insertion into downstream applications**
- **Various downstream applications** such as RNA-Seq library preparation (SENSE Total RNA-Seq, microarrays, SAGE, RACE, cDNA synthesis, cDNA library construction, RT-PCR and others)

Product portfolio in the RNA-Seq workflow






RiBO COP™

Select and Deplete



RiboCop rRNA Depletion Kit V1.2




1 Specific depletion of rRNA sequences through hybridization to affinity probes

2 Magnetic bead-based purification



- **Highly specific elimination (> 98%) of cytoplasmic and mitochondrial rRNA**
- For **human, mouse, and rat** samples
- Suitable for any quality RNA, **intact** and **degraded RNA** (e.g., FFPE samples)
- Standard recovery from **1 - 3% of input RNA** (depending on the input material)
- **Up- and down-scaling** of the protocol (1 ng - 1 µg of total RNA / prep)
- **Convenient** and **fast** protocol – rRNA depleted RNA in only 2 hours
- **Automation**-friendly protocol
- **Simple** and **robust** workflow retaining **full-length transcripts** – no enzymatic reactions or mechanical shearing steps
- **All-in one** solution – no additional purification modules required
- **Seamless integration** in any downstream application (e.g., RNA-Seq library preps including Lexogen's SENSE Total RNA-Seq; Cat. No. 042.04*/08/24/96)

Product portfolio in the RNA-Seq workflow






TELO™ PRIME


Full mRNA - precise ends



TeloPrime Full-Length cDNA Amplification Kit



Lexogen's unique - Cap Dependent Linker Ligation (CDLL)



Poly(A) Tail and cap → ligation → second strand synthesis



- **Exceptional 5'-cap specificity** for full-length cDNA
- **Faithful representation** of the **mRNA transcriptome**
- **Preservation** of **short** and **long RNA** molecules
- Input amounts of **1 ng to 2 µg** of total RNA
- All-in-one protocol – only **60 minutes hands-on** time
- **Extensive tests** on mouse tissue (e.g., liver, lung etc.), Universal Human Reference RNA and Human Brain Reference RNA
- **Target specific application** – exchange Reverse Transcription, PCR Forward, and PCR Reverse Primers with your own primers of interest (e.g., gene-specific)
- **Various downstream applications** such as NGS library generation (incl. PacBio), microarray probe generation, RACE, cloning, and normalization



Global survey of cis-regulation in mammalian translation *Dissertation

Jingyi Hou

[Department of Biology, Chemistry and Pharmacy of Freie Universität Berlin, Germany](#)

Thyroglobulin Represents a Novel Molecular Architecture of Vertebrates

Guillaume Holzer, Yoshiaki Morishita, Jean-Baptiste Fini, Thibault Lorin, Benjamin Gillet, Sandrine Hughes, Marie Tohmé, Gilbert Deléage, Barbara Demeneix, Peter Arvan and Vincent Laudet

[JBC.M116.719047. doi: 10.1074/jbc.M116.719047](#)

cDNA Library Enrichment of Full Length Transcripts for SMRT Long Read Sequencing

Maria Cartolano, Bruno Huettel, Benjamin Hartwig, Richard Reinhardt, Korbinian Schneeberger

[PLoS ONE 11\(6\):e0157779. doi:10.1371/journal.pone.0157779](#)


Pervasive isoform-specific translational regulation via alternative transcription start sites in mammals

Xi Wang, Jingyi Hou, Claudia Quedenau, Wei Chen

[Molecular Systems Biology \(2016\) 12, 875, DOI 10.15252/msb.20166941](#)

Find the latest TeloPrime publications at <https://www.lexogen.com/publications/#teloprimepub>

Product portfolio in the RNA-Seq workflow






MIX²

Accurate Analysis of RNA-Seq Data



Mix² features – Your benefits



Accurate estimation of isoform concentration

Get more relevant cases of differentially expressed isoforms and a reduced number of false positives



Identification of bias types


Detect and classify bias types



Fast run times

Speed up your data analysis and profit from a small memory footprint

Watch Mix² explained by our bioinformatics experts!



<https://www.youtube.com/watch?v=9s5fOWf-s2s>



Andreas Türk, PhD
Head of Bioinformatics, Lexogen



Tristetraprolin binding site atlas in the macrophage transcriptome reveals a switch for inflammation resolution

Vitaly Sedlyarov, Jörg Fallmann, Florian Ebner, Jakob Huemer, Lucy Sneezum, Masa Ivin, Kristina Kreiner, Andrea Tanzer, Claus Vogl, View ORCID ProfileIvo Hofacker, Pavel Kovarik

[doi: 10.15252/msb.20156628](https://doi.org/10.15252/msb.20156628)

Find the latest Mix² publications at <https://www.lexogen.com/publications/#mix2pub>



- **QuantSeq 3' (4)***
- **SENSE mRNA and total (4)***
- **SPLIT (3)***
- **RiboCop (4)***
- **TeloPrime (4)***

- ***Applications support is always available**



Modules and Add-ons

Enabling Complete Transcriptome Sequencing

Questions?

Visit our website: ww.lexogen.com – FAQ
Contact us: +43 1345 1212-41 or info@lexogen.com