



UF-1000i Start-Up

Power On/Log On Information Processing Unit (IPU)


1. Turn on IPU power.
2. When prompted press **[Ctrl]+[Alt]+[Delete]** to logon to Windows®.
3. Input the user name and password, then press **[ENTER]** or click **[OK]**.
4. When IPU logon displays, input user name and password (if required).
5. Press **[ENTER]** or click **[OK]**.



Power Switch on right side



Note : *If a password is defined for the Windows® logon, users have 3 chances to input the correct password. After 3 failed attempts, the account will be disabled for 30 minutes. Unlocking before 30 minutes requires administrative access; call the Sysmex Technical Assistance Center (1-888-879-7639. In Canada, call 1-888-679-7639)*



Power On Main Unit

IPU must display UF main menu screen before powering on main unit.


1. Turn on main unit power by pressing the switch on lower right side of the main unit.
- 2.. Press the Start-Up switch on the right front of the main unit.




UF-1000i Start-Up

Main Unit Start-Up

1. The UF performs a series of self-checks. If any check fails, an error message displays. Refer to the **UF-1000i Instructions for Use**, Troubleshooting chapter.




2. The “Temperature Stability” dialog box displays the temperature of the reaction chambers, the sheath liquid heater and the flowcell module (FCM) until they stabilize. If this box remains on screen for >30 minutes, please turn the main unit power OFF and call the Sysmex Technical Assistance Center.



Background Check

1. After the check of mechanical parts and the temperatures are stabilized, the background check is performed.
 - If the background is within acceptable limits (see chart below), the analyzer icon on the IPU and the “Ready” LED light will be **green** and the instrument is ready to process specimens.

Blank	Limit
RBC	1.0 [/ μ L]
WBC	1.0 [/ μ L]
EC	1.0 [/ μ L]
CAST	0.20 [/ μ L]
BACT	1.0 [/ μ L]
Total Count	300 [count]
BAC Total Count	3000 [count]




“Ready” LED

UF-1000i Start-Up

Background Check

- If any background count is outside the acceptable limits, the Background Check dialog box opens, and parameters exceeding the acceptable limits are displayed in **red**. The “Background Error” message will also display.




The dialog box titled "Background Check - UF-1000i" contains a "Close" button in the top right. Below the title bar, there is a "Background Check" section with a "Sample ID" field containing "BLANKCHECK00002". A table follows with columns for "Blank" and "Limit". The "Total Count" value of 433 is highlighted in red, indicating it is above the limit of 300. The "BAC Total Count" is 660, which is below its limit of 3000. Below the table, a message reads: "Blank Check Please wait for a moment. It will take approx. 3 minutes."

	Blank	Limit
RBC	0.0	1.0 [/uL]
WBC	0.0	1.0 [/uL]
EC	0.0	1.0 [/uL]
CAST	0.00	0.20 [/uL]
BACT	0.0	1.0 [/uL]
Total Count	433	300 [count]
BAC Total Count	660	3000 [count]

Blank Check
Please wait for a moment. It will take approx. 3 minutes.

- When the “Background Error” message displays, click **[OK]** in the dialog box to perform another auto rinse. If the background values are still not within acceptable limits, refer to the **UF-1000i Instructions For Use**, Troubleshooting chapter.



The dialog box titled "Help - UF-1000i" features an "Error List" section with "Background Error" listed. Below this is an "Action" section with the instruction: "1) Press OK button, and execute autorinse sequence." On the right side, there are three buttons: "OK", "Cancel", and "Reset Alarm".

Error List
Background Error

Action
1) Press OK button, and execute autorinse sequence.

UF-1000i Start-Up

Pressure and Vacuum

Pressure and vacuum readings are monitored by the UF. If a pressure or vacuum error occurs, an error message displays. Click **[OK]** in the Help dialog box to display the sensor dialog box.

SV/SNS - UF-1000i							
SV							
1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
SNS							
1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
0.26MPa Plus	0.255MPa	S-ReactorUnit	33.6°C	Air Pressure Temperature	32.6°C		
Sheath Plus	0.123MPa	B-ReactorUnit	40.5°C	Room Temperature	23.7°C		
0.05MPa Plus	0.049MPa	Sheath Fluid Temperature	36.1°C				
0.05MPa minus	-0.054MPa	FCM	34.9°C	LaserElectroCurrent	0mA		
All SV auto ON/OFF		SV No.	1	ON	OFF	OK	

Confirm the pressure and/or vacuum readings are within range, when the UF is in **READY** status, according to the table below. If the error persists, refer to the **UF-1000i Instructions for Use**, Troubleshooting chapter.

UF-1000i	NORMAL RANGE	Location/Action
0.26 MPa pressure	0.22~0.34	Contact Technical Assistance Center
Sheath pressure	0.07~0.15	Contact Technical Assistance Center
0.05 MPa pressure	0.045~0.054	Upper front right interior, Main Unit
-0.053 MPa vacuum	-0.049 ~ -0.057	Upper front right interior, Main Unit

Vacuum (bellows) Unit



Vacuum Trap

Pressure Adjustment Knob

UF-1000i QC Set-Up

QC File Set Up

There are 24 QC files available on the UF. QC file lot information can be entered in the QC dialog box. The lot number, target and limit values of the UFII CONTROL™ product can be automatically entered by using the handheld barcode reader or by manually entering via the keyboard. Target and limit values may vary depending on the production lot.

Enter QC Lot Number—Handheld Barcode Reader Method

1. Select the **[QC Files]** icon from the main menu, or press **[F5]** on the keyboard.
2. Select the appropriate file number (1-24).
3. Select the QC input dialog box by clicking the **[Edit]** icon on the tool bar or **[F9]** on the keyboard.



	Nickname	File No.	Material	Lot No.	Regist. Date	Exp. Day	Analysis Date
	Main1	QC01	UF II CONTROL-H	QC-YA6001	11/14/2006	12/16/2006	11/29/2006 11:13
	Main1	QC02	UF II CONTROL-L	QC-YA6001	11/14/2006	12/16/2006	12/15/2006 11:50
	Main1	QC03	UF II CONTROL-L	QC-YA6002	11/30/2006	02/16/2007	12/13/2006 16:48
	Main1	QC04	UF II CONTROL-H	QC-YA6002	11/30/2006	02/16/2007	12/13/2006 16:46
	Main1	QC05	UF II CONTROL-L	QC-YA6003	11/30/2006	03/13/2007	02/13/2007 10:04
	Main1	QC06	UF II CONTROL-H	QC-YA6003	11/30/2006	03/13/2007	02/13/2007 10:07
	Main1	QC07	UF II CONTROL-L	QC-YA6004	11/30/2006	04/10/2007	02/13/2007 10:10
ERROR	Main1	QC08	UF II CONTROL-H	QC-YA6004	11/30/2006	04/10/2007	02/13/2007 10:13
	Main1	QC09	UF II CONTROL-L	QC-YA6006	01/17/2007	06/05/2007	02/28/2007 11:54
ERROR	Main1	QC10	UF II CONTROL-H	QC-YA6006	01/17/2007	06/05/2007	02/28/2007 11:57
	Main1	QC11					
	Main1	QC12					
	Main1	QC13					

4. Begin with the “UTL/Exp. Date/Lot No” barcode at the top left hand corner of the UFII CONTROL™ assay sheet. This will select the appropriate level (LOW) and automatically enter the 30 day open vial expiration date into the dialog box.
5. Continue scanning the barcodes for each of the low level parameters; the target and limit entry area will automatically populate. Click **[OK]** to save.
6. Select the next sequential file in which to enter the targets for the next (HIGH) control.
7. Select the QC input dialog box by clicking the **[Edit]** icon on the tool bar or by pressing **[F9]** on the keyboard.
8. Scan the “UTH/Exp. Date/Lot No.” barcode in the top left hand corner of the UFII CONTROL™ assay sheet. This will select the appropriate level (HIGH) and automatically enter the 30 day open vial expiration date into the dialog box.
9. Continue scanning the barcodes for each of the high level parameters; the target and limit entry areas will automatically populate. Click **[OK]** to save.

UF-1000i QC Set-Up

Enter QC Lot Number—Manual Entry Method

1. Select the **[QC Files]** icon from the main menu, or press **[F5]** on the keyboard.
2. Select the appropriate file number (1-24).
3. Select the QC input dialog box by clicking the **[Edit]** icon on the tool bar or **[F9]** on the keyboard.
4. Select the appropriate level of control to be entered from the “Material” drop down box.

Input Lot Information

Nickname: UF-1000i File No.: QC15

Material: **UF II CONTROL -H** Lot No. QC- E:

Item	Lower Limit	Target	Upper Limit	Unit
RBC				/uL
WBC				/uL
EC				/uL
CAST				/uL
BACT				/uL
Cond.				mS/cm

5. Using the keyboard, manually enter the 6-digit lot number in the space provided to the right of the “Lot No. QC-” notation on the UFII CONTROL™ assay sheet.
6. Click on the expiration date drop down box to display the calendar. Click on the appropriate expiration date (30 days or date on control sheet or bottle, whichever is earlier)

April, 2007

Unit	Sun	Mon	Tue	Wed	Thu	Fri	Sat
/uL	25	26	27	28	29	30	31
/uL	1	2	3	4	5	6	7
/uL	8	9	10	11	12	13	14
/uL	15	16	17	18	19	20	21
/uL	22	23	24	25	26	27	28
/uL	29	30	1	2	3	4	5
mS/cm	Today: 4/20/2007						

UF-1000i QC Set-Up

Enter QC Lot Number—Manual Entry Method

7. Enter the target and limit values posted on the assay sheet for each measured parameter. Click **[OK]** to save.

Input Lot Information

Nickname File No.

Material Lot No. QC- Exp. Day

Item	Lower Limit	Target	Upper Limit	Unit
RBC	146.1	182.6	219.1	/uL
WBC	608.0	760.0	912.0	/uL
EC	40.3	80.7	121.1	/uL
CAST	9.37	18.74	28.11	/uL
BACT	601.5	802.0	1002.5	/uL
Cond.	33.6	37.3	41.0	mS/cm
S_FSC	99.1	-----	134.1	ch
S_FSCW	22.9	-----	28.1	ch
S_FLH	13.9	-----	32.5	ch
S_FLL	91.1	-----	212.5	ch
S_FLLW		-----		ch
S_SSC		-----		ch
B_FSC		-----		ch
B_FSCW		-----		ch
B_FLH		-----		ch

Manual Setting

Item

Lower Limit Upper Limit

8. For the next control (UFII Control –H), select the next sequential file number. Enter the target and limit values for the measured parameters by following steps 3 through 7 above.
9. Using the keyboard, manually enter the lower and upper limit range of the sensitivity parameters posted on the UFII CONTROL™ assay sheet. Click **[OK]** to save.

UF-1000i QC Analysis


QC Analysis Notes

- ◆ Allow the controls to equilibrate to room temperature for 20-30 minutes before analysis.
- ◆ Once opened, the control product is stable for 30 days if stored at 2-10°C. Do not freeze.
- ◆ After opening the UFII CONTROL™, store in the original box to prevent exposure to direct sunlight.
- ◆ The product expiration date is printed on the box and on the label of each control bottle.

QC Analysis

1. With the UF in “Ready” (LED is **green**), click the **[Manual]** icon from the main menu or **[F2]** on the keyboard.
2. Click on the **[QC]** icon on the right side of the window.
3. Choose the appropriate file from the list displayed.

Note: Only those QC files currently configured with lot number and expiration date will be listed.



File No.	Material	Lot No.	Exp. Day	Analysis Date
QC01	UF II CONTROL -H	QC-YA6006	07/27/2007	
QC02	UF II CONTROL -L	QC-YA6006	07/27/2007	
QC10	UF II CONTROL -L	QC-ya6006	07/20/2007	07/24/2007 15:06:04
QC11	UF II CONTROL -H	QC-ya6006	07/20/2007	07/17/2007 15:45:27
QC14	UF II CONTROL -H	QC-4444	07/06/2007	07/17/2007 15:49:03
QC15	UF II CONTROL -L	QC-5555	07/06/2007	07/24/2007 14:08:57
QC16	UF II CONTROL -H	QC-YA8008	08/24/2007	
QC17	UF II CONTROL -L	QC-YA8008	08/24/2007	
QC23	UF II CONTROL -L	QC-YA7007	08/04/2007	
QC24	UF II CONTROL -H	QC-YA7007	08/04/2007	07/16/2007 15:04:25

4. Click **[OK]**. The “QC Measurement Help” dialog box will display, click **[Close]**.
5. Verify the appropriate lot number and control level is displayed.
6. Mix the UFII Control bottle by shaking until there is no particle sediment at the bottom, then shake vigorously another 20 times to mix.

QC Result L-J - UF-1000i

Nickname: UF-1000i [Accept]

File No.: QC02 [Re-analyze]

Material: UF II CONTROL -L [Cancel]

Lot No.: QC-PD004

Exp. Day: 08/11/2006

QC Measurement Help [Close]

Shake a bottle of UF II CONTROL a few times until there is no particle sediment at the bottom, then shake it vigorously another 20 times.

Immediately (within 10 seconds) after mixing, press the side of the bottle gently to allow 0.9mL (15 to 20 drops) of the reagent to drop from the tip nozzle into a new test tube.


Immediately (within 10 seconds) after distribution, place the test tube into the aspiration pipette and turn the start switch on.

This screen is automatically opened when the QC Measurement.

UF-1000i QC Analysis

QC Analysis

7. Immediately (within 10 seconds) after mixing, dispense 0.9 ml of the UFII CONTROL™ into a new test tube or cup.
8. Immediately (within 10 seconds) after dispensing, place the test tube at the aspiration pipette and hold the tube so the aspiration pipette is well below the meniscus of the control material.
9. Press the manual **[Start]** switch. Hold the tube in place until you hear the last two quick beeps (after 6 slower beeps), which will signal the end of the aspiration. Discard the test tube after measurement; any control material remaining in the tube should not be used.



10. The QC results will display in a dialog box when the sample measurement has completed. Press **[Accept]** for results to plot or **[Cancel]** to abort.
11. To repeat the measurement for the next level of control, close all open windows, click **[Manual]** from the main menu, and continue from step 2.

QC Result L-J - UF-1000i

Nickname	UF-1000i	Accept
File No.	QC17	
Material	<input checked="" type="checkbox"/> UF II CONTROL -L	Re-analyze.
Lot No.	QC-YA8008	Cancel
Exp. Day	08/24/2007	
RBC	36.0 /uL	
WBC	36.2 /uL	
EC	10.5 /uL	
CAST	4.97 /uL	
BACT	211.2 /uL	
Cond.	9.4 mS/cm	
S_FSC		ch
S_FSCW		ch
S_FLH		ch
S_FLL		ch
S_FLLW		ch
S_SSC		ch
B_FSC		ch
B_FSCW		ch
B_FLH		ch
		Scattergram


Press the Accept button to store the QC analysis result, and the QC data will be plotted on the chart.

QC Measurement Help

UF-1000i QC Analysis


QC Analysis

12. Out of range results display with a **red** background. Press **[Accept]** for the results to plot on the graph, or **[Cancel]** to abort. Press **[Re-analyze]** to process a fresh pour of the control, if needed.



13. The L-J graph will automatically appear (if configured to do so in the IPU settings). Confirm the plotted data point is within the ranges previously entered into the QC file.


14. To repeat the measurement for the next level of control, close all open windows, click **[Manual]** from the main menu, and continue from step 2.



UF 1000i QC Analysis

Establishing a New Lot

1. Analyze each control level once and compare to the values entered into the QC file from the assay sheet.
2. Analyze the new control lot at least 10 times before the current lot number expires. For statistical integrity, the accumulation of points should occur over a number of days, for a minimum of five days, two data points per day.
3. To check the mean, SD and %CV of the new QC lot data, click on the dark green cursor on the L-J chart and drag to the left to include all points, or press **[CTRL]** and **[A]** to select all points.
4. The calculated SD, mean and %CV will be displayed in the column to the far right of the chart.




Set the range of data and observe SD, mean and %CV

Comparing New Lot to Existing Lot

1. The results of the new lot number may be compared to an existing lot number by selecting **[F11]** while reviewing the L-J chart of the new QC file, and choosing the appropriate file for the existing QC lot to be compared. The second chart will be superimposed over the first, with the main QC chart data displayed in blue, and the second QC chart data displayed in dark gray.




UF 1000i QC Analysis



Main QC file displays in blue, superimposed QC chart in dark gray

Reviewing QC : Viewing Radar Charts

1. Select [QC Files] from the main menu, or [F5] on the keyboard.
2. Select the appropriate QC file number (1-24). Do not open the file.
3. Once the file is highlighted, the radar chart will display on the right side of the screen. The control data is displayed in blue on the radar chart, and date and time of analysis are displayed on the immediate left of the radar chart in the "Analysis Date" column.
4. Parameters that exceed the acceptable limit will have a red "X" displayed on the radar chart. The parameter name will be backlit in red and the word "ERROR" will be displayed to the left of the file number.




Outside acceptable limits

UF 1000i QC Analysis




Reviewing QC : Viewing Levy—Jennings Charts

1. Select **[QC Files]** from the main menu or **[F5]** on the keyboard.
2. Double-click on the appropriate file number (1-24) to select. The L-J chart will display.
3. Parameters that exceed the acceptable limits will have a red **"X"** displayed as the plotted point. The parameter name as well as the result value will be backlit in red.
4. To view the remaining parameters, use the scroll bar on the right side of the screen or the down arrow on the keyboard.



Sending QC Data to Host

1. Select **[QC Files]** from the main menu or **[F5]** on the keyboard, and double-click the appropriate file number (1-24).
2. Set the range of QC data to be output by clicking on the dark green cursor and dragging it to the left to include all points to be sent. If all data is required, press **[CTRL]** and **[A]** to select all points.
3. Click **[Output]** on the menu bar or press **[F12]**, then select **[Host (HC)]** or press **[F12]** then **[1]**.




F12

UF 1000i QC Analysis


Printing QC Reports

1. Select **[QC Files]** from the main menu or **[F5]** on the keyboard, and double-click on the appropriate file number (1-24) to select.
2. Set the range of QC data to print by clicking on the dark green cursor and dragging it to the left to include the points to be printed. If all data is required, press **[CTRL]** and **[A]** to select all points.
3. Click **[Report]** on the menu bar or press **[F12]**. To print the L-J charts, select **[Report (GP) (R)]**. To print the ledger (list of data points) select **[Report (LP)(L)]**.



Select the range of QC data to print or send to Host

UF 1000i QC Analysis



Deleting Data : Entire QC File

Note: Deleting an entire file will delete all data, targets, limits and lot information.

1. Select **[QC Files]** from the main menu or press **[F5]**, and click on the appropriate file number (1-24) to highlight. Do not open the file.
2. Click the **[DELETE]** icon on the menu bar. A dialog box displays with the following message:
"1 File(s) will be deleted. Are you sure?"
3. Click **[OK]** to delete or **[CANCEL]** to abort.

Nickname	File No.	Material	Lot No.	Regist. Date	Exp. Day	Analysis Date	
Main1	QC01	UF II CONTROL -H	QC-YA6001	11/14/2006	12/16/2006	11/29/2006 11:13:51	
Main1	QC02	UF II CONTROL -L	QC-YA6001	11/14/2006	12/16/2006	12/15/2006 11:50:57	
Main1	QC03	UF II CONTROL -L	QC-YA6002	11/30/2006	02/16/2007	12/13/2006 16:45:20	
Main1	QC04	UF II CONTROL -H	QC-YA6002	11/30/2006	02/16/2007	12/13/2006 16:46:33	
Main1	QC05	UF II CONTROL -L	QC-YA6003	11/30/2006	03/13/2007	02/13/2007 10:04:28	
Main1	QC06	UF II CONTROL -H	QC-YA6003	11/30/2006	03/13/2007	02/13/2007 10:07:18	
Main1	QC07	UF II CONTROL -L	QC-YA6004	11/30/2006	04/10/2007	02/13/2007 10:10:03	
ERROR	Main1	QC08	UF II CONTROL -H	QC-YA6004	11/30/2006	04/10/2007	02/13/2007 10:13:05
ERROR	Main1	QC09	UF II CONTROL -L	QC-YA6006	01/17/2007	06/05/2007	02/28/2007 11:54:00
ERROR	Main1	QC10	UF II CONTROL -H	QC-YA6006	01/17/2007	06/05/2007	02/28/2007 11:57:44
ERROR	Main1	QC11	UF II CONTROL -L	QC-YA6006	05/08/2007	06/05/2007	05/08/2007 14:59:25
ERROR	Main1	QC12	UF II CONTROL -L	QC-YA6006	05/08/2007	06/05/2007	05/08/2007 15:08:15
ERROR	Main1	QC13	UF II CONTROL -L	QC-YA6006	05/08/2007	06/05/2007	05/08/2007 15:57:46
	Main1	QC14				05/08/2007 10:16:46	
	Main1	QC15					
	Main1	QC16					
	Main1	QC17					
	Main1	QC18					
	Main1	QC19					
	Main1	QC20					
	Main1	QC21					
	Main1	QC22					
	Main1	QC23	UF II CONTROL -L	QC-1111	04/20/2007	04/27/2007	
	Main1	QC24	UF II CONTROL -H	QC-2222	04/20/2007	04/27/2007	



Deleting Data : Data Point(s)

1. Select **[QC Files]** from the main menu or press **[F5]**, and click on the appropriate file number (1-24) to highlight. Do not open the file.
2. Select the data for deletion:
 - Single point—click on specific point
 - Multiple points—click on dark green cursor and drag to the left until all desired points are selected

UF 1000i QC Analysis

Deleting Data : Data Point(s)

- Click the **[DELETE]** icon on the menu bar. A dialog box displays with the following message :
“Are you sure you want to delete “X” plot(s) data?”
- Click **[OK]** to delete or **[CANCEL]** to abort.



Item	UL	Target	LL	Data	SD	Mean	CV
WBC	957.0		638.0	-178.5	178.5	---	---
EC	120.3		40.1	-9.0	9.0	---	---
CAST	27.83		9.27	18.94	18.94	---	---
BACT	1046.5		627.9	-107.0	107.0	---	---
Cond.	41.4		33.8	-25.7	25.7	---	---

UF-1000i Sample Processing

Sampler/Auto Mode

1200 μ L aspirated volume, 4 mL minimum tube volume

1. Verify the UF main unit 'Ready' LED is **green**.
2. Place barcoded samples in sample racks with barcodes facing the open side of the rack.
3. Place rack on right side of sampler unit, with the notch to the right.
4. Click **[SAMPLER]** on the IPU main menu screen, or press **[F3]** on the keyboard.
5. Confirm the rack number and tube position in the Sampler Sample No. dialog box.
6. Click **[Sampler Start]** on the IPU main menu screen, or press **[Enter]**.
7. Completed results can be viewed in the IPU Explorer and/or Browser screens.



Sampler Sample No. - UF-1000i

Sample ID
1

Sampler Start
Cancel

Rack No. 1 Tube Position 1

Default Order
SEDIMENT

Confirm the rack number and tube position; click Sampler Start

Out	Sample No.	Date	Time	REVIEW	ERROR	Rack No.	Tube	Discrete	Rackname	Seq
V	DH	0837	B	01/1/2007	14:11:30			SEDIMENT	UF-1000i	197
V	DH	0806	A	01/1/2007	14:11:25			SEDIMENT	UF-1000i	196
V	DH	0805	M	01/1/2007	14:09:28			SEDIMENT	UF-1000i	195
V	DH	0801	M	01/1/2007	13:59:09			SEDIMENT	UF-1000i	194
V	DH	0897	M	01/1/2007	13:57:43			SEDIMENT	UF-1000i	193
V	DH	0893	M	01/1/2007	13:56:27			SEDIMENT	UF-1000i	192
V	DH	0889	M	01/1/2007	13:55:19			SEDIMENT	UF-1000i	191
V	DH	0885	M	01/1/2007	13:53:53			SEDIMENT	UF-1000i	190
V	DH	0881	M	01/1/2007	13:52:32			SEDIMENT	UF-1000i	189
V	DH	0877	M	01/1/2007	13:50:19			SEDIMENT	UF-1000i	188
V	DH	0873	M	01/1/2007	13:48:50			SEDIMENT	UF-1000i	187
V	DH	0869	M	01/1/2007	13:47:31			SEDIMENT	UF-1000i	186
V	DH	0865	M	01/1/2007	13:46:13			SEDIMENT	UF-1000i	185
V	DH	0861	M	01/1/2007	13:44:56			SEDIMENT	UF-1000i	184
V	DH	0857	M	01/1/2007	13:42:35			SEDIMENT	UF-1000i	183
V	DH	0853	M	01/1/2007	13:40:08			SEDIMENT	UF-1000i	182
V	DH	0849	M	01/1/2007	13:38:51			SEDIMENT	UF-1000i	181
V	DH	0845	M	01/1/2007	13:37:23			SEDIMENT	UF-1000i	180
V	DH	0841	M	01/1/2007	13:35:03			SEDIMENT	UF-1000i	179
V	DH	0837	M	01/1/2007	13:33:42			SEDIMENT	UF-1000i	178
V	DH	0833	M	01/1/2007	13:32:23			SEDIMENT	UF-1000i	177
V	DH	0829	M	01/1/2007	13:30:01			SEDIMENT	UF-1000i	176
V	DH	0825	M	01/1/2007	13:28:32			SEDIMENT	UF-1000i	175
V	DH	0821	M	01/1/2007	13:27:13			SEDIMENT	UF-1000i	174
V	DH	0817	M	01/1/2007	13:25:56			SEDIMENT	UF-1000i	173
V	DH	0813	M	01/1/2007	13:24:40			SEDIMENT	UF-1000i	172
V	DH	0809	M	01/1/2007	13:23:20			SPERM	UF-1000i	171

Results can be viewed in Explorer....

.....or in Browser

REVIEW
Validated

05/25/2007 13:45:43

Analysis Parameters

RBC	324.0 /uL	58.3 /HPF
WBC	303.1 /uL	54.6 /HPF
EC	145.8 /uL	26.2 /HPF
CAST	0.74 /uL	2.15 /LPE
BACT	0.0 /uL	0.0x10A0 /mL

Research Parameters

X-TAL	0.0 /uL
YLC	0.0 /uL
SRC	109.5 /uL
Path_CAST	0.74 /uL
MUCUS	0.00 /uL
SPERM	0.0 /uL
Cond.	54.4 mS/cm

Research Information

RBC : Isomorphic?
Cond. : RANK5

REVIEW Comment

Universal
Abn: QC
Sensitivity


0043

UF-1000i Sample Processing

Manual Mode (Open Tube Sampling)

800 µL aspirated volume, 1 mL minimum tube volume

Note: Perform an autorinse after analysis of visibly bloody, pyuric or mucoid samples.



1. Verify UF main unit 'Ready' LED is **green**, and click the **[MANUAL]** icon on the IPU main menu screen, or press **[F2]** on the keyboard.
2. Enter the patient or sample identification number (up to 15 characters) in the "Sample ID" field :
 - Use the handheld barcode reader to scan the barcoded sample ID
 - OR
 - Use the keyboard to manually type in the patient/sample ID number
3. Choose the appropriate collection time, source, color and clarity, if applicable, then click **[OK]**.



Note: Select Discrete testing mode **ONLY** if the Host (LIS) is not providing order information

A screenshot of the UF-1000i software interface. The title bar reads 'Manual Sample No. - UF-1000i'. The main window has a light beige background. At the top, there is a 'Sample ID' field containing the number '0075'. To the right of this field are 'OK' and 'Cancel' buttons. Below the 'Sample ID' field, there are radio buttons for 'Discrete' testing mode: 'SEDIMENT' (selected) and 'ONLY BACTERIA'. Below that is a 'Sample Comment' text box. Further down, there are fields for 'Collection Time' (date: 07/26/2007, time: 14:59:01) and 'Source' (radio buttons for Unknown, OP.CLCT, Morning, Timed, AF.Meal, Cath). Below that are 'Color' options (radio buttons for Unknown, None, LyBrown, Yellow (selected), YBrown, Orange, Red, DBrown, Green, Blue, White) and 'Clarity' options (radio buttons for Unknown, Clear (selected), SlHazy, Hazy, SlClDy, Cloudy). At the bottom, there are fields for 'Patient ID' (containing 569365) and 'Patient Name'.



4. Place a well-mixed patient sample at the aspiration pipette and hold the tube so the pipette is well below the meniscus of the sample.
5. Press the green manual **[Start]** button on the front of the UF.
6. Hold the tube in place until you hear the last two quick beeps (after six slower beeps) which signals the end of the aspiration.
7. When the "ready" LED starts to flash green/orange, the analysis is in progress. Prepare the next sample for analysis, and repeat from step 1.
8. Completed results can be viewed in either Explorer or Browser.



UF-1000i Sample Explorer/Stored Data


Finding A Sample

1. Click on the **[Sample Explorer]** icon, or press **[F7]** on the keyboard. Verify **[Last 20]** icon on the menu bar is de-selected to access all stored sample data.
2. From menu bar, click **[Edit (E)]** then **[Find (F)]**, or press **[Ctrl]+F** on keyboard. The “Find” input dialog box displays.
3. Enter Sample number or Patient ID in the Find dialog window.
 - Use “?” in place of a single digit for an unknown number or letter.
 - Use “*” for multiple unknown letters or numbers and place at the end of a string.
4. Click **[NEXT]** or press **[ENTER]** to locate the sample number entered; the located sample will be backlit in **blue**.
5. Click **[NEXT]** to search below the highlighted sample number or **[PREV]** to search above that sample number.
6. The located sample will be backlit in blue.
7. Select **[CLOSE]** to exit the Find dialogue window.




Editing Sample Number or Patient ID

1. Click on the **[Sample Explorer]** icon, or press **[F7]** on the keyboard. Verify “Last 20” icon is de-selected to access all stored sample data.
2. Click/Highlight the sample number to be changed.
3. Click on the **[Validate]** icon or press **[F11]** to un-validate sample. (“V” disappears in first column of Sample Explorer screen.)
4. From the Menu Bar, click **[Edit]** and **[MODIFY (P)]** or press **[F10]**. The Modify dialog box displays with selected Sample number.



5. In the Modify box, click on the Sample No. or Patient ID to edit. Delete old number and input new ID number.
6. Click **[OK]** and new sample information will be stored. The corrected sample will be backlit in **blue**.
7. Click on the **[Validate]** icon to re-validate the sample.




Note: The newly edited sample ID will NOT automatically transmit to the Host (LIS). Refer to “Reprint to Graphic Printer or Re-transmit Data to Host”.

UF-1000i Sample Explorer/Stored Data

Reprint to Graphic Printer or Retransmit Data To Host

Note: Data that is not validated will not be reprinted or retransmitted.


1. Click **[Explorer]**, or press **[F7]**. Verify the **[Last 20]** icon is de-selected to access all stored sample data.
2. Click on sample(s) you wish to reprint or retransmit to Host :
 - To select any single sample, click on sample ID.
 - To select multiple samples (not necessarily in order), press and hold **[CTRL]** on keyboard and click each sample.



The screenshot shows the IPU Sample Explorer application window. The title bar reads "IPU - [Sample Explorer Filter[ALL] Sort[Date DESC , Time DESC] 809 Record(s)]". The menu bar includes File(F), Edit(E), View(V), Record(R), Action(A), Report(P), Setting(S), Window(W), and Help(H). The toolbar contains icons for Help, MANUAL, SAMPLER, Menu, QC Files, Work list, Explorer, Browser, Modify, Validate, Out, Upper, Lower, Sort, Last20, and Delete. The main data table has the following columns: V, Out, Sample No., Date, Time, REVIEW, ERROR, Rack No., Tube, Discrete, Nickname, and Seq. The data rows are as follows:

V	Out	Sample No.	Date	Time	REVIEW	ERROR	Rack No.	Tube	Discrete	Nickname	Seq.
V	DGH	0075	M 07/26/2007	15:00:48	REVIEW				SEDIMENT	UF-1000i	34
V	DGH	0995	M 07/26/2007	14:56:16	REVIEW				SEDIMENT	UF-1000i	33
V	DGH	0033	B 07/26/2007	14:50:40	REVIEW		3	10	SEDIMENT	UF-1000i	32
V	DGH	0065	B 07/26/2007	14:50:04	REVIEW		3	9	SEDIMENT	UF-1000i	31
V	DGH	0073	B 07/26/2007	14:49:28	REVIEW		3	8	SEDIMENT	UF-1000i	30
V	DGH	0075	B 07/26/2007	14:48:52	REVIEW		3	7	SEDIMENT	UF-1000i	29
V	DGH	0032	M 07/26/2007	14:45:53	REVIEW				SEDIMENT	UF-1000i	28
V	DGH	0077	B 07/26/2007	14:37:29	REVIEW		3	6	SEDIMENT	UF-1000i	27
V	DGH	0072	B 07/26/2007	14:36:53	REVIEW		3	5	SEDIMENT	UF-1000i	26

- To highlight a block of samples, click on first sample to output. Then press and hold **[Shift]** and click on last sample to output.
3. To reprint to graphic printer, click **[Report]** and **[Report(GP)]** or press **[F12]** then **[3]**.
 4. To retransmit to Host (LIS), click **[Report]** and **[Host (HC)]** or press **[F12]** then **[1]**.




UF-1000i Maintenance

UF Daily Shutdown


Perform at the end of each day's analyses, or at least once every 24 hours if instrument is running continuously.

1. Click **[SHUTDOWN]** on the IPU main menu screen; the Shutdown dialog box will appear.



Note: To automatically turn OFF the power to the main unit once the shutdown is complete, select **YES** in the dialog box.


2. Press the manual analysis **[START]** button. The shutdown sequence will begin, with the progress displayed in the dialog box.
3. When the shutdown wash is complete (about 10 minutes) the IPU main menu screen displays "Please power off the Analyzer." To continue analysis, click **[RESTART]**. After auto-rinse and background check is completed, UF is ready.



Weekly: System power off after performing daily shutdown

Note: Once a week, when performing the daily shutdown as outlined above, select **YES** in the shutdown dialog box to turn the main unit OFF, then proceed with the steps below.

1. To completely power off UF, turn off main unit power switch.
2. To power off IPU, click **[FILE]**, **[EXIT]** on screen. Click **[YES]** to exit program.
3. Click **[START]** at bottom of Windows® desktop, then click **[SHUTDOWN]**.
4. Select "Shutdown the Computer" and press **[Enter]** or click **[OK]**.
5. Allow the computer to fully power down, then wait for at least one minute before proceeding.



UF-1000i Maintenance

Check and Drain Vacuum Trap Chamber

At the start of each day's analyses, check for the presence of water in the vacuum trap chamber; discard any liquid that may have accumulated.

Note: *If the trap chamber collects fluid on a daily basis, contact the Sysmex Technical Assistance Center.*


To discard any liquid :

1. Turn OFF the power to the main unit, and wait one minute to let the internal pressure drop.
2. Open the front cover of the analyzer.



3. Discard liquid and rinse the ball and chamber with DI water; replace the ball into the chamber, point side facing up.
4. Replace the chamber and turn clockwise to create a finger-tight seal.

Note: *Pressure and vacuum problems can arise if the chamber is not properly seated or the tubing is not properly connected.*



UF-1000i Maintenance


Clean the Sample Rotor Valve (SRV)

Clean the sample rotor valve (SRV) either monthly or every 9000 cycles. If 9000 or more samples have been analyzed since the previous cleaning, the message “Clean the SRV” will be displayed after powering on the instrument. Always reset the counter after cleaning.


1. Turn OFF the power to the main unit, and wait one minute for the internal pressure to drop.
2. Open the front cover of the instrument.
3. Turn the fixing screw counterclockwise, and remove it from the SRV mounting shaft.
4. Remove the front fixed plate, then the SRV. The valves are held together by suction and are difficult to separate. Pull the valves gently and carefully twist them off. Be careful not to lose the white washer between the valves!



Remove the fixing screw




Remove the front fixed plate



Do not lose the white washer

5. Wash the surfaces of the front and rear fixed valves, the SRV, the shaft and the washer with gauze soaked in distilled water.
6. Reassemble the valve in the reverse order from which it was removed. Turn the fixing screws clockwise to tighten them against the shaft. Do not overtighten.

Note: Do not reverse the front-rear orientation of the SRV, and make sure the metal knob is between the switching plates on the right.



Wash the surfaces with DI water

UF-1000i Maintenance

Clean the Sample Rotor Valve (SRV)

7. Remove, clean and dry the SRV tray, then place it back into position.
8. Close the front cover, and power on the main unit.
9. Verify that the background counts are within limits, and check instrument performance by running QC.
10. Reset the SRV counter according to the directions below.



Resetting the SRV Counter

1. Double-click the **[Controller]** icon from the IPU main menu, then click **[Maintenance]** then **[Counter]**. The counter dialog box will appear.



Counter - UF-1000i		
Counter	<input type="text" value="8890"/>	<input type="button" value="OK"/>
Urine Pump	<input type="text" value="12632"/>	
Sheath Syringe	<input type="text" value="24094"/>	
Charging Pump	<input type="text" value="12157"/>	
SRV	<input type="text" value="0"/>	<input type="button" value="Reset"/>

2. Press the **[Reset]** button displayed next to the SRV cycle counter to reset the count back to zero. The message "Cycle counter of the SRV will be reset"; press **[OK]** to close the message dialog box.
3. Press **[OK]** to save the changes.

UF-1000i Maintenance


Cleaning and/or Replacing Sample Filter

If the sample filter contains a large amount of debris, especially after running samples with high WBC or mucus content, the filter may become clogged. When this occurs, the **Short Sample, Check Sample Filter (Short Sample)** and/or **Check Sample Filter (WBC High)** error messages will appear even though there is sufficient sample in the tube.

Check Sample Filter (Short Sample) Error



Visually check the sample tube that caused the error to confirm adequate sample volume. Remove that sample from the processing queue, and click **[Error Recovery]** to continue.





Restart sample processing when the error recovery has completed.

UF-1000i Maintenance

Check Sample Filter (WBC High) Error

Remove the sample tube that caused the error from the processing queue. The sample filter must be replaced or cleaned to remove possible clogs from the high WBC count in the sample.


1. Click the **[Replace or Clean of Sample Filter]** button in the dialog box.
2. Click **[1. Analysis Stop]**, follow the instructions in the dialog box, then click **[OK]**.




UF-1000i Maintenance

Check Sample Filter (WBC High) Error

3. Click the "Menu" icon on the menu bar, then click **[Shutdown]**. Confirm that the power OFF radio button is set to **NO**.
4. When the Shutdown wash is complete, click **[Restart]** in the dialog box. When the background check has completed and passed, click **[2. Shutdown]**.




5. Open the front cover of the instrument. Remove the sample filter from the clip.
6. Remove the filter fixture (large, center section) from the upper and lower supporting braces (small, threaded pieces with tubing) by turning those braces to loosen.
7. Remove the thinner, top threaded fixture from the filter by turning it counter-clockwise to loosen.
8. Remove the internal (ceramic) filter from the housing and rinse with DI water. Wipe dry with a lint-free cloth or tissue.
9. Check the holes in the filter to confirm there are no clogs. Rinse again and/or replace then put the filter back into the housing.
10. Screw the thin, top threaded fixture back onto the center housing, and re-attach the filter fixture to the supporting braces until finger-tight.



UF-1000i Maintenance

Check Sample Filter (WBC High) Error

11. After replacing the filter assembly back onto the analyzer, click **[3. Replace or clean Sample filter]**, then click **[4. QC Analysis]**.



12. Perform QC analysis to verify accuracy of aspiration, absence of leaks and system performance.
13. If there is a problem with the QC results or instrument function, clean and/or replace the filter again. If the problem persists, contact the Sysmex Technical Assistance Center.

UF-1000i Reagent Replacement

If a reagent runs low during analysis, the instrument will stop automatically after completing the last analysis, and the appropriate message will display in the Help dialog box. Replace the indicated reagent with a new container/bottle.


Note: Perform an autorinse after replacing the reagent, and confirm that the background value is below the allowable limits before (re)starting sample analysis.

Message	Reagent	Abbreviation	Open Stability
UFII PACK-SED Empty Error Air Bubbles in UFII PACK-SED line	UFII PACK™ -SED	UPS	60 Days
UFII PACK-BAC Empty Error Air Bubbles in UFII PACK-BAC line	UFII PACK™-BAC (Buffered Diluent)	UPB	60 Days
UFII SEARCH-SED Empty Error Air Bubbles in UFII SEARCH-SED line	UFII SEARCH™ -SED	USS	60 Days
UFII SEARCH-BAC Empty Error Air Bubbles in UFII SEARCH-BAC line	UFII SEARCH™-BAC (Fluorescent Stain)	USB	60 Days
UFII SHEATH Empty Error	UFII SHEATH™ (Isotonic Sheath Fluid)	UTS	60 Days

Reagent Replacement

Perform the following steps in the order shown :

1. Read the error message displayed. If the Help dialog box does not display, press **[HELP]**. Click **[OK]**, Reagents Replacement Screen is displayed; select the appropriate tab for the reagent to be replaced.



UF-1000i Reagent Replacement

Reagent Replacement

- Using clean technique, remove the cap and tubing from the empty container.
- Check that the expiration date of the new reagent has not passed, and transfer the tubing to the new reagent container.



- Using handheld barcode reader, scan long reagent barcode on reagent box.

If the barcode on the reagent box cannot be read by the wand reader, it may be necessary to manually input the reagent information.

- Highlight the reagent that you are replacing.
 - Click on the LOT field and enter the lot number information.
 - Press **[Tab]** and enter the **on-board** expiration date.
- Updated reagent information is displayed (lot, open exp. date, volume).

Reagents Replacement - UF-1000i

UF II SEARCH -SED | UF II PACK -SED | UF II SEARCH -BAC | UF II PACK -BAC | Reagent Info. | < | >

Reagent Name ■ UF II SEARCH -SED

Lot No.

Exp. Day

Replace UF II SEARCH -SED and press Execute.

RUN

Cancel

- For multiple reagents, repeat steps 2-4.

- Once new reagent is on the instrument click **[RUN]** on the Reagent Replacement Screen and the reagents will automatically be primed.

UF II SHEATH Replace - UF-1000i

Please wait for a moment.
UF II SHEATH is priming.

progress

80%

UF-1000i Reagent Replacement

Reagent Log

Reagent replacement is automatically documented on **Reagent Log**.

- From the IPU main menu screen, click **[Controller]**, then **[Reagent Log]**.
- Comments can be entered for each reagent replaced by double clicking on the line entry and entering the desired comment.





IPU - [Reagent Log Term:All Reagent:All]							
File(F) Edit(E) View(V) Record(R) Action(A) Report(P) Setting(S) Window(W) Help(H) Ver.:00-07_Build51D User Name:sysmex							
F1 Help F2 MANUAL F3 SAMPLER F4 Menu F5 QC Files F6 Work list F7 Explorer F8 Browser F9 Input F12 Out Upper Low							
No.	Date	Time	Logon User	Reagent	Lot No.	Exp. Date	Comments
44	04/20/2007	12:28:06	sysmex	UF II SHEATH	A7001	02/06/2008	
43	04/20/2007	12:25:30	sysmex	UF II PACK -BAC	A6005	01/09/2008	
42	04/20/2007	12:22:54	sysmex	UF II SEARCH -BAC	A6005	01/12/2008	
41	04/20/2007	12:19:48	sysmex	UF II PACK -SED	A6005	01/12/2008	
40	04/20/2007	12:16:38	sysmex	UF II SEARCH -SED	A7001	01/26/2008	
39	02/28/2007	15:29:37	sysmex	UF II SHEATH	A6006	10/06/2007	did not change reagent- prime only
38	02/09/2007	15:25:43		UF II PACK -SED	A6001	04/04/2007	
37	02/05/2007	13:37:33	sysmex	UF II SEARCH -BAC	A6001	03/30/2007	
36	01/18/2007	14:23:24	sysmex	UF II SEARCH -SED	A6001	03/29/2007	
35	01/17/2007	08:52:53	sysmex	UF II SHEATH	A6006	10/06/2007	
34	01/17/2007	08:51:32	sysmex	UF II PACK -BAC	A6002	09/08/2007	
33	01/10/2007	15:47:36	sysmex	UF II PACK -BAC	A6001	04/03/2007	
32	01/10/2007	15:15:55	sysmex	UF II SHEATH	A6001	03/29/2007	
31	01/10/2007	15:09:35	sysmex	UF II SHEATH	A6001	03/29/2007	
30	01/10/2007	11:52:25	sysmex	UF II SHEATH	A6001	03/29/2007	
29	01/09/2007	12:31:35	sysmex	UF II SHEATH	A6002	07/12/2007	
28	01/09/2007	12:30:16	sysmex	UF II SHEATH	A6002	07/12/2007	
27	01/09/2007	12:27:19	sysmex	UF II SHEATH	A6002	07/12/2007	
26	01/09/2007	08:22:45	sysmex	UF II PACK -SED	A6001	04/04/2007	
25	01/08/2007	10:42:45	sysmex	UF II SHEATH	A6002	07/12/2007	
24	01/05/2007	09:51:25	sysmex	UF II SEARCH -BAC	A6001	03/30/2007	
23	01/04/2007	10:20:18	sysmex	UF II PACK -BAC	A6001	04/03/2007	
22	01/04/2007	07:58:44	sysmex	UF II SEARCH -SED	A6001	03/29/2007	
21	01/02/2007	13:23:09	sysmex	UF II SHEATH	A6002	07/12/2007	
20	12/28/2006	09:19:49	sysmex	UF II PACK -BAC	A6001	04/03/2007	

UF-1000i Host Downtime

To Disconnect Host When LIS is Down

1. When the HOST is down, a pop-up error message displays on the IPU. The HOST icon in the lower right area of the IPU screen will be red.
2. Click **[OK]** in the IPU error box to close the message.
3. Go to the **IPU Settings, [HOST(HC)] Setting**.
4. Click the “**Host (HC) Connect**” to remove the check mark. This disconnects the IPU from the LIS.
5. Click **[OK]** to save the changes.



To Re-connect the Host When LIS is Up

1. Go to the **IPU Settings, [HOST (HC)] Setting**.
2. Click the empty box next to “**Host (HC) Connect**” to reconnect IPU and Host.
3. Click **[OK]** to save the changes.
4. The green HOST icon displays. The UF is now reconnected to the HOST.

