
erbe
power your performance.

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1. 一次性软性冷冻探针的开发文件: D144191
2. Schumann C, Hetzel M, Babiak AJ et al. Endobronchial tumor debulking with a flexible cryoprobe for immediate treatment of malignant stenosis. *J Thorac Cardiovasc Surg* 2010; 139: 997–1000
3. Hetzel J, Eberhardt R, Herth FJF et al. Cryobiopsy increases the diagnostic yield of endobronchial biopsy: a multicentre trial. *Eur Respir J* 2012; 39: 685–690
4. Schumann C, Kropf C, Rudiger S et al. Removal of an aspirated foreign body with a flexible cryoprobe. *Respir Care* 2010; 55: 1097–1099
5. Vergnon JM, Mathur PN. Cryotherapy for Endobronchial Disorders. *Interventional Bronchoscopy* 2000; 30: 133–145
6. 内部测试报告: D147560
7. Fruchter O, Kramer MR. Retrieval of various aspirated foreign bodies by flexible cryoprobe: in vitro feasibility study. *Clin Respir J* 2015; 9: 176–179
8. Cryoablation, Mechanisms and influencing factors (Dr. Christiane Nerz, 2018): D135896
9. Tomic R, Podgaetz E, Andrade RS et al. Cryotechnology in diagnosing and treating lung diseases. *J Bronchology Interv Pulmonol* 2015; 22: 76–84
10. Vergnon J-M, Huber RM, Moghissi K. Place of cryotherapy, brachytherapy and photodynamic therapy in therapeutic bronchoscopy of lung cancers. *Eur Respir J* 2006; 28: 200–218
11. Lee S-H, Choi W-J, Sung S-W et al. Endoscopic cryotherapy of lung and bronchial tumors: a systematic review. *Korean J Intern Med* 2011; 26: 137–144
12. 一次性柔性冷冻探针的使用提示: 30402-406
13. Babiak A, Hetzel J, Krishna G et al. Transbronchial cryobiopsy: a new tool for lung biopsies. *Respiration* 2009; 78: 203–208
14. Fruchter O, Fridel L, Rosengarten D et al. Transbronchial cryo-biopsy in lung transplantation patients: first report. *Respirology* 2013; 18: 669–673
15. Yarmus L, Akulian J, Gilbert C et al. Cryoprobe transbronchial lung biopsy in patients after lung transplantation: a pilot safety study. *Chest* 2013; 143: 621–626
16. 一次性柔性冷冻探针的临床评估: D104429
17. Aktas Z, Gunay E, Hoca NT et al. Endobronchial cryobiopsy or forceps biopsy for lung cancer diagnosis. *Ann Thorac Med* 2010; 5: 242–246
18. Chou C-L, Wang C-W, Lin S-M et al. Role of flexible bronchoscopic cryotechnology in diagnosing endobronchial masses. *Ann Thorac Surg* 2013; 95: 982–986
19. Jabari H, Sami R, Fakhri M et al. Different protocols for cryobiopsy versus forceps biopsy in diagnosis of patients with endobronchial tumors. *Pneumologia* 2012; 61: 230–233
20. 带有一次性柔性冷冻探针的 ERBECRYO® 2 的宣传册: 85402-000
21. Hetzel J, Maldonado F, Ravaglia C et al. Transbronchial Cryobiopsies for the Diagnosis of Diffuse Parenchymal Lung Diseases: Expert Statement from the Cryobiopsy Working Group on Safety and Utility and a Call for Standardization of the Procedure. *Respiration* 2018; 95: 188–200
22. Raghu G, Remy-Jardin M, Myers JL et al. Diagnosis of Idiopathic Pulmonary Fibrosis. An Official ATS/ERS/JRS/ALAT Clinical Practice Guideline. *Am J Respir Crit Care Med* 2018; 198: e44–e68
23. Casoni GL, Tomassetti S, Cavazza A et al. Transbronchial lung cryobiopsy in the diagnosis of fibrotic interstitial lung diseases. *PLoS One* 2014; 9: e86716
24. Yilmaz A, Aktas Z, Alici IO et al. Cryorecanalization: keys to success. *Surgical Endoscopy* 2012; 26: 2969–2974
25. Hetzel M, Hetzel J, Schumann C et al. Cryorecanalization: a new approach for the immediate management of acute airway obstruction. *J Thorac Cardiovasc Surg* 2004; 127: 1427–1431
26. Folch E, Mehta AC. Airway interventions in the tracheobronchial tree. *Semin Respir Crit Care Med*. 2008; 29: 441–452
27. Ernst A, Silvestri GA, Johnstone D. Interventional pulmonary procedures: Guidelines from the American College of Chest Physicians. *Chest* 2003; 123: 1693–1717
28. Lee H, Leem CS, Lee JH et al. Successful removal of endobronchial blood clots using bronchoscopic cryotherapy at bedside in the intensive care unit. *Tuberc Respir Dis (Seoul)* 2014; 77: 193–196
29. Sehgal IS, Dhooria S, Agarwal R et al. Use of a Flexible Cryoprobe for Removal of Tracheobronchial Blood Clots. *Respir Care* 2015; 60: e128–31
30. Weerdt SD, Noppen M, Remels L et al. Successful Removal of a Massive Endobronchial Blood Clot by Means of Cryotherapy. *Journal of Bronchology* 2005; 12: 23–24
31. Hammer J, Trachsel D, Nicolai T et al. Caution to use bronchoscopic CO cryotherapy for foreign body removal in children. *Pediatr Pulmonol* 2016; 51: 889–891
32. Ernst A, Feller-Kopman D, Becker HD et al. Central airway obstruction. *Am J Respir Crit Care Med* 2004; 169: 1278–1297
33. Sheski FD, Mathur PN. Endobronchial electrosurgery: argon plasma coagulation and electrocautery. *Semin Respir Crit Care Med* 2004; 25: 367–374
34. Zenker M. Argon plasma coagulation. *GMS Krankenhhyg Interdiszip* 2008; 3: Doc15
35. Tremblay A, Marquette C-H. Endobronchial electrocautery and argon plasma coagulation: a practical approach. *Can Respir J* 2004; 11: 305–310
36. Bolliger CT, Sutedja TG, Strausz J et al. Therapeutic bronchoscopy with immediate effect: laser, electrocautery, argon plasma coagulation and stents. *Eur Respir J* 2006; 27: 1258–1271
37. Sutedja TG, van Boxem TJ, Schramel FM et al. Endobronchial Electrocautery is an Excellent Alternative for Nd: YAG Laser to Treat Airway Tumors. *J Bronchology Interv Pulmonol* 1997; 4
38. Ernst A, Anantham D. Update on interventional bronchoscopy for the thoracic radiologist. *J Thorac Imaging* 2011; 26: 263–277
39. Reichle G, Freitag L, Kullmann HJ et al. Argon plasma coagulation in bronchology: a new method--alternative or complementary? *Pneumologie* 2000; 54: 508–51
40. Morice RC, Ece T, Ece F et al. Endobronchial argon plasma coagulation for treatment of hemoptysis and neoplastic airway obstruction. *Chest* 2001; 119: 781–787
41. 等离子外科基础知识的用户手册: 85800-038
42. VIO 3 使用说明书: D140792
43. Franke K-J, Szyrach M, Nilius G et al. Experimental study on biopsy sampling using new flexible cryoprobes: influence of activation time, probe size, tissue consistency, and contact pressure of the probe on the size of the biopsy specimen. *Lung* 2009; 187: 253–259
44. Amat B, Esselmann A, Reichle G et al. The electrosurgical knife in an optimized intermittent cutting mode for the endoscopic treatment of benign web-like tracheobronchial stenosis. *Arch Bronconeumol* 2012; 48: 14–21
45. APC 探针的临床评估: D099576 v003
46. Jin F, Mu D, Xie Y et al. Application of bronchoscopic argon plasma coagulation in the treatment of tumorous endobronchial tuberculosis: Historical controlled trial. *J Thorac Cardiovasc Surg* 2013; 145: 1650–1653
47. Lentz RJ, Argento AC, Colby TV et al. Transbronchial cryobiopsy for diffuse parenchymal lung disease: A state-of-the-art review of procedural techniques, current evidence, and future challenges. *J Thorac Dis* 2017; 9: 2186–2203

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息肉切除套索, © medwork
Erbe China Ltd.