CellSafe 1 Biopsy Cassette Laser

EBG-04**-03B

EC	REP	EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands
----	-----	---

CellPath 
INNOVATION IN CELLULAR PATHOLOGY

 **CellPath Ltd**
80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com

CellSafe 1 Biopsy Cassette Laser

EBG-04**-03B

A small device intended to be used in the laboratory to support/contain clinical tissue specimens (e.g., biopsy samples) to facilitate specimen processing [e.g., fixation (encapsulation in paraffin wax), dehydration, infiltration] in preparation for subsequent cytological or histological examination, and storage. It is in the form of a synthetic polymer container, with or without a lid, with drainage holes to maximize fluid exposure during tissue submersion in a parent device into which the container is loaded. This is a single-use device.

For laboratory use.

Store in a dry place. Don't store in direct sunlight.

This product does not expire under normal environmental conditions.

Processing Limitations: Suitable for high temperature processing. Suitable for microwave processing.

Marking capabilities: Laser printers, CellMark pens, 2B pencils.

Instructions for Use

STANDARD CASSETTE USE:

1. Insert the thread into the printer. Holding the top cassette, pull the thread up to release.
2. Place tissue into the CellSafe I Cassette. We recommend that only one biopsy of a maximum 4.5mm thick is placed in the cassette. Please ensure tissue is >400µm.
3. Using the lid provided, attach the lid. Close the lid and place the cassette in a fixative ensuring it is fully submerged.

NOTE: Immersion in acidic decalcifying fluids may have an effect upon printed text and barcodes therefore, immersion should be kept to the minimum time possible.

4. Process the cassette through solvents according to your usual protocol until the specimens are ready for embedding. If you are using any unusual reagents, we suggest that you carry out preliminary chemical compatibility tests before using CellSafe I to contain tissue samples.

5. Carefully open the cassette.
6. Embed, cool and section as per your usual protocol.
7. Dispose of the lid in accordance with laboratory guidelines.

DO NOT REUSE THIS PRODUCT.

BIOPSIES PLACED DIRECTLY IN CELLSAFE I:

1. Place tissue into the CellSafe I Cassette. We recommend that only one biopsy of a maximum 4.5mm thick is placed in the cassette. Please ensure tissue is >400µm.

NOTE: To prevent drying of the sample, it is recommended that a tiny amount of fixative is applied to the sample while samples are being collected.

2. Using the lid provided, attach and the lid. Close the lid and place the cassette in a fixative ensuring it is fully submerged.
3. Transfer the pot to the laboratory for processing, if applicable.

NOTE: The CellSafe I is x-ray transparent.

STRAINING TISSUE FRAGMENTS:

1. Presoak the CellSafe I in ethanol or formaldehyde.
2. Ensure the fragments are >400µm.
3. Hold or support the cassette horizontally and slowly pour the suspension through the CellSafe I.
4. The tissue fragments will be retained in the CellSafe I, when all the suspension has been filtered, the CellSafe I can be closed.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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