
# RQ 800 VIPER

6064 McKee Road  
Fitchburg, WI 53719 USA  
Toll Free: 1 (877) 270-1250  
Phone: 1 (608) 270-1230

E-mail: [service@genepro-inc.com](mailto:service@genepro-inc.com)

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---

## TABLA DE REVISIONES

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### PRÓLOGO

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The screenshot shows the 'Batch Log' screen for the RQ800 'VIPER' Packaging System. The title bar indicates the date 'Thu 07/02/15' and the time '12:12:37'. The main content area is titled 'Batch Log' and 'File'. A dropdown menu in the top right corner shows '56'. Below this, a table displays the following information:

Batch ID	56	Qty Requested	0
Boar ID	5	Qty Processed	23
Field 1		Expiration	2
Field 2		Volume	1610
Field 3		Error Cnt	4
Field 4		Calibration Cnt	0
Field 5		Start Time	01/01/1970 00:00:00
Field 6		End Time	07/01/2015 13:32:22
Field 7			

At the bottom of the screen, there are navigation buttons: 'Main', 'Alarm', 'Data', 'Diagnostics', 'Trend', and 'Setup'. A 'Save Log' button is located in the bottom right corner. The bottom right corner also displays 'Logged in as: Administrator' and a temperature indicator '191°C'.

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The screenshot shows the configuration interface for the RQ800 'VIPER' Packaging System. The title bar displays 'Tue 06/30/15', 'RQ800 "VIPER" Packaging System', and '11:51:16'. The interface is divided into three main sections: 'Settings', 'Communication', and 'File'. The 'Settings' section is currently active and displays 'Label Field Mappings'. This section contains two columns of dropdown menus for field assignments (FN0-FN17) and a right-hand column for label format settings. The 'File' section includes a 'Set Logo' button. At the bottom, there are navigation buttons for 'Main', 'Alarm', 'Data', 'Diagnostics', and 'Trend', along with a temperature display showing '199°C' and a 'Setup' button. The user is logged in as 'Supervisor'.

Field	Assignment	Field	Assignment	Setting
FN0:	Boar ID	FN9:	Exp Date	Label Format: Default
FN1:	Field 1	FN10:	Pack Date	Batch Label: Batch0
FN2:	Field 2	FN11:	Exp 1000	User Format 0: user0
FN3:	Not Used	FN12:	Pack 1000	User Format 1: user1
FN4:	Not Used	FN13:	Company	User Format 2: user2
FN5:	Not Used	FN14:	Address 1	User Format 3: user3
FN6:	Not Used	FN15:	Address 2	User Format 4: user4
FN7:	Not Used	FN16:	Telephone	Logo Name: [ ]
FN8:	Batch ID	FN17:	License	Date Format: MM/DD/YY
				Thousand Offset: 648
				Thousand Value: 984

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The screenshot displays the network configuration interface for the RQ800 "VIPER" Packaging System. At the top, the system name and time (13:50:23) are shown. The interface is divided into sections: "Settings", "Network", and "File". The "Data Setup" section contains the following fields:

- IP Address: 0 . 0 . 0 . 0
- Share: [Empty text box]
- User Name: [Empty text box]
- Password: [Empty text box]
- Domain: [Empty text box]
- Host: [Empty text box]
- Protocol: [Empty text box]

An "Enable Network" button is located to the right of the IP address field. The bottom navigation bar includes buttons for "Main", "Alarm", "Data", "Diagnostics", "Trend", and "Setup". A temperature display shows "191°C" and the login status "Logged in as: Administrator".

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Main Alarm Data Diagnostics Trend 190°C Logged in as: Supervisor Setup

Tue 06/16/15 **RQ800 "VIPER" Packaging System** 16:29:40

System	Rack I/O 1	Rack I/O 2	ACOPOS	Remote I/O
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p><b>X20 Digital In 1</b></p> <p>Module ID: 7058</p> <p>Serial Number: 498055</p> <p>Hardware: 5</p> <p>Firmware: 801</p> <p><b>Present</b></p> </div> <div style="width: 45%;"> <p><b>Outside Cutter</b></p> <p>Retract: <input checked="" type="radio"/> Extend: <input type="radio"/></p> <p>Up: <input checked="" type="radio"/> Down: <input type="radio"/></p> </div> </div>				
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p><b>X20 Digital In 2</b></p> <p>Module ID: 7058</p> <p>Serial Number: 498053</p> <p>Hardware: 5</p> <p>Firmware: 801</p> <p><b>Present</b></p> </div> <div style="width: 45%;"> <p>Air: <input checked="" type="radio"/> Vac: <input type="radio"/></p> <p>MICO: <input checked="" type="radio"/> PS: <input checked="" type="radio"/></p> </div> </div>				
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p><b>X20 Digital In 3</b></p> <p>Module ID: 7058</p> <p>Serial Number: 498054</p> <p>Hardware: 5</p> <p>Firmware: 801</p> <p><b>Present</b></p> </div> <div style="width: 45%;"> <p>Fluid: <input checked="" type="radio"/> 2: <input type="radio"/></p> <p>3: <input type="radio"/> 4: <input type="radio"/></p> </div> </div>				

Main Alarm Data Diagnostics Trend 22 °C Logged in as: Administrator Setup

Tue 06/16/15 **RQ800 "VIPER" Packaging System** 16:30:12

System	Rack I/O 1	Rack I/O 2	ACOPOS	Remote I/O
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p><b>X20 Relay</b></p> <p>Module ID: 9851</p> <p>Serial Number: 171491</p> <p>Hardware: 5</p> <p>Firmware: 816</p> <p><b>Present</b></p> <p>Zero Crossing Status: <input type="radio"/></p> </div> <div style="width: 45%;"> <p><b>X20 DO 1</b></p> <p>Module ID: 7064</p> <p>Serial Number: 393692</p> <p>Hardware: 5</p> <p>Firmware: 819</p> <p>Output Status: <input checked="" type="radio"/> 2: <input type="radio"/> Accp: <input checked="" type="radio"/> Air: <input checked="" type="radio"/> MICO: <input type="radio"/> PS: <input type="radio"/></p> </div> </div>				
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p><b>X20 DO 2</b></p> <p>Module ID: 7064</p> <p>Serial Number: 393434</p> <p>Hardware: 5</p> <p>Firmware: 819</p> <p>Output Status: <input type="radio"/> Gr: <input type="radio"/> Or: <input type="radio"/> Red: <input type="radio"/> Buzz: <input type="radio"/> 5: <input type="radio"/> 6: <input type="radio"/></p> </div> <div style="width: 45%;"> <p><b>Bus Transmitter</b></p> <p>Module ID: 41528</p> <p>Serial Number: 173654</p> <p>Hardware: 6</p> <p>Firmware: 784</p> <p>Bus Voltage: 19.8 V</p> <p>Bus Current: 1.1 A</p> <p><b>Present</b> Bus Supply: <b>Supply</b> I/O Supply: <b>Supply</b></p> </div> </div>				

Main Alarm Data Diagnostics Trend 22 °C Logged in as: Administrator Setup

Tue 06/16/15 **RQ800 "VIPER" Packaging System** 16:30:40

System    Rack I/O 1    Rack I/O 2    **ACOPOS**    Remote I/O

**ACOPOS Drive**

Index		Pump	
Velocity	0.0	Velocity	0.0
Position	0	Position	0
Torque	-16.79	Torque	1.85
Initialized	<span style="color: green;">●</span>	Initialized	<span style="color: green;">●</span>
Ready	<span style="color: green;">●</span>	Ready	<span style="color: green;">●</span>
Status	<span style="color: green;">●</span>	Status	<span style="color: green;">●</span>
Homing OK	<span style="color: green;">●</span>	Homing OK	<span style="color: green;">●</span>
Error	<span style="color: green;">●</span>	Error	<span style="color: green;">●</span>
Error	0	Error	0
Enabled	<span style="color: green;">●</span>	Enabled	<span style="color: green;">●</span>
Trigger1	<span style="color: red;">●</span>		

Index Encoder
Pump Encoder

**XV Module** Present

Module ID: 41830 Pwr

Serial Number: 170609

Hardware: 5

Firmware: 771 OvrLd

CtEx 2

3

4

BIEx 6

6

CtDn

CtUp

FIDn

FIUp

Seal

12

LaRo 14

LaEx 16

Vacu 18

Blow 20

21

22

23

24

Tue 06/16/15 **RQ800 "VIPER" Packaging System** 16:30:57

System	Rack I/O 1	Rack I/O 2	ACOPOS	Remote I/O
<b>X67 MM 1</b> <span style="float:right">Present</span> Module ID: 4881 Serial Number: 276500 Hardware: 6 Firmware: 912 Blade R <input checked="" type="checkbox"/> Blade E <input type="checkbox"/> Cut Up <input checked="" type="checkbox"/> Cut Down <input type="checkbox"/> ICut R <input checked="" type="checkbox"/> ICut E <input type="checkbox"/> Bag Tip <input type="checkbox"/> Bag Pos <input type="checkbox"/> Bag Tip Latch <input type="checkbox"/> Ack <input type="checkbox"/>	<b>X67 MM 2</b> <span style="float:right">Present</span> Module ID: 4881 Serial Number: 276534 Hardware: 6 Firmware: 912 Pos Cut <input type="checkbox"/> Pos Fill <input type="checkbox"/> Seal Ret <input checked="" type="checkbox"/> Seal Ext <input type="checkbox"/> Lbl Down <input checked="" type="checkbox"/> Lbl Up <input type="checkbox"/> Lbl Ret <input checked="" type="checkbox"/> Lbl Ext <input type="checkbox"/>	<b>X67 Interface</b> <span style="float:right">Present</span> Module ID: 43279 Serial Number: 169534 Hardware: 1 Firmware: 2 Reset <input type="button"/> Reset <input type="button"/> Reset <input type="button"/> Reset <input type="button"/> 1 <input type="button"/> 2 <input type="button"/> 3 <input type="button"/> 4 <input type="button"/> <input type="button"/> <input type="button"/>		

22 °C
Logged in as: Administrator

## 10. APPENDIX

### 12.1 Axis error Reference Guide

The following table is a list of axis error codes to reference in the case of pump or index errors.

<i>Error ID</i>	<i>Error Text</i>	<i>Notes</i>
1	Invalid parameter ID	
2	Data block for upload is not available	
3	Write access for a read-only parameter	
4	Read access for a write- only parameter	
8	Data block read access already initialized	
9	Data block write access already initialized	
10	Data block read access not initialized	
11	Data block write access not initialized	
16	The data segment is already the last when reading the data block	
17	The data segment is already the last when writing the data block	
18	The data segment is not yet the last when reading the data block	
19	The data segment is not yet the last when writing the data block	
21	Check sum after data block write is invalid	
23	Parameter ID in data block is invalid (data block write)	
25	Burn system module only allowed immediately after download	
27	Operating system not able to be started (operating system is not on the FEPROM)	
40	Value of parameter higher than maximum value	
41	Value of parameter higher than maximum value	
42	Value of parameter higher than maximum value	
52	Value of parameter lower than minimum value	
53	Value of parameter lower than minimum value	
54	Value of parameter lower than minimum value	
64	Hardware ID in BR module is invalid (data block write)	
65	Hardware version in BR module is invalid (data block write)	
66	The operating system on the drive is incompatible to the existing network	
67	Necessary parameter is missing or is invalid	
68	Data block length invalid	
69	Command interface is occupied	
70	Value of a necessary parameter too high	
71	Value of a necessary parameter too low	

1001	Error-FIFO overflow	
1002	Parameter outside the valid range	

<i><b>Error ID</b></i>	<i><b>Error Text</b></i>	<i><b>Notes</b></i>
1003	Parameter cannot be written while loop control is active	
1004	Timeout in network life sign monitor	
1005	Parameter cannot be written while a movement is active	
1006	Invalid parameter for trigger event (digital input + edge)	
1007	Master for network coupling deactivated - one master is already sending	
1008	Master for network coupling deactivated - Encoder error	
1009	Error during memory allocation	
1011	Quickstop input active	
1012	Breakdown of cyclic network communication	
1013	Station is not available for network communication	
1014	Network command interface is occupied	
1016	Maximum cycle time exceeded - CPU load too high	
1017	Invalid parameter ID for cyclic read access	
1018	Invalid parameter ID for cyclic write access	
1021	Parameter cannot be written: Function block active	
1022	Timeout in life sign monitoring of cyclic data to drive	
1023	Network coupling with the cyclic communication mode not allowed	
1024	Cyclic communication mode with current network configuration not possible	
2001	Upload of trace data not allowed: Recording active	
2003	Trace start not allowed: Recording active	
2006	Initialization of trace parameters not allowed: Recording active	
4005	Controller cannot be switched on: Drive in error state	
4007	Lag error stop limit exceeded	Motor slipped during motion, check for drag
4008	Positive limit switch reached	
4009	Negative limit switch reached	
4010	Controller cannot be switched on: Both limit switches are closed	
4011	Controller cannot be switched off: Movement active	
4012	Controller cannot be switched on: Init parameters missing or not valid	
4014	Two encoder control: Stop limit of positions difference exceeded	
5001	Target position exceeds positive SW limit	
5002	Target position exceeds negative SW limit	
5003	Positive SW limit reached	
5004	Negative SW limit reached	

5005	Start of movement not possible: Position controller inactive	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
5006	Start of movement not possible: Axis not referenced	
5010	Move in pos. direction not possible: Pos. limit switch is closed	
5011	Move in neg. direction not possible: Neg. limit switch is closed	
5012	Start of movement not possible: Stop ramp active	
5013	Cyclic set value mode cannot be switched on: Movement active	
5015	Start of movement not possible: Homing procedure active	
5016	Parameter cannot be written: Homing procedure active	
5017	Homing procedure mode not possible: Position controller inactive	
5018	Homing procedure not possible: Movement active	
5019	Homing parameter outside the valid range	
5020	Homing procedure not possible: Both limit switches are closed	
5021	Limit switch closed: No direction change for this homing mode	
5022	Second limit switch signal received: Reference switch not found	
5023	Incorrect limit switch signal received for current movement direction	
5024	Cyclic set value mode aborted: Set positions missing	
5025	Homing offset with counting range correction cannot be set	
5026	Basis movement parameters (with override) exceed speed limit value	
5027	Basis movement parameters (with override) exceed acceleration limit value	
5028	Current movement is not a basis movement	
5029	Trigger ignored - remaining distance exceeds SW limit	
5030	Homing procedure mode not possible: Position controller active	
5031	Homing procedure mode not possible: Cyclic set values mode is off	
5032	Acceleration too low - braking distance exceeds positive SW limit	
5033	Acceleration too low - braking distance exceeds negative SW limit	
5034	Homing procedure not possible: Encoder error	
5035	Reference marks not detected	
5036	Acceleration stop limit exceeded	
5101	Cam profile compensation gears: Limit values exceeded	
5102	Too many changes of cam profile per cycle (master period too short)	

5103	Slave trigger FIFO full	
5104	Slave trigger FIFO empty	
5105	Master trigger FIFO full	
5106	Master trigger FIFO empty	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
5107	Cam coupling cannot be started: Parameter outside the valid range	
5108	Master compensation trigger FIFO full	
5109	Master compensation trigger FIFO empty	
5110	Cam coupling aborted: Cyclic set positions missing	
5111	Cam coupling aborted: Encoder error	
5112	Command not allowed: Cam profile coupling not active	
5113	Command not allowed: Controller is already active	
5114	Parameter cannot be written: Cam coupling active	
5115	Restart command not possible: The cam automat is not active	
5201	Parameter cannot be written: Drumsequencer active	
5202	Cam Control: Switch positions not in ascending order	
5300	Data block for upload is not available	
5301	Start cam automat linkage not possible: Parameter outside the valid range	
5302	Parameter cannot be written: Cam automat active	
5303	Cam profile data not available at index	
5304	Format error in cam profile data	
5311	Cam automat: Event leads to non initialized state	
5315	Download error: Cam profile data in use by cam automat or function block	
5316	Event type is not possible for entry in compensation gears	
5319	Cam profile data not allowed for state 0	
5329	No valid cam profile data or state deactivated	
6000	Master sampling time is not a multiple of position controller sampling time	
6001	Sync controller: Timeout for sync telegram occurred	
6002	Sync controller: Error tolerance of system time difference exceeded	
6008	Controller is already active	
6011	Controller is not in speed mode	
6014	Drive initialization active	
6015	CAN controller: CAN bus disturbance (receive error counter greater 96)	
6016	CAN controller: CAN bus disturbance (transmit error counter greater 96)	
6017	Software: Watchdog active	
6018	Hardware: 15V power supply fail	
6019	ACOPOS: Overcurrent	

6020	Hardware: 24V power supply fail	
6021	Low level at controller enable input	
6023	Voltage sag at controller enable input	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
6024	Current was latched before conversion (OpSys error in ABLs)	
6025	Temperature was latched before conversion (OpSys error in ABLs)	
6026	Holding brake: Stator current limit exceeded during release	
6027	Holding brake: Manual operation not permitted	
6028	Holding brake: Undervoltage/-current (wire breakage, check 24V supply)	
6029	Holding brake: Control signal on and output status off	
6030	Holding brake: Brake output is active, but no brake entered in motor data	
6031	System module already deleted	
6032	Interface: FPGA configuration error	
6033	Type of servo amplifier is not supported by ACOPOS-firmware	
6034	Cyclic set value mode aborted: Set speeds missing	
6036	Motor parameters missing or invalid	
6038	Torque limit higher than peak motor torque	
6040	Operating system version is less than allowed minimum version	
6041	Operating system version is greater than allowed maximum version	
6042	Operating system version is not in the allowed range	
6043	PHASING_MODE is not valid	
6044	Phasing: Rotational direction or position not valid	
6045	Power stage: Connection X5: No current flow	
6046	Phasing: No rotor movement	
6047	Holding brake: Control signal off and output status on	
6048	Motor holding brake movement monitor: Position error too large	
6049	Power stage: Connection X5: Current measurement faulty	
6050	Write parameter not allowed: Set current filter or notch filter active	
6051	Phasing: Speed too high	
6052	Power stage: High-side: Overcurrent	
6053	Power stage: Low-side: Overcurrent	
6054	Power stage: Overcurrent	
6055	Holding brake: Low voltage	
6056	Holding brake: Low current	
6057	Position loop controller: Load encoder error	
6058	Enable1: Voltage sag	

6059	Enable2: Voltage sag	
6060	Power stage: Limit speed exceeded	
6061	CTRL Speed controller: Limit speed exceeded	
6062	CTRL Speed controller: Speed error stop limit exceeded	
6063	Holding brake: External voltage on output over 24V	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
7012	Encoder: Hiperface error bit	
7013	Encoder: Status message	
7014	Encoder: CRC error during parameter transfer	
7015	Encoder: Timeout error during parameter transfer	
7016	Encoder: Busy error during parameter transfer	
7017	Encoder: Error while reading encoder parameter	
7020	OEM data: Data write error	
7021	Encoder: Timeout error while reading absolute position	
7022	Encoder: Initialization is active	
7023	Encoder: Parameter transfer is active	
7029	Encoder: Incremental signal amplitude too small	
7030	Encoder: Incremental signal amplitude too large	
7031	Encoder: Incremental signal amplitude too large (Disturbance)	
7032	Encoder: Incremental signal amplitude too small (Disturbance, no connection)	
7033	Encoder: Incremental position step too large	
7036	Encoder: Interface ID invalid (Check slot and Interface EEPROM data)	
7038	Encoder: Position value not synchronous with absolute value	
7039	Incremental encoder: Cable disturbance track A	
7040	Incremental encoder: Cable disturbance track B	
7041	Incremental encoder: Cable disturbance track R	
7042	Incremental encoder: Edge distance of quadrature signal too small	
7043	Encoder: Cable disturbance track D	
7044	Encoder: Parity	
7045	Resolver: Signal disturbance (plausibility check)	
7046	Resolver: Cable disturbance	
7047	Invalid distance of reference marks	
7048	Error during the reading of encoder memory	
7049	Abnormal encoder current consumption	
7050	Incremental encoder: Illegal AB signal change	
7051	Encoder: Acceleration too large (Disturbance)	
7052	Encoder: Encoder is not Supported	
7100	Parameter function not supported. (Module ?)	

7200	DC bus: Overvoltage	
7210	DC bus: Charging: Voltage unstable	
7211	DC bus: Voltage drop	
7212	DC bus: High voltage drop	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
7214	DC bus: Charging resistor hot (too many power line fails)	
7215	Power mains: At least one phase of the power line failed	
7217	DC bus: Nominal voltage detection: Voltage too high	
7218	DC bus: Nominal voltage detection: Voltage too low	
7219	DC bus: Charging: Voltage too low	
7220	DC bus: Nominal voltage detection: Voltage not allowed	
7221	Mains: Failure	
7222	Power stage: Connection X5: Ground fault	
7223	DC bus: Overvoltage DC-GND	
7224	Connector to back plane: 24V-GND contact monitoring: Voltage too low	
7225	DC bus: Overvoltage	
7226	DC bus: Overcurrent	
7227	Bleeder: Overcurrent	
7300	Digital IO: IO Configuration invalid	
7401	Parameter position exceeds maximum data length	
7402	Processing of parameter sequence aborted: Write error	
7403	Processing of parameter sequence is still active	
7404	Parameter sequence not available at index	
8001	EEPROM select not valid	
8003	Table index not valid	
8004	EEPROM variable type not valid	
8005	EEPROM type not valid	
8006	Value of EEPROM parameter is zero	
8007	Value of EEPROM parameter is not valid	
8020	Invalid switch frequency	
9000	Heatsink temperature sensor: Overtemperature - Movement stopped	
9001	Heatsink temperature sensor: Overtemperature - Limiter active	
9002	Heatsink temperature sensor: Not connected or destroyed	
9003	Heatsink temperature sensor: Not connected or destroyed	
9010	Temperature sensor (Motor Choke External): Overtemperature	
9011	Temperature sensor (Motor Choke External): Overtemperature - Limiter active	

9012	Temperature sensor (Motor Choke External): Not connected or destroyed	
9013	Temperature sensor (Motor Choke External): Short circuit	
9030	Junction temperature model: Overtemperature - Movement stopped	
9031	Junction temperature model: Overtemperature - Limiter active	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
9040	Bleeder temperature model: Overtemperature - Movement stopped	
9041	Bleeder temperature model: Overtemperature - Limiter active	
9050	ACOPOS peak current: Overload - Movement stopped	
9051	ACOPOS peak current: Overload - Limiter active	
9060	ACOPOS continuous current: Overload - Movement stopped	
9061	ACOPOS continuous current: Overload - Limiter active	
9070	Motor temperature model: Overload - Movement stopped	
9071	Motor temperature model: Overload - Limiter active	
9075	ACOPOS continuous power: Overload - Movement stopped	
9076	ACOPOS continuous power: Overload - Limiter active	
9078	ACOPOS temperature sensor: Overtemperature - Movement stopped	
9079	ACOPOS temperature sensor: Overtemperature - Limiter active	
9300	Current controller: Overcurrent	
10000	Identification parameter(s) missing	
10001	Parameter identification: Invalid sub-mode	
10100	Parameter identification: Quality factor not fulfilled	
10101	No ISQ-filter free	
10102	No resonance-frequency for ISQ-filter (band-stop) found	
29200	The axis object is invalid	
29203	Drive is not ready	
29204	Invalid parameter number	
29205	The axis is not homed	
29206	The controller is off	
29207	This movement type is currently not allowed	
29208	The axis object was changed since last FB call	
29209	The drive is in error state	
29210	Parameter initialization (Global-init) failed	
29211	Switching the motor holding brake not possible, Controller is on	
29214	Homing not possible	
29215	Continuous movement not possible	

29217	Invalid input parameter	
29218	Unknown PLCopen axis state	
29219	Invalid PLCopen parameter value	
29221	No cam name	
29222	Error at cam download	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
29225	The target position is outside the axis period	
29226	Error on drive. Use MC_ReadAxisError for details	
29227	No further master position can be sent on the network from this drive	
29228	No further master position can be read from the network on this drive	
29229	Synchronized movement not possible	
29230	Internal error: Error at parameter list transfer	
29231	The master velocity is invalid, 0 or negative	
29232	Internal error: Invalid SPT resource type	
29233	SPT resource of required type not available	
29234	Internal error: Number of requested SPT resources not available	
29235	The functionality is not available for the current axis type	
29237	Error in TriggerInput parameters	
29238	The FB cannot be used in the current state	
29239	This functionality is not available on CAN-Bus	
29241	Wrong data type for specific ParID	
29242	Cyclic read data full	
29246	TouchProbe window invalid	
29247	Master sync position cannot be reached	
29250	CamTableID is invalid	
29251	Error at parameter table download	
29252	Error at parameter list initialization	
29253	Error at parameter sequence download	
29254	Error at parameter sequence initialization	
29255	Initialization not possible, axis coupling is active	
29256	Multiple commands not possible at the same time	
29257	The specified data address is invalid	
29260	No data object name specified	
29261	Invalid data object index	
29263	Slave channel already in use	
29264	Cyclic write data full	
29265	Communication to drive failed	
29266	The MasterParID was changed since last FB call	
29267	Invalid number of cam profile polynomials	

29268	FB aborted by another one	
29269	Error during saving of the NC-INIT parameter module	
29270	Error during loading of the NC-INIT parameter module	
29271	Selected MC_TouchProbe function block is not active	
29272	Error at initialization of Automat	
29273	The given "Subject" is invalid	
29274	Error at initialization of parameters. Use MC-ReadAxisError for details	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
29275	At least one input value was changed while "Enable = TRUE"	
29276	A phase shift is already in progress	
29277	An offset shift is already in progress	
29278	Selected mode not supported without entry of a period	
29279	Output value cannot be calculated	
29280	No valid master defined	
29281	This functionality is not available for ACOPOSmulti drives	
29283	Master or slave position of first cam profile point not equal to 0	
29284	Too few curve points	
29285	Invalid cam profile section type	
29286	Invalid mode for the last cam profile point	
29287	Invalid master or slave position for last cam profile point	
29288	Master positions are not strictly monotonic increasing	
29289	Invalid boundary parameters	
29290	Too many cam profile polynomials	
29291	Turning point outside of cam profile section	
29292	Identical slave positions not permitted	
29293	Specified data length is 0 or too low	
29294	Not able to determine error text. See error text string for details	
29295	An error has occurred. See "ErrorRecord" output for details	
29296	Not enough space in Broadcast channel	
29297	Problem with variable in permanent memory	
29298	Erro in network configuration	
29299	Error occurred during the setup operation	
29300	Invalid number of polynomials in the cam	
29301	No cam value can be calculated	
29302	One instance of the function block is already active	
29303	Specified IntervalTime too small	
29304	This functionality is not available for ACOPOS	

29305	ParID cannot be read in this mode	
29306	Invalid interpolator mode	
29307	Master period is zero	
29308	Internal calculation error	
29309	General internal error	
29489	Internal values of the axis structure are invalid	
29490	Error at internal initialization (global init)	
29491	Error at internal initialization (software limits)	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
29492	Error at internal initialization (homing of virtual axis)	
29498	ACP10_MCLibrary: Initializaiton aborted	
29499	ACP10_MCLibrary: Error with details in "ASCII Data"	
31201	Di/Do Interface: Drive not ready	
31220	Encoder error: Encoder not configured	
31221	Encoder error: Cable disturbance or signal disturbance	
31224	Encoder Interface: HW Module not OK	
31240	Homing procedure mode not allowed with current HW Type	
31247	Drive Interface: DrvOK not set from HW Module	
31248	Trigger Interface: HW Module not OK	
31249	Drive Interface: HW Module not OK	
31250	Di/Do Interface: HW Module not OK	
31260	Current axis configuration only possible in simulation mode	
32001	Error calling CAN_xopen()	
32002	Error defining Write COB for Broadcast Command	
32003	Error defining Write COB for Parameter Read Request	
32004	Error defining Write COB for Parameter Write Request	
32005	Error defining Read COB for Parameter Read Response	
32006	Error defining Read COB for Parameter Write Response	
32007	Error defining Read COB for Monitor Data from the drive	
32008	Error sending Read Request (network error ?)	
32009	Error sending Write Request (network error ?)	
32010	Drive not responding to Read Request (is the drive in the network ?)	
32011	Drive not responding to Write Request (is the drive in the network ?)	
32012	Error reading module description of system module	
32013	No operating system present on the drive	
32014	Operating system version on the drive not compatible with NC manager version	
32015	Error creating message queue	

32016	Error sending an idle time command to the NC Manager Task	
32017	Wrong boot state after start of operating system	
32018	Invalid Parameter ID in system module	
32019	Download of NC system module not allowed (the module is on the PLC)	
32020	System module data could not be read from the drive during NC manager INIT	
32021	System module data could not be read from the drive after download	
32022	Error aborting data block access before download	
32023	Error reading boot state before download	
32025	Wrong boot state after SW Reset before download	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32026	Error during INIT of data block write access for download	
32027	Error sending data segment for download	
32029	Response error after sending data segment for download	
32030	Error at command for system module burn after download	
32031	Error reading status for system module burn after download	
32032	Error while burning system module after download	
32033	Timeout while burning system module after download	
32034	Error at SW Reset before download	
32035	Error at SW Reset after download	
32036	Different system module data after download	
32037	Error message(s) lost because of FIFO overflow (acknowledge errors)	
32038	Error searching for INIT parameter module	
32039	Error reading INIT parameter module	
32040	Version of INIT parameter module is not compatible to NC manager	
32041	The module acp10cfg does not exist	
32042	The module acp10cfg is not an NC data module	
32043	The NC module type of the module acp10cfg is invalid	
32044	The NC module type of the module acp10cfg cannot be read	
32045	The data address in module acp10cfg cannot be read	
32046	The data section of module acp10cfg is empty	
32047	A CAN node number in module acp10cfg is invalid	
32048	A CAN node number in module acp10cfg is used repeatedly	
32049	This NC action is not allowed during Trace is active	
32050	A Trace Data Upload is already active	
32051	Invalid Trace Status for Trace Data Upload	

32053	Error defining Write COB for Parameter Read Request 2	
32054	Error defining Write COB for Parameter Write Request 2	
32055	Error defining Read COB for Parameter Read Response 2	
32056	Error defining Read COB for Parameter Write Response 2	
32057	Error accessing HS task class table	
32058	Error accessing task class table	
32059	Parameter tk_no invalid for access to task class table	
32060	Timeout for cyclic data from drive - Indications invalid (network error ?)	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32061	Timeout sending a Read Request telegram (network error ?)	
32062	Timeout sending a Write Request telegram (network error ?)	
32063	Data address zero (set/read parameter via service interface)	
32064	Convert text into binary data is not possible for this parameter data type	
32065	Convert binary data into text is not possible for this parameter data type	
32066	Parameter ID zero (set/read parameter via service interface)	
32067	Parameter ID invalid (set/read parameter with option ncDATA_TEXT)	
32069	The data address of the ACOPOS parameters in module acp10cfg cannot be read	
32070	Drive for ACOPOS parameters in module acp10cfg not found	
32071	The ACOPOS parameters are invalid (an update of AutomationStudio is necessary)	
32072	Wrong boot state after SW Reset	
32073	Download of NC system module: Error reading NC hardware version of BsLoader	
32074	Incompatible NC hardware version: Download of BsLoader not possible	
32075	Incompatible NC hardware version: Download of operating system not possible	
32076	FIFO for messages with high priority to NC Manager Task is full	
32077	A POWERLINK node number in module acp10cfg is invalid	
32078	A POWERLINK node number in module acp10cfg is used repeatedly	
32079	With this version one CAN interface must be in module acp10cfg	
32080	With this version one POWERLINK interface must be in module acp10cfg	
32081	With this version only one POWERLINK interface is allowed in module acp10cfg	

32082	Module acp10cfg contains a CAN interface without any drive node	
32083	Module acp10cfg contains a POWERLINK interface without any drive node	
32084	The NC configuration does not contain any ACOPOS module	
32085	Module acp10cfg invalid (AutomationStudio V2.2 or higher necessary)	
32086	With this version no CAN interface is allowed in module acp10cfg	
32087	With this version no POWERLINK interface is allowed in module acp10cfg	
32088	The INIT parameter module specified in the NC Mapping Table does not exist	
32089	NC-HW-ID of INIT parameter module is not compatible to NC manager	
32090	NC object type of INIT parameter module is not equal to NC object	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32091	Invalid block data in INIT parameter module (data range exceeded)	
32092	Error sending a command to the NC Idle Task	
32093	NcManCtrl is defined repeatedly with different values	
32094	NetworkInit is defined repeatedly for ncMANAGER with different values	
32095	Value of drive group in CAN-CFG-Module higher than maximum value	
32096	Size of data buffer for trace data upload in module acp10cfg too small	
32097	All counts of used network interfaces in module acp10cfg are zero	
32098	Version of the module acp10cfg is not compatible with NC manager	
32099	Length of data section of module acp10cfg is too small	
32100	Memory for NC error text management cannot be allocated	
32101	Error accessing NC error text module in BR module table	
32102	Version ID of error text module not equal to that of NC manager	
32103	Data section of error text module cannot be read	
32104	Data section of error text module is empty	
32105	Length of data section of error text module is too small	
32106	Error list of error text module not equal with that of NC manager	
32107	Parameter list of error text module not equal with that of NC manager	
32108	The last error number of error text module is not equal to 65535	
32109	The last parameter ID of error text module is not equal to 65535	
32110	Length of data section of CAN-CFG-Module cannot be read	

32111	Length of data section of CAN-CFG-Module is too small	
32112	The data address in the CAN-CFG-Module cannot be read	
32113	The enable code in the CAN-CFG-Module is invalid	
32114	Values not equal to zero in reserved area of CAN-CFG-Module	
32115	The basis CAN ID for WR/RD channel1 in the CAN-CFG-Module is invalid	
32116	The basis CAN ID for WR/RD channel2 in the CAN-CFG-Module is invalid	
32117	The basis CAN ID for WR/RD channel3 in the CAN-CFG-Module is invalid	
32118	The basis CAN ID for monitor data in the CAN-CFG-Module is invalid	
32119	Invalid basis CAN ID for cyclic data to the drive in CAN-CFG-Module	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32120	Invalid basis CAN ID for cyclic data from the drive in CAN-CFG-Module	
32121	The CAN ID for the SYNC telegram in the CAN-CFG-Module is invalid	
32122	The CAN ID for the broadcast command in the CAN-CFG-Module is invalid	
32123	Error defining Read COB for WR2 Request (external set position mode)	
32124	Error defining Read COB for WR2 Response (external set position mode)	
32125	Error defining Read COB for RD2 Request (external set position mode)	
32126	Error defining Read COB for RD2 Response (external set position mode)	
32127	Error deleting Write COB for Broadcast Command (external set position mode)	
32128	Error defining Read COB for Broadcast Command (external set position mode)	
32129	Error defining Read COB for cyclic user data from drive (ext. set pos. mode)	
32130	This external set position mode is only allowed with one CAN interface	
32131	The specified NC data module does not exist	
32132	The specified module is not an NC data module	
32133	The NC module type of the specified NC data module is invalid	
32134	The NC module type of the specified NC data module cannot be read	
32135	The data address of the specified NC data module cannot be read	
32136	The Data section of the specified NC data module is empty	
32137	Data address of structure for a data block operation is zero	
32138	Data address zero (user data for data block operation)	
32139	Data length zero (user data for data block operation)	

32140	Data block operation: Data module name or data address must be zero	
32141	Invalid data format in a parameter sequence	
32142	ID or type of a parameter invalid in parameter sequence with text format	
32143	Data of a parameter in a parameter sequence longer than 6 bytes	
32144	Error for an ACOPOS Parameter Table specified in the NC Mapping Table	
32145	The ACOPOS Parameter Table does not exist	
32146	The ACOPOS Parameter Table is not an NC data module	
32147	The NC module type of the ACOPOS Parameter Table is invalid	
32148	The NC module type of the ACOPOS Parameter Table cannot be read	
32149	The data address in the ACOPOS Parameter Table cannot be read	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32150	The data section of the ACOPOS Parameter Table is empty	
32151	Error initializing memory buffer for XML parser	
32152	No XML elements present in an ACOPOS Parameter Table	
32153	The first XML element is invalid in the ACOPOS Parameter Table	
32154	The ACOPOS Parameter Table does not contain any ACOPOS parameters	
32155	Nesting depth for ACOPOS parameter groups exceeded	
32156	ID or type of an ACOPOS parameter invalid for text conversion	
32157	Length of parameter data too large for ACOPOS parameter in XML data	
32158	ACOPOS parameter: An attribute is not defined (ID)	
32159	ACOPOS parameter: An attribute is not defined (Value)	
32160	Basis movements with mode ncTRG_STOP are not allowed for ncV_AXIS	
32161	ncNC_SYS_RESTART,ncACKNOWLEDGE is not allowed (network.init=ncFALSE)	
32162	Internal task class number wrong (from now on operation is blocked !!!)	
32163	A system module download to all drives is not possible with SwNodeSelect	
32164	The text defined with NetworkInit (global) is invalid	
32165	A CAN node number is equal to NodeNr_SwNodeSelect	
32166	Network initialization during active network initialization not allowed	
32167	The text defined with NetworkInit is invalid	
32168	NodeNr_SwNodeSelect is defined repeatedly with different values	
32169	The node number defined with NodeNr_SwNodeSelect is invalid	

32170	A data module name has to be entered for this data block operation	
32171	Index zero is not allowed (user data for data block operation)	
32172	The specified data module name is not valid for a BR module	
32173	Memory for data module creation cannot be allocated	
32174	Error with installation of data module into BR module table	
32175	Error with installation of data module into BR module table	
32176	Text for parameter data too large for parameter sequence with text format	
32177	Text for parameter data too large for parameter list with text format	
32178	This NC object is not enabled for this ACOPOS (channel number too high)	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32179	ID or type of a parameter invalid in parameter list with text format	
32180	Data address of structure for a parameter list operation is zero	
32181	Data address zero (user data for parameter list operation)	
32182	Data length zero (user data for parameter list operation)	
32183	Data length invalid (user data for parameter list operation)	
32184	Invalid data format in a parameter list	
32185	Data of a parameter in a parameter list longer than 6 bytes	
32186	NetBasisInitNr is defined repeatedly for ncMANAGER with different values	
32187	Error for synchronization of network initialization (details in Logger)	
32188	This NC object is defined in hardware configuration and NC Mapping Table	
32189	Timeout for cyclic data from drive - Indications invalid (network error ?)	
32190	Error defining Write COB for selection of node number via software	
32191	This parameter ID is reserved for the PLCopen MC library	
32192	The specified data module is not an INIT Parameter module	
32193	For this NC object type no INIT parameter module is present	
32194	This function is not implemented for this NC object type	
32195	Error downloading BsLoader to ACOPOS	
32196	Error downloading operating system to ACOPOS	
32197	Error downloading BsLoader to ACOPOS (additional info in Logger)	

32198	Error downloading operating system to ACOPOS (additional info in Logger)	
32200	Error calling plAcycWrite() (read parameter)	
32201	Error calling plAcycWrite() (write parameter)	
32202	Error calling plAcycRead() (read parameter)	
32203	Error calling plAcycRead() (write parameter)	
32204	Timeout while reading par. via acyclic channel (is the drive in the network ?)	
32205	Timeout while writing par. via acyclic channel (is the drive in the network ?)	
32206	Cyclic channel: Read Request in spite of Wait for Response	
32207	Cyclic channel: Write Request in spite of Wait for Response	
32208	Error using plAction(DEVICE_TO_BUS_NR) (additional info in Logger)	
32209	Error using plAction(GET_IDENT) (additional info in Logger)	
32210	Wrong interface ident when calling plState() (additional info in Logger)	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32211	Interface not available when calling plState() (additional info in Logger)	
32212	Fatal interface error when calling plState() (additional info in Logger)	
32213	Timeout for POWERLINK interface (additional info in Logger)	
32214	Error calling plAcycOpen() (additional info in Logger)	
32215	Error calling plCECreate() (additional info in Logger)	
32216	Error using plAction(GET_IF_PAR) (additional info in Logger)	
32217	Broadcast channel: Error calling plAcycWrite() (read parameter)	
32218	Broadcast channel: Error calling plAcycWrite() (write parameter)	
32219	Error using plAction(GET_IF_MUXPRESCALE) (additional info in Logger)	
32220	Error using plAction(GET_IF_CYCLE_TIME) (additional info in Logger)	
32221	Error using plAction(GET_IF_PRESCALE) (additional info in Logger)	
32222	Error using plAction(GET_STATIONFLAG) (additional info in Logger)	
32223	Error calling plGetNodeInfo() (additional info in Logger)	
32224	Error calling plAction(GET_PROTOCOL_VERSION) (additional info in Logger)	
32225	This ACOPOS POWERLINK node does not exist in the AR Configuration	
32226	A SDC node number in module acp10cfg is invalid	
32227	A SDC node number in module acp10cfg is used repeatedly	

32228	There is no network interface (POWERLINK or SDC) contained in acp10cfg	
32229	The object acp10sdc does not exist (necessary for SDC NC objects)	
32230	Error at SDC initialization (see Logger)	
32231	Error at SDC configuration (see Logger)	
32232	The object acp10sim does not exist (necessary for ACOPOS simulation)	
32233	The object acp10bsl contains no NC system module	
32234	The object acp10sys contains no NC system module	
32235	Data address zero for parameter in parameter sequence	
32236	Data address zero for parameter in parameter list	
32237	Channel index for ACOPOS parameters in module acp10cfg is invalid	
32238	This function is not implemented at this time	
32239	Basic memory for INIT parameter modules could not be allocated	
32240	NC object data invalid (PV with INIT value in variable declaration ?)	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32241	Data block operation: Data module name must be zero	
32242	Conditions for cam profile transfer via Systick-Task are not fulfilled	
32243	Error using plAction(GET_PDO_INFO) (additional info in Logger)	
32244	No PDO defined in the cyclic frame for this channel: NC object is disabled	
32245	The Memory for operating the broadcast channel cannot be allocated	
32246	'Direction' of the POWERLINK broadcast channel is not 'Output'	
32251	NcNetCyc: Response timeout	
32252	NcNetCyc: Unexpected Response (invalid counter value)	
32260	Trace: No valid test data point is defined	
32261	Trace trigger: The NC object is invalid	
32262	Trace test data point: The NC object is invalid	
32263	Trace trigger: The NC object is not enabled for this ACOPOS	
32264	Trace test data point: The NC object is not enabled for this ACOPOS	
32265	Trace trigger: Trace for this NC object already active at trace start	
32266	Trace test data point: Trace for this NC object already active at trace start	
32267	Trace trigger: The trace interface for this NC object is occupied	
32268	Trace test data point: The trace interface for this NC object is occupied	
32269	Trace: Too many test data points for one NC object	

32270	ACOPOS coupling: The channel number of send data is invalid	
32271	ACOPOS coupling: The channel number of receive data is invalid	
32272	ACOPOS coupling: Send data with this channel number are not configured	
32273	ACOPOS coupling: The NC object of send data is invalid	
32274	Network coupling: The broadcast channel was not configured	
32275	Network coupling: The broadcast channel was not successfully initialized	
32276	Network coupling: Cycle times not equal (network interfaces, NC task class)	
32277	Network coupling: In the broadcast channel no more data record has place	
32278	ACOPOS coupling: Zero as Parameter ID of send data is not allowed	
32279	Service interface: Data address zero	
32280	Timeout for enable of acyclic network communication	
32294	NC Manager error: Pointer of function nctune_sequ() is zero	
32295	NC Manager error: Data structure address zero when calling nctune_sequ()	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32296	NC Manager error: Data structure size invalid when calling nctune_sequ()	
32297	NC Manager error: NC object address zero when calling nctune_sequ()	
32298	NC Manager error: Unknown error number when calling nctune_sequ()	
32300	For controller setup the library acp10tun must be transferred to the target	
32301	Setup start not allowed: 'count.error' is greater than zero	
32302	Start setup not possible: A setup operation is already active	
32303	Setup start not possible: The axis is not global initialized	
32304	Setup start not possible: The controller is not switched off	
32305	Setup start not possible: The controller is not ready	
32306	Setup start not possible: The trace is already active	
32307	Setup for controller: Percentage for rated current higher than maximum value	
32308	Setup for controller: Percentage for rated current lower than minimum value	
32309	Setup for controller: Percentage for speed higher than maximum value	
32310	Setup for controller: Percentage for speed lower than minimum value	
32311	Setup for controller: Move distance higher than maximum value	

32312	Setup for controller: Move distance lower than minimum value	
32313	Setup start not possible with this mode: The axis is not referenced	
32314	Setup for controller: Mode invalid	
32315	Setup for controller: Mode invalid (Library)	
32316	Setup for controller: The calculated inertia is zero	
32317	Setup for controller: SCTRL_KV is zero	
32318	Setup for controller: POS_CTRL_KV is zero	
32319	Setup for controller: Calculation of result parameters was not possible	
32320	Setup for controller: Speed limit values too high for MOTOR_SPEED_RATED	
32321	Setup for controller: Coefficient a0 is zero ('+ncISQ_F1_NOTCH' not possible)	
32400	Memory cannot be allocated	
32401	No file name specified	
32402	The specified file cannot be opened	
32403	The specified file cannot be created	
32404	Error at writing into a file	
32405	Error at reading from a file	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32406	The specified file cannot be created	
32407	Error at writing into a file	
32408	Error at closing a file	
32409	The specified file cannot be deleted	
32410	The data object type is invalid	
32411	The address of the data object structure is zero	
32412	The name in the data object structure is zero	
32413	The data address in the data object structure is zero	
32414	The data length in the data object structure is zero	
32415	No File Device specified	
32416	The FileIO functions do not exist (is library FileIO existing ?)	
32417	Error at writing into an NC data module	
32418	No XML elements present in MTC data	
32419	MTC data: XML elements invalid or in wrong order	
32420	The MTC data contain an invalid configuration	
32421	The MTC data contain noconfiguration for Multi Axes Trace	
32422	The MTC data contain more than one configuration for Multi Axes Trace	
32423	MTC-Daten: Trigger.NcObject ist invalid	
32424	MTC-Daten: Trigger.Condition ist invalid	
32425	MTC-Daten: Channel.NcObject ist invalid	

32426	Trace trigger: The event is invalid	
32427	Trace trigger: Parameter ID zero not allowed if event unequal to OFF	
32494	AcoposSimulation=Off for one channel although activated for the other channel	
32495	AcoposSimulation: Different values defined for real and virtual axis	
32496	Error creating cyclic task for PLCopen MC (details in Logger)	
32497	Task class for handling of cyclic data with PLCopen in acp10cfg is invalid	
32498	PLCopen_CyclicData_TaskClass is lower than zero or higher than maximum value	
32499	PLCopen_CyclicData_TaskClass: Values for real and virtual axis are not equal	
32500	The Message FIFO already exists	
32501	Error creating Message FIFO	
32502	The Critical Section for Command Semaphore already exists	
32503	Error creating Critical Section for Command Semaphore	
32504	The NC Manager Idle Task already exists	
32505	Error creating NC Manager Idle Task	
32506	Error reading Taskclass Cycle Time	
32507	Error reading Taskclass Tolerance	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
32508	Error sending an idle time command to the NC Manager Task	
32509	The Critical Section for Network Command Trace already exists	
32510	Error creating Critical Section for Network Command Trace	
32511	The Critical Section for messages with high priority already exists	
32512	Error creating Critical Section for messages with high priority	
32513	The Critical Section for global variables already exists	
32514	Error creating Critical Section for global variables	
33000	Master cycle time is not a multiple of communication cycle time	
33002	Floating-Point exception	
33003	Address error exception	
33004	Bus error exception	
33005	Exception	
33006	Access violation Exception	
33007	Violation address	
36001	Parameter limited to valid range	
36002	Total time for the position loop controller limited to prediction time	

37101	Calculated compensation distance on slave axis limited to maximum	
37102	Calculated compensation distance on slave axis limited to minimum	
37103	Slave trigger outside of window	
37104	Slave trigger missing	
37105	Master compensation trigger outside of window	
37107	Displacement actual/set position too high during 'controller switch on'	
37108	Calculated compensation distance of master axis limited to minimum	
37109	Master trigger outside of window	
37110	Master position at start higher than first trigger position	
37111	Cam profile data: Difference between polynomial value $y(x_n)$ and slave period	
37112	Polynomial within cam profile data exceeds limit value	
37113	Cam profile compensation gears: Limit values exceeded	
38000	Current controller: Motor speed is too high during switching on	
38001	Torque limiter: Limit value too high	
38003	Motor holding brake: Test torque was limited	
<b>Error ID</b>	<b>Error Text</b>	<b>Notes</b>
38004	Motor holding brake: Test torque less than load torque	
38005	Motor test: Speed is too high during switch on	
38006	Current controller: Permissible current offset values exceeded	
38007	System administration: BsLoader activ: Jumper plugged	
38008	Bleeder: No current flow	
39001	Encoder: Position correction active	
39002	Resolver: Speed limit for 14 bit resolution exceeded	
39003	EnDat encoder: Alarm bit is set	
39004	EnDat encoder: Alarm bit - Lighting failure	
39005	EnDat encoder: Alarm bit - Signal amplitude too small	
39006	EnDat encoder: Alarm bit - Position value contains an error	
39007	EnDat encoder: Alarm bit - Overvoltage	
39008	EnDat encoder: Alarm bit - Undervoltage	
39009	EnDat encoder: Alarm bit - Overcurrent	
39010	EnDat encoder: Alarm bit - Battery change required	
39011	EnDat encoder: Warning bit - Frequency too high	
39012	EnDat encoder: Warning bit - Temperature too high	
39013	EnDat encoder: Warning bit - Lighting reserve reached	
39014	EnDat encoder: Warning bit - Battery charge too low	

39016	Incremental encoder emulation: Frequency too high	
39017	Encoder: CRC error while reading position	
39018	Reference pulse monitoring: Faulty position, resolution, or reference pulse	
39019	Serial encoder interface: Stop bit error	
39020	Serial encoder interface: Receive data overrun	
39021	Serial encoder interface: Send data error	
39022	EnDat encoder: Warning bit is set	
39023	EnDat encoder: CRC error while reading EnDat2.2 additional information	
39024	EnDat encoder: Alarm bit - Power failure	
39025	EnDat encoder: Alarm bit - multiturn overflow	
39026	EnDat encoder: Type 3 error while reading EnDat2.2 additional information	
39027	Encoder Emulation: 5V power supply fail	
39028	Encoder: Multiturn failure	
39301	Digital IO: 24V power supply fail	
39302	Digital IO 1-4: Diagnose bit active (current, 24V supply)	
39303	Digital IO 5-8: Diagnose bit active (current, 24V supply)	
39305	Digital IO 10: Diagnose bit active (current, temperature)	
39306	Digital IO 9: Diagnose bit active (current, temperature)	
39307	Digital IO: outputs deactivated by output mask after network error	
41001	Heatsink temperature sensor: Overtemperature	
41011	Temperature sensor (Motor Choke External): Overtemperature	
41031	Junction temperature model: Overtemperature	
41041	Bleeder temperature model: Overtemperature	
41051	ACOPOS peak current: Overload	
41061	ACOPOS continuous current: Overload	
41070	Motor temperature model: Overtemperature	
41075	ACOPOS continuous power: Overload	
41078	ACOPOS temperature sensor: Overtemperature	
64001	ncalloc in slower task class than defined for NC Manager Task	
64002	Delay before SW Reset (network with ascending node numbers ?)	
64003	Delay before NC System Start (network with ascending node numbers ?)	
64004	The following boot error could be entered here with a delay	
64005	Timeout for parameter enable after start of operating system	
64006	Drive did not become synchronous with network master	
64007	Timeout for enable of acyclic network communication	

64008	Timeout for enable of cyclic network communication	
64009	Acp10cfg contains a POWERLINK interface, for which no axis is defined	
64010	NC software test version (validity period in Logger)	
64500	Positive speed limiter active	
64501	Negative speed limiter active	
64502	Positive direction acceleration torque limiter active	
64503	Negative direction acceleration torque limiter active	
64504	Positive direction deceleration torque limiter active	
64505	Negative direction deceleration torque limiter active	
64506	Recovery power limiter active (deceleration too high)	
65535	Response error	

## 12.2 Communication Error Reference Guide

<i><b>Error Number</b></i>	<i><b>Description</b></i>	
60	No message received	
8071	No buffer available in the send buffer manager	
8072	Invalid buffer rejected	
8073	Invalid I/O control code	
8078	Send buffer waiting queue full	
8079	Frame buffer with defective character	
8210	Buffer is corrupted	
8251	Device not open	
8252	Device not found	
8253	Syntax error in mode parameter string, or mode parameter not supported by the specified device	
8254	Too many devices opened at once	
8255	/PA, /DB, /SB must always be specified together	
8256	Command not supported by driver (hardware)	
8257	Required resource not available (memory, OS object, etc.)	

## 12.3 Heat Sealer Error Reference Guide

<i>Error ID</i>	<i>Description</i>	<i>Notes</i>
65534	The function block is disabled. enable = FALSE.	
12001	The integral element of the manipulated variable has reached its maximum limit.	
12002	The integral element of the manipulated variable has reached its minimum limit.	
12003	The Y_fbk input exceeds the range defined with Y_max and Y_min and has therefore been limited to this range.	
12004	The cycle time in which the LCRPIDTune() function block is called is too large. Therefore, calculated control parameters could be inaccurate. The cycle time should be less than 1/10 of the smallest time constant in the system (eg dead-time, delay times of PT1 or PT2).	
12005	T1 = 0. The intended PT1 functionality is reduced to purely P functionality. An intended PT2 functionality is reduced to simply PT1 functionality, if T2 > 0. If T2 should also equal 0, then the PT2 element degenerates to a true P-element.	
12007	A time parameter is too small based on the actual cycle time dt of the function block call. Since this is first determined during runtime, the affected time parameter is limited and this warning is output. More detailed information can be found in the description of the status output for the respective function block.	
12008	The A input exceeds the value range defined with Y_max and Y_min while being limited to this range.	
12009	Input Y_man exceeds the value range defined with Y_max and Y_min while being limited to this range.	
12010	T1 value entered is not larger than the corresponding task cycle time.	
12012	T2 = 0. The intended PT2 functionality is reduced to simply PT1 functionality, if T1 > 0. If T1 should also equal 0, then the PT2 element degenerates to a true P-element.	
12013	T2 value entered is not larger than the corresponding task cycle time.	
12015	Tt = 0. The intended dead-time behavior is reduced to an identical image.	
12016	Tt was entered as smaller than the corresponding task cycle time, and was automatically increased to this value.	
12017	Tt was not entered as a whole number multiple of the cycle time. Tt was set to the next lower whole-number multiple of the cycle time.	
12018	An error occurred while allocating memory for the internal ring buffer. This results in a limited dead-time functionality with Tt = 0 (ie only an identical image).	

<b>Error ID</b>	<b>Description</b>	<b>Notes</b>
12019	A change of the tuning mode is not possible and will be ignored.	
12020	Amplitude during oscillation attempt tuning is too low. Therefore, an additional switching hysteresis is used (see Data types: lcrpid_tune_osc_options_typ). This warning no longer occurs after V2.50.	
12021	No turning point detected during heating.	
12022	Values not set or set incorrectly in the lcrtemp_tune_set_typ structure were replaced by default values.	
12023	Not set or incorrectly set filter_base_T rounded up to a multiple of the task cycle time or set to the default value.	
12024	Values not set or set incorrectly in the lcrtemp_pid_set_typ structure were replaced by default values.	
12025	Invalid value for Kfbk	
12026	Invalid value for Kp.	
12027	Invalid value for Kw	
12028	Invalid value for Tf	
12029	Invalid value for Tn.	
12030	Invalid value for Tv	
12031	Invalid value for T	
12032	t_min_pulse is greater than t_period.	
12033	t:period is less than dt (= cycle time of the function block call).	
12034	Invalid value for dy_up or dy_down.	
12035	Invalid value for t_change_up or t_change_down.	
12036	Invalid value for t_impulse (>= t_change_up or t_change_down)	
12037	Y_max is less than or equal to Y_min or Y_max > 100 or Y_min < -100	
12038	Tuning aborted because request was changed before the state output was set to LCRPIDAUTOTUNE_STATE_FINISHED (50).	
12039	Error calculating controller parameters (see addInfo output).	
12040	Tuning aborted because the maximum time addPar.t_autotune_max has been exceeded.	
12041	Tuning aborted because W was changed before the state output was set to LCRPIDAUTOTUNE_STATE_FINISHED (50).	
12042	Invalid value for deadband_mode. Permitted values: 0, 1, 2.	
12043	WX_max is less than or equal to WX_min.	
12044	The dead band is greater than the range from WX_min..WX_max (deadband > (WX_max - WX_min)).	
12045	NoOfPoints is not between 2 and 1000, or ptr_table is zero, or violation of the continuity rule found.	
12046	T1 is negative. y remains 0 until T1 returns to a positive value.	

<b>Error ID</b>	<b>Description</b>	<b>Notes</b>
12048	Tt is negative. y remains 0 until Tt returns to a positive value. A new value is not applied until a positive edge on the enable input.	
31550	Invalid value for Tn.	
31551	Invalid value for base	
31552	ident invalid	
31553	Specified mode is invalid.	
31554	Parameter data indicated as invalid by LCRPIDpara().	
31556	Invalid value for d_mode.	
31557	Invalid value for dy_max.	
31558	Invalid value for fbk_mode.	
31559	Tuning aborted because an invalid request was detected.	
31560	Tuning aborted because the controlled variable X is outside the permissible range ( $X < X_{min}$ or $X > X_{max}$ ).	
31561	The difference between W and X is too minor to start tuning.	
31562	Error determining the slope of the controller value.	
31563	Null pointer specified.	
31564	Y_min or Y_max was changed during tuning.	
31565	Set temperature (Temp_set) is lower than the ambient temperature.	
31566	Parameter structure (Icrtemp_set_typ) for LCRTempTune() is not connected.	
31567	In the Icrtemp_pid_opt_typ structure, an invalid value has been entered for Kp_h or Kp_c.	
31568	Parameter structure (Icrtemp_set_typ) for LCRTempPID() is not connected.	
31569	Specified mode is invalid, see LCRTempPID()	
31570	Task class invalid or not recognized.	
31571	tKnick >= tAlt was selected for aging!	
31572	Invalid internal step causes a calculation error of tAlt. Restart FBK.	
33100	The increase of the current temperature was above the specified increase ( $\Delta T_{sync\_cool}$ ) during the time frame cnt_wait_cool and the temperature has already risen to the upper temperature limit ( $T_{set} + \Delta T_{sync\_stop}$ ). The rdyToCool output is set to TRUE to make synchronization possible. If the okToCool input is not set to TRUE after 5 task cycles (FBK call), the tuning is executed anyways with the warning output LCR_WARN_LCRTEMP Tune_ASYNC.	
33101	The tuning procedure will be executed immediately because the current temperature has exceeded the specified upper limit ( $T_{set} + \Delta T_{sync\_stop}$ ), even though the other zones are not yet ready for the cooling tuning. As a result, tuning of the pure cooling zones is not synchronous to the other zones.	

<i><b>Error ID</b></i>	<i><b>Description</b></i>	<i><b>Notes</b></i>
33102	Warning: The filter time constant Tf entered is not larger than the corresponding task cycle time.	
33150	The combination of enable_cooling and disable_heating would make a tuning procedure unnecessary.	
33151	The combination of enable_cooling and disable_heating would make regulation unnecessary.	
33152	The time constant t_ignore_dT_heat is greater than the heating time.	
33153	Error: The damping constant D is negative. y remains 0 until D returns to a positive value.	
33154	Error: The filter time constant Tf is less than or equal to 0. y remains 0 until Tf returns to a strictly positive value.	

6064 McKee Road

Fitchburg, WI 53719 USA

Toll Free: 1 (877) 270-1250

Phone: 1 (608) 270-1230

E-mail: [service@genepro-inc.com](mailto:service@genepro-inc.com)