pfmmedical

Quality and Experience

Nit-Occlud® PDA

en	Instructions for Use	4
de	Gebrauchsanweisung	12
fr	Mode d'emploi	20
es	Instrucciones de uso	28
it	Istruzioni per l'uso	36
pt	Instruções de utilização	44

 0481

www.pfmmedical.com

Explanation of the symbols on label and packaging // Erklärung der Symbole auf Etikett und Verpackung // Explication des symboles sur l'étiquette et l'emballage // Explicación de los símbolos en la etiqueta y el envoltorio // Illustrazione dei simboli riportati su etichetta e confezione // Explicação dos símbolos nas etiquetas e embalagens.

REF Reference number // Bestellnummer // N° de commande // Número de referencia para pedidos // Numero d'ordine // Número de referência

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
Nit-Occlud® PDA

Product description

Nit-Occlud® PDA is a system for transcatheter occlusion of Patent Ductus Arteriosus (PDA) with spiral coils. The system consists of the following parts:

- Nit-Occlud® PDA

The spiral coil is mounted in a straightened fashion on a flexible delivery system including a disposable handle to which it is connected by means of a patented detachment mechanism.


Nit-Occlud® PDA coils are available in 3 different types: Flexible, Medium, Stiff.

The flexible and medium types are pre-loaded into the transportation sheath. For insertion the transportation sheath must be connected to the implantation catheter. The stiff type is pre-loaded in a 5F implantation catheter. For insertion, an additional long introducer sheath (5 F or larger) is required.

The Nit-Occlud® PDA coil has a cone in cone configuration which results from the fact that the proximal windings of the coil are wound in the reverse direction (see Fig.).

Fig.: Nit-Occlud® PDA spiral coils

(D=Distal diameter, P=Proximal diameter, Lc=Length configured).

- **Implantation catheter**

The implantation catheter is equipped with a marker ring at its distal tip for better orientation during fluoroscopy.

General Information and Warnings

These instructions for use and the information on the packaging should be read carefully before each use. The PDA coil system should be used only by physicians trained in interventional occlusion techniques.

WARNINGS

- Do not use the product if the packaging has been opened, or is damaged, if you are not sure that it is sterile, or if the expiry date has passed.
- Each product is packed separately, and is delivered in an EO-sterilised and pyrogen-free condition. It is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization of single-use devices may result in degraded performance or a loss of functionality. Reuse of single-use devices may result in exposure to pathogens such as viruses, bacteria, fungi, or prions.
- The product must be stored in dry conditions. Do not expose the packaged products to direct sunlight.
- Retrieval devices should be available during implant procedures for interventional retrieval of the coil if required.
- Care must be taken not to damage the coil or to dislodge it from the delivery system while unpacking or inserting it into the implantation catheter.
- Since the delivery system has ferromagnetic properties, implantation must not be carried out in an MR environment.
- The coil should not be removed from the delivery system. It should not be used with another delivery system since this may alter characteristics of configuration and detachability.
- A detached coil should not be remounted on the core wire of the delivery system.
- The configured coil should not be pulled through heart valves or ventricular chambers.

- The implantation catheter is not suitable for application of contrast medium. It must not be connected to high pressure injectors
- The Nit-Occlud® PDA coil consists of a nickel-titanium alloy, which is generally considered safe. In non-clinical testing, nickel has been shown to be released from the device in very small amounts. Nickel is a known sensitizing agent, therefore the medical history of the patient regarding nickel or other metal allergies should be carefully established. Patients should be instructed to report after implantation of the device any symptoms that might indicate an allergic reaction, e.g. difficulty of breathing or inflammation of the face or throat. Non-clinical testing acc. to ISO 10993-10 has shown the device to be non-irritating and non-sensitizing, however there is still a small risk that the nickel released from the device may cause patients to develop an allergy to nickel.
- Nickel is also known to be a genotoxic and carcinogenic agent under certain conditions. Non-clinical testing acc. to ISO 10993-3 has shown the device to be non-mutagenic and not genotoxic. Based on the results of non-clinical testing, the amount of nickel expected to be released from the device after implantation is only a very small fraction of the patient's dietary nickel intake and considerably below the tolerable intake suggested by various Regulatory Agencies/Advisory Bodies.

Product Identification

Each product label has peel-off labels, to allow the product to be identified precisely. These can be used for the patient file and the patient ID card.

MRI Compatibility

The Nit-Occlud® PDA coil was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005. Non-clinical testing demonstrated that the Nit-Occlud® PDA coil is MR conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

- In non-clinical testing, the Nit-Occlud® PDA coil produced the following temperature rise during MRI performed for 15-min in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change +1.6°C.

Therefore, the MRI-related heating experiments for the Nit-Occlud® PDA coil at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

Artifact Information

- MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Nit-Occlud® PDAcoil. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	369 mm ²	118 mm ²	647 mm ²	739 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Indication

The Nit-Occlud® PDA coil for PDA Occlusion is a permanently implanted prosthesis indicated for percutaneous, transcatheter closure of patent ductus arteriosus.

Contraindications

Medical conditions that exclude implantation of a Nit-Occlud® PDA coil include:

- Aneurysmal PDA with risk of rupture
- Endocarditis, endarteritis or active infection at the time of the implantation
- History of contrast medium or nickel allergy
- History of hemorrhagic disease (coagulopathy)
- Patients with a body weight < 5 kg
- Pulmonary hypertension (calculated PVR greater than 5 Wood Units)
- Thrombus in a blood vessel through which access to the PDA must be obtained
- Thrombus in the vicinity of the implantation site.

Possible Complications

- Air embolism / Thromboembolism
- Arrhythmia (implantation of a pacemaker may be required)
- Blood loss requiring transfusion
- Damage to the heart, heart valves or blood vessels
- Death
- Embolization of the occluder requiring percutaneous or surgical intervention
- Endarteritis
- Heart failure
- Hemolysis
- Hypotension or shock
- Infection
- Myocardial infarction
- Occluder fracture or damage
- Stenosis of the left pulmonary artery or descending thoracic aorta


Precautionary Measures

- An angiogram must be performed prior to implantation for measuring length and diameter of the PDA.
- The implantation catheter must be flushed with heparinized saline solution prior to introduction and during the procedure, especially after angiography.
- The pfm medical implantation catheter is specifically designed for the delivery system. Other catheters should not be used to implant the device.
- Contrast media should not be injected through the implantation catheter.
- The coil should not be pulled back into the implantation catheter using strong force.
- Administration of 50 units of heparin per kg body weight is recommended after femoral sheaths are placed.
- Antibiotic coverage before (1 dose) and after implantation (2 doses) is recommended in order to prevent infection during the implant procedure. Antibiotic prophylaxis should be performed to prevent infective endocarditis during first 6 months after coil implantation.


Description of the Procedure

A lateral aortogram should be performed for measurement of PDA dimensions (see Fig. PDA Measurements):

Fig.: PDA Measurements

- D1 = Narrowest diameter
- D2 = Aortic Ampulla diameter
- L3 = PDA length

- D = Distal diameter
- P = Proximal diameter
- Lc = Length configured

Coil Selection

According to the measurements, the ductus type and the following recommendations, an appropriate coil should be selected:

- The distal coil diameter D should be no more than 2 mm larger than D2.
- The distal coil diameter D should be at least 3 to 4 mm larger than D1.
- Length of the configured coil Lc (see product label) should be not longer than L3.

D1	D2	Device
1mm	≤ 3mm	4x4
1mm	4mm	5x4
1mm	≥ 5mm	6x5
2mm	≤ 5mm	6x5
2mm	6-7mm	7x6
2mm	≥ 8mm	9x6

D1	D2	Device
3mm	≤ 7mm	7x6
3mm	8-9mm	9x6
3mm	≥ 9mm	9x6 or 11x6
4 or 5mm	9mm	11x6 or 10x6
4 or 5mm	10-11mm	11x6 or 12x6
4 or 5 mm	≥ 12mm	11x6 or 14x6

Table: Selection of the Nit-Occlud® PDA (according to angiographic PDA dimensions)

Coil Implantation Sequence

The Stiff type of the Nit-Occlud® PDA coil is preloaded in a 5F implantation catheter, therefore steps 01 and 02 are not applicable. After placement of an appropriate long introducer sheath with tip across the PDA into the descending thoracic aorta, the preloaded system is flushed with heparinized saline solution and advanced through the sheath. Starting at step 03, the Stiff type is deployed using the same technique as the Flexible and Medium types.

Step: 01


- Unpack the Nit-Occlud® PDA system consisting of coil with disposable handle and implantation catheter under sterile conditions.

WARNING

Do not pull on the delivery system. If the coil is withdrawn into the Y connector, there is a danger that the system can no longer be loaded. Check all screw connections. Some screw joints may have been loosened by the sterilization process.

- Flush the system carefully through the side access of the Y connector with heparinized saline solution, and ensure that there is no air remaining anywhere in the system.
- Check the coil position inside the transparent transportation sheath. The coil should be inside the sheath. When it is in this position, it is essential not to pull on the delivery system. If the coil is not positioned inside the transportation sheath, or shows visible signs of damage, it must be replaced with a new coil.

Step: 02 Coil Introduction


- Using a soft guide wire, advance the implantation catheter from the right femoral vein through the right heart, across the PDA into the descending thoracic aorta.
- Remove the guide wire from the implantation catheter and flush the catheter with a heparinized saline solution.
- Attach the luer lock connector of the transportation sheath to the implantation catheter.
- Open the hemostatic valve of the Y connector. The coil is now free for advancement.

WARNING

☞ **Do not pull on the delivery system in this position!**

- Advance the coil into the implantation catheter.


Step: 03 **Coil Configuration**

- Under fluoroscopic control advance the coil carefully through the implantation catheter into the aorta. This is done by moving the delivery system forward while arresting the implantation catheter.

WARNING


☞ **Ensure that the pigtail aortography catheter does not become entangled with the loops of the coil.**

- Advance the coil until the first marker M1 is positioned close to the Y connector. At this position, all but one loop is configured outside the implantation catheter.

Step: 04 **Coil Positioning**

Retract the entire system (implantation catheter, delivery system) under fluoroscopic control until the configured coil is positioned in the ampulla of the ductus (close hemostatic valve of Y connector or fix implantation catheter against delivery system).

NOTE: For longer ductus types, coil configuration inside the ductus ampulla is recommended. Here, 2-3 windings of the coil must first be configured in the aorta. Then the entire system is pulled into the ductus ampulla for further configuration of the coil.

Step: 05 **Final Coil Adjustment**

- Open the hemostatic valve of the Y connector.
- Configure the last 1 or 2 loops on the pulmonary side of the ductus by simultaneously pulling back the implantation catheter (with your left hand) and pushing the delivery system (with your right hand). Advance the delivery system until the second marker M2 is close to the Y connector. At this position the coil is outside the catheter.
- Perform an aortogram to confirm that the coil is in the correct position.

NOTE: If the position or size of the coil is not satisfactory, it should be repositioned or exchanged at this point.

Repositioning

- To reposition the coil, pull it back into the implantation catheter by pulling the delivery system.

WARNING

☞ **Close the gap between delivery system and coil before you retrieve the coil into the implantation catheter.**

- To do so, hold the handle with one hand and move the delivery system gently forward while holding it between 2 fingers of the same hand. This movement closes the gap between coil and delivery system and should be done under fluoroscopic control. Once the gap is closed the coil can be pulled smoothly into the implantation catheter.

WARNING

☞ **If a strong resistance is encountered while pulling the delivery system into the catheter, do not pull the system very hard because you risk a premature release of the coil.**


- To reposition the implantation catheter, the coil should be pulled back into the transportation sheath, carefully and under visual control, until the tip of the coil is in line with the marker at the distal end of transportation sheath. Fix the coil position by closing the Y connector.

WARNING

☞ **If the coil is pulled back too far, there is the risk that it may not be possible to reload it into the delivery system.**

- Then flush the implantation catheter with heparinized saline solution and repeat the procedure from Step 02.

Step: 06 Coil Release

- When the coil is properly positioned, it should be released. The rotation screw should lie directly against the pusher ball. If there is any gap between the two, it must be closed.

WARNING

Final release should only be performed if the coil is properly positioned in the PDA. Otherwise, the coil must either be retrieved and repositioned or replaced by an appropriate substitute. Before the coil is finally released, proper position of the coil should be confirmed by angiography.

- Remove the safety clip from the handle.
- Turn the rotation screw under fluoroscopy clockwise until the coil is released. Note that depending on the coil type between 8-15 rotations are needed to release the coil. You will feel an increase in resistance immediately prior to release.
- Remove the delivery system and implantation catheter.
- Perform a final aortogram about 10 minutes later to document position of the coil and PDA occlusion.
- Remove the aortography catheter.

WARNING

Ensure that the catheter does not touch the coil.

Technical Complications and how to avoid them

Complications may be avoided or ameliorated by the following:

- Use of venous access for implants.
- Flushing all components with heparinized saline.
- Keeping delivery system and catheter straight, avoiding loops and curves on the catheterization table.
- Coil recapture in the pulmonary artery or aorta, avoiding pulling the exposed coil across heart valves or through the right ventricle.

Failure of detachment:

Complications may arise if the coil is not released successfully. "Sticking" may occur if positioning of the coil is very time-consuming and/or if the delivery system does not protrude far enough from the end of the implantation catheter. The distal part of the delivery system must be placed outside of the implantation catheter immediately before release. If, however, the coil "sticks", the device must be retrieved and exchanged. As done with all interventional instruments, prior to implantation of the exchanged device the catheter should be flushed thoroughly to prevent coagulation, thus avoiding elevated friction or "sticking" of the system.

Device Embolization/ Premature release

The coil may embolize into the pulmonary artery if the aortic cone is too small or the coil fit is too loose. This can be prevented by accurate measurement of the PDA dimensions and choice of an appropriate coil for the PDA. Correct calibration of the angiographic measurement is a very important factor. In case of coil embolization, interventional retrieval should be performed using a snare or a biopptome.

In the event that the coil embolizes and interventional retrieval is unsuccessful, surgical retrieval should be considered.

Protrusion/ Obstruction

Protrusion of the windings into the aorta and/ or into the pulmonary artery may cause blood flow disturbances or vessel stenosis. This is avoidable by choosing an appropriate coil with a configured length (Lc) equal to or less than the PDA length (L3). The recoil force of the coil tends to return the device to its original configuration whenever possible and will retract the windings into the ampulla and/ or against the vessel wall. Pulmonary artery protrusion may be avoided by correct coil positioning during implant.

Late complications

Delayed complications such as migration or protrusion with a significant blood flow disturbance may require surgical removal of the coil.

Advice to Patients

Vigorous sports and high impact activities (e.g. gymnastics) should be avoided for 4-6 weeks after implantation.

Disposal after Use

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognised medical procedure, and in compliance with the relevant legal regulations and local ordinances.

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