Suction Tube Frazier Electromagnetic

REF 8000-050-005

Suction Tube Eicken Electromagnetic

REF 8000-050-006

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	23
DE Gebrauchsanweisung	45
EL Οδηγίες χρήσης	67
ES Instrucciones de uso	89
FI Käyttöohjeet	111
FR Notice d'utilisation	133
IT Istruzioni per l'uso	155
JA 使用説明書	177
KO 사용 설명서	199
NL Gebruiksaanwijzing	221
NO Bruksanvisning	243
PL Instrukcja użycia	265
PT Instruções de utilização	287
RO Instrucțiuni de utilizare	309
SV Bruksanvisning	331
TR Kullanım talimatları	353
ZH 使用说明	375

1. How to Use this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the user manual of these devices is to be considered as well.

This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient or medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.1. Symbol Definition

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol

Name: Definition

General warning sign: To signify a general warning.

W001

Refer to instruction manual/booklet: To signify that the user instruction manual/booklet must be read.

M002

ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1 General requirements

Symbol/
Number

Name: Definition

5.1.1

Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

5.1.3

Date of manufacture: Indicates the date when the medical device was manufactured.

LOT

5.1.5

Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.

REF

5.1.6

Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.

Symbol/
Number

Name: Definition

5.1.7

Serial number: Indicates the manufacturer's serial number so that the medical device can be identified.

5.2.7

Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.

5.3.2

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

5.3.4

Keep dry: Indicates a medical device that needs to be protected from moisture.

Symbol/
Number

Name: Definition

5.3.7

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

5.3.8

Humidity limitation: indicates the range of humidity to which the medical device can be safely exposed.

5.3.9

Atmospheric pressure limitation: indicates the range of atmospheric pressure to which the medical device can be safely exposed.

5.4.3

Consult instructions for use: Indicates the need for the use to consult the instructions for use.

Symbol/
Number

Name: Definition

5.4.10

Contains hazardous substances: Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties.

5.7.7

Medical device: Indicates the item is a medical device.

IEC 60417 Graphical symbols for use on equipment

Symbol

Name: Definition

5333

Type BF applied part: To identify a type BF applied part complying with IEC 60601-1.

Product-Specific Symbols

Symbol	Name: Definition
	Quantity: Indicates the number of medical devices in the packaging.
	Note symbol: Used to supplement or clarify information.
GTIN	Global Trade Item Number

81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Name: Definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol

Name: Definition

Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

Regulatory marks and logos

Symbol

Definition

Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

2. Safety Information

WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
 - Unauthorized modifications of the product are forbidden for safety reasons.
 - Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
 - The product is a high precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
 - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
 - Do not use the system during the discharge of a defibrillator. All navigated instruments must be removed from the patient and the operating area beforehand.
-

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Intended Use

3.1.1. EU

Suction Tube Frazier Electromagnetic and Suction Tube Eicken Electromagnetic are accessories to the Electromagnetic Navigation Unit and are intended for locating anatomical structures and removing fluids and tissue within ENT procedures.

3.1.2. Outside the EU

Suction Tube Frazier Electromagnetic and Suction Tube Eicken Electromagnetic are accessories to the Electromagnetic Navigation Unit and are intended for pointing at anatomical structures.

3.2. Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

3.3. Contraindications

None known.

3.4. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview


Figure 1: Suction Tubes Electromagnetic

- A Suction Tube Frazier Electromagnetic
- B Suction Tube Eicken Electromagnetic
- C Polymer tubing
- D Sensor tube

Navigation Point

Each navigated instrument has a navigation point. The navigation point is the position that is tracked by the navigation system. The navigation point of the Suction Tube Electromagnetic is visualized in Figure 2.


Figure 2: Visualization of the tracked position

- A Sensor tube
- B Navigation point
- C Tip of the suction tube

Product specifications

	Suction Tube Frazier Electromagnetic	Suction Tube Eicken Electromagnetic
Tube length	139 mm	111 mm
Diameter	4.6 mm	4.8 mm (Olive)
Inner tube diameter	2.5 mm	2.5 mm

5. For Use With

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to compatible software applications, refer to the user manual supplied with the respective software application.

6. Instructions for Application

WARNING

- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
 - Make sure that instruments are not damaged by third party products. Drills and shavers may damage the polymer tubing.
 - Make sure that the polymer tubing is undamaged. If any contamination is detected under the polymer tubing, the product must not be used.
-

When using the product in combination with suction devices, adhere to the following warnings:

WARNING

- Carefully determine the amount of negative pressure that is applied to the patient. Use only the minimum amount of negative pressure necessary to accomplish the suctioning procedure. Failure to comply may result in tissue damage.
 - Check the negative pressure of the vacuum unit before using the suction tubes. Check all connections for any leakages. Too low negative pressure may lead to insufficient suction flow.
 - Do not apply the suction tip too close to blood vessels or vital organs; contact with vital organs or blood vessels may result in serious patient injury.
-

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

7. Reprocessing

For reprocessing instructions, safety directives, and reprocessing equipment, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706).

8. Maintenance


Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition. Sharp or pointed components need to be sent back in a protected state.

9. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

-  • After using the electromagnetic system, make sure that instruments that cannot be reprocessed again are disposed of.
- In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) electronic products should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

10. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature		
Relative humidity		
Atmospheric air pressure		-

Materials

- This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0.
Current scientific evidence supports the position that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- The plug used in the electronic products contains the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
- The silicone glue used in this product contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).

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

Miljøbegrensninger	Drift	Oppbevaring og transport
Temperatur	 A line graph showing a temperature range from 10 °C to 30 °C. The line starts at 10 °C, rises to a peak at 30 °C, and then falls. A small black circle is positioned at the peak of the line.	 A line graph showing a temperature range from -10 °C to 50 °C. The line starts at -10 °C, rises to a peak at 50 °C, and then falls. A small black circle is positioned at the peak of the line.
Relativ luftfuktighet	 A circular icon containing a percentage sign and wavy lines representing water. A line points to the top of the circle labeled 75 %, and another line points to the bottom of the circle labeled 30 %.	 A simple icon of an open umbrella with rain falling from it.
Atmosfærisk trykk	 A circular icon containing a central dot and two horizontal arrows pointing outwards. A line points to the top of the circle labeled 106 kPa, and another line points to the bottom of the circle labeled 80 kPa.	-

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.
Copyright © 2022 Stryker
Stryker Leibinger GmbH & Co. KG