Instrument Clamp, Forceps	REF 8000-060-010
Instrument Clamp, 2-6 mm	REF 8000-060-011
Instrument Clamp, 6-10 mm	REF 8000-060-012
Instrument Clamp, 10-16 mm	REF 8000-060-013

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	18
DE Gebrauchsanweisung	35
EL Οδηγίες χρήσης	52
ES Instrucciones de uso	69
FI Käyttöohjeet	86
FR Notice d'utilisation	103
IT Istruzioni per l'uso	120
JA 使用説明書	137
KO 사용 설명서	154
NL Gebruiksaanwijzing	171
NO Bruksanvisning	188
PL Instrukcja użycia	205
PT Instruções de utilização	222
RO Instrucțiuni de utilizare	239
RU Инструкции по применению	256
SV Bruksanvisning	273
TR Kullanım talimatları	290
ZH 使用说明	307

1. How to Use this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the respective software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well.

This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.1. Symbol Definition

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/number Name: Definition

W001

General warning sign: To signify a general warning.

ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Symbol/Number Name: Definition

5.1.1

Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

5.1.3

Date of manufacture: Indicates the date when the medical device was manufactured.

Symbol/Number Name: Definition

5.1.5

Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.

5.1.6

Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.

5.2.7

Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.

5.3.2

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

5.3.4

Keep dry: Indicates a medical device that needs to be protected from moisture.

Symbol/Number Name: Definition

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

5.3.7

Consult instructions for use: Indicates the need to consult the instructions for use.

5.4.3

Medical device: Indicates the item is a medical device.

5.7.7

Product-Specific Symbols

Symbol Name: Definition

Quantity: Indicates the number of medical devices in the packaging.

Symbol	Name: Definition
--------	------------------

Note: Used to supplement or clarify information

GTIN	Global Trade Item Number
-------------	--------------------------

81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Name: Definition
--------	------------------

Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
----------------	--

Regulatory marks and logos

Symbol	Definition
--------	------------

Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

2. Safety Information

WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
 - Unauthorized modifications of the product are forbidden for safety reasons.
 - Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
 - The product is a high-precision instrument. Avoid subjecting the product to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
 - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
-

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Intended Use

The Instrument Clamp 2-6 mm, Instrument Clamp 6-10 mm, Instrument Clamp 10-16 mm, and Instrument Clamp Forceps are accessories to the Electromagnetic Navigation Unit and is intended for navigating conventional surgical instruments.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical

structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.2. Cranial Guidance System

3.2.1. United States

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview


Figure 1: Instrument Clamps

- A** Instrument Clamp, 2-6 mm
- B** Instrument Clamp, 6-10 mm
- C** Instrument Clamp, 10-16 mm
- D** Instrument Clamp, Forceps

5. Assembly

CAUTION

- Make sure that the attachment of the instrument clamp does not restrict the functionality and moveability of the surgical instrument to be navigated. Otherwise, this instrument must not be navigated via the instrument clamp.
 - If it is not possible to fix the clamp on the instrument without damaging it, the instrument is not suitable for use. Do not use any tools for tightening the fastening nut.
-

The instrument clamp allows you to mount a compatible instrument tracker on surgical instruments that are not originally intended for navigation. To attach an instrument clamp to a surgical instrument, follow steps 1–3.

1. Unscrew the fastening nut counter-clockwise.

2. Press the instrument clamp together and position the instrument clamp on the surgical instrument. To ensure system accuracy, select a position where the instrument tracker can be placed close to the navigated point. The closer the instrument tracker is placed to the navigated point, the more precise the navigation information. For recommended maximum distance values between instrument tracker and navigated point, refer to the instructions for use supplied with the respective software application.

Figure 2: Position of the instrument clamp in relation to the navigated point

3. Tighten the fastening nut clockwise until the clamp is rigidly attached to the surgical instrument.


Figure 3: Tightening of the fastening nut

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

6. For Use with

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to compatible software applications, refer to the user manual supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

Product	For use with
Instrument Clamp, 2-6 mm REF 8000-060-011	Rigid instruments with a round or multi-angular profile. Diameter of instrument: min. 2 mm, max. 6 mm
Instrument Clamp, 6-10 mm REF 8000-060-012	Rigid instruments with a round or multi-angular profile. Diameter of instrument: min. 6 mm, max. 10 mm
Instrument Clamp, 10-16 mm REF 8000-060-013	Rigid instruments with a round or multi-angular profile. Diameter of instrument: min. 10 mm, max. 16 mm
Instrument Clamp, Forceps REF 8000-060-010	Rigid instruments with a round or multi-angular profile. Diameter of instrument: min 5.5 mm, max. 10 mm

7. Disassembly

1. Detach the instrument tracker from the instrument clamp.
2. Detach the surgical instrument from the instrument clamp.

8. Reprocessing

During reprocessing, keep the fastening nut in an open position.

If you use the dedicated Passive Trackers Insert Tray REF 8000-810-000* for cleaning, disinfection and sterilization, refer to the Guide for Cleaning, Disinfection and Steam-Based Sterilization (TD6000005750) for reprocessing instructions, safety directives, and reprocessing equipment. The cleaning group of the instrument clamps is: IV.

If you use sterilization pouches for sterilization, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706) for reprocessing instructions, safety directives, and reprocessing equipment.

* Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.



9. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition.

10. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

11. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature	 10°C — 30°C	 -10°C — 50°C
Relative humidity	 30% — 75%	
Atmospheric air pressure	 80 kPa — 106 kPa	—

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Stryker. All other trademarks are
trademarks of their respective owners
Copyright © 2022 Stryker
Stryker Corporation or its divisions
Stryker Leibinger GmbH & Co. KG