Patient Tracker Tabs

REF 8000-100-001

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	13
DE Gebrauchsanweisung	25
EL Οδηγίες χρήσης	37
ES Instrucciones de uso	49
FI Käyttöohjeet	61
FR Notice d'utilisation	73
IT Istruzioni per l'uso	85
JA 使用説明書	97
KO 사용 설명서	109
NL Gebruiksaanwijzing	121
NO Bruksanvisning	133
PL Instrukcja użycia	145
PT Instruções de utilização	157
RO Instrucțiuni de utilizare	169
RU Инструкции по применению	181
SV Bruksanvisning	193
TR Kullanım talimatları	205
ZH 使用说明	217

1. How to Use this Document

1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the respective software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and/or healthcare staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.2. Symbol Definition

EN ISO 7010 Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/
Number

Name: Definition

W001

General warning sign: To signify a general warning.

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1 General Requirements

Symbol/
Number

Name: Definition

5.1.1

Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

5.1.3

Date of manufacture: Indicates the date when the medical device was manufactured.

Symbol/
Number

Name: Definition

5.1.4

Use-by date: Indicates the date after which the medical device is not to be used.

5.1.5

Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.

5.1.6

Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.

5.2.7

Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.

5.3.2

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

Symbol/
Number

Name: Definition

5.3.4

Keep dry: Indicates a medical device that needs to be protected from moisture.

5.3.7

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

5.4.2

Do not re-use: Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.

5.4.3

Consult instructions for use: Indicates the need to consult the instructions for use.

5.7.7

Medical device: Indicates the item is a medical device.

Product-Specific Symbols

Symbol	Name: Definition
--------	------------------

Quantity: Indicates the number of medical devices in the packaging.

Note: Used to supplement or clarify information.

81 FR 38911 FDA Final rule for the use of symbols in labeling

Symbol	Name: Definition
--------	------------------

Rx Only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Regulatory marks and logos

Symbol	Definition
--------	------------

Indicates that a device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation providing for its affixing.

2. Safety Information

WARNING

- Unauthorized modifications of the product are forbidden for safety reasons.
 - The product is only for single-use and must not be reused.
 - The product is not sterile.
 - Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects are detected.
 - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
-

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Intended Use

The Patient Tracker Tabs are an accessory to the Electromagnetic Navigation Unit and are intended for fixation of the Patient Tracker Electromagnetic to intact skin.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy,

can be identified relative to a CT- or MR-based model of the anatomy.
Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

* Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.2. Cranial Guidance System

3.2.1. United States

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. For Use With

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to compatible software applications, refer to the instructions for use supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

Description	REF
Patient Tracker Electromagnetic	8000-040-001
Patient Tracker Electromagnetic- 10 uses	8000-040-002

5. Instructions for Application

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

For surgical procedures with Cranial Guidance Software application:

WARNING

To avoid the risk of Medical Adhesive Related Skin Injuries (MARSI), consider using alternative adhesives for patients with fragile skin, e.g. pediatric patients. Refer to the user manual of the Cranial Guidance Software for more information on appropriate adhesives.

6. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker. All other trademarks are trademarks of their respective owners or holders.
Copyright © 2022 Stryker
Stryker Leibinger GmbH & Co. KG