



## ABO Discrepancies



Lynsi Rahorst, MHPE, MT(ASCP)SBB<sup>CM</sup>  
Manager, Education and Training IRL/Genomics  
NYBC Enterprise  
lrahorst@cbckc.org



---

---

---

---

---

---

---

---

## *New York Blood Center Enterprises*

EXPANDING OUR ORGANIZATION TO MEET CLINICAL, CELLULAR AND TRANSFUSION PRODUCT AND SERVICE NEEDS FOR PATIENTS. NOW PROVIDING ALMOST ONE MILLION BLOOD PRODUCTS, OVER 450,000 LABORATORY AND MULTI-ASSAY INFECTIOUS DISEASE TESTS AND OVER 12,500 SPECIALTY CLINICAL PROCEDURES ANNUALLY TO HOSPITALS NATIONWIDE.



---

---

---

---

---

---

---

---

## Objectives

- Identify and describe several causes of ABO discrepancies.
- List techniques used to resolve ABO discrepancies.
- Arrive at appropriate ABO interpretations based on laboratory results.



---

---

---

---

---

---

---

---

### ABO Discrepancies

- Weak/missing reactivity
  - Front type
  - Back type
- Extra reactivity
  - Front type
  - Back type

	Front type		Back Type	
	Anti-A	Anti-B	A1 cell	B cell
<b>Weak reactivity</b>	1+	0	0	4+
	0	0	2+	4+
<b>Missing reactivity</b>	0	0	0	4+
	0	4+	0	0
<b>Additional reactivity</b>	1+	1+	4+	4+
	4+	0	2+	4+




---

---

---

---

---

---

---

---

---

---

### ABO Discrepancies: Where to begin?

- Correct sample?? Technical error??
- Patient History
- Very young/old: weak ABO antibodies
  - Alloantibody that might interfere with reverse grouping
  - Strong cold autoantibody
    - May interfere with both forward and reverse grouping
  - Bone Marrow Transplant
  - Recent transfusion
  - Diagnosis
    - Weak antigens in leukemia, pregnancy, cord samples




---

---

---

---

---

---

---

---

---

---

### ABO Discrepancies: How to resolve

- Additional reagents**
- Front type: Anti-A,B (different anti-A/B clones); Anti-A1
  - Back type: A<sub>2</sub>, O, autocontrol
- Wash the cells**
- Removes plasma antibodies/proteins
  - Warm wash cells (37C saline): removes strong cold autos
- Saline replacement**
- Disperses rouleaux
- Pre-warming**
- Circumvent cold autoantibody or cold alloantibody interference
- Increase sensitivity**
- Increase plasma/cell ratio
  - Increase incubation time
  - Decrease temperature (ABO system "cold" reactive)
  - Adsorption/Elution
    - # antigen/cells too low for agglutination
    - Demonstrate by adsorbing and eluting anti-A or B from cells
    - Rarely Used: replaced by molecular testing

Must run appropriate controls for any testing outside limits of package insert.




---

---

---

---

---

---

---


---

---

---

## Case #1

Front Type				Back Type				Auto control
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested




---

---

---

---

---

---

---

---

---

---

## What's anti-A1 lectin?




---

---

---

---

---

---

---

---

---

---

## What's the difference between A<sub>1</sub> and A<sub>2</sub> phenotypes?

- **Frequency in European population**
  - ~80% group A individuals A<sub>1</sub>
  - ~20% group A individuals A<sub>2</sub>
  - Other A subgroups rare (A<sub>3</sub>, A<sub>x</sub>, A<sub>el</sub>, etc.)
- **Antigen differences**
  - Quantity of A antigens on cells
    - A<sub>1</sub> cells have approximately 5 times as many A antigens as A<sub>2</sub> cells
  - Qualitative differences
    - Antigens of A<sub>1</sub> individuals more branched
    - Why A subgroup individuals can make anti-A1.

Fung MK, Eder AF, Spillare SL, Westhoff CM. Technical Manual. 19<sup>th</sup> ed. Bethesda, MD: AABB; 2017: 271-272.




---

---

---

---

---

---

---

---

---

---

## More about anti-A1

- Detected in 1-8% of A<sub>2</sub> individuals, 22-35% of A<sub>2</sub>B individuals
  - Often "naturally occurring"
- Usually IgM, reacts best at room temperature or below.
  - Generally not considered clinically significant
- Reports in literature of hemolytic anti-A1
  - Helmsch F, et al. Acute hemolytic transfusion reaction due to a warm reactive anti-A1. Transfusion. 2018;58:1163-1170.
- Transfusion recommendations: XM compatible units
  - A<sub>2</sub> RBCs
  - O RBCs

Fung MK, Eder AF, Spitalnic SL, Westhoff CM. Technical Manual. 19<sup>th</sup> ed. Bethesda, MD: AABB; 2017: 274.




---

---

---

---

---

---

---


---

---

---

## Case #1

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested




---

---

---

---

---

---

---

---

---

---

## Does the patient's plasma contain cold autoantibody?

	Rh				Kell		Duffy		Kidd		M	MNS			Results	
	D	C	E	c	K	k	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>		N	S	s	5' RT	LISS 37C
1																
2																
3																
4	+	0	0	+	+	0	+	0	+	0	+	+	0	+	1+	0
5	0	+	0	+	+	0	+	0	+	0	+	+	0	0	1+	0
6	0	0	+	+	0	+	+	0	+	0	+	+	0	+	1+	0
7	0	0	0	+	+	+	+	0	+	0	+	+	0	+	1+	0
8	0	0	0	+	+	+	+	0	+	0	+	+	0	+	1+	0
9	0	0	0	+	+	+	+	0	+	0	+	+	0	+	1+	0
10	I-	+	+	0	+	+	+	0	+	+	+	+	0	0	0	0
11	+	0	0	+	+	+	+	0	+	+	+	+	0	+	1+	0
Auto															1+	0

1 nonreactive cell at RT: I-negative cell

Autocontrol reactive at 5' RT




---

---

---

---

---

---

---

---

---

---

# Cold Screen

Cold Screen		
	30' RT	30' 4C
SCI	2+	4+
SCII	2+	4+
I-negative	0	3+
Auto	2+	4+

---

---

---

---

---

---

---

---


---

---



# What is pre-warming?

- Incubate all reagent red cells and patient plasma (and pipet) at 37C prior to testing
- Add patient plasma to reagent cells quickly without changing the temperature of the testing environment from 37C
- Incubate all tubes at 37C (~30 min)
- DO NOT CENTRIFUGE! (centrifugation will quickly cool the sample in the tubes)
- Shake and read settled tubes




---

---

---

---

---

---

---

---

---

---

# One word of caution!

**NEVER utilize pre-warm testing unless you know what you are pre-warming!!**

(demonstrate that the patient has cold autoantibody prior to pre-warming)




---

---

---

---

---

---

---


---

---

---

## Case #1

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested




---

---

---

---

---

---

---

---

---

---

## Scenario #3

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	1+	4+	1+	1+

These same results (same as scenario #2) may be due to rouleaux!!

- In vitro phenomenon due to abnormal patient plasma protein concentration
- Can be seen in any test involving patient plasma (including back type)
- "Stack of coins," refractile aggregation of RBCs
- Can be mistaken for agglutination macroscopically

Fung MK, Eder AF, Spillner SL, Westhoff CM. Technical Manual. 19<sup>th</sup> ed. Bethesda, MD: AABB; 2017. 370-371.




---

---

---

---

---

---


---

---

---

---

## Can you see the difference?



<https://www.atdove.org/video/saline-agglutination-test>




---

---

---

---

---

---


---

---

---

---

### Steps to saline replacement...




---

---

---

---

---

---

---

---

---

---

### Steps to saline replacement...




---

---

---

---

---

---

---


---

---

---

### Case #1

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested




---

---

---

---

---

---


---

---

---

---

### Cells for back type: Pooled



Pooled!




---

---

---

---

---

---

---

---

---

---

---

---

### Identifying the alloantibody

	Rh					Kell		Duffy		Kidd		MNS			Results			
	D	C	E	c	e	K	k	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	M	N	S	s	S <sup>w</sup> RT	LISS 37C	LISS IAT
1	+	+	0	0	+	0	+	+	+	+	+	+	+	+	+	1+	0	0
2	+	+	0	0	+	+	+	0	0	+	0	+	0	+	0	0	0	0
3	+	0	+	+	0	0	+	+	0	+	+	+	0	+	+	2+	1+	1+
4	+	0	0	+	+	0	+	0	0	+	0	+	+	0	+	1+	0	0
5	0	+	0	+	+	0	+	+	0	+	0	+	+	0	0	1+	0	0
6	0	0	+	+	+	0	+	+	+	+	0	+	0	+	+	0	0	0
7	0	0	0	+	+	+	+	0	+	0	+	0	+	+	+	2+	1+	1+
8	0	0	0	+	+	0	+	+	+	0	+	0	+	+	+	0	0	0
9	0	0	0	+	+	0	+	+	+	0	+	+	0	0	+	2+	1+	1+
10	+	+	0	0	+	+	+	0	+	+	+	+	+	0	+	1+	0	0
11	+	0	0	+	+	+	0	0	+	+	0	+	+	+	+	0	0	0
Auto																0	0	0




---

---

---

---

---

---

---

---

---

---

---

---

### Scenario #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	1+	4+	1+	0

Once the antibody is identified, resolve the typing discrepancy by...

- Prewarming the back type
- Using RBCs for your back type testing that don't express the corresponding antigen
  - For example, M-negative A<sub>1</sub> cells & M-negative B cells
  - Enzyme-treated cells




---

---

---

---

---

---

---

---

---

---

---

---



## Review of the 4 scenarios:

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested

← Discrepant Results!

---

---

---

---

---

---

---

---

---

---



## Review of the 4 scenarios:

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested
				1+	0	4+	0	0

← Discrepant Results!

**Possible A subgroup with anti-A1**

- Test patient cells with anti-A1 lectin.

---

---

---

---

---

---

---

---

---

---



## Review of the 4 scenarios:

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested
				1+	0	4+	0	0
				1+	1+	4+	1+	1+

← Discrepant Results!

**Possible rouleaux**  
**Possible cold auto**

- Examine microscopically
- Saline replacement if rouleaux
- Prewarm if cold auto

---

---

---

---

---

---

---

---

---

---



## Review of the 4 scenarios:

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	0	4+	0	1+	Not tested	4+	Not tested	Not tested
				1+	0	4+	0	0
				1+	1+	4+	1+	1+
				1+	1+	4+	1+	0

← Discrepant Results!

**Cold-reacting alloantibody**

- Identify the alloantibody
- Prewarm back type
- Use antigen negative cells to resolve

---

---

---

---

---

---

---

---

---

---

---

---

## Case #2

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
0	0	4+	0	0	0	0	0	0

**Ways to promote/strengthen reactivity of back type:**

- Increase incubation time at 22C
- Decrease temperature\*
- Increase plasma:cell ratio (use 4 drops of plasma & 1 drop of cells in each tube)

\* Be careful: Many individuals have cold autoantibodies!

---

---

---

---

---

---

---

---

---

---

---

---

## Case #2

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
0	0	4+	0	0	0	0	0	0
				2+	1+	2+	0	0

4 drops plasma:1 drop cells, 15 minute 22C

**Patient blood type: O+**

---

---

---

---

---

---

---

---


---

---

---

---

## Wait... What is the Rh control?




---

---

---

---

---

---

---

---

---

---

## Strong cold agglutinins




---

---

---

---

---

---

---

---

---

---

## Case #3

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+	4+	4+	4+	4+	4+	4+	4+	4+
4+	0	4+	0	←Repeat testing following 4X warm wash				
Prewarming the back type →				0	0	2+	0	0

Patient blood type: A+

### Cold adsorption:

Adsorb cold autoantibody at 4C, test adsorbed plasma in the back type

---

---

---

---

---

---

---

---

---

---

# Mixed field




---

---

---

---

---

---

---

---

---

---

---

---

## Case #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+mf	0	4+	0	0	0	3+	0	0

What is the best explanation for these results?

- Patient is group O+ and received A+ blood ←
- Patient is group A+ and received O+ blood
- Patient is group A+ and received O- blood
- Patient is group A+ and received A- blood




---

---

---

---

---

---

---

---

---

---

---

---

## Case #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+mf	0	4+	0	0	0	3+	0	0

Looks like group A




---

---

---

---

---

---

---

---

---

---

---

---

### Case #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+mf	0	4+	0	0	0	3+	0	0



RBCs that express A antigen agglutinate

**Patient cells**



RBCs that don't express A antigen don't agglutinate

**Donor cells**




---

---

---

---

---

---

---

---

---

---

### Case #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+mf	0	4+	0	0	0	3+	0	0



**Patient cells**

All RBCs express D antigen



**Donor cells**




---

---

---

---

---

---

---

---

---

---

### Case #4

Front Type				Back Type				
Anti-A	Anti-B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
4+mf	0	4+	0	0	0	3+	0	0

**Patient is group A+ and received O+ blood**

**Important things to know!!!**

- ✓ Transfusion of non ABO-identical RBCs affects the **FRONT** type
- ✓ Transfusion does not usually interfere with the back type
- ✓ Use the mixed-field reactions to determine what type of blood patient received




---

---

---

---

---

---

---

---

---

---

## Weak A subgroup

Front Type					Back Type				
Anti-A	Anti-B	Anti-A,B	Anti-D	Rh Control	A <sub>1</sub>	A <sub>2</sub>	B	O	Auto control
0	0	1+	3+	0	0	0	4+	0	0

### Weak A subgroups:

- Examples: A<sub>3</sub>, A<sub>x</sub>, A<sub>m</sub>
- Infrequently encountered

### Resolving the discrepancy:

- Anti-A,B reagent
- **Genomic testing**
- Adsorption/elution studies (using polyclonal anti-A)




---

---

---

---

---

---

---

---

---

---

## What do we mean by adsorption/elution?



### Adsorption

Weak A subgroup RBCs




---

---

---

---

---

---


---

---

---

---

## What do we mean by adsorption/elution?




---

---

---

---

---

---

---

---

---

---



## Review: ABO discrepancies

- Weak ( $\leq 2+$ ) reactivity should be investigated
- Decide if missing or additional reactivity present
- Decide if problem is in front or back type
- Use the appropriate tool(s) to investigate
- If testing outside the parameters of the package insert, **RUN APPROPRIATE CONTROLS**

---

---

---

---

---

---

---

---



## Objectives

- Identify and describe several causes of ABO discrepancies.
- List techniques used to resolve ABO discrepancies.
- Arrive at appropriate ABO interpretations based on laboratory results.

---

---

---

---

---

---

---

---



## References

Fung MK, Eder AF, Spitalnic SL, Westhoff CM. Technical Manual. 19<sup>th</sup> ed. Bethesda, MD: AABB; 2017: 274.

Helmich F, et al. Acute hemolytic transfusion reaction due to a warm reactive anti-A1. Transfusion. 2018;58;1163-1170.

---

---

---

---

---

---

---

---







---

---

---

---

---

---

---

---