

Geotrichum **Fusobacterium** Lodderomyces Anaerococcus **Colletotrichum** Edwardsiella Comamonas Pigmentiphaga Moellerella
 Udeniomyces Erysipelothrix **Escherichia** Streptomyces **Bartonella** Hafnia Terrimonas Pseudoxanthomonas **Corynebacterium** Streptococcus
Facklamia Schizosaccharomyces Tetragenococcus Lechevalieria Clavibacter Shmwellia **Brachybacterium** Leclercia Providencia Trablusiella
 Gardnerella Sporothrix Leuconostoc **Pseudoclavibacter** Alkalibacillus Debaryomyces Turicella Roseomonas **Ruminococcus** Scedo
 Staphylococcus Peptostreptococcus Paenibacillus Balneatrix Solibacillus Prototheca Cupriavidus Geobacillus Aspergillus Arthrobacter Mes
 Filobasidium **Propioniferax** Azohydromonas Chromobacterium **Curtobacterium** Kloeckera Austwickia Hyphomicrobium **Cryptococcus** R
 Aquincola Enterobacter Sporobolomyces Brevundimonas Cannocytophaga Tatlockia Neisseria **Salinivibrio** Pullulanibacillus Arc
Eggerthella Methanomonas Mucor **Mobiluncus** Caulobacter Helcococcus **Psychrobacillus** Campylobacter Blastomonas Wohlfahrtiim
 Herminiimonas Tsukamurella Mycobacterium Bordetella Pichia Vibrio Iodobacter Tenacibaculum Listeria **Plesiomonas** Haloarcula S
 Thauera Viridibacillus Yokenella Malassezia **Novosphingobium** Ornithobacterium **Epidermophyton** Oligella Paracoccus Aureobasidi
 Salimicrobium Klebsiella Mycoplasma Variovorax Samsonia Schizophyllum Scopulariopsis Odoribacter Anaerotruncus **Abiotroph**
Empedobacter Sphingopyxis Lactococcus Sphingobium Microsporium Peptoniphilus Vagococcus **Beauveria** Morganella Pasteurella
Micrococcus Propionimicrobium Starkeya Prevotella Histophilus Sphingomonas Acetobacter **Francisella** Photobacterium Propioniba
 Arthrographis Aromatoleum Pedicoccus Phoma Xenorhabdus Methylobacillus Fusarium **Wolinella** Bacteroides Zygosaccharom
Helicobacter Rhizobium Terrabacter Ralstonia **Butyricimonas** Microsporium Castellaniella Borrelia Microbacterium Rheinheimera Wauter
 Rahnella Nocardioides Gluconobacter **Sphingobacterium** Mannheimia Cohnella Aggregatibacter **Cronobacter** Lecythophora Riemerella
 Rhizopus Acidovorax Rothia Kytococcus Chryseobacterium Alishewanella Gemella Methylobacterium Haemophilus Adlercreutzia
 Alloiococcus **Bacillus** Arxiozyma **Halococcus** **Rhodotorula** Pseudochrobactrum Lemnirella Candidatus Xanthomonas Pectoba
 Arthroderma **Slackia** Trueperella Inquilinus Brevibacillus Brachyspira Porphyromonas Aurantimonas **Actinomyces** Eikenella Kitasa
 Psychrobacter Acidiphilium Amycolatopsis Lactobacillus **Marinibacillus** Megamonas Dermatophilus Grimontia Acinetobacter Lysin
 Parvimonas Moesziomyces Legionella Aliivibrio **Dermacoccus** Exiguobacterium Virgibacillus Raoultella Gordonia **Dialister** Parabact
 Stenotrophomonas **Sporobolomyces** Coprobacillus Sporosarcina Brenneria Rathayibacter Arsenophonus Penicillium Pseudomor
 Alloscardovia Nocardia Halomonas Rhodococcus Bergeyella Mallica **Actinocorallia** Aeromonas **Micromonospora** Alcaligenes **Alistipe**
 Kocuria Ochrobactrum Agrococcus Gracilibacillus Chromohalobacter Yersinia Oerskovia Gallibacterium Erwinia Agromyces Filifactor
 Collinsella Finegoldia **Phenylbacterium** Methyloarcula Jonesia **Pantoea** Elizabethkingia Leifsonia Pseudozyma Streptosporangium
 Sporobolomyces **Fusobacterium** Lactococcus **Flavobacterium** **Colletotrichum** **Edwardsiella** **Cytolysis** **Sinomonas** **Caenorhynchus**
 Sphingopyxis **Escherichia** **Bartonella** **Mycobacterium** **Idoneella** **Corynebacterium** **Streptococcus** **Da**
Facklamia **Schizosaccharomyces** **Lechevalieria** **Clavibacter** **Shmwellia** **Brachybacterium** **Leclercia** **Providencia** **Trablusiella** **Geotrichum** **Erysipelothrix** **Escherichia** **Streptomyces** **Bartonella** **Hafnia** **Terrimonas** **Pseudoxanthomonas** **Corynebacterium** **Streptococcus** **Facklamia** **Schizosaccharomyces** **Tetragenococcus** **Lechevalieria** **Clavibacter** **Shmwellia** **Brachybacterium** **Leclercia** **Providencia** **Trablusiella** **Gardnerella** **Sporothrix** **Leuconostoc** **Pseudoclavibacter** **Alkalibacillus** **Debaryomyces** **Turicella** **Roseomonas** **Ruminococcus** **Scedo** **Staphylococcus** **Peptostreptococcus** **Paenibacillus** **Balneatrix** **Solibacillus** **Prototheca** **Cupriavidus** **Geobacillus** **Aspergillus** **Arthrobacter** **Mes** **Filobasidium** **Propioniferax** **Azohydromonas** **Chromobacterium** **Curtobacterium** **Kloeckera** **Austwickia** **Hyphomicrobium** **Cryptococcus** **R** **Aquincola** **Enterobacter** **Sporobolomyces** **Brevundimonas** **Cannocytophaga** **Tatlockia** **Neisseria** **Salinivibrio** **Pullulanibacillus** **Arc** **Eggerthella** **Methanomonas** **Mucor** **Mobiluncus** **Caulobacter** **Helcococcus** **Psychrobacillus** **Campylobacter** **Blastomonas** **Wohlfahrtiim** **Herminiimonas** **Tsukamurella** **Mycobacterium** **Bordetella** **Pichia** **Vibrio** **Iodobacter** **Tenacibaculum** **Listeria** **Plesiomonas** **Haloarcula** **S** **Thauera** **Viridibacillus** **Yokenella** **Malassezia** **Novosphingobium** **Ornithobacterium** **Epidermophyton** **Oligella** **Paracoccus** **Aureobasidi** **Salimicrobium** **Klebsiella** **Mycoplasma** **Variovorax** **Samsonia** **Schizophyllum** **Scopulariopsis** **Odoribacter** **Anaerotruncus** **Abiotroph** **Empedobacter** **Sphingopyxis** **Lactococcus** **Sphingobium** **Microsporium** **Peptoniphilus** **Vagococcus** **Beauveria** **Morganella** **Pasteurella** **Micrococcus** **Propionimicrobium** **Starkeya** **Prevotella** **Histophilus** **Sphingomonas** **Acetobacter** **Francisella** **Photobacterium** **Propioniba** **Arthrographis** **Aromatoleum** **Pedicoccus** **Phoma** **Xenorhabdus** **Methylobacillus** **Fusarium** **Wolinella** **Bacteroides** **Zygosaccharom** **Helicobacter** **Rhizobium** **Terrabacter** **Ralstonia** **Butyricimonas** **Microsporium** **Castellaniella** **Borrelia** **Microbacterium** **Rheinheimera** **Wauter** **Rahnella** **Nocardioides** **Gluconobacter** **Sphingobacterium** **Mannheimia** **Cohnella** **Aggregatibacter** **Cronobacter** **Lecythophora** **Riemerella** **Rhizopus** **Acidovorax** **Rothia** **Kytococcus** **Chryseobacterium** **Alishewanella** **Gemella** **Methylobacterium** **Haemophilus** **Adlercreutzia** **Alloiococcus** **Bacillus** **Arxiozyma** **Halococcus** **Rhodotorula** **Pseudochrobactrum** **Lemnirella** **Candidatus** **Xanthomonas** **Pectoba** **Arthroderma** **Slackia** **Trueperella** **Inquilinus** **Brevibacillus** **Brachyspira** **Porphyromonas** **Aurantimonas** **Actinomyces** **Eikenella** **Kitasa** **Psychrobacter** **Acidiphilium** **Amycolatopsis** **Lactobacillus** **Marinibacillus** **Megamonas** **Dermatophilus** **Grimontia** **Acinetobacter** **Lysin** **Parvimonas** **Moesziomyces** **Legionella** **Aliivibrio** **Dermacoccus** **Exiguobacterium** **Virgibacillus** **Raoultella** **Gordonia** **Dialister** **Parabact**

Industrial Microbiology

MALDI Biotyper

- Fast & Accurate Identification of Microorganisms

High Quality Identifications in just a Couple of Minutes



The MALDI Biotyper System:

- Sequence Quality Data
- Comprehensive Libraries
- Open Microbiology
- 21CFR part 11 support
 - Audit Trail
 - User Management
 - Data Security
 - Electronic Signature
- IQ/OQ/PV
- Robust and Easy to Use
- Compact Bench Top System

A Powerful Technology for Industrial Microbiology Applications

To help solve the challenges of performing microbiology in the industrial market, Bruker has utilized its wealth of experience to create the truly innovative MALDI Biotyper system.

Over the past 5 years the MALDI Biotyper has revolutionised how microorganism identification is performed in more than 700 laboratories around the world.

Identifying Microorganisms by Their Molecular Fingerprint

The MALDI Biotyper identifies micro-organisms using MALDI-TOF (Matrix Assisted Laser Desorption Ionization-Time of Flight) Mass Spectrometry to measure a unique molecular fingerprint of an organism.


This coupled with extensive libraries covering environmental and clinical isolates provides extremely broad coverage of microorganisms found in Food, Pharmaceutical and Water industries.



Installation, Operational and Performance Qualifications

Already in use in a number of Industrial microbiology laboratories, the MALDI Biotyper provides a very rapid and specific microorganism identification with resolution comparable to molecular sequencing techniques with significantly less effort, time and cost.

The implementation and validation of the MALDI Biotyper is assisted by IQ/OQ/PV documentation and 21 CFR Part 11 compliant software.



Quality Assurance

Drawing upon over 20 years of experience in the manufacture of reagents used in Mass Spectrometry, Bruker supplies both certificates of Analysis and Traceability for the Bacterial Test Standard and the HCCA Matrix.

The Bacterial Test Standard is a typical *E. coli* extract containing additional proteins that can be used for instrument mass calibration and as a performance verification standard.



Further reading

McDaniel, A. Validation of an Automated Microbial Identification System, In Microbial Identification: The Keys To A Successful Program. Griffin M. & Reber D. (Eds.) PDA Chapter 5,87-106

Technical Specifications

Dimensions & Operating Parameters

LxWxH:	510 x 680 x 1093mm [20.1" x 26.8" x 43"]
Weight:	84kg (185 lb) net weight
Noise:	<30 dB under normal operating conditions
Temp Range:	10-30°C (50-86°F)
Operating Humidity:	15-85% non-condensing @ 30°C

Instrument: Microflex LT

- Nitrogen Laser with 60Hz repetition rate
- Full Spectrum Resolution (FSR) with broadband focusing mode (PAN™)
- Smart Spectra Acquisition™
- Perpetual Ion Source™ with IR-laser self-cleaning functionality
- FlashDetector™
- Whispermode™
- Oil-free membrane pre-vacuum pump and turbo pump
- Manufactured under QSR regulations

Microbial Identification Applications:

- Gram +/- Bacteria, Yeast, Moulds, Fungi and Mycobacteria
- Direct from Liquid media

Computer, Software & Database:

- Windows 7 operating system with Quad-Core CPU 2.66 GHz, Laser printer and Remote Service Capability via 128-bit SSL
- **MALDI Biotyper** Database
- **MALDI Biotyper** Client Server

Optional System Upgrades and Accessories:

- GPR kit for direct processing of liquid samples
- IQ/OQ/PV
- Security pack

Sample Targets:

- Reusable Polished Stainless Steel Targets: 48 & 96 position with and without barcode
- Disposable 48 position Biotargets with individual barcode
- 24 and 96 position BigAnchorChip™ Targets

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
cam-sales@bruker.com



solarix

- Optimized FTMS Solutions for the Most Challenging Applications

Dedicated to the Most Challenging Applications



solariX, the next-generation line of hybrid Qq-FTMS systems, is the culmination of key technology enhancements that provide unique capabilities in mass spectral performance and versatility.


The technological advances encompassed by solariX provide:

- Superior Sensitivity
- Unmatched Mass Accuracy and Broadband Resolution
- Widest range of structural tools, including Electron Transfer Dissociation (ETD)
- Expansive Mass Range
- Selective Ion Enrichment and Enhanced Dynamic Range
- Application Directed and Optimized Solution Packages

Common applications for solariX

The analytical power and performance of FTMS is well suited to address some of today's most challenging and complex samples. Drawing on years of applications experience, we have combined our unique FTMS instrumentation and comprehensive software tools to provide turn-key solutions for the following areas:

- High End Proteomics Studies (Top-down and Bottom-up workflows)
- Molecular Imaging of Tissue - Distribution of Drugs, Metabolites, and Biomarkers
- Petroleum Product Analysis
- Complex Environmental sample analysis
- Metabolomics Research



● Performance Beyond Compare

At the core of the solariX is dramatically improved sensitivity and dynamic range. This allows researchers to identify and analyze a much wider range of molecules than ever before, and to delve deeper into complex mixtures and analyze lower abundance species.


More power for greater complexity

The broadband, ultra high resolving power (increased 8-fold) and superb mass accuracy of solariX is more powerful than previously possible with any other mass spectrometer. This extraordinary increase in the number of available m/z channels is essential for addressing complex problems such as petroleomics and environmental samples which require resolving powers of greater than 400,000 for effective analysis.


Faster and more advanced LC-MS and LC-MS/MS operation

The new data acquisition functionality enables smarter modes of data dependent operation. Here, spectral acquisition parameters such as data set size or starting mass may be adjusted on-the-fly in a mass dependent mode of operation. Smaller data sets can be selected for MS/MS acquisitions making the overall data acquisition rate faster, while maintaining the high fidelity measurements for the MS precursor acquisitions. Moreover, super stable mass accuracy is maintained throughout the LC-MS analysis using Bruker's proprietary Ion Charge Control (ICC™) mode.

Base Peak Chromatogram for LC-MS and LC-MS/MS of protein mixture



Effective analysis of complex problems



Broad band spectrum of a crude oil measured by APPI in positive ion mode. Insets illustrate the extreme high resolving power (> 550,000) of the solariX-CM.


New Horizons for Biomolecule Analysis

Expanded capabilities for biomolecule analysis


Adding to the existing arsenal of structural fragmentation tools, solarix is fully enabled with Electron Transfer Dissociation (ETD). This exciting new technique is superb for in depth, comprehensive analysis of proteins and peptides and their subtle, posttranslational modifications. For instance, the gentle molecular dissociation chemistry associated with ETD enables researchers to elucidate subtle post-translational modifications such as glycosylation and phosphorylation at levels of accuracy and resolution previously unavailable for such de-novo approaches. Furthermore, ETD can be automated, and combined with LC-MS/MS schemes in combination with quadrupole fragmentation (Q-CID).

Myoglobin, non-denaturing MS conditions

Resolving Power = 450,000 @ m/z 1953



ETD + ECD of Calmodulin



The combination of ECD and ETD performed on Calmodulin. As illustrated above, the expected phosphorylations are not present and the lysine at position 116 is trimethylated.


Definitive molecular identification

With the solariX, high sensitivity, superior quality, exact mass measurements are only a mouse click away. Leveraging sub-ppm levels of mass accuracy for both intact precursor (MS) and product ions (MS/MS) combined with accurate isotope patterns, SmartFormula3D™ provides definitive elemental composition and molecular formula information. This level of confidence is readily achieved without internal standards or recalibration and the high resolution data inherently mitigates the complications resulting from other chemical interferences.



Molecular formulae determination of components in crude oil


High resolution spectrum of crude oil by APPI demonstrating broadband resolving power of 950,000 with average mass measurement error below 100 ppb.



CASI™: For rapid, enhanced selectivity and dynamic range

The unique Qq-FTMS geometry of the solariX can enrich lower abundant or trace species for detection and structural interrogation via MS and MS/MS, respectively. This mode, known as Continuous Accumulation of Selected Ions (CASI), is essential for tissue imaging, analyzing low-copy PTMs, and to extend the general dynamic range of almost any measurement. This can improve signal intensities by as much as an order of magnitude.

CASI Figure





Schematic of CASI workflow where target ions are selected and enriched in the collision cell.



Solution packages optimized for a variety of applications


The solariX product line features a series of tailored solutions around specific applications. Each solution comprises the appropriate hardware, ion source options, and application level software for a turn-key analytical platform. Whether your primary application is complex mixtures (e.g. petroleomics) or high resolution tissue imaging, a solariX package is designed to meet your needs. For laboratories having unique or multiple analysis needs, Bruker will work with you to prepare an appropriate, customized configuration.

● High-End Performance for any Analytical Laboratory




Full mass range MS

$m/z = 313.1485$



CASI boosts sensitivity, enabling the detection of drug and metabolites at lower concentrations



CASI - 20 m/z window

Molecular imaging of olanzapine in kidney Dosage – 5 mg/kg, 6h post dosage analysis.

The solariX FTMS is as intuitive and easy to use as a benchtop instrument. Bruker's Compass™ software enables the researcher to harness the analytical power and versatility with ease so that the focus is on the application and not on the instrumentation.

Weekly cryogen “fills” become annual...

Bruker's patented refrigerated magnet technology means nitrogen-free compact superconducting magnets with very low helium losses to minimize instrument maintenance and service. Patented active shielding technology minimizes the stray magnetic field levels for compact installations and maximum laboratory safety. These magnets are available for the full range of magnetic field strengths (7T, 9.4T, 12T, and 15T).



● Unique Features of solariX

Ultimate versatility in structural analysis


The solariX leverages the complementary nature of traditional quadrupole based Collisional Induced Dissociation (CID), and now Electron Transfer Dissociation (ETD). Additionally, precursors of complex mixtures may be isolated and fragmented using high front-end resolution, in-cell isolation followed by Electron Capture Dissociation (ECD) or Sustained Off-Resonance Irradiation (SORI)-CID. Whether your application is natural product, peptide/protein, carbohydrate, or petroleum product analysis, solariX has an array of fragmentation tools to address practically any compound class.



Ion source flexibility

While FTMS is ideally suited to liquid introduction, Bruker Daltonics offers the unique dual ESI/MALDI source. This is based on patented Ion Funnel technology for maximum sensitivity which allows effortless switching between ESI and MALDI – all at the touch of a button! The intermediate pressure MALDI source offers exquisite sensitivity, and preserves molecular ion fidelity throughout the complete FTMS detection process. The efficiency of the dual ESI/MALDI source combined with the new ion optics of the solariX provides a marriage of ultra-high resolution with high-end “MALDI-TOF sensitivity”.


Along with our array of atmospheric pressure ionization sources (ESI, nano-ESI, APCI, and APPI), Bruker now offers GC-APCI capability.



● Flexible Front-end Solutions

Compass

The powerful Compass software enables the operation of a range of industry standard HPLC and u-HPLC platforms with the solariX. This includes comprehensive system control and seamless integration of mass spectral and UV data.



solariXcontrol™ - the instrument control software for solariX

Technical Specifications



Magnet 7.0T/US/R refrigerated magnet

Field strength 7.0T

Analytical Performance

Mass Range 100 – 10,000 m/z (transmission with RF only)
100 – 6,000 m/z (mass selective)

MS/MS Operation

Isolation efficiency (Qh-Interface) (LHRH, [M+2H]2+) > 60%

MS/MS Efficiency (Qh-Interface) (LHRH, [M+2H]2+) > 50% conversion from isolated precursor

Multistage MS (MS3 guaranteed) (LHRH) LHRH MS/MS (collision cell) -> MS/MS/MS (Infinity Cell™)

Mass Dependent MS/MS Automated isolation and MS/MS of the most intense ions in an LC-MS/MS run (scan ratio set to five)

Sensitivity ECD (Substance P) S/N > 10:1 for 5 fmol (consumed). c5 fragment @ m/z 624

ESI

Mass accuracy (Calibration on any 8 Q-CAD fragments for Angiotensin 1. checked on 4 different masses) <1.0 ppm, m/z range 100 - 1500 (internal)
<1.5 ppm, m/z range 100 – 1500 (external)
*Spec is based on average of errors

Resolving Power @ m/z 400 (lincomycin) > 1,000,000 (FWHH)

Sensitivity (Ubiquitin) S/N > 10:1 for <100 amol (consumed)

High mass (BSA 0,1mg/ml) S/N > 10:1

Negative ions Functionality shown on Fibrinopeptide B

ESI Source Options

Electrospray (standard) (1 µl/min – 1 ml/min)

On-Line Nanospray (standard) (100 nl/min – 500 nl/min)

EZ Nanospray (standard) Zero – Adjust Off line Nanospray

APCI (optional) Optional accessory

APPI (optional) Optional accessory

Off-Line Nanospray (optional) Optional accessory

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

For research use only.
Not for use in diagnostic
procedures.


www.bruker.com



microOTOF-Q III

- The bench-mark in accurate mass LC-MS/MS

The Advantage of Confidence in Routine



The micrOTOF-Q range is widely acknowledged in setting standards in performance and reliability by which all other accurate mass ESI mass spectrometers are judged. micrOTOF-Q III demonstrates the best performance standard in its class:

- 20,000 full sensitivity resolution
- Low picogram sensitivity
- 2 ppm mass accuracy

Each performance parameter whether it is mass accuracy, resolution, or sensitivity is top in its class. Uniquely in the market, micrOTOF-Q III makes no compromise in delivering the best - all parameters are simultaneously delivered for fullest possible insight into your sample.

Complete applications solution software allows your micrOTOF-Q III to become your dedicated partner whatever your challenge in routine formulae confirmation, advanced screening and identification, or intact proteins and biopharmaceutical analysis.



● Your Partner in Innovative Chemistry


Chemical formula generation

Compass OpenAccess™ provides an automated walk-up LC-MS system for chemical formula generation, molecular formula confirmation and generic LC-MS measurements. This client-server based software supports LC-MS workflows especially for chemists in laboratories with various levels of instrumental analysis experience.


GC-APCI

The GC-APCI source enables the combination of LC-QTOF MS with both LC and GC. The micrOTOF-Q III system provides exact mass accuracy and resolving power to expand the horizon for GC-MS based analyses. Data acquisition rates of 20Hz and faster are mandatory for typical GC peak widths of < 2 s (FWHM). The unambiguous formula ID from the GC-TOF-MS run is determined using SmartFormula™.

Compass OpenAccess job results




Compound ID of TNT



Negative ion APCI spectra of the explosive TNT measured on a DirectProbe APCI-TOF system. Here, the accurate mass values allow the direct verification of the compound and its fragment ions.

APCI II ion source with DirectProbe assembly (DIP)



The DirectProbe assembly (left) is an add-on to the Bruker APCI II ion source. Sample preparation involves simply dipping the disposable glass capillary (green) into the solid or liquid sample and sliding it into the APCI II source where vaporization and ionization takes place.

SmartFormula determination

Three dimensions of information simultaneously raise your analytical tasks to unrivaled heights of confidence:

- Measure with unequalled accurate mass
- Validate with True Isotopic Pattern (TIP) analysis
- Also benefit from accurate mass and TIP in analysis of fragments in MS/MS mode



Mass accuracy, chemical knowledge and SmartFormula 3D™ clearly limit the number of possible formulae in molecular formula generation: for confident determination of the elemental composition of a given peak.

This valuable sub-ppm confidence is available for formula determination in pharmaceutical impurity analysis, metabolite identification, pesticide screening and toxicology & doping analysis.

Advanced multi-target screening


Confident quantitative and qualitative multi-target screening for forensics, doping control and residue analysis. Due to a full-scan accurate-mass approach in combination with application specific high-quality screening libraries thousands of compounds are identified and confirmed by SmartFormula3D from 1 single LC-ESI-TOF run.

Retrospective *in-silico* screening for new or unexpected compounds is possible because, unlike in triple-quad based MRM methods, the full molecular information content is retained.




High resolution chromatogram (hrXIC) identifying Pyrrolizidin-Alkaloides 200 ppb each in honey.

Top-Down Analysis of Intact Proteins and Antibodies



Ribonuclease B. A spectrum of the intact protein acquired with the micrOTOF-Q III. The mass difference to the calculated mass is only 1.3 ppm. A resolution of $> 20,000$ FWHM is achieved.




Automated analysis of monoclonal antibodies

Result view of the intact mass workflow for identity testing of discovery mAbs.

The identity of 54 different antibodies coming out of discovery are confirmed here.

Color coding defined by user. Shown here according to average mass accuracy of 50 ppm (left) and an intensity of 3000 counts in the deconvoluted spectrum. For gray, the reference mass was not yet available.



Raw data are available from the Compass OpenAccess server immediately after the analysis, and are accessible to all with a password to the server. In addition, the quantitative data can be exported to Excel.

Data with courtesy of Symphogen A/S, Lyngby, Denmark

Technical Specifications

Best in its class performance

- Proven micrOTOF-Q technology
- Class-leading combination of mass accuracy, resolution and sensitivity without compromise
- SmartFormula3D, the unique combination of accurate mass and true isotopic pattern of both, parent and fragment ions
- Wide dynamic range for ultra-stable accurate mass
- High-transmission dual ion funnel Q-q-front end
- Dimensions 640 x 949 x 1320 mm, weight 160 kg

Source options

- APCI atmospheric pressure chemical ionization source
- APPI atmospheric pressure photo ionization source
- Direct probe option
- Direct GC coupling
- CryoSpray Source
- CaptiveSpray NanoElectrospray source
- CE/MS coupling with grounded ESI needle

Analytical performance

- Mass range 20 – 40,000 m/z
- Mass accuracy 1 - 2 ppm RMS Error
- Mass resolution 20,000 (FWHM) at LC-speed
- Advanced temperature compensation
- Up to 40 Hz Acquisition rate (2GSample/sec sampling rate)

Compass & application software suites

- Integrated LC-MS/MS control and data processing
- Compass OpenAccess: Walk-up LC-MS chemical formula generation
- MetaboliteTools™ – metabolite and impurity identification
- TargetAnalysis™ – multi-target compound screening
- ProfileAnalysis™ – LC-MS based profiling and label-free quantitation
- ProteinScape™ – the bioinformatics platform

Support of – HPLC and sample inlet systems from the following vendors:

Bruker nano-Advance LC, Advion TriVersa, NanoMate, Agilent, Dionex, Shimadzu, VWR/Hitachi, Waters (incl. UPLC), Autosamplers from CTC

For research use only. Not for use in diagnostic procedures.



● Bruker Daltonik GmbH


Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com



solarix XR

- A New Era in Mass Spectrometry



Bruker high resolution FTMS technology enables scientists to see what they have never been able to see before. This technology addresses the needs of many markets (i.e. pharma, petroleum, food safety, etc.) who all want an efficient solution for their analytical problems. The value is in the fine structural elucidation for customer's compounds of interest. This structural info is obtained in the most efficient and cost-effective way using cutting edge technology featuring the enhanced, redesigned Paracell as the key technology that enables **eXtreme Resolution**.

eXtreme Resolution is the ability of the **solariX** to provide "razor thin" peaks in the mass spectrum resulting in significantly greater information content and peak capacity. **eXtreme Resolution** enables interrogation of complex mixtures or compounds very

close in mass without the need for spatial separation providing simpler and more efficient analytical workflows. This is achieved with a combination of technological breakthroughs introduced in the **solariX XR**.



Key Benefits

- ▶ **Analytical power**—unmatched in commercial mass spectrometry provides the capability to create new workflows and explore the chemical landscape in ways not conceivable before.
- ▶ **Flexibility**—combining one of the widest array of sources available with an armada of both traditional and unique dissociation methods providing experimental flexibility that supercharges every application and allowing you to accelerate your workflows.
- ▶ **Speed**—high performance and flexibility translates to faster, streamlined workflows reengineering traditional approaches with tuned methods that save time and money and produce richer datasets in a given time than ever seen before.
- ▶ **Turnkey operation**—advanced software for acquisition, processing and automation combined with a robust source design and fully automated transfer optics provide compelling results with limited user effort.

Key Applications


- ▶ **Small molecule analysis**—elemental compositions are exact and reliable as the **solariX XR** is the first commercial system to routinely provide answers derived from evaluation of fine isotopic structure invisible to most other mass analyzers.
- ▶ **Advanced protein analysis**—**solariX XR** can measure large intact biomolecules with isotopic resolution followed by detailed structural analysis with proven applications in proteomics, biopharmaceutical analysis, and protein science.
- ▶ **Molecular imaging**—leveraging Bruker's unmatched imaging expertise with the power of extreme resolution provides complete competence for spatial localization of small molecules from a variety of samples.
- ▶ **Complex mixtures**—utilize **eXtreme Resolution** to provide selectivity for samples such as petroleum, foods and beverages, environmental, and biological small molecules containing thousands of peaks that cannot be effectively or efficiently separated by conventional chromatographic methods.




eXtreme Resolution for Seeing What Was Missed...

Conventional mass spectrometry only sees nominal mass peaks for the isotopic peak clusters in detected compounds although they are actually comprised of several peaks from heteroatom content. **eXtreme Resolution** allows the user to routinely view this fine isotopic structure to gain powerful insight and unlock the secrets of this previously hidden realm.

$^{13}\text{C}_1$



Automated Fine Structure Interpretation



Automated interpretation of fine isotopic structure

SmartFormula has evolved to handle the increased information provided by **eXtreme Resolution** and allows reading the chemical formula *directly from the mass spectrum*.




3.4.51 ...and Making Quick Work of Complex Mixtures!


Meeting the challenge of analysis for large numbers of samples can quickly outstrip the capability of LC/MS platforms which consume precious time.

solariX XR is optimized for complex samples and can be tasked to handle large sample volumes in a fraction of the time required by traditionally applied methods, conserving time and saving money.

Screening methods can be easily assembled for automated MALDI or ESI that gather complex datasets in less than 1 minute followed by powerful multidimensional statistical analysis to find even the best hidden needle in a haystack.



Detect, Identify, Locate



The **eXtreme Resolution** advantage of **solariX XR** is enhanced by Bruker's industry leading imaging solution to create the ideal environment for high throughput small molecule imaging.


Complex mixtures produced by MALDI imaging are quickly separated in mass space and identified with unmatched specificity allowing seamless workflows that significantly increase chemical information content through spatial localization.

$^{13}C_2$

MALDI MOLECU

Multipole Transfer Optics

Factory-optimized RF ion guides transmit ions 100 – 10,000 m/z based on preloaded methods without the need for specialized tuning.



Powerful qQq Geometry

Enables fast, automated MS/MS and advanced dynamic range MS experiments.



PLANT META

XEN

FOOD TRACEABILITY




LAR IMAGING

UXOMICS

Dual Source Ion Funnel

Orthogonal atmospheric pressure ionization geometry creates a robust, simple yet powerful source with <5 seconds switchover to MALDI.



Forging Productivity from Innovation




Professor Eugene Nikolaev,
ParaCell Inventor,
Russian Academy of Sciences, Moscow

Continuing the tradition of Bruker innovation, the **ParaCell** is a new enabling technology for **solarix XR**. This radical concept is a departure from traditional ICR cell strategies and provides uncommon broadband ion stability resulting in resolution orders of magnitude above other detection schemes.


This power enables the user to effortlessly obtain the extreme resolving power needed to probe isotopic fine structure or highly complex mixtures.




> 10 Million Resolving Power at 7T



ADH Tetramer > 400,000 RP m/z 5460



~700,000 Resolving Power at 1 Hz Duty Cycle



Whether it is increasing the duty cycle or providing additional resolving power, Absorption Mode Processing (AMP) comes with no cost in processing/acquisition time and can provide resolving power in excess of 650,000 at m/z 400 in a 1 second scan.

Simple, Efficient, and Dramatically Robust

The dual stage ion funnel innovates through simplicity. Reducing the number of tunable DC elements means **easy operation** by switching samples and polarity without the need for additional tuning.

The **robust** orthogonal design also eliminates the need for frequent cleaning as contaminants are directed away from the ion optics and removed from the system. Less tuning and **more uptime**

means more concentrating on what matters most.

MALDI has never been easier as the target sits coincidentally in the beam path with the API source, “switching” sources is a thing of the past providing maximum throughput for methods that involve both MALDI and ESI.




Fast polarity switching

Fastest in the industry for high mass accuracy

“Zero Delay Alternating Polarity” adds to the overall efficiency by providing a 4 Hz polarity switching capability at mass accuracies typical of ultra high performance mass spectrometers, providing the only solution to ultra high resolution and high mass accuracy polarity switching.



4 Hz Fast Polarity Switching With Accurate Mass



Flexible Biomolecule Analysis

The new solariX XR is the most **flexible** mass spectrometry platform for the characterization of biomolecules. The solariX XR platform offers:


- ▶ The **largest variety** of collisional- and electron-dissociation techniques available on any mass spectrometry platform.
- ▶ **Highest mass accuracy** and **resolving power** of all mass spectrometers.
- ▶ **Wide m/z range**, from m/z 100 – 10,000.
- ▶ Wide variety of **ion sources**, including MALDI, ESI, nanoESI, APPI, and APCI.



Automated data processing

Automated processing rapidly characterizes the complex, often overlapping dissociation products.

- ▶ Monoisotopic m/z and charge deconvolution.
- ▶ Quickly **combine data sets** from multiple dissociation techniques into an accurate picture of biomolecule sequence, structure, and modifications.




Dynamic Source Configuration

In addition to the included MALDI and ESI sources, **solariX XR** supports a wide range of source options from Bruker and third party vendors, all switchable within seconds.



New easy-to-use software to complete the newest innovations



ftmsControl Features:

- ▶ Absorption mode processing
- ▶ Accumulation during detection
- ▶ Online data reduction
- ▶ Online calibration
- ▶ Online readback traces
- ▶ +/- mode calibration
- ▶ User customizable workspaces
- ▶ Integrated system diagnostics

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH


Bremen · Germany
 Phone +49 (0)421-2205-0
 Fax +49 (0)421-2205-103
 sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
 Phone +1 (978) 663-3660
 Fax +1 (978) 667-5993
 ms-sales@bdal.com

Fremont, CA · USA
 Phone +1 (510) 683-4300
 Fax +1 (510) 490-6586
 ms-sales@bdal.com



solarix XR

- A New Era in Mass Spectrometry



Bruker high resolution FTMS technology enables scientists to see what they have never been able to see before. This technology addresses the needs of many markets (i.e. pharma, petroleum, food safety, etc.) who all want an efficient solution for their analytical problems. The value is in the fine structural elucidation for customer's compounds of interest. This structural info is obtained in the most efficient and cost-effective way using cutting edge technology featuring the enhanced, redesigned Paracell as the key technology that enables **eXtreme Resolution**.

eXtreme Resolution is the ability of the **solariX** to provide "razor thin" peaks in the mass spectrum resulting in significantly greater information content and peak capacity. **eXtreme Resolution** enables interrogation of complex mixtures or compounds very

close in mass without the need for spatial separation providing simpler and more efficient analytical workflows. This is achieved with a combination of technological breakthroughs introduced in the **solariX XR**.



Key Benefits

- ▶ **Analytical power**—unmatched in commercial mass spectrometry provides the capability to create new workflows and explore the chemical landscape in ways not conceivable before.
- ▶ **Flexibility**—combining one of the widest array of sources available with an armada of both traditional and unique dissociation methods providing experimental flexibility that supercharges every application and allowing you to accelerate your workflows.
- ▶ **Speed**—high performance and flexibility translates to faster, streamlined workflows reengineering traditional approaches with tuned methods that save time and money and produce richer datasets in a given time than ever seen before.
- ▶ **Turnkey operation**—advanced software for acquisition, processing and automation combined with a robust source design and fully automated transfer optics provide compelling results with limited user effort.

Key Applications



- ▶ **Small molecule analysis**—elemental compositions are exact and reliable as the **solariX XR** is the first commercial system to routinely provide answers derived from evaluation of fine isotopic structure invisible to most other mass analyzers.
- ▶ **Advanced protein analysis**—**solariX XR** can measure large intact biomolecules with isotopic resolution followed by detailed structural analysis with proven applications in proteomics, biopharmaceutical analysis, and protein science.
- ▶ **Molecular imaging**—leveraging Bruker's unmatched imaging expertise with the power of extreme resolution provides complete competence for spatial localization of small molecules from a variety of samples.
- ▶ **Complex mixtures**—utilize **eXtreme Resolution** to provide selectivity for samples such as petroleum, foods and beverages, environmental, and biological small molecules containing thousands of peaks that cannot be effectively or efficiently separated by conventional chromatographic methods.



eXtreme Resolution for Seeing What Was Missed...


Conventional mass spectrometry only sees nominal mass peaks for the isotopic peak clusters in detected compounds although they are actually comprised of several peaks from heteroatom content. **eXtreme Resolution** allows the user to routinely view this fine isotopic structure to gain powerful insight and unlock the secrets of this previously hidden realm.

$^{13}\text{C}_1$



Automated interpretation of fine isotopic structure

SmartFormula has evolved to handle the increased information provided by **eXtreme Resolution** and allows reading the chemical formula *directly from the mass spectrum*.




3.4.51 ...and Making Quick Work of Complex Mixtures!


Meeting the challenge of analysis for large numbers of samples can quickly outstrip the capability of LC/MS platforms which consume precious time.

solariX XR is optimized for complex samples and can be tasked to handle large sample volumes in a fraction of the time required by traditionally applied methods, conserving time and saving money.

Screening methods can be easily assembled for automated MALDI or ESI that gather complex datasets in less than 1 minute followed by powerful multidimensional statistical analysis to find even the best hidden needle in a haystack.



Detect, Identify, Locate



The **eXtreme Resolution** advantage of **solariX XR** is enhanced by Bruker's industry leading imaging solution to create the ideal environment for high throughput small molecule imaging.


Complex mixtures produced by MALDI imaging are quickly separated in mass space and identified with unmatched specificity allowing seamless workflows that significantly increase chemical information content through spatial localization.

¹³C₂

MALDI MOLECU

Multipole Transfer Optics

Factory-optimized RF ion guides transmit ions 100 – 10,000 m/z based on preloaded methods without the need for specialized tuning.



Powerful qQq Geometry


Enables fast, automated MS/MS and advanced dynamic range MS experiments.



PLANT META

XEN

FOOD TRACEABILITY




LAR IMAGING

UXOMICS

Dual Source Ion Funnel

Orthogonal atmospheric pressure ionization geometry creates a robust, simple yet powerful source with <5 seconds switchover to MALDI.



Forging Productivity from Innovation




Professor Eugene Nikolaev,
ParaCell Inventor,
Russian Academy of Sciences, Moscow

Continuing the tradition of Bruker innovation, the **ParaCell** is a new enabling technology for **solarix XR**. This radical concept is a departure from traditional ICR cell strategies and provides uncommon broadband ion stability resulting in resolution orders of magnitude above other detection schemes.


This power enables the user to effortlessly obtain the extreme resolving power needed to probe isotopic fine structure or highly complex mixtures.




> 10 Million Resolving Power at 7T



ADH Tetramer > 400,000 RP m/z 5460



~700,000 Resolving Power at 1 Hz Duty Cycle



Whether it is increasing the duty cycle or providing additional resolving power, Absorption Mode Processing (AMP) comes with no cost in processing/acquisition time and can provide resolving power in excess of 650,000 at m/z 400 in a 1 second scan.

Simple, Efficient, and Dramatically Robust

The dual stage ion funnel innovates through simplicity. Reducing the number of tunable DC elements means **easy operation** by switching samples and polarity without the need for additional tuning.

The **robust** orthogonal design also eliminates the need for frequent cleaning as contaminants are directed away from the ion optics and removed from the system. Less tuning and **more uptime**

means more concentrating on what matters most.

MALDI has never been easier as the target sits coincidentally in the beam path with the API source, “switching” sources is a thing of the past providing maximum throughput for methods that involve both MALDI and ESI.




Fast polarity switching

Fastest in the industry for high mass accuracy

“Zero Delay Alternating Polarity” adds to the overall efficiency by providing a 4 Hz polarity switching capability at mass accuracies typical of ultra high performance mass spectrometers, providing the only solution to ultra high resolution and high mass accuracy polarity switching.



4 Hz Fast Polarity Switching With Accurate Mass



Flexible Biomolecule Analysis

The new solariX XR is the most **flexible** mass spectrometry platform for the characterization of biomolecules. The solariX XR platform offers:


- ▶ The **largest variety** of collisional- and electron-dissociation techniques available on any mass spectrometry platform.
- ▶ **Highest mass accuracy** and **resolving power** of all mass spectrometers.
- ▶ **Wide m/z range**, from m/z 100 – 10,000.
- ▶ Wide variety of **ion sources**, including MALDI, ESI, nanoESI, APPI, and APCI.



Automated data processing


Automated processing rapidly characterizes the complex, often overlapping dissociation products.

- ▶ Monoisotopic m/z and charge deconvolution.
- ▶ Quickly **combine data sets** from multiple dissociation techniques into an accurate picture of biomolecule sequence, structure, and modifications.




Dynamic Source Configuration

In addition to the included MALDI and ESI sources, **solariX XR** supports a wide range of source options from Bruker and third party vendors, all switchable within seconds.



New easy-to-use software to complete the newest innovations



ftmsControl Features:

- ▶ Absorption mode processing
- ▶ Accumulation during detection
- ▶ Online data reduction
- ▶ Online calibration
- ▶ Online readback traces
- ▶ +/- mode calibration
- ▶ User customizable workspaces
- ▶ Integrated system diagnostics

For research use only. Not for use in diagnostic procedures.

● **Bruker Daltonik GmbH**

Bremen · Germany
 Phone +49 (0)421-2205-0
 Fax +49 (0)421-2205-103
 sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
 Phone +1 (978) 663-3660
 Fax +1 (978) 667-5993
 ms-sales@bdal.com

Fremont, CA · USA
 Phone +1 (510) 683-4300
 Fax +1 (510) 490-6586
 ms-sales@bdal.com

Coming soon!



trichium Fusobacterium Lodderomyces Anaerococcus Colletotrichum Edwardsiella Comamonas Pigmentiphaga Moellerella
ptomyces Bartonella Hafnia Terrimonas Pseudoxanthomonas Corynebacterium Aerococcus Saccharothrix
ibacter Shimwellia Brachy bacterium Leclercia Providencia Trichosporon Gardnerella Sporothrix Leuconostoc
omonas Ruminococcus Scedosporium Dysgonomonas Staphylococcus Bacillus Balneatrix Solibacillus
rhizobium Acholeplasma Filobasidium Propioniferax Azobacterium Kloeckera Austrohalobium
ncola Enterobacter Sporobolomyces Brevundimonas Campylobacter Thermoactinomyces Halobacterium
iluncus Caulobacter Helcococcus Psychrobacillus Campylobacterium Halobacterium
io Iodobacter Tenacibaculum Listeria Plesiomonas Haloquadratum Samsonia Schizophyllum
ella Paracoccus Aureobasidium Eubacterium Dietzia Sarcobacter Vagococcus
holderia Sodalis Empedobacter Sphingopyxis Lactococcus
ionimicrobium Starkeya Prevotella Histophilus Sphingomonas
na Xenorhabdus Methylobacillus Fusarium Wolinella Bacteroides
ellaniella Borrelia Microbacterium Rheinheimera Wautersiella
thophora Riemeirella Chaetomium Atopobium Rhizopus Actinobaculum
ksella Alloicoccus Bacillus Arxiozyma Halobacterium
ilinus Brevibacillus Brachyspira Porphyromonas
amonas Dermatophilus Grimontia
onia Dialister Parabacteroides Clostridium
rivivax Bilophila Alloscardovia Neisseria
obacterium Agrococcus Gracilibacillus
nyloarcula Jonesia Pantoea Elizabethkingia
bacterium Taylorella Delftia Sinorhizobium



● MALDI Biotyper Galaxy

Automated target preparation system for use with the MALDI Biotyper



- **Quality Control**
Using patented light sensor technology the system scans the target positions before and after matrix is deposited, ensuring that each spot is optimally prepared for the MALDI Biotyper workflow.
- **Traceability**
Seamless integration with the MALDI Biotyper server coupled with on-board barcode reading ensures that target plates are matched to their corresponding projects using both the barcode and a date/time stamp.
- **Reproducibility**
A two channel micro volume delivery system not only facilitates support for the extended direct transfer method, but also ensures that the correct volume of matrix is added to each spot thus minimizing the need for further processing of samples.

Completing the workflow ... Pilot your way to the NEW Galaxy!



Multiple Work Benches Supported

Each work bench can access the server using the Satellite software and create projects.

Completely Barcode Driven Workflow

Target barcodes are automatically read by the system and the associated project downloaded from the MALDI Biotyper server.

Traceability Always Assured

Just like the Galaxy, the prepared MALDI target barcode is automatically read by the MALDI Biotyper system and the matching project is selected and processed.



For further information and availability in your country please contact your local Bruker office.

MALDI Biotyper Galaxy: For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

● 450-GC

Specification Sheet

Dimensions and Weights

Size*: Height: 53 cm (20.9 in.),

Width: 66 cm (26.0 in.),

Depth: 56 cm (22.0 in.)

Weight*: 43 kg (95 lb)

* Typical values

Environmental Conditions

Operating temperatures: 10 °C to 40 °C.

Operating humidity (relative): 5 % to 95 %

Line voltage requirements: 101 V, 120 V, 230 V
(±10 % nominal)

Column Oven

Dimensions: 28 cm (w) x 20 cm (d) x 28 cm (h)

Temperature range:

- Ambient +4 °C to 450 °C
- With liquid N₂: -100 °C to 450 °C
- With liquid CO₂: -60 °C to 450 °C

Temperature program ramps/holds: 24/25

Maximum temperature ramp rate: 120 °C/min for all voltages

Cool down rate: 400 °C to 50 °C in 4.5 minutes

Temperature set-point resolution: 1 °C

General Specifications

GC Control:

- External events (digital output):
 - 7 standard
 - 8 optional
- Max number of timed events: 15 x 99
- Heated zones: 7 (including column oven)

Methods:

- Maximum stored methods: 50 (max 30 alphanumeric characters)

Logging:

- Run log file (stored with the chromatogram when using Galaxie™ or MS Workstation)
- Error log file

Local Display:

- TFT full color screen
- VGA resolution (640 x 480)
- Size 8.4" (20 cm)



- Local Control:
 - Touch screen
 - Hard keys
- Languages: English, German, French, Spanish, Italian, Portuguese, Cyrillic, Kanji, Chinese, Thai and Korean (Other languages on request)

Local automation:

- Method lines: 25
- Modes:
 - Infinite looping
 - Dual and duplicate injection

System operational qualities:

- High Inertness: sample path UltiMetal treated, optional
- Low level detection assurance: purged valves, optional

Communication

Ethernet: Protocol: TCP/IP

Data rate: 100 Mbps

Control: GC control and method parameters

Analog output (optional):

- Number of channels: 3
- Time programmable steps: 25
- Output (set individual):
 - 0-1 V (default)
 - 0-10 V

Synchronization signals with other devices and data systems:

- Ready in
- Start out

Data Handling and System control:

- GC: Galaxie™ Chromatography Data System (CDS)
- GC/MS: MSWS (see the GC/MS brochure and datasheet for more information)

Certifications

- CSA:
 - C22.2 61010-1
 - UL 61010-1
- IEC: 61010-1
- EMC:
 - 47 CFR part 15
 - ANSI C63.4
 - EN 61326

Injector Options

Maximum injectors: three, operating concurrently

Pneumatics: Electronic Flow Control (EFC), or manual

Injector types:

- 1177 S/SL Split/Splitless injector
- 1079 PTV Programmable Temperature Vaporizing
- 1093 COC Cold On-Column injector
- 1061 Flash injector
- 1041 PWOC Packed/Wide bore On-Column injector

1177 S/SL Split/Splitless Injector

Pressure range: 0-150 psi
 Total flow: 500 mL/min at 10 psi
 1500 mL/min at 10 psi (He)
 Maximum temperature: 450 °C
 Split range: 1-10,000 (column dependent)

Suited for columns:

- Wide bore: (0.53 mm)
- Narrow bore: (0.05 to 0.32 mm)

1093 COC Cold On-Column Injector

Pressure range: 0-150 psi
 Total Flow: 50 mL/min (Type 23 EFC)
 500 mL/min (Type 24 EFC)

Temperature range:

- Ambient +10 °C to 450 °C using air cooling
- -60 °C to 450 °C using liquid CO₂ cooling
- -160 °C to 450 °C using liquid N₂ cooling

Maximum temperature: 450 °C
 Maximum temperature ramp rate: 200 °C/min
 Temperature ramps/holds: 24/25

Suited for columns:

- Wide bore (0.53 mm)
- Narrow bore (0.32 mm)

Sample Preconcentration Trap (SPT)

Trace level analysis of volatiles in gases
 Fully integrated

Temperature range:

- -60 °C to 450 °C using liquid CO₂ cooling
- -185 °C to 450 °C using liquid N₂ cooling

Temperature rate:

- Ballistic for instant release of adsorbed volatiles

Available traps:

- Two lengths
- A wide range of standard packings and custom packings

Quick-Switch Valve Option


Instantly switch between injectors/columns and detectors
 Configurations: automated or manual, factory or field installed

1079 PTV Programmable Temperature

Vaporizing Injector
 Pressure range: 0-150 psi
 Total flow: 500 mL/min at 10 psi

Temperature range:

- Ambient + 10 °C to 450 °C using air cooling
- -160 °C to 450 °C using liquid N₂ cooling
- -60 °C to 450 °C using liquid CO₂ cooling



1177 S/SL Injector

Maximum temperature ramp rate: 200 °C/min
 Temperature ramps/holds: 24/25
 Split range: 1-10,000 (column dependent)

Operational modes:

- Large volume injection
- Temperature ramped splitless
- Cold on-column
- Split and splitless
- ChromatoProbe solid sample introduction optional

Suited for columns:

- Wide bore (0.53 mm)
- Narrow bore (0.05 to 0.32 mm)

Maximum injection volume: 250 µL (LVI mode)



PFPD Detector

1061 Flash Injector

Pressure range: 0-150 psi

Total flow:

- 50 mL/min (Type 23 EFC)
- 500 mL/min (Type 24 EFC)

Maximum temperature: 450 °C

Suited for columns:

- Wide bore (0.53 mm)
- Packed (1/8 " to 1/4 ")

1041 PWOC Packed/Wide-bore On-Column Injector

Pressure range: 0-150 psi

Total flow: 50 mL/min (Type 23 EFC)

500 mL/min (Type 24 EFC)

Maximum temperature: 450 °C

Suited for columns:

- Wide bore (0.53 mm)
- Packed (1/8 " to 1/4 ")

Electronic Flow Control: Injectors (EFC)

Module types: 4 injector-specific modules

Pressure: 0.1 % Full Scale

Flow: 0.5 % Full Scale and 3% Measured Value

Resolution: 0.1 psi or 0.1 mL/min

Detector Options

Maximum detectors: three: operating concurrently

Pneumatics: Electronic Flow Control (DEFC) or manual

Detector types:

- FID Flame Ionization Detector
- TCD Thermal Conductivity Detector
- ECD Electron Capture Detector
- TSD (NPD) Thermionic Specific Detector
- PFPD Pulsed Flame Photometric Detector
- PDHID Pulsed Discharge Helium Ionization Detector
- MS Mass Spectrometry
(see GC/MS brochure and datasheet)

FID Flame Ionization Detector

Maximum temperature: 450 °C

Detectivity: 2 pg °C/sec

Linear dynamic range: 10⁷

Flame tip type: ceramic (patented)

Operational quality:

- Flame-out detection
- Auto re-ignition

TCD Thermal Conductivity Detector

Maximum temperature: 450 °C

Detectivity: 300 pg/mL (Butane)

Linear dynamic range: 10⁶

Operational quality:

- Filament protection
- Automatic bridge balancing

ECD Electron Capture Detector

Maximum temperature: 450 °C

Detectivity: 7 fg/s Lindane

Linear dynamic range: 10⁴

Radioactive source: 63Ni - 15 mCi (555 Mbq)

TSD Thermionic Specific Detector

Maximum temperature: 450 °C

Detectivity:

- N: 100 fg N/sec (Azobenzene)
- P: 100 fg P/sec (Malathion)

Linear dynamic range:

- N: 10⁵
- P: 10⁴

Operational quality: self-aligning bead

PFPD Pulsed Flame Photometric Detector

Photomultiplier tube:

- S/P
- S/P/N

Maximum temperature: 450 °C

Detectivity:

- S: 1 pg S/sec (S/P tube)
- P: 100 fg P/sec (S/P tube)
- N: 20 pg N/sec (S/P/N tube)

Linear dynamic range:

- S: 10³
- P: 10⁴
- N: 10²

Up to 23 elements can be detected

PDHID Pulsed Discharge Helium Ionization Detector

Detectivity: 50 ppb (Methane)
 Linear dynamic range: 10^4 (Methane)
 Operational quality:

- Gold plated connections
- Welded column connections

Detectors (DEFC)

Module types: 6 detector-specific modules
 Accuracy: $\pm 7\%$ set point flow
 Resolution: 0.1 or 1 mL/min

Automation Options

CP-8410 Auto Injector

Sample capacity:

- 10 x 2 mL vials
- 6 x 5 mL vials
- 5 x 10 mL vials

Large solvent wash vial: 2 x 120 mL*
 Dual and duplicate mode
 Internal standard addition

Modes of operation:

- Liquid
- Ambient headspace*
- SPME (Solid Phase MicroExtraction)*
- Sample heating and cooling*

Pre-programmed modes of injection Syringes:

- 1 μ L, 2 μ L, 5 μ L, 10 μ L, 100 μ L, 250 μ L
 for liquid injection
- SPME

CP-8400 AutoSampler

Sample capacity: 100 x 2 mL vials
 Large solvent wash vial: 2 x 120 mL*
 Dual and duplicate mode
 Internal standard addition

Modes of operation:

- Liquid
- Ambient headspace*
- SPME*
- Sample heating and cooling*

Pre-programmed modes of injection

Syringes:

- 1 μ L, 2 μ L, 5 μ L, 10 μ L, 100 μ L, 250 μ L
 for liquid injection
- SPME

* Optional



Combi PAL AutoSampler

Sample trays: two standard and expandable to four

Tray types:

- 98 x 2 mL vials
- 200 x 1 mL vials
- 32 x 10 mL/20 mL vials
- 96-well plates

Dual and duplicate mode

Internal standard addition

Modes of operation:

- Liquid
- Heated headspace*
- SPME*

Sample heating and cooling

Optional modules: additional sample trays, micro-well plate holders, wash station, SPME fiber bake-out station, dilutor, barcode readers, and flowcell

* Optional

For research use only. Not for use in diagnostic procedures.



For laboratories requiring even greater sample throughput or more extensive sample preparation automation options, Bruker offers the CombiPAL system.



The CP-8400. Automatic access to two injection ports allows you to double your throughput. These can be installed in addition to gas or liquid Sample injection valves for optimum flexibility.

For research use only. Not for use in diagnostic procedure.
Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice.
© Bruker Daltonics 06-2010, #CA-270114
www.bruker.com/chemicalanalysis

Bruker Daltonik GmbH
Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.
Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

● SCION 456-GC

Specification Sheet

Dimensions and Weights

Size*: Height: 57 cm (22.5 in.)

Width: 66 cm (26.0 in.)

Depth: 56 cm (22.0 in.)

Weight*: 43 kg (95 lb)

* Typical values

Environmental Conditions

Operating temperatures: 10 °C to 40 °C.

Operating humidity (relative): 5 % to 95 %

Line voltage requirements: 120 V, 230 V (±10 % nominal)

Column Oven

Dimensions: 28 cm (w) x 20 cm (d) x 28 cm (h)

Temperature range:

- Ambient +4 °C to 450 °C
- Liquid N₂: -100 °C to 450 °C
- Liquid CO₂: -60 °C to 450 °C

Temperature program ramps/holds: 24/25

Maximum temperature ramp rate: 150°C/min for all voltages and 180 °C/min with oven insert accessory

Cool down rate: 400 °C to 50 °C in 4.5 minutes

Temperature set-point resolution: 0.1 °C

Ambient temperature reject <0.01°C change in oven for 1°C change in ambient temp

Retention time Repeatability <0.008% or < 0.0008 min, based on Pentadecane under temperature program conditions

Area repeatability < 1% RSD

General Specifications

Up to 9 EFC modules total, injector, detector and auxiliary

Optional backflush

GC Control:

- External events (digital output):
 - 8 standard
 - 8 optional, total 16
- Max number of timed events: 16
- Heated zones:
 - Standard 5
 - 4 optional, 9 total
- Two power outlets 24V (1A max. each)

Methods:

- Maximum stored internal methods: 50 (max. 30 alphanumeric characters)



Temperature Range (°C)	456-GC Rates (°C/min)
50 -70	150
70 - 115	95
115 - 175	70
175 - 300	45
300 - 450	30

Logging:

- Run log file (stored with the chromatogram when using CompassCDS)
- Error log file

Local Display:

- TFT full color screen
- WVGA resolution (800 x 480)
- Size 23 cm (9")

Local Control:

- Touch screen
- Hard keys

Languages:

- English, German, French, Spanish, Italian, Portuguese, Cyrillic, Kanji, Chinese (standard and traditional), Thai, Korean and Dutch.

Local automation:

- Method lines: 25
- Modes:
 - Infinite looping
 - Dual and duplicate injection

Communication

Ethernet: Protocol: TCP/IP

Data rate: 100 Mbps

Control: GC control and method parameters

Analog output (optional):

- Number of channels: 3
- Time programmable steps: 30
- Output software selectable (set individual):
 - 0-1 V (default)
 - 0-10 V

Synchronization signals with other devices and data systems:

- Ready in and out
- Start in and out

Data Handling and System Control:

- CompassCDS Chromatography Data System

Certifications

- CSA:
 - C22.2 61010-1
 - UL 61010-1
- IEC: 61010-1
- EMC:
 - 47 CFR part 15
 - ANSI C63.4
 - EN 61326


Injector Options

Maximum injectors: three, operating concurrently
 Pneumatics: Electronic Flow Control (EFC), or manual

Injector types:

- S/SL Split/Splitless injector*
- PTV Programmable Temperature Vaporizing*
- COC Cold On-Column injector*
- Flash injector
- PWOC Packed/Wide bore On-Column injector

*Including septum purge



S/SL Split/Splitless Injector

Pressure range: 0-150 psi

Total flow:

- 500 mL/min for N₂/Ar
- 1500 mL/min for He/H₂

Maximum temperature: 450 °C

Split range: 1-10,000 (column dependent)

Suited for columns:

- Wide bore: (0.53 mm)
- Narrow bore: (0.05 to 0.32 mm)

COC Cold On-Column Injector

Pressure range: 0-150 psi

Total Flow: 50 mL/min (Type 23 EFC)

Temperature range:

- Ambient +10 °C to 450 °C using air cooling
- -60 °C to 450 °C using liquid CO₂ cooling
- -160 °C to 450 °C using liquid N₂ cooling

Maximum temperature: 450 °C

Maximum temperature ramp rate: 200 °C/min

Temperature ramps/holds: 24/25

Suited for columns:

- Wide bore (0.53 mm)
- Narrow bore (0.32 mm)

PTV Programmable Temperature

Vaporizing Injector

Pressure range: 0-150 psi

Total flow:

- 500 mL/min for N₂/Ar
- 1500 mL/min for He/H₂

Temperature range:

- Ambient + 10 °C to 450 °C using air cooling
- -160 °C to 450 °C using liquid N₂ cooling
- -60 °C to 450 °C using liquid CO₂ cooling

Maximum temperature ramp rate: 200 °C/min

Temperature ramps/holds: 24/25

Split range: 1-10,000 (column dependent)

Operational capabilities:

- Large volume injection
- Temperature ramped splitless
- Cold on-column
- Split and splitless
- ChromatoProbe solid sample introduction optional

Suited for columns:

- Wide bore (0.53 mm)
- Narrow bore (0.05 to 0.32 mm)

Maximum injection volume: 250 μ L (LVI mode)

Flash Injector

Pressure range: 0-150 psi

Total flow:

- 50 mL/min (Type 23 EFC)

Maximum temperature: 450 °C

Suited for columns:

- Wide bore (0.53 mm)
- Packed (1/8" to 1/4")

PWOC Packed/Wide-bore On-Column Injector

Pressure range: 0-150 psi

Total flow:

- 50 mL/min (Type 23 EFC)

Maximum temperature: 450 °C

Suited for columns:

- Wide bore (0.53 mm)
- Packed (1/8" to 1/4")

Electronic Flow Control: Injectors (EFC)

Module types: 4 injector-specific modules

Pressure: 0.1 % Full Scale

Resolution pressure set points is 0.001psi

Flow sensor accuracy 2% of measured or 0.2% of full scale

Flow sensor repeatability 0.5%

Sample Preconcentration Trap (SPT)

Trace level analysis of volatiles in gases

Fully integrated

Temperature range:

- -60 °C to 450 °C using liquid CO₂ cooling
- -185 °C to 450 °C using liquid N₂ cooling

Temperature rate:

- Ballistic for instant release of adsorbed volatiles

Available traps:

- Two lengths
- A wide range of standard packings and custom packings

Quick-Switch Valve Option

Instantly switch between injectors/columns and detectors
Configurations: automated or manual, factory or field installed


Detector Options

Maximum detectors: four: operating concurrently (one of which is a Single or Triple Quad MS)
Pneumatics: Electronic Flow Control (DEFC) or manual

Detector types:

- FID Flame Ionization Detector
- TCD Thermal Conductivity Detector
- ECD Electron Capture Detector
- NPD (TSD) Nitrogen-Phosphorus Detector
- PFPD Pulsed Flame Photometric Detector
- PDHID Pulsed Discharge Helium Ionization Detector
- MS Mass Spectrometry (see GC/MS brochure and datasheet)

Note: Data Acquisition Rate : 600Hz for all detectors, exception is the PFPD



FID Flame Ionization Detector

Maximum temperature: 450 °C
Detectivity: 2 pg C/sec
Linear dynamic range: 10⁷

Flame tip type: ceramic (patented)

Operational quality:

- Flame-out detection
- Auto re-ignition

TCD Thermal Conductivity Detector

Maximum temperature: 450 °C
Detectivity: 300 pg/mL (Butane)
Linear dynamic range: 10⁶

Operational quality:

- Filament protection
- Automatic bridge balancing

ECD Electron Capture Detector

Maximum temperature: 450 °C
Detectivity: 7 fg/sec Lindane
Linear dynamic range: 10⁴
Radioactive source: 63Ni - 15 mCi (555 Mbq)

NPD (TSD) Nitrogen-Phosphorus Detector

Maximum temperature: 450 °C
Detectivity:
N: 100 fg N/sec (Azobenzene)
P: 100 fg P/sec (Malathion)

Linear dynamic range:

- N: 10⁵
- P: 10⁴

Operational quality: self-aligning bead

PPFD Pulsed Flame Photometric Detector

Photomultiplier tube:

- S/P
- S/P/N

Maximum temperature: 450 °C

Detectivity:

- S: 1 pg S/sec (S/P tube)
- P: 100 fg P/sec (S/P tube)
- N: 20 pg N/sec (S/P/N tube)

Linear dynamic range:

- S: 10^3
- P: 10^4
- N: 10^2

Up to 23 elements can be detected

PDHID Pulsed Discharge Helium Ionization Detector

Detectivity: 50 ppb (Methane)

Linear dynamic range: 10^4 (Methane)

Operational quality:

- Gold plated connections
- Welded column connections

Detectors (DEFC)

Module types: 6 detector-specific modules

Accuracy: $\pm 7\%$ set point flow

Resolution: 0.1 or 1 mL/min

Automation Options

CP-8410 Auto Injector

Sample capacity:

- 10 x 2 mL vials
- 6 x 5 mL vials
- 5 x 10 mL vials

Large solvent wash vial: 2 x 120 mL*

Dual and duplicate mode

Internal standard addition

Modes of operation:

- Liquid
- Ambient headspace*
- SPME (Solid Phase MicroExtraction)*
- Sample heating and cooling*

Pre-programmed modes of injection Syringes:

- 1 μ L, 2 μ L, 5 μ L, 10 μ L, 100 μ L, 250 μ L for liquid injection
- SPME

CP-8400 AutoSampler

Sample capacity: 100 x 2 mL vials
Large solvent wash vial: 2 x 120 mL*
Dual and duplicate mode
Internal standard addition

Modes of operation:

- Liquid
- Ambient headspace*
- SPME*
- Sample heating and cooling*

Pre-programmed modes of injection syringes:

- 1 μ L, 2 μ L, 5 μ L, 10 μ L, 100 μ L, 250 μ L for liquid injection
- SPME

* Optional



PAL Combi-xt AutoSampler

Sample trays: two standard and expandable to four

Tray types:

- 98 x 2 mL vials
- 200 x 1 mL vials
- 32 x 10 mL/20 mL vials
- 96-well plates

Dual and duplicate mode
Internal standard addition

Modes of operation:

- Liquid
- Heated headspace*
- SPME*
- ITEX*


Sample heating and cooling

Additional optional modules: further sample trays, micro-well plate holders, wash station, SPME fiber bake-out station, dilutor, barcode readers, and flowcell


* Optional

SHS-40 Headspace Sampler

Fixed Sample Volume Configuration



Press and Inject Configuration



System Features

Sample Capacity	Up to 40 vials; 20mL or 10 mL w/adaptor Crimp or Screw Top Vials
Injection Modes	Fixed Sample Volume Press and Inject
Injection Valve	6 Way Electric Actuated 1 mL Sample Loop Heat up to 350°C
Incubator	Up to 12 Samples 40-200°C in 1°C increments Integrated Shaker
Instrument Control	Remote with Compass CDS Stand alone
Sample Recognition and Detection	Automated
Carrier Gas Control	Direct from GC
Optional Accessories	Additional Sampling Valve Reactant Gas Catalytic Converter

- Sample Capacity: 40 x 20ml (10ml vials with adaptors).
- Fully integrated
- Sample path inert
- Automated vial queuing system
- Automated Gas Sampling Valve- Electric actuated 6-way injection valve (VICI) with 1 ml sample loop. This entire module may be separately heated up to 350°C.
- Heated incubation oven for 12 samples (40°C - 200°C, in 1°C-steps)
- Integrated Shaker
- Automatic sample detection/recognition
- Compatible with Crimp and Screw Top Vials
- Carrier gas controlled direct from the GC
- AUX gas pre-set internally.
- Up to 9 different parameter sets can be stored.
- Adaptable to different injectors
- Multiple Headspace Extraction mode via single puncture
- Flexible transfer line
- Option: Sample Transfer Line Kit (direct connection to carrier gas, insertion of transfer line via injector no longer required)
- Control: Stand alone or remotely via Compass CDS software
- Voltage: 110V or 230V
- Dimensions: Width: 29cm (12in) Height: 46cm (18in) Depth: 62cm (25in)
- Weight: 21kg (47lb)

For research use only. Not for use in diagnostic procedures.

Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice.

© Bruker Daltonics 04-2012, #CA-270114

www.bruker.com/chemicalanalysis

Description

amaZon SL Dual Funnel Ion Trap Mass Spectrometer System

BDAL #294445

Easy-to-use high performance bench top LCMS system for HPLC detection.

High performance ion trap with dual funnel ion guide.

amaZon SL Ion Trap System

High sensitive small molecule LC/MSⁿ system with fast polarity switching.

The amaZon SL is a robust high performance Electrospray Ion Trap Mass Spectrometer for HPLC detection in routine work.

Compact, small footprint system enclosure for ion trap mass analyzer, ESI source, electronics, and vacuum pumps, containing:

A. Bruker research-grade mass analyzer:

- Advanced multipole High Capacity ion Trap for ultra-fast, high-sensitivity scanning with good resolution
- High sensitivity conversion dynode detector with. For positive and negative ion detection.
- Systems enclosure footprint: width 76 cm x depth 91 cm

B. Robust, computer-controlled APOLLO II electrospray ionization (ESI) source with minimal adjustment and easy maintenance:

- **Apollo II Electrospray Ionization Source with ion funnel**
-
- Highly sensitive ESI Source with proprietary ion funnel guide for gentle mass independent ion focusing and high ion transmission efficiency
- ESI source with grounded needle for safety and easy sample introduction and CE-coupling
- Heated counter current drying gas for gentle and efficient drying
- Pneumatic off-axis nebulization for flow rates up to 1 ml/min., with gradients from 100% aqueous to 100% organic
- SW control of flow and heater of counter current N₂ drying gas
- SW control of high voltage
- High-sensitivity RF ion guide
- Ion lens system including in Source collision induced dissociation possibility (IS-CID)
- Combined Funnel-Octopole-Cartridge with front access for easy maintenance
- Ion lens housing and vacuum system
- Robust setup with easily cleanable glass capillary
- Positive / negative ion operation
- Fast polarity switching for obtaining data of both ion species in a single LC/MS run
- Divert valve option: User exchangeable valve head to suit nano scale or analytical scale applications
- Suitable for HPLC and U-HPLC coupling

C. Vacuum System

- Four-stage pumping system for high ion transmission
- Vacuum system equipped with low maintenance split flow drag pump and single mechanical single stage fore line pump (34 m³/h)
- All gas inlets are fully SW controlled
- Software integrated pressure gauges

D. Research-grade Digital and RF Electronics

- Ultra-stable high-voltage RF generator for high m/z-range
- Low-noise preamplifier for high MS and MSⁿ sensitivity
- Monoisotopic isolation capability in the m/z range 50-2200

High-performance and high-throughput features:**E. Modes of Operation**

- **ultraScan™** Standard Resolution scan from 70-2200 m/z with 32,000 u/sec at 0.6 u FWHM
- **Enhanced Resolution Scan** from 50-2200 m/z with 8,100 u/sec at 0.35 u FWHM for isotopic resolution of up to triply charged ions, optimized for both scan speed and resolution
- **Extended scan range:** 200-4000 m/z with a scan speed of 27,000 u/sec
- **SmartICC™** (ion charge control) unique patented ion charge control for optimal ion trap filling and extended dynamic range for quantification without the need of pre-scans.

Operation modes

- **MSⁿ (n ≤ 11)**
- Selected Reaction Monitoring (**SRM**) and Multiple Reaction Monitoring (**MRM**) for quantitation in complex matrices
- SIM for up to 10 channels
- **Data-dependent Auto-MS/MS** and auto-MSⁿ (n ≤ 5) for alternating MS and MS/MS detection on-the-fly in HPLC runs. MSⁿ detection can be applied for a fixed number of ion signals.

F. Syringe pump

- Low pulsation syringe pump
- Fully integrated in GUI
- User exchangeable syringes
- Flow from 50 nl/min to 1.5 ml/min
- Volume control
- Stall detection
- One 500µL syringe included

G. Data system:

- Workstation: HP z400 or successor
- Windows XP professional pre-installed on image DVD, US English localization
- prepared for installation into LAN
- 3,0 GHz quad core processor, 1333 MHz front side bus
- 4 GB RAM, 1000 GB hard disk
- DVD+/-RW writer 16x DL LightScribe, DVD-ROM 16x
- 2x Ethernet, optical Mouse, Keyboard, Nero 7 Essentials Suite 1
- easy to use emergency backup software Acronis.
- Monitor: HP LP2205wg 22" TFT wide screen display or successor
- Printer: HP LaserJet P3015p USB+parallel or successor

H. Applications software: Compass 1.3 SR2 for amaZon

- Fully integrated software package **Compass 1.3 SR2** for HPLC and Trap control, data acquisition, post processing, and data analysis.

Consisting of:

- **trap control 7.0:**
 - **SmartSuite™** for automated optimization of all instrument parameters without the need of expert skills:
 - **Smart Parameter Setting (SPS):** Auto-adjusting of acquisition parameters to a target mass.
 - **SmartCal:** Auto-Calibration
 - **SmartRamp:** Automated tuning (ramping) of all parameters for best performance, used for easy method development
 - **SmartFrag:** Ramps the collision energy for most efficient and reproducible MS/MS fragmentation
 - **Scheduled Precursor List:** allows defining component specific AutoMS(n) experiments based on known retention times of sample compounds.
 - **Smart Time Segment Editor (STS):** Enables intuitive setup of an LC/MS/MS method with optimized acquisition parameters based on an initial LC/MS run
 - **Expert mode:** extended control over instrument parameters for interactive system optimization of sophisticated MSⁿ methods
- **Data dependent Scans modes:**
 - **ActiveExclusion™**
 - **PassiveExclusion™**
 - Preferred Mass List: contains masses or m/z ranges which should be preferentially selected for AutoMS(n).
 - Preferred Charged State: selection of defined charge states of precursor ions for most efficient AutoMS(n). Avoids as well the fragmentation of the same precursor in various charge states.
- **HyStar 3.2 SR2:** For integrated control of most popular HPLC systems and automation systems. e.g. Bruker EASY-nLC II, Agilent 1100 and 1200 series (also rapid resolution), Dionex/LC Packings Ultimate plus and UltiMate3000, Waters 2795 and Acquity UPLC, LaChrome Elite (Hitachi), Eksigent nanoLC

- **Data Analysis software DA 4.0 SP 4**, including:
 - Advanced data processing with a high degree of automation
 - New QuantAnalysis™ quantitation package
 - LibrarySearch™ module for search of MS/MS and MSⁿ spectra with advanced matching algorithm
 - Charge deconvolution module
 - Neutral loss scans
 - Survey view for density plots of MS and UV-DAD data
 - Export of peak reports to dBase or MS-Excel
 - Export of spectra and ion current profiles as Windows metafiles to word-processing programs
- One SW License ESI Compass 1.3 (incl. esquire control, Data Analysis, LibrarySearch, QuantAnalysis, Dissect)
- One SW License ESI Compass 1.3 charge deconvolution
- One SW License HyStar 3.2 LC/MS (incl. post processing)

- I. Set of manuals and reference CD-ROMs**
- J. Installation**
- K. Familiarization upon installation**
- L. 1 year warranty**
- M. Voucher for a factory-training course - valid for 2 participants.**



Performance Specifications

amaZon SL Dual Funnel Ion Trap Mass Spectrometer System

BDAL #294445

Size	Benchtop: 89 x 76 cm x 51 cm (Height)
Weight	~ 86 kg
Vacuum System	4 stages, 34 m ³ /h rough pump
Apollo II ion funnel electrospray source	Flow rate: 1 µL/min – 1 mL/min
Polarity switching	Fast Polarity switching < 80 ms
Stages of MS(n)	MS(n) for n = 1 through 11
Stages of AutoMS(n)	AutoMS(n) for n = 1 through 5
Mass accuracy in MS and MS/MS	+/- 0.15 u within the calibrated standard mass range at ultraScan resolution in full scan mode, with proper calibration, ICC target and ion statistics, and thermal equilibrium of electronics and ion source
Sensitivity Specification	The following system sensitivity specifications are only applicable when a Agilent 1200 Series HPLC system is purchased together with the amaZon SL System and this LC system is installed in conjunction with the amaZon SL System
Full scan sensitivity in MS	Reserpine 5 pg/µL @ S/N>10: 1 Signal-to-noise ratio of the extracted ion chromatogram of the protonated molecular ion (m/Z 609) as the result of an injection of 1 µL Reserpine (5 pg/µL), measured in positive ion mode, at a flow rate of 200 µL/min when the mass spectrometer is operated in full scan mode using the Enhanced Resolution Mode, scanning from m/z 250 to 750.
Full scan sensitivity in MS/MS	Reserpine 125 fg/µL @ S/N>50: 1 Signal-to-noise ratio of the extracted ion chromatogram of the transition of the protonated molecular ion (m/Z 609) to the most abundant product ion as the result of an injection of 2 µL Reserpine (125 fg/µL), measured in positive ion mode, at a flow rate of 200 µL/min when the mass spectrometer is operated in full scan MS/MS mode scanning the product ion spectrum from m/z 250 to 650 using the Enhanced Resolution Mode.

Scan Mode	Mass Range (m/z)	Resolution FWHM (u)	Scan Speed (u/sec)
ultraScan	70-2200	0.6	32,000
Enhanced Resolution Mode	50-2200	0.35	8,100
Extended Mass Range	200-4000	3	27,000

Optional accessory

Bruker EASY-nLC	Split-free nano-flow HPLC system
Switching valve analytical scale	Optional accessory
APCI	Optional accessory
CE/MS interface	With grounded needle for easy CE-MS set-up (Optional accessory)

• aurora M90

Specification Sheet

Design Overview

The aurora M90 systems is compact, floor-mounted inductively coupled plasma mass spectrometers (ICP-MS) with full PC control of all instrument settings and compatible accessories. It features patented [1] 90 degree reflecting ion optics system for gigahertz sensitivity (1000 Mc/s/mg/L) and low background and interferences. The aurora M90 system includes a sample introduction system, solid state 27 MHz RF generator, and patented [2] Turner Interlaced Coils. Full PC control of plasma positioning, triple stage vacuum system, all plasma gas flows, mass analyzer, and Discrete Dynode Electron Multiplier (DDEM) detector is also included. The vacuum system is fully contained within the instrument for a small instrument footprint. Unique DDEM detector provides nine decades of dynamic range in an all-digital pulse design. Fully web-integrated ICP-MS Expert software uses Bruker's worksheet concept for ease of use and rapid operator training. The aurora M90 system also features a unique and patented [3] Collision Reaction Interface (CRI) providing fast, flexible, interference-free analysis using simple collision and reaction gases.

Instrument Hardware

Sample introduction - ESIA

External Sample Introduction Assembly (ESIA) with Peltier-cooled spraychamber, nebulizer, and peristaltic pump mounted outside the torch box to eliminate temperature-induced drift. Features rapid start-up, easy access for maintenance and cleaning, and rapid switchover to accessories.

- Low-flow (400 μ L/min) glass concentric nebulizer. Optional inert PFA nebulizer and a range of low-flow, inert microconcentric nebulizers
- Peltier-cooled spraychamber with variable temperature control for enhanced stability and reduced oxide ion interferences. PC-controlled temperature from -15 $^{\circ}$ C to room temperature, insulated with inert polypropylene foam
- Fully PC-controlled peristaltic pump, variable speed from 0–50 rpm, three independently pressure adjustable channels for sample, drain, and internal standards or diluent
- Standard one-piece, low flow, ball-and-socket connection torch. Optional semi-demountable torch and inert PFA transfer tubes available

Gas control

- One button automated optimization of all gas flows, including CRI gases
- Standard sheath gas flow allows aerosol dilution of high matrix samples
- Optional external MFC gas control (nitrox 500), for online addition of nitrogen and oxygen to the plasma

RF generator

- 27.12 MHz solid-state, air-cooled, crystal locked RF generator in main instrument housing, 600–1600 W in 10 W increments. Optimum power settings defined and stored within each method for different sample types
- Automatic ignition and shutdown. User-customizable ignition sequence for different accessories and plasma types

Plasma

- Optimize plasma parameters according to specified performance targets (sensitivity, interferences, etc.)
- Full PC control of horizontal, vertical and sampling depth (Z position) of plasma for maximum sensitivity and minimum polyatomic interferences
- Spacious plasma compartment simplifies routine maintenance

Cool plasma

- Patented Interlaced Coils minimize polyatomic interferences without the use of a mechanical torch shield

Instrument Hardware continued

Collision Reaction Interface

- Reduces interferences by injecting simple collision and reaction gases (hydrogen or helium) into the plasma as it passes through the orifice of the cones
- All CRI gas flows controlled by mass flow controllers
- Rapid switchover between gas on and gas off, or between different collision and reaction gases

Plasma interface

- Easy access and removal of sampler and skimmer cone from simple threaded mounts
- One set of high performance nickel cones as standard. Optional high performance platinum cones for corrosive acids and solvents
- Orifice diameters: 1.1 mm sampler, 0.5 mm skimmer
- Water-cooled plasma interface for stability including individual and independent cooling of the cones for faster warm-up, improved stability, and faster cool down

Ion optics

- Patented ion mirror reflects analyte ion beam through 90 degrees while photons and neutrals pass to the vacuum system. Ion mirror creates parabolic electrostatic field to focus the analyte ions with optimum efficiency at the quadrupole entrance aperture. Results in gigahertz sensitivity (1000 Mc/s/mg/L)
- Features easy access to extraction lens 1 and 2 for cleaning without breaking the vacuum
- Auto-optimization of all ion optics settings, including ion mirror, based on selected optimization criteria such as signal and interferences
- Quadrupole mass analyzer is off-axis using patented [4] curved stainless steel entrance rods to ensure low background by further eliminating excited neutrals before the quadrupole

Vacuum system

- Pumps consist of two rotary (SD602 and SD302) and two V-301 turbomolecular pumps for efficient pumping and to eliminate excessive pump wear
- Choice of standard rotary pumps for general applications or inert rotary pumps for more corrosive acids and solvents
- Rotary pumps mounted on forward facing pull-out slides for easy inspection of oil levels and easy access and changing of oil
- All vacuum components located in the main housing of the instrument to reduce noise in the laboratory and eliminate the requirement for extra floor or bench space
- Turbomolecular pumps feature maintenance-free ceramic bearings. First turbopump is positioned immediately behind the ion mirror and skimmer cone for maximum pumping efficiency and removal of unwanted neutrals and particles
- Pneumatic vacuum isolation gate between the first and second vacuum stages. Gate automatically closes in the event of a power failure

Quadrupole

- Precision-machined, stainless steel, round rods manufactured to micrometer tolerances and locked into ceramic mounts for a near-perfect hyperbolic field. Stainless steel construction permits determination of Hg without high memory. Patented curved entrance rods provide a double off-axis design and low background signals
- Easy access to mass analyzer and detector for cleaning or detector replacement. All voltages are fully interlocked and under PC control
- Solid-state air-cooled power supply
- Built-in, on board multi-channel scaler provides up to 40 channels per mass
- Range of 3 to 256 amu with 'zero blast' protection. Resolution adjustable from 0.5 to 1.2 amu
- Mass calibration stability: 0.05 amu per day
- Quadrupole RF frequency: 3.0 MHz
- Scan speed: 2000 amu/s
- Minimum dwell time: 200 μ s

Detector

- All-digital ETP AF250 Discrete Dynode Electron Multiplier (DDEM) as standard provides nine decades of linear dynamic range in an all-pulse-counting detector. No complex, time-consuming analog-to-digital cross calibrations
- Measuring dynodes mounted off-axis for reduced background

Performance

The performance data is typical unless otherwise noted.

High sensitivity mode, Mc/s/mg/L)	⁹ Be	>50
	¹¹⁵ In	>1000
	²³² Th	>500
Precision (10 replicates, 20 minutes)		<3%
Oxide ions	CeO ⁺ /Ce ⁺	<3%
Doubly charged ions	Ba ⁺⁺ /Ba ⁺	<3%
	Ce ⁺⁺ /Ce ⁺	<2%
Background at 5 amu		<5 c/s
Normal sensitivity mode, Mc/s/mg/L)	⁹ Be	>15
	¹¹⁵ In	>120
	²³² Th	>100
Precision (10 replicates, 20 minutes)		<3%
Long term stability (10 μ g/L multi-element standard aspirated for 4 hours)		<4%
Oxides	CeO ⁺ /Ce ⁺	<2%
Doubly charged ions	Ba ⁺⁺ /Ba ⁺	<3%
	Ce ⁺⁺ /Ce ⁺	<2%
Background (at 5 amu)		<2 c/s
Abundance sensitivity	10 ⁻⁶ low, 10 ⁻⁷ high on ²³ Na	
Isotope ratio precision	<0.1% (¹⁰⁷ Ag/ ¹⁰⁹ Ag)	
CRI interference reduction	1 μ g/L As readback in 1% HCl 1 μ g/L \pm 0.1 μ g/L	

For research use only. Not for use in diagnostic procedures.

Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice.

© Bruker Daltonics 03-2011, #CA-270103

www.bruker.com/chemicalanalysis

Instrument Software

Easy-to-use, web-integrated design. Features wizards that guide users through method and sequence development, and method templates for rapid development of commonly used methods.

Features include

- Quantitate analytes on any possible combination of isotopes
- Fully editable interference correction equations
- A range of internal standard assignment options
- Multiple condition sets allowing different element suites to be determined under different conditions with a single sample measurement, including CRI modes, hot plasma, cool plasma, etc.
- Calibration routines for multi-element external calibration, method of standard additions, and isotope ratios
- Automatic method optimization, including ion optics, plasma and CRI gasses and aerosol dilution
- Automatic monitoring and adjustment of nominated elements/isotopes in real time for optimal rinse out
- High speed Time Resolved Signal (TRS) capability for interfacing to chromatographic (such as HPLC) and other separation techniques
- Seamless LC-ICP-MS integration
- Fully editable sample label list of up to 1000 samples for analysis per worksheet
- Autosampler rack and tube positions can be edited for true random access sampling
- Calibrations can be programmed at a user specified rate either amongst the sample tubes or from centralized calibration tubes (rate driven)
- Sequence options include full control of reporting actions at the end of the run, exporting of results at the end of the run, emailing of results, recalibration and resloping error actions, and saving mass scans during the analysis
- Fully automated instrument initialization (start-up) routine, including instrument stabilization time, plasma X/Y position adjustment, mass calibration, and quadrupole resolution
- Simultaneous real-time graphical display of signal as full mass scan, segments of mass scan, and signal response vs time for multiple isotopes or ratios
- Post-run retrospective data editing
- Wide variety of reporting and exporting options
- Comprehensive set of instrument diagnostics and performance tests
- Comprehensive help system

Accessories

Autosamplers

- Compatible with a wide range of commercially available autosamplers and laboratory racks

Productivity Pack

- Optional four-port Switching Valve System (SVS) immediately rinses the sample introduction system while the next sample is being introduced to the instrument
- Reduces carryover, increases sample throughput, and decreases cost per analysis

Integrated speciation (LC-ICP-MS)

- Seamless LC-ICP-MS integration including real-time display of time resolved chromatographic signals
- Full control of LC injections, pumps, autosamplers and data acquisition

High Sensitivity Pack

- Includes non-CRI cones optimized for maximum mid-to-high mass sensitivity when operating the ICP-MS in high sensitivity mode
- Ideal for research applications requiring ultra-low detection limits on non-interfered isotopes

nitrox 500

- Adds oxygen to the plasma to eliminate carbon deposition on the injector tube and sampler cone when analyzing organic solvents.
- Adds nitrogen to the plasma for better detection of elements with high ionization potential (eg. As, Se)

Clean Room Pack

- Includes o-ring free, inert sample introduction system, platinum-tipped cones, and additional exhaust outlet panel
- Removes all exhaust air from the instrument, maintaining the integrity of the highly sensitive clean room environment
- Allows analysis of corrosive samples and acids, such as hydrofluoric acid

Inert sample introduction kits

- Semiconductor kit includes low contamination and inert PFA nebulizer, spraychamber and transfer tubing, platinum-tipped cones and platinum torch injector.
- Chemical kit includes inert PFA nebulizer, spraychamber and transfer tubing, platinum-tipped cones and sapphire torch injector.

Microconcentric nebulizer

- Compatible with a range of microconcentric nebulizers for sample volume limited applications

Laser ablation

- Fully compatible with a range of laser ablation accessories

ICP-MS Installation Requirements

For details of Bruker ICP-MS installation requirements refer to the pre-installation manual, Bruker publication number 8510206700.

Installation qualification (IQ) and operational qualification (OQ) protocols are available for the Bruker ICP-MS.

Bruker Customer Support Policies

Warranty

12 months, though this may vary by location. One-year warranty on DDEM detector with pro-rata warranty on replacements.

Hardware support period

Eight (8) years from date of last unit manufacture. After this time, parts and supplies will be provided if available.

Software support

Software upgrades to fix nonconformances or safety issues will be issued free of charge. Software upgrades for additional functionality will require a fee. Availability of remote diagnostic support may vary according to location.

Further Details

For further details on the following

- PC configurations
 - Installation/Operational Qualification
 - Accessory specifications and application information
 - Part numbers and other ordering information
- Please consult your Bruker office or supplier, or our Web site at www.bruker.com.

For research use only. Not for use in diagnostic procedures.

- [1] US Patent 6,614,021 B1
- [2] US Patent 5,194,731
- [3] US Patent 7,329,863 B2
- [4] US Patent 6,762,407 B2



Fungi Library

● MALDI Biotyper

Expanding MALDI Biotyper's Libraries with the Addition of Filamentous Fungi



The identification of moulds and multicellular fungi has always been one of the most challenging aspects of clinical microbiology. Whilst MALDI-TOF has over recent years revolutionised the identification of bacteria & yeasts, it has had little impact on the identification of filamentous fungi. This has been mainly due to the effect of culture conditions. Bruker has now developed a cultivation method that ensures a stable physiological status and prevents the germination process and the formation of spores.

Standardised liquid cultivation with constant rotation for library construction


In order to reduce the effects of culture conditions and aid in the production of a uniform mycelium, a liquid based cultivation method has been developed which standardises the physiological status. The method has then been used to create library entries. In essence tubes are inoculated with the fungi and placed on a rotator and incubated overnight or until enough biological material is observed.



Rotator SB2


Sample preparation

- Remove cultivation tubes from the rotator and wait for 10 minutes until filamentous fungi sediment to the bottom of the tube
- Harvest up to 1.5 ml from the sediment and centrifuge for 2 min at full speed (e.g. 13,000 upm)
- Carefully remove the supernatant
- Add 1 ml water to the pellet and vortex for one minute, repeat washing and vortexing twice
- Ethanol extract



Daily Routine Workflow – Analysis Possible Direct from Agar

As in this example if a “front mycelium” is clearly visible and can be harvested then it is possible to sample directly from the agar and using the simple ethanol extraction method, good results can usually be obtained without the need for liquid cultivation for most of the samples. In cases where direct harvesting is difficult then the liquid cultivation method should be used.




Top spectrum is achieved after liquid cultivation of *Aspergillus flavus* and the bottom spectrum is achieved by direct harvesting of *Aspergillus flavus* from agar. Note from the image that the “front mycelium is both clearly visible and can be easily harvested.

International fungi consortium

The goal of this consortium was firstly to test the reliability of the liquid based cultivation and sample preparation and then provide securely identified fungi strains for creation of the library.

Currently contributions have been received from over 20 laboratories across 8 countries with an aim of creating an initial library covering >100 species from approximately 40 different genera.



MALDI Biotyper spectra overlay of various species of *Aspergillus*.

Fungi Library

The filamentous fungi library currently includes more than 110 species from approximately 40 different genera. As with other Bruker MALDI Biotyper libraries we will continue to grow and maintain this library.

Material

- Rotator SB2, Order-No. Y549.1, Carl Roth GmbH & Co. KG
- Rotator dish for Rotator SB2, Order-No. Y552.1, Carl Roth GmbH & Co. KG
- Sabouraud Liquid Broth, Modified, 8 ml, Order-No. 221014, Becton Dickinson (BD)



Prof. Dr. med H. Hof • Mycology Lab Laboratory of Limbach, Heidelberg, Germany

“The identification of multicellular fungi to the species level is one of the most challenging tasks of many microbiological laboratories in medicine, hygiene as well as food industries. In cooperation with Bruker’s dedicated microbiology team we worked as part of an international group of fungi experts on the identification of fungi using the MALDI Biotyper approach. Based on Bruker’s existing development on fungi sample preparation procedure, we contributed, established and validated a reference library of a large panel of the most important fungal strains. Our common efforts during the last years’s have shown that MALDI-TOF based molecular fingerprints of fungi provide a high differentiation power both at species and strain level. The analytical performance of the MALDI Biotyper when used with the fungi library is a major technological breakthrough and practical improvement when compared to more conventional approaches and technologies using microscopy and sequencing methods only.”

Order information: Fungi Library 1.0 – #700281

For research use only. Not for use in diagnostic procedures.

Genera list

Acremonium
Alternaria
Arthrinium
Aspergillus
Aureobasidium
Botrytis
Chaetomium
Chrysosporium
Cladosporium
Cunninghamella
Curvularia
Epicoccum
Epidermophyton
Eurotium
Fenellia
Fusarium
Geomyces
Lecytophthora
Lichtheimia
Microsporium
Monilinia
Mucor
Paecilomyces
Penicillium
Phaeoacremonium
Phialemonium
Phialophora
Phoma
Rhizopus
Scedosporium
Schizophyllum
Scopulariopsis
Scytalidium
Thanatephorus
Trichoderma
Trichophyton
Trichurus

● Bruker Daltonik GmbH


Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

www.bruker.com/maldibiotyper



Mycobacteria Library

● MALDI Biotyper

Generation of a *Mycobacterium* Species Library




Next to the clinically important *Mycobacterium tuberculosis complex (MTC)* Nontuberculous Mycobacteria (NTM) are pathogens especially in immunocompromised and senior people. Species identification is important for diagnosis and optimal treatment. Conventional MALDI-TOF identification of mycobacteria using simple extraction methods or direct smearing approaches yield low quality spectra and are unsafe. Hence an inactivation & preparation method for mycobacteria has been optimised using silica beads, which yields high quality spectra and leads to reproducible identifications at the species level. This approach has been used to generate the library entries.

Optimized method using silica beads

- Biomass of mycobacteria in 75% ethanol, centrifugation
- Washing step, 500 μ l water
- Resuspend pellet in 50 μ l of water, 30 min heating 95C
- Addition of 1,2 ml ice cold ethanol
- Centrifugation, discard supernatant
- Suspend dried pellet in acetonitrile
- Addition of silica beads (0.5 mm)
- Vortex for 1 min
- Addition of 70% formic acid
- Vortex for 10 sec
- Centrifugation
- 1 μ l of supernatant on a MALDI target

As can be seen here using two different sample extraction methods, the only extraction that yielded a good quality mass spectrum is the one using silica beads. Hence an optimized method for sample preparation and inactivation using silica beads has been developed to facilitate the generation of good quality mass spectra of every species for library entries.




One Library – Two Sample Cultivation Methods Supported

As can be seen here, it is possible to analyse samples cultured either in liquid media BACTEC™ MGIT™ (Becton Dickinson Company) or from Löwenstein-Jensen media.

Identification of *Mycobacterium* sp. by MALDI-TOF using just 1 ml of BACTEC™ MGIT™ culture is possible.

Figure 1: Spectra of *Mycobacterium peregrinum* DSM 43271^T cultivated on Löwenstein-Jensen medium (top) and in MGIT™ tube (bottom), using “Mycobacterium bead preparation” method.




International mycobacteria consortium

The goal of this consortium was firstly to develop and refine an extraction procedure that works for both liquid and solid media as well as safely inactivating the mycobacterium. And then to provide securely identified mycobacteria strains for creation of the library.

Currently contributions have been received from over 17 laboratories across 9 countries with an aim of achieving both complete species coverage and strain diversity coverage.

Figure 2: The Dendrogram of 19 different species of mycobacterium cultivated on Lowenstein-Jensen medium (blue) and in MGIT™ tubes (red) demonstrates that the same species always cluster together and are well resolved from other species.



Mycobacteria Library



Prof. Dr. Dag Harmsen • Head of R&D Policlinics for Parodontology, WWU Münster

"MALDI-TOF mass spectrometry has been shown to have major clinical impact by a very fast species identification with a very broad species coverage. However, the technology was not readily available for the very demanding analysis of Mycobacteria. In the past years I have been deeply involved in molecular characterization of Mycobacteria, for example by means of sequencing approaches. I am glad to see that now with the new Mycobacteria library of Bruker's MALDI Biotyper, this further dimension of molecular characterization by a proteomic fingerprint is also available for Mycobacteria."

Mycobacteria library

The first version of the mycobacteria library contains more than 90 species and will continue to grow in both species coverage and number of strains per species.

Mycobacteria Library

M. abscessus ssp abscessus	M. colombiense	M. intermedium	M. pseudoshottsii
M. abscessus ssp bolletii	M. conceptionense	M. intracellulare	M. pulveris
M. agri	M. confluentis	M. kansasii	M. rhodesiae
M. alvei	M. conspicuum	M. kumamotoense	M. saskatchewanense
M. arosiense	M. cosmeticum	M. lacus	M. scrofulaceum
M. arupense	M. diernhoferi	M. lentiflavum	M. senegalense
M. asiaticum	M. elephantis	M. mageritense	M. senuense
M. aurum	M. farcinogenes	M. malmoense	M. seoulense
M. austroafricanum	M. florentinum	M. mantenii	M. septicum
M. avium subsp. avium	M. fortuitum subsp. acetamidolyticum	M. marinum	M. setense
M. avium subsp. paratuberculosis	M. fortuitum subsp. fortuitum	M. monacense	M. shimoidei
M. avium subsp. silvaticum	M. gastri	M. montefiorensense	M. shottsii
M. boenickei	M. gilvum	M. mucogenicum	M. simiae
M. bohemicum	M. goodii	M. neoaurum	M. smegmatis
M. botniense	M. gordonae	M. neworleansense	M. szulgai
M. bovis	M. haemophilum	M. nonchromogenicum	M. thermoresistibile
M. branderi	M. hassiacum	M. novocastrense	M. tokaiense
M. brumae	M. heckeshornense	M. parafortuitum	M. triplex
M. brisbanense	M. heidelbergense	M. parascrofulaceum	M. tuberculosis
M. canariense	M. hiberniae	M. paraseoulense	M. vaccae
M. celatum	M. hodleri	M. parmense	M. wolinskyi
M. chelonae subsp. chelonae	M. houstonense	M. peregrinum	M. xenopi
M. chimaera	M. immunogenum	M. phlei	
M. chitae	M. insubricum	M. phocaicum	
M. chlorophenicum	M. interjectum	M. porcinum	

For research use only. Not for use in diagnostic procedures.

Order information: Mycobacterium Library 1.0 – #700279

● Bruker Daltonik GmbH


Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

www.bruker.com/maldibiotyper



SCION™ SQ Series GC-MS

- The Gas Chromatographers' Detector

Innovation with Integrity

GAS CHROMATOGRAPHY

Introducing the SCION SQ GC-MS

Bruker's long tradition of innovation and product reliability have combined to create a new industry standard for gas chromatography single quadrupole mass detection – the SCION SQ series. By designing the GC-MS systems to exceed the most critical performance and reliability needs of GC users, Bruker has delivered systems that are especially for, and all about, the ultimate success of the GC users. The SCION SQ Select, Prime, and Premium models are designed to meet many important user specific requirements – reliable performance, ease-of-use and simple maintenance – all in a small footprint that saves valuable bench space.



SCION SQ GC-MS Benefits

Easy to Use and Maintain

- Simple tuning due to "lens-free" ion-path design
- No multiplier calibration required

Robust

- An inert ion source that requires less frequent cleaning


Sensitive

- Optional active-focusing q0 uses He atoms to focus ions

● Enhanced Library Searching through Cleaner EI Spectra

The SCION SQ models are designed to analyze thousands of samples from complex matrices. With an upper mass limit of m/z 1200, they are exceptionally capable of handling almost any GC application. The innovative lens-free design, combined with the robustness of the axial ion source, delivers unmatched stability and ultra-high sensitivity on a routine basis.

- One-click search of multiple spectral databases
- User-created spectral libraries and full support of spectral libraries such as NIST, Wiley, and Maurer/Pfleger/Weber (MPW)
- Adjustable spectral search parameters to streamline library searches
- Automated workflow to build a SIM method from a full scan data file
- Flexible and easy to use



Flexible and powerful library search showing matching of sample peak to o,p'-DDE in the NIST library




Why Choose the SCION SQ?

The SCION SQ series delivers exceptional performance for a single quadrupole mass spectrometer: a robust axial ion source, ultra-high sensitivity, cleaner spectra, and virtually-zero neutral noise. The series includes the SQ Select, Prime, and Premium GC-MS models.


For routine EI-only applications, the **SCION SQ Select** delivers the best value. It comes with an industry leading 400 L/s high-capacity turbo pump that enables fast pump-down time for quick maintenance, and the use of high carrier gas flow for fast GC separations.

The EI-only **SCION SQ Prime** model comes with the split-flow, ultra-high capacity turbo pump (300-400 L/s) for added robustness and ease of use. It comes with the innovative active-focusing q0 optic that delivers enhanced sensitivity.


The **SCION SQ Premium** model comes CI-ready and therefore is the most versatile GC-MS SQ platform. The active-focusing q0 optic is heated for the demanding analysis that requires the ultimate robustness.



● Linearity, Sensitivity, and Spectral Matching Using the SCION SQ for US EPA Methods



TIC chromatogram of a 5 mL water sample containing 84 VOCs at 10 µg/L (ppb) by US EPA Method 524.3



(A) Excellent linearity of Tetrachloroethylene from 0.1 to 40 ppb with the purge-and-trap, (B) good sensitivity, (C) high quality spectra, (D) down to 0.1 ppb level and showing good match to the NIST library.

Gas Chromatographs

An Infusion of Innovation with a Legacy of Reliability

The GC is a key part to the reliability, robustness, and sensitivity of any GC-MS analysis. Bruker's philosophy of innovation is highlighted by the introduction of two new GCs built to support the ultra-sensitive SCION SQ. The compact SCION 436-GC and the versatile SCION 456-GC can accommodate two columns in the oven and are available with new backflush technology and the innovative ChromatoProbe™. The new GCs are also equipped with the multi-language touchpad display supporting 13 languages and enabling MS control.

SCION 436-GC

Compact design for those focused on routine applications requiring maximum throughput using one or two injectors


- Select up to 2 injectors: Split/Splitless (SSL), Programmable Temperature Vaporization (PTV), Cold-on-Column (COC), Flash and Packed/Wide bore On-Column (PWOC)
- Support one GC detector and the mass spectrometer
- High precision electronic pressure control
- All temperature zones up to 450 °C
- Automated with Model 8400/8410 or CTC liquid/headspace autosampler

SCION 456-GC

Versatile design with additional injector and detector options for laboratories seeking multipurpose analysis using both GC and GC-MS

- Select up to 3 injectors: SSL, PTV, COC and PWOC
- Add up to 3 GC detectors-FID, ECD, TCD, PFPD, NPD (TSD)
- High precision electronic pressure control
- All temperature zones up to 450 °C
- Automated with Model 8400/8410 or CTC liquid/headspace autosampler

Bruker has now set the industry standard for ease-of-use: Basic operation of SCION GC-MS systems can be directly controlled from the multi-language touchpad on the gas chromatograph. Automatic tuning, along with tune-to-target for meeting specific USEPA methods can be done with a touch of a button. And the MS can be vented and pumped down from the interface for the ultimate in easy maintenance.



SCION 436-GC



SCION 456-GC

● Additions to Enhance System Capability and Performance

SHS-40 Automated Headspace Sampler

- Perfect for analyzing VOCs in solid or liquid samples
- 40/125 sample capacity – Crimp cap or screw cap 10 or 20 mL vials
- 12 position oven for increased throughput
- 200 °C sample heating for extended range
- Injection with 1 mL sample loop, designed for EPC GCs MHE mode via single puncture ensuring no leaks

ATOMX® Purge and Trap (P&T) System from Teledyne Tekmar

- Automated VOC Sample Prep System
- Combine an autosampler and purge and trap concentrator into a single platform
- Unique Automated Methanol Extraction (ME) features for high level soil samples
- 80-position carousel design for optimal throughput
- TekLink™ software with fully optimized user interface including diagnostic tools and benchmark tests for instrument validation

The perfect addition for the SCION SQ PTV inlet is the backflush option. Complex sample matrix can quickly ruin the chromatographic performance of your GC column. However, the PTV with Bruker's "backflush" technique can reliably divert the higher boiling sample matrix away from the column. The benefits of this accessory are many:

- Run more samples per day – decrease analysis times as the heavy components are quickly eliminated
- Save time by eliminating column bakeout
- Preserve column performance for extended period of time

ChromatoProbe™


Added versatility for superior analysis of solids, liquids, and slurries


Compounds such as street drugs, industrial solids, synthetic organic products, and plant tissues that normally are not considered amenable to GC-MS analysis can be easily investigated with the ChromatoProbe.

Samples are introduced into the PTV injector via disposable micro-vials. Non-volatile or thermally degraded components from the sample remain in the micro-vial allowing the system to remain clean.

ChromatoProbe benefits:

- Increase uptime
- Minimize system contamination with disposable micro-vials
- Directly desorb samples in the PTV injector without added hardware





SCION
Focused on Results



www.ScionHasArrived.com



www.Bruker.com


For research use only. Not for use in diagnostic procedures.

 www.bruker.com/scion



DHA Analyzer Family

- Optimized Solutions for Detailed Hydrocarbon Analysis



The DHA Analyzer is a complete high resolution gas chromatography solution for the analysis of hydrocarbons in petroleum streams. It is capable of performing all of the standard methods including the analysis of light petroleum streams and crude oil light end.

Key Benefits include:

■ **Compliant with all industry standard methods**

Be confident using Bruker's DHA Analyzers, which are configured in accordance with all the established standard methods including ASTM D6729, D6730, D6733, D5134, IP 344/DHA "Front End" and "Fast DHA."

■ **Complete and fully integrated solution**

DHA Analyzers come complete with everything you need to be up and running quickly.

■ **Powerful and easy-to-use analyzer**

With relatively little training, operators can generate outstanding analysis results day after day.

■ **Save time**

Easily generate reports with a few mouse clicks and reduce analysis time using "Fast DHA," increasing lab productivity.

■ **Single vendor solution**

Bruker's GC analyzers are built and tested at Bruker's factory, as well as installed and performance-verified on-site by Bruker trained and certified engineers. Rest assured that our analyzers can meet or exceed your needs throughout the instrument's lifetime.



● DHA Analyzer Family

Detailed hydrocarbon analysis is often the preferred technique to fully characterize petroleum streams. The technique is based on the identification of individual components using high performance, high resolution capillary gas chromatography.

Software Ensures Accurate Identification

To successfully apply gas chromatography to detailed hydrocarbon analysis (DHA) the analyzer must be able to correctly identify a large number of components (many eluting very closely to one another) in a complex chromatogram. The identification is based on a comparison of their individual retention index values to those in a pre-established database. Therefore, it is extremely important that the analyzer functions in a highly repeatable manner.

Due to the high number of possible components in hydrocarbon streams, even though analyzed on a high resolution column, it cannot prevent coelution of some hydrocarbons. However, the presence of certain hydrocarbons depends on the type of sample stream. For example

the presence of naphthenes is more common in Naphtha than in Reformate where aromatics are more dominant. The DHA analyzer allows presetting a sample type within the DHA method. By setting a component identification preference for a coelution, the most likely main component is assigned and sample characterization is improved. The software also allows defining additional sample presets in case new sample types become available.

KRI or LRI*	Component	Hydrocarbon Type	Component Present in:		
			Reformate	Naphtha	Alkylate
721.4	1,1,3-Trimethyl cyclopentane	Naphthenic C8	-	+	-
	2,2-Dimethylhexane	iso paraffin C8	++	+	-
751.1	Toluene	Aromatic C7	+++	+	-
	2,3,3-Trimethylpentane	Iso paraffin C8	-	+	++
877.9	o-Xylene	Aromatic C8	+++	+	-
	1,1,2-Trimethylcyclohexane	Naphthenic C8	-	+	-

* KRI = Kovats REtention Index. LRI = Linear Retention Index.

Table 1: For some sample streams a single component assignment is possible due to the assumed absence of one of the coeluting components.




Figure 2: Detailed hydrocarbon analysis of a reformate sample showing aromatics identification according ASTM D6730

● Standard Methods

Selecting Individual Peaks and Updating the Database

The DHA software includes a Peak Select and Database Update function to make identification of unknown peaks as straightforward as possible. The system automatically provides the operator with detailed comparative retention index information for each “unknown” peak including a highlighted “best fit” indicator, making it easy for the operator to determine the ID.




Figure 3: Assigning custom peak matching criteria is easy.




Figure 4: DHA provides an easy-to-use graphical means to select peaks and update the database

Integrated Standard Test Methods

Bruker’s DHA analyzers are compliant with the following methods:

- ASTM D6729
- ASTM D5134
- ASTM D6730
- “Fast” DHA
- ASTM D6733
- IP 344 “Front end”




Figure 5: Choosing a preferred standard method is easy with the DHA software




Figure 6: Choosing report options is simple

Although each DHA analyzer is configured, tested and certified at the factory for a standard method specified by the customer, the DHA software permits the operator to utilize any of the other popular standard methods as well. And, because of the outstanding performance and flexibility of the Bruker GC and CompassCDS software design, Bruker is able to quickly modify the existing methods or add new ones if required as a result of the on-going “dynamic” industry standard processes.

Powerful Reporting is Built-in

Bruker’s DHA software includes several report options to accommodate the standard methods and/or to meet the customer’s special needs. These include:

- Carbon number distribution
- PIONA report; (weight and volume percentage by hydrocarbon group)
- Physical properties calculations; specific gravity and molecular weight
- True distillation profile
- RON/MON specification

● Reduce Sample Analysis Time With “Fast DHA”




Figure 7: These chromatograms illustrate the decreased analysis time using the “Fast DHA” method. Chromatogram of a naphtha sample run on a 40 m X 0.10 mm X 0.2 µm film BR-DHA using the “Fast DHA” method (left). Chromatogram of the same sample, but run on a 100 m X 0.25 mm ID X 0.5 µm film BR-DHA column using standard method D6729 (right). Note reduced analysis time from ~80 minutes to ~20 minutes; almost four-fold.

DHA Analyzer Includes These Key Components

- Bruker GC high performance gas chromatograph equipped with:
 - Split/splitless capillary injection port
 - High performance capillary column (dependent on specified method on order)
 - Flame ionization detector (FID)
 - Full electronic flow control (EFC) of all gases
- State of the art backflush capabilities for the IP 344 “Front End” method
- CP-8400 or CP-8410 automatic liquid sampler
- CompassCDS for system control, data acquisition and report generation
- CompassCDS based DHA application software
- Computer/monitor
- Pre-loaded standard methods
- Factory test
- Reference chromatogram
- Reference standard for use in conducting on-site performance verification





Figure 8: Physical properties and detailed hydrocarbon report

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- "Quick connect" fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time



For research use only. Not for use in diagnostic procedures.



www.ScionHasArrived.com



www.Bruker.com



Natural Gas Analyzer

- A Family of optimized GC Solutions

Innovation with Integrity

Gas Chromatography

Bruker Natural Gas Analyzers

Natural gas is bought and sold as a bulk commodity with price based on its energy content. It is very important for all stakeholders in the natural gas supply and consumer chain to accurately determine the heating value of their streams. Bruker offers a full range of GC based solutions for the analysis of natural gas. The analyzer family is designed to offer superior results through the use of industry proven hardware, software, optimized columns and consumables, and is backed by a team of global sales and support specialists.

Key Benefits include:

- A complete range of natural gas analysis (NGA) solutions. Bruker offers many different NGA gas chromatography analyzers to meet the broadest range of stream sample types and throughput needs, whether the analysis is conducted in a laboratory, at-line or in the field.
- Easy to operate, powerful GC solutions. Bruker's GC with CompassCDS Chromatography Software, form a powerful combination and do not require a high degree of skill to be used successfully.
- Flexibility to analyze natural gas, liquified petroleum gas or natural gas liquids (NGL). Bruker's GC based NGA analyzers can be configured to measure the composition of LPG or NGL streams through the use of specialized sample conditioners, ensuring sample integrity is consistently maintained.
- Operational procedures are fully documented. All Bruker NGA analyzers not only incorporate proven GC hardware and software, but arrive with the pre-loaded analysis method(s) and documentation specific to the application.
- Comprehensive single-vendor solution. Bruker is proud to provide complete solutions. The hardware, software, application optimization, documentation, installation and performance verification are all delivered by Bruker.



Figure 1: The Bruker GC based Natural Gas Analyzer.

● Solutions for Natural Gas Analysis

Gas chromatography offers a proven means to determine the composition and heating value of natural gas and related streams quickly and cost effectively. Bruker's natural gas analyzers (NGAs) are standard 'turnkey' systems pre-configured and tuned at the factory to ensure their compliance with standard methods used to determine the heating value of natural gases and related streams. They can also be specially configured to determine other components of interest (eg. sulfur compounds), to ensure suitability for use in downstream processes. The analyzers are based on the Bruker GC gas chromatograph platforms and Bruker's CompassCDS Chromatography Software. All systems employ a proven and optimized multi-channel/multi-dimensional approach to determine the heating/calorific value of natural gas, as well as quantify individual components. Bruker offers several NGA systems to meet the widest range of analysis needs.

Basic NGA (System A)

This is the simplest of all available natural gas analysis systems. As shown in Figure 3, the system employs a single valve column designed for simplicity, a Thermal Conductivity Detector (TCD) and Flame Ionization Detector (FID). The TCD is used for the determination of O₂, N₂, CH₄, CO₂ and Ethane, while the FID, connected in series, determines hydrocarbons in low concentrations, i.e. C3-C5 and C6+ back-flushed grouping peak (late back-flush). A single unheated 4 port Liquid Sampling Valve (LSV) is available for LPG type samples.




Figure 2: Chromatogram of natural gas sample from 'System A'. Up to 10 Min TCD signal, 10 to 31 Min FID signal.




Figure 3: System configuration schematic diagram for Natural Gas Analyzer A.

Natural Gas Analysis

NGA/Natural Gas Liquids (System B)

This system is optimized for the analysis of natural gas or de-methanized hydrocarbon matrices. For natural gas, the components of interest are typically oxygen, nitrogen, carbon dioxide, methane, ethane, propane, butane, isobutane, pentane, hydrogen sulfide, and C6+ as a composite peak. For de-methanized streams (liquid natural gas) the components of interest are typically carbon dioxide, ethane, propane, butane, iso-pentane, hexane, and C7+ as a composite peak.

The system is configured with a 10 and 12 port valve (a third liquid sampling valve is added if liquid streams are to be analyzed) and three analysis columns connected to TCD and FID detectors (Figure 4). The system simultaneously injects the stream onto two column systems, a Molsieve column for the determination of O₂ and N₂ without the use of coolants, and short/long Non-Polar columns for the analysis

of hydrocarbons and CO₂. The Non-Polar columns are set up for early back-flush, which optimizes sensitivity while reducing run time (from 25 minutes using the system 'A' configuration to less than 15 minutes) (Figure 5). The System „B“ is extendable with two standard options. For the analysis of de-methanized liquefied natural gas distillates (i.e. propane, butanes and pentanes) an automated 4 port Liquid Sampling Valve is available to inject the sample as a Liquid. For the analysis of hydrogen and helium in natural gas a He/H₂ channel is available as extra GC channel.




Figure 4: System configuration schematic for Natural Gas Analyzer 'B'.




Figure 5: Chromatogram of natural gas samples from 'System B'. Note separation of oxygen and nitrogen.

● “Rich” Natural Gas Analysis

NGA/Natural Gas Liquids - Extended Analysis (System C)

This system is specifically designed to analyze ‘rich’ natural gas or natural gas liquid streams by separating all hydrocarbon components up to C16. As with System B, it separates and quantifies oxygen and nitrogen, as well as measuring hydrogen sulfide down to ~100 ppm.

The system is configured with a 14 port valve and 6 port valve. The 14 port valve enables the system to introduce the sample stream simultaneously to three independent columns with automated detector switching which provides high sensitivity detection of all components of interest. The valves are installed in the multi-valve oven for flexible operation of the conventional column oven. Two of the sample paths flow onto Molsieve and porous polymer columns to separate oxygen, nitrogen and carbon dioxide, ethane, methane, ethane and H₂S, and the other via a ‘splitter’ onto a high performance non-polar capillary column to separate the hydrogen components up through C16. The 6 port valve is used to direct the separated components fraction




Figure 7: System schematic for Natural Gas Analyzer ‘C’

to the TCD detector while components remaining on the Molsieve column are flushed to vent. If natural gas liquids are to be analyzed, a third valve (liquid sampling) is added to the configuration described above. For natural gas containing hydrogen and helium a He/H₂ channel is configured additionally.




Figure 6: Chromatogram showing natural gas sample from ‘System C’ FID channel. Note individual separation of C6+ components




Figure 8: Chromatogram showing natural gas sample from ‘system C’ TCD channel.

Base Natural Gas Analysis Systems

Analyzer Characteristic	Bruker GC System 'A'	Bruker GC System 'B'	Bruker GC System 'C'
Gas Components Measured			
O ₂ , N ₂ , CO ₂ , CH ₄	YES	YES	YES
Ethane, Propane, Butane, Iso-butane, Neo-pentane, Pentane, Iso-pentane	YES	YES	YES
C6+ as a composite peak	YES	YES	YES
C7+ as a composite peak	YES	YES	YES
O ₂ /N ₂ Separation	YES ⁽¹⁾	YES	YES
He/H ₂ Separation	NO	YES ⁽³⁾	YES ⁽³⁾
Max hydrocarbon number speciated	C6	C7	C16
LPG	YES ^(4 or 5)	YES ^(4 or 5)	YES ^(4 or 5)
De-methanized natural gas/natural gas liquids	NO	YES ⁽⁴⁾	YES ⁽⁴⁾
Typical analysis to analysis repeatability % RSD	<1.0	<1.0	<1.0
Analysis time	30 min	14 min	15-20 min
Standard Gas Methods			
IP-345	YES		
GPA 2172			
GPA 2177		YES	
GPA 2261		YES	
GPA 2286		YES	YES
GPA 2186		YES	YES
ASTM D5504			YES ⁽⁵⁾
ASTM D6228			YES ⁽⁵⁾
Natural Gas Calculation Methods			
ISO 6974		Contact Bruker, several configuration possibilities	
ISO 6976			
GOST-22667			

⁽¹⁾ Requires liquid nitrogen or liquid CO₂ oven cooling

⁽²⁾ Requires 3rd dedicated channel

⁽³⁾ Requires additional channel including valve, columns and TCD

⁽⁴⁾ Requires LSV to be additionally installed

⁽⁵⁾ This method is specifically for sulfur components in natural gas, therefore a sulfur selective detector must be used such as a PFPD

● Sulfur components in Natural Gas

CompassCDS Chromatography Data System
Extended Natural Gas
Analysis Report

Run File NGC.DATA
Method Nat Gas C
Sample Name NGC

Component	Mole %	MW	kJ/Mole (Superior)	kJ/Mole (Inferior)
3,3-Dimethylpentane	0.9586	0.83	40.14	37.17
trans 1,2-Dimethylcyclopentane	0.6972	0.60	29.22	27.06
2,2-Dimethylhexane	0.6100	0.53	25.50	23.60
Nitrogen	1.4815	0.42	0.00	0.00
Methane	65.3595	10.49	582.41	524.61
Carbon Dioxide	0.5229	0.23	0.00	0.00
Ethane	8.2789	2.49	129.52	118.28
Hydrogen Sulphide	4.3573	1.48	24.50	22.57
Propane	3.7037	1.63	82.23	75.68
i-Butane	2.3529	1.37	67.51	62.31
n-Butane	2.4401	1.42	70.24	64.84
neo-Pentane	1.0458	0.75	36.77	33.99
i-Pentane	2.2658	1.63	79.99	73.98
n-Pentane	2.4401	1.76	86.31	79.84
n-Hexane	1.5558	1.43	69.49	64.36
2,2-Dimethylpentane	1.0458	0.90	49.35	40.47
Methylcyclopentane	0.7843	0.68	32.86	30.43
Totals	100.0000	28.63	1,406.03	1,279.18

MJ/kg (Superior)	49.11	Sample Ideal Relative Density	0.9885
MJ/kg (Inferior)	44.68	Sample Real Relative Density	0.9948
MJ/m ³ (Superior)	59.46	Sample Ideal Absolute Density	1.2069 kg/m ³
MJ/m ³ (Inferior)	54.10	Sample Real Absolute Density	1.2151 kg/m ³
Sample Compressibility	.9932	Sample Wobbe Index	59.62

Base Conditions: Temperature = 15 C, Pressure = 101.325 kPa
Reference: ISO 6978: 1996(E)

Figure 9: Natural Gas Report - 'System C' generated using Bruker's CompassCDS Chromatography Data Software

There are several methods specifically used for the analysis of sulfur components in natural gas, e.g. ASTM D5504 and ASTM D6228. Bruker's natural gas analyzers can be modified to measure sulfur components through the addition of an extra, fully inert channel, dedicated for the determination of lowlevel sulfur components only. This sulfur channel will be equipped with the Pulsed Flame Photometric Detector (PFPD), a sulfur specific detector. The systems can also be configured and tested exclusively for the analysis of sulfur components as per standard methods, or enterprise specific requirements.





Figure 10: Detection of trace level sulfur components using the 450-GC, specially treated with Inert Steel surface deactivation and Pulsed Flame Photometric Detector (PFPD).

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- "Quick connect" fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time



For research use only. Not for use in diagnostic procedures.



www.ScionHasArrived.com



www.Bruker.com



Refinery Gas Analyzer

- Optimized GC Analysis Solutions

Bruker Refinery Gas Analyzers

The source and composition of refinery gases varies considerably. Measuring gas composition precisely and accurately is a significant challenge in today's refinery operations. Bruker's refinery gas analyzers are designed to deliver superior, reliable results for a wide range of sources and analysis throughput requirements.

- A range of refinery gas analysis (RGA) solutions. Bruker offers RGA solutions to meet the broadest range of stream sample types and throughput requirements.
- A powerful, easy to use GC solution. Bruker's 456-GC and CompassCDS chromatography software is a very powerful combination designed to achieve the best possible results. In addition, these systems do not require a high degree of operator skills.
- A highly flexible solution for analysis. The Bruker RGA solutions can optional be configured to analyze high pressurized gas and liquefied petroleum gas (LPG) through the use of a fully integrated Micro-Gasifier, giving the flexibility to accommodate a wide range of stream types.
- Operational procedures are fully documented. Bruker RGA analyzers not only incorporate proven GC hardware and software but also arrive pre-loaded with analysis methods, and include documentation specific to the application required.
- A comprehensive, single vendor solution. Bruker provides complete solutions. The hardware, software, application optimization, documentation, installation and performance verification are all provided by Bruker, offering an all inclusive, convenient analysis solution.



Figure 1: The 456-GC RGA has outstanding flexibility, analytical power and robustness.

● Key Benefits

Bruker solutions for refinery gas analysis

Typical sources for refinery gases include atmospheric or FCC overheads, ethylene, propylene production, fuel gas, stack gas and off gas from desulfurization. The physical stream types range from gas to highly pressurized gas or liquefied gases. Bruker's refinery gas analyzers (RGA) are 'turnkey' systems pre-configured and tuned at the factory to conform to industry standard methods including: UOP 539, DIN-51666 and ASTM D2163. The RGA systems are based on the Bruker 456-GC. To perform good analysis, the RGA is optional equipped with an integrated micro-gasifier. This sample conditioning device ensures complete vaporization of LPGs and high pressures samples to prevent any sample discrimination prior injection.

The Analyzers employ a proven and optimized multi-channel approach. They determine the concentration of individual saturated and unsaturated hydrocarbon components up to and including C5 (C6 and higher components as a composite peak) and all permanent gases, including hydrogen and hydrogen sulfide in a single analysis. Included in every system is Bruker's powerful CompassCDS chromatography software to provide complete analyzer control, data acquisition and flexible report generation.




Figure 2: The separation of light hydrocarbons using the Standard RGA.




Figure 3: The analysis of the permanent gases and hydrogen (and helium) using the Standard RGA.



Figure 4: The RGA analyzers are applicable to a variety of different hydrocarbon streams.

● Bruker Refinery Gas Analyzers

Bruker offers two RGA systems to meet the widest range of analysis requirements:

- Standard RGA:** A three channel 456-GC with a multi-valve design using both capillary and packed columns. The first channel is optimized for the analysis of permanent gases, the second is designed for light hydrocarbons, and the third specifically for hydrogen. The system is configured and fully tested in accordance with industry standard methods. Total analysis time for all components is less than 25 minutes.

The standard RGA analyzer is the most powerful tool to analyze the widest range of RGA type streams. This includes sample streams with a high % level of components as in ethylene, propylene and butylene streams.

- Rapid RGA:** A three channel 456-GC that utilizes a multi-valve design in which the packed columns used in the Standard RGA are replaced by micro-packed columns in both the hydrogen and permanent gas channels. Since the micro-packed columns are installed in a separate heated zone, the capillary columns located in the GC oven can be temperature programmed in a more aggressive manner. For high sample analysis demand, the Rapid RGA Analyzer concept provides a substantial reduction in overall analysis time of 5 minutes (7 minutes with H₂S) compared to the 25 minutes with the standard RGA.




Figure 5: The analysis of light hydrocarbons using the Rapid RGA, with complete separation in less than five minutes.




Figure 6: The analysis of permanent gases and hydrogen using the Rapid RGA.




Figure 7: A shows a 'traditional' RGA with all columns mounted in oven. B shows the micro-packed columns mounted in the separate heated zone in the Rapid RGA.

● Bruker Refinery Gas Analyzers

Table 1: RGA Analyzer Characteristics.

Characteristics	Standard RGA	Rapid RGA
No. of Channels/Detectors Used	3	3
No. of Column Ovens	1	2
Analysis Time	25 min	5 min (7 min with H ₂ S)
Repeatability	<1%	<1%
Linear Bench Space Required	66 cm/26 in.	66 cm/26 in.
Minimum Component Detection Level	0.01% all components except H ₂ S = 0.05%	0.01% all components except H ₂ S = 0.05%
Suitability		
Typical Refinery Gas	Excellent	Excellent
Impurities in Bulk Ethylene	Excellent	Excellent
Impurities in Bulk Propylene	Excellent	Good
Impurities in Bulk C4	Good	Good

Table 2: Multiple channels of data are conveniently combined into a single analysis report.


Peak No.	Peak Name	Channel	RT (min.)	Result (g/l)	Norm. (%)	Area (uV/Sec.)
1	Hydrogen	Middle (TCD)	1.6967	36.0300	22.7681	390257
2	Carbon dioxide	Front (TCD)	2.6000	0.1000	0.0632	13376
3	Hydrogen sulfide	Front (TCD)	-	0.0000	0.0000	0
4	Oxygen	Front (TCD)	9.9200	0.0000	0.0000	37325
5	Nitrogen	Front (TCD)	10.3267	1.1990	0.7577	2122071
6	Carbon monoxide	Front (TCD)	-	0.0000	0.0000	0
7a	Methane	Front (TCD)	11.1917	11.9900	7.5767	1394584
7b	Methane	Rear (FID)	3.7350	11.9900	7.5767	1492388
8a	Ethane	Front (TCD)	3.5367	17.9900	11.3682	2867688
8b	Ethane	Rear (FID)	4.1283	17.9900	11.3682	4480322
9a	Ethylene	Front (TCD)	2.9550	29.9800	18.9449	4139442
9b	Ethylene	Rear (FID)	4.7217	29.9800	18.9449	7411134
10	Propane	Rear (FID)	6.1933	0.1990	0.1258	71402
11	Cyclo Propane	Rear (FID)	-	0.0000	0.0000	0
12	Propylene	Rear (FID)	-	0.0000	0.0000	0
13	i-Butane	Rear (FID)	-	0.0000	0.0000	0
14	n-Butane	Rear (FID)	-	0.0000	0.0000	0
15	Propadiene	Rear (FID)	-	0.0000	0.0000	0
16a	Acetylene	Front (TCD)	5.0283	0.5020	0.3172	49786
16b	Acetylene	Rear (FID)	16.4331	0.5020	0.3172	121300
17	t-2-Butene	Rear (FID)	18.5050	0.0990	0.0626	138647
18	1-Butene	Rear (FID)	-	0.0000	0.0000	0
19	i-Butene	Rear (FID)	19.5167	0.0990	0.0626	44492
20	cis-2-Butene	Rear (FID)	-	0.0000	0.0000	0
21	1,3-Butadiene	Rear (FID)	22.1367	0.0000	0.0000	16165
22	Propyne	Rear (FID)	-	0.0000	0.0000	0
23	C5+	Rear (FID)	2.9217	0.1000	0.0632	58164
Totals				158.2480	100.0000	

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- "Quick connect" fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time



For research use only. Not for use in diagnostic procedures.



www.ScionHasArrived.com



www.Bruker.com

Technical Note # CA-270111

Principles and performance of the Collision Reaction Interface for the aurora M90

Bruker's patented 90 degree ion optics provides unsurpassed efficiency of transfer of ions from the interface to the mass analyzer, which enables the aurora M90 ICP-MS to achieve a sensitivity of more than 1000 million c/s per mg/L of analyte, while maintaining oxide ratios (CeO^+/Ce^+) below 3%. The aurora M90 ICP-MS is also equipped with Bruker's unique and patented interference management system, the Collision Reaction Interface (CRI). This CRI technology reduces common polyatomic interferences on elements such as As, Se, Cr, V and Fe, thus achieving lower detection limits in hot plasma, even for samples with complex matrices.

Unlike other interference management systems, the CRI does not use a pressurized multipole prior to the mass analyzer. Instead, reaction and collision gases are injected directly into the plasma through the tips of the interface cones. This innovative approach reduces/removes interferences before the interfering and analyte ions are extracted into the ion optics. Due to highly efficient pumping within the interface region, switching between CRI and normal mode (no CRI gas) is very rapid, enabling multiple condition sets to be run on a single solution.

Basic Principles of Collision Reaction Interface (CRI)

As shown in Figure 1, the CRI works simply by injecting the reactive/collision gases into the plasma through the tips of the sampler and/or skimmer cones to induce collisions and/or ion-molecule reactions with interfering ions. Plasma conditions at the interface cone apertures are ideal for collisions and reactions to occur. The high plasma density and high temperature should lead to a high collision/reaction frequency between the interfering ions and the injected gases. As a result, most argon-based polyatomic interferences are destroyed or removed before they are extracted into the ion optics. Hydrogen and helium are used as CRI gases, as these gases provide efficient interference attenuation, and avoid the need to use expensive or corrosive gases such as methane or ammonia. When a collision/reaction gas is added to the plasma, a number of processes occur including charge transfer,

proton transfer, electron-ion reactions and ion-molecule interactions. For example, when hydrogen is injected, a polyatomic interfering ion such as Ar_2^+ (that interferes with ^{80}Se determinations) collides with a hydrogen molecule. A proton is transferred from the H_2 molecule to the Ar_2^+ ion, forming ArH^+ , a neutral H atom and a neutral Ar atom. The ArH^+ ion then collides with another H_2 molecule and a proton is transferred from the ion to the molecule, forming a neutral Ar atom and an H_3^+ ion. The H_3^+ ion ($m/z = 3$ amu/unit electronic charge) does not interfere with any isotopes of interest in ICP-MS.




Figure 1: Schematic diagram of the Bruker CRI system

When helium is injected, interaction between the electron clouds of the helium atoms and those of interfering polyatomic ions can make a large polyatomic interfering ion such as $^{35}\text{Cl}^{16}\text{O}^+$ (which interferes with ^{51}V) rotationally and vibrationally excited. In subsequent collisions, an excited $^{35}\text{Cl}^{16}\text{O}^+$ ion can receive sufficient energy to bring about its dissociation, removing the $^{35}\text{Cl}^{16}\text{O}^+$ interference from the ^{51}V determination. In general, the analyte ions are not removed by the CRI technique; they simply lose energy as they collide with the CRI gases. Polyatomic ions, in contrast, are destroyed and removed when colliding with the CRI gases, mainly due to the charge transfer reactions, or collisional excitation-dissociation processes.

The Effect of Injecting CRI Gas at the Skimmer and Sampler Cone

Experimental results obtained so far suggest that gas injected into the plasma from the aperture of the sampler cone has less effect on the interferences reduction/removal. Figure 2 is a 3-D graph showing typical results of interference reduction experiments, where the sensitivities for ^{75}As and ^{89}Y were measured at various CRI gas flow rates through both the skimmer and sampler cones. The test solution (Var-IS-1) contains $10\ \mu\text{g/L}$ Y in 0.5% HNO_3 and 0.5% HCl matrix solution, and helium was used as the CRI gas. Since no arsenic was present in the solution, the measured sensitivity for $^{75}\text{As}^+$ is entirely due to the interfering polyatomic ion $^{40}\text{Ar}^{35}\text{Cl}^+$.

It is clear from Figure 2 (A), that the apparent signal for ^{75}As (actually from $^{40}\text{Ar}^{35}\text{Cl}^+$) is decreasing with increasing He flow into the skimmer. At a flow rate around $120\ \text{mL/min}$, the interference from $^{40}\text{Ar}^{35}\text{Cl}^+$ is completely removed. The efficiency of interference removal is demonstrated in Figure 2 (B), where the signal ratio of $^{89}\text{Y}/^{75}\text{As}$ (i.e., ratio of a real analyte ion to the signal from an interfering ion) is improving with increasing He flow rate into the skimmer. In contrast, the He flow at the sampler cone has very little impact on either the $^{40}\text{Ar}^{35}\text{Cl}^+$ signal or the $^{89}\text{Y}/^{75}\text{As}$ ratio, as shown in Figure 2. Most of the results shown in this report were




Figure 2: Typical 3-D graph of isotope sensitivity versus CRI gas flow at both skimmer and sampler apertures

obtained with the CRI gases injected into the plasma at the aperture of the skimmer cone.

Interference Reduction by using the CRI

As a simple example of using the CRI technique to remove the Ar_2^+ interferences, Figure 3 shows a time scan graph for $^{115}\text{In}^+$ and $^{80}\text{Se}^+$ (i.e., Ar_2^+), where the CRI gas flow is increasing stepwise.

Also shown in Figure 3 is the signal ratio of $^{115}\text{In}^+ / ^{80}\text{Se}^+$. The time scan was carried out with a test solution (Var-IS-1) containing $10\ \mu\text{g/L}$ of internal standard elements (Bi, In, Li, Sc, Tb, and Y). During the time scan, the H_2 gas was injected into the plasma through the CRI skimmer cone tip. The H_2 flow rate was stepped from 0, 20, 50, 80, 100, to $120\ \text{mL/min}$ with about 50 seconds between each step. The efficiency of interference reduction/removal using CRI is clearly demonstrated in Figure 3, in which the signal for the interfering species Ar_2^+ (scanned as $^{80}\text{Se}^+$) is progressively decreasing with increasing H_2 gas flow rate, while the signal to interferent level (determined as the ratio $^{115}\text{In}^+ / ^{80}\text{Se}^+$) is continuously improving. As seen in Figure 3, at a H_2 flow rate around $120\ \text{mL/min}$, the interference from Ar_2^+ has been completely removed, and the sensitivity for the analyte $^{115}\text{In}^+$ is still maintained at a high level (i.e., over $50\ 000\ \text{c/s}$ per $1\ \mu\text{g/L}$ of indium).

Another example of using the CRI technique to effectively reduce or remove some common plasma-based polyatomic interferences is shown in Figure 4, which shows mass spectra of high purity water in the m/z range from 42 to 78 amu/unit electronic charge. One spectrum is obtained under the conventional ICP-MS operation mode ('normal mode'), and the other is obtained under the CRI mode where H_2 gas is injected through the skimmer cone into the plasma at a flow rate of $120\ \text{mL/minute}$.




Figure 3: Time scan for $^{115}\text{In}^+$, $^{80}\text{Se}^+$ (i.e., Ar_2^+), and the ratio of $^{115}\text{In}^+ / ^{80}\text{Se}^+$ with increasing H_2 gas flow

Removal of plasma-based polyatomic




Figure 4: Mass spectra of pure water obtained under conventional ICP-MS and CRI-ICP-MS conditions

Calibration curves




Figure 5: Typical calibration curves for ^{56}Fe and ^{80}Se , with hydrogen gas added at the skimmer at 100 mL/min, calibration standards; 0 ng/L, 50 ng/L, 100 ng/L and 1000 ng/L of Fe and Se

Typical detection limits in CRI mode

The principal benefit of using the CRI technique is to reduce common polyatomic interferences on elements such as As, Se, Cr, V and Fe, thereby improving the detection limits for these elements in hot plasma. Table 1 lists the typical detection limit (DL) measured in the CRI mode. All DLs were determined as the concentration corresponding to 3 times the standard deviation of 10 replicates of a blank (i.e., 1% HNO_3). The measurements were made under routine analytical laboratory, not 'clean-room', conditions. Hence, the DL values listed in Table 1 can be routinely achieved outside a 'clean-room' in a clean laboratory.

In addition to the DLs, the background equivalent concentrations (BECs) when the DLs were determined under the CRI mode, are listed in the table. Hydrogen was used as the CRI gas, and it was injected into the plasma through the skimmer cone at a flow rate 80 mL/min for most of the elements listed in the table.

In some cases, the amount of the CRI gas required to remove the interferences can vary depending on the sample. In general (if the interfering ion is formed from the sample matrix), the greater the concentration of sample matrix, the more the CRI gas will be needed.

Determination of Iron using ^{56}Fe , and Selenium using ^{80}Se

In a traditional ICP-MS analysis, both Fe and Se are regarded as difficult elements to determine at trace level. Interferences from ArO^+ and Ar_2^+ make it impossible or extremely difficult to determine trace levels of Fe using the most abundant isotope (^{56}Fe) or trace levels of Se using the most abundant isotope (^{80}Se). In some cases, the cool plasma technique may be used. However, the cool plasma technique is subject to matrix effects that restrict its areas of application. The CRI technique removes polyatomic interferences before the ions enter the ion optics. Hence, even using hot plasma conditions; it is now possible to determine Fe and Se at trace level using their most abundant isotopes. Figure 5 shows typical calibration curves of Fe and Se at sub- $\mu\text{g/L}$ levels.

Determination of Iron in High Calcium Matrices

Another potential interference for ^{56}Fe is $^{40}\text{Ca}^{16}\text{O}^+$ from calcium in the sample matrix. Calcium is a common matrix element in most environmental waters, and its level can vary significantly depending on the source of the water samples. The CRI can effectively remove the CaO^+ interference allowing trace level determinations of Fe using the ^{56}Fe isotope. Figure 6 shows the recovery of 1 $\mu\text{g/L}$ of iron in the presence of various concentrations of calcium.

Table 1: Typical CRI DLs and BECs for selected isotopes

Isotope	CRI (H_2)	Non CRI	CRI (H_2)	Isotope	CRI (H_2)	Non CRI	CRI (H_2)
	DL (ng/L)	DL (ng/L)	BEC (ng/L)		DL (ng/L)	DL (ng/L)	BEC (ng/L)
^9Be	0.5	3	1.1	^{59}Co	0.2	0.2	0.8
^{23}Na	13	200	252	^{60}Ni	10	2	68
^{24}Mg	0.5	2	4.7	^{63}Cu	1	0.3	5.3
^{25}Mg	1	5	5.3	^{65}Cu	1	2	7.0
^{27}Al	0.8	2.0	7.9	^{66}Zn	1.5	5	15
^{39}K	43	500	328	^{68}Zn	1.2	20	17
^{40}Ca	2.5	-	23	^{75}As	0.6	20	2.2
^{44}Ca	6.5	500	81	^{78}Se	1.5	400	26
^{49}Ti	1.3	3	3.2	^{80}Se	8.8	-	193
^{51}V	0.15	3	1.7	^{98}Mo	0.7	0.4	2.7
^{52}Cr	0.6	8	13	^{107}Ag	0.2	0.6	1.2
^{53}Cr	1.5	3	23	^{111}Cd	0.2	0.2	0.2
^{55}Mn	0.4	2	10	$^{206,7,8}\text{Pb}$	0.1	0.3	1.9
^{56}Fe	1.5	4000	167	^{232}Th	0.08	0.04	0.9
^{57}Fe	44	300	2146	^{238}U	0.01	0.06	0.07

Determination of Arsenic in Chloride Matrices

In a chloride matrix, the $^{40}\text{Ar}^{35}\text{Cl}^+$ ion interferes with $^{75}\text{As}^+$ (the only naturally-occurring arsenic isotope). Conventional correction equations may be used, but the resulting detection limits are often inadequate. The CRI efficiently removes the ArCl^+ ion from the plasma before it is extracted into the ion optics, and allows routine detection of As at low ng/L levels. Figure 7 shows the recovery of 1 $\mu\text{g/L}$ As in various HCl concentrations, with and without CRI. Correction equations were not applied in this experiment. Under the CRI mode, H_2 was added through the skimmer cone at a flow rate of 105 mL/min.

Long Term Stability Performance in the CRI Mode

Figure 8 shows the long term signal stability in a high total dissolved solids matrix in the CRI mode. The stability was tested over a period of 5 hours using a solution containing 0.1% w/v (1000 mg/L) NaCl spiked with 1 $\mu\text{g/L}$ of various analytes. The results in Figure 8 show that the aurora M90 ICP-MS, running in the CRI mode, was extremely stable over 5 hours, even in a matrix solution containing 1000 mg/L NaCl (i.e., TDS=0.1%). The relative standard deviations (RSD) of the measured signals for most analytes did not exceed 5% over the 5 hours. Similar results were also obtained using helium as the CRI gas.

Summary

The aurora M90 ICP-MS provides a simple, but very effective approach to removing common polyatomic interferences in ICP-MS analysis. The unique CRI technology used on the aurora M90 ICP-MS reduces common polyatomic interferences on elements such as As, Se, Cr, V and Fe, thus achieving lower detection limits in hot plasma, even for samples with complex matrices.

Keywords

ICP-MS
Collision Reaction Interface
Detection Limits

Instrumentation

aurora M90 ICP-MS




Figure 6: Recovery of 1 $\mu\text{g/L}$ Fe in various Ca concentrations; H_2 added through the skimmer at 120 mL/min




Figure 7: Recovery of 1 $\mu\text{g/L}$ of As in various HCl concentration levels, with and without CRI




Figure 8: Long-term signal stability in 1000 mg/L NaCl; CRI gas H_2 ; flow rate 100 mL/min at the skimmer

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

www.bruker.com/chemicalanalysis

● EVOQ Elite LC Triple Quadrupole Mass Spectrometer

Specification Sheet

Performance Specifications

Mode	Test	Specification
Positive ESI MRM	50 fg of reserpine injected	S/N \geq 2000:1

The Signal-to-Noise ratio S/N values are based on RMS

Dimensions (H x W x D) and Weight

- Mass Spectrometer: 53 cm (H) x 45 cm (W) x 70 cm (D), 193 lbs/68kg
- UHPLC Advance: 71 cm (H) x 76 cm (W), 56 cm (D), 75lb/34kg

Analyzer – EVOQ Elite Specifications

- Scan modes: Full Scan with Q1; Precursor, Product, Neutral Loss/Gain Monitoring, Selected Ion Monitoring (SIM), Multiple Reaction Monitoring (MRM)
- Standard ionization modes: Heated Electrospray Ionization, Atmospheric Pressure Chemical Ionization (APCI)
- Ion source: Fixed 90° spray, single housing for HESI to APCI
- Cone gas orifice temperature up to 400 °C
- HESI and APCI Source temperature: Up to 750 °C
- Mass filters: quadrupole with pre- and post-filters
- Collision cell: 180° curved path with lens-free design
- Collision cell gas: Argon with adjustable pressure up to 2 mTorr
- Collision energy: selectable up to 75 eV
- Mass range (m/z): 10 to 1250 Da
- Scan rate: up to 14,000 Da/sec
- Minimum dwell times: 1 ms
- Maximum acquisition MRM rate: 500 MRM's/sec
- Resolution: adjustable from Unit (0.7 Da) to 4 Da, also with three selectable settings (Unit, Standard, Open) on both Q1 and Q3.
- Mass axis stability: $\leq \pm 0.1$ Da over 24 hours with normal temperature variations (+/- 3°C)
- Manifold temperature: 40-50 °C
- Detector: Electron multiplier with ± 5 kV post acceleration and with an option for on-the-fly multiplier gain optimization for Extended Dynamic Range (EDR); direct ion collection onto multiplier for negative ion detection without dynode loss
- Turbo molecular pump: Three stage, 25/300/400 L/sec

For research use only. Not for use in diagnostic procedures.

Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice.

© Bruker Daltonics 09-2012, #703715

www.bruker.com/chemicalanalysis

- Roughing pump: Two single stage 40 meter³/hour pumps
- Gas Requirements:
 - Nitrogen (nebulizing gas, cone gas and probe gas) 32 L/min @ 80 psi
 - Compressed air 50 L/min @ 80 psi
 - Argon (collision cell) high purity gas @ 20-50 psi
- Power requirements:
 - For mass spectrometer: 16A, 200 - 240V
 - For roughing pump control box, 20A, 200-240 V.
- Operating environment temperature: 15 °C to 30 °C
- Operating environment humidity: 20% to 80% relative humidity (without condensation)
- Syringe Pump: eVol handheld syringe drive
- Divert Valve: six port integrated valve and software controlled

Software

- Bruker MS Workstation equipped with the Compound Base Scanning (CBS) MRM library for data acquisition, data handling and reporting
- PACER data processing with Exception Based Review

Liquid Chromatograph (Bruker Advance)

- Standard Injection Volume Range: 1 - 100 ul
- Injection Precision: <0.25% RSD
- Standard Sample Capacity: 6 deepwell or 6 microplate or 6x 54 2 ml vials
- Sample Temperature Range: 4- 40 °C
- Column Heater: 5 °C above room temperature to 90 °C
- Heater Capacity: Single <250 mm column
- Injector: UHPLC Capable
- Carryover: Less than 0.003% (30 ppm)
- Solvent Channels: 2 binary gradient, with optional third isocratic pump
- Flow Range: 0.005 - 2.5 ml/min
- Max Pressure: HPLC - 8000 PSI, UHPLC - 15,000 PSI
- Dead Volume: <100ul
- pH Range: 2-12
- Degasser: 2 Channel (3 with optional third isocratic pump)
- Flow Rate Precision: 0.06% RSD
- Flow Rate Accuracy: 1%
- Line Voltage: 100 – 240 VAC, ± 10%
- Line Frequency: 50 or 60 Hz, ± 5%
- Power Consumption: 400 VA Maximum
- Communication: USB

• EVOQ Qube LC Triple Quadrupole Mass Spectrometer

Specification Sheet

Performance Specifications

Mode	Test	Specification
Positive ESI MRM	200 fg of reserpine injected	S/N \geq 2000:1

The Signal-to-Noise ratio S/N values are based on RMS

Dimensions (H x W x D) and Weight

- Mass Spectrometer: 53 cm (H) x 45 cm (W) x 70 cm (D), 193 lbs/68kg
- UHPLC Advance: 71 cm (H) x 76 cm (W), 56 cm (D), 75lb/34kg

Analyzer – EVOQ Qube Specifications

- Scan modes: Full Scan with Q1: Precursor, Product, Neutral Loss/Gain Monitoring, Selected Ion Monitoring (SIM), Multiple Reaction Monitoring (MRM)
- Standard ionization mode: Heated Electrospray Ionization, Atmospheric Pressure Chemical Ionization (APCI)
- Ion source: Fixed 90° spray, single housing for from HESI to APCI
- Cone gas orifice temperature up to 400 °C
- HESI and APCI Source temperature: Up to 750 °C
- Mass filters: quadrupole with pre- and post-filters
- Collision cell: 180° curved path with lens-free design
- Collision cell gas: Argon with adjustable pressure up to 2 mTorr
- Collision energy: selectable up to 75 eV
- Mass range (m/z): 10 to 1250 Da
- Scan rate: up to 14,000 Da/sec
- Minimum dwell times: 1 ms
- Maximum acquisition MRM rate: 500 MRM's/sec
- Resolution: adjustable from Unit (0.7 Da) to 4 Da, also with three selectable settings (Unit, Standard, Open) on both Q1 and Q3.
- Mass axis stability: $\leq \pm 0.1$ Da over 24 hours with normal temperature variations (+/- 3°C)
- Manifold temperature: 40-50 °C
- Detector: Electron multiplier with ± 5 kV post acceleration and with an option for on-the-fly multiplier gain optimization for Extended Dynamic Range (EDR); direct ion collection onto multiplier for negative ion detection without dynode loss
- Turbo molecular pump: Three stage, 25/300/400 L/sec

For research use only. Not for use in diagnostic procedures.

Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice.

© Bruker Daltonics 09-2012, #703716

www.bruker.com/chemicalanalysis

- Roughing pump: One single stage 40 meter³/hour pump
- Gas Requirements:
 - Nitrogen (nebulizing gas, cone gas and probe gas) 32 L/min, @ 80 psi
 - Compressed air @ 80 psi
 - Argon (collision cell) high purity gas @ 20-50 psi
- Power requirements:
 - For mass spectrometer: 16A, 200 - 240V
 - For roughing pump control box: 10A, 200-240 V.
- Operating environment temperature: 15 °C to 30 °C
- Operating environment humidity: 20% to 80% relative humidity (without condensation)
- Syringe Pump: eVol handheld syringe drive
- Divert Valve: six port integrated valve and software controlled

Software

- Bruker MS Workstation equipped with the Compound Base Scanning (CBS) MRM library for data acquisition, data handling and reporting
- PACER data processing with Exception Based Review

Liquid Chromatograph (Bruker Advance)

- Standard Injection Volume Range: 1 - 100 ul
- Injection Precision: <0.25% RSD
- Standard Sample Capacity: 6 deepwell or 6 microplate or 6x 54 2 ml vials
- Sample Temperature Range: 4- 40 °C
- Column Heater: 5 °C above room temperature to 90 °C
- Heater Capacity: Single <250 mm column
- Injector: UHPLC Capable
- Carryover: Less than 0.003% (30 ppm)
- Solvent Channels: 2 binary gradient, with optional third isocratic pump
- Flow Range: 0.005 - 2.5 ml/min
- Max Pressure: HPLC - 8000 PSI, UHPLC - 15,000 PSI
- Dead Volume: <100ul
- pH Range: 2-12
- Degasser: 2 Channel (3 with optional third isocratic pump)
- Flow Rate Precision: 0.06% RSD
- Flow Rate Accuracy: 1%
- Line Voltage: 100 – 240 VAC, ± 10%
- Line Frequency: 50 or 60 Hz, ± 5%
- Power Consumption: 400 VA Maximum
- Communication: USB



Application Note # CA-270109

High Matrix Sample Analysis with the aurora M90 ICP-MS using the Collision Reaction Interface (CRI)

Introduction

Inductively Coupled Plasma Mass Spectrometry is a powerful tool for a wide range of applications. The excellent sensitivity, multi-element capability and the wide analytical range give this technique the ability to analyze trace levels in a variety of environmental samples.

The unique design of the aurora M90 ICP-MS features a patented 90-degree ion mirror [1] that significantly reduces the background signals. It allows photons and neutrals to pass directly through the hollow structure of the mirror while analyte ions are directed into the mass analyzer. The amount of analyte passing through the skimmer cone to the quadrupole is over 80%, which provides excellent ion transmission efficiency [2].

The new aurora M90 ICP-MS, in addition to the benefits of the ion mirror, has also been designed to reduce/eliminate interferences. The unique Collision Reaction Interface (CRI) introduces a new era of interference management for ICP-MS systems.

This application note demonstrates the capability of the aurora M90 system to determine sub- $\mu\text{g/L}$ concentrations of a large number of elements in solutions containing high amounts of total dissolved solids (often defined as "high matrix" samples) as waste extracts or digests.

Instrumentation

All measurements were carried out using the aurora M90 ICP-MS equipped with Collision Reaction Interface (CRI) and a Sample Preparation System, SPS3. The system operations of the aurora M90 are fully controlled by Bruker's Quantum software. The software provides one-step instrument setup, optimization, and method development. By using the auto-optimization routine supplied in the Quantum software, the instrument was automatically tuned to chosen CRI conditions.

Materials and reagents

All calibration standards were prepared by diluting multi-element stock standards (Var-Cal-1 and Var-Cal-2, Inorganic Ventures, Inc., Lakewood, NJ, USA) with 1% v/v nitric acid (Ultrapur[®] HNO₃ 60%, Merck, Kilsyth, Victoria, Australia). Working standards below 10 $\mu\text{g/L}$ were prepared immediately before the measurement. An internal standard solution, containing 100 $\mu\text{g/L}$ of ⁶Li, ⁴⁵Sc, ¹¹⁵In, ⁸⁹Y, ¹⁵⁹Tb, and ²⁰⁹Bi, was prepared by diluting a 100 mg/L of internal standard stock (Var-IS-1 Inorganic Ventures, Inc., Lakewood, NJ, USA).

Sample preparation

High matrix samples were prepared to match the EPA 6020 interference check sample matrix (ICS A) [3]. Samples and spikes were prepared by dilution of a stock solution (6020ICS-A, Inorganic Ventures, Inc., Lakewood, NJ, USA) and then spiked to a final concentration of 0.5 $\mu\text{g/L}$ with the multi-element stock standard (Var-Cal-1 and Var-Cal-2, Inorganic Ventures, Inc., Lakewood, NJ, USA). Table 1 lists the elemental concentrations of the high matrix sample (ICS A) and the spiked solution (ICS AB).

The EPA 6020 method requires MDLs to be determined using a spiked, fortified blank solution. That is, a sample at low levels with no matrix. By using ICS-AB to determine the MDLs, the data below shows that even with a highly demanding sample (with a high matrix) the aurora M90 ICP-MS can easily meet the required control limits.

Table 1: Composition of high matrix sample and spike.

Analytes	Concentrations in µg/L	
	ICS A	ICS AB
Al	100,000	100,000
Ca	100,000	100,000
Fe	100,000	100,000
Mg	100,000	100,000
K	100,000	100,000
Na	100,000	100,000
C	200,000	200,000
Cl	1,000,000	1,000,000
P	100,000	100,000
S	100,000	100,000
Mo	2,000	2,000
Ti	2,000	2,000
As		0.5
Cd		0.5
Co		0.5
Cu		0.5
Mn		0.5
Ni		0.5
Ag		0.5
Zn		0.5
Sb		0.5
Ba		0.5
Pb		0.5
Se		0.5
Tl		0.5
V		0.5

Table 2: ICP-MS conditions for analyses of interference check samples.

Parameters	CRI 1	CRI 2
Skimmer gas	H ₂	He
Sampler gas	none	none
Skimmer flow (mL/min)	90	180
Sampler flow	0	0
Outer flow	16.5	16.5
Intermediate flow	1.65	1.65
Sheath gas	0.2	0.2
Nebulizer flow	0.95	0.98
RF power (kW)	1.3	1.3
Sampling depth (mm)	6.5	6.5
Pump rate (rpm)	3	3
Stabilization delay (s)	60	60
Spray chamber (°C)	3	3
First extraction lens (V)	-1	-1
Second extraction lens (V)	-21	-75
Third extraction lens (V)	-195	-235
Corner lens (V)	-177	-197
Mirror lens left (V)	39	38
Mirror lens right (V)	35	39
Mirror lens bottom (V)	23	24
Entrance lens (V)	0	-2
Fringe bias (V)	-2.5	-3.5
Entrance plate (V)	-29	-28
Pole bias (V)	0	0
Scan mode	Peak hopping	Peak hopping
Dwell time (ms)	20	20
Points per peak	1	1
Scans/Replicate	20	20
Replicates/Sample	5	5

Conditions

For the analysis of all samples in this work, normal sensitivity mode was used. Two condition sets were used in this study. Both sets used CRI conditions. The method parameters used for the two conditions sets are summarized in Table 2.

Discussion

In order to verify corrections for elemental and polyatomic isobaric interferences EPA Method 6020 requires the analysis of two interference check samples, ICS A and ICS AB, at the beginning of the analysis run.

In this work, the ICS AB solution (0.5 µg/L) was analyzed to evaluate a lower working range in the presence of a high

matrix. Table 3 summarizes the found results for the ICS A and ICS AB solutions and recoveries of the 0.5 µg/L spike for the different CRI conditions.

This spike level of 0.5 µg/L is 40 times lower than EPA 6020 requires (20 µg/L). This was done to make the analysis more demanding.

The Method Detection Limits (MDLs) and Method 6020 control limits are also listed. Seven subsequent reading of the ICS AB solution were used to calculate MDLs for each of the selected isotopes.

Table 3. Results summary for Interference check samples spike recoveries.

Element	Isotope (m/Z)	Spike level		With CRI		
		ICS AB (µg/L)	Recovery %	MDL* (µg/L)	CRDL (µg/L)	CRI condition
V	51	0.50	94	0.183	5	[2]
Cr	53	0.50	110	0.259	5	[2]
Mn	55	0.50	101	0.103	5	[2]
Co	59	0.50	98	0.032	5	[2]
Ni	60	0.50	103	0.089	5	[2]
Cu	63	0.50	107	0.297	10	[1]
Zn	66	0.50	145	0.253	10	[1]
As	75	0.50	101	0.292	1	[1]
Se	78	0.50	105	0.257	5	[1]
Ag	107	0.50	99	0.050	5	[2]
Cd	111	0.50	102	0.107	1	[2]
Sn	118	0.50	99	0.100	1	[2]
Sb	121	0.50	102	0.027	-	[2]
Ba	137	0.50	102	0.031	5	[1]
Tl	205	0.50	94	0.038	5	[2]
Pb	206+7+8	0.50	97	0.036	10	[1]

[1] Skimmer H₂, [2] Skimmer He

* MDL were calculated using SD of seven subsequent readings of interference check sample spike (ICS AB) multiplied by 3.14

References

- [1] I. Kalinitchenko, Ion Optical system for a Mass Spectrometer, Australian Patent 750860, 14 November 2002
- [2] S. Elliott, M. Knowles and I. Kalinitchenko, "A New Direction in ICP-MS", Spectroscopy, 19(1), 30 (2004)
- [3] EPA Method 6020 Inductively Coupled Plasma – Mass Spectrometry, www.epa.gov/epaoswer/hazwaste/test/pdfs/6020.pdf

Keywords

ICP-MS
Collision Reaction Interface
High Matrix
Multiple Condition Sets

Instrumentation & Software

aurora M90 ICP-MS
Bruker Quantum

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

www.bruker.com/chemicalanalysis



aurora m90

- A new era in ICP-MS

Your Partner in ICP-MS Solutions



Bruker continues to find new and novel ways to meet your changing needs. As a leader in elemental analysis you can be assured that when you buy a Bruker ICP-MS, you're buying more than just an instrument. You're buying a relationship with one of the most respected and experienced instrument companies in the world.

● ICP-MS: It's never been Easier

Bruker innovation, making ICP-MS easier

If you've ever wished that ICP-MS could be simpler, wish no more. The aurora M90 makes light work of it. No matter what your requirements, with a Bruker ICP-MS, you can tackle any application with ease.

Key benefits of aurora M90 include:

- Bruker's patented high-efficiency 90 degree ion optics and double off-axis quadrupole delivers exceptionally low background noise and unmatched sensitivity – at more than 1 million counts per second for 1µg/L.
- Tunable from normal to high sensitivity, the aurora M90 is perfect for both routine and research-grade applications – Flexibility at your fingertips.
- The aurora M90 delivers industry leading detection limit performance. Collision/reaction interface (CRI) technology makes setup of complex cell systems a thing of the past. Simply turn on the gas flow to remove interferences. It's that simple.
- Featuring the only all-digital ICP-MS detector, covering more than nine decades of dynamic range in pulse counting mode, the aurora M90 delivers fast and accurate multi-element analysis from ultra-trace to major levels in a single measurement.

Let Bruker Quantum work for you

If your goal is to spend less time creating methods and optimising conditions, and more time running samples, Bruker Quantum software delivers. Enjoy accurate results in less time with an intuitive yet flexible user interface that takes the hard work out of ICP-MS.


- With auto-tuning of all instrument parameters, you can spend less time on instrument setup and more time on sample analysis. Saving you valuable time and money.
- Fully automated, aerosol dilution extends the high dissolved solids tolerance of your ICP-MS allowing you to directly analyze challenging samples without additional sample preparation.



Innovation you can trust

- Patented 90 degree ion mirror and low noise double off-axis quadrupole provide industry leading sensitivity and background for lowest detection limits.
- New and improved Collision Reaction Interface (CRI II) provides even simpler and more effective removal of troublesome interferences for interference-free analysis of your samples.
- Robust, high-efficiency plasma system and patented Interlaced Coils break down your toughest sample matrices, reduce matrix effects, and minimize ion energy spread for maximum sensitivity and stability.
- All-digital extended range detector – means fewer dilutions, and longer detector lifetime for greater productivity and lower running costs.

Our 90 Degree Reflecting Ion Optics Revolutionizes ICP-MS Performance



● Fast, Flexible, Interference-Free Analysis

Bruker is proud to bring you CRI II, now even simpler to use and more effective at removing troublesome interferences from your sample analysis.

The CRI injects helium (He) and hydrogen (H₂) collision and reaction gasses directly into the plasma as it passes through the orifice of the skimmer cone.

This innovative approach suppresses interferences before the analytes are extracted into the ion optics.

It's that simple!



- No need for expensive or corrosive gasses such as ammonia or methane, so laboratory costs are reduced.
- No additional cleaning as CRI forms part of the cone interface, making this interference management system maintenance-free.

Choose your analysis mode

CRI II universal analysis mode provides fast and accurate results for samples routinely encountered across the wide range of environmental and industrial monitoring processes.

While multi-mode delivers the ultimate in performance and flexibility for any sample type including those less routine in nature.

So if you need absolute confidence in your results, no matter what the sample, CRI II is the answer.



Dramatically reduce or eliminate troublesome plasma and sample matrix based interferences, using CRI II.

The solution to your application needs

With a vast range of options to choose from, Bruker has the solution to your application. Choose from

- CRI II for fast, accurate interference-free analysis of your samples.
- High sensitivity interface. Ideal for research applications on non-interfered isotopes, pushing your detection capability to levels never seen before.
- Make simple work of your most challenging samples. Upgrade to the inert vacuum pump system for low maintenance, high performance ICP-MS operation.
- Clean room package is suited to applications in the semiconductor industry and provides an inert and contamination-free environment for ultra-trace analysis.
- Application-specific sample introduction systems for routine analysis of geochemical and petrochemical samples.
- A range of autosampler and productivity-enhancing accessories provide you with fast, unattended operation of your ICP-MS.
- The nitrox 500 accessory allows online addition of nitrogen or oxygen gas to the plasma. Add nitrogen to lower your detection limits on key elements like As and Se. Add oxygen for routine analysis of organic solvents.
- Fully-integrated speciation options for the analyst wanting to know more about their samples.
- The aurora M90 is compatible with a wide range of laser ablation systems providing you with solution-free analysis.



Clean room package helps semiconductor labs attain clean room conditions for ultra-trace analysis.




All PFA sample introduction kits allow direct analysis of samples containing free HF acid.



Choose from a range of accessories to meet your sample handling requirements.



Boost the performance of your aurora M90 by adding the nitrox 500.



Fully-integrated HPLC interface makes light work of speciation applications.

● The Benchmark in Analytical Performance

Quickly and reproducibly reduce interferences

With CRI II you can quickly switch from CRI to non-CRI, or between different collision and reaction gases.


Multi-mode delivers the ultimate in performance and flexibility for any sample.

Determination of As in Cl containing samples

Use CRI II in H₂ mode to remove the ArCl interference when determining As in high chloride containing samples like blood, serum and urine.



1 µg/L As in HCl matrix




Comparison plots showing 1 µg/L spike recoveries for ⁷⁵As without correction equations. ArCl interferences are removed, allowing accurate trace level quantification of As.

	Certified range µg/L	Measured value µg/L
²⁷ Al	13 – 21.2	20
⁵¹ V	0.27 – 0.37	0.29
⁵² Cr	0.42 – 0.78	0.42
⁵⁶ Fe	404 – 460 mg/L	420 mg/L
⁷⁵ As	1.4 – 2.2	1.8
⁷⁸ Se	74.4 – 85.2	77.2
^{206, 207, 208} Pb	26.2 – 29	27.6
²³⁸ U	0.16 – 0.18	0.17

Obtain accurate results in complex biological matrices. Above, certified and measured values for Reference Whole Blood Seronorm WB1 show that trace and major levels can be determined with accuracy and confidence using CRI II.

● Unrivalled Performance




Maximum dynamic range for food samples

Determine toxic, essential and nutritional elements in a single, all-digital measurement for optimum accuracy and precision.



Unequaled high sensitivity performance

The aurora M90's unique high sensitivity mode opens up a world of new possibilities in research.



Element	Units	Measured	Certified
²⁴ Mg	mg/kg	513	500
³⁹ K	mg/kg	3128	3100
⁴⁴ Ca	mg/kg	422	410
⁵⁶ Fe	mg/kg	39.0	40.7 ± 2.3
⁷⁵ As	mg/kg	0.024	(0.023)
⁷⁸ Se	mg/kg	0.026	(0.025)
¹¹⁴ Cd	mg/kg	0.0270	0.0284 ± 0.0014
²⁰⁶⁻⁸ Pb	mg/kg	0.182	0.187 ± 0.014

Above, certified and measured values for brown bread reference material BCR-191 showing accurate measurement from ultra-trace to major levels.





Achieve levels of detection never seen before in ICP-MS. Above, a typical calibration for ¹⁹⁴Pt in chloroplatinin acid in the ng/L to sub-ng/L range used in the determination of pharmacologically active Pt tracers of anti-cancer drugs.

● Setting the Benchmark for Ease-Of-Use

Bruker Quantum software

Bruker redefines ease-of-use with our Web-integrated ICP-MS worksheet software. Quantum features a range of automated options, including setup and initialization routines, such as plasma alignment, mass calibration and resolution tests. Bruker's AutoMax makes method development easy by automating all ion optics, nebulizer and plasma settings for optimum results. Including auto-optimization for Aerosol Dilution, Bruker Quantum makes light work of your most difficult samples.

The dynamic Instrument Status window provides a quick visual check of the status of all system components. It is an excellent diagnostic tool that maximizes instrument up time.



Each worksheet cell provides all the results you need – including concentrations, intensities, statistics, replicate readings and graphical mass scans.

Bruker Quantum switches automatically between multiple method condition sets within a single sample, giving optimum performance for specific element suites, without having to re-run samples.

Chemical Analysis Solutions

Laboratory gas chromatography systems



The 400 Series consists of two gas chromatographs and an associated range of analyzers and solutions designed for leading applications. These systems allow chemists and engineers to employ standard methods and/or high quality trace sample analysis, in the petrochemical, agrochemical and environmental industries.

The 450-GC is a highly affordable and powerful analytical instrument that offers robust operation in an easy-to-use package. The system gives users a broad choice of injectors, detectors, switching and sampling valves up to three channels. The high resolution color touch screen is intuitive and supports local languages. The Bruker 430-GC offers the same outstanding performance as the 450-GC but in a compact, single channel package that occupies about half the bench space of conventional multi-channel GC.

Triple quadrupole mass spectrometer

The Bruker 320-MS GC/MS stands at the forefront of configurable triple-quadrupole mass spectrometer systems. It offers: femtogram sensitivity, 10 – 2000 Da mass range, and a wide array of chromatographic and ionization configurations to uniquely match your needs - all in less than 72 cm. (28 in.) of linear bench space! In minutes, the 320-MS can be changed from EI to CI modes of operation. Easily, the 320-MS is the most sensitive, robust, and flexible triple-quadrupole MS system currently available.

For research use only.
Not for use in diagnostic procedures.



● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

www.bruker.com/chemicalanalysis

● impact HD LC-MS System

Description

Part Number: # 1819695

Bench-top High-Definition UHR Time-of-Flight Mass Spectrometer



Easy-to-use, instant expertise™ Ultra-High-Resolution electrospray ionization quadrupole time-of-flight LC/MS/MS mass spectrometer designed for exact mass and true isotopic measurements in both MS and MS/MS mode.

Bench-top, small footprint Mass Spectrometer system for exact mass and highest mass resolution at U-HPLC speed in both MS and MS/MS mode:

- Unique FSR technology with Full Sensitivity @ Maximum Resolution achieved without any time constraints, in MS and MS/MS mode
- Outstanding Mass Resolution and Accuracy in both MS and MS/MS
- High-resolution extracted ion chromatograms capabilities
- High performance hyperbolic quadrupole and collision cell for compound fragmentation
- True Isotopic Measurements

A. Apollo II (ESI) Source

- Highly sensitive ESI Source with proprietary dual ion funnel guide for gentle mass independent ion focusing and high ion transmission efficiency
- Combined Funnel-Hexapole-Cartridge with front access for easy maintenance
- Grounded needle for safety and easy sample introduction
- Suitable for U-HPLC, HPLC and CE coupling
- Heated counter current drying gas for gentle and efficient drying
- Ion lens system including in Source collision induced dissociation control (IS-CID)
- Pneumatic off-axis nebulization for flow rates up to 1 ml/min., with gradients from 100% aqueous to 100% organic
- Flow rate 1µl/min to 1ml/min
- Coated glass capillary for physical and electrostatic isolation
- Processor controlled HV and gas controller

B. High Mass Quadrupole Mass Filter:

- Hyperbolic quadrupole mass filter
- Ultra stable monolithic design
- RF-generator for monoisotopic precursor ion selection

C. Novel high-transmission CID Collision cell:

- Hyperbolic quadrupole broad-mass bandwidth design
- Fast radial ion ejection enabling fast MS/MS cycle
- RF-Generator with fast amplitude switching
- Collision gas controller

D. Orthogonal pulsed ion extraction and Time-of-Flight Mass Analyzer

- Interface housing and ion lens system
- In-line detector system for easy maintenance
- Ultra-stable high voltage switches with up to 20 kHz repetition rate and appropriate power supplies.
- TOF analyzer with orthogonal mounted ion source
- In-flight refocusing optics for uncompromised sensitivity
- Dual stage ion reflectron with increased mass resolution and accuracy
- Ultrafast 5 GS/sec 10 bit digitizer
- High-sensitivity and fast ion detector system, mechanical adjustment in micrometer range
- Positive and negative ion modes
- Ultra-stable high voltage power supplies for TOF analyzer and detector

E. Vacuum system

- Analyzer vacuum housing
- Vacuum system with 5 differential pumping stages
- One roughing pump 28 m³/h, tri- and dual-stage turbo pumps for ESI source and UHR-TOF analyzer
- Vacuum measurement and pump control unit

F. Syringe pump**G. Modes of Operation**

- TOF Mass ranges 20-40,000 m/z
- Internal calibration (MS and MS/MS)
- External calibration (MS and MS/MS)
- Exact mass measurements independent from sample concentration over a wide dynamic range without second sprayer.

H. High-resolution-performance and accurate mass features

- Patented ion funnel source
- One-shot acquisition mode, no tuning for mass range optima
- Enhanced low mass sensitivity
- Superior MS/MS sensitivity
- Ultra broad mass-bandwidth

- Long term and ultra stable mass axis stability in MS and MS/MS
- Exact mass independent from sample concentration charge state and collision energy
- Combined calibration for both MS and MS/MS
- Wide dynamic range for quantitation
- Advanced temperature compensated flight tube
- Positive / negative ion operation

I. Data system:

- PC Workstation with 2.66 GHz Single-CPU-Quad-Core-processor, 12 GB RAM, 2 TB hard disk, Ethernet connection for external networks
- DVD-ROM drive
- R/W DVD-ROM drive DL
- ≥ 24 " flat screen colour monitor
- Windows™ 7/32 operating system
- Laser printer
- Remote Service capability via 128-bit SSL-security web connection

J. Applications software

Software package Compass 1.6 for HPLC and MS control, data acquisition, post processing, and data analysis:

- Operation system Windows 7/32
- Compass / HyStar 3.2 for integrated control of most popular U-HPLC and HPLC systems and auto samplers and automation
- Instant Expertise™ features for intelligent autoMS/MS workflows
- Expert mode: extended control over instrument parameters for interactive system optimization of sophisticated exact mass methods
- Compass / Data Analysis software DA 4.1, including:
 - Advanced data processing with a high degree of automation
 - SmartFormula 3D™: Automated sum formula determination using MS and MS/MS data with both, accurate mass and isotopic fit.
 - QuantAnalysis™ quantitation package
 - LibrarySearch™ module for search of MS, MS/MS and MS_n spectra with advanced matching algorithm
 - Charge deconvolution module
 - Export of spectra and ion current profiles as Windows metafiles to word-processing programs
- SW License DataAnalysis 4.1
- SW License Charge Deconvolution for DA 4.1
- SW License MaxEntropy Deconvolution as an option

K. Set of manuals and reference CD-ROMs

L. Installation

M. Familiarization upon installation

N. 1 year warranty

O. Training course - for 2 participants.

● impact HD LC-MS System

Specification Sheet
Part Number: # 1819695

Impact HD UHR Time-of-Flight Mass Spectrometer System



Size	Bench-top: 64cm x 118cm (Footprint) 198cm (Height)
Weight	~ 210 kg
Vacuum System	5 stages, 28 m ³ /h rough pump
Apollo II ion funnel electrospray source	Flow rate: 1 µL/min – 1 mL/min
Mass Range	20 – 40,000 m/z
Quadrupole isolation	Up to 3,000 m/z
Quadrupole Mass Range	Up to 40,000 m/z
Mass accuracy in MS and MS/MS	With internal calibrant: better than 800ppb RMS Error With external calibrant: better than 2 ppm RMS Error
Calibration	ONE calibration valid for MS and MS/MS analysis. Calibration is independent from charge state of calibrant mass
Mass resolution	40,000 FSR (full sensitivity resolution)
Isotopic pattern	The true isotopic pattern is maintained due to TIP™ technology (True Isotopic Pattern) and allows three dimensional chemical characterizations of analytes via SmartFormula™3D algorithm using exact mass, TIP, and MS/MS fragment data.
SmartFormula™3D	Enables unambiguous formula determination at "sub-ppm" confidence level up to 1000 Da.

Mass stability & dynamic range	hrXIC (high resolution Extracted Ion Chromatogram) technology with better than +/- 1.0 mDa stability on centroid data values over an typical LC peak.
Full scan sensitivity in MS	ESI: Reserpine 1 pg S/N>100:1 RMS With Ion-Booster (optional): Reserpine 100 fg S/N>100:1 RMS
Full scan sensitivity in MS/MS	The signal height obtained from a consumption of 2.5 fmol of Glu-Fibrinopeptide B will be better than 100 counts on the most intense y' sequence ion from the MS/MS spectrum of the doubly charged precursor ion. This shall correspond to a signal to noise ratio better than 50:1. The MS/MS sensitivity specification is met while using quadrupole isolation of the precursor ion demonstrating that there is minimal transmission loss through the isolating quad. A solution of 100 fmol/ μ L Glu-Fibrinopeptide B shall be introduced at a flow rate of 3 μ L/min.
TOF repetition rate	Up to 20 kHz
Temperature compensation	Yes
Digitizer	5Gsample/sec ADC with 50 Gbit/sec
Dynamic range	10 bit ADC for high quantitative dynamic range
Acquisition rate	up to 50 Hz MS 50 Hz MS/MS (profile and peak detected spectra to disk)

Optional accessory

IonBooster	Optional ion source
APCI II	Optional ion source
APPI II	Optional ion source
GC-APCI	Allows for direct GC coupling (Optional ion source)
APLI	Optional ion source
CryoSpray	Optional ion source
Bruker CaptiveSpray nanoBooster	Optional ion source
On-/Off-Line Nanospray	Optional ion source
CE/MS interface	With grounded needle for easy CE-TOF set-up (Optional)
DIP	Direct Probe (Optional ion source)

Lock-out specifications for an Ultra High Resolution- Time of flight (UHR-TOF) system maXis HD

System description:

Complete, factory-new, working and ready-to-use LC-MS/MS system including:

1. mass spectrometer with ESI source allowing positive and negative ion detection; the needle of the ESI source has to be on ground potential. (2,3)
2. floorstanding ultra high resolution reflectron Time-of-Flight mass spectrometer with small footprint (maximally 1 by 1,50 m.(2,3)) and near-orthogonal electrospray ion source.
3. The spectrometer has to include the following parts:
 - a. Ion lens system with possibility to perform In Source-CID experiments;
 - b. quadrupole/quadrupolesystem for ion isolation and fragmentation for MS/MS experiments;
 - c. ion transfer line with patented dual ion funnel system for maximum transmission (1,2,3)
 - d. A hyperbolic quadrupole mass filter
 - e. New HDC Collision Cell with DC Gradient (1,2,3,4)
 - f. New 50 Gbit/sec Digitizer(1,2,3,4)
 - g. Ultra high resolution TOF analyzer with dual stage reflectron, in-flight refocusing optics and advanced temperature compensation.
 - h. Fast and robust detector with linear response for high dynamic range
 - i. Five-stage pumping system, incl. >28 m³ rough pump.
 - j. Needs to run with single rough pump (1,3)
 - k. detection system with an New 10 bit analog-to digitizer converter unit allowing recording of the complete analog ion signal (1,2,3,4)
 - l. PC-station with complete software (including Windows operating system and LC-MS control software) that allows acquisition, processing and exporting data + LCD-screen + laser printer;

Required parameters and features:

1. Analyzer resolution in single-reflection mode better than 75'000 (FWHM) at m/z 1222 without the need of multiple ion reflection at 50Hz acquisition rate, available without a loss of sensitivity (compared with the respective instrument specification);^(1, 2, 3, 4)
2. Analyzer resolution in single-reflection mode better than 30'000 (FWHM) at m/z 118 without the need of multiple ion reflection at 50Hz acquisition rate, available without a loss of sensitivity (compared with the respective instrument specification); (1,2,3)
3. MS and MSMS mass accuracy not worse than 600 ppb (calibrated internally) and 2 ppm (calibrated externally and without the use of any additional lock-spray) in a wide dynamic range i.e. independent from sample concentration. (1,2,3,4)
4. Long term and ultra-stable mass axis stability in MS and MS/MS, both with full sensitivity!
5. In scan Dynamic range \geq 4-5 orders of magnitude at 1Hz acquisition (2,3) without any need to split the ion beam dynamically (2;4) and without any need to switch the digitizer into under sampling mode (1).
6. A single calibration must be valid and applicable for both MS and MSMS measurements;
7. Possibility to use any substance for calibration chosen by the user (within available mass range); (3)
8. Mass range not less than 20 – 40'000 m/z; (3,4)
9. Possibility to select MS/MS precursors with m/z values up to 3000 (quadrupole isolation)(4)
10. Scan Speed: MS and MS/MS acquisition rate 50 Hz written to disc in 1 sec (1,2,3,4)
11. Fully automated isotope pattern matching with the generation of a list for the sum formula from both mass accuracy as well as the isotope pattern matching - it has to use both MS and MSMS data (3,4)
12. Possibility to maintain the isotopic pattern with an error of less than 2%.
13. Possibility to analyze intact proteins at m/z values > 4000 to obtain most comprehensive information (3,4)
14. Possibility of remote service diagnostics via secured Internet connection;
15. Possibility to create high resolution Extracted Ion Chromatograms within +/- 0,5 - 1,0 mDa error for screening of complex mixtures (4)
16. Possibility to save the hrEIC to disk (2,3).
17. Flexibility to control HPLC systems of various vendors (Dionex (Ultimate 3000), Waters (UPLC), Agilent, VWR/Hitachi, Shimadzu, Proxeon) within the original MS vendor software offered here (START/STOP signal is insufficient). (1,2,3)

¹ Against Agilent 6538, 6540

² Against Waters Synapt G2 SI

³ Against Thermo OT velos

⁴ Against AB Sciex 6500 TripleTOF

18. Flow rate of 1uL/minute-1mL/minute to facilitate flow injection as well as LC-driven experiments (1[6540])
19. MS sensitivity at reserpine of better than S/N 100:1 RMS at 1 pg Reserpine(at FIA flow rate \geq 200 uL/min) at full instrument resolution.With Ion-Booster (optional): better than S/N 100:1 RMS at 100 fgReserpine (1,2)
20. MS/MS sensitivity equal to or better than a consumption of 2.5 fmol peptide at S/N 100:1 when a solution of 100 fmol/ μ LGlu-Fibrinopeptide B is introduced at a flow rate of 3 μ L/min.
21. Comfortable Charge-state-ruler (1,2,3,4)
22. Possibility to upgrade the system to an LC-NMR-MS setup (4, 2,3)
23. Possibility to upgrade the system with an atmospheric pressure source to couple a GC system (3,4)
24. Changing of API sources from ESI to nano-spray and back without breaking vacuum (4)
25. Possibility to upgrade the system with a CaptiveSpray source with patented nano-booster
26. Possibility to upgrade the system with DIP (direct probe) to analyze insoluble compounds without the need for any sample preparation (1,3,4).
27. The software has to support and allow the quantitative analysis of proteins by the use of all protein labeling techniques, as well as of label-free experiments;
28. The software has to be able to combine results of several proteomics databases like Mascot and Phenyx in order to improve protein identification rates, under strict false positive rates for peptide identifications (1,2,3)
29. package dedicated for small molecules analysis allowing i.a.:
 - statistical comparison analyses of different sample groups
 - prediction of metabolites on the basis of biotransformation rules
 - differential display of LC-MS runs on samples taken from one organism at different times
30. package dedicated for multi-target compound screening allowing i.a.:
 - generation of extracted ion chromatograms (EIC)
 - automated generation of a compound list from those EIC
 - automated comparison of the detected compounds with a predefined target list
 - target list containing drugs, metabolites, natural and toxic products as well as pesticides
31. package for statistical analysis allowing i.a.: (3)
 - PCA analysis;
 - different normalization options;
 - data validation by test set and cross-validation;
 - data recalibration before post-processing;
 - evaluation results visualization;

maXis HD LC-MS System

Specification Sheet
Part Number: # 1820746

maXis HD UHR Time-of-Flight Mass Spectrometer System



Size	Floorstanding: 80 x 132 cm (Footprint) 290 or 320 cm (Height)
Weight	~ 500 kg
Vacuum System	5 stages, 28 m ³ /h rough pump
Apollo II ion funnel electrospray source	Flow rate: 1 µL/min – 1 mL/min
Mass Range	20 – 40,000 m/z
Quadrupole isolation	Up to 3,000 m/z
Quadrupole Mass Range	Up to 40,000 m/z
Mass accuracy in MS and MS/MS	With internal calibrant: better than 600ppb RMS Error With external calibrant: better than 2 ppm RMS Error
Calibration	ONE calibration valid for MS and MS/MS analysis. Calibration is independent from charge state of calibrant mass
Mass resolution	75,000 FSR (full sensitivity resolution)
Isotopic pattern	The true isotopic pattern is maintained due to TIP™ technology (True Isotopic Pattern) and allows three dimensional chemical characterizations of analytes via SmartFormula™3D algorithm using exact mass, TIP, and MS/MS fragment data.
SmartFormula™3D	Enables unambiguous formula determination at “sub-ppm” confidence level up to 1000 Da.

Mass stability & dynamic range	hrXIC (high resolution Extracted Ion Chromatogram) technology with better than +/- 0.5 - 1.0 mDa stability on centroid data values over an typical LC peak.
Full scan sensitivity in MS	ESI: Reserpine 1 pg S/N>100:1 RMS With Ion-Booster (optional): Reserpine 100 fg S/N>100:1 RMS
Full scan sensitivity in MS/MS	The signal height obtained from a consumption of 2.5 fmol of Glu-Fibrinopeptide B will be better than 1000 counts on the most intense y' sequence ion from the MS/MS spectrum of the doubly charged precursor ion. This shall correspond to a signal to noise ratio better than 50:1. The MS/MS sensitivity specification is met while using quadrupole isolation of the precursor ion demonstrating that there is minimal transmission loss through the isolating quad. A solution of 100 fmol/μL Glu-Fibrinopeptide B shall be introduced at a flow rate of 3 μL/min.
TOF repetition rate	Up to 20 kHz
Temperature compensation	Yes
Digitizer	5Gsample/sec ADC with 50 Gbit/sec
Dynamic range	10 bit ADC for high quantitative dynamic range
Acquisition rate	up to 50 Hz MS 50 Hz MS/MS (profile and peak detected spectra to disk)

Optional accessory

IonBooster	Optional ion source
APCI II	Optional ion source
APPI II	Optional ion source
GC-APCI	Allows for direct GC coupling (Optional ion source)
APLI	Optional ion source
CryoSpray	Optional ion source
Bruker CaptiveSpray nanoBooster	Optional ion source
On-/Off-Line Nanospray	Optional ion source
CE/MS interface	With grounded needle for easy CE-TOF set-up (Optional)
DIP	Direct Probe (Optional ion source)



maxis impact

- Maximum Impact – Definitive Answers

Redefining Accurate Mass LC-MS/MS



Until now, mass spectrometry technologies have forced scientists to choose between performance characteristics for a given application. Often, a system is designed for qualitative work, sacrificing quantitation performance in the process.

There is no need to make compromises in mass spectrometry anymore. The maXis impact™ sets a new technology standard where industry leading performance values are all simultaneously available in a single acquisition at full sensitivity.

Powered by a series of patented technology innovations, the maXis impact simply provides the very best results without compromise in a cost effective, benchtop format.


maXis impact: Qualitative and Quantitative results from a single LC/MS analysis

- Sub-ppm mass accuracy and 40,000 full sensitivity resolution
- Extreme sensitivity across the entire mass range from very small molecules to intact proteins
- High speed MS/MS capability with full U-HPLC compatibility
- Simultaneous analysis of major and trace sample components
- Isotopic fidelity for definitive molecular formulae determination
- Robust and easy to operate system

And all in an economical and compact design!



● Unparalleled Performance Now in a Bench-top Package



The only no-compromise full sensitivity and full resolution UHR-TOF on the market

A game changing step forwards in TOF technology featuring nine new patents. The maXis impact is powered by a number of outstanding technology innovations including:

- Unique Full Sensitivity Resolution (FSR) technology in all modes
- Proven Bruker UHR-TOF technology in bench-top format
- New broadband transmission CID cell providing ground-breaking sensitivity
- Fast 4 G sample/sec digitizer enabling data acquisition at up to 50 Hz
- Long life FLASH Detector
- Flexible choice of API sources including new CaptiveSpray™ and solid sample (DIP) and GC inlets

Refining capabilities across the analytical world



With its unique capabilities, the maXis impact liberates your work from restriction and compromise. Whether the need is for small or large molecule analysis, quantitation or molecular identification, the maXis impact is more than capable of providing the results needed to make decisions first time, every time.

Let the maXis impact achieve your goals in:

- Forensics, toxicology & doping control
- Food & environmental testing
- Synthetic chemistry support
- Drug metabolite and impurity identification and quantitation
- Metabolomics
- Intact protein analysis & characterization of biopharmaceuticals
- Biomarker discovery & validation







One-Shot Certainty for Qual/Quant Analysis

maXis impact captures the full picture of your sample in one analysis alleviating the need for multiple runs. More than 4 orders of fundamental dynamic range achieved without artificial compensators allows simultaneous measurement of major and trace sample components in a single analysis - even if coeluent.



MS of Buspiron acquired at 1Hz up to 50 Hz spectra to disk.
Resolution and spectral accuracy maintained even at maximum speed

Qualitative and quantitative results from a single auto MS/MS run



1. TIC from complex sample (60 pesticides spiked into ginger extract (QuEChERS)).
2. High resolution extracted ion chromatogram (hrXIC) of m/z 223.0745 with 5mDa discrimination width (Acetamidrid).
3. Magnified LC peak. High speed acquisition rate in MS and MS/MS mode (MS: 100ms, MS/MS: 20ms spectra time). 100 data points across this LC peak allows precise quantitation.
4. High resolution survey MS with m/z 223.0745 selected as precursor.
5. High spectral accuracy means isotopic fidelity for identification certainty.
6. High sensitivity MS/MS spectrum even at very low MW (56.0498 m/z) for confirmation by fragment information.

● Precisely Targeted Small Molecule Analysis




Ideally suited for broad band multi-target screening

Resolve and identify multiple compounds from extremely complex mixtures across a wide mass range even at low analyte concentrations. With these key performance attributes, the maXis impact is an ideal fit for rapid, high target number screening tasks such as food testing or forensic analysis.


For the laboratory challenged with rapid response to novel analytes, retrospective *in-silico* screening of complete datasets using TargetAnalysis™ software makes maXis impact the most agile screening tool available.

Multi-target screening at 2mDa hrXIC selectivity



2 mDa high resolution XIC from 200 pesticides. Each color represents the hrXIC of one pesticide. Sample at a concentration of 50 ppb each separated with U- HPLC maXis impact.

Triple quad sensitivity and quantitation with high-res, accurate mass



Fluoxastrobin (m/z 459) has been quantified over 4 orders of magnitude in concentration from 1ng down to 50fg. The linearity is $R=0.9996$ with an ultra low LOD coupled with accurate mass certainty.

The maXis impact encapsulates the capabilities of low fg range sensitivity with more than four orders of magnitude dynamic range and quantitative capacity along with the performance of an exceptional high resolution accurate mass instrument. The simultaneous combination of attributes is ideal for many discovery, development, screening and drug metabolism applications which often require both the identification and quantitation of sample components that vary widely in both structure and concentration.

Excellent for use in the analysis of both known and unknown compounds, the maXis impact represents a new level of achievement in the challenge to provide one step qualitative and quantitative data in a single rapid analysis.



Small Molecule ID & Structure Elucidation


Impacting Metabolomics results

A key bottleneck in metabolomics research is structural confirmation and elucidation of metabolites. Diterpene glycosides (DTGs) are abundant plant defence compounds with largely unknown modes of action. Outstanding MS and MS/MS data for *Nicotiana attenuata* (tobacco) DTGs acquired on a maXis impact enabled the rapid identification of the entire compound family. To achieve this, the fragmentation results are visualized and interpreted using Bruker's novel FragmentExplorer™.


The FragmentExplorer has been especially designed for faster interpretation of MS/MS data. It provides an interactive link between SmartFormula 3D™ results, mass spectra and molecular structures.

SmartFormula3D™ co-developed with Pfizer, utilizes maXis impact high isotopic fidelity in MS and MS/MS in a unique relational algorithm to assign molecular formulae of unknown analytes without prior knowledge or assumptions.

Compound identification and structure confirmation




SmartFormula3D




FragmentExplorer with embedded ChemDraw expertly assists fragment assignments.





Biomarker Discovery and Validation


In-depth qualitative analysis of complex proteomics samples



LC auto MS/MS survey. The human cell line HT29 is an established model system to study colon cancer progression. A 1 µg digest of HT29 cells was separated on a RSLCnano system (3 h gradient, 25 cm column) and measured on a maXis impact with a captive spray source.


Get full performance first time, every time

The CaptiveSpray is a revolutionary LC-MS source that combines the sensitivity of nanospray with the ease of use and robustness of electrospray. CaptiveSpray utilizes patented technology to capture and sweep sample ions into the MS independent of LC flow rate.




CaptiveSpray LC-MS source

Identification over a dynamic range of 4 decades of concentration from complex samples



Sigma UPS-2 standard sample


Real-life MS mass accuracy



Plot of all peptides from a digest of 200ng E.coli

BioPharma QC

Isotopic resolution of intact proteins



Quantum leap in data quality for biopharmaceutical characterization

With 30k-40k resolution for intact proteins @ full sensitivity with true isotopic patterns, the maXis impact allows you to simultaneously characterize, impurity detect and quantify your biopharmaceutical products. Combined with MS& MS/MS accuracy < 1ppm the maxis impact ensures certainty first time every time, allowing you to report your results with confidence.

Source Options

Wide choice of ionization and coupling techniques for a broad range of sampling including insoluble compounds:

- GC-MS coupling
- APCI II (atmospheric pressure chemical ionization) source with direct probe
- APPI II (atmospheric pressure photo ionization) source
- CaptiveSpray LC-MS source
- CE-MS coupling with grounded needle

LC Options

Bruker fully supports and integrates a wide range of leading HPLC systems, autosamplers and accessories.



GC-MS coupling



Solid probe mounted on an APCI II source



APPI II source



For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

www.bruker.com/maXis



Description

maXis impact bench-top UHR Time-of-Flight Mass Spectrometer System

BDAL #282000

Easy-to-use, Ultra-High-Resolution electrospray ionization quadrupole time-of-flight LC/MS/MS mass spectrometer designed for exact mass and true isotopic measurements in both MS and MS/MS mode.

Bench-top, small footprint Mass Spectrometer system for exact mass and highest mass resolution at U-HPLC speed in both MS and MS/MS mode:

- Unique FSR technology with Full Sensitivity @ Maximum Resolution achieved without any time constraints, in MS and MS/MS mode
- Outstanding Mass Resolution and Accuracy in both MS and MS/MS
- High-resolution extracted ion chromatograms capabilities
- High performance hyperbolic quadrupole and collision cell for compound fragmentation
- True Isotopic Measurements

A. Apollo II (ESI) Source

- Highly sensitive ESI Source with proprietary ion funnel guide for gentle mass independent ion focusing and high ion transmission efficiency
- Combined Funnel-Octopole-Cartridge with front access for easy maintenance
- Grounded needle for safety and easy sample introduction
- Suitable for U-HPLC, HPLC and CE coupling
- Heated counter current drying gas for gentle and efficient drying
- Ion lens system including in Source collision induced dissociation control (IS-CID)
- Pneumatic off-axis nebulization for flow rates up to 1 ml/min., with gradients from 100% aqueous to 100% organic
- Flow rate 1µl/min to 1ml/min
- Ni coated glass capillary for physical and electrostatic isolation
- Computer controlled HV and gas controller

B. High Mass Quadrupole Mass Filter:

- Hyperbolic quadrupole mass filter
- Ultra stable monolithic design
- RF-generator for monoisotopic precursor ion selection

C. Novel high-transmission CID Collision cell:

- Hexapole broad-mass bandwidth design
- Fast radial ion ejection enabling fast MS/MS cycle
- RF-Generator with fast amplitude switching
- Collision gas controller

D. Orthogonal pulsed ion extraction and UHR Time-of-Flight Mass Analyzer

- Interface housing and ion lens system
- In-line detector system for easy maintenance



- Ultra-stable high voltage switches with up to 20 kHz repetition rate and appropriate power supplies.
- UHR-TOF analyzer with orthogonal mounted ion source
- Novel in-flight refocusing optics for uncompromised sensitivity
- Dual stage ion reflectron with increased mass resolution and accuracy
- High-sensitivity and fast ion detector system, mechanical adjustment in micrometer range
- Positive and negative ion modes
- Ultra-stable high voltage power supplies for TOF analyzer and detector

E. Vacuum system

- Analyzer vacuum housing
- Vacuum system with 5 differential pumping stages
- One roughing pump and quadruple stage turbo-drag pumps for ESI source and UHR-TOF analyzer
- Vacuum measurement and pump control unit

F. Syringe pump

G. Modes of Operation

- TOF Mass ranges 20-40,000 m/z
- Internal calibration (MS and MS/MS)
- External calibration (MS and MS/MS)
- Exact mass measurements independent from sample concentration over a wide dynamic range without second sprayer.

H. High-resolution-performance and accurate mass features

- Patented ion funnel source
- One-shot acquisition mode, no tuning for mass range optima
- Enhanced low mass sensitivity
- Superior MS/MS sensitivity
- Ultra broad mass-bandwidth
- Long term and ultra stable mass axis stability in MS and MS/MS
- Exact mass independent from sample concentration charge state and collision energy
- Combined calibration for both MS and MS/MS
- Wide dynamic range for quantitation
- Advanced temperature compensated flight tube
- Positive / negative ion operation

I. Data system:

- PC Workstation with 2,66 GHz Dual-Pentium Processor, 4 GB RAM, system hard-disk drive plus 1000 GB hard-disk drive for data, 1.44 Mb 3.5" floppy diskette drive
- DVD-ROM drive
- R/W DVD-ROM drive
- 20" flat screen colour monitor
- WinXP OS
- Laser printer



J. Applications software

Software package Compass 1.4 for HPLC and MS control, data acquisition, post processing, and data analysis:

- Operating system Windows XP
- Compass / HyStar 3.2 for integrated control of most popular U-HPLC and HPLC systems and auto samplers and automation
- Expert mode: extended control over instrument parameters for interactive system optimization of sophisticated exact mass methods
- Compass / Data Analysis software DA 4.1, including:
 - Advanced data processing with a high degree of automation
 - SmartFormula 3D™: Automated sum formula determination using MS and MS/MS data with both, accurate mass and isotopic fit.
 - QuantAnalysis™ quantitation package
 - LibrarySearch™ module for search of MS, MS/MS and MS_n spectra with advanced matching algorithm
 - Charge deconvolution module
 - Export of spectra and ion current profiles as Windows metafiles to word-processing programs
- SW License Compass 1.3
- SW License Charge Deconvolution for DA 3.4
- SW License MaxEntropy Deconvolution as an option

K. Set of manuals and reference CD-ROMs

L. Installation

M. Familiarization upon installation

N. 1 year warranty

Voucher for a factory-training course - valid for 2 participants.



Abiotrophia defectiva
Acetobacter aceti
Acetobacter pasteurianus
Acholeplasma laidlawii
Achromobacter denitrificans
Achromobacter insolitus
Achromobacter piechaudii
Achromobacter ruhlandii
Achromobacter sp
Achromobacter spanius
Achromobacter xylooxidans
Acidaminococcus fermentans
Acidaminococcus intestinalis
Acidiphilium acidophilum
Acidovorax avenae
Acidovorax defluvii
Acidovorax delafieldii
Acidovorax facilis
Acidovorax konjaci
Acidovorax temperans
Acinetobacter baumannii
Acinetobacter baylyi
Acinetobacter bouvetii
Acinetobacter calcoaceticus

Acinetobacter gerneri
Acinetobacter guillouiae
Acinetobacter haemolyticus
Acinetobacter johnsonii
Acinetobacter junii
Acinetobacter lwoffii
Acinetobacter nosocomialis
Acinetobacter parvus
Acinetobacter pittii
Acinetobacter radioresistens
Acinetobacter schindleri
Acinetobacter sp
Acinetobacter tandoii
Acinetobacter tjernbergiae
Acinetobacter townneri
Acinetobacter ursingii
Actinobacillus delphinicola
Actinobacillus equuli
Actinobacillus lignieresii
Actinobacillus pleuropneumoniae
Actinobacillus rossii
Actinobacillus suis
Actinobacillus ureae
Actinobaculum schaalii

Actinobaculum suis
Actinobaculum urinale
Actinocorallia libanotica
Actinomyces bovis
Actinomyces bowdenii
Actinomyces canis
Actinomyces cardiffensis
Actinomyces catuli
Actinomyces coleocanis
Actinomyces dentalis
Actinomyces denticolens
Actinomyces europaeus
Actinomyces funkei
Actinomyces gerencseriae
Actinomyces graevenitzi
Actinomyces hordeovulneris
Actinomyces hyovaginalis
Actinomyces israelii
Actinomyces marimammalium
Actinomyces meyeri
Actinomyces naeslundii
Actinomyces nasicola
Actinomyces neuui
Actinomyces odontolyticus



Actinomyces oris
Actinomyces radidentis
Actinomyces radingae
Actinomyces ruminicola
Actinomyces sp
Actinomyces suimastitidis
Actinomyces turicensis
Actinomyces urogenitalis
Actinomyces vaccimaxillae
Actinomyces viscosus
Adlercreutzia equolifaciens
Aerococcus christensenii
Aerococcus sanguinicola
Aerococcus urinae
Aerococcus urinaehominis
Aerococcus viridans
Aeromicrobium flavum
Aeromonas bestiarum
Aeromonas caviae
Aeromonas encheleia
Aeromonas enteropelogenes
Aeromonas eucrenophila
Aeromonas hydrophila
Aeromonas ichthiosmia

Aeromonas jandaei
Aeromonas media
Aeromonas molluscorum
Aeromonas popoffii
Aeromonas punctata
Aeromonas salmonicida
Aeromonas schubertii
Aeromonas simiae
Aeromonas sobria
Aeromonas sp[2]
Aeromonas veronii
Afipia broomeae
Afipia felis
Afipia massiliensis
Aggregatibacter actinomycetemcomitans
Aggregatibacter aphrophilus
Aggregatibacter segnis
Agrococcus jenensis
Agromyces brachium
Agromyces cerinus
Agromyces fucosus
Agromyces hippuratus
Agromyces humatus
Agromyces italicus

Agromyces lapidis
Agromyces mediolanus
Agromyces neolithicus
Agromyces rhizosphaerae
Agromyces salentinus
Agromyces subbeticus
Alcaligenes faecalis
Alcaligenes sp
Alicyclobacillus acidocaldarius
Alicyclobacillus acidoterrestris
Alicyclobacillus cycloheptanicus
Aliivibrio fischeri
Alishewanella fetalis
Alistipes finegoldii
Alistipes onderdonkii
Alistipes shahii
Alkalibacillus haloalkaliphilus
Alloiococcus otitis
Alloscardovia omnicoles
Alternaria alternata
Amycolatopsis alba
Amycolatopsis azurea
Amycolatopsis balhimycina
Amycolatopsis coloradensis



Amycolatopsis fastidiosa
Amycolatopsis japonica
Amycolatopsis kentuckyensis
Amycolatopsis keratiniphila
Amycolatopsis lexingtonensis
Amycolatopsis lurida
Amycolatopsis mediterranei
Amycolatopsis minnesotensis
Amycolatopsis pretoriensis
Amycolatopsis sulphurea
Amycolatopsis tolypomycina
Anaerococcus hydrogenalis
Anaerococcus lactolyticus
Anaerococcus murdochii
Anaerococcus octavius
Anaerococcus prevotii
Anaerococcus sp
Anaerococcus tetradius
Anaerococcus vaginalis
Anaerotruncus colihominis
Aneurinibacillus aneurinilyticus
Aneurinibacillus migulanus
Aquicola tertiaricarbonis
Arcanobacterium haemolyticum

Arcanobacterium hippocoleae
Arcanobacterium phocae
Arcanobacterium pluranimalium
Arcobacter butzleri
Arcobacter cibarius
Arcobacter cryaerophilus
Arcobacter halophilus
Arcobacter nitrofigilis
Arcobacter skirrowii
Aromatoleum alkani
Aromatoleum anaerobicus
Aromatoleum aromaticum
Aromatoleum bremensis
Aromatoleum buckelii
Aromatoleum diolicum
Aromatoleum evansii
Aromatoleum pretroleum
Aromatoleum terpenicum
Aromatoleum tolulyticus
Aromatoleum toluolicum
Aromatoleum toluvorans
Arsenophonus nasoniae
Arthroascus schoenii
Arthrobacter ardeleyensis

Arthrobacter arilaitensis
Arthrobacter aurescens
Arthrobacter bergerei
Arthrobacter castelli
Arthrobacter chlorophenolicus
Arthrobacter citreus
Arthrobacter creatinolyticus
Arthrobacter crystallopoietes
Arthrobacter cumminsii
Arthrobacter flavus
Arthrobacter gandavensis
Arthrobacter gangotriensis
Arthrobacter globiformis
Arthrobacter histidinovorans
Arthrobacter ilicis
Arthrobacter kerguelensis
Arthrobacter koreensis
Arthrobacter luteolus
Arthrobacter monumenti
Arthrobacter mysorens
Arthrobacter nasiphocae
Arthrobacter nicotianae
Arthrobacter nicotinovorans
Arthrobacter nitroguajacolicus



Arthrobacter oxydans
Arthrobacter parietis
Arthrobacter pascens
Arthrobacter pigmenti
Arthrobacter polychromogenes
Arthrobacter protophormiae
Arthrobacter psychrolactophilus
Arthrobacter psychrophenicus
Arthrobacter pyridinolis
Arthrobacter ramosus
Arthrobacter roseus
Arthrobacter russicus
Arthrobacter scleromae
Arthrobacter sp
Arthrobacter stackebrandtii
Arthrobacter sulfonivorans
Arthrobacter sulfureus
Arthrobacter tecti
Arthrobacter tumbae
Arthrobacter uratoxydans
Arthrobacter ureafaciens
Arthrobacter woluwensis
Arthroderma benhamiae
Arthrographis_kalrae[ana]
(Eremomyces_langeronii[teleo])

Aspergillus brasiliensis
Aspergillus flavus
Aspergillus fumigatus
Aspergillus niger
Aspergillus terreus
Aspergillus versicolor
Aspergillus_thermomutatus[ana]
(Neosartorya_pseudofischeri[teleo])
Atopobium minutum
Atopobium parvulum
Atopobium rimae
Atopobium sp
Atopobium vaginae
Aureimonas altamirensis
Aureobasidium pullulans
Austwickia chelonae
Avibacterium avium
Avibacterium endocarditidis
Avibacterium gallinarum
Avibacterium volantium
Azoarcus communis
Azoarcus indigenus
Azoarcus sp
Azohydromonas lata
Bacillus acidicola

Bacillus agaradhaerens
Bacillus akibai
Bacillus alcalophilus
Bacillus algicola
Bacillus altitudinis
Bacillus alveayuensis
Bacillus amyloliquefaciens
Bacillus aquimaris
Bacillus arsenicus
Bacillus asahii
Bacillus atrophaeus
Bacillus azotoformans
Bacillus badius
Bacillus barbaricus
Bacillus bataviensis
Bacillus benzoovorans
Bacillus carboniphilus
Bacillus cellulosilyticus
Bacillus cereus
Bacillus chagannorensis
Bacillus cibi
Bacillus circulans
Bacillus clarkii
Bacillus clausii



Bacillus coagulans	Bacillus infantis	Bacillus patagoniensis
Bacillus cohnii	Bacillus jeotgali	Bacillus pseudocaliphilus
Bacillus decolorationis	Bacillus koreensis	Bacillus pseudofirmus
Bacillus drentensis	Bacillus krulwichiae	Bacillus pseudomycooides
Bacillus endophyticus	Bacillus lentus	Bacillus psychrosaccharolyticus
Bacillus farraginis	Bacillus licheniformis	Bacillus pumilus
Bacillus fastidiosus	Bacillus litoralis	Bacillus ruris
Bacillus firmus	Bacillus luciferensis	Bacillus safensis
Bacillus flexus	Bacillus macauensis	Bacillus salarius
Bacillus fordii	Bacillus mannanyliticus	Bacillus seohaeanensis
Bacillus fortis	Bacillus marisflavi	Bacillus shackletonii
Bacillus funiculus	Bacillus massiliensis	Bacillus simplex
Bacillus galactosidilyticus	Bacillus megaterium	Bacillus siralis
Bacillus gibsonii	Bacillus mojavensis	Bacillus smithii
Bacillus halmapalus	Bacillus muralis	Bacillus soli
Bacillus halodurans	Bacillus mycooides	Bacillus sonorensis
Bacillus hemicellulosilyticus	Bacillus nealsonii	Bacillus sp
Bacillus horikoshii	Bacillus niacini	Bacillus sporothermodurans
Bacillus horneckiae	Bacillus novalis	Bacillus subterraneus
Bacillus horti	Bacillus odysseyi	Bacillus subtilis
Bacillus humi	Bacillus okhensis	Bacillus thermoamylovorans
Bacillus hwajinpoensis	Bacillus okuhidensis	Bacillus thioparans
Bacillus idriensis	Bacillus oleronius	Bacillus thuringiensis
Bacillus indicus	Bacillus oshimensis	Bacillus vallismortis



Bacillus vedderi
Bacillus vietnamensis
Bacillus vireti
Bacillus wakoensis
Bacillus weihenstephanensis
Bacteroides caccae
Bacteroides coagulans
Bacteroides eggerthii
Bacteroides finegoldii
Bacteroides fragilis
Bacteroides gallinarum
Bacteroides intestinalis
Bacteroides massiliensis
Bacteroides nordii
Bacteroides ovatus
Bacteroides pyogenes
Bacteroides salyersiae
Bacteroides stercoris
Bacteroides thetaiotaomicron
Bacteroides uniformis
Bacteroides vulgatus
Balneatrix alpica
Bartonella japonica
Beauveria bassiana

Bergeyella zoohelcum
Bibersteinia trehalosi
Bifidobacterium adolescentis
Bifidobacterium angulatum
Bifidobacterium animalis
Bifidobacterium asteroides
Bifidobacterium bifidum
Bifidobacterium boum
Bifidobacterium breve
Bifidobacterium catenulatum
Bifidobacterium choerinum
Bifidobacterium coryneforme
Bifidobacterium dentium
Bifidobacterium gallicum
Bifidobacterium gallinarum
Bifidobacterium longum
Bifidobacterium magnum
Bifidobacterium merycicum
Bifidobacterium minimum
Bifidobacterium pseudocatenulatum
Bifidobacterium pseudolongum
Bifidobacterium pullorum
Bifidobacterium ruminantium
Bifidobacterium saeculare

Bifidobacterium scardovii
Bifidobacterium thermacidophilum
Bifidobacterium thermophilum
Bilophila sp
Blastomonas natatoria
Blastomonas ursincola
Blautia coccoides
Bordetella avium
Bordetella bronchiseptica
Bordetella hinzii
Bordetella holmesii
Bordetella parapertussis
Bordetella pertussis
Bordetella petrii
Bordetella sp
Bordetella trematum
Borrelia burgdorferi
Borrelia garinii
Borrelia spielmanii
Brachybacterium faecium
Brachybacterium muris
Brachyspira murdochii
Brachyspira pilosicoli
Bradyrhizobium betae



Bradyrhizobium denitrificans

Brenneria alni

Brenneria nigrifluens

Brenneria quercina

Brenneria rubrifaciens

Brenneria salicis

Brevibacillus agri

Brevibacillus borstelensis

Brevibacillus brevis

Brevibacillus centrosporus

Brevibacillus choshinensis

Brevibacillus formosus

Brevibacillus invocatus

Brevibacillus laterosporus

Brevibacillus parabrevis

Brevibacillus reuszeri

Brevibacterium aurantiacum

Brevibacterium casei

Brevibacterium celere

Brevibacterium iodinum

Brevibacterium linens

Brevibacterium marinum

Brevibacterium paucivorans

Brevibacterium picturae

Brevibacterium ravensturgense

Brevibacterium sanguinis

Brevundimonas aurantiaca

Brevundimonas diminuta

Brevundimonas intermedia

Brevundimonas nasdae

Brevundimonas sp

Brevundimonas subvibrioides

Brevundimonas vesicularis

Brochothrix thermosphacta

Budvicia aquatica

Bulleidia extracta

Burkholderia ambifaria

Burkholderia andropogonis

Burkholderia anthina

Burkholderia caledonica

Burkholderia caribensis

Burkholderia cenocepacia

Burkholderia cepacia

Burkholderia diffusa

Burkholderia dolosa

Burkholderia fungorum

Burkholderia gladioli

Burkholderia glathei

Burkholderia glumae

Burkholderia lata

Burkholderia latens

Burkholderia metallica

Burkholderia multivorans

Burkholderia phenazinium

Burkholderia phymatum

Burkholderia plantarii

Burkholderia pyrrocinia

Burkholderia sacchari

Burkholderia seminalis

Burkholderia stabilis

Burkholderia terricola

Burkholderia thailandensis

Burkholderia tropica

Burkholderia tuberum

Burkholderia vietnamiensis

Burkholderia xenovorans

Buttiauxella agrestis

Buttiauxella brennerae

Buttiauxella ferragutiae

Buttiauxella gaviniae

Buttiauxella izardii

Buttiauxella noackiae



Buttiauxella warmboldiae
Butyricimonas virosa
Campylobacter avium
Campylobacter canadensis
Campylobacter coli
Campylobacter concisus
Campylobacter curvus
Campylobacter fetus
Campylobacter gracilis
Campylobacter helveticus
Campylobacter hominis
Campylobacter hyointestinalis
Campylobacter jejuni
Campylobacter lanienae
Campylobacter lari
Campylobacter peloridis
Campylobacter rectus
Campylobacter showae
Campylobacter sputorum
Campylobacter upsaliensis
Campylobacter ureolyticus
Candida alai
Candida albicans
Candida allociferrii

Candida ambrosiae
Candida auris
Candida blattae
Candida bohiensis
Candida boidinii
Candida bracarensis
Candida buenavistaensis
Candida carpophila
Candida castelli
Candida catenulata
Candida cylindracea
Candida dubliniensis
Candida duobushaemulonii
Candida ernobii
Candida freyschussii
Candida friedrichii
Candida frijolesensis
Candida glabrata
Candida haemulonii
Candida inconspicua
Candida infanticola
Candida intermedia
Candida labiduridarum
Candida lactiscondensi

Candida magnoliae
Candida maltosa
Candida membranifaciens
Candida mesenterica
Candida metapsilosis
Candida multigemmis
Candida nemodendra
Candida nitratophila
Candida nivariensis
Candida norvegica
Candida orthopsilosis
Candida palmioleophila
Candida parapsilosis
Candida pararugosa
Candida peltata
Candida pini
Candida pseudohaemulonii
Candida rugosa
Candida saitoana
Candida sake
Candida shehatae
Candida sojae
Candida solani
Candida spandovensis



Candida succiphila
Candida tropicalis
Candida versatilis
Candida vini
Candida viswanathii
Candida zeylanoides
Candida_cacaoi[ana] (Pichia_farinosa[teleo])
Candida_chodatii[ana]
(Hyphopichia_burtonii[teleo])
Candida_ciferrii[ana]
(Stephanoascus_ciferrii[teleo])
Candida_citrea[ana] (Pichia_nakasei[teleo])
Candida_colliculosa[ana]
(Torulaspora_delbrueckii[teleo])
Candida_fabianii[ana] (Pichia_fabianii[teleo])
Candida_guilliermondii[ana]
(Pichia_guilliermondii[teleo])
Candida_guilliermondii_var_membranaefaciens[ana] (Pichia_ohmeri[teleo])
Candida_holmii[ana]
(Kazachstania_exigua[teleo])
Candida_kefyr[ana]
(Kluyveromyces_marxianus[teleo])
Candida_krusei[ana]
(Issatchenkia_orientalis[teleo])
Candida_lambica[ana]
(Pichia_fermentans[teleo])
Candida_lipolytica[ana]
(Yarrowia_lipolytica[teleo])
Candida_lusitaniae[ana]
(Clavispora_lusitaniae[teleo])
Candida_norvegensis[ana]
(Pichia_norvegensis[teleo])

Candida_pelliculosa[ana]
(Pichia_anomala[teleo])
Candida_pintolopesii[ana]
(Kazachstania_pintolopesii[teleo])
Candida_pulcherrima[ana]
(Metschnikowia_pulcherrima[teleo])
Candida_reukaufii[ana]
(Metschnikowia_reukaufii[teleo])
Candida_slooffiae[ana]
(Kazachstania_slooffiae[teleo])
Candida_sorbosa[ana]
(Issatchenkia_occidentalis[teleo])
Candida_sphaerica[ana]
(Kluyveromyces_lactis[teleo])
Candida_thermophila[ana]
(Ogataea_thermophila[teleo])
Candida_utilis[ana] (Pichia_jandinii[teleo])
Candida_valida[ana]
(Pichia_membranifaciens[teleo])#
Candidatus Reyranelia massiliensis
Capnocytophaga canimorsus
Capnocytophaga cynodegmi
Capnocytophaga gingivalis
Capnocytophaga granulosa
Capnocytophaga haemolytica
Capnocytophaga ochracea
Capnocytophaga sp
Capnocytophaga sputigena
Cardiobacterium hominis
Cardiobacterium sp
Cardiobacterium valvarum

Carnobacterium maltaromaticum
Castellaniella defragrans
Caulobacter sp[2]
Caulobacter vibrioides
Cedecea davisae
Cedecea lapagei
Cedecea neteri
Cellulomonas fimi
Cellulomonas flavigena
Cellulomonas gelida
Cellulomonas uda
Cellulosimicrobium cellulans
Chaetomium globosum
Chromobacterium subtsugae
Chromobacterium violaceum
Chromohalobacter salexigenus
Chryseobacterium ginsenosidimutans
Chryseobacterium gleum
Chryseobacterium hagamense
Chryseobacterium hominis
Chryseobacterium indologenes
Chryseobacterium joosteii
Chryseobacterium oranimense
Chryseobacterium scophthalmum



Chryseobacterium sp	Clostridium butyricum	Clostridium halophilum
Citrobacter amalonaticus	Clostridium cadaveris	Clostridium hathewayi
Citrobacter braakii	Clostridium carboxidivorans	Clostridium hiranonis
Citrobacter farmeri	Clostridium carnis	Clostridium histolyticum
Citrobacter freundii	Clostridium celerecrescens	Clostridium homopropionicum
Citrobacter gillenii	Clostridium cellobioparum	Clostridium hylemonae
Citrobacter koseri	Clostridium chauvoei	Clostridium indolis
Citrobacter murliniae	Clostridium citroniae	Clostridium innocuum
Citrobacter rodentium	Clostridium clostridioforme	Clostridium intestinale
Citrobacter sedlakii	Clostridium cochlearium	Clostridium irregulare
Citrobacter youngae	Clostridium colicanis	Clostridium isatidis
Clavibacter michiganensis	Clostridium colinum	Clostridium jejuense
Clostridium acetobutylicum	Clostridium collagenovorans	Clostridium lactatifermentans
Clostridium aciditolerans	Clostridium difficile	Clostridium lentocellum
Clostridium aerotolerans	Clostridium diolis	Clostridium limosum
Clostridium aldenense	Clostridium disporicum	Clostridium lundense
Clostridium algidicarnis	Clostridium drakei	Clostridium magnum
Clostridium algidixylanolyticum	Clostridium fallax	Clostridium malenominatum
Clostridium aminophilum	Clostridium formicaceticum	Clostridium mayombei
Clostridium baratii	Clostridium frigoris	Clostridium novyi
Clostridium bartlettii	Clostridium ghonii	Clostridium papyrosolvens
Clostridium beijerinckii	Clostridium glycolicum	Clostridium paraputrificum
Clostridium bifermentans	Clostridium glycyrrhizinilyticum	Clostridium perfringens
Clostridium bolteae	Clostridium haemolyticum	Clostridium phytofermentans



Clostridium propionicum
Clostridium ramosum
Clostridium saccharobutylicum
Clostridium saccharogumia
Clostridium sardiniense
Clostridium sartagoforme
Clostridium schirmacherense
Clostridium scindens
Clostridium septicum
Clostridium sordellii
Clostridium sp
Clostridium sphenoides
Clostridium spiroforme
Clostridium sporogenes
Clostridium sporosphaeroides
Clostridium subterminale
Clostridium symbiosum
Clostridium tertium
Clostridium tetani
Clostridium thermopalmarium
Clostridium tunisiense
Clostridium xylanovorans
Cohnella fontinalis
Cohnella hongkongensis

Colletotrichum gloeosporioides
Collinsella aerofaciens
Comamonas aquatica
Comamonas kerstersii
Comamonas nitratorans
Comamonas terrigena
Comamonas testosteroni
Coprobacillus cateniformis
Corynebacterium accolens
Corynebacterium afermentans
Corynebacterium ammoniagenes
Corynebacterium amycolatum
Corynebacterium appendicis
Corynebacterium aquilae
Corynebacterium argenteratense
Corynebacterium aurimucosum
Corynebacterium auris
Corynebacterium auriscanis
Corynebacterium bovis
Corynebacterium callunae
Corynebacterium camporealensis
Corynebacterium capitovis
Corynebacterium casei
Corynebacterium ciconiae

Corynebacterium confusum
Corynebacterium coyleae
Corynebacterium cystitidis
Corynebacterium diphtheriae
Corynebacterium durum
Corynebacterium efficiens
Corynebacterium falsenii
Corynebacterium felinum
Corynebacterium flavescens
Corynebacterium freneyi
Corynebacterium glaucum
Corynebacterium glucuronolyticum
Corynebacterium glutamicum
Corynebacterium halotolerans
Corynebacterium hansenii
Corynebacterium imitans
Corynebacterium jeikeium
Corynebacterium kroppenstedtii
Corynebacterium kutscheri
Corynebacterium lipophile_group_F1
Corynebacterium lipophiloflavum
Corynebacterium macginleyi
Corynebacterium mastitidis
Corynebacterium matruchotii



Corynebacterium minutissimum
Corynebacterium mucifaciens
Corynebacterium mycetoides
Corynebacterium phocae
Corynebacterium pilosum
Corynebacterium propinquum
Corynebacterium pseudodiphtheriticum
Corynebacterium pseudotuberculosis
Corynebacterium renale
Corynebacterium resistens
Corynebacterium riegellii
Corynebacterium simulans
Corynebacterium singulare
Corynebacterium sp
Corynebacterium sphenisci
Corynebacterium spheniscorum
Corynebacterium stationis
Corynebacterium striatum
Corynebacterium suicordis
Corynebacterium sundsvallense
Corynebacterium terpenotabidum
Corynebacterium testudinoris
Corynebacterium thomssenii
Corynebacterium tuberculostearicum

Corynebacterium tuscaniense
Corynebacterium ulcerans
Corynebacterium urealyticum
Corynebacterium ureicelerivorans
Corynebacterium variabile
Corynebacterium vitaeruminis
Corynebacterium xerosis
Cronobacter sakazakii
Cryptococcus albidosimilis
Cryptococcus curvatus
Cryptococcus diffluens
Cryptococcus flavescens
Cryptococcus flavus
Cryptococcus gastricus
Cryptococcus humicola
Cryptococcus laurentii
Cryptococcus liquefaciens
Cryptococcus macerans
Cryptococcus magnus
Cryptococcus neoformans
Cryptococcus saitoi
Cryptococcus terreus
Cryptococcus uzbekistanensis
Cryptococcus vishniacii

Cryptococcus_albidus[ana]
Filobasidium_floriforme[teleo])
Cryptococcus_gattii[ana]
(Filobasidiella_bacillispora[teleo])
Cryptococcus_uniguttulatus[ana]
(Filobasidium_uniguttulatum[teleo])
Cryptotrichosporon anacardii
Cupriavidus campinensis
Cupriavidus gilardii
Cupriavidus metallidurans
Cupriavidus necator
Cupriavidus oxalaticus
Cupriavidus pauculus
Cupriavidus respiraculi
Cupriavidus sp[1]
Curtobacterium albidum
Curtobacterium flaccumfaciens
Curtobacterium luteum
Curtobacterium sp
Cyberlindnera mississippiensis
Debaryomyces etchellsii
Debaryomyces hansenii
Deinococcus geothermalis
Delftia acidovorans
Dermabacter hominis
Dermacoccus nishinomiyaensis
Dermatophilus congolensis



Devosia riboflavina
Dialister micraerophilus
Dialister pneumosintes
Dichelobacter nodosus
Dickeya chrysanthemi
Dickeya dadantii
Dickeya dianthicola
Dickeya paradisiaca
Dickeya zeae
Dietzia cinnamea
Dietzia maris
Dietzia natronolimnaea
Dysgonomonas gadei
Edwardsiella hoshinae
Edwardsiella ictaluri
Edwardsiella tarda
Eggerthella lenta
Eggerthia catenaformis
Eikenella corrodens
Elizabethkingia meningoseptica
Elizabethkingia miricola
Emericella nidulans
Empedobacter brevis
Enterobacter aerogenes

Enterobacter amnigenus
Enterobacter asburiae
Enterobacter cancerogenus
Enterobacter cloacae
Enterobacter cowanii
Enterobacter gergoviae
Enterobacter hormaechei
Enterobacter kobei
Enterobacter ludwigii
Enterobacter pyrinus
Enterobacter radicincitans
Enterococcus aquimarinus
Enterococcus asini
Enterococcus avium
Enterococcus caccae
Enterococcus canintestini
Enterococcus canis
Enterococcus casseliflavus
Enterococcus cecorum
Enterococcus columbae
Enterococcus devriesei
Enterococcus dispar
Enterococcus durans
Enterococcus faecalis

Enterococcus faecium
Enterococcus gallinarum
Enterococcus gilvus
Enterococcus haemoperoxidus
Enterococcus hermanniensis
Enterococcus hirae
Enterococcus italicus
Enterococcus malodoratus
Enterococcus moraviensis
Enterococcus mundtii
Enterococcus pallens
Enterococcus phoeniculicola
Enterococcus pseudoavium
Enterococcus raffinosus
Enterococcus ratti
Enterococcus saccharolyticus
Enterococcus silesiacus
Enterococcus sulfureus
Enterococcus termitis
Enterococcus thailandicus
Enterococcus villorum
Epidermophyton floccosum
Erwinia amylovora
Erwinia billiniae



Erwinia mallotivora
Erwinia papayae
Erwinia persicina
Erwinia psidii
Erwinia pyrifoliae
Erwinia rhapontici
Erwinia sp
Erwinia tasmaniensis
Erwinia tracheiphila
Erysipelothrix inopinata
Erysipelothrix rhusiopathiae
Erysipelothrix tonsillarum
Escherichia albertii
Escherichia coli
Escherichia fergusonii
Escherichia hermannii
Escherichia vulneris
Eubacterium brachy
Eubacterium callanderi
Eubacterium limosum
Eubacterium sp[2]
Eubacterium yurii
Ewingella americana
Exiguobacterium aurantiacum

Exiguobacterium sp[4]
Exophiala dermatitidis
Facklamia hominis
Facklamia languida
Filifactor villosus
Filobasidium capsuligenum
Finegoldia magna
Flavobacterium flevense
Flavobacterium gelidilacus
Flavobacterium hibernum
Flavobacterium hydatis
Flavobacterium johnsoniae
Flavobacterium lindanitolerans
Flavobacterium pectinovorum
Flavobacterium saccharophilum
Flavonifractor plautii
Francisella philomiragia
Fusarium poae
Fusarium proliferatum
Fusobacterium canifelinum
Fusobacterium equinum
Fusobacterium gonidiaformans
Fusobacterium mortiferum
Fusobacterium naviforme

Fusobacterium necrophorum
Fusobacterium nucleatum
Fusobacterium periodonticum
Fusobacterium sp
Fusobacterium ulcerans
Fusobacterium varium
Gallibacterium anatis
Gardnerella sp
Gardnerella vaginalis
Gemella bergeri
Gemella haemolysans
Gemella morbillorum
Gemella sanguinis
Geobacillus kaustophilus
Geobacillus stearothermophilus
Geobacillus thermodenitrificans
Geobacillus thermoglucosidasius
Geotrichum silvicola
Geotrichum sp
Geotrichum_candidum[ana]
(Galactomyces_geotrichum[teleo])
Geotrichum_ingens[ana]
(Dipodascus_ingens[teleo])
Globicatella sulfidifaciens
Gluconacetobacter intermedius
Gluconacetobacter liquefaciens



Gluconobacter cerinus
Gluconobacter oxydans
Gordonia aichiensis
Gordonia alkanivorans
Gordonia australis
Gordonia bronchialis
Gordonia rubripertincta
Gordonia sputi
Gordonia terrae
Gracilibacillus dipsosauri
Granulicatella adiacens
Granulicatella balaenopterae
Granulicatella elegans
Grimontia hollisae
Guehomyces pullulans
Haemophilus ducreyi
Haemophilus haemoglobinophilus
Haemophilus haemolyticus
Haemophilus influenzae
Haemophilus parahaemolyticus
Haemophilus parainfluenzae
Haemophilus paraphrohaemolyticus
Haemophilus parasuis
Haemophilus pittmaniae

Hafnia alvei
Haloarcula vallismortis
Halobacterium salinarum
Halobacterium sp
Halococcus morrhuae
Halomonas aquamarina
Halomonas cupida
Halomonas elongata
Halomonas halmophila
Halomonas halodenitrificans
Halomonas pacifica
Halomonas venusta
Halotalea alkalilenta
Hannaella luteola
Hanseniасpora lachancei
Hanseniасpora opuntiae
Helcococcus kunzii
Helcococcus ovis
Helcococcus sueciensis
Helicobacter canadensis
Helicobacter canis
Helicobacter cholecystus
Helicobacter cinaedi
Helicobacter fennelliae

Helicobacter mustelae
Helicobacter pullorum
Helicobacter pylori
Herbaspirillum aquaticum
Herbaspirillum autotrophicum
Herbaspirillum chlorophenolicum
Herbaspirillum frisingense
Herbaspirillum hiltneri
Herbaspirillum huttiense
Herbaspirillum lusitanum
Herbaspirillum rhizosphaerae
Herbaspirillum rubrisubalbicans
Herbaspirillum seropedicae
Herbaspirillum sp
Herminiimonas arsenicoxydans
Herminiimonas fonticola
Histophilus somni
Hydrogenibacillus schlegelii
Hydrogenophaga flava
Hydrogenophaga pseudoflava
Hyphomicrobium sp
Ideonella dechloratans
Inquilinus limosus
Iodobacter fluviatilis



Issatchenkia terricola
Janthinobacterium lividum
Jeotgalicoccus halotolerans
Jonesia denitrificans
Kandleria vitulina
Kazachstania bovina
Kazachstania telluris
Kerstersia gyiorum
Kingella denitrificans
Kingella kingae
Kingella oralis
Kingella potus
Kitasatospora phosalacinea
Klebsiella oxytoca
Klebsiella pneumoniae
Klebsiella variicola
Kloeckera apiculata[ana]
(Hanseniaspora_uvarum[teleo])
Kluyvera ascorbata
Kluyvera cryocrescens
Kluyvera georgiana
Kluyvera intermedia
Kocuria aegyptia
Kocuria carniphila
Kocuria himachalensis

Kocuria kristinae
Kocuria marina
Kocuria palustris
Kocuria polaris
Kocuria rhizophila
Kocuria rosea
Kocuria sp
Kocuria varians
Kytococcus schroeteri
Kytococcus sedentarius
Lachancea fermentati
Lachnoanaerobaculum orale
Lachnoanaerobaculum saburreum
Lachnoanaerobaculum umeaense
Lactobacillus acidifarinae
Lactobacillus acidipiscis
Lactobacillus acidophilus
Lactobacillus agilis
Lactobacillus algidus
Lactobacillus alimentarius
Lactobacillus amylolyticus
Lactobacillus amylophilus
Lactobacillus amylotrophicus
Lactobacillus amylovorus

Lactobacillus antri
Lactobacillus apodemii
Lactobacillus aviarius
Lactobacillus bifermentans
Lactobacillus brevis
Lactobacillus buchneri
Lactobacillus casei
Lactobacillus coleohominis
Lactobacillus collinoides
Lactobacillus concavus
Lactobacillus coryniformis
Lactobacillus crispatus
Lactobacillus curvatus
Lactobacillus delbrueckii
Lactobacillus diolivorans
Lactobacillus equi
Lactobacillus farciminis
Lactobacillus fermentum
Lactobacillus fructivorans
Lactobacillus frumenti
Lactobacillus fuchuensis
Lactobacillus gallinarum
Lactobacillus gasseri
Lactobacillus gastricus



Lactobacillus graminis
Lactobacillus hammesii
Lactobacillus hamsteri
Lactobacillus harbinensis
Lactobacillus helveticus
Lactobacillus hilgardii
Lactobacillus homohiochii
Lactobacillus iners
Lactobacillus ingluviei
Lactobacillus intestinalis
Lactobacillus jensenii
Lactobacillus johnsonii
Lactobacillus kalixensis
Lactobacillus kefiranofaciens
Lactobacillus kefiri
Lactobacillus kimchii
Lactobacillus kitasatonis
Lactobacillus kunkeei
Lactobacillus lindneri
Lactobacillus malefermentans
Lactobacillus mali
Lactobacillus manihotivorans
Lactobacillus mindensis
Lactobacillus mucosae

Lactobacillus murinus
Lactobacillus nagelii
Lactobacillus nantensis
Lactobacillus oligofermentans
Lactobacillus oris
Lactobacillus panis
Lactobacillus pantheris
Lactobacillus parabuchneri
Lactobacillus paracasei
Lactobacillus paracollinoides
Lactobacillus parakefiri
Lactobacillus paralimentarius
Lactobacillus paraplantarum
Lactobacillus pentosus
Lactobacillus perolens
Lactobacillus plantarum
Lactobacillus pontis
Lactobacillus psittaci
Lactobacillus rennini
Lactobacillus reuteri
Lactobacillus rhamnosus
Lactobacillus rossiae
Lactobacillus ruminis
Lactobacillus saerimneri

Lactobacillus sakei
Lactobacillus salivarius
Lactobacillus sanfranciscensis
Lactobacillus satsumensis
Lactobacillus sharpeae
Lactobacillus sp
Lactobacillus spicheri
Lactobacillus suebicus
Lactobacillus ultunensis
Lactobacillus vaccinostercus
Lactobacillus vaginalis
Lactobacillus versmoldensis
Lactobacillus vini
Lactobacillus zeae
Lactobacillus zymae
Lactococcus garvieae
Lactococcus lactis
Lactococcus piscium
Lactococcus plantarum
Lactococcus raffinolactis
Laribacter hongkongensis
Lechevalieria flava
Leclercia adecarboxylata
Lecytophora hoffmannii



Legionella anisa
Legionella beliardensis
Legionella birminghamensis
Legionella bozemanii
Legionella brunensis
Legionella cherrii
Legionella cincinnatiensis
Legionella dresdenensis
Legionella dumoffii
Legionella erythra
Legionella feeleii
Legionella geestiana
Legionella gormanii
Legionella gratiana
Legionella hackeliae
Legionella impletisoli
Legionella israelensis
Legionella jamestowniensis
Legionella jordanis
Legionella lansingensis
Legionella longbeachae
Legionella maceachernii
Legionella micdadei
Legionella moravica

Legionella oakridgensis
Legionella parisiensis
Legionella pneumophila
Legionella rubrilucens
Legionella sainthelensi
Legionella santicrucis
Legionella sp
Legionella tucsonensis
Legionella wadsworthii
Legionella waltersii
Legionella worsleiensis
Legionella yabuuchiae
Leifsonia aquatica
Leminorella grimontii
Leminorella richardii
Leptothrix mobilis
Leptotrichia sp
Leptotrichia trevisanii
Leptotrichia wadei
Leucobacter chironomi
Leucobacter denitrificans
Leuconostoc carnosum
Leuconostoc citreum
Leuconostoc gelidum

Leuconostoc holzapfelii
Leuconostoc inhae
Leuconostoc lactis
Leuconostoc mesenteroides
Leuconostoc palmae
Leuconostoc pseudomesenteroides
Listeria grayi
Listeria innocua
Listeria ivanovii
Listeria monocytogenes
Listeria seeligeri
Listeria welshimeri
Lodderomyces elongisporus
Luteibacter rhizovicinus
Lysinibacillus boronitolerans
Lysinibacillus fusiformis
Lysinibacillus sphaericus
Macrococcus caseolyticus
Magnusiomyces capitatus
Malassezia furfur
Malassezia pachydermatis
Malikia spinosa
Mannheimia glucosida
Mannheimia granulomatis



Mannheimia haemolytica
Mannheimia varigena
Marinibacillus marinus
Marinilactibacillus psychrotolerans
Massilia sp
Massilia timonae
Megamonas sp[2]
Megasphaera micronuciformis
Mesorhizobium loti
Methanomonas methylovora
Methyloarcula marina
Methyloarcula terricola
Methylobacillus glycozenes
Methylobacillus sp
Methylobacterium extorquens
Methylobacterium fujisawaense
Methylobacterium mesophilicum
Methylobacterium organophilum
Methylobacterium radiotolerans
Methylobacterium rhodesianum
Methylobacterium rhodinum
Methylobacterium sp
Methylobacterium zatmanii
Microbacterium aerolatum

Microbacterium arborescens
Microbacterium aurum
Microbacterium barkeri
Microbacterium dextranolyticum
Microbacterium flavescens
Microbacterium flavum
Microbacterium foliorum
Microbacterium halotolerans
Microbacterium hominis
Microbacterium hydrocarbonoxydans
Microbacterium imperiale
Microbacterium keratanolyticum
Microbacterium ketosireducens
Microbacterium koreense
Microbacterium lacticum
Microbacterium laevaniformans
Microbacterium liquefaciens
Microbacterium luteolum
Microbacterium maritypicum
Microbacterium mitrae
Microbacterium natoriense
Microbacterium oleivorans
Microbacterium oxydans
Microbacterium paludicola

Microbacterium phyllosphaerae
Microbacterium resistens
Microbacterium saperdae
Microbacterium schleiferi
Microbacterium sp
Microbacterium terrae
Microbacterium terregens
Microbacterium testaceum
Microbacterium thalassium
Microbacterium trichothecenolyticum
Microbacterium ulmi
Micrococcus flavus
Micrococcus luteus
Micrococcus lylae
Micrococcus terreus
Micromonospora aurantiaca
Micromonospora carbonacea
Micromonospora chalcea
Micromonospora chersina
Micromonospora citrea
Micromonospora coerulea
Micromonospora echinaurantiaca
Micromonospora echinofusca
Micromonospora echinospora



Micromonospora inyonensis	Moraxella_sg_Moraxella lincolnii	Mycobacterium kumamotonense
Micromonospora peucetica	Moraxella_sg_Moraxella nonliquefaciens	Mycobacterium lacus
Micromonospora purpureochromogenes	Moraxella_sg_Moraxella oblonga	Mycobacterium mageritense
Micromonospora sagamiensis	Moraxella_sg_Moraxella osloensis	Mycobacterium malmoense
Micromonospora sp	Moraxella_sg_Moraxella pluranimalium	Mycobacterium marinum
Micromonospora viridifaciens	Morganella morganii	Mycobacterium montefiorensis
Microsporium canis	Mucor circinelloides	Mycobacterium mucogenicum
Microsporium_gypseum[ana] (Arthroderma_gypseum[teleo])	Mycobacterium abscessus	Mycobacterium neoaurum
Mobiluncus curtisii	Mycobacterium agri	Mycobacterium palustre
Mobiluncus sp	Mycobacterium asiaticum	Mycobacterium peregrinum
Moellerella wisconsensis	Mycobacterium avium	Mycobacterium phlei
Moesziomyces bullatus	Mycobacterium boenickei	Mycobacterium pseudoshottsii
Moorella thermoacetica	Mycobacterium bovis	Mycobacterium pulveris
Moraxella sp	Mycobacterium celatum	Mycobacterium rhodesiae
Moraxella_sg_Branhamella catarrhalis	Mycobacterium chelonae	Mycobacterium seoulense
Moraxella_sg_Branhamella ovis	Mycobacterium chlorophenicum	Mycobacterium shottsii
Moraxella_sg_Moraxella atlantae	Mycobacterium conceptionense	Mycobacterium simiae
Moraxella_sg_Moraxella boevrei	Mycobacterium farcinogenes	Mycobacterium smegmatis
Moraxella_sg_Moraxella bovis	Mycobacterium fortuitum	Mycobacterium szulgai
Moraxella_sg_Moraxella bovoculi	Mycobacterium gordonae	Mycobacterium thermoresistibile
Moraxella_sg_Moraxella canis	Mycobacterium heckeshornense	Mycobacterium tokaiense
Moraxella_sg_Moraxella caprae	Mycobacterium hiberniae	Mycobacterium tuberculosis
Moraxella_sg_Moraxella equi	Mycobacterium hodleri	Mycobacterium ulcerans
Moraxella_sg_Moraxella lacunata	Mycobacterium kansasii	Mycobacterium xenopi



Mycoplasma alkalescens
Mycoplasma arginini
Mycoplasma bovirhinis
Mycoplasma bovis
Mycoplasma canis
Mycoplasma gallinaceum
Mycoplasma gallisepticum
Mycoplasma hyorhinis
Mycoplasma ovipneumoniae
Mycoplasma pullorum
Myroides odoratimimus
Myroides odoratus
Neisseria bacilliformis
Neisseria canis
Neisseria cinerea
Neisseria elongata
Neisseria flavescens
Neisseria gonorrhoeae
Neisseria lactamica
Neisseria macacae
Neisseria meningitidis
Neisseria mucosa
Neisseria perflava
Neisseria polysaccharea

Neisseria sicca
Neisseria sp[2]
Neisseria subflava
Neisseria weaveri
Neisseria zoodegmatidis
Nesterenkonia lacusekhoensis
Nocardia abscessus
Nocardia africana
Nocardia anaemiae
Nocardia aobensis
Nocardia araoensis
Nocardia arthritis
Nocardia asiatica
Nocardia asteroides
Nocardia brasiliensis
Nocardia carnea
Nocardia concava
Nocardia cyriacigeorgica
Nocardia elegans
Nocardia exalbida
Nocardia farcinica
Nocardia higoensis
Nocardia ignorata
Nocardia kruzakiae

Nocardia niigatensis
Nocardia nova
Nocardia otitidiscaviarum
Nocardia paucivorans
Nocardia pneumoniae
Nocardia salmonicida
Nocardia seriolae
Nocardia sienata
Nocardia sp
Nocardia testacea
Nocardia thailandica
Nocardia transvalensis
Nocardia vermiculata
Nocardia veterana
Nocardia yamanashiensis
Nocardioides jensenii
Nocardioides simplex
Nocardiopsis alba
Novosphingobium aromaticivorum
Novosphingobium hassiacum
Novosphingobium lentum
Novosphingobium naphthalenivorans
Novosphingobium nitrogenifigens
Novosphingobium pentaromativorans



Novosphingobium resinovorum	Paenibacillus alginolyticus	Paenibacillus gansuensis
Novosphingobium rosa	Paenibacillus alvei	Paenibacillus glucanolyticus
Novosphingobium stygium	Paenibacillus amylolyticus	Paenibacillus glycanilyticus
Novosphingobium subterraneum	Paenibacillus anaericanus	Paenibacillus graminis
Novosphingobium taihuense	Paenibacillus apiarius	Paenibacillus illinoisensis
Novosphingobium tardaugens	Paenibacillus assamensis	Paenibacillus jamilae
Ochrobactrum anthropi	Paenibacillus azoreducens	Paenibacillus kobensis
Ochrobactrum gallinifaecis	Paenibacillus barcinonensis	Paenibacillus lactis
Ochrobactrum grignonense	Paenibacillus barengoltzii	Paenibacillus larvae
Ochrobactrum intermedium	Paenibacillus borealis	Paenibacillus lautus
Ochrobactrum sp[3]	Paenibacillus brasiliensis	Paenibacillus macerans
Ochrobactrum tritici	Paenibacillus chibensis	Paenibacillus macquariensis
Odoribacter splanchnicus	Paenibacillus chinjuensis	Paenibacillus massiliensis
Oerskovia turbata	Paenibacillus chitinolyticus	Paenibacillus mendelii
Ogataea polymorpha	Paenibacillus chondroitinus	Paenibacillus motobuensis
Oligella ureolytica	Paenibacillus cineris	Paenibacillus naphthalenovorans
Oligella urethralis	Paenibacillus cookii	Paenibacillus nematophilus
Olsenella profusa	Paenibacillus curdlanolyticus	Paenibacillus odorifer
Olsenella uli	Paenibacillus daejeonensis	Paenibacillus pabuli
Ornithobacterium rhinotracheale	Paenibacillus dendritiformis	Paenibacillus pasadenensis
Paecilomyces lilanicus	Paenibacillus durus	Paenibacillus peoriae
Paecilomyces variotii	Paenibacillus edaphicus	Paenibacillus phyllosphaerae
Paenibacillus agarexedens	Paenibacillus ehimensis	Paenibacillus polymyxa
Paenibacillus agaridevorans	Paenibacillus favisporus	Paenibacillus rhizosphaerae



Paenibacillus sabiniae	Pantoea dispersa	Pectobacterium carotovorum
Paenibacillus sp	Pantoea gaviniae	Pectobacterium cyripedii
Paenibacillus stellifer	Pantoea septica	Pectobacterium wasabiae
Paenibacillus taiwanensis	Pantoea sp	Pediococcus acidilactici
Paenibacillus terrae	Pantoea stewartii	Pediococcus pentosaceus
Paenibacillus thiaminolyticus	Parabacteroides distasonis	Pelomonas saccharophila
Paenibacillus timonensis	Parabacteroides goldsteinii	Pelomonas sp[3]
Paenibacillus turicensis	Parabacteroides johnsonii	Penicillium camemberti
Paenibacillus urinalis	Paracoccus denitrificans	Penicillium chrysogenum
Paenibacillus validus	Paracoccus versutus	Penicillium sp
Paenibacillus wynnii	Paracoccus yeei	Peptococcus niger
Paenibacillus xinjiangensis	Parascardovia denticolens	Peptoniphilus asaccharolyticus
Paenibacillus xylanilyticus	Parvimonas micra	Peptoniphilus gorbachii
Paenibacillus zanthoxyli	Pasteurella aerogenes	Peptoniphilus harei
Pandoraea apista	Pasteurella bettyae	Peptoniphilus indolicus
Pandoraea norimbergensis	Pasteurella canis	Peptoniphilus ivorii
Pandoraea pnomenusa	Pasteurella dagmatis	Peptoniphilus sp
Pandoraea pulmonicola	Pasteurella mairii	Peptostreptococcus anaerobius
Pandoraea sp[2]	Pasteurella multocida	Peptostreptococcus sp
Pandoraea sputorum	Pasteurella pneumotropica	Phenylobacterium koreense
Pannonibacter phragmitetus	Pasteurella stomatis	Phoma exigua
Pantoea agglomerans	Pectobacterium atrosepticum	Phoma glomerata
Pantoea ananatis	Pectobacterium betavasculorum	Phoma herbarum
Pantoea calida	Pectobacterium cacticida	Photobacterium damsela



Photobacterium iliopiscarium

Photorhabdus asymbiotica

Photorhabdus luminescens

Photorhabdus temperata

Pichia cactophila

Pichia farinosa

Pichia holstii

Pichia kluyveri

Pichia manshurica

Pichia methylivora

Pichia pseudocactophila

Pigmentiphaga daeguensis

Plesiomonas shigelloides

Porphyromonas asaccharolytica

Porphyromonas gingivalis

Porphyromonas gulae

Porphyromonas levii

Porphyromonas macacae

Porphyromonas sp

Pragia fontium

Prevotella amnii

Prevotella baroniae

Prevotella bergensis

Prevotella bivia

Prevotella buccae

Prevotella buccalis

Prevotella corporis

Prevotella denticola

Prevotella disiens

Prevotella histicola

Prevotella intermedia

Prevotella maculosa

Prevotella melaninogenica

Prevotella multisaccharivorax

Prevotella nanceiensis

Prevotella nigrescens

Prevotella oralis

Prevotella oris

Prevotella oulorum

Prevotella pallens

Prevotella salivae

Prevotella shahii

Prevotella sp

Prevotella stercorea

Propionibacterium acidifaciens

Propionibacterium acidipropionici

Propionibacterium acnes

Propionibacterium australiense

Propionibacterium avidum

Propionibacterium cyclohexanicum

Propionibacterium freudenreichii

Propionibacterium granulosum

Propionibacterium jensenii

Propionibacterium microaerophilum

Propionibacterium propionicum

Propionibacterium sp

Propionibacterium thoenii

Propioniferax innocua

Propionimicrobium lymphophilum

Proteus hauseri

Proteus mirabilis

Proteus myxofaciens

Proteus penneri

Proteus vulgaris

Prototheca wickerhamii

Providencia alcalifaciens

Providencia heimbachae

Providencia rettgeri

Providencia rustigianii

Providencia stuartii

Providencia vermicola

Pseudacidovorax intermedius



Pseudochrobactrum asaccharolyticum

Pseudoclavibacter helvolus

Pseudoclavibacter sp

Pseudomonas abietaniphila

Pseudomonas aeruginosa

Pseudomonas agarici

Pseudomonas alcaligenes

Pseudomonas alcaliphila

Pseudomonas anguilliseptica

Pseudomonas antarctica

Pseudomonas asplenii

Pseudomonas avellanae

Pseudomonas azotifigens

Pseudomonas azotoformans

Pseudomonas balearica

Pseudomonas boreopolis

Pseudomonas brassicacearum

Pseudomonas brenneri

Pseudomonas caricapapayae

Pseudomonas cedrina

Pseudomonas chlororaphis

Pseudomonas cichorii

Pseudomonas citronellolis

Pseudomonas congelans

Pseudomonas corrugata

Pseudomonas extremorientalis

Pseudomonas flavescens

Pseudomonas fluorescens

Pseudomonas fragi

Pseudomonas frederiksbergensis

Pseudomonas fulva

Pseudomonas fuscovaginae

Pseudomonas gessardii

Pseudomonas graminis

Pseudomonas grimontii

Pseudomonas indica

Pseudomonas jessenii

Pseudomonas jinjuensis

Pseudomonas kilonensis

Pseudomonas koreensis

Pseudomonas libanensis

Pseudomonas lundensis

Pseudomonas lutea

Pseudomonas luteola

Pseudomonas mandelii

Pseudomonas marginalis

Pseudomonas mendocina

Pseudomonas migulae

Pseudomonas monteilii

Pseudomonas mosselii

Pseudomonas mucidolens

Pseudomonas nitroreducens

Pseudomonas oleovorans

Pseudomonas orientalis

Pseudomonas oryzihabitans

Pseudomonas otitidis

Pseudomonas panipatensis

Pseudomonas pertucinogena

Pseudomonas pictorum

Pseudomonas plecoglossicida

Pseudomonas poae

Pseudomonas pohangensis

Pseudomonas proteolytica

Pseudomonas pseudoalcaligenes

Pseudomonas putida

Pseudomonas resinovorans

Pseudomonas rhizosphaerae

Pseudomonas rhodesiae

Pseudomonas savastanoi

Pseudomonas segetis

Pseudomonas sp

Pseudomonas straminea



Pseudomonas stutzeri
Pseudomonas synxantha
Pseudomonas syringae
Pseudomonas taetrolens
Pseudomonas thermotolerans
Pseudomonas thivervalensis
Pseudomonas tolaasii
Pseudomonas trivialis
Pseudomonas umsongensis
Pseudomonas vancouverensis
Pseudomonas veronii
Pseudomonas viridiflava
Pseudomonas xanthomarina
Pseudonocardia hydrocarbonoxydans
Pseudoxanthomonas kaohsiungensis
Pseudoxanthomonas spadix
Pseudozyma aphidis
Psychrobacillus insolitus
Psychrobacillus psychrodurans
Psychrobacillus psychrotolerans
Psychrobacter lutiphocae
Psychrobacter sp
Pullulanibacillus naganoensis
Rahnella aquatilis

Ralstonia insidiosa
Ralstonia mannitolilytica
Ralstonia picketii
Ralstonia sp
Ralstonia syzygii
Raoultella ornithinolytica
Raoultella planticola
Raoultella terrigena
Rathayibacter rathayi
Rheinheimera soli
Rhizobium radiobacter
Rhizobium rubi
Rhizobium tropici
Rhizopus microsporus
Rhodobacter aestuarii
Rhodobacter veldkampii
Rhodococcus aetherivorans
Rhodococcus baikonurensis
Rhodococcus coprophilus
Rhodococcus corynebacterioides
Rhodococcus equi
Rhodococcus erythropolis
Rhodococcus fascians
Rhodococcus globerulus

Rhodococcus gordoniae
Rhodococcus imtechensis
Rhodococcus jostii
Rhodococcus koreensis
Rhodococcus kroppenstedtii
Rhodococcus kunmingensis
Rhodococcus maanshanensis
Rhodococcus marinonascens
Rhodococcus opacus
Rhodococcus percolatus
Rhodococcus phenolicus
Rhodococcus pyridinivorans
Rhodococcus rhodnii
Rhodococcus rhodochrous
Rhodococcus ruber
Rhodococcus triatomae
Rhodococcus wratislaviensis
Rhodococcus yunnanensis
Rhodococcus zopfii
Rhodospiridium sp
Rhodotorula acheniorum
Rhodotorula bacarum
Rhodotorula bogoriensis
Rhodotorula glutinis



Rhodotorula minuta
Rhodotorula mucilaginosa
Rhodotorula pustula
Riemerella anatipestifer
Riemerella columbina
Roseomonas mucosa
Rothia aeria
Rothia amarae
Rothia dentocariosa
Rothia mucilaginosa
Rothia nasimurium
Rubrivivax gelatinosus
Ruminococcus gnavus
Rummeliibacillus pycnus
Saccharomyces cerevisiae
Saccharopolyspora erythraea
Saccharopolyspora hirsuta
Saccharothrix mutabilis
Salimicrobium halophilum
Salinivibrio costicola
Salmonella sp (bongori)
Salmonella sp (choleraesuis)
Salmonella sp (enterica st Enterica)
Salmonella sp (enteritidis)

Salmonella sp (typhimurium)
Samsonia erythrinae
Saprochaete clavata
Saprochaete suaveolens
Sarocladium strictum
Scedosporium aurantiacum
Schizophyllum commune
Schizosaccharomyces pombe
Scopulariopsis_brevicaulis[ana]
(Microascus_brevicaulis[teleo])
Selenomonas infelix
Selenomonas sputigena
Serratia entomophila
Serratia ficaria
Serratia fonticola
Serratia grimesii
Serratia liquefaciens
Serratia marcescens
Serratia odorifera
Serratia plymuthica
Serratia proteamaculans
Serratia quinivorans
Serratia rubidaea
Serratia ureilytica
Shewanella algae

Shewanella baltica
Shewanella fidelis
Shewanella frigidimarina
Shewanella profunda
Shewanella putrefaciens
Shimwellia blattae
Sinomonas atrocyanea
Slackia exigua
Slackia heliotrinireducens
Sodalis glossinidius
Solibacillus silvestris
Solobacterium moorei
Sphingobacterium faecium
Sphingobacterium mizutaii
Sphingobacterium multivorum
Sphingobacterium spiritivorum
Sphingobacterium thalpophilum
Sphingobium amiense
Sphingobium aromaticiconvertens
Sphingobium chlorophenicum
Sphingobium cloacae
Sphingobium fuliginis
Sphingobium herbicidovorans
Sphingobium indicum



Sphingobium japonicum
Sphingobium olei
Sphingobium xenophagum
Sphingomonas abaci
Sphingomonas adhaesiva
Sphingomonas aerolata
Sphingomonas aquatilis
Sphingomonas asaccharolytica
Sphingomonas aurantiaca
Sphingomonas azotifigens
Sphingomonas desiccabilis
Sphingomonas faeni
Sphingomonas haloaromaticamans
Sphingomonas koreensis
Sphingomonas mali
Sphingomonas melonis
Sphingomonas molluscorum
Sphingomonas mucosissima
Sphingomonas pannii
Sphingomonas parapaucimobilis
Sphingomonas paucimobilis
Sphingomonas phyllosphaerae
Sphingomonas pituitosa
Sphingomonas pruni

Sphingomonas pseudosanguinis
Sphingomonas sanguinis
Sphingomonas soli
Sphingomonas sp
Sphingomonas trueperi
Sphingomonas wittichii
Sphingomonas yabuuchiae
Sphingomonas yunnanensis
Sphingopyxis baekryungensis
Sphingopyxis chilensis
Sphingopyxis macrogoltabida
Sphingopyxis terrae
Sphingopyxis witflariensis
Sporobolomyces roseus
Sporobolomyces_salmonicolor[ana]
(Sporidiobolus_salmonicolor[teleo])
Sporolactobacillus laevolacticus
Sporopachydermia cereana
Sporosarcina globispora
Sporosarcina psychrophila
Sporothrix schenckii
Staphylococcus arlettae
Staphylococcus aureus
Staphylococcus auricularis
Staphylococcus capitis

Staphylococcus caprae
Staphylococcus carnosus
Staphylococcus chromogenes
Staphylococcus cohnii
Staphylococcus condimenti
Staphylococcus delphini
Staphylococcus epidermidis
Staphylococcus equorum
Staphylococcus felis
Staphylococcus fleurettii
Staphylococcus haemolyticus
Staphylococcus hominis
Staphylococcus hyicus
Staphylococcus intermedius
Staphylococcus kloosii
Staphylococcus lentus
Staphylococcus lugdunensis
Staphylococcus lutrae
Staphylococcus muscae
Staphylococcus nepalensis
Staphylococcus pasteurii
Staphylococcus pettenkoferi
Staphylococcus piscifermentans
Staphylococcus pseudintermedius



Staphylococcus saccharolyticus
Staphylococcus saprophyticus
Staphylococcus schleiferi
Staphylococcus sciuri ssp sciuri
Staphylococcus simiae
Staphylococcus simulans
Staphylococcus sp[1]
Staphylococcus succinus
Staphylococcus vitulinus
Staphylococcus warneri
Staphylococcus xylosus

Starkaya novella

Stenotrophomonas acidaminiphila

Stenotrophomonas maltophilia

Stenotrophomonas nitritireducens

Stenotrophomonas rhizophila

Stenotrophomonas sp

Stenotrophomonas_maltophilia
(Pseudomonas_beteli_#)

Stenotrophomonas_maltophilia
(Pseudomonas_geniculata_#)

Stenotrophomonas_maltophilia
(Pseudomonas_hibiscicola_#)

Streptobacillus moniliformis

Streptococcus acidominimus

Streptococcus agalactiae

Streptococcus alactolyticus

Streptococcus anginosus

Streptococcus australis

Streptococcus caballi

Streptococcus canis

Streptococcus castoreus

Streptococcus constellatus

Streptococcus criceti

Streptococcus cristatus

Streptococcus dentirousetti

Streptococcus devriesei

Streptococcus didelphis

Streptococcus downei

Streptococcus dysgalactiae

Streptococcus entericus

Streptococcus equi

Streptococcus equinus

Streptococcus ferus

Streptococcus gallinaceus

Streptococcus gallolyticus

Streptococcus gordonii

Streptococcus halichoeri

Streptococcus henryi

Streptococcus hyointestinalis

Streptococcus hyovaginalis

Streptococcus infantarius

Streptococcus infantis

Streptococcus intermedius

Streptococcus lutetiensis

Streptococcus macacae

Streptococcus marimammalium

Streptococcus massiliensis

Streptococcus merionis

Streptococcus minor

Streptococcus mitis

Streptococcus mutans

Streptococcus oralis

Streptococcus orisratti

Streptococcus orisuis

Streptococcus ovis

Streptococcus parasanguinis

Streptococcus parauberis

Streptococcus peroris

Streptococcus phocae

Streptococcus pleomorphus

Streptococcus pluranimalium

Streptococcus pneumoniae

Streptococcus porcinus

Streptococcus pseudopneumoniae



Streptococcus pseudoporcinus
Streptococcus pyogenes
Streptococcus ratti
Streptococcus salivarius
Streptococcus sanguinis
Streptococcus sinensis
Streptococcus sobrinus
Streptococcus sp
Streptococcus suis
Streptococcus thoralensis
Streptococcus uberis
Streptococcus urinalis
Streptococcus vestibularis
Streptomyces albus
Streptomyces aureofaciens
Streptomyces avidinii
Streptomyces badius
Streptomyces chartreusis
Streptomyces galilaeus
Streptomyces griseus
Streptomyces hirsutus
Streptomyces hygroscopicus
Streptomyces lavendulae
Streptomyces phaeochromogenes

Streptomyces scabiei
Streptomyces sp
Streptomyces violaceoruber
Streptosporangium sibiricum
Suttonella indologenes
Tatumella citrea
Tatumella ptyseos
Tatumella punctata
Tatumella terrea
Taylorella asinigenitalis
Taylorella equigenitalis
Tenacibaculum discolor
Tenacibaculum ovolyticum
Terrabacter tumescens
Terrimonas ferruginea
Tessaracoccus flavescens
Tetragenococcus solitarius
Thauera aminoaromatica
Thauera aromatica
Thauera chlorobenzoica
Thauera linaloolentis
Thauera mechernichensis
Thauera phenylacetica
Thauera terpenica

Thermoactinomyces sp[2]
Thermoanaerobacter thermohydrosulfuricus
Thermoanaerobacterium
thermosaccharolyticum
Thermoanaerobacterium thermosulfurigenes
Tissierella praeacuta
Trabulsiella guamensis
Trichophyton interdigitale
Trichophyton mentagrophytes
Trichophyton rubrum
Trichophyton tonsurans
Trichosporon asahii
Trichosporon coremiiforme
Trichosporon cutaneum
Trichosporon debeurmannianum
Trichosporon dohaense
Trichosporon dulciturum
Trichosporon faecale
Trichosporon gracile
Trichosporon inkin
Trichosporon japonicum
Trichosporon jirovecii
Trichosporon laibachii
Trichosporon loubieri
Trichosporon moniliiforme



Trichosporon montevidense	Veillonella dispar	Vibrio gazogenes
Trichosporon mucoides	Veillonella magna	Vibrio gigantis
Trichosporon mycotoxinivorans	Veillonella montpellierensis	Vibrio harveyi
Trichosporon ovoides	Veillonella parvula	Vibrio hispanicus
Trichosporon sp	Veillonella ratti	Vibrio ichthyenteri
Trichosporon terricola	Veillonella rogosae	Vibrio kanaloae
Trueperella abortisuis	Veillonella sp[3]	Vibrio lentus
Trueperella bernardiae	Vibrio aerogenes	Vibrio mediterranei
Trueperella bialowiezense	Vibrio aestuarianus	Vibrio metschnikovii
Trueperella bonasi	Vibrio agarivorans	Vibrio mimicus
Trueperella pyogenes	Vibrio albensis	Vibrio mytili
Tsukamurella inchonensis	Vibrio alginolyticus	Vibrio natriegens
Tsukamurella paurometabola	Vibrio anguillarum	Vibrio navarrensis
Tsukamurella sp	Vibrio brasiliensis	Vibrio neptunius
Turicella otitidis	Vibrio campbellii	Vibrio nereis
Udeniomyces puniceus	Vibrio chagasii	Vibrio nigripulchritudo
Vagococcus fluvialis	Vibrio cincinnatiensis	Vibrio ordalii
Vagococcus lutrae	Vibrio coralliilyticus	Vibrio orientalis
Varibaculum cambriense	Vibrio cyclitrophicus	Vibrio ostreicida
Variovorax paradoxus	Vibrio diazotrophicus	Vibrio pacinii
Veillonella atypica	Vibrio ezuriae	Vibrio parahaemolyticus
Veillonella caviae	Vibrio fluvialis	Vibrio pectenida
Veillonella criceti	Vibrio fortis	Vibrio pelagius
Veillonella denticariosi	Vibrio furnissii	Vibrio penaeicida



<i>Vibrio pomeroyi</i>	<i>Weissella viridescens</i>	<i>Xenorhabdus ehlersii</i>
<i>Vibrio ponticus</i>	Wohlfahrtiimonas chitiniclastica	<i>Xenorhabdus innexi</i>
<i>Vibrio proteolyticus</i>	Wolinella succinogenes	<i>Xenorhabdus japonica</i>
<i>Vibrio rotiferianus</i>	Xanthobacter autotrophicus	<i>Xenorhabdus nematophila</i>
<i>Vibrio ruber</i>	Xanthomonas arboricola	<i>Xenorhabdus poinarii</i>
<i>Vibrio rumoiensis</i>	<i>Xanthomonas axonopodis</i>	<i>Xenorhabdus szentirmaii</i>
<i>Vibrio scophthalmi</i>	<i>Xanthomonas bromi</i>	Yersinia aldovae
<i>Vibrio shilonii</i>	<i>Xanthomonas campestris</i>	<i>Yersinia aleksiciae</i>
<i>Vibrio splendidus</i>	<i>Xanthomonas cassavae</i>	<i>Yersinia bercovieri</i>
<i>Vibrio superstes</i>	<i>Xanthomonas citri</i>	<i>Yersinia enterocolitica</i>
<i>Vibrio tasmaniensis</i>	<i>Xanthomonas codiae</i>	<i>Yersinia frederiksenii</i>
<i>Vibrio vulnificus</i>	<i>Xanthomonas cucurbitae</i>	<i>Yersinia intermedia</i>
<i>Vibrio xuii</i>	<i>Xanthomonas cynarae</i>	<i>Yersinia kristensenii</i>
Virgibacillus halodenitrificans	<i>Xanthomonas hortorum</i>	<i>Yersinia mollaretii</i>
<i>Virgibacillus pantothenicus</i>	<i>Xanthomonas hyacinthi</i>	<i>Yersinia pseudotuberculosis</i>
<i>Virgibacillus proomii</i>	<i>Xanthomonas melonis</i>	<i>Yersinia rohdei</i>
Viridibacillus arenosi	<i>Xanthomonas perforans</i>	<i>Yersinia ruckeri</i>
<i>Viridibacillus arvi</i>	<i>Xanthomonas pisi</i>	Yokenella regensburgae
<i>Viridibacillus neidei</i>	<i>Xanthomonas theicola</i>	Zygosaccharomyces bailii
Wautersiella falsenii	<i>Xanthomonas translucens</i>	<i>Zygosaccharomyces bisporus</i>
Weeksella virosa	<i>Xanthomonas vasicola</i>	<i>Zygosaccharomyces florentinus</i>
Weissella confusa	Xenorhabdus beddingii	<i>Zygosaccharomyces microellipsoides</i>
<i>Weissella halotolerans</i>	<i>Xenorhabdus bovienii</i>	<i>Zygosaccharomyces rouxii</i>
<i>Weissella minor</i>	<i>Xenorhabdus budapestensis</i>	



~ 2290 species



Release Notes for MBT-BDAL-5627 MSP library

What is new in MBT DB 5627?

The MBT database DB 4613 will be extended by **1020 newly generated reference MSP**, which cover **348 species** overall.

21 new genera and **113 new species** will be implemented and which will be covered by **220 new reference entries**.

Additional 800 MSP will improve the coverage of species already contained in the database.

Main improvements:

- Implementation of *Haemophilus haemolyticus* strains
- Additional strains for improvement of *Streptococcus mitis/oralis/pneumonia* identification
- Many new yeast species and strains
- More *Nocardia* strains to cover the MALDI diversity

The **21 new genera** and **113 new species** cover following groups:

	New genera	New species	Aerobe species	Microaerophile species	Anaerobe species
Gram -	6	23	18	3	2
Gram +	7	27	22	1	4
Yeast	8	63	63	-	-
Fil. Fungi	0	0	-	-	-

Overall improvements:

	New MSP (from 1020)	Species covered	Aerobe strains	Microaerophile strains	Anaerobe strains
Gram -	328	113	265	54	9
Gram +	374	128	185	179	10
Yeast	316	106	106	-	-
Fil. Fungi	2	1	2	-	-

Supported MALDI Biotyper versions (installed on respective operating systems):

- MBT-3.0
- MBT-3.1
- MBT-4.0



Implementation of MSP entries for following 113 new species:

	New genus/species			
1	Actinomyces ruminicola	new species	gram +	anaerobic
2	Aeromicrobium flavum	new genus/species	gram +	aerobic
3	Alicyclobacillus cycloheptanicus	new species	gram +	aerobic
4	Alloiococcus otitis	new genus/species	gram +	aerobic
5	Arthroascus schoenii	new genus/species	yeast	aerobic
6	Arthrobacter flavus	new species	gram +	aerobic
7	Arthrobacter nitroguajacolicus	new species	gram +	aerobic
8	Bacillus altitudinis	new species	gram +	aerobic
9	Bacillus horneckiae	new species	gram +	aerobic
10	Bartonella japonica	new genus/species	gram -	aerobic
11	Brevibacillus invocatus	new species	gram +	aerobic
12	Candida alai	new species	yeast	aerobic
13	Candida ambrosiae	new species	yeast	aerobic
14	Candida auris	new species	yeast	aerobic
15	Candida blattae	new species	yeast	aerobic
16	Candida bohiensis	new species	yeast	aerobic
17	Candida bracarenensis	new species	yeast	aerobic
18	Candida buenavistaensis	new species	yeast	aerobic
19	Candida carpophila	new species	yeast	aerobic
20	Candida castellii	new species	yeast	aerobic
21	Candida duobushaemulonii	new species	yeast	aerobic
22	Candida frijolesensis	new species	yeast	aerobic



	New genus/species			
23	<i>Candida infanticola</i>	new species	yeast	aerobic
24	<i>Candida labiduridarum</i>	new species	yeast	aerobic
25	<i>Candida pseudohaemulonii</i>	new species	yeast	aerobic
26	<i>Candida saitoana</i>	new species	yeast	aerobic
27	<i>Candida shehatae</i> var <i>insectosa</i>	new species	yeast	aerobic
28	<i>Candida sojae</i>	new species	yeast	aerobic
29	<i>Candida viswanathii</i>	new species	yeast	aerobic
30	<i>Candida_citrea</i> [ana] (<i>Pichia_nakasei</i> [teleo]#)	new species	yeast	aerobic
31	<i>Candida_fabianii</i> [ana] (<i>Pichia_fabianii</i> [teleo])	new species	yeast	aerobic
32	<i>Candida_holmii</i> [ana] (<i>Kazachstania_exigua</i> [teleo])	new species	yeast	aerobic
33	<i>Candida_pintolopesii</i> [ana] (<i>Kazachstania_pintolopesii</i> [teleo])	new species	yeast	aerobic
34	<i>Chryseobacterium ginsenosidimutans</i>	new species	gram -	aerobic
35	<i>Chryseobacterium hagamense</i>	new species	gram -	aerobic
36	<i>Cryptococcus albidosimilis</i>	new species	yeast	aerobic
37	<i>Cryptococcus curvatus</i>	new species	yeast	aerobic
38	<i>Cryptococcus diffluens</i>	new species	yeast	aerobic
39	<i>Cryptococcus flavescens</i>	new species	yeast	aerobic
40	<i>Cryptococcus gastricus</i>	new species	yeast	aerobic
41	<i>Cryptococcus humicola</i>	new species	yeast	aerobic
42	<i>Cryptococcus liquefaciens</i>	new species	yeast	aerobic
43	<i>Cryptococcus magnus</i>	new species	yeast	aerobic
44	<i>Cryptococcus saitoi</i>	new species	yeast	aerobic
45	<i>Cryptococcus terreus</i>	new species	yeast	aerobic
46	<i>Cryptococcus uzbekistanensis</i>	new species	yeast	aerobic



	New genus/species			
47	<i>Cryptococcus vishniacii</i>	new species	yeast	aerobic
48	<i>Cryptotrichosporon anacardii</i>	new genus/species	yeast	aerobic
49	<i>Cupriavidus campinensis</i>	new species	gram -	aerobic
50	<i>Cyberlindnera mississippiensis</i>	new genus/species	yeast	aerobic
51	<i>Deinococcus geothermalis</i>	new genus/species	gram +	aerobic
52	<i>Geobacillus thermodenitrificans</i> ssp <i>calidus</i>	new species	gram +	aerobic
53	<i>Geobacillus thermoglucoasidarius</i>	new species	gram +	aerobic
54	<i>Gluconacetobacter intermedius</i>	new genus/species	gram -	aerobic
55	<i>Gluconacetobacter liquefaciens</i>	new species	gram -	aerobic
56	<i>Gluconobacter cerinus</i>	new species	gram -	aerobic
57	<i>Guehomyces pullulans</i>	new genus/species	yeast	aerobic
58	<i>Haemophilus haemoglobinophilus</i>	new species	gram -	microaerophilic
59	<i>Haemophilus haemolyticus</i>	new species	gram -	microaerophilic
60	<i>Haemophilus paraphrohaemolyticus</i>	new species	gram -	microaerophilic
61	<i>Halotalea alkalilenta</i>	new genus/species	gram -	aerobic
62	<i>Hannaella luteola</i>	new genus/species	yeast	aerobic
63	<i>Hanseniaspora lachancei</i>	new species	yeast	aerobic
64	<i>Jeotgalicoccus halotolerans</i>	new genus/species	gram +	aerobic
65	<i>Kazachstania bovina</i>	new species	yeast	aerobic
66	<i>Kytococcus schroeteri</i>	new species	gram +	aerobic
67	<i>Lachancea fermentati</i>	new species	yeast	aerobic
68	<i>Lachnoanaerobaculum orale</i>	new genus/species	gram +	anaerobic
69	<i>Lachnoanaerobaculum saburreum</i>	new species	gram +	anaerobic
70	<i>Lachnoanaerobaculum umeaense</i>	new species	gram +	anaerobic



	New genus/species			
71	Laribacter hongkongensis	new genus/species	gram -	aerobic
72	Legionella dresdenensis	new species	gram -	aerobic
73	Legionella geestiana	new species	gram -	aerobic
74	Legionella gratiana	new species	gram -	aerobic
75	Legionella waltersii	new species	gram -	aerobic
76	Legionella worsleiensis	new species	gram -	aerobic
77	Leucobacter chironomi	new genus/species	gram +	aerobic
78	Leucobacter denitrificans	new species	gram +	aerobic
79	Lysinibacillus boronitolerans	new species	gram +	aerobic
80	Microbacterium mitrae	new species	gram +	aerobic
81	Micrococcus flavus	new species	gram +	aerobic
82	Micrococcus terreus	new species	gram +	aerobic
83	Nocardia asteroides	new species	gram +	aerobic
84	Ogataea polymorpha	new genus/species	yeast	aerobic
85	Paenibacillus barengoltzii	new species	gram +	microaerophilic
86	Pantoea septica	new species	gram -	aerobic
87	Pichia pseudocactophila	new species	yeast	aerobic
88	Porphyromonas levii	new species	gram -	anaerobic
89	Pseudoxanthomonas kaohsiungensis	new genus/species	gram -	aerobic
90	Rhodobacter aestuarii	new genus/species	gram -	aerobic
91	Rhodobacter veldkampii	new species	gram -	aerobic
92	Saprochaete clavata	new genus/species	yeast	aerobic
93	Saprochaete suaveolens	new species	yeast	aerobic
94	Sarocladium strictum	new genus/species	yeast	aerobic



	New genus/species			
95	Tessaracoccus flavescens	new genus/species	gram +	aerobic
96	Trichosporon coremiiforme	new species	yeast	aerobic
97	Trichosporon dohaense	new species	yeast	aerobic
98	Trichosporon dulciturum	new species	yeast	aerobic
99	Trichosporon faecale	new species	yeast	aerobic
100	Trichosporon gracile	new species	yeast	aerobic
101	Trichosporon japonicum	new species	yeast	aerobic
102	Trichosporon jirovecii	new species	yeast	aerobic
103	Trichosporon laibachii	new species	yeast	aerobic
104	Trichosporon loubieri	new species	yeast	aerobic
105	Trichosporon moniliiforme	new species	yeast	aerobic
106	Trichosporon montevideense	new species	yeast	aerobic
107	Trichosporon terricola	new species	yeast	aerobic
108	Veillonella magna	new species	gram -	anaerobic
109	Virgibacillus proomii	new species	gram +	aerobic
110	Zygosaccharomyces bisporus	new species	yeast	aerobic
111	Zygosaccharomyces florentinus	new species	yeast	aerobic
112	Zygosaccharomyces microellipsoides	new species	yeast	aerobic
113	Zygosaccharomyces rouxii	new species	yeast	aerobic



Deletion of MSP entries:

Deletions	Justification
Arthrobacter castelli DSM 16402T DSM	Many polymers. Will be replaced by new measurement.
Arthrobacter chlorophenolicus DSM 12829T DSM	Many polymers. Will be replaced by new measurement.
Arthrobacter pigmenti DSM 16403T DSM	Many polymers. Will be replaced by new measurement.
Mycobacterium avium ssp avium 8671 VAR	match with M. scrofulaceum
Mycobacterium manitobense DSM 44615 DSM	no valid name; match with M. saskatchewanense
Delftia sp[2] 911600013 LBK	Mixed culture with Saccharomyces servazzii. Saccharomyces servazzii is not in the database so far. The mixed culture was detected during project work with Saccharomyces servazzii.



Renaming of MSP entries:

Old MSP name	to	New MSP name	Justification
Anaerococcus vaginalis DSM 7457T DSM		Anaerococcus hydrogenalis DSM 7454T DSM	DB4613 strains A. hydrogenalis and vaginalis are mixed up. Confirmed by new cultivation, measuring and sequencing. Re-naming of old entry necessary.
Anaerococcus hydrogenalis DSM 7454T DSM		Anaerococcus vaginalis DSM 7457T DSM	DB4613 strains A. hydrogenalis and vaginalis are mixed up. Confirmed by new cultivation, measuring and sequencing. Re-naming of old entry necessary.
Aurantimonas altamirensis 284 RLT		Aureimonas altamirensis 284 RLT	Nomenclature change according DSMZ/IJSEM - changes 11/2011: Aurantimonas altamirensis -> Aureimonas
Candida haemulonii MY916_09 ERL		Candida duobushaemulonii MY916_09 ERL	Species was divided into two separate species
Lactobacillus sp M23 101342 CIP		Carnobacterium maltaromaticum CIP 101342 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102035 CIP		Carnobacterium maltaromaticum CIP 102035 CIP	More precise identification of strain by CIP - adaptation of name
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) 29 PSB		Cryptococcus neoformans 29 PSB	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) ATCC 14116 THL		Cryptococcus neoformans ATCC 14116 THL	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) CCM 8312 CCM		Cryptococcus neoformans CCM 8312 CCM	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species



Old MSP name	to	New MSP name	Justification
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) RV07_02 18 VML		Cryptococcus neoformans RV07_02 18 VML	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Dickeya dieffenbachiae DSM 18013T HAM		Dickeya dadantii ssp dieffenbachiae DSM 18013T HAM	Nomenclature change according DSMZ/IJSEM - changes 07/2012: Dickeya dieffenbachiae -> D. dadantii subsp.dieffenbachiae
Lactobacillus catenaformis O7_085 ANA IBS		Eggerthia catenaformis O7_085 ANA IBS	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corr.) -> Eggerthia
Lactobacillus catenaformis CIP 104817T B CIP		Eggerthia catenaformis CIP 104817T B CIP	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corr.) -> Eggerthia
Lactobacillus catenaformis IBS_MS_39 IBS		Eggerthia catenaformis IBS_MS_39 IBS	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corr.) -> Eggerthia
Lactobacillus catenaformis VA12065_1_11 ERL		Eggerthia catenaformis VA12065_1_11 ERL	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corr.) -> Eggerthia
Bacillus schlegelii CIP 106933T CIP		Hydrogenibacillus schlegelii CIP 106933T CIP	Bacillus schlegelii 31:215 (basonym) ≡ Hydrogenibacillus schlegelii
Lactobacillus vitulinus DSM 20405T DSM		Kandleria vitulina DSM 20405T DSM	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus vitulinus -> Kandleria vitulina
Lactobacillus sp 101810 CIP		Lactobacillus delbrueckii ssp lactis CIP 101810 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp 102006 CIP		Lactobacillus fermentum CIP 102006 CIP	More precise identification of strain by CIP - adaptation of name



Old MSP name	to	New MSP name	Justification
Lactobacillus sp 101909 CIP		Lactobacillus gasseri CIP 101909 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102309 CIP		Lactobacillus paracasei ssp paracasei CIP 102309 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102623 CIP		Lactobacillus rhamnosus CIP 102623 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus zymae 108703 CIP		Lactobacillus zymae CIP 108703T CIP	Typo Strain T was missing - Typo
Tatlockia maceachernii 28 RLT		Legionella maceachernii 28 RLT	New nomenclature not very accepted -> http://www.bacterio.net/t/tatlockia.html#maceachernii
Lactobacillus sp CIP 102166 CIP		Marinilactibacillus psychrotolerans CIP 102166 CIP	more precise identification of strain by CIP - adaptation of name
Mycoplasma argini 7SR10 VLW		Mycoplasma arginini 7SR10 VLW	Typo within species name will be corrected.
Mycoplasma argini NCTC 10129T VLW		Mycoplasma arginini NCTC 10129T VLW	Typo within species name will be corrected.
Mycoplasma ovipneumoniae NCTC 10151T VLW		Mycoplasma gallisepticum NCTC 10115T VLW	DB4613 strains Mycoplasma gallisepticum NCTC 10115T VLW and Mycoplasma ovipneumoniae NCTC 10151T VLW are mixed up. Confirmed by new, measuring and sequencing. Re-naming of old entry necessary.



Old MSP name	to	New MSP name	Justification
Mycoplasma gallisepticum NCTC 10115T VLW		Mycoplasma ovipneumoniae NCTC 10151T VLW	DB4613 strains Mycoplasma gallisepticum NCTC 10115T VLW and Mycoplasma ovipneumoniae NCTC 10151T VLW are mixed up. Confirmed by new, measuring and sequencing. Re-naming of old entry necessary.
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) 991400574 LBK		Saccharomyces cerevisiae 991400574 LBK	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) BJ2168 BRL		Saccharomyces cerevisiae BJ2168 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) DTY3 BRL		Saccharomyces cerevisiae DTY3 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) INVSc1 BRL		Saccharomyces cerevisiae INVSc1 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) Isolat LGL Muenchen		Saccharomyces cerevisiae Isolat LGL Muenchen	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) Kontrollstamm Humanmedizin VML		Saccharomyces cerevisiae Kontrollstamm Humanmedizin VML	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) WS LLH		Saccharomyces cerevisiae WS LLH	adaptation to more common nomenclature
Streptococcus thermophilus 37 RLT		Streptococcus salivarius_ssp_thermophilus 37 RLT	adaptation to current nomenclature



Old MSP name	to	New MSP name	Justification
Streptococcus thermophilus 38 RLT		Streptococcus salivarius_ssp_thermophilus 38 RLT	adaptation to current nomenclature
Streptococcus thermophilus 39 RLT		Streptococcus salivarius_ssp_thermophilus 39 RLT	adaptation to current nomenclature
Streptococcus thermophilus DSM 20259 DSM		Streptococcus salivarius_ssp_thermophilus DSM 20259 DSM	adaptation to current nomenclature
Streptococcus thermophilus DSM 20479 DSM		Streptococcus salivarius_ssp_thermophilus DSM 20479 DSM	adaptation to current nomenclature
Streptococcus thermophilus DSM 8713 DSM		Streptococcus salivarius_ssp_thermophilus DSM 8713 DSM	adaptation to current nomenclature
Listonella anguillarum 02 EGS		Vibrio anguillarum 02 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum 03 EGS		Vibrio anguillarum 03 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum DSM 11323 DSM		Vibrio anguillarum DSM 11323 DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum DSM 21597T DSM		Vibrio anguillarum DSM 21597T DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum LMG 4437T HAM		Vibrio anguillarum LMG 4437T HAM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio



Old MSP name	to	New MSP name	Justification
Listonella anguillarum serotype 02 EGS		Vibrio anguillarum serotype 02 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum serotype 03 EGS		Vibrio anguillarum serotype 03 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella pelagia DSM 21205T DSM		Vibrio pelagius DSM 21205T DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio



Matching Hints Changes:

			Justification
Species alkalescens / argini of the genus Mycoplasma have very similar pattern: Therefore distinguishing their species is difficult.	changed to	Species alkalescens / arginini of the genus Mycoplasma have very similar pattern: Therefore distinguishing their species is difficult.	Correction of typo within the species name "arginini"
Corynebacterium amycolatum Corynebacterium durum Corynebacterium minutissimum	deletion of matching hint link	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Matching hint is not necessary for this species, because MALDI can separate this species from others
Acinetobacter baumannii Acinetobacter calcoaceticus Acinetobacter pittii	change to new matching hint	Member of the Acinetobacter baumannii/calcoaceticus complex. Extraction must be performed to permit reliable species identification.	For clarification change from general matching hint (Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.) to this more precise one.
Arthrobacter castelli Arthrobacter chlorophenicus Arthrobacter pigmenti	link with new matching hint	Species tend to produce polymers which can interfere with the identification.	For clarification the matching hint will inform about the production of possibly interfering polymers.
Burkholderia ambifaria, anthina, cenocepacia, cepacia, diffusa, dolosa, lata, latens, metallica, multivorans, pyrrocina, seminalis, stabilis, vietnamensis	link with new matching hint	Is a member of Burkholderia cepacia complex	Additional information. This will help the user to recognize members of the Burkholderia cepacia complex.



Release Notes for MBT-BDAL-5627 MSP library

What is new in MBT DB 5627?

The MBT database DB 4613 will be extended by **1020 newly generated reference MSP**, which cover **348 species** overall.

21 new genera and **113 new species** will be implemented and which will be covered by **220 new reference entries**.

Additional 800 MSP will improve the coverage of species already contained in the database.

Main improvements:

- Implementation of *Haemophilus haemolyticus* strains
- Additional strains for improvement of *Streptococcus mitis/oralis/pneumonia* identification
- Many new yeast species and strains
- More *Nocardia* strains to cover the MALDI diversity

The 21 new genera and 113 new species cover following groups:

	New genera	New species	Aerobe species	Microaerophile species	Anaerobe species
Gram -	6	23	18	3	2
Gram +	7	27	22	1	4
Yeast	8	63	63	-	-
Fil. Fungi	0	0	-	-	-

Overall improvements:

	New MSP (from 1020)	Species covered	Aerobe strains	Microaerophile strains	Anaerobe strains
Gram -	328	113	265	54	9
Gram +	374	128	185	179	10
Yeast	316	106	106	-	-
Fil. Fungi	2	1	2	-	-

Supported MALDI Biotyper versions (installed on respective operating systems):

- MBT-3.0
- MBT-3.1
- MBT-4.0



Implementation of MSP entries for following 113 new species:

	New genus/species			
1	Actinomyces ruminicola	new species	gram +	anaerobic
2	Aeromicrobium flavum	new genus/species	gram +	aerobic
3	Alicyclobacillus cycloheptanicus	new species	gram +	aerobic
4	Alloiococcus otitis	new genus/species	gram +	aerobic
5	Arthroascus schoenii	new genus/species	yeast	aerobic
6	Arthrobacter flavus	new species	gram +	aerobic
7	Arthrobacter nitroguajacolicus	new species	gram +	aerobic
8	Bacillus altitudinis	new species	gram +	aerobic
9	Bacillus horneckiae	new species	gram +	aerobic
10	Bartonella japonica	new genus/species	gram -	aerobic
11	Brevibacillus invocatus	new species	gram +	aerobic
12	Candida alai	new species	yeast	aerobic
13	Candida ambrosiae	new species	yeast	aerobic
14	Candida auris	new species	yeast	aerobic
15	Candida blattae	new species	yeast	aerobic
16	Candida bohiensis	new species	yeast	aerobic
17	Candida bracarenensis	new species	yeast	aerobic
18	Candida buenavistaensis	new species	yeast	aerobic
19	Candida carpophila	new species	yeast	aerobic
20	Candida castellii	new species	yeast	aerobic
21	Candida duobushaemulonii	new species	yeast	aerobic
22	Candida frijolesensis	new species	yeast	aerobic



	New genus/species			
23	<i>Candida infanticola</i>	new species	yeast	aerobic
24	<i>Candida labiduridarum</i>	new species	yeast	aerobic
25	<i>Candida pseudohaemulonii</i>	new species	yeast	aerobic
26	<i>Candida saitoana</i>	new species	yeast	aerobic
27	<i>Candida shehatae</i> var <i>insectosa</i>	new species	yeast	aerobic
28	<i>Candida sojae</i>	new species	yeast	aerobic
29	<i>Candida viswanathii</i>	new species	yeast	aerobic
30	<i>Candida_citrea</i> [ana] (<i>Pichia_nakasei</i> [teleo]#)	new species	yeast	aerobic
31	<i>Candida_fabianii</i> [ana] (<i>Pichia_fabianii</i> [teleo])	new species	yeast	aerobic
32	<i>Candida_holmii</i> [ana] (<i>Kazachstania_exigua</i> [teleo])	new species	yeast	aerobic
33	<i>Candida_pintolopesii</i> [ana] (<i>Kazachstania_pintolopesii</i> [teleo])	new species	yeast	aerobic
34	<i>Chryseobacterium ginsenosidimutans</i>	new species	gram -	aerobic
35	<i>Chryseobacterium hagamense</i>	new species	gram -	aerobic
36	<i>Cryptococcus albidosimilis</i>	new species	yeast	aerobic
37	<i>Cryptococcus curvatus</i>	new species	yeast	aerobic
38	<i>Cryptococcus diffluens</i>	new species	yeast	aerobic
39	<i>Cryptococcus flavescens</i>	new species	yeast	aerobic
40	<i>Cryptococcus gastricus</i>	new species	yeast	aerobic
41	<i>Cryptococcus humicola</i>	new species	yeast	aerobic
42	<i>Cryptococcus liquefaciens</i>	new species	yeast	aerobic
43	<i>Cryptococcus magnus</i>	new species	yeast	aerobic
44	<i>Cryptococcus saitoi</i>	new species	yeast	aerobic
45	<i>Cryptococcus terreus</i>	new species	yeast	aerobic
46	<i>Cryptococcus uzbekistanensis</i>	new species	yeast	aerobic



	New genus/species			
47	<i>Cryptococcus vishniacii</i>	new species	yeast	aerobic
48	<i>Cryptotrichosporon anacardii</i>	new genus/species	yeast	aerobic
49	<i>Cupriavidus campinensis</i>	new species	gram -	aerobic
50	<i>Cyberlindnera mississippiensis</i>	new genus/species	yeast	aerobic
51	<i>Deinococcus geothermalis</i>	new genus/species	gram +	aerobic
52	<i>Geobacillus thermodenitrificans</i> ssp <i>calidus</i>	new species	gram +	aerobic
53	<i>Geobacillus thermoglucoasidius</i>	new species	gram +	aerobic
54	<i>Gluconacetobacter intermedius</i>	new genus/species	gram -	aerobic
55	<i>Gluconacetobacter liquefaciens</i>	new species	gram -	aerobic
56	<i>Gluconobacter cerinus</i>	new species	gram -	aerobic
57	<i>Guehomyces pullulans</i>	new genus/species	yeast	aerobic
58	<i>Haemophilus haemoglobinophilus</i>	new species	gram -	microaerophilic
59	<i>Haemophilus haemolyticus</i>	new species	gram -	microaerophilic
60	<i>Haemophilus paraphrohaemolyticus</i>	new species	gram -	microaerophilic
61	<i>Halotalea alkalilenta</i>	new genus/species	gram -	aerobic
62	<i>Hannaella luteola</i>	new genus/species	yeast	aerobic
63	<i>Hanseniaspora lachancei</i>	new species	yeast	aerobic
64	<i>Jeotgalicoccus halotolerans</i>	new genus/species	gram +	aerobic
65	<i>Kazachstania bovina</i>	new species	yeast	aerobic
66	<i>Kytococcus schroeteri</i>	new species	gram +	aerobic
67	<i>Lachancea fermentati</i>	new species	yeast	aerobic
68	<i>Lachnoanaerobaculum orale</i>	new genus/species	gram +	anaerobic
69	<i>Lachnoanaerobaculum saburreum</i>	new species	gram +	anaerobic
70	<i>Lachnoanaerobaculum umeaense</i>	new species	gram +	anaerobic



	New genus/species			
71	Laribacter hongkongensis	new genus/species	gram -	aerobic
72	Legionella dresdenensis	new species	gram -	aerobic
73	Legionella geestiana	new species	gram -	aerobic
74	Legionella gratiana	new species	gram -	aerobic
75	Legionella waltersii	new species	gram -	aerobic
76	Legionella worsleiensis	new species	gram -	aerobic
77	Leucobacter chironomi	new genus/species	gram +	aerobic
78	Leucobacter denitrificans	new species	gram +	aerobic
79	Lysinibacillus boronitolerans	new species	gram +	aerobic
80	Microbacterium mitrae	new species	gram +	aerobic
81	Micrococcus flavus	new species	gram +	aerobic
82	Micrococcus terreus	new species	gram +	aerobic
83	Nocardia asteroides	new species	gram +	aerobic
84	Ogataea polymorpha	new genus/species	yeast	aerobic
85	Paenibacillus barengoltzii	new species	gram +	microaerophilic
86	Pantoea septica	new species	gram -	aerobic
87	Pichia pseudocactophila	new species	yeast	aerobic
88	Porphyromonas levii	new species	gram -	anaerobic
89	Pseudoxanthomonas kaohsiungensis	new genus/species	gram -	aerobic
90	Rhodobacter aestuarii	new genus/species	gram -	aerobic
91	Rhodobacter veldkampii	new species	gram -	aerobic
92	Saprochaete clavata	new genus/species	yeast	aerobic
93	Saprochaete suaveolens	new species	yeast	aerobic
94	Sarocladium strictum	new genus/species	yeast	aerobic



	New genus/species			
95	Tessaracoccus flavescens	new genus/species	gram +	aerobic
96	Trichosporon coremiiforme	new species	yeast	aerobic
97	Trichosporon dohaense	new species	yeast	aerobic
98	Trichosporon dulcitum	new species	yeast	aerobic
99	Trichosporon faecale	new species	yeast	aerobic
100	Trichosporon gracile	new species	yeast	aerobic
101	Trichosporon japonicum	new species	yeast	aerobic
102	Trichosporon jirovecii	new species	yeast	aerobic
103	Trichosporon laibachii	new species	yeast	aerobic
104	Trichosporon loubieri	new species	yeast	aerobic
105	Trichosporon moniliiforme	new species	yeast	aerobic
106	Trichosporon montevideense	new species	yeast	aerobic
107	Trichosporon terricola	new species	yeast	aerobic
108	Veillonella magna	new species	gram -	anaerobic
109	Virgibacillus proomii	new species	gram +	aerobic
110	Zygosaccharomyces bisporus	new species	yeast	aerobic
111	Zygosaccharomyces florentinus	new species	yeast	aerobic
112	Zygosaccharomyces microellipsoides	new species	yeast	aerobic
113	Zygosaccharomyces rouxii	new species	yeast	aerobic



Deletion of MSP entries:

Deletions	Justification
Arthrobacter castelli DSM 16402T DSM	Many polymers. Will be replaced by new measurement.
Arthrobacter chlorophenolicus DSM 12829T DSM	Many polymers. Will be replaced by new measurement.
Arthrobacter pigmenti DSM 16403T DSM	Many polymers. Will be replaced by new measurement.
Mycobacterium avium ssp avium 8671 VAR	match with M. scrofulaceum
Mycobacterium manitobense DSM 44615 DSM	no valid name; match with M. saskatchewanense
Delftia sp[2] 911600013 LBK	Mixed culture with Saccharomyces servazzii. Saccharomyces servazzii is not in the database so far. The mixed culture was detected during project work with Saccharomyces servazzii.



Renaming of MSP entries:

Old MSP name	to	New MSP name	Justification
Anaerococcus vaginalis DSM 7457T DSM		Anaerococcus hydrogenalis DSM 7454T DSM	DB4613 strains A. hydrogenalis and vaginalis are mixed up. Confirmed by new cultivation, measuring and sequencing. Re-naming of old entry necessary.
Anaerococcus hydrogenalis DSM 7454T DSM		Anaerococcus vaginalis DSM 7457T DSM	DB4613 strains A. hydrogenalis and vaginalis are mixed up. Confirmed by new cultivation, measuring and sequencing. Re-naming of old entry necessary.
Aurantimonas altamirensis 284 RLT		Aureimonas altamirensis 284 RLT	Nomenclature change according DSMZ/IJSEM - changes 11/2011: Aurantimonas altamirensis -> Aureimonas
Candida haemulonii MY916_09 ERL		Candida duobushaemulonii MY916_09 ERL	Species was divided into two separate species
Lactobacillus sp M23 101342 CIP		Carnobacterium maltaromaticum CIP 101342 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102035 CIP		Carnobacterium maltaromaticum CIP 102035 CIP	More precise identification of strain by CIP - adaptation of name
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) 29 PSB		Cryptococcus neoformans 29 PSB	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) ATCC 14116 THL		Cryptococcus neoformans ATCC 14116 THL	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) CCM 8312 CCM		Cryptococcus neoformans CCM 8312 CCM	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species



Old MSP name	to	New MSP name	Justification
Cryptococcus_neoformans[ana]# (Filobasidiella_neoformans[teleo]) RV07_02 18 VML		Cryptococcus neoformans RV07_02 18 VML	Anamorph/telemorph nomenclature is not very accepted - will be changed for this species
Dickeya dieffenbachiae DSM 18013T HAM		Dickeya dadantii ssp dieffenbachiae DSM 18013T HAM	Nomenclature change according DSMZ/IJSEM - changes 07/2012: Dickeya dieffenbachiae -> D. dadantii subsp.dieffenbachiae
Lactobacillus catenaformis O7_085 ANA IBS		Eggerthia catenaformis O7_085 ANA IBS	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corrig.) -> Eggerthia
Lactobacillus catenaformis CIP 104817T B CIP		Eggerthia catenaformis CIP 104817T B CIP	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corrig.) -> Eggerthia
Lactobacillus catenaformis IBS_MS_39 IBS		Eggerthia catenaformis IBS_MS_39 IBS	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corrig.) -> Eggerthia
Lactobacillus catenaformis VA12065_1_11 ERL		Eggerthia catenaformis VA12065_1_11 ERL	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus catenaformis (corrig.) -> Eggerthia
Bacillus schlegelii CIP 106933T CIP		Hydrogenibacillus schlegelii CIP 106933T CIP	Bacillus schlegelii 31:215 (basonym) ≡ Hydrogenibacillus schlegelii
Lactobacillus vitulinus DSM 20405T DSM		Kandleria vitulina DSM 20405T DSM	Nomenclature change according DSMZ/IJSEM - changes 10/2011: Lactobacillus vitulinus -> Kandleria vitulina
Lactobacillus sp 101810 CIP		Lactobacillus delbrueckii ssp lactis CIP 101810 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp 102006 CIP		Lactobacillus fermentum CIP 102006 CIP	More precise identification of strain by CIP - adaptation of name



Old MSP name	to	New MSP name	Justification
Lactobacillus sp 101909 CIP		Lactobacillus gasseri CIP 101909 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102309 CIP		Lactobacillus paracasei ssp paracasei CIP 102309 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus sp CIP 102623 CIP		Lactobacillus rhamnosus CIP 102623 CIP	More precise identification of strain by CIP - adaptation of name
Lactobacillus zymae 108703 CIP		Lactobacillus zymae CIP 108703T CIP	Type Strain T was missing - Typo
Tatlockia maceachernii 28 RLT		Legionella maceachernii 28 RLT	New nomenclature not very accepted -> http://www.bacterio.net/t/tatlockia.html#maceachernii
Lactobacillus sp CIP 102166 CIP		Marinilactibacillus psychrotolerans CIP 102166 CIP	more precise identification of strain by CIP - adaptation of name
Mycoplasma argini 7SR10 VLW		Mycoplasma arginini 7SR10 VLW	Typo within species name will be corrected.
Mycoplasma argini NCTC 10129T VLW		Mycoplasma arginini NCTC 10129T VLW	Typo within species name will be corrected.
Mycoplasma ovipneumoniae NCTC 10151T VLW		Mycoplasma gallisepticum NCTC 10115T VLW	DB4613 strains Mycoplasma gallisepticum NCTC 10115T VLW and Mycoplasma ovipneumoniae NCTC 10151T VLW are mixed up. Confirmed by new, measuring and sequencing. Re-naming of old entry necessary.



Old MSP name	to	New MSP name	Justification
Mycoplasma gallisepticum NCTC 10115T VLW		Mycoplasma ovipneumoniae NCTC 10151T VLW	DB4613 strains Mycoplasma gallisepticum NCTC 10115T VLW and Mycoplasma ovipneumoniae NCTC 10151T VLW are mixed up. Confirmed by new, measuring and sequencing. Re-naming of old entry necessary.
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) 991400574 LBK		Saccharomyces cerevisiae 991400574 LBK	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) BJ2168 BRL		Saccharomyces cerevisiae BJ2168 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) DTY3 BRL		Saccharomyces cerevisiae DTY3 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) INVSc1 BRL		Saccharomyces cerevisiae INVSc1 BRL	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) Isolat LGL Muenchen		Saccharomyces cerevisiae Isolat LGL Muenchen	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) Kontrollstamm Humanmedizin VML		Saccharomyces cerevisiae Kontrollstamm Humanmedizin VML	adaptation to more common nomenclature
Candida_robusta[ana]# (Saccharomyces_cerevisiae[teleo]) WS LLH		Saccharomyces cerevisiae WS LLH	adaptation to more common nomenclature
Streptococcus thermophilus 37 RLT		Streptococcus salivarius_ssp_thermophilus 37 RLT	adaptation to current nomenclature



Old MSP name	to	New MSP name	Justification
Streptococcus thermophilus 38 RLT		Streptococcus salivarius_ssp_thermophilus 38 RLT	adaptation to current nomenclature
Streptococcus thermophilus 39 RLT		Streptococcus salivarius_ssp_thermophilus 39 RLT	adaptation to current nomenclature
Streptococcus thermophilus DSM 20259 DSM		Streptococcus salivarius_ssp_thermophilus DSM 20259 DSM	adaptation to current nomenclature
Streptococcus thermophilus DSM 20479 DSM		Streptococcus salivarius_ssp_thermophilus DSM 20479 DSM	adaptation to current nomenclature
Streptococcus thermophilus DSM 8713 DSM		Streptococcus salivarius_ssp_thermophilus DSM 8713 DSM	adaptation to current nomenclature
Listonella anguillarum 02 EGS		Vibrio anguillarum 02 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum 03 EGS		Vibrio anguillarum 03 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum DSM 11323 DSM		Vibrio anguillarum DSM 11323 DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum DSM 21597T DSM		Vibrio anguillarum DSM 21597T DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum LMG 4437T HAM		Vibrio anguillarum LMG 4437T HAM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio



Old MSP name	to	New MSP name	Justification
Listonella anguillarum serotype 02 EGS		Vibrio anguillarum serotype 02 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella anguillarum serotype 03 EGS		Vibrio anguillarum serotype 03 EGS	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio
Listonella pelagia DSM 21205T DSM		Vibrio pelagius DSM 21205T DSM	nomenclature change according DSMZ/IJSEM - changes 12/2011: Listonella => Vibrio



Matching Hints Changes:

			Justification
Species alkalescens / argini of the genus Mycoplasma have very similar pattern: Therefore distinguishing their species is difficult.	changed to	Species alkalescens / arginini of the genus Mycoplasma have very similar pattern: Therefore distinguishing their species is difficult.	Correction of typo within the species name "arginini"
Corynebacterium amycolatum Corynebacterium durum Corynebacterium minutissimum	deletion of matching hint link	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Matching hint is not necessary for this species, because MALDI can separate this species from others
Acinetobacter baumannii Acinetobacter calcoaceticus Acinetobacter pittii	change to new matching hint	Member of the Acinetobacter baumannii/calcoaceticus complex. Extraction must be performed to permit reliable species identification.	For clarification change from general matching hint (Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.) to this more precise one.
Arthrobacter castelli Arthrobacter chlorophenicus Arthrobacter pigmenti	link with new matching hint	Species tend to produce polymers which can interfere with the identification.	For clarification the matching hint will inform about the production of possibly interfering polymers.
Burkholderia ambifaria, anthina, cenocepacia, cepacia, diffusa, dolosa, lata, latens, metallica, multivorans, pyrrocina, seminalis, stabilis, vietnamensis	link with new matching hint	Is a member of Burkholderia cepacia complex	Additional information. This will help the user to recognize members of the Burkholderia cepacia complex.



Product Description

microflex™ LRF

MALDI Time-of-Flight Mass Spectrometer System

BDAL # 8601800

Small footprint, silently operating MALDI-TOF mass spectrometer with industrial hardened vertical design including self-diagnostics of major components.

A high performance bench-top Time-of-Flight Mass Spectrometer equipped with MALDI Ion Source, Linear/Reflectron mode, high dynamic range FlashDetector™ systems in conjunction with PAN™ wide mass range focusing for unrivaled resolution power and mass accuracy. Additional PSD-MS/MS Capability included.

- Equipped with MALDI-Ion Source
- TOF-Analyzer for Linear and Reflectron mode
- TOF-Analyzer for both positive and negative ion mode
- Resolution performance $RP > 15,000$ FWHM for peptides
- ppm mass accuracy (internal/external calibration) in both modes

A. MALDI-Ion Source

- Highly sensitive microScout™ ion source; medium area target (54mm x 36mm) with exact dimensions of 1/4 microtiter plate
- 2nd generation proprietary PAN™ pulsed ion extraction technology for high mass accuracy and unmatched resolution spectra across an extended mass range
- 60 Hz N₂-Cartridge-Laser including variable power attenuator and UV optics
- microScout™ Target Kit including: two MSP96 ground steel targets, one target each of MSP96 polished steel, MSP AnchorChip 600/96 and NALDI targets (5 plates)
- Manually operated sample inlet
- High resolution magnifying target observation optics with integrated display in Compass™ acquisition software

B. High Performance Time-of-Flight Mass Analyzer

- Ultra-stable electronics for TOF analyser, detector and ion source fully enables a 1-60 Hz data acquisition rate in MS and MS/MS operation
- Integrated pumping system including vacuum measurement and control unit: 70 l/sec turbo-molecular pump including diaphragm-pump
- FlashDetector™ providing unmatched mass resolution and mass accuracy



C. Gridless ion reflectron for increased sensitivity, resolution and accuracy:

- Gridless two stage ion reflectron for superb mass resolution
- Independent reflector power supply with high-precision control
- Includes FlashDetector™ and power supplies

D. FAST™ (Fragmentation Analysis and Structural TOF-MS) accessory for PSD (Post-Source Decay) MS/MS-experiments:

- autoFAST™ software for calibration, pasting of segments, etc.
- FAST-FILTER Pre Cursor Ion Selector for true MS/MS of complex sample mixtures

E. Data system:

- PC Workstation with 2.66 GHz Single-CPU-Quad-Core-Xeon-processor, 12 GB RAM, 2 TB hard disk, Ethernet connection for external networks, ≥ 2 Gs/s Digitizer
- DVD-ROM drive
- R/W DVD-ROM drive DL
- ≥ 24" flat screen colour monitor
- Windows™ WIN7 operating system
- Laser printer
- Remote Service capability via 128-bit SSL-security web connection

F. Applications Software:

Software package Compass 1.4 for flex-series instruments including instrument control, data acquisition, post processing, and data analysis packages:

- Released for Windows™ XP and Windows™ WIN7 operating systems
- Compass / flexControl 3.4 for integrated control of the instrument
- Compass / flexAnalysis 3.4, including:
 - Advanced data processing with a high degree of automation
 - Usage of the new calibration algorithms and High-Precision-Calibration (HPC)
 - Easy export of peak list (e.g. to MS Excel)
 - Interface to bioinformatics software packages as BioTools™ and ProteinScope™
- Compass / AutoXecute™ with fuzzy-logic optimization for automated acquisition



- SW License Compass / flexControl 3.4
- SW Licence Compass / flexAnalysis 3.4
- SW Licence Compass / AutoXecute
- SW License TLC-MALDI 1.0 as an option

G. Set of manuals and reference DVDs

H. Installation

I. Familiarization upon installation

J. 1 year warranty

K. Voucher for a factory-training course - valid for 2 participants.

● **micrOTOF-Q III**

Description

Part Number: # 728889

Benchtop easy-to-use, high-performance electrospray ionization quadrupole time-of-flight LC/MS/MS mass spectrometer designed for exact mass and true isotopic measurements

Small footprint system enclosure for ESI ion funnel source, quadrupole, collision cell, oa-TOF mass analyzer, electronics, and vacuum pumps, containing:

A. Apollo II Electrospray Ionization Source with ion funnel

- ESI source with grounded needle for safety and easy sample introduction
- Heated counter current drying gas for gentle and efficient drying
- Pneumatic off-axis nebulization for flow rates up to 1 ml/min., with gradients from 100% aqueous to 100% organic
- Patented dual ion funnel for mass independent ion transfer
- Ion lens system including in-source collision induced dissociation control (IS-CID)
- Combined Funnel-Hexapole-Cartridge with front access for easy maintenance
- Source HV controller and drying gas controller
- Ion lens housing and vacuum system
- Flow rate with ESI-Source 1µl/min ... 1ml/min
- Suitable for HPLC and CE coupling

B. High Mass Quadrupole Mass Filter:

- Hyperbolic quadrupole mass filter
- Ultra stable monolithic design
- RF-generator for monoisotopic precursor ion selection

C. CID Collision cell:

- Hyperbolic quadrupole design
- RF-Generator with fast amplitude switching
- Collision gas controller

D. Focus™ ion optics

- Ultra precise orthogonal ion beam focusing
- Time controlled ion extraction

E. Orthogonal pulsed ion extraction interface

- Interface housing and ion lens system
- In-line detector system for easy maintenance
- Ultra-stable high voltage switches with up to 20 kHz repetition rate and appropriate power supplies.

F. Time-of-Flight (TOF) Mass Analyzer:

- Reflectron TOF analyzer with orthogonal mounted ion source
- High-sensitivity and fast ion detector system
- Positive and negative ion modes
- Ultra-stable high voltage power supplies for TOF analyzer and detector

G. Electrostatic ion reflectron

- Ion reflectron for increased mass resolution and accuracy
- Includes ion reflectron electronics

H. Vacuum system

- Q-q-TOF analyzer vacuum housing
- Vacuum system with 5 differential pumping stages
- One roughing pump and quadruple stage turbo-drag pumps for ESI source and Q-q-TOF analyzer
- Vacuum measurement and pump control unit

I. Syringe pump

J. Modes of Operation

- TOF Mass ranges 20-80,000 m/z
- FSR resolution 20,000 FWHM
- Internal calibration
- External calibration
- Exact mass measurements independent from sample concentration over a wide dynamic range without second sprayer.
- Dual sprayer (as an option)

K. High-performance and accurate mass features

- Patented ion funnel source
- Superior MS/MS sensitivity
- Long term and ultra stable mass axis stability in MS and MS/MS
- Exact mass independent from sample concentration and collision energy
- Combined calibration for both MS and MS/MS
- TIP™ technology for True-Isotopic-Pattern in MS and MS/MS
- SmartFormula 3D™ for multi-dimensional, unambiguous determination of molecular formula with sub-ppm confidence
- Wide dynamic range for quantitation
- Temperature compensated flight tube
- Max scan rate 40 spectra / sec
- Positive / negative ion operation

L. Data system OTOF:

- 2 GSamples/sec, 16 Gbit/sec sampling rate Digitizer
- 2,66 GHz Quad-Core Processor, 12 GB RAM, system hard-disk drive plus 2 TB hard-disk drive for data
- DVD-ROM drive
- DVD-ROM +/- R/W drive
- 24" LCD flat screen color monitor

- OS WinXP / Win7/32
- Laser printer
- Remote Service capability via 128-bit SSL-security web connection

M. Applications software

Software package Compass 1.5 for HPLC and micrOTOF-Q control, data acquisition, post processing, and data analysis:

- Operating system Windows XP / Win 7/32
- HyStar 3.2 for integrated control of most popular HPLC systems and autosamplers
- OTOFcontrol 3.2 software with smart and expert mode
 - Expert mode: extended control over instrument parameters for interactive system optimization of sophisticated exact mass methods
- Data Analysis software DA 4.1, including:
 - Advanced data processing with a high degree of automation
 - SmartFormula 3D™: Automated sum formula determination using MS and MS/MS data with both, accurate mass and isotopic fit.
 - CompoundCrawler™ for web-based searched of molecular structures
 - FragmentExplorer™ for full annotation of MS/MS spectra
 - QuantAnalysis™ quantitation package
 - LibrarySearch™ module for search of MS, MS/MS and MSn spectra with advanced matching algorithm
 - Charge deconvolution module
 - MaxEntropy Deconvolution as an option
 - Export of spectra and ion current profiles as Windows metafiles to word-processing programs
- SW License Compass 1.5
- SW License Charge Deconvolution for DA4.1

N. Set of manuals and reference CD-ROMs

O. Installation

P. Familiarization upon installation


Q. 1 year warranty

Voucher for a factory-training course - valid for 2 participants.



Petrochemical Gas Chromatographs


● Analyzer Solutions



Bruker has the experience and know-how to provide pre-configured gas chromatographs that are ready at power up to handle your key applications. Our long experience in designing, configuring and manufacturing complete systems, with all their analytical benefits, ensures you get the solution that's right for you. With a host of standard solutions configured to meet the performance specifications outlined in international methods, and the capability to produce unique, tailor-made solutions, we have the answer that you seek.

Bruker Analyzer Solutions

Benefits



Bruker configures and tests GC hardware and software according to widely used industry standard methods (e.g. ASTM, UOP, EN, ISO, GPA,...), to save its clients time and to ensure confidence in results. Solutions are configured to meet the performance specifications outlined in the set method itself.

Included with all Bruker Analyzer solutions:

- All Hardware
- Software (incl. special "plug-ins" where appropriate)
- Pre-Installed methods
- Test Chromatograms
- Installation/ Validation Data
- Trouble Shooting Guide
- User documentation customized for the specific method

Bruker Simulated Distillation Analyzers




Table 1: Overview Simdist methods.

A gas chromatographic (GC) technique, Simulated Distillation (SimDist) reproduces the physical distillation of petroleum materials and products by determining boiling point distribution. Bruker's range of Simulated Distillation Analyzers are designed to meet all industry standard methods, Bruker's analyzer software includes both ASTM D86 and ASTM D1160 correlations. Bruker's highly automated GC, CompassCDS Chromatography Data Handling Software, and integrated SimDist software are also designed to meet worldwide industry standard test methods.



436-GC with Sampler

Key Benefits Include:

- Accurate boiling point distribution up to 750°C
- Integrated standard test methods, applications fully comply with ASTM, IP, DIN and ISO standard test methods
- Complete, single vendor solution
- Complete control from initial setup to final report
- ASTM D86 and ASTM D1160 correlation

Built-in Reports:

Bruker's SimDist software provides a wide variety of report options to meet specific requirements including;

- Chromatogram with merged corrected blank analysis and IBP/FBP marks versus retention time
- Boiling point versus percentage of sample
- Table and plot with retention time versus boiling point
- D86 and D1160 correlations
- DIN Noak and motor oil volatility reports
- Table with cut points and fractions plus residue analysis with recovery calculation up to C120

Hydrocarbon Analysis by Group (PIONA+™)




Figure 1: Traps are easily accessible and do not require any tools to install or replace.

Characterization of Engine Fuels by Hydrocarbon Group Type

Bruker's PIONA+™ Analyzer is a highly flexible GC analysis platform to obtain comprehensive characterization and quantitative information, including hydrocarbon group types, oxygenates and carbon number distribution for spark ignition engine fuels.

Key Analyzer Capabilities:

- Unparalleled operational flexibility
- Compliant with established standards
- A complete and fully integrated solution
- A powerful analyzer, easy to use, generating outstanding analysis results day after day

The PIONA+ Analyzer performs a complete analysis (as described in ASTM D6839 and similar methods) and provides unprecedented analytical flexibility and simplified operation through the use of a novel approach to column/trap heating and exchange (Figure 1).




Figure 2: Chromatogram of a test mix using conventional PIONA mode. (Analysis Time - 180mins)

Reduced Analysis Time and Increased Sample Throughput Efficiency through the Use of “Concurrent” Heating

A unique aspect of the design of the Bruker PIONA+ system is the ability to independently heat the individual traps.

This has 2 major operational benefits:

- Enhanced elution integrity for wide range sample
- Significantly Reduced Analysis Time

The independent and concurrent heating design permits greater trap control and benefits in improved elution integrity of the component groups e.g naphthene, iso-paraffins and n-paraffins even for a wide range sample (C4 - C11). In addition, only a single Molsieve column temperature cycle is employed thus reducing the analysis time by almost half allowing a “fast” PIONA mode of operation (see figures 2 and 3).

By employing this technique, sample throughput can be nearly doubled compared to systems that do not offer this unique capability.




Figure 3: Chromatogram of a test mix using concurrent heating in “fast” PIONA mode. (Analysis Time = 95mins)

Determining Total Olefin Content Is Now Practical

The stability, sample loading and lifetime for all of the critical chromatographic components have been improved and optimized in the Bruker PIONA+ Analyzer. Of special and particular note is the improvement in increasing the sample loading capacity of the "olefin" trap. As a result, it is now possible to analyze streams with olefin content as high as 35-40% or more. This makes it practical to employ a single analytical method to obtain total olefin content (Table 2).

Table 2

Saturates				Unsaturates			Aromatics	Oxygenates	Total
Carbon	Cyclic	Iso	Normal	Cyclic	Iso	Normal			
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.06	0.00	0.03	0.54	0.00	0.00	0.63
5	0.31	11.37	2.98	0.87	9.92	7.58	0.00	0.00	33.03
6	3.19	9.98	1.40	2.40	8.59	4.40	1.72	0.00	31.68
7	4.31	6.77	0.00	2.14	4.76	1.91	7.47	0.00	27.36
8	1.42	3.12	0.00	0.39	2.06	0.00	0.07	0.00	7.06
9	0.01	0.00	0.00	0.04	0.02	0.00	0.03	0.00	0.10
10	0.01	0.00	0.00	0.09	0.00	0.03	0.01	0.00	0.14
11	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Total	9.25	31.24	4.44	5.93	25.38	14.46	9.30	0.00	100.00

Example weight % report for a Naphtha sample with high (46%) olefinic content; (highlighted in blue).

The analyzer design allows the operator to conduct analyses in any one of a number of different operational modes including PNA, PONA, PIONA, O-PONA and O-PIONA in standard and concurrent heating configuration. The system is compliant with established standard methods (see adjacent array).

	PNA	PONA	PIONA	O-PONA	O-PIONA
EN 14517					✓
EN-ISO-22854					✓
ASTM D6839				✓	
DIN 51448-2			✓		
ASTM D1319 (FIA)		✓			
DIN 51448-1	✓				
ASTM D5443	✓				
UOP 870	✓				
IP 382	✓				



Compliant with the method



More information generated than required for the method

Detailed Hydrocarbon Analyzer

The DHA Analyzer is a complete high resolution GC solution for the analysis of hydrocarbons in petroleum streams. It is capable of performing all of the standard methods including ASTM D6729, D6730, D6733, D5134, D6623, IP 344/ DHA "Front End" and "Fast DHA".

Although each DHA analyzer is configured, tested and certified at the factory for a standard method specified by the customer, the DHA software permits the operator to utilize any of the other popular standard methods as well. And, because of the outstanding performance and flexibility of the Bruker GC and CompassCDS software design, Bruker is able to quickly modify the existing methods or add new ones if required as a result of the on-going dynamic industry standard processes.




Figure 4: The analysis of permanent gases and hydrogen using the Rapid RGA.




Figure 5: Report selection output.

Brüker Refinery Gas Analyzers

Peak Identification

1. Hydrogen
2. Carbon Dioxide
3. Hydrogen Sulfide
4. Oxygen
5. Nitrogen
6. Carbon Monoxide
7. Methane
8. Ethane
9. Ethylene
10. Propane
11. Cyclo Propane
12. Propylene
13. i-Butane
14. n-Butane
15. Propadiene
16. Acetylene
17. t-2-Butene
18. i-Butene
19. c-2-Butene
20. i-Pentane
21. n-Pentane
22. 1, 3-Butadiene
23. Propyne
24. Butyne
25. C6+
26. Helium




Figure 6: The separation of light hydrocarbons using the Standard RGA.

The source and composition of refinery gases varies considerably. Measuring gas composition precisely and accurately is a significant challenge in today's refinery operations. Brüker's Refinery Gas Analyzers (RGA) are designed to deliver superior, reliable results for a wide range of sources and analysis throughput requirements.

Key Benefits:

- Pre-configured and tuned
- Standard methods including UOP 539, DIN-51666 and ASTM D2163
- Integrated micro-gasifier ensures complete vaporization of LPGs and high pressure samples to prevent sample discrimination (option)
- Multi-channel approach




Figure 7: The analysis of light hydrocarbons using the Rapid RGA, with complete separation in less than five minutes.




Figure 8: The analysis of permanent gases and hydrogen using the Rapid RGA.

Brüker Offers Two RGA Systems to Meet the Widest Range of Analysis Requirements:

Standard RGA:

A three channel multi-valve design using both capillary and packed columns.

Channel 1 - Analysis of permanent gases

Channel 2 - Light hydrocarbons

Channel 3 - Hydrogen.

Total analysis time for all components in 25 minutes.

Rapid RGA:

The Standard RGA packed columns in the hydrogen and permanent gas channels are replaced by micro packed columns and installed in a separate column oven. Key benefits of this design are:

- Flexibility
- Reduced Analysis Time - 5mins (with H₂S - 7mins)
- Increased Sample Throughput

Characteristics	Standard RGA	Rapid RGA
No. of Channels/Detectors Used	3	3
No. of Column Ovens	1	2
Analysis Time	25 min	5 min (7 min with H ₂ S)
Repeatability	<1%	<1%
Linear Bench Space Required	66 cm/26 in.	66 cm/26 in.
Minimum Component Detection Level	0.01% all components except H ₂ S = 0.05%	0.01% all components except H ₂ S = 0.05%
Suitability		
Typical Refinery Gas	Excellent	Excellent
Impurities in Bulk Ethylene	Excellent	Excellent
Impurities in Bulk Propylene	Excellent	Good
Impurities in Bulk C4	Good	Good

Table 3: Standard RGA vs Rapid RGA.

Low Level Oxygenates Analyzer




Figure 9: Typical chromatogram showing a wide range analysis of a liquid sample stream.

The determination of sub to high ppm levels of ethers, alcohols, aldehydes and ketones in different hydrocarbon matrices is a recurring challenge in the petroleum refining and petrochemical industry. The Bruker Low Level Oxygenates Analyzer is an easy to use solution to meet this challenge and is according ASTM D7423.

The Low Level Oxygenates Analyzer is designed and optimized to quantify ppm and sub ppm levels of ethers (e.g. DME, MTBE, ETBE, DIPE), alcohols (e.g. methanol, ethanol, propanol), ketones (e.g. acetone, MEK) and aldehydes in various hydrocarbon matrices. In general, all oxygenated components with a boiling point of up to 100°C can be analyzed and the sample can be a gas, LPG or liquid under ambient conditions with a final boiling point up to 250°C.

The system is comprised of a Bruker GC configured with gas and liquid sampling valves, two high performance capillary analysis columns, digitally controlled pneumatics including a 'fluidic' switch and Flame Ionization Detector (FID). An optional 'pressure station' can be added to eliminate the possibility of losing sample due to evaporation when analyzing LPG. The GC is controlled via the CompassCDS Chromatography Data Handling Software, which acquires data, processes it and generates analyses reports.

Peak Identification

1. Diethylether
2. Acetaldehyde
3. Ethyl tert. Butyl ether
4. Methyl tert. Butyl ether
5. Diisopropylether
6. Propanal
7. tert amyl methyl ether
8. Propylether
9. Isobutyraldehyde
10. Butyraldehyde
11. Methanol
12. Acetone
13. Isovaleraldehyde
14. Valeraldehyde
15. 2-Butanone
16. Ethanol
17. 1-Propanol
18. tert Butyl alcohol
& Isobutanol
19. 1-Butanol

Bruker 4815 GC Oxygenates Analyzer

The Bruker 4815 GC Analyzer provides a highly cost effective solution for the analysis of oxygenates in gasoline, according to the widely used industry standard method ASTM D4815. The combination of Bruker's reliable GC hardware, powerful software and industry leading pre- and post-sales support teams make this analyzer package the most comprehensive solution available today.

Oxygenated compounds can be present in various hydrocarbon matrices either because they were purposely added (e.g. into gasoline), because they are naturally present, or formed during catalytic processes such as polymer production. In gasoline, oxygenated compounds are added as 'anti-knock' agents to increase the octane number and decrease emissions by replacing organo-lead compounds.

The type and concentration of oxygenated compounds must be measured in reformulated gasolines as part of ongoing product quality assessment, and to confirm the oxygenated components have been added in the correct amounts according to regulatory requirements (e.g. California Air Resources Board).

ASTM D4815 is frequently chosen as the standard method for the determination of oxygenated compounds. Individual ethers and alcohols are quantified in gasoline including: MTBE, ETBE, TAME, DIPE, C1-C4 alcohols and tert-amylalcohol. Individual ether components are measured from 0.1 to 20.0 mass %. The individual alcohols are measured from 0.1 to 12.0 mass %.



456-GC with Sampler




Figure 10 : Typical Chromatogram of the test sample.

Trace Impurity Analyzers

Sulfur Components in LPGs

Low level analysis of sulfur containing components such as H₂S, COS and mercaptanes is extremely challenging and a configured GC offers the solution.

Firstly, the system employs a micro-gasifier enabling the direct coupling of an LPG stream. Secondly, an inert steel sample path ensures a trouble free analysis of sulfur containing components at low concentrations. Finally, a two channel PFPD/ two column approach permits the analysis of all components of interest in one run whatever the LPG matrix. Two differing columns ensures quenching of PFPD signal by the matrix is overcome and full sulfur component analysis is achieved. Figures 11 and 12 show chromatograms obtained in a propane matrix and illustrates the novel benefits of the 2 channel approach.

Permanent Gases in LPGs

Impurities such as CO, CO₂, H₂, O₂ and N₂ need to be determined at low levels in LPGs. Complete separation of these components is done using a two channel single detector (PDHID) system. The GC employs a permanent gas channel for analyzing H₂, O₂, N₂, CH₄ and CO and a specific second channel for the analysis of CO₂. A gasifier is used as a sample introduction device thus giving the capability of handling LPG samples C2 through C4.

Detection limits are at the ppb level (Figure 13), depending on the component of interest.




Figure 11: Sulfur components in propane, BR-1 column.




Figure 12: Sulfur components in propane, BR-Q PLOT.




Figure 13: LPG sample.

Total Characterization of Ethylene Impurities

For a total characterization of impurities in ethylene and also propylene six GC channels are required. By coupling two Bruker GCs with three channels each, a comprehensive solution is available for analyzing these components. The channels used in this analyzer are analytical tools principally developed for the determination of different gases in various hydrocarbon types of gaseous matrices.

GC-1

- H₂ Channel (TCD)
- O₂/N₂ Channel (TCD)
- CO, CO₂ Channel (Methanizer/FID)

GC-2

- Light Hydrocarbon Channel (FID)
- Oxygenates Channel (FID)
- Sulfur Channel (PFPD)

The results (see figures 14 to 19) demonstrate how this 6 channel system is perfectly suited for the total characterization of ethylene and its impurities.




Figure 14: CO, CH₄ and CO₂ on GC-1.




Figure 15: H₂ on GC-1.




Figure 16: O₂ and N₂ on GC-1.




Figure 17: Light hydrocarbons on GC-2.




Figure 18: Oxygenates on GC-2.





Figure 19: Sulfur components on GC-2.

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications (like ASTM, UOP, ISO, GPA and EN) and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- "Quick connect" fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time




For research use only. Not for use in diagnostic procedures.



www.ScionHasArrived.com



www.GlobalEnergyTesting.com




Petrochemical Gas Chromatographs

● Analyzer Solutions

Innovation with Integrity


Solutions



Bruker has the experience and know-how to provide pre-configured gas chromatographs that are ready at power up to handle your key applications. Our long experience in designing, configuring and manufacturing complete systems, with all their analytical benefits, ensures you get the solution that's right for you. With a host of standard solutions configured to meet the performance specifications outlined in international methods, and the capability to produce unique, tailor-made solutions, we have the answer that you seek.

Bruker Analyzer Solutions

Benefits



Bruker configures and tests GC hardware and software according to widely used industry standard methods (e.g. ASTM, UOP, EN, ISO, GPA,...), to save its clients time and to ensure confidence in results. Solutions are configured to meet the performance specifications outlined in the set method itself.

Included with all Bruker Analyzer solutions:

- All Hardware
- Software (incl. special "plug-ins" where appropriate)
- Pre-Installed methods
- Test Chromatograms
- Installation/ Validation Data
- Trouble Shooting Guide
- User documentation customized for the specific method

Bruker Simulated Distillation Analyzers




Table 1: Overview Simdist methods.

A gas chromatographic (GC) technique, Simulated Distillation (SimDist) reproduces the physical distillation of petroleum materials and products by determining boiling point distribution. Bruker's range of Simulated Distillation Analyzers are designed to meet all industry standard methods, Bruker's analyzer software includes both ASTM D86 and ASTM D1160 correlations. Bruker's highly automated GC, CompassCDS Chromatography Data Handling Software, and integrated SimDist software are also designed to meet worldwide industry standard test methods.



436-GC with Sampler

Key Benefits Include:

- Accurate boiling point distribution up to 750°C
- Integrated standard test methods, applications fully comply with ASTM, IP, DIN and ISO standard test methods
- Complete, single vendor solution
- Complete control from initial setup to final report
- ASTM D86 and ASTM D1160 correlation

Built-in Reports:

Bruker's SimDist software provides a wide variety of report options to meet specific requirements including;

- Chromatogram with merged corrected blank analysis and IBP/FBP marks versus retention time
- Boiling point versus percentage of sample
- Table and plot with retention time versus boiling point
- D86 and D1160 correlations
- DIN Noak and motor oil volatility reports
- Table with cut points and fractions plus residue analysis with recovery calculation up to C120

Hydrocarbon Analysis by Group (PIONA+™)




Figure 1: Traps are easily accessible and do not require any tools to install or replace.

Characterization of Engine Fuels by Hydrocarbon Group Type

Bruker's PIONA+™ Analyzer is a highly flexible GC analysis platform to obtain comprehensive characterization and quantitative information, including hydrocarbon group types, oxygenates and carbon number distribution for spark ignition engine fuels.

Key Analyzer Capabilities:

- Unparalleled operational flexibility
- Compliant with established standards
- A complete and fully integrated solution
- A powerful analyzer, easy to use, generating outstanding analysis results day after day

The PIONA+ Analyzer performs a complete analysis (as described in ASTM D6839 and similar methods) and provides unprecedented analytical flexibility and simplified operation through the use of a novel approach to column/trap heating and exchange (Figure 1).




Figure 2: Chromatogram of a test mix in conventional PIONA mode. (Analysis Time - 180mins)

Reduced Analysis Time and Increased Sample Throughput Efficiency through the Use of “Concurrent” Heating

A unique aspect of the design of the Bruker PIONA+ system is the ability to independently heat the individual traps.

This has 2 major operational benefits:

- Enhanced elution integrity for wide range sample
- Significantly Reduced Analysis Time

The independent and concurrent heating design permits greater trap control and benefits in improved elution integrity of the component groups e.g naphthene, iso-paraffins and n-paraffins even for a wide range sample (C4 - C11). In addition, only a single Molsieve column temperature cycle is employed thus reducing the analysis time by almost half allowing a “fast” PIONA mode of operation (see figures 2 and 3).

By employing this technique, sample throughput can be nearly doubled compared to systems that do not offer this unique capability.




Figure 3: Chromatogram of a test mix using concurrent heating in “fast” PIONA mode. (Analysis Time = 95mins)

Determining Total Olefin Content Is Now Practical

The stability, sample loading and lifetime for all of the critical chromatographic components have been improved and optimized in the Bruker PIONA+ Analyzer. Of special and particular note is the improvement in increasing the sample loading capacity of the "olefin" trap. As a result, it is now possible to analyze streams with olefin content as high as 35-40% or more. This makes it practical to employ a single analytical method to obtain total olefin content (Table 2).

Table 2

Saturates				Unsaturates			Aromatics	Oxygenates	Total
Carbon	Cyclic	Iso	Normal	Cyclic	Iso	Normal			
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.06	0.00	0.03	0.54	0.00	0.00	0.63
5	0.31	11.37	2.98	0.87	9.92	7.58	0.00	0.00	33.03
6	3.19	9.98	1.40	2.40	8.59	4.40	1.72	0.00	31.68
7	4.31	6.77	0.00	2.14	4.76	1.91	7.47	0.00	27.36
8	1.42	3.12	0.00	0.39	2.06	0.00	0.07	0.00	7.06
9	0.01	0.00	0.00	0.04	0.02	0.00	0.03	0.00	0.10
10	0.01	0.00	0.00	0.09	0.00	0.03	0.01	0.00	0.14
11	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Total	9.25	31.24	4.44	5.93	25.38	14.46	9.30	0.00	100.00

Example weight % report for a Naphtha sample with high (46%) olefinic content; (highlighted in blue).

The analyzer design allows the operator to conduct analyses in any one of a number of different operational modes including PNA, PONA, PIONA, O-PONA and O-PIONA in standard and concurrent heating configuration. The system is compliant with established standard methods (see adjacent array).

	PNA	PONA	PIONA	O-PONA	O-PIONA
EN 14517					✓
EN-ISO-22854					✓
ASTM D6839				✓	
DIN 51448-2			✓		
ASTM D1319 (FIA)		✓			
DIN 51448-1	✓				
ASTM D5443	✓				
UOP 870	✓				
IP 382	✓				



Compliant with the method



More information generated than required for the method

Detailed Hydrocarbon Analyzer

The DHA Analyzer is a complete high resolution GC solution for the analysis of hydrocarbons in petroleum streams. It is capable of performing all of the standard methods including ASTM D6729, D6730, D6733, D5134, D6623, IP 344/ DHA "Front End" and "Fast DHA".

Although each DHA analyzer is configured, tested and certified at the factory for a standard method specified by the customer, the DHA software permits the operator to utilize any of the other popular standard methods as well. And, because of the outstanding performance and flexibility of the Bruker GC and CompassCDS software design, Bruker is able to quickly modify the existing methods or add new ones if required as a result of the on-going dynamic industry standard processes.




Figure 4: The analysis of permanent gases and hydrogen using the Rapid RGA.




Figure 5: Report selection output.

Brüker Refinery Gas Analyzers

Peak Identification

1. Hydrogen
2. Carbon Dioxide
3. Hydrogen Sulfide
4. Oxygen
5. Nitrogen
6. Carbon Monoxide
7. Methane
8. Ethane
9. Ethylene
10. Propane
11. Cyclo Propane
12. Propylene
13. i-Butane
14. n-Butane
15. Propadiene
16. Acetylene
17. t-2-Butene
18. i-Butene
19. c-2-Butene
20. i-Pentane
21. n-Pentane
22. 1, 3-Butadiene
23. Propyne
24. Butyne
25. C6+
26. Helium




Figure 6: The separation of light hydrocarbons using the Standard RGA.

The source and composition of refinery gases varies considerably. Measuring gas composition precisely and accurately is a significant challenge in today's refinery operations. Brüker's Refinery Gas Analyzers (RGA) are designed to deliver superior, reliable results for a wide range of sources and analysis throughput requirements.

Key Benefits:

- Pre-configured and tuned
- Standard methods including UOP 539, DIN-51666 and ASTM D2163
- Integrated micro-gasifier ensures complete vaporization of LPGs and high pressure samples to prevent sample discrimination (option)
- Multi-channel approach




Figure 7: The analysis of light hydrocarbons using the Rapid RGA, with complete separation in less than five minutes.




Figure 8: The analysis of permanent gases and hydrogen using the Rapid RGA.

Brüker Offers Two RGA Systems to Meet the Widest Range of Analysis Requirements:

Standard RGA:

A three channel multi-valve design using both capillary and packed columns.

Channel 1 - Analysis of permanent gases

Channel 2 - Light hydrocarbons

Channel 3 - Hydrogen.

Total analysis time for all components in 25 minutes.

Rapid RGA:

The Standard RGA packed columns in the hydrogen and permanent gas channels are replaced by micro packed columns and installed in a separate column oven. Key benefits of this design are:

- Flexibility
- Reduced Analysis Time - 5mins (with H₂S - 7mins)
- Increased Sample Throughput

Characteristics	Standard RGA	Rapid RGA
No. of Channels/Detectors Used	3	3
No. of Column Ovens	1	2
Analysis Time	25 min	5 min (7 min with H ₂ S)
Repeatability	<1%	<1%
Linear Bench Space Required	66 cm/26 in.	66 cm/26 in.
Minimum Component Detection Level	0.01% all components except H ₂ S = 0.05%	0.01% all components except H ₂ S = 0.05%
Suitability		
Typical Refinery Gas	Excellent	Excellent
Impurities in Bulk Ethylene	Excellent	Excellent
Impurities in Bulk Propylene	Excellent	Good
Impurities in Bulk C4	Good	Good

Table 3: Standard RGA vs Rapid RGA.

Low Level Oxygenates Analyzer




Figure 9: Typical chromatogram showing a wide range analysis of a liquid sample stream.

The determination of sub to high ppm levels of ethers, alcohols, aldehydes and ketones in different hydrocarbon matrices is a recurring challenge in the petroleum refining and petrochemical industry. The Bruker Low Level Oxygenates Analyzer is an easy to use solution to meet this challenge and is according ASTM D7423.

The Low Level Oxygenates Analyzer is designed and optimized to quantify ppm and sub ppm levels of ethers (e.g. DME, MTBE, ETBE, DIPE), alcohols (e.g. methanol, ethanol, propanol), ketones (e.g. acetone, MEK) and aldehydes in various hydrocarbon matrices. In general, all oxygenated components with a boiling point of up to 100°C can be analyzed and the sample can be a gas, LPG or liquid under ambient conditions with a final boiling point up to 250°C.

The system is comprised of a Bruker GC configured with gas and liquid sampling valves, two high performance capillary analysis columns, digitally controlled pneumatics including a 'fluidic' switch and Flame Ionization Detector (FID). An optional 'pressure station' can be added to eliminate the possibility of losing sample due to evaporation when analyzing LPG. The GC is controlled via the CompassCDS Chromatography Data Handling Software, which acquires data, processes it and generates analyses reports.

Peak Identification

1. Diethylether
2. Acetaldehyde
3. Ethyl tert. Butyl ether
4. Methyl tert. Butyl ether
5. Diisopropylether
6. Propanal
7. tert amyl methyl ether
8. Propylether
9. Isobutyraldehyde
10. Butyraldehyde
11. Methanol
12. Acetone
13. Isovaleraldehyde
14. Valeraldehyde
15. 2-Butanone
16. Ethanol
17. 1-Propanol
18. tert Butyl alcohol
& Isobutanol
19. 1-Butanol

Bruker 4815 GC Oxygenates Analyzer

The Bruker 4815 GC Analyzer provides a highly cost effective solution for the analysis of oxygenates in gasoline, according to the widely used industry standard method ASTM D4815. The combination of Bruker's reliable GC hardware, powerful software and industry leading pre- and post-sales support teams make this analyzer package the most comprehensive solution available today.

Oxygenated compounds can be present in various hydrocarbon matrices either because they were purposely added (e.g. into gasoline), because they are naturally present, or formed during catalytic processes such as polymer production. In gasoline, oxygenated compounds are added as 'anti-knock' agents to increase the octane number and decrease emissions by replacing organo-lead compounds.

The type and concentration of oxygenated compounds must be measured in reformulated gasolines as part of ongoing product quality assessment, and to confirm the oxygenated components have been added in the correct amounts according to regulatory requirements (e.g. California Air Resources Board).

ASTM D4815 is frequently chosen as the standard method for the determination of oxygenated compounds. Individual ethers and alcohols are quantified in gasoline including: MTBE, ETBE, TAME, DIPE, C1-C4 alcohols and tert-amylalcohol. Individual ether components are measured from 0.1 to 20.0 mass %. The individual alcohols are measured from 0.1 to 12.0 mass %.



456-GC with Sampler




Figure 10 : Typical Chromatogram of the test sample.

Trace Impurity Analyzers

Sulfur Components in LPGs

Low level analysis of sulfur containing components such as H₂S, COS and mercaptanes is extremely challenging and a configured GC offers the solution.

Firstly, the system employs a micro-gasifier enabling the direct coupling of an LPG stream. Secondly, an inert steel sample path ensures a trouble free analysis of sulfur containing components at low concentrations. Finally, a two channel PFPD/ two column approach permits the analysis of all components of interest in one run whatever the LPG matrix. Two differing columns ensures quenching of PFPD signal by the matrix is overcome and full sulfur component analysis is achieved. Figures 11 and 12 show chromatograms obtained in a propane matrix and illustrates the novel benefits of the 2 channel approach.

Permanent Gases in LPGs

Impurities such as CO, CO₂, H₂, O₂ and N₂ need to be determined at low levels in LPGs. Complete separation of these components is done using a two channel single detector (PDHID) system. The GC employs a permanent gas channel for analyzing H₂, O₂, N₂, CH₄ and CO and a specific second channel for the analysis of CO₂. A gasifier is used as a sample introduction device thus giving the capability of handling LPG samples C2 through C4.

Detection limits are at the ppb level (Figure 13), depending on the component of interest.




Figure 11: Sulfur components in propane, BR-1 column.




Figure 12: Sulfur components in propane, BR-Q PLOT.




Figure 13: LPG sample.

Total Characterization of Ethylene Impurities

For a total characterization of impurities in ethylene and also propylene six GC channels are required. By coupling two Bruker GCs with three channels each, a comprehensive solution is available for analyzing these components. The channels used in this analyzer are analytical tools principally developed for the determination of different gases in various hydrocarbon types of gaseous matrices.

GC-1

- H₂ Channel (TCD)
- O₂/N₂ Channel (TCD)
- CO, CO₂ Channel (Methanizer/FID)

GC-2

- Light Hydrocarbon Channel (FID)
- Oxygenates Channel (FID)
- Sulfur Channel (PFPD)

The results (see figures 14 to 19) demonstrate how this 6 channel system is perfectly suited for the total characterization of ethylene and its impurities.




Figure 14: CO, CH₄ and CO₂ on GC-1.




Figure 15: H₂ on GC-1.




Figure 16: O₂ and N₂ on GC-1.




Figure 17: Light hydrocarbons on GC-2.




Figure 18: Oxygenates on GC-2.





Figure 19: Sulfur components on GC-2.

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications (like ASTM, UOP, ISO, GPA and EN) and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- "Quick connect" fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time



For research use only. Not for use in diagnostic procedures.



www.ScionHasArrived.com



www.GlobalEnergyTesting.com

The New Precision Gas Generators for GC Instruments

First we made them multi task...then we conquered space.

At Peak Scientific we understand the needs of laboratories and the people who work in them. The lab is a unique environment demanding precision, accuracy, reliability and design steeped in functionality.


PEAK SCIENTIFIC. GO WITH THE FLOW.



Peak expertise defined by our Precision

Leading the way in gas generation for over a decade, we expertly designed a new system solution that can be tailored to any GC laboratory's needs.


Our global presence allows us to support our customers quickly and effectively. No matter where you are, you will benefit from our outstandingly fast response time and a track record of 95% first time fixes!




Zero Air



Nitrogen




Nitrogen Trace




Hydrogen

Why Precision?

- Modular system offering a GC gas supply solution specific to your lab
- Space saving stackable design allowing you to make the most of your lab space
- Combinations available for single GC and multiple GC applications
- Eliminates the inconvenience of gas cylinders, no more changing over, no more supply issues, no safety worries
- Very low maintenance throughout the range
- Long term cost stability
- Precision systems only utilize technology with a proven track record of safety and reliability
- Stylish lighting feature to indicate generator status




Bespoke Combinations



Zero Air + Hydrogen



Zero Air + Nitrogen + Hydrogen Trace



Zero Air + Nitrogen Trace

Which Precision System is right for you?

More often than not, which generator model we recommend comes down to the 'limit of detection' you are looking for in your GC results. The lower the limit of detection, the purer the gases need to be as you will be looking for the lowest possible baseline to achieve the most precise results.

This is why Peak Scientific offers a Standard Analysis Gas Generator as well as a Trace Analysis Gas Generator solution.

Gas Generator recommendation by Detection Limit:

Detection Limit	Gas Generator Purity Recommendation
0 - 1 ppm	Trace Analysis Purity Solution
1 - 1000 ppm	Trace Analysis Purity Solution
1000ppm - 1%	Standard Analysis Solution
1% - 100%	Standard Analysis Solution

Hydrogen - Trace Analysis Solution


Hydrogen Purity	99.9999%
Moisture Content	<1ppm
Flow Rate Options	500cc/min
Delivery Pressure	0- 100psi/ 0- 6.9 bar
Water Purity Requirements	< 1.0 μ S-cm Conductivity / > 1.0M Ω -cm Resistivity
Electrical Requirements	110-230V, 360VA
Dimensions	H = 406mm W = 380mm D = 539.5mm

- Suitable for Carrier Gas and Flame Gas at trace detection limits
- Proven PEM Technology to generate Hydrogen safely and reliably
- Regenerative PSA Dryers to ensure highest level of purity
- Automatic loading pump as standard
- Maintenance limited to replacing de-ionizer cartridge
- Short and easy start-up and shutdown procedures
- Small and stackable
- Creates Hydrogen on demand, minimal storage of Hydrogen in the system
- Internal leak detection with automatic shutdown features
- Series option to combine multiple units for higher flow requirements
- Remote shutdown
- GC In Oven Hydrogen Leak Detector available as an optional extra

Hydrogen – Standard Analysis Solution

Hydrogen Purity	99.9995%
Flow Rate Options	100, 200, 300, 450cc/min
Delivery Pressure	0- 100psi/ 0- 6.9 bar
Water Purity Requirements	< 1.0 $\mu\text{S-cm}$ Conductivity / > 1.0M $\Omega\text{-cm}$ Resistivity
Electrical Requirements	110-230V, 360VA
Dimensions	H = 406mm W = 380mm D = 539.5mm

- Suitable for Carrier Gas and Flame Gas at standard detection limits
- Proven PEM Technology to generate Hydrogen safely and reliably
- Desiccant Dryers to ensure high level of purity
- Automatic loading pump as standard
- Maintenance limited to replacing de-ionizer cartridge and silica gel
- Short and easy start-up and shutdown procedures
- Small and stackable
- Creates Hydrogen on demand, minimal storage of Hydrogen in the system
- Internal leak detection with automatic shutdown features



Nitrogen - Trace Analysis Solution

Nitrogen Purity	99.9995%
Hydrocarbon Content (as methane)	<0.05ppm
Flow Rate Options	250, 600cc/min
Delivery Pressure	0-80psi
Inlet Air Requirements	Minimum 35 l/min at 120-145psi
	Or Peak Precision Compressor
Electrical Requirements	110/230V, 504VA
Dimensions	H = 406mm W = 380mm D = 539.5mm

- Suitable for Carrier Gas and Make Up Gas at trace detection limits
- Generates Zero Nitrogen on demand from compressed air
- Regenerative CMS columns remove Oxygen and moisture
- Catalyst chamber to remove Hydrocarbons (as methane) to <0.05ppm
- Ultra fast start-up time
- Minimum maintenance with an annual filter change
- Small and stackable

Nitrogen - Standard Analysis Solution

Nitrogen Purity	99.9995%
Flow Rate Options	250, 600, 1,000cc/min
Delivery Pressure	0-80psi
Inlet Air Requirements	Minimum 35 l/min at 100-120psi
	Or Peak Precision Compressor
Electrical Requirements	110-230V, 41VA
Dimensions	H = 406mm W = 380mm D = 539.5mm

- Suitable for Carrier Gas and Make Up Gas at standard detection limits
- Generates Nitrogen on demand from compressed air
- Ultra fast start-up time
- Minimum maintenance with an annual filter change
- Small and stackable

Zero Air - Trace and Standard Analysis Solution

Hydrocarbon Content (as methane)	<0.05ppm
Particles	<0.01 μ m
Flow Rate Options	1,500, 3,500cc/min
Delivery Pressure	0-80psi
Inlet Air Requirements	Minimum 1.5 or 3.5 l/min at 90-145psi
	Or Peak Precision Compressor
Electrical Requirements	110/230V, 144VA
Dimensions	H = 156mm W = 380mm D = 539.5mm

- Generates Zero Air on demand from compressed air
- Catalyst chamber to remove Hydrocarbons (as methane) to <0.05ppm
- Minimum maintenance with an annual filter change
- Small and stackable

Compressed Air

Electrical Requirements	110/230, 564VA
Dimensions	H = 406mm W = 380mm D = 539.5mm

- Suitable for a variety of Precision Generator combinations
- Minimal noise emission due to insulated compressor compartment
- Minimal vibration through especially developed compressor anti-vibration mounts
- Compressor service indication
- Serviceable compressor

Model Description	110v	230v	Annual Service Kit
Precision Nitrogen Trace, 250cc	62-0251	62-0250	08-3613
Precision Nitrogen Trace, 600cc	62-0601	62-0600	08-3613
Precision Zero Air, 1.5L	60-1501	60-1500	08-3611
Precision Zero Air, 3.5L	60-3501	60-3500	08-3611
Precision Air Compressor	65-1555	65-0555	08-8343

Model Description	110/230v	Annual Service Kit
Precision Nitrogen, 250cc	61-0250	08-3612
Precision Nitrogen, 600cc	61-0600	08-3612
Precision Nitrogen, 1L	61-1000	08-3612
Precision Hydrogen, 100cc	63-0100	08-3609
Precision Hydrogen, 200cc	63-0200	08-3609
Precision Hydrogen, 300cc	63-0300	08-3609
Precision Hydrogen, 450cc	63-0450	08-3609
Precision Hydrogen Trace, 500cc	64-0500	08-3610

Standard Maintenance Plan 09-3110
Complete Maintenance Plan 09-3010

Website: www.peakscientific.com
Email: precision@peakscientific.com

Product Accreditation:



Our Locations

Peak Scientific UK

Fountain Crescent
Inchinnan Business Park
Inchinnan, PA4 9RE
Scotland, UK

Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8300

Peak Scientific Inc

19 Sterling Road
Suite #1, Billerica
MA 01862, USA

Tel: +1-866-647-1649
Fax: +1-978-608-9503

For a full list of our worldwide
office locations, please visit:
www.peakscientific.com/contact





SCION™ GC Series

- The Gas Chromatographers Choice for Separations

Innovation with Integrity

GAS CHROMATOGRAPHY

Innovation in Gas Chromatography

Bruker's long tradition of innovation and product reliability have combined to create the next generation of Bruker Gas Chromatographs. By understanding and then designing to exceed the most critical performance and reliability needs of GC users, Bruker is delivering systems that are especially for, and all about, the ultimate success of the GC user. The new SCION 436-GC and SCION 456-GC have been designed to meet the most important user specified requirements – reliable performance, ease of use and simple maintenance.


Local User Interface

This large, high resolution display makes all GC functions accessible via touch screen control and “instant access” buttons. Easy to navigate and adapt it comes available in 13 languages for ease of local training and support.

Fast, Flexible Detection

Bruker's comprehensive range of detectors deliver industry leading sensitivity, ease of operation and outstanding reliability. And, now all Bruker detectors feature fast sampling data rate (600Hz) for rapid separations and greater analysis throughput.

Multi-language capability



SCION 436-GC

Choice of multi inlet systems

- 2 injectors from 5 available
- Gas Sampling Valve
- Liquid Sampling Valve

Choice of 7 GC traditional detectors

- Universal
- Specific

Choice of MS detector

- Single Quadrupole (MS)
- Triple Quadrupole (MS/MS)



SCION 456-GC

Choice of multi inlet systems

- 3 injectors from 5 available
- Gas Sampling Valve
- Liquid Sampling Valve

Choice of 7 GC traditional detectors

- Universal
- Specific

Choice of MS detector

- Single Quadrupole (MS)
- Triple Quadrupole (MS/MS)



Enhanced Operator Benefits

Bruker offers a range SCION GCs to meet virtually all application requirements. All SCION GCs are equipped with the convenience of advanced EFC. Whatever the requirement, we have the solution.

GC Control From Anywhere

The unique embedded control architecture incorporated into the GC enables the use of remote user interface software. This offers the user the ability to control the GC in the exact same way and with the same level of functionality as if they were standing at the GC using the User Interface but from a remote location, even from home.


Gas Saving Capability

Essential to any laboratory is the requirement for cost control. Gas consumption is a major source of operational cost and reduction in gas flow when not necessary is a vital function of any modern GC. Bruker GCs have the necessary functionality to save valuable gas and thus costs.

Turnkey Analyzer Solutions

Bruker configures and tests GC hardware and software according to widely used industry standard methods (e.g. ASTM, UOP, EN, ISO, GPA), to save its clients time and to ensure confidence in results. Standard analyzers are configured to meet the performance specifications outlined in the method itself. Included in these analyzer packages:

- All hardware
- Software (including special application "plug-ins" when appropriate)
- Pre-installed methods
- Test chromatograms
- Installation/validation data
- User documentation customized to the specific method



Targeted Solutions for Specific Markets

A series of software customization tools allow users to develop unique calculation modules, that fully integrate with CompassCDS. A large number of standard plug-ins are available that allow special reporting and other post analysis functions. Some examples include:

Simulated Distillation



Provides automated boiling point distributions for a full range of petroleum products for applications that comply with ASTM, IP, DIN and ISO standard test methods.

Detailed Hydrocarbon Analysis

Reports in an automated way the physical properties of gasoline and similar products based on individual components for applications that comply with Bruker developed methods and ASTM, IP and standard test methods.

CompassCDS PeakSynch

Provides rapid visual and quantitative comparison of complex chromatograms and is widely used in the flavor and fragrance industry.





OPERATION FIDsens® Run: 1.95 End: 9.50

End 9.50
1.95

Oven 68.5°C

Status Run

Vial

Run 0 / 0

Force Standby

Signals	
Front	FID 0 µV
Mid	
Rear	

Method Log

BRUKER

SCION T7
455-EC

Clipboard with document and blue pen



Benefits of the SCION-GC

- Multi-language User Interface
- Full EFC Capability
- High Pressure Injection
- IntelliUpdate
- System Suitability Determination
- CompassCDS Software
- 600 Hz Data Sampling Rate On All Detectors
- Inert GC Sample Path
- Constant Linear Velocity Mode
- Fast Cycle Time

SCION 436-GC



- Small foot print
- High performance
- Dual channel architecture

SCION 456-GC



- Solutions platform
- Total flexibility
- Four channel architecture

Increased Productivity


With over 40 years' experience in gas chromatography, we can provide unrivalled expertise, not only in building robust instruments, but also in creating solutions for ensuring productivity. With total control over design and manufacturing, Bruker ensures the quality and technological excellence of its products is complimented and combined with features that deliver the true benefits of productivity.

Speed increased with a factor 6.5

- Small ID from 0.25 to 0.1 mm
- Short column from 15 to 4 mtr
- Increased ramp from 10 to 65 °C/min
- Data rate from 25 to 200 Hz


Fast Cycle Time

The time between injections can considerably improve productivity. The high performance oven incorporating design characteristics that enhances fast heating and cooling ensures maximum productivity. This, in conjunction with high pressure injectors and ultra narrow bore columns will significantly improve and yield fast cycle times without loss of performance (see chromatograms).



IntelliUpdate

In many cases instrument and system effects (column ageing, matrix, etc.) can cause experimental deviations e.g. retention time. CompassCDS IntelliUpdate function can be used automatically to correct and compensate for such deviations. This unique capability is also done without changing fundamental instrument parameters, maintaining accuracy of results and method continuity.







Optimized Switching Valves

Flow splitting, backflushing and Deans switching are valuable techniques in improving cycle times, analytical performance and the robustness of GC methods. Splitting the flow of column effluent into differing detectors can enhance performance, quantitation and confirmation of targeted compounds.

Backflushing is key to reducing analysis time and column protection. It works on the basis of reversing column flow after peaks of interest have been detected. This eliminates the need for time and temperature segments to elute highly retained components injected with compounds of interest. Reversing the flow elutes these materials out through the split vent of the injector with the added benefit of protecting the column from degradation and contamination.

Backflushing capability also allows column changes and injector maintenance without loss of vacuum in the MS detector.




Deans switching enables the use of multiple traps and columns of differing phases in one method/analysis. It is the basis of two dimensional GC and many standard analyzers manufactured by Bruker. Use our expertise to configure the optimum system for you.

Capability and Automation

Bruker offers an injector and detector range to meet virtually all application and market requirements. All are equipped with the convenience of advanced EFC. Whatever the requirement, be it Split/splitless, Cold-on-Column, Packed, Flash or Programmable Temperature Vapourizing injector with a universal or specific detector we have the solution.

	Universal				Specific		
	FID	TCD	PHHID	MS	ECD	NPD (TSD)	PFPD
 Academic	✓	✓		✓		✓	
 Environment			✓	✓	✓	✓	✓
 Food Beverage	✓		✓	✓	✓	✓	✓
 Forensics Toxicology	✓			✓			✓
 Petroleum	✓	✓	✓	✓	✓	✓	✓

Bruker offers a range of differing injector designs for all applications, column dimensions and can be fully automated.



Injector Selection Guide - Sample/Analysis Characteristics or Requirements

Trace Analysis	Separation & Speed	Sample Capacity	Wide Range of Analytes	Preferred Column Type	1st Choice	2nd Choice
✓				Capillary, 0.53 mm ID	Large Volume (LV)	Split/Splitless
	✓			Capillary, 0.1 to 0.53 mm ID	Split/Splitless	Large Volume (LV) SS Mode
	✓	✓		Capillary, 0.53 mm ID	Large Volume (LV)	
	✓			Capillary, retention gap	Cold On-Column	Large Volume PTV Mode
	✓			Capillary, 0.53 mm ID	Packed	Large Volume (LV) PTV Mode
			✓	Capillary, 0.53 mm ID	Cold On-Column	Large Volume (LV) On-Column Mode

Regardless of your sample throughput requirements, Bruker can provide an automated solution to meet your needs. Four samplers are available, the CP-8410, CP-8400, the SHS-40 and the PAL Combi-xt. Each is tailored to meet a differing need and workload.



CP-8400

- High throughput
- 100 x 2ml sample capacity
- Dual/Duplicate Injection
- SPME



CP-8410

- Flexibility
- Accommodates 2, 5, 10 ml vials
- Low cost/high performance
- Ease of use



SHS-40

- Fully automated
- Easily Integrated
- Low maintenance
- Sample Loop or Press and Inject configuration



PAL Combi-xt


- High throughput
- Liquid handling capability
- SPME
- ITEX

Bruker-Certified Consumables for Your SCION GC Series

Bruker GC columns span a broad range of column diameters, stationary phases, and capillary column materials: Fused Silica (FS) and Inert Steel (IS). Ideal for either routine or research type analyses.

Bruker GC column offerings bridge across many important applications and include a number of offerings such as:


- Standard WCOT (Wall Coated Open Tubular)
- Solid Stationary Phase PLOT (Porous Layer Open Tubular)
- Inert Steel Micro-Packed and Packed



Super Clean™ Gas Filters

Bruker Gas Purification Systems have the range to satisfy your needs from individual to combination filters, from Ultra purity combined with Ultra capacity, to all in one solution kits. Innovative features designed into the product yield extensive benefits to the user.

- Ultra-high capacity for long life, less change and improved productivity
- High-purity output ensures 99.9999% Pure Gas
- “Quick connect” fittings for easy, leak-tight filter changes
- Glass internals prevent diffusion; plastic externally for safety
- Easy-to-read indicators for planned maintenance and improved up-time



For research use only. Not for use in diagnostic procedures.





● SCION™ TQ Premium GC Triple Quadrupole Mass Spectrometer

Specification Sheet

The SCION TQ Premium is the chromatographer's choice for triple quadrupole mass detector; it is designed to match your most stringent needs for analytical performance and productivity. The SCION TQ Premium offers superior sensitivity and robustness based on the innovative ion optics, and fast and easy methods development for multi-component quantitation following the unique Compound Based Scanning (CBS) approach with MRM library. The SCION TQ Premium GC-MS/MS system defines a new standard of usability for routine analysis and has the smallest bench footprint in the industry.

Analyzer - MS Specifications

- Scan modes: Full Scan, Precursor, Product, Neutral Loss/Gain Monitoring, Selected Ion Monitoring (SIM), Multiple Reaction Monitoring (MRM) and Result Dependent Scanning
- Standard ionization mode: Electron Ionization (EI)
- Ion source: Auto-aligning EI source constructed of inert materials
- q0 ion guide: 90° curved RF-only entrance quadrupole with active ion beam focusing and heating at 135 °C
- Source temperature: 100 °C to 350 °C
- Filament and emission current: dual filaments; up to 200 µA
- Electron energy: adjustable from 0 to 150 eV
- Mass filters: quadrupole with pre- and post-filters; high ion transmission efficiency lens-less design
- Collision cell: 180° curved path with pre- and post-filter regions
- Collision cell gas: Argon
- Collision energy: selectable up to 75 eV
- Mass range (m/z): 1 to 1200 Da
- Scan rate: up to 14,000 Da/sec
- Minimum scan time (dwell time): 1 ms
- Maximum acquisition MRM rate: 500 MRM's/sec
- Resolution: user-adjustable from 0.7 Da (Unit) to 4 Da, also with three user-selectable settings (Unit, Standard, Open) on both Q1 and Q3.
- Mass axis stability: ± 0.1 Da over 24 hours
- Transfer line temperature: up to 350 °C
- Manifold temperature: 40-50 °C
- Detector: EDR™ Electron multiplier with ± 5 kV post acceleration and with on-the-fly multiplier gain optimization for Extended Dynamic Range (EDR); direct ion collection onto multiplier for negative ion detection without dynode loss
- Turbomolecular pump: dual stage, 310/400 L/sec, air-cooled for helium carrier gas flow up to 25 mL/min.
- Foreline pump: dual-stage rotary vane; voltage 120/230V
- Power requirements: 100-240 Vac, 50/60 Hz ± 3 Hz, 1200 VA
- Operating environment temperature: 15 °C to 33 °C
- Operating environment humidity: 20% to 80% relative humidity (without condensation)



Software

- Bruker MS Workstation equipped with the Compound Base Scanning (CBS) MRM library for data acquisition, data handling, and reporting
- Optional spectral libraries: NIST, Wiley, and Maurer/Pfleger/Weber (PMW) libraries and with user-customizable libraries and automatic searching of multiple libraries

Gas Chromatograph (Bruker 436 and 456 Model GC)

For more specification on GC, refer to the GC Data Sheet

- Injectors: Split/Splitless (SSL), Programmable Temperature Vaporization (PTV), etc. Back-flush option available for all injectors.
- Autosamplers: CP 8400; CP 8410; CTC PAL COMBI-xt
- GC Oven Temperature: Ambient+4 °C to 450 °C
- Temperature Ramps/Holds: 24/25
- Pneumatic: Electronic Flow Control (EFC) or Manual (Model 456)
- ChromatoProbe™: Direct introduction of solids, liquids or slurries (requires PTV injector)
- MS Tuning, tune-to-target, pump-down, and venting controlled by multi-language touchpad on the GC.

Performance Specifications*

Mode	Test (with SSL injector in hot splitless mode)	Specification†
El Full Scan	1 pg Octafluoronaphthalene (OFN) from m/z 50 to 300 for m/z 272	S/N ≥600:1
El SIM	25 fg OFN for m/z 272	S/N ≥50:1
El MRM	100 fg OFN for m/z 272>222	S/N ≥2000:1
El MRM Precision (IDL**)	8 replicate injections of 50 fg OFN in El MRM mode (m/z 272>222)	Peak Area RSD ≤ 6.7% (10 fg)
PCI Full Scan‡	10 pg Benzophenone (BZP) from m/z 80 to 230 for m/z 183	S/N ≥50:1
PCI SIM‡	1 pg BZP for m/z 183.105	S/N≥50:1
PCI MRM‡	100 fg BZP for m/z 183>105	S/N≥100:1
NCI Full Scan‡	1 pg OFN from m/z 200 to 300 for m/z 272	S/N ≥4000:1
NCI SIM‡	10 fg OFN for m/z 272	S/N ≥300:1

* All tests use helium as carrier gas. El MRM sensitivity test will be used as installation checkout specifications; not all other performance tests are confirmed at installation.

† The Signal-to-Noise ratio S/N values are based on RMS

IDL**: Instrument Detection Limit, defined as $IDL=t(0.99, f=7) \times S$, whereas $t(0.99, f=7)$ is the one-sided student's t-distribution value of 2.998 for 99% of confidence and for degree of freedom 7 ($f=n-1$, n the number of injections); S is the peak area standard deviation of 8 replicate injections.

‡ CI tests use methane as reagent gas

Dimensions (H x W x D) and Weight

- SCION TQ: 45 cm (18 in.) x 28 cm (11 in.) x 57 cm (22.5 in.), 40 kg/88 lb
- 436 GC: 57 cm (23 in.) x 32 cm (13 in.) x 61 cm (24 in.); 27 kg/59 lb
- 456 GC: 57 cm (23 in.) x 66 cm (26 in.) x 56 cm (22 in.); 43 kg/95 lb
- CP-8400/8410 Autosamplers: 40 cm (16 in.) x 22 cm (9 in.) x 47 cm (18 in.); 7 kg/15.3 lb



• SCION™ TQ Select GC Triple Quadrupole Mass Spectrometer

Specification Sheet

The SCION TQ Select is the chromatographer's choice for triple quadrupole mass detector; it is designed to match your most stringent needs for analytical performance and productivity. The SCION TQ Select offers superior sensitivity and robustness based on the innovative ion optics, and fast and easy methods development for multi-component quantitation following the unique Compound Based Scanning (CBS) approach with MRM library. The SCION TQ Select GC-MS/MS system defines a new standard of usability for routine analysis and has the smallest bench footprint in the industry.

Analyzer - MS Specifications

- Scan modes: Full Scan, Precursor, Product, Neutral Loss/Gain Monitoring, Selected Ion Monitoring (SIM), Multiple Reaction Monitoring (MRM) and Result Dependent Scanning
- Standard ionization mode: Electron Ionization (EI)
- Ion source: Auto-aligning EI source constructed of inert materials
- q0 ion guide: 90° curved RF-only entrance quadrupole with active ion beam focusing and heating at 135 °C
- Source temperature: 100 °C to 350 °C
- Filament and emission current: dual filaments; up to 200 µA
- Electron energy: adjustable from 0 to 150 eV
- Mass filters: quadrupole with pre- and post-filters; high ion transmission efficiency lens-less design
- Collision cell: 180° curved path with pre- and post-filter regions
- Collision cell gas: Argon
- Collision energy: selectable up to 75 eV
- Mass range (m/z): 1 to 1200 Da
- Scan rate: up to 14,000 Da/sec
- Minimum scan time (dwell time): 1 ms
- Maximum acquisition MRM rate: 500 MRM's/sec
- Resolution: user-adjustable from 0.7 Da (Unit) to 4 Da, also with three user-selectable settings (Unit, Standard, Open) on both Q1 and Q3.
- Mass axis stability: ± 0.1 Da over 24 hours
- Transfer line temperature: up to 350 °C
- Manifold temperature: 40-50 °C
- Detector: EDR™ Electron multiplier with ± 5 kV post acceleration and with on-the-fly multiplier gain optimization for Extended Dynamic Range (EDR); direct ion collection onto multiplier for negative ion detection without dynode loss
- Turbomolecular pump: dual stage, 310/400 L/sec, air-cooled for helium carrier gas flow up to 25 mL/min.
- Foreline pump: dual-stage rotary vane; voltage 120/230V
- Power requirements: 100-240 Vac, 50/60 Hz ± 3 Hz, 1200 VA
- Operating environment temperature: 15 °C to 33 °C
- Operating environment humidity: 20% to 80% relative humidity (without condensation)



Software

- Bruker MS Workstation equipped with the Compound Base Scanning (CBS) MRM library for data acquisition, data handling, and reporting
- Optional spectral libraries: NIST, Wiley, and Maurer/Pfleger/Weber (PMW) libraries and with user-customizable libraries and automatic searching of multiple libraries

Gas Chromatograph (Bruker 436 and 456 Model GC)

For more specification on GC, refer to the GC Data Sheet

- Injectors: Split/Splitless (SSL), Programmable Temperature Vaporization (PTV), etc. Back-flush option available for all injectors.
- Autosamplers: CP 8400; CP 8410; CTC PAL COMBI-xt
- GC Oven Temperature: Ambient+4 °C to 450 °C
- Temperature Ramps/Holds: 24/25
- Pneumatic: Electronic Flow Control (EFC) or Manual (Model 456)
- ChromatoProbe™: Direct introduction of solids, liquids or slurries (requires PTV injector)
- MS Tuning, tune-to-target, pump-down, and venting controlled by multi-language touchpad on the GC.

Performance Specifications*

Mode	Test (with SSL injector in hot splitless mode)	Specification†
EI Full Scan	1 pg Octafluoronaphthalene (OFN) from m/z 50 to 300 for m/z 272	S/N \geq 600:1
EI SIM	25 fg OFN for m/z 272	S/N \geq 50:1
EI MRM	100 fg OFN for m/z 272>222	S/N \geq 2000:1
EI MRM Precision (IDL**)	8 replicate injections of 50 fg OFN in EI MRM mode (m/z 272>222)	Peak Area RSD \leq 6.7% (10 fg)

* All tests use helium as carrier gas. EI MRM sensitivity test will be used as installation checkout specifications; not all other performance tests are confirmed at installation.

† The Signal-to-Noise ratio S/N values are based on RMS

IDL**: Instrument Detection Limit, defined as $IDL = t(0.99, f=7) \times S$, whereas $t(0.99, f=7)$ is the one-sided student's t-distribution value of 2.998 for 99% of confidence and for degree of freedom 7 ($f=n-1$, n the number of injections); S is the peak area standard deviation of 8 replicate injections.

Dimensions (H x W x D) and Weight

- SCION TQ: 45 cm (18 in.) x 28 cm (11 in.) x 57 cm (22.5 in.); 40 kg/88 lb
- 436 GC: 57 cm (23 in.) x 32 cm (13 in.) x 61 cm (24 in.); 27 kg/59 lb
- 456 GC: 57 cm (23 in.) x 66 cm (26 in.) x 56 cm (22 in.); 43 kg/95 lb
- CP-8400/8410 Autosamplers: 40 cm (16 in.) x 22 cm (9 in.) x 47 cm (18 in.); 7 kg/15.3 lb

maXis impact

Specification Sheet
Part Number: # 282000

maXis impact UHR Time-of-Flight Mass Spectrometer System



Size	Bench-top: (Footprint) 198cm (Height)
Weight	~ 210 kg
Vacuum System	5 stages, 28 m ³ /h rough pump
Apollo II ion funnel electrospray source	Flow rate: 1 µL/min – 1 mL/min
Mass Range	20 – 40,000 m/z
Quadrupole isolation	Up to 3,000 m/z
Quadrupole Mass Range	Up to 40,000 m/z
Mass accuracy in MS and MS/MS	With internal calibrant: better than 1ppm RMS Error With external calibrant: better than 2 ppm RMS Error
Calibration	ONE calibration valid for MS and MS/MS analysis. Calibration is independent from charge state of calibrant mass
Mass resolution	40,000 FSR (full sensitivity resolution)
Isotopic pattern	The true isotopic pattern is maintained due to TIP™ technology (True Isotopic Pattern) and allows three dimensional chemical characterizations of analytes via SmartFormula™3D algorithm using exact mass, TIP, and MS/MS fragment data.
SmartFormula™3D	Enables unambiguous formula determination at "sub-ppm" confidence level up to 1000 Da.

Mass stability & dynamic range	hrXIC (high resolution Extracted Ion Chromatogram) technology with better than +/- 1.0 mDa stability on centroid data values over an typical LC peak.
Full scan sensitivity in MS	Reserpine 1 pg S/N>100:1 RMS
Full scan sensitivity in MS/MS	The signal height obtained from a consumption of 2.5 fmol of Glu-Fibrinopeptide B will be better than 100 counts on the most intense y' sequence ion from the MS/MS spectrum of the doubly charged precursor ion. This shall correspond to a signal to noise ratio better than 50:1. The MS/MS sensitivity specification is met while using quadrupole isolation of the precursor ion demonstrating that there is minimal transmission loss through the isolating quad. ? A solution of 100 fmol/μL Glu-Fibrinopeptide B shall be introduced at a flow rate of 3 μL/min.
TOF repetition rate	Up to 20 kHz
Temperature compensation	Yes
Digitizer	4Gsample/sec ADC with 32 Gbit/sec
Acquisition rate	up to 50 Hz MS 50 Hz MS/MS (profile and peak detected spectra to disk)

Optional accessory

APCI II	Optional accessory
APPI II	Optional accessory
GC-APCI	Allows for direct GC coupling (Optional accessory)
Bruker CaptiveSpray	Optional accessory
On-/Off-Line Nanospray	Optional accessory
CE/MS interface	With grounded needle for easy CE-TOF set-up (Optional accessory)
Bruker EASY-nLC	Split-free nano-flow HPLC system

● micrOTOF-Q III

Specification Sheet
Part Number: # 728889

micrOTOF-Q III ESI Quadrupole Time-of-Flight Mass Spectrometer System

Size	Bench-top: 64 x 95 cm (Footprint) 132 cm (Height)
Weight	160 kg
Vacuum System	5 stages, 28 m ³ /h rough pump
ESI dual ion funnel electrospray source	Flow rate: 1 µL/min – 1 mL/min
TOF Mass Range	20 – 40,000 m/z
Quadrupole isolation	Up to 3,000 m/z
Quadrupole Mass Range	Up to 40,000 m/z
Mass accuracy in MS and MS/MS	With internal calibrant: better than 2 ppm RMS Error With external calibrant: better than 5 ppm RMS Error
Calibration	ONE calibration valid for MS and MS/MS analysis. Calibration is independent from charge state of calibrant mass
Mass resolution (FSR: full sensitivity resolution)	20,000 (FWHM) @ 922 m/z at full sensitivity
Isotopic pattern	The true isotopic pattern is maintained due to TIP™ technology (True Isotopic Pattern) and allows three dimensional chemical characterizations of analytes via SmartFormula™3D algorithm using exact mass, TIP, and MS/MS fragment data.
SmartFormula™3D	Enables unambiguous formula determination at "sub-ppm" confidence level up to 1000 Da.
Mass stability & dynamic range	hrEIC (high resolution Extracted Ion Chromatogram) technology with better than 2 mDa stability on centroid data values over an typical LC peak.

ESI sensitivity in MS positive mode	Reserpine 2 pg S/N>100:1 RMS 200fg S/N > 100:1 RMS with IonBooster (Option)
Full scan sensitivity in MS/MS	The signal height obtained from a consumption of 2.5 fmol of Glu-Fibrinopeptide B will be better than 100 counts on the most intense y' sequence ion from the MS/MS spectrum of the doubly charged precursor ion. This shall correspond to a signal to noise ratio better than 50:1. The MS/MS sensitivity specification is met while using quadrupole isolation of the precursor ion demonstrating that there is minimal transmission loss through the isolating quad. ? A solution of 100 fmol/μL Glu-Fibrinopeptide B shall be introduced at a flow rate of 3 μL/min.
TOF repetition rate	Up to 20 kHz
Temperature compensation	Yes
Digitizer	2 GSsamples/sec sampling rate with 16Gbit/sec
Acquisition rate	up to 40 Hz MS 20 Hz MS/MS (profile and peak detected spectra to disk)

Optional accessory

IonBooster	Optional ion source
APCI II	Optional ion source
APPI II	Optional ion source
GC-APCI	Allows for direct GC coupling (Optional ion source)
APLI	Optional ion source
CryoSpray	Optional ion source
Bruker CaptiveSpray	Optional ion source
On-/Off-Line Nanospray	Optional ion source
CE/MS interface	With grounded needle for easy CE-TOF set-up (Optional)
DIP	Direct Probe (Optional ion source)



Application Note # CA-270689

Typical Detection Limits for the aurora M90 ICP-MS in Normal Sensitivity Mode

Introduction

Bruker's aurora M90 is the world's first ICP-MS system with tunable gigahertz sensitivity (over 10^9 cps per mg/L) without compromising oxide and/or doubly charged interference performance, thanks to the patented 90-degree ion mirror [1, 2]. This tunable sensitivity provides the flexibility to choose a suitable sensitivity mode for different sample types or applications, and hence, to achieve the lowest detection limits for a selected application. In this application note, typical instrument detection limits for most elements commonly monitored in environmental analyses are determined using the aurora M90 ICP-MS under normal sensitivity mode.

Basic principles

Limit of detection

According to the recommendation given by the International Union of Pure and Applied Chemistry (IUPAC) [3], the limit of detection expressed as the concentration c_L or the quantity q_L is derived from the smallest measure X_L that can be detected with reasonable certainty for a given analytical procedure. The value of X_L is given by the equation,

$$(1) X_L = X_{bi} + k s_{bi}$$

where X_{bi} is the mean of the blank measures, s_{bi} is the standard deviation of the blank measures, and k is a numerical factor chosen according to the confidence level desired. IUPAC has also recommended a value of 3 for k , which gives a confidence level of about 98%. In ICP-MS, the X_{bi} is the mean cps of the blank measurements, and a typical ICP-MS calibration plot can be expressed by the following equation,

$$(2) X_A = X_{bi} + m C_A$$

where X_A is the cps measured at an analyte concentration of C_A , and m is the sensitivity (i.e. slope of the calibration plot).

Hence in ICP-MS, the value of X_L can also be given as,

$$(3) X_L = X_{bi} + m C_L$$

From equations (1) and (3), the concentration C_L can be calculated by the following equation,

$$(4) C_L = 3 s_{bi} / m$$

Background Equivalent Concentration

In the Bruker ICP-MS Software, the background equivalent concentration (BEC) is calculated by the following equation,

$$(5) BEC = X_{bi} / m$$

Hence, the concentration C_L can be calculated alternatively from the BEC and the % relative standard deviation (RSD) of the blank, that is,

$$(6) C_L = 3 RSD_{blank} BEC / 100$$

where $RSD_{blank} = 100 s_{bi} / X_{bi}$

Experimental

Instrument

Typical instrument detection limits were measured using two aurora M90 ICP-MS instruments, one located in Australia and the other in Japan. Both instruments were installed in non-clean room environment.

Reagents and Samples

A blank solution, 1% v/v HNO_3 , was made using high-purity nitric acid (Ultrapur[®], 60%, Merck) and pure deionized water (18M Ω cm, Millipore Milli-Q, Bedford, USA). Two standards, at 1 μ g/L and 10 μ g/L levels, were prepared by diluting ICP-MS multi-element stocks with the blank. Prior to use, all labwares were thoroughly cleaned (i.e. acid washing and deionized water rinsing), and then the clean containers were left filled with 1% v/v HNO_3 until use.

Results and discussion

Sensitivity Mode

The aurora M90 ICP-MS can be operated in a number of sensitivity modes (eg. normal, high, etc.) without any hardware changes. In general, the “high” sensitivity mode is used for applications requiring highest sensitivity, such as laser ablation and analysis of semiconductor materials; while the “normal” sensitivity mode is used for general chemical analyses, including environmental, agriculture and clinical research applications. Hence the “normal” mode is recommended for most ICP-MS analyses. Typical method parameters used for tuning the instrument to a “normal” sensitivity mode are shown in Table 1.

Table 1: Typical method parameters for the aurora M90 ICP-MS tuned to “normal” mode.

	Parameters	Settings*
Gas flow (L/min)	Plasma Flow	16.5
	Auxiliary Flow	1.65 (1.60)
	Nebulizer Flow	0.95 (1.00)
	Sheath Flow	0.28 (0.45)
RF	RF Power (kW)	1.30 (0.70)
Sample Introduction	Sampling Depth (mm)	5.0 (5.5)
	Pump Rate (rpm)	3 (20)
	Stabilization Time (s)	30
	Spray chamber (°C)	3
Ion Optics (volts)	1 st Extraction Lens	-1
	2 nd Extraction Lens	-140 (-20)
	3 rd Extraction Lens	-200 (-190)
	Corner Lens	-180 (-110)
	Mirror Lens Left	75 (105)
	Mirror Lens Right	5 (15)
	Mirror Lens Bottom	50 (20)
Quadrupole Scan	Entrance Lens	1
	Entrance Plate	-37 (-50)
	Fringe Bias	-5
	Pole Bias	0
	Scan mode	Peak Hopping
Dwell Time (ms)	100	
Points per Peak	1	
Scans/Replicate	30	
Replicates/Sample	10	

* Settings in parentheses are used for “Cool Plasma” conditions.

Settings used for typical “cool” plasma are also listed in Table 1. The “cool” plasma technique is used to minimize polyatomic interferences associated with the plasma gas, such as ArO⁺ and Ar⁺. This technique can improve the detection limits for elements affected by such interferences, including Fe, Ca, Na, K and Mg. More discussions on the use of cool plasma and more detection limit values under cool plasma conditions can be found in other Bruker’s Application Notes (from Bruker’s web site www.bruker.com under the ICP-MS Application Note section).

Better detection limits

It should be noted that tuning an ICP-MS to its highest sensitivity does not necessarily provide the lowest detection limit. To achieve the lowest possible detection limits and accurate analytical results, strict precautions must be taken to eliminate or minimize any potential contamination. Where possible, glassware including volumetric pipettes and flasks should be avoided when preparing and/or storing any solutions (with exception of Hg solutions), because some metals may be leached out from the glass or adsorbed onto the glass surface, which could result in sample contamination or loss of analyte. Prior to use, all the labwares, new or used, should be thoroughly cleaned. A typical cleaning procedure includes acid washing the labwares for at least 24 hours to remove elemental contamination, thoroughly rinsing them with high-purity deionized water, and then leaving clean containers filled with 1% v/v HNO₃ until use.

It is clear from the equations (4) and (6), the lower the standard deviation (or RSD) of the blank, the better the detection limits. A lower RSD can often be obtained by using a relative longer replicate reading and stabilization time. The replicate reading time is dependent on the dwell time and the number of scans per replicate. Typical stabilization and scan settings used in this work are listed in Table 1. Also, the lower the BEC, the better the detection limits. To keep BEC low, high-purity reagents and deionized water should be used in all the samples and standards preparations, and ideally the solutions should be prepared and measured in a class 100 clean room, or at least the sample preparation area should be air-conditioned and dust free. When running the aurora M90 ICP-MS under normal sensitivity mode, the counts for a blank solution should not exceed a few thousand cps per isotope for most isotopes. A higher blank count is often an indication of blank contamination. The lower the contamination (blank counts) the better the detection limits.

Typical detection limits

Table 2 shows the typical detection limits (DLs) for the elements commonly measured in environmental samples. All DLs were calculated using equation (3), i.e. three times the standard deviation of 10 replicates of the blank (i.e. 1% v/v HNO₃). The instrument was tuned to “normal” sensitivity mode under either “hot” or “cool” plasma conditions. All the measurements, however, were made under routine analytical laboratory, not clean-room, conditions. This work indicates typical DL values that can be routinely achieved outside a clean-room in a clean laboratory.

Conclusions

Detection limits are influenced by a number of factors, including the sensitivity of a given isotope, and the presence of background interferences or contamination. It is vitally important to control contamination in the laboratory to achieve the lowest possible detection limits.

Table 2: Typical IDLs for the aurora M90 ICP-MS.

Element	Isotope (m/Z)	Measured in		Back-ground species
		Hot plasma (ng/L)	Cold plasma (ng/L)	
Li	7	1	0.01	
Be	9	3		
B	11	30		
Na	23	200	0.5	
Mg	24	2	0.2	
	25	5	0.08	
Al	27	2	0.2	
Si	28	1000		Co ⁺ , N ₂ ⁺
P	31	700		NOH ⁺
S	34	20000		(OH) ₂ ⁺
K	39	500	0.5	ArH ⁺
Ca	40	500	1	Ar ⁺ , CO ₂ ⁺
	44			
Sc	45	0.8		CO ₂ H ⁺ , N ₂ OH ⁺
Ti	47	3		
V	51	3		ArNH ⁺ , ClO ⁺
Cr	52	8		ArO ⁺ , ArC ⁺
	53	3		
Mn	55	2		ArNH ⁺
Fe	56	4000	0.3	ArO ⁺
	57	300	0.9	ArOH ⁺

Element	Isotope (m/Z)	Measured in		Back-ground species
		Hot plasma (ng/L)	Cold plasma (ng/L)	
Co	59	0.2		ArOH ⁺
Ni	60	2		
Cu	63	0.3		
	65	2		
Zn	66	5		
	68	20		
Ga	69	0.3		
	71	0.2		
Ge	72	4		
As	75	20		
Se	77	30		
	78	400		
	82	300		
Rb	85	1		KrH ⁺
Sr	88	0.7		
Y	89	0.2		
Zr	90	0.4		
Nb	93	0.8		
Mo	98	0.4		
Ru	101	0.4		
Rh	103	0.1		
Pd	108	0.3		
Ag	107	0.6		

Table 2 (cont.): Typical IDLs for the aurora M90 ICP-MS.

Element	Isotope (m/Z)	Measured in		Background species
		Hot plasma (ng/L)	Cold plasma (ng/L)	
Cd	111	0.2		
In	115	0.1		
Sn	118	7		
Sb	121	0.1		
Te	125	4		
Cs	133	0.4		
Ba	138	0.2		
La	139	0.4		
Ce	140	0.06		
Pr	141	0.06		
Nd	146	0.3		
Sm	147	0.2		
Eu	153	0.1		
Gd	157	0.2		
Tb	159	0.09		
Dy	163	0.2		
Ho	165	0.06		
Er	166	0.18		

Element	Isotope (m/Z)	Measured in		Background species
		Hot plasma (ng/L)	Cold plasma (ng/L)	
Tm	169	0.08		
Yb	172	0.3		
Lu	175	0.05		
Hf	178	2		
Ta	181	0.2		
W	182	1.3		
Re	185	0.2		
Ir	193	0.2		
Pt	195	0.3		
Au	197	0.3		
Hg	202	1		
Tl	205	1		
Pb	206+	0.3		
Bi	209	0.3		
Th	232	0.04		
U	238	0.06		

References

- [1] I. Kalinitchenko, Ion Optical System for a Mass Spectrometer, US Patent 6,614,021 B1, 2 September 2003
- [2] S. Elliott, M. Knowles and I. Kalinitchenko, "A New Direction in ICP-MS", Spectroscopy, 19(1), 30 (2004).
- [3] V. Thomsen, D. Schatzlein, and David Mercurio, "Limits of Detection in Spectroscopy", Spectroscopy 18(12), 112 (2003)

Authors

XueDong Wang and Shane Elliott

Keywords

Sensitivity
Detection Limit
Interference
Contamination

Instrumentation

aurora M90 ICP-MS

For research use only. Not for use in diagnostic procedures.


● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com


www.bruker.com/chemicalanalysis



The new **ultraflex**treme

● Beyond Imagination

Explore Your Applications with MALDI-TOF



The new ultrafleXtreme™ revolutionizes what you can achieve with a MALDI-TOF and opens new fields for research and industry.

Applications such as in-depth characterization of Biopharmaceuticals, high speed tissue imaging based biomarker discovery and glycoproteomics extend beyond classical proteomics utilizing a multitude of technology innovations and unique proprietary software.

MALDI Imaging – the complete molecular histology solution from the technology leaders

- Top-Down Biomarker Discovery
- New Bottom-Up ImagerID™ workflow for protein identification and localization – even suitable for FFPE tissue and large proteins
- Highly enhanced data acquisition speed and instrument robustness to perform large scale studies

Proteomics and Glycoproteomics – additional capabilities extending your research scope

- Integrated workflows for LC-MS/MS and 2D gels
- Quantification with any label or label-free Top-Down [1] and Bottom-Up workflow
- Glycoproteomics: Integrated GlycoQuest™ search engine for automated glycan structure elucidation in proteins for high-throughput studies

BioPharma – intact protein MW determination and direct protein sequencing as never before with MALDI-TOF

- Accurate assignment of intact protein molecular weights by novel multiple charge state analysis
- Fast and information-rich protein sequencing
- Side product or terminal variant identification
- Characterization of PEGylated proteins

Glycoprotein Analysis – the next life science frontier at your finger tips

- Full glycan analysis
- N-Glycopeptide targeting in complex mixtures
- Glycopeptide structure elucidation by unique MS/MS patterns in MALDI-TOF and automated glycan database search with GlycoQuest

MALDI Biotyper – microbial identification for the 21st century

- Bacteria ID
- Yeast and fungi ID
- MALDI Sepsityper Kit for bacteria and yeast ID from positive blood culture bottles within 30 min

● More Flexible Than Ever

Accurate Intact Protein Molecular Weight Analysis – completing the picture for biopharmaceutical characterization

- Up to 100x increase of mass accuracy down to the low ppm range for intact protein MW determination due to multiple charged ion signals
- Bruker DHAP matrix proven for Top-Down analysis producing multiple charged ions
- Greater than 30,000 resolution for intact proteins due to FlashDetector™

Gold Standard in Top-Down Protein Sequencing

- Monoisotopic sequence readout of up to 90 residues from N- and C-termini
- Automated, rapid protein QC reports from BioPharma Compass software
- Patented T³-Sequencing for direct N- or C-terminal characterization [2]

High Definition MALDI Imaging

- Continuing innovation leadership with patented enabling technologies from sample preparation to the molecular histology workflow, advanced data treatment and automated tissue classification
- Access to FFPE tissue and large proteins: New proprietary ImageID workflow provides identities and localization of ~ 80% of observed peptides in tissue
- 2 kHz FlatTop smartbeam-II™ laser for non-destructive tissue analysis permitting seamless correlation with high resolution histopathology from the same tissue slice
- High definition MALDI Imaging by industry leading 10 µm laser focus and true edge-to-edge pixel definition
- High spatial resolution for protein images at best image brightness by patented smartbeam laser technology

Targeted Protein Quantification with 1-2% CV [4]

- HPLC-free MALDI-TOF for highest assay speed and accuracy

LC-MALDI – the workhorse for proteomics and glycobiology

- smartbeam-II laser for acquisition of maximum number of MS/MS spectra per fraction
- Highest sensitivity through AnchorChip™ technology – now with up to 1536 spots/plate
- Market leading mass resolution (> 40,000) and accuracy (< 1 ppm)
- ProteinScape™ software covering all proteomics and glycoproteomics workflows
- Direct analysis of glycans and glycopeptides enabled by GlycoQuest

Read more

- [1] Maltmann DJ *et al.* *Proteomics*. 2011;11(20):3992-4006.
- [2] Suckau D, Resemann A. *Anal Chem*. 2003;75(21):5817-24.
- [3] Resemann a *et al.* *Anal Chem*. 2010;82(8):3283-92.
- [4] Anderson NL *et al.* *J Proteome res*. 2012;11(3):1868-78.




Extreme Performance for Bottom-Up Proteomics

The ultrafleXtreme's next-generation MS technologies enable greater proteome coverage than ever before. The ultrafleXtreme offers unparalleled cutting-edge capabilities for Bottom-Up proteomics at your fingertips.

Bruker's 25 years Life Science experience and excellent bioinformatics compliant with HUPO/PSI publication guidelines effortlessly converts unrivalled MALDI-TOF/TOF data into published results.

Identify, quantify and characterize – all your results in the blink of an eye

ProteinScape software makes it possible to access, process, merge, query and report all Bottom-Up proteomics data from high throughput 2D gel identification to label or label-free LC-MALDI quantification.



5 μ g HeLa Cell digest resulting in 2890 protein IDs

First-Class LC-MALDI results from even the most complex Bottom-Up proteomics analysis

5 μ g HeLa Cell digest




- **2890 Protein IDs @ 1% Peptide FDR**
- **11,960 Spectrum matches**
- **38% Spectrum ID rate**





Providing a perfect complement to ESI protein ID workflows

Identification, structural characterization and quantification at a glance



ProteinScape – From data to knowledge

- Data from various workflows (multiple digest enzymes, label or label-free quantification, ESI or MALDI) can be merged into confident, non-redundant protein information
- Just a few mouse clicks to compare results across multiple studies
- Direct access to web resources relates information such as gene ontology or mapping of PTMs on the 3D protein structure




Cross the Next Frontier – Glycoproteomics Expertise Delivered


Get the full picture

The ultrafleXtreme TOF/TOF spectra of N-glycopeptides contain both peptide and glycan information. ProteinScape uses a TOF/TOF specific signature fragmentation pattern generated by intra-ring cleavage of the peptide-binding sugar to assign the aglycone and the glycan molecular weights and triggers their identification in parallel database searches.

The specific signature enables the classification and filtering of glycopeptides from non-glycosylated ones in complex cell lysates. ProteinScape integrates the GlycoQuest search engine for glycan structure DB searches.



LC-MALDI analysis of a human antibody digest, shows at the first glance the separated glycopeptide region on the MS-level.



Automated glycopeptide detection by ProteinScape. Signature ion pattern provides molecular weights for both the peptide and glycan from a single glycopeptide TOF/TOF spectrum. Structures from both moieties are assigned in this spectrum.


Easily and comprehensively ID and assign glycan structures in MS/MS spectra



Representation of an ultraflex glycoprotein analysis in ProteinScape. All glycopeptides were detected in the LC-MS/MS dataset and the glycan were searched using GlycoQuest.

A structure candidate is shown as structure and its fragments annotated to a matching MALDI-TOF/TOF spectrum permitting manual inspection whenever requested.

Analysis of N-glycosylation pattern present on antibodies



Bruker's ProteinScape software detects glycopeptide-specific fragmentation patterns in MALDI-TOF/TOF spectra that allows a separate glycane database search with GlycoQuest to ID the glycane part and a Mascot search to ID the peptide part. ProteinScape can then visualize both results to merge the information on the glycopeptide level.

From Tissue Biomarker Discovery...


MALDI Imaging is emerging into an indispensable discovery tool for molecular tissue and histopathological research, reaching out from high resolution imaging tools to statistical cohort analysis and biomarker discovery and identification.

Bruker's MALDI Molecular Imager solution enables robust and quality-controlled tissue preparations, and links the classical histopathology workflow with the molecular dimension. Unique, patented software and techniques for clinical studies enable powerful biomarker discovery.


The ultrafleXtreme is an essential element in any imaging project.

Rat testis at 20 μm spatial resolution


See how this feature allows interactive exploration of histology and molecular information at <http://www.youtube.com/watch?v=YV2wzVXEvgq?hd=1>



MALDI image of rat testis at 20 μm spatial resolution. Two molecular signals that highlight different features are selected.



Bruker flexImaging software allows cross-fading between the MALDI image and the superimposed virtual slide.




At higher resolution, histological features become visible in the virtual slide image.



Full resolution reveals (1) a capillary vessel (2) nuclei of mature spermatids and demonstrates correlation of histological features and molecular signals.


Virtual microscopy



Digital microscopy enables detailed evaluation of high-resolution histological images on a computer screen.




The alignment of MALDI Imaging data with other images (e.g. histological stains) is crucial for the interpretation of results. Bruker's proprietary solution for virtual microscopy allows co-registration of multiple image modalities with unparalleled ease.



From beautiful pictures to true biomarker discovery – MALDI Imaging is perfectly suited to analyze biomarkers in tissue samples. Bruker’s unique suite of software tools supports the entire workflow, including statistical analysis. Comprehensive bioinformatic tools for statistical analysis of MALDI Imaging data, such as hierarchical clustering, PCA or pLSA allow researchers to mine data efficiently.

... To Identification and Validation

MALDI Imaging covers a multitude of application areas ranging from clinical samples to plant material and animal tissue




A: Human bladder cancer specimen; B: Barley seed; C: Rat testis

Identification of biomarker candidates are essential for validating results and the translation of MALDI Imaging derived biomarkers to subsequent clinical stages, such as treatment prognosis or survival prediction. Top-Down discovery and identification has been demonstrated successfully*.

As well, higher MW protein detection and analysis of formalin-fixed tissues are now enabled by the novel ImageID workflow. In this Bottom-Up workflow highly resolved images from on-tissue digests are merged with classical LC-MS/MS identification.


*References:

Lemaire 2007 Reg- α_2 , Hanrieder 2011 S100-A10, Hardesty 2011 S100-A6, Nipp 2011 S100-A6/10, Meistermann 2006 TrT, Rauser 2010 CRIP1, Lagarrigue 2011 Thym β , LCFA-CoA




Consecutive sections are processed in parallel to provide both protein distribution and identification simultaneously. A tryptic digest is performed maintaining the resolution of protein distributions. Peptide imaging (right) and LC-MALDI-MS/MS (left) yield protein and peptide IDs that is software merged for visualization of protein distribution. 80% of all peptides visible in the images can thus be identified and turned into protein localization of more than 100 proteins per image.

ImageID Results



Increased confidence in peptide identifications: Peptides derived from the same protein show the same spatial distribution.



Bottom-Up Imaging with the ImageID workflow extends the mass range of detectable proteins from approx. 25 kDa to greater than 100 kDa and grants access to proteins cross-linked by formalin fixation.




Biopharmaceutical Characterization

Ultra rapid product confirmation of 150 kDa (MOPC21) antibody

Peptide Mass Fingerprints confirm product identity in 5 seconds without LC-separation

Confident characterization at 1 ppm mass accuracy
 > 40,000 resolution
 Immediate visualization and automated reporting

5 second analysis



- Tryptic digest of MOPC21 IgG
- Ultra fast product confirmation

Rapid combination of data from multiple digests & instruments in 3 clicks in ProteinScope for ultimate sequence coverage. BioPharma Compass software – Fully automated acquisition, processing and reporting of data.

Batch report view for automated QC PMF

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
B	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
C	○	○	○	○	○	○	○	○	○	○	○	○	●	●	●	●	●	●	●	●	●	●	●	●
D	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
E	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
F	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○

Sequence coverage acceptance criteria can be altered by user depending on desired result.


Intact Molecular Weight Analysis

MALDI intact molecular weight analysis

Intact Myoglobin (17 kDa)


- Isotopic resolution up to 30 kDa

5 second acquisition



Intact molecular weight analysis confirms presence of product and/or contaminants in 5 seconds.

High accuracy MW determination of CA based on double charged ions



- < 15 ppm mass accuracy for intact proteins
- Using multiple charged ions
- Generated in DHAP matrix



Full resolution @ full sensitivity with high mass accuracy is unique to Bruker MALDI-TOF.

N/C-Terminal Sequencing

MALDI Top-Down N/C-terminal Sequencing

1 minute acquisition

\$ 0.10 Cost



Top-Down sequencing delivers N/C-terminal sequence information in 1 minute without proteolytic digestion. Terminal modifications can be confirmed with an additional Bruker-patented T³-Sequencing step.

Characterization information generated by TDS for IgG



- Glycan heterogeneity observed at intact level
- > 140 amino acid residues confirmed for light & heavy chain
- N/C-termini identified
- N-terminal pyro-Glu detected
- C-terminal Lys deletion detected
- Isotopically resolved TDS data = characterization certainty
- Automated acquisition & reporting of data

Technical Specifications

Redefining MALDI-TOF/TOF Performance

- 2 kHz smartbeam-II™ laser technology enables ultra-high data acquisition speed
- Laser focus diameters down to 10 µm for high spatial resolution imaging without pixel overlap
- Wide mass range resolution up to 40,000 due to proprietary PAN™ technology delivers:
 - Low ppm accuracy for protein MW determination
 - Long and reliable Top-Down sequence readout, at monoisotopic resolution
- Patented AnchorChip™ technology for unmatched consistency and sensitivity levels
- FlashDetector™ provides 1 ppm mass accuracy for highest confidence
- MALDI Perpetual™ Ion Source with entirely automated self-cleaning in < 15 minutes using a patented IR-laser process.
- Bruker DHAP matrix proven for Top-Down analysis

Optional bioinformatics packages

- BioPharma Compass for comprehensive Biopharmaceutical characterization and QC
- flexImaging™, the leading MALDI Imaging platform – now with enhanced data compression for larger scale studies
- ProteinScape™ for the full execution of proteomics projects featuring GlycoQuest™ search engine for glycoproteomics
- Cutting-edge tools for protein analysis and Top-Down Sequencing
- PolyTools™: Interactive interpretation of MALDI polymer spectra
- Support of a wide variety of external software tools and data export functionality

Support features

- Comprehensive self diagnostics and software update checks
- Remote on-line service and support capability
- Compass Security Pack™: Assisting in 21CFR part 11 compliance
- IQ/OQ/PV procedures for regulated environments
- Maintenance contracts available at various levels

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
ms-sales@bdal.com

www.bruker.com/ultraflexxtreme