Coapt

Manual Surgical Drill Instructions for Use

English	1
Dansk	2
Deutsch	4
Español	6
Français	8
Italiano	10
Nederlands	12
Português	13
Svenska	15
Polski	17
Magyar	19
Ελληνικά	21
Русский	23

Manual Surgical Drill™

Available Models

Part No.	Part Name
CFD-899-4000	Manual Surgical Drill™

INTRODUCTION

This document presents the “Instructions for Use” (IFU) for the Coapt Manual Surgical Drill™. These instructions outline the proper operation, care, and maintenance for this surgical instrument. Following this IFU will insure long life, patient safety, and proper instrument performance.

DESCRIPTION

The Coapt Systems Manual Surgical Drill™ (Drill) is a reusable drill and Chuck Key designed for manual drilling of bone. This drill is for use with the Endotine family of products or compatible drill bits (not included).

INDICATIONS FOR USE

The Coapt Manual Surgical Drill™ is limited to use in conjunction with the ENDOTINE Drill Bits to drill bone holes for placement of the following products; ENDOTINE Forehead™, ENDOTINE Midface™ B, ENDOTINE Triple™ and ENDOTINE TransBleph™. Otherwise, the device is indicated for general use in drilling bone holes.

PRECAUTIONS

When the Coapt Manual Surgical Drill™ is used with the ENDOTINE Forehead™ 3.0/3.5 Devices, please refer to the ENDOTINE Forehead™ Device Instructions for Use for special instructions to be exercised by the practitioner.

OPERATION OF DRILL

The Coapt Systems Manual Surgical Drill™ should be handled and operated by personnel completely familiar with its use and care.

HANDLING

The drill should be inspected to insure proper functionality prior to each use with particular attention paid to all moving parts. If debris emanates from the instrument, is identified in any of the device concavities, or if the device is not functioning smoothly, other alternatives should be employed for operating

the desired drill bit. If any damage or malfunction is detected, the Drill should not be used and should be sent to Coapt Systems or an authorized service center for repair. Failure to make a complete inspection may result in unsatisfactory performance.

Use of the Drill for a task other than that to which it is intended may result in damage, unsatisfactory performance, or patient risk. In order to maintain Coapt Systems warranty, drill **REPAIR MUST BE PERFORMED BY COAPT SYSTEMS** or an authorized service center.

INSTRUCTIONS FOR USE

Operating the Coapt Systems Manual Surgical Drill™ consists of the following steps. Caution: Always perform drill inspection per “Handling” instructions above prior to use.

1. Insert the desired drill bit into the Drill Bit Chuck using the Chuck Key located in the Twisting Grip.

NOTE: Compatible Drill Bit diameters specified on Drill Bit Chuck

2. Grasp the Stationary Grip with one hand.
3. Position the drill bit in the desired location for drilling.
4. Rotate the twisting grip with other hand.

CARE OF DRILL

Prior to each surgical procedure, the Drill and Chuck Key should be cleaned and sterilized as instructed below.

As with any reusable medical device, personnel should follow accepted guidelines such as **ANSI/AAMI ST46:2002 Steam Sterilization and Sterility Assurance in Health Care Facilities and AORN Standards, Recommended Practices, and Guidelines.**

It is important to rinse instruments that have been exposed to blood and saline solution before these substances dry. Blood and saline solutions are corrosive. In addition, blood will also cause a stain that is difficult to remove. The Drill and Chuck Key should, therefore, be kept moist immediately following use, by covering it with a moist towel prior to cleaning. **DO NOT SOAK OR IMMERSE.**

CLEANING

Remove the Chuck Key from the Chuck Key bay. Using a hand-held, soft bristle brush, clean the Drill Bit Chuck, Chuck Key, and Drill body using warm running water and a neutral pH detergent designed for cleaning surgical instruments. Make sure to clean the Chuck Key bay located in the end of the Twisting Grip. Pay close attention to crevices or other areas that could retain contaminants such as fluids or blood. During cleaning, never use steel wool, wire brushes, or highly abrasive detergent cleaners as they will damage the surface finish of the hand drill.

RINSING

Rinse the Drill using warm running water and wipe dry with a clean, soft cloth. Do not immerse the Drill in cleaning solutions, basins, or ultrasonic cleaning devices. For best

results and to avoid residue build up, rinse the Drill with distilled or demineralized water. Be sure to thoroughly rinse the inside of the Chuck Key bay.

INSPECTION

Visually inspect the Drill for cleanliness paying special attention to crevices, recesses, and the Chuck Key bay.

DRYING

Prior to sterilization and storage, the Drill must be thoroughly dry. Remaining moisture may be removed from hard to reach areas by using compressed air.

LUBRICATION

Rotating and sliding surfaces of the Drill Bit Chuck should be lubricated with a water soluble product (such as Instrument Milk or similar lubricant).

PACKAGING FOR STERILIZATION

Prior to sterilization, the Chuck Key should be removed from the Chuck Key bay. The Drill is now ready for sterilization per **ANSI/AAMI ST46:2002 Steam Sterilization and Sterility Assurance in Health Care Facilities**.

STERILIZATION

Independent laboratory testing has validated steam sterilization as an effective sterilization process for the Coapt Systems Manual Surgical Drill™. As recommended by the aforementioned standards the sterilizer manufacturer's written instructions for cycle parameters should be followed in order to achieve adequate instrument exposure.

The following table outlines the minimum cycle temperature and exposure time for sterilization per **ANSI/AAMI ST46:2002 Steam Sterilization and Sterility Assurance in Health Care Facilities**.

Steam Sterilization Method	Temperature	Exposure Time ¹	Dry Time
Vacuum Displacement	270°F / 132°C	4 minutes	30 minutes
Gravity Displacement	250°F / 121°C	30 minutes	45 minutes

MAINTENANCE OF DRILL

THE COAPT MANUAL SURGICAL DRILL™ IS GUARANTEED FOR ONE YEAR. Always inspect the Drill prior to use per "Operation of Drill" section above. If any damage or malfunction is detected, the Drill should not be used and should be sent to Coapt Systems or an authorized service center for repair.

PROTECT INSTRUMENTS

The use of neutral pH detergents is vital to the care of surgical instruments. Contact with acidic or alkaline solutions will remove the instruments' protective barrier and lead to staining or corrosion;

The Manual Surgical Drill™ should never be exposed to: Aqua Regia, Iodine, Hydrochloric Acid, Ferric Chloride, or Sulfuric Acid.

The following substances should be avoided whenever possible: Aluminum Chloride, Barium Chloride, Bichloride of Mercury, Calcium Chloride, Chlorinated Lime, Carbolic Acid, Dakin's Solution, Mercury Chloride, Potassium Permanganate, Potassium Thiocyanate, Saline, Sodium Hypochlorite, or Stannous Chloride.

¹Exposure time: The period during which items are exposed to saturated steam at the specified temperature." (ANSI/AAMI ST46:2002)

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Coapt Systems, Inc.