



Rapid and Cost Effective Immunoassay Development and Validation of a Biomarker Assay Using Spatial Proximity Analyte Reagent Capture Luminescence, SPARCL™

Mark J. Cameron and Wenhua Xie - Lumigen, Southfield Michigan

ABSTRACT

A novel detection technology was used for rapid immunoassay (ligand binding assay) development and validation. SPARCL (Spatial Proximity Analyte Reagent Capture Luminescence) technology is a proximity dependent, homogeneous, chemiluminescent detection method. In a SPARCL assay, a chemiluminescent substrate (acridan) is brought into the proximity of an oxidative enzyme (horseradish peroxidase) through the specific antigen/antibody interaction. A flash of light proportional to the quantity of analyte present in the sample is generated upon addition of a trigger solution containing H₂O₂. Commercially sourced reagents for human IL-8 were purchased and used in assay development and validation. Assay development was completed in less than 3 days. The validated assay features a wide dynamic range (3 – 4,000 pg/mL), low minimum required dilution, short assay run time (60 minutes), and sensitivity (3 pg/mL). Quality control samples (QC's) run in human plasma were used to derive statistics that show the assay is precise and accurate with a total error of all QC levels of under 12%. Furthermore, beyond a luminometer with injectors, the SPARCL assay requires no specialized or capital equipment expense. The SPARCL technology enables rapid immunoassay development, delivers a biomarker assay with desirable performance characteristics and at the same time allows for considerable savings in labor, disposables and capital equipment in development, validation and in sample analysis compared to other immunoassay platforms. Existing software can be used for a SPARCL assay data acquisition, and as a result, there may not be a need to validate new software in existing Laboratory Information Management Systems (LIMS). The SPARCL assay is microtiter plate-based and is suitable for use with robotics and in higher throughput applications and may prove to be a valuable assay development and production tool.




Figure 1. Key components to a SPARCL assay. **A** - HRP labeled Antibody, **B** - Immunoassay target (analyte), **C** - 96 well low binding white plate, **D** - Acridan labeled Antibody, **E** - SPARCL Kit, **F** - Luminometer.




Figure 2. A representative SPARCL assay. Specific antibody and antigen bring acridan and HRP into close proximity. The addition of hydrogen peroxide based trigger solution interacts with acridan and HRP, causing a flash of light.

MATERIALS AND METHODS

Antibodies: A mouse anti-human IL-8 monoclonal antibody was purchased from R&D Systems (Minneapolis MN), catalog number MAB208 and was conjugated with acridan. A HRP-conjugated mouse anti-human IL-8 monoclonal antibody was purchased from Anogen, a division of YES Biotech Laboratories (Mississauga Ontario, Canada), catalog number MO-C40017T.

Calibrator: Recombinant human IL-8 was purchased from R&D Systems (Catalog Number 208-IL).

Matrix: Human plasma was obtained from Bioreclamation. The fetal calf serum and RPMI-1640 was obtained from Sigma.

SPARCL Kit: The SPARCL kit was obtained from Lumigen (Southfield, Michigan). The kit is composed of a hydrogen peroxide based trigger solution, a background reducing agent and acridan.

Acridan Labeling:

Add 500 µL of DMF to Lumigen SPARCL labeling reagent vial.

In a separate 1.5 mL microfuge tube mix the following:

- 708.7 µl of 0.05M sodium borate buffer, pH 8.5
- 41.3 µl of Lumigen SPARCL labeling reagent in DMF (from step 1)
- 250 µl of 1 mg/mL antibody (approximately 160 kDa)

Mix by inverting the tube 4-5 times, cover the tube with aluminum foil and let stand for 30 min. at room temperature. Place the labeling reaction tube on a rocker to mix at 2-8°C overnight.

Assay Buffer: PBS containing 0.1% BSA was used for dilution of the conjugated antibodies.

Assay Development: A fit-for-purpose approach was used for assay development and validation¹.




Figure 3. A typical SPARCL pattern of flash luminescence. Light energy was captured every 0.02 seconds for a total time of 1.0 second. Note peak luminescence at 0.15 seconds after trigger injection.





Figure 4. Screen capture of Softmax Pro 5.4 data acquisition software showing typical standard curve. Note the wide dynamic range of 3 to 4000 pg/mL.

Characteristic	Statistic	QC 200	QC 100	QC 50	QC 25	QC 12.5	QC 6.25
# Results	Mean (N)	6	6	6	6	6	6
Accuracy	Mean Bias (%RE)	0.75	3.5	0.8	3.6	4.8	2.4
Precision	%CV	6.1	3.4	1.9	11.1	4.9	7.7
Total Error	%RE + %CV	6.9	6.9	2.7	14.7	9.7	10

Table 1. Inter assay precision and accuracy. QC's were made in neat matrix. Note that the total error ranged from 2.7 to 14.7 percent.

LUMIGEN SPARCL™ DETECTION KIT



- A: Lumigen SPARCL Labeling Reagent**
- Activated NHS ester of an acridan compound
- B: Lumigen SPARCL BGR**
- Background reducing agent used to enhance S/N ratios
- C: Lumigen SPARCL Trigger Solution**
- Hydrogen peroxide containing solution for signal generation

* Also included in the kit is a Borate Buffer Pack for the labeling reaction

TYPICAL LUMIGEN SPARCL™ ASSAY WORKFLOW

- Step 1:** Label one of the antibody with the activated SPARCL label. No purification required.
- Step 2:** Add antibody solutions and sample; Incubate for 30 minute
- Step 3:** Add SPARCL BGR
- Step 4:** Read the luminescence signal by injecting SPARCL Trigger Solution on a luminescence plate reader

CONCLUSIONS

SPARCL allows for rapid assay development due to short assay run times, a homogenous assay format and flash luminescence.

The SPARCL assay has a wide dynamic range (2.5+ logs), is sensitive (LLOQ of 3 pg/mL) and has a low MRD (1:10).

SPARCL assays are cost effective in terms of capital equipment, labor, reagent usage, disposables, buffers and service contracts on equipment

The SPARCL technology is less complicated than other commonly used platforms while providing the performance of commonly used platforms.

Reference

1. Lee JW, Devanarayan V, Barrett YC, Weiner R, Allinson J, Fountain S, Keller S, Weinryb I, Green M, Duan L, Rogers JA, Millham R, O'Brien PJ, Sailstad J, Khan M, Ray C, Wagner JA. Fit-for-purpose method development and validation for successful biomarker measurement. Pharm Res. 2006 Feb;23(2):312-28

* All trademarks are property of their respective owners. Lumigen, the stylized logo, and SPARCL are registered trademarks of Beckman Coulter, Inc.