


Design, Development and Applications of the MycoSEQ™ Mycoplasma Scan (Myco Scan) Assay



Surekha Karudapuram, Michael T. Brewer, Olga V. Petrauskene, Nikolay V. Sergeev, Somaya Bit, Maxim G. Brevnov, Pius Brzoska, Sueh-Ning Liew, Jen-Kuei Liu and Manohar R. Furtado. Life Technologies Inc., 850 Lincoln Center Drive, Foster City CA

Introduction

Mycoplasma contamination remains a major problem for researchers using mammalian cell culture. Many studies have underlined the need for routine screening of cell lines since Mycoplasma infection is common in laboratories and can affect virtually every cell parameter, subsequently impacting the results of studies performed on those cell lines. Mycoplasma infection have been shown to deplete nutrients, promote metabolic accumulation resulting in pH shifts, induce or suppress cytokine expression, alter metabolism, proliferation characteristics and the morphology of cultured cells¹.


Mycoplasmas pose a challenge to detect and eliminate because of their small cell size, wide diversity of species, limited turbidity produced in culture, requirements for enriched culture media and some species/strains may be cell-invasive. All of these properties have posed a challenge to the development of a reliable rapid molecular test for Mycoplasma detection.

In this poster we present results from the Myco Scan Kit from Applied Biosystems that can be used for rapid detection Mycoplasma contamination using Real-Time PCR. Through intensive bioinformatics and highly optimized multiplexed primer design, the assay allows for highly sensitive, specific and comprehensive Mycoplasma species detection and prevents the detection of other closely related species. Providing results in less than five hours, the kit is a critical tool in testing the authenticity of valuable cell lines used in basic research.

Rapid 4 hour protocol from DNA extraction to presence absence call

PrepSEQ™ 1-2-3 Sample Prep (1 hour)


RNase/Proteinase treatment of sample




Real-Time PCR Assay Set Up (30 minutes)

Component for premix solution	Volume for one 50-µL reaction (µL)	To prepare...	Combine in each tube or well...
Power SYBR® Green PCR Master Mix (2X)	15.0	Negative-control reaction	• 15 µL of Primix Solution • 15 µL of Negative Control (water)
Mycoplasma Real-Time-PCR Primer Mix (10X)	3.0	Your unknown sample reaction	• 15 µL of Primix Solution • 2 to 10 µL of unknown sample • Adjust the final reaction volume to 50 µL with Negative Control (water)
Total premix solution volume	18.0	Inhibition-control reaction	• 15 µL of Primix Solution • 2 to 10 µL of unknown sample • 1 µL of Mycoplasma Real-Time-PCR DNA Control (positive control) • Adjust the final reaction volume to 50 µL with Negative Control (water)
		Positive-control reaction	• 15 µL of Primix Solution • 1 µL of Mycoplasma Real-Time-PCR DNA Control (positive control) • 10 µL of Negative Control (water)

Real-Time PCR Run and Data Review (2.5 hours)




Results from Mycoplasma positive and negative samples



Highly optimized Power-SYBR® Green detection technology allows analysis of multiple parameters (Ct and Tm) for results interpretation


Sensitivity down to 1-10 copies/reaction






10-fold serial dilution of purified genomic DNA was tested using the assay

High sensitivity extraction of M. bovoculi gDNA using PrepSEQ™ 1-2-3 Kit

Commercial Sample Preparation



PrepSEQ™ Sample Preparation



M. Bovoculi DNA was spiked in cells and extracted using the PrepSEQ kit

Broad species detection

Detection of over 90 species

<i>Acholeplasma granularum</i>	<i>Mycoplasma genitalium</i>	<i>Mycoplasma synoviae</i>
<i>Acholeplasma laidlawii</i>	<i>Mycoplasma gypis</i>	<i>Mycoplasma testudinis</i>
<i>Acholeplasma pleciae</i>	<i>Mycoplasma hominis</i>	<i>Mycoplasma timone</i>
<i>Mycoplasma alkalescens</i>	<i>Mycoplasma hyorhinis</i>	<i>Spiroplasma citri</i>
<i>Mycoplasma alvi</i>	<i>Mycoplasma imitans</i>	<i>Spiroplasma endosymbiont</i>
<i>Mycoplasma anseris</i>	<i>Mycoplasma indiensis</i>	<i>Spiroplasma insolitum</i>
<i>Mycoplasma arginini</i>	<i>Mycoplasma lagogenitalium</i>	<i>Spiroplasma kunkelii</i>
<i>Mycoplasma auris</i>	<i>Mycoplasma lipofaciens</i>	<i>Spiroplasma melliferum</i>
<i>Mycoplasma buccale</i>	<i>Mycoplasma mobile</i>	<i>Spiroplasma mirum</i>
<i>Mycoplasma californicum</i>	<i>Mycoplasma molare</i>	<i>Spiroplasma phoeniceum</i>
<i>Mycoplasma canadense</i>	<i>Mycoplasma mycoides</i>	<i>Spiroplasma poulsonii</i>
<i>Mycoplasma capricolum</i>	<i>Mycoplasma neurolyticum</i>	<i>Spiroplasma sp.</i>
<i>Mycoplasma caviae</i>	<i>Mycoplasma orale</i>	<i>Mycoplasma bovirhinis</i>
<i>Mycoplasma collis</i>	<i>Mycoplasma phocidae</i>	<i>Mycoplasma bovis</i>
<i>Mycoplasma cricetuli</i>	<i>Mycoplasma pirum</i>	<i>Mycoplasma bovigenitalium</i>
<i>Mycoplasma equirhinis</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma canis</i>
<i>Mycoplasma fermentans</i>	<i>Mycoplasma salivarium</i>	<i>Mycoplasma felis</i>
<i>Mycoplasma gallinaceum</i>	<i>Mycoplasma simbae</i>	<i>Mycoplasma fastidiosum</i>
<i>Mycoplasma gallisepticum</i>	<i>Mycoplasma sp.</i>	<i>Mycoplasma muris</i>
<i>Mycoplasma gateae</i>	<i>Mycoplasma spumans</i>	<i>Mycoplasma pulmonis</i>

*No cross reactivity with Human, CHO, mouse and different microorganisms such as *E. coli*, *Bacillus cereus*, *Candida albicans*, *Staphylococcus aureus*, *Micrococcus luteus*, *Clostridium perfringens*

Summary

The Applied Biosystems Myco Scan Kit can be used for rapid and reliable detection of Mycoplasma contamination in mammalian cell culture:

- Designed for routine screening of high value cell lines for Mycoplasma contamination using highly sensitive Real-time PCR
- Rapid Time to Results in under five hours
- Detection of greater than 90 species
- Demonstrated sensitivity and specificity
- High DNA recovery sample preparation



For Research Use Only. Not for use in diagnostic procedures.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

© Copyright 2009, 2011. Life Technologies Corporation.

¹ Kagemann et al., Impact of Mycoplasma hyorhinis infection on L-arginine metabolism: differential regulation of the human and murine iNOS gene. Biol. Chem. 2005 Oct; 386(10):1055-63.