

## 1. What is *SDS-eliminant* ?



AE-1390 *SDS-eliminant* is the reagent for removing SDS (Sodium dodecyl sulfate) from samples including SDS by 2 steps. At first, Solution A removes dodecyl sulfate ion to precipitate it by insoluble pellet. The SDS also bound to proteins directly, however, could not be removed by Solution A. Second, the non-ionic surfactant including Solution B makes mixed micelle with SDS. This step additionally removes the SDS direct-binding to the proteins. *SDS-eliminant* is extremely useful method in sample isolation, fractionation and purification.

### ● Purpose

- Removal of SDS from the sample solution
- Removal of SDS from the proteins isolated SDS-PAGE
- Removal of SDS from the cell- or tissue-extracts
- Preparation of Mass Spectrometry, electrophoresis and chromatography
- Re-activation of enzymes\* denatured by SDS \* not all enzymes


### ■ Feature

- Ready to Use
- Over 95% of the ionized SDS is removed briefly
- The recovery of proteins is high
- Available for the samples are small in amount (less than 1 mL)

## 2. Specification

	Solution A	Solution B
Contents	2mL solution / a tube	2mL solution / a tube
Major constituent	KCl, Tris	non-ionic surfactant, Tris
Max sample volume	40 mL	10 mL
Conservation	Cold storage in light-blocking bag, one year guarantee	

## 3. Principle



### A. Sample solution before SDS denaturation

… The protein is not denatured, the construction or activation is kept.

### B. SDS denaturation of proteins

… The protein (represented as pink figure) is denatured\*<sup>1</sup> by SDS (represented as yellow figure, especially yellow clover-shaped shows SDS-micelle), exists in sample solution by two types : dodecyl-sulfate ion or direct-binding type.

### C. Precipitation of dodecyl-sulfate ion

… Add *SDS-eliminant* solution A to Sample solution and mixing gently. Incubate for half an hour at 4 degree C. The potassium ion including solution A (represented as green triangle) binds to dodecyl-sulfate ion and make insoluble pellet. After centrifugation, the supernatant is collected by new tube. The direct-binding SDS, however, is not removed by Solution A.

### D. Addition of non-ionic surfactant

… Add *SDS-eliminant* solution B to collected supernatant and mixing gently. Incubate over 2 hours at 4 degree C. The non-ionic surfactant (represented as red oval) interacts with the direct-binding SDS.

### E. Mixed micellization

… The micelles of non-ionic surfactant incorporate the direct-binding SDS and these forms of mixed micelle. The protein is released and re-activated\*<sup>2</sup> from SDS denaturation.

\*<sup>1</sup> The denaturations of SDS are different from each protein.

\*<sup>2</sup> *SDS-eliminant* can remove the SDS from samples by two steps, but not absolutely.

## 4. Technical data

### I. SDS-rejection ratio & Protein recovery

	BSA	ovalbumin	Casein	Lysozyme	Recombinant Malate Dehydrogenase
M. W.	69.2 kDa	42.9 kDa	25.1 kDa	14.3 kDa	40.0 kDa
Ratio of SDS-rejection	100 %	100 %	99 %	100 %	95 %
Protein recovery	97 %	100 %	100 %	89.5 %	No Data


\* The ratio of SDS-rejection was calculated by the SDS concentration, which measured by the standard curve of Stains-ALL. The plate reader was detected the 450nm absorbance ( $A_{450}$ ) of Stains-ALL which was changed by the SDS amount. It was only the method to measure the amount of dodecyl sulfate ion, not to measure the SDS direct-binding to proteins.

\* Protein recovery was calculated by the Protein concentration, which measured by Pierce® BCA Protein Assay.

### II. Re-activation of enzymes

#### A. Lysozyme


Lysozyme was 92% re-activated comparing to SDS-untreated



These ratio was calculated the 450nm absorbance ( $A_{450}$ ) in bacteriolysis of *Micrococcus lysodeikticus*. The graph bars show the ratio of Lysozyme Re-activation after 60 second of sample adding. These ratios was compared to SDS-untreated sample as 100%.

#### B. Recombinant Malate Dehydrogenase (rMDH)

rMDH was 55.3% re-activated comparing to SDS-untreated




rMDH activity was detected by the 340nm absorbance ( $A_{340}$ ), which was decreased by changing beta-NADH to NAD<sup>+</sup>. The graph bars show the ratio of rMDH Re-activation from 30 to 90 seconds of enzyme adding. These ratios was compared to SDS-untreated sample as 100%.

### III. Competitor's Products


\* The product S. (supplier is T.S. company) was compared with AE-1390 *SDS-eliminant* by 3 points: The ratio of SDS-rejection, Protein Recovery and Re-activation.

\* AE-1390 *SDS-eliminant* showed superior Protein Recovery and Re-activation to the product S.


a. Ratio of SDS-rejection



b. Protein Recovery



c. Re-activation



Each analysis was conducted by the 1% (w/v) SDS-denaturated lysozyme (1mg/mL conc.). Each graph shows the ratios which was compared to the SDS-untreated samples as 100%. The ratio of SDS-rejection was calculated by Stains-ALL. Protein recovery was calculated by Pierce® BCA Protein Assay. Re-activation was calculated by the bacteriolysis of *Micrococcus lysodeikticus*.



**ATTO Corporation**

1-5-32 Yushima, Bunkyo-ku, Tokyo 113-0034, Japan

Tel. +81-3-3814-4824

Fax. +81-3-5684-6646

URL <http://www.atto.co.jp>

Email [eig@atto.co.jp](mailto:eig@atto.co.jp)