



Release Notes

MBT Compass Library Revision L (2020)

covering 3239 species/entries (9607 MSP)

This product has no declared clinical intended purpose and is not for clinical diagnostic use. Any clinical diagnostic use is at the user's own risk and responsibility.

Contents

1	Overview	2
1.1	What is new with the MBT Compass Library?	2
1.2	Main improvements	2
1.3	Additional libraries for Research Use Only workflow	3
1.4	BTyp2.0Sec.Library 1.0 (# 8254705)	3
1.5	Revision status of the MBT Compass Library Release Notes	5
2	New Species	6
3	Deleted MSP entries	17
4	Renaming	19
5	Changes to Matching Hints	29
6	Genera/Species improved by implementation of further reference entries	42
7	Installation Instructions	49
8	Appendix	50
8.1	Matching Hints	50
8.2	Information only required for users of laboratory-specific specimen codes	79
8.3	Additional information for yeast nomenclature update	81
8.4	Additional information for <i>Staphylococcus argenteus</i> / <i>schweitzeri</i>	83



1 Overview

1.1 What is new with the MBT Compass Library?

- **1161 new reference spectra (MSP)** in total.
- **362 MSP** cover **21 new genera** and **278 new species** (see Table 3).
- **799 new MSP** entries **improve the diversity coverage** of genera/species already present in the MBT Compass reference library (see Table 7).
- **9607 MSP** in total

1.2 Main improvements

- New **CLINICAL** species, e.g. *Bordetella* sp., *Clostridium* sp.
- Additional strains of the *Staphylococcus aureus*-related species *Staphylococcus argenteus* and *Staphylococcus schweitzeri* (see also chapter 8.4).
- Extensive increase of species coverage for *Acinetobacter* sp. and *Paenibacillus* sp.
- New **ENVIRONMENTAL** species, e.g. *Aeromicrobium* sp., *Brachybacterium* sp.
- Extensive taxonomical yeast nomenclature update

The 21 new genera and 278 new species cover the following groups:

	New genera	New species	Aerobe species	Microaerophilic species	Anaerobe species
Gram -	15	117	103	8	6
Gram +	6	156	146	-	10
Yeast	-	5	5	-	-

Overall improvements:

	New MSP (from 1161)	Species covered	Aerobe strains	Microaerophilic strains	Anaerobe strains
Gram -	498	269	431	24	43
Gram +	620	436	501	7	112
Yeast	43	25	43	-	-

MBT Compass Library overall numbers:

	MSP	Genus	Species
Gram -	4019	270	1379
Gram +	4672	232	1603
Yeast	851	66	216
Filamentous Fungi	65	25	41
Σ	9607	593	3239



Use the MBT Compass Library revision L with the following MBT RUO software versions:

- MALDI Biotyper 3.1
- MALDI Biotyper 4.0 and 4.0 SR1
- MBT Compass
- MBT Compass HT

Note: The support of software version MALDI Biotyper 3.0 has been discontinued.

1.3 Additional libraries for Research Use Only workflow

New library setups are available and can be installed in parallel with the MBT Compass Library.

Use the following setups after installing MBT Compass Library:

- MBT Mycobacteria Library 6.0 (# 1850766)
- MBT Filamentous Fungi Library 3.0 (# 1829014)

or ask your Bruker Sales contact for further available products.

1.4 BTyp2.0Sec.Library 1.0 (# 8254705)

The BTyp2.0Sec.Library 1.0 contains additional reference spectra of highly pathogenic microorganisms (e.g. *Brucella* sp.), not included in the MBT 8468 MSP Library. It is used for the identification of species listed under the German war weapons control act and therefore is under legal restrictions in relation to its export.

The BTyp2.0Sec.Library 1.0 currently includes 12 species from 10 genera. The 'SR' Library generated after installation will cover the complete MBT SR database content (see Table 1). Due to some limitations of species identification using the MALDI technique, a secondary library ('SR_BBFV') has been created. SR_BBFV is reduced to 4 species/groups (see Table 2).

Table 1: List of species included in 'SR' Library. Bold font indicates the different genera.

No.	Species
1	Bacillus anthracis
2	Brucella melitensis
3	Burkholderia mallei
4	<i>Burkholderia pseudomallei</i>
5	Clostridium botulinum
6	Francisella tularensis
7	Salmonella paratyphi
8	<i>Salmonella typhi</i>
9	Shigella dysenteriae
10	Vibrio cholerae
11	Xanthomonas albilineans
12	Yersinia pestis

Table 2: List of species included in 'SR_BBFV' Library. Bold font indicates the different genera.

No.	Species
1	Brucella melitensis
2	Burkholderia mallei
3	<i>Burkholderia pseudomallei</i>
4	Francisella tularensis
5	Vibrio cholerae

Note: Without installation of BTyp2.0Sec.Library 1.0 the highly pathogenic microorganisms included in this library cannot be correctly identified.

Highly pathogenic microorganisms not included in the BTyp2.0Sec.Library 1.0 (e.g. *Bacillus anthracis*) cannot be identified or will be commented within the Matching Hints of closely related less pathogenic microorganisms included in the MBT 8468 MSP Library (e.g. *Bacillus cereus*).



1.5 Revision status of the MBT Compass Library Release Notes

Previous revisions of the MBT Compass Library Release Notes were handled as individual documents using the MSP DB-XXXX total numbers as a unique identifier.

Starting with the MBT Compass Library DB-7854, this approach will now be altered to better identify the most recent version of the MBT Compass Library Release Notes.

Release Notes	Publication Date	Printed Revision	Corrected Document Revision
MBT Compass Library DB-3995	January 2011	None	A
MBT Compass Library DB-4110	April 2011	None	B
MBT Compass Library DB-4613	July 2012	None	C
MBT Compass Library DB-5627	December 2013	None	D
MBT Compass Library DB-5989	June 2015	A	E
MBT Compass Library DB-6903	April 2016	A	F
MBT Compass Library DB-7311	February 2017	A	G
MBT Compass Library DB-7854	April 2018	H	H
MBT Compass Library DB-7854	November 2018	J	J
MBT Compass Library DB-8468	April 2019	K	K
MBT Compass Library Revision G	November 2020	L	L



2 New Species

Table 3: Implementation of MSP entries for the following 278 new species

	New genus/species				Main relevance
1	<i>Acetobacter orientalis</i>	new species	Gram -	aerobic	FOOD
2	<i>Acidovorax caeni</i>	new species	Gram -	aerobic	ENVIRONMENT
3	<i>Acidovorax cattleyae</i>	new species	Gram -	aerobic	ENVIRONMENT
4	<i>Acidovorax radialis</i>	new species	Gram -	aerobic	ENVIRONMENT
5	<i>Acidovorax soli</i>	new species	Gram -	aerobic	ENVIRONMENT
6	<i>Acidovorax valerianellae</i>	new species	Gram -	aerobic	ENVIRONMENT
7	<i>Acidovorax wautersii</i>	new species	Gram -	aerobic	CLINICAL / ENVIRONMENT
8	<i>Acinetobacter albensis</i>	new species	Gram -	aerobic	ENVIRONMENT
9	<i>Acinetobacter apis</i>	new species	Gram -	aerobic	VETERINARY
10	<i>Acinetobacter boissieri</i>	new species	Gram -	aerobic	ENVIRONMENT
11	<i>Acinetobacter celticus</i>	new species	Gram -	aerobic	ENVIRONMENT
12	<i>Acinetobacter chinensis</i>	new species	Gram -	aerobic	ENVIRONMENT
13	<i>Acinetobacter colistiniresistens</i>	new species	Gram -	aerobic	CLINICAL
14	<i>Acinetobacter cumulans</i>	new species	Gram -	aerobic	ENVIRONMENT
15	<i>Acinetobacter defluvii</i>	new species	Gram -	aerobic	ENVIRONMENT
16	<i>Acinetobacter equi</i>	new species	Gram -	aerobic	VETERINARY
17	<i>Acinetobacter gandensis</i>	new species	Gram -	aerobic	VETERINARY / FOOD
18	<i>Acinetobacter gyllenbergii</i>	new species	Gram -	aerobic	CLINICAL
19	<i>Acinetobacter halotolerans</i>	new species	Gram -	aerobic	ENVIRONMENT
20	<i>Acinetobacter larvae</i>	new species	Gram -	aerobic	VETERINARY
21	<i>Acinetobacter piscicola</i>	new species	Gram -	aerobic	VETERINARY
22	<i>Acinetobacter populi</i>	new species	Gram -	aerobic	ENVIRONMENT
23	<i>Acinetobacter pragensis</i>	new species	Gram -	aerobic	ENVIRONMENT



	New genus/species				Main relevance
24	<i>Acinetobacter pseudolwoffii</i>	new species	Gram -	aerobic	ENVIRONMENT
25	<i>Acinetobacter qingfengensis</i>	new species	Gram -	aerobic	ENVIRONMENT
26	<i>Acinetobacter sichuanensis</i>	new species	Gram -	aerobic	ENVIRONMENT
27	<i>Acinetobacter soli</i>	new species	Gram -	aerobic	CLINICAL / ENVIRONMENT
28	<i>Acinetobacter wuhouensis</i>	new species	Gram -	aerobic	ENVIRONMENT
29	<i>Actinobacillus arthritidis</i>	new species	Gram -	facultative anaerobic	VETERINARY
30	<i>Actinobacillus capsulatus</i>	new species	Gram -	facultative anaerobic	VETERINARY
31	<i>Actinobacillus succinogenes</i>	new species	Gram -	facultative anaerobic	INDUSTRY
32	<i>Actinomyces johnsonii</i>	new species	Gram +	anaerobic	CLINICAL
33	<i>Actinomyces timonensis</i>	new species	Gram +	anaerobic	CLINICAL
34	<i>Aeromicrobium alkaliterrae</i>	new species	Gram +	aerobic	ENVIRONMENT
35	<i>Aeromicrobium camelliae</i>	new species	Gram +	aerobic	ENVIRONMENT
36	<i>Aeromicrobium erythreum</i>	new species	Gram +	aerobic	ENVIRONMENT
37	<i>Aeromicrobium fastidiosum</i>	new species	Gram +	aerobic	ENVIRONMENT
38	<i>Aeromicrobium ginsengisoli</i>	new species	Gram +	aerobic	ENVIRONMENT
39	<i>Aeromicrobium halocynthiae</i>	new species	Gram +	aerobic	ENVIRONMENT
40	<i>Aeromicrobium marinum</i>	new species	Gram +	aerobic	ENVIRONMENT
41	<i>Aeromicrobium ponti</i>	new species	Gram +	aerobic	ENVIRONMENT
42	<i>Aeromicrobium tamlense</i>	new species	Gram +	aerobic	ENVIRONMENT
43	<i>Aliarcobacter thereius</i>	new species	Gram -	microaerophilic	VETERINARY
44	<i>Aliarcobacter trophiarum</i>	new species	Gram -	microaerophilic	VETERINARY
45	<i>Alteromonas australica</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
46	<i>Anoxybacillus amylolyticus</i>	new species	Gram +	aerobic	ENVIRONMENT
47	<i>Anoxybacillus caldiproteolyticus</i>	new species	Gram +	aerobic	FOOD
48	<i>Anoxybacillus contaminans</i>	new species	Gram +	aerobic	FOOD
49	<i>Anoxybacillus rupiensis</i>	new species	Gram +	aerobic	ENVIRONMENT



	New genus/species				Main relevance
50	<i>Anoxybacillus tepidamans</i>	new species	Gram +	aerobic	FOOD
51	<i>Anoxybacillus voinovskiensis</i>	new species	Gram +	aerobic	ENVIRONMENT
52	<i>Bacillus campisalis</i>	new species	Gram +	aerobic	ENVIRONMENT
53	<i>Bacillus chungangensis</i>	new species	Gram +	aerobic	ENVIRONMENT
54	<i>Bacillus foraminis</i>	new species	Gram +	aerobic	ENVIRONMENT
55	<i>Bacillus graminis</i>	new species	Gram +	aerobic	ENVIRONMENT
56	<i>Bacillus kyonggiensis</i>	new species	Gram +	aerobic	ENVIRONMENT
57	<i>Bacillus massiliosenegalensis</i>	new species	Gram +	aerobic	CLINICAL
58	<i>Bacillus niabensis</i>	new species	Gram +	aerobic	ENVIRONMENT
59	<i>Bacillus pocheonensis</i>	new species	Gram +	aerobic	ENVIRONMENT
60	<i>Bacillus timonensis</i>	new species	Gram +	aerobic	CLINICAL
61	<i>Bordetella bronchialis</i>	new species	Gram -	aerobic	CLINICAL
62	<i>Bordetella sputigena</i>	new species	Gram -	aerobic	CLINICAL
63	<i>Brachybacterium aquaticum</i>	new species	Gram +	aerobic	ENVIRONMENT
64	<i>Brachybacterium fresconis</i>	new species	Gram +	aerobic	ENVIRONMENT
65	<i>Brachybacterium ginsengisoli</i>	new species	Gram +	aerobic	ENVIRONMENT
66	<i>Brachybacterium hainanense</i>	new species	Gram +	aerobic	ENVIRONMENT
67	<i>Brachybacterium huguangmaarensis</i>	new species	Gram +	aerobic	ENVIRONMENT
68	<i>Brachybacterium phenoliresistens</i>	new species	Gram +	aerobic	ENVIRONMENT
69	<i>Brachybacterium rhamnosum</i>	new species	Gram +	aerobic	ENVIRONMENT
70	<i>Brachybacterium sacelli</i>	new species	Gram +	aerobic	ENVIRONMENT
71	<i>Brachybacterium saurashtrensis</i>	new species	Gram +	aerobic	ENVIRONMENT
72	<i>Brachybacterium zhongshanense</i>	new species	Gram +	aerobic	ENVIRONMENT
73	<i>Bradyrhizobium</i> sp	new species	Gram -	aerobic	ENVIRONMENT
74	<i>Brevundimonas mediterranea</i>	new species	Gram -	aerobic	ENVIRONMENT
75	<i>Brevundimonas</i> sp[2]	new species	Gram -	aerobic	ENVIRONMENT



	New genus/species				Main relevance
76	<i>Buttiauxella</i> sp	new species	Gram -	microaerophilic	ENVIRONMENT
77	<i>Candida sphagnicola</i>	new species	Yeast	aerobic	ENVIRONMENT
78	<i>Chelatococcus</i> sp[3]	new genus/species	Gram -	aerobic	ENVIRONMENT
79	<i>Chitinophaga jiangningensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
80	<i>Clostridium polynesiense</i>	new species	Gram +	anaerobic	CLINICAL
81	<i>Clostridium senegalense</i>	new species	Gram +	anaerobic	CLINICAL
82	<i>Corynebacterium provencense</i>	new species	Gram +	aerobic	CLINICAL
83	<i>Cupriavidus plantarum</i>	new species	Gram -	aerobic	ENVIRONMENT
84	<i>Cupriavidus taiwanensis</i>	new species	Gram -	aerobic	ENVIRONMENT
85	<i>Debaryomyces maramus</i>	new species	Yeast	aerobic	ENVIRONMENT
86	<i>Dietzia cercidiphylli</i>	new species	Gram +	aerobic	ENVIRONMENT
87	<i>Dietzia lutea</i>	new species	Gram +	aerobic	ENVIRONMENT
88	<i>Eisenbergiella tayi</i>	new species	Gram -	anaerobic	CLINICAL
89	<i>Empedobacter stercoris</i>	new species	Gram -	aerobic	ENVIRONMENT
90	<i>Escherichia marmotae</i>	new species	Gram -	aerobic	ENVIRONMENT
91	<i>Falsiporphyrromonas endometrii</i>	new genus/species	Gram -	anaerobic	VETERINARY
92	<i>Flavobacterium cucumis</i>	new species	Gram -	aerobic	ENVIRONMENT
93	<i>Flavobacterium ginsengisoli</i>	new species	Gram -	aerobic	ENVIRONMENT
94	<i>Flavobacterium jumunjinense</i>	new species	Gram -	aerobic	ENVIRONMENT
95	<i>Fusobacterium simiae</i>	new species	Gram -	anaerobic	VETERINARY
96	<i>Glutamicibacter nicotianae</i>	new species	Gram +	aerobic	INDUSTRY
97	<i>Gordonia otitidis</i>	new species	Gram +	aerobic	CLINICAL
98	<i>Halobacillus alkaliphilus</i>	new species	Gram +	aerobic	ENVIRONMENT
99	<i>Halobacillus faecis</i>	new species	Gram +	aerobic	ENVIRONMENT
100	<i>Halobacillus kuroshimensis</i>	new species	Gram +	aerobic	ENVIRONMENT
101	<i>Halobacillus locisalis</i>	new species	Gram +	aerobic	ENVIRONMENT



	New genus/species				Main relevance
102	<i>Halobacillus mangrovi</i>	new species	Gram +	aerobic	ENVIRONMENT
103	<i>Halobacillus naozhouensi</i>	new species	Gram +	aerobic	ENVIRONMENT
104	<i>Halobacillus salinus</i>	new species	Gram +	aerobic	ENVIRONMENT
105	<i>Halobacillus seohaensis</i>	new species	Gram +	aerobic	ENVIRONMENT
106	<i>Halobacillus yeomjeoni</i>	new species	Gram +	aerobic	ENVIRONMENT
107	<i>Hanseniopsis osmophila</i>	new species	Yeast	aerobic	FOOD
108	<i>Ignatzschineria ureiclastica</i>	new species	Gram -	aerobic	VETERINARY
109	<i>Janthinobacterium agaricidamnosum</i>	new species	Gram -	aerobic	FOOD
110	<i>Jeotgalicoccus coquinae</i>	new species	Gram +	aerobic	ENVIRONMENT
111	<i>Jeotgalicoccus pinnipedialis</i>	new species	Gram +	aerobic	VETERINARY
112	<i>Jeotgalicoccus psychrophilus</i>	new species	Gram +	aerobic	FOOD
113	<i>Kushneria avicenniae</i>	new species	Gram -	aerobic	ENVIRONMENT
114	<i>Kushneria marisflavi</i>	new species	Gram -	aerobic	ENVIRONMENT
115	<i>Lacrimispora saccharolytica</i>	new genus/species	Gram +	anaerobic	ENVIRONMENT
116	<i>Legionella</i> sp	new species	Gram -	aerobic	ENVIRONMENT
117	<i>Leifsonia bigeumensis</i>	new species	Gram +	aerobic	ENVIRONMENT
118	<i>Leifsonia psychrotolerans</i>	new species	Gram +	aerobic	ENVIRONMENT
119	<i>Luteimonas abyssi</i>	new species	Gram -	aerobic	ENVIRONMENT
120	<i>Lysinibacillus massiliensis</i>	new species	Gram +	aerobic	CLINICAL
121	<i>Lysobacter hankyongensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
122	<i>Malaciobacter canalis</i>	new species	Gram -	microaerophilic	ENVIRONMENT
123	<i>Malaciobacter mytili</i>	new species	Gram -	microaerophilic	ENVIRONMENT
124	<i>Mannheimia caviae</i>	new species	Gram -	facultative anaerobic	VETERINARY
125	<i>Mannheimia ruminalis</i>	new species	Gram -	facultative anaerobic	VETERINARY
126	<i>Methylobacterium</i> sp	new species	Gram -	aerobic	ENVIRONMENT
127	<i>Microbacterium gubbeenense</i>	new species	Gram +	aerobic	FOOD



	New genus/species				Main relevance
128	<i>Microbacterium hydrothermale</i>	new species	Gram +	aerobic	ENVIRONMENT
129	<i>Microbacterium invictum</i>	new species	Gram +	aerobic	ENVIRONMENT
130	<i>Microbacterium</i> sp	new species	Gram +	aerobic	CLINICAL
131	<i>Micrococcus endophyticus</i>	new species	Gram +	aerobic	ENVIRONMENT
132	<i>Microvirga flocculans</i>	new species	Gram -	aerobic	ENVIRONMENT
133	<i>Microvirga lotononidis</i>	new species	Gram -	aerobic	ENVIRONMENT
134	<i>Microvirga makkahensis</i>	new species	Gram -	aerobic	ENVIRONMENT
135	<i>Microvirga subterranea</i>	new species	Gram -	aerobic	ENVIRONMENT
136	<i>Mitsuokella multacida</i>	new species	Gram -	anaerobic	CLINICAL
137	<i>Mixta intestinalis</i>	new species	Gram -	aerobic	CLINICAL
138	<i>Muricoccus roseus</i>	new genus/species	Gram +	aerobic	ENVIRONMENT
139	<i>Necropsobacter rosorum</i>	new genus/species	Gram -	facultative anaerobic	CLINICAL
140	<i>Neomicrococcus aestuarii</i>	new species	Gram -	aerobic	ENVIRONMENT
141	<i>Nocardia iowensis</i>	new species	Gram +	aerobic	ENVIRONMENT
142	<i>Nosocomiicoccus ampullae</i>	new species	Gram +	aerobic	ENVIRONMENT
143	<i>Noviherbaspirillum canariense</i>	new species	Gram -	aerobic	ENVIRONMENT
144	<i>Noviherbaspirillum psychrotolerans</i>	new species	Gram -	aerobic	ENVIRONMENT
145	<i>Oceanobacillus locisalsi</i>	new species	Gram +	aerobic	ENVIRONMENT
146	<i>Ornithinibacillus contaminans</i>	new genus/species	Gram +	aerobic	CLINICAL
147	<i>Paenalcaligenes suwonensis</i>	new species	Gram -	aerobic	ENVIRONMENT
148	<i>Paenibacillus abyssii</i>	new species	Gram +	aerobic	ENVIRONMENT
149	<i>Paenibacillus aceris</i>	new species	Gram +	aerobic	ENVIRONMENT
150	<i>Paenibacillus alkaliterrae</i>	new species	Gram +	aerobic	ENVIRONMENT
151	<i>Paenibacillus arcticus</i>	new species	Gram +	aerobic	ENVIRONMENT
152	<i>Paenibacillus azotifigens</i>	new species	Gram +	aerobic	ENVIRONMENT
153	<i>Paenibacillus castaneae</i>	new species	Gram +	aerobic	ENVIRONMENT



	New genus/species				Main relevance
154	<i>Paenibacillus cavernae</i>	new species	Gram +	aerobic	ENVIRONMENT
155	<i>Paenibacillus cellulosityticus</i>	new species	Gram +	aerobic	ENVIRONMENT
156	<i>Paenibacillus chungangensis</i>	new species	Gram +	aerobic	ENVIRONMENT
157	<i>Paenibacillus cisolokensis</i>	new species	Gram +	aerobic	ENVIRONMENT
158	<i>Paenibacillus dakarensis</i>	new species	Gram +	aerobic	CLINICAL
159	<i>Paenibacillus elymi</i>	new species	Gram +	aerobic	ENVIRONMENT
160	<i>Paenibacillus endophyticus</i>	new species	Gram +	aerobic	ENVIRONMENT
161	<i>Paenibacillus etheri</i>	new species	Gram +	aerobic	ENVIRONMENT
162	<i>Paenibacillus eucommiae</i>	new species	Gram +	aerobic	ENVIRONMENT
163	<i>Paenibacillus faecis</i>	new species	Gram +	aerobic	CLINICAL
164	<i>Paenibacillus fonticola</i>	new species	Gram +	aerobic	ENVIRONMENT
165	<i>Paenibacillus forsythiae</i>	new species	Gram +	aerobic	ENVIRONMENT
166	<i>Paenibacillus ginsengarvi</i>	new species	Gram +	aerobic	ENVIRONMENT
167	<i>Paenibacillus ginsengihumi</i>	new species	Gram +	aerobic	ENVIRONMENT
168	<i>Paenibacillus harenae</i>	new species	Gram +	aerobic	ENVIRONMENT
169	<i>Paenibacillus herberti</i>	new species	Gram +	aerobic	ENVIRONMENT
170	<i>Paenibacillus hodogayensis</i>	new species	Gram +	aerobic	ENVIRONMENT
171	<i>Paenibacillus jilunlii</i>	new species	Gram +	aerobic	ENVIRONMENT
172	<i>Paenibacillus konkukensis</i>	new species	Gram +	aerobic	ENVIRONMENT
173	<i>Paenibacillus konsidensis</i>	new species	Gram +	aerobic	CLINICAL
174	<i>Paenibacillus kribbensis</i>	new species	Gram +	aerobic	ENVIRONMENT
175	<i>Paenibacillus lupini</i>	new species	Gram +	aerobic	ENVIRONMENT
176	<i>Paenibacillus marchantiophytorum</i>	new species	Gram +	aerobic	ENVIRONMENT
177	<i>Paenibacillus nanensis</i>	new species	Gram +	aerobic	ENVIRONMENT
178	<i>Paenibacillus panacisoli</i>	new species	Gram +	aerobic	ENVIRONMENT
179	<i>Paenibacillus pectinilyticus</i>	new species	Gram +	aerobic	VETERINARY



	New genus/species				Main relevance
180	<i>Paenibacillus physcomitrellae</i>	new species	Gram +	aerobic	ENVIRONMENT
181	<i>Paenibacillus pinihumi</i>	new species	Gram +	aerobic	ENVIRONMENT
182	<i>Paenibacillus pocheonensis</i>	new species	Gram +	aerobic	ENVIRONMENT
183	<i>Paenibacillus popilliae</i>	new species	Gram +	aerobic	ENVIRONMENT
184	<i>Paenibacillus prosopidis</i>	new species	Gram +	aerobic	ENVIRONMENT
185	<i>Paenibacillus pueri</i>	new species	Gram +	aerobic	ENVIRONMENT / FOOD
186	<i>Paenibacillus qingshengii</i>	new species	Gram +	aerobic	ENVIRONMENT
187	<i>Paenibacillus ripae</i>	new species	Gram +	aerobic	ENVIRONMENT
188	<i>Paenibacillus sacheonensis</i>	new species	Gram +	aerobic	ENVIRONMENT
189	<i>Paenibacillus sediminis</i>	new species	Gram +	aerobic	ENVIRONMENT
190	<i>Paenibacillus segetis</i>	new species	Gram +	aerobic	ENVIRONMENT
191	<i>Paenibacillus sepulcri</i>	new species	Gram +	aerobic	ENVIRONMENT
192	<i>Paenibacillus shenyangensis</i>	new species	Gram +	aerobic	ENVIRONMENT
193	<i>Paenibacillus shirakamiensis</i>	new species	Gram +	aerobic	ENVIRONMENT
194	<i>Paenibacillus silagei</i>	new species	Gram +	aerobic	ENVIRONMENT
195	<i>Paenibacillus silvae</i>	new species	Gram +	aerobic	ENVIRONMENT
196	<i>Paenibacillus sophorae</i>	new species	Gram +	aerobic	ENVIRONMENT
197	<i>Paenibacillus sputi</i>	new species	Gram +	aerobic	CLINICAL
198	<i>Paenibacillus terreus</i>	new species	Gram +	aerobic	ENVIRONMENT
199	<i>Paenibacillus thailandensis</i>	new species	Gram +	aerobic	ENVIRONMENT
200	<i>Paenibacillus thermoaerophilus</i>	new species	Gram +	aerobic	ENVIRONMENT
201	<i>Paenibacillus tibetensis</i>	new species	Gram +	aerobic	ENVIRONMENT
202	<i>Paenibacillus tylopili</i>	new species	Gram +	aerobic	ENVIRONMENT
203	<i>Paenibacillus typhae</i>	new species	Gram +	aerobic	ENVIRONMENT
204	<i>Paenibacillus vulneris</i>	new species	Gram +	aerobic	CLINICAL
205	<i>Paenibacillus wexiniae</i>	new species	Gram +	aerobic	FOOD



	New genus/species				Main relevance
206	<i>Paenisporosarcina macmurdoensis</i>	new genus/species	Gram +	aerobic	ENVIRONMENT
207	<i>Pantoea vagans</i>	new species	Gram -	aerobic	CLINICAL / ENVIRONMENT
208	<i>Paracoccus aminovorans</i>	new species	Gram -	aerobic	ENVIRONMENT
209	<i>Paracoccus</i> sp	new species	Gram -	aerobic	ENVIRONMENT
210	<i>Pasteurella oralis</i>	new species	Gram -	aerobic	CLINICAL / VETERINARY
211	<i>Paucisalibacillus algeriensis</i>	new species	Gram +	aerobic	ENVIRONMENT
212	<i>Phaeobacter italicus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
213	<i>Photobacterium leiognathi</i>	new species	Gram -	aerobic	VETERINARY
214	<i>Phytobacter ursingii</i>	new genus/species	Gram -	facultative anaerobic	CLINICAL
215	<i>Planomicrobium chinense</i>	new species	Gram +	aerobic	ENVIRONMENT
216	<i>Planomicrobium glaciei</i>	new species	Gram +	aerobic	ENVIRONMENT
217	<i>Polynucleobacter rarus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
218	<i>Pontibacter actinarum</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
219	<i>Pontibacter akesuensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
220	<i>Pontibacter amylolyticus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
221	<i>Pontibacter korlensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
222	<i>Pontibacter mucosus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
223	<i>Pontibacter ummariensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
224	<i>Pontibacter virosus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
225	<i>Pseudarcobacter aquimarinus</i>	new genus/species	Gram -	microaerophilic	ENVIRONMENT
226	<i>Pseudarcobacter cloacae</i>	new genus/species	Gram -	microaerophilic	ENVIRONMENT
227	<i>Pseudarcobacter defluvii</i>	new genus/species	Gram -	microaerophilic	ENVIRONMENT
228	<i>Pseudoclavibacter caeni</i>	new species	Gram +	aerobic	ENVIRONMENT
229	<i>Pseudoclavibacter chungangensis</i>	new species	Gram +	aerobic	ENVIRONMENT
230	<i>Pseudoclavibacter endophyticus</i>	new species	Gram +	aerobic	ENVIRONMENT
231	<i>Pseudoclavibacter soli</i>	new species	Gram +	aerobic	ENVIRONMENT



	New genus/species				Main relevance
232	<i>Pseudomonas borbori</i>	new species	Gram -	aerobic	ENVIRONMENT
233	<i>Pseudomonas</i> sp	new species	Gram -	aerobic	CLINICAL
234	<i>Pyramidobacter piscolens</i>	new genus/species	Gram -	anaerobic	CLINICAL
235	<i>Ralstonia</i> sp	new species	Gram -	aerobic	ENVIRONMENT
236	<i>Ramlibacter</i> sp	new genus/species	Gram -	aerobic	ENVIRONMENT
237	<i>Rhizobium selenitireducens</i>	new species	Gram -	aerobic	ENVIRONMENT
238	<i>Rhodococcus defluvii</i>	new species	Gram +	aerobic	CLINICAL / ENVIRONMENT
239	<i>Roseomonas aceris</i>	new species	Gram -	aerobic	ENVIRONMENT
240	<i>Roseomonas aerofrigidensis</i>	new species	Gram -	aerobic	ENVIRONMENT
241	<i>Roseomonas aquatica</i>	new species	Gram -	aerobic	ENVIRONMENT / FOOD
242	<i>Roseomonas lacus</i>	new species	Gram -	aerobic	ENVIRONMENT
243	<i>Roseomonas ludipueritiae</i>	new species	Gram -	aerobic	ENVIRONMENT
244	<i>Roseomonas oryzicola</i>	new species	Gram -	aerobic	ENVIRONMENT
245	<i>Roseomonas vinacea</i>	new species	Gram -	aerobic	ENVIRONMENT
246	<i>Rummeliibacillus stabekisii</i>	new species	Gram +	aerobic	ENVIRONMENT
247	<i>Salinicoccus roseus</i>	new genus/species	Gram +	aerobic	ENVIRONMENT
248	<i>Salinivibrio proteolyticus</i>	new species	Gram -	aerobic	ENVIRONMENT
249	<i>Serratia nematodiphila</i>	new species	Gram -	aerobic	VETERINARY
250	<i>Shewanella xiamenensis</i>	new species	Gram -	aerobic	CLINICAL / ENVIRONMENT
251	<i>Sphingobacterium nematocida</i>	new species	Gram -	aerobic	ENVIRONMENT
252	<i>Sporolactobacillus</i> sp[2]	new species	Gram +	anaerobic	ENVIRONMENT
253	<i>Sporolactobacillus spathodeae</i>	new species	Gram +	anaerobic	ENVIRONMENT
254	<i>Sporosarcina aquimarina</i>	new species	Gram +	aerobic	ENVIRONMENT
255	<i>Sporosarcina pasteurii</i>	new species	Gram +	aerobic	ENVIRONMENT / INDUSTRY
256	<i>Sporosarcina</i> sp	new species	Gram +	aerobic	ENVIRONMENT
257	<i>Staphylococcus massiliensis</i>	new species	Gram +	aerobic	CLINICAL



	New genus/species				Main relevance
258	<i>Starmerella bombicola</i>	new species	Yeast	aerobic	ENVIRONMENT / INDUSTRY
259	<i>Stenotrophomonas koreensis</i>	new species	Gram -	aerobic	ENVIRONMENT
260	<i>Streptococcus penaecida</i>	new species	Gram +	facultative anaerobic	VETERINARY
261	<i>Streptomyces</i> sp[3]	new species	Gram +	aerobic	ENVIRONMENT
262	<i>Terribacillus halophilus</i>	new genus/species	Gram +	aerobic	ENVIRONMENT
263	<i>Tetragenococcus halophilus</i>	new species	Gram +	aerobic	FOOD
264	<i>Undibacterium macrobrachii</i>	new species	Gram -	aerobic	ENVIRONMENT
265	<i>Ureibacillus suwonensis</i>	new genus/species	Gram -	aerobic	ENVIRONMENT
266	<i>Ureibacillus thermosphaericus</i>	new genus/species	Gram -	aerobic	ENVIRONMENT / INDUSTRY
267	<i>Vagococcus fessus</i>	new species	Gram +	anaerobic	VETERINARY
268	<i>Varibaculum anthropi</i>	new species	Gram +	anaerobic	CLINICAL
269	<i>Verticiella sediminum</i>	new species	Gram -	anaerobic	ENVIRONMENT
270	<i>Vibrio europaeus</i>	new species	Gram -	aerobic	VETERINARY
271	<i>Williamsia deligens</i>	new species	Gram +	aerobic	CLINICAL
272	<i>Williamsia faeni</i>	new species	Gram +	aerobic	ENVIRONMENT
273	<i>Williamsia limnetica</i>	new species	Gram +	aerobic	ENVIRONMENT
274	<i>Williamsia maris</i>	new species	Gram +	aerobic	ENVIRONMENT
275	<i>Williamsia phyllosphaerae</i>	new species	Gram +	aerobic	ENVIRONMENT
276	<i>Williamsia serinedens</i>	new species	Gram +	aerobic	CLINICAL / ENVIRONMENT
277	<i>Xylella fastidiosa</i>	new genus/species	Gram -	aerobic	ENVIRONMENT / FOOD
278	<i>Yarrowia galli</i>	new species	Yeast	aerobic	CLINICAL / VETERINARY



3 Deleted MSP entries

Table 4: Deleted MSP entries

Deletions	Reason
Acinetobacter pittii DSM 9320 DSM	Doubtful reference ID
Acinetobacter sp DSM 30009 DSM	
Actinomyces hongkongensis ENR_0065 ENR	Wrong reference ID
Actinomyces hongkongensis ENR_0105 ENR	
Arthrobacter crystallopoietes DSM 20117T DSM	Low spectra quality (polymeric peaks of poly-lysin)
Bacillus sp 57355 RQCL	Unclear and non-solvable taxonomical situation
Dietzia cinnamomea 117 RLT	Contamination of the reference spectra
Dysgonomonas gadei F_8814_1 IMK	Wrong reference ID
Glutamicibacter nicotianae IMET 10353T HKJ	Low spectra quality
Gordonia bronchialis DSM 43247T DSM	
Lactobacillus fructivorans DSM 20353 DSM	Doubtful reference ID
Lactobacillus homohiochii DSM 20354 DSM	
Mycobacterium avium TB_RV422_4_02 UKE	Low spectra quality (polymeric peaks of poly-lysin)
Nocardia sp MB_9090_05 THL	Doubtful reference ID



Deletions	Reason
Paenibacillus lactis CICC 24043 CICC	Doubtful reference ID
Paenibacillus sp 09 CTC	
Pseudarthrobacter polychromogenes DSM 20136T DSM	
Pseudarthrobacter sulfonivorans DSM 14002T DSM	Low spectra quality (polymeric peaks of poly-lysin)
Shewanella putrefaciens CCM 2803 CCM	Wrong reference ID
Sphingomonas aerolata DSM 14746T DSM	Low spectra quality (polymeric peaks of poly-lysin)
Tsukamurella paurometabola DSM 46065 DSM	
Veillonella sp[3] 0807M16032801 IBS	Doubtful reference ID



4 Renaming

Table 5: Renaming of MSP entries

DB-8468	DB-9607	Justification
Acidovorax avenae ssp avenae DSM 7227T HAM	Acidovorax avenae DSM 7227T HAM	According to information from strain collection.
<i>Acidovorax avenae</i> ssp <i>citrulli</i>	<i>Acidovorax citrulli</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Acinetobacter dijkshoorniae</i>	<i>Acinetobacter lactucae</i>	
Actinobacillus equuli ssp equuli DSM 19655T DSM	Actinobacillus equuli DSM 19655T DSM	According to information from strain collection.
<i>Actinomyces canis</i>	<i>Schaalia canis</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Actinomyces cardiffensis</i>	<i>Schaalia cardiffensis</i>	
<i>Actinomyces coleocanis</i>	<i>Gleimia coleocanis</i>	
<i>Actinomyces hominis</i>	<i>Gleimia hominis</i>	
<i>Actinomyces hordeovulneris</i>	<i>Buchananella hordeovulneris</i>	
<i>Actinomyces marimammalium</i>	<i>Boudabousia marimammalium</i>	
<i>Actinomyces nasicola</i>	<i>Bowdeniella nasicola</i>	
Actinomyces sp VA_01434_2_09 ERL	Schaalia sp VA_01434_2_09 ERL	
<i>Actinomyces suimastitidis</i>	<i>Schaalia suimastitidis</i>	
<i>Actinomyces vaccimaxillae</i>	<i>Schaalia vaccimaxillae</i>	
Alcaligenes sp 091029_c SLT	Paenalcaligenes suwonensis 091029_c SLT	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.
Anoxybacillus flavithermus ssp flavithermus DSM 21510 DSM	Anoxybacillus flavithermus DSM 21510 DSM	According to information from strain collection.
Anoxybacillus flavithermus ssp flavithermus DSM 2641T DSM	Anoxybacillus flavithermus DSM 2641T DSM	



DB-8468	DB-9607	Justification
<i>Arcobacter butzleri</i>	<i>Aliarcobacter butzleri</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Arcobacter halophilus</i>	<i>Malaciobacter halophilus</i>	
<i>Arcobacter skirrowii</i>	<i>Aliarcobacter skirrowii</i>	
<i>Arthroascus schoenii</i>	<i>Saccharomycopsis schoenii</i>	
<i>Arthrobacter nasiphocae</i>	<i>Falsarthrobacter nasiphocae</i>	
<i>Atopobium parvulum</i>	<i>Lancefieldella parvula</i>	
<i>Atopobium rimae</i>	<i>Lancefieldella rimae</i>	
<i>Atopobium vaginae</i>	<i>Fannyhessea vaginae</i>	
<i>Bacillus weihenstephanensis</i>	<i>Bacillus mycoides</i>	
Bifidobacterium animalis ssp animalis DSM 20104T DSM	Bifidobacterium animalis DSM 20104T DSM	According to information from strain collection.
<i>Bifidobacterium coryneforme</i>	<i>Bifidobacterium indicum</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Bifidobacterium gallinarum</i>	<i>Bifidobacterium pullorum</i> ssp <i>gallinarum</i>	
Bifidobacterium longum ssp longum DSM 20090 DSM	Bifidobacterium longum ssp infantis DSM 20090 DSM	According to information from strain collection.
Bifidobacterium longum ssp longum DSM 20097 DSM	Bifidobacterium longum ssp suis DSM 20097 DSM	
Bifidobacterium longum ssp longum DSM 20218 DSM	Bifidobacterium longum ssp infantis DSM 20218 DSM	
Bifidobacterium pseudolongum ssp pseudolongum DSM 20094 DSM	Bifidobacterium pseudolongum DSM 20094 DSM	
Bifidobacterium pseudolongum ssp pseudolongum DSM 20099T DSM	Bifidobacterium pseudolongum DSM 20099T DSM	
<i>Bifidobacterium saeculare</i>	<i>Bifidobacterium pullorum</i> ssp <i>saeculare</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.



DB-8468	DB-9607	Justification
<i>Bifidobacterium thermacidophilum</i> ssp <i>porcinum</i>	<i>Bifidobacterium porcinum</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
Bifidobacterium thermacidophilum ssp thermacidophilum DSM 15837T DSM	Bifidobacterium thermacidophilum DSM 15837T DSM	According to information from strain collection.
Campylobacter hyointestinalis CCUG 14169T NVU	Campylobacter hyointestinalis ssp hyointestinalis CCUG 14169T NVU	
<i>Candida carpophila</i>	<i>Meyerozyma carpophila</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Candida catenulata</i>	<i>Diutina catenulata</i>	
<i>Candida davenportii</i>	<i>Starmerella davenportii</i>	
<i>Candida ernobii</i>	<i>Nakazawaea ernobii</i>	
<i>Candida famata</i>	<i>Debaryomyces hansenii</i>	
<i>Candida guilliermondii</i>	<i>Meyerozyma guilliermondii</i>	
<i>Candida inconspicua</i>	<i>Pichia cactophila</i>	
<i>Candida infanticola</i>	<i>Wickerhamiella infanticola</i>	
<i>Candida kefyi</i>	<i>Kluyveromyces marxianus</i>	
<i>Candida krusei</i>	<i>Pichia kudriavzevii</i>	
<i>Candida lactiscondensi</i>	<i>Starmerella lactis-condensi</i>	
<i>Candida lambica</i>	<i>Pichia fermentans</i>	
<i>Candida lusitanae</i>	<i>Clavispora lusitanae</i>	
<i>Candida magnoliae</i>	<i>Starmerella magnoliae</i>	
<i>Candida pararugosa</i>	<i>Wickerhamiella pararugosa</i>	



DB-8468	DB-9607	Justification
<i>Candida peltata</i>	<i>Nakazawaea peltata</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Candida pini</i>	<i>Ogataea pini</i>	
<i>Candida rugosa</i>	<i>Diutina rugosa</i>	
<i>Candida shehatae var insectosa</i>	<i>Scheffersomyces shehatae var insectosa</i>	
<i>Candida spandovensis</i>	<i>Wickerhamiella spandovensis</i>	
<i>Candida tenuis</i>	<i>Yamadazyma tenuis</i>	
<i>Candida utilis</i>	<i>Cyberlindnera jadinii</i>	
<i>Candida valida</i>	<i>Pichia membranifaciens</i>	
<i>Candida versatilis</i>	<i>Wickerhamiella versatilis</i>	
<i>Candida vini</i>	<i>Kregervanrija fluxuum</i>	
<i>Clostridium aldenense</i>	<i>Enterocloster aldensis</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Clostridium difficile</i>	<i>Clostridioides difficile</i>	
Corynebacterium sp 901400365 LBK	Corynebacterium phoceense 901400365 LBK	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.
Corynebacterium sp 901600604 LBK	Corynebacterium phoceense 901600604 LBK	
<i>Cryptococcus albidosimilis</i>	<i>Naganishia albidosimilis</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Cryptococcus albidus</i>	<i>Naganishia albida</i>	
<i>Cryptococcus curvatus</i>	<i>Cutaneotrichosporon curvatum</i>	
<i>Cryptococcus diffluens</i>	<i>Naganishia diffluens</i>	
<i>Cryptococcus flavescens</i>	<i>Papiliotrema flavescens</i>	
<i>Cryptococcus flavus</i>	<i>Saitozyma flava</i>	



DB-8468	DB-9607	Justification
<i>Cryptococcus gastricus</i>	<i>Goffeauzyma gastrica</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Cryptococcus humicola</i>	<i>Vanrija humicola</i>	
<i>Cryptococcus laurentii</i>	<i>Papiliotrema laurentii</i>	
<i>Cryptococcus liquefaciens</i>	<i>Naganishia liquefaciens</i>	
<i>Cryptococcus macerans</i>	<i>Cystofilobasidium macerans</i>	
<i>Cryptococcus magnus</i>	<i>Filobasidium magnum</i>	
<i>Cryptococcus saitoi</i>	<i>Naganishia globosa</i>	
<i>Cryptococcus terreus</i>	<i>Solicoccozyma terrea</i>	
<i>Cryptococcus uzbekistanensis</i>	<i>Naganishia uzbekistanensis</i>	
<i>Cryptococcus vishniacii</i>	<i>Naganishia vishniacii</i>	
<i>Debaryomyces etchellsii</i>	<i>Schwanniomyces etchellsii</i>	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.
Dietzia sp[2] 72094 RQCL	Dietzia papillomatosis 72094 RQCL	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Dipodascus ingens</i>	<i>Magnusiomyces ingens</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Enterobacter xiangfangensis</i>	<i>Enterobacter hormaechei</i> ssp <i>xiangfangensis</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Filobasidium capsuligenum</i>	<i>Piskurozyma capsuligena</i>	



DB-8468	DB-9607	Justification
Flavobacterium aquatile DSM 1132T DSM	Flavobacterium hercynium DSM 18292T DSM	Correction of reference names due to sample mix-up during measurement.
Flavobacterium cutihirudinis DSM 25795T DSM	Flavobacterium frigoris DSM 15719T DSM	
Flavobacterium daejeonense DSM 17708T DSM	Flavobacterium denitrificans DSM 15936T DSM	
Flavobacterium defluvii DSM 17963T DSM	Flavobacterium degerlachei DSM 15718T DSM	
Flavobacterium degerlachei DSM 15718T DSM	Flavobacterium defluvii DSM 17963T DSM	
Flavobacterium flevense DSM 1076T DSM_2	Flavobacterium granuli DSM 19729T DSM	
Flavobacterium frigidarium DSM 17623T DSM	Flavobacterium glaciei DSM 19728T DSM	
Flavobacterium frigoris DSM 15719T DSM	Flavobacterium cutihirudinis DSM 25795T DSM	
Flavobacterium fryxellicola DSM 16209T DSM	Flavobacterium chungbukense DSM 25688 DSM	
Flavobacterium glaciei DSM 19728T DSM	Flavobacterium frigidarium DSM 17623T DSM	
Flavobacterium granuli DSM 19729T DSM	Flavobacterium flevense DSM 1076T DSM_2	
Flavobacterium hercynium DSM 18292T DSM	Flavobacterium aquatile DSM 1132T DSM	
<i>Fluoribacter bozemanae</i>	<i>Legionella bozemanae</i>	Consideration of the more common species name.
<i>Fluoribacter dumoffii</i>	<i>Legionella dumoffii</i>	
<i>Fluoribacter gormanii</i>	<i>Legionella gormanii</i>	
<i>Geotrichum capitatum</i>	<i>Magnusiomyces capitatus</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Geotrichum silvicola</i>	<i>Geotrichum candidum</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.



DB-8468	DB-9607	Justification
<i>Guehomyces pullulans</i>	<i>Tausonia pullulans</i>	Yeast nomenclature was substantially updated according to recent findings. Please consider section 8.3 for detailed explanation.
<i>Issatchenkia terricola</i>	<i>Pichia terricola</i>	
<i>Kocuria koreensis</i>	<i>Rothia koreensis</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
Lachancea cidri CBS 4575 CBS	Lachancea cidri CBS 4575T CBS	Cosmetic change – information “T” of an available Type Strain has been added to strain number within the MSP name.
Lactobacillus aviarius ssp aviarius DSM 20654 DSM	Lactobacillus aviarius DSM 20654 DSM	According to information from strain collection.
Lysinibacillus massiliensis 1112072947_1e MVD	Lysinibacillus halotolerans 1112072947_1e MVD	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.
Millerozyma farinosa DSM 2226T DSM	Millerozyma farinosa DSM 2226 DSM	Cosmetic change – information “T” of an available Type Strain has been removed from strain number within the MSP name.
<i>Mycoplasma bovirhinis</i>	<i>Mycoplasmaopsis bovirhinis</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Mycoplasma bovis</i>	<i>Mycoplasmaopsis bovis</i>	
<i>Mycoplasma gallisepticum</i>	<i>Mycoplasmaoides gallisepticum</i>	
<i>Mycoplasma hominis</i>	<i>Metamycoplasma hominis</i>	
<i>Mycoplasma hyorhinis</i>	<i>Mesomycoplasma hyorhinis</i>	
Peptoniphilus harei 06_622 IBS	Peptoniphilus sp[2] 06_622 IBS	More precise representation of the taxonomical situation.
Peptoniphilus harei 08_570 IBS	Peptoniphilus sp[2] 08_570 IBS	
Peptoniphilus indolicus ENR_0001 ENR	Peptoniphilus sp[2] ENR_0001 ENR	
Peptoniphilus indolicus ENR_0003 ENR	Peptoniphilus sp[2] ENR_0003 ENR	
Peptoniphilus indolicus ENR_0423 ENR	Peptoniphilus sp[2] ENR_0423 ENR	



DB-8468	DB-9607	Justification
Peptoniphilus indolicus ENR_0430 ENR	Peptoniphilus sp[2] ENR_0430 ENR	More precise representation of the taxonomical situation.
<i>Peptoniphilus rhinitidis</i>	<i>Peptoniphilus lacydonensis</i>	Consideration of the validly published species name.
Peptoniphilus sp 110706_K3 LUMC	Peptoniphilus coxii 110706_K3 LUMC	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.
<i>Pichia methylvora</i>	<i>Ogataea methylvora</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Pseudozyma aphidis</i>	<i>Moesziomyces aphidis</i>	
<i>Rhodococcus kunmingensis</i>	<i>Aldersonia kunmingensis</i>	Reclassification of species nomenclature. The former name is still available within the Metadata comment field.
<i>Rhodosporidium</i> sp	<i>Rhodotorula</i> sp	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Rhodosporidium toruloides</i>	<i>Rhodotorula toruloides</i>	
<i>Rhodotorula acheniorum</i>	<i>Farysia acheniorum</i>	
<i>Rhodotorula bacarum</i>	<i>Microstroma album</i>	
<i>Rhodotorula bogoriensis</i>	<i>Pseudohyphozyma bogoriensis</i>	
<i>Rhodotorula minuta</i>	<i>Cystobasidium minutum</i>	
<i>Rhodotorula pustula</i>	<i>Pseudohyphozyma pustula</i>	
<i>Saprochaete clavata</i>	<i>Magnusiomyces clavatus</i>	
<i>Sporidiobolus salmonicolor</i>	<i>Sporobolomyces salmonicolor</i>	
Staphylococcus auricularis DSM 20609 DSM	Staphylococcus auricularis DSM 20609T DSM	



DB-8468	DB-9607	Justification
Staphylococcus equorum ssp equorum DSM 20674T DSM	Staphylococcus equorum DSM 20674T DSM	According to information from strain collection.
Staphylococcus equorum ssp equorum DSM 20675 DSM	Staphylococcus equorum DSM 20675 DSM	
<i>Tatlockia maceachernii</i>	<i>Legionella maceachernii</i>	Consideration of the more common species name.
<i>Tatlockia micdadei</i>	<i>Legionella micdadei</i>	
<i>Trichosporon cutaneum</i>	<i>Cutaneotrichosporon cutaneum</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Trichosporon debeurmannianum</i>	<i>Cutaneotrichosporon debeurmannianum</i>	
<i>Trichosporon dulciturum</i>	<i>Apiotrichum dulciturum</i>	
<i>Trichosporon gracile</i>	<i>Apiotrichum gracile</i>	
<i>Trichosporon jirovecii</i>	<i>Cutaneotrichosporon jirovecii</i>	
<i>Trichosporon laibachii</i>	<i>Apiotrichum laibachii</i>	
<i>Trichosporon loubieri</i>	<i>Apiotrichum loubieri</i>	
<i>Trichosporon moniliiforme</i>	<i>Cutaneotrichosporon moniliiforme</i>	
<i>Trichosporon montevideense</i>	<i>Apiotrichum montevideense</i>	
<i>Trichosporon mucooides</i>	<i>Cutaneotrichosporon mucooides</i>	
<i>Trichosporon mycotoxinivorans</i>	<i>Apiotrichum mycotoxinivorans</i>	
<i>Trichosporon terricola</i>	<i>Cutaneotrichosporon terricola</i>	
<i>Wautersiella falsenii</i>	<i>Empedobacter falsenii</i>	
Wautersiella falsenii 1942_2016 IMHM	Empedobacter tilapiae 1942_2016 IMHM	General library maintenance resulted in improved knowledge about strain ID – renaming to correct species now.



DB-8468	DB-9607	Justification
Xanthomonas campestris DSM 3586T DSM	Xanthomonas campestris DSM 3586T DSM_2	Cosmetic change
Xanthomonas campestris pvar campestris DSM 3586T DSM	Xanthomonas campestris DSM 3586T DSM	According to information from strain collection.
Yarrowia lipolytica DSM 70561 DSM	Yarrowia deformans DSM 70561 DSM	
<i>Zygosaccharomyces florentinus</i>	<i>Zygotorulaspota florentina</i>	Yeast nomenclature was substantially updated according recent findings. Please consider section 8.3 for detailed explanation.
<i>Zygosaccharomyces microellipsoides</i>	<i>Torulaspota microellipsoides</i>	



5 Changes to Matching Hints

Table 6: Matching Hint Changes

Species	Action	DB-9607	Justification
<i>Acidovorax avenae</i> <i>Acidovorax citrulli</i>	Implement new Matching Hint	For the species <i>avenae</i> / <i>cattleyae</i> / <i>citrulli</i> of the genus <i>Acidovorax</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.
<i>Acidovorax delafieldii</i>	Implement new Matching Hint	For the species <i>delafieldii</i> / <i>kalamii</i> of the genus <i>Acidovorax</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Statement regarding separation with MALDI is not possible so far. Therefore, all these species should be taken into consideration as a possible result.
<i>Acinetobacter baumannii</i>	Update wording of Matching Hint	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter baylyi</i>	Update wording of Matching Hint	Species <i>baylyi</i> / <i>soli</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Acinetobacter soli</i> implementation of a more precise Matching Hint is needed.
<i>Acinetobacter bereziniae</i>	Implement new Matching Hint	Species <i>berezinae</i> / <i>guillouiae</i> / <i>wuhouensis</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Acinetobacter bereziniae</i> and <i>Acinetobacter guillouiae</i> manifests their close relationship to each other and to the newly implemented species <i>Acinetobacter wuhouensis</i> . Therefore, a Matching Hint needs to be implemented.



Species	Action	DB-9607	Justification
<i>Acinetobacter bouvetii</i>	Update wording of Matching Hint	Species bouvetii / pragensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Acinetobacter pragensis</i> implementation of a more precise Matching Hint is needed.
<i>Acinetobacter calcoaceticus</i>	Update wording of Matching Hint	Member of the Acinetobacter baumannii complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter courvalinii</i> <i>Acinetobacter dispersus</i>	Implement new Matching Hint	Species colistiniresistens / courvalinii / dispersus / gyllenbergii / proteolyticus / vivianii of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Acinetobacter courvalinii</i> , <i>Acinetobacter dispersus</i> , <i>Acinetobacter proteolyticus</i> and <i>Acinetobacter vivianii</i> manifests their close relationship to each other and to the newly implemented species <i>Acinetobacter colistiniresistens</i> and <i>Acinetobacter gyllenbergii</i> . Therefore, a Matching Hint needs to be implemented.
<i>Acinetobacter gerneri</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.
<i>Acinetobacter guillouiae</i>	Implement new Matching Hint	Species bereziniae / guillouiae / wuhouensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Acinetobacter bereziniae</i> and <i>Acinetobacter guillouiae</i> manifests their close relationship to each other and to the newly implemented species <i>Acinetobacter wuhouensis</i> . Therefore, a Matching Hint needs to be implemented.
<i>Acinetobacter haemolyticus</i> <i>Acinetobacter johnsonii</i> <i>Acinetobacter junii</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.



Species	Action	DB-9607	Justification
<i>Acinetobacter lactucae</i>	Update wording of Matching Hint	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter lwoffii</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.
<i>Acinetobacter nosocomialis</i>	Update wording of Matching Hint	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter parvus</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.
<i>Acinetobacter pittii</i>	Update wording of Matching Hint	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter proteolyticus</i>	Implement new Matching Hint	Species <i>colistiniresistens</i> / <i>courvalinii</i> / <i>dispersus</i> / <i>gyllenbergii</i> / <i>proteolyticus</i> / <i>vivianii</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Acinetobacter courvalinii</i> , <i>Acinetobacter disperses</i> , <i>Acinetobacter proteolyticus</i> and <i>Acinetobacter vivianii</i> manifests their close relationship to each other and to the newly implemented species <i>Acinetobacter colistiniresistens</i> and <i>Acinetobacter gyllenbergii</i> . Therefore, a Matching Hint needs to be implemented.
<i>Acinetobacter puyangensis</i>	Implement new Matching Hint	Species <i>populi</i> / <i>puyangensis</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Acinetobacter populi</i> implementation of a Matching Hint is needed.
<i>Acinetobacter radioresistens</i> <i>Acinetobacter schindleri</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.



Species	Action	DB-9607	Justification
<i>Acinetobacter seifertii</i>	Update wording of Matching Hint	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.	Cosmetical adaptation.
<i>Acinetobacter</i> sp <i>Acinetobacter tandoii</i> <i>Acinetobacter tjernbergiae</i> <i>Acinetobacter townneri</i> <i>Acinetobacter ursingii</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Acinetobacter</i> species in the library the very generally formulated Matching Hint is not needed anymore.
<i>Acinetobacter vivianii</i>	Implement new Matching Hint	Species <i>colistiniresistens</i> / <i>courvalinii</i> / <i>dispersus</i> / <i>gyllenbergii</i> / <i>proteolyticus</i> / <i>vivianii</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Acinetobacter courvalinii</i> , <i>Acinetobacter disperses</i> , <i>Acinetobacter proteolyticus</i> and <i>Acinetobacter vivianii</i> manifests their close relationship to each other and to the newly implemented species <i>Acinetobacter colistiniresistens</i> and <i>Acinetobacter gyllenbergii</i> . Therefore, a Matching Hint needs to be implemented.
<i>Actinobacillus equuli</i> <i>Actinobacillus lignieresii</i> <i>Actinobacillus pleuropneumoniae</i> <i>Actinobacillus suis</i>	Update wording of Matching Hint	For the species <i>arthritidis</i> / <i>capsulatus</i> / <i>equuli</i> / <i>lignieresii</i> / <i>pleuropneumoniae</i> / <i>suis</i> of the genus <i>Actinobacillus</i> the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the newly implemented species <i>Actinobacillus arthritidis</i> and <i>Actinobacillus capsulatus</i> the existing Matching Hint needs to be extended. Recent results have shown that these species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.



Species	Action	DB-9607	Justification
<i>Actinomyces denticolens</i>	Implement new Matching Hint	For the species denticolens / timonensis of the genus Actinomyces the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the newly implemented species <i>Actinomyces timonensis</i> recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing and by MALDI. Therefore, all these species should be taken into consideration as a possible result.
<i>Actinomyces georgiae</i>	Link to Matching Hint deleted	N/A	Due to deletion of the wrong reference entries for <i>Actinomyces hongkongensis</i> the Matching Hint is not needed anymore.
<i>Bacillus algicola</i>	Implement new Matching Hint	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.	Achievement of consistency with all other reference entries for <i>Bacillus</i> species.
<i>Bacillus cereus</i> <i>Bacillus mycoides</i>	Update wording of Matching Hint	<i>Bacillus anthracis</i> , <i>cereus</i> , <i>mycoides</i> , <i>pseudomycoides</i> and <i>thuringiensis</i> are closely related and members of the <i>Bacillus cereus</i> group. In particular <i>Bacillus cereus</i> spectra are very similar to spectra from <i>Bacillus anthracis</i> . <i>Bacillus anthracis</i> is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.	Due to renaming of <i>Bacillus weihenstephanensis</i> to <i>Bacillus mycoides</i> an adaptation of the Matching Hint is needed.
<i>Bacillus oceanisediminis</i> DSM 24771T DSM	Implement new Matching Hint	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.	Achievement of consistency with all other reference entries for <i>Bacillus</i> species.



Species	Action	DB-9607	Justification
<i>Bacillus pseudomycooides</i> <i>Bacillus thuringiensis</i>	Update wording of Matching Hint	Bacillus anthracis, cereus, mycooides, pseudomycooides and thuringiensis are closely related and members of the Bacillus cereus group. In particular Bacillus cereus spectra are very similar to spectra from Bacillus anthracis. Bacillus anthracis is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.	Due to renaming of <i>Bacillus weihenstephanensis</i> to <i>Bacillus mycooides</i> an adaptation of the Matching Hint is needed.
<i>Bifidobacterium porcinum</i> <i>Bifidobacterium thermacidophilum</i> <i>Bifidobacterium thermophilum</i>	Update wording of Matching Hint	Species porcinum / thermacidophilum / thermophilum of the genus Bifidobacterium have very similar patterns: Therefore distinguishing their species is difficult.	Due to renaming of <i>Bifidobacterium thermacidophilum</i> ssp <i>porcinum</i> to <i>Bifidobacterium porcinum</i> the Matching Hint needs to be extended.
<i>Brochothrix campestris</i> <i>Brochothrix thermosphacta</i>	Update wording of Matching Hint	For the species campestris / thermosphacta of the genus Brochothrix the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the improved number of reference entries for <i>Brochothrix thermosphacta</i> the existing Matching Hint needs to be adjusted. Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.
<i>Citrobacter amalonaticus</i>	Update wording of Matching Hint	Is a member of the Citrobacter amalonaticus complex. Species amalonaticus / farmeri of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.	Results from clinical studies have shown that the Matching Hint could be formulated more precisely.
<i>Citrobacter braakii</i>	Update wording of Matching Hint	Is a member of the Citrobacter freundii complex. Species braakii / freundii / gillenii / murliniae / rodentium / sedlakii / werkmannii / youngae of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.	Results from clinical studies have shown that the Matching Hint could be formulated more precisely.



Species	Action	DB-9607	Justification
<i>Citrobacter farmeri</i>	Update wording of Matching Hint	Is a member of the <i>Citrobacter amalonaticus</i> complex. Species <i>amalonaticus</i> / <i>farmeri</i> of the genus <i>Citrobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	Results from clinical studies have shown that the Matching Hint could be formulated more precisely.
<i>Citrobacter freundii</i> <i>Citrobacter gillenii</i>	Update wording of Matching Hint	Is a member of the <i>Citrobacter freundii</i> complex. Species <i>braakii</i> / <i>freundii</i> / <i>gillenii</i> / <i>murlinia</i> / <i>rodentium</i> / <i>sedlakii</i> / <i>werkmannii</i> / <i>youngae</i> of the genus <i>Citrobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	Results from clinical studies have shown that the Matching Hint could be formulated more precisely.
<i>Citrobacter koseri</i>	Link to Matching Hint deleted	N/A	Results from clinical studies have shown that the Matching Hint is not needed anymore.
<i>Citrobacter murlinia</i> <i>Citrobacter rodentium</i> <i>Citrobacter sedlakii</i> <i>Citrobacter youngae</i>	Update wording of Matching Hint	Is a member of the <i>Citrobacter freundii</i> complex. Species <i>braakii</i> / <i>freundii</i> / <i>gillenii</i> / <i>murlinia</i> / <i>rodentium</i> / <i>sedlakii</i> / <i>werkmannii</i> / <i>youngae</i> of the genus <i>Citrobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.	Results from clinical studies have shown that the Matching Hint could be formulated more precisely.
<i>Clostridioides difficile</i>	Implement new Matching Hint	synonym of <i>Clostridium difficile</i>	Displays the former, probably better-known name of the species.
<i>Corynebacterium phoceense</i>	Link to Matching Hint deleted	N/A	Due to the improved number of <i>Corynebacterium</i> species in the library the very generally formulated Matching Hint is not needed anymore.
<i>Dietzia natronolimnaea</i>	Implement new Matching Hint	For the species <i>cercidiphylli</i> / <i>natronolimnaea</i> / <i>psychralcaliphila</i> of the genus <i>Dietzia</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the newly implemented species <i>Dietzia cercidiphylli</i> recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Statement regarding separation with MALDI is not possible so far. Therefore, all these species should be taken into consideration as a possible result.



Species	Action	DB-9607	Justification
<i>Dietzia papillomatosis</i>	Implement new Matching Hint	For the species cinnamea / papillomatosis of the genus Dietzia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the renaming of a reference entry and hereby the implementation of the new species <i>Dietzia papillomatosis</i> recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.
<i>Dysgonomonas capnocytophagoides</i>	Link to Matching Hint deleted	N/A	Due to the improved number of reference entries for <i>Dysgonomonas capnocytophagoides</i> and the replacement of the wrong reference entry for the species <i>Dysgonomonas gadei</i> a clear discrimination of these two species is possible now and Therefore, the Matching Hint not needed anymore.
<i>Elizabethkingia anophelis</i> <i>Elizabethkingia meningoseptica</i> <i>Elizabethkingia miricola</i>	Update wording of Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Due to the improved number of reference entries for the displayed species a more general wording of the Matching Hint is possible now.
<i>Geotrichum candidum</i>	Link to Matching Hint deleted	N/A	Due to the renaming of <i>Geotrichum silvicola</i> to <i>Geotrichum candidum</i> the Matching Hint is not needed anymore.
<i>Globicatella sanguinis</i> <i>Globicatella sulfidifaciens</i>	Implement new Matching Hint	For the species sanguinis / sulfidifaciens of the genus Globicatella the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the improved number of reference entries for the displayed species recent results have shown that these species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.



Species	Action	DB-9607	Justification
<i>Kluyvera ascorbate</i> <i>Kluyvera georgiana</i>	Implement new Matching Hint	Species ascorbata / georgiana of the genus Kluyvera have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for the displayed species manifests their close relationship. Therefore, a Matching Hint needs to be implemented.
<i>Lactobacillus acidophilus</i>	Link to Matching Hint deleted	N/A	Due to the improved number of reference entries for <i>Lactobacillus acidophilus</i> a clear separation is possible now and Therefore, the Matching Hint is not needed anymore.
<i>Lactobacillus amylovorus</i>	Update wording of Matching Hint	Species amylovorus / kitasatonis of the genus Lactobacillus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to an exclusion of the species <i>Lactobacillus acidophilus</i> from the Matching Hint.
<i>Lactobacillus crispatus</i>	Link to Matching Hint deleted	N/A	Due to the improved number of reference entries for <i>Lactobacillus crispatus</i> a clear separation is possible now and Therefore, the Matching Hint is not needed anymore.
<i>Lactobacillus gasseri</i>	Implement new Matching Hint	For the species gasseri / taiwanensis of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Statement regarding separation with MALDI is not possible so far. Therefore, all these species should be taken into consideration as a possible result.
<i>Lactobacillus kitasatonis</i>	Update wording of Matching Hint	Species amylovorus / kitasatonis of the genus Lactobacillus have very similar patterns: Therefore distinguishing their species is difficult.	Due to the improved number of reference entries for <i>Lactobacillus acidophilus</i> a clear separation is possible now and Therefore, the Matching Hint needs to be adjusted.



Species	Action	DB-9607	Justification
<i>Lactobacillus ultunensis</i>	Link to Matching Hint deleted	N/A	Due to the improved number of reference entries for <i>Lactobacillus crispatus</i> a clear separation is possible now and Therefore, the Matching Hint is not needed anymore.
<i>Lactococcus garvieae</i>	Implement new Matching Hint	For the species garvieae / petauri of the genus Lactococcus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Statement regarding separation with MALDI is not possible so far. Therefore, all these species should be taken into consideration as a possible result.
<i>Mannheimia glucosida</i>	Update wording of Matching Hint	For the species glucosida / ruminalis of the genus Mannheimia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the newly implemented species <i>Mannheimia ruminalis</i> the existing Matching Hint needs to be adjusted. Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing Therefore, all these species should be taken into consideration as a possible result.
<i>Mixta calida</i> <i>Mixta gaviniae</i>	Update wording of Matching Hint	Species calida / gaviniae / intestinalis of the genus Mixta have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Mixta intestinalis</i> an extension of the existing Matching Hint is needed.
<i>Moorella thermoacetica</i>	Implement new Matching Hint	For the species thermoacetica / thermoautotrophica of the genus Moorella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the improved number of reference entries for <i>Moorella thermoacetica</i> recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.



Species	Action	DB-9607	Justification
<i>Nocardia abscessus</i> <i>Nocardia arthritidis</i> <i>Nocardia asiatica</i> <i>Nocardia exalbida</i>	Implement new Matching Hint	For the species abscessus / arthritidis / asiatica / exalbida of the genus Nocardia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.
<i>Paenibacillus amylolyticus</i>	Implement new Matching Hint	For the species amylolyticus / cucumis / taichungensis / tundrae / tylopili of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.	Due to the newly implemented species <i>Paenibacillus tylopili</i> recent results have shown that the displayed species could not be distinguished by 16S rRNA gene sequencing. Therefore, all these species should be taken into consideration as a possible result.
<i>Pasteurella canis</i>	Implement new Matching Hint	Species canis / oralis of the genus Pasteurella have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Pasteurella oralis</i> implementation of a Matching Hint is needed.
<i>Peptoniphilus harei</i> <i>Peptoniphilus indolicus</i>	Link to Matching Hint deleted	N/A	Due to renaming of reference entries to <i>Peptoniphilus</i> sp[2] a clear separation is possible now and Therefore, the Matching Hint is not needed anymore.
<i>Peptoniphilus</i> sp[2]	Link to Matching Hint deleted	N/A	Due to renaming of reference entries the Matching Hint is not needed anymore.
<i>Serratia grimesii</i> <i>Serratia liquefaciens</i>	Update wording of Matching Hint	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Serratia proteamaculans</i> manifests their close relationship. Therefore, the existing Matching Hint needs to be adjusted.
<i>Serratia marcescens</i>	Update wording of Matching Hint	Species marcescens / nematodiphila / ureilytica of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Serratia nematodiphila</i> the existing Matching Hint needs to be extended.



Species	Action	DB-9607	Justification
<i>Serratia plymuthica</i>	Implement new Matching Hint	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Serratia proteamaculans</i> manifests their close relationship. Therefore, a Matching Hint needs to be implemented.
<i>Serratia proteamaculans</i>	Update wording of Matching Hint	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Serratia proteamaculans</i> manifests their close relationship. Therefore, the existing Matching Hint needs to be adjusted.
<i>Serratia quinivorans</i>	Implement new Matching Hint	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	The improved number of reference entries for <i>Serratia proteamaculans</i> manifests their close relationship. Therefore, a Matching Hint needs to be implemented.
<i>Serratia ureilytica</i>	Update wording of Matching Hint	Species marcescens / nematodiphila / ureilytica of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.	Due to the newly implemented species <i>Serratia nematodiphila</i> the existing Matching Hint needs to be extended.
<i>Tsukamurella hongkongensis</i>	Implement new Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to the implementation of a Matching Hint.
<i>Tsukamurella inchonensis</i> <i>Tsukamurella paurometabola</i> <i>Tsukamurella pseudospumae</i> <i>Tsukamurella pulmonis</i>	Update wording of Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to a more general wording of the Matching Hint.
<i>Tsukamurella serpentis</i>	Implement new Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to the implementation of a Matching Hint.
<i>Tsukamurella sinensis</i>	Update wording of Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to a more general wording of the Matching Hint.
<i>Tsukamurella soli</i>	Implement new Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to the implementation of a Matching Hint.



Species	Action	DB-9607	Justification
<i>Tsukamurella spumae</i> <i>Tsukamurella strandjordii</i> <i>Tsukamurella tyrosinosolvans</i>	Update wording of Matching Hint	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.	Reassessment of the matching situation leads to a more general wording of the Matching Hint.
Weissella cibaria Weissella confusa	Link to Matching Hint deleted	N/A	Due to the improved number of reference entries for <i>Weissella cibaria</i> a clear discrimination is possible now and Therefore, the Matching Hint is not needed anymore.



6 Genera/Species improved by implementation of further reference entries

Table 7: Genera/Species improved

<i>Achromobacter mucicolens</i>	<i>Acinetobacter lactucae</i>	<i>Acinetobacter vivianii</i>	<i>Aldersonia kunmingensis</i>
<i>Acidovorax delafieldii</i>	<i>Acinetobacter lwoffii</i>	<i>Actinobacillus equuli</i>	<i>Aliarcobacter butzleri</i>
<i>Acinetobacter baylyi</i>	<i>Acinetobacter modestus</i>	<i>Actinobacillus porcinus</i>	<i>Aliarcobacter skirrowii</i>
<i>Acinetobacter beijerinckii</i>	<i>Acinetobacter nectaris</i>	<i>Actinobacillus rossii</i>	<i>Anaerococcus vaginalis</i>
<i>Acinetobacter bereziniae</i>	<i>Acinetobacter parvus</i>	<i>Actinobacillus suis</i>	<i>Anoxybacillus flavithermus</i>
<i>Acinetobacter bouvetii</i>	<i>Acinetobacter pittii</i>	<i>Actinomyces bowdenii</i>	<i>Arcanobacterium hippocoleae</i>
<i>Acinetobacter brisouii</i>	<i>Acinetobacter proteolyticus</i>	<i>Actinomyces denticolens</i>	<i>Arcanobacterium phocisimile</i>
<i>Acinetobacter courvalinii</i>	<i>Acinetobacter puyangensis</i>	<i>Actinomyces naeslundii</i>	<i>Arthrobacter citreus</i>
<i>Acinetobacter dispersus</i>	<i>Acinetobacter radioresistens</i>	<i>Actinomyces oris</i>	<i>Arthrobacter gandavensis</i>
<i>Acinetobacter gernerii</i>	<i>Acinetobacter rudis</i>	<i>Actinomyces radidentis</i>	<i>Arthrobacter koreensis</i>
<i>Acinetobacter guillouiae</i>	<i>Acinetobacter schindleri</i>	<i>Actinomyces viscosus</i>	<i>Arthrobacter monumenti</i>
<i>Acinetobacter haemolyticus</i>	<i>Acinetobacter seifertii</i>	<i>Actinotignum sanguinis</i>	<i>Arthrobacter parietis</i>
<i>Acinetobacter harbinensis</i>	<i>Acinetobacter tandoii</i>	<i>Actinotignum schaalii</i>	<i>Arthrobacter pascens</i>
<i>Acinetobacter indicus</i>	<i>Acinetobacter towneri</i>	<i>Advenella incenata</i>	<i>Arthrobacter psychrolactophilus</i>
<i>Acinetobacter johnsonii</i>	<i>Acinetobacter ursingii</i>	<i>Advenella kashmirensis</i>	<i>Arthrobacter roseus</i>
<i>Acinetobacter junii</i>	<i>Acinetobacter variabilis</i>	<i>Aerococcus urinae</i>	<i>Arthrobacter stackebrandtii</i>
<i>Acinetobacter kookii</i>	<i>Acinetobacter venetianus</i>	<i>Agromyces rhizosphaerae</i>	<i>Arthrobacter tecti</i>



<i>Arthrobacter tumbae</i>	<i>Bacillus mycoides</i>	<i>Bifidobacterium choerinum</i>	<i>Brevibacterium otitidis</i>
<i>Arthrobacter woluwensis</i>	<i>Bacillus okhensis</i>	<i>Bifidobacterium dentium</i>	<i>Brevibacterium paucivorans</i>
<i>Atopobium deltae</i>	<i>Bacillus oshimensis</i>	<i>Bifidobacterium gallicum</i>	<i>Brevibacterium ravensturgense</i>
<i>Atopobium fossor</i>	<i>Bacillus subterraneus</i>	<i>Bifidobacterium indicum</i>	<i>Brevibacterium senegalense</i>
<i>Bacillus algicola</i>	<i>Bacillus thioparans</i>	<i>Bifidobacterium longum</i>	<i>Brochothrix thermosphacta</i>
<i>Bacillus amyloliquefaciens</i>	<i>Bacteroides eggerthii</i>	<i>Bifidobacterium magnum</i>	<i>Buchananella hordeovulneris</i>
<i>Bacillus asahii</i>	<i>Bacteroides finegoldii</i>	<i>Bifidobacterium minimum</i>	<i>Campylobacter fetus</i>
<i>Bacillus badius</i>	<i>Bacteroides helcogenes</i>	<i>Bifidobacterium porcinum</i>	<i>Campylobacter hyointestinalis</i>
<i>Bacillus benzoovorans</i>	<i>Bacteroides massiliensis</i>	<i>Bifidobacterium pseudocatenulatum</i>	<i>Candida nemodendra</i>
<i>Bacillus circulans</i>	<i>Bacteroides stercoris</i>	<i>Bifidobacterium pseudolongum</i>	<i>Candida norvegica</i>
<i>Bacillus clausii</i>	<i>Bergeyella zoohelcum</i>	<i>Bifidobacterium pullorum</i>	<i>Candida parapsilosis</i>
<i>Bacillus coagulans</i>	<i>Bifidobacterium adolescentis</i>	<i>Bifidobacterium ruminantium</i>	<i>Candida tropicalis</i>
<i>Bacillus firmus</i>	<i>Bifidobacterium angulatum</i>	<i>Bifidobacterium thermacidophilum</i>	<i>Carnobacterium gallinarum</i>
<i>Bacillus galliciensis</i>	<i>Bifidobacterium animalis</i>	<i>Bisgaardia hudsonensis</i>	<i>Carnobacterium maltaromaticum</i>
<i>Bacillus halosaccharovorans</i>	<i>Bifidobacterium asteroides</i>	<i>Boudabousia marimammalium</i>	<i>Cedecea davisae</i>
<i>Bacillus idriensis</i>	<i>Bifidobacterium bifidum</i>	<i>Brachybacterium muris</i>	<i>Chryseobacterium pallidum</i>
<i>Bacillus infantis</i>	<i>Bifidobacterium boum</i>	<i>Brevibacterium casei</i>	<i>Chryseobacterium shandongense</i>
<i>Bacillus licheniformis</i>	<i>Bifidobacterium breve</i>	<i>Brevibacterium celere</i>	<i>Clostridium beijerinckii</i>
<i>Bacillus mojavensis</i>	<i>Bifidobacterium catenulatum</i>	<i>Brevibacterium luteolum</i>	<i>Clostridium cadaveris</i>



<i>Clostridium sartagoforme</i>	<i>Deinococcus aerolatus</i>	<i>Escherichia hermannii</i>	<i>Glutamicibacter bergerei</i>
<i>Clostridium saudiense</i>	<i>Deinococcus proteolyticus</i>	<i>Eubacterium tenue</i>	<i>Glutamicibacter creatinolyticus</i>
<i>Clostridium tetani</i>	<i>Deinococcus wulumuqiensis</i>	<i>Falsarthrobacter nasiphocae</i>	<i>Glutamicibacter mysorens</i>
<i>Comamonas aquatica</i>	<i>Dermabacter vaginalis</i>	<i>Fastidiosipila sanguinis</i>	<i>Gordonia bronchialis</i>
<i>Comamonas terrigena</i>	<i>Dermacoccus nishinomiyaensis</i>	<i>Fenollaria massiliensis</i>	<i>Gordonia sputi</i>
<i>Corynebacterium argensoratense</i>	<i>Desemzia incerta</i>	<i>Flavobacterium chungbukense</i>	<i>Hanseniaspora valbyensis</i>
<i>Corynebacterium auriscanis</i>	<i>Desulfovibrio simplex</i>	<i>Flavobacterium columnare</i>	<i>Hanseniaspora vineae</i>
<i>Corynebacterium confusum</i>	<i>Dietzia cinnamea</i>	<i>Flavobacterium daejeonense</i>	<i>Hathewayia histolytica</i>
<i>Corynebacterium durum</i>	<i>Dietzia papillomatosis</i>	<i>Flavobacterium fryxellicola</i>	<i>Histophilus somni</i>
<i>Corynebacterium imitans</i>	<i>Dysgonomonas capnocytophagoideis</i>	<i>Flavobacterium gelidilacus</i>	<i>Hungatella hathewayi</i>
<i>Corynebacterium jeikeium</i>	<i>Dysgonomonas gadei</i>	<i>Flavobacterium limicola</i>	<i>Hydrogenibacillus schlegelii</i>
<i>Corynebacterium kroppenstedtii</i>	<i>Elizabethkingia anophelis</i>	<i>Flavonifractor plautii</i>	<i>Ideonella dechloratans</i>
<i>Corynebacterium phoceense</i>	<i>Elizabethkingia meningoseptica</i>	<i>Frederiksenia canicola</i>	<i>Ignatzschineria indica</i>
<i>Corynebacterium pyruviciproducens</i>	<i>Elizabethkingia miricola</i>	<i>Fusobacterium varium</i>	<i>Inquilinus limosus</i>
<i>Corynebacterium renale</i>	<i>Empedobacter brevis</i>	<i>Geotrichum candidum</i>	<i>Janibacter hoylei</i>
<i>Corynebacterium simulans</i>	<i>Empedobacter falsenii</i>	<i>Globicatella sanguinis</i>	<i>Janibacter indicus</i>
<i>Corynebacterium tuberculostearicum</i>	<i>Enterocloster aldensis</i>	<i>Globicatella sulfidifaciens</i>	<i>Kandleria vitulina</i>
<i>Corynebacterium ulcerans</i>	<i>Enterococcus malodoratus</i>	<i>Glutamicibacter ardleyensis</i>	<i>Kerstersonia gyiorum</i>
<i>Cupriavidus basilensis</i>	<i>Enterococcus raffinosus</i>	<i>Glutamicibacter arilaitensis</i>	<i>Kingella oralis</i>



<i>Kluyvera ascorbata</i>	<i>Lactobacillus collinoides</i>	<i>Lactobacillus vaginalis</i>	<i>Microbacterium imperiale</i>
<i>Kluyvera cryocrescens</i>	<i>Lactobacillus crispatus</i>	<i>Lactobacillus zymae</i>	<i>Microbacterium ketosireducens</i>
<i>Kluyvera georgiana</i>	<i>Lactobacillus delbrueckii</i>	<i>Lactococcus garvieae</i>	<i>Microbacterium koreense</i>
<i>Kocuria palustris</i>	<i>Lactobacillus farciminis</i>	<i>Legionella maceachernii</i>	<i>Microbacterium lacticum</i>
<i>Kodamaea ohmeri</i>	<i>Lactobacillus frumenti</i>	<i>Leptotrichia trevisanii</i>	<i>Microbacterium laevaniformans</i>
<i>Kosakonia cowanii</i>	<i>Lactobacillus fuchuensis</i>	<i>Leuconostoc lactis</i>	<i>Microbacterium liquefaciens</i>
<i>Kurthia gibsonii</i>	<i>Lactobacillus gallinarum</i>	<i>Listeria costaricensis</i>	<i>Microbacterium maritypicum</i>
<i>Kurthia zopfii</i>	<i>Lactobacillus gasseri</i>	<i>Megasphaera elsdenii</i>	<i>Microbacterium oleivorans</i>
<i>Kytococcus schroeteri</i>	<i>Lactobacillus gastricus</i>	<i>Metschnikowia reukaufii</i>	<i>Microbacterium oxydans</i>
<i>Kytococcus sedentarius</i>	<i>Lactobacillus hamsteri</i>	<i>Microbacterium aerolatum</i>	<i>Microbacterium paludicola</i>
<i>Lachancea fermentati</i>	<i>Lactobacillus helveticus</i>	<i>Microbacterium arborescens</i>	<i>Microbacterium paraoxydans</i>
<i>Lachnoanaerobaculum orale</i>	<i>Lactobacillus iners</i>	<i>Microbacterium aurum</i>	<i>Microbacterium phyllosphaerae</i>
<i>Lactobacillus acidifarinae</i>	<i>Lactobacillus johnsonii</i>	<i>Microbacterium barkeri</i>	<i>Microbacterium resistens</i>
<i>Lactobacillus acidipiscis</i>	<i>Lactobacillus kefiri</i>	<i>Microbacterium dextranolyticum</i>	<i>Microbacterium sp</i>
<i>Lactobacillus acidophilus</i>	<i>Lactobacillus paracasei</i>	<i>Microbacterium flavescens</i>	<i>Microbacterium terrae</i>
<i>Lactobacillus amylovorus</i>	<i>Lactobacillus pontis</i>	<i>Microbacterium ginsengisoli</i>	<i>Microbacterium terregens</i>
<i>Lactobacillus apodemi</i>	<i>Lactobacillus rhamnosus</i>	<i>Microbacterium halotolerans</i>	<i>Microbacterium testaceum</i>
<i>Lactobacillus aviarius</i>	<i>Lactobacillus salivarius</i>	<i>Microbacterium hominis</i>	<i>Microbacterium thalassium</i>
<i>Lactobacillus coleohominis</i>	<i>Lactobacillus sanfranciscensis</i>	<i>Microbacterium hydrocarbonoxydans</i>	<i>Microbacterium ulmi</i>



<i>Micrococcus flavus</i>	<i>Nocardiosis dassonvillei</i>	<i>Paeniglutamicibacter sulfureus</i>	<i>Peptoniphilus lacydonensis</i>
<i>Micrococcus terreus</i>	<i>Oceanobacillus caeni</i>	<i>Pantoea septica</i>	<i>Peptoniphilus olsenii</i>
<i>Mixta calida</i>	<i>Oceanobacillus kimchii</i>	<i>Parabacteroides faecis</i>	<i>Peptoniphilus tyrrelliae</i>
<i>Mixta gaviniae</i>	<i>Oceanobacillus oncorhynchi</i>	<i>Parabacteroides johnsonii</i>	<i>Peptostreptococcus canis</i>
<i>Moellerella wisconsensis</i>	<i>Paenicaligenes suwonensis</i>	<i>Parabacteroides merdae</i>	<i>Photobacterium iliopiscarium</i>
<i>Moesziomyces aphidis</i>	<i>Paenarthrobacter histidinovorans</i>	<i>Paracoccus denitrificans</i>	<i>Pichia cactophila</i>
<i>Moorella thermoacetica</i>	<i>Paenarthrobacter ilicis</i>	<i>Paracoccus yeei</i>	<i>Pichia kudriavzevii</i>
<i>Moraxella atlantae</i>	<i>Paenarthrobacter nicotinovorans</i>	<i>Pasteurella aerogenes</i>	<i>Porphyromonas levii</i>
<i>Murdochiella asaccharolytica</i>	<i>Paenibacillus agaridevorans</i>	<i>Pasteurella bettyae</i>	<i>Prevotella baroniae</i>
<i>Nakazawaea ernobii</i>	<i>Paenibacillus anaericus</i>	<i>Pasteurella canis</i>	<i>Prevotella bergensis</i>
<i>Negativicoccus succinicivorans</i>	<i>Paenibacillus chibensis</i>	<i>Pasteurella dagmatis</i>	<i>Prevotella disiens</i>
<i>Neisseria dumasiana</i>	<i>Paenibacillus hunanensis</i>	<i>Pasteurella stomatis</i>	<i>Prevotella nanceiensis</i>
<i>Neisseria oralis</i>	<i>Paenibacillus macquariensis</i>	<i>Pauljensenia hongkongensis</i>	<i>Prevotella oralis</i>
<i>Neisseria zoodegmatis</i>	<i>Paenibacillus pasadenensis</i>	<i>Pediococcus acidilactici</i>	<i>Prevotella oris</i>
<i>Nocardia abscessus</i>	<i>Paenibacillus provencensis</i>	<i>Pediococcus pentosaceus</i>	<i>Prevotella salivae</i>
<i>Nocardia brasiliensis</i>	<i>Paenibacillus residui</i>	<i>Peptococcus niger</i>	<i>Priceomyces carsonii</i>
<i>Nocardia salmonicida</i>	<i>Paenibacillus timonensis</i>	<i>Peptoniphilus coxii</i>	<i>Propionibacterium australiense</i>
<i>Nocardia wallacei</i>	<i>Paenibacillus urinalis</i>	<i>Peptoniphilus harei</i>	<i>Providencia heimbachae</i>
<i>Nocardiosis alba</i>	<i>Paeniglutamicibacter kerguelensis</i>	<i>Peptoniphilus ivorii</i>	<i>Pseudarthrobacter oxydans</i>



<i>Pseudarthrobacter polychromogenes</i>	<i>Rummeliibacillus pycnus</i>	<i>Stenotrophomonas acidaminiphila</i>	<i>Trueperella abortisuis</i>
<i>Pseudarthrobacter scleromae</i>	<i>Saccharomyces paradoxus</i>	<i>Streptococcus dentirousetti</i>	<i>Tsukamurella pulmonis</i>
<i>Pseudarthrobacter sulfonivorans</i>	<i>Salinivibrio costicola</i>	<i>Streptococcus dysgalactiae</i>	<i>Varibaculum cambriense</i>
<i>Pseudescherichia vulneris</i>	<i>Schaalia canis</i>	<i>Streptococcus equi_ssp_equi</i>	<i>Veillonella atypica</i>
<i>Pseudomonas alcaligenes</i>	<i>Serratia proteamaculans</i>	<i>Streptococcus equinus</i>	<i>Vibrio natriegens</i>
<i>Pseudomonas graminis</i>	<i>Sinomonas atrocyanea</i>	<i>Streptococcus gordonii</i>	<i>Virgibacillus halodenitrificans</i>
<i>Pseudomonas guariconensis</i>	<i>Solobacterium moorei</i>	<i>Streptococcus mutans</i>	<i>Virgibacillus pantothenicus</i>
<i>Pseudomonas kuykendallii</i>	<i>Sphingomonas aerolata</i>	<i>Streptococcus pluranimalium</i>	<i>Weissella cibaria</i>
<i>Pseudomonas pohangensis</i>	<i>Sporosarcina luteola</i>	<i>Streptococcus porci</i>	<i>Weissella paramesenteroides</i>
<i>Pseudomonas segetis</i>	<i>Staphylococcus argenteus</i>	<i>Streptococcus pyogenes</i>	<i>Wohlfahrtiimonas chitiniclastica</i>
<i>Pseudopropionibacterium propionicum</i>	<i>Staphylococcus aureus</i>	<i>Streptococcus salivarius</i>	<i>Xanthomonas arboricola</i>
<i>Pseudoxanthomonas spadix</i>	<i>Staphylococcus chromogenes</i>	<i>Streptococcus suis</i>	<i>Xanthomonas bromi</i>
<i>Psychrobacter piechaudii</i>	<i>Staphylococcus cohnii</i>	<i>Streptococcus thoralensis</i>	<i>Xanthomonas campestris</i>
<i>Rahnella aquatilis</i>	<i>Staphylococcus epidermidis</i>	<i>Streptococcus uberis</i>	<i>Xanthomonas codiaei</i>
<i>Ralstonia insidiosa</i>	<i>Staphylococcus equorum</i>	<i>Streptomyces thermoviolaceus</i>	<i>Xanthomonas cucurbitae</i>
<i>Raoultella planticola</i>	<i>Staphylococcus kloosii</i>	<i>Taylorella equigenitalis</i>	<i>Xanthomonas cynarae</i>
<i>Raoultella terrigena</i>	<i>Staphylococcus lentus</i>	<i>Tissierella praeacuta</i>	<i>Xanthomonas hortorum</i>
<i>Roseomonas mucosa</i>	<i>Staphylococcus schweitzeri</i>	<i>Torulaspora delbrueckii</i>	<i>Xanthomonas hyacinthi</i>
<i>Rothia koreensis</i>	<i>Starmerella magnoliae</i>	<i>Trichomonascus ciferrii</i>	<i>Xanthomonas perforans</i>



<i>Xanthomonas pisi</i>	<i>Xanthomonas theicola</i>	<i>Xanthomonas vasicola</i>	<i>Yarrowia deformans</i>
-------------------------	-----------------------------	-----------------------------	---------------------------



7 Installation Instructions

- Make sure that the previously installed library has been updated with the last update: 8468 MSP.

If not, please contact Bruker.

- Run the RUO setup file: *MaldiBiotyperDBUpdate_V10.0.0.0_8468-9607(RUO).exe*.

Note: To run setups requires administrator rights and can take several hours to complete. Please be aware that a once started MSP library update must not be interrupted. Any interruption would leave the installed MSP library in an intermediate and therefore invalid state.

Note: Users of laboratory-specific specimen codes please refer to section 8.2.

- After successful installation, the appropriate MBT-RUO client applications will use the updated MBT Compass Library containing 9607 MSP.
- For **MBT Compass HT** system the library installation is done via installing the appropriate library module. Please do the following:
 - Log into the MBT Compass HT software as **Laboratory Manager**.
 - Visit the menu **Configuration > Modules > Local** and select **Browse for module**.
 - Select the library module archive *MaldiBiotyperLibrary_Module_V10.0.0.0_9607-MSPs_(RUO).zip* from the *MBT Compass Library* installation package (subfolder *MBT-Compass-HT*)
 - Click on **Install**.



8 Appendix

8.1 Matching Hints

The following table lists the Matching Hints included in the MBT Compass Library.

Table 8: Matching Hints

Genus/Species	Matching Hint
<i>Acetobacter aceti</i>	For the species <i>aceti</i> / <i>cerevisiae</i> / <i>farinalis</i> / <i>malorum</i> / <i>orleanensis</i> / <i>persici</i> of the genus <i>Acetobacter</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Acetobacter cerevisiae</i>	
<i>Acetobacter malorum</i>	
<i>Acetobacter persici</i>	
<i>Achromobacter ruhlandii</i>	Species <i>ruhlandii</i> / <i>xylosoxidans</i> of the genus <i>Achromobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.
<i>Achromobacter xylosoxidans</i>	
<i>Acidovorax avenae</i>	For the species <i>avenae</i> / <i>cattleyae</i> / <i>citrulli</i> of the genus <i>Acidovorax</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Acidovorax cattleyae</i>	
<i>Acidovorax citrulli</i>	
<i>Acidovorax delafieldii</i>	For the species <i>de-la-fieldii</i> / <i>kalamii</i> of the genus <i>Acidovorax</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Acidovorax valerianellae</i>	For the species <i>valerianellae</i> / <i>wautersii</i> of the genus <i>Acidovorax</i> the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Acidovorax wautersii</i>	
<i>Acinetobacter baumannii</i>	Member of the <i>Acinetobacter baumannii</i> complex. Extraction must be performed to permit reliable species identification.
<i>Acinetobacter baylyi</i>	Species <i>baylyi</i> / <i>solii</i> of the genus <i>Acinetobacter</i> have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Acinetobacter bereziniae</i>	Species bereziniae / guillouiae / wuhouensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter bouvetii</i>	Species bouvetii / pragensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter calcoaceticus</i>	Member of the Acinetobacter baumannii complex. Extraction must be performed to permit reliable species identification.
<i>Acinetobacter colistiniresistens</i>	Species colistiniresistens / courvalinii / dispersus / gyllenbergii / proteolyticus / vivianii of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter courvalinii</i>	
<i>Acinetobacter dispersus</i>	
<i>Acinetobacter guillouiae</i>	Species bereziniae / guillouiae / wuhouensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter gyllenbergii</i>	Species colistiniresistens / courvalinii / dispersus / gyllenbergii / proteolyticus / vivianii of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter lactucae</i>	Member of the Acinetobacter baumannii complex. Extraction must be performed to permit reliable species identification.
<i>Acinetobacter nosocomialis</i>	
<i>Acinetobacter pittii</i>	
<i>Acinetobacter populi</i>	Species populi / puyangensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter pragensis</i>	Species bouvetii / pragensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter proteolyticus</i>	Species colistiniresistens / courvalinii / dispersus / gyllenbergii / proteolyticus / vivianii of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter puyangensis</i>	Species populi / puyangensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter seifertii</i>	Member of the Acinetobacter baumannii complex. Extraction must be performed to permit reliable species identification.
<i>Acinetobacter soli</i>	Species baylyi / soli of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Acinetobacter vivianii</i>	Species colistiniresistens / courvalinii / dispersus / gyllenbergii / proteolyticus / vivianii of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Acinetobacter wuhouensis</i>	Species bereziniae / guillouiae / wuhouensis of the genus Acinetobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Actinobacillus arthritidis</i>	For the species arthritidis / capsulatus / equuli / lignieresii / pleuropneumoniae / suis of the genus Actinobacillus the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Actinobacillus capsulatus</i>	
<i>Actinobacillus equuli</i>	
<i>Actinobacillus lignieresii</i>	
<i>Actinobacillus pleuropneumoniae</i>	
<i>Actinobacillus suis</i>	
<i>Actinomyces denticolens</i>	For the species denticolens / timonensis of the genus Actinomyces the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Actinomyces timonensis</i>	
<i>Actinotignum sanguinis</i>	Species sanguinis / schaalii of the genus Actinotignum have very similar patterns: Therefore distinguishing their species is difficult.
<i>Actinotignum schaalii</i>	Synonym of Actinobaculum schaalii. Species sanguinis / schaalii of the genus Actinotignum have very similar patterns: Therefore distinguishing their species is difficult.
<i>Advenella incenata</i>	For the species incenata / kashmirensis of the genus Advenella the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Advenella kashmirensis</i>	
<i>Aeribacillus pallidus</i>	For the species composti / pallidus of the genus Aeribacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Aerococcus viridans</i>	Aerococcus viridans is/are closely related to the rare species Aerococcus urinaeequi which is/are not included in the MALDI Biotyper database.
<i>Aeromicrobium ginsengisoli</i>	For the species ginsengisoli / panaciterrae of the genus Aeromicrobium the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Aeromonas bestiarum</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Aeromonas caviae</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Aeromonas encheleia</i>	
<i>Aeromonas enteropelogenes</i>	
<i>Aeromonas eucrenophila</i>	
<i>Aeromonas hydrophila</i>	
<i>Aeromonas ichthiosmia</i>	
<i>Aeromonas jandaei</i>	
<i>Aeromonas media</i>	
<i>Aeromonas molluscorum</i>	
<i>Aeromonas popoffii</i>	
<i>Aeromonas punctata</i>	
<i>Aeromonas salmonicida</i>	
<i>Aeromonas schubertii</i>	
<i>Aeromonas simiae</i>	
<i>Aeromonas sobria</i>	
<i>Aeromonas sp[2]</i>	
<i>Aeromonas veronii</i>	
<i>Aggregatibacter aphrophilus</i>	synonym of <i>Haemophilus aphrophilus</i>
<i>Amycolatopsis alba</i>	For the species <i>alba</i> / <i>azurea</i> / <i>coloradensis</i> / <i>japonica</i> / <i>keratiniphila</i> / <i>lurida</i> / <i>orientalis</i> of the genus <i>Amycolatopsis</i> the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.



Genus/Species	Matching Hint
<i>Amycolatopsis azurea</i>	For the species alba / azurea / coloradensis / japonica / keratiniphila / lurida / orientalis of the genus Amycolatopsis the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Amycolatopsis balhimycina</i>	Species balhimycina / kentuckyensis / lexingtonensis / mediterranei / pretoriensis / sulphurea / tolypomycina of the genus Amycolatopsis have very similar patterns: Therefore distinguishing their species is difficult.
<i>Amycolatopsis coloradensis</i>	For the species alba / azurea / coloradensis / japonica / keratiniphila / lurida / orientalis of the genus Amycolatopsis the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Amycolatopsis japonica</i>	
<i>Amycolatopsis kentuckyensis</i>	Species balhimycina / kentuckyensis / lexingtonensis / mediterranei / pretoriensis / sulphurea / tolypomycina of the genus Amycolatopsis have very similar patterns: Therefore distinguishing their species is difficult.
<i>Amycolatopsis keratiniphila</i>	For the species alba / azurea / coloradensis / japonica / keratiniphila / lurida / orientalis of the genus Amycolatopsis the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Amycolatopsis lexingtonensis</i>	Species balhimycina / kentuckyensis / lexingtonensis / mediterranei / pretoriensis / sulphurea / tolypomycina of the genus Amycolatopsis have very similar patterns: Therefore distinguishing their species is difficult.
<i>Amycolatopsis lurida</i>	For the species alba / azurea / coloradensis / japonica / keratiniphila / lurida / orientalis of the genus Amycolatopsis the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Amycolatopsis mediterranei</i>	Species balhimycina / kentuckyensis / lexingtonensis / mediterranei / pretoriensis / sulphurea / tolypomycina of the genus Amycolatopsis have very similar patterns: Therefore distinguishing their species is difficult.
<i>Amycolatopsis pretoriensis</i>	
<i>Amycolatopsis sulphurea</i>	
<i>Amycolatopsis tolypomycina</i>	
<i>Anoxybacillus amylolyticus</i>	Species amylolyticus / contaminans / voinovskiensis of the genus Anoxybacillus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Anoxybacillus contaminans</i>	
<i>Anoxybacillus flavithermus</i>	For the species eryuanensis / flavithermus / kestanbolensis / mongoliensis of the genus Anoxybacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.



Genus/Species	Matching Hint
<i>Anoxybacillus rupiensis</i>	For the species geothermalis / rupiensis of the genus Anoxybacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Anoxybacillus voinovskiensis</i>	Species amylolyticus / contaminans / voinovskiensis of the genus Anoxybacillus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Arthrobacter castelli</i>	Species tend to produce polymers which can interfere with the identification.
<i>Arthrobacter pigmenti</i>	
<i>Bacillus acidicola</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus agaradhaerens</i>	
<i>Bacillus akibai</i>	
<i>Bacillus alcalophilus</i>	
<i>Bacillus algicola</i>	
<i>Bacillus altitudinis</i>	
<i>Bacillus alveayuensis</i>	
<i>Bacillus amyloliquefaciens</i>	
<i>Bacillus amyloliquefaciens_ssp_plantarum</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus aquimaris</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus asahii</i>	
<i>Bacillus atrophaeus</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus azotoformans</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus badius</i>	
<i>Bacillus bataviensis</i>	



Genus/Species	Matching Hint
<i>Bacillus benzoovorans</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus beringensis</i>	
<i>Bacillus campisalis</i>	
<i>Bacillus carboniphilus</i>	
<i>Bacillus cellulosityticus</i>	
<i>Bacillus cereus</i>	Bacillus anthracis, cereus, mycoides, pseudomycoides and thuringiensis are closely related and members of the Bacillus cereus group. In particular Bacillus cereus spectra are very similar to spectra from Bacillus anthracis. Bacillus anthracis is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus chagannorensis</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus cibi</i>	
<i>Bacillus circulans</i>	
<i>Bacillus clarkii</i>	
<i>Bacillus clausii</i>	
<i>Bacillus coagulans</i>	
<i>Bacillus cohnii</i>	
<i>Bacillus cytotoxicus</i>	
<i>Bacillus decolorationis</i>	
<i>Bacillus drentensis</i>	
<i>Bacillus endophyticus</i>	
<i>Bacillus farraginis</i>	



Genus/Species	Matching Hint
<i>Bacillus fastidiosus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus firmus</i>	
<i>Bacillus flexus</i>	
<i>Bacillus foraminis</i>	
<i>Bacillus fordii</i>	
<i>Bacillus fortis</i>	
<i>Bacillus funiculus</i>	
<i>Bacillus galactosidilyticus</i>	
<i>Bacillus galliciensis</i>	
<i>Bacillus gibsonii</i>	
<i>Bacillus ginsengihumi</i>	
<i>Bacillus gossypii</i>	
<i>Bacillus graminis</i>	
<i>Bacillus halmapalus</i>	
<i>Bacillus halodurans</i>	
<i>Bacillus halosaccharovorans</i>	
<i>Bacillus hemicellulosilyticus</i>	
<i>Bacillus horikoshii</i>	
<i>Bacillus horneckiae</i>	
<i>Bacillus horti</i>	



Genus/Species	Matching Hint
<i>Bacillus humi</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus hwajinpoensis</i>	
<i>Bacillus idriensis</i>	
<i>Bacillus indicus</i>	
<i>Bacillus infantis</i>	
<i>Bacillus jeotgali</i>	
<i>Bacillus kochii</i>	
<i>Bacillus koreensis</i>	
<i>Bacillus krulwichiae</i>	
<i>Bacillus kyonggiensis</i>	
<i>Bacillus lentus</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus licheniformis</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus litoralis</i>	
<i>Bacillus luciferensis</i>	
<i>Bacillus mannanilyticus</i>	
<i>Bacillus marisflavi</i>	
<i>Bacillus massiliosenegalensis</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus megaterium</i>	
<i>Bacillus mojavensis</i>	



Genus/Species	Matching Hint
<i>Bacillus muralis</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus mycooides</i>	Bacillus anthracis, cereus, mycooides, pseudomycooides and thuringiensis are closely related and members of the Bacillus cereus group. In particular Bacillus cereus spectra are very similar to spectra from Bacillus anthracis. Bacillus anthracis is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus nealsonii</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus niabensis</i>	
<i>Bacillus niacini</i>	
<i>Bacillus novalis</i>	
<i>Bacillus oceanisediminis</i>	
<i>Bacillus odysseyi</i>	
<i>Bacillus okhensis</i>	
<i>Bacillus okuhidensis</i>	
<i>Bacillus oleronius</i>	
<i>Bacillus oshimensis</i>	
<i>Bacillus patagoniensis</i>	
<i>Bacillus pocheonensis</i>	
<i>Bacillus pseudalcaliphilus</i>	
<i>Bacillus pseudofirmus</i>	



Genus/Species	Matching Hint
<i>Bacillus pseudomycooides</i>	Bacillus anthracis, cereus, mycooides, pseudomycooides and thuringiensis are closely related and members of the Bacillus cereus group. In particular Bacillus cereus spectra are very similar to spectra from Bacillus anthracis. Bacillus anthracis is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus psychrosaccharolyticus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus pumilus</i>	
<i>Bacillus raris</i>	
<i>Bacillus safensis</i>	
<i>Bacillus salarius</i>	
<i>Bacillus seohaeanensis</i>	
<i>Bacillus shackletonii</i>	
<i>Bacillus simplex</i>	
<i>Bacillus sivalis</i>	
<i>Bacillus smithii</i>	
<i>Bacillus soli</i>	
<i>Bacillus sonorensis</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus sp</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus sp</i>[2]	
<i>Bacillus sp</i>[4]	
<i>Bacillus sporothermodurans</i>	
<i>Bacillus subterraneus</i>	



Genus/Species	Matching Hint
<i>Bacillus subtilis</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus thermoamylovorans</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus thioparans</i>	
<i>Bacillus thuringiensis</i>	Bacillus anthracis, cereus, mycoides, pseudomycoides and thuringiensis are closely related and members of the Bacillus cereus group. In particular Bacillus cereus spectra are very similar to spectra from Bacillus anthracis. Bacillus anthracis is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus timonensis</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus vallismortis</i>	is a member of Bacillus subtilis group. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus vedderi</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Bacillus vietnamensis</i>	
<i>Bacillus vireti</i>	
<i>Bacillus wakoensis</i>	
<i>Bacillus zeae</i>	
<i>Bacteroides eggerthii</i>	Species eggerthii / stercoris of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult.
<i>Bacteroides faecis</i>	Species faecis / thetaiotaomicron of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult.
<i>Bacteroides ovatus</i>	Species ovatus/xylanisolvans of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult. Bacteroides xylanisolvans is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Bacteroides stercoris</i>	Species eggerthii / stercoris of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Bacteroides thetaiotaomicron</i>	Species faecis / thetaiotaomicron of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult.
<i>Bacteroides vulgatus</i>	Species vulgatus/dorei of the genus Bacteroides have very similar patterns: Therefore distinguishing their species is difficult. Bacteroides dorei is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Bifidobacterium porcinum</i>	Species porcinum / thermacidophilum / thermophilum of the genus Bifidobacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Bifidobacterium thermacidophilum</i>	
<i>Bifidobacterium thermophilum</i>	
<i>Bordetella bronchialis</i>	For the species bronchialis / sputigena of the genus Bordetella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Bordetella bronchiseptica</i>	Species bronchiseptica / pertussis / parapertussis of the genus Bordetella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Bordetella parapertussis</i>	
<i>Bordetella pertussis</i>	
<i>Bordetella sputigena</i>	For the species bronchialis / sputigena of the genus Bordetella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Brachybacterium conglomeratum</i>	For the species conglomeratum / paraconglomeratum of the genus Brachybacterium the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Brachybacterium paraconglomeratum</i>	
<i>Brevundimonas abyssalis</i>	Species abyssalis / canariensis of the genus Brevundimonas have very similar patterns: Therefore distinguishing their species is difficult.
<i>Brochothrix campestris</i>	For the species campestris / thermosphacta of the genus Brochothrix the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Brochothrix thermosphacta</i>	
<i>Burkholderia ambifaria</i>	is a member of Burkholderia cepacia complex
<i>Burkholderia anthina</i>	
<i>Burkholderia cenocepacia</i>	



Genus/Species	Matching Hint
<i>Burkholderia cepacia</i>	is a member of Burkholderia cepacia complex
<i>Burkholderia diffusa</i>	
<i>Burkholderia dolosa</i>	
<i>Burkholderia lata</i>	
<i>Burkholderia latens</i>	
<i>Burkholderia metallica</i>	
<i>Burkholderia multivorans</i>	
<i>Burkholderia pyrrocinia</i>	
<i>Burkholderia seminalis</i>	
<i>Burkholderia stabilis</i>	
<i>Burkholderia thailandensis</i>	Burkholderia thailandensis is closely related and shows very similar spectra to the highly pathogenic Burkholderia pseudomallei / mallei which are possibly not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Burkholderia vietnamiensis</i>	is a member of Burkholderia cepacia complex
<i>Campylobacter helveticus</i>	closely related to Campylobacter upsaliensis
<i>Campylobacter upsaliensis</i>	closely related to Campylobacter helveticus
<i>Candida akabanensis</i>	For the species akabanensis / blattae/ dosseyi of the genus Candida the MALDI patterns and the ITS sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Candida blattae</i>	
<i>Citrobacter amalonaticus</i>	is a member of Citrobacter amalonaticus complex. Species amalonaticus / farmeri of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Citrobacter braakii</i>	is a member of Citrobacter freundii complex. Species braakii / freundii / gillenii / murlinae / rodentium / sedlakii / werkmannii / youngae of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Citrobacter farmeri</i>	is a member of Citrobacter amalonaticus complex. Species amalonaticus / farmeri of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Citrobacter freundii</i>	is a member of Citrobacter freundii complex. Species braakii / freundii / gillenii / murlinae / rodentium / sedlakii / werkmannii / youngae of the genus Citrobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Citrobacter gillenii</i>	
<i>Citrobacter murlinae</i>	
<i>Citrobacter rodentium</i>	
<i>Citrobacter sedlakii</i>	
<i>Citrobacter youngae</i>	
<i>Clostridioides difficile</i>	synonym of Clostridium difficile
<i>Clostridium beijerinckii</i>	Species beijerinckii / diolis of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Clostridium bolteae</i>	Species bolteae / clostridioforme of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Clostridium celerecrescens</i>	Species celerecrescens / sphenoides of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Clostridium clostridioforme</i>	Species bolteae / clostridioforme of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Clostridium diolis</i>	Species beijerinckii / diolis of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Clostridium haemolyticum</i>	Clostridium haemolyticum is closely related and shows very similar spectra to the strains of the highly pathogenic Clostridium botulinum groups C and D. Clostridium botulinum is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Clostridium pasteurianum</i>	For the species arbusti / pasteurianum of the genus Clostridium the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Clostridium sphenoides</i>	Species celerecrescens / sphenoides of the genus Clostridium have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Clostridium sporogenes</i>	Clostridium sporogenes is closely related and shows very similar spectra to the strains of the highly pathogenic Clostridium botulinum groups A, B and F. Clostridium botulinum is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Clostridium subterminale</i>	Clostridium subterminale is closely related and shows very similar spectra to the strains of the highly pathogenic Clostridium botulinum group G. Clostridium botulinum is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Corynebacterium glaucum</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium lipophile_group_F1</i>	
<i>Corynebacterium mucifaciens</i>	Species mucifaciens / ureicelerivorans of the genus Corynebacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium propinquum</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium pseudodiphtheriticum</i>	
<i>Corynebacterium pseudotuberculosis</i>	Species pseudotuberculosis / ulcerans of the genus Corynebacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium sp</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium ulcerans</i>	Species pseudotuberculosis / ulcerans of the genus Corynebacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Corynebacterium ureicelerivorans</i>	Species mucifaciens / ureicelerivorans of the genus Corynebacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Cronobacter sp</i>	Cronobacter can only be identified on genus level.
<i>Cryptococcus deneoformans</i>	synonym of Cryptococcus neoformans var. neoformans
<i>Cryptococcus neoformans</i>	synonym of Cryptococcus neoformans var. grubii
<i>Delftia acidovorans</i>	Species acidovorans / lacustris / tsuruhatensis of the genus Delftia are closely related.
<i>Delftia lacustris</i>	



Genus/Species	Matching Hint
<i>Dietzia cercidiphylli</i>	For the species cercidiphylli / natronolimnaea / psychralcaliphila of the genus Dietzia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Dietzia cinnamea</i>	For the species cinnamea / papillomatosis of the genus Dietzia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Dietzia natronolimnaea</i>	For the species cercidiphylli / natronolimnaea / psychralcaliphila of the genus Dietzia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Dietzia papillomatosis</i>	For the species cinnamea / papillomatosis of the genus Dietzia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Elizabethkingia anophelis</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Elizabethkingia meningoseptica</i>	
<i>Elizabethkingia miricola</i>	
<i>Enterobacter asburiae</i>	is a member of Enterobacter cloacae complex
<i>Enterobacter bugandensis</i>	
<i>Enterobacter cancerogenus</i>	
<i>Enterobacter cloacae</i>	
<i>Enterobacter hormaechei</i>	
<i>Enterobacter kobei</i>	
<i>Enterobacter ludwigii</i>	
<i>Escherichia coli</i>	closely related to Shigella / Escherichia fergusonii and not definitely distinguishable at the moment
<i>Escherichia fergusonii</i>	closely related to Shigella / Escherichia coli and not definitely distinguishable at the moment
<i>Fictibacillus arsenicus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Fictibacillus barbaricus</i>	
<i>Fictibacillus macauensis</i>	



Genus/Species	Matching Hint
<i>Fusobacterium gonidiaformans</i>	Species gonidiaformans / necrophorum of the genus Fusobacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Fusobacterium naviforme</i>	Species naviforme / nucleatum of the genus Fusobacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Fusobacterium necrophorum</i>	Species gonidiaformans / necrophorum of the genus Fusobacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Fusobacterium nucleatum</i>	Species naviforme / nucleatum of the genus Fusobacterium have very similar patterns: Therefore distinguishing their species is difficult.
<i>Fusobacterium</i> sp[2]	
<i>Globicatella sanguinis</i>	For the species sanguinis / sulfidifaciens of the genus Globicatella the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Globicatella sulfidifaciens</i>	
<i>Gordonia hongkongensis</i>	For the species hongkongensis / lacunae / terrae of the genus Gordonia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Gordonia terrae</i>	
<i>Halobacillus karajensis</i>	For the species dabanensis / karajensis / litoralis / profundi of the genus Halobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Hungatella effluvii</i>	For the species effluvii / hathewayi of the genus Hungatella the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Hungatella hathewayi</i>	
<i>Isoptricola variabilis</i>	For the species cucumis / muralis / variabilis of the genus Isoptricola the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Klebsiella aerogenes</i>	synonym of Enterobacter aerogenes
<i>Klebsiella oxytoca</i>	Klebsiella oxytoca and species ornithinolytica / planticola / terrigena of the genus Raoultella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Klebsiella pneumoniae</i>	closely related to Klebsiella variicola
<i>Klebsiella variicola</i>	closely related to Klebsiella pneumoniae
<i>Kluyvera ascorbata</i>	Species ascorbata / georgiana of the genus Kluyvera have very similar patterns: Therefore distinguishing their species is difficult.



Genus/Species	Matching Hint
<i>Kluyvera georgiana</i>	Species ascorbata / georgiana of the genus Kluyvera have very similar patterns: Therefore distinguishing their species is difficult.
<i>Kocuria himachalensis</i>	Species himachalensis / polaris / rosea of the genus Kocuria have very similar patterns: Therefore distinguishing their species is difficult.
<i>Kocuria polaris</i>	
<i>Kocuria rosea</i>	
<i>Lactobacillus amylovorus</i>	Species amylovorus / kitasatonis of the genus Lactobacillus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Lactobacillus cerevisiae</i>	For the species cerevisiae / koreensis of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Lactobacillus gallinarum</i>	For the species gallinarum / helveticus of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Lactobacillus gasseri</i>	For the species gasseri / taiwanensis of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Lactobacillus helveticus</i>	For the species gallinarum / helveticus of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Lactobacillus kitasatonis</i>	Species amylovorus / kitasatonis of the genus Lactobacillus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Lactobacillus similis</i>	For the species odoratitofui / similis of the genus Lactobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Lactococcus garvieae</i>	For the species garvieae / petauri of the genus Lactococcus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Leifsonia shinshuensis</i>	For the species poae / shinshuensis / soli of the genus Leifsonia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Leifsonia soli</i>	
<i>Listeria innocua</i>	Extraction must be performed to permit reliable species identification.
<i>Listeria ivanovii</i>	
<i>Listeria monocytogenes</i>	
<i>Listeria seeligeri</i>	



Genus/Species	Matching Hint
<i>Listeria welshimeri</i>	Extraction must be performed to permit reliable species identification.
<i>Lysinibacillus boronitolerans</i>	Species boronitolerans / xylanilyticus of the genus Lysinibacillus have very similar patterns: Therefore distinguishing their species is difficult. The quality of spectra (score) depends on the degree of sporulation: Use fresh material. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Lysinibacillus contaminans</i>	
<i>Lysinibacillus fusiformis</i>	
<i>Lysinibacillus halotolerans</i>	
<i>Lysinibacillus manganicus</i>	
<i>Lysinibacillus massiliensis</i>	
<i>Lysinibacillus meyeri</i>	
<i>Lysinibacillus pakistanensis</i>	
<i>Lysinibacillus sp</i>	
<i>Lysinibacillus sp[2]</i>	
<i>Lysinibacillus sphaericus</i>	
<i>Lysinibacillus xylanilyticus</i>	Species boronitolerans / xylanilyticus of the genus Lysinibacillus have very similar patterns: Therefore distinguishing their species is difficult. The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Macrococcus canis</i>	For the species canis / caseolyticus of the genus Macrocooccus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Macrococcus caseolyticus</i>	
<i>Mannheimia glucosida</i>	For the species glucosida / ruminalis of the genus Mannheimia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Mannheimia granulomatis</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Mannheimia ruminalis</i>	For the species glucosida / ruminalis of the genus Mannheimia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.



Genus/Species	Matching Hint
<i>Mannheimia varigena</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Marinococcus halophilus</i>	For the species halophilus / luteus / tarijensis of the genus Marinococcus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Massilia oculi</i>	Species oculi / timonae / varians of the genus Massilia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Massilia timonae</i>	
<i>Massilia varians</i>	
<i>Methylobacterium hispanicum</i>	For the species gregans / hispanicum of the genus Methylobacterium the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Methylobacterium tardum</i>	For the species longum / phyllostachyos / radiotolerans / tardum of the genus Methylobacterium the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Microbacterium kitamiense</i>	For the species aurantiacum / chocolatatum / kitamiense of the genus Microbacterium the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Mixta calida</i>	Species calida / gaviniae / intestinalis of the genus Mixta have very similar patterns: Therefore distinguishing their species is difficult.
<i>Mixta gaviniae</i>	
<i>Mixta intestinalis</i>	
<i>Moorella thermoacetica</i>	For the species thermoacetica / thermoautotrophica of the genus Moorella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Mycoplasma alkalescens</i>	Species alkalescens / arginini of the genus Mycoplasma have very similar patterns: Therefore distinguishing their species is difficult.
<i>Mycoplasma arginini</i>	
<i>Negativicoccus succinicivorans</i>	For the species massiliensis / succinicivorans of the genus Negativicoccus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Neisseria meningitidis</i>	Non-pathogenic Neisseria species could be misidentified as Neisseria meningitidis. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Nocardia abscessus</i>	For the species abscessus / arthritidis / asiatica / exalbida of the genus Nocardia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Nocardia arthritidis</i>	



Genus/Species	Matching Hint
<i>Nocardia asiatica</i>	For the species abscessus / arthritis / asiatica / exalbida of the genus Nocardia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Nocardia exalbida</i>	
<i>Nocardia farcinica</i>	Nocardia farcinica is/are closely related to the rare species Nocardia kroppenstedtii which is/are not included in the MALDI Biotyper database.
<i>Oceanobacillus oncorhynchi</i>	For the species aidingensis / oncorhynchi of the genus Oceanobacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Paenibacillus amylolyticus</i>	For the species amylolyticus / cucumis / taichungensis / tundrae / tylopili of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Paenibacillus antibioticophila</i>	For the species antibioticophila / apis of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Paenibacillus edaphicus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Paenibacillus provencensis</i>	For the species provencensis / shunpengii / urinalis of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Paenibacillus tylopili</i>	For the species amylolyticus / cucumis / taichungensis / tundrae / tylopili of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Paenibacillus urinalis</i>	For the species provencensis / shunpengii / urinalis of the genus Paenibacillus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Pantoea agglomerans</i>	synonym of Erwinia herbicola
<i>Pantoea anthophila</i>	For the species anthophila / deleyi / eucalypti of the genus Pantoea the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Pasteurella canis</i>	Species canis / oralis of the genus Pasteurella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Pasteurella oralis</i>	
<i>Peptoniphilus senegalensis</i>	For the species senegalensis / tyrrelliae of the genus Peptoniphilus the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Peptoniphilus tyrrelliae</i>	
<i>Proteus hauseri</i>	Species hauseri / penneri / vulgaris of the genus Proteus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Proteus penneri</i>	



Genus/Species	Matching Hint
<i>Proteus vulgaris</i>	Species hauseri / penneri / vulgaris of the genus Proteus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Pseudarcobacter aquimarinus</i>	Species aquimarinus / cloacae / defluvii / ellisii of the genus Pseudarcobacter have very similar patterns: Therefore distinguishing their species is difficult.
<i>Pseudarcobacter cloacae</i>	
<i>Pseudarcobacter defluvii</i>	
<i>Pseudarthrobacter chlorophenolicus</i>	Species tend to produce polymers which can interfere with the identification.
<i>Pseudescherichia vulneris</i>	synonym of Escherichia vulneris
<i>Pseudomonas azotoformans</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas balearica</i>	is a member of Pseudomonas stutzeri group
<i>Pseudomonas brenneri</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas cedrina</i>	
<i>Pseudomonas congelans</i>	
<i>Pseudomonas corrugata</i>	
<i>Pseudomonas extremorientalis</i>	
<i>Pseudomonas fluorescens</i>	
<i>Pseudomonas fluorescens_Group</i>	
<i>Pseudomonas fulva</i>	is a member of Pseudomonas putida group
<i>Pseudomonas gessardii</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas libanensis</i>	
<i>Pseudomonas luteola</i>	is a member of Pseudomonas stutzeri group
<i>Pseudomonas mandelii</i>	is a member of Pseudomonas fluorescens group



Genus/Species	Matching Hint
<i>Pseudomonas marginalis</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas migulae</i>	
<i>Pseudomonas monteilii</i>	is a member of Pseudomonas putida group
<i>Pseudomonas mosselii</i>	
<i>Pseudomonas mucidolens</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas orientalis</i>	
<i>Pseudomonas oryzae</i>	is a member of Pseudomonas putida group
<i>Pseudomonas plecoglossicida</i>	
<i>Pseudomonas poae</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas protegens</i>	
<i>Pseudomonas putida</i>	is a member of Pseudomonas putida group
<i>Pseudomonas putida_Group</i>	
<i>Pseudomonas rhodesiae</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas stutzeri</i>	is a member of Pseudomonas stutzeri group
<i>Pseudomonas synxantha</i>	is a member of Pseudomonas fluorescens group
<i>Pseudomonas tolaasii</i>	
<i>Pseudomonas trivialis</i>	
<i>Pseudomonas veronii</i>	
<i>Pseudomonas xanthomarina</i>	is a member of Pseudomonas stutzeri group
<i>Psychrobacillus insolitus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.



Genus/Species	Matching Hint
<i>Psychrobacillus psychrodurans</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Psychrobacillus psychrotolerans</i>	
<i>Raoultella ornithinolytica</i>	Klebsiella oxytoca and species ornithinolytica / planticola / terrigena of the genus Raoultella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Raoultella planticola</i>	
<i>Raoultella terrigena</i>	
<i>Rhodococcus hoagii</i>	synonym of Rhodococcus equi
<i>Rhodotorula glutinis</i>	For the species araucariae / glutinis / kratochvilovae of the genus Rhodotorula the ITS sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Rouxiella badensis</i>	For the species badensis / chamberiensis / silvae of the genus Rouxiella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Rouxiella chamberiensis</i>	
<i>Rummeliibacillus pycnus</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Salimicrobium halophilum</i>	
<i>Salmonella</i> sp	Salmonella can only be identified on genus level.
<i>Serratia ficaria</i>	For the species ficaria / vespertilionis of the genus Serratia the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Serratia grimesii</i>	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Serratia liquefaciens</i>	
<i>Serratia marcescens</i>	Species marcescens / nematodiphila / ureilytica of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Serratia nematodiphila</i>	
<i>Serratia plymuthica</i>	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Serratia proteamaculans</i>	



Genus/Species	Matching Hint
<i>Serratia quinivorans</i>	Species grimesii / liquefaciens / myotis / plymuthica / proteamaculans / quinivorans of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Serratia ureilytica</i>	Species marcescens / nematodiphila / ureilytica of the genus Serratia have very similar patterns: Therefore distinguishing their species is difficult.
<i>Solibacillus silvestris</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Sphingomonas pseudosanguinis</i>	Species pseudosanguinis / sanguinis of the genus Sphingomonas have very similar patterns: Therefore distinguishing their species is difficult.
<i>Sphingomonas sanguinis</i>	
<i>Staphylococcus carnosus</i>	For the species carnosus / condimenti / piscifermentans of the genus Staphylococcus the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Staphylococcus condimenti</i>	
<i>Staphylococcus piscifermentans</i>	
<i>Stenotrophomonas maltophilia</i>	is a member of Stenotrophomonas maltophilia group or closely related
<i>Streptococcus canis</i>	Species canis / dysgalactiae / pyogenes of the genus Streptococcus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Streptococcus dysgalactiae</i>	
<i>Streptococcus equinus</i>	Species equinus / gallolyticus / infantarius / lutetiensis of the genus Streptococcus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Streptococcus gallolyticus</i>	
<i>Streptococcus infantarius</i>	
<i>Streptococcus lutetiensis</i>	
<i>Streptococcus mitis</i>	Streptococcus mitis / oralis / peroris / pneumoniae / pseudopneumoniae are closely related! The result may be confirmed by a further test, e.g. bile test or optochin test, according to standard clinical microbiological practice.
<i>Streptococcus oralis</i>	
<i>Streptococcus peroris</i>	
<i>Streptococcus pneumoniae</i>	



Genus/Species	Matching Hint
<i>Streptococcus pseudopneumoniae</i>	Streptococcus mitis / oralis / peroris / pneumoniae / pseudopneumoniae are closely related! The result may be confirmed by a further test, e.g. bile test or optochin test, according to standard clinical microbiological practice.
<i>Streptococcus pyogenes</i>	Species canis / dysgalactiae / pyogenes of the genus Streptococcus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Tsukamurella hongkongensis</i>	Species of this genus have very similar patterns: Therefore distinguishing their species is difficult.
<i>Tsukamurella inchonensis</i>	
<i>Tsukamurella paurometabola</i>	
<i>Tsukamurella pseudospumae</i>	
<i>Tsukamurella pulmonis</i>	
<i>Tsukamurella serpentis</i>	
<i>Tsukamurella sinensis</i>	
<i>Tsukamurella soli</i>	
<i>Tsukamurella spumae</i>	
<i>Tsukamurella strandjordii</i>	
<i>Tsukamurella tyrosinosolvens</i>	
<i>Ureibacillus suwonensis</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Ureibacillus thermosphaericus</i>	
<i>Veillonella atypica</i>	Species atypica / caviae / denticariosi / dispar / parvula / rogosae of the genus Veillonella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Veillonella caviae</i>	
<i>Veillonella denticariosi</i>	
<i>Veillonella dispar</i>	



Genus/Species	Matching Hint
<i>Veillonella parvula</i>	Species atypica / caviae / denticariosi / dispar / parvula / rogosae of the genus Veillonella have very similar patterns: Therefore distinguishing their species is difficult.
<i>Veillonella rogosae</i>	
<i>Vibrio albensis</i>	Vibrio albensis (V. cholerae biovar albensis) is closely related and shows very similar spectra to the highly pathogenic Vibrio cholerae which is possibly not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Vibrio alginolyticus</i>	is a member of Vibrio harveyi group
<i>Vibrio campbellii</i>	
<i>Vibrio harveyi</i>	
<i>Vibrio mytili</i>	
<i>Vibrio natriegens</i>	
<i>Vibrio navarrensis</i>	Species navarrensis / vulnificus of the genus Vibrio have very similar patterns: Therefore distinguishing their species is difficult.
<i>Vibrio parahaemolyticus</i>	is a member of Vibrio harveyi group
<i>Vibrio rotiferianus</i>	
<i>Vibrio vulnificus</i>	Species navarrensis / vulnificus of the genus Vibrio have very similar patterns: Therefore distinguishing their species is difficult.
<i>Viridibacillus arenosi</i>	The quality of spectra (score) depends on the degree of sporulation: Use fresh material.
<i>Viridibacillus arvi</i>	
<i>Viridibacillus neidei</i>	
<i>Weissella fabaria</i>	For the species fabalis / fabaria / ghanensis of the genus Weissella the MALDI patterns and the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.
<i>Weissella ghanensis</i>	
<i>Weissella thailandensis</i>	For the species jogaejeotgali / thailandensis of the genus Weissella the 16S rRNA gene sequences are very similar. Therefore distinguishing the mentioned species is difficult.



Genus/Species	Matching Hint
<i>Yersinia frederiksenii</i>	For the species frederiksenii / intermedia / kristensenii of the genus Yersinia it should be considered that they could exhibit pathogenic characteristics.
<i>Yersinia intermedia</i>	
<i>Yersinia kristensenii</i>	
<i>Yersinia pseudotuberculosis</i>	Yersinia pseudotuberculosis is closely related and shows very similar spectra to the highly pathogenic Yersinia pestis which is not included in the MALDI Biotyper database. These species cannot be distinguished by the MALDI Biotyper method. For differentiation an adequate identification method has to be selected by an experienced professional.
<i>Zygoascus meyeriae</i>	For the species hellenicus / meyeriae of the genus Zygoascus the ITS sequences are very similar. Therefore distinguishing the mentioned species is difficult.



8.2 Information only required for users of laboratory-specific specimen codes

Please refer to the tables provided above regarding the genus/species renaming of individual MSP or the general genus/species renaming.

Note: Users of laboratory-specific specimen codes may use this information to update their code lists. These lists can then be used as basis to update dedicated code list used, for example, in LIMS integration. For further details, please contact your LIMS provider or contact MALDI Biotyper Software support at the following Email address: biotyper.sw.support@bruker.com.

MBT 3.1/4.0 software

The update of appropriate code lists can be performed using the Biotyper Abbreviation Editor. You will find *MspAbbreviationEditor.exe* on the hard disk in the same folder as the MALDI Biotyper RTC program. Default path:

C:\Program Files (x86)\Bruker Daltonik\BiotyperRTC

Run the *MspAbbreviationEditor.exe* and a window will open. Select **Group by Species** and adapt the code list accordingly.

Note: The process of validation of these codes is the user's responsibility. It is recommended that validation must be repeated at reasonable intervals.

MBT Compass software

The update of code lists can be performed by selecting the sprocket-icon in the top right corner. Click the **Advanced** button at the bottom of the configuration window, enter the credentials of a Windows-administrator and select **MSP Abbr.** On the next open page click **Start edit** to run the Abbreviation Editor. A window will open. Select **Group by Species** and adapt the code list accordingly.

Note: The process of validation of these codes is the user's responsibility. It is recommended that validation must be repeated at reasonable intervals.

MBT Compass HT software

The MBT Compass HT software allows to import, export and display laboratory specific species code list:

- The appropriate functionality is available to the Laboratory Manager at this menu **Configuration > Msp abbreviations**.
- More than one code lists are allowed – new (initially empty) code lists are created (and named) using the **Add** button.
- New code lists are filled from a simple CSV file using the **Import** button (CSV row: *Laboratory,Abbreviation,Name*).

Release Notes

MBT Compass Library

Revision L (2020)



- Present code lists may be exported into a CSV file using the **Export** button (same CSV format).
- The menu displays a sortable list of all currently present lab specific code tabled (Species | Abbreviation).

Note: The process of validation of these codes is the user's responsibility. It is recommended that validation must be repeated at reasonable intervals.

8.3 Additional information for yeast nomenclature update

Specialist statement by Dr. Andrey Yurkov (Curator Fungi and Yeasts at Leibniz Institute – DSMZ; Braunschweig, Germany)

Reclassification of yeasts in the Bruker Yeast Library

Recent changes in the nomenclature of fungi (10, 13) discontinued the use of dual nomenclature, different names for same species of sexual and asexual fungi, and adopted the 'One fungus, one name' principle (10). These new rules have a major impact on the classification of many asexual species in large polyphyletic genera *Candida*, *Cryptococcus*, *Pseudozyma*, *Rhodotorula*, *Tilletiopsis* and *Trichosporon*. Subsequent analyses showed that certain species may be assigned to already existing sexual genera with high confidence using multigene phylogenies (e.g. 3, 6, 9, 14, Species that form well-circumscribed phylogenetic clades without any teleomorph member justified the creation of new genera (e.g. 2, 4, 7).

Genera of basidiomycetous yeasts, including genera *Cryptococcus*, *Pseudozyma*, *Rhodotorula*, *Tilletiopsis* and *Trichosporon*, have been reclassified recently resulting in 46 new genera and over 300 new combinations (8, 14, 16). Reclassification of the genus *Candida* is underway will continue in the future so that the genus *Candida* will be restricted to the monophyletic group that includes clinically relevant species *Candida albicans* and *Candida tropicalis* (1). As the result of the ongoing classification, many yeast names have been already changed, and these changes are now implemented in the database. Although the introduction of many new and unfamiliar names may be irritating for users of the database, these changes are extremely important. Whether identified from clinical or environmental sources, classification of yeasts in large genera *Candida* may lead to a wrong assumption that all these species are closely related (1, 11). It has been demonstrated that antifungal resistance profiles of yeasts are correlated with their phylogenetic position in clades and families Debaryomycetaceae (e.g. *Candida albicans*), Metschnikowiaceae (e.g. *Candida auris*), Pichiaceae (e.g. *Candida krusei*), and Saccharomycetaceae (e.g. *Candida glabrata*). Thus, separation distantly related species in new genera will give non-specialists a clear signal that these species have different properties, e.g. antifungal susceptibility (12).

The reclassification of yeast names in the database was necessary because of many reclassifications made by taxonomists during the last 10 years. The nomenclature is following the latest published edition of the taxonomic compendium *The Yeasts: A Taxonomic Study* (5) and several later publications on yeast taxonomy. High-ranking classification of yeasts reflects phylogenetic relationships inferred from multi-gene and whole-genome sequence analyses (1, 3, 9, 11, 14, 15).

Out of 212 species in the Bruker database, 109 names remained unchanged. A total of 81 species were transferred into a new genus; 21 species received new genus and new species names; one species was assigned to another existing species in the same genus. In the genus *Candida*, 32 species that await reclassification have been provisionally assigned to a genus or family.



References

1. Daniel, H. M., Lachance, M. A., & Kurtzman, C. P. (2014). On the reclassification of species assigned to *Candida* and other anamorphic ascomycetous yeast genera based on phylogenetic circumscription. *Antonie van Leeuwenhoek*, 106(1), 67-84.
2. Khunnamwong, P., Lertwattanasakul, N., Jindamorakot, S., Limtong, S., & Lachance, M. A. (2015). Description of *Diutina* gen. nov., *Diutina siamensis*, f.a. sp. nov., and reassignment of *Candida catenulata*, *Candida mesorugosa*, *Candida neorugosa*, *Candida pseudorugosa*, *Candida ranongensis*, *Candida rugosa* and *Candida scorzettiae* to the genus *Diutina*. *International Journal of Systematic and Evolutionary Microbiology*, 65(12), 4701-4709.
3. Kurtzman, C. P., & Robnett, C. J. (2013). Relationships among genera of the Saccharomycotina (Ascomycota) from multigene phylogenetic analysis of type species. *FEMS Yeast Research*, 13(1), 23-33.
4. Kurtzman, C. P., & Robnett, C. J. (2014) Three new anascosporic genera of the Saccharomycotina: *Danielozyma* gen. nov., *Deakozyma* gen. nov. and *Middelhovenomyces* gen. nov. *Antonie van Leeuwenhoek*, 105(5), 933-942.
5. Kurtzman, C. P., Fell, J. W., & Boekhout, T. (Eds.). (2011). *The Yeasts: A Taxonomic Study*. Elsevier, Amsterdam, The Netherlands.
6. Kurtzman, C. P., Robnett, C. J., & Basehoar-Powers, E. (2008). Phylogenetic relationships among species of *Pichia*, *Issatchenkia* and *Williopsis* determined from multigene sequence analysis, and the proposal of *Barnettozyma* gen. nov., *Lindnera* gen. nov. and *Wickerhamomyces* gen. nov. *FEMS Yeast Research*, 8(6), 939-954.
7. Lachance, M. A., & Kurtzman, C. P. (2013). The yeast genus *Tortispora* gen. nov., description of *Tortispora ganteri* sp. nov., *Tortispora mauiana* f.a., sp. nov., *Tortispora agaves* f.a., sp. nov., *Tortispora sangerardonensis* f.a., sp. nov., *Tortispora cuajiniquilana* f.a., sp. nov., *Tortispora starmeri* f.a., sp. nov. and *Tortispora phaffii* f.a., sp. nov., reassignment of *Candida caseinolytica* to *Tortispora caseinolytica* f.a., comb. nov., emendation of *Botryozyma*, and assignment of *Botryozyma*, *Tortispora* gen. nov. and *Trigonopsis* to the family Trigonopsidaceae fam. nov. *International Journal of Systematic and Evolutionary Microbiology*, 63(8), 3104-3114.
8. Liu, X. Z., Wang, Q. M., Göker, M., Groenewald, M., Kachalkin, A. V., Lumbsch, H. T., Millanes, A. M., Wedin, M., Yurkov, A. M., Boekhout, T., & Bai, F. Y. (2015). Towards an integrated phylogenetic classification of the Tremellomycetes. *Studies in Mycology*, 81, 85-147.
9. Liu, X. Z., Wang, Q. M., Theelen, B., Groenewald, M., Bai, F. Y., & Boekhout, T. (2015). Phylogeny of tremellomycetous yeasts and related dimorphic and filamentous basidiomycetes reconstructed from multiple gene sequence analyses. *Studies in Mycology*, 81, 1-26.
10. McNeill, J., Barrie, F. R., Buck, W. R., Demoulin, V., Greuter, W., Hawksworth, D. L., Herendeen, P. S., Knapp, S., Marhold, K., Prado, J., Prud'homme Van Reine, W. F., Smith, G. F., & Wiersema, J. H. (Eds.) (2012). *International Code of Nomenclature for Algae, Fungi, and Plants (Melbourne Code)*. Regnum Vegetabile 154. Koelz Scientific Books, Koenigstein, Germany.
11. Shen, X.X., Oplente, D.A., Kominek, J., Zhou, X., Steenwyk, J.L., Buh, K.V., Haase, M.A., Wisecaver, J.H., Wang, M., Doering, D.T., & Boudouris, J. T. (2018). Tempo and mode of genome evolution in the budding yeast subphylum. *Cell*, 175(6), 1533-1545.



12. Stavrou, A. A., Lackner, M., Lass-Flörl, C., & Boekhout, T. (2019). The changing spectrum of Saccharomycotina yeasts causing candidemia: phylogeny mirrors antifungal susceptibility patterns for azole drugs and amphotericin B. *FEMS Yeast Research*, 19(4), foz037.
13. Turland, N. J., Wiersema, J. H., Barrie, F. R., Greuter, W., Hawksworth, D. L., Herendeen, P. S., Knapp, S., Kusber, W.-H., Li, D.-Z., Marhold, K., May, T. W., McNeill, J., Monro, A. M., Prado, J., Price, M. J. & Smith, G. F. (Eds.) (2018). International Code of Nomenclature for algae, fungi, and plants (Shenzhen Code) adopted by the Nineteenth International Botanical Congress Shenzhen, China, July 2017. *Regnum Vegetabile* 159. Koeltz Botanical Books, Glashütten, Germany.
14. Wang, Q. M., Begerow, D., Groenewald, M., Liu, X. Z., Theelen, B., Bai, F. Y., & Boekhout, T. (2015). Multigene phylogeny and taxonomic revision of yeasts and related fungi in the Ustilaginomycotina. *Studies in Mycology*, 81, 55-83.
15. Wang, Q. M., Groenewald, M., Takashima, M., Theelen, B., Han, P. J., Liu, X. Z., Boekhout, T., & Bai, F. Y. (2015). Phylogeny of yeasts and related filamentous fungi within Pucciniomycotina determined from multigene sequence analyses. *Studies in Mycology*, 81, 27-53.
16. Wang, Q. M., Yurkov, A. M., Göker, M., Lumbsch, H. T., Leavitt, S. D., Groenewald, M., Theelen, B., Liu, X. Z., Boekhout, T., & Bai, F. Y. (2015). Phylogenetic classification of yeasts and related taxa within Pucciniomycotina. *Studies in Mycology*, 81, 149-189.

8.4 Additional information for *Staphylococcus argenteus* / *schweitzeri*

The species *Staphylococcus argenteus* and *Staphylococcus schweitzeri* form parts of a "*Staphylococcus aureus*-related complex" (19, 20) and show similar virulence potential compared to *Staphylococcus aureus* (17, 18). Therefore, they should be further evaluated.


References

17. Hansen TA, Bartels MD, Høgh SV, Dons LE, Pedersen M, et al. Whole Genome Sequencing of Danish *Staphylococcus argenteus* Reveals a Genetically Diverse Collection with Clear Separation from *Staphylococcus aureus*. *Frontiers in Microbiology*;8. Epub ahead of print 9 August 2017. DOI: 10.3389/fmicb.2017.01512.
18. Kaden R, Engstrand L, Rautelin H, Johansson C. Which methods are appropriate for the detection of *Staphylococcus argenteus* and is it worthwhile to distinguish *S. argenteus* from *S. aureus*? *Infection and Drug Resistance* 2018;Volume 11:2335–2344.
19. Schuster D, Rickmeyer J, Gajdiss M, Thye T, Lorenzen S, et al. Differentiation of *Staphylococcus argenteus* (formerly: *Staphylococcus aureus* clonal complex 75) by mass spectrometry from *S. aureus* using the first strain isolated from a wild African great ape. *International Journal of Medical Microbiology* 2017;307:57–63.
20. Tong SYC, Schaumburg F, Ellington MJ, Corander J, Pichon B, et al. Novel staphylococcal species that form part of a *Staphylococcus aureus*-related complex: the non-pigmented *Staphylococcus argenteus* sp. nov. and the non-human primate-associated *Staphylococcus schweitzeri* sp. nov. *International Journal of Systematic and Evolutionary Microbiology* 2015;65:15–22.



Bruker Daltonik GmbH
Fahrenheitstraße 4
28359 Bremen
Germany

Email: biotyper.support@bruker.com
Phone: +49 (421) 2205-1401
Fax: +49 (421) 2205-106

Symbol	Description
	Manufacturer
REF	Catalogue number

*Descriptions and specifications supersede all previous information and are subject to change without notice.
Copyright © 2020, Bruker Daltonik GmbH, November 2020*