

BactoReal[®] Kit Brachyspira hyodysenteriae



For veterinary use only

BactoReal® Kit <i>Brachyspira hyodysenteriae</i>			
Order no.	Reactions	Pathogen	Internal positive control
DVEB01113	100	FAM channel	Cy5 channel
DVEB01153	50	FAM channel	Cy5 channel
DVEB01111	100	FAM channel	VIC/HEX channel
DVEB01151	50	FAM channel	VIC/HEX channel

Kit contents:

- Detection assay for B. hyodysenteriae
- Detection assay for internal positive control (control of amplification)
- DNA reaction mix (contains uracil-N glycosylase, UNG)
- Positive control for B. hyodysenteriae
- Water



Background: Brachyspira hyodysenteriae (previously called Serpulina hyodysenteriae or Treponema hyodysenteriae) is an anaerobic intestinal spirochaete and the causative agent of swine dysentery. Swine dysentery is a disease characterized by mucohaemorrhagic diarrhoea with lesions confined to the large intestine of pigs. Its incubation time is 10-14 days. The faecal excretion of B. hyodysenteriae starts 2 days after infection. Brachyspira hyodysenteriae is an anaerobic bacterium but is aerotolerant due, at least in part, to high NADH oxidase activity.

Description: BactoReal® Kit *Brachyspira hyodysenteriae* is based on the amplification and detection of the nox gene of *B. hyodysenteriae* using real-time PCR. It allows the rapid and sensitive detection of the nox gene of *B. hyodysenteriae* from DNA samples purified from faecal samples, or biopsies of the intestinal mucosa (e.g. with the QIAamp DNA Stool Mini Kit or the QIAamp DNA Mini Kit, respectively).


PCR-platforms: BactoReal® Kit *Brachyspira hyodysenteriae* is developed and validated for the ABI PRISM® 7500 instrument (Life Technologies), LightCycler® 480 (Roche) and Mx3005P® QPCR System (Agilent), but is also suitable for other real-time PCR instruments.

Sensitivity and specificity: BactoReal® Kit *Brachyspira hyodysenteriae* has an analytical sensitivity of 10 target copies/reaction. The limit of detection (LoD95 = smallest number of copies of target DNA which can be detected in 95% of cases) of 18 target copies/reaction was determined by several replicates around the detection limit. The kit is specific for *B. hyodysenteriae*. Specificity was tested on isolates of *B. pilosicoli, B. innocens, B. murdochii, E. coli, H. parasuis, L. intracellularis, L. innocua, L. monocytogenes, P. multocida, S. aureus, S. agalactiae and S. pyogenes.* No cross reactions were observed. A total of 31 field samples were tested and correctly analysed.

References: Atyeo, R. F., T. B. Stanton, N. S. Jensen, D. S. Suriyaarachichi, and D. J. Hampson. 1999. Differentiation of *Serpulina* species by NADH oxidase gene (*nox*) sequence comparisons and nox-based polymerase chain reaction tests. Vet. Microbiol. 67:47–60.

Product Description





BactoReal®, MycoReal, ParoReal and ViroReal® Kits run with the same thermal cycling conditions.

RNA and DNA material can be analysed in one PCR run.

For further information on our products please visit our homepage (www.ingenetix.com)