Mouse anti-CD30 (Ki-1), Clone Ber-H2

For In Vitro Diagnostic Use

Lot No.

- [] 08-0155 6.0 mL Predilute Antibody, Ready-To-Use
- [] 08-1155 6.0 mL 2ndGen Predilute Antibody, Ready-To-Use
- [] 18-0155 1.0 mL Concentrate Antibody
- [X] 18-7155 0.5 mL Concentrate Antibody

INTENDED USE

For In Vitro Diagnostic Use

Invitrogen's monoclonal Mouse anti-CD30 antibody (clone: Ber-H2) is intended to qualitatively stain a subset of B cells and T cells or certain abnormal lymphocytes in frozen tissue or formalin-fixed, paraffin-embedded tissue sections. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

BACKGROUND

CD30 (Ki-1)^{1,2} is a transmembrane glycoprotein (120 kDa non-reduced, 105 kDa reduced). CD30 is a cytokine receptor which belongs to the tumor necrosis factor receptor/nerve growth factor receptor (TNFR/NGFR) superfamily.³ A soluble form of CD30 (85 kDa) is produced by metalloproteinase.

In normal tissue, CD30 is found in activated B and T lymphocytes, activated NK cells, and monocytes. Reactivity in normal lymphoid tissues is localized to scattered, large lymphoid cells around lymph follicles and at the margin of germinal centers. Ber-H2 does not label tissues of non-lymphoid origin, except for some diffuse staining of the cytoplasm of exocrine pancreatic cells and a portion of cerebral cortical neurons and Purkinje cells of the cerebellum.^{2,4}

In abnormal tissue, CD30 is expressed by Hodgkins and Reed-Sternberg cells in Hodgkins disease, anaplastic large cell lymphomas (ALCL), embryonic carcinomas, and mixed germ cell tumors.^{2,4}

REAGENT PROVIDED

Mouse anti-CD30 (Clone: Ber-H2) is purified from mouse ascites and diluted in phosphate buffered saline (PBS), pH 7.4, and 1% bovine serum albumin (BSA) with 0.1% sodium azide (NaN₃) as a preservative.

<u>Immunogen:</u> CO cell line cells	<u>Total protein concentration:</u>	g/L
<u>Clone:</u> Ber-H2 ⁴	<u>Isotype:</u> IgG ₁	<u>Antibody concentration:</u>
		mg/L

STORAGE: 2-8°C

POSITIVE CONTROL TISSUE: Tonsil or lymph node

EXPECTED STAINING PATTERN: Primarily membrane with limited cytoplasmic

INSTRUCTIONS FOR USE

PRETREATMENT REQUIREMENTS:

- Epitope Retrieval: Required (EDTA, pH 8) (See page 2 for protocol)
- Enzyme Digestion: Not required

Mouse anti-CD30, (Clone: Ber-H2) may be diluted according to Table 1 when using the Invitrogen detection systems below.

Table 1. Dilution Table

Invitrogen Kit	Predilute Antibody	Dilution for Concentrate	Incubation Time
Histostain-SP or SAP kits*	Ready-To-Use	1: 50 - 1: 100	60 min.
Histostain [®] -Plus Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.
SuperPicTure [™] Polymer Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.
Cap-Plus [™] Kits	Ready-To-Use	1: 100 - 1: 200	30-60 min.

* Use Histostain-SP or SAP kits only for Cat. No. 08-0XXX and 18-X001 to 18-X200 primary antibodies.

This is a guideline only. Optimal antibody concentrations may vary based on specimen and preparation method used, and should be determined by each individual laboratory.

SPECIMEN PREPARATION

1. Use tissue fixed in 10% neutral buffered formalin or other fixative on regular basis, or frozen tissue sections.
2. Cut 3-4 μm sections and place on positively charged slides.
3. Dry overnight at 37° C or for 2-4 hrs at 58°C.

PRETREATMENT

Heat Induced Epitope Retrieval (HIER), if required

1. Deparaffinize slides. Tissue sections should be mounted on silane, poly-L-Lysine, or HistoGrip (Cat. No. 00-8050) coated slides.
2. Wash slides with distilled water 3 times for 2 min each.
3. Place a 1L glass (Pyrex) beaker containing 500 mL of 0.01 M citrate buffer (Cat. No. 00-5000) or EDTA solution (Cat. No. 00-5500) on a hot plate. Heat the buffer solution until it boils. (*This step may be prepared before slide deparaffinization, as the buffer may take several minutes to boil*).
4. Put the slides in a slide rack and place in the beaker with boiling buffer. Keep it boiling for 15 minutes.
5. After heating, remove beaker from the hot plate and allow it to cool down for at least 15-20 minutes at room temperature.
6. Rinse slides with PBS (Cat. No. 00-3000) and begin the immunostaining protocol.

Enzyme Digestion, if required

1. Prewarm the enzyme of choice at 37°C for 10 min.
2. Add the prewarmed enzyme to a tissue section and incubate at 37°C for 10 min.
3. Wash in several changes of PBS (Cat. No. 00-3000) and begin the immunostaining protocol.

RECOMMENDED MANUAL STAINING PROCEDURE

1. Submerge slides in peroxidase quenching solution and rinse with PBS.
2. Apply serum blocking solution.
3. Apply primary antibody and incubate for 30-60 min at room temperature; rinse with PBS.
4. Apply secondary antibody and incubate for 10 min at room temperature; rinse with PBS.
5. Apply enzyme conjugate and incubate for 10 min at room temperature; rinse with PBS.
6. Apply chromogen and incubate for 5-10 min at room temperature; rinse with PBS.

MATERIALS REQUIRED BUT NOT PROVIDED

Reagent	Catalog No.
1. HistoGrip™	00-8050
2. Super PAP Pen	00-8899
3. Isotype Control for Rabbit or Mouse Primary Antibody	08-6199 or 08-6599
4. Antibody Diluent	00-3118
5. PBS (0.01 M PBS)	00-3000
6. Mayer's Hematoxylin	00-8011
7. Citrate Buffer pH 6.0 (if required for HIER)	00-5000
8. EDTA Solution (if required for HIER)	00-5500
9. Digest-All™ 1, Digest-All™ 2, or Digest-All™ 3 (if required for Enzyme Digestion)	00-3007 or 00-3008 or 00-3009
10. SuperPicTure™ polymer kit, or LAB-SA kit (Histostain®-Plus, and Cap-Plus™) (see Table 2).	
11. Chromogen/substrate (if not included in detection kit): Single Solution AEC (Cat. No. 00-1111), or DAB (Cat. No. 00-2014), or Fast-Red (Cat. No. 00-2234).	
12. Mounting solution: Histomount™ (for DAB: Cat. No.00-8030), GVA (for AEC, or Fast-Red: Cat. No. 00-8000), or Clearmount™ (for AEC, DAB, or Fast-Red: Cat. No. 00-8010).	

REFERENCES

1. Falini B, et al. *Blood* 85(1): 1-14, 1995.
2. Chiarle R, et al. *Clin Immunol* 90(2): 157-64, 1999.
3. Durkop H. *Cell* 68(3): 421-7, 1992.
4. Schwarting R, et al. *Blood* 74(5): 1678-89, 1989.

TRADEMARKS

Cap-Plus™, Clearmount™, Digest-All™, HistoGrip™, Histomount™, and Histostain®, SuperPicTure™, and Zymed® are trademarks of Zymed Laboratories, Inc.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PI 18-7155

(Rev 08/08) DCC-08-1089

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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