VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

H. ducreyi + C. trachomatis (LGV)

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For Research Use Only (RUO)

This product has no declared clinical intended purpose and is not for clinical diagnostic use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

ENGLISH

1. Intended use

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Haemophilus ducreyi* and/or LGV-associated strains of *Chlamydia trachomatis* in urogenital specimens, urine samples and rectal samples from patients with signs and symptoms of *Haemophilus ducreyi* infection and/or lymphogranuloma venereum (LGV). This test is intended to be used for research purposes, without any medical objective is not regarded as devices for performance evaluation. DNA is extracted from specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific for *Haemophilus ducreyi* and *Chlamydia trachomatis* (LGV strains).

2. Summary and Explanation

Haemophilus ducreyi, a gram-negative bacterium, is the causative agent of chancroid, a genital ulcer disease (GUD). The organism is usually spread during sexual intercourse through microabrasions, and the disease usually manifests as multiple painful superficial ulcers associated with inguinal lymphadenitis. As a result of the painful nature of the lesions, patients usually seek immediate treatment, and asymptomatic carriage is therefore uncommon. In addition to causing GUD, *H. ducreyi* has been found in several recent studies to be a major cause of chronic skin ulceration in children from developing countries.

Chancroid or soft chancre is a sexually transmitted disease (STD) of humans, characterized by painful genital ulcers. The primary ulcer is often followed by multiple lesions, which can extend the duration of the illness to 1-3 months if appropriate treatment is not given. Traditionally, chancroid has been difficult to diagnose clinically, because similar symptoms can occur in infections with *Treponema pallidum* (syphilis) and Herpes simplex virus. Laboratory diagnoses can be equally frustrating, as *Haemophilus ducreyi*, the causative agent of chancroid, does not survive well in most transport media, has complex nutritional requirements for growth, and once isolated is poorly reactive biochemically.

Chlamydia trachomatis (Ct) is a gram-negative obligate intracellular pathogen responsible for over 110 million global sexually transmitted disease (STD) cases annually, with sequelae including pelvic inflammatory disease (PID), ectopic pregnancy, and infertility. Ct is also responsible for outbreaks of lymphogranuloma venereum (LGV) and trachoma, a chronic ocular disease that can lead to blindness.

LGV is caused by specific serovars of *Chlamydia trachomatis* (L1, L2 and L3) and these strains are associated with a more chronic and invasive infection than other serovars. Although LGV symptoms can vary according to site of entry and stage of infection, genital ulceration and inguinal lymphadenopathy are the classical presentations of this disease. The laboratory identification of LGV can be problematic, as routine culture of *C. trachomatis* for diagnostic purposes has largely been replaced with nucleic acid amplification techniques.

The detection method applied at the beginning was based on the culture, which takes a long time to complete and many difficulties. To shorten detection time and improve sensitivity, real-time PCR assays have proven to be a tool for the detection of *Haemophilus ducreyi* and *Chlamydia trachomatis* (LGV).

3. Principle of the procedure

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Haemophilus ducreyi* and/or LGV-associated strains of *Chlamydia trachomatis* in urogenital specimens, urine samples and rectal samples from patients with signs and symptoms of *Haemophilus ducreyi* infection and/or lymphogranuloma venereum (LGV). After DNA isolation, the identification of *Haemophilus ducreyi* and LGV-associated strains of *Chlamydia trachomatis* is performed by the amplification of a conserved region of the 16S rRNA gene for *Haemophilus ducreyi* and *pmpH* gene for LGV-associated *C. trachomatis* serovars, using specific primers and a fluorescent-labelled probe.

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in an stabilized format, as well as an internal control to monitor PCR inhibition. *Chlamydia trachomatis* (LGV strains) DNA targets are amplified and detected in the FAM channel, *Haemophilus ducreyi* DNA targets are amplified and detected in the ROX channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1 and Table 2:

Reference	Reagent/Material	Description	Colour	Amount
VS-HCT1SL/ VS-HCT1SH	<i>H. ducreyi</i> + <i>C. trachomatis</i> (LGV) 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-HCT1C	<i>H. ducreyi</i> + <i>C. trachomatis</i> (LGV) Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit with Ref. VS-HCT106L, VS-HCT106H, VS-HCT112L and VS-HCT112H.

Reference	Reagent/Material	Description	Color	Amount
VS-HCT1PL/ VS-HCT1PH	<i>H. ducreyi</i> + <i>C. trachomatis</i> (LGV) 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-HCT1C	<i>H. ducreyi</i> + <i>C. trachomatis</i> (LGV) Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNAse/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit with Ref VS-HCT113L and VS-HCT113H.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DPrime Real-time Detection Thermal Cycler, DNA-Technology DLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.

- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- This VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection kit is for Research Use Only.
- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-HCT113L and VS-HCT113H). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.

8. Test procedure

8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from urogenital specimens, urine samples and rectal samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.

- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- US121 NX-48 Urine/Swab DNA Kit using the Nexttractor® NX-48 (Genolution Inc.).
- MagDEA® Dx SV, using magLEAD® 6gC (Precision System Science)).
- Invisorb® Spin Universal Kit (Strattec).

8.2. Lyophilized positive control

H. ducreyi + *C. trachomatis* (LGV) Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *H. ducreyi* + *C. trachomatis* (LGV) Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *H. ducreyi* + *C. trachomatis* (LGV) Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up your thermocycler (to check compatibility see Annex 1).

Program your thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*C. trachomatis* (LGV)), ROX (*H. ducreyi*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection

channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System, StepOne Plus™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *H. ducreyi* + *C. trachomatis* (LGV) in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>C. trachomatis</i> (LGV) (FAM)	<i>H. ducreyi</i> (ROX)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	<i>C. trachomatis</i> (LGV) and <i>H. ducreyi</i> Positives
-	-	+	-	+	<i>C. trachomatis</i> (LGV) and <i>H. ducreyi</i> Negative
+	-	+/-	-	+	<i>C. trachomatis</i> (LGV) Positive and <i>H. ducreyi</i> Negative
-	+	+/-	-	+	<i>H. ducreyi</i> Positive and <i>C. trachomatis</i> (LGV) Negative
-	-	-	-	+	Experiment fail
+	+	+	+	-	Experiment fail

Table 4. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

10. Limitations of the test

- The test must be used for RUO purposes (Exclusive Use in Research), without any medical objective.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from urogenital specimens (included endocervical specimens and urethral samples), urine and rectal samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *H. ducreyi* and *C. trachomatis* (LGV) either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Preliminary results of the test characteristics

12.1. Clinical sensitivity and specificity

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit was evaluated using 4 different QCMD panels (2 from the QCMD 2017 *Chlamydia trachomatis* and *Neisseria gonorrhoeae* DNA EQA Programme and 2 other from the QCMD 2017 *Chlamydia trachomatis* DNA EQA Programme). These panels consist of 19 clinical specimens dissolved in transport medium. VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit assay found 10/19 positive samples for *C. trachomatis* (LGV) (6 Simulated Swabs and 4 urines). The results were compared with the QCMD 2017 *Chlamydia trachomatis* and *Neisseria gonorrhoeae* DNA EQA Programme and the QCMD 2017 *Chlamydia trachomatis* DNA EQA Programme final report. All *Chlamydia trachomatis* (LGV strains) samples could be detected. Non-LGV *Chlamydia trachomatis* were not detected.

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit was evaluated in a Multicenter Evaluation conducted through collaboration with the National Microbiology Departments from different entities. VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit was evaluated using in total 115 clinical specimens (69 endocervical specimens, 4 urethral samples, 21 urogenital specimens, 9 rectal sample, 11 urine samples and 1 serum sample). VIASURE assay found 16/115 positive samples for *C. trachomatis* (LGV) infection (7 endocervical biological matrix, 4 urogenital specimens, 4 rectal samples, and 1 urine sample). These samples were analysed for *H. ducreyi* as well, proving to all be negative.

Besides, the clinical performance of VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit for *H. ducreyi* was testing using urine samples from symptomatic patients. Positive results were obtained for 10 samples which were spiked with synthetic *H. ducreyi* DNA a prior to nucleic acid extraction.

In conclusion, the results show a high sensitivity and specificity to detect LGV-associated strains of *Chlamydia trachomatis* and *H. ducreyi* using VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for *H. ducreyi* + *C. trachomatis* (LGV) (Figure 2 and 3).

Figure 2. Dilution series of *C. trachomatis* (LGV) (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).



Figure 3. Dilution series of *H. ducreyi* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).

12.3. Analytical specificity

The specificity of the *H. ducreyi* + *C. trachomatis* (LGV) assay was confirmed by testing a panel consisting of different microorganisms representing the most common sexually transmitted diseases pathogens. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Candida parapsilosis</i>	-	<i>Enterococcus faecium</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Candida tropicalis</i>	-	<i>Escherichia coli</i> O.1285; O18:H7:K1	-	<i>Serratia marcescens</i>	-
<i>Candida glabrata</i>	-	<i>Gardnerella vaginalis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
<i>Candida krusei</i>	-	<i>Haemophilus influenza</i> MinnA	-	<i>Stenotrophomonas maltophilia</i>	-
<i>Candida dubliniensis</i>	-	<i>Haemophilus ducreyi</i> class 1	-/+	<i>Streptococcus agalactiae</i>	-
<i>Candida albicans</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Streptococcus pneumoniae</i>	-
<i>Aspergillus fumigatus</i>	-	<i>Klebsiella pneumoniae</i>	-	<i>Treponema pallidum</i>	-
<i>Bacteroides fragilis</i>	-	<i>Listeria ivanovii</i>	-	<i>Trichomonas vaginalis</i>	-
<i>Acinetobacter baumannii</i>	-	<i>Listeria monocytogenes</i>	-	<i>Ureaplasma parvum</i>	-
<i>Chlamydia trachomatis</i> (Genovar F)	-	<i>Listeria innocua</i>	-	<i>Ureaplasma urealyticum</i>	-
<i>Chlamydia trachomatis</i> (LGV)	-/+	<i>Mycoplasma genitalium</i>	-	Cytomegalovirus	-
<i>Chlamydia trachomatis</i> (SW)	-	<i>Mycoplasma hominis</i>	-	Herpes simplex virus 1	-
<i>Enterobacter cloacae</i>	-	<i>Neisseria gonorrhoeae</i>	-	Herpes simplex virus 2	-
<i>Enterobacter aerogenes</i>	-	<i>Neisseria meningitidis</i>	-	Human papillomavirus 16	-
<i>Enterococcus faecalis</i>	-	<i>Proteus mirabilis</i>	-	Human papillomavirus 18	-

Table 5. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit for *C. trachomatis* (LGV strains) was evaluated against *C. trachomatis* (LGV) showing positive result.

The reactivity of VIASURE *H. ducreyi* + *C. trachomatis* (LGV) Real Time PCR Detection Kit for *H. ducreyi* was evaluated against *Haemophilus ducreyi* class 1 showing positive result.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure, section 8.3) and transferred into the specific tubes designed to perform on Rotor-Gene® Q or SmartCycler® instruments.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". Set exposition values as follow:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -150, HEX channel – 3000, ROX channel – 2000 and Cy5 channel - 1500.
- DTLite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 150, HEX channel - 500, ROX channel – 500 and Cy5 channel – 100.

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3. R.S Turingan *et al.* Rapid detection and strain typing of *Chlamydia trachomatis* using a highly multiplexed microfluidic PCR assay. *PLoS ONE* 2017; 12(5): e0178653.
4. S. Alexander *et al.* A comparison of two methods for the diagnosis of lymphogranuloma venereum. *Journal of Medical Microbiology* 2008; 57, 962–965.
5. M. Rodríguez-Domínguez *et al.* Clinical and epidemiological characterization of a lymphogranuloma venereum outbreak in Madrid, Spain: co-circulation of two variants. *Clinical Microbiology and Infection* 2014; 20(3): 219-225.
6. M. Glatz *et al.* A multicenter prospective trial to asses a new real-time polymerase chain reaction for detection of *Treponema pallidum*, herpes simplex-1/2 and *Haemophilus ducreyi* in genital, anal and oropharyngeal ulcers. *Clinical Microbiology and Infection* 2014; 20(12): O1020-7.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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For information purposes only

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: October 2018

For information purposes only

For information purposes only

For information purposes only

For information purposes only

CerTest Biotec, S.L.
Pol. Industrial Río Gállego II - Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-349 rev00

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC