



| | | |
|---|---|--------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 1 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |
| Prepared by QA Committee | | |
| Issued by: Laboratory Manager | Revision Date: 8/14/2023 | |
| Approved by Laboratory Director: Microbiologist-in-Chief | Next Review Date: 8/14/2025 | |

Uncontrolled When Printed

Table of Contents



| | |
|--|----|
| Table of Contents | 1 |
| Introduction..... | 6 |
| Secondary Specimen Labelling | 6 |
| General Processing Instructions: | 6 |
| Aliquoting..... | 7 |
| Planting Specimens in ES swabs..... | 8 |
| Planting by WASP..... | 8 |
| Manual Planting..... | 8 |
| Broth Inoculation for Infection Control Outbreak Investigation..... | 8 |
| BLOOD CULTURE | 9 |
| Operation of the Virtuo System :..... | 9 |
| Accessioning BC Bottles:..... | 10 |
| Loading Bottles: | 10 |
| Unloading Negative Bottles: | 11 |
| Processing macroscopically positive bottles not loaded in Virtuo:..... | 11 |
| Processing BC Bottles Flagged Positive: | 12 |
| Subculturing: bottles for WASP Streaking..... | 12 |
| Subculturing: if WASP rejects the BC+ tube or the second bottle | 13 |
| Manual Back-up Protocol..... | 15 |
| SPS Blood Culture Protocol for Dimorphic Fungi | 16 |
| Handling of Body Fluids in Blood Culture Bottles | 17 |
| Enterics | 19 |
| Faeces/Rectal Swabs | 19 |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 2 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| | |
|--|----|
| Rectal/Large Bowel (Colon) Biopsies | 21 |
| Genital Specimens | 22 |
| Cervical Swabs | 22 |
| Group B Streptococcus Screen | 22 |
| Vaginal Swab for Vaginitis/Vaginosis Screen | 25 |
| Vaginal Swab for Culture | 26 |
| Urethral Swab | 26 |
| Penis Swab..... | 26 |
| Seminal Fluid..... | 27 |
| Endometrial Swabs, Biopsies and Curettings, Placenta Swab/Tissue,..... | 28 |
| Products of Conception, Endometrial/Uterine, Cul de sac/Transvaginal,..... | 28 |
| Fallopian Tube, Tubo-Ovarian Swabs or Aspirates | 28 |
| Genital Ulcer Swab..... | 28 |
| Intra-Uterine Device | 29 |
| Infection Control Screening | 30 |
| Methicillin Resistant Staphylococcus aureus (MRSA) Screen | 30 |
| Vancomycin Resistant Enterococcus (VRE) Screen | 30 |
| Resistant Gram Negative Bacilli Screen | 31 |
| ESBL & CRE Screen..... | 32 |
| CRE Screen | 32 |
| RESISTANT Pseudomonas aeruginosa SCREEN | 32 |
| Mycology | 34 |
| Respiratory Specimens..... | 37 |
| Bronchoalveolar Lavage (BAL) | 37 |
| CMV Surveillance Bronchoscopy Specimens..... | 39 |
| Bronchial Brush Specimens..... | 39 |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 3 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Epiglottal Swabs 40

Open Lung/Transthoracic Needle/Transbronchial Lung Biopsies/Lung Aspirates 40

Mouth Swabs 41

Nasal Swabs..... 42

Nasopharyngeal Swabs/Auger Suctions for Bordetella pertussis 42

Oral Abscess Swabs 42

Sinus/Antral Specimens..... 43

Sputum including Endotracheal tube and Tracheotomy Specimens, Bronchoscopy Aspirates/Washings 43

Throat Swabs 45

Gastric Aspirates/Biopsies (for Helicobacter pylori) 46

Gastric Aspirates/Swabs from Neonates or Stillborn 47

Respiratory Specimen Treatment Solutions 48

Sterile Fluids 49

 Cerebrospinal Fluids..... 49

 Other Sterile Fluids – Amniotic, Pleural (Thoracentesis/Empyema), Peritoneal (Ascites), Synovial (Joint), Pericardial, Tympanocentesis, Intraocular, Hydrocele Fluids etc 52

 Peritoneal Dialysis Effluent..... 54

 Predialysis Fluid 55

 Bone Marrow (Aspirates or biopsies)..... 56

 Blood, Platelets, & Other Transfusion Products 57

 Handling of Body Fluids in Blood Culture Bottles 58

Urines 59

 Urine Hold Specimens:..... 61



Wounds/Tissues/Aspirates 65

 Duodenal or Small Bowel Aspirate/Swab/Biopsy 65

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 4 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| | |
|--|----|
| Intraoperative/Interventional Swabs..... | 66 |
| Wound/Abscess Swabs and Drain | 67 |
| Bite Wound Swabs | 67 |
| Intravenous & Central Line Catheter Exit Site Swabs | 68 |
| Intraoperative/Interventional Abscess (Pus, Cyst Fluid or Aspirate) | 68 |
| Pus & Abscess Material (other than Intraoperative/Interventional, Rectal or Bartholin) | 70 |
| Rectal Abscess | 71 |
| Bartholin’s Abscess Swab/Aspirate | 71 |
| TISSUES, PROSTHETIC DEVICES, AND AUTOPSY SPECIMENS | 72 |
| Tissues/Biopsies (other than skin or transplant tissues)..... | 72 |
| Skin Biopsies | 73 |
| Transplant Specimens - Bone Graft & Cadaver Fascia/Tissue/ Swab Specimens/Donor Amniotic Fluid/Membrane; Donor Corneal Ring Material | 75 |
| Prosthetic Devices (e.g. Pacemaker Wire, Dacron Graft, Prosthetic Valve, Coronary/Vascular Stents) | 75 |
| Autopsy Specimens | 76 |
| CATHETER..... | 77 |
| Intravascular Catheter Tips..... | 77 |
| Peritoneal Dialysis Catheter/Canula..... | 77 |
| BILE SPECIMENS..... | 78 |
| Bile and Bile Stents (including urinary, urethral, prostatic stents)..... | 78 |
| MISCELLANEOUS FLUID SPECIMENS | 78 |
| Breast Milk – Clinical samples..... | 78 |
| Total Parental Nutrition (TPN)..... | 79 |
| EAR SPECIMENS..... | 80 |
| Ear Swab..... | 80 |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 5 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Tympanocentesis Fluid..... 80

EYE SPECIMENS..... 81

 Eye / Conjunctival / Lid Swabs 81

 Eye / Corneal Scrapings 81



 Intraocular Aspirates..... 83

 Lacrimal (Tear Duct) Stone / Secretions 83

FACIAL SPECIMENS 84

 Facial Swabs 84

Record of Edited Revisions 85

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 6 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Introduction

Specimen processing is an essential for accurate testing results. Many microbiological tests rely heavily on specified processing instructions and is the most important step in the recovery of pathogenic organisms. Treatment and handling of any specimen must be done in a safe and ethical manner following the procedures described below. Specimens are to be used tested only as requested and according to laboratory protocols. Proper approval and consent must be received for samples to be used for any other purpose

Secondary Specimen Labelling

Adequate specimen information must be provided on any secondary items related to a sample throughout all phases of testing. Secondary items include aliquots, tubes, slides, culture plates, extraction vessels, images etc created during testing.

- A single unique identifier (e.g. specimen number - LIS or auxiliary number) is sufficient to labels materials derived from the primary specimen for use in testing and must be able to be linked to the full patient and specimen specific information. This includes barcodes.
- The identifier must be legible and clear throughout all processing and conditions of storage.

General Processing Instructions:



- Print required labels one sample at a time, removing all labels from the printer before proceeding with subsequent samples.
- Process one specimen in the Biological Safety Cabinet (BSC).
- Refer to the [Laboratory Safety Manual](#) for best working practices in a BSC.
- Change gloves and perform hand hygiene after manipulating biological samples.
- Follow specific processing procedures per testing method used as described in each testing manual. (ie. PCR, Serology etc)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 7 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Aliquoting



Aliquoting or transferring primary specimens may be necessary for the testing process. Aliquots are considered a secondary specimen container and must follow secondary specimen labeling requirements.

Potential risk exists for errors when aliquoting including cross contamination.

When aliquoting:

- Work with one specimen at a time.
- Compare sample identifiers of primary and secondary containers.
- Use a new pipette for each sample.

Never transfer an aliquot back into its original container.

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 8 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Planting Specimens in ESwabs

Planting by WASP

Use the WASP routinely to plate all specimens collected with ESwab.
 WASP will inoculate 30µL of the ESwab transport medium to each media plate.
 For WASP operation instructions, please see

Manual Planting

For ESwabs with fluid:

1. Vortex the ESwab.
2. Use a transfer pipette and inoculate 1 drop of ESwab transport medium onto each media.

For ESwabs with NO fluid:

1. Use swab to inoculate each media
2. In the culture “Test Comment” add comment }NFLD found in the rejection keypad



Broth Inoculation for Infection Control Outbreak Investigation

1. Place 500uL (0.5 mL mark of transfer pipette) of the ESwab transport medium into a 3 mL tube of Brain Heart Infusion broth (BHIB), incubate overnight on the shaker at 37°C in the walin incubator.
2. Place 30uL (1 drop from transfer pipette) of the ESwab transport medium onto Specific Screening Agar, incubate agar as per routine screening protocol.
3. After overnight incubation, subculture the BHIB onto Specific Screening Agar, incubate agar as per routine screening protocol.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|-----------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 9 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

BLOOD CULTURE

Operation of the Virtuo System :

Blood culture Bottles are tracked by the Virtuo system and communicate directly with the LIS. All information about the Virtuo and BC bottles can be accessed on the large or small display screen.

The Virtuo manual is located beside the system for reference.

The large display screen



This screen monitors a number of functions including positive and negative bottles, warnings and errors, and is used to generate reports, searches and edit bottle data. These are controlled through 5 icons:

1. Home: Select this icon anytime to display the main large display screen.
2. Search: Use this icon to search for an order/ bottle to view or edit any related information (Order number, incubation time etc)
3. Display Bank: Select this icon to display information specific to each connected unit.
4. Reports: Select this icon to generate different reports.
5. Configure: Select this icon to view Virtuo configuration.

The small display screen

The small display screen is located above the bottle return chute and display 2 icons.

1. Unit temperature: Displays unit temperature in degrees Celsius.
2. Waste level: Displays waste level used in percentage.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 10 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Accessioning BC Bottles:

- i. For Isolator tubes follow [Isolator 10 protocol](#)
- ii. For C&S order on BC bottles or Fluids received in BC bottles refer to section: [Fluids received in BC bottles](#)
- iii. Enter any specimen comments or BCs with **special incubation times** (see table below) in “Order Comment” in the LIS

| Blood Culture | 5 days | 7 days | 14 days | 21 days |
|--|--------|--------|---------|---------|
| Routine (Blood, Bone marrow, Sterile fluids, Blood products) | X | | | |
| Fungus/Yeast | X | | | |
| Dimorphic Fungi & Cryptococcus* (Histoplasma, Blastomyces) <i>*notify ward / physician must also collect Isolator collection tube if not received. Process BC Bottles as per normal protocol</i> | X | | | |
| Bone Bank Blood | | X | | |
| ENACT Sterilities | | | X | |
| Brucella | | | | X |
| PD Effluent | | | | X |
| SBE/IE, PUO/FUO (Subacute bacterial endocarditis/Infective endocarditis & Pyrexia of unknown origin/Fever of Unknown origin | | | | X |

- iv. Place barcodes vertically without covering the bottle barcode.

Loading Bottles:



DO NOT load macroscopically positive bottles (yellow septum). Follow protocol for [processing positive BC bottles](#)

- i. Ensure each set of BC bottles are matching (same LIS/aux #). If mismatched bottles arrive, wait for other set or phone ward to correct issue.
- ii. Ensure labels are affixed vertically and bottle barcode is visible. If bottle barcode is covered, see [Bottles with Unreadable or Generic Labels](#)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 11 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

- iii. Clean bottles with 3% hydrogen peroxide or Virox wipes
- iv. Place bottles upright on the conveyor belt which automatically load into the Virtuo.
DO NOT LOAD MORE THAN 40 BOTTLES AT A TIME
 - Before walking away, verify the instrument is operating properly and bottles remain upright and are being properly loaded into the bottle indexer. Monitor instrument for any flashing indicatory lights.

For bottles with unreadable or generic labels:

If a bottle barcode or LIS label is unreadable/covered or a generic label was used it will be flagged as anonymous bottle and an alert will appear.

- i. Log in using your user name and password
- ii. Select the anonymous bottle alert.
 - If multiple anonymous bottles exist, they appear one record at a time for editing
- iii. If visible from the bottle image, enter bottle ID and/or type and save record.
- iv. If the bottle ID cannot be seen / generic label not applied, unload bottle.
 - An alert indicating a label is required will be displayed
 - Place a generic label on the bottle and reload bottle.
 - You must specify the bottle type in Virtuo

Unloading Negative Bottles:

- Negative bottles automatically unload and drop into the garbage container.
- Dump garbage into double-bagged, cardboard waste box. Label outside of the box with date and controller the bottles were removed from.
- A negative result is automatically released and finalized upon unloading.
- The waste capacity is noted on the top right corner of the small display.
 - An alert will display on the large display if the capacity is full or almost full.

Processing macroscopically positive bottles not loaded in Virtuo:

Macroscopically positive bottles have a lemon-yellow sensor disc or a bulging septum.



- Process **ONLY** the suspected positive bottle (do not sub-culture the matching bottle unless it is also macroscopically positive or flagged positive).

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 12 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

- Follow the appropriate positive [subculturing procedure](#) for the positive bottle.
- Do not load macroscopically positive bottles into the Virtuo. They may not be flagged as positive and result in false negative results.

Processing BC Bottles Flagged Positive:

Positive bottles will alert on the Virtuo and show a count of positive bottles on the screen.

- Touch “+” sign to unload positive bottles to the bottle retrieval area.
- Only one bottle can unload to each chute in the retrieval area at one time. Remove these in order for further positive bottles to be unloaded.
- Write controller number on the bottle when you remove it.
 - Bottles which need to be reloaded, must be returned to the controller where they originated.
- When bottle is removed the positive bottle count updates and the unload event is recorded.

*Note: If a culture bottle marked as “**Brucella**” is flagged positive, remove a small amount of the culture for a Gram smear ONLY. Notify technologist. DO NOT SUBCULTURE until further instructions from technologist.*

*Note: If a culture bottle is marked with ESO flag “**VDE alert for micro lab**”, please include subculture to BVRE in addition to routine culture media to aid in isolation of vancomycin dependent enterococci..*



Subculturing: bottles for WASP Streaking

- In softmic open **result entry**
- Scan the bottle LIS barcode
- Go to the media section, the appropriate bottle (FO2V, FNV or PEDV)
- From the keypad, pick **1. {BC+}**
 - If there is already an anaerobic bottle subbed, [sub the bottle manually](#).
- Click the Order/Entry icon
- Click “Yes” to save changes to Result entry
- From the order entry keypad, Click **Blood Culture** |>BLD
- Pick “**Aerobic BC subculture**” or “**Anaerobic BC subculture**” or “**Paediatric BC subculture**”
- Click the “Save” icon
- Print **one** specimen label for subculture tube.
- Put the large barcode label, the FO2 or FN or PED bottle and one new red capped BC+ subculture tube into a cup in the rack.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 13 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

- xii. Placing only one cup in the hood at a time, affix the large barcode label to a red capped BC+ tube for the WASP
- xiii. Clean the top of the blood culture bottle and the top of the BC+ tube with an alcohol swab.
 - Inspect bottle septum. If obviously bulging, use a venting unit to safely offload the gas.
- xiv. Place a needle adaptor into the blood culture bottle.
- xv. Turn the bottle upside down (45degree angle), with the adaptor in place. Making sure the needle is submerged in the blood.
- xvi. Allow the resin to settle (5-10 seconds)
- xvii. Insert the BC+ tube on the other end of the needle adaptor; allow 1.5 mL to 2 mL of blood to transfer from the bottle into the tube. There is an indicator (red zone) on the BC+ tube label. Ensure the blood is within the red zone.
- xviii. Place the blood culture in the numbered storage rack. Date the first positive bottle with today's date.
- xix. Place the tube in a WASP palette. Load no more than 3 BC+ blood tubes per pallet.
 - Run wash loop after every 3 BC+ tubes loaded
- xx. Place wasp specimen rack on Wasp 3. It is programmed to subculture and do a gram stain.
- xxi. Ensure that plates are streaked and put into appropriate incubators STAT
- xxii. Stain gram smear as per [Gram Stain protocol](#) STAT
- xxiii. Give the stained smear to the STAT technologist (or designate) STAT
- xxiv. Store red-capped BC tube from the wasp in storage boxes (as e-swabs are saved) at room temperature.

Subculturing: if WASP rejects the BC+ tube or the second bottle

- i. In softmic open **result entry**
- ii. Scan the bottle LIS barcode
- iii. Go to the media section, the appropriate bottle (FNV, FO2V, SAV/SNV, or PEDV)
- iv. From the keypad, pick **1. {BC+**
- v. Press the "Save" icon
- vi. Print media labels
- vii. From labels printed, on one small barcode label, write the bottle code of the positive bottle (e.g. FN/FA/P/SAV/SNV)
 - Place one label on a slide for the gram stain with date i.e. FN6
- viii. On the large Media labels, write the bottle code and date bottle turned positive (e.g. FN 21)
- ix. Label media and sub-culture 2 drops onto the following (whole) plates:



| Media | Incubation | |
|------------------------|-----------------|-----------------|
| Blood Agar (SUBBA) | CO ₂ | 35°C x 48 hours |
| Chocolate Agar (SUBCH) | CO ₂ | 35°C x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 14 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |



| | | |
|-----------------------------------|------------------|-----------------|
| MacConkey Agar (SUBMC) | CO ₂ | 35°C x 48 hours |
| Fastidious Anaerobic Agar (SUBBR) | AnO ₂ | 35°C x 48 hours |

- x. Inoculate slide with one drop of specimen for gram stain and spread drop with a clean slide.
 - Place slide on heading block to dry.
 - Stain as per [Gram Stain protocol](#)
- xi. Put plates on ISOplater for streaking
- xii. Place ANO₂ plates in ANO₂ holding jar.
- xiii. Give stained slide to technologist
- xiv. Place CO₂ plates in rack labeled “New Positive. No Maldi Yet” in CO₂ incubator in planting area
- xv. Keep bottles in the Positive bottles tray until the isolate has been frozen and the final report has been issued.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 15 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Manual Back-up Protocol

If the automated system should fail for >24 hours or if expired bottles, non-Virtuo bottles are received, the following procedures will be used.

1. Examine all bottles for macroscopic growth twice daily for day 1 and day 2, record results in the LIS.
2. Ensure BCBC media is checkmarked under column R, this prevents the sample from automatically finalizing.
3. For the remainder of the incubation period check bottles macroscopically at least daily and record in the LIS
4. For macroscopically positive bottles, process as indicated under "Positive Cultures".
5. If macroscopically negative, sub-culture as follows:

Direct Examination: Not required.



Sub-culture:

| Incubation Day | Bottle Type | Medium | Incubation |
|----------------|------------------------------|-------------------|---|
| 1 | Aerobic bottle | CHOC | CO ₂ , 35°C, x 48 hours |
| 2 & 5 | Aerobic and Anaerobic bottle | CHOC & BRUC | CO ₂ , 35°C, x 48 hours AnO ₂ , 35°C, x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 16 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

SPS Blood Culture Protocol for Dimorphic Fungi

Isolate 10 system collection tubes are currently unavailable. Follow this protocol utilizing SPS vacutainer tubes in lieu as a secondary protocol only to used until isolate tubes again become available.

Order Entry in Site Field: Add “SPS TUBE”

1. Centrifuge specimen at 3000g for 30 minutes

Working in a biosafety cabinet:

2. Clean vacutainer stopped and top of tube with virox or equivalent. Remove vacutainer stopper
3. Using an individually wrapped sterile pipette, remove and discard supernatant, not disrupting the buffy coat.
4. Place the original stop back onto the tube and vortex tube for at least 10 seconds on the highest setting.
5. Using a new individually wrapped sterile pipette, slowly aspirate concentrate left in tube.
6. Dispense concentrate in a straight line along the surface of the agar. Keep inoculum away from the edge of the plate

| Media | Incubation |
|---|---------------------------------|
| Inhibitory Mold Agar (IMA) | O ₂ , 28°C x 4 weeks |
| Esculin Base Medium (EBM) | O ₂ , 28°C x 4 weeks |
| Two Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) | O ₂ , 28°C x 4 weeks |



7. Using the tip of the pipetting, streak the plate. Use 15-20 passes perpendicular to the original inoculum line.
8. Dry plates. Seal with parafilm and forward plates to Mycology for incubation and processing.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 17 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Handling of Body Fluids in Blood Culture Bottles

Order/Entry:

1. PD effluent fluids arrive into microbiology in Blood Culture Bottles

- i. At Order/Entry Screen, call up order by Auxiliary number.
- ii. Add the test "Fluid in BC Bottle" (?BTLE)
- iii. Receive and plate as usual.
- iv. Add bottles received under media for C&S test
- v. Apply the specimen labels on the bottles (the labels with the lab number that contains the extension 93). **DO NOT** use the media labels or the specimen label for the C&S or FLDM test.
- vi. Load the bottles to Virtuo as usual by scanning the bottle barcode and then the specimen barcode.
- vii. Ensure the incubation time is changed 21 days for both bottles.



2. Sterile Body fluids, Blood bank products and Bone Marrow put in a Blood Culture Bottle ONLY.

- i. At Order/Entry Screen, call up order by Auxiliary number.
- ii. Add the test "Fluid in BC Bottle" (?BTLE)
- iii. Receive and plate as usual.
- iv. Before saving, to go to media screen. Cancel the media that are attached to the C&S or FLDM test and add appropriate bottle media codes (ex FO2, FN, ANAO2...)
- v. Save the order and print the labels
- vi. Apply the specimen labels on the bottles (the labels with the lab number that contains the extension 93). **DO NOT** use the media labels or the specimen label for the C&S or FLDM test.
- vii. In result entry, result Gram as "Not Applicable" and finalize.
- viii. Load the bottles to Virtuo as usual by scanning the bottle barcode and then the specimen barcode.
- ix. These all have an incubation time of 5 days for the bottles.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 18 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

3. Sterile Body fluids, Blood bank products and Bone Marrow that are sent to the lab in a sterile container AND a set of blood culture bottles.

- i. At Order/Entry Screen, call up order by Auxiliary number.
- ii. Add the test "Fluid in BC Bottle" (option 1)
- iii. Receive and plate as usual.
- iv. Do not cancel the media that are attached to the C&S / FLDM test as they will be processed. add appropriate bottle media codes (ex FO2, FN, ANAO2...)
- v. Save the order and print the labels
- vi. Apply the specimen labels on the bottles (the labels with the lab number that contains the extension 93). DO NOT use the media labels or the specimen label for the C&S or FLDM test.
- vii. Load the bottles to Virtuo as usual by scanning the bottle barcode and then the specimen barcode.
- viii. These all have an incubation time of 5 days for the bottles.

4. Aspirates from non-sterile sites e.g. wound drainage, pus and abscess material that are normally sent to the lab in a sterile container but arrive into the lab in a Blood Culture Bottle ONLY.

- i. Treat the bottle as a specimen transport medium. Planting bench will inoculate the assigned media from aspiration of the bottle contents.

5. Aspirates from non-sterile sites e.g. wound drainage, pus and abscess material that are sent to the lab in a sterile container AND a Blood Culture Bottle.

- i. Plant the specimen from the sterile container as usual and discard the blood culture bottle.

6. Ascitis Fluid in Blood Culture Bottles ordered in UHN HIS as “ASCBT” test.



- i. Receive the order in LIS as in routine Bacteriology orders.
- ii. Save the order and print the labels
- iii. Apply the specimen labels on the bottles (the labels with the lab number that contains the extension 91).
- iv. Load the bottles to Virtuo as usual by scanning the bottle barcode and then the specimen barcode.
- v. (NO NEED to order “?BTLE – Fluid in BC Bottle”)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 19 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Enterics

Faeces/Rectal Swabs

- a) Direct Examination: Not routinely performed.
For special request, prepare a Gram stain for faecal leukocytes.
- b) Culture:

| Media | Incubation | |
|-----------------------------------|------------------|----------------------|
| MacConkey Agar (MAC) | O ₂ , | 35°C x 18 - 24 hours |
| Hektoen Agar (HEK) | O ₂ , | 35°C x 18 - 24 hours |
| MacConkey Sorbitol Agar (SMAC) | O ₂ , | 35°C x 18 - 24 hours |
| Campylobacter Agar (CAMPY) | Campy Jar | 42°C x 48 hours |
| Selenite Broth (SEL) ¹ | O ₂ , | 35°C x 12-18 hours |

If *Yersinia* is requested or patient is >1 month–12 years old (except for NICU), **add:**

| | | |
|--|------------------|-----------------|
| Cefsulodin Irgasan Novobiocin Agar (CIN) | O ₂ , | 30°C x 48 hours |
|--|------------------|-----------------|

If *Vibrio* is requested, **add:**

| | | |
|---|------------------|----------------------|
| Thiosulphate Citrate Bile Salt Sucrose Agar (TCBS) | O ₂ , | 35°C x 18 - 24 hours |
| Alkaline Phosphate Broth (APB) ² | O ₂ , | 35°C x 5 - 8 hours |

If *Pleiomonas* is requested, **add:**

| | | |
|-----------------|------------------|-----------------|
| Blood Agar (BA) | O ₂ , | 35°C x 24 hours |
|-----------------|------------------|-----------------|

If *Aeromonas* is requested, **add:**

| | | |
|--|------------------|-----------------|
| Blood Agar (BA) | O ₂ , | 35°C x 24 hours |
| Cefsulodin Irgasan Novobiocin Agar (CIN) | O ₂ , | 35°C x 24 hours |

If *Neisseria gonorrhoeae* (GC) is requested (rectal swab only), **inoculate only:**



| | | |
|------------------------|-------------------|-----------------|
| Martin-Lewis Agar (ML) | CO ₂ , | 35°C x 72 hours |
|------------------------|-------------------|-----------------|

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 20 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

If *C. difficile* toxin assay is requested set up:

[Cepheid GeneXpert PCR](#)

For Mount Sinai Hospital patients: Specimen must also set up for Vancomycin Resistant Enterococcus (VRE) screen:

Brilliance VRE Agar (BVRE)

O₂,

37°C x 36hrs in the dark

- Notes:**
1. Subculture Selenite broth following overnight incubation onto HEK. Incubate the sub-cultured HEK at 35°C in O₂ for 12 - 18 hours.
 2. Subculture APB to TCBS Agar after 5-8 hours incubation. Planter has to notify the technologist on the Enteric bench at the time of processing. Incubate the TCBS Agar at 35°C in O₂ for 18 – 24 hours.

WASPLab Stool Processing:

Stool cultures received in Copan FecalSwab will be processed by the WASP and incubated within the

Processing Stool:



1. Upon receipt of specimens, put EPT specimens in a WASP palette marked for Stool Culture. The “test ordered” should only be for Stool C&S.
2. Palettes will be loaded onto the WASP at 9am, 2pm, 9 pm and 2 am. Fill the silos of the WASP2 with the appropriate plates prior to loading the palette.
3. Let WASP inoculate and streak the media and transfer the media plates to WASPLab O2 incubator. On the WASP, flip the orientation of the loop after stools are complete. Perform a wasp cycle on the loop.
4. The CAMPY plates will be sent to stacker 202 of the WASPLab line. Put the CAMPY plates into an offline incubation box, set up the box up with CAMPY gas pack. Incubate the box in the 42oC incubator.
5. Unload the Selenite broths, put them in a rack by their incubation time and incubate them in the 37oC.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 21 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Subculturing selenite broths:

1. Remove Selenite broths to subculture by the “load time” in the selenite rack (12 hrs after incubated). Eg. At 9am, load selenites that were set up at 9pm.
2. In the LIS:
 - a. Open “Enterics” **Resulting Worklist**
 - b. At the bottom left, click “Mark”. Scan all selenite LIS barcodes
 - c. On the right panel, click “Define MC (^M)”. Type “SEL” and hit OK.
 - d. From the keypad, select “^SFHEK”. Select “OK”
 - e. On the right panel, click “Add Result (^F7)”. Exit worklist.
3. Load the selenites onto a palette and onto WASP2. Store selenites in a selenite done box. Keep 2 boxes worth. Discard older box when both are full.

Rectal/Large Bowel (Colon) Biopsies



- a) Direct Examination: Gram stain not indicated.
- b) Culture:

| Media | Incubation |
|--------------------------------|-------------------------------------|
| MacConkey Agar (MAC) | O ₂ 35°C x 18 – 24 hours |
| Hekton Agar (HEK) | O ₂ 35°C x 18 – 24 hours |
| MacConkey Sorbitol Agar (SMAC) | O ₂ 35°C x 18 – 24 hours |
| Camyplobacter Agar (CAMPY) | Campy Jar 42°C x 48 hours |
| Selenite Broth (SEL) | O ₂ 35°C x 18 – 24 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 22 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Genital Specimens

Cervical Swabs

- a) Direct Examination: Not indicated.
- b) Culture:

| Media | Incubation |
|------------------------|-----------------------------------|
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 72 hours |
| Martin-Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |

If Group B streptococcus is requested, refer to the Group B streptococcus screen section.

Group B Streptococcus Screen

- a) Direct Examination: Not indicated.
- b) Culture:

| Media | Incubation |
|--|----------------------------------|
| Carrot Broth for Group B Strep (CAROT) | O ₂ , 35°C x 24 hours |

Processing on the WASP:

1. Load wasp with CAROT broths
2. Load GBS specimens into the WASP
3. When WASP processing is complete, in the hood place one full tile into each CAROT broth.
4. Incubate CAROT broths.



Manual Processing: (when WASP is unavailable)

1. Inoculate 2-3 drops of each GBS specimen into each CAROT broth
2. Add one full tile into each CAROT broth.
3. Incubate CAROT broths.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 23 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

c) Subculture:

| Media | Incubation |
|------------------------------------|----------------------------------|
| Colistin Nalidixic Acid Agar (CNA) | O ₂ , 35°C x 24 hours |



To Subculture Carrot Broths by Wasp:

1. Gynae bench technologist will bring colourless Carrot broths to the planting area and process specimens on WASP themselves once WASP is unoccupied.
 - a. First thing in the morning after cleaning is an ideal time
2. Ensure broth has WASP label. Otherwise, generate a large LIS barcode label and affix to swab.
3. Place broths in WASP **BLUE** palette(s) labelled 'BLU URISTRAW'.
4. Obtain adequate number of CNA plates for the broths being subcultured
5. Ensure that the WASP is not processing specimens before proceeding.
6. Take plates and palettes to the WASP:
 - Check status of WASP, by clicking on the "MAIN" icon on bottom toolbar
 - If "Paused", click "STOP", then "RESET".
 - Load palettes onto conveyor belt
 - If "READY", and there are plates to remove from the WASP, click "OPERATIONS", then click "UNLOAD CAROUSEL". Four stacks will appear on the screen. Click on each stack that has plates in it. Unload plates to a trolley. Click "OK". Otherwise, proceed to the next step
 - Click "CLOSE" (bottom right corner)
 - Click "SWITCH MODE"
 - Click "TECHNOLOGIST"
 - Click "SETTINGS" (icon on bottom toolbar)
 - Click "CAROUSEL" (icon in the middle of the screen)
 - Click "LOAD CONFIGURATION"
 - Click "GBS-CNA"
 - Click on the number on the carousel designated (6)

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 24 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

- Load CNA plates into the carousel (If you need to load more than the tower will hold, fill up to the max and the WASP will alert you when to add the remainder.)
- Click “OK”
- Click “Start” icon
- Click on the icon “COPAN 12”
- Using the mouse, slide the cursor down to you see the icon and highlight the icon of a wasp with “GBS-CNA”)
- Click “OK”
- Load palettes
- WASP will begin processing

7. After subculturing is complete:

- Click “MAIN”
- Click “Stop”
- Click “Operations” ”, then click “UNLOAD CAROUSEL”. Four stacks will appear on the screen. Click on each stack that has plates in it. Unload plates to a trolley. Click “OK”.
- Click “MAIN”
- Click “RESET”
- Click “SETTINGS” (icon on bottom toolbar)
- Click “CAROUSEL” (icon in the middle of the screen)
- Click “Load Configuration”
- Choose IC-UR-WND-BRUC-KV
- Click “OK”
- Click “Close”
- Click “SWITCH MODE”
- Click “LIS”
- Check bin for rejects. Sub those manually.
- Unload palettes



8. Return Carrot broths and SUBCN plates to the Gynae bench.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 25 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Vaginal Swab for Vaginitis/Vaginosis Screen



a) Direct Examination:

- i. Wet mount: To be set up immediately. Gently press the swab into a drop of sterile saline on a slide. For Eswabs, vortex swab and add 2 drops on a slide. Place a cover slip on the slide and examine under the microscope using the 40X objective. Examine for the presence of *Trichomonas vaginalis*.
- ii. Gram stain: Examine for the presence of yeast, clue cells and organisms associated with bacterial vaginosis.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 26 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Vaginal Swab for Culture

- a) Direct Examination: not required
- b) Culture:

| Media | Incubation |
|---|-----------------------------------|
| Colistin Nalidixic Acid Agar (CNA) | CO ₂ , 35°C x 48 hours |
| Carrot Broth for Group B Strep (CAROT) | O ₂ , 35°C x 24 hours |
| Chocolate Agar (CHOC) (if GC is requested) | CO ₂ , 35°C x 72 hours |
| Martin-Lewis Agar (ML) (if GC is requested) | CO ₂ , 35°C x 72 hours |
| MacConkey (MAC) (for <12 years old) | CO ₂ , 35°C x 48 hours |

Urethral Swab

- a) Direct Examination: Gram stain - Quantitate the presence of pus cells and intracellular gram negative diplococci.
- b) Culture:



| Media | Incubation |
|------------------------|-----------------------------------|
| Martin-Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |

Penis Swab

- a) Direct Examination: Gram stain.
- b) Culture:

| Media | Incubation |
|-------|------------------------------------|
| BA | CO ₂ , 35°C, x 48 hours |
| MAC | O ₂ , 35°C, x 48 hours |
| CNA | O ₂ , 35°C, x 48 hours |
| CHOC* | CO ₂ , 35°C, x 72 hours |
| ML* | CO ₂ , 35°C, x 72 hours |

*Occasionally, urethral swabs may be labelled as penile swabs. *Neisseria gonorrhoeae* culture is set up for this reason.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 27 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

The ML should be inoculated by rotating the swab in a Z streak manner. The inoculum is then streaked by the ISOplater to obtain discrete colonies. Examine ML plates at 48 and 72 hours.

Seminal Fluid

Set up all specimens for culture and sensitivity testing as outlined below. After inoculating the culture media below, forward the remainder of the specimen to the Virology lab for molecular detection of GC and CT

a) Direct Examination: Not required



b) Culture:

| Media | Incubation |
|--------------------------------------|-----------------------------------|
| Blood Agar (BA)* | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) ^{ab} | CO ₂ , 35°C x 72 hours |
| Martin-Lewis Agar (ML) ^{ab} | CO ₂ , 35°C x 72 hours |
| MacConkey Agar (MAC)* | O ₂ , 35°C x 48 hours |

*Use a 10 µl disposable culture loop to inoculate media by using urine colonies count culture technique.

^aDilute specimen 1:2 in a 15mL-conical tube with sterile saline before inoculating CHOC and ML agar. Store the original sample and the labelled conical tube with miscellaneous specimens.

^bUse a sterile pipette apply one drop to media and streak

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 28 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

**Endometrial Swabs, Biopsies and Curettings, Placenta Swab/Tissue,
 Products of Conception, Endometrial/Uterine, Cul de sac/Transvaginal,
 Fallopian Tube, Tubo-Ovarian Swabs or Aspirates**

- a) Preparation of specimen for culture:
1. If tissue is received, aseptically macerate the tissue with the use of a tissue grinder or stomacher.
 2. If STAT TB is requested and approved by the Microbiologist, prepare a slide for Acid Fast Bacilli ([FLUOROCHROME STAIN for MYCOBACTERIUM](#)). A portion of the specimen should be forwarded to the Public Health Laboratory (PHL) for processing.
 3. If TB culture is requested, send half of the specimen to PHL.

b) Direct examination: Gram stain.



c) Culture:

| Media | Incubation | |
|------------------------------------|--------------------|-----------------|
| Blood Agar (BA) | CO ₂ , | 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , | 35°C x 48 hours |
| Martin-Lewis Agar (ML) | CO ₂ , | 35°C x 72 hours |
| MacConkey Agar (MAC) | O ₂ , | 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC)* | AnO ₂ , | 35°C x 48 hours |
| Kanamycin-Vancomycin Agar (KV)* | AnO ₂ , | 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO)* | O ₂ , | 35°C x 48 hours |

* If tissue/aspirate or anaerobic swab is received, add BRUC, KV and THIO.
 For tissues and biopsies, freeze remaining tissue in the -70°C freezer for minimum of 3 months.



Genital Ulcer Swab

[Refer to Ministry of Health Specimen Collection Guide](#)

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 29 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Intra-Uterine Device

- a) Direct Examination: Gram stain of secretions.
Examine for the presence of branching gram positive bacilli suggestive of *Actinomyces* species.
- b) Culture: Not indicated.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 30 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Infection Control Screening

Methicillin Resistant Staphylococcus aureus (MRSA) Screen

- a) Direct Examination: Not indicated
- b) Culture:
Routine Screening:

| Media | Incubation |
|--------------------------------|---|
| OXOID Denim Blue Agar (DBLUE)* | O ₂ , 37°C x 24 h -in the dark |

If multiple swabs from a single patient are received individually, then process as separate specimens. If multiple swabs from a single patient are received as a "bundle" with a single label and order number, then process all swabs in the bundle on a single "DBLUE" plate.

On Fridays and Saturdays, specimens will not be planted past 3 pm. Any specimens received after this time will be held and planted the following morning.

These will be stored in a basket labelled for this purpose in the planting refrigerator.

Vancomycin Resistant Enterococcus (VRE) Screen

- a) Direct Examination: Not indicated
- b) Culture in non-outbreak setting:

| Media | Incubation |
|----------------------------|---|
| Brilliance VRE Agar (BVRE) | O ₂ , 37°C x 48h in the dark |



Inoculate the BVRE Agar by rotating the swab on the primary inoculum area to the size of 2.5 cm in diameter (size of a Loonie).

Put inoculated/streaked plates into the VRE bin in the incubator.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 31 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

On Fridays and Saturdays, specimens will not be planted past 3 pm. Any specimens received after this time will be held and planted the following morning.

NB: In the event of an outbreak investigation, the above protocol (b) may not apply. VRE PCR maybe ordered on prior arrangement by Infection Control or Microbiologist. All swabs for VRE PCR are to be received by 8 am so they may be processed as a batch. VRE PCR positive specimens will be processed as per protocol (c) below.

c) Culture for **VRE PCR positive samples** in outbreak setting:

| Media | Incubation time (all O₂ at 37°C) |
|---|--|
| i) Place 500uL (0.5 mL mark of transfer pipette) of the eSwab transport medium into: 3 mL Brain Heart Infusion broth (BHIB) Place 30uL (1 drop from transfer pipette) of the eSwab transport medium onto: Brilliance VRE Agar (BVRE) | overnight on shaker 48h in the dark |
| ii) If BVRE is no growth after 24 hours of incubation, subculture swab from BHIB to: Brilliance VRE Agar (BVRE) | 48h in the dark |

Resistant Gram Negative Bacilli Screen

a) Direct Examination: Not indicated



b) Culture:

| Media | Incubation |
|---|------------------------------|
| For <i>Enterobacteriaceae</i> with fluoroquinolone and/or aminoglycoside resistance but susceptibility to cefpodoxime: MacConkey Agar (Mac) –no antibiotic | O ₂ , 35°C x 18 h |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 32 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

For *Serratia marcescens* outbreaks, use MacConkey Agar + 7.5 mg/L colistin; in-house media preparation to be cleared by Microbiologist.

ESBL & CRE Screen

- a) Direct Examination: Not indicated
- b) Culture:

| Media | Incubation |
|---|-------------------------------------|
| ESBL Isolation Agar - MacConkey agar with 2 µg/ml cefpodoxime (Media code: MCPOD) | O ₂ , 37°C x 18-24 hours |

CRE Screen

- a) Direct Examination: Not indicated
- b) Culture:

| Media | Incubation |
|---|-------------------------------------|
| ESBL Isolation Agar - MacConkey agar with 2 µg/ml cefpodoxime (Media code: MCPOD) | O ₂ , 37°C x 18-24 hours |



RESISTANT *Pseudomonas aeruginosa* SCREEN

| Specimen | Processing | Media | Incubation |
|----------|--|-----------|----------------------------------|
| Water | Centrifuge the entire sample at 3500 rpm for 20 minutes. Pour off all supernatant. Transfer the contents of a 2 mL tube of BHI | BHI Broth | O ₂ at 35°C overnight |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 33 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| Specimen | Processing | Media | Incubation |
|--|--|------------------------------------|--|
| | broth into in the falcon tube containing the sediment Subculture after overnight incubation to MCPOD by the IC Bench technologist | MCPOD | O ₂ at 35°C for 24 hours |
| Environmental swabs | Incubate the BHI Broth Subculture BHI broth after overnight incubation to MCPOD by the IC Bench tech using a new sterile swab | MCPOD | O ₂ at 35°C overnight O ₂ at 35°C for 24 hours |
| Patient pharmaceutical infusates/injectables | ≤1 mL >1 mL | TH14 SD14 ETH14 ESD14 | O ₂ at 35°C for 14 days O ₂ at RT ^o for 4 days O ₂ at 35°C for 14 days O ₂ at RT ^o for 4 days |
| Swabs from patients | Directly inoculate MCPOD plate with specimen | MCPOD | O ₂ at 35°C for 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 34 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Mycology

1. All specimens, except for skin scrapings, nail and hair (see below) are processed initially in the general planting section of the Microbiology Laboratory. Requests for skin scrapings, nails and hair should be forwarded to the Public Health Lab (PHL) for processing. Fungal smears and inoculated culture plates should be forwarded to the Mycology section for staining, interpretation, incubation and work-up.

2. **Fungal smears:**
 - i) At least one, and where applicable, two smears from each specimen should be prepared and forwarded to the mycology section. One smear is for fungal staining and the second is held in reserve for special stains (if indicated).
 - ii) In the planting area, the smears should be fixed in 99.9% methanol for 5 minutes before staining.
 - iii) Fungal smears should be stained using Fungi-Fluor™ (Refer to Staining Methods).

3. **Culture:**
 - i) Refer to Bacteriology manual and/or Table 1 below, for appropriate media to be inoculated depending on specimen type. Between 0.1 and 0.5 mL of fluid specimens is recommended for each media.
 - ii) Tissue specimens should be mixed with filtered sterile saline. Gently cut the tissue with a sterile scalpel and use sterile forceps to inoculate the media
NB: DO NOT grind the specimen or put it into the Stomacher to maintain viability of possible fungi including Zygomycete .
 - iii) Sterile body fluids and urine specimens should be centrifuged at 3500 rpm x 20 min prior to inoculating them. The supernate should be removed leaving sufficient volume to resuspend the sediment and inoculate the appropriate number of culture media.
 - iv) Inoculated culture plates should be held unsealed overnight at room temperature to dry the plates. The following morning, all plates should be placed in a plastic bag and sent to Mycology to be sealed with Parafilm to prevent contamination and incubated at 28°C.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 35 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Table 1 - PROCESSING AND INCUBATION OF SPECIMENS FOR MYCOLOGY

| SPECIMEN | FS | IMA | EBM | BHIM | SAB G | MYCOSEL | Incubation ⁴ O ₂ , 28°C |
|---|----|-----|-----|------|-------|---------|--|
| Bone Marrow ⁵ | X | X | X | X | | | 4 weeks |
| Brain Biopsy/Abscess | X | X | X | X | | | 4 weeks |
| Sterile Fluids: | | | | | | | |
| CSF ¹ | | X | X | X | | | 4 weeks |
| Pleural Fluid | X | X | X | X | | | 4 weeks |
| Synovial Fluid | X | X | | X | | | 4 weeks |
| Pericardial/Ascitic Fluid | X | X | | X | | | 4 weeks |
| Peritoneal/Dialysis Effluent ² | X | X | | X | | | 4 weeks |
| Aqueous/Vitreous | X | X | | X | | | 4 weeks |
| Eye/Corneal Ulcer | X | X | | | | | 3 weeks |
| Blood Culture ^{3,5} | | X | X | X | | | 4 weeks |
| Lung Biopsy | X | X | X | X | | | 4 weeks |
| Tissue Biopsy | X | X | | X | | | 4 weeks |
| Lymph Node | X | X | X | X | | | 4 weeks |
| Skin Biopsy | X | X | | X | | | 4 weeks |
| Clot/Thrombus | X | X | | X | | | 4 weeks |
| Bone/Joint | X | X | | X | | | 4 weeks |
| Heart Valves | X | X | X | X | | | 4 weeks |
| Intravascular Grafts | X | X | X | X | | | 4 weeks |
| Abscess/Pus/Wound ² | X | X | | X | | | 4 weeks |
| Autopsy | X | X | | | | | 3 weeks |
| Respiratory Secretions: | | | | | | | |
| Sputum ² | X | X | X | X | | | 4 weeks |
| Tracheal Aspirate ² | X | X | X | X | | | 4 weeks |
| Bronchial Washing | X | X | X | X | | | 4 weeks |
| BAL (Routine) | X | X | X | X | | | 4 weeks |
| BAL(Lung Transplant Screen) | X | X | X | X | | | 4 weeks |
| Sinus/Antral Wash | X | X | | | | | 3 weeks |
| Nose ² | | X | | | | | 3 weeks |
| Mastoid | X | X | | | | | 3 weeks |
| Ear | X | X | | | | | 3 weeks |
| Urine ² | | X | X | | | | 3 weeks |



1. Cryptococcal Antigen in the Planting Area
2. Only if Fungal Culture is specifically requested
3. Collected in "Isolator 10" Tube
4. If dimorphic fungus is requested, regardless of site, media should be incubated for 6 weeks.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 36 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |



5. For Malassezia request, incubate IMA And IMAO for one week Only

| | |
|---------|---|
| FS | Fungal Stain |
| IMA | Inhibitory Mould Agar |
| IMAO | Inhibitory Mould Agar with Chloramphenicol, sterile olive oil to be overlaid at Mycology bench. |
| BHIM | Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide |
| EBM | Aesculin agar with chloramphenicol, gentamycin |
| SAB G | Sabouraud Dextrose Agar |
| MYCOSEL | contains Chloramphenicol and Cyclohexamide |
| PDA | Potato Dextrose Agar |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 37 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Respiratory Specimens

Bronchoalveolar Lavage (BAL)

a) Specimen Processing

Bronchoalveolar Lavage specimens from the same patient should be pooled together and processed as a single sample EXCEPT specimen from Trillium Gift of Life Donor (identified by MRN prefixed with “DONOR”)

BAL Test specimen requirement if requested*:

- 1mL for bacteriology culture
- 1 mL for fungus (if requested OR specimen is from lung transplant patients)
- 5 mL for TB (if requests, send to PHL)
- 2 mL for Virology (if requested) - Supernatant
- 1 mL for Legionella (if requested, send to PHL)
- 1 mL for *Chlamydia pneumoniae* (if requested, send to PHL)
- 1 mL for Galactomannan - Supernatant

C&S

1. Prepare a bacterial C&S gram smear with unspun specimen using a Cytospin funnel.
2. For bacterial culture requests, process specimens with [liquefying solution](#) before placing on the WASP for processing using 1mL loop.

During WASP downtime: process [liquefied specimen](#) and inoculate media manually using a calibrated 1 uL loop to inoculate unspun sample onto each routine bacteriology plates below. Streak the loop down the center of each of the first set of plates and then cross-streaked at a 90° angle to the inoculum.

Other requested tests processing



If the remaining specimen volume is less than the amount the requested tests volumes, call the physician to determine the priority of test(s) to be performed. Result the test(s) that cannot be done as “Insufficient sample for _____.” Document in the LIS call window the name of the physician contacted.

1. Centrifuge the entire specimen at 3500 rpm for 20 minutes.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*
 Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 38 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

2. Remove all but 0.5 mL of the supernatant to a new sterile container labelled as “Supernatant”.
3. Vortex the sediment to mix well.
4. Prepare smears (for Fungifluor stain and PCP IFA if tests are requested, respectively). [See Direct Examination Preparation below.](#)
5. Split/aliquot specimen for remaining tests requested requiring supernatant. Reconstitute the rest of the sediment with 1-10 mL of reserved supernatant for remaining requested tests. [See test requirements above.](#)

b) Direct Examination Preparation:

Stain the prepared smears as follows:

- Cytospin smear for Gram stain. Following Gram Stain procedure and deliver stained slide to the Gram Stain Bench.
- Concentrate smears (for Fungifluor stain and PCP IFA if C&S, fungus culture and PCP are requested respectively):
 - Forward slide to Mycology Section for Fungifluor staining and interpretation (if fungal culture is requested).
 - Forward slide to Virology Section for PCP IFA staining and interpretation (if PCP is requested).

c) Culture Preparation:

Inoculate the specimen onto the following media:



| Media | Incubation |
|--|-----------------------------------|
| Inoculate with unspun specimen using 1 uL loop : | |
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | CO ₂ , 35°C x 48 hours |
| If <i>B. cepacia</i> is requested or specimen is from a patient with Cystic Fibrosis, add : | |
| OF Base, Colistin, Bacitracin & Lactose Agar (OCBL) | O ₂ , 35°C x 5 days |
| Keep the BA, HI and MAC plates | CO ₂ , 35°C x 5 days |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 39 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Inoculate with the **resuspended** specimen:

If Fungus is requested OR specimen is from lung transplant patients, **add:**

| | | |
|---|------------------|----------------|
| Inhibitory Mold Agar (IMA) * | O ₂ , | 28°C x 4 weeks |
| Esculin Base Medium (EBM)* | O ₂ , | 28°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM)* | O ₂ , | 28°C x 4 weeks |

If Nocardia is requested, **add:**

| | | |
|-----------------------------|------------------|----------------|
| Sodium Pyruvate Agar (PYRA) | O ₂ , | 35°C x 4 weeks |
|-----------------------------|------------------|----------------|

* Forward inoculated fungal media to Mycology Section for incubation and work-up.

CMV Surveillance Bronchoscopy Specimens

Specimens for **routine CMV surveillance** require CMV culture only. Refer specimens to Virology immediately upon receipt.

Bronchial Brush Specimens

All brushes received in small cap (bijou) bottles in 1 mL Ringer's lactate solution should be processed "STAT" because colony counts may change if processing is delayed.

If the brush is received in <1 mL of fluid, add sufficient Ringer's lactate solution to make up 1 mL. Document in the "Test Comment" field of the LIS as "Brush received in wrong volume of fluid".

If a dry brush is received, add 1 mL of Ringer's lactate solution into the vial. Document in the "Test Comment" field of the LIS as "Dry brush received".

Mix the vial well by vortexing vigorously.

Measure 0.01 mL of the sample (using a calibrated **10 uL loop**) and inoculate onto routine media for culture of aerobic bacteria. Streak the loop down the center of the plate and then cross-streaked at a 90° angle to the inoculum top to bottom.



a) Direct Examination: Not indicated.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 40 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

b) Culture:

| Media | Incubation |
|------------------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilias Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | CO ₂ , 35°C x 48 hours |

If *B. cepacia* is required or specimen is from a patient with Cystic Fibrosis, **add:**

| | |
|---|---------------------------------|
| Of Base, Colistin, Bacitracin & Lactose Agar (OCBL) | O ₂ , 35°C x 5 day |
| Keep the BA, HI and MAC plates | CO ₂ , 35°C x 5 days |

Epiglottal Swabs

- a) Direct examination: Not indicated
- b) Culture:

| Medium | Incubation |
|-----------------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |



Open Lung/Transthoracic Needle/Transbronchial Lung Biopsies/Lung Aspirates

- a) Direct examination: Using clean, sterile slides make 3 touch preparation smears from a cut surface of the biopsy for:
- i) Gram stain: Examine for presence of organisms and cells.
 - ii) Fungifluor stain: Forward to Mycology Section for staining and interpretation.
 - iii) Extra smear held in Mycology Section for special stains.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 41 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Macerate (using a hand grinder) the remaining biopsy tissue using sterile, saline as the diluent.

b) Culture:



| <u>Media</u> | <u>Incubation</u> |
|---|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Fastidious Anaerobe Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |
| Inhibitory Mold Agar (IMA) * | O ₂ , 28°C x 4 weeks |
| Esculin Base Medium (EBM) * | O ₂ , 28°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM)* | O ₂ , 28°C x 4 weeks |
| ----- | |
| If <i>B. cepacia</i> is requested or the specimen is from a patient with Cystic Fibrosis, add: | |
| OF Base, Colistin, Bacitracin & Lactose Agar (OCBL) | O ₂ , 35°C x 5 days |
| Keep the BA, HI and MAC plates for 5 days | |
| ----- | |
| If <i>Nocardia</i> is requested, add: | |
| Sodium Pyruvate Agar (PYRA) | O ₂ , 35°C x 4 weeks |

* Forward inoculated fungal cultures to Mycology for incubation and work-up.

Mouth Swabs

Applies to both C&S and Fungus test requests:

- a) Direct Examination: Gram stain
- b) Culture: Not indicated.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 42 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Nasal Swabs

- a) Direct Examination: Not indicated.
- b) Culture:

| Media | Incubation |
|--|---|
| If C&S is requested: Colistin-Nalidixic Agar (CNA) | CO ₂ , 35°C x 48 hours |
| If <i>N. meningitidis</i> is requested, add : Martin-Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |
| If MRSA is requested: OXOID Denim Blue Agar (DBLUE)* | O ₂ , 37°C x 24 h -in the dark |

Nasopharyngeal Swabs/Auger Suctions for *Bordetella pertussis*

The specimen should be forwarded to the Public Health Laboratory (PHL) for processing.

Oral Abscess Swabs

If an anaerobic swab is received and *Actinomyces* is not requested, process for aerobic culture only. Do not process for anaerobic culture. Document in the “TEST Comment” field in the LIS as “Inappropriate specimen for anaerobic culture, therefore not processed”



- a) Direct Examination: Gram stain
- b) Culture:

| Media | Incubation |
|-----------------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
 Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 43 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

.....

If *Actinomyces* requested or suggested on Gram stain, **add:**

| | |
|----------------------------------|----------------------------------|
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 7 days |
| Kanamycin / Vancomycin Agar (KV) | AnO ₂ , 35°C x 7 days |

Sinus/Antral Specimens

- a) Direct examination:
- i) Gram stain - Examine for and quantitate the presence of pus cells and organisms.
 - ii) Fungifluor - Forward to Mycology Section for staining and interpretation.
- b) Culture:

| Media | Incubation |
|-----------------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | CO ₂ , 35°C x 48 hours |
| Inhibitory Mold agar (IMA)* | O ₂ , 28°C x 3 weeks |

.....

If anaerobic culture requested, **add:**

| | |
|-----------------------------------|------------------------------------|
| Kanamycin Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 48 hours |

*Forward inoculated fungal media to Mycology section for incubation and work-up.

Sputum including Endotracheal tube and Tracheotomy Specimens, Bronchoscopy Aspirates/Washings



- a) Specimen preparation
- For sputum bacterial culture requests, process specimens with [liquefying solution](#) before placing on the WASP for processing. During WASP downtime, use [liquefied specimen](#) to inoculate manually.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 44 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

- Flush small volumes of endotracheal tube (ETT) secretions from suction tubings of neonates with a small volume of WASP solution and transfer to a WASP tube for processing on the WASP
- An ETT tip from an adult is an unacceptable specimen and will not be processed. Notify the ward immediately and issue a report stating “ETT tips are not acceptable for culture”.
- Bronchoscopy Aspirates/Washings from the same patient should be pooled together and processed as a single sample EXCEPT specimen from Trillium Gift of Life Donor (identified by MRN prefixed with “DONOR”). These samples are processed in the same manner as sputum
- For requests for viruses, forward sample to the Virology section for processing.
- Requests for “**STAT**” acid-fast bacilli (AFB) staining must be approved by a medical microbiologist or infection control practitioner before processing in-house. A portion of the sample should also always be forwarded to the Public Health Laboratory for processing.
- If mycobacterium (TB) is requested, forward a portion of the specimen to the Public Health Laboratory (PHL) for processing.
- If Legionella is requested, forward a portion of the specimen to the Public Health Laboratory (PHL) for processing.
- For Bronchoscopy Aspirates/Washings from **Donors or ICU**, process as BAL. Add Test Comment in both Gram and Culture results:
This specimen }BRON
“This specimen was received labelled as bronchial washing but was worked up as BAL”

b) Direct Examination:

Gram smears:

Prepared by the WASP except during downtime.

During WASP downtime select the most purulent portion of the specimen for gram staining. Follow Gram Stain procedure and deliver stained slide to the Gram stain bench.

Approved requests for **STAT** acid fast/TB requests:



- See *Mycobacterium tuberculosis* & Rifampicin resistance PCR by Cepheid GeneXpert.
- For smears forward to PHOL.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 45 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Fungus requests:

Prepare smear for Fungifluor and forward to Mycology Section for staining and interpretation.

c) Culture:

Prepared by the WASP except during downtime.

During WASP downtime, select the most purulent portion of the specimen for plate inoculation.

Specimens will be inoculated onto the following media:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | CO ₂ , 35°C x 48 hours |
| If <i>B. cepacia</i> is requested or specimen is from a patient with Cystic Fibrosis, add: | |
| OF base, colistin, bacitracin & lactose Agar (OCBL) | O ₂ , 35°C x 5 days |
| Keep the BA, HI and MAC plates | CO ₂ , 35°C x 5 days |
| If Nocardia culture is requested, add: | |
| Sodium Pyruvate Agar (PYRA) | O ₂ , 35°C x 4 weeks |
| If fungal culture is requested, add: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 28°C x 4 weeks |
| Esculin Base Medium (EBM)* | O ₂ , 28°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM)* | O ₂ , 28°C x 4 weeks |

* Forward inoculated fungal media to Mycology section for incubation and work-up.

Throat Swabs



- a) Direct Examination: Not indicated for Group A streptococcus, *Neisseria gonorrhoeae* or *Neisseria meningitidis*.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is *Uncontrolled When Printed*.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 46 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

If yeast (thrush) is suspected / requested: Gram stain. Examine for presence of pseudohyphae and/or budding yeast.

If Vincent's angina is suspected / requested: Gram stain. Examine for presence of spirochetes and/or fusiform bacilli and pus cells.

b) Culture:

| Media | Incubation |
|---|---------------------------------------|
| Blood Agar (BA) | AnO ₂ , 35°C x 18-24 hours |
| If <i>Neisseria gonorrhoeae</i> requested, add: | |
| Martin-Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |
| If <i>Corynebacterium diphtheriae</i> requested, forward swab to Public Health Laboratory (PHL) for processing. | |

Note: The ML plate is inoculated by rolling the swab in a "Z" pattern over the medium followed by cross streaking with a sterile loop over the entire plate.

Set up anaerobic jar at 1600 hours and at the end of the evening shift.

Gastric Aspirates/Biopsies (for *Helicobacter pylori*)

a) Direct Examination:

- i) Use one half of the tissue to prepare a smear for Gram stain and to directly inoculate a urea slant.
- ii) Incubate the urea slant at 35°C for 1 – 4 hours.
- iii) Macerate the remaining tissue for culture:

b) Culture:



| Media | Incubation |
|----------------------------|---------------------------------|
| Blood Agar (BA) | Microaerophilic, 35°C x 7 days |
| Campylobacter Agar (CAMPY) | Microaerophilic, 35°C x 7 days |
| Urea (Rapid) | O ₂ , 35°C x 4 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 47 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Gastric Aspirates/Swabs from Neonates or Stillborn

a) Direct Examination: Gram Stain



b) Culture:

| Media | Incubation |
|-----------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 48 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Respiratory Specimen Treatment Solutions:

Liquefying solutions are ready-to-use mucolytic agents for liquefying respiratory specimens prior to planting and streaking for the isolation of bacteria and fungi that cause respiratory tract infections.



Copan Sputum Liquefying (SL) Solution Procedure

1. Take one tube containing 1 ml SL solution and label with a LIS collection label.
2. Using the Copan Transfer Device, transfer an equal amount of specimen (~1mL) to SL solution into the tube.
Note: best results can be achieved by transferring the sample into the solution in a 1:1
3. Dip and rotate the device in the medium 5-10 times.
4. Break off the tip end of the dipper and re-cap tube.
5. Vortex in the test tubes for 30 seconds at 2000/2500 rpm.
6. Leave the test tube at room temperature for at least 15 minutes. A prolonged contact time up to 6 hours doesn't impair the survival of bacteria or fungi contained in the sample and may be necessary for especially viscous samples.
7. Place tube onto the WASP. During WASP downtime for mucoid or aseptically remove aliquots of the liquefied sample and inoculate onto appropriate bacteriology culture medium.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 49 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Sterile Fluids

Cerebrospinal Fluids

NB: If only VDRL is requested, the specimen will **NOT** be cultured. All other requests require C&S. These specimens must be processed immediately.

Unspun sample is required for Biofire meningoencephalitis panel

1. Examine the fluid for colour, xanthochromia, turbidity and the presence of visible blood or clots. Record the observations in the LIS under COMMENT in the SOURCE screen.
2. If any of the following is requested:
 - Cryptococcus or cryptococcal antigen,
 - India ink preparation or
 - Fungal culture
 Perform a cryptococcal antigen test on the supernatant or from the original specimen (see Technical Manual)
3. If bacterial antigens are requested, this test is no longer available.



If >1 mL of specimen is received:

4. Make a direct cytospin smear using 2-3 drops of unspun CSF. Set cytospin to centrifuge for 900 rpm for 3 minutes.
5. Centrifuge the remaining specimen at 3500 rpm x 20 minutes.
6. Transfer the supernatant to a sterile tube which can be used for biochemical, viral, molecular, Cryptococcal antigen or serological tests if required. Label clearly remaining supernatant and refrigerate.
7. Use the sediment to inoculate the culture media listed below. Label clearly any remaining sediment and refrigerate.
8. If Mycobacteria (TB) is requested, reconstitute sediment with 1 mL of supernatant and label as such. Forward the specimen to the Public Health Laboratory (PHL) for processing.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY



*NOTE: This document is **Uncontrolled When Printed**.*
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 50 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

If <1 mL of specimen is received:

4. Make a direct cytospin smear using 2-3 drops of unspun CSF. Set cytospin to centrifuge for 900 rpm for 3 minutes.
5. **DO NOT** centrifuge the specimen. Consider calling other departments for more specimen. Consult a Microbiologist if the sterility of the extra specimen is questionable.
6. If more than one test is ordered, call the physician to prioritize the test ordered.
7. For routine C&S, prepare a Gram smear using a drop of CSF and then directly inoculate the culture media listed below.
8. If Mycobacteria (TB) is requested, forward a portion of the specimen to the Public Health Laboratory (PHL) for processing.
9. If virology is requested, forward a portion of the specimen to the Virology Section for processing.
10. If a swab is received,
 - a. in the Order/Entry Specimen QA Comment field of the LIS as “UNST”
 - b. in the Result Comment field enter “Swab received – insufficient specimen. A negative result may not rule out infection.”
Process the swab as per routine.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 51 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

A. Direct Examination:



Gram stain

Make 2 smears if the specimen is grossly bloody.

B. Culture:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar | CO ₂ , 35°C x 48 hours |
| Chocolate Agar | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth | O ₂ , 35°C x 5 days |
| ----- | |
| If fungus or cryptococcus is requested, add: | |
| Inhibitory Mould Agar ¹ | O ₂ , 28°C x 4 weeks |
| Esculin Base Medium ¹ | O ₂ , 28°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) ¹ | O ₂ , 28°C x 4 weeks |

¹Forward inoculated fungal media to Mycology section for incubation and work-up.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 52 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Other Sterile Fluids – Amniotic, Pleural (Thoracentesis/Empyema), Peritoneal (Ascites), Synovial (Joint), Pericardial, Tympanocentesis, Intraocular, Hydrocele Fluids etc

NB: Chest tube drainage will be treated as a wound specimen (Refer to Wound/Tissues/Aspirates Culture Manual).

1. If a clotted specimen is received, emulsify the coagulated material using a sterile swab or loop.
2. If >1 mL of fluid is received, centrifuge the specimen (3500 rpm x 20 minutes). If <1 mL of specimen is received, do **NOT** centrifuge the specimen.
3. Transfer the supernatant into a sterile container and refrigerate until the final report is issued. Label the container as supernatant.
4. Use the sediment to prepare the smears and inoculate the culture media. If <1 mL of specimen is received, prepare smears and inoculate culture media directly.
5. If viral culture is requested, forward a portion of the specimen to the Virology Section for processing.
6. If Mycobacteria (TB) is requested, forward a portion of the specimen to the Public Health Laboratory (PHL) for processing.
7. If a swab is received,
 - a. In the Order/Entry Specimen QA Comment field of the LIS as “UNST”
 - b. In the Result Comment field enter “Swab received – insufficient specimen. A negative result may not rule out infection.”
 Process the swab as per routine.

A. Direct Examination: Gram stain
Calcofluor white stain (If fungus is requested).



B. Culture:

| Media | Incubation |
|-----------------------|---------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 2 days |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 2 days |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 53 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| | |
|---|------------------------------------|
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) ² | AnO ₂ , 35°C x 48 hours |



For sterile fluids other than Peritoneal (Ascites) fluid:
 Fastidious Anaerobic Broth (THIO) O₂, 35°C x 5 days

For Peritoneal (Ascites) fluid:
 Blood Culture bottles (FA and FN) BacT/Alert 35°C x 5 days

If fungus is requested, **add:**
 Inhibitory Mould Agar (IMA)¹ O₂, 28°C x 3 weeks
 Esculin Base Medium (EBM)¹ O₂, 28°C x 3 weeks

¹ Forward inoculated fungal media to Mycology section for incubation and work-up.

² Not needed for **any** Eye fluids or Tympanocentesis specimens

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 54 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Peritoneal Dialysis Effluent

NB: No more than one dialysis fluid per patient should be processed every other day. If a bag of cloudy fluid is received after a clean one is processed, culture and sensitivity is always done. A portion of the fluid should be sent to Hematology for cell count if requested.

- 1) Examine the bags noting the colour and turbidity of the fluid. Record the observations in the LIS under COMMENT in the SOURCE screen.
- 2) Shake bag well to mix. Working at the sink, disinfect the port with 70% ethanol and let stand for several minutes. Open clamp and run off some of the dialysate.
- 3) Fill a red top vacutainer tube with approximately 5 mL of dialysate and send to Hematology for cell count.
- 4) Fill three 50 mL sterile plastic centrifuge tubes with dialysate. Centrifuge two of the tubes (3500 rpm x 20 minutes). Store the third centrifuge tube at 4°C for at least 7 days. Discard the bag after collecting dialysate.
- 5) Decant the supernatant leaving approximately 10 mL in each tube. Vortex well to resuspend the sediment and pool the two tubes into one. Prepare Gram smear and inoculate culture media as outlined below.

a) **Direct Examination:**

- i) If specimen is cloudy: Gram stain
- ii) If specimen is clear: Gram stain is not needed.

b) **Culture:**



| Media | Incubation |
|------------------------|---|
| If specimen is clear: | |
| Bact/Alert bottles* | Processed as per blood culture protocol |
| ----- | |
| If specimen is cloudy: | |
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| BacT/Alert bottles* | Processed as per blood culture protocol |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 55 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

*Inoculate both an aerobic and anaerobic bottle with 8 to 9 mL of fluid each.

Predialysis Fluid



1. Examine the fluid noting the colour and turbidity. Record the observations in the LIS under COMMENT in the SOURCE screen.
2. Centrifuge specimens (3500 rpm x 20 minutes).
3. Decant supernatant leaving approximately 2-3 mL. Vortex well to resuspend the sediment.
 - a) Direct Examination:
 - i) If specimen is cloudy: Gram stain
 - ii) If specimen is clear: No Gram stain is needed.
 - b) Culture:

| Media | Incubation |
|-----------------------------------|-----------------------------------|
| If specimen is clear: | |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 4 days |
| If specimen is cloudy: | |
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 56 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Bone Marrow (Aspirates or biopsies)

1. If SPS tube received, inoculate to media listed below

For bone marrow specimens not received in SPS

- If <1 ml is received, sterile saline may be added to a final volume of 1 ml so that the media listed below can be inoculated.
- If a clotted specimen is received, emulsify the coagulated material using a sterile swab or loop.
- If only heparinized sample is received and bacterial or fungal C&S is requested, we can process it with the comment:

“This specimen was received in a heparin tube which may suppress bacterial growth leading to a false negative result. Please take this into consideration when interpreting results”

- A portion of **ALL** specimens should be forwarded to the Public Health Laboratory (PHL) for Mycobacteria (TB) culture. (There is no limitation of using heparin for TB per the PHOL website and Canadian Tuberculosis Standards)
- a) Direct Examination: Gram stain.
Calcofluor white stain.

- b) Culture:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin / Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |
| Inhibitory Mould Agar (IMA) ¹ | O ₂ , 28°C x 4 weeks |
| Esculin Base Medium (EBM) ¹ | O ₂ , 28°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) ¹ | O ₂ , 28°C x 4 weeks |



¹Forward inoculated fungal media to the Mycology Section for incubation and work-up.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 57 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Blood, Platelets, & Other Transfusion Products

A sample of the transfusion product from the bag should be processed. Order entry into the LIS

- Open the seal of an unused port and insert a sterile spike
- Remove cap on spike outside the bag and attached a 10mL syringe
- Aspirate sample for processing to direct agar/gram smear and a maximum of 10mL in each blood culture bottle as below:
 - For volumes greater than 20mL
 - Aspirate 10mL of sample, attached needle to syringe, inject 10mL of sample into blood culture bottle. Discard syringe
 - With a new sterile syringe, aspirate another 10mL of sample, attached needle to syringe, inject 10mL of sample into blood culture bottle
 - With a new sterile syringe, aspirate 1mL of sample for inoculation of gram and direct agar
 - For volumes 1ml-20mL
 - Follow the procedure for >20mL using 1mL to inoculate the gram/direct agar and 1mL split into both blood culture bottles
 - For volumes <1mL
 - Inoculate gram / direct agar. If any remaining volume, inoculate a THIO broth.

a) Direct Examination: Gram Stain.

b) Culture:

| Media | Incubation |
|--|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| FAN Aerobic Blood Culture bottle (FO2)* | in BacT/Alert 35°C x 5 days |
| FAN Anaerobic Blood Culture bottle (FN)* | in BacT/Alert 35°C x 5 days |

*Inoculate 10 ml of product from blood bag into each blood culture bottles



*If <10 ml of product received, aseptically inject 10 to 20 cc of Tryptone Soya broth into the blood bag, the bag shaken and the broth re-aspirated and aseptically inoculated into the blood culture bottles

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 58 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Handling of Body Fluids in Blood Culture Bottles

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 59 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Urines

a) Direct Examination:

Gram stain: Not routinely performed. If specifically requested, perform Gram stain directly on unspun specimen.


Fungal stain: Not routinely performed. If dimorphic fungus or cryptococcus specifically requested, prepare 2 slides with spun urine (3500 rpm x 15 min.). Send unstained slides to mycology for staining and interpretation.

Eosinophil stain: Not routinely performed. If requested, prepare smear by cytospin method and fix in absolute alcohol. Stain slide as per Technical manual.

b) Culture:

Automatic WASP culture method: Place Urine ChromID CPS4/CNA bi-plate in WASP culture machine and perform urine culture in WASP.

WASP mix the urine specimen and dip a sterile calibrated 0.001 mL or 0.01 mL loop vertically into the sample. Streak the loop down the center of the plate and then cross-streaked at a 90° angle to the inoculum. (See diagram below).




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.*

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.


Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|  Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 60 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |



Manual method: Mix the urine specimen and dip a sterile calibrated 0.001 mL or 0.01 mL loop vertically into the sample. Streak the loop down the centre of the plate and then cross-streaked at a 90° angle to the inoculum. (See diagram below). If using the Autostreaker, inoculate the plate with a small streak using a sterile calibrated loop (See diagram below). The amount of urine inoculated is dependent on the type of urine specimen that is processed. (See Table below)



i) Autostreaker



ii) Manual inoculation

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 61 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Urine Hold Specimens:

Mount Sinai Hospital urine samples from 11S, 11N, 14N, 14S and 14 Step Down will have urines ordered from Cerner as group test “UCHX1 - Urine Hold”.

- Urine Hold specimens seen or loaded in to the WASP and rejected.
- Go to Order Entry and Receive the specimen.
- Leave Plated blank
- In 2 hours the LIS will automatically add the canned message }UCH1 and give the Test a Preliminary status.

“The majority of positive urine cultures from inpatients without an indwelling urinary catheter represent asymptomatic bacteriuria. If you strongly suspect that your patient has developed a urinary tract infection, please call the Microbiology laboratory within the next 48 hours to have this specimen processed. If a phone call is not received, this specimen will not be processed.”

A. If a request is made to process a urine sample:

- Document call
- In SoftMIC Result Entry, Remove Urine Hold comment “The majority of positive urine...”
- Result test code UCH1 with canned message from Keypad }req - **Request for culture received, results to follow.**
- Finalize Urine Hold test UCHI.
- In order entry use the same LIS order as the Urine Hold order.
- Select “MSU” from order entry keypad and Save.
- Load specimen onto the WASP for routine urine processing.
- The WASP will add the Plated time.

B. If NO request is made, urine hold test will automatically be Finalized.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.**
Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 62 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

TABLE 1. Specimen type, inoculum, media and incubation conditions.



| Media | Incubation |
|---|--|
| <u>Voided urine:</u> use 0.001 mL loop | |
| ChromID BUTI/CNA bi-plate | O ₂ 35°C x 18 – 24 hours |
| <u>In and Out catheter/catheter insertion urine:</u> use 0.001 mL loop | |
| ChromID BUTI/CNA bi-plate | O ₂ 35°C x 18 – 24 hours |
| <u>Aseptically collected urine:</u> use both 0.01 mL and 0.001 mL loop | |
| ChromID BUTI/CNA bi-plate Label the plates with the inoculum volume | O ₂ 35°C x 18 – 24 hours |
| <u>Suprapubic Bladder Aspirate:</u> use both 0.01 mL and 0.001 mL loop | |
| ChromID BUTI/CNA bi-plate Label the plates with the inoculum volume | O ₂ 35°C x 48 hours O ₂ 35°C x 48 hours AnO ₂ 35°C x 48 hours |
| <u>Segmented urines</u> labeled as VB1, VB2, VB3: use both 0.01 mL and 0.001 mL loops (2 sets) | |
| ChromID BUTI/CNA bi-plate Label the plates with the inoculum volume. | O ₂ 35°C x 18 – 24 hours |
| <u>Prostatic Fluid</u> labeled as EPS sent together with VB1, 2 or 3 or labeled as Seminal Fluid sent together with VB1, 2 or 3: use both 0.01 mL and 0.001 mL loops (2 sets) | |
| ChromID BUTI/CNA bi-plate Label the plates with the inoculum volume. | O ₂ 35°C x 18 – 24 hours |
| <u>Urine for routine fungus:</u> use 0.001 mL loop | |
| ChromID BUTI/CNA bi-plate Label the plates with the inoculum volume. | Room temperature x 18-24 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 63 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| Media | Incubation |
|--|--------------------------------|
| Urine for <u>Cryptococcus or specific filamentous fungi</u> : centrifuge urine at 3500 rpm, for 15 minutes | O ₂ 28°C x 3 weeks. |
| EBM | O ₂ 28°C x 3 weeks. |
| IMA | |
| Inoculate 1 drop of sediment per plate. Forward plates to Mycology | |

Note: For specimens processed after 1600 hours, place the plates in a separate basket in the incubator.

Special Requests:

Eosinophil Stain –

Not routinely performed. If requested, prepare smear by cytospin method and fix in absolute alcohol. Stain slide as per Technical manual.

Bacterial Latex Agglutination –

Not done due to poor sensitivity and specificity.

Anaerobes –

Appropriate specimen is bladder suprapubic aspirate – see above Table. Process these specimens immediately. If requested on other urines, consult the Medical Microbiologist.

Chlamydia Detection –

Send the specimen to Virology immediately.

Legionella Antigen Detection –

Send the specimen to the Provincial Health Laboratory.

Leptospira Detection –

Send the specimen to the Provincial Health Laboratory.

Cryptococcus/Systemic Fungi –

See above Table.

TB Culture –

Send the specimen to the Provincial Health Laboratory.

Viral Culture –

Send the specimen to Virology immediately.



Parasitology – *Schistosomiasis* –

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 64 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Send the specimen to Parasitology immediately.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 65 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Wounds/Tissues/Aspirates

Duodenal or Small Bowel Aspirate/Swab/Biopsy



- a) Direct Examination: Gram stain not performed.
- b) Culture:
 - i) Duodenal or Small Bowel Aspirates
 Process Duodenal and Small Bowel aspirates for O&P only. Forward the specimen to the Parasitology section immediately. If a delay is anticipated and the specimen is not received in SAF, transfer the specimen to SAF. DO NOT process these specimens routinely for bacterial culture unless requested.
 - ii) Duodenal or Small Bowel Swab

| Media | Incubation |
|----------------------------------|---------------------------------------|
| Blood Agar (BA) | O ₂ , 35°C x 18 - 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 18 - 24 hours |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin / Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 66 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Intraoperative/Interventional Swabs

a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.

If fungus is requested, **add:**



Calcofluor white stain - Refer to Mycology Manual.

b) Culture:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) ¹ | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) ¹ | O ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) ¹ | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) ¹ | O ₂ , 35°C x 5 days |
| If anaerobic swab is received, add: | |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| If fungus is requested, add: | |
| Inhibitory Mold Agar (IMA) [*] | O ₂ , 30°C x 4 weeks |
| Esculin Base Medium (EBM) [*] | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) [*] | O ₂ , 30°C x 4 weeks |

¹ If organisms were seen in direct Gram stain and cultures yield no corresponding growth after 48 hours of incubation, check direct Gram stain (if discrepant compared to original report, check with the Charge technologist), and re-incubate all aerobic plates and broth for 7 days. If there is no evidence of corresponding growth after 7 days, subculture the THIO to CHOC and BRUC.

* Forward fungus culture media to Mycology section for incubation and processing.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 67 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Wound/Abscess Swabs and Drain

- a) Direct Examination: Gram stain - Quantitate the presence of pus cells, squamous epithelial cells, and organisms.
- b) Culture:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Colistin Nalidixic Acid Agar (CNA) | O ₂ , 35°C x 48 hours |
| For chest tube drainage and tracheal swabs, add: | |
| Haemophilus Isolation Medium (HI) | CO ₂ , 35°C x 48 hours |
| If anaerobic transport media is received, add: | |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin / Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |

Bite Wound Swabs



- a) Direct Examination: Gram stain – Quantitate the presence of pus cells, squamous epithelial cells, and organisms.
- b) Culture:

| Media | Incubation |
|--|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| If anaerobic transport media received, add: | |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin / Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 68 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Intravenous & Central Line Catheter Exit Site Swabs

a) Direct Examination: Not indicated.

b) Culture:

| Media | Incubation |
|----------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |

Intraoperative/Interventional Abscess (Pus, Cyst Fluid or Aspirate)

1. If >1 ml of thin fluid is received, centrifuge specimen at 3500 rpm for 20 minutes. For purulent and thick specimens or if <1 ml is received, centrifugation is not required.
2. If parasitology is requested, send a portion of the fresh specimen to the Public Health Laboratories (PHL) on Mondays to Fridays. On weekends and holidays, mix an equal volume of specimen with SAF and send to the Public Health Laboratories (PHL). Note on the specimen label that it has been mixed with SAF.
3. If TB culture is requested, send a portion of the specimen to the Public Health Laboratory (PHL) for processing.



a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.

If fungus is requested, **add:**
Calcofluor white stain - Refer to Mycology Manual.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 69 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

b) Culture:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) ¹ | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) ¹ | O ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) ¹ | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) ¹ | O ₂ , 35°C x 5 days |
| If fungus is requested, add: | |
| Inhibitory Mold Agar (IMA) [*] | O ₂ , 30°C x 4 weeks |
| Esculin Base Medium (EBM) [*] | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) [*] | O ₂ , 30°C x 4 weeks |

¹ If organisms were seen in direct Gram stain and cultures yield no corresponding growth after 48 hours of incubation, check direct Gram stain (if discrepant compared to original report, check with the Charge technologist), and re-incubate all aerobic plates and broth for 7 days. If there is no evidence of corresponding growth after 7 days, subculture the THIO to CHOC and BRUC.

* Forward fungus culture media to Mycology section for incubation and processing.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 70 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Pus & Abscess Material (other than Intraoperative/Interventional, Rectal or Bartholin)

- a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.
 Modified acid fast stain - If Actinomyces or Nocardia is requested or suggested on Gram stain.
 Calcofluor white stain - If fungus is requested. (Refer to Mycology Manual).
- b) Culture:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) ¹ | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) ¹ | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC) ² | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) ² | AnO ₂ , 35°C x 48 hours |
| <hr/> | |
| If Nocardia is requested, add: | |
| Na Pyruvate Agar (PYRU)* | O ₂ , 35°C x 4 weeks |
| AND fungus media below | |
| <hr/> | |
| If fungus culture is requested, add: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM)* | O ₂ , 30°C x 4 weeks |



*Forward the fungus culture media and PYRU to the Mycology section for incubation and work-up.

- NOTE:**
1. If Nocardia is requested, send the BA and CHOC plates to mycology after 48 hours incubation. The plates will be incubated in mycology for 4 weeks.
 2. If Actinomyces is requested, set up a second set of anaerobic media to be incubated for 7 days before opening jar.
 3. If Nocardia or Actinomyces is suggested on Gram stain, set up a second set of anaerobic media to be incubated for 7 days before opening jar and send BA and CHOC plates to Mycology after 48 hours incubation.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 71 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Rectal Abscess



- a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.
- b) Culture:

| Media | Incubation |
|------------------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Colistin Nalidixic Acid Agar (CNA) | O ₂ , 35°C x 48 hours |

Bartholin's Abscess Swab/Aspirate

- a) Direct Examination: Gram stain. - Quantitate the presence of pus cells, squamous epithelial cells, and organisms.
- b) Culture:

| Media | Incubation |
|-------------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| Martin –Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 72 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

TISSUES, PROSTHETIC DEVICES, AND AUTOPSY SPECIMENS

Tissues/Biopsies (other than skin or transplant tissues)

Processing of Specimen:

1. When sample size allows and a fungal culture is requested, macerate the tissue using scalpels before adding to stomacher. Small pieces of tissue should be inserted into each of the fungal plates. If small sample is received then macerate the tissue using only a stomacher. Tissue around bone specimen should be removed and macerate. Dry bone can be inoculated directly into Fastidious Anaerobic Broth and not macerated. If it is a mixture of bone and tissue, macerate the tissue as usual and submerge the bone into Fastidious Anaerobic Broth
2. If fungal stains are requested, prepare 3 slides from the macerated material: one for Gram stain, one for Calcofluor white stain and one unstained (Stored in the “extra smear” slide box).
3. Fungus culture is **NOT** set up for wound debridement tissue, joint capsules, gas gangrene tissue and necrotizing fasciitis tissue unless specifically requested.
4. Send a portion of **ALL** macerated specimens or tissues to the Public Health Laboratory (PHL) for TB **except** for wound debridement tissue, joint capsules, gas gangrene tissue, and necrotizing fasciitis tissue unless specifically requested. If viral isolation is requested, send a portion of the specimen to the Virology section for processing.
5. Inoculate the following media with the remaining sample:



| Media | Incubation |
|---|--|
| Blood Agar (BA) ¹ | CO ₂ , 35°C x 48 hours ¹ |
| MacConkey Agar (MAC) ¹ | O ₂ , 35°C x 48 hours ¹ |
| Chocolate Agar (CHOC) ¹ | CO ₂ , 35°C x 48 hours ¹ |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) ¹ | O ₂ , 35°C x 5 days ¹ |
| If Fungus is requested, add: | |
| Inhibitory Mold Agar (IMA) [*] | O ₂ , 30°C x 4 weeks |
| Esculin Base Medium (EBM) [*] | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM) [*] | O ₂ , 30°C x 4 weeks |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 73 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

¹ If organisms were seen in direct Gram stain and cultures yield no corresponding growth after 48 hours of incubation, check direct Gram stain (if discrepant compared to original report, check with the Charge technologist), and re-incubate all aerobic plates and broth for 7 days. If there is no evidence of corresponding growth after 7 days, subculture the THIO to CHOC and BRUC.

* Forward the fungal culture media to the Mycology section for incubation and work-up.

Skin Biopsies

Processing of Specimen:

1. Macerate the entire specimen with a small volume of sterile saline from the original container using scalpels.
2. Prepare 3 slides from the ground material: one for Gram stain, one for Calcofluor white stain and one unstained (stored in the “extra smear” slide box).
3. If TB has been requested, forward half the specimen to the Public Health Laboratory (PHL). If viral isolation is requested, a portion of the ground material should be given to the Virology section for processing. If parasitology is requested, forward a portion of the specimen (before macerating) to the Public Health Laboratory (PHL).

a) Direct Examination: Gram stain: Quantitate the presence of pus cells and organisms.

If fungus is requested, **add:**

Calcofluor white stain - Refer to Mycology Manual.

b) Culture:



| Media | Incubation |
|--|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| If fungus is requested, add:: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide | O ₂ , 30°C x 4 weeks |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 74 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |



(BHIM)*

* Forward the fungus culture media to the Mycology section for incubation and work-up.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 75 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Transplant Specimens - Bone Graft & Cadaver Fascia/Tissue/ Swab Specimens/Donor Amniotic Fluid/Membrane; Donor Corneal Ring Material

a) Direct Examination: Not indicated

b) Culture:

| Media | Incubation |
|------------------------------------|--------------------------------|
| Fastidious Anaerobic Broth (THIO)* | O ₂ , 35°C x 7 days |

* A separate THIO should be inoculated for each specimen / swab received.



Prosthetic Devices (e.g. Pacemaker Wire, Dacron Graft, Prosthetic Valve, Coronary/Vascular Stents)

Please note: for Urinary, urethral, prostatic and Bile duct stents, see [Bile and Bile Stent section](#)

a) Direct Examination: Not indicated.

b) Culture:

| Media | Incubation |
|-----------------------------------|--------------------------------|
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 76 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Autopsy Specimens

- a) Direct Examination: Gram stain - Quantitate the presence of pus cells, squamous epithelial cells, and organisms.
- b) Culture:
 1. Macerate the tissue using scalpels. Tissue around bone specimen should be removed and macerate. Dry bone can be inoculated directly into Fastidious Anaerobic Broth and not macerated.
 2. Prepare 2 slides from the macerated material: one for a Gram stain and a second unstained. (Stored in the “extra smear” slide box).
 3. A portion of all autopsy lung tissue (except newborns) is to be sent to the Public Health Lab (PHL) for *Legionella* detection and TB culture.
 4. Inoculate the following media with the remaining sample:

| Media | Incubation |
|--|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| Colistin Nalidixic Acid Agar (CAN) | O ₂ , 35°C x 48 hours |
| For all lung tissue or if fungal culture is requested, add: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 3 weeks |

* Forward the fungus culture media to the Mycology section for incubation and work-up.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 77 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

CATHETER

Intravascular Catheter Tips

- a) Direct Examination: Not indicated.
- b) Culture:



| Media | Incubation |
|-----------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |

Roll the segment back and forth 4 times across the surface of the BA using sterile forceps.

Peritoneal Dialysis Catheter/Canula

- a) Direct Examination: Not indicated.
- b) Culture:

| Media | Incubation |
|-----------------------------------|--------------------------------|
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 78 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

BILE SPECIMENS

Bile and Bile Stents (including urinary, urethral, prostatic stents)

Inoculate one drop of specimen onto each media and smear. If stent is received, “wash” the specimen in 0.5 ml of sterile saline and inoculate the media with the solution.

a) Direct Examination: Gram stain – Examine for the presence of pus cells and organisms.

b) Culture:

| Media | Incubation |
|--|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| If anaerobic culture is requested or bile is collected by PTC, add: | |
| Fastidious Anaerobic Agar (BRUC) | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |

MISCELLANEOUS FLUID SPECIMENS

Breast Milk – Clinical samples

a) Direct Examination: Not required

b) Culture:

| Media | Incubation |
|----------------------|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |


UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

Total Parental Nutrition (TPN)



1. Determine if tubings are clipped at positions A, B & C at bedside (See diagram below).
If not, clip at these positions and note on work card.
2. Using aseptic technique, remove ~ 1 mL each from **1** & **2** and culture as outlined below.
3. Save the entire TPN set up for Infection Control.



b) Culture:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 3 weeks |
| IMA with sterile olive oil overlay (olive oil is stored in media room)* | O ₂ , 30°C x 1 week |

*Forward these plates to the Mycology section for incubation and work-up.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 80 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

EAR SPECIMENS

Ear Swab



- a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms. Calcofluor white stain (If fungus is requested). - Refer to Mycology Manual.
- b) Culture:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| Colistin Nalidixic Acid Agar (CNA) | O ₂ , 35°C x 48 hours |
| If fungus culture is requested, add: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 3 weeks |

* Forward the fungal culture media to the Mycology section for incubation and work-up.

Tympanocentesis Fluid

See [Other Sterile Fluids – Amniotic, Pleural \(Thoracentesis/Empyema\), Peritoneal \(Ascites\), Synovial \(Joint\), Pericardial, Tympanocentesis, Intraocular, Hydrocele Fluids etc](#)

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 81 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

EYE SPECIMENS

Eye / Conjunctival / Lid Swabs

Note: If pre-inoculated culture plates are received, these should be incubated as listed below. No Gram stain will be performed.

- a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.
- b) Culture:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| For all neonates ≤ 1 week of age, or if <i>N. gonorrhoeae</i> is requested, add: | |
| Martin-Lewis Agar (ML) | CO ₂ , 35°C x 72 hours |

Eye / Corneal Scrapings

Note: If pre-inoculated plates are received and no smear or additional specimen is received, direct smear stains will not be performed. Result in the “TEST” field in the LIS as “No smear received, test not performed.”



- a) Direct Examination: Gram stain – Examine for the presence of pus cells and organisms. Calcofluor white stain (if two smears are provided). Refer to Mycology Manual.
- b) Culture:

| Media | Incubation |
|-----------------------------------|---------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 4 days |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 4 days |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 3 weeks |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.
 Any documents appearing in paper form that do not state “CONTROLLED COPY” in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\



| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 82 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

*Forward the fungal culture media to the Mycology section for incubation and workup.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 83 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Intraocular Aspirates



See [Other Sterile Fluids – Amniotic, Pleural \(Thoracentesis/Emphyema\), Peritoneal \(Ascites\), Synovial \(Joint\), Pericardial, Tympanocentesis, Intraocular, Hydrocele Fluids etc](#)

Lacrimal (Tear Duct) Stone / Secretions

- a) Direct examination: Crush specimen on a glass slide to obtain a thin smear for Gram stain.
- a) Culture: Crush specimen using a sterile wood applicator stick or urine loop before planting onto the following media:

| Media | Incubation |
|---|------------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Agar (BRUC) ¹ | AnO ₂ , 35°C x 48 hours |
| Kanamycin/Vancomycin Agar (KV) ¹ | AnO ₂ , 35°C x 48 hours |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 5 days |

¹If Actinomyces is suggested on direct Gram stain, set up a second set of anaerobic media to be incubated for 7 days before opening jar.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 84 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

FACIAL SPECIMENS

Facial Swabs



a) Direct Examination: Gram stain – Quantitate the presence of pus cells and organisms.
Calcofluor white stain (If fungus is requested).

b) Culture:

| Media | Incubation |
|---|-----------------------------------|
| Blood Agar (BA) | CO ₂ , 35°C x 48 hours |
| Chocolate Agar (CHOC) | CO ₂ , 35°C x 48 hours |
| MacConkey Agar (MAC) | O ₂ , 35°C x 48 hours |
| If Actinomyces is requested or suggested on Gram stain or an anaerobic swab collected or thick pus is received, add: | |
| Fastidious Anaerobic Agar (BRUC) ¹ | AnO ₂ , 35°C x 7 days |
| Kanamycin/Vancomycin (KV) ¹ | AnO ₂ , 35°C x 7 days |
| Fastidious Anaerobic Broth (THIO) | O ₂ , 35°C x 7 days |
| If fungus culture is requested, add: | |
| Inhibitory Mold Agar (IMA)* | O ₂ , 30°C x 4 weeks |
| Brain Heart Infusion Agar with 5% Sheep Blood, Gentamicin, Chloramphenicol, Cyclohexamide (BHIM)* | O ₂ , 30°C x 4 weeks |

¹If Actinomyces is requested or suggested on direct Gram stain, set up a second set of anaerobic media to be incubated for 7 days before opening jar.

*Forward the fungal culture media to the Mycology section for incubation and work-up.

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 85 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

Record of Edited Revisions

Manual Section Name: **Specimen Processing Manual**



| Page Number / Item | Date of Revision | Signature of Approval |
|--|--------------------|-----------------------|
| Faces/Rectal Swab C&S – <i>Pleisiomonas</i> and <i>Aeromonas</i> added | April 10, 2006 | Dr. T. Mazzulli |
| CIN agar for <i>Yersinia</i> – changed to 48 hour incubation | April 10, 2006 | Dr. T. Mazzulli |
| Annual Review | April 10, 2006 | Dr. T. Mazzulli |
| CSF processing – added cytospin rpm and time | August 2, 2006 | Dr. T. Mazzulli |
| Annual Review | April 10, 2007 | Dr. T. Mazzulli |
| BAL smear – change smear made from spun sediment before reconstituting specimen for planting | June 29, 2007 | Dr. T. Mazzulli |
| New Mycology planting media table | June 29, 2007 | Dr. T. Mazzulli |
| MRSA screening media change | June 29, 2007 | Dr. T. Mazzulli |
| Annual Review | June 23, 2008 | Dr. T. Mazzulli |
| Sterile fluids incubation time – changed to 5 days | July 27, 2009 | Dr. T. Mazzulli |
| Miscellaneous Bench specimens added | July 27, 2009 | Dr. T. Mazzulli |
| Annual Review | July 27, 2009 | Dr. T. Mazzulli |
| Added Chocolate Agar to GC cultures | April 5, 2010 | Dr. T. Mazzulli |
| Annual Review | April 5, 2010 | Dr. T. Mazzulli |
| Added 0.01uL loop to suprapubic urine and removed BRUC | June 18, 2010 | Dr. T. Mazzulli |
| Tissues, macerate with scalpels only | June 18, 2010 | Dr. T. Mazzulli |
| Annual Review | May 31, 2011 | Dr. T. Mazzulli |
| BAL change to planting bacteriology culture with calibrated loop | November 7, 2011 | Dr. T. Mazzulli |
| Moved Bronchoscopy Aspirates/Washings to Sputum section | November 7, 2011 | Dr. T. Mazzulli |
| Annual Review | May 31, 2012 | Dr. T. Mazzulli |
| Modified planting volume into BHI broth for VRE/MRSA | August 28, 2012 | Dr. T. Mazzulli |
| Modified Sterile Specimens incubation to 48 hours for aerobic plates | September 24, 2012 | Dr. T. Mazzulli |
| Added ESwab handling | September 24, 2012 | Dr. T. Mazzulli |
| Annual Review | May 31, 2013 | Dr. T. Mazzulli |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 86 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |



| Page Number / Item | Date of Revision | Signature of Approval |
|--|-------------------|-----------------------|
| Changed Seminal fluid inoculation to CHOC and ML by drops | April 16, 2014 | Dr. T. Mazzulli |
| Replace THIO with FA and FN bottles for Peritoneal/Ascites Fluid | April 16, 2014 | Dr. T. Mazzulli |
| Annual Review | April 16, 2014 | Dr. T. Mazzulli |
| Enterococcus agar changed to BVRE in dark in C&S Added <i>CRE Screen</i> Added SL solution procedure under BAL/Sputum section | November 12, 2014 | Dr. T. Mazzulli |
| ChromID CPS4/CNA bi-plate for urine culture Tissue planting procedure | March 8, 2015 | Dr. T. Mazzulli |
| Addition of link in urine section to Urine Hold procedure | April 30, 2015 | Dr. T. Mazzulli |
| Updated Bloods, Platelets, & Other Transfusion Products processing procedure for direct examination and culture , (page 39) | May 26, 2015 | Dr. T. Mazzulli |
| GBS section inserted instructions for processing subcultures. Fixed formatting and TOC headers for eye specimens | August 18, 2015 | Dr. T. Mazzulli |
| GBS updated instructions: -specified to load in BLUE palette(s) labelled 'BLU URISTRAW' -removed step to ensure tiles are at the bottom of tubes -wasp program changed to all capital letter (GBS-CNA) -Added GBS processing instructions | September 4, 2015 | Dr. T. Mazzulli |
| Updated GBS procedure Technologist loads GBS subs | September 9, 2015 | Dr. T. Mazzulli |
| Added STAT AFB by Fluorochrome stain method p.31 | October 20, 2015 | Dr. T. Mazzulli |
| Annual Reivew | May 26, 2016 | Dr. T. Mazzulli |
| Updated MSH logo in header Removed from fungal smears/culture processing: For special requests, send skin, hair and nail specimens to mycology for processing. – all sent to PHOL | July 4, 2016 | Dr. T. Mazzulli |
| Streaking pattern for bronchial brush added. | July 16, 2016 | Dr. T. Mazzulli |
| CIN added to <i>Plesiomonas</i> set up | January 4, 2017 | Dr. T. Mazzulli |
| Asymptomatic urine hold study 68 link removed. Study 68 removed and incorporated into procedure manual. | January 31, 2017 | Dr. T. Mazzulli |
| Reviewed BAL and Sputum sections. Updated SL Liquefying solutions to use transfer device. Removed | July 19, 2017 | Dr. T. Mazzulli |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

NOTE: This document is Uncontrolled When Printed.

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 87 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| Page Number / Item | Date of Revision | Signature of Approval |
|---|--------------------|-----------------------|
| EDM solution. | | |
| Updated positive Blood Culture subculture by WASP | December 5, 2017 | Dr. T. Mazzulli |
| Updated urine hold procedure. | January 25, 2018 | Dr. T. Mazzulli |
| Updated urine hold procedure. | March 6, 2018 | Dr. T. Mazzulli |
| Addition of secondary specimen labeling. All aliquoting procedures transferred and summarized in introduction. | June 13, 2018 | Dr. T. Mazzulli |
| Annual Review Addition of specimen safe/ethical handling in introduction. | September 13, 2018 | Dr. T. Mazzulli |
| Enterics update to WASPlab procedure | February 5, 2019 | Dr. T. Mazzulli |
| General Processing Instructions added. | March 15, 2019 | Dr. T. Mazzulli |

Full document review included in all updates. Bi-annual review conducted when no revision had been made within 2 years.



| Page Number / Item | Date of Revision | Edited by: |
|--|------------------|----------------|
| p.34 Mycology – removed grinding and stomacher use for tissues processed for fungi. Tissue to be (only) gently minced to maintain fungi viability. | January 6, 2020 | Jessica Bourke |
| p.56 Blood, Platelets, & Other Transfusion Products – addition of detailed procedure for processing bags. | January 6, 2020 | Jessica Bourke |
| Included Stents as example in prosthetic device, added note to process bile stents per bile procedure | March 19, 2021 | Wayne Chiu |
| Minor formatting changes | April 11, 2021 | Jessica Bourke |
| Added info regarding urinary/urethral/prostatic vs coronary/vascular stents | April 19, 2021 | Wayne Chiu |
| Added comment for heparinized bone marrow | May 14, 2021 | Wayne Chiu |
| Added note regarding BVRE for flag VDE | June 7, 2021 | Wayne Chiu |
| Fixed subculture volume for BC+ tube | July 8, 2021 | Wayne Chiu |
| Clarified agar name for ESBL/CPE/Resistant pseudo screen culture | Nov 16, 2021 | Wayne Chiu |
| For manual processing of blood culture bottles, ensure BCBC media is check\marked under R | May 13, 2022 | Wayne Chiu |
| Updated bone marrow section with SPS tube | June 15, 2022 | Wayne Chiu |
| Added Alternate SPS Blood Culture Protocol for Dimorphic Fungi | July 15, 2022 | Qin Liu |
| Added note to CSF section – unspun sample required for | Oct 4, 2022 | Wayne Chiu |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\

| | | |
|--|---|------------------|
|   Department of Microbiology Quality Manual | Policy # MI_SM_PROC | Page 88 of 88 |
| | Version: 3.10 CURRENT | |
| Section: Specimen Processing Procedures | Subject Title: Specimen Processing Procedure | |

| | | |
|--|----------------|---------|
| biofire Removed isolator 10 tube for blood, discontinued by vendor | | |
| Pg.44 process bronchial washing from ICU/Donor as BAL and add the following comment: }BRONC | August 9, 2023 | Qin Liu |

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Specimen Processing Procedures\