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
Agilent xCELLigence RTCA 手册

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图片摘自 *Journal of Clinical Microbiology*, volume 48, Ryder, A. B. et al. "Clostridium difficile Infections by Quantitative Detection of tcdB Toxin by Use of a Real-Time Cell Analysis System," pages 4129–4134. 版权所有 2010，经美国微生物学会授权。

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图片摘自 *Journal of Clinical Microbiology*, volume 48, Ryder, A. B. et al. "Clostridium difficile Infections by Quantitative Detection of tcdB Toxin by Use of a Real-Time Cell Analysis System," pages 4129–4134. 版权所有 2010, 经美国微生物学会授权。

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1. Jin, D. *et al.* Quantitative Detection of *Vibrio cholerae* Toxin by Real-Time And Dynamic Cytotoxicity Monitoring. *J. Clin. Microbiol.* **2013** Dec, 51(12), 3968–74
2. Real-Time Cell-Based Toxicology Testing Might Replace Animal Testing for Product Release and Drug Safety. *Biochemica* **2008**, 11–13
3. Ryder, A. B. *et al.* Assessment of *Clostridium difficile* Infections by Quantitative Detection of tcdB Toxin by Use of a Real-Time Cell Analysis System. *J. Clin. Microbiol.* **2010**, 48, 4129–4134
4. Huang, B. *et al.* Real-Time Cellular Analysis Coupled with a Specimen Enrichment Accurately Detects and Quantifies *Clostridium difficile* Toxins in Stool. *J. Clin. Microbiol.* **2014**, 52, 1105–1111
5. Pollock, N. R. Ultrasensitive Detection and Quantification of Toxins for Optimized Diagnosis of *Clostridium difficile* Infection. *J. Clin. Microbiol.* **2016**, 54, 259–264
6. Planche, T. *et al.* Diagnosis of *Clostridium difficile* Infection by Toxin Detection Kits: A Systematic Review. *Lancet Infect. Dis.* **2008**, 8, 777–784
7. Polage, C. R. *et al.* Overdiagnosis of *Clostridium difficile* Infection in the Molecular Test Era. *JAMA Intern. Med.* **2015**, 175, 1792–1801

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转载自 *Journal of Microbiological Methods*, volume 84 (1), Slaninak, H. et al. "Real-Time Impedance Analysis of Host Cell Response to Meningococcal Infection," pages 101–8. 版权所有 2011, 经 Elsevier 授权。

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1. Coureuil, M. *et al.* Mechanism of Meningeal Invasion by *Neisseria meningitidis*. *Virulence* **2012**, Mar-Apr, 3(2), 164–72
2. Slanina, H. *et al.* Real-Time Impedance Analysis of Host Cell Response to Meningococcal Infection. *J. Microbiol. Methods* **2011** Jan, 84(1), 101–8
3. Mou, X. *et al.* Phenotypic Pattern-Based Assay for Dynamically Monitoring Host Cellular Responses to *Salmonella* Infections. *PLoS One* **2011**, 6(11), e26544
4. Miranda-CasoLuengo, A. A. *et al.* A Real-Time Impedance Based Method to Assess *Rhodococcus equi* Virulence. *PLoS One* **2013**, 8(3), e60612
5. Hidalgo-Cantabrana, C. *et al.* Effect of Bacteria Used in Food Industry on the Proliferation and Cytokine Production of Epithelial Intestinal Cellular Lines. *J. Functional Foods* **2014**, 6, 348–355

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1. Junka, A. F. *et al.* Use of the Real Time xCelligence System for Purposes of Medical Microbiology. *Polish Journal of Microbiol.* **2012**, 61(3), 191–197
2. Cihalova, K. *et al.* *Staphylococcus aureus* and MRSA Growth and Biofilm Formation after Treatment with Antibiotics and SeNPs. *Int. J. Mol. Sci.* **2015** Oct 16, 16(10), 24656–72
3. Gutiérrez, D. *et al.* Monitoring in Real Time the Formation and Removal of Biofilms from Clinical Related Pathogens Using an Impedance-Based Technology. *PLoS One* **2016**, Oct 3, 11(10), e0163966
4. Ferrer, M. D. *et al.* Effect of Antibiotics on Biofilm Inhibition and Induction measured by Real-Time Cell Analysis. *J. Appl. Microbiol.* **2016** Dec 8
5. Wang, T.; Su, J. *Bacillus subtilis* from Soybean Food Shows Antimicrobial Activity for Multidrug-Resistant *Acinetobacter baumannii* by Affecting the *adeS* Gene. *J. Microbiol. Biotechnol.* **2016** Dec 28, 26(12), 2043–2050
6. van Duuren, J. B. J. H. *et al.* Use of Single-Frequency Impedance Spectroscopy to Characterize the Growth Dynamics of Biofilm Formation in *Pseudomonas aeruginosa*. *Sci. Rep.* **2017** Jul 12, 7(1), 5223
7. Gutiérrez, D. *et al.* Real-Time Assessment of *Staphylococcus aureus* Biofilm Disruption by Phage-Derived Proteins. *Front Microbiol.* **2017** Aug 24, 8, 1632

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1. Fang, Y. *et al.* Real-Time Monitoring of Flavivirus Induced Cytopathogenesis Using Cell Electric Impedance Technology. *J. Virol. Methods*. **2011** May, 173(2), 251–8
2. Witkowski, P. T. *et al.* Cellular Impedance Measurement as a New Tool for Poxvirus Titration, Antibody Neutralization Testing and Evaluation of Antiviral Substances. *Biochem. Biophys. Res. Commun.* **2010** Oct 8, 401(1), 37–41
3. Charretier, C. *et al.* Robust Real-Time Cell Analysis Method for Determining Viral Infectious Titers During Development of a Viral Vaccine Production Process. *J. Virol. Methods* **2017** Nov 14, 252, 57–64
4. Charretier, C. *et al.* Robust Real-Time Cell Analysis Method for Determining Viral Infectious Titers During Development of a Viral Vaccine Production Process. *J. Virol. Methods* **2018** Feb, 252, 57–64 (republished)
5. Burmakina, G. *et al.* Real-Time Analysis of the Cytopathic Effect of African Swine Fever Virus. *J. Virol. Methods* **2018** Jul, 257, 58–61
6. Lebourgeois, S. *et al.* Development of a Real-Time Cell Analysis (RTCA) Method as a Fast and Accurate Method for Detecting Infectious Particles of the Adapted Strain of Hepatitis A Virus. *Front. Cell Infect. Microbiol.* **2018** Sep 25, 8, 335

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图片摘自 *Asia Pacific Biotech News*, volume 14(10), Lu, H. et al. "Label-free Real-time Cell Based Assay System for Evaluating H1N1 Vaccination Success," pages 31–32. 版权所有 2010, 经 Asia Pacific Biotech News 授权。

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1. Fang, Y. et al. Real-time Monitoring of Flavivirus Induced Cytopathogenesis Using Cell Electric Impedance Technology. *J. Virol. Methods*. **2011** May, 173(2), 251–8
2. Lu, H.; Xu, X. Label-Free Real-Time Cell Based Assay System for Evaluating H1N1 Vaccination Success. *Asiabiochem*. **2010**, 14(10), 31–32
3. Teng, Z. et al. Real-Time Cell Analysis—A New Method For Dynamic, Quantitative Measurement of Infectious Viruses and Antiserum Neutralizing Activity. *J. Virol. Methods*. **2013** Nov, 193(2), 364–70
4. Burmakina, G. et al. Real-Time Analysis of the Cytopathic Effect of African Swine Fever Virus. *J. Virol. Methods* **2018** Jul, 257, 58–61
5. Tian, D. et al. Novel, Real-Time Cell Analysis for Measuring Viral Cytopathogenesis and the Efficacy of Neutralizing Antibodies to the 2009 Influenza A (H1N1) virus. *PLoS One* **2012**, 7(2), e31965

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图片经授权摘自 *ACS Infectious Diseases Volume 3 (6)*, Prasad, V. et al. "Cell Cycle-Dependent Kinase Cdk9 Is a Postexposure Drug Target Against Human Adenoviruses," pages 398–405. 版权所有 2017, 美国化学学会

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数据引自 *Journal of Clinical Microbiology*, volume 54 (8), Piret, J. et al. "Novel Method Based on Real-Time Cell Analysis for Drug Susceptibility Testing of Herpes Simplex Virus and Human Cytomegalovirus." 版权所有 2016，经美国微生物学会授权。

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1. Watterson, D. *et al.* A Generic Screening Platform for Inhibitors of Virus Induced Cell Fusion Using Cellular Electrical Impedance. *Sci. Rep.* **2016** Mar 15, 6, 22791
2. Zandi, K. A Real-Time Cell Analyzing Assay for Identification of Novel Antiviral Compounds against Chikungunya Virus. *Methods Mol. Biol.* **2016**, 1426, 255–62
3. Sharma, B. N. *et al.* Antiviral Effects of Artesunate on Polyomavirus BK Replication in Primary Human Kidney Cell. *Antimicrob. Agents Chemother.* **2014**, 58(1), 279–89
4. Witkowski, P. T. *et al.* Cellular Impedance Measurement as a New Tool for Poxvirus Titration, Antibody Neutralization Testing and Evaluation of Antiviral Substances. *Biochem. Biophys. Res. Commun.* **2010** Oct 8, 401(1), 37–41
5. Piret, J. *et al.* Novel Method Based on Real-Time Cell Analysis for Drug Susceptibility Testing of Herpes Simplex Virus and Human Cytomegalovirus. *J. Clin. Microbiol.* **2016** Aug, 54(8), 2120–7
6. Cymerys, J. *et al.* Primary Cultures of Murine Neurons for Studying Herpes Simplex Virus 1 Infection and Its Inhibition by Antivirals. *Acta Virol.* **2013**, 57(3), 339–45
7. Prasad, V. *et al.* Cell Cycle-Dependent Kinase Cdk9 Is a Postexposure Drug Target Against Human Adenoviruses. *ACS Infect. Dis.* **2017** Jun 9, 3(6), 398–405

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图片摘自 *Public Library of Science Neglected Tropical Diseases*, volume 4 (11), Smout, M. J. et al. "Novel High Throughput Assay for Anthelmintic Drug Screening and Resistance Diagnosis by Real-Time Monitoring of Parasite Motility" 2010.

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图片摘自 *Public Library of Science Neglected Tropical Diseases*, volume 4 (11), Smout, M. J. et al. "A Novel High Throughput Assay for Anthelmintic Drug Screening and Resistance Diagnosis by Real-Time Monitoring of Parasite Motility" 2010.

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1. Smout, M. J. *et al.* A Novel High Throughput Assay for Anthelmintic Drug Screening and Resistance Diagnosis by Real-Time Monitoring of Parasite Motility. *PLoS Negl. Trop. Dis.* **2010** Nov 16, 4(11), e885
2. Rinaldi, G. *et al.* Viability of Developmental Stages of *Schistosoma mansoni* Quantified with xCELLigence Worm Real-Time Motility Assay (xWORM). *Int. J. Parasitol. Drugs Drug Resist.* **2015** Aug 6, 5(3), 141–8
3. Zeraik, A. E. *et al.* Reversible Paralysis of *Schistosoma mansoni* by Forchlorfenuron, a Phenylurea Cytokinin that Affects Septins. *Int. J. Parasitol.* **2014** Jul, 44(8), 523–31
4. Hong, Y. *et al.* Transcriptional Responses of *In Vivo* Praziquantel Exposure in Schistosomes Identifies a Functional Role for Calcium Signaling Pathway Member CamKII. *PLoS Pathog.* **2013** Mar, 9(3), e1003254
5. Silbereisen, A. *et al.* Exploration of Novel *In Vitro* Assays to Study Drugs Against *Trichuris* spp. *J. Microbiol. Methods.* **2011** Nov, 87(2), 169–75
6. Tritten, L. *et al.* Comparison of Novel and Existing Tools for Studying Drug Sensitivity Against the Hookworm *Ancylostoma ceylanicum* *In Vitro*. *Parasitology* **2012** Mar, 139(3), 348–57

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© 安捷伦科技（中国）有限公司，2019