


Microbiology catalogue



in partnership with



Excellence through measurement

Contents

- 3** About ATCC
- 4** About the ATCC / LGC Standards partnership
- 5** Quality control strains
- 7** Pharmaceutical applications
- 12** Food microbiology
- 20** Food microbiology, Genomic DNA
- 21** Water microbiology
- 24** Catalogue listing
- 38** Nucleic Acids prepared from ATCC Genuine Cultures®
- 39** Fully sequenced microbes
- 45** Microbial technical resource
- 61** Instructions for rehydrating freeze-dried cultures
- 63** ATCC MTA

All care has been taken in the compilation of the information contained in this catalogue. However, any applications for products are suggestions only, and LGC Standards makes no express or implied representations or warranties regarding the accuracy, content, completeness, or reliability of the information or any suggestions provided, and specifically disclaims any and all implied warranties with respect to the same, including without limitation any warranties of merchantability, fitness or suitability for particular purpose.

About ATCC

ATCC is an independent, private, nonprofit biological resource centre (BRC) and research organisation.

As a biological resource centre, ATCC authenticates microorganisms and cell lines and manages logistics of long-term preservation and distribution of cultures for the scientific community. ATCC supports the cultures it acquires and authenticates with expert technical support, intellectual property management and characterisation data.

As a research organisation, ATCC works to generate new knowledge and technology, as well as to continuously improve its function as a BRC. ATCC scientists develop new *in vitro* model systems, describe new species, search for new disease biomarkers and build bodies of characterisation data for valuable biological materials. In addition, they study improved methods for characterisation, long-term preservation and optimal authentication of biological materials.

Mission

ATCC is a global nonprofit bioresource centre and research organisation that provides biological products, technical services and educational programmes to private industry, government and academic organisations. The mission of ATCC is to acquire, authenticate, preserve, develop and distribute biological materials, information, technology, intellectual property and standards for the advancement and application of scientific knowledge.

Vision

The ATCC vision is to use its resources and experience as a BRC to become the world leader in standard biological reference materials management, intellectual property resource management and translational research as applied to biomaterial development, standardisation and certification.

Status

ATCC is a nonprofit 501(c)(3) organisation. The culture fees paid by purchasers support the functions of the mission.

About the ATCC / LGC Standards partnership

LGC Standards' partnership with ATCC facilitates the distribution of ATCC cultures and bioproducts to life science researchers throughout Europe and India. Our specific aims are to make access to the important resources of ATCC more easily accessible to the European scientific community through a local stock holding of more than 5,000 individual culture items supported by our local office network delivering the highest levels of customer service and technical support.

For materials which are part of the European stock holding, delivery times are typically 3-5 working days from order receipt, and we are able to assist with permit requirements as necessary for many of the "controlled" materials in the ATCC collection. For ATCC products which are not able to be included in the European stock, delivery times are typically 2-3 weeks.

ATCC Genuine Cultures®: Direct from the source

The ATCC Genuine Cultures emblem was designed to help scientists identify authentic ATCC microbial strains that come directly from ATCC. ATCC Genuine Cultures are backed by meticulous laboratory procedures and 80 years of experience and represent ATCC's high standards of quality for microbial strains:

- Full characterisation of each strain to establish identity
- Utilisation of a seed stock system to minimise subculturing
- Careful preservation and storage protocols to maintain the culture safely and effectively

As always, a customer who purchases a genuine ATCC culture from LGC Standards is receiving a culture that is a direct, minimal-passage descendant of the original material deposited with ATCC and has been cultured only by ATCC. These products are backed by our warranty and covered by our expert technical support.

Don't take chances on the quality of your cultures. Insist on products that meet ATCC's high standards of full characterisation and low passage number. ATCC Genuine Cultures are available in Europe only through LGC Standards.




Quality control strains

Commercial firms specify ATCC® strains as controls for rapid identification, minimum inhibitory concentration of antibiotics and antibiotic susceptibility panels. Species and strains specified or recommended for use are listed in this brochure. Repeated subculturing can be detrimental to quality control strains.

Every passage carries the potential for contamination, genetic drift and mutation, taking your culture farther and farther from the original. Most professional organisations that draft standards require that cultures not be used beyond five passages from the original ATCC strain. ATCC defines the first passage to be the broth or agar culture prepared from the vial supplied by ATCC. Each subsequent transfer to a new broth or agar is counted as an additional passage.


Come to the source for ATCC Genuine Cultures® for your laboratory quality control protocols.

ATCC Genuine Cultures® are available throughout Europe exclusively from LGC Standards.



Microbial applications directory

- Pharmaceutical applications
- Food microbiology
- Water microbiology



Pharmaceutical applications

Pharmaceutical companies are well versed in the important role of microbiological testing in its basic functions – product research and development, process validation, manufacturing and quality control.

ATCC cultures are cited in many international standards including most of the major pharmacopoeia, and in many laboratories. ATCC Genuine Cultures® have always been an integral part of pharmaceutical products performance testing, both as positive controls, negative controls and as identification standards.

QC testing guidelines recommend that a strain may not be used in testing which is more than 5 passages from the original reference strain – defined as the original strain and provided by a recognised culture collection such as ATCC. Organisms supplied for QC testing by ATCC are considered passage zero. LGC Standards is the

only source for ATCC cultures in Europe. A key objective of the ATCC partnership is to make these important ATCC materials available to the pharmaceutical quality control community, efficiently and cost effectively, and to support laboratory quality programmes across the pharmaceutical industry. All ATCC cultures sold by LGC Standards are the original strain provided by ATCC. They are never regrown or repackaged and are fully supported by the ATCC comprehensive quality programme which utilises a polyphasic approach whenever possible to authenticate and identify all ATCC Genuine Cultures.

The following lists provide an easy to use reference for all ATCC microbial strains quoted in the USP and EP, including easy to use cross-referenced tables that list all of the microbial strains specified for each of the pharmacopoeial QC tests. If you are not able to find what you are looking for, please contact us.

With a full range of ATCC cultures, LGC Standards is uniquely positioned to meet all your biological standards needs.



Pharmaceutical application

	ATCC® No.
EP 6.5 2.6.1 Sterility	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Candida albicans</i>	10231™
<i>Clostridium sporogenes</i>	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
EP 6.5 2.6.12 Microbial examination of non-sterile products: Microbial enumeration tests	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Candida albicans</i>	10231™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
EP 6.5 2.6.13 Microbial examination of non-sterile products: Test for specified microorganisms	
<i>Candida albicans</i>	10231™
<i>Clostridium sporogenes</i>	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Escherichia coli</i>	8739™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>	14028™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
EP 6.5 2.6.27 Microbiological control of cellular products	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Bacteroides fragilis</i>	25285™
<i>Candida albicans</i>	10231™
<i>Clostridium sporogenes</i>	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Propionibacterium acnes</i>	11827™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
<i>Streptococcus pyogenes</i>	19615™
<i>Yersinia enterocolitica</i>	9610™
EP 6.5 2.6.7 Mycoplasmas	
<i>Acholeplasma laidlawii</i>	23206™
<i>Mycoplasma fermentans</i>	19989™
<i>Mycoplasma gallisepticum</i>	19610™
<i>Mycoplasma hyorhinis</i>	17981™
<i>Mycoplasma hyorhinis</i>	29052™
<i>Mycoplasma orale</i>	23714™
<i>Mycoplasma pneumoniae</i>	15531™
<i>Mycoplasma synoviae</i>	25204™

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Pharmaceutical application

ATCC® No.

EP 6.5 2.7.2 Microbiological assays of antibiotics

<i>Bacillus pumilus</i>	14884™
<i>Bacillus subtilis</i>	11774™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Bordetella bronchiseptica</i>	4617™
<i>Enterococcus hirae</i>	10541™
<i>Escherichia coli</i>	10536™
<i>Escherichia coli</i>	9637™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	10031™
<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	9341™
<i>Micrococcus luteus</i>	10240™
<i>Mycobacterium smegmatis</i>	607™
<i>Saccharomyces cerevisiae</i>	9763™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	9144™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538P™
<i>Staphylococcus epidermidis</i>	12228™

EP 6.5 5.1.1 Methods of preparation of sterile products

<i>Brevundimonas diminuta</i> (formerly <i>Pseudomonas</i> sp.)	19146™
---	--------

EP 6.5 5.1.2 Biological indicators of sterilisation

<i>Bacillus atrophaeus</i> (formerly <i>Bacillus subtilis</i> var. <i>niger</i>)	9372™
<i>Bacillus pumilus</i> (strain recommended by ATCC)	27142™
<i>Geobacillus stearothermophilus</i> (formerly <i>Bacillus stearothermophilus</i>)	7953™

EP 6.5 5.1.3 Efficacy of antimicrobial preservation

<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Candida albicans</i>	10231™
<i>Escherichia coli</i>	8739™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™

EP 6.5 Monographs – Penicillamine

<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	9341™
---	-------

EP 6.5 Monographs – Water for injections

<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Pseudomonas aeruginosa</i>	9027™

EP 6.5 Monographs – Water, highly purified

<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Pseudomonas aeruginosa</i>	9027™

EP 6.5 Monographs – Water, purified

<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Pseudomonas aeruginosa</i>	9027™

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Pharmaceutical application

	ATCC® No.
USP32, Antimicrobial effectiveness testing <51>	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Candida albicans</i>	10231™
<i>Escherichia coli</i>	8739™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Microbiological examination of nonsterile products: Microbial enumeration tests <61>	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Candida albicans</i>	10231™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Microbial examination of nonsterile products: Tests for specified microorganisms <62>	
<i>Candida albicans</i>	10231™
<i>Clostridium sporogenes</i>	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Escherichia coli</i>	8739™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>	14028™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Sterility tests <71>	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Bacteroides vulgatus</i>	8482™
<i>Candida albicans</i>	10231™
<i>Clostridium sporogenes</i>	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	9341™
<i>Pseudomonas aeruginosa</i>	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Antibiotics – Microbial assays <81>	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Bordetella bronchiseptica</i>	4617™
<i>Enterococcus hirae</i>	10541™
<i>Escherichia coli</i>	10536™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	10031™
<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	9341™
<i>Micrococcus luteus</i>	10240™
<i>Mycobacterium smegmatis</i>	607™
<i>Pseudomonas aeruginosa</i>	25619™
<i>Saccharomyces cerevisiae</i>	2601™
<i>Saccharomyces cerevisiae</i>	9763™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29737™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	9144™
<i>Staphylococcus epidermidis</i>	12228™
USP32, Biological indicators, <1035> (Official Monographs)	
<i>Bacillus atrophaeus</i> (formerly <i>Bacillus subtilis</i> var. <i>niger</i>)	9372™
<i>Geobacillus stearothermophilus</i> (formerly <i>Bacillus stearothermophilus</i>)	12980™
<i>Geobacillus stearothermophilus</i> (formerly <i>Bacillus stearothermophilus</i>)	7953™

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Pharmaceutical application

	ATCC® No.
USP32 Disinfectants and antiseptics <1072> (according to AOAC)	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i>	19659™
<i>Candida albicans</i>	10231™
<i>Candida albicans</i>	2091™
<i>Escherichia coli</i>	11229™
<i>Penicillium chrysogenum</i>	11709™
<i>Pseudomonas aeruginosa</i>	15442™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Sterilization and sterility assurance of compendial articles <1211>	
<i>Bacillus pumilus</i> (strain recommended by ATCC)	27142™
USP32 Sterilization and sterility assurance of compendial articles <1211>	
<i>Bacillus atrophaeus</i> (formerly <i>Bacillus subtilis</i> var. <i>niger</i>) (strain recommended by ATCC)	9372™
<i>Geobacillus stearothermophilus</i> (formerly <i>Bacillus stearothermophilus</i>) (strain recommended by ATCC)	7953™
USP32 Sterilization and sterility assurance of compendial articles <1211>	
<i>Brevundimonas diminuta</i> (formerly <i>Pseudomonas</i> sp.)	19146™
<i>Serratia marcescens</i>	14756™
USP32 Microbial limit tests – Dietary supplements – Preparatory testing, growth promotion testing <2021>	
<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	16404™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	6633™
<i>Candida albicans</i>	10231™
<i>Escherichia coli</i>	8739™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>	13311™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	6538™
USP32, Oil and water-soluble vitamins with mineral tablets (biotin, cyanocobalamin) (Official Monographs)	
<i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> (formerly <i>Lactobacillus leichmannii</i>)	7830™
<i>Lactobacillus plantarum</i>	8014™
USP32, Penicillamine (Official Monographs)	
<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	9341™
USP32, Penicillin G procaine and novobiocin sodium intramammary infusion (Official Monographs)	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	12692™
USP32, Tetracycline hydrochloride and novobiocin sodium tablets (Official Monographs)	
<i>Escherichia coli</i>	10536™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29737™
USP32, Water-soluble vitamins capsules (dexpanthenol, panthenol) (Official Monographs)	
<i>Pediococcus acidilactici</i>	8042™

European Pharmacopoeia 6th Edition Supplement 6.5
U.S. Pharmacopoeia 32nd rev.; 2009

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Food microbiology

Top performing food processors realise the importance of effective microbiological testing. Safety, reputation and business performance depend on it. Microbial strains with confirmed identity, viability and purity, backed by meticulous laboratory procedures that minimise sub-culturing are important to product safety, and they are important to ATCC. Food processors look to ATCC to provide the top quality microbial strains needed to maintain outstanding quality control (QC) microbiology programs. Whether in quality control testing, process validation or research and development, microbial strains should be considered vital factors in generating valid, accurate results. To identify microbial strains that are handled, stored and tested with the unparalleled expertise that comes from over 80 years of experience, look for the ATCC symbols of quality.

Some standards organisations have recommended, and ATCC concurs, that strains being used in standard test methods should not be used beyond 5 passages from the original ATCC strain. ATCC defines the first passage as the first broth or agar culture started from the culture vial supplied by ATCC. Each subsequent transfer to new broth or agar is counted as an additional passage. Since 1925, ATCC has set the standard for authenticating and distributing biological materials for research and testing in food and other industries. The following lists summarise the cultures specified or recommended for use in standards (including specifications, assays, and other tests) written by government agencies and professional organisations. Most standards name specific strains, often by ATCC number. Species are listed here by the name currently used at ATCC. If the name used in the standard is different, it has most likely been reclassified based on recent molecular phylogeny and genotyping methods following the Rules of the Bacteriological Code and current taxonomic interpretation.



Food microbiology standards

International Organisation for Standardization (ISO)

ATCC® No.

ISO 6888-3:2003 – Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of coagulase – positive staphylococci (*Staphylococcus aureus* and other species) – Part 3: Detection and MPN technique for low numbers.

<i>Escherichia coli</i> FDA strain Seattle 1946	25922™
<i>Escherichia coli</i> Crooks	8739™
<i>Penicillium aurantiogriseum</i>	8732™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Seattle 1945	25923™

ISO 11290-2:1998 – Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of *Listeria monocytogenes* – Part 2: Enumeration method.

<i>Rhodococcus equi</i>	6939™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Seattle 1945	25923™

ISO 11133-2:2003 – Microbiology of food and animal feeding stuffs – Guidelines on preparation and production of culture media – Part 2: Practical guidelines on performance testing of culture media.

<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Bacillus cereus</i> FDA strain PCI 213	11778™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Candida albicans</i> 3147	10231™
<i>Citrobacter freundii</i> LRA 117.03.76	43864™
<i>Clostridium perfringens</i>	13124™
<i>Enterococcus faecalis</i>	19433™
<i>Enterococcus faecalis</i> Portland	29212™
<i>Escherichia coli</i> FDA strain Seattle 1946	25922™
<i>Escherichia coli</i> Crooks	8739™
<i>Escherichia coli</i>	11775™
<i>Lactobacillus sakei</i> subsp. <i>sakei</i> T.S.	15521™
<i>Lactococcus lactis</i> subsp. <i>lactis</i>	19435™
<i>Listeria monocytogenes</i> Li 20	19111™
<i>Listeria monocytogenes</i> 1071/53	13932™
<i>Pediococcus damnosus</i> NCDO 1832	29358™
<i>Penicillium aurantiogriseum</i> IMI 19759	16025™
<i>Proteus mirabilis</i> CDC PR 14	29906™
<i>Pseudomonas aeruginosa</i> Boston 41501	27853™
<i>Saccharomyces cerevisiae</i>	9763™
<i>Salmonella enterica</i> subsp. <i>enterica</i> CDC K-1891	13076™
<i>Salmonella enterica</i> subsp. <i>enterica</i> CDC 6516-60	14028™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Seattle 1945	25923™
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
<i>Yersinia enterocolitica</i> 33114	9610™
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> Billups-1803-68	23715™

ISO 21527-1 – Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of yeasts and moulds – Part 1: Colony count technique in products with water activity greater than 0,95.

<i>Saccharomyces cerevisiae</i>	9763™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Mucor racemosus</i>	42647™
<i>Escherichia coli</i> FDA Strain Seattle 1946	25922™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™

Food microbiology standards

International Organisation for Standardization (ISO)

ATCC® No.

ISO 21527-2 – Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of yeasts and moulds – Part 2: Colony count technique in products with water activity less than or equal to 0,95.

Saccharomyces cerevisiae

9763™

Wallemia sebi

42694™

Aspergillus restrictus

42693™

Eurotium rubrum

42690™

Escherichia coli FDA Strain Seattle 1946

25922™

Bacillus subtilis subsp. *spizizenii* NRS 231

6633™

British Standards Institution (BSI)

ATCC® No.

BS EN 1104:2005 – Paper and board intended to come into contact with foodstuffs – Determination of the transfer of antimicrobial constituents.

Aspergillus niger 4247

6275™

BS EN 13697:2001 – Chemical disinfectants and antiseptics – Quantitative non-porous surface test for the evaluation of bactericidal and/or fungicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas – Test method and requirements without mechanical action (phase 2/step 2).

Aspergillus brasiliensis WLRI 034(120)

16404™

Candida albicans 3147

10231™

Enterococcus hirae FDA M19

10541™

Escherichia coli MacLeod

10536™

Pseudomonas aeruginosa PRD-10

15442™

Saccharomyces cerevisiae

9763™

Salmonella enterica subsp. *enterica*

13311™

Staphylococcus aureus subsp. *aureus* FDA 209

6538™

BS EN 13704:2002 – Chemical disinfectants – Quantitative suspension test for the evaluation of sporicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas – Test method and requirements (phase 2, step 1).

Bacillus cereus Type Strain A, variant IV

12826™

Bacillus subtilis subsp. *spizizenii* NRS 231

6633™

BS EN 14131:2003 – Foodstuffs – Determination of folate by microbiological assay.

Lactobacillus rhamnosus

7469™

BS EN 1650:1998 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal activity of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas – Test method and requirements (phase 2, step 1).

Aspergillus brasiliensis WLRI 034(120)

16404™

Candida albicans 3147

10231™

Saccharomyces cerevisiae

9763™

BS EN ISO 11290-1:1997 – Microbiology of food and animal feeding stuffs - Horizontal method for the detection and enumeration of *Listeria monocytogenes* - Part 1: Detection method, Annex B.

Enterococcus faecalis Portland

29212™

Escherichia coli FDA strain Seattle 1946

25922™

Listeria innocua SLCC 3379

33090™

Listeria monocytogenes 1071/53

13932™

Listeria monocytogenes Li20

19111™

Rhodococcus equi

6939™

Staphylococcus aureus subsp. *aureus* Seattle 1945

25923™

Food microbiology standards

British Standards Institution (BSI)	ATCC® No.
BS EN ISO 11290-2:1998 – Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of <i>Listeria monocytogenes</i> – Part 2: Enumeration method, Annex B.	
<i>Rhodococcus equi</i>	6939™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Seattle 1945	25923™
BS EN ISO 21871:2006 – Microbiology of food and animal feeding stuffs.	
<i>Bacillus cereus</i> FDA strain PCI 213	11778™
<i>Escherichia coli</i> FDA strain Seattle 1946	25922™
<i>Escherichia coli</i> Crooks	8739™
BS EN ISO 6888-3:2003 – Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coagulase-positive staphylococci (<i>Staphylococcus aureus</i> and other species) – Part 3: Detection and MPN technique for low numbers.	
<i>Penicillium aurantiogriseum</i>	8732™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
DD ENV 14166:2001 – Foodstuffs – Determination of vitamin B6 by microbiological assay.	
<i>Saccharomyces cerevisiae</i> 4228	9080™
Association of Official Analytical Chemists (AOAC)	ATCC® No.
AOAC 955.12 – Testing disinfectants against <i>Staphylococcus aureus</i>, phenol coefficient method.	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
AOAC 955.13 – Testing disinfectants against <i>Pseudomonas aeruginosa</i>, phenol coefficient method.	
<i>Pseudomonas aeruginosa</i> PRD-10	15442™
AOAC 955.14 – Testing disinfectants against <i>Salmonella choleraesuis</i>, use-dilution methods.	
<i>Salmonella enterica</i> subsp. <i>enterica</i> ETS 34	10708™
AOAC 955.15 – Testing disinfectants against <i>Staphylococcus aureus</i>, use-dilution methods.	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
AOAC 955.16 – Chlorine (available) in disinfectants, germicidal equivalent concentration.	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
AOAC 955.17 – Fungicidal activity of disinfectants.	
<i>Trichophyton mentagrophytes</i> 640	9533™
AOAC 957.23 – Antibiotics in feeds, microbiological methods.	
<i>Bacillus cereus</i> FDA strain PCI 213	11778™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Escherichia coli</i> UC 527	29998™
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
<i>Micrococcus luteus</i> Mercedita	7468™
<i>Micrococcus luteus</i> 130.21	10240™
<i>Saccharomyces cerevisiae</i>	9763™
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
AOAC 960.09 – Germicidal and detergent sanitising action of disinfectants.	
<i>Escherichia coli</i> AMC 198	11229™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™

Food microbiology standards

Association of Official Analytical Chemists (AOAC)	ATCC® No.
AOAC 960.46 – Vitamin assays, microbiological method. <i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> 313 <i>Lactobacillus rhamnosus</i>	7830™ 7469™
AOAC 960.47 – Amino acids in vitamin preparations. <i>Enterococcus hirae</i> R <i>Lactobacillus plantarum</i> 17-5 <i>Pediococcus acidilactici</i>	9790™ 8014™ 8042™
AOAC 960.67 – Hygromycin B in feeds, microbiological method. <i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
AOAC 961.02 – Germicidal spray products as disinfectants. <i>Pseudomonas aeruginosa</i> PRD-10 <i>Salmonella enterica</i> subsp. <i>enterica</i> ETS 34 <i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209 <i>Trichophyton mentagrophytes</i> 640	15442™ 10708™ 6538™ 9533™
AOAC 961.15 – Vitamin B6 (pyridoxine, pyridoxal, pyridoxamine) in food extracts, microbiological method. <i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
AOAC 962.14 – Beta-lactam antibiotics in milk, qualitative field disk assay. <i>Saccharomyces cerevisiae</i> 4228	9080™
AOAC 964.02 – Testing disinfectants against <i>Pseudomonas aeruginosa</i>, use-dilution method. <i>Pseudomonas aeruginosa</i> PRD-10	15442™
AOAC 972.56 – Monensin in feeds, microbiological method. <i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
AOAC 975.56 – Virus in beef (ground), microbiological method. <i>Cercopithecus aethiops</i> Vero	CCL-81™
AOAC 976.37 – Monensin in feeds, turbidimetric method. <i>Enterococcus hirae</i> R	8043™
AOAC 977.37 – Chlortetracycline HCl in feeds, turbidimetric method. <i>Staphylococcus aureus</i> subsp. <i>aureus</i> 3R7089 strain Oxford	9144™
AOAC 979.14 – Beta-lactam antibiotics, qualitative disc method I. <i>Geobacillus stearothermophilus</i> NRS T15	10149™
AOAC 982.16 – Beta-lactam antibiotics in milk, quantitative disc method. <i>Geobacillus stearothermophilus</i> NRS T15	10149™
AOAC 982.17 – Beta-lactam antibiotics in milk, qualitative disc method II. <i>Geobacillus stearothermophilus</i> NRS T15	10149™
AOAC 982.36 – Invasiveness by <i>Escherichia coli</i> of mammalian cells, microbiological method. <i>Homo sapiens</i> (human) HeLa	CCL-2™
AOAC 982.43 – Bacitracin in premix feeds. <i>Micrococcus luteus</i> 130.21	10240™

Food microbiology standards

Association of Official Analytical Chemists (AOAC)	ATCC® No.
AOAC 984.34 – Detection of <i>Escherichia coli</i> producing heat labile enterotoxin, DNA colony hybridisation method. <i>Escherichia coli</i> H10407 <i>Escherichia coli</i> pBR313	35401™ 37018™
AOAC 985.32 – Vitamin B6 in ready-to-feed milk based infant formula, microbiological method. <i>Saccharomyces cerevisiae</i> 4228	9080™
AOAC 986.23 – Vitamin B12 activity in milk based infant formula, turbidimetric method. <i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> 313 <i>Weissella confusa</i> 548-D	7830™ 10881™
AOAC 991.38 – Salmonella in foods. <i>Escherichia coli</i> FDA strain Seattle 1946	25922™
AOAC 991.47 – Testing disinfectants against <i>Salmonella choleraesuis</i>, hard surface carrier test method. <i>Salmonella enterica</i> subsp. <i>enterica</i> ETS 34	10708™
AOAC 991.48 – Testing disinfectants against <i>Staphylococcus aureus</i>, hard surface carrier test method. <i>Saccharomyces cerevisiae</i> <i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	9763™ 6538™
AOAC 991.49 – Testing disinfectants against <i>Pseudomonas aeruginosa</i>, hard surface carrier test method. <i>Pseudomonas aeruginosa</i> PRD-10	15442™
AOAC 992.05 – Folic acid (pteroylglutamic acid) in infant formula, microbiological methods. <i>Lactobacillus rhamnosus</i>	7469™
AOAC 992.18 – Listeria species – Biochemical identification method (MICRO-ID) Listeria. <i>Lactococcus lactis</i> subsp. <i>cremoris</i> NCDO 607 <i>Listeria grayi</i> V-1 <i>Listeria monocytogenes</i> Li 20 <i>Listeria seeligeri</i> CIP 100100 <i>Streptococcus mitis</i>	19257™ 25400™ 19111™ 35967™ 6249™
AOAC 992.19 – Listeria species – Biochemical identification method (Vitek GPI and GNI+). <i>Acinetobacter baumannii</i> 2208 <i>Bordetella bronchiseptica</i> 3127 <i>Enterococcus durans</i> 23C2 <i>Enterococcus faecalis</i> Portland <i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> <i>Proteus mirabilis</i> <i>Pseudomonas aeruginosa</i> Boston 41501 <i>Serratia odorifera</i> 1073 <i>Shigella sonnei</i> <i>Staphylococcus xylosum</i> KL 162 <i>Streptococcus equi</i> subsp. <i>equi</i> 2-1-23 <i>Streptococcus gallolyticus</i> 38 <i>Streptococcus pneumoniae</i> R36a rough phase <i>Streptococcus pyogenes</i> Bruno	19606™ 10580™ 6056™ 29212™ 13883™ 7002™ 27853™ 33077™ 25931™ 29971™ 9528™ 9809™ 27336™ 19615™

Food microbiology standards

Association of Official Analytical Chemists (AOAC)	ATCC® No.
AOAC 993.29 – Bacitracin-MD (BMD) in complete feed, microbiological plate assay method. <i>Micrococcus luteus</i> 130.21	10240™
AOAC 997.17 – Microbial ranking of porous packaging materials (Exposure Chamber Method). <i>Bacillus atrophaeus</i> NRS 1221A	9372™
AOAC 998.02 – Neomycin in feeds - stahl microbiological agar diffusion assay. <i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
AOAC 2004.04 – Identification of presumptive isolates of <i>Bacillus anthracis</i>. <i>Bacillus cereus</i>	14579™
AOAC 2004.05 – Total folates in cereals and cereal Foods. <i>Lactobacillus rhamnosus</i>	7469™
AOAC 2004.11 – Identification of <i>Bacillus anthracis</i> from culture - Gas chromatographic analysis of Fatty Acid Methyl Esters (FAMES). <i>Bacillus cereus</i>	14679™
<hr/>	
Bacteriological Analytical Manual, U.S. Food and Drug Administration (BAM)	ATCC® No.
BAM 10.F – Detection and enumeration of <i>Listeria monocytogenes</i> in foods, the CAMP Test. <i>Rhodococcus equi</i> <i>Staphylococcus aureus</i> subsp. <i>aureus</i> Seattle 1945 <i>Staphylococcus pseudintermedius</i>	6939™ 25923™ 49444™
BAM 13b – Electrophoretic and immunoblot analysis of Staphylococcal Enterotoxins in food. <i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 196E	13565™
BAM 20A – Inhibitory substances in milk. <i>Geobacillus stearothermophilus</i> NRS T15 <i>Kocuria rhizophila</i> FDA strain PCI 1001	10149™ 9341™
BAM 24 – Identification of foodborne bacterial pathogens by gene probes: Enterotoxigenic <i>Escherichia coli</i>: heat stable Enterotoxin (Human), heat stable Enterotoxin (Porcine), and heat labile Enterotoxin. <i>Escherichia coli</i> FDA strain Seattle 1946	25922™
BAM 24 – Identification of foodborne bacterial pathogens by gene probes, <i>Listeria monocytogenes</i>: combination of invasion associated protein (iap) and Hemolysin (hly) Gene Probes – AD713. <i>Listeria innocua</i> SLCC 3379	33090™
BAM 24 – Identification of foodborne bacterial pathogens by gene probes, <i>Vibrio vulnificus</i> VV6. <i>Vibrio vulnificus</i> 324	27562™
BAM 24 – Identification of foodborne bacterial pathogens by gene probes, <i>Vibrio parahaemolyticus</i> tdh3. <i>Vibrio parahaemolyticus</i> EB 101	17802™
BAM 4.II.3 – Enumeration of <i>Escherichia coli</i> and the Coliform bacteria: LST-MUG method for detecting <i>Escherichia coli</i> in chilled or frozen foods exclusive of bivalve molluscan shellfish. <i>Enterobacter aerogenes</i> NCDC 819-56 <i>Escherichia coli</i> FDA strain Seattle 1946	13048™ 25922™
BAM 5.D.7 – Salmonella: Isolation of Salmonella <i>Salmonella enterica</i> subsp. <i>diarizonae</i> 62	29934™

Food microbiology standards

Other important food microbiology organisms

ATCC® No.

<i>Salmonella enterica</i> subsp. <i>enterica</i>	BAA-215™
<i>Salmonella enterica</i> subsp. <i>diarizonae</i>	BAA-216™
<i>Cronobacter muytjensii</i>	51329™
<i>Escherichia coli</i>	11303™
<i>Campylobacter coli</i>	33559™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	7080™
<i>Campylobacter lari</i>	35221™
<i>Campylobacter coli</i>	43478™

AOAC International = Association of Official Analytical Chemists

ISO = International Organization for Standardization

DD = Draft for Development

BSI = British Standards Institution

BAM = Bacteriological Analytical Manual, U.S. Food and Drug Administration

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Food microbiology, Genomic DNA

ATCC® No.	Strain	Product related to standard/application
16404D-2	<i>Aspergillus brasiliensis</i> (formerly <i>Aspergillus niger</i>)	ISO 11133-2:2003; ISO 21527-1
6633D-5	<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	ISO 11133-2:2003; ISO 21527-1
33560D-5	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	(ISO 10272)
10231D-5	<i>Candida albicans</i> , serotype A	ISO 11133-2:2003; ISO 21527-1
13124D-5	<i>Clostridium perfringens</i>	ISO 11133-2:2003
8739D-5	<i>Escherichia coli</i>	ISO 11133-2:2003, ISO 6888-3:2003, ISO 11290-1/2
25922D-5	<i>Escherichia coli</i>	ISO 11133-2:2003, ISO 6888-3:2003, ISO 11290-1/2; ISO 21527-1
9341D-5	<i>Kocuria rhizophila</i> (formerly <i>Micrococcus luteus</i>)	Antibacterial substances in animal feeding stuffs
19435D-5	<i>Lactococcus lactis</i> subsp. <i>lactis</i>	ISO 11133-2:2003
27853D-5	<i>Pseudomonas aeruginosa</i>	ISO 11133-2:2003
9763D-5	<i>Saccharomyces cerevisiae</i>	ISO 11133-2:2003
25923D-5	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	ISO 6888-3:2003, CAMP Listeria (FDA)
6538D-5	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	ISO 11133-2:2003
12228D-5	<i>Staphylococcus epidermidis</i>	ISO 11133-2:2003
23715D-5	<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i>	ISO 11133-2:2003
9610D	<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i>	ISO 11133-2:2003

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Water microbiology


Microbial analyses are a vital part of many companies' environmental laboratory testing programmes. The control strains used to support and measure the results of these analyses are an important laboratory tool. The analysis of water and environmental samples can have a number of different objectives including, regulatory compliance, resource allocation, safety as well as protection of the environment and your business.

ATCC cultures are specified by standard-setting organisations and regulatory bodies for use in their official assays and tests. Manufacturers of microbial identification systems also specify ATCC microorganisms for the quality testing of instrumentation, kits or reagents where the use of authenticated control cultures is an essential element to ensuring the quality and comparability of laboratory results. One of the primary objectives of the ATCC/LGC partnership is to help make the important resources of the ATCC collection

available to European quality control analysts across a range of industries including environmental and water testing. The following lists describe a number of the key tests used in environmental and water sample analysis laboratories. All ATCC cultures sold by LGC standards are the original strain provided by ATCC. They are never regrown or repackaged and are fully supported by the ATCC comprehensive quality programme which utilises a polyphasic approach whenever possible to authenticate and identify all ATCC Genuine Cultures®.

The following lists are cross-referenced to provide an alphabetical listing of relevant microbial strains and a listing by test method with the recommended control strains from the ATCC. The catalogue tables provide an easy to use guide to accessing the right materials.

Used in conjunction with a comprehensive Proficiency Testing programme; the use of authenticated, highly characterised microbial controls such as those available from ATCC provides the most effective tools for the ongoing monitoring and improvement of analytical results.



Water microbiology

ATCC® No.	Strain	Standard/Application
9144™	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	AS/NZS 4276.16:1999, AS/NZS 4276.20:2003
10145™	<i>Pseudomonas aeruginosa</i>	IDEXX (Colilert, Colisure), AS 4276.4-1995, AS 4276.5-1995, AS 4276.10-1995, AS 4276.11-1995, AS 4276.12-1995, AS 4276.13-1995, AS/NZS 4276.16:1999, AS/NZS 4276.18:2001, EN ISO 16266:2008
10400™	<i>Aerococcus viridans</i>	IDEXX Enterolert-E
11775™	<i>Escherichia coli</i>	EPA 1103.1 (2005), EPA 1106.1 (2005), EPA 1600 (2006), EPA 1603 (2005), EPA 1680 (2005), AS 4276.6-1995, AS 4276.7-1995, AS 4276.12-1995, AS 4276.13-1995, AS/NZS 4276.16:1999, EN ISO 16266:2008
11778™	<i>Bacillus cereus</i>	EPA 1605 (2005)
13048™	<i>Enterobacter aerogenes</i>	EPA 1103.1 (2005), EPA 1603 (2005), EPA 1680 (2005), EPA 1681 (2005), AS 4276.4-1995, AS 4276.5-1995, AS 4276.6-1995, AS 4276.7-1995, AS 4276.10-1995, AS 4276.11-1995, AS/NZS 4276.18:2001
13124™	<i>Clostridium perfringens</i>	AS/NZS 4276.17.1:2000, AS/NZS 4276.17.2:2000, AS/NZS 4276.18:2001
14028™	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>	EPA 1682 (2005)
15597™	<i>Escherichia coli</i>	ISO 10705-1:1995, EPA 1601 (2001), EPA 1601 (2001)
15597-B1™	<i>Escherichia coli</i> , MS2	ISO 10705-1:1995, EPA 1601 (2001), EPA 1601 (2001)
19433™	<i>Enterococcus faecalis</i>	EPA 1103.1 (2005), EPA 1106.1 (2005), EPA 1600 (2006), EPA 1603 (2005), AS 4276.8-1995, AS 4276.9-1995
23631™	<i>Escherichia coli</i>	ISO 10705-1:1995
25922™	<i>Escherichia coli</i>	IDEXX (Colilert, Colisure), EPA 1605 (2005), EPA 1681 (2005), EPA 1682 (2005)
25923™	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	AS 4276.8-1995, AS 4276.9-1995
27853™	<i>Pseudomonas aeruginosa</i>	EPA 1605 (2005), EPA 1680 (2005)

Water microbiology

ATCC® No.	Strain	Standard/Application
31488™	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	IDEXX (Colilert, Colisure)
35667™	<i>Enterococcus faecium</i>	IDEXX Enterolert-E
43862™	<i>Serratia marcescens</i>	IDEXX Enterolert-E
700078™	<i>Escherichia coli</i>	ISO 10705-2:2000
700609™	<i>Escherichia coli</i>	EPA 1601 (2001), EPA 1602 (2001)
700786™	<i>Bacteroides fragilis</i>	ISO 10705-4:2001
700786-B1™	<i>Bacteroides fragilis</i> , B56-3	ISO 10705-4:2001
700891™	<i>Escherichia coli</i>	EPA 1601 (2001), EPA 1601 (2001)
7966™	<i>Aeromonas hydrophila</i>	EPA 1605 (2005), AS/NZS 4276.18:2001
8043™	<i>Enterococcus hirae</i>	ISO 7899-1:1998

ISO = International Organization for Standardization

AS/NZS = Australian/New Zealand standards published by Standards Australia and Standards New Zealand

EPA = United States Environmental Protection Agency

Colilert, Colisure and Enterolert are trademarks or registered trademarks of IDEXX Laboratories, Inc.

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Catalogue listing

While this list of organisms represents some of the most important microbes used in industrial testing, it is only a small section of the ATCC® microbial collection.

The most complete listing can be found by visiting the LGC Standards ATCC website:
www.lgcstandards-atcc.org

Still can't find what you are looking for?
Contact your local LGC Standards office.

LGC Standards offices

France

Tel: +33 (0)3 88 04 82 82
Fax: +33 (0)3 88 04 82 90
fr@lgcstandards.com

Countries served: Belgium, France, Luxembourg, Monaco, Switzerland.

Germany

Tel: +49 (0)281 9887 230
Fax: +49 (0)281 9887 239
atcc.de@lgcstandards.com

Countries served: Austria, Croatia, Germany, Greece, Hungary, Kosovo, Moldova, Mongolia, Montenegro, Netherlands, Romania, Serbia, Slovenia, Turkey.

India

Tel: +91 80 6701 2000
Fax: +91 80 6701 2046
in@lgcpromochem.com

Italy

Tel: +39 02 24126 830
Fax: +39 02 24126 831
it@lgcstandards.com

Poland

Tel: +48 (0)22 751 31 40
Fax: +48 (0)22 751 58 45
pl@lgcstandards.com

Countries served: Armenia, Belarus, Bulgaria, Czech Republic, Lithuania, Macedonia, Poland, Russia, Slovakia, Ukraine.

Spain

Tel: +34 93 308 41 81
Fax: +34 93 307 36 12
es@lgcstandards.com

Countries served: Andorra, Portugal, Spain.

Sweden

Tel: +46 (0)33 20 90 60
Fax: +46 (0)33 20 90 79
atcc.se@lgcstandards.com

Countries served: Denmark, Estonia, Finland, Iceland, Latvia, Norway, Sweden.

United Kingdom

Tel: +44 (0)20 8943 8489
Fax: +44 (0)20 8943 8405
atcc@lgcstandards.com

Countries served: Albania, Bosnia Herzegovina, Channel Islands, Cyprus, Ireland, Liechtenstein, Malta, United Kingdom.

Catalogue

Organism	Designation	ATCC® No.
<i>Acetobacter aceti</i>	NCIB 8621	15973™
<i>Acholeplasma laidlawii</i> deposited as <i>Mycoplasma laidlawii</i>	PG8	23206™
<i>Achromobacter xylosoxidans</i> deposited as <i>Achromobacter xylosoxidans</i>	KM 543	27061™
<i>Acinetobacter baumannii</i>	Vitek #109234	BAA-747™
<i>Acinetobacter baumannii</i> deposited as <i>Bacterium anitratum</i>	2208	19606™
<i>Acinetobacter baumannii</i> deposited as <i>Moraxella glucidolytica</i> subsp. <i>nonliquefaciens</i>	5377	17978™
<i>Acinetobacter baylyi</i> deposited as <i>Acinetobacter calcoaceticus</i>	BD413	33305™
<i>Acinetobacter calcoaceticus</i> deposited as <i>Moraxella calcoacetica</i>	46	23055™
<i>Acinetobacter haemolyticus</i> deposited as <i>Herellea caseolytica</i>	57.073.192	19002™
<i>Acinetobacter lwoffii</i>	NCTC 5866	15309™
<i>Acinetobacter lwoffii</i> deposited as <i>Achromobacter lwoffii</i>	NCIB 9020	17925™
<i>Acinetobacter</i> sp. deposited as <i>Acinetobacter anitratus</i>	AmMS 203	49137™
<i>Acinetobacter</i> sp. deposited as <i>Acinetobacter anitratus</i>	AmMS 202	49139™
<i>Acinetobacter</i> sp. deposited as <i>Acinetobacter calcoaceticus</i>	AmMS 243	49466™
<i>Actinobacillus pleuropneumoniae</i> deposited as <i>Haemophilus parahaemolyticus</i>	4074	27088™
<i>Actinobacillus pleuropneumoniae</i> deposited as <i>Haemophilus parahaemolyticus</i>	S 1421	27090™
<i>Actinomyces odontolyticus</i>	CDC X363	17929™
<i>Actinomyces viscosus</i>	MG-1	43146™
<i>Actinomyces viscosus</i> deposited as <i>Odontomyces viscosus</i>	T-6	15987™
<i>Aerococcus viridans</i>	NCTC 8251	11563™
<i>Aerococcus viridans</i>	API 78-12-095	700406™
<i>Aerococcus viridans</i> deposited as <i>Gaffkya homari</i>	[ICPB 4308]	10400™
<i>Aeromonas caviae</i> deposited as <i>Aeromonas punctata</i> subsp. <i>caviae</i>	[NRRL B-968]	15468™
<i>Aeromonas hydrophila</i>	LRA 3300 776	35654™
<i>Aeromonas hydrophila</i>	AmMS 199	49140™
<i>Aeromonas hydrophila</i>	NCTC 7812	7965™
<i>Aeromonas hydrophila</i> deposited as <i>Proteus ichthyosmius</i>	[CDC 359-60, IAM 12460, NCIB 9240, NCMB 86, NCTC 8049, RH 250]	7966™
<i>Aeromonas salmonicida</i> subsp. <i>salmonicida</i>	NCMB 1102	33658™
<i>Aggregatibacter aphrophilus</i> deposited as <i>Haemophilus aphrophilus</i>	NCTC 5906	33389™
<i>Alcaligenes faecalis</i> subsp. <i>faecalis</i>	LRA 41 02 82	35655™
<i>Alcaligenes faecalis</i> subsp. <i>faecalis</i>	16, 104-1A	8750™
<i>Alicyclobacillus acidocaldarius</i> subsp. <i>acidocaldarius</i> deposited as <i>Bacillus acidocaldarius</i>	[HAMBI 2073, IFO 15652, JCM 5260, KCTC 1825, LMG 7119, NCCB 89167, NCIMB 11725, NRRL B-14509]	27009™
<i>Alicyclobacillus acidoterrestris</i> deposited as <i>Bacillus acidoterrestris</i>	GD3B	49025™
<i>Alloicoccus otitis</i>	NCFB 2890	51267™
<i>Alternaria alternata</i> deposited as <i>Alternaria tenuis</i>		6663™
<i>Anabaena variabilis</i>	[IUCC 1444, MSU A-37]	29413-U™
<i>Aneurinibacillus aneurinolyticus</i> deposited as <i>Bacillus thiaminolyticus</i>	[CIP 104323, JCM 9023, LMG 12387, Vitek #202274]	11376™
<i>Aquaspirillum</i> sp. deposited as <i>Spirillum</i> sp.	NOX	49643™
<i>Arcanobacterium pyogenes</i> deposited as <i>Actinomyces pyogenes</i>	LRA 212.10.89	49698™
<i>Arcanobacterium pyogenes</i> deposited as <i>Corynebacterium pyogenes</i>	NCTC 5224	19411™
<i>Arthrobacter psychrolactophilus</i>	B7	700733™
<i>Aspergillus brasiliensis</i> deposited as <i>Aspergillus niger</i>	WLRI 034(120)	16404™
<i>Aspergillus brasiliensis</i> deposited as <i>Aspergillus niger</i>	SN 26	9642™
<i>Aspergillus brasiliensis</i> deposited as <i>Aspergillus niger</i>	SN 26	9642-U™
<i>Aspergillus flavus</i> var. <i>flavus</i>	WB 1957	16883™
<i>Aspergillus flavus</i>	MCV-C#1	204304™
<i>Aspergillus fumigatus</i> deposited as <i>Aspergillus fumigatus</i>	FG 1432	MYA-3627™
<i>Aspergillus niger</i> 16404™ renamed (2008) as <i>Aspergillus brasiliensis</i>		16404™
<i>Aspergillus niger</i>		10535™
<i>Aspergillus niger</i>	WB 326	16888™
<i>Aspergillus niger</i>	4247	6275™
<i>Aspergillus niger</i>		64958™
<i>Aspergillus oryzae</i> deposited as <i>Aspergillus flavus</i>	NRRL 484	10124™
<i>Aspergillus restrictus</i>	FRR 2176	42693™
<i>Aspergillus ustus</i>	NRRL A-310	10760™
<i>Bacillus atrophaeus</i> deposited as <i>Bacillus subtilis</i>	NCTC 10073	51189™
<i>Bacillus atrophaeus</i> deposited as <i>Bacillus subtilis</i> var. <i>niger</i>	NRS 1221A	9372™
<i>Bacillus cereus</i>	NRRL B-569	10876™
<i>Bacillus cereus</i>	PCI 246	13061™
<i>Bacillus cereus</i>	[BCRC 10603, CCM 2010, CCUG 7414, CIP 66.24, DSM 31, HAMBI 1887, HAMBI 1905, IAM 12605, JCM 2152, LMG 6923, NBRC 15305, NCCB 75008, NCIMB 9373, NCTC 2599, NRRL B-3711, VKM B-504]	14579™
<i>Bacillus cereus</i>	ENSP6	33019™

Catalogue

Organism	Designation	ATCC® No.
<i>Bacillus cereus</i> deposited as <i>Bacillus agri</i>	AMC 800	2™
<i>Bacillus cereus</i> deposited as <i>Bacillus mycoides</i>	FDA strain PCI 213	11778™
<i>Bacillus circulans</i>		4513™
<i>Bacillus circulans</i>	7	4516™
<i>Bacillus circulans</i> deposited as <i>Bacillus fusiformis</i>	AMC 732	61™
<i>Bacillus coagulans</i>	NRS 609	7050™
<i>Bacillus licheniformis</i>	[ATCC 11560, Damodaron P-8, LMG 7560, NRS 1415]	12759™
<i>Bacillus licheniformis</i>	[46, NCIB 9375, NCTC 10341, NRS 1264]	14580™
<i>Bacillus megaterium</i>	[BCRC 10608, CCM 2007, CCUG 1817, CIP 66.20, DSM 32, HAMB I 2018, IAM 13418, JCM 2506, KCTC 3007, LMG 7127, NBRC 15308, NCCB 75016, NCIMB 9376, NCTC 10342, NRIC 1710, NRRL B-14308, VKM B-512]	14581™
<i>Bacillus megaterium</i>	4R6259	9885™
<i>Bacillus pumilus</i>	NCTC 8241	14884™
<i>Bacillus pumilus</i>	E601	27142™
<i>Bacillus pumilus</i>	GB34	700814™
<i>Bacillus pumilus</i>	200211	BAA-1434™
<i>Bacillus subtilis</i>	NCTC 8236	11774™
<i>Bacillus subtilis</i>	PRD 66	19659™
<i>Bacillus subtilis</i>	168	23857™
<i>Bacillus subtilis</i>	BD170	33608™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> deposited as <i>Bacillus subtilis</i>	NRS 231	6633™
<i>Bacillus subtilis</i> subsp. <i>subtilis</i> deposited as <i>Bacillus subtilis</i>	Marburg strain	6051™
<i>Bacillus subtilis</i> subsp. <i>subtilis</i> deposited as <i>Bacillus subtilis</i>	Marburg strain	6051-U™
<i>Bacillus thuringiensis</i> deposited as <i>Bacillus cereus</i> subsp. <i>thuringiensis</i>	[CCUG 7429, CIP 53.137, DSM 2046, HAMB I 478, LMG 7138, NCAIM B.01292, NCCB 70008, NRRL HD-735, VKM B-1544]	10792™
<i>Bacillus thuringiensis</i> deposited as <i>Bacillus thuringiensis</i> subsp. <i>kurstaki</i>	NRRL B-3792	33679™
<i>Bacteroides fragilis</i>	VPI 2553	25285™
<i>Bacteroides fragilis</i>	RYC2056	700786™
<i>Bacteroides fragilis</i> ATCC 700786 bacteriophage	B56-3	700786-B1™
<i>Bacteroides fragilis</i> deposited as " <i>Sphaerophorus intermedium</i> "	[ICPB 3498, NCTC I0581]	23745™
<i>Bacteroides ovatus</i>	[NCTC 11153]	8483™
<i>Bacteroides ovatus</i>	Vitek 400963; NSB 24257	BAA-1296™
<i>Bacteroides ovatus</i>	Vitek 400956	BAA-1304™
<i>Bacteroides thetaiotaomicron</i>	WAL 2926	29741™
<i>Bacteroides uniformis</i> deposited as <i>Bacteroides thetaiotaomicron</i>		8492™
<i>Bacteroides ureolyticus</i>	NCTC 10941	33387™
<i>Bacteroides vulgatus</i>	WAL 1887	29327™
<i>Bacteroides vulgatus</i>	[NCTC 11154]	8482™
<i>Bartonella henselae</i> deposited as <i>Rochalimaea henselae</i>	Houston-1	49882™
<i>Bifidobacterium breve</i>	S1 (Variant a)	15700™
<i>Bordetella bronchiseptica</i>	3127	10580™
<i>Bordetella bronchiseptica</i>	NRRL B-140	4617™
<i>Bordetella parapertussis</i> deposited as <i>Acinetobacter parapertussis</i>	NCTC 5952	15311™
<i>Bordetella pertussis</i>	F	8467™
<i>Bordetella pertussis</i>	5	9340™
<i>Bordetella pertussis</i>	18323	9797™
<i>Brevibacillus agri</i> deposited as <i>Bacillus agri</i>	NRS 1219	51663™
<i>Brevibacillus brevis</i> deposited as <i>Bacillus brevis</i>	NRS 604	8246™
<i>Brevibacillus laterosporus</i> deposited as <i>Bacillus laterosporus</i>	AMC 797	64™
<i>Brevundimonas diminuta</i> deposited as <i>Pseudomonas diminuta</i>	DSM 7234	11568™
<i>Brevundimonas diminuta</i> deposited as <i>Pseudomonas</i> sp.	FDA strain PCI 818	19146™
<i>Brochothrix thermosphacta</i> deposited as <i>Microbacterium thermosphactum</i>	SW 26	11509™
<i>Budvicia aquatica</i>	85-01-010	51341™
<i>Burkholderia cepacia</i> deposited as <i>Pseudomonas cepacia</i>	UCB 717	25416™
<i>Burkholderia cepacia</i> deposited as <i>Pseudomonas kingii</i>	[NCDC 7119 (EO-1 group)]	25608™
<i>Campylobacter coli</i>	CIP 7080	33559™
<i>Campylobacter coli</i>	76-GA2	43478™
<i>Campylobacter fetus</i> subsp. <i>fetus</i> deposited as <i>Vibrio fetus</i>	NADL 1083-2255	25936™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	CIP 702	33560™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	Strain LRA 094.06.89	49943™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	Vitek 109169	BAA-1153™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	CJC-25	43430™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> deposited as <i>Campylobacter fetus</i> subsp. <i>jejuni</i>	VPI H840	29428™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> deposited as <i>Campylobacter fetus</i> subsp. <i>jejuni</i>	AS-83-79	33291™

Catalogue

Organism	Designation	ATCC® No.
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> deposited as <i>Campylobacter fetus</i> subsp. <i>jejuni</i>	AS-84-79	33292 TM
<i>Campylobacter lari</i>	NCTC 11352	35221 TM
<i>Candida albicans</i>	3147	10231 TM
<i>Candida albicans</i>	NIH 3172	14053 TM
<i>Candida albicans</i>	CBS 562	18804 TM
<i>Candida albicans</i>	Wasson	24433 TM
<i>Candida albicans</i>	H-29	26790 TM
<i>Candida albicans</i>	S-24	60193 TM
<i>Candida albicans</i>	FC18	62376 TM
<i>Candida albicans</i>	AmMS 225	66027 TM
<i>Candida albicans</i>	NCCLS 11	90028 TM
<i>Candida albicans</i>	NCCLS 67	90029 TM
<i>Candida albicans</i> deposited as <i>Candida stellatoidea</i>	[ATCC 42266, CBS 1905, CCRC 20518, CECT 1439, DBVPG 6157, IFO 1397, JCM 1543]	11006 TM
<i>Candida albicans</i> deposited as <i>Candida stellatoidea</i>	MCV52.10	36232 TM
<i>Candida albicans</i> deposited as <i>Endomyces albicans</i>	132	2091 TM
<i>Candida catenulata</i>	NRRL Y-1508	10565 TM
<i>Candida geochares</i> deposited as <i>Torulopsis geochares</i>	CBS 6870	36852 TM
<i>Candida glabrata</i>	LRA 85.10.75	64677 TM
<i>Candida glabrata</i>	303542	MYA-2950 TM
<i>Candida glabrata</i> deposited as <i>Torulopsis glabrata</i>	mutant TMAGR-23	15126 TM
<i>Candida glabrata</i> deposited as <i>Torulopsis glabrata</i>	CBS 138	2001 TM
<i>Candida glabrata</i> deposited as <i>Torulopsis glabrata</i>	AmMS 231	66032 TM
<i>Candida glabrata</i> deposited as <i>Torulopsis glabrata</i>	NCCLS 84	90030 TM
<i>Candida guilliermondii</i> deposited as <i>Monilia guilliermondii</i>	Yamadazyma guilliermondii (Wickerham) Billon-Grand, teleomorph	6260 TM
<i>Candida kefyr</i>	Vitek 302993	204093 TM
<i>Candida kefyr</i>	CBS 607	4135 TM
<i>Candida kefyr</i> deposited as <i>Candida pseudotropicalis</i>	AmMS 226	66028 TM
<i>Candida kefyr</i> deposited as <i>Torula cremoris</i>	[BI CZAS 652/6, CCY 29-8-9]	2512 TM
<i>Candida krusei</i>		14243 TM
<i>Candida krusei</i>	ST-112	34135 TM
<i>Candida lusitanae</i>	IFO 1019	34449 TM
<i>Candida lusitanae</i>	AmMS 233	66035 TM
<i>Candida magnoliae</i>	API 85 02 008	201379 TM
<i>Candida membranaefaciens</i>	API 85 03 203	201377 TM
<i>Candida parapsilosis</i>	CBS 604	22019 TM
<i>Candida parapsilosis</i>	ST-89	34136 TM
<i>Candida parapsilosis</i>	7312/52.493	90018 TM
<i>Candida tropicalis</i>	API 90 01 105	201380 TM
<i>Candida tropicalis</i>	FDA PCI M-59	13803 TM
<i>Candida tropicalis</i>	API 90 06 019	201381 TM
<i>Candida tropicalis</i>	AmMS 227	66029 TM
<i>Candida tropicalis</i> deposited as <i>Monilia murmannica</i>	[CBS 2317, CCRC 21436, NRRL Y-1378, NRRL Y-1767, VTT C-78087]	9968 TM
<i>Candida tropicalis</i> deposited as <i>Monilia tropicalis</i>	[1909, ATCC 4563, ATCC 7349, CBS 94, CCRC 20520, DSY 59, IFO 1070, IFO 1400, JCM 1541, NRRL Y-12699, NRRL Y-12968, NRRL Y-607]	750 TM
<i>Candida utilis</i> deposited as <i>Torulopsis utilis</i>	<i>Pichia jadinii</i> (Sartory et al.) Kurtzman, teleomorph	9950 TM
<i>Capnocytophaga sputigena</i>		33612 TM
<i>Cedecea neteri</i>	CDC 0621-75	33855 TM
<i>Cellulosimicrobium cellulans</i> deposited as <i>Oerskovia xanthineolytica</i>	LL G62	27402 TM
<i>Chlamydomyces pneumoniae</i> deposited as <i>Chlamydia pneumoniae</i>	TWAR strain AR-39	53592 TM
<i>Citrobacter braakii</i>	CDC 80-58	51113 TM
<i>Citrobacter braakii</i> deposited as <i>Citrobacter freundii</i>	MMS 4477	43162 TM
<i>Citrobacter freundii</i>	LRA 117.03.76	43864 TM
<i>Citrobacter freundii</i>	[ATCC 13316, NCTC 9750]	8090 TM
<i>Citrobacter koseri</i> deposited as <i>Citrobacter diversus</i>	CDC 3613-63	27156 TM
<i>Cladosporium cladosporioides</i> deposited as <i>Cladosporium herbarum</i>	IMI 45534	16022 TM
<i>Clostridium beijerinckii</i> deposited as <i>Clostridium butyricum</i>	CP Hegarty	8260 TM
<i>Clostridium bifermentans</i> deposited as <i>Bacillus centrosporogenes</i>	76	638 TM
<i>Clostridium difficile</i>	VPI 10463	43255 TM
<i>Clostridium difficile</i>	1351	43593 TM
<i>Clostridium difficile</i>	VPI 11186	700057 TM
<i>Clostridium difficile</i>	90556-M6S	9689 TM
<i>Clostridium difficile</i>	630	BAA-1382 TM
<i>Clostridium histolyticum</i>	NCTC 503	19401 TM

Catalogue

Organism	Designation	ATCC® No.
<i>Clostridium novyi</i>	1323	7659™
<i>Clostridium perfringens</i>		13124™
<i>Clostridium septicum</i>	NCTC 547	12464™
<i>Clostridium sordellii</i>	211	9714™
<i>Clostridium sporogenes</i>		19404™
<i>Clostridium sporogenes</i>	L.S. McClung 2006	11437™
<i>Clostridium sporogenes</i>	388	3584™
<i>Clostridium tetani</i>	43415	10779™
<i>Corynebacterium diphtheriae</i>	48255	11913™
<i>Corynebacterium diphtheriae</i>	5159	13812™
<i>Corynebacterium minutissimum</i>	[IFO 15361, NCTC 10288]	23348™
<i>Corynebacterium pseudodiphtheriticum</i>	153	10700™
<i>Corynebacterium pseudodiphtheriticum</i>	39555	10701™
<i>Corynebacterium renale</i>	NCTC 7448	19412™
<i>Corynebacterium striatum</i>	NSB 23348	BAA-1293™
<i>Corynebacterium striatum</i> deposited as <i>Corynebacterium flavidum</i>	NCTC 764	6940™
<i>Corynebacterium xerosis</i>		373™
<i>Cronobacter muytjensii</i> deposited as <i>Enterobacter sakazakii</i>	83-07-023	51329™
<i>Cryptocodinium cohnii</i>	PGM-1	30772™
<i>Cryptococcus albidus</i> var. <i>albidus</i>	NRRL Y-1400	10666™
<i>Cryptococcus albidus</i> var. <i>albidus</i>	CR-180	34140™
<i>Cryptococcus albidus</i>	AmMS 228	66030™
<i>Cryptococcus humicola</i> deposited as <i>Candida humicola</i>	LRA 229.03.85	64676™
<i>Cryptococcus humicola</i> deposited as <i>Candida humicola</i>	NRRL Y-1266	9949™
<i>Cryptococcus laurentii</i>	CBS 139	18803™
<i>Cryptococcus laurentii</i>	AmMS 230	66036™
<i>Cryptococcus laurentii</i>	YBC 505	76483™
<i>Cryptococcus luteolus</i>	NRRL Y-986	10671™
<i>Cryptococcus neoformans</i>	3756	14116™
<i>Cryptococcus neoformans</i>	Vitek 302898	204092™
<i>Cryptococcus neoformans</i>	CBS 132	32045™
<i>Cryptococcus neoformans</i>	NIH 76	34877™
<i>Cryptococcus neoformans</i>	AmMS 229	66031™
<i>Cryptococcus neoformans</i>	YBC 81	76484™
<i>Cryptococcus uniguttulatus</i>	AmMS 234	66033™
<i>Deinococcus radiodurans</i> deposited as <i>Micrococcus radiodurans</i>	R1 (smooth)	13939™
<i>Dekkera bruxellensis</i>	KT412	200341™
<i>Delftia acidovorans</i> deposited as <i>Pseudomonas acidovorans</i>	LRA 16.01.76	43868™
<i>Desulfovibrio vulgaris</i> deposited as <i>Spirillum desulfuricans</i> previously named	C-6	7757™
<i>Desulfovibrio desulfuricans</i> subsp. <i>desulfuricans</i>		
<i>Edwardsiella tarda</i>	CDC 1483-59	15947™
<i>Eggerthella lenta</i> deposited as <i>Eubacterium lentum</i>	MSMC 77-67	43055™
<i>Eikenella corrodens</i>	Vitek #101889	BAA-1152™
<i>Eikenella corrodens</i> deposited as <i>Bacteroides corrodens</i>	333/54-55	23834™
<i>Elizabethkingia meningoseptica</i> deposited as <i>Flavobacterium meningosepticum</i>	14	13253™
<i>Enterobacter aerogenes</i>	413	35028™
<i>Enterobacter aerogenes</i>	203	35029™
<i>Enterobacter aerogenes</i>	AmMS 264	51697™
<i>Enterobacter aerogenes</i> deposited as <i>Aerobacter aerogenes</i>	NCDC 819-56	13048™
<i>Enterobacter amnigenus</i> deposited as <i>Enterobacter cloacae</i>	C3	51816™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> subsp. <i>nov.</i> deposited as <i>Aerobacter cloacae</i>	NCDC 279-56	13047™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> subsp. <i>nov.</i> deposited as <i>Aerobacter cloacae</i>	A-8	23355™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> subsp. <i>nov.</i> deposited as <i>Enterobacter cloacae</i>	2581	35030™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> subsp. <i>nov.</i> deposited as <i>Enterobacter cloacae</i>	AmMS 204	49141™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> subsp. <i>nov.</i> deposited as <i>Enterobacter cloacae</i>	LBM 93.03.067	700323™
<i>Enterobacter gergoviae</i>	CDC 604-77	33028™
<i>Enterobacter sakazakii</i> (51329™) renamed <i>Cronobacter muytjensii</i>		51329™
<i>Enterococcus avium</i> deposited as <i>Streptococcus</i> sp.	NCTC 9938	14025™
<i>Enterococcus casseliflavus</i>	9199	700327™
<i>Enterococcus casseliflavus</i> deposited as <i>Streptococcus faecium</i> var. <i>casseliflavus</i>	20	25788™
<i>Enterococcus durans</i>	AmMS 254	49479™
<i>Enterococcus durans</i>	AmMS 206	49135™
<i>Enterococcus durans</i> deposited as <i>Streptococcus durans</i>	SD-A	11576™
<i>Enterococcus durans</i> deposited as <i>Streptococcus durans</i>	NCTC 8307	19432™
<i>Enterococcus durans</i> deposited as <i>Streptococcus durans</i>	23C2	6056™
<i>Enterococcus faecalis</i>	AmMS 161	49149™
<i>Enterococcus faecalis</i>	API 639	49452™
<i>Enterococcus faecalis</i>	UWH 1921	49532™
<i>Enterococcus faecalis</i>	UWH 1936	49533™
<i>Enterococcus faecalis</i>	AmMS 260	51188™

Catalogue

Organism	Designation	ATCC® No.
<i>Enterococcus faecalis</i>	NJ-3	51299™
<i>Enterococcus faecalis</i>	V583	700802™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	PCI 1325	14506™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	NCTC 775	19433™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	[Portland]	29212™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	CN 478	33186™
<i>Enterococcus faecalis</i> deposited as <i>Streptococcus faecalis</i>	110	7080™
<i>Enterococcus faecium</i>	MMC4	51559™
<i>Enterococcus faecium</i>	VRE	700221™
<i>Enterococcus faecium</i> deposited as <i>Streptococcus faecalis</i>	24	6057™
<i>Enterococcus faecium</i> deposited as <i>Streptococcus faecalis</i>	PRD	6569™
<i>Enterococcus faecium</i> deposited as <i>Streptococcus faecium</i>	X3	27270™
<i>Enterococcus faecium</i> deposited as <i>Streptococcus faecium</i>	LRA 55 03 77	35667™
<i>Enterococcus faecium</i> deposited as <i>Streptococcus</i> sp.	NCTC 7171	19434™
<i>Enterococcus gallinarum</i>	NCDO 2313	49573™
<i>Enterococcus gallinarum</i>	API 84-10-088	700425™
<i>Enterococcus hirae</i> deposited as <i>Streptococcus faecalis</i>	FDA M19	10541™
<i>Enterococcus hirae</i> deposited as <i>Streptococcus lactis</i>	R	8043™
<i>Enterococcus raffinosus</i> deposited as <i>Enterococcus avium</i>	AmMS 239	49464™
<i>Enterococcus saccharolyticus</i> deposited as <i>Streptococcus saccharolyticus</i>	NCDO 2594	43076™
<i>Epidermophyton floccosum</i>	ES 3115	52066™
<i>Erysipelothrix rhusiopathiae</i> deposited as <i>Erysipelothrix insidiosa</i>	NCTC 8163	19414™
<i>Escherichia coli</i>	C600	23724
<i>Escherichia coli</i>	HB101	33694
<i>Escherichia coli</i>	MacLeod	10536™
<i>Escherichia coli</i>	K-12	10798™
<i>Escherichia coli</i>		11105™
<i>Escherichia coli</i>	AMC 198	11229™
<i>Escherichia coli</i>	B	11303™
<i>Escherichia coli</i>	B	11303-U™
<i>Escherichia coli</i>	NCTC 9001	11775™
<i>Escherichia coli</i>	CDC 5624-50	12014™
<i>Escherichia coli</i>	W1485	12435™
<i>Escherichia coli</i>	C	13706™
<i>Escherichia coli</i>	W-mutant 99-1	13762™
<i>Escherichia coli</i>	NCIB 9270	14169™
<i>Escherichia coli</i>	W3100	14948™
<i>Escherichia coli</i>	FDA strain PCI 1657	15222™
<i>Escherichia coli</i>	ML308	15224™
<i>Escherichia coli</i>	C-3000	15597™
<i>Escherichia coli</i>		23631™
<i>Escherichia coli</i>	[EMG 2: K (lambda)]	23716™
<i>Escherichia coli</i>	[B, EMG 31]	23848™
<i>Escherichia coli</i>	K-12	25404™
<i>Escherichia coli</i>	FDA strain Seattle 1946	25922™
<i>Escherichia coli</i>	K-12 RV308	31608™
<i>Escherichia coli</i>	K380, 81E1301	33605™
<i>Escherichia coli</i>	1532	35218™
<i>Escherichia coli</i>	API 1157	35421™
<i>Escherichia coli</i>	NCTC 86 (original Escherich strain)	4157™
<i>Escherichia coli</i>	V1076	51446™
<i>Escherichia coli</i>	DUP-101	51739™
<i>Escherichia coli</i>	BDMS 605	51755™
<i>Escherichia coli</i>	DG1H9	51813™
<i>Escherichia coli</i>	27	53498™
<i>Escherichia coli</i>	C Na1(r)	700078™
<i>Escherichia coli</i>	CN13	700609™
<i>Escherichia coli</i>	HS(pFamp)R	700891™
<i>Escherichia coli</i>	MG1655	700926™
<i>Escherichia coli</i>	Crooks	8739™
<i>Escherichia coli</i>	[397E, CCM 2024, DSM 1116, IFO 13500, NCIB 8666, NRRL B-766, W]	9637™
<i>Escherichia coli</i> bacteriophage	T2	11303-B2™
<i>Escherichia coli</i> bacteriophage	T4	11303-B4™
<i>Escherichia coli</i> bacteriophage	phi X174	13706-B1™
<i>Escherichia coli</i> bacteriophage	MS2	15597-B1™
<i>Escherichia coli</i> bacteriophage	P1	25404-B1™
<i>Eurotium rubrum</i>	FRR 1968	42690™
<i>Exiguobacterium aurantiacum</i> deposited as <i>Corynebacterium</i> sp.	LRA 143.06.89	49676™
<i>Finegoldia magna</i> deposited as <i>Peptococcus magnus</i>	WAL2508	29328™

Catalogue

Organism	Designation	ATCC® No.
<i>Flavobacterium</i> sp.	HC6	51823™
<i>Fluoribacter bozemanii</i> deposited as <i>Legionella bozemanii</i>	WIGA	33217™
<i>Fusarium solani</i>	FIV/74	36031™
<i>Fusobacterium necrophorum</i> subsp. <i>necrophorum</i>	VPI 2891	25286™
<i>Fusobacterium nucleatum</i> subsp. <i>nucleatum</i> Knorr	VPI 4355	25586™
<i>Fusobacterium nucleatum</i> subsp. <i>polymorphum</i> deposited as <i>Fusobacterium polymorphum</i>	[NCTC 10562]	10953™
<i>Fusobacterium varium</i>		27725™
<i>Gardnerella vaginalis</i>	AmMS 117	49145™
<i>Gardnerella vaginalis</i> deposited as <i>Haemophilus vaginalis</i>	594	14018™
<i>Gardnerella vaginalis</i> deposited as <i>Haemophilus vaginalis</i>	317	14019™
<i>Gemella morbillorum</i> deposited as <i>Streptococcus morbillorum</i>	VPI 5424	27824™
<i>Geobacillus stearothermophilus</i> deposited as <i>Bacillus calidolactis</i>	NRS T15	10149™
<i>Geobacillus stearothermophilus</i> deposited as <i>Bacillus stearothermophilus</i>	NCA 1805	12978™
<i>Geobacillus stearothermophilus</i> deposited as <i>Bacillus stearothermophilus</i>	NCA 26	12980™
<i>Geobacillus stearothermophilus</i> deposited as <i>Bacillus stearothermophilus</i>		7953™
<i>Geotrichum candidum</i>		34614™
<i>Geotrichum capitatum</i> deposited as <i>Trichosporon capitatum</i>	NRRL Y-1487	10663™
<i>Geotrichum capitatum</i> deposited as <i>Trichosporon capitatum</i>	CBS 5882	28576™
<i>Haemophilus actinomycetemcomitans</i> deposited as <i>Actinobacillus actinomycetemcomitans</i>	CDC A1916	29523™
<i>Haemophilus aphrophilus</i>	NCTC 5886	19415™
<i>Haemophilus haemolyticus</i>	NCTC 10659	33390™
<i>Haemophilus influenzae</i>	AMC 36-A-1	10211™
<i>Haemophilus influenzae</i>	NCTC 4560	19418™
<i>Haemophilus influenzae</i>		33533™
<i>Haemophilus influenzae</i>	R387	33930™
<i>Haemophilus influenzae</i>	[CIP 103777]	35056™
<i>Haemophilus influenzae</i>	MMS 34D-F	43163™
<i>Haemophilus influenzae</i>	AmMS 120	49144™
<i>Haemophilus influenzae</i>	TD-4	49247™
<i>Haemophilus influenzae</i>	L-378	49766™
<i>Haemophilus influenzae</i>	NCTC 8143	33391™
<i>Haemophilus influenzae</i>	3591	43065™
<i>Haemophilus influenzae</i>	AMC 36-A-3	9006™
<i>Haemophilus influenzae</i>	AMC 36-A-5	9007™
<i>Haemophilus parahaemolyticus</i> deposited as <i>Haemophilus haemolyticus</i>	536	10014™
<i>Haemophilus parainfluenzae</i>	H30	51505™
<i>Haemophilus parainfluenzae</i>	429	7901™
<i>Haemophilus paraphrophilus</i>	AmMS 115	49146™
<i>Haemophilus paraphrophilus</i>	LRA 116.07.89	49917™
<i>Haemophilus paraphrophilus</i>	NCTC 10557	29241™
<i>Haemophilus somnus</i>	8917974	700025™
<i>Hafnia alvei</i> deposited as <i>Enterobacter</i> sp.	C2	51815™
<i>Halobacterium salinarum</i>	1	19700™
<i>Helicobacter pylori</i> deposited as <i>Campylobacter pyloridis</i>	NCTC 11637	43504™
<i>Issatchenkia orientalis</i> deposited as <i>Candida krusei</i>	[ATCC 749, CBS 573, CCRC 20514, IFO 1064, IFO 1395, JCM 1609, NRRL Y-413, NRRL Y-7179]	6258™
<i>Kingella denitrificans</i>	NCTC 10995	33394™
<i>Klebsiella oxytoca</i>	Pasco 201	43086™
<i>Klebsiella oxytoca</i>	MMS 4377	43165™
<i>Klebsiella oxytoca</i>	LRA 06.01.73	43863™
<i>Klebsiella oxytoca</i>	AmMS 101	49131™
<i>Klebsiella oxytoca</i>	C4	51817™
<i>Klebsiella oxytoca</i>	LBM 90.11.033	700324™
<i>Klebsiella oxytoca</i> deposited as <i>Aerobacter aerogenes</i>	NRRL B-199	8724™
<i>Klebsiella oxytoca</i> deposited as <i>Klebsiella pneumoniae</i>	479-2	13182™
<i>Klebsiella pneumoniae</i>	ART 2008133	BAA-1705™
<i>Klebsiella pneumoniae</i>	AIS 2007023	BAA-1706™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>		10031™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	NCTC 9633	13883™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	S 389	27736™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	PS-53	31488™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	C122	33495™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	LRA 24 11 73	35657™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	[CIP 104034, NCDC 400-68]	9997™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> deposited as <i>Aerobacter aerogenes</i>	[NCDC 379-68, NCIB 8152]	8308™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> deposited as <i>Escherichia coli</i> var. <i>communior</i>	[CIP 104216, NCIB 10341]	4352™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> deposited as <i>Klebsiella aerogenes</i>	NCTC 8172	13882™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> deposited as <i>Klebsiella pneumoniae</i>	K6	700603™
<i>Kloeckera apiculata</i> var. <i>apis</i>	CBS 2591	32857™

Catalogue

Organism	Designation	ATCC® No.
<i>Kocuria kristinae</i>	Vitek #12917	BAA-752™
<i>Kocuria rhizophila</i> deposited as <i>Sarcina lutea</i>	FDA strain PCI 1001/E	15957™
<i>Kocuria rhizophila</i> deposited as <i>Sarcina lutea</i>		533™
<i>Kocuria rhizophila</i> deposited as <i>Sarcina lutea</i>	FDA strain PCI 1001/D	9341a™
<i>Kocuria rhizophila</i> deposited as <i>Sarcina lutea</i>	FDA strain PCI 1001	9341™
<i>Kocuria rosea</i> deposited as <i>Micrococcus rubens</i>	[CCM 679, IAM 1315, IFO 3768, NCIB 11696, NCTC 7523]	186™
<i>Kocuria varians</i> deposited as <i>Micrococcus varians</i>	D757	51820™
<i>Lactobacillus acidophilus</i>	Scav	4356™
<i>Lactobacillus acidophilus</i> deposited as <i>Bacillus acidophilus</i>	[43]	314™
<i>Lactobacillus brevis</i>	Bb14	14869™
<i>Lactobacillus casei</i> deposited as <i>Lactobacillus casei</i> subsp. <i>casei</i>		334™
<i>Lactobacillus casei</i> deposited as <i>Lactobacillus casei</i> subsp. <i>casei</i>	03	393™
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> deposited as <i>Lactobacillus bulgaricus</i>	Lb14	11842™
<i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> deposited as <i>Lactobacillus leichmannii</i>	326	4797™
<i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> deposited as <i>Lactobacillus leichmannii</i>	313	7830™
<i>Lactobacillus fermentum</i>	36	9338™
<i>Lactobacillus gasserii</i> deposited as <i>Lactobacillus acidophilus</i>	1SL4	19992™
<i>Lactobacillus paracasei</i> subsp. <i>paracasei</i>	NCDO 206	BAA-52™
<i>Lactobacillus plantarum</i>	Lp 39	14917™
<i>Lactobacillus plantarum</i> deposited as <i>Lactobacillus arabinosus</i>	17-5	8014™
<i>Lactobacillus rhamnosus</i> deposited as <i>Lactobacillus acidophilus</i>	GG	53103™
<i>Lactobacillus rhamnosus</i> deposited as <i>Lactobacillus casei</i>	4R2127	9595™
<i>Lactobacillus rhamnosus</i> deposited as <i>Lactobacillus casei</i>	[BUCSAV 227, M. Rogosa, V300, M.E. Sharpe H2, NCDO 243, NCIB 6375, NCIB 8010, NCTC 6375, NRC 488, P.A. Hansen 300, R.P. Tittler 300]	7469™
<i>Lactobacillus sakei</i> subsp. <i>sakei</i> Katagiri deposited as <i>Lactobacillus sakei</i>	T.S	15521™
<i>Lactococcus lactis</i> deposited as <i>Enterococcus faecium</i>	9155	49032™
<i>Lactococcus lactis</i> subsp. <i>cremoris</i> deposited as <i>Streptococcus cremoris</i>	NCDO 607	19257™
<i>Lactococcus lactis</i> subsp. <i>lactis</i> deposited as <i>Streptococcus lactis</i>	Berridge X 13	11454™
<i>Lactococcus lactis</i> subsp. <i>lactis</i> deposited as <i>Streptococcus</i> sp.	NCTC 6681	19435™
<i>Leclercia adecarboxylata</i>	LBM 84.05.644	700325™
<i>Leclercia adecarboxylata</i> deposited as <i>Escherichia adecarboxylata</i>	1783	23216™
<i>Legionella longbeachae</i>	Long Beach 4	33462™
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i>	Philadelphia-1	33152™
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i>	Bloomington-2	33155™
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i> deposited as <i>Legionella pneumophila</i>	Knoxville-1	33153™
<i>Leishmania braziliensis</i>	MHOM/BR/75/M2903	50135™
<i>Listeria grayi</i>	API 88-05-013	700545™
<i>Listeria grayi</i> deposited as <i>Listeria murrayi</i>	F-9	25401™
<i>Listeria innocua</i>	SLCC 3379	33090™
<i>Listeria innocua</i>	DUP-104	51742™
<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> deposited as <i>Listeria ivanovii</i>	BE3229	BAA-139™
<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> deposited as <i>Listeria monocytogenes</i>	Li 1979	19119™
<i>Listeria monocytogenes</i>	1071/53	13932™
<i>Listeria monocytogenes</i>	53 XXIII	15313™
<i>Listeria monocytogenes</i>	Li 20	19111™
<i>Listeria monocytogenes</i>	Li 21	19112™
<i>Listeria monocytogenes</i>	Li 23	19114™
<i>Listeria monocytogenes</i>	Li 2	19115™
<i>Listeria monocytogenes</i>	Li 2109	19118™
<i>Listeria monocytogenes</i>	Gibson	7644™
<i>Listeria monocytogenes</i>	[API 87-10-028, NSB 22072]	BAA-751™
<i>Listeria seeligeri</i>	CIP 100100	35967™
<i>Lysinibacillus sphaericus</i> deposited as <i>Bacillus sphaericus</i>	25	4525™
<i>Macrococcus caseolyticus</i> deposited as <i>Staphylococcus cohnii</i>	LRA 041 575	35662™
<i>Magnetospirillum</i> sp.	AMB-1	700264™
<i>Mannheimia haemolytica</i> deposited as <i>Pasteurella haemolytica</i>	NCTC 9380	33396™
<i>Methylobacterium mesophilicum</i> deposited as <i>Pseudomonas mesophila</i>	A47	29983™
<i>Microbacterium esteraromaticum</i>	MC3	51822™
<i>Microbacterium testaceum</i> deposited as <i>Brevibacterium testaceum</i>	Rp-3	15829™
<i>Micrococcus luteus</i>	AmMS 533	49732™
<i>Micrococcus luteus</i> deposited as <i>Micrococcus flavus</i>	130.21	10240™
<i>Micrococcus luteus</i> deposited as <i>Micrococcus lysodeikticus</i>	[ATCC 15307, CCM 169, CIP A-270, IAM 1056, IFO 333, NCIB 9278, NCTC 2665, NRRL B-287]	4698™
<i>Micrococcus luteus</i> deposited as <i>Sarcina subflava</i>	Mercedita	7468™
<i>Micrococcus luteus</i> ; renamed <i>Kocuria rhizophila</i>	refer to strains 9341™, 9341a™, 15957™	
<i>Micrococcus</i> sp.	API 82-06-057	700405™

Catalogue

Organism	Designation	ATCC® No.
<i>Micrococcus</i> sp.	MC7	51819™
<i>Microsporium canis</i>	Vasquez	11621™
<i>Microsporium canis</i>	A 3697 (2)	36299™
<i>Microsporium gypseum</i>	VH/3032	24102™
<i>Moraxella (Branhamella) catarrhalis</i> deposited as <i>Branhamella catarrhalis</i>	Ne 11	25238™
<i>Moraxella (Branhamella) catarrhalis</i> deposited as <i>Branhamella catarrhalis</i>	N9	25240™
<i>Moraxella (Branhamella) catarrhalis</i> deposited as <i>Branhamella catarrhalis</i>	AmMS 116	49143™
<i>Moraxella (Branhamella) catarrhalis</i> deposited as <i>Neisseria catarrhalis</i>	NCTC 4103	23246™
<i>Moraxella (Branhamella) catarrhalis</i> deposited as <i>Neisseria catarrhalis</i>	20	8176™
<i>Moraxella (Moraxella) nonliquefaciens</i> deposited as <i>Moraxella nonliquefaciens</i>	A170	17953™
<i>Moraxella (Moraxella) osloensis</i> deposited as <i>Mima polymorpha</i>	CDC	10973™
<i>Morganella morganii</i> subsp. <i>morganii</i> deposited as <i>Proteus morganii</i>	M4	25829™
<i>Morganella morganii</i> subsp. <i>morganii</i> deposited as <i>Proteus morganii</i>	M11	25830™
<i>Mucor racemosus</i> Fresenius	NRRL 6341	42647™
<i>Mycobacterium avium</i> subsp. <i>avium</i>	1982	15769™
<i>Mycobacterium avium</i> subsp. <i>avium</i>	Vet. 1387	25291™
<i>Mycobacterium fortuitum</i> subsp. <i>fortuitum</i>	[TMC 1529]	6841™
<i>Mycobacterium gordonae</i>	[L. Wayne W-1609, P-15, TMC 1324]	14470™
<i>Mycobacterium intracellulare</i>	3600	13950™
<i>Mycobacterium kansasii</i>	G133 Bostrom	12478™
<i>Mycobacterium peregrinum</i>	MF1353	700686™
<i>Mycobacterium phlei</i>	Timoteo	11758™
<i>Mycobacterium smegmatis</i>	W-113	14468™
<i>Mycobacterium smegmatis</i>	NCTC 8159	19420™
<i>Mycobacterium smegmatis</i>		607™
<i>Mycobacterium smegmatis</i>	mc(2)155	700084™
<i>Mycobacterium terrae</i>	W-45	15755™
<i>Mycobacterium tuberculosis</i>	H37Ra	25177™
<i>Mycoplasma arginini</i> deposited as <i>Mycoplasma arginini</i>	G230	23838™
<i>Mycoplasma fermentans</i> deposited as <i>Mycoplasma fermentans</i>	PG18	19989™
<i>Mycoplasma gallisepticum</i> deposited as <i>Mycoplasma gallisepticum</i>	[NCTC 10115, PG 31, X95]	19610™
<i>Mycoplasma genitalium</i>	G37	33530™
<i>Mycoplasma hyorhinis</i> deposited as <i>Mycoplasma hominis</i>	[H27]	15488™
<i>Mycoplasma hyorhinis</i> deposited as <i>Mycoplasma hominis</i>	PG21	23114™
<i>Mycoplasma hyorhinis</i> deposited as <i>Mycoplasma hyorhinis</i>	BTS-7	17981™
<i>Mycoplasma hyorhinis</i> deposited as <i>Mycoplasma hyorhinis</i>	DBS 1050	29052™
<i>Mycoplasma orale</i> deposited as <i>Mycoplasma orale</i>	CH 19299	23714™
<i>Mycoplasma pneumoniae</i> deposited as <i>Mycoplasma pneumoniae</i>	FH strain of Eaton Agent	15531™
<i>Mycoplasma synoviae</i>	WVU 1853	25204™
<i>Mycoides odoratus</i> deposited as <i>Flavobacterium odoratum</i>	[CCUG 7321, CIP 103105, DSM 2801, IFO 14945, JCM 7458, LMG 1233, NCTC 11036]	4651™
<i>Neisseria gonorrhoeae</i>	NCTC 8375	19424™
<i>Neisseria gonorrhoeae</i>		31426™
<i>Neisseria gonorrhoeae</i>	CDC Ng-116	43069™
<i>Neisseria gonorrhoeae</i>	CDC Ng-98	43070™
<i>Neisseria gonorrhoeae</i>	F-18	49226™
<i>Neisseria gonorrhoeae</i>	WHO V	49498™
<i>Neisseria gonorrhoeae</i>	NHI 1	49981™
<i>Neisseria lactamica</i>	NCDC A5906	23971™
<i>Neisseria lactamica</i>	NCDC A7515	23970™
<i>Neisseria lactamica</i>	AmMS 118	49142™
<i>Neisseria meningitidis</i>	M1027	13077™
<i>Neisseria meningitidis</i>	M2092	13090™
<i>Neisseria meningitidis</i>	M1628	13102™
<i>Neisseria meningitidis</i>	M-112	35561™
<i>Neisseria mucosa</i>	X71	19695™
<i>Neisseria mucosa</i>	AmMS 138	49233™
<i>Neisseria perflava</i> deposited as <i>Neisseria subflava</i>	28	14799™
<i>Neisseria sicca</i>	CN	29193™
<i>Neisseria sicca</i>	AMC 14-D-1	9913™
<i>Nitrosomonas europaea</i>	EL Schmidt	19718™
<i>Nocardia brasiliensis</i>	IMRU 845	19296™
<i>Novosphingobium capsulatum</i> deposited as <i>Flavobacterium capsulatum</i>	28	14666™
<i>Ochrobactrum anthropi</i>	Biolog 1185	49687™
<i>Ochrobactrum anthropi</i>	Vitek 109243	BAA-749™
<i>Octosporomyces octosporus</i> deposited as <i>Schizosaccharomyces octosporus</i>	CBS 371	4206™
<i>Oligella ureolytica</i>	CCUG 1465	43534™
<i>Oligella ureolytica</i>	CCUG 17791	43535™
<i>Oligella urethralis</i> deposited as <i>Mima polymorpha</i> var. <i>oxidans</i>	CDC 7603	17960™

Catalogue

Organism	Designation	ATCC® No.
<i>Paecilomyces marquandii</i>	NRRL 901	10525™
<i>Paenibacillus gordonae</i> deposited as <i>Bacillus gordonae</i>	Q1	29948™
<i>Paenibacillus larvae</i> subsp. <i>larvae</i> deposited as <i>Bacillus larvae</i>	846	9545™
<i>Paenibacillus macerans</i> deposited as <i>Bacillus macerans</i>	18	8509™
<i>Paenibacillus polymyxa</i> deposited as <i>Bacillus polymyxa</i>	LRA 88.01.80	43865™
<i>Paenibacillus polymyxa</i> deposited as <i>Bacillus polymyxa</i>	[Difco Labs. 8277, N.R. Smith NRS 279, Vitek #200115]	7070™
<i>Paenibacillus polymyxa</i> deposited as <i>Bacillus polymyxa</i>	[BUCSAV 162, CCM 1459, NCIB 8158, NCTC 10343, NRS 1105]	842™
<i>Parabacteroides distasonis</i> deposited as <i>Bacteroides distasonis</i>	[NCTC 11152]	8503™
<i>Parabacteroides distasonis</i> deposited as <i>Bacteroides distasonis</i>	Vitek 400127	BAA-1295™
<i>Parvimonas micra</i> formerly <i>Micromonas micros</i> and <i>Peptostreptococcus micros</i> deposited as <i>Peptostreptococcus micros</i>	VPI 5464	33270™
<i>Pasteurella aerogenes</i>	P-172-71	27883™
<i>Pasteurella multocida</i> subsp. <i>multocida</i>	P-931	12945™
<i>Pasteurella multocida</i> subsp. <i>multocida</i>	NCTC 10322	43137™
<i>Pediococcus acidilactici</i> deposited as <i>Leuconostoc mesenteroides</i>		8042™
<i>Pediococcus damnosus</i>	NCDO 1832	29358™
<i>Penicillium aurantiogriseum</i> deposited as <i>Penicillium cyclopium</i>	IMI 19759	16025™
<i>Penicillium aurantiogriseum</i> deposited as <i>Penicillium puberulum</i>	NRRL 1889	8732™
<i>Penicillium chrysogenum</i>	[CBS 277.47, IFO 4626, IMI 37767, PRL 1910, QM 943, Wis. Q-176]	10002™
<i>Penicillium chrysogenum</i>	NRRL 807	10106™
<i>Penicillium chrysogenum</i>	Wis. 49-133	11709™
<i>Penicillium chrysogenum</i> deposited as <i>Penicillium notatum</i>	NRRL 832	9179™
<i>Penicillium chrysogenum</i> var. <i>dipodomys</i>	NRRL 13485	64185™
<i>Penicillium variabile</i>	NRRL 3765	32333™
<i>Peptoniphilus asaccharolyticus</i> deposited as <i>Peptococcus aerogenes</i>	UW 228	14963™
<i>Peptoniphilus asaccharolyticus</i> deposited as <i>Peptococcus asaccharolyticus</i>	WAL 3218	29743™
<i>Peptostreptococcus anaerobius</i>	VPI 4330	27337™
<i>Pityrosporum ovale</i>		12078™
<i>Plesiomonas shigelloides</i>	GNI 14	51903™
<i>Plesiomonas shigelloides</i> deposited as <i>Aeromonas shigelloides</i>	CDC 3085-55	14029™
<i>Porphyromonas gingivalis</i> deposited as <i>Bacteroides gingivalis</i>	2561	33277™
<i>Porphyromonas gingivalis</i> deposited as <i>Bacteroides gingivalis</i>	W50	53978™
<i>Porphyromonas levii</i> deposited as <i>Bacteroides melaninogenicus</i>	1	29147™
<i>Prevotella loescheii</i> deposited as <i>Bacteroides oratus</i>	8B	15930™
<i>Prevotella melaninogenica</i> deposited as <i>Bacteroides melaninogenicus</i> subsp. <i>melaninogenicus</i>	VPI 2381	25845™
<i>Prevotella nigrescens</i> deposited as <i>Bacteroides intermedius</i>	VPI 8944	33563™
<i>Propionibacterium acnes</i> deposited as <i>Corynebacterium acnes</i>	Gerath	11827™
<i>Propionibacterium acnes</i> deposited as <i>Corynebacterium acnes</i>	NCTC 737	6919™
<i>Propionibacterium granulosum</i>	VPI 0507	25564™
<i>Proteus hauseri</i> deposited as <i>Proteus vulgaris</i>	NCTC 4125 strain Lehmann	13315™
<i>Proteus mirabilis</i>	D1	12453™
<i>Proteus mirabilis</i>	FDA strain PCI 765	14153™
<i>Proteus mirabilis</i>	NCDC 2059-70	25933™
<i>Proteus mirabilis</i>	CDC 73-57740	29245™
<i>Proteus mirabilis</i>	CDC PR 14	29906™
<i>Proteus mirabilis</i>	571101	33583™
<i>Proteus mirabilis</i>	LRA 08 01 73	35659™
<i>Proteus mirabilis</i>	CDC S-17	43071™
<i>Proteus mirabilis</i> deposited as <i>Proteus ammoniae</i>		7002™
<i>Proteus vulgaris</i>		33420™
<i>Proteus vulgaris</i>	AmMS 105	49132™
<i>Proteus vulgaris</i>		6380™
<i>Proteus vulgaris</i>		8427™
<i>Proteus vulgaris</i> deposited as <i>Bacterium proteus</i>	NCTC 4636	6896™
<i>Prototheca wickerhamii</i>	NRRL YB-4330	16529™
<i>Providencia alcalifaciens</i>	GNI 3	51902™
<i>Providencia alcalifaciens</i>		9886™
<i>Providencia rettgeri</i> deposited as <i>Shigella rettgeri</i>	[ATCC 6907, IFO 13501, NCTC 1501]	9250™
<i>Providencia stuartii</i>	495	33672™
<i>Providencia stuartii</i>	AmMS 256	49809™
<i>Pseudomonas aeruginosa</i>	[CCEB 481, MDB strain BU 277, NCIB 8295, NCPPB 1965, NCTC 10332, NRRL B-771, R.Hugh 815]	10145™
<i>Pseudomonas aeruginosa</i>	[CCEB 481, MDB strain BU 277, NCIB 8295, NCPPB 1965, NCTC 10332, NRRL B-771, R.Hugh 815]	10145-U™

Catalogue

Organism	Designation	ATCC® No.
<i>Pseudomonas aeruginosa</i>	PRD-10	15442™
<i>Pseudomonas aeruginosa</i>	1C	15692™
<i>Pseudomonas aeruginosa</i>	NCTC 6750	19429™
<i>Pseudomonas aeruginosa</i>		25619™
<i>Pseudomonas aeruginosa</i>	NCTC 10662	25668™
<i>Pseudomonas aeruginosa</i>	Boston 41501	27853™
<i>Pseudomonas aeruginosa</i>	C118	35032™
<i>Pseudomonas aeruginosa</i>	API 1099	35422™
<i>Pseudomonas aeruginosa</i>	API 1143	35554™
<i>Pseudomonas aeruginosa</i>		9027™
<i>Pseudomonas aeruginosa</i>	GMC #6	BAA-427™
<i>Pseudomonas fluorescens</i>	NCTC 10038	13525™
<i>Pseudomonas fluorescens</i>	1	17386™
<i>Pseudomonas fluorescens</i>	P17	49642™
<i>Pseudomonas fluorescens</i>	AmMS 257	49838™
<i>Pseudomonas putida</i>	KT2440	47054™
<i>Pseudomonas fluorescens</i>	AmMS 106	49128™
<i>Pseudomonas putida</i> deposited as <i>Pseudomonas fluorescens</i>	[A.3.12, ATCC 23467, NCIB 9494, NCTC 10936, R.Y. Stanier 90]	12633™
<i>Pseudomonas</i> sp.	PS8	51821™
<i>Pseudomonas stutzeri</i>	221	17588™
<i>Ralstonia pickettii</i> deposited as <i>Pseudomonas cepacia</i>	AmMS 155	49129™
<i>Raoultella terrigena</i> deposited as <i>Klebsiella terrigena</i>	CUETM 77-176	33257™
<i>Rhizobium rhizogenes</i> deposited as <i>Agrobacterium rhizogenes</i>	[CIP 104786]	15834™
<i>Rhizopus stolonifer</i>	UNB-1	14037™
<i>Rhizopus stolonifer</i> deposited as <i>Rhizopus nigricans</i>	C-25A	6227a™
<i>Rhizopus stolonifer</i> deposited as <i>Rhizopus nigricans</i>	C-25B	6227b™
<i>Rhodococcus equi</i> deposited as <i>Corynebacterium equi</i>	NCTC 1621	6939™
<i>Rhodotorula mucilaginosa</i> var. <i>mucilaginosa</i> , deposited as <i>Rhodotorula rubra</i>	AmMS 232	66034™
<i>Rhodotorula mucilaginosa</i> var. <i>mucilaginosa</i> , deposited as <i>Rhodotorula rubra</i>	NRRL Y-1592	9449-U™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces carlsbergensis</i>	4228	9080™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces cerevisiae</i>	S288C	204508™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces cerevisiae</i>	NRRL Y-53	2601™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces cerevisiae</i>	[CBS 1368, CCRC 21727, DSM 2155, IFO 1346, NCYC 79, NRRL Y-977]	7754™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces cerevisiae</i>		9763™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces ellipsoideus</i>	[657, IFO 1661]	4098™
<i>Saccharomyces cerevisiae</i> deposited as <i>Saccharomyces sake</i>	B73.P	32701™
<i>Salmonella bongori</i> deposited as <i>Salmonella enterica</i> subsp. <i>bongori</i>	CIP 82.33	43975™
<i>Salmonella enterica</i> subsp. <i>arizonae</i> deposited as <i>Arizona arizonae</i>	NCTC 8297	13314™
<i>Salmonella enterica</i> subsp. <i>diarizonae</i>	0563-95	BAA-216™
<i>Salmonella enterica</i> subsp. <i>diarizonae</i> deposited as <i>Arizona hinshawii</i>	CDC 656/75	29226™
<i>Salmonella enterica</i> subsp. <i>enterica</i>	0267-95	BAA-215™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Abaetetuba</i> deposited as <i>Salmonella</i> sp.		35640™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Anatum</i> deposited as <i>Salmonella</i> sp.	5101	9270™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i> deposited as <i>Salmonella choleraesuis</i> subsp. <i>choleraesuis</i> serovar <i>Choleraesuis</i>	ETS 34	10708™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i> deposited as <i>Salmonella choleraesuis</i> subsp. <i>choleraesuis</i> serovar <i>Choleraesuis</i>	NCTC 5735	13312™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i> deposited as <i>Salmonella choleraesuis</i> subsp. <i>choleraesuis</i> serovar <i>Choleraesuis</i>		7001™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i> var. <i>Kunzendorf</i> deposited as <i>Salmonella choleraesuis</i> subsp. <i>Kunzendorf</i>	CDC (5210) (37)	12011™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Enteritidis</i> deposited as <i>Salmonella enteritidis</i>	CDC K-1891	13076™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Infantis</i> deposited as <i>Salmonella infantis</i>	DUP-103	51741™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Paratyphi A</i> deposited as <i>Salmonella paratyphi-A</i>		9150™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Pullorum</i> deposited as <i>Aeromonas</i> sp.	NRRL B-663	13036™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Pullorum</i> deposited as <i>Salmonella</i> sp.	X-12	19945™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Tallahassee</i> deposited as <i>Salmonella</i> sp.	CDC	12002™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> bacteriophage	P22	19585-B1™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella choleraesuis</i> subsp. <i>choleraesuis</i> serotype <i>Typhimurium</i>	U24	51812™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella enterica</i> subsp. I serotype <i>Typhimurium</i>	0267-95	BAA-215™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	NCTC 74	13311™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	CDC 6516-60	14028™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	LT2	19585™

Catalogue

Organism	Designation	ATCC® No.
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	met-A-15	25241™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	TA 1535	29629™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	TA 1537	29630™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	TA 1538	29631™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> deposited as <i>Salmonella typhimurium</i>	Leu 130	49416™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Vellore</i> deposited as <i>Salmonella</i> sp.	V1796	15611™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Paratyphi B</i> deposited as <i>Salmonella schottmuelleri</i>		8759™
<i>Schizochytrium</i> sp.	S31	20888™
<i>Scopulariopsis acremonium</i>	UAMH 541	58636™
<i>Serratia liquefaciens</i>	CDC 1284-57	27592™
<i>Serratia marcescens</i>	BS 303	13880™
<i>Serratia marcescens</i>	NRRL B-1481	14041™
<i>Serratia marcescens</i>	PCI 1107	14756™
<i>Serratia marcescens</i>	LRA 13.05.73	43861™
<i>Serratia marcescens</i>	LRA 04.03.73	43862™
<i>Serratia marcescens</i>	NRS 116-175	8100™
<i>Serratia odorifera</i>	1073	33077™
<i>Serratia rubidaea</i>	E1344	33670™
<i>Shewanella putrefaciens</i> deposited as <i>Pseudomonas putrefaciens</i>	AmMS 201	49138™
<i>Shewanella putrefaciens</i> deposited as <i>Pseudomonas putrefaciens</i>	95	8071™
<i>Shigella boydii</i> deposited as <i>Shigella paradysenteriae</i>	AMC 43-G-58	9207™
<i>Shigella flexneri</i>	CDC 3591-52	12022™
<i>Shigella flexneri</i> deposited as <i>Shigella paradysenteriae</i>	AMC 43-G-68	9199™
<i>Shigella sonnei</i>	NCDC 1120-66	25931™
<i>Shigella sonnei</i>	WRAIR I virulent	29930™
<i>Shigella sonnei</i>	AMC 43-GG9	9290™
<i>Sordaria fimicola</i>	A-1, wild type g+	14517™
<i>Sordaria fimicola</i>	t, tan-spored	16096™
<i>Sphingobacterium multivorum</i> deposited as <i>Flavobacterium multivorum</i>	LRA 26 07 76	35656™
<i>Sphingobacterium spiritivorum</i> deposited as <i>Flavobacterium spiritivorum</i>	CDC E7288	33861™
<i>Sphingomonas paucimobilis</i> deposited as <i>Pseudomonas paucimobilis</i>	NCTC 11030	29837™
<i>Sphingomonas trueperi</i> deposited as <i>Pseudomonas azotocolligans</i>	DSM 7225	12417™
<i>Spiroplasma citri</i>	AZ-1729	33274™
<i>Staphylococcus aureus</i>	Vitek #8753	BAA-1026™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Wood 46	10832™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	S13	11632™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	NCTC 8530	12598™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	NCTC 8532	12600™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	NCTC 8532	12600-U™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Newman D2C	25904™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	655	27660™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Wichita	29213™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	FDA	29737™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	I.J.7	31153™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	328	33591™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	1063	33592™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	SS697	33862™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	F-182	43300™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	AmMS 241	49476™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	V8	49775™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	91089	51153™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	DUP-102	51740™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	FDA 209P	6538P™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	FDA 209	6538™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	3R7089 strain Oxford	9144™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Micrococcus pyogenes</i>	PCI 1209/N	12692™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus aureus</i>	Seattle 1945	25923™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus aureus</i>	Mu3	700698™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus aureus</i>	Mu50	700699™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus aureus</i>	UT 25	BAA-976™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus aureus</i>	UT 32	BAA-977™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> deposited as <i>Staphylococcus</i> sp.	FDA 196E	13565™
<i>Staphylococcus capitis</i> deposited as <i>Sarcina aurantiaca</i>	[NRRL B-2616]	146™
<i>Staphylococcus capitis</i> subsp. <i>capitis</i>	LRA 360 677	35661™

Catalogue

Organism	Designation	ATCC® No.
<i>Staphylococcus epidermidis</i>	FDA strain PCI 1200	12228™
<i>Staphylococcus epidermidis</i>	Fussel	14990™
<i>Staphylococcus epidermidis</i>	255-01B	29887™
<i>Staphylococcus epidermidis</i>	RP62A	35984™
<i>Staphylococcus epidermidis</i>	AmMS 205	49134™
<i>Staphylococcus epidermidis</i>	AmMS 242	49461™
<i>Staphylococcus epidermidis</i>	CCF 15990	51625™
<i>Staphylococcus epidermidis</i>	C078	700296™
<i>Staphylococcus epidermidis</i> deposited as <i>Staphylococcus saprophyticus</i>	MDB strain BS 747	13518™
<i>Staphylococcus gallinarum</i>	AP 91-07-073	700401™
<i>Staphylococcus haemolyticus</i>	SM 131	29970™
<i>Staphylococcus hominis</i> subsp. <i>hominis</i>	DM 122	27844™
<i>Staphylococcus intermedius</i>	H11/68	29663™
<i>Staphylococcus lentus</i>	API 86-01-034	700403™
<i>Staphylococcus lugdunensis</i>	LRA 260.05.79	49576™
<i>Staphylococcus lugdunensis</i>	6733	700328™
<i>Staphylococcus pseudintermedius</i> deposited as <i>Staphylococcus aureus</i>		49444™
<i>Staphylococcus saprophyticus</i>	NCTC 7292	15305™
<i>Staphylococcus saprophyticus</i>	API 1101	35552™
<i>Staphylococcus saprophyticus</i>	LRA 27.02.80	43867™
<i>Staphylococcus saprophyticus</i>	API 222	49453™
<i>Staphylococcus saprophyticus</i>	AmMS 143	49907™
<i>Staphylococcus saprophyticus</i>	Vitek #8935	BAA-750™
<i>Staphylococcus sciuri</i> subsp. <i>sciuri</i>	GH9	29060™
<i>Staphylococcus sciuri</i> subsp. <i>sciuri</i>	GV252	29061™
<i>Staphylococcus simulans</i>	KH 1	27851™
<i>Staphylococcus</i> sp. deposited as <i>Staphylococcus epidermidis</i>	AMC 263	155™
<i>Staphylococcus</i> sp. deposited as <i>Staphylococcus epidermidis</i>	AMC 263	155-U™
<i>Staphylococcus warneri</i>	API 379	49454™
<i>Staphylococcus xylosus</i>	SL 8	29967™
<i>Staphylococcus xylosus</i>	KL 162	29971™
<i>Staphylococcus xylosus</i>	9280	35033™
<i>Staphylococcus xylosus</i>	LRA 1641 575	35663™
<i>Staphylococcus xylosus</i>	AmMS 151	49148™
<i>Staphylococcus xylosus</i>	API 85-12-228	700404™
<i>Stenotrophomonas maltophilia</i> deposited as <i>Pseudomonas maltophilia</i>	810-2	13637™
<i>Stenotrophomonas maltophilia</i> deposited as <i>Pseudomonas maltophilia</i>	300	17666™
<i>Stenotrophomonas maltophilia</i> deposited as <i>Pseudomonas maltophilia</i>	AmMS 194	49130™
<i>Stenotrophomonas maltophilia</i> deposited as <i>Xanthomonas maltophilia</i>	89-02-019	51331™
<i>Streptococcus agalactiae</i>	NCTC 8181	13813™
<i>Streptococcus agalactiae</i>	NADC 44	27956™
<i>Streptococcus agalactiae</i>	2603 V/R	BAA-611™
<i>Streptococcus agalactiae</i> deposited as <i>Streptococcus</i> sp.	grouping strain O90R	12386™
<i>Streptococcus agalactiae</i> deposited as <i>Streptococcus</i> sp.	typing strain D136C(3)	12403™
<i>Streptococcus anginosus</i> deposited as <i>Streptococcus anginosus</i>	NCTC 10713	33397™
<i>Streptococcus bovis</i>	NCDO 597	33317™
<i>Streptococcus bovis</i>	9145	35034™
<i>Streptococcus criceti</i> deposited as <i>Streptococcus cricetus</i>	HS-6	19642™
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i> deposited as <i>Streptococcus equisimilis</i>	LRA 06 11 76	35666™
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i> deposited as <i>Streptococcus equisimilis</i>	1180	9542™
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i> deposited as <i>Streptococcus</i> sp.	grouping strain C74	12388™
<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i> deposited as <i>Streptococcus</i> sp.	grouping strain D166B	12394™
<i>Streptococcus equi</i> subsp. <i>equi</i>	NCTC 9682	33398™
<i>Streptococcus equi</i> subsp. <i>equi</i>	2-1-23	9528™
<i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>	NCDO 1358	43079™
<i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>	API 77-01-036	700400™
<i>Streptococcus gallolyticus</i> deposited as <i>Streptococcus bovis</i>	AmMS 119	49147™
<i>Streptococcus gallolyticus</i> deposited as <i>Streptococcus bovis</i>	38	9809™
<i>Streptococcus gallolyticus</i> deposited as <i>Streptococcus bovis</i>	AmMS 236	49475™
<i>Streptococcus mitis</i> deposited as <i>Streptococcus viridans</i>		6249™
<i>Streptococcus mutans</i>	NCTC 10449	25175™
<i>Streptococcus mutans</i>	LRA 28 02 81	35668™
<i>Streptococcus oralis</i> deposited as <i>Streptococcus mitis</i>	M7A	9811™
<i>Streptococcus parasanguinis</i> deposited as <i>Streptococcus mitis</i>	[SS 895]	15909™
<i>Streptococcus pasteurianus</i> deposited as <i>Streptococcus bovis</i>	AmMS 207	49133™
<i>Streptococcus pneumoniae</i>	AmMS 208	49136™
<i>Streptococcus pneumoniae</i>	AmMS 149	49150™
<i>Streptococcus pneumoniae</i>	262	49619™
<i>Streptococcus pneumoniae</i>	CDC CS111	51916™
<i>Streptococcus pneumoniae</i>	Slovakia 14-10	700677™

Catalogue

Organism	Designation	ATCC® No.
<i>Streptococcus pneumoniae</i>	TIGR4	BAA-334™
<i>Streptococcus pneumoniae</i> deposited as <i>Diplococcus pneumoniae</i>	SVI	10015™
<i>Streptococcus pneumoniae</i> deposited as <i>Diplococcus pneumoniae</i>	R36a rough phase	27336™
<i>Streptococcus pneumoniae</i> deposited as <i>Diplococcus pneumoniae</i>		6301™
<i>Streptococcus pneumoniae</i> deposited as <i>Diplococcus pneumoniae</i>	[CIP 104225]	6303™
<i>Streptococcus pneumoniae</i> deposited as <i>Diplococcus pneumoniae</i>		6305™
<i>Streptococcus pyogenes</i>	typing strain T1	12344™
<i>Streptococcus pyogenes</i>	typing strain C203	12384™
<i>Streptococcus pyogenes</i>	C203 S	14289™
<i>Streptococcus pyogenes</i>	Bruno	19615™
<i>Streptococcus pyogenes</i>	QC A62	49399™
<i>Streptococcus pyogenes</i>	SF370; M1 GAS	700294™
<i>Streptococcus salivarius</i>	C699	13419™
<i>Streptococcus sanguinis</i>	SK36	BAA-1455™
<i>Streptococcus sanguinis</i> deposited as <i>Streptococcus sanguis</i>	DSS-10	10556™
<i>Streptococcus</i> sp.	grouping strain H60R	12392™
<i>Streptococcus</i> sp.	NCTC 4725	27284™
<i>Streptococcus</i> sp.	typing strain H36B	12401™
<i>Streptococcus thermophilus</i>	NCDO 573	19258™
<i>Streptococcus uberis</i>	API 78-11-025	700407™
<i>Streptococcus uberis</i>	[CIP 105801, NCDO 2055]	9927™
<i>Streptococcus vestibularis</i>	MM1	49124™
<i>Streptomyces griseus</i> subsp. <i>griseus</i>	4	10137™
<i>Streptomyces somaliensis</i>	IMRU 1274	33201™
<i>Tannerella forsythensis</i> deposited as <i>Bacteroides forsythus</i>	FDC 338	43037™
<i>Thermoanaerobacterium thermosaccharolyticum</i>	NCA 3814 (thermophilic anaerobe)	7956™
<i>Treponema phagedenis</i> deposited as <i>Treponema pallidum</i>	Kazan 8	27087™
<i>Trichomonas vaginalis</i>	C-1:NIH	30001™
<i>Trichomonas vaginalis</i>	123414	30184™
<i>Trichophyton mentagrophytes</i> deposited as <i>Trichophyton interdigitale</i>	640	9533™
<i>Trichophyton rubrum</i>	379	28188™
<i>Trichophyton tonsurans</i>		28942™
<i>Trichosporon mucoides</i>	API 85 02 177	201382™
<i>Trichosporon mucoides</i>	API 94 09 047	201383™
<i>Trichosporon mucoides</i>	Vitek 303483	204094™
<i>Ureaplasma urealyticum</i> deposited as <i>Ureaplasma urealyticum</i>		27619™
<i>Ureaplasma urealyticum</i> deposited as <i>Ureaplasma urealyticum</i>	7	27813™
<i>Veillonella parvula</i>	[ATCC 17742, Te 3]	10790™
<i>Vibrio alginolyticus</i> deposited as <i>Oceanomonas alginolytica</i>	XII-53	17749™
<i>Vibrio fischeri</i> deposited as <i>Achromobacter fischeri</i>	[NCMB 1281]	7744™
<i>Vibrio harveyi</i>	BB120	BAA-1116™
<i>Vibrio natriegens</i> deposited as <i>Pseudomonas natriegens</i>	[P Baumann 111]	14048™
<i>Vibrio parahaemolyticus</i>	EB 101	17802™
<i>Virgibacillus pantothenicus</i> deposited as <i>Bacillus pantothenicus</i>	NRS 1321	14576™
<i>Walleimia sebi</i>	FRR 1471	42694™
<i>Yarrowia lipolytica</i> deposited as <i>Mycoderma lipolytica</i>	251	9773™
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> deposited as <i>Bacterium enterocoliticum</i>	33114	9610™
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> deposited as <i>Yersinia enterocolitica</i>	Billups-1803-68	23715™
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> deposited as <i>Yersinia enterocolitica</i>	WA	27729™
<i>Yersinia kristensenii</i>	CDC 1459-81	33639™
<i>Yersinia ruckeri</i>	CDC 2396-61	29473™
<i>Zalerion arboricola</i>	MF5533	74030™
<i>Zygosaccharomyces rouxii</i> deposited as <i>Saccharomyces bisporus</i> var. <i>mellis</i>	NARD 3344	34890™
<i>Zygosaccharomyces rouxii</i> deposited as <i>Saccharomyces rouxii</i>	59-4	28253™

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Nucleic Acids prepared from ATCC Genuine Cultures®

For a complete listing see our online catalogue: www.lgcstandards-atcc.org

ATCC offers Nucleic Acids from well characterised and authenticated microbial strains, viruses and molecularly cloned viruses. Nucleic Acids from ATCC Genuine Cultures® fall into the following categories:

- Genomic DNA from microbial strains, namely bacteria, archaea, mycoplasma, protists, and fungi/yeast
- Viral DNA or RNA including viral genomic material, from infected cells or allantoic fluid
- Plasmid DNA from molecularly cloned viruses.

Nucleic Acids from ATCC Genuine Cultures can save you the time and expense of isolating DNA and RNA yourself.

Applications include:

- Construction of genomic libraries
- Southern hybridisation
- PCR
- Method development
- Positive controls for PCR/RT-PCR and other detection methods
- Identification and comparison at the genus, species, and strain levels.

This high quality DNA and RNA have been isolated under aseptic conditions to prevent cross-contamination. Batches are evaluated for integrity, purity and quality by several methods, including:


- Agarose gel electrophoresis
- Spectrophotometry
- Suitability for amplification by PCR
- Sequencing of the 16S ribosomal RNA gene (first ~ 500 base pairs) consistent with the sequence of the source organism for bacteria
- ITS sequencing consistent with the sequence of the source organism for fungi or yeast
- Sequence of PCR amplicon consistent with the sequence of the infecting agent for viruses.

Genomic DNA is available from bacteria and most yeasts in a 5 µg package size and from most filamentous fungi in a 2 µg package size. Most ATCC catalogue numbers have changed to add either '-5' or '-2' as a suffix. For example: ATCC catalogue number 10556D has been changed to 10556D-5 for the new 5 µg size. We've expanded our quality control testing for the new 5 µg and 2 µg sizes to include the determination of the total amount of nucleic acids by PicoGreen® measurement. Original packaging as indicated by a 'D' suffix was tested by an alternate method.

The package size of DNA and RNA prepared from Mycoplasma remains unchanged at either 10 ng or 50 ng.

The package size of viral DNA is 100 µl per vial, dilutable ten-fold or more for amplification. All viral nucleic acids have been tested to determine the total amount of the appropriate nucleic acid by PicoGreen® or RiboGreen® measurement, and have been tested for viral inactivation.

For more information, please refer to the individual product sheets.



Fully Sequenced Microbes

Bacteria and Archea

Name and strain	Culture ATCC® No.	Genomic DNA ATCC® No.	GenBank® No.
<i>Acidithiobacillus ferrooxidans</i> ATCC 23270™	23270™	23270D	CP001219
<i>Acidithiobacillus ferrooxidans</i> BNL-5-31	53993™	53993D-5	CP001132
<i>Acidothermus cellulolyticus</i> 11B	43068™	43068D-5	CP000481
<i>Acidovorax</i> sp. JS42	700441™		CP000539
<i>Acinetobacter baumannii</i> 5377	17978™	17978D-5	CP000521
<i>Actinobacillus succinogenes</i> 130Z	55618™	55618D-5	CP000746
<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i> RH 250	7966™	7966D-5	CP000462
<i>Aeropyrum pernix</i> K1	700893™		BA000002
<i>Akkermansia muciniphila</i> Muc	BAA-835™	BAA-835D-5	CP001071
<i>Anabaena variabilis</i> IUCC 1444	29413™	29413D-5	CP000117
<i>Anaeromyxobacter dehalogenans</i> 2CP-1	BAA-258™		CP001359
<i>Anaeromyxobacter dehalogenans</i> 2CP-C	BAA-259™	BAA-259D-5	CP000251
<i>Archaeoglobus fulgidus</i> VC16	49558™	49558D-5	AE000782
<i>Arthrobacter aurescens</i> TC1	BAA-1386™	BAA-1386D-5	CP000474
<i>Arthrobacter chlorophenolicus</i> A6	700700™	700700D-5	CP001341
<i>Bacillus cereus</i> ATCC 10987™	10987™	10987D-5	AE017194
<i>Bacillus cereus</i> ATCC 14579™	14579™	14579D-5	AE016877
<i>Bacillus halodurans</i> C-125	BAA-125™	BAA-125D-5	BA000004
<i>Bacillus licheniformis</i> 46	14580™	14580D-5	CP000002
<i>Bacteroides fragilis</i> EN-2	25285™	25285D	CR626927
<i>Bacteroides thetaiotaomicron</i> E50	29148™	29148D-5	AE015928
<i>Bacteroides vulgatus</i> ATCC 8482™	8482™	8482D-5	CP000139
<i>Bartonella bacilliformis</i> KC583	35685™	35685D-5	CP000524
<i>Bartonella henselae</i> Houston-1	49882™	49882D-5	BX897699
<i>Bdellovibrio bacteriovorus</i> HD100	15356™	1536D-5	BX842601
<i>Beijerinckia indica</i> subsp. <i>indica</i> 11	9039™	9039D-5	CP001016
<i>Bifidobacterium adolescentis</i> E194a	15703™	15703D	AP009256
<i>Bifidobacterium infantis</i> S12	15697™	15697D	CP001095
<i>Bordetella bronchiseptica</i> RB50	BAA-588™	BAA-588D-5	BX470250
<i>Bordetella parapertussis</i> 12822	BAA-587™	BAA-587D-5	BX470249
<i>Bordetella pertussis</i> Tohama I	BAA-589™	BAA-589D-5	BX470248
<i>Borrelia burgdorferi</i> B31	35210™	35210D-5	AE001115
<i>Bradyrhizobium</i> sp. BTAi1	BAA-1182™		CP000494
<i>Brucella ovis</i> 63/290	25840™	25840D-5	CP000709
<i>Burkholderia ambifaria</i> AMMD	BAA-244™	BAA-244D-5	CP000440
<i>Burkholderia cepacia</i> 249	17616™	17616D-5	CP000869
<i>Burkholderia thailandensis</i> E264	700388™	700388D-5	CP000086
<i>Caldicellulosiruptor saccharolyticus</i> Tp8T.6331	43494™	43494D-5	CP000679
<i>Campylobacter concisus</i> 13826	BAA-1457™	BAA-1457D-5	CP000792
<i>Campylobacter hominis</i> CH001A	BAA-381™	BAA-381D-5	CP000776
<i>Campylobacter jejuni</i> RM1221	BAA-1062™	BAA-1062D-5	CP000025
<i>Campylobacter jejuni</i> subsp. <i>doylei</i> 269.97	BAA-1458™	BAA-1458D-5	CP000768
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> NCTC 11168	700819™	700819D-5	AL111168
<i>Campylobacter lari</i> RM2100	BAA-1060™	BAA-1060D-5	CP000932
<i>Caulobacter vibrioides</i> CB 15	19089™	19089D-5	AE005673
<i>Chlorobium tepidum</i> TLS	49652™	49652D-5	AE006470
<i>Chloroflexus aurantiacus</i> J-10-fl	29366™	29366D-5	CP000909
<i>Chloroflexus aurantiacus</i> Y-400-fl	29364™		CP001364
<i>Chloroherpeton thalassium</i> GB-78	35110™		CP001100
<i>Chromobacterium violaceum</i> MK	12472™	12472D	AE016825
<i>Chromohalobacter salexigens</i> 1H11	BAA-138™	BAA-138D-5	CP000285
<i>Citrobacter koseri</i> 4225-83	BAA-895™	BAA-895D-5	CP000822
<i>Clostridium acetobutylicum</i> ATCC 824™	824™	824D-5	AE001437
<i>Clostridium difficile</i> 630	BAA-1382™	BAA-1382D-5	AM180355
<i>Clostridium perfringens</i> ATCC 13124™	13124™	13124D-5	CP000246
<i>Clostridium phytofermentans</i> ISDg	700394™	700394D-5	CP000885
<i>Clostridium thermocellum</i> 157	27405™	27405D-5	CP000568
<i>Colwellia psychrerythraea</i> 34H	BAA-681™	BAA-681D	CP000083
<i>Corynebacterium diphtheriae</i> ATCC 700971™	700971™	700971D-5	BX248353
<i>Corynebacterium glutamicum</i> 534	13032™	13032D-5	BA000036
<i>Cupriavidus metallidurans</i> CH34	43123™	43123D-5	CP000352
<i>Cupriavidus necator</i> H16	17699™	17699D-5	AM260479
<i>Cyanothece</i> sp. BH68	51142™	51142D-5	CP000806
<i>Cyanothece</i> sp. PCC 7425	29141™		CP001344
<i>Cytophaga hutchinsonii</i> ATCC 33406™	33406™	33406D-5	CP000383
<i>Deinococcus radiodurans</i> R1	BAA-816™	BAA-816D	AE000513
<i>Desulfovibrio desulfuricans</i> subsp. <i>desulfuricans</i> MB	27774™	27774D-5	CP001358
<i>Desulfovibrio vulgaris</i> subsp. <i>vulgaris</i> Hildenborough	29579™	29579D-5	AE017285
<i>Dictyoglomus thermophilum</i> H-6-12	35947™		CP001146
<i>Ensifer meliloti</i> Rm1021	51124™	51124D-5	AL591688

Fully Sequenced Microbes

Bacteria and Archea

Name and strain	Culture ATCC® No.	Genomic DNA ATCC® No.	GenBank® No.
<i>Enterobacter sakazakii</i> 2001-10-01	BAA-894™	BAA-894D-5	CP000783
<i>Enterococcus faecalis</i> V583	700802™	700802D-5	AE016830
<i>Escherichia coli</i> Crooks	8739™	8739D-5	CP000946
<i>Escherichia coli</i> K12-MG1655	700926™	700926D-5	U00096
<i>Escherichia coli</i> UPEC-CFT073	700928™	700928D-5	AE014075
<i>Escherichia coli</i> RIMD 0509952	BAA-460™	BAA-460D-5	BA000007
<i>Escherichia fergusonii</i> CDC 0568-73	35469™	35469D-5	CU928158
<i>Fingoldia magna</i> WAL2508	29328™	29328D-5	AP008971
<i>Flavobacterium johnsoniae</i> MYX.1.1.1	17061™	17061D-5	CP000685
<i>Flavobacterium psychrophilum</i> TG 02/86	49511™	49511D-5	AM398681
<i>Francisella philomiragia</i> O#319-036	25017™	25017D-5	CP000937
<i>Fusobacterium nucleatum</i> subsp. <i>nucleatum</i> 1612A	25586™	25586D-5	AE009951
<i>Geobacter lovleyi</i> SZ	BAA-1151™		CP001089
<i>Geobacter metallireducens</i> GS-15	53774™	53774D-5	CP000148
<i>Geobacter sulfurreducens</i> PCA	51573™	51573D-5	AE017180
<i>Gloeobacter violaceus</i> PCC 7421	29082™	29082D-5	BA000045
<i>Gluconacetobacter diazotrophicus</i> PAI 5	49037™	49037D-5	AM889285
<i>Gluconobacter oxydans</i> ATCC 621H™	621H™	621HD-5	CP000009
<i>Granulobacter thebesensis</i> CGDNIH1	BAA-1260™	BAA-1260D-5	CP000394
<i>Haemophilus ducreyi</i> 35000HP	700724™	700724D-5	AE017143
<i>Haemophilus influenzae</i> KW20 Rd	51907™	51907D	L42023
<i>Haloarcula marismortui</i> ATCC 43049™	43049™	43049D-5	AY596297
<i>Halobacterium salinarum</i> NRC-1	700922™	700922D	AE004437
<i>Halorubrum lacusprofundi</i> ACAM34	49239™		CP001365
<i>Helicobacter pylori</i> 26695	700392™	700392D-5	AE000511
<i>Helicobacter pylori</i> J99	700824™	700824D-5	AE001439
<i>Heliobacterium modesticaldum</i> Ice1	51547™	51547D-5	CP000930
<i>Herpetosiphon aurantiacus</i> 114-95	23779™		CP000875
<i>Hyphomonas neptunium</i> 14-15	15444™	15444D-5	CP000158
<i>Idiomarina loihiensis</i> L2-TR	BAA-735™	BAA-735D-5	AE017340
<i>Kineococcus radiotolerans</i> SRS30216	BAA-149™	BAA-149D-5	CP000750
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> MGH 78578	700721™	700721D-5	CP000647
<i>Kocuria rhizophila</i> DC2201	9341™	9341D-5	AP009152
<i>Lactobacillus brevis</i> 118-8	367™	367D-5	CP000416
<i>Lactobacillus casei</i> ATCC 334™	334™	334D-5	CP000423
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> Lb14	11842™	11842D-5	CR954253
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> Lb-18	BAA-365™	BAA-365D	CP000412
<i>Lactobacillus gasseri</i> 63 AM	33323™	33323D-5	CP000413
<i>Lactobacillus plantarum</i> WCFS1	BAA-793™	BAA-793D	AL935263
<i>Lactobacillus reuteri</i> F275	23272™	23272D-5	CP000705
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i> Philadelphia 1	33152™	33152D-5	AE017354
<i>Leptospira interrogans</i> serovar <i>Copenhageni</i> , Fiocruz L1-130	BAA-1198™	BAA-1198D-5	AE016823
<i>Leuconostoc mesenteroides</i> subsp. <i>mesenteroides</i> 37Y	8293™	8293D-5	CP000414
<i>Listeria innocua</i> CLIP 11262	BAA-680™	BAA-680D	AL592022
<i>Listeria monocytogenes</i> EGDe	BAA-679™	BAA-679D-5	AL591824
<i>Listeria welshimeri</i> serovar 6b, SLCC5334	35897™	35897D-5	AM263198
<i>Magnetospirillum magneticum</i> AMB-1	700264™	700264D-5	CP007255
<i>Marinobacter aquaeolei</i> VT8	700491™	700491D-5	CP000514
<i>Mesoplasma florum</i> L1	33453™	33453D-5	AE017263
<i>Metallosphaera sedula</i> TH2	51363™	51363D-5	CP000682
<i>Methanobrevibacter smithii</i> P5	35061™	35061D-5	CP000678
<i>Methanocaldococcus jannaschii</i> JAL-1	43067™	43067D-5	L77117
<i>Methanococcus aeolicus</i> Nankai-3	BAA-1280™	BAA-1280D-5	CP000743
<i>Methanococcus maripaludis</i> C5	BAA-1333™	BAA-1333D-5	CP000609
<i>Methanococcus maripaludis</i> C6	BAA-1332™	BAA-1332D-5	CP000867
<i>Methanococcus maripaludis</i> C7	BAA-1331™	BAA-1331D-5	CP000745
<i>Methanococcus vannielii</i> SB	35089™	35089D-5	CP000742
<i>Methanocorpusculum labreanum</i> Z	43576™	43576D-5	CP000559
<i>Methanoculleus marisnigri</i> JR1	35101™	35101D-5	CP000562
<i>Methanosarcina acetivorans</i> C2A	35395™	35395D-5	AE010299
<i>Methanosarcina mazei</i> Go1	BAA-159™	BAA-159D-5	AE008384
<i>Methanospirillum hungatei</i> JF-1	27890™	27890D-5	CP000254
<i>Methanothermobacter thermautotrophicus</i> delta H	29096™	29096D-5	AE000666
<i>Methylobium petroleiphilum</i> PM1	BAA-1232™	BAA-1232D-5	CP000555
<i>Methylobacillus flagellatus</i> KT	51484™	51484D	CP000284
<i>Methylobacterium radiotolerans</i> O-1	27329™	27329D-5	CP001001
<i>Methylococcus capsulatus</i> Bath	33009™	33009D-5	AE017282
<i>Moorella thermoacetica</i> ATCC 39073™	39073™	39073D-5	CP000232
<i>Mycobacterium abscessus</i> L948	19977™		CU458896
<i>Mycobacterium avium</i> subsp. <i>paratuberculosis</i> K-10	BAA-968™	BAA-968D-5	AE016958

Fully Sequenced Microbes

Bacteria and Archea

Name and strain	Culture ATCC® No.	Genomic DNA ATCC® No.	GenBank® No.
<i>Mycobacterium bovis</i> AF2122/97	BAA-935™		BX248333
<i>Mycobacterium marinum</i> M	BAA-535™	BAA-535D-5	CP000854
<i>Mycobacterium smegmatis</i> MC2 155	700084™	700084D-5	CP000480
<i>Mycobacterium tuberculosis</i> H37Ra	25177™	25177D-5	CP000611
<i>Mycobacterium tuberculosis</i> H37Rv	25618™		AL123456
<i>Mycoplasma genitalium</i> G-37	33530™	33530D	L43967
<i>Mycoplasma hyopneumoniae</i> J	25934™	25934D	AE017243
<i>Mycoplasma mobile</i> 163K	43663™	43663D	AE017308
<i>Mycoplasma pneumoniae</i> M129	29342™	29342D	U00089
<i>Natronomonas pharaonis</i> Gabara	35678™	35678D-5	CR936257
<i>Nautilia profundicola</i> AmH	BAA-1463™		CP001279
<i>Neisseria gonorrhoeae</i> FA 1090	700825™	700825D-5	AE004969
<i>Neisseria meningitidis</i> FAM18	700532™	700532D-5	AM421808
<i>Neisseria meningitidis</i> MC58	BAA-335™	BAA-335D-5	AE002098
<i>Nitrobacter winogradskyi</i> Nb-255	25391™		CP000115
<i>Nitrosococcus oceani</i> C-107	19707™	19707D-5	CP000127
<i>Nitrosomonas europaea</i> ATCC 19718™	19718™	19718D	AL954747
<i>Nitrospira multififormis</i> C 71	25196™	25196D-5	CP000103
<i>Nocardioides</i> sp. JS614	BAA-499™	BAA-499D-5	CP000509
<i>Nostoc punctiforme</i> PCC 73102	29133™	29133D	CP001037
<i>Nostoc</i> sp. PCC 7120	27893™	27893D-5	BA000019
<i>Novosphingobium aromaticivorans</i> SMCC F199	700278™	700278D-5	CP000248
<i>Ochrobactrum anthropi</i> CL350/82	49188™	49188D-5	CP000758
<i>Oenococcus oeni</i> PSU-1	BAA-331™	BAA-331D-5	CP000411
<i>Oligotropha carboxidovorans</i> OM5	49405™		CP001196
<i>Parabacteroides distasonis</i> ATCC 8503™	8503™	8503D-5	CP000140
<i>Pectobacterium atrosepticum</i> SCRI 1043	BAA-672™	BAA-672D	BX950851
<i>Pediococcus pentosaceus</i> 183-1w	25745™	25745D-5	CP000422
<i>Photobacterium profundum</i> SS9	BAA-1253™	BAA-1253D-5	CR354531
<i>Polaromonas naphthalenivorans</i> CJ2	BAA-779™		CP000529
<i>Polaromonas</i> sp. JS666	BAA-500™	BAA-500D-5	CP000316
<i>Porphyromonas gingivalis</i> 2561	33277™	33277D-5	AP009380
<i>Porphyromonas gingivalis</i> W83	BAA-308™	BAA-308D-5	AE015924
<i>Pseudoalteromonas atlantica</i> T6c	BAA-1087™	BAA-1087D-5	CP000388
<i>Pseudomonas fluorescens</i> Pf-5	BAA-477™	BAA-477D-5	CP000076
<i>Pseudomonas putida</i> F1	700007™	700007D-5	CP000712
<i>Pseudomonas putida</i> KT2440	47054™	47054D-5	AE015451
<i>Pseudomonas syringae</i> pathovar <i>Phaseolicola</i> , 1448A	BAA-978™	BAA-978D	CP000058
<i>Pseudomonas syringae</i> pathovar <i>Tomato</i> , DC3000	BAA-871™	BAA-871D-5	AE016853
<i>Psychrobacter cryohalolentis</i> K5	BAA-1226™	BAA-1226D-5	CP000323
<i>Pyrococcus furiosus</i> Vc1	43587™	43587D-5	AE009950
<i>Pyrococcus horikoshii</i> OT3	700860™	700860D-5	BA000001
<i>Ralstonia solanacearum</i> GMI1000	BAA-1114™	BAA-1114D-5	AL646052
<i>Renibacterium salmoninarum</i> Lea-1-74	33209™	33209D-5	CP000910
<i>Rhizobium etli</i> CFN 42	51251™	51251D-5	CP000133
<i>Rhizobium radiobacter</i> C58	33970™	33970D	AE008688
<i>Rhizobium radiobacter</i> K84	49644™		CP000628
<i>Rhizobium vitis</i> S4	BAA-846™	BAA-846D	CP000633
<i>Rhodobacter sphaeroides</i> 2.4.1	BAA-808™	BAA-808D	CP000143
<i>Rhodobacter sphaeroides</i> ATH 2.4.9	17029™	17029D-5	CP000577
<i>Rhodoferrax ferrireducens</i> T118	BAA-621™	BAA-621D-5	CP000267
<i>Rhodopseudomonas palustris</i> BisA53	BAA-1125™	BAA-1125D-5	CP000463
<i>Rhodopseudomonas palustris</i> BisB18	BAA-1124™	BAA-1124D-5	CP000301
<i>Rhodopseudomonas palustris</i> BisB5	BAA-1123™	BAA-1123D-5	CP000283
<i>Rhodopseudomonas palustris</i> CGA009	BAA-98™	BAA-98D-5	BX571963
<i>Rhodopseudomonas palustris</i> HaA2	BAA-1122™	BAA-1122D-5	CP000250
<i>Rhodospirillum centenum</i> SW	51521™		CP000613
<i>Rhodospirillum rubrum</i> ATCC 11170™	11170™	11170D-5	CP000230
<i>Roseobacter denitrificans</i> OCh 114	33942™	33942D-5	CP000362
<i>Saccharophagus degradans</i> 2-40	43961™	43961D-5	CP000282
<i>Saccharopolyspora erythraea</i> M5-12259	11635™	11635D-5	AM420293
<i>Salinibacter ruber</i> M31	BAA-605™	BAA-605D	CP000159
<i>Salinispora tropica</i> CNB-440	BAA-916™	BAA-916D-5	CP000667
<i>Salmonella enterica</i> subsp. <i>arizonae</i> RSK2980	BAA-731™	BAA-731D-5	CP000880
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Paratyphi</i> A, ATCC 9150™	9150™	9150D-5	CP000026
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Paratyphi</i> B, SPB7	BAA-1250™	BAA-1250D-5	CP000886
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> , LT2	700720™	700720D-5	AE006468
<i>Shewanella amazonensis</i> SB2B	700329™	700329D-5	CP000507
<i>Shewanella baltica</i> OS155	BAA-1091™	BAA-1091D-5	CP000563
<i>Shewanella denitrificans</i> OS217	BAA-1090™	BAA-1090D-5	CP000302

Fully Sequenced Microbes

Bacteria and Archea

Name and strain	Culture ATCC® No.	Genomic DNA ATCC® No.	GenBank® No.
<i>Shewanella frigidimarina</i> NCIMB 400	BAA-1089™	BAA-1089D-5	CP000447
<i>Shewanella oneidensis</i> MR-1	700550™	700550D	AE014299
<i>Shewanella pealeana</i> ANG-SQ1	700345™	700345D-5	CP000851
<i>Shewanella putrefaciens</i> CN-32	BAA-453™	BAA-453D-5	CP000681
<i>Shewanella woodyi</i> MS32	51908™	51908D-5	CP000961
<i>Shigella flexneri</i> 2457T	700930™	700930D-5	AE014073
<i>Silicibacter pomeroyi</i> DSS-3	700808™	700808D	CP000031
<i>Staphylococcus aureus</i> MW2 USA 400	BAA-1707™		BA000033
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Mu3	700698™	700698D-5	AP009324
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Mu50	700699™	700699D-5	BA000017
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> USA 300, FPR3757	BAA-1556™	BAA-1556D-5	CP000255
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™	12228D-5	AE015929
<i>Staphylococcus epidermidis</i> RP62A	35984™	35984D-5	CP000029
<i>Staphylococcus saprophyticus</i> S-41	15305™	15305D-5	AP008934
<i>Staphylothermus marinus</i> F1	43588™	43588D-5	CP000575
<i>Streptococcus agalactiae</i> 2603V/R	BAA-611™	BAA-611D-5	AE009948
<i>Streptococcus agalactiae</i> A909	BAA-1138™	BAA-1138D-5	CP000114
<i>Streptococcus gordonii</i> Challis, CH1	35105™	35105D-5	CP000725
<i>Streptococcus mutans</i> UA159	700610™	700610D-5	AE014133
<i>Streptococcus pneumoniae</i> Mu50	700669™	700669D-5	FM211187
<i>Streptococcus pneumoniae</i> R6	BAA-255™	BAA-255D-5	AE007317
<i>Streptococcus pneumoniae</i> TIGR4	BAA-334™	BAA-334D-5	AE005672
<i>Streptococcus pyogenes</i> M1 GAS, SF370	700294™	700294D-5	AE004092
<i>Streptococcus pyogenes</i> MGAS 8232	BAA-572™	BAA-572D-5	AE009949
<i>Streptococcus pyogenes</i> MGAS 315	BAA-595™	BAA-595D-5	AE014074
<i>Streptococcus pyogenes</i> MGAS 10394	BAA-946™	BAA-946D-5	CP000003
<i>Streptococcus pyogenes</i> MGAS 5005	BAA-947™	BAA-947D-5	CP000017
<i>Streptococcus pyogenes</i> MGAS 10270	BAA-1063™	BAA-1063D-5	CP000260
<i>Streptococcus pyogenes</i> MGAS 6180	BAA-1064™	BAA-1064D-5	CP000056
<i>Streptococcus pyogenes</i> MGAS 2096	BAA-1065™	BAA-1065D-5	CP000261
<i>Streptococcus pyogenes</i> MGAS 10750	BAA-1066™	BAA-1066D-5	CP000262
<i>Streptococcus pyogenes</i> MGAS 9429	BAA-1315™	BAA-1315D-5	CP000259
<i>Streptococcus sanguinis</i> SK36	BAA-1455™	BAA-1455D-5	CP000387
<i>Streptococcus thermophilus</i> LMD-9	BAA-491™	BAA-491D-5	CP000419
<i>Streptococcus thermophilus</i> LMG 18311	BAA-250™	BAA-250D-5	CP000023
<i>Streptococcus uberis</i> O140J	BAA-854™	BAA-854D-5	AM946015
<i>Streptomyces avermitilis</i> MA-4680	31267™	31267D-5	BA000030
<i>Streptomyces violaceoruber</i> M145	BAA-471™	BAA-471D-5	AL645882
<i>Sulfolobus acidocaldarius</i> 98-3	33909™	33909D-5	CP000077
<i>Sulfolobus solfataricus</i> P2	35092™	35092D-5	AE006641
<i>Sulfurimonas denitrificans</i> ATCC 33889™	33889™	33889D-5	CP000153
<i>Synechococcus elongatus</i> PCC 7942	33912™	33912D-5	CP000100
<i>Synechococcus</i> sp. PR-6	27264™	27264D-5	CP000951
<i>Synechocystis</i> sp. N-1	27184™	27184D-5	BA000022
<i>Thermoanaerobacter pseudethanolicus</i> 39E	33223™	33223D-5	CP000924
<i>Thermoanaerobacter</i> sp. X514	BAA-938™	BAA-938D-5	CP000923
<i>Thermobifida fusca</i> YX	BAA-629™	BAA-629D-5	CP000088
<i>Thermodesulfovibrio yellowstonii</i> YP87	51303™		CP001147
<i>Thermomicrobium roseum</i> P-2	27502™		CP001275
<i>Thermoplasma acidophilum</i> 122-1B2	25905™	25905D-5	AL139299
<i>Thermoplasma volcanium</i> GSS1	51530™	51530D-5	BA000011
<i>Thermotoga lettingae</i> TMO	BAA-301™	BAA-301D-5	CP000812
<i>Thermotoga maritima</i> MSB8	43589™	43589D-5	AE000512
<i>Thermotoga neapolitana</i> NS-E	49049™	49049D	CP000916
<i>Thermotoga petrophila</i> RKU-1	BAA-488™	BAA-488D-5	CP000702
<i>Thermus thermophilus</i> HB27	BAA-163™	BAA-163D-5	AE017221
<i>Thermus thermophilus</i> HB8	27634™	27634D-5	AP008226
<i>Thiobacillus denitrificans</i> T1	25259™	25259D-5	CP000116
<i>Treponema denticola</i> a	35405™	35405D-5	AE017226
<i>Ureaplasma parvum</i> ATCC 700970™	700970™		AF222894
<i>Ureaplasma parvum</i> serotype III, 27	27815™	27815D	CP000942
<i>Ureaplasma urealyticum</i> Western	33699™		CP001184
<i>Verminephrobacter eiseniae</i> EF01-2	BAA-1489™	BAA-1489D-5	CP000542
<i>Vibrio fischeri</i> ES114	700601™	700601D-5	CP000020
<i>Vibrio Harveyi</i> BB120	BAA-1116™	BAA-1116D-5	CP000789
<i>Wolinella succinogenes</i> FDC 602W	29543™	29543D-5	BX571656
<i>Xanthobacter autotrophicus</i> Py2	BAA-1158™	BAA-1158D-5	CP000781
<i>Xanthomonas campestris</i> pathovar <i>Campestris</i> , ATCC 33913™	33913™	33913D-5	AE008922
<i>Xylella fastidiosa</i> Temecula	700964™	700964D-5	AE009442
<i>Zymomonas mobilis</i> ZM4	31821™	31821D-5	AE008692

Fully Sequenced Microbes

Eukaryotes

Name and strain - Fungi	Culture ATCC® No.	Genomic DNA ATCC® No.
<i>Ashbya gossypii</i> NRRL Y-1056	10895™	10895D-5
<i>Candida glabrata</i> CBS 138	2001™	2001D-5
<i>Cryptococcus neoformans</i> JEC21	MYA-565™	MYA-565D-5
<i>Debaryomyces hansenii</i> CBS 767	36239™	36239D-5
<i>Kluyveromyces lactis</i> 61	8585™	8585D-5
<i>Kluyveromyces waltii</i> UCD 72-13	56500™	56500D-5
<i>Pichia stipitis</i> CBS 6054	58785™	58785D-2
<i>Saccharomyces cerevisiae</i> S288C	204508™	204508D-5
<i>Saccharomyces mikatae</i> CBS 8839	MYA-4448™	MYA-4448D-5
<i>Schizosaccharomyces pombe</i> 972 h-	24843™	24843D-5
<i>Schizosaccharomyces pombe</i> 972H-	26189™	26189D-5
<i>Vanderwaltozyma polyspora</i> CBS 2163	22028™	22028D-5
Name and strain - Protozoa	Culture ATCC® No.	Genomic DNA ATCC® No.
<i>Cryptosporidium parvum</i> Iowa		PRA-67D
<i>Entamoeba histolytica</i> HM1:IMSS	30459™	

A complete listing of ATCC Genuine Cultures® can be found at: www.lgcstandards-atcc.org

Chill Out!

Ready-to-use Microbial Freeze Media



Safely cryopreserve non-fastidious bacterial
and transformed hosts!

ATCC offers two long-term microbial freeze media products. Just add actively growing cells to the freeze medium, mix well, and store at -80°C . When needed, thaw the suspension and inoculate fresh culture medium.

- ✓ Ready to use
- ✓ Packaged in 12 convenient single-use cryo-vials
- ✓ Space on the label to add culture information

TSB with 10% glycerol 12 x 1.0 ml.....Catalogue no. 20-2200

Tryptic soy broth with glycerol is used for cryopreservation of non-fastidious bacteria.

LB with 30% glycerol 12 x 0.5 ml.....Catalogue no. 60-2200

Luria-Bertani broth (Miller formulation) with glycerol is used for cryopreservation of bacteria transformed with plasmids.

Discounts are available for large orders. Please enquire:

atcc@lgcstandards.com



in partnership with

Excellence through measurement

Microbial technical resource



Technical resource

Preparation of biological materials

How materials are processed in preparation for freezing can have an affect on the outcome of the preservation process. For non-replicable materials such as tissues, nucleic acids and proteins, the preparation process consists of ensuring that the materials are in the proper solution or freezing medium in order to maximize the intended use of the materials when recovered. However the stability and recoverability of living cells and organisms is affected by the growth conditions and pre-freezing processing. Several factors must be considered when preparing cells for cryopreservation. These include the type of cell, cell viability, growth conditions, physiological state of the cells, the number of cells, and how the cells are handled.

When preparing the initial seed stock of a new isolate or cell line, the culture should be examined for identity and contaminating microorganisms at a minimum. This examination should be repeated after preservation and each time a new lot of the culture is prepared.

Microorganisms

Microbial cells, particularly bacteria and yeast, grown under aerated conditions demonstrate a greater resistance to the detrimental effects of cooling and freezing than non-aerated cells.¹ T. Nei et.al.² have demonstrated that cell permeability is greater in aerated cultures, and that the aerated cells dehydrate faster during cooling than non-aerated cells. Microbial cells harvested from late log or early stationary cultures also demonstrate greater resistance to the freezing process than younger or older cells.¹

Generally, the greater the number of cells present initially, the greater the recovery. For most bacteria and yeast, approximately 10^7 /ml cells are required to ensure adequate recovery.³ These can be conveniently harvested from agar slants or plates, or when greater quantities are required, grown in broth culture and harvested by centrifugation. In either case, cells are generally suspended in fresh growth medium containing the cryoprotective agent. Protists can also be concentrated by centrifugation, but are often suspended in the used medium and then diluted by adding an equal volume of fresh growth medium containing the cryoprotective agent.³ Spore forming fungi require harvesting of spores and suspension of the spores in fresh growth medium containing the cryoprotective agent. When freezing fungal spores, care must be taken not to delay the freezing process too long to ensure that germination does not occur prior to freezing. For fungi that do not form spores, special procedures for harvesting mycelia prior to freezing must be utilised. For fungi with tough mycelia, the culture is harvested from agar growth by cutting and removing agar plugs containing the mycelia and placing the plugs into fresh growth medium containing the cryoprotective agent. Tough mycelia that do not adhere well to agar cultures are grown in broth culture and the mycelial mass is blended prior to freezing.³

The viability and an estimate of recovery should be determined both before and after freezing the culture. Viability is a measure of the culture's ability to grow and reproduce. For some material, such as protozoan cultures, this should include several passages to ensure stability. An estimate of the number of cells recovered can be made by several means including serial dilution, plate counts, or direct cell counting. A comparison of the counts prior to and after freezing gives an indication of the degree of recovery or the success of the preservation procedure.

Genetically modified materials

Genetically modified cells and organisms can be cryopreserved in a manner similar to the unmodified host cell.^{4,5}

Inventory control

Appropriate record keeping is important in any laboratory and there are a number of methods available for keeping records on cryopreserved materials.⁶ When establishing your own method, keep in mind that there is key information which will be important for future use: (a) the preservation methodology used; (b) the location and identification of the stored material; (c) preservation date; and (d) number of passages for replicable material. The item number should be linked to associated data for that material, and for some purposes each container may require a unique identifying code linked to specific information for that particular aliquot.

Identification begins with proper labeling of the storage container. The label information should include a name or identification code for the frozen material, as well as a lot number. The information on the label should be kept with the inventory records that include the location code for each vial. These records can be maintained as paper documents, or preferably as electronic files. Duplicate inventory records should be maintained in a location separate from working records. Locator codes should be specific enough to allow rapid and easy retrieval for a specific lot and should include freezer unit number, a code for a freezer section or inventory rack, a box or canister number, and possibly even a grid spot within the box or a cane number when canes are used. Detailed locator codes minimise hunting for material which risks warming the freezer unit, exposure of other materials to warmer temperatures, and prolonged exposure of laboratory personnel to extremely cold temperatures.

References

1. Simone, F.P. 1992. Key issues relating to the genetic stability and preservation of cells and cell banks. *J. Parent. Science and Technology* 46: 226-232.
2. Nei, T., T. Araki and T. Matsusaka. 1969. Freezing injury to aerated and non-aerated cultures of *Escherichia coli*. In T. Nei, Ed. *Freezing and Drying of Microorganisms*. University of Tokyo Press, Tokyo Japan.
3. Ed. F.P. Simone and E.M. Brown. 1991. *ATCC Preservation Methods: Freezing and Freeze Drying*. American Type Culture Collection, Manassas, VA.
4. Nierman, W.C. and T. Feldblyum. 1985. Cryopreservation of cultures that contain plasmids. *Dev. Ind. Microbiol.* 26: 423-434.
5. Nierman, W.C., C. Trypus and L.L. Deaven. 1987. Preservation and stability of bacteriophage lambda libraries by freezing in liquid nitrogen. *Biotechniques* 5: 724-727.
6. Simone, F.P. Cryopreservation: Storage and Documentation Systems, In: *Biotechnology: Quality Assurance and Validation, Drug Manufacturing Technologies Series, Vol. 4*, Interpharm Press, Buffalo Grove, Illinois, 1999, pgs 7-31.

Come directly to the source

Useful terms for ATCC Microbiology

Organisations and agencies such as AOAC, ASTM, FDA, CLSI (formerly NCCLS) and USP specify ATCC cultures in their standards and specifications, then describe the propagation and growth conditions to be used. As the collection responsible for these cultures, ATCC uses the following working definitions derived from the Bacteriological Code,¹ Brown and Gilbert's Microbiological Quality Assurance² and NCCLS (now CLSI) M22-A2.³

Authentication – the process of comparing the characteristics of a strain to the description of the species to determine or verify the strain's identity.

Characterisation – the process of subjecting a strain to a battery of morphological, physiological, molecular and/or other tests. The data derived can be used to describe the strain and objectively compare it to other strains.

Culture – a population of cells in a given place at a given time, e.g., in a test tube or on an agar plate.

Distribution stock – cultures derived from seed stock and intended for distribution. At least one vial of each batch of distribution stock is checked for purity and viability, then characterised and compared with the data for the seed stock.

Master culture – a culture derived directly from the reference culture (i.e., the vial obtained from ATCC or another recognised culture collection). The master culture can be stored as freeze-dried ampoules, low temperature suspensions, low temperature suspension on beads, or through suitably validated alternatives. The master culture serves as the in-house reference standard. The master culture is used to prepare stock cultures.

Passage – transfer of an inoculum of cells from an existing culture to fresh growth medium in another vessel. Each transfer (i.e., subculture) is counted as an additional passage. Some organisations have recommended, and ATCC concurs, that strains being used in standard test methods should not be used beyond five

passages from the original ATCC vial. ATCC defines the first passage as the first broth or agar culture started from the culture vial supplied by ATCC. Actively growing cultures are subject to genetic drift, accidental loss, contamination at time of subculture, human error in transfer and labelling, and selection of mutant lines.

Pure culture – a culture of only one organism (i.e., not contaminated, not a mixed culture).

Quality control strain – a strain previously determined to meet the performance criteria of a particular test. Seed stocks are maintained at low passage numbers. Each passage is re-authenticated and characterised.

Seed stock – cultures stored under optimal conditions and used to prepare distribution stock. Ideally, seed stock is from very low-passage material. At ATCC, every effort is made to minimise the number of transfers. Contamination and genetic drift are minimised. One vial of seed stock is used to test for viability and purity, as well as to authenticate and characterise the strain.

Stock culture – a culture derived from a master culture or a reference culture and used as a reagent on a regular basis. Official standards and specifications mandate propagation procedures, growth conditions and limits on passage numbers for stock cultures used in individual tests.

Strain – cells descended from a single isolation in pure culture, usually from a single colony, though not necessarily from a single cell (i.e., not necessarily a clone).

Subculture – (see Passage)

Test organism – an organism used in a particular assay or other test.

References

1. International Code of Nomenclature of Bacteria. Washington DC: American Society for Microbiology, 1992.
2. Brown MRW; Gilbert P. Microbiological Quality Assurance: a guide towards relevance and reproducibility. Boca Raton, FL: CRC Press, 1994.
3. National Committee for Clinical Laboratory Standard (now CLSI) M22-A2. Quality assurance for commercially prepared microbiological culture media. Second edition. Approved standard. Wayne, PA: NCCLS, 1996.

Reprinted from ATCC® Technical Bulletin no. 1.

Preservation and recovery of filamentous fungi

Preservation methods for filamentous fungi vary depending on the type and degree of sporulation. Spore forming strains (with the exception of zoosporic fungi) can usually be freeze dried successfully. Similar success with nonsporulating strains is far less likely. Both types can be frozen and stored for long periods in liquid nitrogen or liquid nitrogen vapour. All plasmid containing or mutant strains should be frozen directly from the original material if possible to prevent any alteration or loss of desired characteristics. The following overview discusses preservation methods used at ATCC and is presented for descriptive purposes only; it is not intended as a laboratory protocol. Anyone planning to preserve cultures by these methods is strongly advised to study detailed, published protocols before proceeding.

Freeze drying sporeformers; preparing the cultures

1. Grow fungi under conditions that will induce maximum sporulation so that sufficient spores will survive the freezing and drying process. Optimum media and growth conditions are listed in the strain descriptions in ATCC's online catalogue at www.lgcstandards-atcc.org. Literature cited there may also give further guidance on appropriate cultivation procedures.
2. Prepare a 20% solution of skim milk and autoclave at 116°C for 20 minutes in 10 ml tubes. One tube is usually more than enough for 10 freeze dried vials unless the culture is very mycelial. Store the milk at 2-8°C until needed so that it will be cold when used.
3. Prepare the spore suspension by slowly introducing about 2 ml of milk into the culture tube or plate while gently scraping the surface of the culture with a pipette. Take care to avoid raising a cloud of spores, especially with *Aspergillus*, *Penicillium*, and other fungi that produce large amounts of dry spores. For example, *Neurospora* spores are very difficult to contain; cultures

for freeze drying are usually grown on agar in 250 ml Erlenmeyer flasks.

4. Transfer the suspension back into the tube containing the remainder of the skim milk and mix thoroughly. If more than one plate or tube is used, repeat the procedure for each and pool the suspensions into one tube. A concentration of at least 10⁶ spores per ml of milk is needed.
5. Dispense 0.2 ml of the suspension into each vial for freeze drying. Many spores will begin to germinate when suspended in liquid, so timing is critical when filling vials. Spores should not be in the skim milk for more than two hours before being processed. Refrigerate filled vials while waiting for further processing.

Freeze drying methods

Although the process is somewhat labour intensive, freeze drying spore forming fungi greatly facilitates their distribution and storage. The following four methods are described in detail in Simione and Brown (1991).

1. **Component freeze dryer.** Samples are freeze dried in cotton plugged glass inner vials which are then sealed inside glass outer ampoules under vacuum. The components of the system (vacuum pump and condenser) are assembled on a bench top.
2. **Commercial freeze dryer.** Samples are prepared in glass vials in ampoules (as in no.1 above), and then freeze dried in a commercial freeze dryer.
3. **Serum vial.** Samples are processed in glass serum vials sealed with rubbers stoppers and metal caps. A commercial freeze dryer is used.
4. **Manifold.** Samples are processed in bulb shaped or tubular glass ampoules attached with latex tubing to a manifold. (This method is not used at ATCC, but is relatively inexpensive and uses equipment that a lab may already own.)

Recovery

1. If the culture to be recovered was obtained directly from ATCC, check it thoroughly upon receipt. If it is found to be unsatisfactory in any respect, notify ATCC so that the strain in question can be investigated.

2. Open the ampoule as directed and, using sterile distilled water and a sterile pipette, transfer the contents of the preparation to approximately 5 ml of sterile distilled water in a test tube.

3. Allow the contents to rehydrate for at least one hour (two hours is better, and overnight is not too long) before transfer of a few drops to broth or agar. Use the media and growth conditions specified in the strain descriptions when first subculturing to ensure optimal recovery.

4. Incubate at the appropriate temperature. The remainder of the suspension may be stored for a few days if refrigerated, allowing for another recovery attempt if the first should fail. Given proper treatment and conditions, most cultures will grow in a few days. However, some may exhibit a prolonged lag period and should be given twice the normal incubation time before being discarded as nonviable.

5. For special media, growth conditions, and tips on maintenance, carefully read the literature cited in the strain descriptions. This is especially important for producers of secondary metabolites and quality control strains.

Freezing filamentous fungi; preparing the culture

Grow sporulating strains on solid media as for freeze-drying; grow nonsporulating strains on either solid or liquid medium. If the mycelium is easily broken, grow the culture on agar in test tubes, scrape with a pipette, and suspend the fragments in sterile 10% glycerol. Dispense 0.5 ml into each plastic vial. If the mycelium is sticky, will not break up, or grows embedded in the agar, grow the culture on agar in plates, cut out plugs containing new growth (hyphal tips) with a sterile cork borer, and place three or four plugs into

each plastic freezing vial with approximately 0.4 ml of 10% glycerol. If the mycelium does not break up or does not adhere to the agar well enough to cut with a cork borer, grow the culture in liquid medium in flasks and use a blender to prepare the culture for freezing.

1. Before blending a culture, sterilise the blender jar and lid by wrapping them separately in foil and autoclaving at 121°C for 15 minutes. Unwrap the jar and check the bolt to ensure that it is tight.

2. Pour the contents of the flask (growing culture and liquid medium) down the side of the blender jar, being careful to minimise splatter. Replace the lid and set the blender jar firmly on its base.

3. Blend on medium or high for 2 or 3 seconds. Observe the appearance of the contents directly or draw up some of the suspension into a wide mouthed pipette to ensure that small segments of mycelia are homogeneously suspended throughout the medium. If the culture is especially tough and does not blend well, blend for 2 or 3 seconds longer.

4. Pour the contents back into the flask, again being careful to minimise splatter. Incubate slow growing strains for several days before freezing to allow for healing. Fast-growing strains can be frozen one day after blending.

5. Add a portion of this suspension to an equal amount of sterile 20% glycerol, for a final concentration of 10% glycerol. Mix well and dispense approximately 0.5 ml into each freezing vial.

Warning: Pathogens should be handled under a hood and should not be blended. *Histoplasma*, *Paracoccidioides*, and *Blastomyces* should be frozen in the yeast phase and *Coccidioides* in the young mycelial stage to avoid contamination from airborne spores and to minimise exposure of laboratory personnel.

Freezing methods

To ensure long-term viability of fungal cultures, ATCC recommends freezing and storage at liquid nitrogen temperatures. Storage in the liquid itself is not always convenient or safe. Storage in liquid nitrogen vapour is a more practical alternative. It is critical, however, to constantly monitor the liquid level of the liquid nitrogen freezer to ensure that material is maintained below -130°C . Storage at warmer temperatures can compromise the stability of many strains.

1. Select a container appropriate for the material to be preserved. Plastic screw-capped vials with internal tube threads (1.0 to 2.0 ml) are appropriate for most fungi. They are sterilised by the manufacturer and, when properly handled, remain sterile throughout labelling and dispensing.
2. Label each vial clearly and accurately. Whatever labelling method is chosen, labels must be able to withstand subsequent freezing and thawing procedures.
3. Prepare the cultures for freezing as described above. Seal the plastic vials as tightly as possible with the screw cap.
4. Load the vials onto aluminium canes and record the location of each culture.
5. Cool the chamber of a controlled-rate cooling apparatus to 4°C and place the canes into the unit. Insert the thermocouple probe into one sample containing only sterile medium. Prior to starting the cooling program, allow the material to cool to within 2°C of the chamber temperature.
6. Cool the material at a rate of 1°C per minute to -40°C , then cool 10°C per minute to -90°C .
7. When the program is complete, transfer the canes to boxes in the vapour phase of a liquid nitrogen unit for storage. If the unit is more than a few feet from the programmable freezer, transport the canes in an insulated container with liquid nitrogen. Be careful not to store the

vials directly in liquid nitrogen. See notes under "Safety" below. Liquid nitrogen has a temperature of -196°C . The vapour will have a temperature gradient which is near -196°C at the level of the liquid and which gradually becomes warmer near the top of the freezer. For best long-term storage, keep frozen materials below -130°C . To ensure that materials are stored at proper temperatures, liquid nitrogen freezers must first be validated by placing a thermometer at the top of the unit and adding liquid nitrogen until a working temperature of at least -130°C is maintained at the top of the freezer. This level is then continuously monitored and an alarm system activated if the levels fluctuate above or below predetermined limits. Occasionally, strains may require special handling. If a strain does not survive freezing in glycerol, try 5% DMSO. Some cultures, such as some *Agaricus* strains, may be grown on sterile seeds, grains or pollen, and may be frozen without a cryoprotectant.

Recovery

Thaw frozen cultures quickly in a 37°C water bath, transfer immediately to appropriate growth media, and incubate at an appropriate temperature. Given proper treatment and conditions, most cultures will grow in a few days. However, some may exhibit a prolonged lag period and should be given twice the normal incubation time before being discarded as non-viable. For special media, growth conditions, and tips on maintenance, carefully read the literature cited in the strain descriptions. This is especially important for transformation hosts, genetic mutants, producers of secondary metabolites, and quality control strains.

Safety

Safety precautions must be considered when preserving living cells and microorganisms by freeze drying, freezing, and storing at cryogenic temperatures.

Cryogenic storage

Because of its extremely cold temperature, liquid nitrogen can be hazardous if improperly

used. When handling liquid nitrogen, take precautions to protect your face and exposed skin from exposure to the liquid. Wear protective clothing, including a laboratory coat, gloves designed for handling material at cryogenic temperatures, and a face shield. To reduce your exposure to cryogenic temperatures, design inventory systems for storing frozen specimens that allow for easy retrieval and that minimise the time required to look for specimens. Prolonged exposure to cryogenic temperatures can lead to a loss of sensation in the hands that can only be recovered after warming. This loss of sensation can lead to a false sense of security regarding damage to tissues by the low temperatures. When the temperature in a liquid nitrogen unit becomes tolerable and working in the unit is no longer uncomfortable, the operator has reached a point where damage from the cryogenic temperatures is likely. When liquid nitrogen is used in confined and inadequately ventilated areas, the nitrogen can quickly displace the room air. Liquid nitrogen freezers should be located in well-ventilated areas, and special precautions should be taken during fill operations. In facilities with several liquid nitrogen freezers, an oxygen monitor should be installed to warn occupants of any deterioration in the air quality due to the nitrogen gas. Plastic screw-capped vials can present a hazard if stored directly in liquid nitrogen. Vials with an inadequate seal between the cap and the vial can fill with liquid nitrogen. Upon retrieval to warmer temperatures the vials may explode violently or may spray liquid, potentially disseminating the contents of the vial. Likewise when opening plastic vials after thawing some dissemination of the contents may occur. Material in plastic ampoules is therefore stored in the vapour above the liquid nitrogen.

Freeze drying

When freeze drying microorganisms in vials or ampoules without cotton plugs or other bacteriological filters, the microorganisms can be carried from the container and contaminate the freeze drying system. Microbial contamination can be found on the outside of the vial or

ampoule, and on parts of the freeze drying system such as the condenser. A system should be designed to monitor the contamination level, and decontamination procedures should be implemented if necessary. Take care to properly treat freeze dried cultures prior to disposal. To autoclave freeze dried cultures, open the vial or ampoule to allow penetration of the steam. An alternative to autoclaving is to heat the preparations in a hot air oven at 180°C for four hours.

Culture handling

When opening frozen or freeze dried cultures take care to prevent dispersion of the ampoule contents. Open these preparations in a biological safety cabinet if possible, and perform all work with hazardous cultures in a biological safety cabinet. There are varying degrees of pathogenicity among microorganisms. All laboratory personnel should be aware of the hazards posed by the cultures they are handling. Detailed discussions of laboratory safety procedures are provided in the U.S. Dept. of Health and Human Services / CDC publication *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington D.C. U.S. Government Printing Office; 2007. This publication is available in its entirety on the CDC Office of Health and Safety website at www.cdc.gov/od/ohs.

Cultivation and preservation literature

For any strain listed in ATCC's electronic catalogue, note the recommended media and incubation conditions. The literature cited here is useful for general knowledge on the cultivation and preservation of a wide variety of fungi.

Atlas RM. *Handbook of Microbiological Media*. 2nd ed. Boca Raton, Fla.: CRC Press; 1996.

Atlas RM. *Handbook of Media for Environmental Microbiology*. Boca Raton, Fla.: CRC Press; 1995.

Brown RW, Gilbert P. *Microbiological Quality Assurance: A Guide Towards Relevance and Reproducibility of Inocula*. Boca Raton, Fla.: CRC Press; 1995.

Booth C (ed.). *Methods in Microbiology*. Volume 4. London: Academic Press; 1971.

Chang ST, Hayes WA. *The Biology and Cultivation of Edible Mushrooms*. New York: Academic Press; 1978.

Chang ST, Miles PG. *Edible Mushrooms and Their Cultivation*. Boca Raton, Fla.: CRC Press; 1989.

Chang ST, Quimio TH. *Tropical Mushrooms: Biological Nature and Cultivation Methods*. Hong Kong: The Chinese University Press; 1982.

Dhingra OD, Sinclair JB. *Basic Plant Pathology Methods*. 2nd ed. Boca Raton, Fla.: CRC Press; 1995.

Elliott TJ (ed). *Science and Cultivation of Edible Fungi*. Volumes 1 and 2. (Mushroom Science XIV) Rotterdam: A.A. Balkema; 1995.

Fassatiava O. *Moulds and Filamentous Fungi in Technical Microbiology*. Amsterdam: Elsevier; 1986.

Flegg PB, Spencer DB, Wood DA (eds.). *The Biology and Technology of the Cultivated Mushroom*. New York: John Wiley; 1985.

Fletcher JT, White PF, Gaze RH. *Mushrooms: Pest and Disease Control*. 2nd ed. Andover, Hants, England: Intercept; 1989.

Fuller MS, Jaworski A (eds). *Zoospore Fungi in Teaching and Research*. Athens, GA: Southeastern Publishing Corp; 1987.

Hunter-Cevera JC, Belt A. *Maintaining Cultures for Biotechnology and Industry*. New York: Academic Press; 1996.

Johnston A, Booth C. *Plant Pathologist's Pocketbook*. 2nd ed. Farnham Royal, Slough, England: CAB International; 1983.

King AD Jr., Pitt JL, Beuchat LR, Corry JEL (eds). *Methods for the Mycological Examination of Food*. New York: Plenum Press; 1986.

Simione F, Brown EM. *ATCC Preservation Methods: Freezing and Freeze-drying*. 2nd ed. Rockville, Md.: ATCC; 1991.

Smith D, Onions AHS. *The Preservation and Maintenance of Living Fungi*. 2nd ed. Wallingford, Oxon, England: CAB International; 1994.

Stamets P. 1993. *Growing Gourmet & Medicinal Mushrooms*. Berkeley, Calif.: Ten Speed Press.

Tuite J. *Plant Pathological Methods: Fungi and Bacteria*. Minneapolis, Minn.: Burgess Publishing Co; 1969.

van Griensven LJLD. *The Cultivation of Mushrooms*. Rustington, Sussex, England: Darlington Mushroom Laboratories; 1988.

Wuest PJ, Royse DJ, Beelman RB (eds.). *Cultivating Edible Fungi*. Amsterdam: Elsevier; 1987.

ATCC® and Preceptrol® are trademarks of the American Type Culture Collection.

Reprinted from ATCC® Technical Bulletin no. 2

Authentication of prokaryotes at ATCC

ATCC's goal is to fully characterise every prokaryote that is accepted into the collection. Considering the diversity of organisms housed at ATCC it is a constant challenge to streamline the authentication process so that microbes can be characterised with the greatest efficiency. By utilising several diverse identification strategies at the phenotypic and genotypic levels ATCC has developed protocols that ensure thorough characterisation of every strain.

Quality control (QC) begins when a new strain arrives at ATCC and continues through each step of the manufacturing process (Figure 1). Testing protocols and results become part of the laboratory record of the strain. This technical bulletin provides an overview of prokaryotic characterisation at ATCC.

Basic testing

The first step is to check the growth, purity, and cell and colony morphologies of each culture that arrives at ATCC for deposit. Cultures then undergo biochemical testing if appropriate. For many of the more common bacteria, ATCC has developed a standard set of growth and biochemical methods based on the phenotypic traits of these organisms. These biochemical tests, which include the use of API strips (bioMérieux, Inc.) and other commercial rapid tests, have been developed for different bacterial groups and are summarised in Table 1. ATCC has refined this schema over the years to minimise the number of tests yet still provide robust identification. Even with these refinements there are 23 different characterisation schemes encompassing nearly 350 individual tests that are performed on a routine basis on newly deposited organisms.

The biochemical tests are instrumental for authenticating many important phenotypic properties of ATCC microbes. They are time consuming, however, and the repertoire

of biochemical tests is very limited for some groups, such as extremophiles and fastidious organisms like mycobacteria or actinomycetes. Therefore whenever possible we take a polyphasic approach to authentication that elucidates both phenotypic and genotypic traits.

Automated phenotypic testing

In the last few years instrumentation has been developed that allows standardised testing of phenotypic and genotypic traits across a wide range of organisms. ATCC is constantly evaluating and utilising these methods to balance selectivity, throughput, cost, and effort to ensure the best quality for our authentication procedures.

Phenotypic methods currently used at ATCC include fatty acid methyl ester analysis (FAME), the Biolog system, and matrix-assisted laser desorption ionisation time-of-flight mass spectrometry (MALDI-TOF-MS or MALDI). FAME analysis provides a fatty acid profile using a standardised procedure for deriving the methyl esters of bacterial fatty acids and then analysing the products on a gas chromatograph. The resulting chromatogram provides a quantitative profile of an organism's fatty acids that can be used for identification¹. ATCC uses a commercial FAME system produced by MIDI (Newark, Delaware). The MIDI system includes databases to which new organisms can be compared to obtain presumptive identification.




Figure 1. Processing of new accessions at ATCC.

FAME analysis is useful in this respect, since there is a correspondence between phylogenetic groups of organisms and the dominant fatty acids in their cell membranes. FAME analysis has limitations, however it does not work with archaea, which do not have conventional fatty acids in their cells walls. It is also quite sensitive to both the growth medium and incubation time and temperature, since bacteria modulate the fatty acid contents of their membranes in response to environmental factors.

The Biolog system (Biolog, Inc., Hayward, California) tests the metabolic phenotype of an organism using a 96-well plate format. Each well contains a different carbon source. The test organism is inoculated into the plate and incubated, usually for 24 hours. The pattern of wells that are positive or negative for substrate utilisation is then analysed on a plate reader and matched to a database which provides a presumptive identification for the organism. This identification has no phylogenetic context;

however, it does provide useful information about the substrates an organism has the potential to metabolise. There are a variety of plates available for Gram-negative and Gram-positive organisms as well as some other types.

MALDI-TOF-MS is a technique that was first introduced a decade ago, and only in the last few years has it began being used as a tool for bacterial typing^{2,3}. MALDI uses a laser to irradiate intact cells, which are placed in a special matrix that ionises cell surface components, principally proteins. The masses of these ionised molecules are registered and a spectral profile or fingerprint is produced for an organism.

The potential advantages that MALDI offers over other techniques include minimal sample preparation, rapid results, and very low reagent costs. Perhaps most important is that this technique is dependent on a trait common to all prokaryotes: a complex macromolecular cell wall. Therefore it is a single authentication technique that covers the entire prokaryotic world.

ATCC evaluated a MALDI system produced by Micromass (Manchester, UK) that is specifically designed for bacterial identification including development of databases. We have evaluated over 50 prokaryotic genera and have found that MALDI produced consistent patterns for all of them. An example of organisms that have been analysed is shown in Figure 2.

Genotypic characterisation

A major thrust of microbial identification research has been the development of molecular methods for genotyping microorganisms. Genotypic methods can be highly specific and sensitive and are largely independent of the physiological or growth state of the organism.

For several years, ATCC has been using ribotyping as part of its QC process for bacteria. Ribotyping is done using the RiboPrinter microbial characterisation system (Qualicon, Wilmington, Delaware). The RiboPrinter is a

molecular workstation that performs a restriction digest (using EcoRI or other restriction enzymes) of the chromosomal DNA, separates the restriction fragments by gel electrophoresis, and simultaneously blots the DNA fragments to a membrane which is used for Southern blot analysis⁴. Restriction digest fragments are hybridised to a bacterial probe that is based on the conserved regions of the genes for the ribosomal DNA operon. The result is a DNA fingerprint which is strain specific. Each fingerprint is stored in a database so it can be accessed for future comparisons and identifications^{5,6}.




Figure 2. A sample of MALDI-TOF profiles. The two strains of *T. fusca*, a thermophilic actinomycete, show similar but distinct profiles. *P. furiosus* is a thermophilic archaeon that can grow at temperatures of 100°C. The three *Mycobacterium* strains are typically difficult to differentiate; their mass spectral profiles are similar but unique.




Figure 3. Riboprints of two colony variants of *Bacillus subtilis* ATCC® 6051™.

ATCC now ribotypes over 250 genera on the RiboPrinter and found that over 80% yield reproducible riboprints. As a result ATCC now ribotypes the original isolate that is received, the frozen seed stock, and each batch of distribution stock as part of our standard QC process for many bacteria. In addition the RiboPrinter can be useful for resolving phenotypic abnormalities among strains and ruling out contamination. For example, in the case of the *Bacillus subtilis* shown in Figure 3, ribotyping can show that two colony types from the same strain are genotypically identical.

The RiboPrinter system allows the flexibility in using other restriction enzymes beside the standard EcoRI. Some genera yield better patterns with PvuII or PstI. Other restriction enzymes can also be substituted in the process. Other advantages include the ease of operation and data analysis. The system requires only a single colony as inoculum and there are no restrictions on media and growth conditions.

Rep-PCR analysis

Ribotyping does not work for all organisms, including archaea and other challenging environmental microbes and some fastidious bacteria like Mycobacteria. In some cases ATCC uses sequencing of the 16S rRNA gene to confirm an organism's identity. ATCC uses rep-PCR as another genotyping tool. This technique takes advantage of specific repetitive DNA sequences that are randomly interspersed around the chromosome. By using PCR primers specific for these elements, it is possible to amplify the intervening sequences and separate them by gel electrophoresis, thus creating a genotypic fingerprint that is strain specific⁷. What is remarkable about this technique is that primer sets are available that allow the genotyping of archaea and other difficult organisms (Figure 4). Furthermore, rep-PCR appears to work well for fungi and protists, giving it the potential to be a universal typing technique.

ATCC incorporates the bioMérieux Diversilab system that uses lab-on-a-chip technology to enhance the speed, precision and reproducibility of the entire process.

Authentication of all our microorganisms is a task ATCC takes seriously. ATCC will continue to use a polyphasic approach that balances traditional phenotypic methods with the newest technologies to ensure the delivery of high quality microbial strains.




Figure 4. Fingerprinting of different archaea using rep-PCR. Lane A: MW standards; B: *Methanococcus maripaludis* ATCC® 43000™; C: *Methanocaldococcus jannaschii* ATCC® 43067™; D: *Haloferax volcanii* ATCC® 29605™; E: *Thermococcus gorgonarius* ATCC® 700654™

References

1. Suzuki K et al. Cell envelopes and classification. In: O'Donnell AG and Goodfellow M, eds. Handbook of new bacterial systematics. London: Academic Press; 1993.
2. Holland RD et al. Rapid Commun. Mass Spectrom. 10(10): 1227-32, 1996.
3. Lay JO Jr. and Holland RD. Methods Mol. Biol. 146: 461-87, 2000.
4. Bruce J. Food Technol. 50: 77-81, 1996.
5. Tang J et al. ATCC Connection 23 (1): 1 and 4-5, 2003.
6. Cordevant C et al. J. Clin. Microbio. 41(1): 34-43, 2003.
7. Versalovic J et al. Meth. Cell. Mol. Biol. 5: 25-40, 1994.

Table 1. Schema of biochemical characterisation tests performed at ATCC.

Phenotype	Number of tests	Types of tests
Actinobacilli not requiring X or V	37	Acid production from 17 sugars; enzymatic tests; misc. biochemical tests
Agromyces	35	Utilisation of 13 substrates; temperature range; enzymatic tests; misc. biochemical tests
Anaerobic cocci	38	Utilisation of 10 substrates; enzymatic tests; misc. biochemical tests
Bacteroides	45	Utilisation of 16 substrates; enzymatic tests; misc. biochemical tests
Bifidobacterium	48	Acid production from 19 substrates; misc. biochemical tests; enzymatic tests
Campylobacter/Helicobacter	27	Temperature range; sensitivity to 3 antibiotics; misc. biochemical tests
Clostridium	40	Utilisation of 12 substrates; endospore formation and location; enzymatic tests; misc. biochemical tests
Corynebacterium	34	Acid production; enzymatic tests; misc. biochemical tests
Erysipelothrix	44	Misc. biochemical and hydrolytic tests; acid production from 21 substrates; hemolytic activity
Eubacterium	43	Utilisation of 15 substrates; enzymatic tests; misc. biochemical tests
Gram-negative aerobic bacteria	120	Temperature range; NaCl tolerance; O ₂ tolerance; NO ₃ utilisation; acid production from 16 substrates; 54 sole carbon source tests
Gram-negative facultative bacteria	79	NO ₃ utilisation; acid production from 16-26 substrates; gas production; misc. biochemical and hydrolytic tests
Gram-positive aerobic cocci	48	NaCl tolerance; temperature range; sugar utilisation; misc. hydrolytic and biochemical tests
Gram-positive aerobic sporeformers	66	Substrates utilisation; NaCl tolerance; pH, and temperature ranges; O ₂ tolerance and NO ₃ utilisation; misc. biochemical and hydrolytic tests
Gram-negative fastidious bacteria	61	Growth on 7 complex substrates; acid production from 8 substrates; NaCl tolerance; misc. biochemical and hydrolytic tests
Lactic-acid bacteria	60	Acid production from 24 substrates; temperature, NaCl and pH optima; misc. biochemical and hydrolytic tests
Mycobacteria	25	Temperature range and growth rate; pigment production; misc. biochemical tests; misc. biochemical tests
Neisseria	30	Pigmentation; penicillin sensitivity; misc. biochemical tests
Organisms requiring X and/or V factors	30	Utilisation of 10 sugars; enzymatic tests; misc. biochemical tests
Pasteurellaceae not requiring X or V	56	Acid production from 28 sugars; enzymatic tests; misc. biochemical tests
Propionibacterium	38	Acid production; enzymatic tests; misc. biochemical tests
Streptomycetes/Actinomycetes	28	Utilisation of 11 carbon sources; morphology of conidia and sporangia; pigmentation; enzymatic tests
Treponema	37	Utilisation of 8 sugars; enzymatic tests; misc. biochemical tests

Reprinted from ATCC® Technical Bulletin no. 5

Reference strains: How many passages are too many?

ATCC cultures have been cited in national and international standards for many years. Examples include the standards promulgated by the Clinical and Laboratory Standards Institute (formerly NCCLS) for the healthcare community and by the United States Pharmacopeia (USP) for the pharmaceutical and biopharmaceutical industries. Cultures are used in performance testing of products, as positive and negative controls, as indicator organisms and as identification standards.

Though the use of microbiological standards is widely accepted, there is still some confusion as to specific laboratory guidelines, especially when determining the number of subcultures allowed beyond the reference strain. As recently as 2009 discussions about passages took place on the Pharmaceutical Microbiological Mail List (PMFLIST)1-3. This article will attempt to clear up some of this confusion and provide some definitions and recommendations.

Strain definitions

The confusion starts with the different names that are ascribed to these cultures. In various CLSI and USP publications these cultures are called control strains, standard cultures, reference strains, test strains, and quality control strains. These terms can generally be used interchangeably, though the preference seems to be reference strain or reference culture. Both the CLSI and USP agree that reference strains should come from a reliable source; both organisations cite ATCC. There is agreement that the reference strains from ATCC are subcultured to make “stock cultures,” which are subcultured weekly or monthly to make the “working cultures” used daily. Working cultures are often kept as slants, and it is these subcultures that raise the questions about passages from the original reference strain.

A subculture is a passage. The USP26 Antimicrobial Effectiveness Testing, section 51 states: “For the purposes of the test, one passage is defined as the transfer of organisms from an established culture to fresh medium. All transfers are counted⁴.”

This definition was updated in the Pharmacopeial Previews to read: “One passage is defined as the transfer of organisms from a viable culture to fresh medium with growth of the microorganisms. Any form of subculturing is considered to be a transfer/passage⁵.”

This updated definition is preferable. The earlier definition left questions about the meaning of “an established culture.” There were several questions raised on the PMFLIST as to whether the frozen or freeze dried vial from ATCC was an established culture. To some, the phrase “established culture” implied a growing culture. It is clear, however, that these frozen or freeze dried vials of reference strains from ATCC are indeed “viable” cultures. A passage involves growing the microorganism with fresh medium, either on solid agar or in broth. Resuscitating frozen or freeze dried cultures by thawing or rehydrating is not by itself considered a passage. Thus subculturing the reference strain from ATCC’s frozen or freeze dried vial to the stock culture is the first passage. Subculturing stock cultures to working cultures is the second passage from the original reference strain. Any subsequent subculture is another cumulative passage.

Authenticity of cultures

Another term used by commercial sources other than ATCC is “ATCC derived.” A culture is derived from the ATCC reference strain by subculturing in other words, one or more passages. In most cases it is not known how many passages are used by the commercial source to produce their cultures. Some commercial sources charge a premium for cultures certified to be “only” two passages

from the ATCC. CLSI recognises reliable commercial sources for reference strains. However, you must take into account that these cultures are at least one or two passages from the ATCC reference strains.

Maintaining cultures

The USP recommends a seed lot system to maintain reference strains in laboratories. In this system the ATCC reference strain is subcultured to several replicates at one time, all of which are within one passage. These replicates of the stock culture are the seed vials for the laboratory. The seed stock is subcultured, the second passage, to make replicates of working cultures. This is the same system that the ATCC uses to minimise passages in its culture collection (Figure 1).

But how many passages are recommended or acceptable in the laboratory? There is agreement that the number of passages should be minimised to reduce the possibility of phenotypic variations, genetic drift, and contamination as much as possible, but standards organisations differ as to how many passages are acceptable.




Figure 1. Seed lot system

Over the years the recommendations from CLSI have varied. The least specific recommendations call for subculturing stock cultures weekly with new working cultures subcultured monthly, with no maximum number of passages noted. Another CLSI standard recommends up to three subcultures of the stock cultures and up to three subcultures of the working cultures⁶. This would add up to up to seven passages from the original ATCC reference culture (one to make the first stock culture, plus three subcultures and three additional subcultures).

USP standards have been more specific. USP clearly states that the working cultures used for testing should not be more than five passages from the ATCC reference culture. The USP26 states: “The viable microorganisms used in this test must not be more than five passages from the original ATCC culture⁴.”

The USP26 also contains the following definition: “Microbial Strains – Where a microbial strain is cited and identified by its ATCC catalogue number, the specified strain shall be used directly or, if subcultured, shall be used not more than five passages removed from the original strain⁷.”

The recommendation of five passages or less from the ATCC reference culture has been broadly accepted in the healthcare community and the pharmaceutical and biopharmaceutical industries. ATCC agrees with this recommendation.

Utilising cold storage

Storage temperature of stock and working cultures can affect growth characteristics and viability. The CLSI recommendations include storage at -50 to -70°C for one year or below -70°C indefinitely⁵, or -20°C or below (preferably below -70°C) for “prolonged” storage⁸. Storage of slants is recommended at 2 to 8°C for either one week^{8,9} or two weeks¹⁰. The USP26 recommends storage in liquid nitrogen or a mechanical freezer below -50°C .

For long-term storage of frozen cultures, ATCC recommends the vapour phase of liquid nitrogen or a mechanical freezer at -80°C . Immersion in liquid nitrogen is not recommended. Frozen cultures may be kept at -20°C for short-term storage (less than one month). Do not store frozen cultures in a freezer with a defrost cycle; this will expose the cultures to higher temperatures. Freeze dried cultures should be stored at 2 to 8°C . Slants can be kept at 2 to 8°C for up to a week.

Conclusion: ATCC recommendations:

1. ATCC reference strains should be subcultured to replicate stock cultures in the laboratory. Stock cultures can be subcultured for working cultures weekly, typically kept as slants. A seed lot system is recommended.

2. A passage is defined as a subculture involving growth of the viable microorganism with fresh medium. Thawing or rehydrating ATCC reference cultures is not a passage.

3. Microorganisms for standard protocols should be used within five passages of the ATCC reference culture.



4. Frozen cultures should be stored in the vapour phase of liquid nitrogen or in a mechanical freezer at -80°C or below. Freeze dried cultures should be stored at 2 to 8°C . Slants may be stored at 2 to 8°C for up to a week.

References

1. Friedel B. Culture passage guidelines. In: PMFLIST Archives [Internet]. [The Pharmaceutical Microbiology Mail List] 2003 Jan 2. Available from: <http://peach.ease.lsoft.com/archives/pmflist.html>
2. Kramer K. Seed stock culture. In: PMFLIST Archives [Internet]. [The Pharmaceutical Microbiology Mail List] 2003 June 10. Available from: <http://peach.ease.lsoft.com/archives/pmflist.html>
3. Anger C. Seed stock culture. In: PMFLIST Archives [Internet]. [The Pharmaceutical Microbiology Mail List] 2003 June 11. Available from: <http://peach.ease.lsoft.com/archives/pmflist.html>
4. U.S. Pharmacopeia. Antimicrobial Effectiveness Testing. In: U.S. Pharmacopeia 26th rev. Rockville, MD: U.S. Pharmacopeia; 2003, 2nd Supplement <51>.
5. U.S. Pharmacopeia. Microbiological Good Laboratory Practices, <1117>. Pharmacopeial Forum 29(3): Pharmacopeial Previews, 2003.
6. NCCLS. Quality assurance for commercially prepared microbiological culture media. 2nd ed. Wayne, PA: NCCLS; M22-A2, 1996.
7. U.S. Pharmacopeia. General Notices. In: U.S. Pharmacopeia 26th rev. Rockville, MD: U.S. Pharmacopeia; 2003, 2nd Supplement.
8. NCCLS. Performance standards for antimicrobial disk susceptibility tests. 8th ed. Wayne, PA: NCCLS; M2-A8, 2003.
9. NCCLS. Performance standards for antimicrobial disk and dilution susceptibility tests for bacteria isolated from animals; approved standard. 2nd ed. Wayne, PA: NCCLS; 2002.
10. NCCLS. Reference method for broth dilution antifungal susceptibility testing of yeasts; approved standard. Wayne, PA: NCCLS; M27-A, 1997.

Reprinted from ATCC Connection 23(2): 6-7, 2003, with 2010 updates.

- All cultures should be regarded as potentially hazardous and should be opened by persons trained in microbiological techniques working in facilities with containment requirements appropriate for the biosafety level of the cultures.
- Work in a biological safety cabinet. If this is not possible, wear suitable eye protection. Hold vials away from face and over a can or tray.
- Wear gloves.
- Sterilise all empty vials and fragments before disposal.





Rehydrating Freeze-Dried Cultures Instructional Guide

Check each culture thoroughly upon receipt. If a culture is unsatisfactory, notify ATCC so that the strain in question can be investigated. Use the medium and incubation conditions specified on the product information sheet when first reviving strains to ensure optimal conditions for recovery.

If cultures are not revived immediately, store vials at 5°C or colder (plant viruses at -20°C or colder).

Bacteria and algae

Aseptically add 0.5 mL of liquid medium to the freeze-dried material with a sterile Pasteur pipette and mix well. For bacteria, transfer the total mixture to a test tube containing 5 to 6 mL of the recommended broth medium. The last few drops of this suspension also may be transferred to an agar slant. Algae cultures must be initiated on agar plates.

Incubate cultures under the appropriate conditions. Given proper treatment and conditions, most freeze-dried cultures will grow out in a few days. However, some may exhibit a prolonged lag period and should be given twice the normal incubation time before discarding as nonviable.

Bacteriophages

1. Prepare an actively growing broth culture of the host before opening the phage specimen. The host should be 18 to 24 hours old prior to inoculation.
2. To propagate phage: Aseptically rehydrate the freeze-dried phage specimen with 1.0 mL of appropriate broth (see maintenance conditions for host) and mix well. Add 0.1 mL of this mixture to 0.9 mL of broth medium and serially dilute as desired. One drop of each dilution can be spotted onto prepared plates (three or four dilutions per plate). Lysis should be visible after overnight incubation.
3. To prepare a high-titre phage stock: Refer to the product sheet for instructions.

Filamentous fungi and yeast

Use a Pasteur pipette to add 0.5 to 1 mL sterile water to the freeze-dried pellet then draw up the entire contents into the pipette and transfer to a test tube with about 5 mL sterile water.

Let the yeast or fungus rehydrate for a minimum of 2 hours (overnight is not too long) before transferring to broth or solid agar. Incubate at the recommended temperature. Keep in mind that some cultures may exhibit a prolonged lag period and should be given twice the normal incubation time before discarding as nonviable.

Save the mixture of lyophilised material and water until you know you have growth. If not contaminated, it will keep for several days in a refrigerator.

If you subsequently contaminate your culture, you can recover the desired microorganism by serial dilution and picking single colonies.

Plant viruses

Plant viruses are usually distributed in freeze-dried plant tissues within single, sealed vials. Open the vials carefully, removing the metal retaining cap and the rubber stopper. If the vial has been flame-sealed, open according to directions on the previous page.

For tissue reconstitution and inoculum preparation, place the contents of the vial in a precooled (4°C) mortar with 2 to 3 mL of an inoculation buffer (e.g., 0.05 M sodium phosphate buffer, pH 7.0, with 10 mM sodium sulfite). Triturate the tissue thoroughly with a pestle to prepare the inoculum.

Rub the inoculum onto host plants using a sterile cotton swab and a fine abrasive, such as 500- to 600-mesh Carborundum® (silicon carbide) or Celite® (diatomaceous earth). The abrasive may be added to the inoculum (50 to 100 mg/mL) or dusted onto the plants prior to inoculation. After inoculation, spray the plants with water to remove buffer salts and abrasive.

Cultures in stoppered serum vials

Open serum vials carefully by aseptically removing the metal retaining cap and the rubber stopper, then follow instructions for the appropriate type of organism.

If you have any questions, please contact a technical service representative or your local distributor. Details regarding our warranty can be found in the Material Transfer Agreement packed in your shipment or available at www.lgcstandards-atcc.org.

MATERIAL TRANSFER AGREEMENT ("MTA")

Last updated February 1, 2010



IMPORTANT! PLEASE READ CAREFULLY BEFORE SUBMITTING AN ORDER. THIS IS A CONTRACT.

This Material Transfer Agreement ("MTA") is between the _____, a _____ type of organization, having its principal place of business at _____ ("Purchaser") and the American Type Culture Collection, a not-for-profit organization, having its principal place of business at 10801 University Boulevard, Manassas, VA 20110-2209 ("ATCC"). Purchaser must have an approved, current ATCC account to place an order. This MTA is effective as of the last date of execution by the parties and governs the purchase and use of all ATCC Materials under the terms and conditions set forth below.

TERMS AND CONDITIONS

Definitions

"**ATCC Material(s)**" means materials acquired from ATCC as documented on an ATCC Sales Order.

"**ATCC Sales Order**" means an order submitted for ATCC Materials in a form and format as determined by ATCC from time to time.

"**Biological Material(s)**" means ATCC Materials, Progeny, Unmodified Derivatives and any Unmodified Derivatives within Modifications, either individually or jointly.

"**Commercial Use**" means the sale, license, lease, export, transfer or other distribution of the Biological Materials to a person or entity not party to this MTA for financial gain or other commercial purposes and/or the use of the Biological Material: (a) to provide a service to a person or entity not party to this MTA for financial gain; (b) to produce or manufacture products for general sale or products for use in the manufacture of products ultimately intended for general sale (c) in connection with ADME (Absorption, Distribution, Metabolism and Excretion) testing; (d) in connection with drug potency or toxicity testing which does not include either screening multiple cell lines for potential inclusion in a screening assay system or screening multiple compounds in a system for internal research purposes only; (e) in connection with proficiency testing service(s), including but not limited to, providing the service of determining laboratory performance by means of comparing and evaluating calibrations or tests on the same or similar items or materials in accordance with predetermined conditions; or (f) for research conducted under an agreement wherein a for-profit entity receives a right whether actual or contingent to the results of the research. Commercial Use specifically does not include Industry Sponsored Academic Research.

"**Contributor(s)**" means an organization(s) and/or individual(s) providing original material to ATCC for deposit.

"**Industry Sponsored Academic Research**" means research sponsored by a for-profit organization carried out at a non-profit organization and by the non-profit organization's employees.

"**Investigator**" means the Purchaser's principal scientist or researcher using the Biological Material(s).

"**Modification(s)**" mean substances created by Purchaser which contain and/or incorporate a significant or substantial portion of ATCC Material.

"**Progeny**" means an unmodified descendant from the ATCC Materials, such as virus from virus, cell from cell, or organism from organism.

"**Purchaser(s)**" means the organization purchasing and receiving ATCC Material pursuant to this MTA.

"**Unmodified Derivative(s)**" mean substances created by Purchaser that constitute an unmodified functional sub-unit or product not changed in form or character and expressed by the ATCC Material provided by ATCC. Unmodified Derivatives include, but are not limited to, subclones of unmodified cell lines, purified or fractionated subsets of materials provided by ATCC, proteins expressed by DNA/RNA supplied by ATCC, or monoclonal antibodies secreted by a hybridoma cell line.

Scope of Use

Subject to the terms of this MTA, Purchaser's Investigator may make and use the Biological Materials provided to Purchaser by ATCC for research purposes only in Purchaser's Investigator's laboratory only. The Biological Materials are not intended for use in humans. Purchaser agrees that Biological Materials designated as biosafety level 2 or 3 constitute known pathogens and that other Biological Materials not so designated may be pathogenic under certain conditions. Purchaser assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, transfer and Purchaser's use of the Biological Materials including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. Purchaser agrees that any activity undertaken with the Biological Materials will be conducted in compliance with all applicable guidelines, laws and regulations, and that Purchaser will obtain all permits, licenses or other approvals required by any governmental authority in connection with Purchaser's receipt, handling, storage, disposal, transfer and use of the Biological Materials.

Purchaser shall not distribute, sell, lend or otherwise transfer, to a person other than the Purchaser's Investigator, or entity not party to this MTA, the Biological Material, as defined above, for any reason, without ATCC's prior written agreement.

Any Commercial Use of the Biological Material is strictly prohibited without ATCC's prior written consent. Purchaser acknowledges and agrees that Purchaser's use of certain Biological Material may require a license from a person or entity not party to this MTA, or be subject to restrictions that may be imposed by a person or entity not party to this MTA ("Third Party Terms"). To the

extent of ATCC's knowledge of the existence of any such applicable rights or restrictions, ATCC will take reasonable steps to identify the same, either in ATCC's catalog of ATCC Materials and/or through ATCC's customer service representatives, and to the extent they are in the possession of ATCC, ATCC shall make information regarding such Third Party Terms reasonably available for review by Purchaser upon request. Purchaser expressly acknowledges that if there is a conflict between this MTA and the Third Party Terms, the Third Party Terms shall govern. Use of the Biological Materials may be subject to the intellectual property rights of a person or entity not party to this MTA, the existence of which rights may or may not be identified in the ATCC catalog or website, and ATCC makes no representation or warranty regarding the existence or the validity of such rights. Purchaser shall have the sole responsibility for obtaining any intellectual property licenses necessitated by its possession and use of the Biological Materials.

The use permitted under this MTA for Industry Sponsored Academic Research extends only to the academic research carried out at the non-profit organization and the non-profit organization's employees. Any non-profit Purchaser using the Biological Materials in connection with Industry Sponsored Academic Research agrees to notify the industrial sponsor that any use of the Biological Materials by the industry sponsor will require a separate license from ATCC and/or its Contributors and that ATCC and/or its Contributors are under no obligation hereunder to license any Biological Materials to any such industry sponsor.

Warranty; Warranty Disclaimer

ATCC warrants that (a) cells and microorganisms included in the ATCC Material shall be viable upon initiation of culture for a period of thirty (30) days after shipment thereof from ATCC and (b) any ATCC Material other than cells and microorganisms shall meet the specifications on the applicable ATCC Material product information sheet, certificate of analysis, and/or catalog description until the expiration date on the applicable ATCC Material's product label (such thirty (30) day period, or period until the expiration date, referred to herein as the "Warranty Period"). Purchaser's exclusive remedy, and ATCC's sole liability, for breach of the warranties set forth in this paragraph is for ATCC to, at ATCC's sole option, either (i) refund the fee paid to ATCC for such ATCC Material (exclusive of shipping and handling charges), or (ii) replace the ATCC Material. The warranties set forth in this paragraph apply only if Purchaser handles and stores the ATCC Material as described in the applicable ATCC Material product information sheet. To obtain the exclusive remedy, Purchaser must report the lack of viability or non-conformation to specifications to ATCC's Technical Service Department within the applicable Warranty Period. Any expiration date specified on the ATCC Material shipment documentation states the expected remaining useful life, but does not constitute a warranty or extend any applicable Warranty Period. **Except as expressly provided above, the ATCC Material and any technical information and assistance provided by ATCC are provided as-is, without warranties of any kind, express or implied, including but not limited to any implied warranties of merchantability, fitness for a particular purpose, typicity, safety, accuracy and/or non-infringement.**

Compliance with Laws

Purchaser is solely responsible for compliance with all foreign and domestic, federal, state and local statutes, ordinances and regulations applicable to use of the Biological Material. Without limiting the generality of the foregoing, any shipment of Biological Materials to countries outside the United States must comply with all applicable foreign and U.S. laws, including the U.S. export control laws and related regulations. Distribution by ATCC of Budapest Treaty patent deposits are made pursuant to, and in compliance with, all applicable laws and regulations, including the Budapest Treaty and related 37 C.F.R. provisions. If there is any conflict between the terms of this MTA and any applicable law or regulation with respect to Materials that are supplied hereunder by ATCC from the stock of a Budapest Treaty deposit, then the terms of the applicable law or regulation shall govern.

Indemnification

If Purchaser is a for-profit or private non-profit organization:

Purchaser hereby agrees to indemnify, defend and hold harmless ATCC and its Contributors against all third party claims, losses, expenses and damages, including reasonable attorneys' fees (collectively "Claims") arising out of or relating to Purchaser's use, receipt, handling, storage, transfer, disposal and other activities relating to Biological Materials, provided that Purchaser's liability shall be limited to the extent that any such Claim arises out of ATCC's gross negligence or willful misconduct. All non-monetary settlements of any such Claims are subject to ATCC's prior written consent, such consent not to be unreasonably withheld.

If Purchaser is a Federal or State non-profit organization or a foreign organization that is prohibited by law from entering into the indemnification obligation set forth in the above paragraph:

Purchaser assumes all liability for any and all third party claims, losses, expenses and damages, including reasonable attorneys' fees (collectively "Claims") arising out of or relating to Purchaser's use, receipt, handling, storage, transfer, disposal and other activities relating to Biological Materials, provided that Purchaser's liability shall be limited to the extent that any such Claim arises out of ATCC's gross negligence or willful misconduct, and provided further that if the Purchaser is the U.S. federal government or a state institution such Purchaser assumes such liability only to the extent provided under the Federal Tort Claims Act, 28 U.S.C. §§ 2671 et seq. or under equivalent applicable State or foreign law.

Limitation of Liability

In no event will ATCC or its Contributors be liable for any indirect, special, incidental or consequential damages of any kind in connection with or arising out of the MTA or Biological Materials (whether in contract, tort, negligence, strict liability, statute or otherwise) even if ATCC has been advised of the possibility of such damages. In no event shall ATCC's cumulative liability to the Purchaser exceed the fees paid by Purchaser under this MTA for the twelve (12) month period preceding the date of the event giving

rise to the claim. Purchaser agrees that the limitations of liability set forth in this MTA shall apply even if a limited remedy provided hereunder fails of its essential purpose.

Intellectual Property; Identification

As between the parties, ATCC and/or its Contributors shall retain ownership of all right, title and interest in the ATCC Materials, Progeny, Unmodified Derivatives and Biological Materials contained or incorporated in Modifications. Purchaser retains ownership of: (a) Modifications (except that, as between the parties, ATCC retains ownership rights to Biological Material included therein) and (b) those substances created through the use of Biological Material, but which do not contain Biological Material. Notwithstanding the foregoing, Purchaser acknowledges and agrees that the Biological Materials are subject to the restrictions noted in the "Scope of Use" section above. Purchaser agrees to acknowledge ATCC and any Contributor indicated by ATCC as the source of the Biological Material in all research, academic or scholarly publications and in patent applications that reference the Biological Material. If required by the Contributor of the ATCC Material, ATCC may inform the Contributor of Purchaser's identity. Purchaser explicitly acknowledges that ATCC retains all right, title and interest in the ATCC trademarks, trade-names, logos, ATCC catalog numbers and ATCC specific designations of ATCC Materials sold by ATCC (including but not limited to ATCC[®], UNIPLUS[™], YOUR DISCOVERIES BEGIN WITH US[®], THE GLOBAL BIORESOURCE CENTER[™], Authenticult[™], SafeTsource[™], ATCC CULTURES[™], ATCC BIOPRODUCTS[™], ATCC SPECIAL COLLECTIONS[™], ATCC SERVICES[™], ATCC Genuine Cultures[®], ATCC Licensed Derivative[®], BioEscrow[®], ATCC Standards Resource[®], ATCC Proficiency Standard[®], ATCC Standard Reference Material[™]). Purchaser expressly agrees not to use the ATCC trademarks, trade-names, logos, ATCC catalog numbers or ATCC specific designations of ATCC Materials sold by ATCC in any way without ATCC's prior written agreement.



Miscellaneous

Any disputes arising under this Agreement shall be tried exclusively in the United States District Court for the Eastern District of Virginia or if subject matter jurisdiction does not exist in that court, then in the state courts of Virginia for Prince William County, and Purchaser hereby expressly consents to, submits to and waives any objection to the jurisdiction of such courts; provided however, if Purchaser is a US Federal or State non-profit organization; then any disputes arising under this Agreement shall be tried in exclusively in a court of competent jurisdiction.

Purchaser agrees that any breach of this Agreement, including but not limited to any breach of the scope of use provisions of this Agreement, will entitle ATCC to immediately cease without notice to Purchaser further shipments of Biological Materials and may create such irreparable injury as to entitle ATCC to seek temporary restraining orders and other preliminary or permanent injunctive relief in addition to all other equitable and legal remedies that may be afforded under US or foreign laws.

Purchaser may not assign or otherwise transfer this MTA or any rights or obligations under this MTA, whether by operation of law or otherwise. Any such attempted assignment or transfer will be void and of no force or effect. This MTA, including all documents incorporated herein by reference, constitutes the entire agreement between ATCC and Purchaser with respect to the Biological Material and supersedes all previous agreements or representations (whether written or oral) between ATCC and Purchaser relating to the same subject matter. This MTA may not be modified, waived or terminated except in writing and signed by the parties hereto. No term or provision contained herein shall be deemed waived and no breach excused unless such waiver or consent shall be in writing and signed by the parties. If any provision of this MTA is for any reason found to be unenforceable, the remainder of this Agreement will continue in full force and effect. None of the provisions of this MTA are intended to create, nor shall be deemed or construed to create, any relationship between ATCC or Purchaser other than that of independent entities contracting with each other hereunder solely for the purpose of effecting the provisions of this MTA.

Any correspondence concerning the ATCC Material Transfer Agreement should be addressed to ATCC, Attention: Office of IP, Licensing and Services, P.O. Box 1549, Manassas, VA 20108, Phone: (703) 365-2700 or contact us at licensing@atcc.org



CERTIFIED REFERENCE MATERIAL

BIOLOGICAL STANDARDS WITH UNPRECEDENTED RELIABILITY ONLY FROM ATCC

ATCC is pleased to introduce Certified Reference Material (CRM) produced under ISO Guide 34: 2000 and ISO/IEC 17025: 2005 accredited and ISO 9001: 2008 certified laboratories.

Why use certified reference material?

The inherent variability of biological materials brings unique challenges to establishing standards for in vitro model systems and to establishing ISO compliant processes to produce them. Biological reference materials produced under an ISO Guide 34 accredited process have confirmed identity, well-defined characteristics and an established chain of custody. These qualities make CRMs ideal as biological standards for research and development purposes.

What are the features of ATCC Certified Reference Materials?

- Optimised for homogeneity and consistency from lot to lot
- Tested to a stated level of confidence, where appropriate
- Traceable pedigree to the ATCC seed culture
- Verified using polyphasic (genotypic and phenotypic) testing to confirm identity
- Unique level of accuracy and lot specific results
- Comprehensive documentation, including property values and their uncertainty calculations, expiration dating and serial numbering of vials

What are the benefits of Certified Reference Materials from ATCC?

Certified Reference Material offers the highest level of quality assurance, accuracy and traceability. They give you complete confidence that your results are reliable and reproducible.

When are ATCC Certified Reference Materials recommended for use?

- For establishing sensitivity, linearity and specificity during assay validation or implementation
- For challenging assay performance

- For validating or comparing test methods
- For testing and calibration in ISO 17025 labs which stipulate the use of reference materials
- For benchmarking critical assay performance for regulatory submissions and production lot release
- For Pharmacopeia compendial tests



CURRENTLY AVAILABLE

ATCC® No.	Item Description	Designation
CELL LINES		
CRM-CCL-2™	HeLa cells	
CRM-CRL-1550™	CaSki cells	
CRM-HTB-31™	C-33 A	
BACTERIA		
CRM-6633™	<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	NRS-231
CRM-11437™	<i>Clostridium sporogenes</i>	L.S. McClung 2006
CRM-8739™	<i>Escherichia coli</i>	Crooks
CRM-11229™	<i>Escherichia coli</i>	AMC 198
CRM-9341™	<i>Kocuria rhizophila</i>	FDA strain PCI 1001
CRM-9027™	<i>Pseudomonas aeruginosa</i>	R. Hugh 813
CRM-6538™	<i>Staphylococcus aureus</i>	FDA 209
CRM-12228™	<i>Staphylococcus epidermidis</i>	FDA strain PCI 1200
YEAST AND FUNGI		
CRM-16404™	<i>Aspergillus brasiliensis</i>	WLRI 034 (120)
CRM-10231™	<i>Candida albicans</i>	3147
CRM-9763™	<i>Saccharomyces cerevisiae</i>	NRRL Y-567

LGC Standards offices

France

LGC Standards S.a.r.l.
6, rue Alfred Kastler
BP 83076
67123 MOLSHEIM CEDEX
FRANCE

Tel: +33 (0)3 88 04 82 82
Fax: +33 (0)3 88 04 82 90
fr@lgcstandards.com

Countries served: Belgium,
France, Luxembourg, Monaco,
Switzerland.

Italy

LGC Standards S.r.l.
Via Venezia, 23
20099 Sesto San Giovanni (MI)
ITALY

Tel: +39 02 24126 830
Fax: +39 02 24126 831
it@lgcstandards.com

Countries served: Italy.

Sweden

LGC Standards AB
PO Box 1737
SE-501 17 BORÅS
SWEDEN

Tel: +46 (0)33 20 90 60
Fax: +46 (0)33 20 90 79
atcc.se@lgcstandards.com

Countries served: Denmark,
Estonia, Finland, Iceland, Latvia,
Norway, Sweden.

Germany

LGC Standards GmbH
Mercatorstrasse 51
46485 Wesel
GERMANY

Tel: +49 (0)281 9887 230
Fax: +49 (0)281 9887 239
atcc.de@lgcstandards.com

Countries served: Austria, Croatia,
Germany, Greece, Hungary, Kosovo,
Moldova, Mongolia, Montenegro,
Netherlands, Romania, Serbia,
Slovenia, Turkey.

Poland

LGC Standards Sp. z o.o.
ul. M. Konopnickiej 1
Dziekanow Lesny
05-092 Lomianki
POLAND

Tel: +48 (0)22 751 31 40
Fax: +48 (0)22 751 58 45
pl@lgcstandards.com

Countries served: Armenia,
Azerbaijan, Belarus, Bulgaria, Czech
Republic, Georgia, Kazakhstan,
Kyrgyzstan, Lithuania, Macedonia,
Poland, Slovak Republic, Tajikistan,
Turkmenistan, Ukraine, Uzbekistan.

United Kingdom

LGC Standards
Queens Road
Teddington
Middlesex TW11 0LY
UK

Tel: +44 (0)20 8943 8489
Fax: +44 (0)20 8943 8405
atcc@lgcstandards.com

Countries served: Albania, Bosnia
Herzegovina, Channel Islands,
Cyprus, Ireland, Liechtenstein, Malta,
United Kingdom.

India

LGC Promochem India Private Limited
488B, Fourth Floor
VITC Model Export Bhavan
14th Cross, IV Phase
Peenya Industrial Area
Bangalore 560 058
INDIA

Tel: +91 80 6701 2000
Fax: +91 80 6701 2046
in@lgcpromochem.com

Countries served: India.

Spain

LGC Standards S.L.U
C/ Salvador Espriu, N° 59 2° planta
08005 Barcelona
SPAIN

Tel: +34 93 308 41 81
Fax: +34 93 307 36 12
es@lgcstandards.com

Countries served: Andorra,
Portugal, Spain.

LGC Standards, Queens Road, Teddington, Middlesex TW11 0LY, UK

www.lgcstandards.com



Excellence through measurement

No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2011. All rights reserved.

© 2011 American Type Culture Collection. The ATCC trademark and trade name, any and all ATCC catalogue numbers and any other trademarks listed in this publication are trademarks of the American Type Culture Collection unless indicated otherwise. Picogreen® and Ribogreen® are registered trademarks of Invitrogen. GenBank® is a trademark of the United States Department of Health and Human Sciences.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services or to manufacture commercial products without prior ATCC written approval.

2339/CF/1211