
BinaxNOW[®]

Staphylococcus aureus Card

Product Instructions / fr-Instructions produit / de-Gebrauchsanweisung / it-Istruzioni per il prodotto / pt-Instruções do produto / es-Instrucciones del producto / da-Produktvejledning / nl-Productinstructies / sv-Produktanvisningar / no-Produktinstruksjoner / el-Οδηγίες προϊόντος

Materials Provided / fr-Matériel Fourni / de-Im Lieferumfang enthaltene Materialien / it-Materiali Forniti /
pt-Materiais Fornecidos / es-Materiales Suministrados / da-Medfølgende materialer / nl-Meegeleverde materialen /
sv-Tillhandahållet Material / no-Materialer Som Følger Med / el-Παρεχόμενα υλικά

1

2


3

4

5

6

Intended Use

The Alere BinaxNOW® *Staphylococcus aureus* Card is a qualitative, *in vitro* immunochromatographic assay for the presumptive identification of *Staphylococcus aureus*. The test is performed directly on blood culture samples positive for Gram-positive cocci in clusters. The Alere BinaxNOW® *Staphylococcus aureus* Card is not intended to diagnose *Staphylococcus aureus* nor to guide or monitor treatment for *Staphylococcus aureus* infections. Subculturing positive blood cultures is necessary to recover organisms for susceptibility testing and/or differentiation of mixed growth.

Summary and Explanation of the Test

Staphylococcus aureus, a Gram-positive, coagulase-positive bacteria, is a major cause of infection in both healthcare and community settings. *Staphylococcus aureus* is the leading cause of nosocomial bloodstream infections¹ and over half of all sepsis cases are estimated to be caused by Gram-positive cocci in clusters (GPCC) bacteria.²

Hospitalized patients with a *Staphylococcus aureus* infection have five times the risk of in-hospital mortality compared with inpatients without this infection.³ Rapid identification of the pathogen associated with bloodstream infection is crucial to improving patient outcomes, length of hospitalization, and health care costs.⁴ Alere BinaxNOW® *Staphylococcus aureus* Card provides a simple and rapid method to detect *Staphylococcus aureus* directly from blood culture bottles positive for Gram-positive cocci in clusters.

Principles of the Procedure

Alere BinaxNOW® *Staphylococcus aureus* Card is a rapid immunochromatographic membrane assay that uses highly sensitive polyclonal antibodies to detect a *Staphylococcus aureus*-specific protein directly from blood cultures which have been identified as being positive for Gram-positive cocci in clusters. The protein-specific antibodies and a control antibody are immobilized onto a test strip as two distinct lines and combined with other reagents/pads. This test strip is mounted inside a cardboard, book-shaped hinged test card.

Specimens are aliquots from blood cultures which have been identified as positive for Gram-positive cocci in clusters. After the sample is prepared, it is added to the sample pad at the top of the test strip, and the card is closed. Results are read at 10 minutes.

Reagents and Materials

See pull-out flap for illustrations.

Materials Provided in the Alere BinaxNOW® *Staphylococcus aureus* Card Kit

- 1 **Test Cards:** A cardboard, book-shaped, hinged test card containing the test strip.
- 2 **Reagent A:** Clear, alkaline solution that lyses red blood cell membranes and exposes the protein target.
- 3 **Reagent B:** Clear, slightly acidic solution used to neutralize the sample.
- 4 **Reagent C:** Lyophilized enzyme which lyses the bacterial membrane and further exposes the target.
- 5 **Reagent C Diluent:** Clear, acidic solution used to reconstitute Reagent C.
- 6 **Reagent D:** Clear, slightly alkaline solution that clears particles and reduces background on the membrane.

Materials Required, But Not Provided:

Clock, timer or stopwatch; a centrifuge capable of 1,500 x g or equivalent speed in rpm, sample tubes suitable for processing and centrifuging samples, calibrated pipettes, vortex mixer.

Precautions

1. For *in vitro* diagnostic use.
2. Allow all kit components to equilibrate to room temperature (15-30°C) before use.
3. Leave test card sealed in its foil pouch until just before use.
4. Avoid skin and eye contact with reagents and card strip.
5. Do not use kit past its expiration date.
6. Do not interchange or mix components from different kit lots.
7. As blood culture samples, controls and test cards may contain pathogenic organisms, handle with appropriate precautions and dispose of materials safely in biohazard waste receptacles.
8. Alere BinaxNOW® *Staphylococcus aureus* Card should only be performed on blood culture samples known to be positive for Gram-positive cocci in clusters.
9. Reagent C must be added using a calibrated pipette; failure to add the correct volume of reagent to the sample may result in erroneous results.
10. The processed sample must be added to the card using a calibrated pipette; failure to add the correct volume of sample to the card may result in erroneous results.
11. Once reconstituted, Reagent C must be stored at 2-8°C. After use, return Reagent C or the entire kit to the refrigerator (2-8°C).

12. Centrifuge speed must be set precisely at 1,500 x g. For centrifuges set in rpm speed, it is important to convert the 1,500 x g into equivalent rpm according to the centrifuge's manual. Failure to use the correct centrifuge speed may result in erroneous results.

Storage and Stability

Store test cards, Reagent A, Reagent B, un-reconstituted Reagent C, Reagent C diluent and Reagent D at 2-30°C.

Once reconstituted with Reagent C diluent, Reagent C must be refrigerated (2-8°C). For convenience, the entire kit may be refrigerated (2-8°C), or reconstituted Reagent C may be removed from the kit box and refrigerated at 2-8°C. However, if Reagent C is removed from the test kit for refrigeration, do not interchange or mix components from different kit lots.

Alere BinaxNOW® *Staphylococcus aureus* Card and reagents are stable until the expiration dates marked on their outer packaging and containers.

Quality Control

Daily Quality Control:

Alere BinaxNOW® *Staphylococcus aureus* Card has built-in positive and negative procedural controls. For quality control, Alere suggests that you record these controls for each test run.

Procedural Controls:

- A. The appearance of a pink-purple line at the “control line” position can be considered an internal positive procedural control. If capillary flow has occurred, this line will always appear.
- B. The clearing of the background color from the result window is a negative background control. In comparison to the color of the control line, the background color in the window should be white within 10 minutes.

External Positive and Negative Controls:

Good laboratory practice suggests the use of positive and negative controls to ensure that:

- test reagents are working, and
- the test is correctly performed.

For each new lot, external positive and negative controls should be tested. These will monitor the entire assay. The recommended controls for Alere BinaxNOW® *Staphylococcus aureus* Card are as listed below. The results of these control samples should be recorded. Please refer to the Control Sample Preparation Procedure section for further instructions.

Positive control: *Staphylococcus aureus* ATCC#12598

Negative control: *Staphylococcus epidermidis* ATCC#14990

Other controls may be tested in order to conform with:

- local, state and/or federal regulations,
- accrediting groups, and/or
- your laboratory's standard Quality Control procedures.

If the correct control results are not obtained, do not report results. Contact Alere™ Technical Service during normal business hours.

Specimen Collection and Handling

Specimens are aliquots from blood cultures which have been identified as containing Gram-positive cocci in clusters. Test as soon as possible after the blood culture instrument signals the culture is positive. If it is necessary to delay testing, blood cultures can be held in the blood culture instrument or at room temperature (15 - 30°C) for up to 24 hours after signaling positive.

Within the 24 hours after the sample has signaled positive, it is acceptable to freeze the sample at -20°C or -80°C for later testing in Alere BinaxNOW® *Staphylococcus aureus* Card. If frozen blood culture samples are to be tested, allow samples to warm to room temperature before testing in Alere BinaxNOW® *Staphylococcus aureus* Card.

Blood Culture Systems:

The performance of Alere BinaxNOW® *Staphylococcus aureus* Card was validated in Alere's clinical study using BacT/ALERT® Automated Microbial Detection System and the BACTEC® Instrumented Blood Culture System. Use of any other blood culture systems would require independent validation.

Reagent Preparation Procedure

1. Allow reagents and cards to equilibrate to room temperature (15-30°C) before testing.
2. Reconstitute **Reagent C**.
3. Remove the aluminum seal and the rubber stopper from the **Reagent C** vial. Twist the tip off the **Reagent C diluent**.
4. Empty contents of **Reagent C diluent** into the **Reagent C** vial. Place rubber stopper firmly into the vial.
5. Thoroughly mix by vortexing.

Control Sample Preparation:

1. Prepare a 4.0 McFarland Standard suspension in sterile saline using a fresh isolate of the recommended staphylococcal strains.
2. Transfer 0.5 mL of the suspension to a tube appropriate for centrifuging.

- Obtain a patient blood culture bottle that has recently been determined to be negative for microbial growth by an automated blood culture system. Gently mix by inverting and remove a 0.5 mL aliquot from the blood culture bottle and add it to the sample tube.
- Add 1 mL of **Reagent A** to the sample tube.
- Thoroughly mix by vortexing.
- Centrifuge the sample tube for 5 minutes at 1500 x g (or equivalent speed in rpm).
- Carefully aspirate and discard as much of the supernatant as possible. Do not disturb the pellet. **Note:** *It is important to remove all of the supernatant without disturbing the pellet. If the pellet is pulled into the pipette, aspirate the sample back into the tube and repeat procedure starting at Step 5.*
- Add 1 mL **Reagent A** to the pellet.
- Thoroughly resuspend the pellet by vortexing.
- Centrifuge the sample tube for 5 minutes at 1500 x g (or equivalent speed in rpm).
- Carefully aspirate and discard as much of the supernatant as possible. Do not disturb the pellet. **Note:** *It is important to remove all of the supernatant without disturbing the pellet. If the pellet is pulled into the pipette, aspirate the sample back into the tube and repeat procedure starting at Step 9.*
- Add 5 drops **Reagent B** to the pellet.
- Thoroughly resuspend the pellet by vortexing.
- Using a calibrated pipette, add 25 μ L of **Reagent C** to the sample tube and vortex thoroughly.
- Remove card from foil pouch just prior to testing. Open card and lay it flat on work surface. Ensure that the test card is aligned so the adhesive strip is on the right hand side.

For charcoal blood bottles ONLY, skip to step 17. For other blood culture bottles, proceed to step 16.

Non-charcoal blood bottles

- Using a calibrated pipette, carefully withdraw 50 μ L of sample for testing and go to **ASSAY PROCEDURE** below.

Charcoal blood bottles

- Centrifuge the sample for 5 minutes at 1500 x g (or equivalent speed in rpm).
- Using a calibrated pipette, carefully withdraw 50 μ L of supernatant for testing without disturbing the pellet, and go to **ASSAY PROCEDURE** below.

Patient Sample Preparation:

- Gently mix the blood bottle by inversion.
- Transfer 1 mL of blood culture bottle aliquot to a tube appropriate for centrifuging.
- Add 1 mL of **Reagent A** to the sample tube.
- Thoroughly mix by vortexing.
- Centrifuge the sample tube for 5 minutes at 1500 x g (or equivalent speed in rpm).
- Carefully aspirate and discard the supernatant. **Note:** It is important to remove all of the supernatant without disturbing the pellet. If the pellet is pulled into the pipette, aspirate the sample back into the tube and repeat procedure starting at Step 4.
- Add 1 mL of **Reagent A** to the pellet.
- Thoroughly resuspend the pellet by vortexing.
- Centrifuge the sample tube for 5 minutes at 1500 x g (or equivalent speed in rpm).

- Carefully aspirate and discard the supernatant. **Note:** It is important to remove all of the supernatant without disturbing the pellet. If the pellet is pulled into the pipette, aspirate the sample back into the tube and repeat procedure starting at Step 8.
- Add 5 drops of **Reagent B** to the pellet.
- Thoroughly resuspend the pellet by vortexing.
- Using a calibrated pipette, add 25 μL of **Reagent C** to the sample tube and vortex thoroughly.
- Remove test card from foil pouch. Open card and lay it flat on work surface. Ensure that the test card is aligned so that the adhesive strip is on the right hand side.

For charcoal blood bottles **ONLY**, skip to step 16. For other blood culture bottles, proceed to step 15.

Non-charcoal blood bottles

- Using a calibrated pipette, carefully withdraw 50 μL of sample for testing and go to **ASSAY PROCEDURE** below.

Charcoal blood bottles

- Centrifuge the sample for 5 minutes at 1500 x g (or equivalent speed in rpm).
- Using a calibrated pipette, carefully withdraw 50 μL of supernatant for testing without disturbing the pellet, and go to **ASSAY PROCEDURE** below.

Assay Procedure:

- Slowly add 50 μL of sample to the middle of the **WHITE** sample pad, located on the upper right hand side of the card.
- Add 4 drops of **Reagent D** to the **WHITE** wash pad, located on the top left hand side of the card.

- Slowly add one drop of **Reagent D** to the **TOP** of the **WHITE** sample pad located on the upper right hand side of the test card. Hold dropper bottle vertically while dispensing reagent.
- Immediately peel off adhesive liner from the test card. Close and securely seal the card by pressing firmly along the area to the right of the window. Read the result through the viewing window at 10 minutes. Results read before or after 10 minutes may be inaccurate. **Note:** When reading test results, tilt the card to reduce glare on the result window, if necessary.

Result Interpretation

For a **NEGATIVE SAMPLE**, the blue Control Line turns a pink-purple color. No other line appears.

For a **POSITIVE SAMPLE**, the blue Control Line turns a pink to purple color and a second pink to purple Sample Line appears above it in the top half of the window. Any Sample Line, even when very faint, is positive. **Note:** Strong positive results may cause a Control Line that is lighter in color than the Sample Line. This is still a valid result.

A test is **INVALID** if the blue Control Line does not turn a pink-purple color, or does not appear at all, whether a Sample Line is present or not. Repeat invalid tests with a new test card. Call Alere™ Technical Service if the problem persists.

	Positive	Negative	Invalid	Invalid
Sample				
Control				

Reporting of Result

Result Suggested Report

Positive Positive for *Staphylococcus aureus*

Negative Negative for *Staphylococcus aureus*

Invalid Indeterminate; do not report. Re-test sample and contact Alere™ Technical Service if problem persists.

Limitations

In analytical testing, the *Staphylococcus aureus* strain ATCC 14993 produced a negative result on Alere BinaxNOW® *Staphylococcus aureus* Card.

In the event of co-infection, the accuracy of this card for detecting *S. aureus* at the Limit of Detection (5.42×10^8 cells/mL) concentration in the presence of other bacteria at the concentration higher than 1×10^9 cells/mL has not been established and is therefore unknown.

Certain non-*S.aureus* strains may cause false positive results on Alere BinaxNOW® *Staphylococcus aureus* Card. These include strains of *Clostridium bifermentans*, *Clostridium perfringens*, and *Clostridium histolyticum*. Therefore, the use of samples with a mixed culture of Gram-positive rods and Gram-positive cocci in clusters on Alere BinaxNOW® *Staphylococcus aureus* Card will require further testing for verification.

Additionally, *Staphylococcus schleiferi* subsp. *coagulans* produced positive results on Alere BinaxNOW® *Staphylococcus aureus* Card in analytical testing.

The performance of Alere BinaxNOW® *Staphylococcus aureus* Card in blood culture systems other than the BacT/ALERT® Automated Microbial Detection System and the BACTEC® Instrumented Blood Culture System have not been evaluated and are therefore unknown.

Alere BinaxNOW® *Staphylococcus aureus* Card has not been validated with samples other than blood culture.

Performance Characteristics

Clinical Performance

The clinical performance of Alere BinaxNOW® *Staphylococcus aureus* Card was established in a multi-center clinical study conducted in 2008-9 at five geographically diverse hospital laboratories within the US.

A total of 623 blood culture samples with Gram-positive cocci in clusters were evaluated in Alere BinaxNOW® *Staphylococcus aureus* Card and compared to standard methods used routinely by the testing laboratories. Alere BinaxNOW® *Staphylococcus aureus* Card identified 93.3% of the specimens positive for *Staphylococcus aureus* and 99.6% of the specimens negative for *Staphylococcus aureus* relative to the reference method.

Alere BinaxNOW® *Staphylococcus aureus* Card Compared to Reference Methods

Alere BinaxNOW® <i>Staphylococcus aureus</i>	Reference Method	
	Positive	Negative
Positive	154	2
Negative	11	456

95% C.I.

Sensitivity: 93.3% (88.5-96.2%)

Specificity: 99.6% (98.4-99.9%)

Analytical Reactivity

The 54 human pathogenic Network on Antimicrobial Resistance in *Staphylococcus aureus* (NARSA) and American Type Culture Collection (ATCC) *Staphylococcus aureus* strains (see Table A on page 81) tested positive on Alere BinaxNOW® *Staphylococcus aureus* Card.

Analytical Specificity (Cross-Reactivity)

To determine the analytical specificity of Alere BinaxNOW® *Staphylococcus aureus* Card, coagulase-negative staphylococcal strains, yeasts and non-staphylococcal strains (see Tables B and C on pages 81 to 83) were tested with all strains testing negative.

Interfering Substances

None of the 20 potentially interfering substances (see Table D on pages 83 to 84) produced false results in Alere BinaxNOW® *Staphylococcus aureus* Card.

Reproducibility Study:

A study of Alere BinaxNOW® *Staphylococcus aureus* Card was conducted at 3 separate sites using panels of blind coded specimens containing negative and positive samples. Participants tested each sample twice on 5 different days. There was 98% (588/600) agreement with expected test results, with no significant differences within run (replicates tested by one operator), between run (5 different days), between sites (3 sites), or between operators (6 operators).

Ordering and Contact Information

Reorder numbers:

800-000: Alere BinaxNOW® *Staphylococcus aureus* 22-Card Kit

US: +1-877-441-7440

OUS: +1-321-441-7200

Technical Support

Advice Line

Further information can be obtained from your distributor, or by contacting Alere™ Technical Support on:

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Table B – Coagulase Negative Staphylococcal Strains

Bacterium
<i>Staphylococcus auricularis</i> ATCC 33753
<i>Staphylococcus capitis</i> ATCC 35661
<i>Staphylococcus caprae</i> ATCC 51548
<i>Staphylococcus cohnii</i> ATCC 29972
<i>Staphylococcus delphini</i> ATCC 49171
<i>Staphylococcus epidermidis</i> ATCC 700579
<i>Staphylococcus haemolyticus</i> ATCC 29970
<i>Staphylococcus hominis</i> ATCC 27844
<i>Staphylococcus hyicus</i> ATCC 11249
<i>Staphylococcus intermedius</i> ATCC 29663
<i>Staphylococcus kloosii</i> ATCC 43959
<i>Staphylococcus lentus</i> ATCC 700403

Bacterium
<i>Staphylococcus lugdunensis</i> ATCC 43809
<i>Staphylococcus lutrae</i> ATCC 700373
<i>Staphylococcus pasteurii</i> ATCC 51128
<i>Staphylococcus pseudintermedius</i> ATCC 49444
<i>Staphylococcus pulvereri</i> ATCC 51698
<i>Staphylococcus saprophyticus</i> ATCC 35552
<i>Staphylococcus schleiferi</i> ATCC 43808
<i>Staphylococcus sciuri</i> ATCC 49575
<i>Staphylococcus stimulans</i> ATCC 27851
<i>Staphylococcus vitulinus</i> ATCC 51162
<i>Staphylococcus warneri</i> ATCC 49454
<i>Staphylococcus xylosus</i> ATCC 49148

Table C – Non-Staphylococcal Strains

Bacterium	Bacterium
<i>Acinetobacter calcoaceticus</i> ATCC 51432	<i>Pasteurella multocida</i> ATCC 51687
<i>Aerococcus urinae</i> ATCC 700306	<i>Pediococcus acidilactici</i> ATCC 12697
<i>Aerococcus viridans</i> ATCC 10400	<i>Peptostreptococcus anaerobius</i> ATCC 27337

Bacterium	Bacterium
<i>Aeromonas hydrophila</i> ATCC 35654	<i>Planococcus citreus</i> ATCC 14404
<i>Bacillus cereus</i> ATCC 11778	<i>Proteus mirabilis</i> ATCC 7002
<i>Bacillus subtilis</i> ATCC 6633	<i>Proteus vulgaris</i> ATCC 33420
<i>Bacteroides fragilis</i> ATCC 23745	<i>Providencia stuartii</i> ATCC 49809
Beta Strep Group F ATCC 12392	<i>Pseudomonas aeruginosa</i> ATCC 15442
<i>Burkholderia cepacia</i> ATCC 25416-T	<i>Pseudomonas fluorescens</i> ATCC 49271
<i>Citrobacter freundii</i> ATCC 8090	<i>Pseudomonas putida</i> ATCC 49128
<i>Clostridium septicum</i> ATCC 12646	<i>Rhodococcus equi</i> ATCC 10146
<i>Clostridium sordelli</i> ATCC 9714	<i>Salmonella adelaide</i> ATCC 10718
<i>Corynebacterium amycolatum</i> ATCC 49368	<i>Serratia marcescens</i> ATCC 13880
<i>Corynebacterium diphtheriae</i> ATCC 13812	<i>Stenotrophomonas maltophilia</i> ATCC 13637-T
<i>Corynebacterium glutamicum</i> ATCC 13869	<i>Stomatococcus (Rothia mucilaginosa)</i> ATCC 25296
<i>Corynebacterium jeikeium</i> ATCC 43734	<i>Stomatococcus (Rothia mucilaginosa)</i> ATCC 49040
<i>Corynebacterium pseudodiphtheriticum</i> ATCC 10700-T	<i>Stomatococcus (Rothia mucilaginosa)</i> ATCC 49041
<i>Corynebacterium urealyticum</i> ATCC 43042	<i>Stomatococcus (Rothia mucilaginosa)</i> ATCC 49042
<i>Corynebacterium xerosis</i> ATCC 7711	<i>Streptococcus agalactiae</i> (Beta Strep Group B) ATCC 13813
<i>Enterobacter aerogenes</i> ATCC 35029	<i>Streptococcus anginosus (milleri)</i> ATCC 33397

Bacterium	Bacterium
<i>Enterobacter cloacae</i> ATCC 49141	<i>Streptococcus dysgalactiae</i> (Group C) ATCC 12388
<i>Enterococcus avium</i> ATCC 49462	<i>Streptococcus dysgalactiae</i> (Group G) ATCC 12394
<i>Enterococcus casseliflavus</i> ATCC 12817	<i>Streptococcus intermedius</i> (milleri) ATCC 27355
<i>Enterococcus durans</i> ATCC 49135	<i>Streptococcus mitis</i> ATCC 49456
<i>Enterococcus faecalis</i> ATCC 49474	<i>Streptococcus mutans</i> ATCC 25175
<i>Enterococcus faecium</i> ATCC 12952	<i>Streptococcus pasteurans</i> (bovis) ATCC 49133
<i>Enterococcus gallinarum</i> ATCC 49608	<i>Streptococcus pneumoniae</i> ATCC 33400
<i>Enterococcus hirae</i> ATCC 10541	<i>Streptococcus pneumoniae</i> ATCC 39938
<i>Enterococcus mundtii</i> ATCC 43187	<i>Streptococcus pneumoniae</i> ATCC 49136
<i>Enterococcus raffinosus</i> ATCC 49464	<i>Streptococcus pneumoniae</i> ATCC 49619
<i>Escherichia coli</i> ATCC 10798	<i>Streptococcus pneumoniae</i> ATCC 51937
<i>Gemella spp. bergeri</i> ATCC 700627	<i>Streptococcus pneumoniae</i> ATCC 51938
<i>Haemophilus influenzae</i> ATCC 49247	<i>Streptococcus pneumoniae</i> ATCC 6301
<i>Klebsiella oxytoca</i> ATCC 49131	<i>Streptococcus pneumoniae</i> ATCC SSI-1
<i>Klebsiella pneumoniae</i> ATCC 49472	<i>Streptococcus pneumoniae</i> ATCC SSI-10A
<i>Lactobacillus casei</i> ATCC 393	<i>Streptococcus pneumoniae</i> ATCC SSI-14
<i>Lactococcus spp. garvieae</i> ATCC 49157	<i>Streptococcus pneumoniae</i> ATCC SSI-7F
<i>Leuconostoc mesenteroides</i> ATCC 10877	<i>Streptococcus pyogenes</i> , Group A ATCC 12384

Bacterium	Bacterium
<i>Listeria monocytogenes</i> ATCC 19115	<i>Streptococcus salivarius</i> ATCC 13419
<i>Macrococcus caseolyticus</i> (formerly <i>Staph. cohnii</i> subsp. <i>cohnii</i>) ATCC 35662	<i>Streptococcus salivarius</i> ATCC 13419
<i>Micrococcus luteus</i> ATCC 27141	Yeasts
<i>Moraxella catarrhalis</i> ATCC 25238	<i>Candida albicans</i> ATCC 60193
<i>Morganella morganii</i> ATCC 25830-T	<i>Candida glabrata</i> ATCC 66032
<i>Neisseria meningitidis</i> (Serogroup A) ATCC 13077	<i>Candida tropicalis</i> ATCC 750
<i>Neisseria sicca</i> ATCC 9913	

Table D – Interfering Substances

Anti-Inflammatory Drugs	Test Concentration
Acetaminophen	1324 µmol/L
Acetylsalicylic acid	3.62 mmol/L
Ibuprofen	2425 µmol/L
Antibiotics	Test Concentration
Amoxicillin	206 µmol/L
Cephalexin	337 µmol/L
Chloramphenicol	155 µmol/L

Endogenous Blood Components	Test Concentration
Hemoglobin	2 g/L
Triglyceride sera	37 mmol/L
Conjugated bilirubin	342 $\mu\text{mol/L}$
Unconjugated bilirubin	342 $\mu\text{mol/L}$
γ - globulin	120g/L
Anti-coagulant	Test Concentration
Sodium Polyanetholesulfonate (SPS)	1%

References / fr -Références / de-Literaturhinweise / it - Bibliografia / pt - Referências / es-Referencias / da-Referencer / nl-Literatuur / sv-Källförteckning / no-Referanser / el-Βιβλιογραφία

1. Noskin G, Rubin R, Schentag J, Kluytmans J, Hedblom E, Smulders M, Lapetina, E, Gemmen E. The Burden of *Staphylococcus aureus* Infections on Hospitals in the United States: An Analysis of the 2000 and 2001 Nationwide Inpatient Sample Database. Arch Intern Med. 2005;165:1756-1761.
2. Wisplinghoff H, Bischoff T, Tallent SM, Seifert H, Wenzel RP, Edmond MB. Nosocomial bloodstream infections in US hospitals: analysis of 24,179 cases from a prospective nationwide surveillance study. Clin Inf Dis. 2004;39:309-317.
3. Kuehnert M, Hill H, Kurponis B, Tokars J, Solomon S, Jernigan D, CDC, Atlanta, Georgia, "Methicillin-Resistant *Staphylococcus aureus* Hospitalizations, United States", Emerging Infectious Diseases, CDC Website, June 2005.
4. Beekmann SE, Diekema DJ, Chapin KC, Doern GV. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. J Clin Microbiol. 2003 Jul;41(7):3119-25.

Alere Scarborough, Inc.
10 Southgate Road
Scarborough, Maine 04074 USA
www.alere.com

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands

© 2012 Alere. All rights reserved.
The Alere Logo, Alere and BinaxNOW are trademarks of the Alere group of companies.
BACTEC is a registered trademark of Becton, Dickinson and Company.
BacT/ALERT is a registered trademark of Biomerieux, Inc.
All trademarks referenced are trademarks of their respective owners.