

ELDOR CSEN

The Global Regional Anesthesia Distribution

ELDOR CSEN

cse_international@cse.com

[Mobile/Whatsapp/Wechat : +972 050 5582284](tel:+9720505582284)

Innovation - Quality – Price – Service

ELDOR CSEN Catalogue

CSEN 1: - Eldor Spinal Needle 26G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.cse.com/Eldor-dispersion.jpg>

http://cse.com/eldor_needle-2.jpg

<http://www.kumj.com.np/issue/4/263-266.pdf>

<https://www.yumpu.com/en/document/view/26271693/eldor-spinal-needle-versus-pencan-spinal-needle>

ELDOR SPINAL NEEDLE (Double Hole Pencil Point Spinal Needle – DHPP)- has ZERO Post Dural Puncture Headache (PDPH)

<http://www.csen.com/ESN-PDPH.pdf>

<https://www.youtube.com/watch?v=3IJNP-7hI9M>

<http://www.csen.com/DHPP.pdf>

<http://www.csen.com/TNS.pdf>

<http://www.csen.com/newf.pdf>

<http://www.csen.com/eldvid.pdf>

<http://www.csen.com/DMT.pdf>

<http://www.csen.com/OPUS%20MEDICAL%20-%20ELDOR%20Products.pdf>

www.csen.com/Epidural.pps

<http://www.csen.com/510K.pdf>

<http://www.csen.com/fda/>

<http://www.csen.com/cat12.pdf>

<http://www.csen.com/SPINAL.pdf>

(ELDOR spinal needle is fig. 23)

https://industrytoday.co.uk/health_and_safety/spinal-needle-market-2018--major-developments-in-healthcare-sector---transform-the-competitive-landscape-by-2023-with-6-8--of-cagr

A New Frontier in Spinal Anesthesia : ELDOR Spinal Needle

<http://www.csen.com/newf.pdf>

<https://encrypted-tbn2.gstatic.com/images?q=tbn:ANd9GcSPAasy7obfDbPO7eRvRg9fI9vE6jVYzsKJoHa72KZT0u3FLUR6>

https://www.google.co.il/imgres?imgurl=http%3A%2F%2Fpatentimages.storage.googleapis.com%2FEP0727187B1%2F00180001.png&imgrefurl=http%3A%2F%2Fwww.csen.com%2Fpic1.htm&docid=vY5pUKnUwxy_FM&tbnid=FajOYvyzzBzqZM%3A&vet=10ahUKEwit6fe7y-LbAhWhAsAKHdxjCBUQMwg-KAlwAg..i&w=2368&h=1412&bih=662&biw=1366&q=eldor%20spinal%20needle&ved=0ahUKEwit6fe7y-LbAhWhAsAKHdxjCBUQMwg-KAlwAg&iact=mrc&uact=8

<https://www.google.co.il/imgres?imgurl=http%3A%2F%2Fwww.csen.com%2Fanesthesia%2Fbook%2F12.jpg&imgrefurl=http%3A%2F%2Fwww.csen.com%2Fanesthesia%2Fproducts.htm&docid=RFh-Nr0-PQq7WM&tbnid=MtOE0C8DRoTehM%3A&vet=10ahUKEwit6fe7y-LbAhWhAsAKHdxjCBUQMwg8KAAwAA..i&w=498&h=305&bih=662&biw=1366&q=eldor%20spinal%20needle&ved=0ahUKEwit6fe7y-LbAhWhAsAKHdxjCBUQMwg8KAAwAA&iact=mrc&uact=8>

<http://www.csen.com/pic1.htm>

<http://www.csen.com/TSEN.pdf>
(Fig 3.D is the Eldor spinal needle)

<http://www.csen.com/pics/eldor-ad.jpg>

ELDOR Spinal Needle (CSEN International Ltd.) - The Perfect Needle for Spinal Anesthesia

<https://www.prnewswire.com/news-releases/eldor-spinal-needle-csen-international-ltd---the-perfect-needle-for-spinal-anesthesia-277579591.html>

A comparison of 25 gauge Quincke spinal needle with 26 gauge Eldor spinal needle for the elective Caesarian sections: Insertion characteristics and complications
<http://www.kumj.com.np/issue/4/263-266.pdf>

<http://www.csen.com/Eldor%20spinal%20needle%20-%20Nepal.pdf>

<http://www.csen.com/fda>

<http://www.csen.com/pencan.pdf>

<http://www.csen.com/benefits.htm>

<http://www.csen.com/anesthesia/newsletter.htm>

<http://www.csen.com/newf.pdf>

<http://www.csen.com/DMT.pdf>

<http://csen.com/GETTER%20-%20ELDOR%20products.pdf>

<http://www.csen.com/Focused.pdf>

<http://www.csen.com/tray.pdf>

<http://www.csen.com/anesthesia/journal.htm>

<http://www.csen.com/TNS.htm>

<http://www.csen.com/pic1.htm>

<http://www.csen.com/510K.pdf>

CSEN 2: - Eldor Spinal Needle 26G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 2.1: - Eldor Spinal Needle 25G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 2.2: - Eldor Spinal Needle 24G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 2.3: - Eldor Spinal Needle 22G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 2.4: - Eldor Spinal Needle 20G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 2.5: - Eldor Spinal Needle 25G, 125 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 2.6: - Eldor Spinal Needle 24G, 125 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 2.7: - Eldor Spinal Needle 22G, 125 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 2.8: - Eldor Spinal Needle 20G, 125 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 3: - Eldor Spinal Needle 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 4: - Eldor Spinal Needle 24G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 5: - Eldor Spinal Needle 25G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 5a: Eldor spinal needle, 25G, 105 mm, with an introducer, packed, sterilized, CE marked

CSEN 6: - Eldor Spinal Needle 26G, 125 mm, packed, sterilized, CE marked.

CSEN 7: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

<https://encrypted->

[tbn3.gstatic.com/images?q=tbn:ANd9GcTMgaZLuYGyJJIGeA1etyPS0JspcywzS_iUz9Sf_uelH1m5STz-odQ](https://encrypted-tbn3.gstatic.com/images?q=tbn:ANd9GcTMgaZLuYGyJJIGeA1etyPS0JspcywzS_iUz9Sf_uelH1m5STz-odQ)

<http://www.csen.com/csen7-1.JPG>

<http://www.csen.com/anesthesia/label-eldor-needle.htm>

<http://www.csen.com/anesthesia/sides.htm>

<http://www.csen.com/anesthesia/commadd.htm>

<http://www.csen.com/CSEK.pdf>

<http://www.csen.com/CSEA-pain.htm>

<http://www.csen.com/anesthesia/CSEA.htm>

<http://www.csen.com/anesthesia/pintore.htm>

<http://www.csen.com/comment2.pdf>

<http://www.csen.com/delayed.pdf>

CSEN 7a: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 90 mm, Eldor Spinal Needle 26G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

<http://csen.com/CSEN.png>

http://images.slideplayer.com/7/1728521/slides/slide_59.jpg

<http://www.csen.com/anesthesia/book/3.jpg>

CSEN 7A: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Assist Guide, LOR Syringe. Packed, sterilized, CE marked.

<http://www.csen.com/csen7A=1.JPG>

<http://www.csen.com/csen7A=.JPG>

CSEN 7B: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, CSEN Spinal

Needle 25G, 123 mm, Eldor Epidural Catheter, Adapter, Filter,
Thread Assist Guide, LOR Syringe. Packed, sterilized, CE marked.

CSEN 7C: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm; Eldor Spinal Needle 25G, 125 mm; Eldor Epidural Catheter; Adapter; Filter; Thread Assist Guide; LOR Syringe. Packed, sterilized, CE marked.

CSEN 7D: - Eldor Combined Spinal Epidural Kit: - Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit 120 mm; -CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 125 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2 µm, - Thread Assist Guide, - LOR (Loss of Resistance) Syringe 6 ml, packed, sterilized, CE marked.

<http://www.csen.com/CSEN%207D-1.jpg>

<http://www.csen.com/CSEN%207D-2.jpg>

CSEN 8: - Eldor Epidural Kit: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 9: - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Needle 18G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 10: - Eldor Insulated Plexus Block Needle 24G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 11: - Eldor Insulated Plexus Block Needle 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 12: - Eldor Continuous Plexus Block Kit: Tuohy Insulated Needle 18G, 90 mm, Eldor Plexus Block Catheter, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 13: - Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, packed, sterilized, CE marked.

<http://www.csen.com/csen7-1.JPG>

<http://www.csen.com/anesthesia/book/15.jpg>

<http://www.csen.com/pics/Aspirating%20Needles.jpg>

CSEN 14: - Eldor Spinal Needle 26G, 105 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 15: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 27G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 16 : - Eldor Epidural Catheter Kit: Eldor Epidural Catheter, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 17: - Eldor Spinal Needle Kit: Eldor Spinal Needle 26G, 105 mm, with an Introducer, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 18 : - Eldor Epidural Catheter Kit : Eldor Epidural Catheter, Adapter, Filter, packed, sterilized, CE marked.

CSEN 19: - Eldor CSENBR Kit: Tuohy Needle 18G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 19a: - Eldor CSENBR Kit: Tuohy Needle 18G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Epidural Catheter, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 20: - Eldor Spinal Needle 27G, 125 mm, packed, sterilized, CE marked.

CSEN 21: - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 22: - Eldor Spinal Needle 27G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 23: - Eldor Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

<http://csen.com/16G.JPG>

CSEN 24: - Tuohy epidural needle 18G, 90 mm, packed, sterilized, CE marked.

CSEN 25: - Tuohy epidural needle 16G, 90 mm, packed, sterilized, CE marked.

CSEN 26: - Tuohy epidural needle 18G, 120 mm, packed, sterilized, CE marked.

CSEN 26A: - Tuohy needle 18G, 120 mm.

CSEN 27: - Epidural catheter filter 0.2 micron, packed, sterilized, CE marked.

CSEN 28: - Loss of Resistance (LOR) syringe 10 cc, packed, sterilized, CE marked.

<http://www.csen.com/pics/Aspirating%20Needles.jpg>

CSEN 29: - Eldor "Hanging Drop Tube" 5 cm, L-L male, packed, sterilized, CE marked.

CSEN 30: - Epidural catheter thread assist guide, packed, sterilized, CE marked.

CSEN 31: - Epidural catheter adapter (Tuohy Borst), packed, sterilized, CE marked.

CSEN 32: Eldor Epidural Kit: Tuohy Needle 18G, 120 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor "Hanging Drop Tube", packed, sterilized, CE marked.

CSEN 33: Eldor Epidural Kit: Tuohy Needle 18G, 140 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor "Hanging Drop Tube", packed, sterilized, CE marked.

CSEN 34: Eldor Epidural Kit: Tuohy Needle 18G, 160 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor "Hanging Drop Tube", packed, sterilized, CE marked.

CSEN 35: - Eldor Spinal Needle 22G, 200 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 35.1: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 22G, 200 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

[/http://www.csen.com/eldor-amniocentesis-needle](http://www.csen.com/eldor-amniocentesis-needle)

<http://www.csen.com/eldor-amniocentesis-needle/eldor-directions2.htm>

CSEN 35.2: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 22G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.3: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 22G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.4: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 22G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/amnio.pdf>

CSEN 35.5: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 21G, 200 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.6: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 21G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.7: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 21G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.8: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 21G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.9: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 20G, 200 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.10: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 20G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.11: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 20G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 35.12: - Eldor Amniocentesis Needle (Double Hole Pencil Point Needle), Sharp tip, 20G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 36: - Eldor Spinal Needle 22G, 250 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 36.1: - Eldor Spinal Needle 22G, 200 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 36.2 : - Eldor Spinal Needle 22G, 150 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 37: - Eldor Spinal Needle 18G, 200 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 38: - Quincke Spinal Needle 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 38.1 (Z41643) : CSENCAN- Quincke Spinal Needle 22G, 90 mm, packed, sterilized, CE marked.

<http://www.csen.com/CSEN38.1.JPG>

<http://www.csen.com/38.1.JPG>

CSEN 38.2: - Quincke Peripheral Block Needle 22G, 50 mm. Packed, sterilized, CE marked.

CSEN 38.3 (Z41C43) : - Quincke Spinal Needle 18G, 90 mm, packed, sterilized, CE marked.

<http://www.csen.com/38.3.JPG>

CSEN 38.4 (Z41743) : - Quincke Spinal Needle 20G, 90 mm, packed, sterilized, CE marked.

<http://www.csen.com/38.4.JPG>

CSEN 38.5 (Z41613) : - Quincke Spinal Needle 22G, 38 mm, packed, sterilized, CE marked.

<http://www.csen.com/38.5.JPG>

CSEN 38.6 (Z41343) : - Quincke Spinal Needle 25G, 90 mm, with an introducer, packed, sterilized, CE marked.

CSEN 38.7 (Z41363) : - Quincke Spinal Needle 25G, 110 mm, with an introducer, packed, sterilized, CE marked.

<http://www.csen.com/38.7.JPG>

CSEN 38.8 (Z41373) : - Quincke Spinal Needle 25G, 120 mm, with an introducer, packed, sterilized, CE marked.

CSEN 38.9 (Z41143) : - Quincke Spinal Needle 27G, 90 mm, with an introducer, packed, sterilized, CE marked.

CSEN 39: - Quincke Spinal Needle 26G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/csen39.JPG>

CSEN 39a: Quincke spinal needle, 26G, 105 mm, with an introducer, packed, sterilized, CE marked.

CSEN 39b: Quincke spinal needle, 26G, 90 mm, packed, sterilized, CE marked.

CSEN 39c: CSEN Quincke spinal needle, 26G, 120 mm, packed, sterilized, CE marked.

CSEN 39d: CSEN Quincke spinal needle, 26G, 120 mm, with an introducer, packed, sterilized, CE marked.

CSEN 40: - Quincke Spinal Needle 27G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/csen40.JPG>

CSEN 40.1: - Quincke Spinal Needle 27G, 90 mm, packed, sterilized, CE marked.

CSEN 40a: Quincke spinal needle, 27G, 105 mm, with an introducer, packed, sterilized, CE marked.

CSEN 41: - Quincke Spinal Needle 29G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 42: - Tuohy Epidural Needle 18G, 110 mm, packed, sterilized, CE marked.

CSEN 43: - Tuohy Epidural Needle 18G, 140 mm, packed, sterilized, CE marked.

CSEN 44: - Epidural Catheter Kit : Epidural Catheter, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 44.1: CSEN Epidural Catheter Kit - CSEN Epidural Catheter 20G 3 lateral holes closed end , Tuohy Borst Adapter, Thread Assist Guide.

CSEN 45: - Tuohy Epidural Needle 18G, 160 mm, packed, sterilized, CE marked.

CSEN 45a: Tuohy epidural needle, 18G, 140 mm, packed, sterilized, CE marked.

CSEN 46: - Tuohy Epidural Needle 18G, 40 mm, packed, sterilized, CE marked.

CSEN 47: - Tuohy Epidural Needle 18G, 60 mm, packed, sterilized, CE marked.

CSEN 47.1: - **CSENPED** Tuohy Epidural Needle 18G, 64 mm, packed, sterilized, CE marked.

CSEN 48: - Tuohy Epidural Needle 18G, 80 mm, packed, sterilized, CE marked.

<http://www.csen.com/csen48=.jpg>

CSEN 49: - Eldor Spinal Needle Introducer 20G, 30 mm, packed, sterilized, CE marked.

CSEN 50:-Epidural Kit: Tuohy Needle 18G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 51:-Epidural Kit: Tuohy Needle 16G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 52:-Epidural Kit: Tuohy Needle 18G, 80 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 53:-Epidural Kit: Tuohy Needle 18G, 60 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 54:-Epidural Kit: Tuohy Needle 18G, 40 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 55: - Eldor Spinal Needle 26G, 115 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 56: - Combined Spinal Epidural Kit: Tuohy Needle 18G, 90 mm, Spinal Needle 26G, 125 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 56.1: - Combined Spinal Epidural Kit: Tuohy Needle 18G, 90 mm, CSEN Pencil Point Spinal Needle 26G, 125 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 57: - Pencil Point Spinal Needle 26G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN57-2.jpg>

<http://www.csen.com/CSEN57-2.jpg>

<http://www.csen.com/csen57-3.JPG>

CSEN 58: - Pencil Point Spinal Needle 27G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN58-1.jpg>

<http://www.csen.com/CSEN58-2.jpg>

<http://www.csen.com/CSEN58-3.jpg>

<http://www.csen.com/CSEN58.JPG>

<http://www.csen.com/CSEN58.pdf>

<http://www.csen.com/58-J.png>

<http://www.csen.com/58-J2.pdf>

<http://www.csen.com/58-J.png>

<http://www.csen.com/58-J1.png>

<http://www.csen.com/58.jpg>

<http://www.csen.com/csen58A.pdf>

<http://www.csen.com/CSEN58-3.jpg>

<http://www.csen.com/CSEN58-2.jpg>

<http://www.csen.com/CSEN58-1.jpg>

CSEN 58.1: - Pencil Point Spinal Needle 25G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 59 : - Epidural Catheter Kit : Epidural Catheter, Adapter, Filter, packed, sterilized, CE marked.

CSEN 60 : - Epidural Catheter Kit : Epidural Catheter, Adapter, packed, sterilized, CE marked.

CSEN 61: - Quincke Spinal Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 61.1: Quincke Spinal Needle, 25G, 90 mm, with an introducer, packed, sterilized, CE marked.

<http://www.csen.com/csen61-1.JPG>

CSEN 62: - Quincke Spinal Needle 25G, 90 mm, packed, sterilized, CE marked.

CSEN 62.1: Hypodermic needle, 25G, 90 mm, with an introducer, packed, sterilized, CE marked.

CSEN 62a: Quincke spinal needle, 25G, 105 mm, with an introducer, packed, sterilized, CE marked.

CSEN 62.2 : - Quincke Spinal Needle 27G, 90 mm, packed, sterilized, CE marked

CSEN 63: - Tuohy Epidural Needle 16G, 110 mm, packed, sterilized, CE marked.

CSEN 64: - Eldor Spinal Needle 26G, 130 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 65: - Eldor Spinal Needle 26G, 135 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 66: - Eldor Spinal Needle 26G, 140 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 67: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 68 : - **Eldor Epidural Kit:** Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter (6 lateral holes), Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

<http://www.csen.com/pics/DSC03560.JPG>

<http://www.csen.com/csen68.jpg>

<http://csen.com/F.png>

<http://www.csen.com/pics/Catheter%20Needles.JPG>

<http://www.csen.com/csen68=.jpg>

<http://www.csen.com/pics/DSC03553.JPG>

<http://www.csen.com/chir-csen68.pdf>

<http://www.csen.com/CSEN68.pdf>

<http://www.csen.com/68-J.png>

<http://www.csen.com/68-J2.pdf>

<http://www.csen.com/68-J.png>

<http://www.csen.com/68.jpg>

<http://www.csen.com/epidB.jpg>

<http://www.csen.com/eld6.pdf>

<http://www.csen.com/epidural-hematoma>

<http://www.csen.com/68-J2.pdf>

CSEN 68.1NRFit : Epidural Kit :--Epidural Tuohy needle 18G (OD: 1.2mm), effective length:90mm, overall length :118mm

The needle tube is marked with distance bands which are 10mm apart and originate from the hub.

--Epidural catheter 3 lateral holes OD: 0.8mm , \geq 900mm, RADIO-OPAQUE. The catheter is marked with graduations at 10mm intervals from 10mm to 250mm from the tip of patient end . Two bands indicate 100mm, three bands 150mm, four bands 200mm, five bands 250mm.

--Thread guide

--Catheter connector

--Liquid filter 0.22 μ m

--Lor syringe 10 ml

Fitted with connectors that comply with ISO 80369-6

CSEN 68.1 : Epidural Kit : --Epidural Tuohy needle 18G (OD: 1.2mm), effective length:90mm, overall length :118mm. The needle tube is marked with distance bands which are 10mm apart and originate from the hub.

--Epidural catheter 3 lateral holes OD: 0.8mm , \geq 900mm, RADIO-OPAQUE. The catheter is marked with graduations at 10mm intervals from 10mm to 250mm from the tip of patient end . Two bands indicate 100mm, three bands 150mm, four bands 200mm, five bands 250mm.

--Thread guide

--Catheter connector

--Liquid filter 0.22 μ m

--Lor syringe 10 ml

<http://www.csen.com/csen681A.pdf>

CSEN 68.1SNAP : Epidural kit

a. Epidural needle 18GX80MM

b. 0.8mm*900 Epidural anesthesia catheter

c. LOR syring (plastic) 10ml

d. 0.22um epidural filter (flat)

e. Catheter snap connector

1pc per hard plastic with tyvek paper, 10 pcs/clinical box, 120pcs/carton

G.W.: 11KGS

Master carton size: 650*408*420mm

<http://www.csen.com/SNAP.jpg>

CSEN 68.1METAL : Epidural kit

- a. Epidural needle 16GX80MM**
- b. Epidural anesthesia catheter with metal inside**
- c. LOR syring (plastic) 10ml**
- d. 0.22um epidural filter (flat)**
- e. Catheter connector**

1pc per hard plastic with tyvek paper, 10 pcs/clinical box, 120pcs/carton

G.W.: 11KGS

Master carton size: 650*408*420mm

<http://www.csen.com/metal1.jpg>

<http://www.csen.com/metal2.png>

<http://www.csen.com/metal3.png>

CSEN 68.1A: - CSEN Epidural Kit: CSEN Tuohy Epidural Needle 18G, 80 mm, CSEN Epidural Catheter 3 lateral holes , Tuohy Borst Adapter, Filter, Thread Assist Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 68.2 : - ELDOR CSEN Extra Epidural Kit : Eldor Tuohy Needle 16G, 200 mm, Eldor Epidural Catheter (6 lateral holes), Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 69: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 90 mm, - Eldor Spinal Needle 26G, 125 mm,- Eldor Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.

CSEN 69a: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2 µm, - CSEN Epidural Catheter Fixation Device, - Thread Assist Guide,- LOR (Loss of Resistance) Syringe 10 ml, packed, sterilized, CE marked.

<http://csen.com/K.png>

<http://www.csen.com/pics/Eldor%20CSEN%2069.jpg>

CSEN 69b: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2 µm, - Thread Assist Guide, - LOR (Loss of Resistance) Syringe 6 ml. Packed, sterilized, CE marked.

CSEN 69c: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, - Eldor Epidural Catheter 20G, 6 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2 µm, - Thread Assist Guide, - LOR (Loss of Resistance) Syringe 6 ml, packed, sterilized, CE marked.

CSEN 69.1: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - Eldor Spinal Needle 25G, 123 mm, - Eldor Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, - Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; - Syringes 2.5 ml/ 10 ml/ 20 ml; - Swab (3x); Clamp; - Compress 10 x 10 cm, 8 layers (3x); - Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm;- Tray for the products; - Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 69.2: CSEN Combined Spinal Epidural Kit: - Tuohy Epidural Needle 18G, 80 mm, metal stylet, - CSEN Pencil Point Spinal Needle 27G, 125 mm, - Epidural Catheter 20G, 3 lateral holes, closed end, - Tuohy-Borst Adapter, - Epidural Catheter Flat Filter, 0.2 µm, - Thread Assist Guide, - LOR Syringe 6 ml. Packed, sterilized, CE marked.

<http://csen.com/CSEN69=1.jpg>

<http://csen.com/CSEN69=1.jpg>

<http://csen.com/CSEN69=2.jpg>

<http://www.csen.com/ATA-csen%2069.2.pdf>

<http://www.csen.com/ATA-csen%2069.2.pdf>

<http://www.csen.com/csen692.pdf>

<http://www.csen.com/csen692ZR.jpg>

<http://www.csen.com/epsp.jpg>

<http://www.csen.com/ATA-csen%2069.2.pdf>

CSEN 69.3: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - Eldor Pencil Point Spinal Needle 25G, 125 mm, - Eldor six holes Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.

CSEN 69.4: - Eldor/Tuohy Combined Spinal Epidural Kit: - Tuohy Needle 18G, 80 mm, - CSEN Pencil Point Spinal Needle 25G, 125 mm, - CSEN three holes Epidural Catheter, - Adapter, - Filter, - Thread Guide, - LOR Syringe, packed, sterilized, CE marked.

CSEN 69.5: CSEN Combined Spinal Epidural Kit: Tuohy Epidural Needle 18G, 80 mm, metal stylet, CSEN Pencil Point Spinal Needle 25G, 110 mm, Epidural Catheter 20G, 3 lateral holes, closed end, Tuohy-Borst Adapter, Epidural Catheter Flat Filter, 0.2 µm, Thread Assist Guide, LOR Syringe 10 ml. Packed, sterilized, CE marked.

<http://www.csen.com/CSEN69-5.pdf>

<http://www.csen.com/69.5-J.png>

<http://www.csen.com/69.5-J2.pdf>

<http://www.csen.com/69.5-J.png>

<http://www.csen.com/695.jpg>

<http://www.csen.com/CSEN69-5.pdf>

CSEN 69.5SNAP : Combined spinal & epidural kit

- a. Epidural needle 18GX80MM
- b. 0.8mm*900 Epidural anesthesia catheter
- c. LOR syring (plastic) 10ml
- d. 0.22um epidural filter (flat)
- e. Catheter snap connector
- f. spinal needle with pencil point 25GX110MM

1pc per hard plastic with tyvek paper, 10 pcs/clinical box, 120pcs/carton

G.W.: 11KGS

Master carton size: 650*408*420mm

<http://www.csen.com/SNAP.jpg>

CSEN 69.5METAL : Combined spinal & epidural kit

- a. Epidural needle 16GX80MM
- b. Epidural anesthesia catheter with metal inside
- c. LOR syring (plastic) 10ml
- d. 0.22um epidural filter (flat)
- e. Catheter connector
- f. spinal needle with pencil point 25GX110MM

1pc per hard plastic with tyvek paper, 10 pcs/clinical box, 120pcs/carton

G.W.: 11KGS

Master carton size: 650*408*420mm

<http://www.csen.com/metal1.jpg>

<http://www.csen.com/metal2.png>

<http://www.csen.com/metal3.png>

CSEN 70: - Eldor Combined Needles: Eldor Combined Spinal Epidural

Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle

26G, 125 mm, packed, sterilized, CE marked.

CSEN 71: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal

Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal

Needle 27G, 125 mm, Eldor Epidural Catheter, Adapter, Filter,

Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 72 : - Eldor Epidural Catheter Kit: Eldor Epidural Catheter, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 73: - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe,packed, sterilized, CE marked.

CSEN 73.1: - **Eldor/Tuohy Combined Spinal Epidural Kit**: Tuohy Needle 16G, 90 mm, CSEN Pencil Point Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe,packed, sterilized, CE marked.

CSEN 74: - Eldor Epidural Kit: Tuohy Needle 16G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 75: - Eldor Epidural Kit: Tuohy Needle 17G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

<http://www.csen.com/csen75.JPG>

CSEN 76: - Eldor Epidural Kit: Tuohy Needle 17G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 77: - Tuohy epidural needle 17G, 90 mm, packed, sterilized, CE marked.

CSEN 77.1 : - CSEN Tuohy epidural needle 17G, 90 mm, removable wing.

CSEN 78: - Tuohy epidural needle 17G, 120 mm, packed, sterilized, CE marked.

CSEN 78.1 : - Tuohy epidural needle 17G, 127 mm, packed, sterilized, CE marked.

CSEN 79: Eldor Epidural Kit: Tuohy Needle 18G, 120 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 80: Eldor Epidural Kit: Tuohy Needle 17G, 120 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 81: Eldor Epidural Kit: Tuohy Needle 17G, 120 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 82: Eldor Epidural Kit: Tuohy Needle 18G, 140 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 83: Eldor Epidural Kit: Tuohy Needle 17G, 140 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”, packed, sterilized, CE marked.

CSEN 84: Eldor Epidural Kit: Tuohy Needle 17G, 140 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 85: Eldor Epidural Kit: Tuohy Needle 18G, 160 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 86: - Quincke Spinal Needle 32G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 87: - Tuohy Epidural Needle 18G, 40 mm, Insulated, packed, sterilized, CE marked.

CSEN 88: - Tuohy Epidural Needle 18G, 60 mm, Insulated, packed, sterilized, CE marked.

CSEN 89: - Tuohy Epidural Needle 18G, 80 mm, Insulated, packed, sterilized, CE marked.

CSEN 90:-Epidural Kit: Tuohy Needle 17G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 90.1: CSEN Epidural Anesthesia Kit - Tuohy Epidural Needle 17G 80 mm, CSEN Epidural Catheter 20G 3 lateral holes closed end, Tuohy Borst Adapter, Epidural Catheter Filter, Thread Assist Guide, Loss of Resistance (LOR) Syringe 10 ml.

CSEN 91: - Epidural Catheter Kit : Epidural Catheter, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 92: - Eldor Spinal Needle 27G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 93: - Pencil Point Spinal Needle 25G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/csen101-14.jpg>

<http://www.csen.com/csen93=3.JPG>

<http://www.csen.com/csen93=.JPG>

<http://www.csen.com/CSEN93-1.jpg>

<http://www.csen.com/CSEN93-2.jpg>

<http://www.csen.com/CSEN93.JPG>

<http://www.csen.com/Spin25.jpg>

<http://www.csen.com/CSEN93-2.jpg>

<http://www.csen.com/CSEN93-1.jpg>

CSEN 94: - Pencil Point Spinal Needle 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN94-2.jpg>

<http://www.csen.com/CSEN94-3.jpg>

<http://www.csen.com/CSEN94-3.jpg>

<http://www.csen.com/CSEN94-2.jpg>

CSEN 94.1: - CSEN Pencil Point Spinal Needle 22G, 120 mm, Sharp tip with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN94.1.pdf>

CSEN 94.2: - CSEN Single Hole Pencil Point Spinal Needle 22G, 103 mm, with an Introducer, packed, sterilized.

CSEN 95:-Epidural Kit: Tuohy Needle 19G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 96:-Epidural Kit: Tuohy Needle 22G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 97: - Tuohy epidural needle 19G, 90 mm, packed, sterilized, CE marked.

CSEN 98: - Tuohy epidural needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 99: - Quincke Spinal Needle 20G, 90 mm, packed, sterilized, CE marked.

CSEN 99.1 (Z41740) : - Quincke Spinal Needle 20G, 90 mm.

CSEN 100: - Quincke Spinal Needle 18G, 90 mm, packed, sterilized, CE marked.

CSEN 101: - Quincke Spinal Needle 19G, 90 mm, packed, sterilized, CE marked.

<http://www.csen.com/csen101-14.jpg>

<http://www.csen.com/CSEN101.pdf>

CSEN 102: - Quincke Spinal Needle 22G, 38 mm, packed, sterilized, CE marked.

CSEN 103:-Pediatric Epidural Kit: Tuohy Needle 19G, 50 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 103.1 :- ELDOR CSEN Pediatric Epidural Kit: Tuohy Needle 18G, 50 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe. Packed, sterilized, CE marked.

CSEN 103.2 :- ELDOR CSEN Pediatric Epidural Kit: Tuohy Needle 20G, 50 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe. Packed, sterilized, CE marked.

CSEN 103.3:- ELDOR CSEN Pediatric Epidural Kit: Tuohy Needle 22G, 50 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe. Packed, sterilized, CE marked.

CSEN 104: - Pencil Point Spinal Needle 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 105:-Epidural Kit: Tuohy Needle 18G, 90 mm, Epidural Catheter, Adapter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 106 : - Epidural Catheter Kit: Epidural Catheter 24G, with a Wire, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 107 : - Epidural Catheter Kit: Epidural Catheter 24G, with a Wire, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 108:-Plexus Kit: Tuohy Needle 20G, 90 mm, Insulated, Epidural Catheter 24G, with a Wire, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 109:-Plexus Kit: Tuohy Needle 20G, 90 mm, Insulated, Epidural Catheter 24G, with a Wire, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 110:-Plexus Kit: Tuohy Needle 22G, 90 mm, Insulated, Epidural Catheter 26G, with a Wire, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 111:-Plexus Kit: Tuohy Needle 22G, 90 mm, Insulated, Epidural Catheter 26G, with a Wire, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 112 : - Epidural Catheter Kit: Epidural Catheter 26G, with a Wire, Adapter, Thread Guide, packed, sterilized, CE marked.

CSEN 113: - Tuohy epidural needle 24G, 90 mm, packed, sterilized, CE marked.

CSEN 114: - Tuohy needle 24G, 90 mm, Insulated, packed, sterilized, CE marked.

CSEN 115:-Epidural Kit: Tuohy Needle 20G, 150 mm, Epidural Catheter 24G, with a Wire, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 116:-Plexus Kit: Tuohy Needle 20G, 150 mm, Insulated, Epidural Catheter 24G, with a Wire, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 117: - Tuohy needle 20G, 150 mm, Insulated, packed, sterilized, CE marked.

CSEN 118: - Tuohy needle 20G, 150 mm, packed, sterilized, CE marked.

CSEN 119: - Polycarbonate Loss of Resistance (LOR) syringe, 10 ml, packed, sterilized, CE marked.

CSEN 120: - CSEN Spinal/Epidural Tray - Polycarbonate LOR syringe 10 ml; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x);Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 121: - CSEN Spinal/Epidural Tray - LOR syringe 10 ml; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x);Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 122: - CSEN Spinal/Epidural Tray - Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x);Cannulas OD 0.5

mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 123: - Eldor Spinal Tray: Eldor Spinal Needle 26G, 90 mm, with an Introducer; Eldor "Hanging Drop Tube" 5 cm, L-L male; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm;

Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 124: - Eldor Spinal Tray: Eldor Spinal Needle 26G, 90 mm, with an Introducer, Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 124.1: - Eldor Spinal Tray: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with an Introducer, Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 124.2: - Eldor Spinal Tray: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer, Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD

1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 125: - Eldor Epidural Tray : Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 126: - Eldor Epidural Tray: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 127: - Eldor Combined Spinal Epidural Tray: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, Eldor “Hanging Drop Tube”; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 128: - Eldor Combined Spinal Epidural Tray: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 26G, 125 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 129: - Eldor Spinal Tray: Eldor Spinal Needle 26G, 90 mm; Introducer 20G, 30 mm; Syringe 5 ml; Swab (3x); Cup OD 11cm x 2.5 cm; Cannulas 25G OD 0.5 mm x L 16 mm/ 21G OD 0.8 mm x L 40 mm/ 18G OD 1.2 mm x L 40 mm; Tray 135 x 260 cm for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 130: - Eldor Spinal Tray: Syringe 5 ml; Swab (3x); Cup OD 11cm x 2.5 cm; Cannulas 25G OD 0.5 mm x L 16 mm/ 21G OD 0.8 mm x L 40 mm/ 18G OD 1.2 mm x L 40 mm; Tray 135 x 260 cm for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 131: Nerve Stimulator : Complete stimulator with accessories (battery, cable for connection to needles and monopolar catheters, user's manual), CE marked.

CSEN 132: Cable for connection to needles and monopolar catheters, packed, sterilized, CE marked.

CSEN 133: Re-sterilisable metal cape, packed, sterilized, CE marked.

CSEN 134: Cable for connection to bipolar catheters, packed, sterilized, CE marked.

CSEN 135: Clip for stimulating catheters, packed, sterilized, CE marked.

CSEN 136: Atraumatic Pencil Point stimulation needle 24G x 25 mm, packed, sterilized, CE marked.

CSEN 137: Atraumatic Pencil Point stimulation needle 22G x 50 mm, packed, sterilized, CE marked.

CSEN 138: Atraumatic Pencil Point stimulation needle 22G x 70 mm, packed, sterilized, CE marked.

CSEN 139: Atraumatic Pencil Point stimulation needle 21G x 90 mm, packed, sterilized, CE marked.

CSEN 140: Atraumatic Pencil Point stimulation needle 21G x 150 mm, packed, sterilized, CE marked.

CSEN 141: Stimulation needles with long insulated bevel 20° 24G x 25 mm, packed, sterilized, CE marked.

CSEN 142: Stimulation needles with long insulated bevel 20° 24G x 35 mm, packed, sterilized, CE marked.

CSEN 143: Stimulation needles with long insulated bevel 20° 22G x 50 mm, packed, sterilized, CE marked.

CSEN 144: Stimulation needles with long insulated bevel 20° 22G x 70 mm, packed, sterilized, CE marked.

CSEN 145: Stimulation needles with long insulated bevel 20° 21G x 90 mm, packed, sterilized, CE marked.

CSEN 146: Stimulation needles with long insulated bevel 20° 21G x 120 mm, packed, sterilized, CE marked.

CSEN 147: Stimulation needles with long insulated bevel 20° 21G x 150 mm, packed, sterilized, CE marked.

CSEN 148: Stimulation needles with short bevel 45° 24G x 25 mm, packed, sterilized, CE marked.

CSEN 149: Stimulation needles with short bevel 45° 22G x 50 mm, packed, sterilized, CE marked.

CSEN 150: Stimulation needles with short bevel 45° 22G x 70 mm, packed, sterilized, CE marked.

CSEN 151: Stimulation needles with short bevel 45° 21G x 90 mm, packed, sterilized, CE marked.

CSEN 152: Stimulation needles with short bevel 45° 21G x 150 mm, packed, sterilized, CE marked.

CSEN 153: CSEN Set with Pencil Point stimulation needle N30, packed, sterilized, CE marked.

CSEN 154: CSEN Set with Pencil Point stimulation needle N50, packed, sterilized, CE marked.

CSEN 155: CSEN Set with Pencil Point stimulation needle N70, packed, sterilized, CE marked.

CSEN 156: CSEN Set with Pencil Point stimulation needle N90, packed, sterilized, CE marked.

CSEN 157: CSEN N30-T with IV canula, packed, sterilized, CE marked.

CSEN 158: CSEN N50-T with IV canula, packed, sterilized, CE marked.

CSEN 159: CSEN N50 plus with stimulating catheter, packed, sterilized, CE marked.

CSEN 160: CSEN N90 plus with stimulating catheter, packed, sterilized, CE marked.

CSEN 161: CSEN N50 biplus with bipolar stimulating catheter, packed, sterilized, CE marked.

CSEN 162: CSEN C with long bevel 20° C30, packed, sterilized, CE marked.

CSEN 163: CSEN C with long bevel 20° C50, packed, sterilized, CE marked.

CSEN 164: CSEN C with long bevel 20° C90, packed, sterilized, CE marked.

CSEN 165: CSEN C50 plus with stimulating catheter, packed, sterilized, CE marked.

CSEN 166: CSEN C90 plus with stimulating catheter, packed, sterilized, CE marked.

CSEN 167: CSEN C50 biplus with bipolar stimulating catheter, packed, sterilized, CE marked.

CSEN 168: CSEN C90 biplus with bipolar stimulating catheter, packed, sterilized, CE marked.

CSEN 169: CSEN C.K with long bevel 20° Catheter trough needle technic / without IV canula C30K, packed, sterilized, CE marked.

CSEN 170: CSEN C.K with long bevel 20° Catheter trough needle technic / Without IV canula C50K, packed, sterilized, CE marked.

CSEN 171: CSEN C.K with long bevel 20° Catheter trough needle technic / Without IV canula C70K, packed, sterilized, CE marked.

CSEN 172: CSEN C.K with long bevel 20° Catheter trough needle technic / Without IV canula C120K, packed, sterilized, CE marked.

CSEN 173: CSEN C.K with long bevel 20° Catheter trough needle technic /

Without IV canula C150K, packed, sterilized, CE marked.

CSEN 174: CSEN C.K with long bevel 20° Catheter trough needle technic /

Without IV canula C30K plus, packed, sterilized, CE marked.

CSEN 175: CSEN C.K with long bevel 20° Catheter trough needle technic /

Without IV canula C50K plus, packed, sterilized, CE marked.

CSEN 176: CSEN C150K plus with stimulating catheter, packed, sterilized, CE marked.

CSEN 177: CSEN C150K biplus with bipolar stimulating catheter, packed, sterilized, CE marked.

CSEN 178: CSEN M with short bevel 45° M50, packed, sterilized, CE marked.

CSEN 179: CSEN T with stimulating Tuohy needle T50, packed, sterilized, CE marked.

CSEN 180: CSEN T with stimulating Tuohy needle T90, packed, sterilized, CE marked.

CSEN 181: CSEN T with stimulating Tuohy needle T120, packed, sterilized, CE marked.

CSEN 182: CSEN T with stimulating Tuohy needle T150, packed, sterilized, CE marked.

CSEN 183: Set with 21G x 50 cm catheter only + filter + connector, packed, sterilized, CE marked.

CSEN 184: Catheter fixation, packed, sterilized, CE marked.

CSEN 185: - Eldor Epidural Tray: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 5 ml/ 10 ml; Swab

(3x); Clamp; Compress 10 x 10 cm, 8 layers ; Cannulas OD 0.4 mm x L 40 mm/ OD 1.2 mm x L 50 mm; Cup for NaCl; Tray for the products; Cover 70 x 90 cm, packed, sterilized, CE marked.

CSEN 186: - Epidural Tray: Tuohy Needle 18G, 90 mm, Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 5 ml/ 10 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers ; Cannulas OD 0.4 mm x L 40 mm/ OD 1.2 mm x L 50 mm; Cup for NaCl; Tray for the products; Cover 70 x 90 cm, packed, sterilized, CE marked.

CSEN 187: - Epidural Tray: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe; Syringes 2 ml 20 ml; Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm, packed, sterilized, CE marked.

CSEN 188: - Epidural Kit: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter, Adapter, Thread Guide; Cover 70 x 90 cm, packed, sterilized, CE marked.

CSEN 189: - Quincke Spinal Needle 25G, 120 mm, packed, sterilized, CE marked.

CSEN 189.1: - Quincke Spinal Needle 25G, 120 mm, with an introducer, packed, sterilized, CE marked.

CSEN 189.2: - Quincke Spinal Needle 25G, 125 mm, with an introducer, packed, sterilized, CE marked.

CSEN 189.3: - Quincke Spinal Needle 22G, 125 mm, packed, sterilized, CE marked.

CSEN 190: - Whitacre Spinal Needle 25G, 120 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 191: - Eldor Spinal Needle 25G, 120 mm, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN191.pdf>

CSEN 192: - Eldor Spinal Needle 27G, 120 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 193: - Eldor Spinal Tray: Eldor Spinal Needle 26G, 90 mm, with an Introducer, Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 5 ml 2-component / 5 ml 3-component; Filtercannulla, Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (2x); Cannulas OD 0.4 mm x L 40 mm/ OD 1.2 mm x L 50 mm; Cup for NaCl; Tray for the products; Cover 70 x 90 cm, packed, sterilized, CE marked.

CSEN 194: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 80 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

<http://csen.com/Postoperative%20Wound%20Analgesia.jpg>

<http://www.admedus.com/au/wp-content/uploads/sites/4/2016/04/Eldor-Brochure-160421.pdf>

<http://www.csen.com/csen194.jpg>

<http://www.csen.com/CSEN194-2017.jpg>

<http://www.csen.com/ch88.jpg>

<http://www.csen.com/194b.jpg>

<http://www.csen.com/194c.jpg>

<http://www.csen.com/194a.jpg>

<http://www.csen.com/aus6.JPG>

<http://www.csen.com/aus7.JPG>

<http://www.csen.com/aus8.JPG>

<http://www.csen.com/aus10.JPG>

<http://www.csen.com/aus11.JPG>

<http://www.csen.com/194-J.jpg>

<http://www.csen.com/194.jpg>

<http://www.csen.com/194IFU.pdf>

<http://www.csen.com/194LTR1.JPG>

<http://www.csen.com/194LTR2.JPG>

<http://www.csen.com/BTC194.pdf>

<http://www.csen.com/BTC.pdf>

<http://www.csen.com/194ZS.pdf>

<http://www.csen.com/ZS194A.pdf>

<http://www.csen.com/ZS194C.pdf>

<http://www.csen.com/ZS194.pdf>

<http://www.csen.com/194IFU.pdf>

<http://www.csen.com/Admedus1.pdf>

CSEN 194.1: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 20 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 80 mm; Wound needle filter (blue); Wound catheter filter (white).). Packed, Sterilized, CE marked.

<http://www.csen.com/csen194.1.JPG>

CSEN 194.2: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

<https://www.youtube.com/watch?v=P1Z3fd4NpN4>

[http://www.csen.com/194.2\(8\).jpg](http://www.csen.com/194.2(8).jpg)

[http://www.csen.com/194.2\(5\).jpg](http://www.csen.com/194.2(5).jpg)

<http://www.csen.com/194.2.jpg>

<http://www.csen.com/wound3.JPG>

<http://www.csen.com/CSEN194-2-2017.jpg>

[http://www.csen.com/194.2\(4\).jpg](http://www.csen.com/194.2(4).jpg)

<http://www.csen.com/194.pdf>

<http://www.csen.com/aus1.JPG>

<http://www.csen.com/aus2.JPG>

<http://www.csen.com/aus3.JPG>

<http://www.csen.com/aus4.JPG>

<http://www.csen.com/aus5.JPG>

<http://www.csen.com/EL1.jpg>

<http://www.csen.com/EL2.jpg>

<http://www.csen.com/EL3.jpg>

<http://www.csen.com/EL4.jpg>

<http://www.csen.com/EL5.jpg>

<http://www.csen.com/EL6.jpg>

<http://www.csen.com/EL7.jpg>

<http://www.csen.com/EL8.jpg>

<http://www.csen.com/EL9.jpg>

<http://www.csen.com/led1.jpg>

<http://www.csen.com/led2.jpg>

<http://www.csen.com/led3.jpg>

<http://www.csen.com/led4.jpg>

<http://www.csen.com/OE.jpg>

<http://www.csen.com/1942LTR1.JPG>

<http://www.csen.com/1942LTR2.JPG>

<http://www.csen.com/1942BTC.pdf>

<http://www.csen.com/1942BTC2.pdf>

<http://www.csen.com/1942BTCADM.jpg>

<http://www.csen.com/BTC.pdf>

<http://www.csen.com/1942F1.jpg>

<http://www.csen.com/1942F2.jpg>

<http://www.csen.com/1942F3.jpg>

<http://www.csen.com/1942F4.jpg>

<http://www.csen.com/ZS1942.pdf>

<http://www.csen.com/ZS1942A.pdf>

<http://www.csen.com/ZS1942B.pdf>

<http://www.csen.com/1942S.pdf>

<http://www.csen.com/BTC.pdf>

CSEN 194.2A: - ELDOR Peripheral Nerve Block kit : Eldor Peripheral Nerve block Catheter 20G, 3 lateral holes, closed end, OD 1 mm x L 1000 mm; Tuohy-Borst Adapter ; Thread Assist Guide; Peripheral nerve block catheter Flat Filter, 0.2 µm; Eldor Peripheral Nerve Block needle - 16G x 200 mm, Echogenic ; Packed, Sterilized, CE marked.

CSEN 194.2B: - ELDOR Peripheral Nerve Block kit : Eldor Peripheral Nerve block Catheter 20G, end hole, OD 1 mm x L 1000 mm; Tuohy-Borst Adapter ; Thread Assist Guide; Peripheral nerve block catheter Flat Filter, 0.2 µm; Eldor Peripheral Nerve Block needle - 16G x 200 mm, Echogenic ; Packed, Sterilized, CE marked.

CSEN 194.2C: - ELDOR Ultrasound Needle : Eldor Ultrasound needle - 16G x 200 mm, Echogenic ; Flat Filter, 0.2 µm; Extension tube; Packed, Sterilized, CE marked.

CSEN 194.2D: - ELDOR Wound Catheter kit - OD 1 mm x L 900 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

<http://www.csen.com/194-3.jpg>

CSEN 194.2 - Cat: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; In Bulk.

CSEN 194.3: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.4: - CSEN Wound needle - 16G x 200 mm. Packed, Sterilized, CE marked.

CSEN 194.5: - ELDOR Wound Catheter - OD 1 mm x L 1000 mm, 20 lateral holes, 10 mm apart, closed tip. Packed, Sterilized, CE marked.

CSEN 194.6: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 20 lateral holes, 5 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.6A: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 30 lateral holes, 5 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.7: - ELDOR Wound Catheter kit - OD 1 mm x L 1000 mm, 40 lateral holes, 5 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Wound needle filter (blue); Wound catheter filter (white). Packed, Sterilized, CE marked.

CSEN 194.8 : - ELDOR Wound Catheter kit - OD 1 mm x L 900 mm, 9 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 80 mm; Wound catheter filter. Packed, Sterilized, CE marked.

<http://www.csen.com/CSEN194-8.png>

CSEN 194.9 : - ELDOR Wound Catheter kit – OD 1 mm x L 1000 mm, 9 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle -16G x 80 mm; Wound needle flat filter (white).Packed, Sterilized, CE marked.

<http://www.csen.com/194-9.JPG>

<http://www.csen.com/194--9.jpg>

<http://www.csen.com/194-9F.JPG>

<http://www.csen.com/194-9FF.JPG>

CSEN 195: CSEN Disposable Infusion Pump 200 ml, 5 ml/hr. Packed,
Sterilized, CE marked.

<http://www.csen.com/pump1.jpg>

<http://www.csen.com/pics/pump-1.JPG>

<http://www.csen.com/pump.pdf>

CSEN 195.1: CSEN Disposable Infusion Pump 200 ml, 3 ml/hr. Packed,
Sterilized, CE marked.

CSEN 195.1A: CSEN Disposable Infusion Pump 200 ml, 4 ml/hr. Packed,
Sterilized, CE marked.

<http://www.csen.com/195-1A.JPG>

CSEN 195.2: CSEN Disposable Infusion Pump 200 ml, 2 ml/hr. Packed,
Sterilized, CE marked.

CSEN 195.3: CSEN Disposable Infusion Pump 100 ml, 5 ml/hr. Packed,
Sterilized, CE marked.

CSEN 195.4: CSEN Disposable Infusion Pump 100 ml, 3 ml/hr. Packed,
Sterilized, CE marked.

CSEN 195.5: CSEN Disposable Infusion Pump 100 ml, 2 ml/hr. Packed,
Sterilized, CE marked.

<http://www.csen.com/195-5.JPG>

CSEN 195.6: CSEN Disposable Infusion Pump 275 ml, 5 ml/hr. Packed,
Sterilized, CE marked.

<http://csen.com/CDIP-.jpg>

[http://www.csen.com/195.6%20\(2\).png](http://www.csen.com/195.6%20(2).png)

<http://www.csen.com/195.6.png>

<http://www.csen.com/DIP.jpg>

<http://www.csen.com/195.pdf>

<http://www.csen.com/195-6.JPG>

<http://www.csen.com/195--6.png>

<http://www.csen.com/195-6F.JPG>

<http://www.csen.com/195-6FF.JPG>

CSEN 195.7: CSEN Disposable Infusion Pump 275 ml, 3 ml/hr. Packed, Sterilized, CE marked.

CSEN 195.8: CSEN Disposable Infusion Pump 275 ml, 2 ml/hr. Packed, Sterilized, CE marked.

<http://www.csen.com/Eldor%20Pump.pdf>

CSEN 195.9: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 1 ml/hr, PCA Bolus 0.5 ml/30 min. Packed, Sterilized, CE marked.

CSEN 195.10: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 1 ml/hr, PCA Bolus 0.5 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.11: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 2 ml/hr, PCA Bolus 0.5 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.12: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 0.5 ml/hr, PCA Bolus 0.5 ml/8 min. Packed, Sterilized, CE marked.

CSEN 195.13: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 1 ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.14: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 2 ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.15: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 3 ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.16: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 4ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.17: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 5 ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.18: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 2 ml/hr, PCA Bolus 2 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.19: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 5 ml/hr, PCA Bolus 2 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.20: CSEN Disposable Infusion Pump and PCA Bolus 100 ml, 1 ml/hr, PCA Bolus 0.5 ml/30 min. Packed, Sterilized, CE marked.

CSEN 195.21: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 1 ml/hr, PCA Bolus 0.5 ml/30 min. Packed, Sterilized, CE marked.

CSEN 195.22: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 2 ml/hr, PCA Bolus 0.5 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.23: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 3 ml/hr, PCA Bolus 0.5 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.24: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 5 ml/hr, PCA Bolus 1 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.25: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 5 ml/hr, PCA Bolus 2 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.26: CSEN Disposable Infusion Pump and PCA Bolus 275 ml, 8 ml/hr, PCA Bolus 2 ml/15 min. Packed, Sterilized, CE marked.

CSEN 195.27 : CSEN Disposable Infusion Pump 275 ml, 4 ml/hr. Packed, Sterilized, CE marked.

CSEN 196: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 196.1: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 24G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 196.2: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 22G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 196.3: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 20G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 197: CSEN Epidural Catheter kit: CSEN Epidural Catheter 20G, L 90 cm, 3 lateral holes, closed end; Thread assist guide; Tuohy-Borst Adapter; packed, sterilized, CE marked.

CSEN 198: ELDOR Epidural Catheter kit: ELDOR Epidural Catheter 20G, L 90 cm, 6 lateral holes, closed end; Thread assist guide; Tuohy-Borst Adapter; packed, sterilized, CE marked.

<http://www.csen.com/pics/CSEN198a.JPG>

CSEN 199: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 200: CSEN Disposable Infusion Pump 100 ml, 2ml/hr.

CSEN 201 : Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with an Introducer, with Markings, packed, sterilized, CE marked.

CSEN 202: - Eldor Combined Spinal Epidural Kit: Eldor Combined Spinal Epidural Needle 18G/20G laser welded spinal conduit, 120 mm, Eldor Spinal Needle 25G, 123 mm, Eldor Epidural Catheter, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 203: - Eldor/Tuohy Combined Spinal Epidural Kit: Tuohy Epidural Needle 18G, 80 mm, metal stylet, CSEN Pencil Point Spinal Needle , Sharp tip, 25G, 123 mm, Eldor Epidural Catheter 20G, 6 lateral holes, closed end, Tuohy-Borst Adapter, Epidural Catheter Flat Filter, 0.2 µm, CSEN Epidural Catheter Fixation Device, Thread Assist Guide, LOR (Loss of Resistance) Syringe 10 ml, packed, sterilized, CE marked.

CSEN 204: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/CSEN204-3.jpg>

<http://csen.com/CSEN%20204.jpg>

<http://www.csen.com/IND4.jpg>

<http://www.csen.com/IND1.jpg>

<http://www.csen.com/CSEN204-1.jpg>

<http://www.csen.com/CSEN204-2.jpg>

<http://www.csen.com/CSEN204-3.jpg>

<http://www.csen.com/CSEN204-4.jpg>

https://youtu.be/_sjZalNdILc

<http://www.csen.com/EN2.jpg>

<http://www.csen.com/cardinal204.JPG>

<http://www.csen.com/CSEN204-4.jpg>

<http://www.csen.com/CSEN204-3.jpg>

<http://www.csen.com/CSEN204-2.jpg>

<http://www.csen.com/CSEN204-1.jpg>

CSEN 204.1: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 204.2: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized.

CSEN 204.3: Eldor Double Hole Pencil Point Stem Cell Needle, Sharp tip, 18G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 204.4: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 24G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized.

CSEN 204.5: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 20G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 205: - CSEN Single Hole Pencil Point Spinal Needle 25G, 190 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 206: Eldor Double Hole Pencil Point Spinal Needle 25G, 190 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 207: - CSEN Single Hole Pencil Point Spinal Needle 26G, 190 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 208: - CSEN Single Hole Pencil Point Spinal Needle 27G, 103 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 209: - CSEN Single Hole Pencil Point Spinal Needle 25G, 103 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 210: - Tuohy epidural needle 18G, 150 mm, metal stylet, packed, sterilized, CE marked.

CSEN 211: - Tuohy epidural needle 18G, 200 mm, metal stylet, packed, sterilized, CE marked.

CSEN 212: - Tuohy epidural needle 16G, 150 mm, metal stylet, packed, sterilized, CE marked.

CSEN 213: - Tuohy epidural needle 16G, 200 mm, metal stylet, packed, sterilized, CE marked.

CSEN 214: - Quincke Spinal Needle 25G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 214.1: - **Quincke Pediatric Spinal Needle** 25G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 214.2: - **Quincke Peripheral Block Needle** 25G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 215: - **Quincke Spinal Needle** 26G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 216: - **Quincke Spinal Needle** 27G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 217: - **CSEN Single Hole Pencil Point Spinal Needle** 24G, 103 mm, Sharp tip, with an Introducer, packed, sterilized, CE marked.

CSEN 217.1 : - **CSEN Pencil Point Spinal Needle** 24G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 218: - **CSEN Single Hole Pencil Point Spinal Needle** 25G, 103 mm, Sharp tip, with an Introducer, packed, sterilized, CE marked.

CSEN 219: - **CSEN Single Hole Pencil Point Spinal Needle** 26G, 103 mm, Sharp tip, with an Introducer, packed, sterilized, CE marked.

CSEN 219.1: - **CSEN Single Hole Pencil Point Spinal Needle** 26G, 103 mm, with an Introducer, packed, sterilized.

CSEN 220: - **CSEN Single Hole Pencil Point Spinal Needle** 27G, 103 mm, Sharp tip, with an Introducer, packed, sterilized, CE marked.

CSEN 221: - **CSEN Single Hole Pencil Point Spinal Needle** 25G, 90 mm, Sharp tip, with an Introducer 20 mm, packed, sterilized.

CSEN 222: - **CSEN Single Hole Pencil Point Spinal Needle** 27G, 90 mm, Sharp tip, with an Introducer 20 mm, packed, sterilized.

<http://www.csen.com/Spin27.jpg>

CSEN 223: - **CSEN Single Hole Pencil Point Spinal Needle** 25G, 90 mm, Sharp tip, packed, sterilized.

<http://www.csen.com/Spin25.jpg>

CSEN 224: - Tuohy epidural needle 20G, 90 mm, packed, sterilized, CE marked.

CSEN 224.1 : ELDOR CSEN Epidural Kit : - Tuohy Epidural Needle 20G, 90 mm, - Eldor End Hole Epidural Catheter 24G, 100 cm, - Tuohy-Borst Adapter, - Epidural Catheter Filter, - Thread Assist Guide, - LOR Syringe 6 ml. Packed, sterilized, CE marked.

CSEN 224.2: - **CSENEPI** -Tuohy epidural needle 20G, 100 mm, packed, sterilized, CE marked.

CSEN 225: ELDOR Continuous Peripheral Nerve Block kit for

Ultrasound - Eldor Catheter 20G, 6 lateral holes, closed end ; Catheter adapter (Tuohy Borst); Catheter thread assist guide; CSEN needle (Tuohy Huber) - 16G x 200 mm; Needle filter (blue); Catheter filter (white). Packed, Sterilized, CE marked.

CSEN 226: ELDOR CSEN Needle for Peripheral Nerve Block with

Ultrasound - Eldor CSEN needle (Tuohy Huber) - 16G x 200 mm; Needle filter. Packed, Sterilized, CE marked.

CSEN 227: CSEN Intraosseous Needle 16G. Packed, Sterilized, CE marked.

CSEN 227.1: CSEN Intraosseous Needle 18G. Packed, Sterilized, CE marked.

CSEN 227.2: CSEN Intraosseous Needle 14G. Packed, Sterilized, CE marked.

CSEN 227.3: CSEN Intraosseous Needle 12G. Packed, Sterilized, CE marked.

CSEN 227.4: CSEN Intraosseous Needle 11G. Packed, Sterilized, CE marked.

CSEN 228: CSEN Epidural Catheter Flat Filter, 0.2 µm. Packed, Sterilized, CE marked.

CSEN 229: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 229.1: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 24G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 229.2: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 22G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 229.3: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 20G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 229.4: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 25 mm, packed, sterilized, CE marked.

CSEN 230 : Eldor Atraumatic Double Hole Pencil Point Ultrasound Needle, Sharp tip, 16G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/USN.htm>

CSEN 230.1 : Eldor Atraumatic Double Hole Pencil Point Ultrasound Needle, Sharp tip, 18G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 230.2 : Eldor Atraumatic Double Hole Pencil Point Ultrasound Needle, Sharp tip, 20G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 230.3 : Eldor Atraumatic Double Hole Pencil Point Ultrasound

Needle, Sharp tip, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 231: Eldor Double Hole Pencil Point Needle, Sharp tip, 20G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 231.1: Eldor Double Hole Pencil Point Needle, Sharp tip, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 231.2: Eldor Double Hole Pencil Point Needle, Sharp tip, 24G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 232: Eldor Transforaminal Double Hole Pencil Point Needle, Sharp tip, 20G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 232.1: Eldor Transforaminal Double Hole Pencil Point Needle, Sharp tip, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 232.2: Eldor Transforaminal Double Hole Pencil Point Needle, Sharp tip, 24G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 233: - CSEN Biopsy needle 18g, packed, sterilized, CE marked.

CSEN 234: - CSEN Tracheal tube, ID 7 mm, with cuff, packed, sterilized, CE marked.

CSEN 235: - CSEN Tracheostomy tube, ID 6 mm, packed, sterilized, CE marked.

CSEN 236: - CSEN Reinforced Tracheal tube, ID 7.5 mm, packed, sterilized, CE marked.

CSEN 237: - CSEN Tracheal tube, ID 3 mm, without cuff, packed, sterilized, CE marked.

CSEN 238: - CSEN Ultrasound Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 239: - CSEN Ultrasound Needle 22G, 50 mm, packed, sterilized, CE marked.

CSEN 240: - CSEN Ultrasound Needle 18G, 90 mm, packed, sterilized, CE marked.

CSEN 241: - CSEN Ultrasound Needle 18G, 50 mm, packed, sterilized, CE marked.

CSEN 241.1: - CSEN Ultrasound Needle 18G, 60 mm, packed, sterilized, CE marked.

CSEN 242: - CSEN Ultrasound Kit: Tuohy Needle 18G, 90 mm, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

<http://www.admedus.com/au/wp-content/uploads/sites/4/2016/04/Eldor-Brochure-160421.pdf>

<http://www.csen.com/242-2.jpg>

<http://www.csen.com/242b.jpg>

<http://www.csen.com/242a.jpg>

<http://www.csen.com/242-J1.jpg>

<http://www.csen.com/242.jpg>

<http://www.csen.com/242IFU.pdf>

CSEN 242.1: - ELDOR CSEN Ultrasound Kit: Tuohy Needle 18G, 90 mm, ELDOR CSEN nerve block Catheter (6 lateral holes) , Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 243: - CSEN Ultrasound Kit: Tuohy Needle 18G, 90 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

<http://www.csen.com/243-2.jpg>

CSEN 244: - CSEN Ultrasound Kit: Tuohy Needle 18G, 50 mm, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 244.1: - CSEN Ultrasound Kit: Tuohy Needle 18G, 60 mm, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

<http://www.csen.com/244-2.jpg>

<http://www.csen.com/244-1a.jpg>

<http://www.csen.com/244-1b.jpg>

<http://www.csen.com/244.1-J.jpg>

CSEN 244.2: - ELDOR CSEN Ultrasound Kit: Tuohy Needle 18G, 60 mm, ELDOR CSEN nerve block Catheter (6 lateral holes), Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 245: - CSEN Ultrasound Kit: Tuohy Needle 18G, 50 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 245.1: - CSEN Ultrasound Kit: Tuohy Needle 18G, 60 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

<http://www.csen.com/245-2.jpg>

CSEN 246: - ELDOR CSEN Ultrasound Needle 16G, 90 mm, packed, sterilized, CE marked.

CSEN 247: - ELDOR CSEN Ultrasound Needle 16G, 50 mm, packed, sterilized, CE marked.

CSEN 248: - ELDOR CSEN Ultrasound Needle 16G, 150 mm, packed, sterilized, CE marked.

CSEN 249: - ELDOR CSEN Ultrasound Needle 16G, 200 mm, packed, sterilized, CE marked.

CSEN 250 : - Eldor Spinal Tray: - Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, - Introducer, - Lidocaine 1%, 5 ml, - Sensorcaine 0.75% w/ Dextrose 8.25%, 2 ml, - Syringe 3 cc Luer Slip, - Syringe 5 cc Luer Slip, - Hypo Needle 20g x 1.5", - Hypo Needle 25g x 1.5", - Clear 4" Fenestrated Drape, - Medicine Cup, 20z., - 4x4 Gauze, 8ply, - Povidone iodine, 1oz, - Povidone iodine swab stick triples. FDA approved.

CSEN 250.1 : - Eldor Spinal Tray: - Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, - Introducer, - Syringe 3 cc Luer Slip, - Syringe 5 cc Luer Slip, - Hypo Needle 20g x 1.5", - Hypo Needle 25g x 1.5", -

Clear 4" Fenestrated Drape, - Medicine Cup, 20z., - 4x4 Gauze, 8ply. FDA approved.

CSEN 251: - ELDOR CSEN Ultrasound Needle 22G, 70 mm, packed, sterilized, CE marked.

CSEN 252: - ELDOR CSEN Ultrasound Needle 20G, 70 mm, packed, sterilized, CE marked.

CSEN 253: - ELDOR CSEN Ultrasound Needle 18G, 70 mm, packed, sterilized, CE marked.

CSEN 254: - ELDOR CSEN Ultrasound Needle 16G, 70 mm, packed, sterilized, CE marked.

CSEN 255: - ELDOR CSEN / Tuohy Ultrasound Needle 22G, 70 mm, packed, sterilized, CE marked.

CSEN 256: - ELDOR CSEN / Tuohy Ultrasound Needle 20G, 70 mm, packed, sterilized, CE marked.

CSEN 257: - ELDOR CSEN / Tuohy Ultrasound Needle 18G, 70 mm, packed, sterilized, CE marked.

CSEN 258: - ELDOR CSEN / Tuohy Ultrasound Needle 16G, 70 mm, packed, sterilized, CE marked.

CSEN 259: - CSEN Ultrasound Kit: Quincke Needle 18G, 90 m, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 260: - CSEN Ultrasound Kit: Quincke Needle 18G, 90 m, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 261: - CSEN Ultrasound Kit: Quincke Needle 18G, 70 m, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 262: - CSEN Ultrasound Kit: Quincke Needle 18G, 70 m, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 263: - CSEN Ultrasound Kit: Quincke Needle 18G, 50 m, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 264: - CSEN Ultrasound Kit: Quincke Needle 18G, 50 m, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 265: - ELDOR CSEN / Quincke Ultrasound Needle 18G, 90 mm, packed, sterilized, CE marked.

CSEN 266: - ELDOR CSEN / Quincke Ultrasound Needle 18G, 70 mm, packed, sterilized, CE marked.

CSEN 267: - ELDOR CSEN / Quincke Ultrasound Needle 18G, 50 mm, packed, sterilized, CE marked.

CSEN 268: - ELDOR CSEN / Quincke Ultrasound Needle 16G, 90 mm, packed, sterilized, CE marked.

CSEN 269: - ELDOR CSEN / Quincke Ultrasound Needle 16G, 70 mm, packed, sterilized, CE marked.

CSEN 270: - ELDOR CSEN / Quincke Ultrasound Needle 16G, 50 mm, packed, sterilized, CE marked.

CSEN 271: - ELDOR CSEN / Quincke Ultrasound Needle 16G, 150 mm, packed, sterilized, CE marked.

CSEN 272: - ELDOR CSEN / Quincke Ultrasound Needle 16G, 200 mm, packed, sterilized, CE marked.

CSEN 273: ELDOR Amniocentesis kit : Tuohy Needle 18G, 90 mm; Eldor Amniocentesis Catheter 20G, 6 lateral holes, closed end; Tuohy-Borst Adapter; Thread Assist Guide; packed, sterilized, CE marked.

CSEN 274: ELDOR Amniocentesis kit : Tuohy Needle 18G, 120 mm; Eldor Amniocentesis Catheter 20G, 6 lateral holes, closed end; Tuohy-Borst Adapter; Thread Assist Guide; packed, sterilized, CE marked.

CSEN 275: ELDOR Amniocentesis kit : Tuohy Needle 18G, 150 mm; Eldor Amniocentesis Catheter 20G, 6 lateral holes, closed end; Tuohy-Borst Adapter; Thread Assist Guide; packed, sterilized, CE marked.

CSEN 276: ELDOR Amniocentesis kit : Tuohy Needle 18G, 200 mm; Eldor Amniocentesis Catheter 20G, 6 lateral holes, closed end; Tuohy-Borst Adapter; Thread Assist Guide; packed, sterilized, CE marked.

CSEN 277: ELDOR CSEN - Quincke Ultrasound Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 278: ELDOR CSEN - Quincke Ultrasound Needle 22G, 60 mm, packed, sterilized, CE marked.

CSEN 279: ELDOR CSEN - Pencil Point Single Hole Ultrasound Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 280: ELDOR CSEN - Pencil Point Single Hole Ultrasound Needle 22G, 60 mm, packed, sterilized, CE marked.

CSEN 281: ELDOR CSEN - Eldor Pencil Point Double Hole

Ultrasound Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 282: ELDOR CSEN - Eldor Pencil Point Double Hole Ultrasound

Needle 22G, 60 mm, packed, sterilized, CE marked.

CSEN 283: - CSEN Pencil Point Spinal Needle Single Hole 26G, 90 mm,
with an Introducer, packed, sterilized.

<http://www.csen.com/csen283-13.jpg>

CSEN 283.1: - CSEN Pencil Point Spinal Needle Single Hole 26G, 110
mm, with an Introducer, packed, sterilized.

CSEN 283.2: - CSEN Pencil Point Spinal Needle Single Hole 26G, 115
mm, with an Introducer, packed, sterilized.

CSEN 284: - CSEN Pencil Point Spinal Needle Single Hole 25G, 120 mm,
with an Introducer, packed, sterilized.

<http://www.csen.com/csen284-13.jpg>

CSEN 284.1: - CSEN Pencil Point Spinal Needle Single Hole 26G, 120 mm,
with an Introducer, packed, sterilized.

CSEN 285: - CSEN Pencil Point Spinal Needle Single Hole 19G, 90 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 286: - CSEN Pencil Point Spinal Needle Single Hole 19G, 100 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 287: - CSEN Pencil Point Spinal Needle Single Hole 19G, 110 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 288: - CSEN Pencil Point Spinal Needle Single Hole 21G, 90 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 289: - CSEN Pencil Point Spinal Needle Single Hole 21G, 100 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 290: - CSEN Pencil Point Spinal Needle Single Hole 21G, 110 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 291: - CSEN Pencil Point Spinal Needle Single Hole 22G, 90 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 292: - CSEN Pencil Point Spinal Needle Single Hole 22G, 100 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 293: - CSEN Pencil Point Spinal Needle Single Hole 22G, 110 mm,
with an Introducer, packed, sterilized, CE marked.

CSEN 294: - CSEN Quincke Spinal Needle 19G, 90 mm, packed,
sterilized, CE marked.

CSEN 295: - CSEN Quincke Spinal Needle 19G, 100 mm, packed,
sterilized, CE marked.

CSEN 296: - CSEN Quincke Spinal Needle 19G, 110 mm, packed,
sterilized, CE marked.

CSEN 297: - CSEN Quincke Spinal Needle 21G, 90 mm, packed,
sterilized, CE marked.

CSEN 298: - CSEN Quincke Spinal Needle 21G, 100 mm, packed,
sterilized, CE marked.

CSEN 299: - CSEN Quincke Spinal Needle 21G, 110 mm, packed,
sterilized, CE marked.

CSEN 300: - CSEN Quincke Spinal Needle 22G, 90 mm, packed,
sterilized, CE marked.

<http://www.csen.com/Eldor-Quincke.jpg>

CSEN 301: - CSEN Quincke Spinal Needle 22G, 100 mm, packed, sterilized, CE marked.

CSEN 302: - CSEN Quincke Spinal Needle 22G, 110 mm, packed, sterilized, CE marked.

CSEN 302.1 : - CSEN Quincke Spinal Needle 22G, 103 mm, packed, sterilized, CE marked.

CSEN 302.2 : - CSEN Quincke Spinal Needle 22G, 120 mm, packed, sterilized, CE marked.

CSEN 302.3 (Z41670) : - CSEN Quincke Spinal Needle 22G, 120 mm.

CSEN 303 : - ELDOR CSEN Continuous Spinal Kit for Elderly Patients :
CSEN Tuohy Epidural Needle 18G, 90 mm, CSEN Spinal Catheter 20G, 3 lateral holes , Tuohy Borst Adapter, Spinal Catheter Flat Filter, 0.2 µm, Thread Assist Guide, LOR (Loss of Resistance) Syringe 6 ml., packed, sterilized, CE marked.

CSEN 304 : - ELDOR CSEN Continuous Spinal Kit for Elderly Patients :
CSEN Tuohy Epidural Needle 18G, 90 mm, ELDOR Spinal Catheter 20G, 6 lateral holes , Tuohy Borst Adapter, Spinal Catheter Flat Filter, 0.2 µm, Thread Assist Guide, LOR (Loss of Resistance) Syringe 6 ml., packed, sterilized, CE marked.

CSEN 305 : - ELDOR CSEN Continuous Spinal Kit for Elderly Patients :
CSEN Tuohy Epidural Needle 18G, 90 mm, CSEN Spinal Catheter 20G, end hole, Tuohy Borst Adapter, Spinal Catheter Flat Filter, 0.2 µm, Thread Assist Guide, LOR (Loss of Resistance) Syringe 6 ml., packed, sterilized, CE marked.

CSEN 306: - CSEN Quincke Amniocentesis Needle 22G, 90 mm, packed, sterilized, CE marked.

CSEN 307: ELDOR Pencil Point Amniocentesis Needle, Double Hole, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 308: - CSEN Pencil Point Amniocentesis, Single Hole 22G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 309: - CSEN Pencil Point Amniocentesis Needle, Single Hole 22G, 120 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 310: - CSEN Quincke Ultrasound Needle 19G, 100 mm, packed, sterilized, CE marked.

CSEN 311: - CSEN Quincke Ultrasound Needle 19G, 50 mm, packed, sterilized, CE marked.

CSEN 312 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

<http://www.csen.com/ECBN.pdf>

CSEN 313 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 314 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 315 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 316 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 317 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 318 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 319 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 25G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 320 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 321 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 322 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 323 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 324 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 325 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 326 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 327 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 23G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 328 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 329 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 330 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 331 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 332 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 333 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 334 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 335 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 22G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 336 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 337 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 338 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 339 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 340 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 341 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 342 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 343 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 21G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 344 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 345 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 346 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 347 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 348 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 349 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 350 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 351 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 19G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 352 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 353 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 354 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 355 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 356 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 357 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 358 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 359 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 18G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 360 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 361 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 362 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 363 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 364 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 365 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 366 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 367 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 17G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 368 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 369 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 370 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 371 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 372 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 373 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 374 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 375 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 16G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 376 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 40 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 377 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 50 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 378 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 70 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 379 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 90 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 380 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 100 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 381 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 382 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 150 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 383 : Eldor Biopsy Needle (Double Hole Pencil Point Aspiration Needle) Sharp tip, 15G, 230 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 384: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 385: - Eldor DHPP Needle (Double Hole Pencil Point) 24G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 386: - Eldor DHPP Needle (Double Hole Pencil Point) 22G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 387: - Eldor DHPP Needle (Double Hole Pencil Point) 20G, 125 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 388: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 125
mm, with markings every 10 mm from tip, with an Introducer, packed,
sterilized, CE marked.

CSEN 389: - Eldor DHPP Needle (Double Hole Pencil Point) 24G, 125
mm, with markings every 10 mm from tip, with an Introducer, packed,
sterilized, CE marked.

CSEN 390: - Eldor DHPP Needle (Double Hole Pencil Point) 22G, 125
mm, with markings every 10 mm from tip, with an Introducer, packed,
sterilized, CE marked.

CSEN 391: - Eldor DHPP Needle (Double Hole Pencil Point) 20G, 125
mm, with markings every 10 mm from tip, with an Introducer, packed,
sterilized, CE marked.

CSEN 392 : - Eldor DHPP Needle (Double Hole Pencil Point) 22G, 90
mm, with an Introducer, packed, sterilized, CE marked.

CSEN 393: - Eldor DHPP Needle (Double Hole Pencil Point) 24G, 90
mm, with an Introducer, packed, sterilized, CE marked.

CSEN 394: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 90
mm, with an Introducer, packed, sterilized, CE marked.

CSEN 395: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 105
mm, with an introducer, packed, sterilized, CE marked

CSEN 396: - Eldor DHPP Needle (Double Hole Pencil Point) 22G, 200
mm, with an Introducer, packed, sterilized, CE marked.

CSEN 397: - Eldor DHPP Needle (Double Hole Pencil Point) 22G, 150
mm, with an Introducer, packed, sterilized, CE marked.

CSEN 398: - Eldor DHPP Needle (Double Hole Pencil Point) 18G, 200 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 399: - Eldor DHPP Needle (Double Hole Pencil Point) Introducer 20G, 30 mm, packed, sterilized, CE marked.

CSEN 400: - Eldor DHPP Needle (Double Hole Pencil Point) Tray: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with an Introducer, Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 401: - Eldor DHPP Needle (Double Hole Pencil Point) Tray: Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer , Fenestrated Cover 45 x 70 cm with hole 10 cm, edge adhesive tape; Syringes 2.5 ml/ 10 ml/ 20 ml; Swab (3x); Clamp; Compress 10 x 10 cm, 8 layers (3x); Cannulas OD 0.5 mm x L 16 mm/ OD 0.8 mm x L 40 mm/ OD 1.2 mm x L 40 mm; Tray for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 402: - Eldor DHPP Needle (Double Hole Pencil Point) Tray: Syringe 5 ml; Swab (3x); Cup OD 11cm x 2.5 cm; Cannulas 25G OD 0.5 mm x L 16 mm/ 21G OD 0.8 mm x L 40 mm/ 18G OD 1.2 mm x L 40 mm; Tray 135 x 260 cm for the products; Cover 70 x 90 cm; packed, sterilized, CE marked.

CSEN 403: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 120 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 404:- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 405: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 24G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 406: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 22G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 407: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 20G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 408:- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 103 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 409 :- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 103 mm, with an Introducer, with Markings, packed, sterilized, CE marked.

CSEN 410 :- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 411:- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 412:- Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized.

CSEN 413: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 24G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized.

CSEN 414: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 20G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 415: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 190 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 416: - Eldor DHPP Needle (Double Hole Pencil Point) 25G, 190 mm, Sharp tip, with an Introducer, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 417: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 418: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 24G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 419: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 22G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 420: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 20G, 50 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 421: - Eldor DHPP Needle (Double Hole Pencil Point) , Sharp tip, 25G, 25 mm, packed, sterilized, CE marked.

CSEN 422 : - Eldor DHPP Needle (Double Hole Pencil Point) Tray: -

Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, -

Introducer, - Lidocaine 1%, 5 ml, - Sensorcaine 0.75% w/ Dextrose 8.25%, 2 ml, - Syringe 3 cc Luer Slip, - Syringe 5 cc Luer Slip, - Hypo Needle 20g x 1.5", - Hypo Needle 25g x 1.5", - Clear 4" Fenestrated Drape, - Medicine Cup,

20z., - 4x4 Gauze, 8ply, - Povidone iodine, 1oz, - Povidone iodine swab stick triples. FDA approved.

CSEN 423: - Eldor DHPP Needle (Double Hole Pencil Point) Tray: - Eldor Double Hole Pencil Point Spinal Needle, Sharp tip, 25G, 103 mm, - Introducer, - Syringe 3 cc Luer Slip, - Syringe 5 cc Luer Slip, - Hypo Needle 20g x 1.5", - Hypo Needle 25g x 1.5", - Clear 4" Fenestrated Drape, - Medicine Cup, 20z., - 4x4 Gauze, 8ply. FDA approved.

CSEN 424 : Eldor Subcutaneous Continuous Drug Infusion Kit: Tuohy Needle 18G, 50 mm, Eldor subcutaneous infusion Catheter end hole, 20G 100 mm, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 425 : Eldor Subcutaneous Continuous Drug Infusion Kit: Tuohy Needle 18G, 50 mm, Eldor subcutaneous infusion Catheter 3 lateral holes, 20G 100 mm, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 426 : Eldor Subcutaneous Continuous Drug Infusion Kit: Tuohy Needle 18G, 50 mm, Eldor subcutaneous infusion Catheter 6 lateral holes, 20G 100 mm, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 427: - Eldor Epidural Hematoma Immediate Evacuation Kit: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter 6 lateral holes, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 428: - Eldor Epidural Hematoma Immediate Evacuation Kit: Tuohy Needle 18G, 90 mm, CSEN Epidural Catheter 3 lateral holes, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 429: - Eldor Epidural Hematoma Immediate Evacuation Kit: Tuohy Needle 18G, 90 mm, CSEN Epidural Catheter end hole, Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 431 : - Eldor Pericardiocentesis kit: CSEN Tuohy Needle 16G x 200 mm, packed, sterilized, CE marked.

CSEN 432 : - Eldor Pericardiocentesis kit: CSEN Tuohy Needle 17G x 200 mm, packed, sterilized, CE marked.

CSEN 433 : - Eldor Pericardiocentesis kit: CSEN Tuohy Needle 18G x 200 mm, packed, sterilized, CE marked.

CSEN 434 : - Eldor Pericardiocentesis kit: CSEN Tuohy Needle 19G x 200 mm, packed, sterilized, CE marked.

CSEN 435 : - Eldor Pericardiocentesis kit: CSEN Tuohy Needle 20G x 200 mm, packed, sterilized, CE marked.

CSEN 436 : - Eldor Epicardial ablation kit: CSEN Tuohy Needle 16G x 200 mm, packed, sterilized, CE marked.

CSEN 437 : - Eldor Epicardial ablation kit: CSEN Tuohy Needle 17G x 200 mm, packed, sterilized, CE marked.

CSEN 438 : - Eldor Epicardial ablation kit: CSEN Tuohy Needle 18G x 200 mm, packed, sterilized, CE marked.

CSEN 439 : - Eldor Epicardial ablation kit: CSEN Tuohy Needle 19G x 200 mm, packed, sterilized, CE marked.

CSEN 440 : - Eldor Epicardial ablation kit: CSEN Tuohy Needle 20G x 200 mm, packed, sterilized, CE marked.

CSEN 441 : - Eldor Epicardial pacing kit: CSEN Tuohy Needle 16G x 200 mm, packed, sterilized, CE marked

CSEN 442 : - Eldor Epicardial pacing kit: CSEN Tuohy Needle 17G x 200 mm, packed, sterilized, CE marked

CSEN 443 : - Eldor Epicardial pacing kit: CSEN Tuohy Needle 18G x 200 mm, packed, sterilized, CE marked

CSEN 444 : - Eldor Epicardial pacing kit: CSEN Tuohy Needle 19G x 200 mm, packed, sterilized, CE marked

CSEN 445 : - Eldor Epicardial pacing kit: CSEN Tuohy Needle 20G x 200 mm, packed, sterilized, CE marked

CSEN 446 : ELDOR Atraumatic Intracardiac Needle for immediate Adrenaline injection in Cardiac Arrest - Eldor Double Hole Pencil Point Needle 22G, 103 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.
<http://www.csen.com/intracardiac1.pdf>

CSEN 447. CSEN Epidural Catheter 19G, 90 cm, 3 lateral holes, closed tip, Thread Assist Guide, Tuohy Borst Adapter.
<http://www.csen.com/CSEN447.JPG>
<http://www.csen.com/CSEN447-2.JPG>
<http://www.csen.com/CSEN447A.JPG>

CSEN 448. CSEN Epidural Catheter 19G, 90 cm, end hole, Thread Assist Guide, Tuohy Borst Adapter.
<http://www.csen.com/CSEN448.JPG>

CSEN 449. CSEN Epidural Catheter 19G, 90 cm, 6 lateral holes, closed tip, Thread Assist Guide, Tuohy Borst Adapter.
<http://www.csen.com/CSEN449.JPG>
<http://www.csen.com/CSEN449A.JPG>

CSEN 450. CSEN Epidural kit - CSEN Tuohy epidural needle 17G, 90 mm, CSEN epidural catheter 19G, 90 cm, 3 lateral holes, closed tip, Thread Assist Guide, Tuohy Borst Adapter, LOR syringe, Filter 0.2 micron, packed, sterilized, CE marked .

CSEN 451. CSEN Epidural kit - CSEN Tuohy epidural needle 17G, 90 mm, CSEN epidural catheter 19G, 90 cm, end hole, Thread Assist Guide, Tuohy Borst Adapter, LOR syringe, Filter 0.2 micron, packed, sterilized, CE marked.

CSEN 452. CSEN Epidural kit - CSEN Tuohy epidural needle 17G, 90 mm, CSEN epidural catheter 19G, 90 cm, 6 lateral holes, closed tip, Thread Assist Guide, Tuohy Borst Adapter , LOR syringe, Filter 0.2 micron, packed, sterilized, CE marked.

CSEN 453 : CSEN Epidural Anaesthesia Kit – 17G x 3.5" Tuohy needle with Huber point, 18G x 1.5" Safety needle, 21G x 1.5" Safety needle, 25G x 5/8" Safety

needle, 3 cc Luer slip LOR glass syringe, 5 cc Luer slip plastic syringe, 50 ml Drug cup, 3" x 3" Gauze sponges x 6, 0.22 micron Epidural filter, Adhesive bandages, Needle stick pad, Towel/Underpad, Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

<http://csen.com/s1.jpg>

CSEN 454 : CSEN Epidural Anaesthesia Kit – 17G x 3.5" Tuohy needle with Huber point, 18G x 1.5" Safety needle, 21G x 1.5" Safety needle, 25G x 5/8" Safety needle, 5 cc Luer slip LOR glass syringe, 10 cc Luer slip plastic syringe, 20 cc Luer slip LOR glass syringe, Prep tray, 50 ml Drug cup, 3" x 3" Gauze sponges x 3, 4" x 4" Gauze sponges x 2, Protective shield, Towel/Underpad, Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

CSEN 455 : CSEN Epidural Anaesthesia Kit – 18G x 1.5" Short bevel needle, 25G x 1.5" needle, 3 cc Luer lock plastic syringe, 5 cc Luer slip LOR glass syringe with glass tip, 5 cc Luer lock plastic syringe, 20 cc Luer lock plastic syringe, Prep tray, 0.22 micron Epidural filter, 3" x 3" Gauze sponges x 7, Foam prep stick x 3, Towel/Underpad, Clear Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

CSEN 456 : CSEN Epidural catheter kit – 19G Epidural catheter closed end, 20G Tuohy Borst adapter, Thread assist guide, Epidural catheter label, packed, sterilized, CE marked.

CSEN 457 : CSEN Epidural Anaesthesia Kit – 17G x 3.5" Tuohy needle with Huber point, 18G x 1.5" short bevel needle, 20G x 1.5" Filter needle, 21G x 1.5" Needle, 25G x 1.5" Needle, 3cc Luer slip plastic syringe, 7 cc Luer slip LOR plastic syringe, 20 cc Luer slip plastic syringe, Prep tray, Foam prep stick, 3" x 3" Gauze sponges x 3, 4" x 4" Gauze sponges x 3, Towel/Underpad, Clear Fenestrated drape with patch, CSR wrap, packed, sterilized, CE marked.

CSEN 458 : CSEN Epidural Anaesthesia Kit – 22G x 3.5" Tuohy needle with Huber point, 18G x 1.5" short bevel needle, 25G x 1.5" Needle, 7 cc Luer slip LOR plastic syringe, 12 cc Luer lock plastic syringe, 12 cc Luer slip plastic syringe, Prep tray, Foam prep stick x 3, 4" x 4" Gauze sponges x 8, Towel/Underpad x 2, Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

CSEN 459 : CSEN Epidural Anaesthesia Kit – 22G x 3.5" Tuohy needle with Huber point, 18G x 1.5" short bevel needle, 20G x 1.5" Filter needle, 25G x 5/8" Needle, 3 cc Luer slip plastic syringe, 5 cc Luer slip plastic syringe, 5 cc Luer slip LOR glass syringe, 0.22 micron Epidural filter, 3" x 3" Gauze sponges x 3, 4" x 4" Gauze sponges x 2, Needle stick pad, Towel/Underpad x 2, Clear Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

CSEN 460 : CSEN Epidural Anaesthesia Kit – 17G x 3.5" Tuohy – Schliff needle, 18G x 1.5" short bevel needle, 25G x 1.5" Needle, 3 cc Luer slip plastic syringe, 5 cc Luer slip plastic syringe, 10 cc Luer slip LOR plastic syringe, 20 cc Luer slip plastic syringe, 3" x 3" Gauze sponges x 3, 19G Epidural catheter with closed end, Thread assist guide, Epidural catheter label, 0.22 micron Epidural filter, Catheter pad, Needle stick pad, Towel/Underpad, Clear Fenestrated drape, CSR wrap, packed, sterilized, CE marked.

CSEN 461 : CSEN Epidural Anaesthesia Kit – 16G x 3.5" Tuohy needle with Huber point, 18G x 1.5" Safety needle, 21G x 1.5" Safety needle, 25G x 1.5" Safety needle, 5 cc Luer lock plastic syringe, 10 cc Luer lock LOR glass syringe, 20 cc Luer lock plastic syringe, 19G Epidural catheter closed tip multi port, 18G / 19G catheter connector clamp, Epidural catheter label, Prep tray, Foam prep stick x 3, 3" x 3" Gauze sponges x5, Towel/Underpad, Fenestrated drape , CSR wrap, packed, sterilized, CE marked.

CSEN 462 : CSEN Epidural Anaesthesia Kit – 18GA X 1.5" SHORT BEVEL NEEDLE X2, 18GA X 1.5" BLUNT ASPIRATION FILTER NEEDLE X1, 25GA X 1.5" NEEDLE X1, 1CC LUER SLIP PLASTIC SYRINGE X1, 3CC LUER SLIP PLASTIC SYRINGE X1, 5CC LUER LOCK PLASTIC SYRINGE X1, 20CC LUER LOCK LOR GLASS SYRINGE X1, PREP TRAY X1, FOAM PREP STICK X3, 3" X 3 GAUZE SPONGES X8, TOWEL/UNDERPAD X2, CLEAR FENESTRATED DRAPE X1, CSR WRAP X1, PACKED, STERILIZED, CE MARKED.

CSEN 463 : CSEN Epidural Anaesthesia Kit –
18GA X 1.5" SHORT BEVEL NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
7 CC LUER SLIP LOR PLASTIC SYRINGE X1
12CC LUER LOCK PLASTIC SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
NEEDLE STICK PAD X1
0.22 UM EPIDURAL FILTER X1
EPIDURAL CATHETER LABEL X1
4" X 4" GAUZE SPONGES X5
FOAM PREP STICK X3
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 464 : CSEN Epidural Anaesthesia Kit –
17 GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1
18 GA X 1.5" BLUNT FILL NEEDLE X1
21 GA X 1.5" NEEDLE X1
25 GA X 1.5" NEEDLE X1
30 GA X 0.5" NEEDLE X1
5 CC LUER SLIP LOR GLASS SYRINGE X1
5 CC LUER LOCK PLASTIC SYRINGE X1
10 CC LUER SLIP LOR PLASTIC SYRINGE X1
20 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
3" X 3" GAUZE SPONGES X3
FOAM PREP STICK X3
TRANSPARENT DRESSING WITH LABEL X1
0.22 UM EPIDURAL FILTER X1

SAFETY NEEDLE X1
NEEDLE STICK PAD X1
CLEAR FENESTRATED DRAPE X1
TOWEL/UNDERPAD X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 465 : CSEN Epidural Anaesthesia Kit –
17GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1
18 GA X 1.5" BLUNT FILL NEEDLE X1
25 GA X 1.5" NEEDLE X1
30 GA X 0.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER SLIP LOR GLASS SYRINGE X1
20 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
3" X 3" GAUZE SPONGES X3
SAFETY NEEDLE X1
NEEDLE STICK PAD X1
TOWEL/UNDERPAD X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 466 : CSEN Epidural Anaesthesia Kit –
18GA X 3.5" TUOHY-SCHLIFF NEEDLE X1
22 GA X 1.5" SAFETY NEEDLE X1
25 GA X 1" SAFETY NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
10 CC LUER LOCK LOR PLASTIC SYRINGE X1
PREP TRAY X1
4" X 4" GAUZE SPONGES X4
FOAM PREP STICK X3
TOWEL/UNDERPAD X1
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 467 : CSEN Epidural Anaesthesia Kit –
17 GA X 3.5" TUOHY NEEDLE WITH HUSTEAD POINT X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
18 GA X 1.5" FILTER NEEDLE X1
22GA X 1.5" NEEDLE X1
25 GA X 1" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER SLIP LOR GLASS SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X1
3" X 3" GAUZE SPONGES X3
NEEDLE STICK PAD X1

TOWEL/UNDERPAD X1
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1

PACKED, STERILIZED, CE MARKED.

[http://intranet.tdmu.edu.ua/data/kafedra/internal/xirtop/classes_stud/en/stomat/p
tn/operative%20surgery/2/Employment%201.%20General%20and%20special%20sur
gical%20instruments,%20suture%20material.files/image004.jpg](http://intranet.tdmu.edu.ua/data/kafedra/internal/xirtop/classes_stud/en/stomat/p
tn/operative%20surgery/2/Employment%201.%20General%20and%20special%20sur
gical%20instruments,%20suture%20material.files/image004.jpg)

CSEN 468 : CSEN Epidural Anaesthesia Kit –

17 GA X 3.5" TUOHY NEEDLE X1

18 GA X 1.5" SHORT BEVEL NEEDLE X1

22 GA X 1.5" NEEDLE X1

25GA X 1.5" NEEDLE X1

3 CC LUER SLIP PLASTIC SYRINGE X1

5 CC LUER LOCK LOR GLASS SYRINGE X1

20 CC LUER LOCK PLASTIC SYRINGE X1

PREP TRAY X1

FOAM PREP STICK X3

4" X 4" GAUZE SPONGES X4

19 GA EPIDURAL CATHETER X1

18GA/19GA CATHETER CONNECTOR CLAMPING X1

EPIDURAL CATHETER LABEL X1

FILTER STRAW X1

SAFETY NEEDLE X1

CATHETER PAD X1

PROTECTIVE SHIELD X1

TRANSPARENT DRESSING WITH LABEL X1

CLEAR FENESTRATED DRAPE X1

ABSORBENT TOWEL X1

CSR WRAP X1

PACKED, STERILIZED, CE MARKED.

CSEN 469 : CSEN Spinal Epidural Anaesthesia kit –

17 GA X 3.5" TUOHY NEEDLE X1

27 GA X 5.5" SPINAL NEEDLE X1

18 GA X 1.5" FILTER NEEDLE X1

22 GA X 1.5" NEEDLE X1

25 GA X 5/8" NEEDLE X1

25 GA X 1.5" NEEDLE X1

1 CC LUER SLIP PLASTIC SYRINGE X1

3 CC LUER LOCK PLASTIC SYRINGE X1

3CC LUER SLIP PLASTIC SYRINGE X1

5 CC LUER LOCK LOR GLASS SYRINGE X1

PREP TRAY X1

4"X 4" GAUZE SPONGES X8

19 GA EPIDURAL CATHETER CLOSED TIP MULIT PORT X1

18GA/19GA CATHETER CONNECTOR CLAMPING X1

EPIDURAL CATHETER LABEL X1

CATHETER PAD X1

DOCKING SYSTEM X1

PIN PAD X1
0.22 UM EPIDURAL FILTER X1
CLEAR FENESTRATED DRAPE WITH PATCH X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 470 : CSEN Epidural Anaesthesia Kit –
18 GA X 3.5" TUOHY EPIDURAL NEEDLE X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
22 GA X 1.5" NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER SLIP PLASTIC SYRINGE X1
8 CC LUER SLIP LOR PLASTIC SRYINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
3" X 3" GAUZE SPONGES X4
20 GA EPIDURAL CATHETER CLOSED END X1
20 GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
NEEDLE STICK PAD IN SAFETY CUP X1
0.22 UM EPIDURAL FILTER X1
CLEAR FENESTRATED DRAPE WITH PATCH X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 471 : CSEN Epidural Anaesthesia Kit –
17 GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
25 GA X 5/8" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER LOCK PLASTIC SYRINGE X1
8 CC LUER SLIP LOR PLASTIC SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
3" X 3" GAUZE SPONGES X4
19 GA EPIDURAL CATHETER WITH CLOSED END X1
18GA/19GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
FILTER STRAW X1
CLEAR FENESTRATED DRAPE X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 472 : CSEN Epidural Anaesthesia Kit –
16 GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1

16 GA X 1.5" NEEDLE X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER LOCK LOR GLASS SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
2" X 2" GAUZE SPONGES X4
19 GA EPIDURAL CATHETER WITH CLOSED TIP MULTI PORT X1
18GA/19GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
CLEAR FENESTRATED DRAPE X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 473 : CSEN Epidural Anaesthesia Kit –
17 GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
22 GA X 1.5" NEEDLE X1
25 GA X 1.5" NEEDLE X1
27 GA X 1.25" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER LOCK LOR GLASS SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
4" X 4" GAUZE SPONGES X4
19 GA EPIDURAL CATHETER CLOSED TIP MULTI-PORT X1
18GA/19GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
0.22 UM EPIDURAL FILTER X1
FILTER STRAW X1
FEMALE LUER LOCK CONNECTOR X1
NEEDLE STICK PAD IN SAFETY CUP X1
CATHETER PAD X1
CLEAR FENESTRATED DRAPE X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 474 : CSEN Epidural Anaesthesia Kit –
18 GA X 3.5" TUOHY EPIDURAL NEEDLE X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
22 GA X 1.5" NEEDLE X1
25 GA X 5/8" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
8 CC LUER LOCK LOR PLASTIC SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
EPIDURAL TRAY X1

PREP TRAY X1
FOAM PREP STICK X3
4" X 4" GAUZE SPONGES X4
20 GA EPIDURAL CATHETER CLOSED END X1
20 GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
FILTER STRAW X1
CATHETERED PAD X1
CLEAR FENESTRATED DRAPE WITH PATCH X1
ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 475 : CSEN Epidural Anaesthesia Kit –
17 GA X 3.5" TUOHY NEEDLE WITH HUBER POINT X1
18 GA X 1.5" SHORT BEVEL NEEDLE X1
22 GA X 1.5" NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER SLIP PLASTIC SYRINGE X1
8 CC LUER SLIP LOR PLASTIC SYRINGE X1
20 CC LUER LOCK PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
3" X 3" GAUZE SPONGES X4
19 GA EPIDURAL CATHETER CLOSED TIP MULTI-PORT X1
18GA/19GA CATHETER CONNECTOR CLAMPING X1
EPIDURAL CATHETER LABEL X1
0.22 UM EPIDURAL FILTER X1
NEEDLE STICK PAD IN SAFETY CUP X1
CLEAR FENESTRATED DRAPE WITH PATCH X1
ABSORBENT TOWEL X1
TOWEL/UNDERPAD X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 476 : CSEN Spinal Anaesthesia Kit –
18 GA X 1.5" SHORT BEVEL NEEDLE X1
20 GA X 1.25" INTRODUCER NEEDLE X1
20 GA X 1.5" FILTER NEEDLE X1
25 GA X 1.5" NEEDLE X1
27 GA X 0.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
4" X 4" GAUZE SPONGES X5
NEEDLE STICK PAD X1
TOWEL/UNDERPAD X1
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1

PACKED, STERILIZED, CE MARKED.

CSEN 477 : CSEN Spinal Anaesthesia Kit –
18GA X 1.5" SHORT BEVEL NEEDLE X1
20 GA X 1.25" INTRODUCER NEEDLE X1
20 GA X 1.5" FILTER NEEDLE X1
25 GA X 5/8" NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
3 CC LUER SLIP PLASTIC SYRINGE X1
5 CC LUER LOCK PLASTIC SYRINGE X1
5 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
4" X 4" GAUZE SPONGES X4
NEEDLE STICK PAD X1
SAFETY NEEDLE X1
TOWEL/UNDERPAD X1
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 478 : CSEN Spinal Anaesthesia Kit –
25 GA X 3.5" SPINAL NEEDLE X1
20 GA X 1.25" INTRODUCER NEEDLE X1
20 GA X 1.5" FILTER NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER SLIP PLASTIC SYRINGE X1
5 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
4" X 4" GAUZE SPONGES X3
TOWEL/UNDERPAD X1
CLEAR FENESTRATED DRAPE X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 479 : CSEN Spinal Anaesthesia Kit –
18 GA X 1.5" SHORT BEVEL NEEDLE X1
20 GA X 1.25" INTRODUCER NEEDLE X1
18 GA X 1.5" FILTER NEEDLE X1
25 GA X 1.5" NEEDLE X1
3 CC LUER LOCK PLASTIC SYRINGE X1
3 CC LUER SLIP PLASTIC SYRINGE X1
5 CC LUER SLIP PLASTIC SYRINGE X1
PREP TRAY X1
FOAM PREP STICK X3
3" X 3" GAUZE SPONGES X4
NEEDLE STICK PAD IN SAFETY CUP X1
CLEAR FENESTRATED DRAPE X1

ABSORBENT TOWEL X1
CSR WRAP X1
PACKED, STERILIZED, CE MARKED.

CSEN 480 : CSEN Spinal Anaesthesia Kit –
25 GA X 3.5" WHITACRE SPINAL NEEDLE X1
20 GA X 1.25" INTRODUCER NEEDLE X1
19 GA X 1.5" FILTER NEEDLE X1
19 GA X 1.5" NEEDLE X1
5 ML LUER LOCK PLASTIC SYRINGE X1
3 ML LUER LOCK PLASTIC SYRINGE X1
1 ML LUER SLIP PLASTIC SYRINGE WITH 25GA X 5/8" NEEDLE X1
NEEDLE STICK PAD X1
CLEAR FENESTRATED DRAPE X1
SPONGE APPLICATORS X3
TOWEL X1
GAUZE SPONGES X3
PACKED, STERILIZED, CE MARKED.

CSEN 481 : CSEN Epidural Anaesthesia Kit –
19 GA X 35.5" EPIDURAL CATHETER WITH THREAD ASSISTING
DEVICE X1
17 GA X 3.5" EPIDURAL NEEDLE X1
18 GA X 1.5" FILTER NEEDLE X1
18 GA X 1.5" INJECTION NEEDLE X1
22 GA X 1.5" INJECTION NEEDLE X1
25 GA X 1.5" INJECTION NEEDLE X1
10 ML LUER SLIP LOR SYRINGE X1
20 ML LUER SLIP SYRINGE X1
3 ML LUER SLIP SYRINGE X1
5 ML LUER LOCK LOR SYRINGE X1
25 GA X 5/8" INJECTION NEEDLE WITH 3ML LUER LOCK SYRINGE X1
CATHETER SYRINGE ADAPTER X1
CSR WRAP X1
24" X 36" DRAPE X1
24" X 36" DRAPE WITH 4" FENESTRATION AND ADHESIVE X1
0.2 MICRON FLAT FILTER X1
EPIDURAL PAD X1
4" X 4" GAUZE PADS X5
PREP SPONGE SWABS X3
MEDICINE CUP X1
PREP TRAY X1
PACKED, STERILIZED, CE MARKED.

CSEN 482: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block
Needle 18G, 90 mm, packed , sterilized, CE marked.

CSEN 483: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block
Needle 18G, 50 mm, packed , sterilized, CE marked.

CSEN 484: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 18G, 38 mm, packed , sterilized, CE marked.

CSEN 485: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 17G, 90 mm, packed , sterilized, CE marked.

CSEN 486: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 17G, 50 mm, packed , sterilized, CE marked.

CSEN 487: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 17G, 38 mm, packed , sterilized, CE marked.

CSEN 488: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 16G, 90 mm, packed , sterilized, CE marked.

CSEN 489: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 16G, 50 mm, packed , sterilized, CE marked.

CSEN 490: ELDOR CSEN Quincke Ultrasound Guided Peripheral Block Needle 16G, 38 mm, packed , sterilized, CE marked.

CSEN 491 (B5163-1) : 0.22 micron Epidural filter.

CSEN 492 (F361) : Loss of Resistance (LOR), Luer slip, plastic syringe, 10 cc.
<http://www.csen.com/CSEN492.JPG>

CSEN 493 : Plastic syringe , Luer lock, 12 cc.

CSEN 494 : Plastic syringe, Luer slip, 12 cc.

CSEN 495 (Z11840) : Tuohy needle with Huber point, 16G, 3.5".
<http://csen.com/G.JPG>
<http://www.csen.com/Z11840.JPG>

CSEN 496 (Z11B40) : Tuohy needle with Huber point, 17G, 3.5".

CSEN 497 (F141-2) : Epidural catheter, 19G, closed end.

CSEN 498 : Epidural catheter, 19G, closed tip, multi port (6 holes).

CSEN 499 : Loss of Resistance (LOR), Luer lock, Glass syringe, 20 cc.

CSEN 500 (F401-1) : Tuohy Borst Adapter (TBA), 20G.
<http://csen.com/D.JPG>
<http://www.csen.com/CSEN500.png>

CSEN 501 (Z11640) : Tuohy needle with Huber point, 22G, 3.5".

CSEN 502.1 (D214-2) : Loss of Resistance (LOR), Luer slip, Glass syringe, 2 ml.

CSEN 502.2 (D214-1) : Loss of Resistance (LOR), Luer slip, Glass syringe, 5 ml.

CSEN 503 (D214-1) : Loss of Resistance (LOR), Luer slip, Glass syringe, with Glass tip, 5 ml.

CSEN 504 (F361-1) : Loss of Resistance (LOR), Luer slip, Plastic syringe, 7 cc.
<http://www.csen.com/F361-1.png>

CSEN 505 (D307) : Adhesive Bandage.
<http://www.csen.com/D307.JPG>

CSEN 506 : Clear Fenestrated Drape with Patch.

CSEN 507 (F410-2) : Epidural catheter label.

CSEN 508 (G105) : Thread Assist Guide (TAG).

CSEN 509 (J170W-D201) : Loss of Resistance (LOR), Luer lock, Plastic syringe, 10 cc.

CSEN 510 (F106) : Epidural catheter, 20G, closed end,
<http://www.csen.com/F106.JPG>

CSEN 511 (D213) : Loss of Resistance (LOR), Luer lock, Plastic syringe, 5 cc.

CSEN 512 (J170W-D301) : Pin pad

CSEN 513 : Epidural catheter, 20G, single open end hole.

CSEN 514 : Tuohy Epidural needle, 18G, 1-3/4".

CSEN 515 : Crawford Epidural needle, 18G, 1-3/4".

CSEN 516 : Quincke spinal needle , 16G~27G, 38~200mm.

CSEN 517 : Pencil point spinal needle with single hole, 18G~27G, 90~120mm.

CSEN 518 : Epidural needle (15G/16G/17G/18G; 20G/22G).

CSEN 519 : Safety Vacuum Blood Collection Needle, pen like, without holder, 21G, 22G, 23G, visible type.

CSEN 520 : Safety Vacuum Blood Collection Needle, pen like, with holder, 21G, 22G, 23G, visible type.

CSEN 521 : Safety Vacuum Holder, bulking packing

CSEN 522 : Safety Holder with Vacuum Blood Collection Needle, pen like, blister packing

CSEN 523 : Tegaderm Transparent Dressing

CSEN 524 : Insulin Syringe 1ml, 30G.

CSEN 525 : Insulin Syringe 0.5ml, 30G.

CSEN 526 : Insulin Syringe 0.3ml, 30G.

CSEN 527 : Insulin Pen Needle 32G.

CSEN 528 : Safety Scalp Vein Set, DEHP-free material, 17" needle bevel,19G - 27G, 1/2" - 3/4" length.

CSEN 529 : Normal Scalp Vein Set 19G-27G, 3/4".

CSEN 530 : Three-way Stopcock without extension tube.

CSEN 531 : Three-way Stopcock with extension tube, 10cm-50cm.

CSEN 532 : Three-way Stopcock with extension tube, 80CM-100CM.

CSEN 533 : Three-way Stopcock with extension tube, 120CM-150CM.

CSEN 534 : Extension Tube 50CM-80CM.

CSEN 535 : Extension Tube 80CM-100CM.

CSEN 536 : Extension Tube 100CM-150CM.

CSEN 537 : Normal Hypodermic Needle 16G-20G.

CSEN 538 : Normal Hypodermic Needle 21G-25G.

CSEN 539 : Normal Hypodermic Needle 26G-27G.

CSEN 540 : Normal Hypodermic Needle 29G-30G.

CSEN 541 : Surgical Brush.

CSEN 542 : - ELDOR CSEN Epidural Kit : Tuohy Epidural needle 18G (0.8 mm) 80 mm, Eldor epidural catheter 6 lateral holes (0.8 mm), Epidural catheter adapter, LOR syringe 10 ml, Filter 0.22 micron, packed, sterilized , CE marked.

CSEN 543 : ELDOR CSEN Combined Spinal Epidural Kit : Tuohy Epidural needle 18G (0.8 mm) 80 mm, Pencil point spinal needle 25G 110 mm, Eldor epidural catheter 6 lateral holes (0.8 mm), Epidural catheter adapter, LOR syringe 10 ml, Filter 0.22 micron, packed, sterilized , CE marked.

CSEN 544 : ELDOR CSEN Quincke Spinal needle 26G 90 mm , packed, sterilized , CE marked.

<http://www.csen.com/544.jpg>

CSEN 545: - CSEN Pencil Point Spinal Needle Single Hole 27G, 110 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 546: – CSEN Quincke Spinal Needle 27G, 110 mm, with an Introducer, packed, sterilized , CE marked.

CSEN 547: - (Z11B40) :17GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.

<http://www.csen.com/Z11B40.png>

<http://www.csen.com/CSEN547.JPG>

<http://www.csen.com/CSEN547—.png>

CSEN 548: - (B5163-1) : 0.22UM EPIDURAL FILTER.

<http://www.csen.com/B5163-1.png>

<http://csen.com/J.png>

<http://www.csen.com/B5163-1.JPG>

<http://www.csen.com/CSEN548.JPG>

CSEN 549: - (F401-1) : Universal Tuohy Borst Adapter.

<http://www.csen.com/F401-1.JPG>

<http://www.csen.com/CSEN549.JPG>

<http://www.csen.com/CSEN549A.JPG>

CSEN 550: - (M6/019) : Epidural catheter, 19G closed end multi port (6 holes).

CSEN 551: - (T0/030) : ELDOR FILTER Tube , 5 micron filter, 10 cm tube.

<http://www.csen.com/ELDOR%20Filter%20Tube.png>

<http://www.csen.com/CSEN551.JPG>

<http://www.csen.com/CSEN551—.png>

CSEN 552: - (F141-2) : Epidural catheter, 19G closed end (3 holes).

<http://www.csen.com/CSEN552.JPG>

CSEN 553: - (Z11840) : 16GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.

<http://www.csen.com/Z11840.png>

<http://www.csen.com/CSEN553.JPG>

CSEN 554: - (F361-1) : 7CC LUER SLIP LOR PLASTIC SYRINGE.

<http://www.csen.com/F361.png>

<http://csen.com/C.png>

<http://www.csen.com/CSEN%20554.JPG>

<http://www.csen.com/CSEN554.png>

CSEN 555: - (G105) : Thread assist guide.

<http://www.csen.com/G105.png>

<http://www.csen.com/G105.JPG>

<http://www.csen.com/CSEN555.JPG>

CSEN 556: - (F361) : 10CC LUER SLIP LOR PLASTIC SYRINGE.

<http://www.csen.com/F361.png>

<http://www.csen.com/F361.JPG>

<http://www.csen.com/CSEN556.JPG>

CSEN 557: - (Z11640) : 22GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.

<http://www.csen.com/Z11640.png>

CSEN 558: - (M6/019+G105+F401-1) : Epidural Catheter 19G x 90 cm closed tip (6 holes), Thread Assist Guide , Tuohy Borst Adapter .

<http://www.csen.com/M6019+F401-1+G105.png>

CSEN 559: - (D307): ADHESIVE BANDAGES.

<http://www.csen.com/D307.png>

CSEN 560: - (J170W-D201): 10CC LUER LOCK LOR PLASTIC SYRINGE

<http://www.csen.com/J170W-D201.png>

CSEN 561: - (F106) : 20GA EPIDURAL CATHETER CLOSED END (3 holes).

<http://www.csen.com/F106%206PCS.png>

<http://www.csen.com/F106%2010PCS.png>

<http://www.csen.com/CSEN561.JPG>

CSEN 562: - (M6/020) : Epidural catheter, 19G open tip.

<http://www.csen.com/M6-020.png>

CSEN 563: - (J170W-D301): PIN PAD.

<http://www.csen.com/J170W-D301.png>

CSEN 564: - (Z11C40) : Epidural needle Tuohy 18G x 3.5".

<http://www.csen.com/Z11C40.png>

<http://csen.com/I.JPG>

CSEN 565: - (M6/021) : 20GA EPIDURAL CATHETER SINGLE OPEN END HOLE.

<http://www.csen.com/M6021.png>

CSEN 566: - (M6/021+F401-1+G105): 20GA EPIDURAL CATHETER SINGLE OPEN END HOLE , TUOHY BORST ADAPTER, THREAD ASSISTING DEVICE.

CSEN 567: - (M62120 + F401-1 + G105) : 19GA EPIDURAL CATHETER 2 SIDE HOLES WITH STYLET, OD 1 mm, TUOHY BORST ADAPTER, THREAD ASSISTING GUIDE.

<http://www.csen.com/M62120+F401-1+G105.png>

<http://csen.com/E.png>

<http://www.csen.com/CSEN567-3.JPG>

<http://www.csen.com/CSEN567.JPG>

<http://www.csen.com/CSEN567—.png>

CSEN 568. ELDOR Quincke Safety Spinal needle : 25G, 90 mm, 20% smaller hub opening

CSEN 569. ELDOR Safety Pencil point spinal needle (single hole): 25G, 90 mm, 20% smaller hub opening.

CSEN 570. ELDOR Safety spinal needle (double hole) : 25G, 90 mm, 20% smaller hub opening.

CSEN 571. ELDOR Safety Epidural needle : Tuohy needle 18G, 90 mm, Huber tip, 20% smaller hub opening.

CSEN 572. ELDOR Safety Peripheral Block needle : Tuohy needle 18G, 90 mm, Huber tip, 20% smaller hub opening.

CSEN 573. ELDOR Safety Wound Analgesia needle : Tuohy needle 18G, 90 mm, Huber tip, 20% smaller hub opening.

CSEN 574. ELDOR CSEN Safety Neuroaxial Syringe 10 ml, 20% smaller hub opening , Luer Lock.

CSEN 575. ELDOR CSEN Safety Neuroaxial Syringe 10 ml, 20% smaller hub opening , Non Luer Lock.

CSEN 576. ELDOR CSEN Safety Neuroaxial Syringe 5 ml, 20% smaller hub opening , Luer Lock.

CSEN 577. ELDOR CSEN Safety Neuroaxial Syringe 5 ml, 20% smaller hub opening , Non Luer Lock.

CSEN 578. ELDOR CSEN Safety Neuroaxial Syringe 3 ml, 20% smaller hub opening , Luer Lock.

CSEN 579. ELDOR CSEN Safety Neuroaxial Syringe 3 ml, 20% smaller hub opening , Non Luer Lock.

CSEN 580: ELDOR CSEN Safety LOR (Loss of Resistance) Syringe 10 ml, 20% smaller hub opening.

CSEN 581: ELDOR CSEN Safety Epidural Filter, 0.22 micron, 20% smaller hub openings

CSEN 582: ELDOR CSEN Safety Tuohy Borst Adapter, 20% smaller hub opening.

CSEN 583: ELDOR CSEN Safety Neuroaxial Infusion Set, 20% smaller hub opening.

CSEN 584: - ELDOR CSEN Ultrasound Kit: ECHOGENIC Tuohy Needle 16G, 90 mm, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 585: - ELDOR CSEN Ultrasound Kit: ECHOGENIC Tuohy Needle 16G, 90 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 586: - ELDOR CSEN Ultrasound Kit: ECHOGENIC Tuohy Needle 16G, 60 mm, CSEN nerve block Catheter, 3 lateral holes, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 587: - ELDOR CSEN Ultrasound Kit: ECHOGENIC Tuohy Needle 16G, 60 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 588: - ELDOR CSEN Ultrasound Needle: ECHOGENIC Tuohy Needle 16G, 90 mm, packed, sterilized, CE marked.

CSEN 589: - ELDOR CSEN Ultrasound Needle: ECHOGENIC Tuohy Needle 16G, 60 mm, packed, sterilized, CE marked.

CSEN 590 : - ELDOR Safety Ampoule Syringe Band.

CSEN 591 : (Z41640): Quincke spinal needle, 22 G x 90 mm. Sterile, single use.

<http://www.csen.com/Z41640.png>

<http://www.csen.com/CSEN591.JPG>

<http://www.csen.com/591B.JPG>

<http://www.csen.com/591.JPG>

<http://www.csen.com/CSEN591--.JPG>

<http://www.csen.com/591-J.png>

CSEN 592 : (Z41670): Quincke spinal needle, 22 G x 5 in. Sterile, single use.

<http://www.csen.com/CSEN592.JPG>

<http://www.csen.com/592B.JPG>

<http://www.csen.com/592-J.png>

CSEN 593 : (Z41620): Quincke spinal needle, 22G x 2.5 in. Sterile, single use.

CSEN 594 : (Z41340): Quincke spinal needle, 25G x 3.5 in. Sterile, single use.

<http://www.csen.com/CSEN594.JPG>

<http://www.csen.com/CSEN594A.JPG>
<http://www.csen.com/594B.JPG>
<http://www.csen.com/594.JPG>
<http://www.csen.com/CSEN594--.JPG>

CSEN 595 : (Z41740): Quincke spinal needle, 20 G x 3.5 in. Sterile, single use.
<http://www.csen.com/CSEN595.JPG>
<http://www.csen.com/595B.JPG>
<http://www.csen.com/595.JPG>
<http://www.csen.com/CSEN595-.JPG>

CSEN 596 : (Z41C40): Quincke spinal needle, 18 G x 3.5 in. Sterile, single use.
<http://www.csen.com/CSEN596.JPG>
<http://www.csen.com/596B.jpg>
<http://www.csen.com/596.JPG>
<http://www.csen.com/CSEN596-.JPG>

CSEN 597 : (Z4/023): Quincke spinal needle, 25G x 2 in. Sterile, single use.

CSEN 598 : (Z4/024): Quincke spinal needle, 22 G x 7 in. Sterile, single use.
<http://www.csen.com/CSEN598.JPG>
<http://www.csen.com/598B.jpg>
<http://www.csen.com/598.JPG>
<http://www.csen.com/CSEN598-.JPG>
<http://www.csen.com/598-J.png>

CSEN 599 : Tuohy 20G x 3.5in. Needle, Weiss style winged.
<http://www.csen.com/csen284-13.jpg>

CSEN 600 : (Z41540): Whitacre Spinal, 23 G x 3 1/2 in., Blue Hub. High-flow Pencil Point.

CSEN 601 : Quincke spinal introducer, yellow hub, 20 G x 1 1/4 in.
<http://www.csen.com/CSEN601.JPG>

CSEN 602 : (Z4/025): Quincke spinal needle, 18 G x 3 in. Sterile, single use.

CSEN 603 :
(Z8/012): Quincke Needle, 25G x 6in. with Plastic Hub, Calibrations and Needle Stop, Sterile.

CSEN 604 : (Z11840):
Needle, 16G x 3.5" Tuohy, Graduated with plastic hub and wings.
<http://www.csen.com/Z11840.png>

CSEN 605 : Tuohy needle 20G x 6in. Needle, Weiss style winged.

CSEN 606 : Quincke Needle, 22 G x 6 in. with Metal Hub and Stylet, Sterile.

CSEN 607: (F111-1) Epidural catheter with 3 holes JPG

<http://csen.com/A.JPG>

CSEN 608 : (F111-1+G105) Epidural catheter+Thread assist guide

<http://csen.com/B.png>

CSEN 609 : (Z51640) Spinal needle pencil point

<http://csen.com/H.JPG>

CSEN 610 : - Eldor Spinal Needle 22G, 120 mm, with markings every 10 mm from tip, packed, sterilized, CE marked.

CSEN 611 : - Eldor Spinal Needle 22G, 120 mm, with markings every 10 mm from tip, with an Introducer, packed, sterilized, CE marked.

CSEN 612 (M64120) : Epidural catheter with Soft tip, 1.0 x 900 mm.

<http://www.csen.com/CSEN612.png>

CSEN 613 (M64220) : Epidural catheter with Soft tip, 0.8 x 900 mm.

CSEN 614 : ELDOR CSEN Safety Adapter (Luer opening - 20% smaller hub opening) - Red.

CSEN 615 : ELDOR CSEN Safety Adapter (20% smaller hub opening - Luer opening) – Blue.

CSEN 616 (ZC1710): CSEN GWI (Guide Wire Insertion) Needle 20G x 1-1/2".

<http://www.csen.com/GWI1.JPG>

<http://www.csen.com/GWI2.JPG>

CSEN 617 (PAMA – NG – 25315) : Quincke spinal needle, 25G x 3.5 in. Sterile, single use.

<http://www.csen.com/NG25315.jpg>

CSEN 618 (ECK-1904) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, open tip, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK1904.pdf>

CSEN 619 (ECK-1906) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, closed tip, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK1906.pdf>

CSEN 620 (ECK-1908) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, closed tip multi, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK1908.pdf>

CSEN 621 (ECK-2004) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, open tip , Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2004.pdf>

CSEN 622 (ECK-2006) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, closed tip , Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2006.pdf>

CSEN 623 (ECK-2008) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, closed tip multi, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2008.pdf>

CSEN 624 : - NeedleEye : Epidural needle 18G, 120 mm combined by laser welding to ultrasound/fiberoptic Eye conduit 20G, 80 mm, packed, sterilized, CE marked.

<http://www.csen.com/csen7-1.JPG>

<http://www.csen.com/NEEDLEYE%20-6.pdf>

CSEN 625: – Aiguille Epidurale Tuohy 18G x 3.5" / Tuohy Epidural Needle 18G x 3.5".

<http://www.csen.com/CESEN625-2.JPG>

<http://www.csen.com/CESEN625.JPG>

CSEN 626: – Mini Epidural Kit : Tuohy needle 18G x 8 cm with separate wings, Epidural catheter open tip, Thread Assist Device, Tuohy Borst Adapter, Epidural Filter 0.22 um, LOR syringe 10 ml Luer slip, Packed.

<http://www.csen.com/CESEN626.JPG>

CSEN 627: – Mini Epidural Kit : Tuohy needle 18G x 8 cm with separate wings, Epidural catheter open tip, Thread Assist Device, Tuohy Borst Adapter, Epidural Filter 0.22 um, LOR syringe 10 ml Luer slip.

<http://www.csen.com/CSEN627.JPG>

CSEN 628 : - Epidural catheter 21G, closed tip, 3 lateral holes.

<http://www.csen.com/CSEN628.JPG>

CSEN 629 : - Quincke Spinal Needle 20G, 150 mm.

<http://www.csen.com/CSEN629.JPG>

CSEN 630 : - Removable Epidural Needle Wings.

<http://www.csen.com/RENEW.jpg>

CSEN 631 : CSEN Epidural Mini Kit : Aiguille Tuohy à pointe Huber 17G x 3½" / Tuohy needle with Huber 17G x 3½", Epidural Catheter 19G x 90 cm closed tip multi (6 lateral holes), Thread Assist Guide, Tuohy Borst Adapter, Seringue LOR 10ml Luer Slip / 10ml LOR Syringe Luer Slip, Filtre épidural 0.22µm / Epidural filter 0.22µm. Packed.

<http://www.csen.com/CSEN%20631.JPG>

CSEN 632 : CSEN Epidural Mini Kit : Aiguille Tuohy à pointe Huber 17G x 3½" / Tuohy needle with Huber 17G x 3½", Epidural Catheter 19G x 90 cm closed tip multi (6 lateral holes), Thread Assist Guide, Tuohy Borst Adapter, Seringue LOR 10ml Luer Slip / 10ml LOR Syringe Luer Slip, Filtre épidural 0.22µm / Epidural filter 0.22µm, Removable Epidural Needle Wings. Packed.

CSEN 633 : CSEN Epidural Mini Kit : Aiguille Tuohy à pointe Huber 17G x 3½" / Tuohy needle with Huber 17G x 3½", Adaptateur Tuohy Borst universel / Universal Tuohy Borst adapter, Seringue LOR 7ml Luer Slip / 7ml LOR Syringe Luer Slip, 21G EPIDURAL CATHETER CLOSED END 3 lateral holes, Guide à cathéter / Thread assist guide , Filtre épidural 0.22µm / Epidural filter 0.22µm, Epidural label (Yellow label written "Epidural Catheter" on it). Packed.

<http://www.csen.com/CSEN%20633.JPG>

CSEN 634 : -Introducer needle 20G x 1.25".

CSEN 635 : CSEN Acupuncture Needle : Size 0.25 x 40 mm , packed, sterilized, CE marked.

CSEN 636 : CSEN Acupuncture Needle : Size 0.25 x 25 mm , packed, sterilized, CE marked.

CSEN 637.1 : CE Sterile Non-absorbable Monofilament Nylon Surgical Sutures
<http://www.csen.com/SS4.pdf>

CSEN 638.1 (BN-Z) : Disposable biopsy needle intending for Bard Magnum Gun and Palium Gun
<http://www.csen.com/ECBN.pdf>

CSEN 639 : Sternal metal cable (suture with needle)
<http://www.csen.com/sternal1.jpg>
<http://www.csen.com/sternal2.jpg>

CSEN 640.1 (Barcode : 5070809050303) : CSEN Disposable Infusion Pump
100 ml, 2 ml/hr. Packed, Sterilized, CE marked

CSEN 640.2 (Barcode : 5090807050303) : CSEN Disposable Infusion Pump
100 ml, 3 ml/hr. Packed, Sterilized, CE marked

CSEN 640.3 (Barcode : 5070809070103) : CSEN Disposable Infusion Pump
100 ml, 4 ml/hr. Packed, Sterilized, CE marked

CSEN 640.4 (Barcode : 5070809040403) : CSEN Disposable Infusion Pump
100 ml, 5 ml/hr. Packed, Sterilized, CE marked

CSEN 640.5 (Barcode : 5070908050303) : CSEN Disposable Infusion Pump
100 ml, 6 ml/hr. Packed, Sterilized, CE marked

CSEN 640.6 (Barcode : 6080609050303) : CSEN Disposable Infusion Pump
100 ml, 7 ml/hr. Packed, Sterilized, CE marked

CSEN 640.7 (Barcode : 8040809050303) : CSEN Disposable Infusion Pump
100 ml, 8 ml/hr. Packed, Sterilized, CE marked

CSEN 640.8 (Barcode : 5020809050803) : CSEN Disposable Infusion Pump
100 ml, 9 ml/hr. Packed, Sterilized, CE marked

CSEN 640.9 (Barcode : 5070909050203) : CSEN Disposable Infusion Pump
100 ml, 10 ml/hr. Packed, Sterilized, CE marked

CSEN 640.10 (Barcode : 5060909050303) : CSEN Disposable Infusion
Pump 100 ml, 11 ml/hr. Packed, Sterilized, CE marked

CSEN 640.11 (Barcode : 5070409050703) : CSEN Disposable Infusion
Pump 100 ml, 12 ml/hr. Packed, Sterilized, CE marked

CSEN 640.12 (Barcode : 5070807050503) : CSEN Disposable Infusion
Pump 200 ml, 2 ml/hr. Packed, Sterilized, CE marked

CSEN 640.13 (Barcode : 5070609050503) : CSEN Disposable Infusion
Pump 200 ml, 3 ml/hr. Packed, Sterilized, CE marked

CSEN 640.14 (Barcode : 8070809020303) : CSEN Disposable Infusion
Pump 200 ml, 4 ml/hr. Packed, Sterilized, CE marked

CSEN 640.15 (Barcode : 8070806050303) : CSEN Disposable Infusion
Pump 200 ml, 5 ml/hr. Packed, Sterilized, CE marked

CSEN 640.16 (Barcode : 5050809070303) : CSEN Disposable Infusion
Pump 200 ml, 6 ml/hr. Packed, Sterilized, CE marked

CSEN 640.17 (Barcode : 5040809050603) : CSEN Disposable Infusion Pump
200 ml, 7 ml/hr. Packed, Sterilized, CE marked

CSEN 640.18 (Barcode : 9070809010303) : CSEN Disposable Infusion
Pump 200 ml, 8 ml/hr. Packed, Sterilized, CE marked

CSEN 640.19 (Barcode : 5090809050103) : CSEN Disposable Infusion
Pump 200 ml, 9 ml/hr. Packed, Sterilized, CE marked

CSEN 640.20 (Barcode : 6060809050303) : CSEN Disposable Infusion
Pump 200 ml, 10 ml/hr. Packed, Sterilized, CE marked

CSEN 640.21 (Barcode : 8070809020303) : CSEN Disposable Infusion
Pump 200 ml, 11 ml/hr. Packed, Sterilized, CE marked

CSEN 640.22 (Barcode : 5050809070303) : CSEN Disposable Infusion
Pump 200 ml, 12 ml/hr. Packed, Sterilized, CE marked

CSEN 640.23 (Barcode : 8070809020303) : CSEN Disposable Infusion
Pump 275 ml, 2 ml/hr. Packed, Sterilized, CE marked

CSEN 640.24 (Barcode : 6060809050303) : CSEN Disposable Infusion
Pump 275 ml, 3 ml/hr. Packed, Sterilized, CE marked

CSEN 640.25 (Barcode : 5070209050903) : CSEN Disposable Infusion
Pump 275 ml, 4 ml/hr. Packed, Sterilized, CE marked

CSEN 640.26 (Barcode : 5070808060303) : CSEN Disposable Infusion
Pump 275 ml, 5 ml/hr. Packed, Sterilized, CE marked

CSEN 640.27 (Barcode : 507080 9020603) : CSEN Disposable Infusion
Pump 275 ml, 6 ml/hr. Packed, Sterilized, CE marked

CSEN 640.28 (Barcode : 5080709050303) : CSEN Disposable Infusion
Pump 275 ml, 7 ml/hr. Packed, Sterilized, CE marked

CSEN 640.29 (Barcode : 5070609070303) : CSEN Disposable Infusion
Pump 275 ml, 8 ml/hr. Packed, Sterilized, CE marked

CSEN 640.30 (Barcode : 5090807050303) : CSEN Disposable Infusion
Pump 275 ml, 9 ml/hr. Packed, Sterilized, CE marked

CSEN 640.31 (Barcode : 7070609050303) : CSEN Disposable Infusion
Pump 275 ml, 10 ml/hr. Packed, Sterilized, CE marked

CSEN 640.32 (Barcode : 9030809050303) : CSEN Disposable Infusion
Pump 275 ml, 11 ml/hr. Packed, Sterilized, CE marked

CSEN 640.33 (Barcode : 3090809050303) : CSEN Disposable Infusion
Pump 275 ml, 12 ml/hr. Packed, Sterilized, CE marked

CSEN 641.1 : ELDOR CSEN Disposable Blood Pressure Cuff - New born
Cuff 25 x 5.5 cm

CSEN 641.2 : ELDOR CSEN Disposable Blood Pressure Cuff – Infant Cuff
30 x 7 cm

CSEN 641.3 : ELDOR CSEN Disposable Blood Pressure Cuff – Child Cuff 38
x 11 cm

CSEN 641.4 : ELDOR CSEN Disposable Blood Pressure Cuff – Standard
Cuff 52 x 14.5 cm

CSEN 641.5 : ELDOR CSEN Disposable Blood Pressure Cuff – Adult Cuff 56
x 14.5 cm

CSEN 641.6 : ELDOR CSEN Disposable Blood Pressure Cuff – Adult Cuff 54
x 14.5 cm

CSEN 641.7 : ELDOR CSEN Disposable Blood Pressure Cuff – Adult Cuff 52
x 14.5 cm

CSEN 641.8 : ELDOR CSEN Disposable Blood Pressure Cuff – Thigh Cuff
78 x 22 cm

CSEN 641.9 : ELDOR CSEN Disposable Blood Pressure Cuff – Cuff with metal hook 56 x 14.5 cm

CSEN 641.10 : ELDOR CSEN Disposable Blood Pressure Cuff – Cuff with metal chip 56 x 14.5 cm

CSEN 641.11 : ELDOR CSEN Disposable Blood Pressure Cuff – Nylon cuff with D-ring 56 x 14.5 cm

CSEN 641.12 : ELDOR CSEN Disposable Blood Pressure Cuff – Cotton cuff with D-ring 56 x 14.5 cm

CSEN 641.13 : ELDOR CSEN Disposable Blood Pressure Cuff – Large adult Cuff 66 x 17.5 cm

CSEN 641.14 : ELDOR CSEN Disposable Blood Pressure Cuff – Lengthen Cuff 100 x 14.5 cm

CSEN 642 : ELDOR CSEN Epidural Tray

<http://www.csen.com/TRAY1.JPG>

CSEN 643 : ELDOR CSEN Combined Spinal Epidural Tray

<http://www.csen.com/TRAY2.JPG>

CSEN 644: - Disposable Central Venous Catheter

CSEN 645: - Anti Infection Central Venous Catheter

CSEN 646: - Disposable Hemodialysis Catheter

CSEN 647: - Disposable Blood Pressure Transducer

CSEN 648: - Peripherally Inserted Central Catheter (PICC)

CSEN 649: - Artery Compression Device Sampler

CSEN 650: - Disposable General Anesthesia Catheter Kit

CSEN 651 - Disposable Infusion Pump - Introducer Sheath Kit

CSEN 652: - PTFE Coated Guide Wire

CSEN 653: - Y Valve set

CSEN 654: - Cardio Inflation Device

CSEN 655: - Angiography Control Syringe

CSEN 656: - Foley Catheter

CSEN 657: - Stomach Tube

CSEN 658: - HSG Catheter

CSEN 659: - Temperature Sensor Urethral Catheterization Kit

CSEN 660: - Cervical Ripening Balloon

CSEN 58NRFit: - Pencil Point Spinal Needle 27G, 90 mm, with an Introducer, packed, sterilized, CE marked.

CSEN 68NRFit: - Eldor Epidural Kit: Tuohy Needle 18G, 90 mm, Eldor Epidural Catheter (6 lateral holes), Adapter, Filter, Thread Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 68.1NRFit : - CSEN Epidural Kit: CSEN Tuohy Epidural Needle 18G, 80 mm, CSEN Epidural Catheter 3 lateral holes , Tuohy Borst Adapter, Filter, Thread Assist Guide, LOR Syringe, packed, sterilized, CE marked.

CSEN 543NRFit : ELDOR CSEN Combined Spinal Epidural Kit : Tuohy Epidural needle 18G (0.8 mm) 80 mm, Pencil point spinal needle 25G 110 mm, Eldor epidural catheter 6 lateral holes (0.8 mm), Epidural catheter adapter, LOR syringe 10 ml, Filter 0.22 micron, packed, sterilized , CE marked.

CSEN 69.5NRFit: CSEN Combined Spinal Epidural Kit: Tuohy Epidural Needle 18G, 80 mm, metal stylet, CSEN Pencil Point Spinal Needle 25G, 110 mm, Epidural Catheter 20G, 3 lateral holes, closed end, Tuohy-Borst Adapter, Epidural Catheter Flat Filter, 0.2 µm, Thread Assist Guide, LOR Syringe 10 ml. Packed, sterilized, CE marked.

CSEN 194.2NRFit: - ELDOR Wound Catheter kit : OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip; Wound catheter adapter (Tuohy Borst); Wound catheter thread assist guide; CSEN Wound needle - 16G x 200 mm; Filter 0.22 micron. Packed, Sterilized, CE marked.

CSEN 243NRFit: - CSEN Ultrasound Kit: Tuohy Needle 18G, 90 mm, CSEN nerve block Catheter, end hole, Adapter, Filter, Thread Guide, packed, sterilized, CE marked.

CSEN 194NRFit : - ELDOR Wound Catheter kit -

OD 1 mm x L 1000 mm, 10 lateral holes, 10 mm apart, closed tip;
Wound catheter adapter (Tuohy Borst);
Wound catheter thread assist guide;
CSEN Wound needle - 16G x 80 mm;
Wound catheter filter .
Packed, Sterilized, CE marked.

CSEN 242NRFit: - CSEN Ultrasound Kit:

Tuohy Needle 18G, 90 mm,
CSEN nerve block Catheter, 3 lateral holes,
Adapter,
Filter,
Thread Guide,
packed, sterilized, CE marked.

CSEN 243NRFit: - CSEN Ultrasound Kit:

Tuohy Needle 18G, 90 mm,
CSEN nerve block Catheter end hole with metal inside ,
Adapter,
Filter,
Thread Guide,
packed, sterilized, CE marked.

CSEN 661 : Mouth Mask

CSEN 662 : Face mask KN95 with EN149:2001. CE and FDA approved.

CSEN 662.1 : Face Mask Filter Respirator KN95 – GB2626-2006 CE
marked.

CSEN 663 : Face mask KN95 with GB (China standard). CE and FDA
approved.

CSEN 664 : Medical face mask . CE and FDA approved.

CSEN 665 : CSEN Epidural Kit :

-- CSEN Tuohy Epidural needle 18G (OD: 1.2mm), effective length:90mm, overall length :118mm. The needle tube is marked with distance bands which are 10mm apart and originate from the hub.

--Epidural catheter OD: 0.8mm , ≥900mm, RADIO-OPAQUE

The catheter is marked with graduations at 10mm intervals from 10mm to 250mm from the tip of patient end . Two bands indicate 100mm, three bands 150mm, four bands 200mm, five bands 250mm.

--Thread guide

--Catheter connector

--Liquid filter 0.22μ m

--Lor syringe 10 m

CSEN 666NRFit : CSEN Epidural Kit :

-- CSEN Tuohy Epidural needle 18G(OD: 1.2mm), effective length:90mm, overall length :118mm. The needle tube is marked with distance bands which are 10mm apart and originate from the hub.

--Epidural catheter OD: 0.8mm , ≥900mm, RADIO-OPAQUE

The catheter is marked with graduation at 10mm intervals from 10mm to 250mm from the tip of patient end . Two bands indicated 100mm, three bands 150mm, four bands 200mm, five bands 250mm.

--Thread guide

--Catheter connector

--Liquid filter 0.22µ m

--Lor syringe 10 ml

CSEN 667 : Blue Liquid Filter, packed, sterilized, CE marked.

<http://www.csen.com/BLF.jpg>

CSEN 668 : ELDOR PRODUCTS - TESTIMONIALS

<http://www.csen.com/testimonials.htm>

CSEN 669 : INTRALIPID IS A MAGIC BULLET

<https://www.amazon.com/dp/B07DX5LWGY>

<https://www.amazon.com/dp/B07B67TPH2>

<https://www.amazon.com/dp/B07FQPV55T>

<https://www.amazon.com/dp/B07CQ3LDCC>

<https://www.amazon.com/dp/B08PC8DQN2>

<https://www.amazon.com/dp/B07R4BNSS5>

<https://www.amazon.com/dp/B07B66FS74>

<https://www.amazon.com/dp/B078R1JB5B>

<https://www.amazon.com/dp/B079PYTLM8>

<https://www.amazon.com/dp/B0763S3CF5>

<https://www.amazon.com/dp/B078WCKJRP>

CSEN 670 : Browning Surgical, South Africa

<http://www.csen.com/CSEN38.1.pdf>

<http://www.csen.com/CSEN100.pdf>

<http://www.csen.com/CSEN101.pdf>

<http://www.csen.com/CSEN283.pdf>

<http://www.csen.com/CSEN284.pdf>

CSEN 671 : Chirmedical SNC, Messina, Italy

<http://www.csen.com/chir-csen1.pdf>

<http://www.csen.com/chir-csen2.pdf>

<http://www.csen.com/chir-csen3.pdf>

<http://www.csen.com/chir-csen4.pdf>

<http://www.csen.com/chir-csen5.pdf>

<http://www.csen.com/chir-csen7A.pdf>

<http://www.csen.com/chir-csen7B.pdf>

<http://www.csen.com/chir-csen57.pdf>

<http://www.csen.com/chir-csen58.pdf>

<http://www.csen.com/chir-csen68.pdf>

<http://www.csen.com/chir-csen69.pdf>

<http://www.csen.com/chir-csen74.pdf>

<http://www.csen.com/chir-csen75.pdf>

<http://www.csen.com/chir-csen93.pdf>

<http://www.csen.com/chir-csen94.pdf>

<http://www.csen.com/chir-csen94.1.pdf>

<http://www.csen.com/chir-csen191.pdf>

CSEN 672 : Orsan Medikal, Turkey

<http://www.csen.com/194.pdf>

<http://www.csen.com/195.pdf>

<http://www.csen.com/194.2.jpg>

[http://www.csen.com/194.2\(4\).jpg](http://www.csen.com/194.2(4).jpg)

[http://www.csen.com/194.2\(5\).jpg](http://www.csen.com/194.2(5).jpg)

<http://www.csen.com/195.6.png>

[http://www.csen.com/195.6%20\(2\).png](http://www.csen.com/195.6%20(2).png)

CSEN 673 : BTC Health, Australia

<https://www.btchealth.com.au/wp-content/uploads/2019/10/Eldor-CSEN-product-range.pdf>

<https://www.btchealth.com.au/wp-content/uploads/2019/05/INF0028-Eldor-Brochure-141217-1.pdf>

CSEN 674 : LTR Medical, Australia

<https://www.ltrmedical.com/eldor-wound-catheter-kits/>

CSEN 675 : PAMA Manufacturing, Canada

<http://www.csen.com/NG25315.jpg>

CSEN 492 (F361) : Loss of Resistance (LOR), Luer slip, plastic syringe, 10 cc.
<http://www.csen.com/CSEN492.JPG>

CSEN 495 (Z11840) : Tuohy needle with Huber point, 16G, 3.5".
<http://csen.com/G.JPG>
<http://www.csen.com/Z11840.JPG>

CSEN 500 (F401-1) : Tuohy Borst Adapter (TBA), 20G.
<http://csen.com/D.JPG>
<http://www.csen.com/CSEN500.png>

CSEN 504 (F361-1) : Loss of Resistance (LOR), Luer slip, Plastic syringe, 7 cc.
<http://www.csen.com/F361-1.png>

CSEN 505 (D307) : Adhesive Bandage.
<http://www.csen.com/D307.JPG>

CSEN 510 (F106) : Epidural catheter, 20G, closed end,
<http://www.csen.com/F106.JPG>

CSEN 544 : ELDOR CSEN Quincke Spinal needle 26G 90 mm , packed, sterilized ,
CE marked.
<http://www.csen.com/544.jpg>

CSEN 547: - (Z11B40) :17GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.
<http://www.csen.com/Z11B40.png>
<http://www.csen.com/CSEN547.JPG>
<http://www.csen.com/CSEN547—.png>

CSEN 548: - (B5163-1) : 0.22UM EPIDURAL FILTER.
<http://www.csen.com/B5163-1.png>
<http://csen.com/J.png>
<http://www.csen.com/B5163-1.JPG>
<http://www.csen.com/CSEN548.JPG>

CSEN 549: - (F401-1) : Universal Tuohy Borst Adapter.
<http://www.csen.com/F401-1.JPG>
<http://www.csen.com/CSEN549.JPG>
<http://www.csen.com/CSEN549A.JPG>

CSEN 551: - (T0/030) : ELDOR FILTER Tube , 5 micron filter, 10 cm tube.
<http://www.csen.com/ELDOR%20Filter%20Tube.png>
<http://www.csen.com/CSEN551.JPG>
<http://www.csen.com/CSEN551—.png>

CSEN 552: - (F141-2) : Epidural catheter, 19G closed end (3 holes).
<http://www.csen.com/CSEN552.JPG>

CSEN 553: - (Z11840) : 16GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.
<http://www.csen.com/Z11840.png>
<http://www.csen.com/CSEN553.JPG>

CSEN 554: - (F361-1) : 7CC LUER SLIP LOR PLASTIC SYRINGE.
<http://www.csen.com/F361.png>
<http://csen.com/C.png>
<http://www.csen.com/CSEN%20554.JPG>
<http://www.csen.com/CSEN554.png>

CSEN 555: - (G105) : Thread assist guide.
<http://www.csen.com/G105.png>
<http://www.csen.com/G105.JPG>
<http://www.csen.com/CSEN555.JPG>

CSEN 556: - (F361) : 10CC LUER SLIP LOR PLASTIC SYRINGE.
<http://www.csen.com/F361.png>
<http://www.csen.com/F361.JPG>
<http://www.csen.com/CSEN556.JPG>

CSEN 557: - (Z11640) : 22GA X 3.5" TUOHY NEEDLE WITH HUBER POINT.
<http://www.csen.com/Z11640.png>

CSEN 558: - (M6/019+G105+F401-1) : Epidural Catheter 19G x 90 cm closed tip (6 holes), Thread Assist Guide , Tuohy Borst Adapter .
<http://www.csen.com/M6019+F401-1+G105.png>

CSEN 559: - (D307): ADHESIVE BANDAGES.
<http://www.csen.com/D307.png>

CSEN 560: - (J170W-D201): 10CC LUER LOCK LOR PLASTIC SYRINGE
<http://www.csen.com/J170W-D201.png>.

CSEN 561: - (F106) : 20GA EPIDURAL CATHETER CLOSED END (3 holes).
<http://www.csen.com/F106%206PCS.png>
<http://www.csen.com/F106%2010PCS.png>
<http://www.csen.com/CSEN561.JPG>

CSEN 562: - (M6/020) : Epidural catheter, 19G open tip.
<http://www.csen.com/M6-020.png>

CSEN 563: - (J170W-D301): PIN PAD.
<http://www.csen.com/J170W-D301.png>

CSEN 564: - (Z11C40) : Epidural needle Tuohy 18G x 3.5".
<http://www.csen.com/Z11C40.png>

<http://csen.com/I.JPG>

CSEN 565: - (M6/021) : 20GA EPIDURAL CATHETER SINGLE OPEN END HOLE.

<http://www.csen.com/M6021.png>

CSEN 567: - (M62120 + F401-1 + G105) : 19GA EPIDURAL CATHETER 2 SIDE HOLES WITH STYLET, OD 1 mm, TUOHY BORST ADAPTER, THREAD ASSISTING GUIDE. <http://www.csen.com/M62120+F401-1+G105.png>

<http://csen.com/E.png>

<http://www.csen.com/CSEN567-3.JPG>

<http://www.csen.com/CSEN567.JPG>

<http://www.csen.com/CSEN567—.png>

CSEN 591 : (Z41640): Quincke spinal needle, 22 G x 90 mm. Sterile, single use.

<http://www.csen.com/Z41640.png>

<http://www.csen.com/CSEN591.JPG>

<http://www.csen.com/591B.JPG>

<http://www.csen.com/591.JPG>

<http://www.csen.com/CSEN591—.JPG>

<http://www.csen.com/591-J.png>

CSEN 592 : (Z41670): Quincke spinal needle, 22 G x 5 in. Sterile, single use.

<http://www.csen.com/CSEN592.JPG>

<http://www.csen.com/592B.JPG>

<http://www.csen.com/592-J.png>

CSEN 594 : (Z41340): Quincke spinal needle, 25G x 3.5 in. Sterile, single use.

<http://www.csen.com/CSEN594.JPG>

<http://www.csen.com/CSEN594A.JPG>

<http://www.csen.com/594B.JPG>

<http://www.csen.com/594.JPG>

<http://www.csen.com/CSEN594—.JPG>

CSEN 595 : (Z41740): Quincke spinal needle, 20 G x 3.5 in. Sterile, single use.

<http://www.csen.com/CSEN595.JPG>

<http://www.csen.com/595B.JPG>

<http://www.csen.com/595.JPG>

<http://www.csen.com/CSEN595-.JPG>

CSEN 596 : (Z41C40): Quincke spinal needle, 18 G x 3.5 in. Sterile, single use.

<http://www.csen.com/CSEN596.JPG>

<http://www.csen.com/596B.jpg>

<http://www.csen.com/596.JPG>

<http://www.csen.com/CSEN596-.JPG>

CSEN 598 : (Z4/024): Quincke spinal needle, 22 G x 7 in. Sterile, single use.

<http://www.csen.com/CSEN598.JPG>

<http://www.csen.com/598B.jpg>

<http://www.csen.com/598.JPG>

<http://www.csen.com/CSEN598-.JPG>

<http://www.csen.com/598-J.png>

CSEN 599 : Tuohy 20G x 3.5in. Needle, Weiss style winged.

<http://www.csen.com/csen284-13.jpg>

CSEN 601 : Quincke spinal introducer, yellow hub, 20 G x 1 1/4 in.

<http://www.csen.com/CSEN601.JPG>

CSEN 604 : (Z11840): Needle, 16G x 3.5" Tuohy, Graduated with plastic hub and wings.

<http://www.csen.com/Z11840.png>

CSEN 607: (F111-1) Epidural catheter with 3 holes JPG

<http://csen.com/A.JPG>

CSEN 608 : (F111-1+G105) Epidural catheter+Thread assist guide

<http://csen.com/B.png>

CSEN 609 : (Z51640) Spinal needle pencil point

<http://csen.com/H.JPG>

CSEN 612 (M64120) : Epidural catheter with Soft tip, 1.0 x 900 mm.

<http://www.csen.com/CSEN612.png>

CSEN 616 (ZC1710): CSEN GWI (Guide Wire Insertion) Needle 20G x 1-1/2".

<http://www.csen.com/GWI1.JPG>

<http://www.csen.com/GWI2.JPG>

CSEN 617 (PAMA – NG – 25315) : Quincke spinal needle, 25G x 3.5 in. Sterile, single use.

<http://www.csen.com/NG25315.jpg>

CSEN 618 (ECK-1904) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, open tip, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock. <http://www.csen.com/ECK1904.pdf>

CSEN 619 (ECK-1906) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, closed tip, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK1906.pdf>

CSEN 620 (ECK-1908) : Tuohy epidural needle 17G x 3.5", Epidural catheter 19G x 90 cm, closed tip multi, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK1908.pdf>

CSEN 621 (ECK-2004) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, open tip , Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2004.pdf>

CSEN 622 (ECK-2006) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, closed tip , Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2006.pdf>

CSEN 623 (ECK-2008) : Tuohy epidural needle 17G x 3.5", Epidural catheter 20G x 90 cm, closed tip multi, Tuohy Borst adapter, 0.22 um Epidural filter, Thread assist device, 10 ml LOR syringe, Luer Lock.

<http://www.csen.com/ECK2008.pdf>

CSEN 625: – Aiguille Epidurale Tuohy 18G x 3.5" / Tuohy Epidural Needle 18G x 3.5".

<http://www.csen.com/CSEN625-2.JPG>

<http://www.csen.com/CSEN625.JPG>

CSEN 626: – Mini Epidural Kit : Tuohy needle 18G x 8 cm with separate wings, Epidural catheter open tip, Thread Assist Device, Tuohy Borst Adapter, Epidural Filter 0.22 um, LOR syringe 10 ml Luer slip, Packed.

<http://www.csen.com/CSEN626.JPG>

CSEN 627: – Mini Epidural Kit : Tuohy needle 18G x 8 cm with separate wings, Epidural catheter open tip, Thread Assist Device, Tuohy Borst Adapter, Epidural Filter 0.22 um, LOR syringe 10 ml Luer slip.

<http://www.csen.com/CSEN627.JPG>

CSEN 628 : - Epidural catheter 21G, closed tip, 3 lateral holes.

<http://www.csen.com/CSEN628.JPG>

CSEN 629 : - Quincke Spinal Needle 20G, 150 mm.

<http://www.csen.com/CSEN629.JPG>

CSEN 630 : - Removable Epidural Needle Wings.

<http://www.csen.com/RENEW.jpg>

CSEN 631 : CSEN Epidural Mini Kit : Aiguille Tuohy à pointe Huber 17G x 3½" / Tuohy needle with Huber 17G x 3½", Epidural Catheter 19G x 90 cm closed tip multi (6 lateral holes), Thread Assist Guide, Tuohy Borst Adapter, Seringue LOR 10ml Luer Slip / 10ml LOR Syringe Luer Slip, Filtre épidural 0.22µm / Epidural filter 0.22µm. Packed.

<http://www.csen.com/CSEN%20631.JPG>

CSEN 633 : CSEN Epidural Mini Kit : Aiguille Tuohy à pointe Huber 17G x 3½" / Tuohy needle with Huber 17G x 3½", Adaptateur Tuohy Borst universel / Universal Tuohy Borst adapter, Seringue LOR 7ml Luer Slip / 7ml LOR Syringe Luer Slip, 21G

EPIDURAL CATHETER CLOSED END 3 lateral holes, Guide à cathéter / Thread assist guide , Filtre épidural 0.22µm / Epidural filter 0.22µm, Epidural label (Yellow label written “Epidural Catheter” on it). Packed.

<http://www.csen.com/CSEN%20633.JPG>

CSEN 676 : Disposable Surgical Gown X , packed, sterilized.

<http://www.csen.com/gown1.jpg>

<http://www.csen.com/gown2.jpg>

CSEN 677 : Disposable Surgical Gown XL , packed, sterilized.

CSEN 678 : Disposable Surgical Gown XXL , packed, sterilized.

CSEN 679 : Machine needle 18G, 120 mm.

CSEN 680 : DISPOSABLE INFUSION PUMP

<https://www.smashwords.com/books/view/1070971>

CSEN 681 : INTRALIPID : A Magic Bullet in COVID 19 ?

<https://www.smashwords.com/books/view/1056670>

CSEN 682 : Quantum Consciousness

<https://www.smashwords.com/books/view/986828>

CSEN 683 : Intralipid is a Magic Bullet in Cancer treatment?

<https://www.smashwords.com/books/view/936029>

CSEN 684 : Intralipid-Iodine for Imaging

<https://www.smashwords.com/books/view/880608>

CSEN 685 : Intralipid Is a Magic Bullet

<https://www.smashwords.com/books/view/843448>

CSEN 686 : Mitochondrial theory of aging : to be or not to be ? (The role of Linoleic acid)

<https://www.smashwords.com/books/view/821212>

CSEN 687 : Intralipid and other Lipid Emulsions

<https://www.smashwords.com/books/view/798442>

CSEN 688 : Intralipid Fatty Acids : The Code of Life ?

<https://www.smashwords.com/books/view/798434>

CSEN 689 : Intralipid infusion for Myelin Sheath Repair in Multiple Sclerosis and Trigeminal Neuralgia?

<https://www.smashwords.com/books/view/790008>

CSEN 690 : Intralipid

<https://www.smashwords.com/books/view/777382>

CSEN 691 : The Spectrum of Amniotic Fluid Embolism - Is Intralipid the solution ?

<https://www.smashwords.com/books/view/774210>

CSEN 692 : Match Cognitive Brain Test

<https://www.smashwords.com/books/view/768631>

CSEN 693 : Serendipity Based Medicine (SBM) : To Infinity and Beyond

<https://www.smashwords.com/books/view/760572>

CSEN 694 : Combined Spinal-Epidural Anesthesia

<https://www.smashwords.com/books/view/760568>

CSEN 695 : Bread Fat

<https://www.smashwords.com/books/view/758878>

CSEN 696 : [Can intralipid infusion open a coronary occlusion causing acute myocardial infarction?](https://www.oatext.com/can-intralipid-infusion-open-a-coronary-occlusion-causing-acute-myocardial-infarction.php)

<https://www.oatext.com/can-intralipid-infusion-open-a-coronary-occlusion-causing-acute-myocardial-infarction.php>

CSEN 697 : Intralipid treatment: Is it only the tip of an iceberg? A new suggestion: Bone cement implantation syndrome (BCIS)

<https://www.hoajonline.com/jacs/2049-9752/1/12>

CSEN 698 : [Intralipid treatment for Alzheimer Disease](https://www.oatext.com/intralipid-treatment-for-alzheimer-disease.php)

<https://www.oatext.com/intralipid-treatment-for-alzheimer-disease.php>

CSEN 699 : Intralipid Treatment for Post Operative Delirium

<http://medcraveonline.com/JACCOA/JACCOA-07-00273.php>

CSEN 700 : Combined Spinal Epidural Catheters for Epidural Cooling, Cerebrospinal Fluid Aspiration and Spinal Intralipid Infusion for Treatment of Spinal and Brain Injuries, Diseases and Protection

<https://www.scirp.org/journal/paperinformation.aspx?paperid=42375>

CSEN 701 : NESA Intralipid Project The New Game Player

<https://slidetodoc.com/nesa-intralipid-project-the-new-game-player-joseph/>

CSEN 702 : Intralipid Infusion for Myelin Sheath Repair in Multiple Sclerosis and Trigeminal Neuralgia?

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Intralipid-Infusion-for-Myelin-Sheath-Repair-in-Multiple-Sclerosis-and-Trigeminal-Neuralgia-#fulltext>

CSEN 703 : Lipid Emulsion for Local Anesthesia Reversal (LAR) after Prolonged Spinal/Epidural Anesthesia

<http://innovationinfo.org/articles/JHSD-106.pdf>

CSEN 704 : Intralipid Treatment of Preeclampsia/Eclampsia?

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Intralipid-Treatment-of-Preeclampsia-Eclampsia-#fulltext>

CSEN 705 : Intralipid and the Immune system

<http://www.csen.com/immune.pdf>

CSEN 706 : Local Anesthesia Reversal (LAR) of Total Spinal Anesthesia (TSA) by Lipofundin (Lipid Emulsion)

<https://innovationinfo.org/journal-of-health-science-and-development/article/Local-Anesthesia-Reversal-LAR-of-Total-Spinal-Anesthesia-TSA-by-Lipofundin-Lipid-Emulsion-#fulltext>

CSEN 707 : Intralipid for Amniotic Fluid Embolism (AFE)?

<http://www.csen.com/oja.pdf>

CSEN 708 : INTRALIPID FOR THERAPY

<http://www.csen.com/Pat1.pdf>

CSEN 709 : Intralipid-A Magic Bullet?

<http://www.csen.com/IMB2.pdf>

CSEN 710 : Intralipid Treatment for New-borns with Meconium Stained Amniotic Fluid (MSAF)

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Intralipid-Treatment-for-New-borns-with-Meconium-Stained-Amniotic-Fluid-MSAF-#fulltext>

CSEN 711 : Lipid Emulsion Treatment (LET) of Post-Operative Cognitive Dysfunction (POCD)

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Lipid-Emulsion-Treatment-LET-of-Post-Operative-Cognitive-Dysfunction-POCD-#fulltext>

CSEN 712 : Warburg effect and Intralipid effect on the mitochondria of cancer cells

<http://www.csen.com/Warburg.pdf>

CSEN 713 : Joseph Eldor (Editor in Chief) : Frontiers Journal of Pediatrics and Neonatology (FJPN)

<https://fmpublishers.org/journals/Frontiers-Journal-of-Pediatrics-and-Neonatology-FJPN-editorial-board>

CSEN 714 : Intralipid-Iodine for Hysterosalpingography in infertile women

<http://www.csen.com/ILH.pdf>

CSEN 715 : Editor in Chief: [Joseph Eldor](#) : Journal of Health Science and Development

<https://www.innovationinfo.org/journal-of-health-science-and-development>

CSEN 716 : Intralipid Eye Drops for Cataract Treatment?

<http://www.csen.com/cataract.pdf>

CSEN 717 : Postoperative Shivering after Spinal Anesthesia Treated by Lipid Emulsion

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Postoperative-Shivering-after-Spinal-Anesthesia-Treated-by-Lipid-Emulsion#fulltext>

CSEN 718 : Local Anesthesia Reversal (LAR) of Spinal Anesthesia by Lipid Emulsion (Lipofundin 20%) for Day Case Surgery

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Local-Anesthesia-Reversal-LAR-of-Spinal-Anesthesia-by-Lipid-Emulsion-Lipofundin-20-for-Day-Case-Surgery#fulltext>

CSEN 719 : Upper Arm Local Anesthesia Reversal (LAR) Using Lipofundin (3 Case Reports)

<http://innovationinfo.org/articles/JHSD-110.pdf>

CSEN 720 : General Anesthesia Reversal (GAR) by Lipofundin

<https://innovationinfo.org/articles/JHSD/JHSD-113.pdf>

CSEN 721 : Lipid Emulsion for Local Anesthesia Reversal (LAR) after Prolonged Spinal/Epidural Anesthesia

<https://www.innovationinfo.org/articles/JHSD/JHSD-106.pdf>

CSEN 722 : Lipid Emulsion Treatment for Post Spinal Anesthesia Myoclonus

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Lipid-Emulsion-Treatment-for-Post-Spinal-Anesthesia-Myoclonus#fulltext>

CSEN 723 : The First Case Report of Local Anesthesia Reversal (LAR) of the Upper Arm Brachial Plexus Block by Lipid Emulsion

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/The-First-Case-Report-of-Local-Anesthesia-Reversal-LAR-of-the-Upper-Arm-Brachial-Plexus-Block-by-Lipid-Emulsion#fulltext>

CSEN 724 : INTRALIPID Rescue of Amniotic Fluid Embolism: from Theory to Existence

<https://www.innovationinfo.org/articles/JHSD/JHSD-118.pdf>

CSEN 725 : Linoleic Acid: The Code of Life?

<https://www.innovationinfo.org/journal-of-health-science-and-development/article/Linoleic-Acid-The-Code-of-Life-#fulltext>

CSEN 726 : ELDOR CSEN patents in needles and catheter

CSEN 726.1 : Multiple hole epidural catheter

<https://patents.google.com/patent/US5800407A/en?inventor=%22joseph+eldor%22&dq=%22joseph+eldor%22>

CSEN 726.2 : Multiple hole spinal needle

<https://patents.google.com/patent/EP0727187B1/en?inventor=%22joseph+eldor%22&dq=%22joseph+eldor%22>

CSEN 726.3 : Device for combined spinal and epidural anesthesia

<https://patents.google.com/patent/US5163901A/en?inventor=%22joseph+eldor%22&dq=%22joseph+eldor%22>

