

Neutron Spin Echo Spectroscopy



Piotr A. Żołnierczuk

Jülich Centre for Neutron Science, Forschungszentrum Jülich

1. Introduction
2. The Principles of NSE
3. The SNS-NSE Spectrometer
4. Examples
5. Summary

Neutron Spin Echo in A Nutshell


- NSE measures very small velocity changes by using neutron magnetic moment as a “clock”
- NSE “end-product” is the Intermediate Scattering Function: $I(Q, \tau)$
- NSE can use broad wavelength band while maintaining very good resolution
- NSE is a complementary technique to SANS in that it provides the dynamic information about the system
- NSE is “counting intensive”: long counting times, large samples



Neutron and X-Ray Instruments Zoology

With NSE one can study:

- Coherent Dynamics
 - Diffusion
 - Shape fluctuations
 - Polymer dynamics
 - Glassy systems
- Incoherent Dynamics
 - Hydrogen
- Magnetic Dynamics
 - Spin Glasses



Source: ESS

NSE Spectrometers in the World



Classic “IN11-Type”

- IN11 - Institute Laue-Langevin, Grenoble, France
- IN15 - Institute Laue-Langevin, Grenoble, France
- J-NSE – JCNS hosted by FRM-II, Garching, Germany
- NG-A NSE – NIST, Gaithersburg, USA
- C2-3-1 iNSE, ISSP JRR-3M, Tokai, Japan
- BL-15 SNS-NSE – JCNS/ORNL, Oak Ridge, USA

- MUSES [NSRE] – Laboratoire Léon Brillouin, Saclay, France
- RESEDA [NSRE] – FRM-II, Garching, Germany
- VIN-ROSE [NSRE] – J-PARC, Tokai, Japan
- MIRA [TAS+MIEZE] – FRM-II, Garching, Germany
- FLEXX [TAS+NSRE] – HZB, Berlin, Germany
- SESANS [SE+SANS] – TU Delft, Holland

The Principles of Neutron Spin Echo

Inelastic Neutron Scattering

Dynamic Structure Factor

$$\frac{d\sigma}{d\Omega dE} \sim S(Q, \omega)$$


Neutron scattering kinematics:

$$\Delta E = \hbar\omega = (mv_1^2 - mv_2^2)/2$$

$$\vec{Q} = \vec{k}_1 - \vec{k}_2$$

$$|\vec{k}| = \frac{2\pi}{\lambda}$$

$$\lambda = \frac{h}{mv} \quad (\text{Louis de Broglie, 1924})$$



Bloch Equation:


$$\frac{d\vec{s}}{dt} = \gamma \vec{s} \times \vec{B}$$

Larmor Frequency

$$\omega_L = |\gamma B|$$



Neutron Gyromagnetic Ratio

$$\gamma/2\pi \approx 30 \text{ MHz/T}$$



Source: Wikipedia

Larmor Precession (II)




$$B=1.6 \text{ mT} \quad \lambda=1.6 \text{ \AA}$$


Neutron Spin Manipulation

Mezei Flippers

π -flipper



$\pi/2$ -flipper



Larmor Precession (III)

Accumulated phase:

$$\varphi = \omega_L t = \gamma B l \frac{1}{v}$$


or

$$\varphi = \gamma \frac{m}{h} J \lambda$$

for example:

$$\lambda = 8\text{\AA}, B = 0.4\text{T}, l = 1.2\text{m}$$

$$\rightarrow \varphi / 2\pi \cong 3 \times 10^4 \text{ [Number of Turns]}$$





Notes:

$$J \stackrel{\text{def}}{=} B l \quad \text{or more precisely} \quad J \stackrel{\text{def}}{=} \int \vec{B} \cdot \vec{dl}$$

Neutron Spin Echo (I)

neutron spin


spin rotation



Neutron Spin Echo (II)

neutron spin



spin rotation



Neutron Spin Echo (Inelastic)

neutron spin

spin rotation



Neutron Spin Echo

Energy change (quasi-elastic scattering)

$$\Delta E = \hbar\omega = \frac{1}{2}mv_1^2 - \frac{1}{2}mv_2^2 \cong m v dv$$

Accumulated phase:

$$\phi_2 - \phi_1 = \gamma \left(\frac{B_2 l_2}{v_2} - \frac{B_1 l_1}{v_1} \right) = \gamma B l \left(\frac{1}{v_2} - \frac{1}{v_1} \right) \quad \leftarrow B_1 l_1 = B_2 l_2$$

$$\Delta\phi \cong \gamma B l \frac{\hbar\omega}{mv^3} = \frac{\gamma}{2\pi} \left(\frac{m}{h} \right)^2 J \lambda^3 \omega$$

or

$$\Delta\phi = \tau\omega$$

where we define Fourier time $\tau \stackrel{\text{def}}{=} \frac{\gamma}{2\pi} \left(\frac{m}{h} \right)^2 J \lambda^3$

For quick computation $\tau \cong 0.186 J \lambda^3$

where $[\tau] = \text{ns}$, $[J] = \text{Tm}$ and $[\lambda] = \text{Å}$

Example: $\lambda = 8\text{Å}$, $J = 0.56 \text{Tm} \rightarrow \tau \cong 50\text{ns} \rightarrow \Delta E \cong 0.01 \mu\text{eV}$

Signal: $I \sim \langle \cos \phi \rangle = \langle \cos \omega \tau \rangle$


$$I \sim \frac{1}{2} \left[S(Q) \pm \int S(Q, \omega) \cos(\omega \tau) d\omega \right]$$

Fourier transform
(Real part)

$$B \sim S(Q)$$


$$A \sim I(Q, \tau) = \mathcal{F}[S(Q, \omega)]$$

$$\frac{I(Q, \tau)}{I(Q, 0)} = \frac{2 A}{U - D}$$



Standard 4-point echo evaluation

$\Delta\lambda/\lambda_0 \sim 15\text{-}20\%$



90-degree steps

$$E_1 = B + A \sin(\varphi_0)$$

$$E_2 = B + A \cos(\varphi_0)$$

$$E_3 = B - A \sin(\varphi_0)$$

$$E_4 = B - A \cos(\varphi_0)$$




$$B = \frac{1}{4} \sum_{i=1..4} E_i$$

$$A = \sqrt{\frac{1}{2} \sum_{i=1..4} (E_i - B)^2}$$

$\Delta\lambda/\lambda_0$ up to ~50%

Echo Signals





- 1) F. Mezei, Nucl. Inst. Methods **164**, 153-156 (1979)
- 2) B. Farago, Time-of-Flight Neutron Spin Echo: Present Status in F.Mezei, C.Pappas, T.Gutberlet Neutron Spin Echo Spectroscopy, Springer (2003)

Energy and Time Domain (Backscattering ↔ NSE)

Dynamic Structure Factor
 $S(Q, \omega)$

Intermediate Scattering Function
 $I(Q, \tau)$

F.T.






$$I_D(Q, \omega) = S(Q, \omega) * R(Q, \omega)$$

$$I_D(Q, \tau) = I(Q, \tau) R(Q, \tau)$$

To get $S(Q, \omega)$ we have to
de-convolute instrument resolution

To get $I(Q, \tau)$ we just
divide out instrument resolution

SamSDSMicelle (./fits_8327 echo[10])



1. Resolution

- symmetry phase (echo)
- get $R(Q, \tau; \text{pixel})$

2. Sample: $I_{raw}(Q, \tau; \text{pixel})$


3. Background:

- $I_{bgr}(Q, \tau; \text{pixel})$
- correction $\rightarrow I_{sig}(Q, \tau; \text{pixel})$

4. Compute $I(Q, \tau; \text{pixel}) = \frac{I_{sig}(Q, \tau)}{R(Q, \tau)}$


5. Group results

Complications




$$I_{\text{up}} = I_{\text{coh}} + \frac{1}{3} I_{\text{inc}}$$

$$I_{\text{down}} = \frac{2}{3} I_{\text{inc}}$$




$$I_{\text{amp}} = I_{\text{coh}} f_{\text{coh}}(\tau) - \frac{1}{3} I_{\text{inc}} f_{\text{inc}}(\tau)$$

Correction Coils






- Problem:
 - Field integral $J = Bl$ the same for all trajectories
 - Solenoid field $B(r) \sim r^2$
- Solution:
 - correction coils with current distribution that varies as r^2




Theme Variations

Paramagnetic NSE




“Non-standard” NSE Instruments

NSRE




RF field instead of solenoid

SESANS



Spin Echo encoded SANS

SNS-NSE: NSE Spectrometer at SNS



Backscattering Spectrometer (BASIS) • BL - 2

Dynamics of macromolecules, constrained molecular systems, polymers, biology, chemistry, materials science

Eugene Mamontov • 865.574.5109 • mamontove@ornl.gov

Nanoscale-Ordered Materials Diffractometer (NOMAD) • BL - 1B

Liquids, solutions, glasses, polymers, nanocrystalline and partially ordered complex materials

Joerg Neufeind • 865.241.1635 • neufeindjc@ornl.gov

Wide Angular-Range Chopper Spectrometer (ARCS) • BL - 18

Atomic-level dynamics in materials science, chemistry, condensed matter sciences

Doug Abernathy • 865.617.5624 • abernathyd@ornl.gov

Fine-Resolution Fermi Chopper Spectrometer (SEQUOIA) • BL - 17

Dynamics of complex fluids, quantum fluids, magnetism, condensed matter, materials science

Matthew Stone • 865.202.6896 • stonemb@ornl.gov

Spallation Neutrons and Pressure Diffractometer (SNAP) • BL - 3

Materials science, geology, earth and environmental sciences

Chris Tulk • 865.574.5764 • tulkca@ornl.gov

Ultra-Small-Angle Neutron Scattering Instrument (USANS) • BL - 1A

Life sciences, polymers, materials science, earth and environmental sciences

Michael Agamalian • 865.382.6367 • magamalian@ornl.gov

Vibrational Spectrometer (VISION) • BL - 16B

Vibrational dynamics in molecular systems, chemistry

Yongqiang Cheng • 865.5576.9064 • chengy@ornl.gov

BL - 16A

Neutron Spin Echo Spectrometer (NSE) • BL - 15

High-resolution dynamics of slow processes, polymers, biological macromolecules

Michael Ohl • 865.574.8426 • ohlme@ornl.gov

Hybrid Spectrometer (HYSPEC) • BL - 14B

Atomic-level dynamics in single crystals, magnetism, condensed matter sciences

Barry Winn • 865.805.6819 • winnbl@ornl.gov

BL - 14A

Magnetism Reflectometer • BL - 4A

Chemistry, magnetism of layered systems and interfaces

Valeria Lauter • 865.387.5389 • lauterv@ornl.gov

Liquids Reflectometer • BL - 4B

Interfaces in complex fluids, polymers, chemistry

John Anknier • 865.377.4523 • anknierj@ornl.gov

Cold Neutron Chopper Spectrometer (CNCS) • BL - 5

Condensed matter physics, materials science, chemistry, biology, environmental science

Georg Ehlers • 865.206.0824 • ehlersg@ornl.gov

Fundamental Neutron Physics Beam Line • BL - 13

Fundamental properties of neutrons

Geoffrey Greene • 865.576.0023 • greenejl@ornl.gov

Extended Q-Range Small-Angle Neutron Scattering Diffractometer (EQ-SANS) • BL - 6

Life science, polymer and colloidal systems, materials science, earth and environmental sciences

William Heller • 865.241.0093 • hellerwt@ornl.gov

BL - 8A
BL - 8B

Elastic Diffuse Scattering Spectrometer (CORELLI) • BL - 9

Detailed studies of disorder in crystalline materials

Feng Ye • 865.576.0931 • yef1@ornl.gov

Versatile Neutron Imaging Instrument at SNS (VENUS) • BL - 10

Energy selective imaging in materials science, engineering, materials processing, environmental sciences and biology

Hassina Bilheux • 865.384.9630 • bilheuxhn@ornl.gov

Macromolecular Neutron Diffractometer (MaNDI) • BL - 11B

Atomic-level structures of membrane proteins, drug complexes, DNA

Leighton Coates • 865.964.8651 • coatesl@ornl.gov

Single-Crystal Diffractometer (TOPAZ) • BL - 12

Atomic-level structures in chemistry, biology, earth science, materials science, condensed matter physics

Christina Hoffmann • 865.576.5127 • hoffmann@ornl.gov

Powder Diffractometer (POWGEN) • BL - 11A

Atomic-level structures in chemistry, materials science, and condensed matter physics including magnetic spin structures

Ashfia Huq • 630.986.7321 • huqa@ornl.gov

Engineering Materials Diffractometer (VULCAN) • BL - 7


Mechanical behaviors, materials science, materials processing

Ke An • 865.919.5226 • kean@ornl.gov


***Scheduled commissioning date**

LEGEND

- Operating instrument in user program
- In design or construction
- Under consideration





SNS-NSE: IN11-Type NSE



Typical neutron wavelengths:
 $\lambda = 2 - 14 \text{ \AA}$, $\Delta\lambda = 2.5, 3.6 \text{ \AA}$



$J_{\max} \sim 0.56 \text{ Tm}$
 $2\theta = 3.5^\circ - 79.5^\circ$

SNS-NSE: the first NSE with superconducting main coils



“Pythagorean” correction
coils: $r^2 = y^2 + z^2$

Magnetic Shielding






$\Delta\phi \rightarrow 1 \text{ deg}$

Two "historical" examples:

1. SDS Micelles (1981)
2. Reptation Model (2002)

Several recent results from SNS-NSE




$$\frac{I(Q, t)}{I(Q, 0)} = e^{-D_{\text{eff}}(Q)Q^2 t}$$

$$D_{\text{eff}} = D_0 H(Q) / S(Q)$$


$$D_0 = \frac{kT}{6\pi\eta R}$$




Reptation Model (de Gennes/Edwards)



Coherent scattering, labeled chain

Melt of long chain linear PE









$$\frac{I(Q, t)}{I(Q, 0)} = \exp(-(\Gamma_{\text{Bend}} t)^{2/3})$$


$$\Gamma_{\text{Bend}} = 0.025 \gamma_k \left(\frac{k_B T}{\kappa}\right)^{1/2} \left(\frac{k_B T}{\eta}\right) Q^3$$

Mechanical Properties of Nanoscopic Lipid Domains



$$\frac{S(Q, t)}{S(Q, 0)} = \exp(-(\Gamma(Q)t)^{2/3})$$








$$F(Q, t) = F_{\text{trans}}(Q, t) \cdot F_{\text{rot}}(Q, t) \cdot F_{\text{int}}(Q, t)$$

$$F_{\text{trans}}(Q, t) = \exp(-Q^2 D_T t)$$

$$F_{\text{int}}(Q, t) = \left\langle \sum_{\alpha, \beta} b_{\alpha} b_{\beta} \exp(iQR_{\beta}^{eq}) \exp(-iQR_{\alpha}^{eq}) \cdot f_{\alpha\beta}(Q, \infty) \cdot f'_{\alpha\beta}(Q, t) \right\rangle$$



Example of Paramagnetic NSE



- NSE
 - a high resolution neutron spectroscopy
 - measures the intermediate scattering function $S(Q, \tau)$
 - need large samples, scattering intensity
- SNS-NSE (BL-15)
 - the first NSE at a Spallation Source
 - the first one with superconducting coils
 - the only one with magnetic shielding
 - available in user program – write a proposal!!

Acknowledgements: M. Monkenbusch, D. Richter, L. Stingaciu, M. Cochran, G. Ehlers, W-R Chen, T. Kozielski, M. Ohl, N. Arend and M. Sharp

1. R. Pynn, Neutron Scattering, Neutron Spin Echo
http://www.iub.edu/~neutron/notes/20061204_Pynn.pdf
2. F. Mezei, Z. Physik (1972) 255: 146.
3. F. Mezei (Ed.): Neutron Spin-Echo, Lecture Notes in Physics 128, Springer, 1980.
4. F. Mezei, C. Pappas, T. Gutberlet (Eds.): Neutron Spin-Echo Spectroscopy, Lecture Notes in Physics 601, Springer, 2003.
5. D. Richter, M. Monkenbusch, A. Arbe, J. Colmenero, Neutron Spin Echo in Polymer Systems, Adv. in Polymer Science 174, Springer, 2005
6. M. Monkenbusch, D. Richter, High Resolution Neutron Spectroscopy
<http://doi.org/10.1016/j.crhy.2007.10.001>
7. C.Pappas, G. Ehlers, F. Mezei, Neutron Spin Echo and Magnetism in Neutron Scattering from Magnetic Materials, ed. T. Chatterji, Spinger 2006
8. B.Farago, Basics of Neutron Spin-Echo, ILL Neutron Data Booklet,
https://www.ill.eu/fileadmin/users_files/documents/links/documentation/NeutronDataBooklet.pdf
9. G.L. Squires, Introduction to the Theory of Thermal Neutron Scattering, Cambridge Univ. Press, 3rd edition (2012)