

# Aesculap® Acculan 3Ti Dermatome

Look. Feel. Experience.



Aesculap Power Systems

# Acculan 3Ti dermatome

## Look. Feel. Experience.

Aesculap is one of the leading suppliers of motor systems. Our objective is to provide you with instruments upon which you can depend. Driven by this objective, we rely on the most innovative technologies and the highest-quality materials in the development of our systems.

With the new Acculan 3Ti dermatome, Aesculap has set a new milestone – the most modern drive technology in a robust titanium housing combined with extensively proven non-sterile battery technology – a dermatome you can rely on.





The difference  
is in the detail.

The difference is in the detail – Acculan 3Ti avoids sterilization not only of the temperature-sensitive battery cells, but also of the equally vulnerable control electronics by integrating it in the battery unit.

A concept which offers you the reliability that you expect and also spares your budget.

# Acculan 3Ti dermatome at a glance



## ■ Acculan 3Ti dermatome GA670

The best dermatome ever from Aesculap. More performance and reliability thanks to the patented motor technology – an instrument you can rely on.

Due to its titanium housing, this dermatome is suitable for use in conjunction with the Eccos storage system for machine reprocessing using alkaline cleaning agents and disinfectants.

- Quick and simple blade changes
- Cutting width max. 78 mm, adjustable by 10 flaps of width 8 mm
- 0.2 – 1.2 mm cutting thickness in 1/10 mm steps, keyless adjustment



### ■ A clever combination – GA677 charger and GA666 battery

- Intuitive, clear indicator elements
- Automatic battery check
- Charge status display
- Electronically controlled charging to avoid memory effect
- Compatible with predecessor Acculan II through simple-to-replace charging bays
- Innovative battery exchange concept spares the environment and your budget

# Aesculap skin graft mesher



## ■ BA720R skin graft mesher

The ideal supplement to the Acculan 3Ti dermatome for creating meshed split-thickness skin grafts.

- Precise guidance at the skin carrier for best results
- Keyless disassembly for cleaning
- Machine cleanable with Eccos

## ■ Skin carriers



- BA721  
1:1.5
- BA722  
1:3
- BA723  
1:6

# Ordering information

GA670 Acculan 3Ti dermatome  
incl. battery, lid and sterile funnel



## Accessories







GA666 Small NiMH battery for GA670  
 GA675 Lid  
 GA677 Charger  
 GA678 Sterile funnel  
 GA679 Battery removal device  
 GB228R Dermatome blades   10  
 GB256R Eccos set for GA670  
 GB600 Sterilit Power Systems oil spray  
 GA059 Sterilit Power Systems drip lubricator




BA720R Skin graft mesher incl. BA726R ratchet wrench and  
BA725R cutting roller



## Accessories

BA721 Skin carrier 1:1.5   10  
 BA722 Skin carrier 1:3   10  
 BA723 Skin carrier 1:6   10  
 BA725R Spare cutting roller  
 GB688R Eccos fixation for BA720R  
 JC254R Tray suitable for GB688R



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

The main product mark 'Aesculap' and the product marks '.....', '...' and '...' are registered marks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. 035902

0910/1/2