

rbmedical

Providing the instruments for good health

rbmedical



Ophthalmic Instruments

www.rbmedical.co.uk

Contents

Reusable

Speculums	Page 3-8
Retractors	Page 9-10
Knives & Spuds	Page 11
Punches	Page 12
Handles & Blade Breakers	Page 13
Curettes & Spatulas.....	Page 14
Lasik, Lasek & Dalk	Page 15-19
Spatulas, Chopper & Manipulators	Page 20-28
Lens Loops & Spoons	Page 29-30
Spoons, Hooks, Manipulators & Fixation Rings	Page 31-34
Muscle Hooks	Page 35-36
Corneal Markers	Page 37-39
Calipers	Page 40
Lens Injector	Page 41
Scissors	Page 42-52
Forceps.....	Page 53-93
Needle Holders	Page 94-98
Lacrimal Instruments	Page 99-103
Cannulas	Page 104-107
I/A Bimanual Handles.....	Page 108

Single Use


Speculums	Page 110
Chopper & Manipulators	Page 111-112
Scissors	Page 113
Forceps.....	Page 114-117
Forceps, Calipers & Gauge	Page 118
Fixation Rings, Lasik & Sub Tenon Set.....	Page 119
Needle Holders & Cannulas	Page 120
Lacrimal.....	Page 121
Gauges & Scissors.....	Page 122
Foceps & Holders.....	Page 123
Index	Page 124-127
Contacts	Page 128

Speculums

Barraquer Wire Speculums

standard wire 0.8mm, 10mm blades
heavy wire 1mm, 15mm blades
standard wire 0.8mm, 15mm blades

OP001




Barraquer Wire Speculums

Wire Speculum
10mm blades (infant)
15mm blades

Grandon-Barraquer
15mm blades

OP002



Temporal Approach

Wire Speculum
with spring on nasal side
15mm closed blades
15mm open blades


OP003



Kratz-Barraquer

Wire Speculum
with open blades
For unobstructed access to the
cataract incision with phaco or I/A tip.
heavy wire 1mm, 15 mm blades
standard wire 0.8mm, 15mm blades
standard wire 0.8mm, 13mm blades
standard wire 0.8mm, 15mm blades
heavy wire 1mm, 15 mm blades

OP004



Barraquer Solid Blade

Wire Speculum
Barraquer, 10mm blades
Barraquer, 15mm blades
Feaster, 20mm blades

OP005



Sauer

Lid Speculum
for infants, 11mm blades
for premature infants, 2mm blades

OP006



Frankel

Infant Lid Speculum adjustable
for precise control of
infant eyelid retraction.
10mm blades

OP007



Infant Wire Speculum

Blades size 5mm, Open
Blades size 5mm, Close
Blades size 6mm, Solid


OP008



Alphonso

Infant Lid Speculum finger grips and
interlacing blades to facilitate insertion
one blade 5mm, one blade 7mm


OP009



Infant Lid Speculum

Lid Speculum closed wire blades
gently curved spring
4mm blades

OP010



Cook Eye Speculum

Eye Speculum with locking screw
#1, 8mm blades
#2, 11mm blades
#3, 15mm blades

OP011



Speculums

Murdoch

Eye Speculum with locking screw
#2, 12mm blades
#3, 15mm blades

OP012



Mellinger

Eye Speculum self retaining
for unobstructed access to the
cataract incision with phaco or I/A tip.
15mm blades

OP013



Williams

Eye Speculum with locking screw
11mm blades
15mm blades

OP014



Knapp

Eye Speculum
13mm blades

OP015



Lancaster Speculum

Eye Speculum
adjustable mechanism
with locking nut
16mm blades

OP016



Speculums

Lester - Burch Speculum

Eye Speculum with flanges on blades for retaining eye lashes
adjustable mechanism
with locking screw

OP017









Lieberman Speculum

adjustable mechanism

- Lieberman, K-wire, 5mm blades, infant
- Lieberman, K-wire, 10mm blades, small
- Lieberman, K-wire, 15mm blades, standard
- Lieberman, V-shaped, 15mm blades
- Tanna-Lieberman, 15mm blades
- Thorlakson, rounded, 13mm blades
- Lieberman, triple post, 13mm blades
- Titanium 15mm K-wire, blades
- Titanium 15mm W-wire, blades
- Titanium rounded, 13mm blades

OP018



Speculums

Lieberman Speculum

Eye Speculum
Adjustable mechanism
non-reflective finish
5mm blades, Infant
10mm blades, Small
15mm blades, Standard

OP019



Jaffe

Wire Lid Retractor
pair, 15mm blades (large)
pair, 13mm blades (small)

OP020



Castroviejo Speculum

Eye Speculum
adjustable mechanism
5mm blades
10mm blades
14mm blades
16mm blades

OP021



Bangerter Speculum

15mm blades
10mm blades

OP022



Stevenson

Lacrimal Sac Retractor
3x3 blunt prongs, adjustable

OP023



Stevenson

Lacrimal Sac Retractor
solid blades with serrated
edge adjustable
13mm blades





Weiss
EYE SPECULUM Fenestrated

OP024




Lieberman Style, aspirating adjustable mechanism
Features aspirating ports to remove excess fluid from the surgical field. Supplied with silicone tubing and luer-lock adaptor.
Lieberman, Kratz style, 15mm blades
Thorlakson, rounded, 11mm blades
Lieberman, V-shaped, 15mm blades
Chu, solid, reversible, 15mm blades
Replacement silicone tubing, 100cm (39")

OP025



Flieringa
Scleral Fixation Ring
set of eight, 15mm-22mm

OP026



McNeill-Goldman
Scleral Fixation Ring and Blepharostat
Small, 15mm I.D.
Medium, 16mm I.D.
Large, 17mm I.D.

OP027



Retractors

Schocket

Scleral Depressor double-end
with pocket clip 6mm wide bar
and 2.5mm tear drop

OP028



Josephberg-Besser

OP029



Schepens Orbital Retractor

OP030



Senn

Senn Miller Lid Retractor
Sharp, 7x8mm, 16 cm 6¼"
Blunt, 7x8mm, 16 cm 6¼"

OP031



Senn

Senn Miller Lid Retractor,
Double Ended
Sharp
Blunt

OP032



Desmarres

Lid Retractor
thin solid blades

size 0, 11mm wide

size 1, 13mm wide

size 2, 15mm wide

size 3, 17mm wide

OP033



Jaeger

Lid Plate
stainless steel, Small
stainless steel, Large

OP034



Wright

Fascia Needle
1mm x 6mm oval hole
140mm long

OP035



Knapp

Lacrimal Retractor
4 prongs
blunt
sharp

OP036



Rollet

Rollet retractor
4-blunt
3-blunt

OP037



Helveston

Tissue Retractor
thin curved blade
"Barbie" 7mm wide
"Big Barbie" 9mm wide
"Great Big Barbie"
11mm wide concave blade

OP038



Knives & Spuds

Tooke

Corneal Knife
straight blade
with curved cutting edge

OP039



Gill

Corneal Knife
curved cutting edge
straight
light curve

OP040



Fukasaku

"Hockey" Knife
sharp 3mm front edge
semi-sharp 6mm
bottom edge

OP041



Foreign body spuds

Golf Club, angled, lightly rounded tip
Francis, angled, flag shaped tip
Ellis, lightly curved and rounded tip


OP042



Spud & Gouge

in reversible handle for
removing foreign bodies.

OP043



Gass

Corneoscleral Punch
1.5mm diameter head
0.75mm deep bite
rotatable tip for
convenient positioning


OP044



Kelly

Descemet's Membrane Punch
1.0mm diameter head
0.75mm deep bite

OP045



Fukasaku

Micro Trabeculectomy Punch
0.6mm diameter head
with bullet shaped tip
0.3mm x 0.6mm deep bite


OP046



Christensen

Punctum Punch
1mm wide bite
for opening stenotic punctums.

OP047



Handles/Blade Breakers

Castroviejo
Blade breaker and holder
with sliding lock.

OP048



Swiss Model
Blade breaker and holder
spherical tip.

OP049




Chuck Handle
round, serrated
10cm
7.5cm

OP050



Bard-Parker Handle
for disposable slot blades
#3 with gauge
#3

OP051



Barron Handle # 3

OP052



Curettes/Spatulas

Meyerhoefer

Chalazion Curette
 size 0, 1.5mm dia.
 size 1, 1.75mm dia.
 size 2, 2.25mm dia.
 size 3, 2.5mm dia.
 size 4, 3.5mm dia.

OP053



Double End Curette

size 1.5mm / 2.5mm dia.
 size 2mm / 3mm dia.
 size 2.5mm / 3.5mm dia.

OP054



Capsule

Polishing Curette
 1.75mm diameter open ring
 Nightingale,
 semi-sharp, posterior edge
 Shepherd,
 semi-sharp, anterior and
 posterior edges

OP055



Culler

Iris Spatula
 1mm wide
 2mm wide
 nettleship iris reposer

OP056



Wheeler

Double-End Spatula
 lightly curved

OP057



Barraquer

Iris Spatula
 round, smooth
 0.5mm diameter
 0.25mm diameter

OP058



Buratto

LASIK Flap Protector
For retracting and
protecting nasally hinged
corneal flaps.

OP059



Buratto

LASIK Flap Protector
blades angled 15° For retracting
and protecting superiorly
hinged corneal flaps.

OP060



Mannis - Buratto

LASIK Flap Protector
blades angled 45°
For retracting and protecting
superiorly hinged corneal flaps.

OP061



Rowen

Combination LASIK Flap
Protector and Spatula

OP062




Dalk

Corneal Dissector

12mm long, flattened,
vaulted spatula with
beveled tip

Rosenwasser DALK Trisector
with triangular shaft and
spatulated tip


OP063



Endothelial Stripper

1mm wide, thin beveled tip
smooth surface
sosenwasser with textured surface

OP064



Sloane

LASEK Flap Repositor
gently curved spatulated tip
with beveled smooth edges
for repositing and smoothing
the epithelial flap after laser
ablation

OP065



Camellin style

LASEK Trephine &
alcohol retaining well
8mm trephine and 8.5mm well

OP066



Camellin Style

LASEK Alcohol Well
For retaining alcohol solution
on the corneal epithelium.
8.5mm diameter S.S
8.5mm diameter titanium
9.5mm diameter S.S
9.5mm diameter titanium


OP067



Anis

Capsule Scrubber flexible
wire textured disc


OP068



Moran

LASIK Enhancement
Spatula Fukasaku spatula
and a modified paton spatula


OP069



Corneal Dissector

for intrastromal dissection.
angled, straight blade
angled, curved blade

OP070



Fogla Pointed Dissector

Fogla Pointed Dissector

OP071



Carlson DSEK Smoother

Carlson DSEK Smoother

OP072



Fogla Trifacet Spatula

Fogla Trifacet Spatula

OP073




Lasik/Lasek/Dalk

Hersh

Flap spatula beveled tips
a combination of a 6mm
spatula and 10mm spatula

OP074



Binder

Flap elevator 1mm wide 12mm
semi-sharp edges

OP075



Fukasaku

LASIK Spatula rounded anterior
and flat posterior surface with
beveled tip for lifting and
dissecting the corneal flap.

OP076



Guimaraes

LASIK spatula
gently curved, bi-convex
semi-sharp edges


OP077



Katzen

LASIK Flap "Unzipper"
combines a modified
sinskey hook and a
bi-convex spatula

OP078



Maddox

LASIK Spatula
combines a flat spatula and
a cylindrical spatula for
elevating and repositing the flap


OP079



Thorlakson

Flap Lifter combines a modified Sinskey hook and a bi-convex spatula with beveled notches.
long handle
short handle


OP080



Lasik Flap Lifter

Combines a 45° angled Sinskey hook and a thin hook with a tip
S.S
Titanium

OP081



Spatulas/Choppers/Manipulators

Hirschman

Lens spatula flat
0.5mm wide with
multiple notches

OP082



Jaffe

Lens spatula flat
0.5mm wide

OP083



Lester

IOL Manipulator
0.25mm diameter tip ideal
for multi-directional manipulation
straight, with guard

OP084



IOL Manipulator
0.25mm diameter tip
vaulted shank
with guard



Spatulas/Choppers

Phaco chopper

Chopper 10mm angle to tip
blunt cutting tip

S.S

Titanium

Chopper 10mm angle to tip
rounded front, 90° cutting edge.

S.S

Titanium

Delicate 10mm angle to tip
short tip, 90° cutting edge.

S.S

Titanium

OP085



Mics Phaco Pre-chopper

For bimanual nucleus pre-chopping

Alio-Rosen,
with two defined edges.

Alio-Scimitar, curved
with blunt tip and sharp inferior edge.

OP086



Milder

Phaco Chopper

60° offset inferior wedge
for horizontal chopping and a
beveled distal tip for vertical chopping

S.S for use in left hand

Titanium for use in left hand

Brierley

Nucleus Splitter

with incision sealing shaft

S.S

Titanium

OP087



Dodick - Kammann

Nucleus Chopper

for bimanual nucleus division.

1.5mm long tip, str

1.5mm long tip, angled

2mm long tip, straight

2mm long tip, angled

OP088



Escaf

Nucleus Divider

1.25mm long wedge shaped tip
for vertical chopping.

S.S

Titanium



Spatulas/Nucleus Manipulators


Gimbel
Nucleus Spatula
round, smooth
angled 30° at
10mm from the tip
0.25mm diameter

OP089




Knolle
Nucleus Spatula
malleable stainless steel,
flat 0.5mm wide

OP090




Koch
Nucleus Spatula
gently curved duckbill
shaped tip
smooth
with notches

OP091




Escaf
Nucleus Manipulator
flattened olive-shaped tip angled at 45°
S S
Titanium

OP092



Akahoshi
0.5mm diameter bulbous tip,
tip length 1.5mm 45° angled
shaft, tip to angle length 10.0mm
S S
Titanium



Escaf
Nucleus Sustainer
smooth ball tip textured inferior surface

Lieberman
MicroFinger Nucleus Manipulator
for use in right hand
for use in left hand



Spatulas/Choppers Manipulators

Neuhann
Nucleus Divider
strongly curved spatulated tip

OP093




Tennant
Nucleus Rotator
with 0.6mm diameter
ball shaped tip
S.S
titanium

OP094



Hunkeler
IOL Manipulator
with 0.3mm diameter
hemi-ball tip,
vaulted shank

Mushroom
IOL manipulator angled
with 0.65mm diameter ball tip
S.S
Titanium



Nucleus sustainer
Nucleus splitter and sustainer,
curved, blunt tip to be used
with the prechoppers

OP095



Scoring reverse hook
"Scoring" reverse hook for DSEK,
sharp, angled
DSEK stripper,
"T" shaped tips

OP096



Spatulas/Choppers Manipulators

OP097

Masket

a combination of the Rosen Phaco Splitter and a 12mm long round Phaco Spatula with disc shaped tip
S.S
Titanium

Chang

Phaco Chopper a combination of a modified MicroFinger and a 0.9mm long "Quick Chop" tip.
for use in right hand
for use in left hand

Milder

Phaco Chopper, double-end 60° for a beveled distal tip for vertical chopping and 12mm long spatula disc-shaped tip
for use in right hand
for use in left hand

Hwang

Phaco Chopper 1mm long vertical chopper and a chisel tipped spatula
S.S
Titanium

Warren

Combination of a 3.5mm long MicroFinger and 2mm long pointed "Quick Chop" tip for veterinary use only.

Rowen

Phaco Spatula a combination of the Rosen Phaco Splitter and the Bechert Nucleus Rotator
S.S
Titanium

Chang - Seibel

Phaco Chopper, a combination chopper with the Chang MicroFinger end and Seibel Vertical Chopper
for use in right hand
for use in left hand

Sharp Chopper & Blunt Chopper


A combination sharp chopper & blunt chopper
S.S
Titanium

Mushroom & Chopper

A combination of the mushroom and nagahara chopper
S.S
Titanium

Mushroom & Y Pusher

A combination of mushroom and Y pusher bechert
S.S
Titanium



Spatulas/Choppers Manipulators

OP098

Double Ended

Combination of 0.65mm mushroom tip and Nucleus Splitter
S.S
Titanium

Combination of nagahara blunt chopper and sinskey hook
S.S
Titanium

Combination of 0.65mm mushroom and sinskey hook
S.S
Titanium

Combination of nagahara blunt chopper and spatula
S.S
Titanium

Combination of sinskey hook and spatula
S.S
Titanium

Combination of 0.65mm mushroom tip and spatula
S.S
Titanium

Combination of 0.65mm mushroom tip and Y bechurt
S.S
Titanium


Combination of sinskey hook and Y bechurt
S.S
Titanium

Combination of nagahara blunt chopper and Y bechurt
S.S
Titanium

Combination of 0.65mm mushroom tip and sharp chopper
S.S
Titanium

Combination of sinskey hook and sharp chopper
S.S
Titanium

Combination of nagahara blunt chopper and sharp chopper
S.S
Titanium



Manipulators

OP099

Beichert

Nucleus Rotator blunt forked tip
 angled at 7mm S.S
 angled at 7mm titanium
 angled at 9mm S.S
 angled at 9mm titanium

Beichert

Lens Pusher vertical Y-shaped tip angled.

Drysdale

Nucleus Manipulator paddle shaped tip
 rounded edges S.S
 rounded edges titanium
 defined anterior S.S
 edge for capsule polishing

Pisacano

Nucleus Rotator flat disc shaped tip
 horizontal S.S
 horizontal titanium
 vertical S.S
 Vertical titanium

Siepser

Endocapsular Controller trumpet shaped tip

Anis

Nucleus Manipulator
 with 0.5mm diameter ball tip
 straight
 angled


Grayson

Nucleus Manipulator flattened,
 smoothly polished 1.1mm long
 with "tear-drop" shaped tip.

Lester

IOL Manipulator 0.25mm diameter
 tip Ideal for
 multi-directional manipulation.
 S.S
 Titanium

IOL Manipulator 0.25mm
 diameter tip vaulted shank
 For engaging the inferior
 haptic without touching the
 dome of the IOL.
 S.S
 Titanium



Hooks/Manipulators


Reversed Sinskey Hook

0.2mm diameter blunt tip for scoring the recipient endothelium and positioning the donor lamella.

OP100

Sinskey

Lens Manipulating Hook round handle
0.25mm diameter blunt tip
S.S straight, with guard
S.S angled
Titanium angled



Sinskey II

Lens Manipulating Hook
0.2mm diameter blunt tip round handle
S.S straight, with guard
S.S angled
Titanium angled



Lewicky

Lens Manipulating Hook
0.15mm dia blunt tip ang vaulted shank
S.S
Titanium



Hirschman

Iris Hook blunt, angled

Bonn

Micro Iris Hook blunt, flat handle

Shepard

Iris Hook blunt

Lens Manipulating Hook

0.15mm diameter blunt tip
Jaffe, straight
Clayman-Jaffe, angled

OP101



Sinskey

Lens Manipulating Hook
0.25mm diameter blunt tip straight



Harris - Sinskey

Lens Manipulating Hook
0.25mm diameter blunt tip angled



Sinskey II

Lens Manipulating Hook
0.2mm diameter blunt tip
straight
angled

Spatulas/ Manipulators

Lippman
Nucleus Spatula
0.5mm diameter round
shaft with blunt forked tip

OP102



Girard
Synechia Spatula
0.5mm wide with
semi-sharp tip

OP103



Obstbaum
Synechia Spatula
0.25mm wide with
semi-sharp tip

OP104



Benedetto
Spatula 10mm long,
smoothly polished
flattened spatula
with beveled notches

OP105



Elsching
Cyclodialysis Spatula
flat, 1mm wide.

OP106



Spatulas/ Manipulators

Castroviejo
Synechiaie Spatula
0.5mm wide with
semi-sharp tips.

OP107



Castroviejo
Cyclodialysis Spatula
0.25mm wide
0.5mm wide


Kimura
Platinum Spatula

OP108



Nucleus Dividing and Separating
size: 0.6x1.2mm tip oval
size: 0.4x0.65mm tip oval
size: 0.4x0.65mm tip oval

OP109



Straight Ring Manipulator

OP110

Straight shaft
0.5mm inside diameter, 0.8mm outside diameter, rounded tip
Round handle, length 121mm

0.5mm inside diameter, 0.8mm outside diameter, rounded tip
Angled shaft, tip to angle length 10mm
Round handle, length 118mm

0.8mm diameter, rounded tip
Angled shaft, tip to angle length 10mm
Round handle, length 118mm



Lens Loops

Keener
Lens Divider
with three preformed snare
wire loops
complete
wire loops

OP111



Lens Loop

OP112



Kirby
Lens Loop and
Muscle Hook

OP113



Morrison
Lens Loop and
Probe

OP114



Lens Loop
serrated, 3mm wide
Kansas, light curve
Keener-Arlt, strong curve

OP115



Lens Loops/Spoons

Alfonso

Cutting Platform
Used as a cutting
block for nucleus
trisection.
2mm wide
3mm wide

OP116



Giannetti

Nucleus Capture Spoon

OP117



Culler

Lens Spoon
sterling silver

OP118



Paton

Spatula & Spoon
For corneal
button transfer


OP119



Bunge

Evisceration Spoon
small
large

OP120



Spoons/Hooks

Wells
Enucleation
Spoon

OP121




Graether
Collar Button
micro iris retractor
straight, with guard
angled

OP122



Keuch
Small Pupil Dilator
with guard for enlarging
a small pupil through the
sideport and primary incision.

OP123



ICL Manipulator

OP124



Bonn
Micro Iris Hook
extra delicate blunt,
round handle with guard.

OP125



Hooks/Manipulators

OP126

Kuglen Iris & Hook Lens Manipulator


- Push-pull model, flat handle, straight
- Push-pull model, flat handle, angled
- Clover leaf model, flat handle, straight
- Push-pull model, round handle, straight w/guard
- Push-pull model, round handle, angled



Lens Manipulator

- straight, with guard
- Rentsch, "Boat Hook"
- Osher, Y-hook


OP127



Lens Manipulating Hook

0.15mm diameter tip


OP128



Deltz

0.15mm diameter tip

OP129



Alio - Kelman

IOL Hook 0.2mm diameter tip



Hooks/Fixation Rings

Fixation Hook
double, sharp

OP130



Twist Hook
for scleral fixation
right
left

OP131




Silverman
Fixation Pick
blunt cone shaped tips
spaced 3mm apart
for atraumatic scleral fixation

OP132



Globe Ring
Globe Fixation Ring
with concentric grooves
for refractive procedures
for cataract procedures

OP133



Lu - Mendez
LRI Guide and Fixation Ring
with delicate blunt teeth marks 10°
increments 11mm I.D.
S.S
Titanium

OP134



Gimbel

Stabilization Ring
with swivel handle
13mm diameter with blunt teeth
on both sides
For stabilizing either the right or
left eye during LASIK.

OP135



Fine - Thornton

Swivel Fixation Ring
blunt teeth features marks on anterior
surface in 15° increments as a guide
for placement of limbal relaxing incisions.
DT30-6161 13mm diameter titanium
13mm diameter S.S
16mm diameter S.S

OP136



Thornton

Globe Fixation Ring
multiple blunt teeth
16mm diameter ring S.S
16mm diameter ring titanium

OP137



Hofmann - Thornton

Globe Fixation Ring
with swivel handle
16mm diameter ring
multiple blunt teeth

OP138



Muscle Hooks

Stevens
Tenotomy Hook
blunt, small

OP139




Helveston
"Teaser" Hook
extra delicate, blunt

OP140



Graefe
Muscle Hook
size 1
size 2
size 3

OP141



Gass
Retinal Detachment
Hook with oval hole

OP142



Helveston
Muscle Hook
extra delicate,
with blunt spatulated tip
large
small

OP143



Muscle Hooks

Helveston

Muscle Hook
extra delicate, with blunt
spatulated tip
8mm
10mm
12mm

OP144



Jameson

Muscle Hook
small
large

OP145



Du Plessis

Muscle Hook
double arm

OP146



Guyton

Small Incision
Muscle Hook

OP147



Green

Muscle Hook

OP148



Moorfields

Moorfields
Strabismus Hook

OP149



Corneal Markers

Mendez

Degree Gauge 12mm I.D.
marked in 10° increments
for accurate axial
alignment on the cornea
S.S
Titanium

OP150



Bores

Axis Marker
S.S
Titanium


OP151



Henderson

Degree Gauge 11.7mm I.D.
Beveled face with markings
in 5° increments 4 notches
on internal edge to align
with cardinal positions
Titanium

OP152



Henderson

Toric IOL Marker
two 2.5mm radial
marks extending to 11mm
(used with ST30-7904)
Titanium

OP153



Henderson

Henderson Alignment Marker
with atraumatic disc-shaped
marking surfaces
at 3, 6 & 9 o'clock
Titanium

OP154



Optic Zone Markers

Bores	Hoffer, with crosshairs
3.0mm	3.0mm
3.25mm	3.25mm
3.5mm	3.75mm
4.0mm	4.0mm
4.25mm	4.25mm
4.5mm	4.5mm
4.75mm	4.75mm
5.0mm	5.0mm
6.0mm	5.5mm
7.0mm	6.0mm
8.0mm	7.0mm
9.0mm	8.0mm

OP155



Anis

Radial Marker
4mm diameter center opening
12.5mm O.D.
8 blades
12 blades
16 blades

OP156



Anis

Suture Placement Marker
8 radial blades
8mm diameter center opening
for corneal grafts.

OP157



Fechtner

Trabeculectomy Marker

OP158



MARKING PATTERN



Corneal Markers

McDonald

Optic Zone Marker
double end, with cross hairs
6.5 and 7.0mm diameters

OP159



Osher Neumann

Osher-Neumann
Corneal Marker
4 Line
8 Line
12 Line

OP160



PRK Scraper

combination of
PRK Scraper and Marker
6.5mm marker
total length 112mm

OP161



Double Ended Marker

5mm diameter on one side
5.5mm opposite

OP162



Maloney

Intra-Operative
Keratometer titanium
for qualitative measurement
of astigmatism.
S.S
Titanium

OP163



Castroviejo

Caliper
20mm spread
straight
angled
curved


OP164



Helveston

Scleral Ruler
15mm long, with notches
in 5mm increments single
marking tooth at tip for
measuring and marking the
distance from the limbus to
the location of the muscle
reinsertion.

OP165



Jameson

Caliper
chrome plated brass


OP166



Ruler

Stainless Steel
inches and mm

OP167



Braunstein

Caliper
3.5/4.0mm
Stainless steel

OP168



OP169

PRECISE, CONTROLLED DELIVERY, TIME AFTER TIME



- Rounded, Capsular-friendly dimensions assist lens manipulation post-implantation
- Distinctive blue color designed for easy visibility in the cartridge and reduced microscope reflections

CARTRIDGE LOCK

Quick and easy snap-in design, for secure cartridge loading, assures correct position of cartridge every time.



SCREW-STYLE

Consistent insertion forceps result in predictable lens delivery from start to finish.

High quality titanium material and durable construction for safe, reliable implantation

Vitreotomy Scissors

"Squeeze Handle" Scissors
these "squeeze handle"
instruments feature titanium
handles to reduce overall
weight, and front ends
of hardened stainless steel
for secure cutting

OP170



Scissors
horizontal, straight
20G 0.8 S.S
20G 0.8 titanium
23G 0.6 S.S
23G 0.6 titanium



Scissors
horizontal, curved
20G 0.8 S.S
20G 0.8 titanium
23G 0.6 S.S
23G 0.6 titanium
25G 0.50 S.S
25G 0.50 titanium



Scissors
horizontal, angled 45°
20G 0.8 S.S
20G 0.8 titanium
23G 0.6 S.S
23G 0.6 titanium



Scissors
vertical opening straight
20G 0.8 S.S
20G 0.8 titanium
23G 0.6 S.S
23G 0.6 titanium



Scissors
vertical opening angled 60°
20G 0.8 S.S
20G 0.8 titanium



Scissors
vertical opening angled 80°
20G 0.8 S.S
20G 0.8 titanium

Alio

MICS scissors
23 gauge curved shaft with horizontal
micro blades for cutting delicate membranes
Titanium

OP171



Horizontal blade



Vertical blade

YAE

Multi - Purpose 23G
vertical blade
3mm blade 24mm long shaft 112mm titanium
2.5mm blade, 15mm long curved shaft 100mm titanium

OP172



Scissors

Castroviejo
Corneoscleral Scissors
small blades
right
left

OP173



Castroviejo
Corneoscleral Scissors
medium blades
right
left

OP174



Troutman - Castroviejo
Corneoscleral Scissors
with stop, small blades
right
left

OP175



Troutman - Castroviejo
Corneoscleral Scissors
with stop, medium blades
right
left

OP176



Castroviejo
Corneal Transplant Scissors
with stop
strong curve, small blades
right
left

OP177



Retinal Scissors

Rich

Corneal Transplant Scissors
extremely delicate, strongly
curved blades for trimming the
recipient corneal bed and donor
button. Can also be used for open
sky anterior capsulotomy.
right
left

OP178



Terry

Micro Corneal Scissors
round handle for
meticulous control.
right
left

OP179



Katzin

Corneal Transplant Scissors
strong curve, medium blades
right
left


OP180



Castroviejo

Universal Corneal Scissors
curved, blunt tips
small blades
medium blades
large blades
strongly curved
large blades

OP181



Castroviejo

Keratoplasty Scissors
angled to side
blunt tips, medium blades

OP182



Scissors

Corneal scissors
curved, 10mm long blades,
pointed tips curved
blunt tips



OP183




Conjunctival scissors
curved, 10mm long blades,
pointed tips curved
reverse blades
blunt tips



OP184



Conjunctival scissors
curved, 11mm long blades,
blunt tips



OP185



Conjunctival scissors
curved, 8mm long blades,
blunt tips



OP186



IOL removal scissors
Cutting scissors for removal of
IOL straight, blunt tips, 8mm long
with serrated lower blade

OP187



IOL removal scissors
Cutting scissors for removal of
IOL angled, blunt tips, 8mm long
with serrated lower blade

OP188



McPherson - Westcott OP189

Conjunctival Scissors curved
blunt tips, small blades.



Westcott OP190

Tenotomy Scissors
blunt tips, medium blades
curved.



Westcott OP191

Tenotomy Scissors
blunt tips, standard blades curved
right
left



Tibolt OP192

Punctal Scissors
one conically shaped blade
with micro grooves and
textured surface for simple
"Two-snip Punctoplasty"



Mini Westcott OP193

Tenotomy Scissors
4mm long curved blades
blunt tips.



Scissors

McPherson - Westcott

Stitch Scissors
very sharp pointed
tips small blades curved.

OP194



Westcott

Stitch Scissors
very sharp pointed tips
standard blades curved

OP195



Jaffe

Stitch Scissors
very sharp pointed tips
medium blades curved

OP196



Iris

Scissors
sharp tips, small blades
straight
curved

OP197



Fine

Stitch Scissors
sharp tips, medium blades
straight
curved

OP198



Uthoff

Capsulotomy scissors,
convex curved,
7 mm long, thin blades,
pointed tips, right

OP199



Fogla Scissors Right

The lower blade is thinner and
1mm longer than the upper blade

OP200



Multipurpose scissors Noyes

Multipurpose scissors Noyes,
straight line, blunt tips

OP201



Corneal scissors

Corneal scissors, strong, angled,
8 mm long blades, blunt tips

OP202



Scissors

Barraquer

Iris Scissors
blunt tips, 7mm long blades.


OP203



Capsulotomy Scissors Vannas

5mm long blades, sharp tips
straight
curved
angled forward
angled to side

OP204



Gills - Vannas

7mm long blades, sharp tips
straight
curved
angled forward




Mcperson - Vannas

8mm long blades, sharp tips.
curved

Capsulotomy Scissors Vannas

7mm extra thin long blades, sharp tips
straight
curved
angled forward
10mm extra thin long blades, sharp tips
straight
curved
angled forward
curved, angled, forward
10mm long angled blades, curved tips
sharp points
blunt points

OP205



Eye Scissors

pointed tips small
straight
curved

OP206



Eye Scissors

pointed tips
standard
straight S.S
curved S.S
Tungsten Carbide Str
Tungsten Carbide Cvd
Super Cut Str
Super Cut Cvd

OP207



Eye Scissors

large rings
pointed tips
standard
straight
curved

OP208



Stitch Scissors

ribbon type
needle points curved

OP209



Stitch Scissors

needle points
light curve

OP210



Scissors

Stevens

Tenotomy Scissors
blunt tips
standard
straight S.S
curved S.S
Tungsten Carbide Str
Tungsten Carbide Cvd
Super Cut Str
Super Cut Cvd


OP211



Stevens

Tenotomy Scissors
ribbon type blunt tips
straight
curved

OP212



Knapp

Strabismus Scissors
blunt tips
straight S.S
curved S.S
Tungsten Carbide Str
Tungsten Carbide Cvd
Super Cut Str
Super Cut Cvd


OP213



Strabismus Scissors

ribbon type
blunt tips
straight
curved

OP214



Metzenbaum

Dissecting Scissors
gently curved blades
blunt tips
5 3/4" (14.5cm)
S.S
Tungsten Carbide
Super Cut

OP215



Utility Scissors

one pointed and one
blunt tip
5 3/4" (14.5cm) long
straight
S.S
Tungsten Carbide
Super Cut


OP216



Enucleation Scissors

blunt tips
medium curve
strong curve

OP217



Black Coated Super Cut Scissors

black coated super cut Iris scissors
Str
Cvd

OP218



Forceps

ICL forceps
ICL forceps for
loading, angled tips
S.S
Titanium

OP219



Foldable Lens Cutter Squeeze
Titanium Handle 112mm

OP220



IOL Lens Cutter
IOL Lends cutter/Holder with one serrated
blade stainless steel
S.S

OP221



Foldable Lens Removing Forceps
Titanium Handle 110mm

OP222



Vitreo - Retinal Forceps

Stainless Steel Autoclavable

OP223




Stainless Steel Autoclavable



Kuhn Forceps 45°
Power Gripping
23 gauge



Asymmetrical
End-Gripping
20 gauge



Asymmetrical
End-Gripping, 1x2 teeth
20 gauge, straight



Subretinal Forceps
Angled 135°,
4.5mm from the tip
20G/0.9



Serrated jaws
straight
20G/0.9
23G/0.6



ILM Forceps
2.2mm long with
shaft size 27.5mm
20G/0.9



End gripping forceps
straight
20G/0.9
23G/0.6



Serrated ILM
Forceps, 1.25mm
Long – with fine pick.
20G/0.9
23G/0.6



Spaide Pic forceps,
20 G/0.9
23 G/0.6




Asymmetrical
End-gripping Forceps,
20 G/0.9
23 G/0.6



Eckardt End-gripping
Forceps,
20 G/0.9
23 G/0.6



ILM forceps,
20 G/0.9
23 G/0.6



Foreign Body Forceps
Basket style with three
gripping arms.
20 G/0.9

Titanium Retinal Forceps

Vitreo - Retinal Forceps

OP224




Titanium



Kuhn Forceps 45°
Power Gripping
23 gauge titanium



Asymmetrical
End-Gripping
20 gauge titanium



Asymmetrical
End-Gripping, 1x2 teeth
20 gauge, straight
titanium



Subretinal Forceps
Angled 135°,
4.5mm from the tip
titanium



Serrated jaws
straight titanium



ILM Forceps
2.2mm long with
shaft size 27.5mm
20G/0.9 titanium
23G/0.6 titanium



End gripping forceps
straight

20G/0.9 titanium
23G/0.6 titanium



Serrated ILM
Forceps, 1.25mm
Long – with fine pick.
20G/0.9 titanium
23G/0.6 titanium



Spaide Pic forceps,
20 G/0.9 titanium
23 G/0.6 titanium



Asymmetrical
End-gripping Forceps,
20 G/0.9 titanium
23 G/0.6 titanium



Eckardt End-gripping
Forceps,
20 G/0.9 titanium
23 G/0.6 titanium



ILM forceps,
20 G/0.9 titanium
23 G/0.6 titanium



Foreign Body Forceps
Basket style with three
gripping arms.
20 G/0.9 titanium

Titanium Retinal Forceps

OP225

Micro Capsulorhexis Forceps

Capsulorhexis forceps 4mm long jaws with angled, grasping platform. Curved shaft 24mm long.
 21G/0.8 S.S
 21G/0.8 titanium
 23G/0.6 S.S
 23G/0.6 titanium



Ikeda Capsulorhexis forceps
 0.6mm (23G) Delicate 1.2mm long jaws with angled, grasping platform.
 Curved shaft 24mm
 S.S
 Titanium



Ikeda angled micro Capsulorhexis
 0.9mm long jaw with 90° angled, grasping platform.
 Strongly curved shaft 24mm
 S.S
 Titanium

23 G, MICS/MCCS Coaxial Capsulorhexis Forceps

MICS/MCCS Capsulorhexis 23G
 1.3 mm jaws, for 1 mm incisions
 Titanium

OP226



MICS capsulorhexis forceps

23G Ikeda capsulorhexis forceps 0.65 mm delicate 1.2 mm long jaws with angled, grasping platform, curved shaft
 titanium

OP227

23G komaki micro capsulorhexis forceps universal capsulorhexis forceps 0.7 mm long jaws with 90° angled grasping platform
 titanium



Botvin

Iris Forceps
For peripheral iridectomies.
serrated tips
1x2 teeth

OP228



Colibri Style

very delicate, 1x2 teeth, 0.12mm
Colibri S.S
Harms-Colibri, with tying platform S.S
Harms-Colibri, titanium, with tying platform titanium
Dyson-Colibri, short tips, extremely thin shanks,
with tying platform Standard, 1x2 teeth with tying platform S.S
Colibri, 0.4mm Corneal Suturing Forceps delicate,
1mm spread, 1x2 teeth, 0.12mm S.S
Polack Corneal Forceps
very delicate, 1x1 teeth, 0.1mm S.S
Pierse-Colibri S.S

OP229



Gills - Colibri

Tissue Forceps
very fine pointed tips
with tying platform

OP230



Colibri Forceps

Colibri corneal forceps,
delicate, 1x2 straight,
with tying platform
0.12 mm
0.40 mm

OP231



Harms

Hoskin

Tissue Forceps, straight
#18, fine grooved tips
#19, micro grooved tips
#28, very delicate, micro
grooved tips


OP232



Hoskin

Tissue Forceps,
Colibri style
#22, micro grooved tips
#30, very delicate micro
grooved tips
#34, pointed tips


OP233



Pierse

Corneal Forceps
very delicate 1x1 teeth,
0.1mm straight


OP234



Fechtner


Conjunctiva Forceps
very delicate ring shaped
tips with tying platform

OP235



Castroviejo - Colibri

OP236



Hoskin

Tissue Forceps, curved
#16, fine grooved tips

OP237



Forceps

Bonn

French pattern bonn forceps
0.1 mm with tying S.S
11cm long handle S.S
7.5cm short handle Titanium

OP238



Bonn

Forceps
very delicate 1x2 teeth, 0.12mm
with tying

OP239



Bonn

Forceps
1x2 teeth, 0.12mm
with tying platform
S.S
Titanium

OP240



Thorpe

Corneal Forceps
2x3 teeth, 0.12mm
with tying platform


OP241



Colibri Corneal Forceps

Colibri Corneal Forcep
with tying platform titanium
0.10 1x2 titanium
0.30 1x2 titanium

OP242



Bonn Corneal Forceps

Bonn Corneal Forcep
delicate, 1x2 straight 0.12mm teeth
with tying platform

OP243



Corneal Forceps

Corneal Forcep
delicate, 1x2 straight 0.12mm teeth
with tying platform, ergonomic handle

OP244



Corneal Forceps

Corneal Forcep 1x2 delicate w/tying
platform titanium ergonomic handle
0.10 titanium
0.30 titanium
0.10 oblique titanium
0.30 oblique titanium

OP245



Corneal Forceps

Corneal Forceps 90mm titanium
0.12 1x2 titanium
0.3 1x2 titanium
0.12 1x2 oblique titanium
0.3 notched titanium
Plain tying forceps
5mm titanium

OP246



Kelmann McPherson

Kelmann Forcep 90mm titanium
tying Ang 9mm titanium
tying Ang 7mm titanium
tying Ang 11mm titanium

OP247



Forceps

Iris forceps

Iris forceps, 1x2 straight,
0.20 mm teeth, with tying platform

OP248



Corneal forceps

Corneal forceps, delicate, 1x2 straight
0.12 mm teeth, with tying platform


OP249



Conjunctival and Fixation Forceps

Conjunctival and fixation, forceps
1x2 straight, 0.50 mm teeth

OP250



Scleral forceps, straight,
1x2 straight 0.70 mm teeth, with tying platform

Conjunctival and Fixation Forceps

Conjunctival and fixation, forceps
1x2 straight, 0.50 mm teeth

OP251



Suturing forceps

Suturing forceps, straight,
1x2 oblique teeth, with tying platform
0.12mm
0.35mm
0.60mm
0.90mm

OP252



Jayle Corneal Tying Forceps

0.15mm teeth
0.3mm teeth




OP253



Jayle Corneal Tying Forceps

0.12mm teeth titanium
0.3mm teeth titanium




OP254




St.Martin Corneal Tying Forceps

0.15mm teeth S.S
0.3mm teeth S.S
0.15mm teeth titanium
0.3mm teeth titanium



OP255



Bonn Micro Tying Forceps

corneal 0.12mm teeth
corneal 0.12mm teeth titanium
iris 0.12mm teeth S.S
iris
0.12mm teeth titanium



OP256



Mcpherson

Tissue Forceps

1x2 teeth, 0.4mm with tying platform


straight

angled

8mm long 1x2 teeth, 0.3mm

angled

OP257



Hoffer - Mcpherson

Forceps very delicate

11mm long jaws

1x2 teeth, 0.12mm

OP258



Castroviejo Suturing Forceps

1x2 teeth with tying platform

0.12mm S.S

0.3mm S.S

0.5mm S.S

0.9mm S.S

0.12mm, titanium

OP259



Castroviejo Suturing Forceps

1x2 teeth with tying platform

0.12mm S.S

0.3mm S.S

0.5mm S.S

0.9mm S.S

0.12mm, titanium

OP260



Castroviejo Suturing Forceps

1x2 teeth with tying platform

0.12mm S.S

0.3mm S.S

0.5mm S.S

0.9mm S.S

OP261



Moody

Fixation Forceps
1x2 teeth, 0.5mm
with lock
left
right

OP262



Castroviejo-Colibri

Castroviejo-Colibri
type corneal forceps 0.12mm, 1x2

OP263



Girard

Corneoscleral Forceps
1x2 teeth, 0.12mm
with tying platform
right hand
left hand

OP264



Stern-Castroviejo

Fixation Forceps
1x2 teeth, 0.5mm with lock

OP265



Bores

Corneal Fixation Forceps
U - Shaped, 3mm spread
1x2 teeth at each tip, 0.12mm

OP266



Bishop-Harmon

Tissue Forceps
1x2 teeth
delicate
standard
heavy

OP267



Bishop-Harmon

Dressing Forceps
criss-cross serrated tips
serrated 0.5
serrated 0.8

OP268



Fixation Forceps

1x2 teeth, 0.3mm
sauer
manhattan Eye & Ear
with tying platform

OP269



Dressing Forceps

Serrated, delicate
straight
curved

OP270



Tissue Forceps

1x2 teeth
delicate, straight
delicate, curved
standard, straight
standard, curved


OP271



Adson

straight, 5" (125mm) long
 serrated jaws, delicate S.S
 serrated jaws, standard S.S
 1x2 teeth, delicate S.S
 1x2 teeth, standard S.S
 serrated titanium
 1x2 delicate titanium

OP272



Bayonet

Dressing Forceps serrated jaws
 5 3/4" (146mm) long

OP273



Lester

Fixation Forceps straight
 1x2 teeth
 2x3 teeth

OP274



Fixation Forceps

1x2 teeth
 Bracken delicate teeth
 Elschnig, standard teeth

OP275



Troutman

Superior Rectus
 Forceps 1x2 teeth angled

OP276



Tying Forceps

For 8-0 to 11-0 sutures forceps

McPherson

5mm long tying straight S.S

5mm long tying straight, titanium

5mm long tying angled S.S

angled, 5mm long tying titanium

Kelman - McPherson forceps

7.5mm long tying angled S.S

7.5mm long tying angled, titanium

Bechert Forceps


10mm long tying angled S.S

Gaskin Forceps

11mm long tying angled S.S

11mm long tying angled Titanium

OP277



Haberle-McPerson

Lens Holding Forceps

7.5mm long jaws

with groove, angled

For grasping IOL haptic

OP278



Kurstin


Flap Stretching Forceps

disc shaped

serrated jaws

Designed for stretching the flap to eliminate striae.

OP279



LASIK Flap Forceps

disc shaped

serrated jaws

Inamura

Capsulorrhexis Forceps
delicate tapered shanks
to seal the incision during
continuous curvilinear
capsulorrhexis
Excellent for MICS

OP280



Utrata

Capsulorrhexis Forceps
very delicate, short handle
triangular grasping tips
extremely thin
with landmarks at 2.5mm and
5.0mm
straight

OP281



Masket

Capsulorrhexis Forceps
very delicate, short handle
triangular grasping tips
extremely thin
with landmarks at 2.5mm and
5.0mm
vaulted

OP282



Kraff

Capsule Tag Forceps
11mm long shanks
1.5mm long grasping tips

OP283



Utrata Capsulorrhexis Forceps
Sharp Tips, 80mm

OP284



Utrata

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin
11mm long shanks
S.S
titanium

OP285



Utrata

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin
11mm long shanks
round handle with guide pin
for alignment

OP286



Nevyas

Capsulorrhexis Forceps
cystotome shaped tips
extremely thin
11mm long shanks
For tearing and grasping
the anterior capsule.
S.S
titanium

OP287



Masket

Capsulorrhexis Forceps
very delicate
triangular grasping tips
extremely thin 11mm long
vaulted shanks
S.S
titanium

OP288



Inamura

Capsulorrhexis Forceps
round handle cross action
2 mm incision

OP289



Masket

Capsulorhexis Forceps
 very deliate triangular
 grasping tips extremely
 thin 11mm long vaulted
 shanks round handle with
 guide pin for alignment
 S.S
 titanium



OP290



Capsulorhexis Forceps

Corydon forceps
 ergonomic curved, cystotome tips,
 11.5 cm long



OP291



Capsulorhexis Forceps

Capsulorhexis Forceps
 angled, 12 mm graduated slightly curved
 arms, sharp tips, ergonomic handle
 11mm angled
 12mm angled



OP292



Capsulorhexis Forceps

Round handle
 11mm long shanks straight



OP293



Capsulorhexis Forceps

Round handle, curved,
 cystotome tips, 11.5cm long



OP294



Forceps

Inamura Capsulorhexis Forceps

Sharp pointed serrated tips angled 90° from shaft
 Curved shaft, tip to pivot point 8.5mm
 Cross action tips, 1.2mm width at pivot box
 Round handle, overall length 111mm desired size of capsulorhexis
 Titanium
 marks on shaft at 2.5mm and 5.0mm denote Titanium

OP295



Capsulorhexis forceps


Capsulorhexis forceps, lightweight, angled, 12 mm, slightly curved arms, sharp tips, titanium

Capsulorhexis forceps, lightweight, angled, 12 mm, slightly curved arms, strongly angled, sharp tips, titanium

Capsulorhexis forceps, lightweight, angled, 11 mm, slightly curved arms, notched tips, titanium

Capsulorhexis forceps, lightweight, angled, 11 mm, slightly curved arms, sharp, strongly angled tips, titanium

OP296



Capsulorhexis forceps

Capsulorhexis forceps, angled, 11 mm slightly curved arms, sharp, notched tips

OP297



Capsulorhexis forceps

Capsulorhexis forceps, angled, 11 mm slightly curved arms, sharp tips, ergonomic handle, titanium

OP298



Capsulorhexis forceps

Capsulorhexis forceps
flat handle, 11 mm sharp tips,
micro incision
titanium

OP299



MICS capsulorrhexis forceps

MICS capsulorrhexis forceps,
controlled opening for 2 mm incisions
titanium

OP300



Capsulorhexis forceps for micro coaxial phaco

Capsulorhexis forceps with 2mm
controlled opening for 2mm incisions
titanium

OP301



MICS/MCCS Capsulorhexis Forceps

MICS/MCCS Capsulorhexis
with 2 mm controlled opening,
for 2 mm incisions

OP302



Capsulorhexis Forceps

extremely delicate

OP303



Prechopper

Combo Prechopper,
for grades 1-4 nuclei,
2 mm opening, sharp top side


OP304




Prechopper

Akahoshi Super Combo III Prechopper
For nucleus grades 1 to 3.
is suitable for the vertical Karate
Prechop technique within a 2mm incision.
Prechopping can be completed without
nuclear support for nuclei grade 1 to 3


OP305



Akahoshi Hybrid Combo prechopper
For Grades 1-4 Nuclei



Akahoshi Universal Nucleus Splitter



Akahoshi

Phaco PreChopper
4.5mm long delicate jaws
with sharp tips
For cracking, dividing and rotating
the nucleus without sculpting

OP306



Akahoshi

Combo PreChopper
For cracking, dividing and rotating
both hard and soft nuclei without
sculpting

OP307



Anis

Tying Forceps
extra delicate
smooth jaws
For 9-0 to 11-0 sutures.
straight
curved

OP308



Freeman - Jaffe

Lens Holding Forceps
with longitudinal groove
For grasping IOL haptic.


OP309



Jaffe

Tying Forceps
very delicate
smooth jaws
6mm long platform
For 8-0 to 11-0 sutures.
straight S.S
curved S.S
straight titanium
curved titanium

OP310



Gerl

Tying Forceps
curved shanks
6mm long platform
For 8-0 to 11-0 sutures.
titanium


OP311



Tennant

Tying Forceps
extra delicate smooth jaws
with 6mm long platform
round handle with guide pin
For 9-0 to 11-0 sutures.
straight S.S
straight, titanium
curved S.S
curved, titanium

OP312



Maumenee

Tying Forceps
 concave-convex jaws
 round handle with guide pin
 For 8-0 to 11-0 sutures.
 straight titanium
 angled titanium
 straight S.S
 angled S.S

OP313




Jaws Cross Section



Harms

Tying Forceps
 standard, smooth jaws
 for 7-0 to 9-0 sutures.
 straight S.S
 curved S.S
 straight titanium
 curved titanium

OP314



Tuebingen Model

Tying forceps straight
 S.S
 titanium

OP315



Tuebingen Model

Tying forceps curved
 S.S
 titanium


OP316



McPherson Tying Forcep

Recommended for 8-0 to 10-0 sutures
 4mm tying straight
 4mm tying curved
 4mm tying angled
 8mm tying angled


OP317



Kelmann forcep

Kelmann forcep angled
ergonomic handle
11 mm tying
7 mm tying


OP318



Troutmann forcep

Troutmann forcep very delicate
ergonomic handle
7 mm tying str
7 mm tying cvd

OP319



Troutman forceps

Straight, delicate, bullet shape.
Curved, delicate, bullet shape.

OP320



IOL Loading forceps

Loading forceps for foldable lenses
into the cartridge, 11 mm angle to tip,
cylindrical arms
titanium


OP321



Golf Ball Tying Forcep

Golf ball tying forcep titanium
angled 8mm titanium
straight 5mm titanium
curved 5mm titanium
delicate str 5mm titanium

OP322




Forceps

Tying Forceps

Troutman forceps 6mm
straight
curved


OP323



Kelman McPherson

French pattern tying forceps
11mm angled


OP324



McPherson

Tying Forceps 5mm str
Tying Forceps 5mm angled
Tying Forceps 10mm angled


OP325



Fixation forceps

Fixation forceps, straight,
3mm
4mm

OP326



Green Fixation Forceps

Green fixation forceps

OP327



Helveston

Tying Forceps
standard, smooth jaws
For 4-0 to 6-0 sutures.

OP328



Fine

Tying Forceps
standard, smooth jaws
For 7-0 to 9-0 sutures.

OP329



Arita

Meibomian Gland Compressor
broad, smoothly polished jaws

OP330



Cilia Forceps with spring tension
S.S
titanium

OP331



Barraquer Cilia Forceps

OP332



Beer Forceps Cilia

OP333



Barraquer

Cilia Forceps lightly angled
Textured jaws provide
superior grasping capabilities.

OP334



Cilia Forceps

Cilia Forceps
slanted jaws gold plated

OP335



Gradle

Cilia Forceps smooth jaws


OP336



Bonaccolto

Conjunctiva Forceps
longitudinally serrated jaws
with cross serrations at tips
1.2mm wide
1.7mm wide

OP337



Giannetti

MICS Capsulorrhexis Forceps
ultra delicate shanks and tips
maximum jaw opening 1.75mm
For performing capsulorrhexis
through a micro incision.


OP338



Moorfields Conjunctival Forceps

standard S.S
fine S.S
standard titanium
fine titanium

OP339



Stallard Tying Forceps

OP340



Rycroft Corneal Tying Forceps

plain
0.3mm teeth

OP341



Scleral Plug Forceps

OP342



Vessel Dilator

11 cm long, flat handle 9 mm wide,
angulated 10°

OP343



Micro Forcep

15 cm long, round handle
6mm tying 1x2, 0.30 straight
6mm tying straight
6mm tying curved
vessel dilator

OP344



Nugent

Utility Forceps
smooth jaws
serrated jaws

OP345



Foreign

Body Forceps straight

OP346



Utility & Splinter

Forceps pointed, short jaws
Suitable for epilation.
#H

OP347



Jewelers

Forceps
standard, straight
#1

OP348



Jewelers

Forceps
delicate, straight
#3-C

OP349



Jewelers

Forceps
delicate, straight
#3

OP350



Jewelers

Forceps
delicate, straight
#4

OP351



Jewelers

Forceps
extra delicate, straight
#5

OP352



Forceps

Osher
Iris Forceps
very fine textured jaws
angled

OP353



Jewelers
Forceps
extra delicate, angled
#5/45

OP354



Jewelers
Forceps
curved
#7

OP355



Suture
Removal Forceps
very fine pointed tips with textured
jaws Ideal for grasping wet sutures

OP356



Jaffe
Utility Forceps
very delicate smooth jaws curved

OP357



Arruga
Capsule Forceps

OP358



Rigid IOL implantation forceps
Rigid IOL implantation forceps, angled

OP359



Implantation forceps
Implantation forceps
for 5.5 and 6.0 mm acrylic IOLs

OP360



Ernest
Nucleus Cracker
extra delicate cross-action
serrated paddle tips
For cracking the nucleus
through a 2.5mm incision

OP361



Ernest
Micro Nucleus Cracker
extra delicate cross-action
serrated paddle tips
For cracking the nucleus
through a 2mm incision.

OP362



Dodick

Nucleus Cracker
cross-action
longitudinally serrated
paddle shaped jaws
For cracking the nucleus
through a 3mm incision

OP363



Hu Femto Forceps

Blunted 0.3mm tips
cross-action mechanism
For opening incisions
created with a
femtosecond laser

OP364



Alfonso

Nucleus Grasping Forceps
two short rows of
delicate interlocking teeth
For grasping and extracting
nucleus fragments

OP365



Coaxial forceps

Coaxial forceps for endothelial
lenticule insertion in DSEK technique,
atraumatic

OP366



Corneal-scleral punch

Corneal-scleral punch,
straight, \varnothing 1.5 mm \varnothing 1.7 mm


OP367



Stripping forceps

Stripping forceps for the Schlemm's canal, straight

OP368



Stripping forceps

Stripping forceps for the Schlemm canal, curved tip

OP369



Dodick

IOL folding forceps for soft IOLs

OP370



Livernois

IOL folding forceps concave jaw surfaces for soft IOLs

OP371



Faulkner

Lens Holding Forceps
round, smooth jaws
For holding soft IOLs

OP372



Buratto

Soft IOL Inserting Forceps
extra delicate highly polished
convex jaws for foldable acrylic and
silicone IOLs

OP373



Spaleck

Soft IOL Inserting Forceps
delicate highly polished convex jaws
For foldable acrylic and silicone IOLs

OP374



Ernest-McDonald

IOL Inserting Forceps
polished bi-convex jaws
patented cross-action for
acrylic IOLs

OP375



JAWS CROSS SECTION



Bechert

Lens Holding Forceps
very delicate flat jaws
gently curved

OP376



Blaydes

Lens Holding Forceps
very delicate narrow jaws
angled 45°

OP377



Shepard

Lens Holding Forceps
duckbill shaped jaws
gently curved

OP378



Clayman

Lens Holding Forceps
delicate angled jaws
gently curved

OP379



Jameson

Muscle Forceps
with slide lock
4 teeth, child size, right
4 teeth, child size, left
6 teeth, adult size, right
6 teeth, adult size, left

OP380



Berke

Ptosis Forceps
with slide lock
longitudinally grooved jaws
20mm long
27mm long

OP381



Steinert

IOL Inserting Forceps
conically shaped jaws
with flat inside surfaces
Used for folding and
inserting soft IOLs.

OP382



Livernois

IOL Holding Forceps
highly polished round jaws
For holding a soft IOL
while folding with an
inserting forceps.
straight
angled

OP383



Livernois-McDonald

IOL Inserting Forceps
5mm long concave jaws
patented cross-action
For inserting silicone IOLs

OP384



Ernest-McDonald

IOL Inserting Forceps
5mm long straight jaws
low profile
patented cross-action
For inserting silicone IOLs

OP385



Ernest-McDonald II

Lens Cartridge Loading Forceps
straight, 20 gauge jaws feature a
series of concave and convex surfaces


OP386



Chalazion Forceps

round
Lambert, 8mm I.D.
Hunt, 12mm I.D.


OP387



Chalazion Forceps

oval
Baird, 11 x 8mm I.D. small
Francis, 14 x 11mm I.D. medium
Desmarres, 20 x 13mm I.D. large

OP388



Desmarres Forceps

Oval
25mm

OP389



Watzke

Sleeve Spreading Forceps

OP390



Forceps

Hartmann
Mosquito Forceps
extra delicate
3 3/4" (95mm) long
straight
curved

OP391



Halsted
Mosquito Forceps
delicate
5" (125mm) long
straight
curved

OP392



Halsted
Mosquito Forceps
delicate 1x2
5" (125mm) long
straight
curved

OP393



Backhaus

Towel Clamp
3½" (90mm) S.S
5¼" (135mm) S.S
3½" (90mm) titanium
5¼" (135mm) titanium

OP394



Baby Jones

Towel Clamp
Cross Action
(55mm) long

OP395



Schaedel

Towel Clamp
Cross Action
(90mm) long

OP396



Serrefine

serrated jaws
small (38mm)

OP397



Dieffenbach Serrefine

Straight, 50 mm
curved, 50 mm

OP398



Forceps


Pean Forceps
Mosquito forceps,
straight, serrated jaws

OP399



Ramplays Forceps
Ramplay sponge forceps
18cm
25cm

OP400



Frazier
Suction and wash tube,
with hole on/off, 12cm long

OP401

Suction and wash tube,
with hole on/off, 17cm long



Needle Holder

McPherson
Needle Holder
very delicate jaws
with lock straight
without lock straight
with lock curved
without lock curved

OP402



Jaffe
Needle Holder
delicate conical jaws
with lock straight
with lock curved
without lock curved

OP403



Castroviejo
Needle Holder
delicate jaws straight
with lock
without lock

OP404



Castroviejo
with lock
without lock

OP405



Silcock Needle Holder, catch

OP406



Needle Holder

Castroviejo

Needle Holder
delicate jaws
curved
with lock
without lock

OP407



Castroviejo

Needle Holder
standard jaws
straight, with lock
S.S
Tungsten Tip TC

OP408



Castroviejo

Needle Holder
standard jaws curved
with lock
without lock
with lock
without lock

OP409



Kalt

Needle Holder
delicate smooth jaws
straight, with lock

OP410



Drews

Needle Holder
extra delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock

OP411



Needle Holder

Castroviejo

Needle Holder
delicate jaws titanium
with lock str
with lock cvd
without lock cvd
with lock standard

OP412



Barraquer

Needle Holder
short model
very delicate jaws
curved
with lock
without lock

OP413



Troutman

Needle Holder
delicate jaws
strongly curved
without lock
with lock

OP414



Barraquer

Needle Holder
very delicate jaws
curved
S.S with lock
S.S without lock
titanium, with lock
titanium, without lock

OP415



Barraquer

Needle Holder
delicate jaws
straight
with lock
without lock

OP416



Needle Holder

Barraquer

Needle Holder
delicate jaws
curved
with lock
without lock

OP417

Needle Holder
delicate jaws
curved large
with lock
without lock



Barraquer

Needle Holder
standard jaws
curved
with lock S.S
without lock S.S
Tungsten Carbide with lock
Tungsten Carbide without lock
with lock titanium
without lock titanium

OP418



Anis

Needle Holder
extra delicate jaws
with lock straight
without lock straight
with lock curved
without lock curved

OP419



Micro Needle Holder

Needle Holder
extra delicate jaws
straight
with lock S.S
without lock S.S
titanium, with lock
titanium, without lock

OP420



Needle Holder

Barraquer (Golf Pattern)

Needle Holder
short model
very delicate jaws
Curved
with lock
without lock

OP421



Barraquer

Needle Holder
short model
very delicate jaws
curved
with lock
without lock

OP422



Micro Needle Holder

Needle Holder
delicate jaws straight
without lock
with lock

OP423



Needle Holder

delicate,
curved, without lock,
overall length: 115 mm, titanium
Ideal for 9-0 to 10-0 sutures
S.S with lock
S.S without lock
titanium, with lock
titanium, without lock

OP424



Barraquer

Needle Holder
delicate jaws
straight
S.S with lock
S.S without lock
with lock, titanium
without lock, titanium

OP425



Lacrimal Instruments

Kerrison

Bone Nibbling Rongeur
9mm opening
size 0, 3.0mm wide
size 1, 3.5mm wide
size 2, 4.0mm wide
size 3, 4.5mm wide

OP426



Citelli

Bone Nibbling Rongeur
7mm opening
2.0mm wide

OP427



Beyer

Rongeur
single action

OP428



Lacrimal Instruments

Nasal
Spaculum
infant
child
adult

OP429



2/3

Freer
Periosteal Elevator
one end blunt
other end semi-sharp
4.5mm wide blades

OP430



2/3

Joseph
Periosteal Elevator
gently curved tip

OP431



3/4

West
Lacrimal Chisel
straight

OP432



3/4

Freer
Lacrimal Chisel
curved

OP433



Lacrimal Instruments

Mallet
lead filled
stainless steel head

OP434



Harms
Trabeculotomy Probe
right
left

OP435



Infant
Lacrimal Dilator
very fine blunt tip

OP436



Muldoon
Lacrimal Dilator

OP437



Wilder
Lacrimal Dilator
size 1, fine taper
size 2, medium taper
size 3, heavy taper

OP438



Lacrimal Instruments


Castroviejo
Lacrimal Dilator
double-end
fine and
medium taper

OP439



Bowman
Lacrimal Probe
malleable stainless steel
set of 6
size 0000 & 000
size 00 & 0
size 1 & 2
size 3 & 4
size 5 & 6
size 7 & 8

OP440



3/4

Catalano
Lacrimal Intubation Set
includes 2 blunt
0.7mm diameter probes
attached to 30.5cm long
silicone tubing
.012" I.D., .025" O.D.
set
silicone tubing only

OP441



3/4

Crawford
Lacrimal Intubation Set
includes 2 olive tip
0.6mm diameter probes
attached to 30.5cm long
silicone tubing
.012" I.D., .025" O.D.
set
silicone tubing only

OP442



3/4

Lacrimal Instruments

Crawford
Retrieving Hook
for Crawford lacrimal
intubation probes

OP443



Pigtail
Probe
with suture holes

OP444



Rollet Ruginé

OP445



Nettleship Lacrimal Dilator
Nettleship lacrimal dilator
single end

OP446

Nettleship lacrimal dilator
double end



Howarth Periosteal Elevator
Howarth Periosteal Elevator,
for DCR, 215mm length

OP447



Fison Retractor
Fison Retractor, for detachment

OP448



Cannulas

OP449



Lacrimal
Cannula
malleable tip
23 gauge

straight
curved

OP450




Air Injection
Cannula
angled at 7mm
30 gauge
27 gauge



Bottom View

OP451



Dalk
Air Injection Cannula
flattened, 27 gauge
vaulted 0.2mm
posterior port near tip
27 gauge



Top View

Bottom View

OP452



Dalk
DALK Cannula
flattened "cobra-shaped"
tip 0.2mm posterior port
27 gauge




OP453



Lasek
Alcohol Cannula
delicate olive shaped tip
angled at 10mm
25 gauge



OP454




Stevens
Sub-Tenon's Anesthesia
Cannula gently curved,
flat tip 19 gauge

OP455



Gess
Irrigating-Aspirating
Cannula U-shaped,
with end opening
aspiration port,
23 gauge Supplied
with silicone tubing
and adaptor
curved right
curved left


OP456



Simcoe
Irrigating-Aspirating
Cannula original model,
15mm long,
0.3mm diameter
aspiration port




OP457



Corbin
Sub-Tenon's Anesthesia
Cannula smooth rounded
tip 0.3mm side port
23 gauge

Cannulas

OP458




Reversed Simcoe
Irrigating-Aspirating
Cannula 15mm long
0.3mm diameter
aspiration port front
opening irrigation port
23 gauge

OP459



Simcoe
Irrigating-Aspirating Cannula
U-shaped, with 0.3mm anterior
aspiration port, 23 gauge
curved right
curved left

OP460



Strauss-Simcoe
MSICS I/A Cannula
23 gauge, 0.3mm port
elongated and gently curved
with 15° angled tip.



Kelman
Cystotome
sharp point and
sharp cutting edge
22 gauge



McIntyre
Cystotome
angled at 10mm
sharp point
23 gauge




Lieppman
Cystotome
extra small sharp point
and sharp cutting edge
22 gauge



Shammas
Posterior Capsule
Polishing Curette
1mm diameter cup
with textured base
22 gauge

OP465



Shapiro
Posterior Capsule
Polisher carbide
impregnated tip dual
0.3mm side ports
curved, 23 gauge

OP466



Kratz
Cystotome
sharp point and
sharp cutting edge
formed, 22 gauge



Jensen
Posterior Capsule
Polisher carbide
impregnated on posterior
surface of olive tip
front opening
25 gauge
23 gauge

Bishop-Harmon

Anterior Chamber Irrigator includes 19 gauge angled cannula, adaptor and blue silicone bulb complete silicone bulb and adaptor
silicone bulb and adaptor
19 gauge cannula
21 gauge cannula


OP468



OP469

Lu


Posterior Capsule Polisher flattened olive shaped tip carbide impregnated on the posterior surface front opening, 23 gauge



OP470

Dahan


Nucleus Prolapsing Cannula angled at 4mm miniature olive tip with 0.15mm end opening 25 gauge



OP471

Hydrodissection


Cannula angled at 8mm flat tip, front opening 25 gauge



OP472

Irrigating Cannula


dual 0.4mm side ports 23 gauge thin wall



OP473

Fukasaku


LASIK Cannula low profile tip, end opening 23 gauge



OP474

Banaji

LASIK Cannula bullet shaped tip 6 ports for multi-directional irrigation 25 gauge smooth textured



OP475

Buratto


LASIK Cannula slightly flattened one front and two side ports 25 gauge 23 gauge



OP476

McIntyre

Anterior Chamber Cannula smooth blunt tip, front opening 23 gauge, straight 23 gauge, angled 26 gauge, straight 26 gauge, angled 30 gauge, straight 30 gauge, angled




OP477

Gimbel


Cortex Aspirating Cannula end opening 30 gauge

Scissors Cannulas




OP478

Vidaurri
LASIK Cannula
double-armed
8 irrigating ports
25 gauge
vertical
horizontal




OP479

Simcoe
Cortex Extractor
0.3mm side port
13mm long, 23 gauge
smooth tip
textured tip




OP480

Simcoe
Cortex Extractor
0.3mm side port
21mm long, 23 gauge
smooth tip
textured tip




OP481

McIntyre-Binkhorst
Cannula, "J" shape
smooth blunt tip
end opening, 26 gauge
angled right
straight
angled left




OP482

Chang
Hydrodissection Cannula
flat tip angled 90°
with beveled opening
27 gauge
standard tip
long tip



OP483

Knolle-Pearce
Irrigating Vectis
fluid flow through
three front opening
ports 23 gauge loop



OP484

Simcoe
Irrigating Lens Nucleus Loop
two front opening ports
finely serrated, strong curve
25 gauge
right
left



OP485

Simcoe
Cortex Extractor
0.4mm side port
21mm long
23 gauge thin wall



OP486

Alfonso-McIntyre
Irrigating Nucleus Spoon
single 0.3mm
front opening port
2.5mm wide
4.0mm wide

I/A Bimanual Handles

Infusion Handle

Infusion Handle
male/male connectors S.S
male/female connectors S.S

OP487



Aspiration Handpieces

with male connector
side port facing into curve
magenta color titanium
titanium
23 gauge, 0.3mm port
23 gauge, 0.3mm port, textured tip
21 gauge, 0.35mm port
21 gauge, 0.35mm port, textured tip
21 gauge, 0.4mm port
21 gauge, 0.4mm port, textured tip
22 gauge, 0.3mm port, ova
22 gauge, 0.3mm port, oval, textured tip
20 gauge, 0.3mm port, oval

OP488



Irrigation Handpieces

with female connector, blue color
23 gauge, dual 0.3mm side ports
22 gauge, dual 0.4mm side ports, oval
21 gauge, dual 0.5mm side ports

OP489



Aspiration Handpieces Ø21G

Curved sandblasted tip
titanium handle

OP490



Irrigation Handpieces Ø21G

Curved tip, dual sideports,
titanium handle



OP491



Scoring Hook and Aspirating Handpiece

Scoring hook and aspirating handpiece
for DSEK procedure

OP492



2



Single Use
Ophthalmic
Instruments



Barraquer Speculum

SUOP01


Curved 10mm blade
Curved 15mm blade



Kratz Speculum

SUOP02


6.5mmx14mm, 1mm Dia Wire
6.5mmx14mm, 0.9mm Dia Wire



Barraquer Speculum

SUOP03

Flat blade 10mm
Flat blade 15mm



Barraquer Speculum

SUOP04

solid blade 10mm
solid blade 15mm



Infant Speculum

SUOP05

4.3mm x 4.3mm Closed Blades,
Ø0.7mm Wire. 25mm Wide



Neonatal Speculum

SUOP06

4.0mm x 10mm Closed Blades;
Ø0.8mm Wire. 25mm Wide
3.0mm x 7.0mm Closed Blades



Kratz Speculum

SUOP07


Kratz-Barraquer Speculum, Curved
6.5mm x 14.5mm Open Curved Blades,
Ø0.9mm Wire. 30mm Wide



Jaffe Lid Retractors

SUOP08

Jaffe Style wire speculum (Pair)
15mm blades (large)



Liberman Speculum

SUOP09

K Wire 10mm
K Wire 15mm
Round Wire 14mm for lasik
V Wire 15mm
Solid 10mm
Solid 15mm



Lancaster Speculum

Eye Speculum
adjustable mechanism
with locking nut
16mm blades

SUOP10



Bangerter Speculum

SUOP11

15mm blades
10mm blades



Bard-Parker Handle
 Chuck Handle round, serrated
 10cm
 7.5cm
 # 3
 for disposable slot blades

SUOP12



Nucleus Manipulators & Hooks



Sinskey Hook Kuglen Hook Rotator Y Shaped Musroom Manipulator Drysdale

SUOP13



Choppers & Manipulators



Spatula Nucleus Spatula round smooth Lester Manipulator iol Nagahara Chopper Blunt Rosen Splitter Claw Chopper

SUOP14



Double Ended



Phaco Chopper/Mushroom Sinskey Hook / Phaco Chopper

SUOP15



Kirby Lens Loop and Muscle Hook

SUOP16



SUOP17

Jameson Muscle Hooks
Small
Large
Desmarres Lid Retractor
Size 1, 11mm wide
Size 2, 15mm wide



Chalazion Forceps
round
Lambert, 8mm I.D.
Hunt, 12mm I.D.

SUOP18



Chalazion Forceps
oval
Francis,
14 x 11mm I.D. medium
Desmarres,
20 x 13mm I.D. large

SUOP19



Westcott Scissor
Tenotomy Scissors
blunt tips, curved

SUOP20



Westcott Scissor
Tenotomy Scissors
sharp tips, curved

SUOP21



Westcott Scissor
Tenotomy Scissors
blunt tips, medium blades
curved

SUOP22



Corneal Scissor
Corneal Scissors
curved, blunt tips

SUOP23



Corneal Scissor
Corneal Scissors
curved, blunt tips

SUOP24



Vannas Scissor
Str 5mm Blade
Cvd 5mm Blade
Str 7mm Blade
Cvd 7mm Blade
Str 11mm Blade
Ang 11mm Blade

SUOP25



Iris Scissor
Str 9cm
Cvd 9cm
Str 11cm
Cvd 11cm

SUOP26



Moorfield Forceps

Standard
Micro

SUOP27



Moorfield Dressing Forceps

Standard
Micro

SUOP28



Iris Forceps

Str 0.7mm
Cvd 0.7mm
Str Teeth 1x2
Cvd Teeth 1x2

SUOP29



St Martin Forceps

with tying 1x2
Plain tying

SUOP30



Jewelers Forceps

delicate, straight # 1
standard, straight # 3
standard, straight # 4
standard, straight # 5

SUOP31



Adson Forceps

Serrated 1mm
Micro Srt 0.5mm
1x2 Tissue
1x2 Micro 0.5mm

SUOP32



Bishop Harmon Forceps
Serrated
1x2

SUOP33




TUEBINGEN Tying Forceps
Str
Cvd

SUOP34



Jaffe Tying Forceps
Str 6mm
Cvd 6mm

SUOP35



Mc-Pherson Forceps
Str 5mm tying
Ang 5mm tying

SUOP36



Kelman Mcpherson Forceps
Ang 8mm tying
Ang 10mm tying

SUOP37



Bonn Forceps
0.12 1x2
5mm tying

SUOP38



Bonn Forceps
0.12 1x2



Castroviejo Suturing Forceps
0.3mm 1x2
0.8mm 1x2

SUOP39



Colibri Forceps
Tying 1x2, 0.12mm
Tying

SUOP40



Colibri Forceps Notched
Fine Groove



Corneal Fine Notched Forceps
0.2mm notch

SUOP41



Jayle Corneal Tying Forceps
1x2 0.3mm
1x2 0.5mm

SUOP42



Utrata Capsulorhexis Forceps

Capsulorhexis Forceps
very delicate
11mm long shanks



SUOP43

Curved Capsulorhexis Forceps

Capsulorhexis Forceps
very delicate tips
11mm long



SUOP44

Capsulorhexis Forceps Round Bodied

straight shafts.
45° angled shafts.
Tip to angled length 14 mm



SUOP45

Lens Insertion Forceps



SUOP46

Lasik Flap Forcep



SUOP47

Buratto
LASIK Flap Forceps
disc shaped
serrated jaws
Str
Cvd



SUOP48

Barraquer Cilia Forceps



SUOP49

BEER Forceps Cilia



SUOP50

Castroveijo Caliper
Marking 20mm Str



SUOP51

Braunstein Gauge
3,5/4.0mm



SUOP52



Fixation Rings & Lasik

Fine - Thornton
Swivel Fixation Ring blunt teeth
13mm diameter
16mm diameter

SUOP53



Thornton
Glove Fixation Ring
multiple blunt teeth
16mm diameter

SUOP54



Lasek Alcohol Well
Liquid holding well
Well 7 mm diameter
Well 8 mm diameter
Well 8.5 mm diameter
Well 9 mm diameter

SUOP55



Lasik Flap Lifter
Combines a 45° angled
sinskey hook and a thin
hook with a tip

SUOP56



Sub Tenon Set

SUOP57



Barraquer Needle Holder

With lock
Without lock

SUOP58



Castroviejo Needle Holder

Str with lock
Str without lock
Cvd with lock
Cvd without lock

SUOP59



Jaffe Needle Holder

Without lock

SUOP60



Barraquer Needle Holder

Without lock cvd

SUOP61



Simcoe I/A Cannula
Simcoe Cannula



SUOP62

Air Injection
Cannula angled at 5mm
30Gauge
27Gauge



SUOP63

Hydrodissection
Cannula angled at 8mm
flat tip, front opening
25Gauge



SUOP64

SUOP65

Lacrimonal Dilator And Probe
Double end



Double Curette
1.5/2mm
2.5/3mm

SUOP66

Double Ended
Punctal Dilator
Double - Ended

SUOP67



Punctal Dilator
Sharp end
Blunt end

SUOP68



Schocket
Scleral Depressor

SUOP69



Foriegn body spuds

Golf Club, angled,
lightly rounded tip

Francis, angled,
flag shaped tip

Ellis, lightly curved and
rounded tip



SUOP70



Epi - Peeler **SUOP71**



Flap Elevator

Blunt tipped,
1 mm spatula

SUOP72



Mendez

Degree Gauge 12mm I.D.

SUOP73



Stevens Scissors

Str 11cm
Cvd 11cm

SUOP74



Strabismus Scissor

Str b/b 11cm
Cvd b/b 11cm

SUOP75



Dressing Scissor

Str 14cm S/S
Str 14cm B/B
Cvd 14cm B/B
Cvd 14cm S/S

SUOP76



Stitch Scissor
9cm

SUOP77



Mosquito Forceps
Str 9cm
Cvd 9cm
Str 12.5cm
Cvd 12.5cm

SUOP78



Kelly Forceps
14cm
18cm

SUOP79



Sponge Forceps
18cm
25cm

SUOP80



Mayo Needle Holder
13cm
16cm
18cm

SUOP81



Adson-Forceps	66, 114	Braunstein-Marker	40
Air Injection-Cannula	104	Brierley-Splitter	20
Akahoshi Combo Prechopper	73	Bunge-Spoon	30
Akahoshi Prechopper	73	Buratto-Cannula	106
Akahoshi-Manipulator	21	Buratto-Forceps	87, 118
Alfonso	30	Buratto-Lasik Flap Protector	15
Alfonso-Forceps	85	Buratto-Mannis Lasik Flap Protector	15
Alfonso-McIntyre Cannula	107	Camellin Style-Lasek	16
Alio-Kelman-Hook	32	Capsule-Curette	14
Alio-Scissors	42	Capsulorhexis Forceps	70, 71, 72, 117
Alphonso-Speculum	4	Capsulorhexis Forceps For Micro Coaxial Phaco ...	72
Anis	16	Capsulotomy Scissors Vannas	49
Anis-Forceps	74	Carlson DSEK Smoother	17
Anis-Manipulator	25	Castoviejo-Caliper	40, 118
Anis-Marker	38	Castoviejo-Spatula	28
Anis-Needle Holder	97	Castoviejo-Colibri Forceps	64
Arita-Forceps	78	Castoviejo Suturing Forceps	63, 116
Arruga-Forceps	84	Castoviejo-Blade Breaker & Holder	13
Aspiration Handpieces	108	Castoviejo-Colibri Forceps	58
Baby Jones-Forceps	92	Castoviejo-Dilator	102
Backhaus-Forceps	92	Castoviejo-Needle Holder	94, 95, 96, 120
Banaji-Cannula	106	Castoviejo-Scissors	43, 44
Bangerter-Speculum	7, 110	Castoviejo-Speculum	7
Bard-Parker Handle	13, 111	Catalano-Intubation Set	102
Barraquer Cilia Forcep	78, 118	Chalazion Forceps	90, 112
Barraquer Golf Pattern-Needle Holder	98	Chang-Cannula	107
Barraquer-Forceps	79	Chang-Chopper	23
Barraquer-Kratz Speculum	3	Chang-Seibel Chopper	23
Barraquer-Needle Holder	96, 97, 98, 120	Choppers & Manipulators	111
Barraquer-Scissors	49	Christensen-Punch	12
Barraquer-Solid Blade Speculum	3	Chuck Handle	13
Barraquer-Spatula	14	Cilia Forceps	79
Barraquer-Speculum	3, 110	Cilia Forceps-with spring tension	78
Barron Handle # 3	13	Citelli-Rongeur	99
Bayonet-Forceps	66	Clayman-Forceps	88
Bechert	25	Coaxial Forceps	85
Bechert-Forceps	87	Colibri Corneal Forceps	59
Beer Cilia Forceps	78, 118	Colibri Forceps	57, 116
Benedetto-Spatula	27	Colibri Style-Forceps	57
Berke-Forceps	88	Conjunctival Scleral Fixation Forceps	61
Beyer-Rongeur	99	Conjunctival-Scissors	45
Binder	18	Cook Eye-Speculum	4
Bishop-Harmon Cannula	106	Corbin-Cannula	104
Bishop-Harmon Forceps	65, 115	Corneal Dissector	17
Black coated Supercut-Scissors	52	Corneal Forceps	60, 61, 116
Blaydes-Forceps	88	Corneal Scissors	45, 48, 113
Bonaccolto-Forceps	79	Corneal Scleral Punch	85
Bonn Corneal Forceps	60	Crawford-Hook	103
Bonn Micro Tying Forceps	62	Crawford-Intubation Set	102
Bonn-Forceps	59, 116	Culler-Spatula	14
Bonn-Hook	26, 31	Culler-Spoon	30
Bores-Axis Marker	37	Dahan-Cannula	106
Bores-Forceps	64	Dalk-Cannula	104
Botvin-Forceps	57	Dalk-Corneal Dissector	15
Bowman-Probe	102	Deltz-Hook	32
Braunstein-Gauge	118	Desmarres Forceps	90

Alphabetical Index

Desmarres-Retractor	9	Giannetti-Forceps	79
Diffenbach Serrefine-Forceps	92	Giannetti-Spoon	30
Dodick-Forceps	86	Gill-Knife	11
Dodick-Kammann Chopper	20	Gills-Colibri Forceps	57
Dodik-Forceps	85	Gills-Vannas Scissors	49
Double Ended	24, 111, 121	Gimbel-Cannula	106
Double Ended Marker	39	Gimbel-Spatula	21
Double Ended-Curette	14, 121	Gimbel-Stabilization Ring	34
Dressing Forceps	65	Girard-Forceps	64
Dressing Scissors	122	Girard-Spatula	27
Drews-Needle Holder	95	Globe Ring-Fixation Ring	33
Drysdale-Manipulator	25	Golf Ball Tying Forcep	76
Du Plessis-Hook	36	Gradle-Forceps	79
Elsching-Spatula	27	Graefe-Hook	35
Endothelial Stripper	16	Graether	31
Enucleation-Scissors	52	Grayson-Manipulator	25
Epi-Peeler	122	Green Fixation Forcep	77
Ernest McDonald II-Forceps	89	Green-Hook	36
Ernest McDonald-Forceps	87, 89	Guimares-Lasik	18
Ernest-Forceps	84	Guyton-Hook	36
Escaf-Divider	20	Habele-McPherson Forceps	67
Escaf-Manipulator	21	Halsted-Forceps	91
Eye Scissors	49, 50	Harms-Forceps	75
Eye Scissors	49, 50	Harms-Probe	101
Faulkner-Forceps	87	Harris-Sinskey Hook	26
Fechtner-Forceps	58	Hartmann-Forceps	91
Fechtner-Marker	38	Helveston-Forceps	78
Fine-Forceps	78	Helveston-Hook	35, 36
Fine-Scissors	47	Helveston-Retractor	10
Fine-Thornton Fixation Ring	34	Helveston-Ruler	40
Fison Retractor	103	Henderson-Gauge	37
Fixation Forcep	77	Hersh	18
Fixation Forceps	65, 66	Hirschman-Hook	26
Fixation Hook	33	Hirschman-Spatula	19
Flap Elevator	122	Hoffer-McPherson Forceps	63
Fleirina-Fixation Ring	8	Hofmann-Thornton Fixation Ring	34
Fogla Pointed Dissector	17	Hoskin-Forceps	57, 58
Fogla Scissors Right	48	Howarth Periosteal Elevator	103
Fogla Trifacet Spatula	17	Hu Femto Forceps	85
Foldable Lens Cutter Squeeze	53	Hunkeler-Manipulator	22
Foldable Lens Removing Forceps	53	Hwang-Chopper	23
Foreign body-spud	11, 121	Hydrodissection-Cannula	106
Foreign-Forceps	81	ICL Forceps	53
Frankel-Speculum	4	ICL Manipulator	31
Frazier Suction Tube	93	Ikeda Capsulorhexis Forceps	56
Freeman Jaffe-Forceps	74	Implantation Forceps	84
Freer-Chisel	100	Inamura Capsulorhexis-Forceps	69, 71
Freer-Elevator	100	Inamura-Forceps	68, 69
Fukasaku-Cannula	106	Infant Lid-Speculum	4, 110
Fukasaku-Knife	11	Infant-Dilator	101
Fukasaku-Lasik	18	Infusion Handle	108
Fukasaku-Punch	12	IOL Lens Cutter	53
Gass-Hook	35	IOL Loading Forcep	76
Gass-Punch	12	IOL Removal Scissors	45
Gerl-Forceps	74	Iris Forceps	61, 114
Gess-Cannula	104	Iris-Scissors	47, 113

Irrigating Cannula.....	106	Lester-Burch Speculum.....	6
Irrigation Handpieces.....	108	Lester-Forceps.....	66
Jaeger-Lid Plate.....	10	Lester-Manipulator.....	19, 30
Jaffe-Forceps.....	74, 83, 115	Lewicky-Hook.....	26
Jaffe Lid Retractor.....	110	Lieberman-Spatula.....	21
Jaffe-Needle Holder.....	94, 120	Lieberman-Speculum.....	6, 110
Jaffe-Scissors.....	47	Lieppman-Cannula.....	105
Jaffe-Spatula.....	19	Lippman-Spatula.....	27
Jaffe-Speculum.....	7	Livernois-Forceps.....	86, 89
Jameson-Caliper.....	40	Livernois-McDonald Forceps.....	89
Jameson-Forceps.....	88	Lu-Cannula.....	106
Jameson-Hook.....	36, 112	Lu-Mendez-Fixation Ring.....	33
Jayle Corneal Tying Forceps.....	62, 116	Maddox-Lasik.....	18
Jensen-Cannula.....	105	Mallet.....	101
Jewelers-Forceps.....	82, 83, 114	Maloney.....	39
Josephberg-Besser Retractor.....	9	Masket-Forceps.....	68, 69, 70
Joseph-Elevator.....	100	Masket-Splitter.....	23
Kalt-Needle Holder.....	95	Maumenee-Forceps.....	75
Katzen-Lasik.....	18	Mayo-Needle Holder.....	123
Katzin-Scissors.....	44	McDonald-Marker.....	39
Keener-Divider.....	29	McIntyre-Binkhorst Cannula.....	107
Kelly-Forceps.....	123	McIntyre-Cannula.....	105, 106
Kelly-Punch.....	12	McNeill Goldman-Fixation Ring & Blepharostat.....	8
Kelman McPherson-Forceps.....	77, 115	McPherson Tying Forceps.....	75
Kelman-Cannula.....	105	McPherson-Forceps.....	63, 77, 115
Kelmann Forcep.....	76	McPherson-Needle Holder.....	94
Kelmann McPherson-Forceps.....	60	McPherson-Vannas Scissors.....	49
Kerrison-Rongeur.....	99	McPherson-Westcott Scissors.....	46
Keuch-Dilator.....	31	Mellinger-Speculum.....	5
Kimkura-Spatula.....	28	Mendez-Gauge.....	37, 122
Kirby-Loop.....	29	Metzenbaum-Scissors.....	52
Kirby-Lens Loop & Muscle Hook.....	111	Meyerhoefer-Curette.....	14
Knapp-Retractor.....	10	Micro Capsulorhexis Forceps.....	56
Knapp-Scissors.....	51	Micro Forcep.....	81
Knapp-Speculum.....	5	Micro Needle Holder.....	97, 98
Knolle-Pearce Cannula.....	107	MICS Capsulorhexis Forceps.....	56, 72
Knolle-Spatula.....	21	Mics Phaco Pre Chopper.....	20
Koch-Spatula.....	21	MICS/MCCS Capsulorhexis Forceps.....	72
Kraff-Forceps.....	68	MICS/MCCS Coaxial Capsulorhexis Forceps.....	56
Kratz-Cannula.....	105	Milder-Chopper.....	20, 23
Kratz-Speculum.....	110	Mini Westcott-Scissors.....	46
Kuglen Iris & Hook Lens Manipulator.....	32	Moody-Forceps.....	64
Kurstin-Forceps.....	67	Moorfields-Forceps.....	80, 114
Lacrimal-Cannula.....	104	Moorfields-Hook.....	36
Lacrimal Dilator and Probe.....	121	Moran-Lasik.....	17
Lancaster-Speculum.....	5, 110	Morrison-Loop.....	29
Lasek Alcohol Well.....	119	Mosquito Forceps.....	123
Lasek-Cannula.....	104	Muldoon-Dilator.....	101
Lasik Flap Lifter/Forceps.....	19, 117, 119	Multipurpose Scissors Noyes.....	48
Lieberman Style-Aspirating Speculum.....	8	Murdoch-Speculum.....	5
Lieberman-Speculum.....	6, 7	Mushroom and Chopper.....	23
Lens-Forceps.....	117	Mushroom and Y Pusher.....	23
Lens Injector.....	41	Mushroom-Manipulator.....	22
Lens Loop.....	29	Nasal-Speculum.....	100
Lens Manipulating Hook.....	26, 32	Neonatal Speculum.....	110
Lens Manipulator.....	32	Nettleship Lacrimal Dilator.....	103

Alphabetical Index

Neucleus Sustainer	22	St Martin Tying Forceps	62, 114
Neuhann-Divider	22	Stallard Tying Forceps	80
Neyvas-Forceps	69	Stern-Castroviejo Forceps	64
Nucleus Dividing and Spreading	28	Stevens-Cannula	104
Nucleus Manipulators & Hooks	111	Stevens-Hook	35
Nugent-Forceps	81	Stevenson-Speculum	7
Obstbaum-Spatula	27	Stevens-Scissors	51, 122
Optic Zone Markers	38	Stitch Scissors	50, 123
Osher Neumann-Marker	39	Strabismus-Scissors	51, 122
Osher-Forceps	83	Straight Ring-Manipulator	28
Paton-Spoon	30	Strauss-Simcoe Cannula	105
Pean Forceps	93	Stripping Forceps	86
Phaco Chopper	20	Sub Tenon Set	119
Pierse-Forceps	58	Suture-Forceps	61, 83
Pigtail-Probe	103	Swiss Model-Blade Breaker & Holder	13
Pisacano-Rotator	25	Temporal Approach-Speculum	3
Prechopper	73	Tennant-Forceps	74
PRK Scraper--Marker	39	Tennant-Rotator	22
Punctual Dilator	121	Terry-Scissors	44
Ramplays Forceps	93	Thoralakson-Flap Lifter	19
Reversed Simcoe-Cannula	105	Thornton-Fixation Ring	34
Reversed Sinskey-Hook	26	Thornton-Forceps	119
Rich-Scissors	44	Thorpe-Forceps	59
Rigid IOL Implantation Forceps	84	Tibolt-Scissors	46
Rollet Rugine	103	Tissue Forceps	65
Rollet-Retractor	10	Tooke-Knife	11
Rowen-Lasik Flap Protector & Spatula	15	Troutman-Castroviejo Scissors	43
Rowen-Spatula	23	Troutman-Forceps	66, 76
Ruler	40	Troutman-Needle Holder	96
Rycroft Corneal Tying Forceps	80	Tubingen Model-Forceps	75
Saur-Speculum	4	Tuebingen Tying Forceps	115
Schaedel-Forceps	92	Twist Hook	33
Schepens-Retractor	9	Tying Forceps	67, 77
Schocket-Retractor	9	Ultrata Capsulorhexis Forcep	68, 69, 117
Schocket/Scleral Depressor	121	Ultrata-Forceps	68, 69
Scleral Plug Forceps	80	Utility and Splinter-Forceps	81
Scoring Hook & Irrigation Handpiece	108	Utility-Scissors	52
Scoring reverse hook	22	Utoff-Scissors	48
Seinert-Forceps	89	Vannas-Scissors	113
Seipser	25	Vessel Dialator	81
Senn-Retractor	9	Vidaurri-Cannula	107
Serrefine-Forceps	92	Vitrectomy-Scissors	42
Shammas-Cannula	105	Vitreo-Retinal Forceps	54, 55
Shapiro-Cannula	105	Warren	23
Sharp Chopper and Blunt Chopper	23	Watzke-Forceps	90
Shepard-Hook	26	Weiss-Speculum	8
Sheppard-Forceps	88	Wells-Spoon	31
Silcock, Catch-Needle Holder	94	West-Chisel	100
Silverman-Pick	33	Westcott-Scissors	46, 47, 112, 113
Simcoe-Cannula	104, 105, 107	Wheeler-Spatula	14
Sinskey II-Hook	26	Wilder-Dilator	101
Sinskey-Hook	26	Williams-Speculum	5
Slone-Lasek Flap Repositor	16	Wire Speculum-Infant Speculum	4
Spaleck-Forceps	87	Wright-Fascia Needle	10
Sponge Forceps	123	YAE-Scissors	42
Spud and Gouge	11		

Contacts

North Wales, North Midlands & North West England

Ian Simpson

07872 512601

i.simpson@rbmedical.co.uk

North Thames, East Midlands & East England

Jason Miller

07500 537153

j.miller@rbmedical.co.uk

South Thames & South England

Adam Knight

07471 865733

a.knight@rbmedical.co.uk

South Wales, South Midlands & South West England

Ian Simpson

07872 512601

i.simpson@rbmedical.co.uk

stephen@plsmedical.co.uk

Stephen Sewell

07484 332995

stephen@plsmedical.co.uk

rbmedical

Providing the instruments for good health



rbmedical

RB Medical Engineering Ltd
Unit 2, Alton Road Industrial Estate
Ross-on-Wye, Herefordshire
HR9 5NSt

T: +44 (0)1989 563958

F: +44 (0)1989 768267

eMail: info@rbmedical.co.uk

