

Article



In Vitro Screening of Chicken-Derived *Lactobacillus* Strains that Effectively Inhibit *Salmonella* Colonization and Adhesion

Dan Hai¹, Zhaoxin Lu¹, Xianqing Huang², Fengxia Lv¹ and Xiaomei Bie^{1,*}

- ¹ College of Food Science and Technology, Nanjing Agricultural University, Nanjing 210095, China; 2017208019@njau.edu.cn (D.H.); fmb@njau.edu.cn (Z.L.); lvfengxia@njau.edu.cn (F.L.)
- ² College of Food Science and Technology, Henan Agricultural University, Zhengzhou 450000, China; hxq8210@126.com
- * Correspondence: bxm43@njau.edu.cn; Tel.: +86-25-8439-6570

Abstract: Inhibition of Salmonella by Lactobacillus has been a popular research topic for decades; however, the inhibition potential of chicken-derived Salmonella by chicken-derived Lactobacillus has not yet been studied. In this study, 89 strains of Lactobacillus from chicken intestines were isolated by national standard method, Gram staining, physiological, and biochemical experiments and molecular sequencing; The inhibition characteristics of 89 strains of chicken derived Lactobacillus against 10 strains Salmonella (S. Enteritidis SE05, SC31, SC21, SC72 SC74, SC79, SC83, SC87; S. bongori SE47; S. Typhimurium, SC85) were detected by agar inhibition zone, The results showed that the inhibition zone of 24 strains of chicken derived Lactobacillus was more than 10 mm, which indicated that the isolated chicken derived Lactobacillus could effectively inhibit the growth of Salmonella; The drug resistance and bile salt tolerance of these 24 strains were analyzed, The results showed that the standard strains LG and L76 were not resistant, and the other 22 Lactobacillus strains showed different degrees of resistance. The strains LAB24, LAB26, LAB53, LAB69, and L76 showed good tolerance at the concentration of 3 g/L bile salt; Caco-2 cell experiment and flow cytometry were used to analyze the inhibitory effect of chicken derived Lactobacillus on the adhesion of Salmonella to Caco-2 cells, The results showed that 16 probiotics could effectively inhibit the adhesion of Salmonella to Caco-2 cells. Twelve probiotics were identified by molecular biology. The results showed that L76 was Enterococcus faecalis, and the other 11 strains were Lactobacillus.

Keywords: chicken-derived Lactobacillus; chicken-derived Salmonella; bacteriostatic; Caco-2 cells

1. Introduction

Salmonellosis is a zoonotic food borne disease which causes outbreaks and sporadic cases of gastroenteritis in human worldwide [1]. Chickens have been known as the major source of *Salmonella* contaminated food products such as chicken eggs and meat that cause human salmonellosis in many countries [2,3]. Although poultry are asymptomatic *Salmonella* carriers, and their production performance is unaffected, *Salmonella* can continuously colonize the cecum of broilers [4,5]. During slaughtering and processing, the infected broiler contaminates the meat products, thereby causing food poisoning to humans through the food chain. Infected laying hens can contaminate their eggs, thereby vertically transmitting the infection to the offspring. Salmonellosis may also be caused by ingesting egg products contaminated by *Salmonella*. Notably, in recent years, multiple resistant *Salmonella* strains have been detected. Moreover, the continued use of antibiotics as growth promoters and to control *Salmonella* colonization of broilers may result in more resistant strains. In many Chinese regions, meat products and poultry have been found positive for a variety of foodborne *Salmonella* strains [6,7].

Lactic acid bacteria (LAB) are widely used probiotic organisms, and their strains usually occur in human and animal intestinal microbiota. LAB promote the development



Citation: Hai, D.; Lu, Z.; Huang, X.; Lv, F.; Bie, X. In Vitro Screening of Chicken-Derived *Lactobacillus* Strains that Effectively Inhibit *Salmonella* Colonization and Adhesion. *Foods* **2021**, *10*, 569. https://doi.org/ 10.3390/foods10030569

Academic Editor: Dario De Medici

Received: 19 February 2021 Accepted: 4 March 2021 Published: 9 March 2021

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2021 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). of host-favorable microbiota reduce or prevent the colonization of harmful pathogens, enhance mucosal immunity, improve the digestibility of the gastrointestinal tract and reduce its pH, and enhance the maturity and integrity of the intestinal tissue [4]. Additionally, some LAB strains are tolerant to the effects of digestive products, such as gastric acid and bile salts [8], and they can adhere to the host intestinal epithelium [9], thereby gaining a competitive advantage and are important for bacterial maintenance [10]. In vitro models with intestinal cell lines are widely used. Laboratory models using human intestinal cell lines, such as Caco-2 [11], have been developed to study the adhesion of probiotic LAB and the competitive exclusion of pathogenic bacteria.

The aim of this study is to study the efficacy of chicken-derived probiotics on *Salmonella* (*S.* Enteritidis SE05, SC31, SC21, SC72 SC74, SC79, SC83, SC87; *S.* bongori SE47; *S.* Typhimurium, SC85) by preliminary screening. The isolated strains of chicken-derived probiotics with better antibacterial effect were selected, their surface characteristics, bile salt tolerance and drug resistance of lactic acid bacteria from chicken were studied; The possible reasons of lactic acid bacteria inhibiting *Salmonella* were further analyzed: Study on the adhesion of lactic acid bacteria and the effect of protecting cell membrane.

2. Materials and Methods

2.1. Strain Isolation, Culture Media, and Used Cell Lines

120 samples of chicken intestines from different slaughterhouses were collected, and *Lactobacillus* were isolated by national standard method [12]. Added 25 g mixed sample to 225 mL normal saline. Mixed at 37 °C for 30 min and serially diluted, coated on de Man, Rogosa, and Sharpe (MRS) agar solid medium, then the coated plates were incubated at 37 °C for 48 h. Typical samples with obvious sour taste were transferred to MRS liquid medium, incubated at 37 °C for 48 h, and then coated again. Then, 107 suspected *Lactobacillus* were selected for Gram staining, and 89 of them turned purple after adding Gram staining agent, microscopic observation showed spherical, short chain, rod-shaped, etc. [13], and then combined with the physiological and biochemical identification of *Lactobacillus*, finally the 89 strains of *Lactobacillus* were isolated. The result show in the Tables S1 and S2 and Figure S1. *S. Enteritidis* SE05, SC31, SC21, SC72 SC74, SC79, SC83, SC87; *S. bongori* SE47; *S. Typhimurium*, SC85 were isolated previously [14] and preserved in our laboratory, Caco-2 cells (Shanghai Bogu Biotechnology Co., Ltd., Shanghai, China) were purchased by our laboratory.

2.2. Inhibitory Effect of Lactobacillus Strains on Salmonella

The 89 *Lactobacillus* strains were cultured in MRS liquid medium at 37 °C for 48 h, centrifuged at 5000 rpm/min for 5 min, and the supernatant was taken for use. Next, we added 100 μ L of *Salmonella* suspension (SE05, SE47, SC21, SC31, SC72, SC74, SC79, SC83, SC85, and SC87) to 100 mL of LB liquid medium, and cultured the mixture at 37 °C and 180 rpm/min for 24 h. Then, 100 mL bottles of solid LB culture medium melted completely and cooled to 40–50 °C, we added 100 μ L of culture liquid that contained 10 *Salmonella* isolates, which were cultured overnight; then, the mixture was shaken and poured on the agar plates, which were marked in advance. After the plates cooled, we used a 5 mm punch to drill holes on them. The number of holes drilled depended on the plate size; the discarded agar pieces were picked out with sterilized toothpicks. We added the samples according to the marks on the backs of the plates, as well as 100 μ L of LAB supernatant to each hole. Following this, we incubated the plates at 37 °C for 24 h [15]. Twenty-four *Lactobacillus* strains that had the best bacteriostatic effect were selected for subsequent experiments. Three parallel experiments were performed in each group

2.3. Antibiotics Resistance of Lactobacillus Strains

Antimicrobial susceptibility of the 24 *Lactobacillus* isolates was evaluated by the disk diffusion test, using Mueller–Hinton agar, according to the guidelines of the Clinical and Laboratory Standards Institute (CLSI, 2012) [16]. The isolates were screened for resistance

to the following antibiotics: tetracycline (TET, 30 mg); amoxicillin (AMX, 30 mg); ceftriaxone (CRO, 30 μ g); chloramphenicol (CHL, 30 μ g); gentamicin (GEN, 10 μ g); trimethoprim/sulfamethoxazole (SXT, 1.25/23.75 μ g); kanamycin (KAN, 30 μ g); erythromycin (ERY, 15 μ g). Three parallel experiments were performed in each group. The results were interpreted according to the established CLSI guidelines (CLSI, 2012).

2.4. Lactobacillus Bile Salt Tolerance Test

The concentration of bile salt in most animals ranged from 0.3 g/L to 3.0 g/L. in this experiment, the concentration of bile salt was selected in the range of 1.0 g/L–3.0 g/L [17]. The activated bacterial suspension was centrifuged at 4000 rpm for 25 min, and the supernatant was discarded. The precipitated bacteria were washed three times with sterile MRS liquid medium. The bacteria were suspended in the medium, and the bacterial suspension concentration was adjusted to 3.0×10^8 CFU/mL. The reconstituted bacterial suspension was divided into different groups. Next, by adding 1.0, 2.0, and 3.0 g/L of pig and bovine bile salt, the bacterial suspensions were incubated at a constant temperature of 37 °C for 3 h and then removed the supernatant. Gradient dilution was performed immediately by PBS, 1 mL mixture was used to coat the MRS plate. After 48 h of culture at 37 °C, we determined the colony count. The survival rate was calculated as follows: survival rate (%) = number of "2 H" colonies/number of "0 h" colonies × 100% [18]. Three parallel experiments were performed in each group.

2.5. Adhesion of Different Bacteria to Caco-2 Cells

Firstly, 1 mL of each *Lactobacillus* or *Salmonella* suspension (2×10^8 CFU/mL) was added to the corresponding 6-well plate containing Caco-2 cells, and cultured at 37 °C and 5% CO₂, Standard *Lactobacillus* strain LG (*L. rhamnosus* GG, FSMM22) as control. Each treatment was repeated three times. After 1 h culture, the culture medium was taken out, the supernatant was discarded, and Caco-2 cells were washed with PBS (pH 7.4) for three times to remove the non-adherent bacteria. The sterile 1% Triton X-100 PBS (1 mL, pH 7.4) solution was added into each well, and it was allowed to stand for 10 min. After the Caco-2 cells were completely removed, the supernatant were transferred into the sterile 1.5 mL centrifuge tube and mixed well. A total of 100 µL of the prepared mixture was diluted at a gradient of 10^{-1} – 10^{-6} . Then *Lactobacillus* and *Salmonella* were counted on the medium plate [19,20]. The adhesion number of bacteria was calculated using Equation (1):

A dhesion index = number of bacteria/cells adhered(1)

2.6. Inhibition of Salmonella Adhesion to Caco-2 Cells by Lactobacillus

Caco-2 cells were passaged for 3-4 generations in DMEM complete medium, and then lactic acid bacteria adhesion test was started. The density of Caco-2 cell suspension was adjusted to 2×10^5 cells/mL. A total of 16 strains of *Lactobacillus* were isolated and cultured for 24 h. After centrifugation, the supernatant was removed and the OD value of the suspension was adjusted to 0.6 with PBS. Three Salmonella strains were labeled with 1 mL FITC fluorescent solution and incubated in dark for 30 min. The supernatant was removed by centrifugation and resuspended with PBS. A 48 well plate, 200 µL/well was used. First, 48 wells were wetted with 100 μ L DMEM solution, then 100 μ L cell suspension was added, and 300 μ L DMEM incomplete culture medium was added to 500 μ L. In the competitive adhesion experiment group, 500 μ L LAB + 500 μ L fluorescent labeled Salmonella were added to each well at the same time; In the rejection adhesion experiment group, 500 μ L LAB was added for incubation at 37 °C for 1 h, PBS was washed for 3 times, and then 500 µL fluorescent labeled Salmonella for 1 h was added for incubation at 37 °C for 1 h; In the replacement adhesion experiment group, 500 µL fluorescent labeled Salmonella was added for 1 h at 37 $^\circ C$ for 1 h, washed with PBS for 3 times, and then 500 μL LAB was added for 1 h at 37 °C for labeling. Finally, all treatments were washed with PBS three times, digested with 0.2 mL trypsin for 2 min, and then added with 0.4 mL DMEM complete culture medium to terminate the reaction. The liquid was collected, and the

fluorescence values of each group were detected by the multi-function fluorescent enzyme reader (Infinite M200 Pro multi-function enzyme reader, TECAN, Männedorf, Switzerland). Three parallel experiments were performed in each group.

2.7. Effect of Lactobacillus Adhesion to Caco-2 Cells on Their Physiological Metabolism

A total of 1 mL Caco-2 cell suspension (2×10^5 cells/mL) was transferred to a 6-well cell culture plate, and 1 mL DMEM complete culture medium was added and incubated at 37 °C and 5% CO₂. After the second day of culture, the new DMEM solution was replaced. After the Caco-2 cells completely adhered to the 6-well plate, the used DMEM medium was discarded, and the Caco-2 cells were washed twice with PBS. In each well, 1 mL DMEM medium containing 10% fetal bovine serum was added to a 6-well cell culture plate. The components are as follows: The suspension of CT group (1 mL sterile PBS), LAB group (1 mL 2 \times 10⁸ CFU/mL *Lactobacillus*), SAL group (1 mL 2 \times 10⁸ CFU/mL Salmonella), and rejection group (1 mL of Lactobacillus was cultured for 1 h, then 1 mL of Salmonella was cultured for 1 h), each group was cultured for 2 h. After culture, the supernatant was collected, and the Caco-2 cells were washed three times with sterile PBS. The Caco-2 cells were dissolved with 1% Triton X-100 (1 mL per well), and the collected lysate stored at -20 °C. Then the alkaline phosphatase (AKP, BC2145, Solarbio, BeiJing, China) [21] and lactate dehydrogenase (LDH, BC0685, Solarbio, BeiJing, China) activities were tested according to the instructions of the sword test box. Three parallel experiments were performed in each group.

2.8. Adhesion of Lactobacillus and Salmonella Isolates to Caco-2 Cells and Apoptosis Test

Firstly, 5 mL Cultured Caco-2 cells ((1 – 5) \times 10⁶/mL per sample) were directly collected into 10 mL centrifuge tubes and centrifuged at 500–1000 r/min for 5 min. Then, the culture medium was discarded, and the cells were washed with the incubation buffer. Next, they were centrifuged again at 1000 r/min for 5 min. Then, they were labeled, re-suspended in a 100 μ L labeling solution, incubated in the dark for 10–15 min at room temperature, centrifuged at 1000 r/min for 5 min, and washed with the incubation buffer. SA-FLOUS (Fluorescent solution) was added, and the mixture was incubated at 4 °C for 20 min, in the dark, with no shaking. Flow cytometry analysis (CytoFLEX LX, BECKMAN, Bria, Florida, USA) showed that the fluorochromes were excited with a 488 nm wavelength. Fluorescein isothiocyanate (FITC) fluorescence was detected with a passband filter at 515 nm wavelength, while PI was detected with a filter at a wavelength greater than 560 nm. Specific parameters of the instrument: Flow pressure: 0–100 psi, step-less adjustable. The liquid flow system can effectively filter 99.999% of the impurities $>0.1 \mu m$ in the sheath. We used the automatic sample loading system with 1.5 mL sample tubes, which offers the advantage of automatic cleaning, backwashing, bubble removal, and temperature control. Separation speed: 100,000 cells/s, separation purity: >99% (using 70 µm nozzle, at the speed of 70,000 s, 60 psi pressure, the sample with 1% target cell content was separated under the separation mode). The recovery rate: >90% of the theoretically predicted cells. Effective acquisition speed: >100,000 cells/s. Effective sorting speed: >70,000 cells/s. Three parallel experiments were performed in each group.

2.9. Purification and Sequencing of PCR Products

Firstly, 12 strains of *Lactobacillus* with typical physiological and biochemical characteristics were screened out by acid-base, bile salt, drug resistance and cell adhesion tests, and then the molecular sequencing test was carried out. 27F(5'-AGAGTTTGATCCTGGCTCAG-3') and 1492R(5'-GGTTACCTTGTTACGACTT-3') primers (Primer production, Henan ShangYa Biotechnology Co., Ltd., Zhengzhou, China) were used for PCR amplification and sequencing identification All the PCR products were purified with a High Pure PCR purification kit (Roche, Mannheim, Germany) and sequenced by Nanjing GenScript Biotech Co., Ltd., China. The resultant DNA sequences were analysed using BLAST (http://www.ncbi.nih.gov/BLAST/. accessed date: 19 May 2020).

2.10. Statistical Analysis

The experiment was performed triplicate as a replicated test. Data were analyzed using one-way analysis of variance using the SPSS 16.0 software (SPSS Inc., Chicago, Illinois, USA) for Windows. The results are expressed as mean (M) \pm standard deviation (SD). Mean separation was performed via Duncan's multiple range tests (p < 0.05).

3. Results

3.1. Inhibition of Salmonella by Lactobacillus

Figure 1 shows the inhibition results obtained for 10 different *Salmonella* strains by using the bacteriostatic circle method. Different *Lactobacillus* strains had different bacteriostatic effects on different *Salmonella* strains. The diameter of inhibition circle is less than 5 mm, which indicates that the bacteriostatic effect of lactic acid bacteria is not obvious. The diameter of inhibition circle is more than 8 mm, which indicates that lactic acid bacteria have bacteriostatic effect. The diameter of inhibition circle is more than 10 mm, which indicates that the bacteriostatic effect of lactic acid bacteria is obvious. From the Figure 1, the average inhibition zone diameter of the strains LAB1, LAB4, LAB7, LAB16, LAB22, LAB24, LAB35, LAB44, LAB53, LAB54, LAB56, LAB60, LAB64, LAB69, LAB70, LAB73, LAB76, and LG were more than 8 mm and were selected. The inhibition zone of LAB69 (17.2 mm) and LAB76 (17.5 mm) is the largest, indicating that the antibacterial effect is the best.

SC87	13	12.5	12	11.5	12	11.5	12.5	11.5	6.5	6.5	7	10.5	SC87	13	9	10	8.5	12.5	13.8	8.9	16	15	12.7	12.7	15	13
SC85	11	13	10.5	12	12	12	13	12	12	12	6.5	12.5	SC85	13	17.5	13	12	14	12	9.5	14.5	13	12	13.5	16	12.8
SC83	12	12.5	10	10.5	13	10.5	12.5	10.5	10.5	10.5	10.5	15	SC83	16.5	14	14.5	11	11.5	9.5	8.5	15	12.8	11	15	15	13
SC79	11	11	11	11	11	11	11.5	11	11	11	5.5	10.5	SC79	10	11	9	10	12	8	9	11.5	11	8.5	10	11.5	10
SC74	10.5	14	12	11	11	11	14	11	11	11	5	12	SC74	13	12.5	10.5	12.5	11.5	13	13.5	14.5	15	13.5	15.5	17.5	12
SC72	12	13	12	10	10	10	13	10	10	10	10	13.5	SC72	11.5	13	13.5	13	12	11.5	13	13	13	11.5	10	16.5	11.5
SC31	11.5	10.5	15	12	11.5	12	10.5	12	12	7.5	5	12	SC31	12.5	14	11	12	13	13	12	17.2	14	11	13	13	12
SC21	12	16	11.5	6	7	8	16	11	11	11	12	11	SC21	14	13.5	10	10	9.5	11	12.5	13	12	13	11.5	15	13.5
SE47	5.5	12.5	10.5	11	11	6	11	12	12	16	10.5	13.5	SE47	11.5	12	12	11	11.5	12.2	11	16.8	12.5	12	11	15.6	11.5
SE05	10.5	10	11	5	6.5	10	10	10	7	12	11	8.5	SE05	10	12	10	9	12	11	10	12	13	6	7	12	10
	LAB1		LAB7		LAB14	4 ່ ເ	AB22	2 ' 1	LAB26	5 1	LAB34	4		LAB4	4 L	AB54	1 L	AB60	ו כ	AB6	5 1	AB7	ו כ	AB73	3	LG
		LAB4	L	AB12		LAB16	5 L	AB24	1	LAB3	LI	LAB35	5	L	AB53	ι Ι	AB56	5 L	AB6	1 L	AB69	L	AB72	: I	LAB7	6
							1																			
					5		1	6										5				17.	5			

Figure 1. Experimental results of inhibitory zone of 24 *Lactobacillus* isolates to 10 *Salmonella* isolates. Note: The vertical axis is 10 strains of *Salmonella*, the horizontal axis is 24 strains of *Lactobacillus*, and the intersection of the horizontal axis and the vertical axis is the diameter of inhibition zone of *Lactobacillus* inhibiting *Salmonella*. (The diameter of inhibition zone: mm).

3.2. Antibiotic Resistance of Lactobacillus Strains

The resistance of 24 strains of *Lactobacillus* to antibiotics TET, AMX, CRO, CHL, GEN, SXT, KAN, and ERY were shown in Table 1. The standard strains LG and LAB76 were not drug-resistant. The other 22 strains of LAB showed different degrees of drug resistance. LAB12, LAB14, LAB22, and LAB34 were resistant to 3 antibiotics, 5 strains of lactic acid bacteria were resistant to 2 kinds of antibiotics, and 14 strains were resistant to 1 kind of antibiotics. This may be related to the type and frequency of antibiotics used in poultry breeding. *Lactobacillus* resistant to three kinds of antibiotics were not used in subsequent experiments. Among the 24 strains of *Lactobacillus*, 25% (6/24) were resistant to CRO, 29.2% (7/24) to SXT, 20.8% (5/24) to CHL, and 25% (6/24) to TET. These indicators provide a data base for screening target lactic acid bacteria.

3.3. Lactobacillus Bile Salt Tolerance Test

Figure 2 shows the growth of 24 strains of LAB within 8 h under three different bile salt concentrations. It can be seen from Figure 2 that the growth of 24 *Lactobacillus strains* is not consistent with the increase of bile salt concentration. The OD₆₀₀ value of LAB4, LAB7, LAB24, LAB26, LAB31, LAB35, LAB44, LAB56, LAB65, LAB69, LAB70, L76, and LG were more than 0.20, show that most of the lactic acid bacteria detected in this paper

can maintain a high survival rate under the condition of 3 g/L bile salt, and the strains still show good tolerance to the bile salt. At present, many studies have proved that most of the lactic acid bacteria can survive well in the intestine (trypsin, pH 8.0, and 0.3% bile salt), and can return to the initial number of bacteria when passing through the small intestine [22,23].

Strain Antib		LAB 4	LAB 7	LAB 12	LAB 14	LAB 16	LAB 22	LAB 24	LAB 26	LAB 31	LAB 34	LAB 35	
GEN	S	S	S	S	R	S	R	R	S	S	S	S	
KAN	S	S	S	R	S	S	S	S	S	S	R	S	
AMX	S	S	S	R	S	R	S	S	S	S	S	S	
CRO	R	R	S	S	R	S	R	S	R	S	S	S	
SXT	S	S	R	S	S	R	S	S	S	S	R	S	
CHL	R	S	S	S	R	S	S	S	S	R	S	S	
TET	S	R	S	R	S	S	R	S	S	S	S	R	
ERY	S	S	S	S	S	S	S	S	R	S	R	S	
Total													
num-	2	2	1	3	3	2	3	1	2	1	3	1	
ber													
	LAB	LAB	LAB	LAB	LAB	LABV	LAB	LAB	LAB	LAB	LAB	LAB	
GEN	44	53	54	56	60	64	65	69	70	72	73	76	LG
KAN	S	S	S	S	S	R	S	S	R	S	S	S	S
AMX	S	S	S	S	S	S	S	S	S	S	S	S	S
CRO	S	S	S	S	S	S	S	S	S	S	S	S	S
SXT	S	S	S	S	S	S	S	S	S	S	R	S	S
CHL	S	R	R	S	R	S	S	S	S	R	S	S	S
TET	S	S	S	R	S	S	R	S	S	S	S	S	S
ERY	R	S	S	S	S	S	S	R	S	S	S	S	S
Total													
num-	S	S	S	S	S	S	S	S	S	S	R	S	S
ber													
	1	1	1	1	1	1	1	1	1	1	2	0	0

 Table 1. Drug resistance of 24 strains of lactic acid bacteria.

Note: S (susceptible), R (Resistant).




Figure 2. The bile salt tolerance of 24 strains of *Lactobacillus*. Note: (**a**)-The growth of 8 strains of *Lactobacillus* (LAB1, LAB4, LAB7, LAB12, LAB14, LAB16, LAB22, LAB24) under the condition of 1 g/L, 2 g/L, 3 g/L bile salt concentration for 8 h; (**b**)-The growth of 9 strains of *Lactobacillus* (LAB26, LAB31, LAB34, LAB35, LAB44, LAB53, LAB54, LAB56, LAB60) under the condition of 1 g/L, 2 g/L, 3 g/L bile salt concentration for 8 h; (**c**)-The growth of 8 strains of *Lactobacillus* (LAB64, LAB65, LAB69, LAB69, LAB70, LAB72, LAB73, LAB76, LG) under the condition of 1 g/L, 2 g/L, 3 g/L bile salt concentration for 8 h.

3.4. Bacterial Adhesion to Caco-2 Cells

Adhesion is an important basis for the colonization of *Lactobacillus* in the intestinal tract, and it is necessary to play the role of biological barrier, inhibit the growth and reproduction of pathogenic bacteria, and improve the structure of intestinal flora [24,25]. We screened the *Lactobacillus* with strong adhesion ability from 15 strains of chicken derived *Lactobacillus*. Figure 3 is the result about the vitro evaluation of the adhesion ability of 15 *Lactobacillus* isolates and 3 *Salmonella* strains using Caco-2 cells as a model. Adhesion index = number of bacteria/cells adhered, the higher the ratio of adhesion index, the stronger the adhesion ability. Fifteen *Lactobacillus* strains had significant differences in their ability to adhere to Caco-2 cells. L76 had the strongest adhesion ability of three *Salmonella* strains was 4.86, 6.03, and 7.57 bacteria/cell). The adhesion ability of three *Salmonella* strains with higher adhesion ability than *Salmonella* as the target strains of subsequent experiments.




Figure 3. Adhesion index of different bacteria to CaCO-2 cells.

3.5. Inhibition of Salmonella Adhesion to Caco-2 Cells by Lactobacillus

The treatment of Caco-2 cells by different bacterial combinations is shown in Table 2. *Salmonella* were incubated with fluorescence for 30 min, and *Lactobacillus* were not treated. The experimental group was incubated with Caco-2 cells according to different combinations. The stronger the fluorescence intensity was, the redder the color was, and the more *Salmonella* adhered to Caco-2 cells.

Single Bacteria Group	Competitive Treatment Group	Rejection Treatment Group	Replacement Treatment Group
SC79 (<i>n</i> = 3)	LAB44 + SC79 $(n = 3)$	LAB44—SC79 $(n = 3)$	SC79—LAB44 (<i>n</i> = 3)
LAB35 (<i>n</i> = 3)	LAB53 + SC79 $(n = 3)$	LAB53—SC79 $(n = 3)$	SC79—LAB53 (<i>n</i> = 3)
LAB44 (<i>n</i> = 3)	LAB35 + SC79 (n = 3)	LAB35—SC79 (<i>n</i> = 3)	SC79—LAB35 (<i>n</i> = 3)
LAB53 $(n = 3)$	LAB69 + SC79 $(n = 3)$	LAB69—SC79 $(n = 3)$	SC79—LAB69 (<i>n</i> = 3)
LAB69 $(n = 3)$	LAB76 + SC79 $(n = 3)$	LAB76—SC79 $(n = 3)$	SC79—LAB76 (<i>n</i> = 3)
LAB76 (<i>n</i> = 3)	LG + SC79 (n = 3)	LG—SC79 (<i>n</i> = 3)	SC79—LG (<i>n</i> = 3)
LG $(n = 3)$			

Table 2. Treatment of Caco-2 cells with different combinations of bacteria.

On the contrary, the fluorescence intensity was weaker, and the color was blue, indicating that the number of *Salmonella* adhered to Caco-2 cells was less. Adding *Lactobacillus* and *Salmonella* labeled with fluorescence to Caco-2 cells can show the change of fluorescence intensity. The changes of cell adhesion of *Lactobacillus* and *Salmonella* in different treatments were studied by color change. In the competition experiment, *Lactobacillus* and *Salmonella* were added to Caco-2 cells at the same time. The results showed that *Salmonella* in SC31-LAB84, SC79-LAB34 and SC79-LAB72 groups adhered more to Caco-2 cells than other groups (Figure 4a). In the rejection experiment, Caco-2 cells were incubated with *Lactobacillus* for 1 h, and then incubated with *Salmonella* (Figure 4b). The results showed that *Salmonella* had more adhesion to Caco-2 cells in LAB54-SC31 group. In other treatments, *Lactobacillus* adhered more to Caco-2 cells. In the substitution experiment (Figure 4c), Caco-2 cells were incubated with *Salmonella* for 1 h, and then incubated with *Lactobacillus*. The results showed that seven groups of treatments (SC31-LAB35, SC31-LAB44, SC31-LAB53, SC31-LAB54, SC79-LAB24, SC79-LAB34, and SC79-LAB44) had stronger fluorescence intensity, and indicated more *Salmonella* adhesion to Caco-2 cells. While the adhesion rate in other groups was inhibited by *Lactobacillus*, based on the color change.




Figure 4. (a) Competitive adhesion test of Caco-2 cells incubated with *Lactobacillus* and *Salmonella*. In Figure 4a, The target *Salmonella* are SC31, SC85, and SE05. In this experiment, *Lactobacillus* and *Salmonella*

were added at the same time and incubated for 1 h, then the fluorescence intensity was measured. The stronger the fluorescence intensity was, the more red the color was. It showed that the number of Salmonella incubated on Caco-2 cells was more, while the fluorescence intensity was weak, the color was blue, indicating that the number of Salmonella incubated on Caco-2 cells was small. (b) Rejection and adhesion test of Caco-2 cells incubated by lactic acid bacteria first and then by Salmonella. In Figure 4b, the Lactobacillus isolated strains were added for incubation for 1 h, then the target Salmonella SC21, SC85 and SE05 were added for incubation for 1 h, and the results were measured by multi-functional enzyme labeling instrument. The stronger the fluorescence intensity was, the more red the color was. It showed that the number of Salmonella incubated on Caco-2 cells was more, while the fluorescence intensity was weak, the color was blue, indicating that the number of Salmonella incubated on Caco-2 cells was small. (c) Displacement adhesion test of Caco-2 cells incubated by Salmonella and then by lactic acid bacteria. In Figure 4c, the target Salmonella SC21, SC85 and SE05 were added to incubate for 1 h, and then the Lactobacillus isolate was added to incubate for 1 h, and then the results were measured by multi-functional enzyme labeling instrument. The stronger the fluorescence intensity was, the more red the color was. It showed that the number of Salmonella incubated on Caco-2 cells was more, while the fluorescence intensity was weak, the color was blue, indicating that the number of Salmonella incubated on Caco-2 cells was small.

3.6. Lactobacilli Adherence Effect on the Physiological Metabolism of Caco-2 Cells

Caco-2 cells incubated with PBS was the control group. LAB group were incubated with Lactobacillus, and the SAL group were incubated with Salmonella. The AKP activity of Caco-2 cell lysate and supernatant in the control group were 1378.21 U/gprot and 79.95 U/gprot, respectively (Table 3). LAB group and SAL group (SC31, SC79 and SE05) had different effects on the metabolism of Caco-2 cells. Compared with the control group, the AKP activity of Caco-2 cell supernatant in SAL group was significantly higher than that of Caco-2 cell lysate. SAL significantly inhibited the activity of AKP in Caco-2 cells, destroyed the integrity of cell membrane, and resulted in a large amount of AKP leakage into the cell culture supernatant. Compared with the control group, the AKP activity in the supernatant of Caco-2 cells treated with LAB decreased by 15–59.6% (p < 0.01), and the AKP activity in Caco-2 cells treated with LAB also decreased significantly, but it was still 3.45–4.92 times higher than that in SAL group (p < 0.01). The AKP activity of Caco-2 cells treated with LAB-SAL was slightly higher than that of SAL, while the AKP content in the supernatant of Caco-2 cells treated with LAB-SAL was significantly lower than that of control group. These results indicate that *Lactobacillus* can repel *Salmonella* and protect Caco-2 cell membrane.

Treatment Group	Cell Culture Supernatant AKP (U/100 mL)	Cell Lysate AKP (U/gprot)	Treatment Group	Cell Culture Supernatant AKP (U/100 mL)	Cell Lysate AKP (U/gprot)
CT(PBS)	1378.21 ± 1.36	79.95 ± 1.36	LAB76-SC31	28.56 ± 1.36	178.25 ± 1.36
LAB24	321.97 ± 1.36	29.66 ± 1.36	LG-SC31	29.37 ± 1.36	167.91 ± 1.36
LAB31	391.54 ± 1.36	18.44 ± 1.36	LAB24-SC79	39.97 ± 1.36	152.09 ± 1.36
LAB35	401.62 ± 1.36	14.25 ± 1.36	LAB31-SC79	33.56 ± 1.36	161.21 ± 1.36
LAB44	339.61 ± 1.36	47.67 ± 1.36	LAB35-SC79	29.77 ± 1.36	179.74 ± 1.36
LAB53	413.39 ± 1.36	14.09 ± 1.36	LAB44-SC79	40.92 ± 1.36	146.34 ± 1.36
LAB60	400.83 ± 1.36	12.93 ± 1.36	LAB53-SC79	27.42 ± 1.36	173.99 ± 1.36
LAB65	319.57 ± 1.36	34.72 ± 1.36	LAB60-SC79	33.88 ± 1.36	166.83 ± 1.36
LAB69	414.06 ± 1.36	13.55 ± 1.36	LAB65-SC79	49.37 ± 1.36	130.65 ± 1.36
LAB70	399.51 ± 1.36	15.78 ± 1.36	LAB69-SC79	21.89 ± 1.36	174.64 ± 1.36
LAB72	344.21 ± 1.36	30.19 ± 1.36	LAB70-SC79	39.77 ± 1.36	159.82 ± 1.36
LAB76	423.01 ± 1.36	10.68 ± 1.36	LAB72-SC79	40.03 ± 1.36	153.53 ± 1.36
LG	405.3 ± 1.36	17.83 ± 1.36	LAB76-SC79	23.56 ± 1.36	175.11 ± 1.36
SC31	92.04 ± 1.36	133.77 ± 1.36	LG-SC79	28.64 ± 1.36	169.48 ± 1.36
SC79	87.21 ± 1.36	156.22 ± 1.36	LAB24-SE05	51.77 ± 1.36	138.79 ± 1.36
SE05	86.49 ± 1.36	163.06 ± 1.36	LAB31-SE05	47.32 ± 1.36	156.05 ± 1.36

Table 3. Activity of alkaline phosphatase (AKP) in lysate and supernatant of Caco-2 cells incubated with 12 strains of *Lactobacillus* and 3 strains of *Salmonella*.

Treatment Group	Cell Culture Supernatant AKP (U/100 mL)	Cell Lysate AKP (U/gprot)	Treatment Group	Cell Culture Supernatant AKP (U/100 mL)	Cell Lysate AKP (U/gprot)
LAB24-SC31	47.42 ± 1.36	143.35 ± 1.36	LAB35-SE05	18.88 ± 1.36	176.43 ± 1.36
LAB31-SC31	41.32 ± 1.36	150.65 ± 1.36	LAB44-SE05	52.73 ± 1.36	139.03 ± 1.36
LAB35-SC31	39.88 ± 1.36	169.44 ± 1.36	LAB53-SE05	37.11 ± 1.36	159.04 ± 1.36
LAB44-SC31	40.77 ± 1.36	146.59 ± 1.36	LAB60-SE05	37.28 ± 1.36	164.03 ± 1.36
LAB53-SC31	37.97 ± 1.36	153.29 ± 1.36	LAB65-SE05	55.42 ± 1.36	133.57 ± 1.36
LAB60-SC31	40.32 ± 1.36	156.35 ± 1.36	LAB69-SE05	21.32 ± 1.36	170.65 ± 1.36
LAB65-SC31	59.77 ± 1.36	139.74 ± 1.36	LAB70-SE05	29.48 ± 1.36	163.48 ± 1.36
LAB69-SC31	30.56 ± 1.36	171.21 ± 1.36	LAB72-SE05	46.83 ± 1.36	159.52 ± 1.36
LAB70-SC31	35.88 ± 1.36	166.43 ± 1.36	LAB76-SE05	20.36 ± 1.36	173.05 ± 1.36
LAB72-SC31	43.73 ± 1.36	149.53 ± 1.36	LG-SE05	22.63 ± 1.36	168.57 ± 1.36

Table 3. Cont.

The effects of different treatments on LDH metabolism of Caco-2 cells were significantly different (Table 4). Compared with the control group, LDH activity was not detected in Caco-2 cells incubated with SAL, while LDH activity in the supernatant of Caco-2 cells incubated with SAL increased by 64.73% (p < 0.01). It may be that SAL treatment destroys the integrity of Caco-2 cell membrane and releases LDH into cell culture supernatant. The LDH activity in Caco-2 cells incubated with LAB increased significantly (p < 0.05), but decreased by 98.49% (p < 0.01). The LDH activity of Caco-2 cells incubated with LAB-SAL was significantly higher than that of Caco-2 cells incubated with SAL (p < 0.05); the LDH activity content in the supernatant of Caco-2 cells incubated with LAB-SAL increased by 35.26% (p < 0.01), but it was still lower than that of SAL. These results suggest that LAB can enhance the LDH metabolic activity of Caco-2 cells, which further confirms the protective effect of SAL on the integrity of Caco-2 cell membrane by preventing the destruction of Caco-2 cell membrane.

Table 4. Lactate dehydrogenase (LDH) activity in lysate and supernatant of Caco-2 cells incubated with 12 strains of *Lactobacillus* and 3 strains of *Salmonella*.

Treatment Group	Cell Culture Supernatant LDH (U/L)	Cell Lysate LDH (U/gprot)	Treatment Group	Cell Culture Supernatant LDH (U/L)	Cell Lysate LDH (U/gprot)
CT(PBS)	0 ± 0.22	0 ± 0.22	LAB24-SC79	577.32 ± 0.22	30.66 ± 0.22
LAB24	23.1 ± 0.22	18.76 ± 0.22	LAB31-SC79	596.38 ± 0.22	32.54 ± 0.22
LAB31	24.58 ± 0.22	11.27 ± 0.22	LAB35-SC79	561.88 ± 0.22	40.33 ± 0.22
LAB35	39.87 ± 0.22	30.5 ± 0.22	LAB44-SC79	592.37 ± 0.22	21.46 ± 0.22
LAB44	46.25 ± 0.22	31.84 ± 0.22	LAB53-SC79	560.34 ± 0.22	39.27 ± 0.22
LAB53	35.57 ± 0.22	27.09 ± 0.22	LAB60-SC79	603.66 ± 0.22	21.88 ± 0.22
LAB60	40.09 ± 0.22	31.59 ± 0.22	LAB65-SC79	705.39 ± 0.22	18.41 ± 0.22
LAB65	27.89 ± 0.22	17.51 ± 0.22	LAB69-SC79	548.17 ± 0.22	43.67 ± 0.22
LAB69	42.55 ± 0.22	34.27 ± 0.22	LAB70-SC79	613.05 ± 0.22	23.49 ± 0.22
LAB70	33.59 ± 0.22	29.66 ± 0.22	LAB72-SC79	578.02 ± 0.22	28.44 ± 0.22
LAB72	31.67 ± 0.22	21.55 ± 0.22	LAB76-SC79	559.55 ± 0.22	41.81 ± 0.22
LAB76	47.66 ± 0.22	38.78 ± 0.22	LG-SC79	573.22 ± 0.22	36.19 ± 0.22
LG	41.92 ± 0.22	33.76 ± 0.22	LAB24-SE05	651.44 ± 0.22	19.31 ± 0.22
SC31	913.5 ± 0.22	0 ± 0.22	LAB31-SE05	600.23 ± 0.22	23.37 ± 0.22
SC79	885.76 ± 0.22	0 ± 0.22	LAB35-SE05	552.19 ± 0.22	40.55 ± 0.22
SE05	848.27 ± 0.22	0 ± 0.22	LAB44-SE05	636.15 ± 0.22	20.94 ± 0.22
LAB24-SC31	611.22 ± 0.22	19.36 ± 0.22	LAB53-SE05	549.92 ± 0.22	44.52 ± 0.22
LAB31-SC31	605 ± 0.22	20.41 ± 0.22	LAB60-SE05	579.11 ± 0.22	24.68 ± 0.22
LAB35-SC31	598.17 ± 0.22	22.57 ± 0.22	LAB65-SE05	644.24 ± 0.22	19.35 ± 0.22
LAB44-SC31	573.05 ± 0.22	27.49 ± 0.22	LAB69-SE05	590.03 ± 0.22	29.44 ± 0.22
LAB53-SC31	571.55 ± 0.22	33.88 ± 0.22	LAB70-SE05	588.64 ± 0.22	30.39 ± 0.22
LAB60-SC31	608.02 ± 0.22	21.46 ± 0.22	LAB72-SE05	599.65 ± 0.22	27.56 ± 0.22
LAB65-SC31	631.42 ± 0.22	20.31 ± 0.22	LAB76-SE05	553.09 ± 0.22	46.43 ± 0.22
LAB69-SC31	561.09 ± 0.22	32.65 ± 0.22	LG-SE05	564.57 ± 0.22	41.82 ± 0.22
LAB70-SC31	599.14 ± 0.22	21.59 ± 0.22	LAB76-SE05	553.09 ± 0.22	46.43 ± 0.22
LAB72-SC31	578.06 ± 0.22	20.58 ± 0.22	LG-SE05	564.57 ± 0.22	41.82 ± 0.22

3.7. Lactobacillus and Salmonella Isolate Adhesion to Caco-2 Cells and Apoptosis Test Results

Apoptosis was detected using flow cytometry (CytoFLEX LX, BECKMAN, USA). Propidium iodide (PI) is a type of nucleic acid dye. Under normal circumstances, PI cannot penetrate the entire cell membrane; this is not true for mid- and late apoptotic as well as dead cells, the cell membrane of which PI penetrates and dyes the nucleus red. Caco-2 cells treated with PBS were used as blank control group, multidrug resistant Salmonella SC79 was used as negative control group, and standard Lactobacillus strain LG was used as positive control group. Figure 5 shows the four quadrant method to observe the apoptosis of Caco-2 cells, that is, the scatter plot of bivariate flow cytometry. In the Figure 5, P1 gate showed the ratio of living cells, which shows the number of living cells detected by flow cytometry after Caco-2 cells were incubated with different bacteria. Compared with the blank control group, the number of Caco-2 cells in LAB69 group and LG group were significantly increased, while the number of Caco-2 cells in LAB69 group was higher than that in the positive control group; the number of Caco-2 cells in the negative control group was significantly decreased. P2 gate and P3 gate show the apoptotic and dead cells detected by flow cytometry after Caco-2 cells were incubated with different bacteria. Compared with the blank control group, the number of apoptotic cells and dead cells of Caco-2 cells in LAB69 group and LG group decreased, and the number of apoptotic cells and dead cells in LAB69 group was lower than that in the positive control group; while the number of apoptotic cells and dead cells in Caco-2 cells in SC79 group increased significantly. As can be seen from Figure 5, compared with the blank control group, the results of Caco-2 cells incubated with different bacteria were significantly different. The invasion and destruction of Caco-2 cells in Salmonella negative control group were more serious, Compared with positive control group LG and LAB69, the number of live cells increased and the number of dead cells decreased. These results indicate that LAB69 has a certain protective effect on Caco-2.




Figure 5. Four quadrants results of Caco-2 cells treated with different bacteria in flow cytometry. Note: P1 gate, the red part of P1 gate is the ratio of the number of living cells. The higher the red ratio, the more living Caco-2 cells; P2 gate, the green part of P2 gate is the ratio of apoptotic cells. The larger the green area, the more apoptotic cells; P3 gate, The blue part of P3 gate is the ratio of death cells. The larger the blue area, the more death cells.

3.8. Molecular Identification of Lactobacillus Strains

After PCR amplification, electrophoresis results of Lactobacillus amplification products was show in the Figure S2. The sequencing results in Figure 6 and Table 5. In this experiment, 12 strains of chicken derived lactic acid bacteria with good inhibitory effect on *Salmonella*, weak drug resistance, strong bile salt resistance, and strong adhesion to Caco-2 cells were selected for molecular sequencing analysis, according to the experimental results in Table 5, LAB53, LAB60, and LAB64 are the same strain; that is, *Lactobacillus salinus* strain. Except *Lactobacillus reuteri* strain LAB4 and *Enterococcus faecium* strain IAH L76, the rest were *Lactobacillus plantarum* strains.




Figure 6. Neighbor-joining phylogenetic tree based on 16 S rRNA gene sequences showing relationship between *Lactobacillus* isolates.

Strains	Gen Bank Homologous Sequence	The Highest Homologous Strain	Maximum Homology (%)
LAB4	NR_113820.1	Lactobacillus reuteri strain NBRC 15892	99.5%
LAB31	NR_112690.1	Lactobacillus plantarum strain NBRC 15891	99.9%
LAB35	NR_028725.	Lactobacillus salinus strain HO 66	99.8%
LAB44	MF108394.1	Lactobacillus plantarum strain cau6549	99.4%
LAB53	NR_112759.1	Lactobacillus salinus strain JCM 1231	98%
LAB60	NR_112759.1	Lactobacillus salinus strain JCM 1231	99%
LAB64	NR_112759.1	Lactobacillus salinus strain JCM 1231	98%
LAB65	KP317684.1	Lactobacillus salinus strain L13	99.5%.
LAB69	JQ046410.1	Lactobacillus sp. wx213	98%
LAB73	MK311261.1	Lactobacillus plantarum strain AR503	99.1%
L76	MK990071.1	Enterococcus faecium strain IAH_28	98.20%
LAB79	MN389601.1	Lactobacillus plantarum strain 42	98.5%

Table 5. 16S rDNA sequence analysis results of BLAST.

4. Discussion

In this experiment, 89 strains of lactic acid bacteria were isolated from Chicken Intestines of different slaughterhouses according to the national standard method and Gram staining. Using *Salmonella* as indicator bacteria and agar diffusion inhibition zone method, 24 strains of chicken-derived lactic acid bacteria with an average inhibition diameter of more than 10 mm were preliminarily screened. Due to the cost, all 24 strains of chickenderived lactic acid bacteria were not identified by molecular sequencing. The drug resistance of these 24 strains of lactic acid bacteria was analyzed. Except for LAB 76 and LG, the drug resistance rate of 24 strains of lactic acid bacteria was 33%, gentamicin 21%, chloramphenicol 21%, ceftriaxone sodium 25% and tetracycline 25%. LAB12, LAB22, and LAB34 were resistant to three kinds of antibiotics, so the selected lactic acid bacteria were selected from the remaining 21 strains. At the same time, the bile salt tolerance of 24 lactic acid bacteria isolates was detected, the OD600 value of LAB4, LAB7, LAB24, LAB26, LAB31, LAB35, LAB44, LAB56, LAB65, LAB69, LAB70, LAB76, and LG were more than 0.20, and the strains still show good tolerance to the bile salt. The cell adhesion ability of 15 LAB strains differed depending on the species and strain. Murphy (2009) [26] also confirmed great adhesion differences among different strains and even within the same strain. Thus, adhesion ability is suggested to have an inevitable relationship with LAB characteristics. Adhesion can be non-specific (related to surface properties) or specific (related to adhesives) [27]. The existing adhesives include surface proteins, pili, peptidoglycans, and lipopolysaccharides. Studies have shown that the S-layer protein on the surface of Lactobacillus strains may be the active adhesion site [28]. Based on the aforementioned considerations, the adhesion of Lactobacilli to Caco-2 cells may result from the orchestrated action of various adhesives. The mechanism of *Lactobacillus* antagonizing the adhesion of *Salmonella* to Caco-2 cells was explored by three ways of exclusion, competition and substitution; the protective effect of *Lactobacillus* on intestinal cells was preliminarily analyzed by flow cytometry to detect the permeability of cell membrane; and the biological function of Lactobacillus SLP was preliminarily explored These results suggest that LAB group bacteria protected the AKP metabolism activity and Caco-2 cell membrane integrity by reducing the damage caused by the SAL group to the Caco-2 cell membrane.

Our results showed that the bacteriostatic zones of LAB69 and L76 were the widest (17.2 mm and 17.5 mm, respectively). Additionally, most of the 24 *Lactobacillus* strains inhibited the adhesion of *Salmonella* to the Caco-2 cell surface bacteria effectively. Moreover, after attaching to the cell surface, *Lactobacillus* effectively controlled the reattachment of *Salmonella* bacteria to cells. Further, 16 *Lactobacillus* strains could replace *Salmonella* strains effectively and adhere to the surface of Caco-2 cells. In general, incubation with *Lactobacillus* had protective effects on Caco-2 cell membranes, and the detection results were close to that of the control group.

5. Conclusions

In this experiment, we isolated chicken-derived *Lactobacillus*, which can effectively inhibit the colonization and adhesion of *Salmonella*, from chicken intestines. The LAB53, LAB60, LAB69, LAB72, and LAB76 *Lactobacillus* isolates had slightly higher tolerance and adhesion capacity, as well as probiotic potential than the standard LG strain; however, the probiotic effect of *Lactobacillus* isolates on the host after adhering to the intestinal tract should be studied further. We speculate that the incorporation of chicken-derived LAB during the later stage of chicken feeding might curtail *Salmonella* outbreaks in the breeding and production processes of chicken eggs, ultimately improving food safety, Additionally, the LAB69 and LAB76 isolates had good adhesion ability and may be rich in surface proteins. Their adhesion properties and mechanism need to be studied in detail in the future.

Supplementary Materials: The following are available online at https://www.mdpi.com/2304-8 158/10/3/569/s1, Figure S1: Gram staining results of strains isolated from *Lactobacillus*., Figure S2: Electrophoresis results of *Lactobacillus* amplification products., Table S1: Colony characteristics and Bacterial forms of different strains., Table S2: Test results of physiological and biochemical characteristics of *Lactobacillus*.

Author Contributions: Conceptualization and study design: D.H. and X.B.; Performing experiments: D.H., Z.L. and X.H.; Review and final edit: F.L. All authors have read and agreed to the published version of the manuscript.

Funding: This work was supported by the National Key R&D Program of China, Grant No. 2018YFC1602500, NSFC (Grant No. 31972174), the Agricultural Independent Innovation Fund of Jiangsu Province (Grant number CX (18)3053), and the Priority Academic Program Development of Jiangsu Higher Education Institutions (PAPD).

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Conflicts of Interest: The authors declare no competing financial interest.

References

- Nagappa, K.; Tamuly, S.; Saxena, M.; Singh, S. Isolation of salmonella typhimurium from poultry eggs and meat of tarai region of uttaranchal. *Indian J. Biotechnol.* 2007, *6*, 407–409.
- Jalali, M.; Abedi, D.; Pourbakhsh, S.A.; Ghoukasin, K. Prevalence of Salmonella spp. in raw and cooked foods in Isfahan-Iran. J. Food Saf. 2008, 28, 442–452. [CrossRef]
- 3. Akiba, M.; Kusumoto, M.; Iwata, T. Rapid identification of Salmonella enterica serovars, typhimurium, choleraesuis, infantis, hadar, enteritidis, dublin and gallinarum, by multiplex PCR. *J. Microbiol. Methods* **2011**, *85*, 9–15. [CrossRef] [PubMed]
- 4. Tellez, G.; Pixley, C.; Wolfenden, R.E.; Layton, S.L.; Hargis, B.M. Probiotics/direct fed microbials for Salmonella control in poultry. *Food Res. Int.* **2012**, *45*, 628–633. [CrossRef]
- Chambers, J.R.; Gong, J. The intestinal microbiota and its modulation for Salmonella control in chickens. *Food Res. Int.* 2011, 44, 3159. [CrossRef]
- Revolledo, L.; Ferreira, A.J.P.; Mead, G.C. Prospects in salmonella control: Competitive exclusion, probiotics, and enhancement of avian intestinal immunity. J. Appl. Poult. Res. 2006, 15, 341–351. [CrossRef]
- Sadler, W.W.; Brownell, J.R.; Fanelli, M.J. Influence of age and inoculum level on shed pattern of Salmonella typhimurium in chickens. *Avian Dis.* 1969, 13, 793–803. [CrossRef]
- 8. Gast, R.K.; Holt, P.S. Persistence of Salmonella enteritidis from one day of age until maturity in experimentally infected layer chickens. *Poult. Sci.* **1998**, 77, 1759–1762. [CrossRef]
- 9. Yang, B.; Qu, D.; Zhang, X.; Shen, J.; Cui, S.; Shi, Y.; Xia, M.; Shenga, M.; Zhi, S.; Mengad, J. Prevalence and characterization of Salmonella serovars in retail meats of marketplace in Shaanxi, China. *Int. J. Food Microbiol.* **2010**, *141*, 63–72. [CrossRef]
- 10. Yin, M.; Yang, B.; Wu, Y.; Wang, L.; Wu, H.; Zhang, T.; Tuohetaribayi, G. Prevalence and characterization of Salmonella enterica serovar in retail meats in market place in Uighur, Xinjiang, China. *Food Control* **2016**, *64*, 165–172. [CrossRef]
- 11. Chou, S.L.; Weimer, B. Isolation and characterization of acid- and bile-tolerant isolates from strains of lactobacillus aci-dophilus. *J. Dairy Sci.* **1999**, *82*, 23–31. [CrossRef]
- 12. Xu, Z. Update of national inspection standard for lactic acid bacteria. Mod. Food 2018, 4, 47–49.
- 13. Yin, Q.; Zheng, Q. Isolation and identification of the dominant lactobacillus in gut and faeces of pigs using carbohydrate fermentation and 16S rDNA analysis (microbial physiology and biotechnology). *J. Biosci. Bioeng.* **2005**, *99*, 68–71. [CrossRef]
- 14. Hai, D.; Yin, X.; Lu, Z.; Lv, F.; Zhao, H.; Bie, X. Occurrence, drug resistance, and virulence genes of Salmonella isolated from chicken and eggs. *Food Control* **2020**, *113*, 107109. [CrossRef]
- 15. Aween, M.M.; Hassan, Z.; Muhialdin, B.J.; Noor, H.M.; Eljamel, Y.A. Evaluation on antibacterial activity of lactobacillus acidophilus strains isolated from honey. *Am. J. Appl. Sci.* **2012**, *9*, 807–817. [CrossRef]
- 16. Fothergill, A.W. Antifungal susceptibility testing: Clinical Laboratory and Standards Institute (CLSI) methods. In *Interactions of Yeasts, Moulds, and Antifungal Agents;* Humana Press: Totowa, NJ, USA, 2012.
- 17. Lei, R.; Wu, H.; Yongxing, H. Screening of lactic acid bacteria with high antibacterial activity from Campylobacter. *Chin. Anim. Husb. Vet.* **2013**, *40*, 129–134.
- 18. Stuart, M.R.; Chou, L.S.; Weimer, B.C. Influence of carbohydrate starvation and arginine on culturability and amino acid utilization of lactococcus lactis subsp. lactis. *Appl. Environ. Microbiol.* **1999**, *65*, 665–673. [CrossRef] [PubMed]
- 19. Jankowska, A.; Laubitz, D.; Antushevich, H.; Zabielski, R.; Grzesiuk, E. Competition of Lactobacillus paracasei with Salmonella enterica for Adhesion to Caco-2 Cells. *J. Biomed. Biotechnol.* **2008**, 2008, 357964. [CrossRef] [PubMed]
- Adlerberth, I.; Ahrne, S.I.V.; Johansson, M.L.; Molin, G.; Hanson, L.A.; Wold, A.E. A mannose-specific adherence mechanism in Lactobacillus plantarum con-ferring binding to the human colonic cell line HT-29. *Appl. Environ. Microbiol.* 1996, 62, 2244–2251. [CrossRef]
- 21. Chu, Y.H.; Yu, X.X.; Jin, X.; Wang, Y.T.; Zhao, D.J.; Zhang, P.; Sun, G.-M.; Zhang, Y.-H. Purification and characterization of alkaline phosphatase from lactic acid bacteria. *RSC Adv.* **2018**, *9*, 354–360. [CrossRef]
- 22. Jensen, H.; Grimmer, S.; Naterstad, K.; Axelsson, L. In vitro testing of commercial and potential probiotic lactic acid bacteria. *Int. J. Food Microbiol.* **2012**, 153, 216–220. [CrossRef] [PubMed]
- 23. Fernández, M.F.; Boris, S.; Barbés, C. Probiotic properties of human lactobacilli strains to be used in the gastrointestinal tract. *J. Appl. Microbiol.* 2003, *94*, 449–455. [CrossRef] [PubMed]
- Tuo, Y.; Yu, H.; Ai, L.; Wu, Z.; Guo, B.; Chen, W. Aggregation and adhesion properties of 22 Lactobacillus strains. *J. Dairy Sci.* 2013, 96, 4252–4257. [CrossRef]
- 25. Sidira, M.; Kourkoutas, Y.; Kanellaki, M.; Charalampopoulos, D. In vitro study on the cell adhesion ability of immobilized lactobacilli on natural supports. *Food Res. Int.* 2015, *76*, 532–539. [CrossRef] [PubMed]
- Murphy, L. Assessment of potential probiotic strains: Evaluation of their establishment, persistence, and localisation in the mu-rine gastrointestinal tract. *Microb. Ecol. Health Dis.* 2009, 11, 149–157.
- 27. García-Cayuela, T.; Korany, A.M.; Bustos, I.; Gómez de Cadiñanos, L.P. Adhesion abilities of dairy Lactobacillus plantarum, strains showing an aggre-gation phenotype. *Food Res. Int.* **2014**, *57*, 44–50. [CrossRef]
- 28. Taverniti, V.; Stuknyte, M.; Minuzzo, M.; Taverniti, V.; Stuknyte, M.; Minuzzo, M.; Arioli, S.; De Noni, I.; Scabiosi, C.; Martinez Cordova, Z.; et al. MIMLh5 on innate immunity. *Appl. Environ. Microbiol.* **2013**, *79*, 1221–1231. [CrossRef] [PubMed]