



# Internal Validation of the GlobalFiler STR Kit

24th International Symposium


On

**Human Identification** 

Atlanta, Georgia 2013

Jason Kokoszka, Ph.D.
Forensic Biology Section Chief
Alabama Department of Forensic Sciences
Mobile Regional Laboratory

## The GlobalFiler™ Kits: 6-Dye Configuration



GlobalFiler™ Kit 10 Mini-STRs PI\* = 9.2 x 10<sup>-12</sup> MiniFiler ™ Kit 8 Mini-STRs PI\* = 1.05 x 10<sup>-10</sup> 97% of alleles in allelic ladder <400bp

Life Technologies\*\* | 7

\* Preliminary calculations based on a US Caucasian Database

life

# FBI's Quality Assurance Standards for Forensic DNA Testing Laboratories

- 8.3.1 For Internal Validation Studies:
  - a. Have internal validation studies been documented and summarized?
  - Have all internal validation studies conducted on or after July 1, 2009, included, as applicable:



N/A



1.	Known and non pa			lence	samp	les or	
	mock evidence sa	Yes	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	No	П	N/A	
2.	Reproducibility an						
		Yes		No		N/A	
3.	Sensitivity and sto	chasti	c stud	ies?			
		Yes		No		N/A	
4.	Mixture studies?						
		Yes		No		N/A	
5.	Contamination as	sessm	ent?				

Yes

## **Internal Validation Studies Background**

**Current Test Kit:** 

**Detection Platform:** 

3130 Analytical Threshold:

**Identifiler STR Kit** 

3130 Genetic Analyzer

75 rfu



**Test Kit:** 

**Detection Platform:** 

GlobalFiler STR Kit

3500 Genetic Analyzer

(8 capillary)

3500 Analytical Threshold: 175 rfu – starting point



## **Known and Non-Probative Evidence Samples**

Previously extracted non-probative evidence samples, which mimic other sample types routinely encountered in forensic casework, were tested with the GlobalFiler STR Kit. The results obtained using GlobalFiler were compared to those previously obtained with the Identifiler STR Kit.

#### The sample types included in this study included:

- -sexual assault samples (VS, condom, bedding, and panties)
- -wearer items
- -cigarette butts
- -bottle/can swabs
- -blood swabs
- -degraded/inhibited samples
- -reference samples (buccal swabs and blood cards)



## Known and Non-Probative Evidence Samples

#### **GlobalFiler:**

- tests 24 loci
- 21 loci are autosomal STRs
- -12 loci < 250 bp

#### **Identifiler:**

- tests 16 loci
- -15 loci are autosomal STRs
- -9 loci < 250 bp

#### variables:

- input volume differences (10ul vs 15ul)
- chemistry differences between the GF and ID kits?
- 3500 and 3130 differences?
  - injection parameters
  - analysis algorithms




- Good concordance between GlobalFiler and Identifiler
- More information from GlobalFiler due to the additional dye and the number of loci < 250 bp
- GlobalFiler more effective in overcoming inhibition and degradation

#### 1) Semen stain from a quilt

- 0.041 ng/μl
- 0.41 ng amplified with Identifiler
- maximum injection parameters on the 3130 (2 μl/ 9 secs)
- partial male DNA profile
   (10 loci of data with 7 of the CODIS core)

## GlobalFiler



- 0.62 ng amplified with GlobalFiler
- standard injection parameters on the 3500 (1 μl/15 secs)
- full male DNA profile (13 CODIS core)

#### GATICGA AGCATICA AGCATICA


### **Known and Non-Probative Evidence Samples**

#### 2) Semen stain from a pair of panties

- 0.0091 ng/μl
- 0.09 ng amplified with Identifiler
- maximum injection parameters on the 3130 (2  $\mu$ l/ 9 secs)
- partial male DNA profile
- 12 loci of data with 9 of the CODIS core



#### GlobalFiler




- 0.14 ng amplified with GlobalFiler
- standard injection parameters on the 3500 (1 μl/15 secs)
- partial male DNA profile (23 loci of data and 12 CODIS core)


## **Known and Non-Probative Evidence Samples**

- 3) Swabbing of a concrete block
  - $0.17 \text{ ng/}\mu\text{l}$
  - 1 ng amplified with Identifiler
  - standard injection on the 3500 (1 μl/15 secs)
  - no results

inhibition suspected due to extract color



#### GlobalFiler




- 1 ng amplified with GlobalFiler
- standard injection parameters on the 3500 (1 μl/15 secs)
- partial DNA profile (19 loci of data; 10 CODIS core)

## Reproducibility, Precision, and Accuracy


**Precision:** characterizes the amount of agreement among a series of individual (or <u>repeated</u>) measurements

-Precision does not relate to the accuracy of the measurements

Accuracy: describes how close the measured values are to the true value








## Reproducibility and Precision Study Design

GlobalFiler allelic ladder was prepared in one column and injected eight (8) times on the 3500 Genetic Analyzer

-size determinations were conducted on all alleles of the ladder for each injection


-standard deviations in base pair sizing were calculated for each allele in the ladder



GlobalFiler allelic ladder – 343 alleles across
24 loci
2,744 alleles per 8-cap injection

21,952 alleles per run of 8 injections

#### Reproducibility and Precision of GlobalFiler Allelic Ladder




-standard deviation was averaged across all the alleles at each locus

#### **Accuracy Study Design**


One (1) ng of the 007 Human DNA control was amplified concurrently four times with each amplified product prepared for electrophoresis in triplicate for a total of twelve (12) samples. The twelve (12) samples were injected on three separate days generating data for thirty-six (36) samples across three separate runs.

**Accuracy:** describes how close the measured values are to the true value


The measurement of each allele in the 007 sample was compared to the measurement of the corresponding allele in the GlobalFiler allelic ladder




## **Accuracy Study Day 1**



## **Accuracy Study Day 2**



## **Accuracy Study Day 3**



### Sensitivity and Stochastic Study Design



The following amounts of the positive control 007 were amplified, *in quadruplicate*, using the Globalfiler STR Kit to assess the sensitivity and stochastic effects of the GlobalFiler STR kit on the 3500 Genetic Analyzer:

15.625 pg, 31.25 pg, 62.5 pg, 125 pg, 250 pg, 500 pg, and 1 ng

The above samples were injected on the 3500 Genetic Analyzer using the standard injection parameters (1  $\mu$ l/15 secs).




## **Sensitivity and Stochastic Study**

#### 007 dilution series (43 total alleles)

ng/μl	pg/µl	input template	average # of alleles detected	
0.1	100	1 ng	43	
0.05	50	500 pg	43	
0.025	25	250 pg	43	
0.0125	12.5	125 pg	42	
0.00625	6.25	62.5 pg	32.3	
0.003125 3.125		31.25 pg	6.3	
0.0015625	1.5625	15.625 pg	1	

## **Sensitivity and Stochastic Study**



#### **Mixture Study Design**





The positive controls 9947A and 007 were mixed in the following ratios:

19:1, 9:1, 4:1, 3:2, 1:1, 1:0, 0:1, 1:1, 2:3, 1:4, 1:9, and 1:19

Approximately one nanogram of each mixture was amplified with the GlobalFiler STR Kit, in triplicate, and then injected on the 3500 Genetic Analyzer (1  $\mu$ l/15 secs and 2  $\mu$ l/20 secs).

#### **Mixture Studies**


#### Standard injection parameters (1 µl\15 secs)

	007 (26 unique alleles)		9947A (20 unique alleles)		
007:9947A	# 007 alleles detected	% complete 007	# 9947A alleles detected	% complete 9947A	
19:1	26	100 %	13	65 %	
9:1	26	100 %	19.7	98.5 %	
4:1	26	100 %	20	100 %	
3:2	26	100 %	20	100 %	
1:1	26	100 %	20	100 %	
2:3	26	100 %	20	100 %	
1:4	26	100 %	20	100 %	
1:9	24	92.3 %	20	100 %	
1:19	16.3	62.7 %	20	100 %	

#### Maximum injection parameters (2 μl\20 secs)

	007 (26 unic	ue alleles)	9947A (20 unique alleles)		
007:9947A	# 007 alleles detected	% complete 007	# 9947A alleles detected	% complete 9947A	
19:1	26	100 %	17	85 %	
9:1	26	100 %	20	100 %	
1:9	24.7	95 %	20	100 %	
1:19	20	76.9 %	20	100 %	

## The GlobalFiler™ Kits: Configuration



Life Te cha diagies™ | 12

Applied Pringerens\*

(Filtra)

Malacalan Property


Hesta"

Tagitter




#### What is an Indel?

- Indels are short length polymorphisms, consisting of the presence (INsertion)
  or absence (DELetion) of a short sequence (1-50 nucleotides)
- In the GlobalFiler™ Kits, an Indel marker on the Y chromosome was incorporated into the multiplex




#### **Mixture Studies**



### **Contamination Assessment Study Design**

To evaluate the possibility of carry-over signal from injection to injection, a variety of configurations of excessive GS600-LIZ size standard (4X) and blank formamide were processed on the 3500



Injections occurred from left to right with single injections for columns 1-4.

Columns 5 and 6 injected three (3) times each consecutively.


Single injections for columns 7-12.

Form only


Form + ESS

### **Contamination Assessment**

**GS 600 LIZ** 




Blank




#### **Contamination Assessment**









## **Final Thoughts**

- 1) Reproducibility and precision of the GlobalFiler STR Kit on the ABI 3500: the standard deviation across multiple allelic ladders injections was  $\leq 0.05$  bp.
- 2) Accuracy of the GlobalFiler STR Kit on the ABI 3500: with the exception of D8S1179, all alleles detected from the 007 sample across three separate 3500 runs were within ± 0.25 bp of the allelic ladder measurement. The alleles detected at D8S1179 were within ± 0.30 bp.
- 3) <u>Contamination assessment on the ABI 3500</u>: unexpected signal was detected in the formamide-only blanks no carry-over detected in any amplification blanks employed in the rest of the validation studies.
- 4) <u>Mixture studies</u>: significant minor component could be detected in the 1:19 and 19:1 mixture samples.
- 5) Sensitivity and stochastic studies: full DNA profiles were obtained from 125 pg of template and significant data were obtained from 62.5 pg of template.
- 6) <u>Casework samples</u>: concordance between DNA profiles previously obtained with Identifiler. The GlobalFiler STR Kit was successful in overcoming both inhibition and degradation.

## Acknowledgements

Jacquelyn M. Jenkins, PhD, MSFS Forensic Scientist Birmingham Regional Laboratory Carl Mauterer, MSFS
Forensic Biology Section Chief
Montgomery Regional Laboratory



#### **Alabama Department of Forensic Sciences**

Angelo Della Manna, MSFS Assistant Director Chief of Forensic Biology & DNA A.Dellamanna@adfs.alabama.gov Jason E. Kokoszka, PhD
Forensic Biology Section Chief
Mobile Regional Laboratory
Jason.Kokoszka@adfs.alabama.gov





## Questions?



life technologies™



Go faster. Go further. Go global.

## **GlobalFiler**

Powered by 6-Dye™