



1 Dihydroergocornine mesylate

2 Dihydro- $\alpha$ -ergocryptine mesylate

3 Dihydroergocristine mesylate

4 Dihydro- $\beta$ -ergocryptine mesylate

#### HPLC Conditions

Column CAPCELL PAK C<sub>18</sub> UG120 S5 : 4.6 mm i.d.×250 mm

Mobile phase H<sub>2</sub>O / CH<sub>3</sub>CN / Triethylamine = 70 / 30 / 3

Flow rate 1.0 mL/min

Temperature 30 °C

Detect UV 280 nm